

# LOGIQ™ P10 XDclear™ Performance Series


Make it easy. Make it your own.



GE HealthCare

LOGIQ and XDclear are trademarks of GE HealthCare.

© 2023 GE HealthCare. GE is a trademark of General Electric Company used under trademark license.



# Challenging ultrasound environment

## Time & quality challenges

- Growing patient volume
- Need for fast, accurate diagnosis

## Operator essentials

- User-friendly equipment
- Avoid work related injuries

## Economic pressures

- Budget constraints
- Declining reimbursement
- ROI and life cycle costs

## Physical environment

- Easily disinfected and cleanable
- Cramped exam rooms
- Increasing need for portable exams

## Patient expectations

- Well-informed patients
- Comfort and safety



# GE HealthCare Ultrasound innovation

Designed with your patients in mind

Radiology/  
Vascular



LOGIQ™ Family

Cardiovascular



Vivid™ Family

Women's Health



Voluson™ Family

Breast  
Screening



Invenia™ ABUS 2.0

Point of Care



Venue™ Family

Hand-held



Vscan™ Family

# GI Ultrasound

Continued innovation



LOGIQ™ P10 XDclear/  
P9 XDclear/P8



LOGIQ S8  
XDclear™ 2.0



LOGIQ E10s



LOGIQ E10



# LOGIQ™ P10 XDclear™

Make it easy. Make it your own.

## Personalized workflow tools and automation...for efficiency

The new LOGIQ P10 XDclear ultrasound system brings you advanced ways to increase everyday efficiency—tools and functions you can customize to your own preferences, so exams flow easily, your way.



## Expanded patient-centric diagnostic capabilities...for flexibility

Equipped with a robust set of advanced features, the LOGIQ P10 XDclear system gives you the flexibility to manage simple to complex cases.



## Powerful support...for long-term value

The LOGIQ P10 XDclear system gives you access to industry-leading security, data management, and educational resources to help you optimize system uptime and utilization, while enhancing user skills.



Advanced features. Everyday affordability.

# Organizing schema for LOGIQ™ P10 XDclear™ (R4)

Make it easy. Make it your own.

## Personalized workflow tools and automation...for efficiency

- Touch Control, UI – Joystick, LOGIQ brand usability
- Ergonomics, transportable, button probe
- LOGIQ apps: Photo Assistant & Remote Control
- My Page, Start Assistant, My Preset
- Power Assistant
- Vnav Import
- Auto tools: Auto TGC, AutoEF, Auto IMT, Measure Assistant, Compare Assistant, Scan Assistant, Measurement enhancements
- CTO
- SonoNT, SonoIT
- AI tools: Auto Lesion Segmentation, OB Measurement Assistant, Auto Doppler Assistant, Breast Assistant, powered by Koios DS™

## Expanded patient-centric diagnostic capabilities...for flexibility

- Excellent IQ, LCD display, HDU display
- Advanced features: SWE enhancements & Strain, Raw Data, CEUS, B-Flow, B-Flow Color, HD Color, HDlive™, Stress Echo, STIC/Omniview, TVI/TVD, Cardiac Strain, B-Steer +, urology dual display, UGAP, SRI, 3D/4D enhancements with SonoRenderlive
- XDclear™ probes
- Sono AVC
- Urology enhancement
- CF/PW Enhancement
- Multi-purpose – Abdomen, Liver, Cardiac, OB/GYN, Vascular, Breast, Thyroid, MSK, Urology, Pediatrics

## Powerful support...for long-term value

- My Trainer
- Data Management
- Tricefy™
- SonoDefense enhancements
- Life Cycle Solutions

# Personalized workflow



# LOGIQ™ P10 XDclear™

Personalized workflow tools and automation for efficiency

## Personalized set-up tools

- Users can customize their own workflow preferences and launch in seconds
  - Start Assistant
  - My Page
  - My Preset

## AI-based tools

- Harnesses artificial intelligence for imaging standardization and speed
  - Auto Lesion Segmentation
  - OB Measurement Assistant
  - Auto Doppler Assistant
  - Breast Assistant powered by Koios DS™\*

## Automated scanning tools

- Help reduce exam time and increase user efficiency
  - Continuous Tissue Optimization (CTO)
  - Auto IMT
  - AutoEF
  - Measure Assistant
  - Compare Assistant
  - Scan Assistant

## LOGIQ brand usability

- Redesigned user interface — with customizable keys — offers the simplicity and ease of operation you've come to expect





# Patient-centric diagnosis



# LOGIQ™ P10 XDclear™

Helps you improve patient care

## Superb image quality with XDclear probes

- Powerful high fidelity and broad bandwidth produce high resolution images
- Scan superficial or deep targets – or any point in between
- XDclear probes include C1-6-D, C3-10-D and M5Sc-RS
- New probes include C2-7-D, 10C-D

## Advanced imaging and visualization tools<sup>1</sup>

- Ultrasound-Guided Attenuation Parameter (UGAP) a real-time, image-guided method to measure fatty content in the liver tissue and assess non-alcoholic fatty liver disease (NAFLD)
- 2D Shear Wave Elastography with Quality Indicator and Anatomical Site Reporting
- 3D/4D with SonoRender/live

## Multi-purpose capabilities for a diverse caseload

- Including liver, cardiac, OB/GYN, breast, thyroid, musculoskeletal, urologic, vascular and pediatric studies




# Spotlight on liver

## 2D Shear Wave Elastography with Quality Indicator, C1-6-D

LOGIQ™ P10

Quantitative measurement of tissue elasticity displayed in color-coded elastograms. Useful in liver, breast and small parts applications.




# Spotlight on MSK

## 2D Shear Wave Elastography with Quality Indicator, patellar tendon, ML6-15-RS

LOGIQ™ P10

Quantitative measurement of tissue elasticity displayed in color-coded elastograms. Useful in liver, breast and small parts applications.






# Spotlight on GYN

## 2D Shear Wave Elastography with Quality Indicator, cervix, IC9-RS

LOGIQ™ P10


Quantitative measurement of tissue elasticity displayed in color-coded elastograms.



# Spotlight on 3D rendering

## Gynecological assessment tools

Advanced Volume Contrast Imaging (VCI) with OmniView – Diagnostic confidence in sonography requires the ability to differentiate irregular shapes with precision. This tool can help improve contrast resolution and visualization of the rendered anatomy with clarity in any image plane.




# Spotlight on pediatrics

Neonatal head with XDclear™ transducer technology, C3-10-D

LOGIQ™ P10


Superb image quality and transducer technology combine to help clinicians quickly diagnose a wide range of pediatric conditions.




# Abdominal



Liver B-Mode, C1-5-RS




B-Flow dual liver, C2-7-D




Appendix, L3-12-RS




# Liver




Liver SWE with Quality Indicator, C1-6 -D



Liver CEUS, 9L-RS




Liver UGAP, C1-6-D



Liver contour using Virtual Convex, L3-12-RS


# Vascular



Auto IMT measurement of CCA using L3-12-RS




Carotid B-Flow HD Color, 12L-RS



Carotid hemodynamics with PW using L4-12t-RS


# Small parts



Thyroid PDI, ML6-15-RS




Testis Virtual Convex, L3-12-RS




Thyroid Strain Elasto, L3-12-RS


# Musculoskeletal



Supraspinatus tendon, ML6-15-RS




2D Shear Wave Elastography with Quality Indicator,  
patellar tendon, ML6-15-RS




Flexor tendon of finger, ML6-15-RS



# Musculoskeletal




CF of wrist with Virtual Convex, 12L-RS




Extensor tendon of wrist, ML6-15-RS


# Cardiology



TVI TVD, M5Sc-RS



Strain image of 4CH, 3Sc-RS




Mitral valve, 6Tc-RS

# OB/GYN



3D HDlive™ imaging of the fetus with RAB2-6-RS



2D Shear Wave Elastography with Quality Indicator, cervix, IC9-RS




4CH heart, C1-6-D


# Pediatrics



Neonatal head, C3-10-D




Neonatal spine B-Mode, 12L-RS




Neonatal abdomen, ML6-15-RS

# Breast



Breast Strain Elastography, ML6-15-RS



Breast, ML6-15-RS




# Interventional




Liver, C1-5-RS and CT abdominal scan



Shoulder bursa injection, ML6-15-RS




MSK B-Steer+, 12L-RS



Wrist injection, L8-18i-RS

# B-Flow imaging



B-Flow HD Color thyroid L3-12 -RS



B-Flow dual LEV, L3-12-RS



Neonatal head B-Flow, 9L-RS

# Powerful support



# LOGIQ™ P10 XDclear™

Advanced data security and educational resources

Optimize system uptime and utilization,  
while enhancing user skills

## SonoDefense

Built on the Windows® 10 Internet of Things (IoT) operating system, SonoDefense provides multi-layer security to protect system integrity and patient data privacy.

## My Trainer

Self-training modules — available on-board, online or in the app — help new users become proficient quickly.

## Data management tools

Easy data migration supported by DICOM® and cloud-based image collection with anonymization.

LOGIQ and XDclear are trademarks of GE HealthCare. Windows is a registered trademark of Microsoft Corporation. DICOM is a registered trademark of the National Electrical Manufacturers Association.



LOGIQ P10 XDclear Customer Presentation



© 2023 GE HealthCare. GE is a trademark of General Electric Company used under trademark license.

29

# Customer care





# LOGIQ™ P10 XDclear™

Make it easy. Make it your own.

Tailored solutions<sup>1</sup> to help you address your daily challenges

## Education

- My Trainer – On-board training modules
- LOGIQ Club community and product educational services
- Access to additional purchasable online and in-person education

## Service and care

- Three years of service included with purchase
- InSite™ remote support
- Dedicated financial solutions
- Probe care
- <http://www.gehealthcare.com/transducers>



1. Refer to the LOGIQ P10 XDclear Product Data Sheet for a list of purchasable options. LOGIQ, XDclear, and InSite are trademarks of GE HealthCare.

© 2023 GE HealthCare. GE is a trademark of General Electric Company used under trademark license.

# LOGIQ™ P10 remote support: InSite™

**Designed for:** Clinical and biomed staff

**Problem:** Waiting for help on technical or applications problems impacts patients and operations

**Solution:** Fast live response from OEM experts with InSite connection

InSite remote connectivity can help maximize your ultrasound uptime and investment. When you need help, our most experienced expert sonographers and engineers can electronically link to your system and quickly get you back up and running. **InSite experts have an average of 15-20 years experience.\*** They resolve most issues remotely and reduce onsite time by proactively diagnosing problems.



Access GE HealthCare ultrasound experts directly with **one-touch** on your system



60% of ultrasound systems are repaired remotely, and often in **less than 20 minutes**



If your ultrasound cannot be repaired remotely, remote diagnostics can reduce on site repair time by **up to 70%**



**Console**  
with “Contact  
GE” icon



**Phone**  
In the USA,  
call 800-437-1171  
and enter your  
system ID



**Web**  
iCenter users can  
create a service  
request online

Fast live response from GE HealthCare clinical and service experts




# LOGIQ™ P10 productivity: Performance analytics

**Designed for:** Ultrasound department manager with multiple assets

**Problem:** Need visibility to exam volumes and throughput to optimize utilization of assets, staffing, and patient scheduling

**Solution:** iCenter™ and Imaging Insights† uncover opportunities to improve use of assets and personnel



## iCenter\*

Built-in, web-based

Included w/ contract

GE HealthCare product utilization & service data

More info [here](#)

## Imaging Insights†

On-premise app

Install w/ subscription

GE HealthCare and multi-vendor asset & exam data

More info [here](#)



**Example:** Renal exams take longer than scheduled time. Exam history shows most longer exams occur Tuesday mornings. Action is taken to review of Renal exam training and procedures with Tuesday AM staff.

**Example:** Renal exams take longer than scheduled time. Analysis tools show scan time is 10 min longer than planned and longer exams do not use scan assistant. Actions are taken to extend renal exam schedule slots with targeted training on using scan assistant to improve workflow & protocol.

Optimize assets, capacity and scheduling with iCenter\* or Imaging Insights†

# LOGIQ™ P10 biomed support: System Health Dashboard\*


**Designed for:** Biomed and equipment maintenance personnel

**Problem:** Diagnosing problems on a complex ultrasound device can be time consuming

**Solution:** Advanced diagnostics with a simple visual interface can help reduce repair time and improve asset availability

Visual status of key subsystems, intelligent alerts, and click-through details make it easier to manage preventative and corrective maintenance for qualified users. Includes new diagnostics for connected probes that run-in air (no probe manipulation).

Dashboard features are enabled remotely and are included with proprietary service license (offered with in house support contracts).



Intelligent visual diagnostics simplify maintenance and support biomed productivity

# LOGIQ™ P10 XDclear™ Performance Series

Make it easy. Make it your own.

Its sleek, lightweight design incorporates powerful technologies while enhancing workflow and user experience.

All at an affordable price that helps meet budget expectations.







# Personalized workflow



# My Page


## Configurable on each My Page

- GE HealthCare innovative personalization feature
- Preset depended (Abd/Carotid/...has My Page)
- Each use case (preset) has My Page


### Abd




### Carotid



### OB



### Breast



# Start Assistant




## Current workflow

- 1 Select Patient from Worklist  
*Auto Populate Exam Description*
- 2 Select Exam Category Tab
- 3 Select Scan Assistant
- 4 Exit Patient Screen
- 5 Select Probe
- 6 Select Application

## With Start Assistant


- 1 Select Patient from Worklist  
*Auto Populate Exam Description*  
*Auto Select Exam Category Tab*  
*Auto Select Scan Assistant*
- 2 Exit Patient Screen  
*Auto Select Probe*  
*Auto Select Application*




# My Preset

With My Preset, users can customize their own workflow preferences and use case presets, and then launch these settings in seconds.


Utility page



My Preset off



My Preset on






# Additional configurable capability



1: 7 physical user defined keys

2: Freeze & P1  
*Exchangeable*

5 digital user defined keys



Configurable smart keys

Imaging mode	Function
Live B, Con, B-Flow	Focal Zone Up/Down
	Frequency Up/Down
CF, TVI, B-Flow Color, PDI	Box Steer
	Scale (PRF) Up/Down Auto Doppler Assistant
PWD, CW, TVD	Baseline Up/Down
	Scale (PRF) Up/Down Auto Doppler Assistant
M/D cursor	Sample volume size

# Digital TGC control assigned to rotary knobs



# Ergonomic adaptable system<sup>1</sup>



# Ergonomic adaptable system<sup>1</sup>




# Ergonomic adaptable system<sup>1</sup>



**Rotate**  
+/- 30 degrees




**Up/down**  
810 – 910 mm





# Ergonomic adaptable system<sup>1</sup>



200 mm

200 mm



1. Refer to the LOGIQ™ P10 XDclear Product Data Sheet for a list of purchasable options. LOGIQ is a trademark of GE HealthCare.

© 2023 GE HealthCare. GE is a trademark of General Electric Company used under trademark license.



# LOGIQ™ brand usability

The redesigned user interface — with customizable keys — delivers the simplicity and ease of operation you've come to expect from LOGIQ systems.




# Additional configurable capability



1: 7 physical user defined keys

2: Freeze & P1  
Exchangeable

4 digital user defined keys



Configurable smart keys

Imaging mode	Function
Live B, Con, B-Flow	Focal Zone Up/Down
	Frequency Up/Down
CF, TVI, B-Flow Color, PDI	Box Steer
	Scale (PRF) Up/Down Auto Doppler Assistant
PWD, CW, TVD	Baseline Up/Down
	Scale (PRF) Up/Down Auto Doppler Assistant
M/D cursor	Sample volume size

# Intuitive interface

- Touch control for easy imaging parameter adjustment, even with gloves on
- 10.4" touch screen with simplified user interface
- User-friendly keys and backlight design
- Joystick



# System adjustability

Up and down



## Compact system – Easy to move



External battery for  
off-line scanning,  
up to one hour






# Easy access to PACS

## Highlights


- Enables the console to receive exams from PACs in the background without interrupting workflow
- Easily pull images via Query Retrieve to console
- Operates both wired and wirelessly





# Fetal assessment tools


SonoNT and SonoIT semi-automated measurement tools support exam consistency and reproducibility, while SonoRender/live automates render line placement in 3D/4D imaging.



# AI-based workflow and decision support tools

LOGIQ™ P10 XDclear™ harnesses artificial intelligence for imaging standardization, speed and decision support:

- Auto Lesion Segmentation
- OB Measurement Assistant
- Auto Doppler Assistant
- Breast Assistant, powered by Koios DS™\*



\*Not all products or features are available in all geographies. Check with your local GE HealthCare representative for availability in your country.  
LOGIQ and XDclear are trademarks of GE HealthCare.  
Koios DS is a trademark of Koios Medical.

# AI-based workflow tools

LOGIQ™ P10 XDclear™ harnesses artificial intelligence for imaging standardization, speed and decision support:

- Auto Lesion Segmentation
- OB Measurement Assistant
- Auto Doppler Assistant






# LOGIQ™ P10: Harnessing the power of AI

## AI-based Auto Lesion Segmentation




User identifies a breast, thyroid or liver lesion




**Benefit:** Users scan as they normally would

User clicks on the lesion and simply expands a graphical circle to encompass it



**Benefit:** User identifies a lesion with a single click

The algorithm segments the lesion, providing a trace and extents of the lesion



**Benefit:** Calipers are automatically placed to measure the lesion, saving keystrokes and providing consistency

Manual view classification

Semi-automated structure detection


# LOGIQ™ P10: Enhancing user productivity and experience over time | The journey to AI: Auto OB Assistant




User identifies a view appropriate for measuring the BPD, HC, AC, FL and HL and initiates a measurement



Auto OB Assistant automatically segments out the appropriate structure from the image



The measurement result and associated fetal age is presented to the user



**Benefit:** Users scan and initiates a measurement as they normally would

**Benefit:** Measurement is automatically performed, saving keystrokes and providing consistency

**Benefit:** Results presented to the user without extra steps

Manual view classification


Automated structure detection and results

# Harnessing the power of AI

AI-based Auto Doppler Assistant on the LOGIQ™ P10




User guides system to the vasculature of interest and pushes a button




**Benefit:** Users scan as they normally would and control initiation of the algorithm

Doppler Assistant determines the location and direction of vessels



**Benefit:** Keystrokes are reduced as the color ROI and Doppler gate are automatically placed

Doppler Assistant determines the location and direction of vessels



**Benefit:** Results are appropriate for venous versus atrial exams even when both types of vessels are present

Manual view classification

Automated structure detection and classification






# Breast Assistant, powered by Koios DS™★



Automatically provides an AI-based quantitative risk assessment that aligns to a BI-RADS® category

## In Breast Assistant:

- Measure the lesion (either manually or with Auto Contour)
- Select the “Koios” button to analyze the image
- In 2 seconds or less, a quantitative color-coded confidence scale is generated that includes the Likelihood of Malignancy (LoM)



Benign breast mass, L3-12-RS



Designed to increase clinical confidence and reduce patient anxiety

# Breast Assistant, powered by Koios DS™\*



Automatically provides an AI-based quantitative risk assessment that aligns to a BI-RADS® category:

- Based on machine learning
- Uses a proprietary AI algorithm that includes >400,000 clinical breast lesion images
- One button click to analyze
- A color-coded confidence scale generates a likelihood of malignancy (LoM) aligned to the appropriate BI-RADS category
- Results on scanner
- Select the “Koios” button to analyze the image
- Available in two seconds, or less
- Based on machine learning



\*Not all products or features are available in all geographies. Check with your local GE HealthCare representative for availability in your country.  
Koios DS is a trademark of Koios Medical. BI-RADS is a trademark of the American College of Radiology.  
© 2023 GE HealthCare. GE is a trademark of General Electric Company used under trademark license.

# Breast Assistant, powered by Koios DS™★

## Analysis return examples

B = Benign

P = Probable Benign

S = Suspicious

M = Malignant

**B** Risk Alignment: BI-RADS 2

B

P

S

M

BI-RADS DESCRIPTORS

shape **Oval**

orientation **Parallel**

Accept Reject

**P** Risk Alignment: BI-RADS 3

B

P

S

M

BI-RADS DESCRIPTORS

shape **Round**

orientation **Not parallel**

Accept Reject

**S** Risk Alignment: BI-RADS 4A-4B

B

P

S

M

BI-RADS DESCRIPTORS

shape **Irregular**

orientation **Parallel**

Accept Reject

**M** Risk Alignment: BI-RADS 4C+

B

P

S

M

BI-RADS DESCRIPTORS

shape **Irregular**

orientation **Not parallel**

Accept Reject


\*Not all products or features are available in all geographies. Check with your local GE HealthCare representative for availability in your country.  
Koios DS is a trademark of Koios Medical.

© 2023 GE HealthCare. GE is a trademark of General Electric Company used under trademark license.



# Continuous Tissue Optimization (CTO)

Delivers continuously optimized images with incredible uniformity from near to far field



# Patient-centric diagnosis



# Traditional ultrasound

Rigid assumptions about how sounds interact with the body

- Speed of sound is constant (1,540 m/s)
- Sound attenuation is constant
- Ultrasound beam is pencil shaped
- Vessels are straight
- Blood flow is laminar



Looking at the body as a phantom



## Reality is not so simple:

- Sound attenuation and velocity vary
- Ultrasound beam is a volume
- Vessels are tortuous
- Blood flow can be turbulent

## Simplifying assumptions can lead to:


- Poor image quality
- Distortion of image geometry
- Lots of adjustments to optimize images






# Agile Acoustic Architecture


- Flexible clinically-based mathematical models of the body
- Dynamically optimizes image acquisition for every body type
- Including the reality of body types
- Speed of sound is variable: 1450-1560 m/s
- Sound attenuation varies based on tissue type
- Ultrasound beam is dynamic



# P-Agile Acoustic Architecture



Agile System



P-Agile System




# Agile Acoustic Architecture

Helping improve your diagnostic confidence

## Clinical benefits

- Image uniformity
- Spatial resolution
- High frequency imaging at depth
- Reduced keystrokes to help optimize the image




# LOGIQ™ P Series probes

Innovative technologies

## ML6-15-RS


High frequency linear probe

- Matrix arrays provide multiple rows of crystals, allow focusing in the near, mid and far field
- Great spatial resolution and image uniformity from near to far field
- Footprint: 50 mm
- Ultra-high Doppler & color frequency for excellent slow flow sensitivity





# LOGIQ™ P Series probes

## Acoustic Amplifier



# XDclear™

Extended to more probes





# XDclear™ probes


## Cool stack

Penetration,  
sensitivity,  
bandwidth

## Single crystal

Axial resolution,  
bandwidth


## Acoustic Amplifier



**Single crystal:** Advanced piezoelectric material that delivers high quality acoustic signal

**Acoustic Amplifier:** Innovative design that captures and redirects the unused energy that passes through the crystal

**Cool stack:** Patented technology integrated into the transducer's internal architecture that relieves inherent heat generation that can otherwise reduce sensitivity and penetration





# LOGIQ™ P10 XDclear™ – Probe line up

Addressing a wide range of patients

## Abdominal

C1-6-D, C2-7-D, C1-5-RS, C3-10-D, M5Sc-RS, 3Sc-RS, RAB2-6-RS, 9L-RS, L3-12-RS, L3-9i-RS

## Small Parts

ML6-15-RS, L4-12t-RS, 12L-RS, L8-18i-RS, L3-12-RS, 9L-RS, L6-12-RS, L10-22-RS, C1-6-D, C1-5-RS

## Vascular

C1-6-D, C1-5-RS, C3-10-D, C2-7-D, 10C-D, ML6-15-RS, 9L-RS, L4-12t-RS, 12L-RS, L3-12-RS, L8-18i-RS, L6-12-RS, P2D, P6D, P8D

## OB/GYN

RAB2-6-RS, C1-6-D, C1-5-RS, IC9-RS, E8C-RS, E8CS-RS, RIC5-9A-RS, 9L-RS, L3-12-RS

## Cardiac

M5Sc-RS, 3Sc-RS, 6S-RS, 12S-RS, 6Tc-RS, P2D

## Breast

ML6-15-RS, L4-12t-RS, 12L-RS, L3-12RS, 9L-RS, L6-12-RS, L8-18i-RS



## Pediatrics/Neonatal

C3-10-D, 10C-D, C2-7-D, 8C-RS, 6S-RS, 12S-RS, ML6-15-RS, L8-18i-RS, L4-12t-RS, 12L-RS, L3-12RS, 9L-RS, C1-6-D, C1-5-RS

## Urology

C1-6-D, C1-5-RS, E8CS-RS, RIC5-9A-RS, BE9CS-RS, IC9-RS, L3-12-RS, ML6-15-RS, 12L-RS

# CrossXBeam™

## Spatial compounding imaging

- Provides 3, 5, 7 of spatial compounding
- Live side-by-side dual view display
- Compatible with:
  - Color Mode
  - PW
  - SRI-HD
  - Coded Harmonic Imaging
  - Virtual Convex




Help increase clinical confidence in all imaging modes

# B-Steer+

B-Steer+<sup>1</sup> enables enhanced visualization of the needles structure during interventional procedures, helping improve user confidence and exam accuracy.


## Highlights:

- Up to 8 selectable steering angles available (4 each direction)
- Separate gain control for needle reflection
- Available on all linear and convex probes
- Quick one-button operation
- Dynamic FOV
- Needle thickness




# Virtual Convex

- Provides a convex field of view
- 20% increase in image size
- Compatible with CrossXBeam™
- Available on linear and phased array probes



No VC thyroid



VC thyroid

View large anatomy in greater detail

# Contrast Imaging<sup>1,2</sup>

## Amplitude Modulation Technique

- Great penetration and contrast sensitivity
- Excellent tissue suppression
- High image uniformity

## Hi-Res Technique

- Superb spatial and temporal resolution

## Imaging modes


- Dual or single display
- Hybrid contrast
- Accumulation
- 'Easy 3D' CEUS<sup>1</sup> imaging

## Features


- Dual caliper, dual timer
- TIC Q-analysis package
- Retrospective/prospective storage
- One button background storage

## Supported probes

- BE9CS-RS, C1-5-RS, C1-6-D, C2-7-D, 9L-RS, 3Sc-RS, IC9-RS, M5Sc-RS



C1-5 CEUS



1. The LOGIQ™ P10 XDclear™ has been designed for compatibility with most commercially available ultrasound contrast agents. Availability of these agents is subject to government regulation and approval. Contrast imaging should be performed within the approved indications for use of the contrast agent used in the exam.

2. Refer to the LOGIQ P10 XDclear Product Data Sheet for a list of purchasable options. LOGIQ and XDclear are trademarks of GE HealthCare.


© 2023 GE HealthCare. GE is a trademark of General Electric Company used under trademark license.



# Contrast Imaging<sup>1,2</sup>

## Time Intensity Curve (TIC)

- Raw data processing for contrast uptake
- Q-Analysis of both compressed and uncompressed data
- Up to 8 selectable ROI's
- Up to 10 parameters
- Ellipsoid or manual ROI tracing
- Anchor tracking function
- Automated motion tracking
- Automatic enabling/disabling of frames
- Trace export in ASCII format



TIC/Q-Analysis

1. The LOGIQ™ P10 XDclear™ has been designed for compatibility with most commercially available ultrasound contrast agents. Availability of these agents is subject to government regulation and approval. Contrast imaging should be performed within the approved indications for use of the contrast agent used in the exam.

2. Refer to the LOGIQ P10 XDclear Product Data Sheet for a list of purchasable options. LOGIQ and XDclear are trademarks of GE HealthCare.

© 2023 GE HealthCare. GE is a trademark of General Electric Company used under trademark license.


# B-Flow/B-Flow Color

## A GE HealthCare innovation<sup>1</sup>

- Does not use Doppler processing
- Based on GE HealthCare patented coded technology
- Display real hemodynamics
- Direct visualization of blood reflectors

## True hemodynamics and anatomy

- Dynamic appearance of flow
- Minimal tissue overwrite compared to Color Flow
- Excellent control of flash artifacts compared to Color Flow
- Clearly visible background image
- Independent mode control



1. Refer to the LOGIQ™ P10 XDclear™ Product Data Sheet for a list of purchasable options. LOGIQ and XDclear are trademarks of GE HealthCare.


# B-Flow<sup>1</sup>/B-Flow Color

## Benefits over Doppler

- No tissue overwrite
- No impact on frame rate
- Less angle dependency
- High resolution Imaging
- Background image clearly visible

## Available on probes

- |             |             |
|-------------|-------------|
| • 9L-RS     | • L6-12-RS  |
| • 12L-RS    | • L4-12t-RS |
| • ML6-15-RS | • L10-22-RS |
| • L8-18i-RS | • L3-9i-RS  |
| • C1-5-RS   | • L3-12-RS  |
| • C1-6-D    | • E8CS-RS   |
| • C2-7-D    | • BE9CS-RS  |
| • C3-10-D   | • 10C-D     |
| • 8C-RS     | • IC9-RS    |



Neonatal head B-Flow, 9L-RS



B-Flow Color ICA, L6-12-RS


# Hybrid B-Flow

## Hybrid B-Flow<sup>1</sup>

- Visualization Modes – Toggle between three Modes: B-Flow; Dual B-Flow (simultaneous, side-by-side displays of B-Mode reference image and B-Flow); and Hybrid B-Flow (reference B-Mode image overlaid on the B-Flow for improved detail and less noise)
- Available with Easy 3D, PW, Accumulation Mode, and ON/OFF tissue background information



B-Flow Hybrid Liver, C1-5-RS



B-Flow Hybrid Jugular Vein, L3-12-RS

# HD Color

## HD Color

New functionality within B-Flow Color – offers sensitivity in visualizing small vessels and slow flow to enhance diagnostic confidence in:


- Venous, liver, and renal imaging
- Musculoskeletal imaging
- Pediatric imaging
- Small parts imaging

## Probe availability:

- C1-5-RS
- C1-6-D
- ML6-15-RS
- 12L-RS
- L3-12-RS
- L4-12t-RS



B-Flow HD Color obstetric, C1-5-RS



B-Flow HD Color thyroid, L3-12-RS

# Elastography – Clinical applications

As a non-invasive technique to help clinicians characterize abnormal tissue, ultrasound elastography offers exceptional versatility:

## Chronic liver disease

- Evaluation of liver fibrosis
- Monitor response to therapy
- Assist in patient management decisions

## Oncology

- Evaluation of soft tissue lesions, including breast, small parts and urology
- Soft tissue lesion classification
- Surgical planning

## Musculoskeletal

- Provide additional information for diagnosis of tendinopathy, tendinosis, synovial hypertrophy, tears, and other conditions
- Assist in injury management from early diagnosis (e.g., acute muscle injury) to treatment response

## Gynecology

- Additional information to assess cervical strength to help identify women at risk of preterm birth

## Urology


- Assists in prostate cancer identification by improving localization of abnormal foci thus enabling more targeted biopsies

## Breast

- Useful in the detection and characterization of breast disease
- Well-suited for dense breast patients
- Provides a quantitative measurement of tissue elasticity

## Thyroid

- Useful in nodules with indeterminate US or cytologic characteristics
- Well-suited to follow up on lesions







# 2D Shear Wave Elastography

Provides a quantitative estimate of tissue elasticity and displays this measurement in color-coded elastograms. Ultrasound imaging and 2D Shear Wave Elastography can be combined in a single exam.

- Assists in a growing range of applications, including chronic liver disease, and diseases of the breast, musculoskeletal (MSK) and small parts
- Semi-automatic measurements provide quick, quantifiable information
- Multiple measurement regions of interest (ROIs) help increase exam speed
- Probe availability:
  - C1-6-D
  - C1-5-RS
  - L3-12-RS
  - ML6-15-RS
  - 12L-RS
  - IC9-RS



2D Shear Wave Elastography with Quality Indicator, liver, C1-6-D



2D Shear Wave Elastography with Quality Indicator, cervix, IC9-RS

# Strain Elastography<sup>1</sup>

Strain imaging technology requiring a light manual compression or patient breath to perform tissue deformation. A qualitative and Semi-Quantitative<sup>2</sup> solution.

## Highlights:


- High sensitivity and persistence
- Consistent pattern
- User selectable color maps
- Dual measurements
- User support by pressure quality bar and graph

## Shear Wave elastography


- Focused burst of acoustic energy to perform tissue deformation through a comb-push excitation
- Color coded elastogram and quantitative measurements
- Focus: Chronic liver disease, oncology

## Available on the following probes:

- |             |             |            |
|-------------|-------------|------------|
| • C1-5-RS   | • 12L-RS    | • E8CS-RS  |
| • C1-6-D    |             |            |
| • ML6-15-RS | • L6-12-RS  | • BE9CS-RS |
| • 9L-RS     | • L4-12t-RS | • IC9-RS   |
|             | • L3-12-RS  |            |




This soft lesion deforms under pressure.



This hard lesion does not deform under pressure.

$$\text{Elasticity} = \frac{\text{Stress}}{\text{Strain}}$$

**Stress** = axial force applied to lesion  
**Strain** = tissue deformation due to applied stress



1. Refer to the LOGIQ™ P10 XDclear™ Product Data Sheet for a list of purchasable options.

2. Not cleared by the U.S. FDA and not for sale in the USA.

LOGIQ and XDclear are trademarks of GE HealthCare.

© 2023 GE HealthCare. GE is a trademark of General Electric Company used under trademark license.

# Strain Elastography<sup>1</sup>

## Semi-Quantification<sup>2</sup>

### E-Index



- Selected ROI's elasticity value
- Round or manual tracing of ROI
- E-index range from 0 – 6
- Based on GE HealthCare Raw Data processing

### E-Ratio

- Calculated E-Index between ROI's (up to 8)
- Representing relative stiffness

### Q-Analysis over multi-frame acquisition

- Automatic skip of low-quality frames
- Anchor function
- Trace export in ASCII format



1. Refer to the LOGIQ™ P10 XDclear™ Product Data Sheet for a list of purchasable options.

2. Not cleared by the U.S. FDA and not for sale in the USA.


LOGIQ and XDclear are trademarks of GE HealthCare.

© 2023 GE HealthCare. GE is a trademark of General Electric Company used under trademark license.

# LOGIQ™ View

**LOGIQ View<sup>1</sup>** allows a real-time representation of long anatomical areas (e.g., Patellar Tendon)

- Available on all probes
- Combine with CrossXBeam™ for linear probes and convex probes
- Auto detection of scan direction
- Up to 60 cm scan length





# 3D/4D ultrasound<sup>1</sup>

GE HealthCare Volume Imaging with 3D/4D dedicated probe

## Highlights:

- Easy, quick, reproducible
- High volume data accuracy
- Tint Render Map
- Comprehensive settings

## Volume Modes:

- Multi-planar Imaging
- Surface rendering
- TUI – Tomographic Ultrasound Imaging
- VCI – Volume Contrast Imaging
- Vocal– Volume Calculation
- STIC
- OmniView
- SonoRenderlive
- HDlive™

## Probe:

- RAB2-6-RS
- RIC5-9A-RS



1. Refer to the LOGIQ™ P10 XDclear™ Product Data Sheet for a list of purchasable options. HDlive, LOGIQ, and XDclear are trademarks of GE HealthCare.


# 3D/4D ultrasound<sup>1</sup>

## TUI – Tomographic Ultrasound Imaging

Visualization mode that presents data as parallel slices (planes) through acquired Volume dataset. It works with CFM/PDI and SRI-HD.

### Highlights:

- Information consistent to CT & MR format
- 3D static with color
- Up to 9 slices, with user selectable distance (min 0.5 mm, step by 0.1 mm) and angle
- Top left held as reference image
- Works with SRI



Pelvic TUI



# 3D/4D ultrasound<sup>1</sup>

## VCI – Volume Contrast Imaging

VCI is a volume acquisition technique enhancing B-Mode contrast resolution and speckle suppression.

### Highlights:

- Help improve assessment of lesions size, margins and internal structures for comprehensive patient management
- 3D manipulation in A, B, and C planes



VCI uterus, RIC5-9A-RS


# 3D/4D ultrasound<sup>1</sup>

## Vocal

Vocal is a 'Volume Computer-aided analysis' based on a volume acquisition enabling fast and accurate volume calculations.

### Highlights:

- Manual, semi-automatic or automatic borders definition
- Basic measurements as length, angle and area
- Easy corrections and contour modifications
- Can be used with any lesion or volume to measure




# 3D/4D ultrasound<sup>1</sup>

## STIC

Spatio-Temporal Image Correlation (STIC) captures one fetal heart cycle in 3D cine.

### Highlights:

- Adjustable acquisition time
- Use with Color Doppler or Power Doppler modes
- 3D manipulation in A, B, and C planes




# 3D/4D ultrasound<sup>1</sup>

## OmniView

OmniView is the “any-plane” function for 3D and 4D data

### Highlights:

- Gives the ability to trace along any shape or structure
- Can be started from A, B, or C. Up to 3 planes can be displayed simultaneously




VCI OmniView, RIC5-9A-RS

# 3D/4D ultrasound<sup>1</sup>

## SonoRender<sup>live</sup>

- Automates render line placement in 3D/4D imaging
- Helps remove artifacts, improves rendered image without having to use the scalpel tool



SonoRender<sup>live</sup> off



SonoRender<sup>live</sup> on

# 3D ultrasound<sup>1</sup>

HDlive™

- Reveal a clinical perspective of fetal and female anatomy that brings anatomical realism to surface structures
- Utilizes a combination of advanced skin illuminating and shadowing techniques to expose possible hidden details for a deeper understanding of relational anatomy and developing structures



1. Refer to the LOGIQ™ P10 XDclear™ Product Data Sheet for a list of purchasable options. HDlive, LOGIQ, and XDclear are trademarks of GE HealthCare.

© 2023 GE HealthCare. GE is a trademark of General Electric Company used under trademark license.





# Flow Quantification<sup>1</sup>


2D CFM/PDI Quantitative assessment of vascular feeding in a selected ROI

## Highlights:

- Consistent, repeatable and objective measurement
- Can help in treatment planning and monitoring protocols
- Provides data to support outcome measurements

## Features:


- Up to 8 selectable ROI's
- Analysis over 4/5 heartbeat cycles
- Automatic or manual ROI tracing
- “Save ROI” feature for monitoring
- Manual disabling & enabling of frames
- Export traces in ASCII format



QA analysis

# AMM – Anatomical M-Mode<sup>1</sup>

- M-Mode cursor adjustable at any plane to be vertical to myocardium
- Compatible with:
  - Live image
  - Stored image – Raw Data
  - Color Flow Mode



# Connectivity



## Output ports

- USB
- Composite
- HDMI
- S-Video
- Ethernet

## Network storage

- DICOM®
- SaveAs
- MPEGView
- Report Save As




# Auto Optimization

## One button press




## ATO (Auto Tissue Optimization)




Optimize B-Mode image to help improve contrast resolution

## Auto TGC



Continuous 2D auto TGC control

## ASO (Auto Spectral Optimization)



Baseline and scale (PRF)

# Barcode reader

Help streamline exams

Barcode reader helps to:

Simplify patient entry



0001LastNameFirstNameMiddleName19901010F



MeasureReportScan Assistant ManagerImaging Preset ManagerLDAPDisk EncryptAudit ReportBarcode

Barcode

Input Mode

Patient

○ Patient ID

○ Complexatio

○ Off

Dicom Worklist

○ Patient ID

○ Accession #

○ Off

Barcode Scanner

Xenon 1900 Area-Imaging Scanner

Input Data

0001LastNameFirstNameMiddleName19901010F

Complexation

Patient ID

1

2

00

Other ID

3

4

01

Last Name

5

12

LastName

First Name

13

21

FirstName

Middle Name

22

31

MiddleName

Birth Year

32

35

1990

Birth Month

36

37

10

Birth Day

38

39

10

Gender

40

40

F

Male

M

Female

F

Save

Exit

Cancel

# Scan Assistant<sup>1</sup>

- The personal assistant to user's exam
- Up to 63% time reduction<sup>2</sup>
- Up to 87% keystroke reduction<sup>2</sup>
- Help increase exam consistency
- Compatible with LOGIQ™ E9, LOGIQ S8 and LOGIQ S7 workflows

## Feature highlights


Initiates and completes user selected required measurements

Automatically reorders images to reader's preference

Automatically steers Color Doppler

Automatically inserts comments

Automatically sets up Imaging controls and modes



1. Refer to the LOGIQ P10 XDclear™ Product Data Sheet for a list of purchasable options.  
2. Internal GE HealthCare engineering study using standardized protocols for an abdominal exam compared with prior version GE HealthCare LOGIQ P6 ultrasound system.  
LOGIQ and XDclear are trademarks of GE HealthCare.



# Compare Assistant<sup>1</sup>

Help streamline comparison to prior exams

## Opportunity

Many ultrasound exams are follow ups to prior ultrasound or other modality exams

## Goal

Drive productivity for acquiring and reading the exam by designing a workflow that uses prior exam data

## Result

A quick image comparison or a replicated prior exam

1. Refer to the LOGIQ™ P10 XDclear™ Product Data Sheet for a list of purchasable options. LOGIQ and XDclear are trademarks of GE HealthCare.



# Compare Assistant<sup>1</sup>

Help streamline comparison to prior exams



## At the scanner...

Easy access to past exam data  
on the scanner



Side-by-side compare and store past exam  
image to today's image



Set image setting of live scanning to match  
past exam image<sup>2</sup>



Create entire new exam  
to match old exam

## At the reading station...

Help reduce the time spent to find, open, sort, compare  
to prior exams



1. Refer to the LOGIQ P10 XDclear™ Product Data Sheet for a list of purchasable options.

2. Automated when LOGIQ P10 XDclear data is used.

LOGIQ and XDclear are trademarks of GE HealthCare.


© 2023 GE HealthCare. GE is a trademark of General Electric Company used under trademark license.

# Breast Productivity<sup>1</sup>

## Measurement package

A dedicated breast-specific measurement package that allows user to:

- Make labeling, measuring and describing lesion easy
- Leverage the BI-RADS® lexicon criteria/assessment
- Organizes multiple measurements into a convenient worksheet
- Send results via DICOM® SR



Directly from BI-RADS lexicon

Rt	
Lesion	
Shape	Irregular
Orientation	Parallel
Margin*	Circumscribed
Echo Pattern	Complex cystic and solid
Posterior Features	Enhancement
Associated features*	Absent
Calcifications*	Calcifications in a mass
Special Cases*	Lymph nodes-axillary
BI-RADS Assessment	
<input checked="" type="radio"/> 0 <input type="radio"/> 1 <input type="radio"/> 2 <input type="radio"/> 3 <input type="radio"/> 4a <input type="radio"/> 4b <input type="radio"/> 4c <input type="radio"/> 5 <input type="radio"/> 6	
Koios Assessment	

# Thyroid Productivity<sup>1</sup>

## Measurement package

### Thyroid-specific measurement package



- Enables labeling, measuring and describing nodules, lymph nodes and parathyroids
- Multiple measurements can be organized into a convenient worksheet and sends results via DICOM® SR
- TI-RADS® compatible

Overall  
Thyroid

---

Resected -  
Appearance -

Comment



### Show features


1. Refer to the LOGIQ™ P10 XDclear™ Product Data Sheet for a list of purchasable options. LOGIQ and XDclear are trademarks of GE HealthCare. DICOM is a registered trademark of the National Electrical Manufacturers Association. TI-RADS is a registered trademark of the American College of Radiology.

© 2023 GE HealthCare. GE is a trademark of General Electric Company used under trademark license.

# Auto IMT<sup>1</sup>

**Auto IMT** is an automated method of measuring the intima media thickness of the CCA or ICA from multiple samples across a user defined length.

- Simple and easy to operate
- Direct export of measurements to a worksheet and report page
- Including ECG trigger to help increase consistency and reliability
- Save offset distance and IMT measurement lengths to help increase reproducibility





An efficient, reproducible method of carotid artery analysis

# Auto EF<sup>1</sup>

**Automated Ejection Fraction (AutoEF)** is a semi-automatic measurement tool used for measurement of the global EF (Ejection Fraction).

- The AutoEF tool tracks and calculates the myocardial tissue deformation based on feature tracking on
- B-Mode cine loops
- AutoEF is performed on either one or both apical
- 4-chamber or 2-chamber views, in any order
- Result is presented as Ejection Fraction value for each view and average Ejection Fraction for the whole LV. All values are stored to the worksheet after the results are approved






# Cardiac Strain<sup>1</sup>


**Cardiac Strain** is the support tool for evaluation of the wall motion of the entire left ventricle or local region.

- Reduce exam time with easy workflow
- Wall motion tracing in 2D view from Apex (APLAX, 4-Ch, 2-Ch)
- Display strain graph, Parametric imaging and Anatomical M mode for each plane
- Bull's eye is provided from 3 plane
- Display Peak Systolic Strain and Post Systolic Index




# Ultrasound-Guided Attenuation Parameter (UGAP)<sup>1</sup>

- A real-time, image-guided method to measure fatty content in the liver tissue
- A valuable metric in assessing non-alcoholic fatty liver disease (NAFLD)




# Speckle Reduction Imaging (SRI)

- New algorithms to enhance border differentiation, contrast resolution, improve tendon visualization, removes unwanted speckle and improves tissue homogeneity
- Available on all LOGIQ™ P10 transducers and models



# SonoAVC

Automatic Volume Calculation – helps to identify and measure follicles and cystic areas within a 3D volume. It standardizes the process of follicular assessment and decreases inter-observer and intra-observer variability.




# Urology enhancement

## IC9-RS features

- 2D Shear Wave Elastography
- Strain Elastography
- CEUS

## BE9CS-RS features

- Strain Elastography
- CEUS





# Multi-purpose capabilities for a diverse workload

Designed with your patients in mind



Liver



Cardiac



OB/GYN



Breast



Thyroid



Musculoskeletal



Urology



Vascular



Pediatrics



# Power Assistant<sup>1</sup>

**Power Assistant** is an innovative solution that provides the system battery power during transport to help decrease system shut-down and reboot time – helping achieve excellent productivity for excellent portable exams.

## Highlights

- Prompt in & out of battery operation mode to help improve system's portability
- Simple plug in and out operation
- System safely shuts down automatically before battery runs out
- Wireless LAN capability
- Always ready to scan with Power Assist battery operation
- External battery enables up to one hour of offline scanning




1. Refer to the LOGIQ™ P10 XDclear™ Product Data Sheet for a list of purchasable options. LOGIQ and XDclear are trademarks of GE HealthCare.

© 2023 GE HealthCare. GE is a trademark of General Electric Company used under trademark license.

# Raw data

The foundation for simplified workflow

**Raw data** capture enables to build a thorough exam while helping reduce scan time. This proprietary raw data format from GE HealthCare captures data earlier in the image processing chain enabling users to make changes to the data during or even after the exam has ended.



# Powerful support



# LOGIQ™ apps

## Remote Control

- Remotely operate the system from tablet or phone that has LOGIQ Smart App loaded
- Focused on ergonomics
- Includes:
  - Major modes
  - Freeze/print
  - Depth
  - Gain
  - ROI placement
  - Dual Image





## LOGIQ™ apps

### Photo Assistant

#### A picture is worth a 1,000 words

- Photograph relevant anatomy and include photos with the clinical images
- Provides value context for documentation and comparison after a procedure
- Utilizes Android™ tablet or phone



# LOGIQ™ P10 XDclear™ and Trice® – Sharing ultrasound images




## Clinical communication

- Ultrasound lab sends images, reports for consultation before referring the patient
- Ultrasound lab shares examination results for second opinion
- Ultrasound lab shares examination results with surgeon for treatment and surgery planning
- Sharing cases with medical community, accreditation support
- Routing of studies from remote sites

## Doctor-patient communication

- Sharing exam results with patients






# Data privacy & cybersecurity



## Purpose

- Keep the machine functional in the face of cyber threats
- Protect the patient data on the machine from unauthorized access
- Enable you to successfully implement their HIPPA and security policies and still manage your daily workflow

## SonoDefense solution




# Defense-in-depth strategy

SonoDefense is designed for maximum security protection with a defense-in-depth strategy that incorporates security controls deployed in multiple layers. This approach enhances security by protecting the system against any particular attack using several independent methods.

Limits what can be run on the LOGIQ™ P10 XDclear™


- Customizable, role-based access
- Federated Identity Management
- Session management
- Auditing
- Secure remote access

- Customizable patient data encryption
- Enterprise wireless encryption
- IPv6 Internet Protocol address standard



# Raw data

TruScan™ architecture | Capturing raw data early in the image chain





Raw data processing with TruScan architecture enabling “virtual rescanning”

# Raw Data Processing

Original Acoustic Data are stored before scan converting in a GE HealthCare “raw” format to be easily accessed and re-processed any time after the exam completion.

## Highlights:


- Sub-optimal studies can be optimized
- Measurements can be re-done and reports regenerated
- All Imaging control parameters can be changed as:
  - B-Mode: Gain, DR, AO, Zoom, SRI...
  - CFM: Gain, Threshold, DualView...
  - PW: Baseline, Invert, Angle, Gain...




# Multi-modality Query Retrieve<sup>1</sup>

## Volume Navigation Import


- View real-time ultrasound with retrieved Ultrasound, Mammo, CT or MR datasets
- Potential time saver
- Simplified workflow
- Complementarity of information
  - Volume data set retrieval
  - Select desired image plane by scrolling
- Helps improve diagnostic confidence



# My Trainer on-board training modules



Expandable  
menus



Content  
description

Digitally encoded technology; using digital codes to enhance weak signals from small particulate reflectors and suppress signals from strong reflectors (tissue). Coded excitation and tissue equalization are key to B-Flow.

## Activation

Created by the B-Mode processor.

B-Flow is located on the keyboard next to the B-Mode button.

B-Flow gain is controlled by turning the B-Mode gain button.

**Smaller:** Reduces the size of the My Trainer display  
**Edit:** Allows you to remove content that is not required  
**Exit (Esc):** Close My Trainer



# LOGIQ™ Club users' community

[www.LOGIQClub.net](http://www.LOGIQClub.net)

## Access to the Club websites

- Local contents, events
- Clinical courses in your country
- Application tips & tracks
- Publications and cases & papers
- Clinical cases & technical presentations
- Downloadable educational materials and DVD
- Online education

## Personalized mailings & newsletters

- Be first to learn about new ultrasound products and software upgrades

## User days and VIP lounges

- Learn about best practices from specialists around the globe
- Discuss and exchange information with ultrasound users worldwide





# LOGIQ™ Club users' community

## Education highlights

### Online education

- GE HealthCare white papers
- Application videos
- Clinical tutorials
- Video tutorials on advanced features
- Access to peer-to-peer reviewed papers
- International Academy of Medical Education (IAMU) at <http://iamu.logiqclub.net>

### Residential education

- Volume Navigation training schools in:
  - Abdominal, Interventional
  - Musculoskeletal
  - Breast
  - Urology

### Mobile app – now available

- Watch course videos and product tutorials directly on your smartphone

# Welcome to the LOGIQ Club



# Clinical skills training: Digital Expert\*

**Designed for:** Ultrasound department managers and console users

**Problem:** It is difficult to maintain the skill level of all users for all features on advanced ultrasound systems, especially for new staff or infrequent exam types

**Solution:** Digital Expert provides easy on-demand access to live expert applications training, delivered on your own device

GE HealthCare Digital Expert delivers one-on-one personalized and hands-on learning from your ultrasound console to provide ongoing education to fully utilize the functions and features of your system. See more at <https://youtu.be/TZjISiOkPFM>

## Hands-on learning, personalized experience

- Users learn on their own equipment
- Instructors see the users primary monitor allowing for them to understand what the student is seeing

## Flexible and convenient training to fit your schedule

- Users schedule classes that fit their work and personal schedules
- 30-60 minutes courses with convenient online booking and flexible cancellation ensure that learning does not disrupt patient scanning



## Digital technology connecting you to clinical experts

- Face-to-face technology provides a more personal experience vs phone calls
- Our trainers have many years of experience and are extensively trained on the content and technology

Live one-on-one expert training from your own system on your own schedule

