ABBOTT MECHANICAL HEART VALVES

Regent, Masters HP and Masters

AORTIC VALVE

MITRAL VALVE

IMPLANT
CONFIDENTLY.

Information contained herein for DISTRIBUTION in Europe, Middle East and Africa ONLY.
Check the regulatory status of the device in areas where CE marking is not the regulation in force.
©2021 Abbott. All rights reserved. 9-EH-1-12392-01 06-2021

**CONFIDENTLY
IMPLANT
THE MOST
TRUSTED
MECHANICAL
VALVES
IN THE WORLD**

**3 MILLION PATIENTS
TREATED WITH ABBOTT
MECHANICAL HEART VALVES**

**MORE THAN 1,000
PEER-REVIEWED
PUBLICATIONS PROVIDE
EVIDENCE FOR ABBOTT
MECHANICAL HEART VALVES**

**LOW THROMBOGENICITY AND
EXCELLENT PATIENT OUTCOMES**

PROVEN DESIGN TO RESTORE NATIVE VALVE HEMODYNAMICS

AN ABBOTT HALLMARK

Featured in all Abbott Mechanical Heart Valves the unique Pivot Guard Design offers benefits both during implant and post-implant.

Shields pivot mechanism from pannus ingrowth

Minimizes interaction with sub-annular native valve apparatus in the mitral position and ensure coronary ostia clearance in the aortic position

Enables for an 85° leaflet opening angle, minimizing leaflet flutter and leading to smoother laminar flow through the orifice*


Can lessen thrombus formation by minimizing carbon surface area and thanks to the washout flow through the hinges*

Information contained herein for DISTRIBUTION in Europe, Middle East and Africa ONLY.
Check the regulatory status of the device in areas where CE marking is not the regulation in force.
©2021 Abbott. All rights reserved. 9-EH-1-12392-01 06-2021

**FIND OUT MORE
ON ABBOTT
MHV DESIGN**

2017 ESC/EACTS GUIDELINES¹

Target INR for mechanical prostheses

Prosthesis thrombogenicity	Patient-related factors ^a	
	None	≥1 risk factor
Low ^b	2.5	3.0
Medium ^c	3.0	3.5
High ^d	3.5	4.0

INR = international normalized ratio; LVEF = left ventricular ejection fraction.

^a Mitral or tricuspid valve replacement; previous thromboembolism; atrial fibrillation; mitral stenosis of any degree; LVEF <35%.

^b Carbomedics, Medtronic Hall, ATS, Medtronic Open-Pivot, St Jude Medical, OnX, Sorin Bicarbon.

^c Other bileaflet valves with insufficient data.

^d Lillehei-Kaster, Omniscience, Starr-Edwards (ball-cage), Björk-Shiley and other tilting-disc valves.

OPERATE WITH THE FACTS

ABBOTT MECHANICAL HEART VALVES SHOW LOWER THROMBOEMBOLISM, THROMBOSIS AND BLEEDING EVEN AT A LOW INR RANGE

INR 1.5 ————— 2.0 ————— 2.5

LOWERING-IT² (target INR 1.5–2.5)

Randomized Study	197 patients implanted* (44 with Abbott valves, 153 with LivaNova valves), 5 years
Thromboembolism	0.09%/pt-year
Thrombosis	0%/pt-year
Bleeding Events	0.56%/pt-year

ESCAT III³ (target INR 1.6–2.1)

Randomized Study	1137 patients** (all Abbott valves), 2 years
Thromboembolism	0%/pt-year, 0.58%/pt-year***
Bleeding Events	0.58%/pt-year, 1.07%/pt-year†

PROACT⁴ (target INR 1.5–2.0)

Randomized Study	375 patients (all On-X Valves), 3 years
Thromboembolism	2.67%/pt-year
Bleeding Events	2.67%/pt-year

*44/197 patients in the Lowering-IT study were implanted with Abbott Valves. **This was further stratified into a control group, a very low INR (monitored 1x weekly), and a very low INR (monitored 2x weekly) group. ***Thromboembolic events for VL1 and VL2 groups are listed together, respectively. †Bleeding events for VL1 and VL2 groups are listed together, respectively.

CAUTION: This product is intended for use by or under the direction of a physician. Prior to use, reference the Instructions for Use provided inside the product carton (when available), at eifu.abbottvascular.com or at medical.abbott/manuals for more detailed information on Indications, Contraindications, Warnings, Precautions and Adverse Events. Information contained herein for Europe, Middle East and Africa ONLY. Check the regulatory status of the device in areas where CE marking is not the regulation in force. Tests performed by and data on file at Abbott unless otherwise noted. Illustrations are artist's representations only and should not be considered as engineering drawings or photographs. Photos on file at Abbott.


Abbott Vascular International BVBA
Park Lane, Culliganlaan 2B, B-1831 Diegem, Belgium, Tel: +32 2 714 14 11 www.Cardiovascular.Abbott
©2019 Abbott. All rights reserved 9-EH-8-9151-01 03-2019

**FIND OUT MORE
ON LOW INR
CLINICAL EVIDENCE**

A BROAD RANGE OF SOLUTIONS TO TAILOR THE IMPLANTATION TO EVERY PATIENT

AORTIC VALVE

Regent
EXCEPTIONAL HEMODYNAMICS
IN THE AORTIC POSITION IN SMALL
AORTIC ROOT PATIENTS

Masters HP
HEMODYNAMICS MEETS
IMPLANTABILITY IN BOTH THE
AORTIC AND MITRAL POSITION,
NOW AVAILABLE ALSO FOR BABIES
AND NEWBORNS

Masters
OPTIMAL IMPLANTABILITY IN BOTH
THE AORTIC AND MITRAL POSITION

FIND OUT HOW TO MINIMIZE AV BLOCK
WHEN IMPLANTING
A MECHANICAL HEART VALVE

TREATING THE TINIEST PATIENTS INCLUDING NEWBORNS AND BABIES IS NOW POSSIBLE...

with the world's smallest
mechanical heart valves

15mm MASTERS HP
MITRAL AND AORTIC VALVE

WATCH
SADIE'S STORY

VALVE ORDERING GUIDE

STANDARD CUFF

Compact, double
velour Dacron

FLEXCUFF

Flanged and more
conformable than the standard
cuff, to accommodate variable
anatomy

EXPANDED CUFF

25% more cuff material
than the standard cuff,
for even more anatomic
accommodation

PTFE

Known for its low
friction, very easy
to suture

EXPANDED PTFE

Easy to suture with
16% more material,
for extra anatomical
conformability

AORTIC VALVE

	REGENT™ VALVE		MASTERS HP Series		MASTERS Series		
SIZE (MM)	STANDARD CUFF	FLEXCUFF™ 4	STANDARD CUFF	EXPANDED CUFF 1,3	STANDARD CUFF	EXPANDED CUFF 1,3	PTFE
15			15AHPJ-505				
17	17AGN-751	17AGFN-756	17AHPJ-505	17AEHPJ-505			
19	19AGN-751	19AGFN-756	19AHPJ-505	19AEHPJ-505	19AJ-501	19AECJ-502	19ATJ-503
21	21AGN-751	21AGFN-756	21AHPJ-505	21AEHPJ-505	21AJ-501	21AECJ-502	21ATJ-503
23	23AGN-751	23AGFN-756	23AHPJ-505	23AEHPJ-505	23AJ-501	23AECJ-502	23ATJ-503
25	25AGN-751	25AGFN-756	25AHPJ-505	25AEHPJ-505	25AJ-501	25AECJ-502	25ATJ-503
27	27AGN-751	27AGFN-756	27AHPJ-505	27AEHPJ-505	27AJ-501	27AECJ-502	27ATJ-503
29	29AGN-751	29AGFN-756			29AJ-501	29AECJ-502	29ATJ-503
31					31AJ-501	31AECJ-502	31ATJ-503

MITRAL VALVE

	MASTERS HP Series	MASTERS Series			
SIZE (MM)	STANDARD CUFF	STANDARD CUFF	EXPANDED CUFF 1,3	PTFE	EXPANDED PTFE 2
15	15MHPJ-505				
17	17MHPJ-505				
19	19MHPJ-505	19MJ-501	19MECJ-502	19MTJ-503	19METJ-504
21	21MHPJ-505	21MJ-501	21MECJ-502	21MTJ-503	21METJ-504
23	23MHPJ-505	23MJ-501	23MECJ-502	23MTJ-503	23METJ-504
25	25MHPJ-505	25MJ-501	25MECJ-502	25MTJ-503	25METJ-504
27	27MHPJ-505	27MJ-501	27MECJ-502	27MTJ-503	27METJ-504
29		29MJ-501	29MECJ-502	29MTJ-503	29METJ-504
31		31MJ-501	31MECJ-502	31MTJ-503	31METJ-504
33		33MJ-501	33MECJ-502	33MTJ-503	33METJ-504
35		35MJ-501			
37		37MJ-501			

1. The Expanded Aortic and Mitral Cuff has approximately 25% more cuff material than the standard cuff, for even more anatomic accommodation.
2. The Expanded PTFE Cuff easy to suture with 16% more material, for extra anatomical conformability
3. The Expanded HP Cuff has approximately 15% more cuff than the HP Series cuff.
4. The FlexCuff™ is flanged and more conformable than the standard cuff, to accommodate variable anatomy.

ACCESSORIES

ORDERING GUIDE

MECHANICAL VALVE SIZER SETS AND ACCESSORIES

Model No.	DESCRIPTION
905	Universal Sizer Set contains 17-33 mm valve sizers and valve holder handle model 905-HH.
905-15	Mitral and Aortic Double-Ended Sizer for Masters HP 15 mm.
905-35	35mm Masters Series valve sizer.
905-37	37mm Masters Series valve sizer.
905-MHH	Mitral valve holder handle.
905-RHH	Rigid valve holder handle.
907	Regent Sizer Set contains 17-29mm valve sizers and valve holder handle model 905-HH.
A-RHR	Contains sizes 19mm-31mm aortic Masters Series holder/rotators.
M-RHR	Contains sizes 19mm-37mm mitral Masters Series holder/rotators.
AHP-RHR	Contains sizes 17mm-27mm Masters Series Hemodynamic Plus holder/rotators.
AG-RHR	Contains sizes 17mm-29mm Regent™ holder/rotators.
LT100	Mechanical valve leaflet tester.

REFERENCES

1. European Heart Journal, Volume 38, Issue 36, 21 September 2017, Pages 2739–2791, <https://doi.org/10.1093/eurheartj/ehx391>
2. Torella, Michele, et al. "LOWERING the INTensity of oral anticoagulant Therapy in patients with bileaflet mechanical aortic valve replacement: results from the "LOWERING-IT" Trial." American heart journal 160.1 (2010): 171-178.
3. Koertke, Heinrich, et al. "Efficacy and safety of very low-dose self-management of oral anticoagulation in patients with mechanical heart valve replacement." The Annals of thoracic surgery 90.5 (2010): 1487-1493.
4. Puskas JD et al. Reduced anticoagulation after mechanical aortic valve replacement: interim results from the prospective randomized on-X valve anticoagulation clinical trial randomized Food and Drug Administration investigational device exemption trial. J Thorac Cardiovasc Surg 2014;147:1202-11

CAUTION: This product is intended for use by or under the direction of a physician. Prior to use, reference the Instructions for Use, inside the product carton (when available) or at eifu.abbottvascular.com or at medical.abbott.com for more detailed information on Indications, Contraindications, Warnings, Precautions and Adverse Events.

This Material is intended for use with healthcare professionals only. Data on File at Abbott. **Information contained herein for DISTRIBUTION in Europe, Middle East and Africa ONLY.** Check the regulatory status of the device in areas where CE marking is not the regulation in force.

Abbott International BVBA

Park Lane, Culliganlaan 2B, 1831 Diegem, Belgium, Tel: 32.2.714.14.11

™ Indicates a trademark of the Abbott Group of Companies.

‡ Indicates a third party trademark, which is property of its respective owner.

www.cardiovascular.abbott


©2021 Abbott. All rights reserved. 9-EH-1-12392-01 06-2021

Regent™ Mechanical Heart Valve

Standard Cuff

Aortic Valves

Product Highlights

- Designed to deliver exceptional hemodynamics and performance while maintaining low complication rates, structural integrity and durability
- Single-digit gradients, even in sizes as small as 19 mm¹
- Up to 84% orifice to annulus ratio
- 85 degree leaflet opening angle offers improved laminar flow and reduces turbulence²⁻⁴
- Supra-annular placement and low implant height
- The FlexCuff™ is flanged and more pliable than the standard cuff
- St. Jude Medical heart valves are MR Conditional⁵

Ordering Information

Contents: Aortic Valve (1 unit per box)

Model/Reorder* Number	Tissue Annulus Diameter (mm)	Valve Orifice Inner Diameter (mm)	Geometric Orifice Area (cm ²)	Effective Orifice Area ⁶ (cm ²)	Cuff Style
17AGN-751	17	15.9	1.87	1.42	Standard
19AGN-751	19	17.8	2.39	1.84	Standard
21AGN-751	21	19.6	2.90	2.47	Standard
23AGN-751	23	21.4	3.45	2.91	Standard
25AGN-751	25	23.0	4.02	3.34	Standard
27AGN-751	27	24.9	4.69	4.28	Standard
29AGN-751	29	26.8	5.44	Not reported	Standard

*All sizes currently not available in all markets.

1. St. Jude Medical, Inc. Pre-Market Approval Application Supplement. Regent Heart Valve. Washington, D.C.; U.S. Food and Drug Administration. 2002. P810002/S57.
2. King, M., David T, Fisher J. An initial parametric study on fluid flow through bileaflet mechanical heart valves using computational fluid dynamics. *J Eng Med.* 1994;208:63-72.
3. Feng Z, Nakamura T, Fujimoto T, et al. In vitro investigation of opening behavior and hydrodynamics of bileaflet valves in the mitral position. *Artificial Organs.* 2002;26(1):32-9.
4. Shipkowitz T, Ambrus J, Kurk J, et al. Evaluation technique for bileaflet mechanical valves. *J Heart Valve Dis.* 2002;11(2):275-82.
5. St. Jude Medical MRI Safety Document.
6. Walker DK, Brendzel AM, Scotten LN. The new St. Jude Medical Regent mechanical heart valve: Laboratory measurements of hydrodynamic performance. *J Heart Valve Dis.* 1999;8(6):687-96.

Please review the Instructions for Use prior to using these devices for a complete listing of indications, contraindications, warnings, precautions, potential adverse events and directions for use.

Product referenced is approved for CE Mark.

Regent, FlexCuff, ST. JUDE MEDICAL, the nine-squares symbol and MORE CONTROL, LESS RISK. are registered and unregistered trademarks and service marks of St. Jude Medical, Inc. and its related companies.

SJM Regent[®]

Heart Valve

REDEFINING HEMODYNAMIC PERFORMANCE.[™]

THE NEWEST MEMBER OF THE ST. JUDE MEDICAL
MECHANICAL HEART VALVES FAMILY

REACHING NEW LEVELS.


The SJM Regent® heart valve represents a significant step forward in prosthetic valve design. It provides outstanding hemodynamics while maintaining the traditional quality and proven features that have established St. Jude Medical® mechanical valves as the gold standard.

Regent™ Valve

Evolutionary improvement

Development of the St. Jude Medical® mechanical heart valve has led to a progressively greater geometric orifice area within a given tissue annulus dimension.

SJM® Masters Series mechanical heart valve

- Intraannular cuff.
- Intraannular carbon rim.


SJM® Masters Series Hemodynamic Plus valve

- Supraannular cuff.
- Carbon rim remains intraannular.

SJM Regent® heart valve

- Supraannular cuff.
- Carbon rim shifts to supraannular position.
- New rotation mechanism is completely housed within the carbon rim.

Orifice to annulus ration — 19mm valve²

"The energy loss results show an approximate one-size increment improvement of the SJM Regent™ heart valve over the St. Jude Medical® HP Series valve, which is equivalent to a two-size increment improvement over the standard valve."

(Walker et al, 1999)¹





SJM REGENT® HEART VALVE

Redefining Hemodynamic Performance.™

Unprecedented hemodynamics

The SJM Regent® heart valve delivers exceptional hemodynamics and performance while maintaining the design features that have set the standard for low complication rates, structural integrity, and durability in a mechanical valve. The SJM Regent® heart valve provides:

- Single digit in vivo pressure gradients—even in valve sizes as small as 19mm.³
- Significantly larger EOAs.³
- Excellent patient-prosthesis match⁴—even in small sizes.
- Significant reduction of LV Mass.³ Numerous studies have shown a direct correlation between LV hypertrophy and morbidity and mortality.⁵⁻¹⁶ Even moderate LV hypertrophy can result in: congestive heart failure,¹⁷ arrhythmias,⁵ myocardial infarction,¹⁷ and sudden death.^{5,17}

"Significant and rapid left ventricular mass regression with all valve sizes confirms excellent hemodynamics associated with the valve. . . ."

(Bach et al, 2001)¹⁸

Traditional reliability

The SJM Regent® heart valve retains the same design features that have made St. Jude Medical® mechanical heart valves the standard for reliability and proven performance for more than 25 years.

- All blood-contact surfaces remain unchanged.
- In seven broad categories of structural integrity and durability tests, the SJM Regent® heart valve met all the demanding standards set by previous St. Jude Medical® mechanical heart valves.²
- The SJM Regent® heart valve is made of the same pyrolytic carbon used in more than one million St. Jude Medical® mechanical heart valve implants.

FAMILIAR TECHNIQUE

The technique for implanting the SJM Regent® heart valve remains the same as for other St. Jude Medical® mechanical heart valves. For added convenience, it's available in two cuff configurations to suit your implant preferences.

FlexCuff™ sewing ring

- Provides enhanced conformability and maximum suture target area while offering lateral flange to accommodate the tissue annulus.

Standard cuff sewing ring

- Rounded cuff designed for excellent conformability and suturability.


Product Specifications*

Valve Size (mm)	Tissue Annulus Diameter (mm)	Overall Height Open (mm)	Implant Height (mm)	Geometric Orifice Area (cm ²)
17	17.0	10.6	5.3	1.87
19	19.0	11.5	5.9	2.39
21	21.0	12.5	6.7	2.90
23	23.0	13.7	7.3	3.45
25	25.0	13.9	7.6	4.02
27	27.0	14.9	8.5	4.69
29	29.0	16.1	9.1	5.44

*From manufacturer's data

Sizes 17 mm and 29 mm are not currently approved in the United States

References:

- Walker DK, Brendzel AM, Scotten LN. The new St. Jude Medical Regent mechanical heart valve: laboratory measurements of hydrodynamic performance. *J Heart Valve Dis.* 1999 Nov;8(6):687-696.
- Product and Design Evaluation Center, St. Jude Medical, Inc.
- St. Jude Medical, Inc., Pre-Market Approval Application—Summary of Safety and Effectiveness. SJM Regent heart valve, P810002/S57.
- Pibarot P, Dumesnil JG. Hemodynamic and clinical impact of prosthesis-patient mismatch in the aortic valve position and its prevention. *J Am Coll Cardiol.* 2000 Oct;36(4):1131-1141.
- De Paulis R, Sommariva L, Colagrande L, et al. Regression of left ventricular hypertrophy after aortic valve replacement for aortic stenosis with different valve substitutes. *J Thorac Cardiovasc Surg.* 1998 Oct;116(4):590-598.
- Jin XY, Gibson DG, Yacoub MH, Pepper JR. Perioperative assessment of aortic homograft, Toronto stentless valve, and stented valve in the aortic position. *Ann Thorac Surg.* 1995 Aug;60(2 Suppl):S395-401.
- Jin XY, Zhang ZM, Gibson DG, Yacoub MH, Pepper JR. Effects of valve substitute on changes in left ventricular function and hypertrophy after aortic valve replacement. *Ann Thorac Surg.* 1996 Sep;62(3):683-690.
- Del Rizzo DF, Goldman BS, Christakis GT, David TE. Hemodynamic benefits of the Toronto Stentless Valve. *J Thorac Cardiovasc Surg.* 1996 Dec;112(6):1431-1445; discussion 1445-1446.
- Lund O. Valve replacement for aortic stenosis: the curative potential of early operation. *Scand J Thorac Cardiovasc Surg.* 1993;40(Suppl):1-137.
- He GW, Grunkemeier GL, Gately HL, Furnary AP, Starr A. Up to thirty-year survival after aortic valve replacement in the small aortic root. *Ann Thorac Surg.* 1995 May;59(5):1056-1062.
- Kannel WB, Cobb J. Left ventricular hypertrophy and mortality—results from the Framingham Study. *Cardiology.* 1992;81(4-5):291-298.
- Levy D, Garrison RJ, Savage DD, Kannel WB, Castelli WP. Prognostic implications of echocardiographically determined left ventricular mass in the Framingham Heart Study. *N Engl J Med.* 1990 Dec 13;323(24):1706-1707.
- Dunn FG, McLenachan J, Isles CG, et al. Left ventricular hypertrophy and mortality in hypertension: an analysis of data from the Glasgow Blood Pressure Clinic. *J Hypertens.* 1990 Aug;8(8):775-782.
- Krayerbuehl HP, Hess OM, Monrad ES, Schneider J, Mall G, Turina M. Left ventricular myocardial structure in aortic valve disease before, intermediate, and late after aortic valve replacement. *Circulation.* 1989 Apr;79(4):744-755.
- Orsinelli DA, Aurigemma GP, Battista S, Krendel S, Gaasch WH. Left ventricular hypertrophy and mortality after aortic valve replacement for aortic stenosis. A high risk subgroup identified by preoperative relative wall thickness. *J Am Coll Cardiol.* 1993 Nov 15;22(6):1679-1683.
- Villari B, Vassalli G, Monrad ES, Chiariello M, Turina M, Hess OM. Normalization of diastolic dysfunction in aortic stenosis late after valve replacement. *Circulation.* 1995 May 1;91(9):2353-2358.
- Christakis GT, Joyner CD, Morgan CD, et al. Left ventricular mass regression early after aortic valve replacement. *Ann Thorac Surg.* 1996 Oct;62(4):1084-1089.
- Bach DS, Goldbach M, Sakwa MP, Petracek M, Errett L, Mohr F. Hemodynamics and early performance of the St. Jude Medical Regent aortic valve prosthesis. *J Heart Valve Dis.* 2001 Jul;10(4):436-442.

Visit our website: www.sjm.com

St. Jude Medical
One Lillehei Plaza
St. Paul, MN 55117, USA
651 483 2000
651 482 8318 Fax

St. Jude Medical Europe, Inc.
The Corporate Village
Avenue Da Vinci laan, 11 - Box F1
B-1935 - Zaventem
Belgium
+32 2 774 68 11
+32 2 772 83 84 Fax

St. Jude Medical Brasil, Ltda.
Rua Frei, Caneca 1380-9° A-CJ91/92
Sao Paulo - SP - Brasil
CEP 01307-002
+55 11 5080 5400
+55 11 5080 5423 Fax

St. Jude Medical (Hong Kong) Limited
Unit 2701-07, COSCO Tower
Grand Millenium Plaza
183 Queen's Road
Central, Hong Kong
+852 2996 7688
+852 2956 0622 Fax

CAUTION: FEDERAL LAW RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN OR PROPERLY LICENSED PRACTITIONER.

St. Jude Medical prosthetic heart valves are indicated for use as replacement valves in patients with a diseased, damaged, or malfunctioning native or prosthetic valve. Possible side effects for all valvular implants include, but are not limited to, regurgitation, thromboembolic phenomena, resistance to flow, infection, hemolysis, dysrhythmias, and prosthetic dehiscence or failure. Anticoagulation is recommended for patients with mechanical valve implants. See the physician's manual for a full description of indications, contraindications, side effects, precautions, warnings and instructions for use.

SJM and St. Jude Medical are registered trademarks of St. Jude Medical, Inc. FlexCuff and The Global Leader in Heart Valve Devices are trademarks of St. Jude Medical, Inc. US Pat # 6,007,577 and 5,876,436; other US and foreign patents pending. St. Jude Medical Cardiac Surgery. ©2002, 2004 St. Jude Medical, Inc. All rights reserved. Printed in the USA. ITEM 1420B/0604/3.5M/EN/GRG

The Global Leader in Heart Valve Devices.™

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