**Generic Name: Bilirubin Direct Kit (DSA Method)****Abbreviated name: Bil-D (DSA)****Order Information**

| Cat. No. | Package size              |
|----------|---------------------------|
| DBI0102  | R1 4×20 mL + R2 1×20 mL   |
| DBI0103  | R1 4×32 mL + R2 4×8 mL    |
| DBI0104  | R1 4×48 mL + R2 4×12 mL   |
| DBI0105  | R1 2×250 mL + R2 1×125 mL |

**Intended use**


In vitro test for the quantitative determination of Bil-D concentration in serum and plasma on photometric systems.

**Summary**<sup>1</sup>

Bilirubin consists of conjugated and unconjugated forms. Direct bilirubin is mainly the former and exists in the liver. It is conjugated with glucuronic acid to convert into soluble substance and excreted via the bile ducts. Elevated levels of direct bilirubin are usually related with hepatitis, occlusion of bile ducts and hepatatrophia caused by acute hepatitis. It is reported that when the total bilirubin  $\leq 3$  mg/dL, the ratio of Bil-D/Bil-T can be used for differentiating the haemolytic jaundice with liver or gall jaundice. The critical ratio is 3.3. In the haemolytic icterus, the ratio is lower than 3.3.

**Method**

Diazotized Sulfanilic Acid (DSA) Method

**Reaction Principle**

Direct bilirubin couples with diazo salt at an acid condition to form a red product of azobilirubin. The absorbency increase is directly proportional to the concentration of direct bilirubin in the sample.

**Reagents****Components and concentrations**

|            |                   |            |
|------------|-------------------|------------|
| <b>R1:</b> | Hydrochloric acid | 170 mmol/L |
|            | Sulfanilic acid   | 29 mmol/L  |
| <b>R2:</b> | Sodium nitrite    | 72 mmol/L  |

**Warnings and precautions**

1. For in vitro diagnostic use only.
2. Take the necessary precautions for the use of laboratory reagents.
3. Preservative contained. Do not swallow. Avoid contact with skin and mucous membranes.
4. Disposal of all waste material should be in accordance with local guidelines.

5. Material safety data sheet is available for professional user on request.

**Reagent Preparation**

R1 and R2 are ready to use.

**Storage and stability**

Up to expiration date indicated on the label, when stored unopened at 2-8°C and protected from light.

Once opened, the reagents are stable for 40 days at 2-8°C when refrigerated on the analyzer or refrigerator.

Contamination of the reagents must be avoided.

Do not freeze the reagents.

**Reagent blank absorbency**

The absorbance of reagent blank at 546 nm should be <0.2 A.

**Materials required but not provided**

1. Calibrator and controls as indicated below.
2. NaCl solution 9 g/L.
3. General laboratory equipments.

**Specimen collection and preparation <sup>2</sup>**

1. Serum, EDTA plasma is suitable for samples. Whole blood and hemolysis are not recommended for use as a sample. Freshly drawn serum is the preferred specimen.
2. Use the suitable tubes or collection containers and follow the instruction of the manufacturer; avoid effect of the materials of the tubes or other collection containers.
3. Centrifuge samples containing precipitate before performing the assay.
4. Stability:

Store samples protected from light. Perform the assay immediately.

**Assay procedure**

|                                                                             | Blank  | Sample |
|-----------------------------------------------------------------------------|--------|--------|
| <b>R1</b>                                                                   | 160 µL | 160 µL |
| <b>Dist. water</b>                                                          | 16 µL  | —      |
| <b>Sample</b>                                                               | —      | 16 µL  |
| Mix, incubate for 5 min. at 37°C, read the absorbance, and then add:        |        |        |
| <b>R2</b>                                                                   | 40 µL  | 40 µL  |
| Mix thoroughly, incubate for 5 min. at 37°C, and read the absorbance again. |        |        |
| $\Delta A = [\Delta A \text{ sample}] - [\Delta A \text{ blank}]$           |        |        |

Application sheets for BS series analyzers are available in this document. Refer to the appropriate operator manual for the analyzer-specific assay instructions.

**Calibration**

1. It is recommended to use the Human multi-calibrator from Mindray and 9

g/L NaCl for two-point calibration. Traceability of the multi-calibrator can refer to the calibrator instructions for use of Mindray Company.

**2. Calibration frequency:**

After reagent lot changed.

As required following quality control procedures.

**Quality control**

At least two levels of control material should be analyzed with each batch of samples. In addition, these controls should be run with each new calibration, with each new reagent cartridge, and after specific maintenance or troubleshooting procedures as detailed in the appropriate system manual. We recommend using the Human Assayed Control made by Mindray to verify the performance of the measurement procedure; other suitable control material can be used in addition.

Each laboratory should establish its own internal quality control scheme and procedures for corrective action if controls do not recover within the acceptable tolerances.

**Calculation**

The analyzer calculates the Bil-D concentration of each sample automatically after calibration.

Conversion factor:  $\text{mg/dL} \times 17.1 = \mu\text{mol/L}$

Or:  $C_{\text{sample}} = (\Delta A_{\text{sample}} / \Delta A_{\text{calibration}}) \times C_{\text{calibration}}$

**Reference Intervals<sup>1</sup>**

Each laboratory should establish its own reference intervals based upon its patient population. The reference intervals measured at 37°C listed below were taken from literature:

| <b>Sample Type</b> | <b>Conventional Units</b> | <b>S.I. Units</b> |
|--------------------|---------------------------|-------------------|
| Serum / Plasma     | ≤0.3 mg/dL                | ≤5.13 μmol/L      |

**Performance Characteristics**

Representative performance data obtained from Mindray system (Mindray BS series analyzers / Mindray Bil-D Reagent) is given below. Results may vary if a different instrument, an individual laboratory or a manual procedure is used.

**Limitations-interference**

The following substance was tested for interference with this methodology. Criterion: Recovery within ±10 % of initial value.

| <b>Substance</b> | <b>Level Tested</b> | <b>Observed Effect</b> |
|------------------|---------------------|------------------------|
| Ascorbic acid    | 20 mg/dL            | NSI*                   |

\* NSI: No Significant Interference (within ± 10 %)

**Linearity range**

The Mindray System provides the following linearity range:

| Sample Type    | S.I. Units              |
|----------------|-------------------------|
| Serum / Plasma | 1-260 $\mu\text{mol/L}$ |

If the value of sample exceeds 260  $\mu\text{mol/L}$ , the sample should be diluted with 9 g/L NaCl solution (e.g. 1+ 4) and rerun; the result should be multiplied by 5.

## Analytic Sensitivity/Limit of Detection

The lowest measurable Bil-D concentration that can be distinguished from zero is 1  $\mu\text{mol/L}$  with 99.7% confidence.

## Precision

Precision performance using the CLSI Approved Guideline EP5-A2 to assay serum control appears in the table below<sup>3</sup>. U:  $\mu\text{mol/L}$

| Type of Imprecision | Level I |      |      | Level II |      |      |
|---------------------|---------|------|------|----------|------|------|
|                     | Mean    | SD   | CV % | Mean     | SD   | CV % |
| Within-run          | 13.08   | 0.01 | 0.75 | 59.04    | 0.45 | 0.77 |
| Between-run         |         | 0.15 | 1.14 |          | 0.64 | 1.09 |
| Between-day         |         | 0.16 | 1.23 |          | 0.75 | 1.26 |
| Within-device       |         | 0.24 | 1.84 |          | 1.08 | 1.83 |

## Method Comparison

A comparison between Mindray System (Mindray BS series analyzers /Mindray Bil-D Reagent) (y) and Hitachi/Roche System (Hitachi /Roche Bil-D) (x) using 40 samples gave following correlation ( $\mu\text{mol/L}$ ):  $y=0.9342x + 2.437$ ,  $R^2 = 0.9984$ .

Details of the comparison experiments are available on request.

## References

- 1.Thomas L ed. Clinical Laboratory Diagnostics. 1st ed. Frankfurt: TH-Books Verlags - gesellschaft; 1998. p. 192-202.
- 2.Tietz NW. Clinical Guide to Laboratory Tests.3rd ed. Philadelphia, Pa: WB Saunders Co; 1995:88
- 3.CLSI. Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline-Second Edition. CLSI document EP5-A2 [ISBN 1-56238-542-9. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087 USA, 2008.

## Graphical symbols

In Vitro Diagnostic  
medical device

Batch Code

European  
Conformity

Authorized representative in  
the European Community

Use By

Consult  
Instructions for useTemperature  
Limit

Manufacturer

Catalogue  
number

© 2024 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights Reserved

**Manufacturer:** Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

**Address:** Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P.R.China

**E-mail Address:** service@mindray.com

**Website:** www.mindray.com

**Tel:** +86-755-81888998

**Fax:** +86-755-26582680

**EC-Representative:** Shanghai International Holding Corp.  
GmbH(Europe)

**Address:** Eiffeustraße 80, Hamburg 20537, Germany

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**Referências**

1. Thomas L ed. Clinical Laboratory Diagnostics. 1st ed. Frankfurt: TH-Books Verlags - gesellschaft; 1998. p. 192-202.
2. Tietz NW. Clinical Guide to Laboratory Tests. 3rd ed. Philadelphia, Pa: WB Saunders Co; 1995:88
3. CLSI. Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline-Second Edition. CLSI document EP5-A2 [ISBN 1-56238-542-9. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087 USA, 2008.

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**Fabricante:** Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

**Endereço:** Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P.R.China

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© 2024 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. Tüm hakları saklıdır
**Üretici:**Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
**AT Temsilcisi:**Shanghai International Holding Corp.GmbH(Avrupa)