Familia AW

AW VolumeShare 7

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SECRET DE
AFACERI

Volume Viewer

GE Healthcare's premium 3D Advanced Visualization and Image Processing platform brings you powerful capabilities to help you visualize and interpret your multimodality imaging data with confidence and ease.

AutoLaunch

Preprocessing enabler

Imaging modalities provide a wealth of diagnostic information, but also present a challenge to your Radiology workflow as image volumes continue to increase while your time for reading and reporting decreases.

Volume Viewer provides you with a rich 3D image processing toolset aimed at creating and displaying the views you need with little user input and streamlining interpretation and reporting by providing the visualization tools you need with minimum clicks.

With an intuitive, modern user interface, Volume Viewer helps a wide diversity of users to learn and master the expanding portfolio of sophisticated tools and applications that it provides. This customizable user interface maximizes the real estate allocated for displaying images and provides access to the users' favorite tools directly from the image viewport. In addition, a customizable toolbar provides one click access to commonly used tools to facilitate a productive review.

Overview

Volume Viewer provides excellent 3D visualization and processing capabilities for reading and comparing CT, MR, 3D X-ray, PET, PET/MR and PET/CT datasets. Volume Viewer also features a broad portfolio of high-performance analysis tools, automating routine tasks and helping to make 3D image processing a stress-free component of your routine workflow.

Volume Viewer is available on VolumeShare 7, a multi-modality advanced visualization workflow solution that helps to enhance diagnostic precision and productivity.

Highlights

- **Customizable protocol-driven workflows** supporting CT, MR, PET and Innova 3D XA modalities, designed to help enhance your clinical review productivity
- **Compare multiple exams from different modalities**, at loading² and on the fly²
- **Adaptive layouts** that take advantage of the dual monitor capabilities in landscape and portrait orientations
- **Easy point-and-click segmentation tools** for segmenting structures of interest
- **CT Dual Energy and GSI⁴ support** with dedicated workflows and tools
- **Advanced 4D PET support** with Q Freeze⁵ tool and 4D PET Dynamic workflow
- **Interactive Summary Table** collecting measurements as they are deposited on the images.
- **Rich set of export tools** including export for 3D printing³.

General Features

- **Unified user interface**, with one click access to your favorite tools and more space for displaying images.
- **Customizable protocol-driven workflows** designed to help enhance your clinical review productivity for a wide range of care areas, with 3D and 4D real time interaction.
- **"Smart layout"** protocol, which adapts the screen layout to the type of loaded exams.
- **High resolution, real-time rendering modes:** MIP/Min IP/Average; standard volume rendering and enhanced volume rendering with Volume Illumination¹ on the fly; Navigator view for endoluminal fly-through; lumen and curved views.
- **Compare mode** for reviewing multiple exams and modalities, simplified with **Dynamic Load to load and register²** new series into the current review session.
- Rich set of **2D/3D ROI tools**, including **AutoContour** for 3D semi-automatic contouring of structures of interest (CT, MR, PET), and **ROI color coded** for tissue classification based on voxel values.
- **Advanced tools**, to take advantage of full 3D capabilities: **AutoSelect**, for easy point-and-click segmentation; One or two clicks **Quick Vessel Trace** to analyze all vessels, in curved reformat, lumen, or MPR view; Advanced 3D processing tools.
- **Summary Table** extends reporting capabilities, while collecting measurements as they are deposited on the images.

- **Multiple Export capabilities**, in different DICOM and non-DICOM formats, including **Key Image Notes**.
- **3D Suite³:** STL, OBJ, VRML and 3MF file formats to export single or multi mesh objects, for 3D printing. 3D Suite is not intended for clinical use.
- **Save state**, which allows to save and restore to pause and resume the review at a later time or to facilitate collaboration workflows.

Modality Features

- **Dual energy and GSI⁴ GE Healthcare CT images** supported with dedicated review protocols.
- **Bone VCAR⁵** combines spine labelling on the fly for all CT review and a dedicated layout to improve spine review.
- **MR review support:** Recognition of PSD name, and specific MR parameters at image loading, which simplifies protocol customization.
- **Dedicated review workflows for PET/CT and PET/MR images**, including **PET SUV measurements**.
- **Q.Freeze⁶** combines the quantitative benefits of **4D PET gated** imaging into an image that provides both frozen patient motion and reduced image noise.
- **Dedicated protocol-driven workflows for Innova 3D XA images⁹.**

Platform features

- **"Smart Compression"** technology for fast and responsive client performance on AW Server.
- **Fast automated Pre-processing⁷**

- Multi-sessions switch on AW workstation with **AutoLaunch⁸**.
- Support of Chinese, Japanese, Korean and Russian languages.

System Requirements

- AW Server 3.2 and above, with recommended monitor resolution up to dual 2MP or single 3MP.
- AW 4.7 Workstation and above (some performances limitations can occur with Z800)
- CentricityTM Universal Viewer

Intended use / Indications for use

Volume Viewer is a medical diagnostic software that allows the processing, review, analysis and communication of 3D reconstructed images and their relationship to originally acquired images from CT, MR, X-Ray Angiography and PET Scanning devices. The combination of acquired images, reconstructed images, annotations and measurements performed by the clinician are intended to provide to the referring physician clinically relevant information for diagnosis, surgery and treatment planning.

Regulatory compliance

This product complies with the European Council Directive 93/42/EEC Medical Device Directive as amended by European Council Directive 2007/47/EC. This product or its feature may not be available in some other countries or regions. Please contact your sale associate.

Rx Only

¹Requires Volume Illumination license

Volume Illumination may not be available in all countries or regions.

Please contact your sale associate.

²Requires Integrated Registration license

³Requires 3D Suite license

⁴Requires Volume Viewer GSI license

⁵Requires Bone VCAR license

⁶Requires Q.Freeze license

⁷Preprocessing is compatible with the following purchasable options: AutoBoneTM Xpress, CardIQ Xpress Reveal, CardIQ Xpress Function, CardIQ Xpress Process, Advanced CTC Pro 3D EC, Colon VCAR, and CT Perfusion 4D and FastStroke. ColonVCAR is not for sale in the US.

⁸Requires AutoLaunch license for AW workstation

⁹Requires Innova Volume Viewer license

CONFIDENTIAL

**SECRET DE
AFACERI**

Volume Viewer


General Features

Unified user interface

- The user interface provides one click access to favorite tools and more space for clinical images.
- The page layout is organized to facilitate navigation through the review steps. Guided protocols walk the unfamiliar user through a clear workflow of creating and manipulating volumetric data, helping to reduce the learning curve by providing a consistent workflow.
- The customizable toolbar can be docked at the top, bottom, left, or right, by a simple drag & drop. This preferred location is saved for future use for each user.
- An extensive set of advanced tools are grouped into categories on the toolbar to facilitate quick retrieval when needed.
- Each category and tool are easily configurable by a simple drag & drop in the Toolbar so that they are only one click away for your future reviews
- To reduce long mouse trips, up to 6 favorite tools can be accessible from each image via a right mouse menu.
- Customization of the toolbar and tool behaviors with a user's preferences is provided in a single intuitive Preferences menu.
- Keyboard shortcuts and quick access tools for on-viewport controls

The following mouse modes are available to help you manipulate the images:

- Interactive VR adjustment lets you adjust VR opacity interactively with the mouse as Window Level/Window Width
- Direct Paging (continuous or contiguous)
- Free hand rotation on 3D and oblique views
- Percentage fusion, to easily change the transparency of objects on fused images. There are several fusion modes available for your use.

- Mag glass displays zoomed-up imagery around the cursor.

Protocol management and loading

- Customizable protocol-driven workflows for a wide range of clinical applications, including support for multiple review steps. You are guided through the review steps by a Navigation bar at the top of the User Interface.
 - Review protocols created by user with Volume Viewer can be saved as favorite protocols and accessed easily in routine.
 - Progressive Load enhances image loading performance. This feature has two modes:
 - Images are loaded in sequential mode (e.g., Reformat)
 - Images are loaded in interlaced mode (e.g., 3D/VR)
- In either case, image review can begin as soon as the first image is displayed.
- Navigation through series and exams can be performed without exiting the patient list.

Smart layout


- Volume Viewer General Review protocol presents an adaptive layout for single and multi-modality PACS-like reading of data. This protocol distributes and displays all loaded series of a single exam in an intelligent layout to maximize use of available monitor space. With multiple exams, the protocol automatically detects and distributes "like" series for optimal comparison.

Advanced visualization capabilities

Volume Viewer offers various capabilities to display advanced rendering modes:

- 3D Volume Rendering / Volume Illumination¹ images
- MPR views with different slice thickness. The following rendering modes are available for all Thick Slab: MIP, Min IP, Volume Rendering, Average.
- Navigator views which display interactive endoluminal views exportable as movies.
- Curvilinear reformatting allowing display of curved, lumen and cross section views for various structures (vessels, spine, etc).
- Merge multiple volume rendered models into a single view.

Volume Rendering

Volume Illumination

Predefined cut planes

Volume Viewer lets the user define cut planes to isolate specific structures in the VR model, in 16 configurations, including Left / Right / Inferior / Superior / Anterior / Posterior / Front. You can then display the voxels values on each plane.

Multiple Volume Rendering and Volume Illumination¹ objects

- Volume Viewer allows you to merge up to 8 Volume Rendered models from the same series into a single 3D view.
- You can adjust independently the threshold, colors, and transparency of each Volume Rendering or Volume Illumination


Auto-center fly through with smart cursor

- In fly-through studies (airways, colon, angiography) navigation along the centerline of the structure of interest can be recorded step by step and can be exported as a movie.
- In addition to the regular 180° or less Navigator view, Volume Viewer enables a Fish-Eye View with any wide camera angle value from 180° to 360°. This provides a view of structures both in front of and behind the users' virtual location on the same image.

Lumen view

- Lumen View provides an unfolded 3D view around a user defined centerline. The lumen view can be interactively adjusted (rotation around the centerline, width, field of view).

3D and 4D capabilities

- Volume Viewer allows 3D cursor synchronization in any orientation on the fly.

- Real time reformat in oblique planes is available for your exams, as well as simple and double oblique interactive modes.
- Quick access to cross reference display gives you the ability to display the current slice location or all slices, or other planes.
- Real time interaction in 4D mode, by using Cine tool on all planes and 3D views for CT, MR, and PET multiphase data.

Compare mode


- With Volume Viewer several volumes can be loaded either from the same exam or from different exams. This is useful for exam comparison (Compare Mode) or for Multiphase examinations (of the liver for example). The Compare Mode is also available for PET/CT examinations and takes advantage of a dual screen configuration.
- Multi-phase CT and MR data are listed as separate phases in the Series Selection screen, allowing selection of all phases or a subset of phase data for review. The Series Selection Panel displays the multi-phase data split into arterial and portal phases, for example.
- Dynamic Load, compatible with both Volume Viewer and Fusion protocols, lets you drag & drop 3D volumes from CT, MR, PET, and 3D XA modalities into a desired viewport. Together with the Integrated Registration² option, Dynamic Load allows new volumes to be registered and loaded on the fly. Save State series can be restored as a separate session. For dual monitor configurations, a Save State series can be displayed on the right monitor together with a current session on the left monitor. For single monitor configurations, a user may simply switch between the two sessions.
- The zoom and pan functions are propagated to all images displayed in the same orientation.

2D/3D ROI for quantitative measurements


- Volume Viewer contains a set of basic 2D/3D tools: distance, angle, report cursor, arrow tool, annotation, free hand ROI tool, that can be placed on fused images as well.
- Easy deposition and labelling of measurements.
- All 3D ROIs can be customized, and color coded to display statistics computed on different intensity range and help assessing tissue classification based on voxel values.
- All measurements are considered bookmarks and can be accessed by navigating through the list of measurements.

AutoContour tool

- AutoContour provides consistent tools for 3D contouring ROI, with the same workflow used for CT, MR, and PET images.
- For CT and MR images, the contouring tool will suggest contour boundaries after defining 1 point or a diameter on a structure of interest, which can be edited prior to validation.

- For PET images, the contouring tool has been adapted for standard SUV calculations.

Advanced 3D tools

Segmentation tools

- AutoSelect segmentation tool helps adding or removing any continuous structures of interest by simple point and click.
- You can also access advanced segmentation tools such as Scalpel, Threshold, Dilate, Erode, Subtraction, Intersection, Addition, Filter Floaters.

Quick Vessel Trace (CT and MR images)

- You can perform vessel analysis with one or two points clicks from any CT or MR Angiography images.
- Quick Vessel Trace automatically extracts the vessel after user clicks and launches vessel tracking for fast review in curved reformat, cross-sectional, best L-section, lumen, and MRP view.

- Tracking can be extended proximally or distally for a full view of the vessel.

Summary table

- The Summary Table collects measurements and organizes them according to finding and exam date.
- The Summary Table allows interactive navigation through findings in the image viewports and can be docked or expanded.
- The Summary Table can be copied into the client clipboard and then pasted in email, reports, office software etc. (on AWS only)

Summary Table		Findings	
Findings	Date	Measurements	Measurements
1 - Finding 1	Sep 11 2007	CT 1	CT 2
		ROI 1	ROI 2
		2Dmax	2Dmax
		Short Axis	Short Axis
		Vol	Vol
2 - Finding 2	Oct 17 2007	CT 3	CT 4
		ROI 3	ROI 4
		2Dmax	2Dmax
		Short Axis	Short Axis
		Vol	Vol

Export capabilities

Volume Viewer contains multiple standard options for exporting the results of a review session:

- Save images to the database, as new DICOM series.
- Save the significant images as Key Image Note objects in the database. Key Image Notes and End Review allow you to flag images of interest as Key Image Notes (IHE profile) and push them to the archiving system when you exit the application.
- Save State object are used to save the work in a new series containing all the post processing data (3D Model, displays, measurements, annotations, etc), for future review.
- Save STL, OBJ, VRML or 3MF file from 3D Volume Rendering / Volume Illumination¹ viewport for 3D printing purposes through 3D Suite³

Additional export tools are available to create specific sequences of images to be exported:

- The Batch tool allows creating a sequence of reformatted images or a sequence of rotating 3D views.
- The Movie tool creates a comprehensive movie including different rotations, zooms, and pan of the image, which can be exported as DICOM series or mpeg file.
- Quick Export: Exports in a single click a full batch of contiguous images at the displayed thickness for 2D images, or a batch of rotations of a 3D View.
- Cardiac Review and Export: Processing and reviewing cardiac exams for CT, MR and PET with manual oblique


reformatted protocols can be exported into a multi-phase Cine movie that allows the referring physician to review the exam in a dynamic mode.

- Capture and send images from the viewport to your computer clipboard on the fly (Ctrl+C/Ctrl+V on AW Server only).

Modality Features

CT image quality


- CT filters to enhance or smooth CT images to help improve image quality on the fly.
- The lung filter helps enhance contours of images reconstructed in standard mode for excellent visualization of lung structures

Standard recon w/o filter

Standard recon with Lung filter

- A smoothing filter reduces noise while maintaining high image contrast.

No filter

Smooth filter

No filter

Smooth 3D+ filter

CT Dual energy protocols

- Three protocols are available to enable review of dual energy images acquired on compatible GE Healthcare scanners and enhancing pixels with specific ratio: DE Calcium ($HU_{80}/HU_{140} \geq 1.25$), DE Uric Acid ($HU_{80}/HU_{140} \leq 1.25$) and DE Custom (custom threshold on HU_{Low}/HU_{High}).

GE imagination at work

- In addition, all Volume Viewer standard protocols are compatible with GSI⁴ GE Healthcare CT images.

Bone VCAR⁵


- Bone VCAR is available in Volume Viewer toolbar and can be launched on the fly
- Vertebrae are automatically recognized and labelled (deep learning algorithm)
- A dedicated protocol provides a specific layout of the spine including curved views

Review of MR studies

Volume Viewer offers support for review of MR studies, through dedicated protocols.

- The MR-specific parameters (PSD name, image weighting, scan plane, etc) are recognized at loading and are used to customize the layout display.
- Anatomy-based protocols for standard review (e.g. spine, brain stroke...), with dedicated and guided review steps, to help review standard MR exams.
- Dedicated MR Cardiac Viewer for single and multiple Cine Review. Enables comparison between cardiac series such as Time Course and Myocardial Delayed Enhancement (MDE) images.
- Support of MR Multi Echo, Multiphase, Diffusion series.
- Support of parametric series (e.g., ADC, MTT, and other parametric maps) created in READY View, allowing direct functional measurements.


- 2D and 3D ROI propagation to other phases / series / exams allowing you to easily correlate information from multiple sources.
- Improved MR annotation consistency.
- In combination with Integrated Registration², enables direct access to MR image registration (inter/intra exam), image fusion and Whole Body MR Review protocols
- Ability to automatically bind several MR axial series corresponding to sequential axial locations into a single series. This is useful for consolidating multi-stage acquisitions for displaying Whole Body MR scanning.

SUV measurements on PET images

- Volume Viewer supports Standardized Uptake Values (SUV) for image display and measurements. Several SUV scales are available like SUVbw, SUV lbm, SUVbsa, as well as SUV Peak.
- Window/level presets may be user-defined for PET images based on SUV values.
- The SUV values are available in all the basic 2D/3D ROI tools of Volume Viewer and in the AutoContour tool.

PET/CT and PET/MR dedicated protocols

- Volume Viewer contains a list of predefined protocols allowing the review of PET/CT and PET/MR data. The screen layouts contain fused views between the morphological and functional images, as well as 3DMIP rendering of PET images.
- Factory protocols are customizable for your own review.
- PET 4D protocol allows users to load and display gated and dynamic PET series, as well as summing or reframing these series.

PET IQ improvements

- Enhanced 3D visualization algorithm to fully support Q.Clear PET images.
- Q.Freeze⁶ option provides a respiratory motion-corrected Gated PET volume, statistically similar to conventional static PET with significantly reduced or eliminated blurring effects due to patient respiration. It provides also a corrected Gated PET series where counts are summed back to each respiratory phase.

Dedicated features for 3D XA images⁹

The optional Innova 3D XR application produces XA images (512x512 images or 256x256 images depending on your settings). This format opens the Volume Viewer world to Innova 3D XA images:

- Single exam: Innova Navigator, Curved Reformat, Multi Oblique, MPVR 3mm Average, Segment Structure
- Compare exam: Compare Axial, Compare Coronal, Compare Sagittal, Reformat
- Layouts and protocols customization
- The Innova options provide access to Workflow booster tools, easy multi-modality review "on-the-fly" dedicated Innova review protocols with dedicated Review Steps.

Smart compression

- The Smart Compression technology automatically displays images with full fidelity when the image is still, then uses the selected compression level for increased interaction speed during user interaction. This allows for diagnostic reads on full fidelity static images with responsive dynamic display even at low bandwidth. On-image visual indicators notify the user when compression is in effect.

Pre-processing⁷ features on AW Server

- Preprocessing is a Volume Viewer feature that automatically performs routine processing tasks and saves the results so they are ready for you. When new exams are transferred to the database, Volume Viewer recognizes user-defined keywords in the Series Description and launches the appropriate pre-processing. The results are stored in a Save State object with the original study so you can load it when you are ready for review.

AutoLaunch⁸ features on AW workstation

- AutoLaunch is a Volume Viewer feature that automatically preloads exams, so they are ready for your review.

- When you are processing a study and a new exam is transferred to the AW Workstation, this feature automatically launches the new exam with an appropriate Volume Viewer protocol in the background. When you are ready, one click in the AutoLaunch window instantly brings up the exam in the Volume Viewer application, eliminating waiting time and extra steps to load the exam into computer memory for processing.
- AutoLaunch is compatible with CT, MR and PET single volume protocols of Volume Viewer.
- When combined with optional applications⁸, this feature gives access to data already preprocessed, that you can review by just clicking on AutoLaunch window.

Summary of Operation

Volumetric models are loaded by selecting the exam or series. The user can select a protocol category from an anatomical selector or go directly to a Review Layout. In either case, images are loaded progressively in the background; this gives control to the user in just a few seconds after selecting the images. Selecting a Review Layout launches a volumetric display protocol with predefined layout preferences. Review Layouts may be customized and combined to suit the user's workflow. Selecting a protocol category unlocks a variety of visual protocols that include the layout, threshold, rendering mode and filming formats. Some of these protocols direct the user through the process providing capabilities to interactively view and manipulate the model, increasing productivity and consistency for all modalities.

References

¹Requires Volume Illumination license

Volume Illumination may not be available in some other countries or regions.

Please contact your sale associate

²Requires Integrated Registration license

³Requires 3D Suite license

⁴Requires Volume Viewer GSI license

⁵Requires Bone VCAR license

⁶Requires Q.Freeze license

⁷Preprocessing is compatible with the following purchasable options: AutoBone™ Xpress, CardIQ Xpress Reveal, CardIQ Xpress Function, CardIQ Xpress Process, Advanced CTC Pro 3D EC, Colon VCAR, and CT Perfusion 4D and FastStroke. ColonVCAR is not for sale in the US.

⁸Requires AutoLaunch license for AW workstation

⁹Requires Innova Volume Viewer license

CONFIDENTIAL

SECRET DE
AFACERI

GE Medical Systems SCS
283 rue de la Miniere
78530 Buc France

DOC1460639 rev9

General Electric reserves the right to make changes in specifications and features, or discontinue the product or service described at any time, without notice or obligation. This does not constitute a representation or warranty or documentation regarding the product or service featured. Illustrations are provided for informational purposes, and your configuration may differ.

This information does not constitute legal, financial, coding, or regulatory advice in connection with your use of the product or service. Please consult your professional advisors for any such advice. Operation of GE Healthcare products should neither circumvent nor take precedence over required patient care, including human intervention of healthcare providers. GE Healthcare products and services do not code medical procedures. Accurate coding is the responsibility of the provider or billing professional.

GE, the GE Monogram, imagination at work, Centricity and AutoBone are trademarks of General Electric company.

General Electric Company, by and through its GE Healthcare division.

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
General Electric Company, prin intermediul diviziei sale GE Healthcare.

©2019 General Electric Company

[Scrieți aici]


GE imagination at work

CONFIDENTIAL

SECRET DE
AFACERI

VessellQ™ Xpress with AutoBone™ Xpress

Fast and Easy to Use Vessel Analysis with Automated Real-Time Tracking, Auto-Labeling and Simplified User Interface.

Peripheral vascular disease (PVD) is a condition similar to coronary artery disease and carotid artery disease. In PVD, fatty deposits build up in the inner linings of the artery walls. These blockages restrict blood circulation, mainly in arteries leading to the kidneys, stomach, arms, legs and feet.

There are different clinical techniques used to diagnose PVD. Computed Tomography Angiography (CTA) can potentially provide sufficient information for therapeutic decision-making and reduced imaging costs for the management of PVD. One of the advantages of computed tomography (CT) is that it shows both vascular & nonvascular anatomy and also provides 3D information for best view. Also, CT may be more useful for patients with limited or no peripheral access or with stents, who may not be candidates for magnetic resonance angiography (MRA). Complementary to the ability to perform CTA is to perform clinical assessment by a dedicated post-processing tool.

Overview

VessellQ™ Xpress with Autobone Xpress is a post-processing software package for the Advantage Workstation (AW) platform, AW Server, CT scanners and PACS reading stations. **It is a tool to analyze CT angiographic data including stenosis analysis, thrombus, pre and post stent planning procedures, and directional vessel tortuosity visualization.** Autobone Xpress provides zero-click bone segmentation for head, neck and other anatomy. The automation and ease of use streamlines workflow.

Highlights

- Automated detection of aorta and iliacs with auto-labeling of vessels
- Automated real-time fast tracking for all vasculature
- New contextual wheel menu including easy access to advanced options
- New centerline edition capabilities
- Improved predefined measurements and edit capabilities
- New reporting interface which includes measurements in summary table
- New user interface

CONFIDENTIAL

SECRET DE
AFACERI

Visit us:

<https://www.gehealthcare.com/products/advanced-visualization/all-applications/autobone-vessellq-xpress>

Features

- Anatomy based protocols for guided workflow in the analysis for your peripheral vascular CTA acquisitions.
- Automatic bone removal for efficient review of vasculature.
- Automated abdominal aorta CTA protocol with vessel tracking and vessel labeling.
- Fast tracking detects entire length of vessel and displays the centerline based on a single start point.
- Advanced editing tools bringing new flexibility for faster editing based on targeted anatomy.
- Thrombus detection with measurements and comparison over time.
- Dynamic AVA provides access to vessel tracking at any time, from any protocol, on any vessel.
- **Advanced measurement tools for stenosis detection** and color identification of HU value for quantification of different densities.
- Create a queue of cases that are pre-loaded and pre-processed in the background while you read another exam.
- Measurements are organized in the summary table allowing for easy viewing and exporting.
- Fully customizable workflows with adjustable layouts, personalized parameter and settings, custom review steps.

- Save State lets you save and restore the state of the processed images at any stage

System Requirements

- AW Workstation AW 4.7 ext 14 or higher
- AW Server 3.2 Ext. 3.2 or higher

Indications for Use

VessellQ Xpress is an optional, non-invasive, optimized, post-processing application intended to provide images and tools to analyze vascular anatomy and pathology, aiding physicians in diagnosis and determination of treatment paths, from a set of Computed Tomography (CT) Angiographic images. VessellQ Xpress is an option for the Advantage Workstation (AW) platform, CT Scanner, and/or PACS, which can be used in the analysis of 2D and 3D CT Angiography images/data for the purpose of cardiovascular and vascular disease assessment. This software-only device is designed to support physician assessment for a wide variety of clinical uses such as stenosis analysis, pre/post stent planning, pre/post valve replacement planning, and directional vessel tortuosity visualization. VessellQ Xpress' automatic visualization tools provide users the capability to segment bony structures for accurate identification of the vessels. Additional tools enable analysis of the vascular anatomy including the aorta, valves, and branching vessels for: anatomical sizing;

density and volume analysis of segmented vasculature and calcified / non-calcified plaque; and measurements of abnormalities. The TAVI Analysis option for VessellQ Xpress is a planning tool used for Trans Aortic Valve Implantation (TAVI) procedures. It automatically segments the aorta and displays the aortic valve in multiple views for measurements of anatomic structures commonly needed for aortic annulus replacement planning. TAVI Analysis provides guided workflow and semi-automated tools to aid in evaluation of appropriate access pathways for interventional procedure planning.

AutoBone Xpress option is a software package that is intended to facilitate segmentation of bony structures and calcifications from abdominal and extremity CT Angiography data.

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SECRET DE
AFACERI

GE imagination at work

© 2021 General Electric Company.
All rights reserved. Data subject to change.
GE, GE Monogram, imagination at work, AutoBone Xpress, VessellQ Xpress are
trademarks of General Electric Company.
All other third-party trademarks are the property of their respective owners
Doc1471990

Features Detail

Anatomy-based protocols

- The program offers you predefined protocols for fast, easy review of all vascular structures:
 - Intra-cerebral
 - Carotids
 - Thoracic Aorta
 - Abdominal Aorta
 - Runoff
 - Spinal
 - Extremity
- Protocols may be customized to meet your requirements and saved for quick access.

AutoBone™ Xpress

- Anatomical driven protocols with automatic bone removal.
- Improve vessel visualization by removing obstructive bony detail


- Automatic calcium removal for vessel lumen visualization

Auto Aorta

- With a single click, the program automatically removes bones and performs centerline tracking of the aorta and iliacs and labels vessels.

- Simply confirm the automated tracking or modify it by adding additional branch points, deleting a single point, or deleting multiple points from any view.


Fast Tracking

- Semi-automatic detection of complete vessel from a single starting point

- Advanced editing with variable interpolation tool for quick centerline edits adjusted for length of vessel.

- Bulls eye adjustment capability on x-section views to perfect tracking of the smallest most tortuous vessels.
- Volume measurements of the thrombus allow you to track change over time. Additional measurements include :
 - Minimum diameter
 - Maximum diameter
 - Mean diameter
 - Cross-sectional area

Vessel Label Data Base

- Each protocol is associated with a database of vessel names
- Customize vessel label database to meet your requirements for every anatomical region.
- A right mouse click provides access to the list of vessel names per anatomy.


Dynamic AVA

- Real-time tracking of the vessel center line and instant visualization of an unfolded view for quick vessel inspection.
- Immediate access to vessel analysis tools.
- Edit contours with a selection of drawing tools with variable interpolation for fast editing.


CONFIDENTIAL

SECRET DE
AFACERI

Advanced Vessel Analysis Tools

- Semi-automated vessel analysis lets you track, extract, visualize, and size vessels.
- Track a vessel and its branches for quick review in curved, oblique longitudinal or cross-sectional views.
- Multi-point deposition capabilities for improved center line tracking of challenging cases.
- When you edit a vessel contour, the program automatically interpolates the edit to neighboring slices.
- Generate and save rotating movies
 - Curved planar reformation and cross-sectional (lumen) views
 - Best L-section reformatted images
- Select vessel segments for quantitative or qualitative analysis.


- Available measurements include:
 - Distance
 - Cross-sectional area

- Single or dual reference point comparison
- Relative percent stenosis / aneurysm
- Color Identification
- 3D Angle
- Thrombus Extraction
- Volume
- Mean, min, max diameter
-
- Move freely between review steps and protocols without leaving the exam.

Stenosis Detection & ColorID

- With the customizable one- or two-point deposit technique, you can automatically quantify degree of stenosis and segment length.
- By assigning customizable color mapping to CT number ranges within an image, you can identify plaque on an axial image, curved reformat, or MPVR image. Up to ten levels of color coding are available. Colors transition smoothly from one density to another.
- Make volume and area measurements to track plaque size over time.

- Custom sizeable 3D ROI lets you apply color over all anatomy.
- Generate multiple ROIs on one vessel or multiple vessels.
- Automatically include each defined plaque, its volume and area into the report.

Quick AVA

- Dynamic AVA, provides access to vessel analysis from any protocol.
- Launch the Quick AVA feature from the main tool bar.
- Perform vessel analysis from any 3D or reformatted image with one- or two-point clicks.
- If you detect a lesion, you can analyze it by dropping a point above and below the section. The vessel analysis tool launches automatically for fast review in curved reformat, cross-sectional, lumen, and MRP view.
- Extend tracking proximally or distally for extended views at any time.

Summary Table & Report

- Measurements are automatically captured in summary table for easy review.


- Simply click on any measurement in the summary table and you will be taken to that measurement location.

CONFIDENTIAL

SECRET DE AFACERI

- Prepare report with a single click or generate an advanced report with customized options
- Send the report to database, filer or PDF.

Save State & Protocol

- Save the current state of tracking and measurements for easy reloading with a single mouse click.
Save state captures:
 - 3D images
 - Vessel tracking
 - Measurements
 - Vessel contours.
- Transfer the Save State to the PACS or other AW workstations for future review.
- Save custom protocols including anatomical labels and measurements.

Summary

VessellQ Xpress provides user- friendly tools an intuitive interface to help streamline your vessel analysis workflow and help you generate more comprehensive reports for referring physicians.

Multiple Study Comparison

- With VessellQ Xpress, compare a previous study with the current study. You also have the flexibility to launch comparison within Quick AVA.

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
**SECRET DE
AFACERI**

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CT Perfusion 4D

Comprehensive perfusion analysis available for neurology, oncology and cardiology.

Easy-to-use, protocol-driven workflow for all organs.

Perfusion 4D has been expanded beyond stroke and oncology to now include dynamic myocardial perfusion. For any organ where blood perfusion is of interest Perfusion 4D is a complete package to quantitatively provide the necessary information to guide treatment decisions.

The protocol-driven design leads the user step-by-step through the process, reducing keystrokes and improving repeatability. Get the information you want quickly and reliably.

CONFIDENTIAL

SECRET DE
AFACERI


Overview

GE's Perfusion 4D is a fast, easy-to-use automated software for analyzing CT Perfusion images related to stroke, tumor angiogenesis and dynamic myocardial perfusion. Its simple user interface and automated perfusion post-processing make it easy for you to diagnose quickly and accurately – and make treatment decisions more confidently.

GE leverages its deconvolution expertise, with its innovative delay correction algorithm. Perfusion 4D takes perfusion assessment to the next level, with the addition of CT Dynamic Myocardial Perfusion. Providing the classic color overlays with quantitative measurements to evaluate myocardial blood perfusion.

Highlights

- Improved IQ of functional maps in presence of noise
- Improved processing times
- Improved ROI mirroring for brain protocols
- Deep Learning algorithm for ventricle removal¹
- Automated midline detection for brain symmetry
- Absolute and relative values with customizable thresholds for tissue classification¹
- Permits injection rates of 4cc/sec
- Dynamic registration
- Dynamic myocardial perfusion maps

Visit us:

<https://www.gehealthcare.com/en/products/advanced-visualization/all-applications/ct-perfusion-4d-neuro>

Visit us:

<https://www.gehealthcare.com/en/products/advanced-visualization/all-applications/ct-perfusion-4d-multi>

Features

- Visualize all information in true volumetric form with the ability to employ all the volumetric-based image analysis tools in the AW Volume Viewer.
- Deep Learning algorithm for brain ventricles removal¹
- Automated midline detection for brain symmetry
- Absolute and relative values with customizable thresholds for tissue classification¹
- Background processing will send processed images and functional maps in email² format to stroke team²
- Whole heart perfusion with semi-automatic guided workflow to aid in assessment of myocardial ischemia
- Smart Maps using the intelligent 4D noise suppression algorithm improves functional map image quality in the presence of noise.
- Dynamic non-rigid registration protocols for dynamic acquisitions within the body.
- Productivity enhancements include:
 - Ability to select automated or manual vessel detection for Brain Perfusion
 - Faster processing times for functional map computations
 - Volumetric visualization of functional maps
 - Access from anywhere using AW Server.
- Time stamps on functional maps and Regions of interest (ROIs).
- Interactive large vessels exclusion.
- Many basic features, including the ability to:
 - Process double phase acquisition data.
 - Create 4D Regions of Interest (ROIs).
 - Display ROIs statistical information.
 - Display averaged time intensity information of the ROIs.
 - Simultaneously review any number of functional maps.
 - Customize protocols
 - Save parameters
- **Includes CT perfusion protocols for stroke, brain & body tumor, myocardium as well as dynamic registration**
- Functional maps:
 - Regional Blood Volume (BV)
 - Regional Blood Flow (BF)
 - Regional Mean Transit Time (MTT)
 - Contrast arrival delay (IRF T0)
 - Transit time to peak of impulse residue function (TMax)
 - Mean Slope of Increase
 - Capillary permeability surface area for lengthened acquisition protocols (PS)

- Average image
- Base image

Tissue Classification

System Requirements³

- AW Workstations AW 4.7 Ext. 14 or higher
- AW Server 3.2 Ext. 3.2 or higher

Regulatory Compliance

This product complies with Regulation (EU) 2017/745 of the European Parliament and of the Council on medical devices (MDR). This product or its features may not be available in some other countries or regions. Please contact your sales associate.

Rx Only

Indications for Use

CT Perfusion 4D is an image analysis software package that allows the user to produce dynamic image data and to generate information with regard to changes in image intensity over time. It supports the analysis of CT Perfusion images (in the head and body) after the intravenous injection of contrast, in calculation of the various perfusion-related parameters (i.e. regional blood flow, regional blood volume, mean transit time and capillary permeability). The results are displayed in a user-friendly graphic format as parametric images.

This software will aid in the assessment of the extent and type of perfusion, blood volume and capillary permeability changes, which may be related to stroke or tumor angiogenesis and the treatment thereof.

CONFIDENTIAL

SECRET DE AFACERI

1. Requires CT Perfusion 4D Neuro Enhanced license

CT Perfusion 4D Neuro Enhanced may not be available in all countries or regions. Please contact your sales associate.

2. Requires Send by Email license and AW Server version with email service, not compatible with stand-alone AW systems

Send by Email may not be available in all countries or regions. Please contact your sales associate.

3. Send by email functionality requires an AW Server version with email service available

GE imagination at work

Features Detail

Basic Features

- Process double phase acquisition data (data with minimal temporal sampling and data with longer temporal sampling in one acquisition to cover entire functional perfusion mechanism).
- Create 4D ROIs (multi-locations ROI with time-varying position).
- Display ROI statistical information (Area, Min., Max, Std Dev.).
- Display input ROI from several locations on the same graph.
- Display the averaged time-intensity information from an ROI (VOI).
- Simultaneously review any number of functional maps.
- Customize any viewport to series, graph or functional view.
- View series or functional data in axial, sagittal, coronal or oblique view.
- View series or functional data in 3D or VR mode.
- Customize protocols (advanced settings, color ramp and/or W/L used for each generated functional map) and to save these user preferences.
- Save the state of a protocol including parameters, ROIs, functional maps and pre/post-processing results for further reloading in CT Perfusion.
- Save all images in the database.
- Save ROIs statistical information using screensave.
- Print to film or color paper printer.
- Save parameters, used to generate functional maps, as DICOM objects and to restore them including ROIs with the corresponding maps in CT Perfusion.

Protocols

CT Perfusion 4D has fully integrated protocols that assist the user in obtaining results efficiently for neurology, oncology and cardiology care areas.

CT Perfusion 4D package protocols are:

CT Brain Stroke Protocol. The CT brain stroke protocol will open a loading panel containing information about the selected perfusion series in a tabulated format. It is possible to remove unwanted time points from the loading sequence.


Once the exam is loaded the SW will perform the following algorithmic steps:

- **Image Registration** to minimize the effects of patient movement.
- **Bone Removal** to remove the cranium automatically from the dataset.
- **Arterial Input** used as the input for the deconvolution algorithm. Artery is automatically selected and can be easily adjusted if needed.
- **Venous output** for partial volume correction. Vein is automatically selected and can be easily adjusted if needed.
- **Last pre-enhancement / First post-enhancement** image are automatically selected to ascertain correct input for functional maps calculations.
- **Final Settings and Compute.** The software will compute the functional maps with the saved settings, but can be made to stop and recompute the functional maps if changes to the algorithm settings were implemented. The changes may include any of the following:
 - Algorithm Resolution, Smart Smoothing size.
 - Maximum Blood Flow, Brain Density, Hematocrit Ratio.

The software will use deconvolution to generate all functional maps, *Blood Flow, Blood Volume, Mean Transit Time, and Transit Time to IRF Peak (Tmax)*.

The Brain Stroke protocol has multiple review steps to guide the user through the functional information.


- **Review settings:** When loading the dynamic exam to process, this first review step is initiated. The general purpose of this step is to confirm the automatic selection of the input vessels for the brain stroke or tumor protocols, or to conduct a general review of the time resolved perfusion data and place ROI's to select vessels.


- **Set Symmetry:** The software will automatically define the symmetry plane to be used for mirroring ROI's and relative thresholds. For the brain stroke protocol, the application provides a "Check view" image for a fast, intuitive way to indicate the plane and reference hemisphere identified.


If edition is needed a dedicated review step allows users to quickly adjust the symmetry axis.

- **Review Functional Maps** enables the viewing and assessment of functional maps such as Blood Flow, Blood Volume, Mean Transit Time, and Tmax. ROI mirroring functionality makes it easier to compare hemispheres.


- **Tissue Classification,** specific to the brain stroke protocol, enables the visualization of regions that are segmented from absolute or relative values, customizable thresholds and user selectable input maps.

The Brain Stroke protocol is typically used for the assessment of cerebral hemodynamic disturbances e.g. stroke and traumatic brain injury.

- **Email,** when working on AW Server and invoking the preprocessing for full automation of the Brain Stroke protocol, the processed images can be configured to auto send to the stroke team in email format.

CT Brain Tumor Protocol. Essentially identical to the brain stroke protocol but without the Tissue Classification step.

GE imagination at work

CONFIDENTIAL

SECRET DE AFACERI

CT Body Tumor Protocol. The Body tumor protocol produces functional maps automatically once the user selects an artery using a simple and intuitive workflow.

The Body Tumor protocol also leverages most of the advanced algorithm settings that are part of the Brain Stroke protocol. The protocol has the following algorithms

- **Noise Threshold.** Thresholds can be customized to exclude air and bone from the calculations.
- **Artery Input.** System allows for an intuitive selection of the artery.
- **Last pre-enhancement / Last post-enhancement image / Last Second Phase Image** (last image of the phase 2 data). These can be selected to ascertain correct input into the parameter calculations.
- **Final Settings and Compute.** The user can compute immediately or customize the following advanced settings:
 - Algorithm Resolution, Smart Smoothing.
 - Maximum Blood Flow, Brain Density, Hematocrit Ratio.

The Compute command then engages the exclusive de-convolution algorithm to produce the parametric images.

CT Pancreas Perfusion Protocol. Identical to the Body Tumor Protocol but optimized for pancreas tumors.

CT Prostate Perfusion Protocol. Identical to the Body Tumor Protocol but optimized for prostate tumors.

CT Kidney Perfusion Protocol. Identical to the Body Tumor Protocol but optimized for kidney tumors.

CT Soft Tissue Perfusion Protocol. Identical to the Body Tumor Protocol, with an additional step to define the vein output for partial volume correction (automatic or manual venous selection).


CT Spleen Perfusion Protocol. Identical to the Body Tumor Protocol but optimized for spleen tumors.

CT Bone Perfusion Protocol. Identical to the Body Tumor Protocol but optimized for bone tumors.


CT Non-Rigid Registration. Fully automated protocols used for dynamic body acquisitions to register data for use in perfusion analysis. Provides simultaneous review of original data with the registered results for full confidence in the newly created series. Saved registered results are to be loaded into the appropriate body or cardiac perfusion protocol for functional analysis

Two registration protocols are provided. CT Cardiac Dynamic Registration and CT Liver Dynamic Registration.

CT Liver Perfusion Protocol. Identical to the Body Tumor Protocol, with an additional step to define the portal vein input (automatic or manual venous selection).

Processing times have improved exponentially, now with an average processing time of one to three minutes depending on slice thickness. **Hepatic Arterial Fraction maps have been improved and can be an efficient tool for tumor detection. A new map, hepatic arterial blood flow (HABF), has been added to assist in tumor assessment.**

GE imagination at work


CONFIDENTIAL

SECRET DE
AFACERI


Dynamic CT Myocardial Perfusion Protocol

Complete guided workflow for assessing dynamic myocardial perfusion. The dynamic myocardial perfusion protocol is laid out in a series of three easy to follow steps. Start by selecting the series that has been generated from the non-rigid registration protocol and then follow the guided workflow:

- 1) **Cardiac Reformat** – images are automatically orientated into short and long axis views,
- 2) **Segmentation** – define the valve plane and the apex of the heart. The software will then calculate the endocardial and epicardial borders and provide a 17 segmented map of the myocardium along with a graphed view for each segment

- 3) **Functional Maps** – the third and final step produces the quantitative data with the color overlay and bullseye maps for blood flow, blood volume, mean slope of increase and mean transit time.

Quantitative data can be exported as a .csv file for a complete segmental and regional quantitative evaluation.

Functional Maps Available

The key perfusion parameters that CT Perfusion 4D package offers (measured in their appropriate physiological units) are:

1. **Regional Blood Volume (BV; ml/100g).**
Characterizes the amount of blood present in a local region; this is affected by the size of the blood vessels and the number of open capillaries.
2. **Regional Blood Flow (BF; ml/min/100 g).**
It is a measure of the volume of blood that flows through a local tissue region within a period of time.
3. **Regional Mean Transit Time (MTT; sec).**
Mean Transit time is the time that blood requires to traverse the tissue circulation from arterial inlet(s) to the venous outlet(s). Mathematically:
 $rMTT = rCBV / rCBF$ (Cerebral), or
 $rMTT = rTBV / rTBF$ (Tissue)
4. **Contrast Arrival Delay (IRF T0)**
It is the delay of arrival of contrast to any given location from its arrival at the artery ROI.
5. **Transit Time to IRF Peak (Tmax; sec)**
It is the transit time to the peak of the Impulse Residue Function (IRF).
Mathematically:
 $Tmax = IRF\ T0 + MTT/2$
6. **Capillary Permeability Surface Area Product (PS; ml/min/100 g)**
It is the unidirectional rate of transfer of blood solutes (contrast media) through the capillary endothelium into the interstitial space. PS is given in the same units as blood flow, ml/min/100g. Permeability Surface maps are intended to be used with lengthened acquisition times, typically 60 seconds are greater.
7. **Base Image** corresponds to an average of all pre-enhancement images. It is used in baseline subtraction and may also be as a reference image.
8. **Average Image** corresponds to the as the average of all images at a given slice location.


CONFIDENTIAL

SECRET DE
AFACERI

GE imagination at work

© 2021 General Electric Company.
All rights reserved. Data subject to change.
GE and GE Monogram are trademarks of General Electric Company.
* Trademark of General Electric Company
DOC1996217

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SECRET DE AFACERI

Visit us:

http://www3.gehealthcare.com/en/products/categories/advanced_visualization/applications/ct_perfusion_4d_multi-organs

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SECRET DE
AFACERI

GE imagination at work

© 2012 General Electric Company.
All rights reserved. Design subject to change.
GE and GE Monogram are trademarks of General Electric Company.
© 2012 General Electric Company.
* Trademark of General Electric Company. All rights reserved. Data subject to change.

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SECRET DE AFACERI

GE imagination at work

© 2017 General Electric Company.
All rights reserved. Date subject to change.
GE and GE Monogram are trademarks of General Electric Company.
* Trademark of General Electric Company
Doc1996217.

CardIQ Xpress 2.0

Bringing integration and automation to your CT Cardiac post processing needs.

Cardiac disease is one of the leading health concerns worldwide. Successful treatment of the many conditions that cause and perpetuate heart disease requires that physician's approach specific cardiac problems with as much information as today's technology can provide. To supply your referring physicians with this information you need software that gives insight into cardiac anatomy and tools for thorough evaluation into the extent of disease.

Overview

CardIQ Xpress 2.0 is an integrated post processing imaging analysis application dedicated to cardiovascular imaging on GE Advantage Workstation (AW) and AW Server Enterprise System (AWS). The CardIQ Xpress 2.0 software option can be used to display, reformat, and analyze 2D or 3D cardiac CT images for qualitative or quantitative assessment of heart anatomy and coronary artery vessels from a single or multiple cardiac phase image data sets. Cardiac motion is a very real challenge that can occur at any heart rate. CardIQ Xpress 2.0 is designed to work with SnapShot Freeze* images to automatically process and display images generated with reduced motion blur artifact.

Highlights

- Automatically segment coronary tree across phases.
- Automatically tracks and labels coronary arteries.
- Improved centerline editing tools for faster edits.
- Right mouse wheel menu for quick access to renaming, deleting, and editing centerlines.
- Plaque ID tool assists in visualizing and quantifying plaque burden
- Relative perfusion highlights and quantifies hypo-dense areas of myocardium

CONFIDENTIAL

SECRET DE
AFACERI

Visit us:

<https://www.gehealthcare.com/products/advanced-visualization/all-applications/cardiq-xpress-2.0-manual>

Features

- Pre-processing automatically recognizes cardiac datasets and performs all segmentations as data arrives on the system saving processing time
- Auto Coronary analysis automatically segments the coronary tree, tracks and labels the coronary arteries
- Three pre-defined orientations provide quick access to planes for best visualization of the coronaries in 2D
- Quick AVA allows access to vessel tracking at any time
- Relative perfusion color codes and quantifies percent of hypo-dense areas of myocardium with four selectable color maps and fusion overlays with the coronary tree
- Plaque ID provides volume measurements for four distinct Hounsfield ranges to aid with identification and visualization of coronary plaque
- Automatically display SnapShot Freeze¹ processed images for reduced motion blur
- Robust, automatic calculation of ejection fraction and stroke volumes from the 3D endocardium volumes
- Stenosis measurements, IVUS views and preset volume rendering models assist in communication of specific findings
- Measure ES and ED for ejection fraction & volume with automatic extraction of the left ventricle

- Create multiphase beating hearts
- Select oblique reformat views in the standard cath angles for easy analysis of the coronary vessels
- Display 4D valve views with a single click

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Rx Only

Indications for Use

CardIQ Xpress 2.0 is intended to provide an optimized non-invasive application to analyze cardiovascular anatomy and pathology and aid in determining treatment paths from a set of Computed Tomography (CT) Angiographic images. CardIQ Xpress 2.0 is a CT, non-invasive, image analysis software package, which aids in diagnosing of cardiovascular disease to include, coronary artery disease, functional parameters of the heart,

heart structures and follow-up for stent placement, bypasses and plaque imaging. CardIQ Xpress 2.0 offers unique tools such as automatic tracking, which will pre-process the CT data into multiple viewing ports to allow for an expedited read time improving workflow. With CardIQ Xpress 2.0, the user can color code the myocardial tissue to show hypo/hyper-dense areas in the myocardial tissue of the heart. With the IVUS-like view the user can color code the HU units of the plaque to better visualize the difference between calcified and non-calcified plaque in the wall of the vessel and the lumen to determine the amount of atherosclerosis. The user can see the different valve planes along with a variety of new layouts to align the heart. The IVUS-like view is created by applying GE's Volume Rendering on a cross-section perpendicular to the detected centerline. This view merely displays a cross section as in IVUS imaging and color codes like IVUS images. No new or additional diagnostic information is added. CardIQ Xpress 2.0 is for use on the Advantage Workstation (AW) platform, CT Scanner, PAC or Centricity stations, which can be used in the analysis of 2D or 3D CT angiography images/data derived from DICOM 3.0 CT scans.

CONFIDENTIAL

SECRET DE
AFACERI

GE imagination at work

© 2021 General Electric Company.
All rights reserved. Data subject to change.
GE and GE Monogram are trademarks of General Electric Company.
* Trademark of General Electric Company.
DOC1036477

Features Detail

Pre-Processing


- Automated processing of single or multiphase data
- Customizable menu to choose which protocols are pre-processed
- Pre-processed protocols include: Auto Coronary Analysis, Labeling of vessels, VR heart, Angiographic View, and Auto EF

Auto Launch

- Multiple exams can be loaded into the auto launch to ready them for review (AW only)
- Color-code display of ready-to-read exam listing (AW only)
- The ability to switch between exams without having to quit out of the application and reload


Automated Coronary Vessel Analysis (CVA)

- CVA allows the users to track, extract, visualize and measure coronary arteries from either single or multiple cardiac phase images
- Auto launch case selection and display of automatically tracked and labeled coronary vessels
- Simple, corrections to vessel branch tracking
- Curved, oblique, longitudinal, and cross-sectional reformatted views are automatically generated in the software for quick review

- Current-state tracking points within vessel analysis may be saved for future review and/or manipulation
- Color-coded plaque analysis with volume measurements for evaluating the change in plaque size over time
- Generate and save rotational movies from curved planar reformation, best section and cross-sectional (lumen) views

- Quantitative or qualitative coronary vessel assessment on user-selected vessel segments.
- Vessel measurements including:
 - Distance and volume
 - Cross-sectional area and mean diameter
 - Single or Dual reference point comparison
 - Relative percent stenosis
- Single-screen filming capabilities with multi-views within the screen to show the entire picture of the vessel
- Automatically display of the coronary vessel tree
- Using pre-set protocols, 2D or 3D coronary vessel tree models can be generated and displayed in an automatic fashion to give users a qualitative overview of coronary vessel structures

- New right mouse wheel menu for quick access to renaming, deleting and editing centerlines
- Improved centerline editing tools for faster edits
- The vessel tree models can be in the form of either 3D volume rendered or Maximum Intensity Projection (MIP)

2D Review

- Three pre-defined orientation protocols for easy review of the coronary vessels and chambers
- Dual reformat review allows automatic review of the coronaries from an axial image while linking to oblique views in longitudinal and cross-sectional planes
- Reformat review predefines workflow steps to automatically present thick, multi-planer views of cardiac anatomy
- One-click access to quick vessel analysis

Quick/Dynamic AVA

- Single or dual click vessel analysis from any protocol; 3D, reformat or vessel tree
- Real-time tracking of the center line with instant visualization of a vessel for quick inspection
- Ability to extend tracking proximally or distally for full view of the vessel
- Ability to select a point above and below potential lesions for automated vessel analysis

Vessel Data Base

- A vessel name data base for quick access to vessel labels
- Anatomically relevant listings
- Ability to add new vessel names into the data base
- Automated processing of single or multiphase data

Relative Perfusion

- Hypo-dense areas of myocardium can be high-lighted with color to aid in the visualization of ischemic heart disease
- Pre-defined layouts for density enhancement
- A hybrid display view to show the relationship of the vessel tree with the perfusion defect
- Quantification of a hypo-perfused area as it relates to myocardial defects
- Four selectable color maps to help display the hypo-dense areas
- Exportable statistics via one touch

IVUS Views

- Interactive volume rendered images to better define the calcium, lumen and non-calcified plaque related to lumen narrowing or wall abnormalities
- Can be applied to any 2D MIP image to include; cross-section, MPR, best L-section images

- Display cross-sections perpendicular to the vessel centerline to create an IVUS-like image. No new diagnostic information is added

PlaqID

- Customizable color mapping to HU ranges for easy identification of plaque in axial, curved reformat or MPVR images
- Four distinct colors to aid in differentiating vessel lumen, non-calcified plaque and calcified plaque
- Easy modification of the color and transparency of plaque
- Volume, area to track size of plaque over time
- Smooth transition of color from one density to another
- Plaque volume and area automatically captured in summary table.

VR Heart

- One-touch automatic segmentation algorithm designed to extract the cardiac anatomy from within the chest image
- Optimize algorithms for the easy review of heart and bypass grafts
- Multiple optimized VR curve settings to enhance different structures within the heart.
- Automatic segmentation of both single or multi-phase cardiac image data sets
- 4D beating heart images with the ability to page through the phases, rotate the image with real-time functionality

Multi-Phase Image Review

- Multiphase images can be reviewed in any protocol
- 10 phase images show the heart throughout the complete cycle
- Quickly edit phase images to keep only the phase or phases needed for analysis of the coronary vessels once multi-phase review is complete

GE imagination at work

CONFIDENTIAL

SECRET DE
AFACERI

Transparency View

- 3D vessel tree overlaid on the heart chambers with the ability to adjust opacity of chambers of heart, myocardium vessels tree and bones

3D Angiographic View

- One-touch automatic visualization of the coronary arteries in an x-ray angiographic view
- Single click cardiac cath views to orientate the angiographic view into RAO/LAO, or CRANIAL/CAUDAL orientation
- Added or removed vessels are easily archived with the Auto Select feature
- Toggle between gray scale or inverse video views

One Touch Cath Views

- Default cardiac cath views are provided with CardIQ Xpress 2.0 software to provide 3D VR heart, 3D vessel tree or reformatted image in the orientation of a standard cardiac catheterization
- Single click on the angulation area to change between one-touch views
- Interactive RAO/LAO & CRANIAL/CAUDAL orientations allow for manually positioning the view
- Customizable orientations can be easily generated and saved for future use

3D Ejection Fraction

- Automatic extraction of the left ventricle across all phases and the automatic detection of systole and diastole delivers calculated ejection fraction and stroke volumes
- Robust, automatic calculation of ejection fraction and stroke volumes from the 3D endocardium volumes

Layouts

- Configurable layouts
- Dual monitor support
- Configurable default layout for use at AVA launch
- Easily capture the screen layout for identical visualization on PACs or film

4D Movie Mode

- Movie mode allows the user to view and save multiple views of the heart and display as a movie sequence
- Load multiple phases from systole to diastole to create a beating heart movie
- Save movies as a DICOM image set or export in JPEG/MPEG format

SnapShot Freeze (SSF)

- Automatically display motion blur-reduced SnapShot Freeze processed images
- Reprocess SnapShot Freeze images for artifacts caused by motion blurring after manual vessel editing
- Motion blurring correction requires the use images produced by a CT using the optional SnapShot Freeze feature

Protocols

- CardIQ Xpress 2.0 is supplied with a set of pre-defined protocols that are easily adaptable to customized protocols
- The protocols have comprehensive instructions and tools which minimize the need to memorize the procedures or to refer continuously to user documentation.

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5818595-1EN
Revision 7

Clinical Benefits

CardIQ Xpress 2.0 is a software application that provides anatomical and functional information from a gated cardiac CT dataset. It provides automatic segmentation of the heart and coronary arteries along with labeling, centerline tracking and tools for manual editing to assist physicians in visualization of the cardiac anatomy and disease assessment.

Performance Characteristics

- Automatically segments coronary tree across phases
- Automatically tracks and labels coronary arteries
- Guided workflow including pre-defined orientations
- Unrestricted access to vessel tracking
- Plaque ID tool assists in visualizing and quantifying plaque burden
- Relative perfusion highlights and quantifies relative hypo-perfused areas of myocardium
- Automatic segmentation of Left Ventricle across multi-phase data presenting volume and ejection fraction quantitative results
- Guided workflow with quantitative measurement tools and predefined visual displays
- Automatic segmentation of cardiac structures along with Volume Viewer cine capability

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
SECRET DE
AFACERI

CardIQ Function Xpress

Real time cardiac review the instant you're ready to read.

Cardiac disease is one of the leading health concerns world wide. Successful treatment of the many conditions that cause and perpetuate heart disease requires that physicians approach specific cardiac problems with as much information as today's technology can provide. To supply your referring physicians with this information you need software that gives you insight into cardiac anatomy and function supported by quantifiable data.

Overview


Cardiac IQ Function Xpress post-processing software helps you evaluate cardiac function and diagnose cardiovascular disease with a high degree of confidence. Providing accurate and reproducible quantification of left and right ventricular volumes, ejection fractions, and myocardial mass, Cardiac IQ Function Xpress is optimized to assess cardiac function using multiphase, multi-slice cardiac CT images. The CardIQ Function Xpress option provides an easy-to-use and time-effective way for you to perform cardiovascular functional analysis.

CONFIDENTIAL

SECRET DE AFACERI

What's new

- Automatically pre-processes and loads exams.
- Automatically detects all chambers in all phases.
- Automatically performs ejection fraction and chamber volume analysis with 91% reliability on LV and RV automatic ejection fraction.
- Provides single click myocardial analysis with bull's eye images.
- Performs simultaneous case review in real time.
- Lets you load multiphase data with no load time.
- Lets you review multiple cases in real time.
- Lets you generate reports with imported analytics.

Visit us:

www.gehealthcare.com/aw/applications/cardiac-function-xpress/

Features

- Apply the CardIQ Function Xpress application to standard axial CT images acquired on a GE scanner using the cardiac CT SnapShot imaging acquisition option.
- Automatically selects each heart chamber for individual volume analysis.
- Performs behind-the-scene function data process for real time ejection fraction review, volume analysis, and myocardial analysis.
- Enables you to extract, render, and display 3D volumetric models of the endocardium for ejection fraction calculation.
- Automatic detection of the endocardial and epicardial walls for wall motion, wall thickness, wall thickening and myocardial mass analysis. Lets you perform volume analysis of all chambers.
- You can automatically calculate left atrium volume while excluding the pulmonary vein.
- With a single click you can activate visual wall motion with short axis images in basal, mid, and distal orientations along with a two-chamber long axis view.
- You can perform myocardial analysis with wall motion, wall thickness, and mass calculations.
- The application's flexible reporting tool lets you include graphical representations.
- Display on screen a table of key functional parameters for instant visualization.

Intended Use

CardIQ Function Xpress is intended to provide an optimized non-invasive application to analyze cardiovascular anatomy and pathology and aid in determining treatment paths from a set of Computed Tomography (CT) Angiographic images.

Image Requirements

CardIQ Function Xpress accepts standard CT image sets acquired on qualified GE CT scanners using the appropriate cardiac imaging software. These images must meet the same image requirements as those for the basic Volume Share 2 application.

System Requirements

- AW Volume Share Workstation with Volume Share 2 running 4 8400 hardware or higher (running 64 bit).
- Auto launch and preprocessing are available only on 8400 workstations with 16 GB RAM.
- 2 monitor configuration.
- Color landscape monitor.

Recommended Printers

- Codonics NP-1660M
- Kodak 3600 DMI
- Codonics 1660M, 1660MD, or Horizon
- Lexmark Optra 1650N, 1855N, SC1275N, C710N, C72N, T612, or T614
- Seiko 1720D
- Quantum GL2101HD with film/thick paper
- Quantum GL2101HD with plain paper (see PI-008)
- Tally T8106
- HP LaserJet
- Xerox Phaser

Regulatory Compliance

This product complies with the European CE Marking regulation for Medical Devices Directive: Directive 93/42/EEC, dated 14 June 1993.

CONFIDENTIAL

SECRET DE
AFACERI


GE imagination at work

© 2012 General Electric Company.
All rights reserved. Data subject to change.
GE and GE Monogram are trademarks of General Electric Company.


Overview

CardIQ Function Xpress allows you to accurately and reproducibly quantify myocardial mass, and left and right ventricular volumes and ejection fractions. The software automatically detects left ventricular endocardial and epicardial contours so that you can more easily assess left ventricular (LV) and right ventricular (RV) functional parameters. The program is optimized to perform assessment of cardiac function using GE multiphase, multi-slice cardiac CT images.

The CardIQ Function Xpress menu page consists of multiple image processing tools for single or multiphase cardiac studies. These protocols can be accessed either from the AW Volume Viewer application or directly from the application menu page on the browser.

CardIQ Function Xpress benefits to cardiac patients include:

- Minimally invasive alternative to help diagnose heart function.
- Data is acquired during a coronary CTA exam, so no additional dose or contrast is required.
- Offers a diagnostic option to patient who may not be able to undergo other non-invasive exams.

Feature Details

CardIQ Function Xpress uses single or multiphase image data sets to measure, graphically plot, and tabulate a number of ventricular function and wall motion parameters.

Pre-processing

- The program automatically processes multiphase data for functional analysis.
- You can select configurations to process only the chambers you require.
- For more streamlined reading, ejection fraction, volume analysis, and myocardial analysis are completed before you view the case.

Auto Launch


- Load multiple cases into auto launch so they are ready to review when you are.
- Switch between exams without having to quit the application and reload.
- Exams ready for review display in green.

Automatic Chamber Detection

- The program automatically loads single or multiphase data sets and selects the heart's chambers for volume analysis.
- Analyze individual chambers by automatically contouring any of the heart's four chambers.
- The system automatically detects the left and right ventricles and the left atrium and contours the endocardium for ejection fraction calculation.
- Choose single or multi-phase chamber detection.

Automatic Ejection Fraction

- The system features a protocol that, when chosen, automatically detects the end diastolic and end systolic phases and calculates an accurate ejection fraction of the left ventricle, right ventricle, and left atrium.

- Ejection fraction measurements are available for


the right and left ventricles.

- Stroke volume, end diastolic and end systolic volumes are automatically calculated and shown in a graphical representation for assessment of heart failure.
- You can easily send ejection fractions and other quantification results to your report.
- With dual screens, you can display quantification assessments on the left monitor and a qualitative overview of the full chamber from apex to base on the right monitor.
- Customize layout designs to accommodate your viewing preferences:
 - Base to apex
 - Basal, mid, and distal along with a two-chamber long axis.
- You can change the selection of end diastolic and end systolic phase.
- While contour detection is automatic, scalpel and auto select tools are available with which you can edit the contours.

Display Table

- A display table shows volume and provides a graphical update of end diastolic and end systolic locations.

- Click on available tabs to select different chambers for analysis.
- Left and right ventricle measurements include:
 - Stroke volume
 - Ejection fraction
 - Cardiac output
 - Myocardial mass
 - End systolic volume
 - End diastolic volume
 - Pulmonary vascular resistance
 - Systemic vascular resistance
 - Body surface area
 - Heart rate

- Left atrium measurements include:
 - Min and max volume
 - Ejection fraction
 - Left atrium index
 - Body surface
 - Heart rate
- Right atrium measurements include:
 - Min and max volume
 - Heart rate
- You can film, save, and export the volume graphs.

Myocardial Analysis

- CardIQ Function Xpress automatically detects the endocardial and epicardial walls to analyze wall motion, wall thickness, wall thickening, and myocardial mass analysis.
- You can edit the contours with a freehand trace and click to spline tools.
- The program displays multiple bull's eye graphs to graphically represent these cardiac characteristics:
 - Wall motion
 - Wall thickness
 - Myocardial mass
 - Cardiac output
- Choose from four bull's eye graph overlays to show 17, 20, and 49 segment model, or a coronary artery territory model.
- Manual scoring is available for wall motion and thickening in 17 and 20 segment.
- Cord display is available for visual assessment of wall abnormalities.

Reporting

- The application provides you with automatic report generation capabilities.
- You can import hospital or clinic logo and address to make your reports site specific.
- You can capture 3D volume rendering endocardium images with ejection fraction calculations displayed for inclusion in your report.
- You can export the volume graph of all heart chambers to your report.
- The program automatically fills these report fields:
 - Ejection fraction
 - Stroke volume
 - End diastolic volume
 - End systolic volume
 - Pulmonary vascular resistance
 - Systemic vascular resistance
 - Myocardial mass
 - Cardiac output
- You can create predefined conclusion templates for various types of findings.
 - Report templates are fully customizable, so you

can include logos, diagrams, signatures, patient history, and messages.

- Generate reports in PDF, HTML, or plain text.
- You can combine CardIQ Function Xpress reports with CardIQ Xpress 2.0 Elite reports.

Simplified Workflow

CardIQ Function Xpress simplifies your workflow with:

- Automatic ejection fraction protocol
- Streamlined multiphase beating images for easy visualization of cardiac wall motion.
- Single click myocardial analysis with bull's eye images.
- Simplified reporting.

Summary

CardIQ Function Xpress provides you with the tools you need to accurately analyze, evaluate, and report on cardiac anatomy and function. Its ease of use will streamline your workflow and its accuracy will give you greater diagnostic confidence.

Indication for Use

CardIQ Function Xpress is used in conjunction with CT cardiac images to automatically calculate and display various left ventricular and right ventricular functional parameters such as ejection fraction, end systolic and end diastolic volumes, stroke volumes, wall motion, wall thickening, cardiac output, myocardial mass, and systemic and pulmonary vascular resistance. Volume measurement of each chamber of the heart is also available. With CardIQ Function Xpress atrium volumes may be used to determine volume assessment of atrial disease to including, but not limited to, atrial fibrillation. CardIQ Function Xpress is a CT, non-invasive image analysis software package, which aids in the assessment of cardiac function and in determination of cardiovascular disease diagnosis and management. CardIQ Function Xpress is available on the AW Workstation, CT Scanner, PAC, or Centricity stations, which can be used in the analysis of 2D or 3D CT angiography images/data derived from DICOM 3.0 CT scans.

CardIQ Function Xpress is an extension of Volume Viewer 4 and is dedicated to cardiovascular imaging. It provides an optimized non-invasive application to analyze cardiovascular anatomy and pathology and aid in determining treatment paths from a set of CT Angiographic images.


CONFIDENTIAL

SECRET DE
AFACERI

GE imagination at work

© 2012 General Electric Company.
All rights reserved. Data subject to change.
GE and GE Monogram are trademarks of General Electric Company.
* Trademark of General Electric Company.

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
*Marcă comercială a General Electric Company

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
SmartScore 4.0

Advanced imaging software that detects, quantifies and scores cardiac calcium plaque burden. Instantly.

Cardiovascular disease remains one of the most common health issues in the world today. As with many conditions, early detection and patient risk assessment are vital to preventing or minimizing long-term negative effects. But many cases of this disease aren't diagnosed until the patient presents with symptoms. Conventional procedures to assess risk can be expensive, time consuming, and uncomfortable. An accessible, patient-friendly risk assessment method could help clinicians devise a regimen for their patients that might lessen the chance for serious cardiac events.

What's new

- Non-invasive alternative to conventional assessment procedures.
- Score can be correlated to age group cohort to determine patient's risk per population
- Provides information on coronary artery wall calcium plaque buildup.
- Automatically detects calcium and highlights it in green.
- Streamlines workflow by networking patient demographics from CT scanner directly to SmartScore program.
- Free Hand Trace lets you outline specific ROIs.

Overview

SmartScore 4.0 is designed to identify the presence of regional and global coronary artery calcification from a CT scan, then measure and score the results. Scores can be calculated using a standard Agatston/Janowitz (AJ) method. When correlated with a patient's personal information, the score can yield an estimation of a patient's risk for coronary artery disease.

CONFIDENTIAL

SECRET DE
AFACERI

Visit us:

www.gehealthcare.com/aw/applications/smartscore-4/

Features

- AJ 130 scoring method uses conventional Agatston/Janowitz technique with a threshold of 130 HU which is adjusted to the appropriate image slice thickness.
- Volume scoring calculates volumes (mm^3) of calcified plaque above the 130 HU threshold.
- Mass scoring calculates mass (mg) of calcified plaque above the 130 HU threshold.
- Individual and aggregate scores are computed for each artery type.
- Provides two methods of calcium scoring.
- Customizable parameters can be predefined based on your site's preferences.
- User interface streamlines your workflow and reduces reading time.
- Report options allow you to custom tailor your reports and distribute them in a variety of formats.

System Requirements

- EKG monitor with recording device and x-ray translucent lead cable.
- Gantry hardware upgrade kit for those scanner systems already in operation.

Minimum platform release:

- Software for Advantage Workstation 4.2P or higher.
- AW Server 2.

Recommended Options

Postscript Printers:

- Codonics: NP-1660M
- Kodak: 3600 DMI
- Codonics 1660M, 1660MD or Horizon
- Lexmark Optra 1650N, 1855N, SC1275N, C710N, C72N, T612, or T614
- Seiko 1720D
- Kodak Dmi3600
- Quantum GL2101HD with film/thick paper
- Quantum GL2101HD with plain paper (see PI-008)
- Tally T8106
- HP LaserJet
- Xerox Phaser

Intended Use

SmartScore is a non-invasive software option that can be used to evaluate calcified plaques in the coronary arteries, which may be a risk factor for coronary artery disease. SmartScore may be used to monitor the progression/regression of calcium in coronary arteries over time, which may aid in the prognosis of cardiac disease.

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SECRET DE
AFACERI

GE imagination at work

© 2012 General Electric Company.
All rights reserved. Data subject to change.

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Vizitați-ne:

www.gehealthcare.com/aw/application/gsi-viewer/

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
RADIATOR AUTORIZAT ENCLEZĂ-FRANCEZ
FODOR
ALEXANDRA
IOANA
2810318125819

Bone VCAR

Automated spine labeling with curved reformat generation to improve ease of reading and reporting

Background


Approximately six percent of the total CT exam procedure volume is dedicated spine imaging. Add to that the number of other routine procedures that include the spine and are referenced within an imaging report you will find a significant number of imaging series generated every day for every scanner where identification and review of the vertebral bodies and their discs are necessary. Bone VCAR will help simplify the reading experience and improve reporting efficiency by identifying and displaying the anatomy in an automated yet flexible implementation.

Overview

Bone VCAR provides an automated spine labeling application based on a deep learning model. It can identify and label segments of the spine or the entire spine in a matter of seconds. Additionally, let Bone VCAR generate the curved reformats of the spine along with the perpendicular oblique views to easily see the true cross section of the vertebral bodies and disc spaces. Utilize these time saving features to assist in your reading experience. Bone VCAR can also improve the reporting efficiency for identifying and remembering key landmarks that are included in the report dictation.

Highlights

- Automated spine labeling
- Automated curved reformats with oblique views parallel to the bodies and disc space
- Accessible for any exam type; trauma, oncology, dedicated spine, general imaging
- Accessible with all workflows

CONFIDENTIAL

SECRET DE
AFACERI

Features

- Automated vertebrae labeling in 5 seconds or less for volumes of 300mm or less *
- Deep learning application trained on global datasets acquired with a broad range of acquisition parameters achieving >90% labeling accuracy
- Automated generation of a 3D trace to generate curved reformat views
- Automated generation of oblique views perpendicular to vertebral bodies and disc spaces
- Compatible with axial source images as well as previously generated reformatted series
- Manually insert, edit or delete vertebral body labels as necessary
- Compatible with data acquisitions for dedicated spine work, trauma, oncology and routine imaging
- Load multiple series and vertebrae labeling will propagate to the different volumes
- Easily review multiple reconstructions with simultaneous display of different kernels
- Access from the Volume Viewer tool bar to easily apply during your standard read at any time
- Compatible with all DICOM CT data that meets image input requirements as described in the user guide

System Requirements

Minimum platform release:

- AW Workstations AW 4.7 Ext. 12 with VV14.0 Ext 6 or higher **Timing specifications with 2440 platform configuration and normal anatomy*
- AW Server 3.2 or higher
- Color Monitors
- Single or Dual Display Monitors for AW Server in landscape or portrait orientations

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Regulatory Compliance

This product complies with the European Council Directive 93/42/EEC Medical Device Directive as amended by European Council Directive 2007/47/EC. This product or its features may not be available in some other countries or regions. Please contact your sales associate.

Rx Only

Indications for Use

Bone VCAR is a post processing application for use in the analysis of CT images. The software is intended to support clinicians in the review of images that include the spine by providing tools to label the spine and optimize the display of anatomy within the CT image.

Bone VCAR is designed to support the clinician in visualization of the spine, by providing initial identification of vertebrae to assist in report dictation.

The software also assists the user by providing optimized display settings for easier identification of anatomy to facilitate fast image review and reporting of findings. Bone VCAR may be used for multiple care areas and is not specific to any disease state. It can be utilized during the review of various types of exams including trauma, oncology, and routine body.


CONFIDENTIAL

**SECRET DE
AFACERI**

GE imagination at work

© 2019 General Electric Company.
All rights reserved. Data subject to change.
GE and GE Monogram are trademarks of General Electric Company.
* Trademark of General Electric Company.
DOC2185114

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Colon VCAR EC

A comprehensive reading workflow solution for colonic lesion detection.


Colorectal cancer is a major health concern for both men and women. Recent studies conclude that early detection of colonic polyps can reduce the mortality rate of colon cancer. You need a protocol that lets you quickly and accurately detect and report results of CT colon scans; a program that optimizes your workflow and improves your reading productivity.

Overview

Colon VCAR EC gives you a comprehensive reading workflow solution for detecting colonic lesions with flexibility and efficient performance. With it you can visualize anatomy that would otherwise be obscured by tagged fluid. Plus digital contrast agent filter highlights potentially suspicious regions.

What's new

- The program allows you to read and problem solve in 2D, 3D, or 360 degree dissection views.
- Prone/Supine polyp linking provides quick navigation between datasets.
- Digital Contrast Agent (DCA) provides automatic visualization of shapes characteristic of polyps.
- Electronic Cleansing (EC) subtracts tagged stool and fluid, making it easier to identify lesions.

CONFIDENTIAL

SECRET DE
AFACERI


Visit us:

www.gehealthcare.com/aw/applications/


Features

- Prone and supine synchronized reviewing for problem solving.
- You synchronize prone and supine navigation based on polyp bookmark location and distance from center line.
- Provides you with quick, easy to use polyp and center line measurement tools for complete reporting.
- You can leverage any combination of synchronized 2D, 3D, and 360 views.
- Dissection views make primary reading and problem solving more efficient.

- You can vary navigation speed.
- Lets you track entire colon from rectum to cecum for thorough, precise exams.
- You are provided with unrestricted viewing angles of all colon surfaces.
- You can apply color to bookmarked lesions for better visualization.
- The program allows you to assess luminal, mural and extracolonic anatomy and pelvic structures.

- A selection of application-specific tools ensures ease of use and thorough exams.
- Your workflow benefits from fast, interactive electronic patient reporting.

System Requirements

- AW Workstation
- AW Server

Indications for Use

Colon VCAR EC is a CT, non-invasive, image analysis software package that allows the visualization of 2D, 3D and dissected medical image data of the colon derived from DICOM 3.0 compliant CT scans. ColonVCAR is designed to aid the physician in evaluating the lumen and internal wall of the colon to confirm the presence or absence of colonic lesions (e.g. polyps). It provides functionality for 2D/3D rendering, bookmarking of suspected lesions, synchronized viewing of the 2D, 3D and 360 dissection views for data sets acquired in any position, and an object oriented endoluminal display. In comparison to Colonoscopy, this tool has an advantage of depth penetration due to its 3D presentation capability.

It is intended for use by clinicians to process, render, review, archive, print and distribute colon image studies.

The Colon VCAR DCA (Digital Contrast Agent) module is an automated highlight feature for the visual identification of spherical structures within the colon and is intended to be used as concurrent reading device. Digital Contrast Agent is a 3D filter that produces images that highlight spherical anatomical regions, such as polyps, and/or stool. Colon VCAR uses color to display these highlighted spheres. Images are made available to the physician to aid in characterization of potential polyps and thus, the patient management care decision process.

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This product has not been approved for sale in the United States by the United States Food and Drug Administration (FDA).


CONFIDENTIAL

SECRET DE
AFACERI

GE imagination at work

© 2015 General Electric Company.
All rights reserved. Data subject to change.
GE and GE Monogram are trademarks of General Electric Company.

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
Lung VCAR

For lung nodule analysis - automatic visualization, measurement, reporting and follow-up

Clinical and Technical Background

CT scanning is utilized as the exam of choice to visualize and analyze complex lung pathology. The detection of pulmonary nodules and assessment of their evolution with CT are of major importance in chest imaging.

CONFIDENTIAL


SECRET DE
AFACERI

Overview

Lung VCAR brings efficient CT pulmonary nodule assessment and diagnosis. The innovative Digital Contrast Agent (DCA) feature automatically visualizes lung nodules to help you confirm the presence or absence of suspicious lesions from 2 to 12 mm in size. Lung VCAR allows automated follow-up for lesion matching by the registration of two or more datasets, automatic lesion classification.

Highlights

- Synchronized 2D, Digital Contrast Agent (DCA) and segmentation analysis.
- Automatic nodule visualization.
- Automatic nodule analysis (volume, doubling time, % growth).
- Automatic follow-up
- Reporting workflow

Visit us:

<https://www.gehealthcare.com/products/advanced-visualization/all->

Features

Review

- Ability to synchronize multiple images for nodule comparison.
- Ability to review single or multiple exams and compare axial, sagittal, oblique, coronal, and volume-rendered images.
- Automatically propagates previous exam bookmarks to current exam.
- Automatically segments both right and left lung to reduce visual distraction.
- Digital Contrast Agent (DCA) automatically highlights spherical shapes to enhance visualization of suspicious nodules.

Analysis

- Performs automatic segmentation of all nodule types.
- Provides automatic nodule analysis, including % growth, doubling time
- Automatic Nodule Contour to verify pixels within the volume

Indications for Use

Lung VCAR/AdvantageALA is intended to provide an optimized non-invasive application to measure abnormalities in the lung (for example, nodules, lesions, etc.) from a set of computed tomography (CT) images.

The software is designed to support the physician in confirming the presence or absence of physician identified lung lesions (e.g. nodules). The software allows measurement of volume over time using a consistent standardized measurement protocol, thus providing an estimation of the volume doubling time. The Lung VCAR/AdvantageALA software allows analysis and displays statistics for nodule characterization all the different nodule types. Lung VCAR optional Digital Contrast Agent (DCA) module is an automated highlight feature for the visual identification of possible lesions. Digital Contrast Agent (DCA) is a 3D filter that produces images that highlight spherical (S) or cylindrical (C) anatomical regions, such as nodules, cysts, scars, and vessels. Images are made available to the physician to aid in characterization of suspicious nodules and thus, the patient management care decision process. Lung VCAR/AdvantageALA provides the physician with additional information, meant to complement diagnosis based on classical techniques.

System Requirements

Lung VCAR/AdvantageALA can be installed on AW Workstations AW 4.7 Ext. 14 or later and AW Server 3.2 Ext 3.2 or higher.

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SECRET DE
AFACERI

GE imagination at work

© 2021 General Electric Company
All rights reserved. Data subject to change
GE and GE Monogram are trademarks of General Electric Company
* Trademark of General Electric Company
DOC0751472

Features Detail

Scanning

Lung VCAR sensitivity and segmentation performance is a function of slice thickness and dose. Lung VCAR was designed to be compatible with slice thicknesses ranging from 0.625 –1.25 mm with a dose range of >40.

Review Correlated Workflow

The Lung VCAR user interface is designed to organize routine chest radiology reading and make it more efficient. Viewports and displayed anatomy are synchronized and correlated, facilitating reading by making review and problem solving more transparent.


Lung VCAR gives you the ability to:

- Synchronize multiple images for nodule comparison.
- Review single or multiple exams and compare axial, sagittal, coronal, oblique, and volume rendered images.
- Automatically propagate and display previous exam bookmarks on the current exam when two exams are loaded.

Lung VCAR provides you with tools to help improve your workflow and analysis precision.


Lung Segmentation

Lung VCAR automatically segments both the right and left lungs. Segmentation reduces visual distractions associated with irrelevant anatomy during lung nodule evaluation. Total lung volume is captured in the patient report.


Automatic Follow-up synchronizes and displays current and prior patient exams and their corresponding bookmarks for quick and easy comparison review.

Digital Contrast Agent (DCA) automatically highlights spherical shapes to enhance visualization of suspicious nodules.

DCA Color Mapping highlights nodules in your choice of red, or for a colorblind-friendly display, orange.

DCA Object Size (0-10mm) gives you full control over visualization sensitivity to a desired nodule size and the number of false positives displayed.

DCA Shape Choice (Sphere / Cylindrical) lets you select from two different filters to help visualize specific shapes and assist in problem solving:


- Cylindrical for vessels
- Spheres for nodules

DCA Toggle enables you to toggle the DCA-highlighted nodules on and off.

X-Reference / Correlation Bar provides you with a quick reference for localizing a nodule's global location when reading 2D axial images. The x-reference bar is synchronized with both review and analysis layouts for immediate screen refresh when deposited at any location or on any anatomy.

CONFIDENTIAL

SECRET DE
AFACERI

Bookmark Management

Lung VCAR provides you with custom-designed (select dot or small box bookmark format) bookmark management capability linked to review, paging and Digital Contrast Agent for a more organized nodule tracking and analysis workflow.

Interface lets you dynamically and easily view, page and deposit bookmarks within an exam. It provides:

- On/Off toggle capability for DCA-highlighted nodules
- Smooth, contiguous paging through large data sets
- Thick slab MIP paging
- Bookmark paging (next / prior and add / delete)

Automatic Bookmark Propagation lets you automatically propagate bookmarks from previous to current exam or current to previous exam.

Deferred bookmarking lets you queue bookmarked nodules for later segmentation (operator initiated segmentation).


Immediate bookmarking automatically segments the nodule and displays the analysis review layouts immediately after you place a bookmark

Layouts

Layout Preferences let you create and save a layout to match your reading style. Layout elements such as viewport, image type in each viewport, slice thickness, window level/width, magnification, rendering mode, can be easily customized.


Analysis

Analysis mode allows you to maximize use of a second monitor by displaying a segmented nodule of interest to an orientation useful for quick problem solving and analysis.

Automatic segmentation of all nodule types provides automated segmentation and sizing of all nodules types. This feature can:

- Display 3D volume rendered nodule view.
- On follow-up exams, display doubling time and percentage of nodule growth.
- Compare nodules from multiple exams

Shutter view provides you with a quick Volume Rendered framed view of a point of interest in relationship to surrounding structures for fast visualization and problem solving.

CONFIDENTIAL

SECRET DE
AFACERI


Manual or automatic nodule characterization can automatically characterize and label nodule type and apply the necessary segmentation algorithm for accurate isolation. You can at any time manually override this feature and apply a different segmentation algorithm and labeling for any nodule.

Automatic Nodule Analysis provides the following nodule information:

- % growth
- Doubling time
- Volumes (solid component, non-solid component, total volume)
- Maximum distance in 3D (long axis)
- 2D measurements
- Hounsfield units (Min, Max, Avg)
- Scan interval

Paging enhancements provide a special scrolling tool that lets you page through large data sets while the image displays seamlessly.

Automatic nodule contour capability lets you verify pixels within the volume. It displays segmentation nodule contour on the 2D image for quick segmentation validation.

Patient Report

DICOM SR allows a reporting workflow to present and archive images and measurements.

Report feature only available on AW Workstation.

Summary

With automated nodule visualization and localization, increased detection sensitivity, accurate and repeatable segmentations analysis, Lung VCAR could streamline your reading workflow, enhance your efficiency, and increase your diagnostic confidence.


CONFIDENTIAL

SECRET DE
AFACERI

GE imagination at work

© 2021 General Electric Company
All rights reserved. Data subject to change
GE and GE Monogram are trademarks of General Electric Company
* Trademark of General Electric Company
DOC0751472

Stroke VCAR

Part of a comprehensive acute Neuro workup—Semi-automated segmentation and assessment of Hematoma and Aneurysms in CT images.

Clinical and Technical Background

Stroke and trauma to the head may result in intra cerebral or subdural hematomas. Patients under suspicion of such internal bleeding are triaged using CT imaging. It is important to determine both the presence as well as the amount of such bleeding to properly triage a patient for treatment or follow-up. Bleeding in the brain may also be caused by the rupture of aneurysms. It is important to detect the presence of the aneurysm and assess its potential hazard.

Overview

Stroke VCAR (Volume Computer Assisted Reading) gives you a complete reading workflow solution for a comprehensive and robust analysis of hematoma and aneurysms. Hematoma assessment is done using semi-automated segmentation tools married with innovative interactive editing capability in the form of SmartMesh. Aneurysm assessment is done through an innovative user guided aneurysm segmentation and visualization. The program lets you generate a clear, concise clinical summary for sharing with referring physicians.

Highlights

- The program allows you to quickly segment the hematomas in a non-contrast CT.
- Introducing SmartMesh: an intuitive volumetric editing tool.
- Intelligent user guided aneurysm segmentation.
- Fully integrated with Spectral CT when used with GSI datasets.

CONFIDENTIAL

SECRET DE
AFACERI

Visit us:

www.gehealthcare.com/en/products/ct/ct-scans/advanced_visualization/computed_tomography_imaging_software

Features

- Semi-automated segmentation and sizing of hematoma.
- Introducing SmartMesh—an interactive volumetric editing tool.
- Comprehensive reporting of segmented hematomas.
- Fully integrated with Spectral CT when using GSI data.
- Track hematoma changes both visually and quantitatively with longitudinal exam workflow.
- Simple user driven aneurysm detection workflow.
- Easy and intuitive workflow for aneurysm segmentation. Initial segmentation done with four guided clicks.
- A selection of application specific tools ensures ease of use.
- Create layouts to match your reading style.

System Requirements

Minimum platform release:

- AW VolumeShare 7 and AW Server 3.2 or later.
- Color monitors
- Single or Dual Display monitors for AW Server in landscape or portrait orientations.

Regulatory Compliance

This product complies with the European CE marking regulation following Medical Devices Directive: Directive 93/42/EEC.

Indications for Use

Stroke VCAR is a CT image analysis software package that allows for the analysis and visualization of Brain CT data derived from DICOM 3.0 compliant CT scans. Stroke VCAR is designed for the purpose of segmenting and assessing hemorrhages in the brain using automated tools on non contrast CT exams. Additionally Stroke VCAR provides a set of workflow tools for the segmentation and visualization of aneurysms in the brain from contrast enhanced CT exams. It is intended for use by clinicians to process, review, archive, print and distribute CT studies.

This software will assist the user by providing initial 3D segmentation, measurements and visualization of hemorrhages and aneurysm in the brain. The user has the ability to adjust, review and has to confirm the final segmentation.

CONFIDENTIAL

SECRET DE
AFACERI

GE imagination at work


© 2016 General Electric Company.
All rights reserved. Data subject to change.
GE and GE Monogram are trademarks of General Electric Company.
* Trademark of General Electric Company.

Features Detail

Hematoma


Segmentation

The Stroke VCAR hematoma tool uses an intelligent user guided algorithm for semi-automated segmentation of all types of hematomas in the brain. These segmentations can be converted into volume of interest to measure and assess their density.

SmartMesh editing

In cases when the semi-automated segmentation does not fully meet a user's expectation, SmartMesh allows for innovative editing capability to quickly complete the segmentation. A novel feature of the editing tool is to provide users with immediate visual feedback on what is included or excluded from the edition.

CONFIDENTIAL

SECRET DE
AFACERI

Reporting

Simple and intuitive reporting enables users to confidently manage and share hematoma findings in a consistent and standardized manner.

Follow up

Load repeat scans in Stroke VCAR Hematoma protocol and use its longitudinal analysis capability to track changes, both visually and quantitatively, due to progression or therapy response.

Jan 10 2014 11:11:00		Jan 12 2014 14:51:00	
CT - ARNDT/MZ/D		CT - ARNDT/MZ/D	
FACIL 1 - Hematoma		FACIL 2 - Hematoma	
Vol	94.91 cm3	Vol	90.40 cm3 (+5%)
Short Axis	30.0 mm	Short Axis	30.0 mm (+0%)
Long Axis	55.7 mm	Long Axis	54.4 mm (+1%)
Cumulative Volume for selected volumes 91.81 cm3		Cumulative Volume for selected volumes 90.40 cm3 (+5%)	

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
SECRET DE
AFACERI

Aneurysm


Detection Workflow

A simple and intuitive workflow that helps users with a fast and easy assessment for the detection and treatment planning of aneurysms.


Segmentation

A simple guided workflow for the segmentation of an aneurysm. Fast initial segmentation is done with four guided clicks.

Define visualization planes to better see how the aneurysm is connected to the feeding vessel. Once the segmentation is complete, the user is able to interactively review 3D images of the segmented aneurysm and measurements.

Reporting

Simple and intuitive reporting enables users to share images and measurements of aneurysms in a consistent and standardized manner.


Summary

Stroke VCAR provides workflow benefits using semi-automated hematoma and aneurysm segmentation tools enabling a standardized approach to your facilities review & analysis of these exams.

GE imagination at work

© 2016 General Electric Company.
All rights reserved. Data subject to change.
GE and GE Monogram are trademarks of General Electric Company.
* Trademark of General Electric Company.

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
CardEP

Streamlined post-processing for enhanced electrophysiology procedures.


CT cardiac exams provide critical information to practitioners of cardiac therapy. Precise information pertaining to the left atrium's complex anatomy, the pulmonary veins, the coronary sinus, or the cardiac veins has a major impact on the efficacy of subsequent cardiac therapy: It can speed procedures and facilitate treatment. But analyzing and reporting the wealth of information CT cardiac studies provides can be time consuming. A program that automates many routine functions and gives you tools to easily quantify and qualify many aspects of cardiac function would streamline your workflow and give you greater diagnostic confidence.

Overview

CardEP is an integrated post processing image analysis software for the application of cardiovascular and electrophysiology imaging on the AW Workstation (AW) and AW Server platforms. With it, you can process display, reformat, and analyze 2D or 3D cardiac images for qualitative or quantitative assessment of heart anatomy and pulmonary veins from single or multi-phase cardiac image data sets.

Highlights

- Visualize the anatomical layout of the left atrium and pulmonary vessels.
- Visualize the origins of the pulmonary vessels.
- Automatic heart segmentation.
- Helps you visualize and quantify cardiac venous pathways.
- Offers you a variety of 2D, 3D, or reformatted protocols with which to perform image analysis.

CONFIDENTIAL

SECRET DE
AFACERI

Features

CardEP offers a variety of features that give you image analysis flexibility:

- Volume render 3D heart and chamber views
- Obtain segmented 3D views of the left atrium and pulmonary veins
- View hollow 3D VR heart models
- Measure the pulmonary veins and left atrial appendage with predefined protocols
- Render and display 2D/3D views of the left chambers in single or multiple phases of the cardiac cycle with one-touch automatic extraction
- Reformat standard axial CT images of single or multiple cardiac phases into short or long axis and save the image series for further analysis.
- Automatically track, extract, and display pulmonary veins with the pulmonary vessel analysis tool
- Generate pulmonary vein measurements
- Automatically track, extract, and display coronary sinus and cardiac veins with the coronary sinus analysis tool
- Perform one-touch standard 2D reformation and 3D volume rendering, including hollow view, of cardiac anatomy with optimized Volume Viewer and 3D rendering tools

- Use predefined navigator views of the pulmonary veins to get a fly-through perspective of the vessels.

Image Requirements

CardEP accepts standard cardiac gated CT image sets acquired on CT scanners equipped with DICOM 3.0 standards. Images must meet the same image requirements as those for the basic Volume Viewer application.

System Requirements

The CardEP option can be installed on Advantage Workstation (AW) or AW Server with VolumeShare 7 or later software.

Compatible hardware:

- Z800, Z820, Z440 (and later)
- AW Server 3.2 and above

(recommended monitor resolution is up to dual 2MP (1600 x 1200) or a single 3MP (1536 x 2048))

Indications for Use

CardEP is a post-processing software option for the Advantage Workstation (AW) Platform. This product can be used for the analysis of CT angiographic images for the assessment of the heart to include the atria, pulmonary veins and coronary sinus. It provides quantitative analysis tools which include a number of display, measurement and model export capabilities. This product can be used to aid trained physicians in the visualization and assessment of cardiac anatomy.

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SECRET DE
AFACERI

GE imagination at work

© 2016 General Electric Company.
All rights reserved. Data subject to change.
GE and GE Monogram are trademarks of General Electric Company.
* Trademark of General Electric Company
DOC1117927

Feature Details

Operation

- On the menu page, multiple image processing tools for both single and multiple cardiac phase studies are provided. Access each pre-set protocol from the AW Volume Viewer application or from the application pull-down menu page on the browser.

- CardEP images present you with a detailed display of patient heart rate, cardiac scanning parameters, and each scan's temporal resolution.

3D VR Atrium

- Pre-set protocols enable you to semi-automatically generate and display 2D or 3D left atrial/ventricular views for a qualitative overview of the heart's chambers.

- The atrial/ventricular model can be either 3D volume rendered or Maximum Intensity Projection (MIP).

- With the endocardial view feature, you can view inside the atrium to analyze the relationship of the pulmonary vein ostia from the chamber's endocardial surface.
- With a one-touch protocol you can quickly change from a 3D model to an endoscopic view of the atrium and pulmonary veins.
- Anatomical landmark selection lets you deposit points on the 3D model, embed them, and then export the model for future registration use with a mapping system.
- Save predefined landmarks and models and export them to CD.
- Save current model state to the AW browser and reload for future analysis.
- Optimized VR and curve setting provides clearer contrast enhanced cardiac display.
- Choose from single or multi cardiac phase sets.
- Add/remove structure (Autoselect), scalping and brush threshold tools are all available for quick optimization of a 3D model.

3D VR Heart

- Automatically extracts cardiac anatomy from within the chest cavity.


- Optimized VR and curve setting

provides clearer contrast enhanced cardiac display.

- Automatic segmentation for single or multi-phase cardiac data.
- The program supports multiphase thresholding, filtering, and scalping techniques.
- Create movies of rotating and beating hearts with easy export capabilities.
- Export enhanced 3D models, including anatomical markers.

Image Reformation

- With batch processing tools, you can reformat either single or multi-phase cardiac images in various non-axial orientations, e.g. short axis, long axis, chamber views. Save reformatted images in separate image series for further analysis (e.g. from PACs).

- Offers predefined protocols to generate automatically short and long axis views.
- Multi-phase image capability lets you quickly and efficiently assess phases or function.
- View the cross-sectional diameter of the pulmonary vein and left atrial appendage with specifically designed protocols.
- Measurement tools integrated


CONFIDENTIAL

SECRET DE
AFACERI

with protocols enable you to quickly ascertain length and cross-sectional diameter measurements.

Coronary Sinus and Pulmonary Vein Analysis

- Load single or multi-phase cardiac data sets in the vessel analysis tool.

- With the semi-automated pulmonary vein/coronary sinus analysis tool (CSA), based on the Advanced Vessel Analysis (AVA) technique, you can track, extract, visualize, and measure the veins and sinus (physical dimension and stenosis sizing) from either single or multiple cardiac phase image data sets.
- Generate curved, oblique, longitudinal, and cross-sectional reformatted views of individual vessels.
- Generate 3D models of the atrium and pulmonary veins during your vessel analysis procedure.
- Save current state tracking points within your vessel analysis for future image review or manipulation.
- Generate and save rotational movies from curved planar reformation and cross-sectional (lumen) views.
- Generate and save movies from Best L Section reformatted images.
- Perform quantitative and qualitative coronary sinus assessments on segment you

choose:

- Distance measurement
- Cross-sectional area
- Single or dual reference point comparison
- Relative percentage stenosis
- Volume measurement
- Mean diameter
- Generate beating, rotating movies for visualization and communicating findings.

Filming

CardEP offers you filming flexibility. Batch filming capability is provided, and any image can be queued for filming.

One Touch EP Views

- Default views provide you with visualization of the atria and pulmonary veins in correct anatomical orientation. You are provided the capability to create your own presets to save and apply to future analyses.
- Use the interactive RAO/LAO and CRA/CAU annotation feature to manually orient the view.

Pulmonary Vein (PV) and Left Atrial Appendage Reformation

- Predefined protocols allow you to generate reformation of the pulmonary veins and left atrial appendage in the long axis.
- Generate cross-sectional views of the veins or appendage to determine the ostium's diameter.

EP Navigator Views

- With one-touch, the navigator displays endoscopic views of the left atrium.
- View the pulmonary vein and atrial appendage ostia from either the left or right side using

a predefined orientation layout preset.

- Optimize interior atrium views with automated threshold settings.
- Global virtual endoscopic view provides you ostial displays of the right and left pulmonary veins and the left atrial appendage.

Batch Movie Mode

- Batch movie mode lets you define first and last images for movie creation.
- Create rotating and beating movies
- Save automatically generated movies to the viewer as DICOM image sets or send them to data export in JPEG/MPEG/AVI formats.
- Use the movie mode with MIP, 3D VR, and reformatted images.
- Generate and save movies from Best L Section Reformatted images, curved planar reformatted, cross-sectional, and lumen views.

Protocols

CardEP contains predefined protocols that you can edit or modify to adapt them to your specific customized protocols and save under a new name. The CardEP protocols comprise comprehensive instructions and the tools you need to perform analyses. Tool icons are displayed in the protocol panels for fast and easy access.

Summary

CardEP gives you the tools and flexibility to quickly analyze cardiac studies. The information and statistics you can generate with this program may help referring physicians with their assessment of therapy options for their patients in the visualization and assessment of cardiac anatomy.

Indications for Use

CardEP is a post-processing software option for the Advantage Workstation (AW) Platform. This product can be used for the analysis of CT angiographic images for the assessment of the heart to include the atria, pulmonary veins and coronary sinus. It provides quantitative analysis tools which include a number of display, measurement and model export capabilities. This product can be used to aid trained physicians in the visualization and assessment of cardiac anatomy.

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SECRET DE
AFACERI

GE imagination at work

© 2016 General Electric Company.
All rights reserved. Data subject to change.
GE and GE Monogram are trademarks of General Electric Company.
* Trademark of General Electric Company.
D0C1117327

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