EIAgen

HCV Ab (v.4) Kit

REF 071067

Σ 96

REF 071064

Σ 192

REF 071068

Σ 480

IVD

CE 0459

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BIBLIOGRAFIE

1. CDC. Public Health Service inter-agency guidelines for screening donors of blood, plasma, organs, tissues, and semen for evidence of hepatitis B and hepatitis C. *MMWR* 1991;40(No. RR-4):1-17.
2. Alter MJ. Epidemiology of hepatitis C. *Hepatology* 1997;26:62S-5S.
3. McQuillan GM, Alter MJ, Moyer LA, Lambert SB, Margolis HS. A population based serologic study of hepatitis C virus infection in the United States. In Rizzetto M, Purcell RH, Gerin JL, Verme G, eds. *Viral Hepatitis and Liver Disease*, Edizioni Minerva Medica, Turin, 1997, 267-70.
4. Dufour MC. Chronic liver disease and cirrhosis. In Everhart JE, ed. *Digestive diseases in the United States: epidemiology and impact*. US Department of Health and Human Services, Public Health Service, National Institutes of Health, National Institute of Diabetes and Digestive and Kidney Diseases. Washington, DC: US Government Printing Office, 1994; NIH publication no. 94-1447, 615-45.
5. Alter MJ, Hadler SC, Judson FN, et al. Risk factors for acute non-A, non-B hepatitis in the United States and association with hepatitis C virus infection. *JAMA* 1990;264:2231-35.
6. Alter HJ, Holland PV, Purcell RH, et al. Posttransfusion hepatitis after exclusion of commercial and hepatitis-B antigen-positive donors. *Ann Intern Med* 1972;77:691-9.
7. Alter HJ, Purcell RH, Holland PV, Feinstone SM, Morrow AG, Moritsugu Y. Clinical and serological analysis of transfusion-associated hepatitis. *Lancet* 1975;2:838-41.
8. Seeff LB, Wright EC, Zimmerman HJ, McCollum RW, VA Cooperative Studies Group. VA cooperative study of post-transfusion hepatitis and responsible risk factors. *Am J Med Sci* 1975;270:355-62.
9. Feinstone SM, Kapikian AZ, Purcell RH, Alter HJ, Holland PV. Transfusion-associated hepatitis not due to viral hepatitis type A or B. *N Engl J Med* 1975;292:767-70.
10. Choo QL, Kuo G, Weiner AJ, Overby LR, Bradley DW. Isolation of a cDNA clone derived from a blood-borne non-A, non-B viral hepatitis genome. *Science* 1989;244:359-62.
11. Kuo G, Choo QL, Alter HJ, et al. An assay for circulating antibodies to a major etiologic virus of human non-A, non-B hepatitis. *Science* 1989;244:362-4.
12. Alter HJ, Purcell RH, Shih JW, et al. Detection of antibody to hepatitis C virus in prospectively followed transfusion recipients with acute and chronic non-A, non-B hepatitis. *N Engl J Med* 1989;321:1494-1500.
13. Aach RD, Stevens CE, Hollinger FB, et al. Hepatitis C virus infection in post-transfusion hepatitis. An analysis with first- and second-generation assays. *N Engl J Med* 1991;325:1325-9.
14. Alter MJ, Margolis HS, Krawczynski K, Judson, FN, Mares A, Alexander WJ, et al. The natural history of community-acquired hepatitis C in the United States. *N Engl J Med* 1992;327:1899-1905.
15. Alter, MJ. Epidemiology of hepatitis C in the west. *Semin Liver Dis* 1995;15:5-14.
16. Donahue JG, Nelson KE, Muñoz A, et al. Antibody to hepatitis C virus among cardiac surgery patients, homosexual men, and intravenous drug users in Baltimore, Maryland. *Am J Epidemiol* 1991;134:1206-11.
17. Zeldis JB, Jain S, Kuramoto IK, et al. Seroepidemiology of viral infections among intravenous drug users in northern California. *West J Med* 1992;156:30-5.
18. Fingerhood MI, Jasinski DR, Sullivan JT. Prevalence of hepatitis C in a chemically dependent population. *Arch Intern Med* 1993;153:2025-30.
19. Garfein RS, Vlahov D, Galai N, Doherty, MC, Nelson, KE. Viral infections in short-term injection drug users: the prevalence of the hepatitis C, hepatitis B, human immunodeficiency, and human Tlymphotropic viruses. *Am J Pub Health* 1996;86:655-61.
20. Brettler DB, Alter HJ, Deinstag JL, Forsberg AD, Levine PH. Prevalence of hepatitis C virus antibody in a cohort of hemophilia patients. *Blood* 1990;76:254-6.
21. Troisi CL, Hollinger FB, Hoots WK, et al. A multicenter study of viral hepatitis in a United States hemophilic population. *Blood* 1993;81:412-8.
22. Kumar A, Kulkarni R, Murray DL, et al. Serologic markers of viral hepatitis A, B, C, and D in patients with hemophilia. *J Med Virology* 1993;41:205-9.
23. Tokars JI, Miller ER, Alter MJ, Arduino MJ. National surveillance of dialysis associated diseases in the United States, 1995. *ASAIO Journal* 1998;44:98-107.
24. Osmond DH, Charlebois E, Sheppard HW, et al. Comparison of risk factors for hepatitis C and hepatitis B virus infection in homosexual men. *J Infect Dis* 1993;167:66-71.
25. Weinstock HS, Bolan G, Reingold AL, Polish LB: Hepatitis C virus infection among patients attending a clinic for sexually transmitted diseases. *JAMA* 1993;269:392-4.
26. Thomas DL, Cannon RO, Shapiro CN, Hook EW III, Alter MJ. Hepatitis C, hepatitis B, and human immunodeficiency virus infections among non-intravenous drug-using patients attending clinics for sexually transmitted diseases. *J Infect Dis* 1994;169:990-5.
27. Buchbinder SP, Katz MH, Hessel NA, Liu J, O'Malley PM, Alter, MJ. Hepatitis C virus infection in sexually active homosexual men. *J Infect* 1994;29:263-9.
28. Thomas DL, Zenilman JM, Alter HJ, et al. Sexual transmission of hepatitis C virus among patients attending sexually transmitted diseases clinics in Baltimore--an analysis of 309 sex partnerships. *J Infect Dis* 1995;171:768-75.
29. Thomas DL, Factor SH, Kelen GD, Washington AS, Taylor E Jr, Quinn TC. Viral hepatitis in health care personnel at The Johns Hopkins Hospital. *Arch Intern Med* 1993;153:1705-12.
30. Cooper BW, Krusell A, Tilton RC, Goodwin R, Levitz RE. Seroprevalence of antibodies to hepatitis C virus in high-risk hospital personnel. *Infect Control Hosp Epidemiol* 1992;13:82-5.
31. Abdel-Hamid, M., M. El-Day, S. El-Kafrawy, N. Mikhail, G.T. Strickland, and A.D. Fix. 2002. Comparison of second- and third-generation enzyme immunoassays for detecting antibodies to hepatitis C virus. *J. Clin. Microbiol.* 40:1656-1659.
32. Dusheiko, G., H. Schlimovitz-Weiss, D. Brown, F. McOmish, P.-L. Yap, S. Sherlock, N. McIntyre, and P. Simmonds. 1994. Hepatitis C virus genotypes: an investigation of type-specific differences in geographic origin and disease. *Hepatology* 19:13-18.
33. Gretch, D. Diagnostic tests for hepatitis C. The article can be found at: <http://www.hepnet.com/nih/gretch.html>. It was written as part of a National Institute of Health Conference on Hepatitis C, held from March 24-26, 1997 in Bethesda Maryland
34. Mondelli, M.U., A. Cerino, F. Bono, A. Cividini, A. Maccabruni, M. Aricò, A. Malfitano, G. Barbarini, V. Piazza, L. Minoli, and E. Silini. 1994. Hepatitis C virus (HCV) core serotype in chronic HCV infection. *J. Clin. Microbiol.* 32:2523-2527.
35. Ohno, T., M. Mizokami, R.-R. Wu, M.G. Saleh, K.-I. Ohba, E. Orito, M. Mukaide, R. Williams, and J.Y.N. Lau. 1997. New hepatitis C virus (HCV) genotyping system that allows for the identification of HCV genotype 1a, 1b, 2a, 2b, 3a, 3b, 4, 5a, and 6a. *J. Clin. Microbiol.* 35:201-207.
36. Takada, N., S. Takase, N. Enomoto, A. Takada, and T. Date. 1992. Clinical backgrounds of the patients having different types of hepatitis C virus genome. *J. Hepatol.* 14:35-40.
37. Yoshioka, K., S. Kakumu, T. Wakita, T. Ishikawa, Y. Itoh, M. Takayanagi, Y. Higashi, M. Shibata, and T. Morishima. 1992. Detection of hepatitis C virus by polymerase chain reaction and response to interferon- α therapy: relationships to genotypes of hepatitis C virus. *Hepatology* 16:293-299.

ATLAS SLE SLIDE TEST

IVD For *in vitro* diagnostic and professional use only

2°C  8°C
Store at 2°-8°C

INTENDED USE

Atlas SLE Slide Test is a slide agglutination assay for the qualitative and semi quantitative detection of anti-deoxyribonucleoprotein (anti-DNP) in human serum. No initial dilution of patient samples is required for this test. These materials are intended to be acquired, possessed and used only by health professionals.

INTRODUCTION

The detection of antinuclear antibodies, by such laboratory methods as immunofluorescence, LE cell test, and agglutination of coated particles, can aid in the diagnosis of such autoimmune diseases as systemic lupus erythematosus (SLE). The antibodies most associated with SLE are those directed against DNP. These antibodies are believed to cause the formation of the LE cell *in vitro*, occurring in 75-80% of patients diagnosed as having SLE. Given that 20-25% of SLE patients do not exhibit the formation of LE cells, other methods can be used to detect antinuclear antibodies.

PRINCIPLE

Atlas SLE Slide Test provides a means of detecting anti-DNP in human serum. SLE Slide reagent is a stabilized buffered suspension of polystyrene latex particles that have been coated with DNP. When the latex reagent is mixed with the serum containing antibodies to DNP, agglutination occurs. Using dilutions of a reactive patient sample, the anti-DNP titer can be determined.

MATERIALS

MATERIALS PROVIDED

- SLE Latex Reagent: Suspended inert latex particles coated with DNP, with 0.1% sodium azide as preservative.
- SLE Positive Human serum or defibrinated plasma (liquid), with 0.1% sodium azide as preservative.

- SLE Negative Control: Non-reactive buffer containing BSA and 0.1% sodium azide.
- Stirring sticks.
- Glass slide.
- Package insert.

MATERIALS NEEDED BUT NOT PROVIDED

- Timing device.
- 13 x 75 mm test tubes
- Volumetric pipet to deliver 0.25 ml
- Saline (0.9% NaCl solution)
- Mechanical rotator (optional)

PACKAGING CONTENTS

REF 8.00.11.0.0025 (1x1 mL Latex, 1x0.5 mL Positive Control, 1x0.5 mL Negative Control)

REF 8.00.11.0.0050 (1x2 mL Latex, 1x0.5 mL Positive Control, 1x0.5 mL Negative Control)

REF 8.00.11.0.0100 (1x4 mL Latex, 1x1 mL Positive Control, 1x1 mL Negative Control)

PRECAUTIONS

- For *in vitro* diagnostic use.
- Latex reagent and controls contain sodium azide. Azides in contact with lead and copper plumbing may react to form highly explosive metal azides. When disposing of reagents containing azide, flush down the drain with large quantities of water to prevent azide build-up.
- The controls contain human serum or plasma which has been tested at the donor level for HBsAg and for HIV-1, HIV-2 and HCV antibodies and found to be nonreactive. As no known test offers complete assurance that infectious agents are absent, the controls should be considered potentially infectious and universal precautions should be used.
- Do not pipet by mouth.
- Do not smoke, eat, drink or apply cosmetics in areas where plasma/serum samples are handled.
- Any cuts, abrasions or other skin lesions should be suitably protected.
- In order to obtain reliable and consistent results, the instructions in the package insert must be strictly followed.
- Do not use past the expiration date indicated on the kit.
- Do not interchange components of one kit with those of another kit.

- Turbidity or precipitation in controls is indicative of deterioration and the component should not be used.
- Bacterial contamination of reagents or specimens may cause false positive results.

STORAGE & STABILITY

- Store all reagents at 2-8°C in an upright position when not in use.
- Do not freeze reagents.

SPECIMEN COLLECTION and STORAGE

- Use only serum that is free from contamination. Test samples should not be heat-inactivated.
- It is preferable to test samples on the day of their collection. If samples cannot be tested immediately, maintain them in their original tubes at 2-8°C and test within 48 hours.
- Serum samples stored longer than 48 hours should be stored at -20°C or below until testing. Avoid repeated freezing and thawing of specimens.
- If necessary before testing, centrifuge the specimens at a force sufficient to sediment cellular components.
- Samples to be sent out for testing should be placed on ice packs and packaged like any other biohazardous material that could potentially transmit infection.

REAGENT PREPARATION

- Allow all reagents and samples to warm to room temperature (20-30°C) before use. Do not heat reagents in a water bath.
- All reagents are ready for use as supplied. Gently mix the reagents before use; avoid foaming.
- Gently mix the latex reagent before each use to ensure homogeneity.

PROCEDURES

A. Method I (Qualitative)

1. Dispense (35 µL) of each serum sample onto a separate circle on the test slide. Add one drop of Positive and negative controls from the dropper vials supplied onto a separate circle on the test slide.
2. Dispense one drop of latex reagent (35 µL) to each serum specimen and to each control.
3. Using the flat end of the stirring sticks, mix each specimen and control serum with the latex reagent, in a circular manner, over the entire area in the circles of the card.

- Gently tilt and rotate the card for one (1) minute and observe for agglutination. All test results should be compared to both positive and negative controls.

INTERPRETATION OF RESULTS (QUALITATIVE)

Agglutination indicates a reactive SLE sample. Sera that elicit a reactive result should be retested and tittered using the "Semi quantitative Assay Protocol".

B. Method II (Semi-Quantitative)

- Prepare serial dilutions of patient serum, in saline, in test tubes as follows:

Tube	Dilution	Composition
1	1:2	0.25 ml of serum + 0.25 ml saline.
2	1:4	0.25 ml from tube 1 + 0.25 ml saline.
3	1:8	0.25 ml from tube 2 + 0.25 ml saline.
4	1:16	0.25 ml from tube 3 + 0.25 ml saline.
5	1:32	0.25 ml from tube 4 + 0.25 ml saline.
6	1:64	0.25 ml from tube 5 + 0.25 ml saline.

Note: Testing on additional dilutions should be performed as needed.

- Using each dilution as a separate test specimen, apply the samples to the slide as described in Step 1 of the "Qualitative method" and proceed with Steps 2 through 4 of the "Qualitative method". Include undiluted sample if not tested previously on that day with the same lot of latex reagent.

INTERPRETATION OF RESULTS (SEMI-QUANTITATIVE)

The highest dilution in which visible agglutination occurs is considered the endpoint titer.

QUALITY CONTROL

Quality Control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control Procedures. Controls with graded reactivity should be included. If control samples do not yield the expected response, the assay should be considered invalid and the assay repeated. If the repeat assay does not elicit the expected results for the control samples, discontinue use of the kit and contact your local distributor.

EXPECTED VALUES

Serum samples from 155 individuals were tested using the **SLE Slide Test**. Of the 155 individuals, 29 had active SLE, 23 had clinically inactive SLE, 8 had connective tissue diseases and the remaining 95 were either clinically normal or had some nonrelated disease (including anemia, infectious mononucleosis and rheumatic heart disease) and were used

as controls. Results from testing with the **SLE Slide Test** were compared with the results from testing of the samples using a standard LE cell preparation assay and a fluorescent ANA assay.

Of the 29 active SLE patients, 82% were positive using the SLE Slide Test, 86% were positive by the LE cell prep, and 82% positive by the ANA test. For the 23 clinically inactive SLE patients, 19% were positive by both the SLE Slide Test and the LE cell prep; and 71% were positive by the ANA test. None of the 8 patients having connective tissue disease tested positive with the SLE Latex Test, whereas 17% and 50% tested positive by the LE cell prep and the ANA procedures, respectively. Of the controls, 1% tested positive by both the SLE Latex Test and the LE cell prep, while 6% tested positive by the ANA assay.

LIMITATION

- Serum from patients with scleroderma, rheumatoid arthritis, dermatomyositis, and a variety of connective tissue diseases may elicit agglutination in the SLE slide test.
- Because extremely high levels of antibodies might affect the degree of agglutination, positive samples should be reassayed using the semi quantitative procedure.
- Contaminated, lipemic, or grossly hemolyzed sera should not be used because of the possibility of nonspecific results.
- Plasma samples should not be used because of the possibility of nonspecific results.
- Samples yielding indeterminate results may be resolved by repeating the test utilizing a two (2) minute slide rotation period. Reaction times longer than two minutes might cause false positive results due to a drying effect.
- Drugs such as hydralazine, isoniazid, procainamide and a number of anticonvulsant drugs can induce an SLE syndrome.
- In accord with all diagnostic methods, a final diagnosis should not be made on the result of a single test, but should be based on a correlation of test results with other clinical findings.

REFERENCES


- Christian CL, Mendez-Bryan R, Larson DL. 1958. *Proc Exp Biol Med*, **98**:820-823.
- Friou GJ, Finch SC, Detre KD. 1958. *J Immunol*, **80**:324-329.
- Hargraves MM, Richmond H, Morton R. 1948. *Proc Mayo Clin*, **23**:25-28.
- Holman HR, Kunkel HG. 1957. *Science*, **126**:163.
- Miescher PA, Strassel R. 1957. *Vox Sang*, **2**:283-287.
- Miescher PA, Rothfield N, Miescher A. 1966. *Lupus Erythematosus*, EL Dubois, Ed., Blakiston Co., New York.

- Rothfield NF, Phythyon JJ, McEwan C., Miescher P. 1961. *Arth Rheum*, **4**:223-229.

 **ATLAS Medical GmbH**
Ludwig-Erhard Ring 3
15827 Blankenfelde-Mahlow
Germany
Tel: +49 - 33708 – 3550 30
Email: Info@atlas-medical.com
Website: www.atlas-medical.com

PP12339A01

Rev B (22.06.2023)

	Catalogue Number		Temperature limit
	In Vitro diagnostic medical device		Caution
	Contains sufficient for <n> tests and Relative size		Consult instructions for use (IFU)
	Batch code		Manufacturer
	Fragile, handle with care		Use-by date
	Manufacturer fax number		Do not use if package is damaged
	Manufacturer telephone number		Date of Manufacture
	Keep away from sunlight		Keep dry
	Positive control		Negative control

RPR SYPHILIS CARD TEST

IVD For In-Vitro diagnostic and professional use only

8°C
2°C 1 Store at 2 to 8 °C

INTRODUCTION

Syphilis is a disease caused by infection with the spirochete *Treponema pallidum*. The infection is systemic and the disease is characterized by periods of latency. These features, together with the fact that *T pallidum* cannot be isolated in culture, mean that serologic techniques play a major role in the diagnosis and follow-up of treatment for syphilis. Syphilis is categorized by an early primary infection in which patients may have non-specific symptoms, and potentially, genital lesions. Patients tested by serology during the primary phase may be negative for antibodies, especially if testing is performed during the first 1 to 2 weeks after symptom onset. As the disease progresses into the secondary phase, antibodies to *T pallidum* reach peak titers, and may persist indefinitely regardless of the disease state or prior therapy. Therefore, detection of antibodies to nontreponemal antigens, such as cardiolipin (a lipoidal antigen released by host cells damaged by *T pallidum*) may help to differentiate between active and past syphilis infection. Nontreponemal antibodies are detected by the rapid plasma reagin (RPR) assay, which is typically positive during current infection and negative following treatment or during late/latent forms of syphilis.

PRINCIPLE

RPR utilises carbon particles coated with cardiolipin antigen to detect reagin antibodies present in serum or plasma of syphilitic persons. Specimens that contain reagin cause aggregation of the carbon particles which appear as dark clumps against a white background. The aggregation can be read macroscopically. Non-reactive samples typically appear as a smooth non-aggregated pattern which may form buttons in the centre of the test area.

MATERIALS

MATERIALS PROVIDED

- **RPR carbon antigen reagent:** Contains less than 0.1% sodium azide.
- **Positive Control:** Contains less than 0.1% sodium azide.
- **Negative control:** Contains less than 0.1% sodium azide.

- RPR test cards (Optional).
- Plastic sticks.
- Package insert.

NOTE: This package insert is also used for Individually MATERIALS NEEDED BUT NOT PROVIDED

- Rotator (100rpm).
- Timer.
- Pipettes.

PACKAGING CONTENT

- REF 8.00.18.0.0100** (2mL Latex, 1x0.5mL Positive Control, 1x0.5mL Negative Control)
- REF 8.00.18.0.0500** (10mL Latex, 1x1mL Positive Control, 1x1mL Negative Control)

SAMPLES

Fresh serum or plasma. The samples with presence of fibrin should be centrifuged before testing. Do not use highly hemolized or lipemic samples.

PRECAUTIONS

- For professional in vitro diagnostic use only. Do not use after expiration date.
- Do not eat, drink or smoke in the area where the specimens or kits are handled.
- Always use a fresh pipette tip for every test.
- Handle all negative and positive in the manner as patient specimens.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are assayed.
- The used test should be discarded according to local regulations.
- Components of different human origin have been tested and found to be negative for the presence of antibodies anti- HIV 1+2 and anti-HCV, as well as for HBsAg. However, the controls should be handled cautiously as potentially infectious.

STORAGE AND STABILITY

All components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C.

PROCEDURES

QUALITATIVE PROCEDURE

- **Mix well the RPR reagent before use.**
- 1. Bring the reagents and samples to room temperature.
- 2. Dispense 50 µL of each sample into a separate circle on the card. Use a separate tip for each sample.
- 3. Dispense 1 drop of each of positive and negative controls into two additional circles.

4. Gently shake the dispensing vial and slightly press to remove air bubbles from the needle and the drop obtained is correct.
5. Dispense 1 drop (17.5 µl) of RPR antigen to each circle next to the sample to be tested.
6. Place the card on a mechanical rotator and rotate at 100 r.p.m. for 8 minutes.
7. Observe macroscopically for agglutination within a minute after removing the card from the rotator.

SEMI-QUANTITATIVE PROCEDURE

- **Mix well the RPR reagent before use.**
- 1. Make doubling dilutions from Undiluted to 1:16 normal saline.
- 2. Place 50 µl of each dilution in to a separate circle on the test card.
- 3. Spread each dilution evenly over the test circle.
- 4. Continue as from Qualitative procedure. The titer of the sample is expressed as the final dilution which shows aggregation of the carbon particles.

PERFORMANCE CHARACTERISTICS

1. Sensitivity: 100%.
2. Specificity: 100%.

INTERPRETATION OF TEST RESULTS

1. **Strong Reactive:** Large clumps of carbon particles with a clear background.

2. **Reactive:** Large clumps of carbon particles somewhat more disperse than Strong Reactive pattern.

3. **Weak Reactive:** Small clumps of carbon particles with light grey background.

4. **Trace Reactive:** Slight clumping of carbon particles typically seen as a button of aggregates in the centre of the test circle or dispersed around the edge of the test circle.

5. **Non-Reactive:** Typically a smooth grey pattern or a button of non-aggregated carbon particles in the centre of the test circle.

REFERENCES

- Falcone V.H., Stout G.W. and Moore M.B. Jr., PHR 79: 491-495, 1964.

ATLAS Medical GmbH
 Ludwig-Erhard Ring 3
 15827 Blankenfelde-Mahlow
 Germany
 Tel: +49 - 33708 - 3550 30
 Email: info@atlas-medical.com
 Website: www.atlas-medical.com

PPI2280A01
 Rev B (06.05.2023)


REF	Catalogue Number	Temperature limit
IVD	In Vitro diagnostic medical device	Caution
Biohazard	Contains sufficient for C_{50} tests and Relative site	Consult instructions for use (IFU)
LOT	Batch code	Manufacturer
Fragile, handle with care		Use-by date
Manufacturer fax number		Do not use if package is damaged
Manufacturer telephone number		Date of Manufacture
Keep away from sunlight		Keep dry

HAV IgM

**ELISA kit for the qualitative detection of IgM
antibodies to hepatitis A virus**

Instruction for use

IVD

REF
EI-031

Σ 96
tests

CE

EQUI HAV IgM

ELISA kit for the qualitative detection
of IgM antibodies to hepatitis A virus

1. INTENDED USE

The «EQUI HAV IgM» is ELISA kit intended for the qualitative detection of IgM antibodies to hepatitis A virus in human serum or plasma by enzyme-linked immunosorbent assay (ELISA) to diagnose acute hepatitis A. The testing procedure is designed for both manual arrangement with automatic pipettes and standard equipment, and for automated «open» immunoassay analysers.

Target group: blood or organ donors; pregnant women and children born to infected mothers; patients with symptoms of liver disease.

Usage: ELISA kit is used in clinical diagnostic laboratories, blood transfusion stations, as well as in other institutions working in the field of *in vitro* diagnostics.

2. CLINICAL SIGNIFICANCE

One of the most common foodborne infections is hepatitis A. The hepatitis A virus (HAV) causes acute liver disease, which can be mild or severe. Unlike hepatitis B and C, this hepatitis does not become chronic, but can lead to acute liver failure, which is characterized by high mortality.

Hepatitis A virus is a small shellless RNA virus from the *Picornaviridae* family. It is characterized by being highly stable in different environments and can be stored at + 4 ° C for several months, but becomes inactive after 5 minutes at 100 ° C. The virus replicates in liver cells, and then is released through the bile into the environment with the fecal masses of the infected person. The cellular immune response to HAV infection leads to the destruction of hepatocytes, liver dysfunction and the development of symptoms, typical for other types of hepatitis.

Acute hepatitis A, even with the clinical manifestations, does not differ from other viral hepatitis. Therefore, serological markers of infection are used for diagnosis, namely the detection of specific antibodies to HAV antigens. IgM antibodies are detected in the serum 1-2 weeks after infection, with the onset of symptoms or a few days before. In the maximum titer of anti-CAA IgM are detected in the jaundice period, after that their level gradually decreases. In most patients, specific IgM ceases to be detected after 6 months, and may occasionally circulate in the blood for more than a year. IgG antibodies to HAV antigens begin to be released shortly after IgM antibodies and stay in the blood throughout life. Also, specific IgG is produced after a vaccination. Detection of IgG antibodies to CAA indicates the formation of a stable immunity due to infection or immunization.

3. ANALYSIS PRINCIPLE

Detection of specific IgM antibodies to hepatitis A virus in the «EQUI HAV IgM» ELISA kit is based on the principle of «IgM capture» of solid-phase ELISA in a two-stage incubation. Monoclonal antibodies specific for human IgM immunoglobulins are adsorbed into the wells of the plate. During the first step of incubation of the

test samples in the wells of the ELISA plate, IgM immunoglobulins, if present in the samples, bind to monoclonal antibodies in the solid phase. The wells are washed to remove unbound components, leaving only specific antibody-antibody complexes. A mixture of HAV antigen and peroxidase conjugate of HAV-specific antibodies that bind to solid-phase immune complexes is then added. Unbound components are removed during washing. Immune complexes are detected by adding a solution of chromogen 3,3', 5,5'-tetramethylbenzidine (TMB) with hydrogen peroxide. After a 30-minute incubation, the reaction is halted by adding a stop solution. Optical density (OG) in the wells is determined on a spectrophotometer at a wavelength of 450 / 620-695 nm. The intensity of the yellow color is proportional to the number of antibodies in the sample.

4. MATERIALS AND EQUIPMENT

4.1. Contents of the ELISA kit

STRIPS	1 x 96 wells	<p>Microplate Each plate well is coated with monoclonal antibodies specific for human IgM immunoglobulins. The wells are detachable. After the first opening, store unused strips in the package at 2-8 °C for a maximum of 6 months</p>
CONTROL +	1 x 0,25 ml	<p>Positive control The solution of human IgM immunoglobulins crosslinked with monoclonal antibodies specific for horseradish peroxidase, with a preservative (pink). Store at 2-8 °C</p>
CONTROL -	1 x 0,6 ml	<p>Negative control Negative human serum with a preservative (yellow). Store at 2-8 °C</p>
DIL SAMPLE	1 x 13 ml	<p>Serum dilution solution Buffer solution with monoclonal antibodies to human IgG, milk extract, detergent and preservative (brown). Store at 2-8 °C</p>
CONJ 11x	1 x 1,3 ml	<p>Conjugate (11x concentrated) 11-fold concentrate of conjugate of antibodies to hepatitis A virus with horseradish peroxidase in buffer solution with stabilizers (purple). Dilute the conjugate (11x) 1:11 with the conjugate dilution solution before use (eg 100 µl concentrate + 1 ml conjugate dilution solution, enough for 8 wells). Diluted solution should be stored at 2-8 °C for no more than 1 day</p>
DIL CONJ	1 x 13 ml	<p>Conjugate dilution solution Buffer solution of inactivated hepatitis A virus antigen with detergent and preservative (yellow). Store at 2-8 °C</p>
SOLN TMB	1 x 13 ml	<p>TMB solution (ready to use) TMB solution, H₂O₂, a stabilizer, a preservative (colourless). Store at 2-8 °C</p>

TWEEN WASH 20x	1 x 50 ml	Washing solution TWEEN (20x concentrated) 20-fold phosphate buffer concentrate with Tween-20 (colourless). Dilute TWEEN detergent (20x) at 1:20 with distilled or deionized water (e. g., 5 mL of concentrate + 95 mL of water for 8 wells) before use. Store the diluted solution at 2-8 °C for a maximum of 7 days
SOLN STOP	1 x 13 ml	Stop Solution (ready to use) 0.5 mol H ₂ SO ₄ solution (colourless). Store at 2-8 °C

The ELISA kit also includes adhesive films (2 items), sample application plan (1 item), checklist, and instruction for use.

4.2. Optional reagents, materials and equipment

Automatic single and multichannel pipettes 10–1000 µL, tips, volumetric laboratory glassware (10–1,000 mL), deionized or distilled water, thermostat at 37 °C, automatic or semi-automatic plate washer, spectrophotometer (reader) for microplates at 450/620-695 nm, appropriate containers for potentially contaminated waste, timer, filter paper, disposable powder-free gloves, disinfectants.

5. PRECAUTIONS AND SAFETY

5.1. Precautions

Be sure to read the instructions for use carefully before the test. The validity of the test results depends on strict following of the test procedure.

- do not use the ELISA kit components after the expiry date;
- do not use for analysis or mix components of different batches, components of kits for different nosologies, or reagents from other manufacturers with the «EQUI HAV IgM» ELISA kit;
- do not freeze the ELISA kit or its contents;
- after using a reagent, close each vial with its cap;
- when washing, control filling and complete aspiration of solution from the wells;
- use a new pipette tip each time you add samples or reagents;
- prevent direct sunlight from reaching the reagents from the ELISA kit;
- [SOLN|TMB] solution must be colourless before use. Do not use the solution if its colour is blue or yellow. Avoid contact of [SOLN|TMB] with metals or metal ions. Use only clean glassware thoroughly rinsed with distilled water;
- do not use reagents with colour not in line with para. 4.1;
- under no circumstances should the same glassware be used for conjugate solution and [SOLN|TMB];
- do not evaluate the test results visually (without a reader);
- any optional equipment that is in direct contact with biological material or kit components should be considered contaminated and requires cleaning and decontamination;
- the ELISA kit includes materials for 96 tests. Dispose of the used components as well as any remaining unused components.

5.2. Safety requirements

- all reagents in the ELISA kit are for laboratory professional use for *in vitro* diagnosis only and may only be used by qualified personnel;
- conduct the tests in disposable powder-free gloves and goggles only;
- do not eat, drink, smoke, or apply make-up in the test room;
- do not mouth-pipette the solutions;
- controls from the «EQUI HAV IgM» ELISA kit have been tested and found to be for anti-HIV1/2, anti-HCV and anti-*Treponema pallidum* antibodies and HBsAg negative; however, controls and test samples should be handled as potentially hazardous infectious materials;
- some of the kit components contain low concentrations of harmful substances and can damage skin or mucosa. In case of contact of [SOLN|TMB], [SOLN|STOP] and conjugate solution with mucous membranes or skin, immediately wash the affected area with plenty of water;
- in case of spillage of acid-free solutions, e. g. sera, treat the surface with a disinfectant solution and then wipe dry with filter paper. Otherwise first neutralize acid with sodium bicarbonate solution and then wipe the surface dry as described above.

5.3. Waste inactivation and disposal

- the liquid waste must be inactivated, for example, with hydrogen peroxide solution at the final concentration of 6% for 3 hours at room temperature, or with sodium hypochlorite at the final concentration of 5% for 30 minutes, or with other approved disinfectants;
- the solid waste must be inactivated by autoclaving at a temperature not less than 132°C;
- do not autoclave the solutions that contain sodium azide or sodium hypochlorite;
- disposal of inactivated waste must be conducted due to national laws and regulations.

6. STORAGE AND STABILITY

ELISA kit is stable up to the expiry date stated on the label when stored at 2-8°C. The kit should be transported at 2-8°C. Single transportation at a temperature up to 23°C for two days is possible.

7. SAMPLE COLLECTION, TRANSPORTATION AND STORAGE GUIDELINES

Collect blood from the vein into the sterile test tube. Test tube must be marked with patient ID and date of sample collecting. Blood before serum separation can be stored at 2-8 °C for 24 hours, avoiding freezing.

Serum or plasma can be stored at 2-8 °C for maximum 3 days. Frozen serum can be stored for longer periods of time at -20 °C or -70 °C. Thaw frozen samples and keep them at room temperature for 30 minutes before use. After thawing, the stir samples to achieve homogeneity. Avoid repeated freezing-thawing cycles for test samples. If serum (or plasma) is turbid, remove insoluble inclusions by

centrifugation at 3000 rpm for 10-15 minutes. Do not use serum samples with hyperlipidemia, hemolysis, and bacterial growth.

Transport serum samples in insulated containers. To do that, put closed labelled tubes in a plastic bag, tightly seal it and place in the centre of an insulated container. Put the frozen cold packs on the bottom, along the side walls of the insulated container and on top of the serum samples.

8. REAGENT PREPARATION

NOTE! It is very important to keep all ELISA kit components for at least 30 min at room temperature 18-25 °C before the assay!

8.1. Microplate preparation

To prevent water condensation in the wells, keep the **STRIPS** for 30 minutes at a room temperature before opening. Open the vacuum pack, detach the appropriate number of wells, and carefully pack the remaining wells with a desiccant and store tightly zip-locked at 2-8 °C. Storing the packed plate this way ensures its stability for 6 months.

8.2. Washing solution preparation

To prepare detergent, dilute **TWEEN|WASH|20x** at 1:20 (1+19) with distilled or deionized water and stir. E. g., 5 mL of concentrate + 95 mL of water, which is enough for 8 wells. If there are crystals present in the detergent concentrate, heat the vial at 37 °C until the crystals dissolve completely (15–20 minutes). Store the diluted solution at 2-8 °C for a maximum of 7 days.

8.3. Conjugate solution preparation

Working dilution of the conjugate is prepared as follows: dilute **CONJ|11x** (purple) in a clean vial of solution **DIL|CONJ** (yellow) in the ratio 1:11 (ie, 1 + 10), the solution turns green. For example, for 8 well analysis add to 1 ml **DIL|CONJ** 100 µl **CONJ|11x**. The solution of the conjugate in the working dilution is stable during the day when stored at 2-8 °C.

9. ASSAY PROCEDURE

9.1. Prepare the necessary number of wells (four wells for controls and a necessary number of wells for test samples) and insert them into the ELISA plate frame.

Be sure to add control wells in every test run.

9.2. Fill in the sample application plan.

9.3. Prepare the detergent as per para. 8.2.

9.4. Add 90 µL of **DIL|SAMPLE** into each plate well.

9.5. Add 10 µL of controls and test samples into the wells:

CONTROL|+ – into well A1,

CONTROL|- – into wells B1, C1 and D1,

and test samples into the remaining wells.

At the time of adding, the solution changes its colour from brown to blue.

Pipette the mix in the wells carefully to avoid foaming.

9.6. Cover the strips up with adhesive film and incubate for **30 minutes at 37 °C**.

- 9.7. Remove and discard the adhesive film and wash all wells 5 times with automatic washer or 8-channel pipette as follows:
- aspirate the content of all wells into a liquid waste container;
 - add a minimum of 300 µl of diluted washing solution to each well, soak each well for 30 seconds;
 - aspirate the content of all wells again. The residual volume after every aspiration should be less than 5 µl;
 - repeat the washing step 4 more times;
 - after the final aspiration, eliminate extra moisture by tapping the plate against a piece of filter paper.
- 9.8. Prepare conjugate solution as per para. 8.3.
- 9.9. Add 100 µL of conjugate solution into each well. Cover the strips with a new piece of adhesive film and incubate for **60 minutes at 37 °C**.
- 9.10. Following incubation, remove the film carefully and wash the wells five times as described in para. 9.7.
- 9.11. Add 100 µL of [SOLN|TMB] into the wells; do not touch the bottom and the walls of the plate wells.
- 9.12. Incubate the strips for **30 minutes** in a dark place at a room temperature of 18-25 °C. Do not use adhesive film at this stage.
- 9.13. Add 100 µL of [SOLN|STOP] into each strip well to stop the enzymatic reaction; adhere to the same sequence of actions as when adding [SOLN|TMB]. At the time of adding, the solution colour changes from blue to yellow, and clear solution slightly changes its shade.
- 9.14. Measure the optical density (OD) of the wells at 450/620-695 nm wavelength using an ELISA microplate reader within 5 minutes after stopping the reaction. Pay attention to the cleanness of the plate bottom and the absence of bubbles in the wells before reading.

Measurement at the single wavelength of 450 nm is possible, in that case, it is needed to leave one well for blank (only [SOLN|TMB] and [SOLN|STOP] must be added in blank well).

10. CALCULATION AND INTERPRETATION OF RESULTS

10.1. Calculation of results

Calculate the average OD of the negative control (\bar{Nc}), Cut off (CO) and a sample positivity index (IP_{sample}).

$$\bar{Nc} = (Nc1 + Nc2 + Nc3)/3; \quad CO = \bar{Nc} + 0,3$$

$$IP_{\text{sample}} = OD_{\text{sample}} / CO, \text{ where } OD_{\text{sample}} \text{ is the OD sample.}$$

10.2. Quality control (assay validation)

The test results are considered valid if they meet the following requirements:

CONTROL +

OD ≥ 1,2

CONTROL -

OD ≤ 0,150

CONTROL -

$\bar{N}_c \times 0,5 \leq N_{cn} \leq \bar{N}_c \times 2,0$ where N_{cn} is the OD for each N_c run

If any of the OD values for the negative control is beyond the above interval, it should be discarded, and \bar{N}_c is calculated based on the remaining OD values for the negative control. If several OD values for the negative control fail to meet the above requirements, the test is considered invalid and requires a new run.

10.3. Interpretation of results

$IP_{\text{sample}} > 1,1$	POSITIVE
$0,9 \leq IP_{\text{sample}} \leq 1,1$	BORDERLINE*
$IP_{\text{sample}} < 0,9$	NEGATIVE

* Uncertain samples are recommended to be re-examined in two wells of the ELISA kit. If the results are again uncertain, a new sample should be selected and analyzed in 2-4 weeks. In case of repeated indeterminate results, such samples shall be considered negative.

11. PERFORMANCE CHARACTERISTICS

11.1. Analytical performance characteristics

Precision of measurement

Intra assay repeatability

The coefficient of variation (CV) for two sera with different levels of specific antibodies was evaluated in 32 replicates on one series of ELISA kits.

Sample No.	OD _{av}	IP _{av}	CV, %
37s/2	1,576	4,85	5,3
24s	2,462	7,57	4,3

Inter assay reproducibility

The coefficient of variation (CV) for two sera with different levels of specific antibodies was evaluated for 4 days in 4 sets of analysis, 8 replicates in each analysis.

Sample No.	OD _{av}	IP _{av}	CV, %
37s/2	1,600	4,78	6,3
24s	2,463	7,36	7,1

Analytical specificity

The test results are not affected by bilirubin at up to 0.21 mg/mL (361.8 μmol/L), haemoglobin at up to 10 mg/mL and triglycerides at up to 10 mg/mL (11.3 mmol/l) present in the sample.

11.2. Diagnostic characteristics

The diagnostic characteristics of the ELISA kit were evaluated by examining a set of 30 samples containing IgM antibodies to hepatitis A virus, a set of donor serum samples (188 samples) and samples of PHT202 Anti-Hepatitis A Virus (HAVed) panel. Performance Panel (contains 8 positive and 13 negative samples) - a total of 239 samples - was compared to similar commercial kits. For this set (239 samples) the relative sensitivity of the «EQUI HAV IgM» ELISA kit is 100%, the relative specificity - 100%, the percentage of coincidence - 100%.

12. LIMITATIONS OF ASSAY

A positive result in the «EQUI HAV IgM» ELISA kit is the evidence that the patient has IgM antibodies specific for hepatitis A virus. Anti-CAA specific IgM antibodies are usually markers of active replication of hepatitis A virus.

In order to counteract the false-positive results caused by the presence of autoantibodies specific for class G immunoglobulins (rheumatoid factor) in human serum samples, the kit uses a special block component that prevents the formation of immune complexes with anti-human antibodies in the solid phase.

The final diagnosis cannot be established solely on the basis of serological test results. When making a diagnosis the results of a set of laboratory and instrumental studies, as well as clinical manifestations of the disease should all be taken into account.

13. DIFFICULTIES THAT CAN OCCUR DURING THE ASSAY PROCEDURE

Possible reasons	Solution
<i>High background in all wells</i>	
Contaminated washer	Clean the washer head and rinse according to the instructions for use
Poor quality or contaminated water	Use purified water with specific resistance $\geq 10 \text{ M}\Omega \cdot \text{cm}$
Use of poorly washed glassware	Use chemically clean utensils
Use of chlorinated disinfectants	Do not use chlorine disinfectants
Use of contaminated tips	Use new tips
Increased incubation times or change in the temperature conditions	Adhere to the incubation regime according to the instructions for use
<i>High background in a row of wells</i>	
Repeat application of TMB solution	TMB solution should be applied once
Contamination of the automatic pipette nozzle with conjugate solution	Clean the pipette and dial carefully liquid
Contamination of one of the washer's channel	Clean the flush channel, rinse washer

<i>Received OD of the positive control is below the border value</i>	
One of the reagents (conjugate solution or TMB solution) was not prepared in a correct way or was not added	Re-conduct ELISA, pay attention to the correctness of the introduction of these reagents
Reduced incubation times at any stage	Incubate according to instructions for use
<i>The colour density of the wells fails to meet the obtained optical density value</i>	
This may suggest that the optical beam has been displaced	Check the correct operation of the reader

14. TECHNICAL ASSISTANCE AND CUSTOMER SERVICE

In case of technical problems, you can obtain assistance by contacting the manufacturer.

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Manufacturer

Authorized Representative in the European Community

In vitro diagnostic medical device

Catalogue number

Date of manufacture

Use by date

Batch code

Temperature limit

Contains sufficient for <n> tests

Caution

Non-Sterile

Consult instructions for use

Keep away from sunlight

Keep dry

Compliance with EU safety requirements

Edition 6, 01.11.2021

For questions and suggestions regarding the ELISA kit contact:

Obelis s.a.
Bd Général Wahis 53
1030 Brussels
Belgium
Tel: +(32)2 732-59-54
Fax: +(32)2 732-60-03
mail@obelis.net

Ekvitestlab LLC
Velyka Vasylykivska St. 114, Kyiv, Ukraine, 03150
Tel: 0(800)31-89-87, +38 (044)334-89-87,
e-mail: info@equitest.com.ua, www.equitest.com.ua

ASSAY PROCEDURE SCHEME

Keep all reagents for 30 min at temperature 18-25°C before use

Dispense 90 µl [DIL|SAMPLE] into the wells
(brown)

Add to 10 µl of controls and samples into the wells:
A1 – [CONTROL|+], B1, C1, D1 – [CONTROL|-],
other wells – examined samples
(change of colour from brown to blue)

Cover strips with an adhesive film, incubate for **30 min at 37°C**

Rinse the wells 5 times with prepared 1:20 (1+19) washing solution TWEEN
(300 µl per well)

Add 100 µl of prepared 1:11 (1+10) conjugate solution into all wells
(green)

Cover strips with an adhesive film, incubate for **60 min at 37°C**

Rinse the wells 5 times with prepared 1:20 (1+19) washing solution TWEEN
(300 µl per well)

Add 100 µl of [SOLN|TMB] into all wells

Incubate for **30 min in the dark at 18-25°C**

Add 100 µl of [SOLN|STOP] into all wells
(change of colour from blue to yellow)

Measure the optical density (OD) with an ELISA microplate reader at
450/620-695 nm

CALCULATION OF RESULTS

$$\bar{Nc} = (Nc1 + Nc2 + Nc3)/3;$$

$$CO = \bar{Nc} + 0,3;$$

$$IP_{\text{sample}} = OD_{\text{sample}}/CO$$

\bar{Nc} - the average value of OD 3-x [CONTROL|-]

CO - Cut off

IP_{sample} - sample positivity index

INTERPRETATION OF RESULTS


$IP_{\text{sample}} > 1,1$	POSITIVE
$0,9 \leq IP_{\text{sample}} \leq 1,1$	BORDERLINE
$IP_{\text{sample}} < 0,9$	NEGATIVE

Ascaris lumbricoides IgG

ELISA kit for the qualitative detection of IgG
antibodies to *Ascaris lumbricoides*

Instructions for use

IVD

REF
EI-601

Σ 96
tests

CE

EQUI *Ascaris lumbricoides* IgG

ELISA kit for the qualitative detection of
IgG antibodies to *Ascaris lumbricoides*

1. INTENDED USE

The «EQUI *Ascaris lumbricoides* IgG» is ELISA kit intended to qualitatively detect anti-*Ascaris lumbricoides* IgG in human serum or plasma by enzyme-linked immunosorbent assay (ELISA) in order to diagnose lumbricosis. The testing procedure is designed for both manual arrangement with automatic pipettes and standard equipment, and for automated «open» immunoassay analysers.

Target group: children, rural people, summer visitors.

Usage: ELISA kit is used in clinical diagnostic laboratories and other institutions engaged in *in vitro* diagnostics.

2. CLINICAL SIGNIFICANCE

Ascaris lumbricoides is a human parasite resulting in lumbricosis — one of the most common helminthiases in the world. By some estimates, over a milliard of people infested with acaricides are on earth.

Human ascaris belongs to *Nematoda* roundworms infesting the small intestine of a man who is its exclusive host. *Ascaris lumbricoides* eggs are excreted in the environment with faeces of the infested man. In a warm, wet soil, ascaris larvae develops in the eggs, therefore eggs become invasive only after a maturation period (2 to 3 weeks at 25–30 °C, lower temperatures require longer term). After infestation, larvae leave eggs in the human intestine, penetrates blood circulation and migrate to the liver and lungs with blood flow. The larvae move to the pharynx from the lungs, and here they are re-ingested and further enter the small intestine. In 2 to 3 months, adult ascaris able to propagate develops from larvae in the small intestine.

The helminths are transferred by faecal-oral route upon injection of mature eggs of *Ascaris lumbricoides* with soil-contaminated vegetables, fruits, water, as well as through dirty hands after contact with soil. Lumbricosis is conditionally divided into the early stage (migration of larvae) and late stage (parasitism of adults in the intestine). Invasion is asymptomatic in most cases. Primary feeling of being unwell occurs as early as several days after infestation and is accompanied by weakness, abdominal pain, nausea. Migration of larvae to the lungs may manifest as rales and cough. In some cases, intense invasion may result in pneumonia and liver damage. However, the most common symptom of early lumbricosis are allergic reactions due to hypersensitivity to metabolic products of larvae.

Late stage manifests as decreases appetite, abdominal pain, vomiting, diarrhoea, constipation. Massive ascaris invasion may result in the intestinal obstruction with a lump of helminths or rupture of the walls with peritonitis. When ascarides penetrate other organs, complications may develop such as

hepatitis, cholangitis, pancreatitis and even asphyxia. Cases of neurological disorders sometimes develop in lumbricosis, namely: headache, irritability, sleep impairment, inattention, etc. If no timely treatment is started for intense invasion, it may lead to death, especially in younger children.

Strong immune response to *Ascaris lumbricoides* invasion develops as early as at the early stage. It includes cellular and humoral immunity. Antigens of ascaris larvae stimulate secretion of all-class specific immunoglobulins, however, the level of specific and total IgE antibodies is the highest. The intensity of the immune response (including increased IgG titres) correlates with the massiveness of the invasion.

For diagnosis of lumbricosis, parasitologic stool test for presence of ascaris larvae and eggs is the most common. X-ray imaging of the lungs is additionally applied at the early stage of invasion. Complete blood count (eosinophilia develops in lumbricosis) and detection of serum anti-*Ascaris lumbricoides* antibodies also is included in the set of exams. The presence of specific anti-ascaris antibodies may suggest asymptomatic invasion, and allows initiation of treatment before complications develop in conjunction with other diagnostic instruments.

3. ANALYSIS PRINCIPLE

The procedure of testing for anti-*Ascaris lumbricoides* IgG in «EQUI *Ascaris lumbricoides* IgG» ELISA kit is based on «indirect» solid-phase ELISA with a two-stage incubation. Antigens of *Ascaris lumbricoides* larvae are entrapped in the wells. During the first step of incubation of ELISA plate wells with test samples, specific anti-*Ascaris lumbricoides* antibodies (if present in the samples) bind to the solid-phase antigens. The wells are washed to remove unbound antibodies and have only specific antigen-antibody complexes left. Then, a conjugate of anti-species IgG monoclonal antibodies with horseradish peroxidase is added, which binds to solid-phase immune complexes. Unbound components are removed by washing. Antigen-antibody complexes are detected by adding a solution of chromogen 3,3',5,5'-tetramethylbenzidine (TMB) with hydrogen peroxide. After 30-minute incubation, the reaction is stopped by adding the stop solution. The optical density (OD) in the wells is determined using a spectrophotometer at 450/620-695 nm. The intensity of the yellow colour is proportional to the level of antibodies in the sample.

4. MATERIALS AND EQUIPMENT

4.1. Contents of the ELISA kit

STRIPS

1 x 96
wells

Microplate

Each plate well is coated with *Ascaris lumbricoides* antigen. The wells are detachable. After the first opening, store unused strips in the package at 2-8 °C for a maximum of 6 months

CONTROL +	1 x 0,25 ml	Positive control Conjugated specific monoclonal antibody solution with preservative (pink). Store at 2-8 °C
CONTROL -	1 x 0,6 ml	Negative control Negative human serum with a preservative (yellow). Store at 2-8 °C
DIL SAMPLE	1 x 13 ml	Serum dilution solution Buffer solution with a milk extract, a detergent and a preservative (brown). Store at 2-8 °C
SOLN CONJ	1 x 13 ml	Conjugate solution (ready to use) Buffer solution of monoclonal antibodies to human IgG, conjugated with horseradish peroxidase, with stabilizers and preservative (green). Store at 2-8 °C
SOLN TMB	1 x 13 ml	TMB solution (ready to use) TMB solution, H ₂ O ₂ , a stabilizer, a preservative (colourless). Store at 2-8 °C
TWEEN WASH 20x	1 x 50 ml	Washing solution TWEEN (20x concentrated) 20-fold phosphate buffer concentrate with Tween-20 (colourless). Dilute TWEEN detergent (20x) at 1:20 with distilled or deionized water (e. g., 5 mL of concentrate + 95 mL of water for 8 wells) before use. Store the diluted solution at 2-8 °C for a maximum of 7 days
SOLN STOP	1 x 13 ml	Stop Solution (ready to use) 0.5 mol H ₂ SO ₄ solution (colourless). Store at 2-8 °C

The ELISA kit also includes adhesive films (2 items), sample application plan (1 item), checklist, and instruction for use.

4.2. Optional reagents, materials and equipment

Automatic single and multichannel pipettes 10–1000 µL, tips, volumetric laboratory glassware (10–1,000 mL), deionized or distilled water, thermostat at 37 °C, automatic or semi-automatic plate washer, spectrophotometer (reader) for microplates at 450/620-695 nm, appropriate containers for potentially contaminated waste, timer, filter paper, disposable powder-free gloves, disinfectants.

5. PRECAUTIONS AND SAFETY

5.1. Precautions

Be sure to read the instructions for use carefully before the test. The validity of the test results depends on strict following of the test procedure.

- do not use the ELISA kit components after the expiry date;
- do not use for analysis or mix components of different batches, components of kits for different nosologies, or reagents from other manufacturers with the «EQUI Ascaris lumbricoides IgG» ELISA kit;
- do not freeze the ELISA kit or its contents;
- after using a reagent, close each vial with its cap;

- when washing, control filling and complete aspiration of solution from the wells;
- use a new pipette tip each time you add samples or reagents;
- prevent direct sunlight from reaching the reagents from the ELISA kit;
- [SOLN|TMB] solution must be colourless before use. Do not use the solution if its colour is blue or yellow. Avoid contact of [SOLN|TMB] with metals or metal ions. Use only clean glassware thoroughly rinsed with distilled water;
- do not use reagents with colour not in line with para. 4.1;
- under no circumstances should the same glassware be used for [SOLN|CONJ] and [SOLN|TMB];
- do not evaluate the test results visually (without a reader);
- any optional equipment that is in direct contact with biological material or kit components should be considered contaminated and requires cleaning and decontamination;
- the ELISA kit includes materials for 96 tests. Dispose of the used components as well as any remaining unused components.

5.2. Safety requirements

- all reagents in the ELISA kit are for laboratory professional use for *in vitro* diagnosis only and may only be used by qualified personnel;
- conduct the tests in disposable powder-free gloves and goggles only;
- do not eat, drink, smoke, or apply make-up in the test room;
- do not mouth-pipette the solutions;
- controls from the «EQUI *Ascaris lumbricoides* IgG» ELISA kit have been tested and found to be for anti-HIV1/2, anti-HCV and anti-*Treponema pallidum* antibodies and HBsAg negative; however, controls and test samples should be handled as potentially hazardous infectious materials;
- some of the kit components contain low concentrations of harmful substances and can damage skin or mucosa. In case of contact of [SOLN|TMB], [SOLN|STOP] and [SOLN|CONJ] with mucous membranes or skin, immediately wash the affected area with plenty of water;
- in case of spillage of acid-free solutions, e. g. sera, treat the surface with a disinfectant solution and then wipe dry with filter paper. Otherwise first neutralize acid with sodium bicarbonate solution and then wipe the surface dry as described above.

5.3. Waste inactivation and disposal

- the liquid waste must be inactivated, for example, with hydrogen peroxide solution at the final concentration of 6% for 3 hours at room temperature, or with sodium hypochlorite at the final concentration of 5% for 30 minutes, or with other approved disinfectants;
- the solid waste must be inactivated by autoclaving at a temperature not less than 132°C;

- do not autoclave the solutions that contain sodium azide or sodium hypochlorite;
- disposal of inactivated waste must be conducted due to national laws and regulations.

6. STORAGE AND STABILITY

ELISA kit is stable up to the expiry date stated on the label when stored at 2-8°C. The kit should be transported at 2-8°C. Single transportation at a temperature up to 23°C for two days is possible.

7. SAMPLE COLLECTION, TRANSPORTATION AND STORAGE GUIDELINES

Collect blood from the vein into the sterile test tube. Test tube must be marked with patient ID and date of sample collecting. Blood before serum separation can be stored at 2-8 °C for 24 hours, avoiding freezing.

Serum or plasma can be stored at 2-8 °C for maximum 3 days. Frozen serum can be stored for longer periods of time at -20 °C or -70 °C. Thaw frozen samples and keep them at room temperature for 30 minutes before use. After thawing, the stir samples to achieve homogeneity. Avoid repeated freezing-thawing cycles for test samples. If serum (or plasma) is turbid, remove insoluble inclusions by centrifugation at 3000 rpm for 10-15 minutes. Do not use serum samples with hyperlipidemia, hemolysis, and bacterial growth.

Transport serum samples in insulated containers. To do that, put closed labelled tubes in a plastic bag, tightly seal it and place in the centre of an insulated container. Put the frozen cold packs on the bottom, along the side walls of the insulated container and on top of the serum samples.

8. REAGENT PREPARATION

NOTE! It is very important to keep all ELISA kit components for at least 30 min at room temperature 18-25 °C before the assay!

8.1. Microplate preparation

To prevent water condensation in the wells, keep the **STRIPS** for 30 minutes at a room temperature before opening. Open the vacuum pack, detach the appropriate number of wells, and carefully pack the remaining wells with a desiccant and store tightly zip-locked at 2-8 °C. Storing the packed plate this way ensures its stability for 6 months.

8.2. Washing solution preparation

To prepare detergent, dilute **TWEEN|WASH|20x** at 1:20 (1+19) with distilled or deionized water and stir. E. g., 5 mL of concentrate + 95 mL of water, which is enough for 8 wells. If there are crystals present in the detergent concentrate, heat the vial at 37 °C until the crystals dissolve completely (15–20 minutes). Store the diluted solution at 2-8 °C for a maximum of 7 days.

9. ASSAY PROCEDURE

- 9.1. Prepare the necessary number of wells (four wells for controls and a necessary number of wells for test samples) and insert them into the ELISA plate frame. Be sure to add control wells in every test run.
- 9.2. Fill in the sample application plan.
- 9.3. Prepare the detergent as per para. 8.2.
- 9.4. Add 90 µL of [DIL|SAMPLE] into each plate well.
- 9.5. Add 10 µL of controls and test samples into the wells:
[CONTROL|+] – into well A1,
[CONTROL|-] – into wells B1, C1 and D1,
and test samples into the remaining wells.
At the time of adding, the solution changes its colour from brown to blue. Pipette the mix in the wells carefully to avoid foaming.
- 9.6. Cover the strips up with adhesive film and incubate for **30 minutes at 37 °C**.
- 9.7. Remove and discard the adhesive film and wash all wells 5 times with automatic washer or 8-channel pipette as follows:
 - aspirate the content of all wells into a liquid waste container;
 - add a minimum of 300 µl of diluted washing solution to each well, soak each well for 30 seconds;
 - aspirate the content of all wells again. The residual volume after every aspiration should be less than 5 µl;
 - repeat the washing step 4 more times;
 - after the final aspiration, eliminate extra moisture by tapping the plate against a piece of filter paper.
- 9.8. Add 100 µL of [SOLN|CONJ] into each well. Cover the strips with a new piece of adhesive film and incubate for **30 minutes at 37 °C**.
- 9.9. Following incubation, remove the film carefully and wash the wells five times as described in para. 9.7.
- 9.10. Add 100 µL of [SOLN|TMB] into the wells; do not touch the bottom and the walls of the plate wells.
- 9.11. Incubate the strips for **30 minutes** in a dark place at a room temperature of 18-25 °C. Do not use adhesive film at this stage.
- 9.12. Add 100 µL of [SOLN|STOP] into each strip well to stop the enzymatic reaction; adhere to the same sequence of actions as when adding [SOLN|TMB]. At the time of adding, the solution colour changes from blue to yellow, and clear solution slightly changes its shade.
- 9.13. Measure the optical density (OD) of the wells at 450/620-695 nm wavelength using an ELISA microplate reader within 5 minutes after stopping the reaction. Pay attention to the cleanness of the plate bottom and the absence of bubbles in the wells before reading.

Measurement at the single wavelength of 450 nm is possible, in that case, it is needed to leave one well for blank (only [SOLN|TMB] and [SOLN|STOP] must be added

in blank well).

10. CALCULATION AND INTERPRETATION OF RESULTS

10.1. Calculation of results

Calculate the average OD for the negative control (\bar{Nc}), Cut off (CO) and a sample positivity index (IP_{sample}).

$$\bar{Nc} = (Nc1 + Nc2 + Nc3)/3; \quad CO = \bar{Nc} + 0,3$$

$$IP_{\text{sample}} = OD_{\text{sample}}/CO, \text{ where } OD_{\text{sample}} \text{ is the OD sample.}$$

10.2. Quality control (assay validation)

The test results are considered valid if they meet the following requirements:

$$\boxed{\text{CONTROL}+} \quad OD \geq 1,0$$

$$\boxed{\text{CONTROL}-} \quad OD \leq 0,150$$

$$\boxed{\text{CONTROL}-} \quad \bar{Nc} \times 0,5 \leq Ncn \leq \bar{Nc} \times 2,0 \quad \text{where } Ncn \text{ is the OD for each } Nc \text{ run}$$

If any of the OD values for the negative control is beyond the above interval, it should be discarded, and \bar{Nc} is calculated based on the remaining OD values for the negative control. If several OD values for the negative control fail to meet the above requirements, the test is considered invalid and requires a new run.

10.3. Interpretation of results

$$\begin{array}{ll} IP_{\text{sample}} > 1,1 & \text{POSITIVE} \\ 0,9 \leq IP_{\text{sample}} \leq 1,1 & \text{BORDERLINE*} \\ IP_{\text{sample}} < 0,9 & \text{NEGATIVE} \end{array}$$

* Uncertain samples are recommended to be re-examined in two wells of the ELISA kit. If the results are again uncertain, a new sample should be selected and analyzed in 2-4 weeks. In case of repeated indeterminate results, such samples shall be considered negative.

11. PERFORMANCE CHARACTERISTICS

11.1. Analytical performance characteristics

Precision of measurement

Intra assay repeatability

The coefficient of variation (CV) for three sera with different levels of specific antibodies was evaluated in 24 replicates on one series of ELISA kits.

Sample No.	OD _{av}	IP _{av}	CV, %
547	0,504	1,43	2,9
671	0,753	2,13	3,6
413	1,165	3,30	3,1

Inter assay reproducibility

The coefficient of variation (CV) for three sera with different levels of specific antibodies was evaluated for 4 days in 4 sets of analysis, 8 replicates in each analysis.

Sample No.	OD _{av}	IP _{av}	CV, %
547	0,534	1,55	5,0
671	0,750	2,17	4,6
413	1,159	3,36	3,6

Analytical specificity

The test results are not affected by bilirubin at up to 0.21 mg/mL (361.8 µmol/L), haemoglobin at up to 10 mg/mL and triglycerides at up to 10 mg/mL (11.3 mmol/l) present in the sample.

11.2. Diagnostic characteristics

To evaluate clinical sensitivity and specificity of «EQUI *Ascaris lumbricoides* IgG» ELISA kits, 55 serum samples from patients with clinical symptoms typical for lumbricosis and 60 serum samples from patients without clinical manifestations (seronegative in terms of *Ascaris lumbricoides*) were used. Clinical sensitivity of «EQUI *Ascaris lumbricoides* IgG» ELISA kits was 94.55 % and clinical specificity — 93.3 %.

Method characteristics in comparison with equal commercial ELISA kit was studied in target paediatric population (160 samples) and population of donors (346 samples). For paediatric population serum, relative specificity of «EQUI *Ascaris lumbricoides* IgG» ELISA kits was established at the level of 97.92 % and percent agreement was 95.51 %. For donor population serum, relative specificity of was 89.74 %, relative specificity — 96.30 % and percent agreement was 95.47 %.

12. LIMITATIONS OF ASSAY

Positive result in «EQUI *Ascaris lumbricoides* IgG» ELISA kit supports presence of anti-*Ascaris lumbricoides* specific IgG antibodies. Presence of this class antibodies in newborns is not an evidence of *Ascaris lumbricoides* invasion.

Inconclusive results may suggest a history of *Ascaris lumbricoides* invasion.

Negative result of «EQUI *Ascaris lumbricoides* IgG» ELISA kit supports the absence of anti-*Ascaris lumbricoides* IgG specific antibodies in the test sample or concentration of specific antibodies is below the sensitivity limit of the assay.

The results of serological test only are not the basis for final diagnosis. When establishing the diagnosis, the results of complex laboratory and instrumental tests, as well as clinical manifestations should be considered. Cross-reactions with antibodies to antigens of other helminths cannot be fully ruled out.

13. DIFFICULTIES THAT CAN OCCUR DURING THE ASSAY PROCEDURE

Possible reasons	Solution
<i>High background in all wells</i>	
Contaminated washer	Clean the washer head and rinse according to the instructions for use
Poor quality or contaminated water	Use purified water with specific resistance $\geq 10 \text{ M}\Omega \cdot \text{cm}$
Use of poorly washed glassware	Use chemically clean utensils
Use of chlorinated disinfectants	Do not use chlorine disinfectants
Use of contaminated tips	Use new tips
Increased incubation times or change in the temperature conditions	Adhere to the incubation regime according to the instructions for use
<i>High background in a row of wells</i>	
Repeat application of TMB solution	TMB solution should be applied once
Contamination of the automatic pipette nozzle with conjugate solution	Clean the pipette and dial carefully liquid
Contamination of one of the washer's channel	Clean the flush channel, rinse washer
<i>Received OD of the positive control is below the border value</i>	
One of the reagents (conjugate solution or TMB solution) was not prepared in a correct way or was not added	Re-conduct ELISA, pay attention to the correctness of the introduction of these reagents
Reduced incubation times at any stage	Incubate according to instructions for use
<i>The colour density of the wells fails to meet the obtained optical density value</i>	
This may suggest that the optical beam has been displaced	Check the correct operation of the reader

14. TECHNICAL ASSISTANCE AND CUSTOMER SERVICE

In case of technical problems, you can obtain assistance by contacting the manufacturer.

REFERENCES

1. CDC - Ascariasis // <https://www.cdc.gov/parasites/ascariasis/index.html>.
2. Cooper P. J., Chico M. E. et al. Human Infection with *Ascaris lumbricoides* Is Associated with a Polarized Cytokine Response // *The Journal of Infectious Diseases*. - 2003. - Vol. 182 (4). - P. 1207–1213.
3. Gupta S., Kumar S. et al. *Ascaris lumbricoides*: an unusual aetiology of gastric perforation // *Journal of Surgical Case Reports*. - 2012. - Vol. 2012. - rjs008.
4. Li Q., Zhao D. et al. Life-threatening complications of ascariasis in trauma patients: a review of the literature // *World Journal of Emergency Medicine*. - 2014. - Vol. 5 (3). - P. 165–170.
5. McSharry C., Xia Y. et al. Natural Immunity to *Ascaris lumbricoides* Associated with Immunoglobulin E Antibody to ABA-1 Allergen and Inflammation Indicators in Children // *Infection and Immunity*. - 1999. - Vol. 67(2). - P. 484–489.
6. Palmer L. J., Celedón J. C. et al. *Ascaris lumbricoides* Infection Is Associated with Increased Risk of Childhood Asthma and Atopy in Rural China // *American Journal of Respiratory and Critical Care Medicine*. - 2002. - Vol. 165, No. 11. - P. 1489–1493.
7. Shalaby N. Effect of *Ascaris lumbricoides* infection on T helper cell type 2 in rural Egyptian children // *Therapeutics and Clinical Risk Management*. - 2016. - Vol. 12. - P. 379–385.
8. WHO. Water related diseases: ascariasis. 2013 // http://www.who.int/water_sanitation_health/diseases/ascariasis/en/.
9. Regulation (EU) 2017/746 of the European Parliament and of the Council of 5 April 2017 on in vitro diagnostic medical devices and repealing Directive 98/79/EC and Commission Decision 2010/227/EU.
10. Закон України «Про відходи» // Відомості Верховної Ради України. - 1998. - №36-37.
11. Наказ МОЗ України №325 від 08.06.2015 «Про затвердження Державних санітарно-протиепідемічних правил і норм щодо поводження з медичними відходами».
12. Постанова КМУ від 02 жовтня 2013р. №754 «Про затвердження технічного регламенту щодо медичних виробів для діагностики in vitro».
13. Hanna Tolonen, Kari Kuulasmaa, Tiina Laatikainen, Hermann Wolf and the European Health Risk Monitoring Project. Recommendation for indicators, international collaboration, protocol and manual of operations for chronic disease risk factor surveys Part 4.Storage and transfer of serum/plasma samples// Finnish National Public Health Institute 2002// https://thl.fi/publications/ehrm/product2/part_iii4.htm
14. Surveillance Guidelines for Measles, Rubella and Congenital Rubella Syndrome in the WHO European Region. Annex 3.Collection, storage and shipment of specimens for laboratory diagnosis and interpretation of results// Geneva: World Health Organization; 2012 Dec.

Manufacturer

Authorized Representative in the European Community

In vitro diagnostic medical device

Catalogue number

Date of manufacture

Use by date

Batch code

Temperature limit

Contains sufficient for <n> tests

Caution

Non-Sterile

Consult instructions for use

Keep away from sunlight

Keep dry

Compliance with EU safety requirements

Edition 8, 10.02.2022

For questions and suggestions regarding the ELISA kit contact:

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ASSAY PROCEDURE SCHEME

Keep all reagents for 30 min at temperature 18-25°C before use

Dispense 90 µl [DIL|SAMPLE] into the wells (brown)

Add to 10 µl of controls and samples into the wells:
A1 – [CONTROL|+], B1, C1, D1 – [CONTROL|-],
other wells – examined samples
(change of colour from brown to blue)

Cover strips with an adhesive film, incubate for **30 min at 37°C**

Rinse the wells 5 times with prepared 1:20 (1+19) washing solution TWEEN (300 µl per well)

Add 100 µl of [SOLN|CONJ] into all wells (green)

Cover strips with an adhesive film, incubate for **30 min at 37°C**

Rinse the wells 5 times with prepared 1:20 (1+19) washing solution TWEEN (300 µl per well)

Add 100 µl of [SOLN|TMB] into all wells

Incubate for **30 min** in the dark at **18-25°C**

Add 100 µl of [SOLN|STOP] into all wells
(change of colour from blue to yellow)

Measure the optical density (OD) with an ELISA microplate reader at 450/620-695 nm

CALCULATION OF RESULTS

$$\bar{Nc} = (Nc1 + Nc2 + Nc3)/3;$$

$$CO = \bar{Nc} + 0,3;$$

$$IP_{\text{sample}} = OD_{\text{sample}} / CO$$

\bar{Nc} - the average value of OD 3-x [CONTROL|-]

CO - Cut off

IP_{sample} - sample positivity index

INTERPRETATION OF RESULTS


$IP_{\text{sample}} > 1,1$	POSITIVE
$0,9 \leq IP_{\text{sample}} \leq 1,1$	BORDERLINE
$IP_{\text{sample}} < 0,9$	NEGATIVE

Toxocara canis IgG

ELISA kit for the qualitative detection of IgG
antibodies to *Toxocara canis*

Instructions for use

IVD

REF
EI-603

Σ 96
tests

CE

EQUI *Toxocara canis* IgG

ELISA kit for the qualitative detection of
IgG antibodies to *Toxocara canis*

1. INTENDED USE

The «EQUI *Toxocara canis* IgG» is ELISA kit intended to qualitatively detect anti-*Toxocara canis* IgG in human serum or plasma by enzyme-linked immunosorbent assay (ELISA) in order to diagnose toxocariasis. The testing procedure is designed for both manual arrangement with automatic pipettes and standard equipment, and for automated «open» immunoassay analysers.

Target group: children, pet owners, rural people, summer visitors, forest guards, veterinarians.

Usage: ELISA kit is used in clinical diagnostic laboratories and other institutions engaged in *in vitro* diagnostics.

2. CLINICAL SIGNIFICANCE

Toxocariasis is a common disease induced by *Toxocara* helminth which is transmitted from animals to human. Toxocariasis is spread throughout the world, however, it is more common in depressed areas with poor hygienic conditions. In some regions, up to 90 % of puppies and up to 10 % of adult domesticated dogs are infested with toxocara. The risk of infestation is higher for owners of cats and dogs and for children due to playing in the sandpits and on the playgrounds contaminated with animal faeces.

Toxocara are threadworms belonging to *Nematoda*. Human conditions are mostly caused by *Toxocara canis*, which infested canids, rare - *Toxocara cati*, which is more common in felids. Adult toxocara in the body of infested animals reaches 5–15 cm in length; their propagation takes place here. Female helminths lay about 200 thous eggs daily, which are excreted in the environment with faeces. If conditions are favourable, following several weeks of maturation in the soil they become invasive — a larva is developed in the eggs. In the paratenic host (mice, poultry, cows, pigs, etc.). larva develops without propagation. If the conditions are unfavourable, larvae are encapsulated and may maintain viability for a long time (up to 10 years). They may also be the source of invasion.

People are infested through faecal-oral route when ingesting *Toxocara canis* mature eggs with soil-contaminated vegetables, fruits, berries, via dirty hands or when consuming meat of paratenic hosts. In the small intestine, larvae leave their cover and penetrates blood circulation through the intestinal walls. The larvae migrate to other organs and tissues with blood, namely: liver, lungs, muscles, eyes, CNS, etc. In the most of the infested, toxocariasis is asymptomatic. Clinical manifestations of this disease are associated with the site of larvae migration and depend on the intensity of invasion and age of the host. Visceral syndrome larva migrans is typical after infestation of the internal organs with *Toxocara canis* and ocular

toxocariasis, when eye and optic nerve are involved. Symptoms of visceral toxocariasis: fever, fatigue, abdominal pain, anorexia, hepatomegaly, cough and others. Heart and respiratory failure may develop in severe cases. Due to a strong immune response to larvae antigens, immediate and delayed hypersensitivity reactions develop. Granulomatosis in ocular toxocariasis may result in retinal detachment and loss of vision.

Diagnosis of toxocariasis is complicated due to the lack of specific manifestations of the disease, even upon intense invasion. Furthermore, a man is an intermediate host of *Toxocara canis* and does not excrete parasites in the environment, whereas it is difficult to localise larvae in certain organs via non-invasive methods. Eosinophilia may appear in blood tests, however, serological tests are more common to detect toxocariasis (immunofluorescence reaction, ELISA and immunoblotting). Detection of specific anti-*Toxocara canis* IgG to larvae antigens may suggest current or previous invasion. High titer of IgE antibodies is also typical for active invasion. However, the combination of clinical manifestations and laboratory findings are necessary for diagnosis.

3. ANALYSIS PRINCIPLE

The procedure of testing for anti-*Toxocara canis* IgG in «EQUI *Toxocara canis* IgG» ELISA kit is based on «indirect» solid-phase ELISA with a two-stage incubation. Antigens of *Toxocara canis* larvae are entrapped in the wells. During the first step of incubation of ELISA plate wells with test samples, specific anti-*Toxocara canis* antibodies (if present in the samples) bind to the solid-phase antigens. The wells are washed to remove unbound antibodies and have only specific antigen-antibody complexes left. Then, a conjugate of anti-species IgG monoclonal antibodies with horseradish peroxidase is added, which binds to solid-phase immune complexes. Unbound components are removed by washing. Antigen-antibody complexes are detected by adding a solution of chromogen 3,3',5,5'-tetramethylbenzidine (TMB) with hydrogen peroxide. After 30-minute incubation, the reaction is stopped by adding the stop solution. The optical density (OD) in the wells is determined using a spectrophotometer at 450/620-695 nm. The intensity of the yellow colour is proportional to the level of antibodies in the sample.

4. MATERIALS AND EQUIPMENT

4.1. Contents of the ELISA kit

STRIPS

1 x 96 wells

Microplate

Each plate well is coated with *Toxocara canis* larval antigens. The wells are detachable. After the first opening, store unused strips in the package at 2-8 °C for a maximum of 6 months

CONTROL +	1 x 0,25 ml	Positive control Conjugated specific monoclonal antibody solution with preservative (pink). Store at 2-8 °C
CONTROL -	1 x 0,6 ml	Negative control Negative human serum with a preservative (yellow). Store at 2-8 °C
DIL SAMPLE	1 x 13 ml	Serum dilution solution Buffer solution with a milk extract, a detergent and a preservative (brown). Store at 2-8 °C
SOLN CONJ	1 x 13 ml	Conjugate solution (ready to use) Buffer solution of monoclonal antibodies to human IgG, conjugated with horseradish peroxidase, with stabilizers and preservative (green). Store at 2-8 °C
SOLN TMB	1 x 13 ml	TMB solution (ready to use) TMB solution, H ₂ O ₂ , a stabilizer, a preservative (colourless). Store at 2-8 °C
TWEEN WASH 20x	1 x 50 ml	Washing solution TWEEN (20x concentrated) 20-fold phosphate buffer concentrate with Tween-20 (colourless). Dilute TWEEN detergent (20x) at 1:20 with distilled or deionized water (e. g., 5 mL of concentrate + 95 mL of water for 8 wells) before use. Store the diluted solution at 2-8 °C for a maximum of 7 days
SOLN STOP	1 x 13 ml	Stop Solution (ready to use) 0.5 mol H ₂ SO ₄ solution (colourless). Store at 2-8 °C

The ELISA kit also includes adhesive films (2 items), sample application plan (1 item), checklist, and instruction for use.

4.2. Optional reagents, materials and equipment

Automatic single and multichannel pipettes 10–1000 µL, tips, volumetric laboratory glassware (10–1,000 mL), deionized or distilled water, thermostat at 37 °C, automatic or semi-automatic plate washer, spectrophotometer (reader) for microplates at 450/620–695 nm, appropriate containers for potentially contaminated waste, timer, filter paper, disposable powder-free gloves, disinfectants.

5. PRECAUTIONS AND SAFETY

5.1. Precautions

Be sure to read the instructions for use carefully before the test. The validity of the test results depends on strict following of the test procedure.

- do not use the ELISA kit components after the expiry date;
- do not use for analysis or mix components of different batches, components of kits for different nosologies, or reagents from other manufacturers with the «EQUI Toxocara canis IgG» ELISA kit;

- do not freeze the ELISA kit or its contents;
- after using a reagent, close each vial with its cap;
- when washing, control filling and complete aspiration of solution from the wells;
- use a new pipette tip each time you add samples or reagents;
- prevent direct sunlight from reaching the reagents from the ELISA kit;
- [SOLN|TMB] solution must be colourless before use. Do not use the solution if its colour is blue or yellow. Avoid contact of [SOLN|TMB] with metals or metal ions. Use only clean glassware thoroughly rinsed with distilled water;
- do not use reagents with colour not in line with para. 4.1;
- under no circumstances should the same glassware be used for [SOLN|CONJ] and [SOLN|TMB];
- do not evaluate the test results visually (without a reader);
- any optional equipment that is in direct contact with biological material or kit components should be considered contaminated and requires cleaning and decontamination;
- the ELISA kit includes materials for 96 tests. Dispose of the used components as well as any remaining unused components.

5.2. Safety requirements

- all reagents in the ELISA kit are for laboratory professional use for *in vitro* diagnosis only and may only be used by qualified personnel;
- conduct the tests in disposable powder-free gloves and goggles only;
- do not eat, drink, smoke, or apply make-up in the test room;
- do not mouth-pipette the solutions;
- controls from the «EQUI *Toxocara canis* IgG» ELISA kit have been tested and found to be for anti-HIV1/2, anti-HCV and anti-*Treponema pallidum* antibodies and HBsAg negative; however, controls and test samples should be handled as potentially hazardous infectious materials;
- some of the kit components contain low concentrations of harmful substances and can damage skin or mucosa. In case of contact of [SOLN|TMB], [SOLN|STOP] and [SOLN|CONJ] with mucous membranes or skin, immediately wash the affected area with plenty of water;
- in case of spillage of acid-free solutions, e. g. sera, treat the surface with a disinfectant solution and then wipe dry with filter paper. Otherwise first neutralize acid with sodium bicarbonate solution and then wipe the surface dry as described above.

5.3. Waste inactivation and disposal

- the liquid waste must be inactivated, for example, with hydrogen peroxide solution at the final concentration of 6% for 3 hours at room temperature, or with sodium hypochlorite at the final concentration of 5% for 30 minutes, or with other approved disinfectants;

- the solid waste must be inactivated by autoclaving at a temperature not less than 132°C;
- do not autoclave the solutions that contain sodium azide or sodium hypochlorite;
- disposal of inactivated waste must be conducted due to national laws and regulations.

6. STORAGE AND STABILITY

ELISA kit is stable up to the expiry date stated on the label when stored at 2-8°C. The kit should be transported at 2-8°C. Single transportation at a temperature up to 23°C for two days is possible.

7. SAMPLE COLLECTION, TRANSPORTATION AND STORAGE GUIDELINES

Collect blood from the vein into the sterile test tube. Test tube must be marked with patient ID and date of sample collecting. Blood before serum separation can be stored at 2-8 °C for 24 hours, avoiding freezing.

Serum or plasma can be stored at 2-8 °C for maximum 3 days. Frozen serum can be stored for longer periods of time at -20 °C or -70 °C. Thaw frozen samples and keep them at room temperature for 30 minutes before use. After thawing, the stir samples to achieve homogeneity. Avoid repeated freezing-thawing cycles for test samples. If serum (or plasma) is turbid, remove insoluble inclusions by centrifugation at 3000 rpm for 10-15 minutes. Do not use serum samples with hyperlipidemia, hemolysis, and bacterial growth.

Transport serum samples in insulated containers. To do that, put closed labelled tubes in a plastic bag, tightly seal it and place in the centre of an insulated container. Put the frozen cold packs on the bottom, along the side walls of the insulated container and on top of the serum samples.

8. REAGENT PREPARATION

NOTE! It is very important to keep all ELISA kit components for at least 30 min at room temperature 18-25 °C before the assay!

8.1. Microplate preparation

To prevent water condensation in the wells, keep the **STRIPS** for 30 minutes at a room temperature before opening. Open the vacuum pack, detach the appropriate number of wells, and carefully pack the remaining wells with a desiccant and store tightly zip-locked at 2-8 °C. Storing the packed plate this way ensures its stability for 6 months.

8.2. Washing solution preparation

To prepare detergent, dilute **TWEEN|WASH|20x** at 1:20 (1+19) with distilled or deionized water and stir. E. g., 5 mL of concentrate + 95 mL of water, which is enough for 8 wells. If there are crystals present in the detergent concentrate, heat the vial at 37 °C until the crystals dissolve completely (15–20 minutes). Store the diluted solution at 2-8 °C for a maximum of 7 days.

9. ASSAY PROCEDURE

- 9.1. Prepare the necessary number of wells (four wells for controls and a necessary number of wells for test samples) and insert them into the ELISA plate frame. Be sure to add control wells in every test run.
- 9.2. Fill in the sample application plan.
- 9.3. Prepare the detergent as per para. 8.2.
- 9.4. Add 90 µL of [DIL|SAMPLE] into each plate well.
- 9.5. Add 10 µL of controls and test samples into the wells:
[CONTROL|+] – into well A1,
[CONTROL|-] – into wells B1, C1 and D1,
and test samples into the remaining wells.
At the time of adding, the solution changes its colour from brown to blue. Pipette the mix in the wells carefully to avoid foaming.
- 9.6. Cover the strips up with adhesive film and incubate for **30 minutes at 37 °C**.
- 9.7. Remove and discard the adhesive film and wash all wells 5 times with automatic washer or 8-channel pipette as follows:
 - aspirate the content of all wells into a liquid waste container;
 - add a minimum of 300 µl of diluted washing solution to each well, soak each well for 30 seconds;
 - aspirate the content of all wells again. The residual volume after every aspiration should be less than 5 µl;
 - repeat the washing step 4 more times;
 - after the final aspiration, eliminate extra moisture by tapping the plate against a piece of filter paper.
- 9.8. Add 100 µL of [SOLN|CONJ] into each well. Cover the strips with a new piece of adhesive film and incubate for **30 minutes at 37 °C**.
- 9.9. Following incubation, remove the film carefully and wash the wells five times as described in para. 9.7.
- 9.10. Add 100 µL of [SOLN|TMB] into the wells; do not touch the bottom and the walls of the plate wells.
- 9.11. Incubate the strips for **30 minutes** in a dark place at a room temperature of 18-25 °C. Do not use adhesive film at this stage.
- 9.12. Add 100 µL of [SOLN|STOP] into each strip well to stop the enzymatic reaction; adhere to the same sequence of actions as when adding [SOLN|TMB]. At the time of adding, the solution colour changes from blue to yellow, and clear solution slightly changes its shade.
- 9.13. Measure the optical density (OD) of the wells at 450/620-695 nm wavelength using an ELISA microplate reader within 5 minutes after stopping the reaction. Pay attention to the cleanness of the plate bottom and the absence of bubbles in the wells before reading.

Measurement at the single wavelength of 450 nm is possible, in that case, it is needed to leave one well for blank (only [SOLN|TMB] and [SOLN|STOP] must be added

in blank well).

10. CALCULATION AND INTERPRETATION OF RESULTS

10.1. Calculation of results

Calculate the average OD for the negative control (\bar{Nc}), Cut off (CO) and a sample positivity index (IP_{sample}).

$$\bar{Nc} = (Nc1 + Nc2 + Nc3)/3; \quad CO = \bar{Nc} + 0,3$$

$$IP_{\text{sample}} = OD_{\text{sample}}/CO, \text{ where: } OD_{\text{sample}} \text{ is the OD sample.}$$

10.2. Quality control (assay validation)

The test results are considered valid if they meet the following requirements:

CONTROL +	$OD \geq 1,0$
-----------	---------------

CONTROL -	$OD \leq 0,150$
-----------	-----------------

CONTROL -	$\bar{Nc} \times 0,5 \leq Ncn \leq \bar{Nc} \times 2,0$	where Ncn is the OD for each Nc run
-----------	---	---

If any of the OD values for the negative control is beyond the above interval, it should be discarded, and Nc is calculated based on the remaining OD values for the negative control. If several OD values for the negative control fail to meet the above requirements, the test is considered invalid and requires a new run.

10.3. Interpretation of results

$IP_{\text{sample}} > 1,1$	POSITIVE
$0,9 \leq IP_{\text{sample}} \leq 1,1$	BORDERLINE*
$IP_{\text{sample}} < 0,9$	NEGATIVE

* Uncertain samples are recommended to be re-examined in two wells of the ELISA kit. If the results are again uncertain, a new sample should be selected and analyzed in 2-4 weeks. In case of repeated indeterminate results, such samples shall be considered negative.

11. PERFORMANCE CHARACTERISTICS

11.1. Analytical performance characteristics

Precision of measurement

Intra assay repeatability

The coefficient of variation (CV) for three sera with different levels of specific antibodies was evaluated in 24 replicates on one series of ELISA kits.

Sample No.	OD_{av}	IP_{av}	CV, %
669	0,927	2,81	4,8
544	1,503	4,56	1,4
666	1,694	5,14	4,5

Inter assay reproducibility

The coefficient of variation (CV) for three sera with different levels of specific antibodies was evaluated for 4 days in 4 sets of analysis, 8 replicates in each analysis.

Sample No.	OD _{av}	IP _{av}	CV, %
669	1,016	3,04	4,7
544	1,516	4,54	1,9
666	1,683	5,04	4,1

Analytical specificity

The test results are not affected by bilirubin at up to 0.21 mg/mL (361.8 µmol/L), haemoglobin at up to 10 mg/mL and triglycerides at up to 10 mg/mL (11.3 mmol/l) present in the sample.

11.2. Diagnostic characteristics

To evaluate diagnostic characteristics of «EQUI *Toxocara canis* IgG» ELISA kits, 78 serum samples from patients with clinical symptoms typical for toxocariasis and 60 serum samples from patients without clinical manifestations (seronegative in terms of *Toxocara canis*) were used. Clinical sensitivity of «EQUI *Toxocara canis* IgG» ELISA kits was 98.7 %, clinical specificity — 96.7 %.

Method characteristics in comparison with equal commercial ELISA kit was studied in target paediatric population (160 samples) and population of donors (298 samples). For paediatric population serum, relative specificity of «EQUI *Toxocara canis* IgG» ELISA kits was established at the level of 99.28 % and percent agreement was 97.45 %. For donor population serum, relative specificity of was 89.19 %, relative specificity — 93.55 % and percent agreement was 91.73 %.

12. LIMITATIONS OF ASSAY

Positive result in «EQUI *Toxocara canis* IgG» ELISA kit supports presence of anti-*Toxocara canis* specific IgG antibodies. Presence of this class antibodies in newborns is not an evidence of *Toxocara canis* invasion.

Inconclusive results may suggest a history of *Toxocara canis* invasion.

Negative result of «EQUI *Toxocara canis* IgG» ELISA kit supports the absence of anti-*Toxocara canis* specific IgG antibodies in the test sample or concentration of specific antibodies is below the sensitivity limit of the assay.

The results of serological test only are not the basis for final diagnosis. When establishing the diagnosis, the results of complex laboratory and instrumental tests, as well as clinical manifestations should be considered. Cross-reactions with antibodies to antigens of other helminths cannot be fully ruled out.

13. DIFFICULTIES THAT CAN OCCUR DURING THE ASSAY PROCEDURE

Possible reasons	Solution
<i>High background in all wells</i>	
Contaminated washer	Clean the washer head and rinse according to the instructions for use
Poor quality or contaminated water	Use purified water with specific resistance $\geq 10 \text{ M}\Omega \cdot \text{cm}$
Use of poorly washed glassware	Use chemically clean utensils
Use of chlorinated disinfectants	Do not use chlorine disinfectants
Use of contaminated tips	Use new tips
Increased incubation times or change in the temperature conditions	Adhere to the incubation regime according to the instructions for use
<i>High background in a row of wells</i>	
Repeat application of TMB solution	TMB solution should be applied once
Contamination of the automatic pipette nozzle with conjugate solution	Clean the pipette and dial carefully liquid
Contamination of one of the washer's channel	Clean the flush channel, rinse washer
<i>Received OD of the positive control is below the border value</i>	
One of the reagents (conjugate solution or TMB solution) was not prepared in a correct way or was not added	Re-conduct ELISA, pay attention to the correctness of the introduction of these reagents
Reduced incubation times at any stage	Incubate according to instructions for use
<i>The colour density of the wells fails to meet the obtained optical density value</i>	
This may suggest that the optical beam has been displaced	Check the correct operation of the reader

14. TECHNICAL ASSISTANCE AND CUSTOMER SERVICE

In case of technical problems, you can obtain assistance by contacting the manufacturer.

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ASSAY PROCEDURE SCHEME

Keep all reagents for 30 min at temperature 18-25°C before use

Dispense 90 µl [DIL|SAMPLE] into the wells
(brown)

Add to 10 µl of controls and samples into the wells:
A1 – [CONTROL|+], B1, C1, D1 – [CONTROL|-],
other wells – examined samples
(change of colour from brown to blue)

Cover strips with an adhesive film, incubate for **30 min at 37°C**

Rinse the wells 5 times with prepared 1:20 (1+19) washing solution TWEEN
(300 µl per well)

Add 100 µl of [SOLN|CONJ] into all wells
(green)

Cover strips with an adhesive film, incubate for **30 min at 37°C**

Rinse the wells 5 times with prepared 1:20 (1+19) washing solution TWEEN
(300 µl per well)

Add 100 µl of [SOLN|TMB] into all wells

Incubate for **30 min** in the dark at **18-25°C**

Add 100 µl of [SOLN|STOP] into all wells
(change of colour from blue to yellow)

Measure the optical density (OD) with an ELISA microplate reader at
450/620-695 nm

CALCULATION OF RESULTS

$$\bar{N}_c = (N_{c1} + N_{c2} + N_{c3})/3;$$

$$CO = \bar{N}_c + 0,3;$$

$$IP_{\text{sample}} = OD_{\text{sample}}/CO$$

\bar{N}_c - the average value of OD 3-x [CONTROL|-]

CO - Cut off

IP_{sample} - sample positivity index

INTERPRETATION OF RESULTS


$IP_{\text{sample}} > 1,1$	POSITIVE
$0,9 \leq IP_{\text{sample}} \leq 1,1$	BORDERLINE
$IP_{\text{sample}} < 0,9$	NEGATIVE

anti-*Trichinella* *spiralis*

ELISA kit for the qualitative detection of antibodies to
Trichinella spiralis

Instructions for use

IVD

REF
EI-605

Σ 96
tests

CE

EQUI anti-*Trichinella spiralis*

ELISA kit for the qualitative detection of antibodies to *Trichinella spiralis*

1. INTENDED USE

The «EQUI anti-*Trichinella spiralis*» is ELISA kit intended to qualitatively detect anti-*Trichinella spiralis* antibodies in human serum or plasma by enzyme-linked immunosorbent assay (ELISA) in order to diagnose trichinosis. The testing procedure is designed for both manual arrangement with automatic pipettes and standard equipment, and for automated «open» immunoassay analysers.

Target group: villagers, hunters and their families, people whose diet includes meat, particularly pork.

Usage: ELISA kit is used in clinical diagnostic laboratories and other institutions engaged in *in vitro* diagnostics.

2. CLINICAL SIGNIFICANCE

Among helminthiases, which are often caused by meat consumption, trichinosis has a special place. The causative agent of this disease is the *Trichinella spiralis* roundworm, the larvae of which were found even in the mummies of Ancient Egypt. The greatest danger of this helminthiasis is the possibility of severe damage to the central nervous system, which leads to death.

Adult *Trichinella spiralis* have a filamentous body up to 0.5 mm in diameter and 2 mm in length. *Trichinella* is characterized by live birth – while in the human body, one female can give birth to about one and a half thousand rod-shaped larvae. The larvae actively migrate to other organs, but only settle massively in the striated muscles. Connective tissue capsule with feeding vessels gets formed around the larva after 1-1.5 months. About a year later, the capsule begins to calcify, but the larva inside remains alive and invasive for many years. Most mammals can become infected with trichinosis in the natural environment by eating fresh or decomposed meat that contains *Trichinella* larvae.

Trichinella enters the human body through the meat infested with encapsulated larvae (pork, horse meat, nutria's meat, wild boar, bear, etc.). Acute (intestinal) trichinosis may be asymptomatic or show symptoms 1-2 days after eating raw or undercooked meat. This stage of the disease is characterized by abdominal pain, nausea or diarrhea. However, the classic symptoms of trichinosis appear at the stage of larval migration about 2 weeks after the invasion and can last up to two months. These include fever, fatigue, muscle aches and pain, swelling of the face and eyes, skin rashes, and more. *Trichinella spiralis* affects the most active human muscles: diaphragm, intercostal and abdominal muscles, tongue, as well as facial and skeletal muscles. The overall extent of clinical manifestations is directly correlated with the quantity of consumed parasites.

During larval migration, the organism is hyper-sensitized by highly immunogenic products of *Trichinella* metabolism. Patients have eosinophilia and an increase in IgE antibodies. However, a strong immune response does not lead to the elimination of parasites, although it may be crucial for protection from re-infection.

Thus, recurring cases of trichinosis are much easier to tolerate than the primary ones.

To diagnose this helminthiasis, specific tests are conducted on meat that could be a source of invasion, to test for the presence of encapsulated larvae of *Trichinella spiralis*. During the early stages of infection, PCR is used to detect *Trichinella* DNA. A few weeks after the probable invasion, a biopsy of the affected patient's muscle tissue may be performed to detect the parasite. But from the second week of the invasion, serodiagnosis with the detection of specific antibodies to antigens of larvae of the pathogen trichinosis and an increase in the titer of these antibodies in paired sera is particularly useful. 4-5 months after the invasion, the number of specific antibodies decreases, but may remain at a level sufficient for detection for several years. Immunoblotting is used to confirm the positive result of the detection of specific IgG antibodies. However, for a comprehensive diagnosis of trichinosis we should use not only laboratory-generated data but also clinical and epidemiological studies.

3. ANALYSIS PRINCIPLE

The procedure of testing for anti-*Trichinella spiralis* antibodies in «EQUI anti-*Trichinella spiralis*» ELISA kit is based on «indirect» solid-phase ELISA with a two-stage incubation. Antigens of *Trichinella spiralis* larvae are entrapped in the wells. During the first step of incubation of ELISA plate wells with test samples, specific anti-*Trichinella spiralis* antibodies (if present in the samples) bind to the solid-phase antigens. The wells are washed to remove unbound antibodies and have only specific antigen-antibody complexes left. Then, a conjugate of anti-species (anti-IgG and anti-IgA) monoclonal antibodies with horseradish peroxidase is added, which binds to solid-phase immune complexes. Unbound components are removed by washing. Antigen-antibody complexes are detected by adding a solution of chromogen 3,3',5,5'-tetramethylbenzidine (TMB) with hydrogen peroxide. After 30-minute incubation, the reaction is stopped by adding the stop solution. The optical density (OD) in the wells is determined using a spectrophotometer at 450/620-695 nm. The intensity of the yellow colour is proportional to the level of antibodies in the sample.

4. MATERIALS AND EQUIPMENT

4.1. Contents of the ELISA kit

STRIPS	1 x 96 wells	Microplate Each plate well is coated with <i>Trichinella spiralis</i> larvae antigens. The wells are detachable. After the first opening, store unused strips in the package at 2-8 °C for a maximum of 6 months
CONTROL +	1 x 0,25 ml	Positive control Conjugated specific monoclonal antibody solution with preservative (pink). Store at 2-8 °C

CONTROL -	1 x 0,6 ml	Negative control Negative human serum with a preservative (yellow). Store at 2-8 °C
DIL SAMPLE	1 x 13 ml	Serum dilution solution Buffer solution with a milk extract, a detergent and a preservative (purple). Store at 2-8 °C
SOLN CONJ	1 x 13 ml	Conjugate solution (ready to use) Buffer solution of monoclonal antibodies to human IgG and IgA, conjugated with horseradish peroxidase, with stabilizers and preservative (green). Store at 2-8 °C
SOLN TMB	1 x 13 ml	TMB solution (ready to use) TMB solution, H ₂ O ₂ , a stabilizer, a preservative (colourless). Store at 2-8 °C
TWEEN WASH 20x	1 x 50 ml	Washing solution TWEEN (20x concentrated) 20-fold phosphate buffer concentrate with Tween-20 (colourless). Dilute TWEEN detergent (20x) at 1:20 with distilled or deionized water (e. g., 5 mL of concentrate + 95 mL of water for 8 wells) before use. Store the diluted solution at 2-8 °C for a maximum of 7 days
SOLN STOP	1 x 13 ml	Stop Solution (ready to use) 0.5 mol H ₂ SO ₄ solution (colourless). Store at 2-8 °C

The ELISA kit also includes adhesive films (2 items), sample application plan (1 item), checklist, and instruction for use.

4.2. Optional reagents, materials and equipment

Automatic single and multichannel pipettes 10–1000 µL, tips, volumetric laboratory glassware (10–1,000 mL), deionized or distilled water, thermostat at 37 °C, automatic or semi-automatic plate washer, spectrophotometer (reader) for microplates at 450/620-695 nm, appropriate containers for potentially contaminated waste, timer, filter paper, disposable powder-free gloves, disinfectants.

5. PRECAUTIONS AND SAFETY

5.1. Precautions

Be sure to read the instructions for use carefully before the test. The validity of the test results depends on strict following of the test procedure.

- do not use the ELISA kit components after the expiry date;
- do not use for analysis or mix components of different batches, components of kits for different nosologies, or reagents from other manufacturers with the «EQUI anti-Trichinella spiralis» ELISA kit;
- do not freeze the ELISA kit or its contents;
- after using a reagent, close each vial with its cap;
- when washing, control filling and complete aspiration of solution from the wells;
- use a new pipette tip each time you add samples or reagents;
- prevent direct sunlight from reaching the reagents from the ELISA kit;

- [SOLN|TMB] solution must be colourless before use. Do not use the solution if its colour is blue or yellow. Avoid contact of [SOLN|TMB] with metals or metal ions. Use only clean glassware thoroughly rinsed with distilled water;
- do not use reagents with colour not in line with para. 4.1;
- under no circumstances should the same glassware be used for [SOLN|CONJ] and [SOLN|TMB];
- do not evaluate the test results visually (without a reader);
- any optional equipment that is in direct contact with biological material or kit components should be considered contaminated and requires cleaning and decontamination;
- the ELISA kit includes materials for 96 tests. Dispose of the used components as well as any remaining unused components.

5.2. Safety requirements

- all reagents in the ELISA kit are for laboratory professional use for *in vitro* diagnosis only and may only be used by qualified personnel;
- conduct the tests in disposable powder-free gloves and goggles only;
- do not eat, drink, smoke, or apply make-up in the test room;
- do not mouth-pipette the solutions;
- controls from the «EQUI anti-Trichinella spiralis» kit have been tested and found to be for anti-HIV1/2, anti-HCV and anti-*Treponema pallidum* antibodies and HBsAg negative; however, controls and test samples should be handled as potentially hazardous infectious materials;
- some of the kit components contain low concentrations of harmful substances and can damage skin or mucosa. In case of contact of [SOLN|TMB], [SOLN|STOP] and [SOLN|CONJ] with mucous membranes or skin, immediately wash the affected area with plenty of water;
- in case of spillage of acid-free solutions, e. g. sera, treat the surface with a disinfectant solution and then wipe dry with filter paper. Otherwise first neutralize acid with sodium bicarbonate solution and then wipe the surface dry as described above.

5.3. Waste inactivation and disposal

- the liquid waste must be inactivated, for example, with hydrogen peroxide solution at the final concentration of 6% for 3 hours at room temperature, or with sodium hypochlorite at the final concentration of 5% for 30 minutes, or with other approved disinfectants;
- the solid waste must be inactivated by autoclaving at a temperature not less than 132°C;
- do not autoclave the solutions that contain sodium azide or sodium hypochlorite;
- disposal of inactivated waste must be conducted due to national laws and regulations.

6. STORAGE AND STABILITY

ELISA kit is stable up to the expiry date stated on the label when stored at 2-8°C. The kit should be transported at 2-8°C. Single transportation at a temperature up to 23°C for two days is possible.

7. SAMPLE COLLECTION, TRANSPORTATION AND STORAGE GUIDELINES

Collect blood from the vein into the sterile test tube. Test tube must be marked with patient ID and date of sample collecting. Blood before serum separation can be stored at 2-8 °C for 24 hours, avoiding freezing.

Serum or plasma can be stored at 2-8 °C for maximum 3 days. Frozen serum can be stored for longer periods of time at -20 °C or -70 °C. Thaw frozen samples and keep them at room temperature for 30 minutes before use. After thawing, the stir samples to achieve homogeneity. Avoid repeated freezing-thawing cycles for test samples. If serum (or plasma) is turbid, remove insoluble inclusions by centrifugation at 3000 rpm for 10-15 minutes. Do not use serum samples with hyperlipidemia, hemolysis, and bacterial growth.

Transport serum samples in insulated containers. To do that, put closed labelled tubes in a plastic bag, tightly seal it and place in the centre of an insulated container. Put the frozen cold packs on the bottom, along the side walls of the insulated container and on top of the serum samples.

8. REAGENT PREPARATION

NOTE! It is very important to keep all ELISA kit components for at least 30 min at room temperature 18-25 °C before the assay!

8.1. Microplate preparation

To prevent water condensation in the wells, keep the **STRIPS** for 30 minutes at a room temperature before opening. Open the vacuum pack, detach the appropriate number of wells, and carefully pack the remaining wells with a desiccant and store tightly zip-locked at 2-8 °C. Storing the packed plate this way ensures its stability for 6 months.

8.2. Washing solution preparation

To prepare detergent, dilute **TWEEN|WASH|20x** at 1:20 (1+19) with distilled or deionized water and stir. E. g., 5 mL of concentrate + 95 mL of water, which is enough for 8 wells. If there are crystals present in the detergent concentrate, heat the vial at 37 °C until the crystals dissolve completely (15–20 minutes). Store the diluted solution at 2-8 °C for a maximum of 7 days.

9. ASSAY PROCEDURE

- 9.1. Prepare the necessary number of wells (four wells for controls and a necessary number of wells for test samples) and insert them into the ELISA plate frame. Be sure to add control wells in every test run.
- 9.2. Fill in the sample application plan.
- 9.3. Prepare the detergent as per para. 8.2.

- 9.4. Add 90 μL of [DIL|SAMPLE] into each plate well.
- 9.5. Add 10 μL of controls and test samples into the wells:
[CONTROL +] – into well A1,
[CONTROL -] – into wells B1, C1 and D1,
and test samples into the remaining wells.
At the time of adding, the solution changes its colour from purple to blue.
Pipette the mix in the wells carefully to avoid foaming.
- 9.6. Cover the strips up with adhesive film and incubate for **30 minutes at 37 °C**.
- 9.7. Remove and discard the adhesive film and wash all wells 5 times with automatic washer or 8-channel pipette as follows:
– aspirate the content of all wells into a liquid waste container;
– add a minimum of 300 μl of diluted washing solution to each well, soak each well for 30 seconds;
– aspirate the content of all wells again. The residual volume after every aspiration should be less than 5 μl ;
– repeat the washing step 4 more times;
– after the final aspiration, eliminate extra moisture by tapping the plate against a piece of filter paper.
- 9.8. Add 100 μL of [SOLN|CONJ] into each well. Cover the strips with a new piece of adhesive film and incubate for **30 minutes at 37 °C**.
- 9.9. Following incubation, remove the film carefully and wash the wells five times as described in para. 9.7.
- 9.10. Add 100 μL of [SOLN|TMB] into the wells; do not touch the bottom and the walls of the plate wells.
- 9.11. Incubate the strips for **30 minutes** in a dark place at a room temperature of 18-25 °C. Do not use adhesive film at this stage.
- 9.12. Add 100 μL of [SOLN|STOP] into each strip well to stop the enzymatic reaction; adhere to the same sequence of actions as when adding [SOLN|TMB]. At the time of adding, the solution colour changes from blue to yellow, and clear solution slightly changes its shade.
- 9.13. Measure the optical density (OD) of the wells at 450/620-695 nm wavelength using an ELISA microplate reader within 5 minutes after stopping the reaction. Pay attention to the cleanness of the plate bottom and the absence of bubbles in the wells before reading.

Measurement at the single wavelength of 450 nm is possible, in that case, it is needed to leave one well for blank (only [SOLN|TMB] and [SOLN|STOP] must be added in blank well).

10. CALCULATION AND INTERPRETATION OF RESULTS

10.1. Calculation of results

Calculate the average OD for the negative control (\bar{Nc}), Cut off (CO) and a sample positivity index (IP_{sample}).

$$\bar{Nc} = (Nc1 + Nc2 + Nc3)/3; \quad CO = \bar{Nc} + 0,3$$

$$IP_{\text{sample}} = OD_{\text{sample}}/CO, \text{ where } OD_{\text{sample}} \text{ is the OD sample.}$$

10.2. Quality control (assay validation)

The test results are considered valid if they meet the following requirements:

$$\boxed{\text{CONTROL} +} \quad OD \geq 1,0$$

$$\boxed{\text{CONTROL} -} \quad OD \leq 0,150$$

$$\boxed{\text{CONTROL} -} \quad \bar{Nc} \times 0,5 \leq Ncn \leq \bar{Nc} \times 2,0 \quad \text{where } Ncn \text{ is the OD for each } Nc \text{ run}$$

If any of the OD values for the negative control is beyond the above interval, it should be discarded, and Nc is calculated based on the remaining OD values for the negative control. If several OD values for the negative control fail to meet the above requirements, the test is considered invalid and requires a new run.

10.3. Interpretation of results

$$IP_{\text{sample}} > 1,1 \quad \text{POSITIVE}$$

$$0,9 \leq IP_{\text{sample}} \leq 1,1 \quad \text{BORDERLINE*}$$

$$IP_{\text{sample}} < 0,9 \quad \text{NEGATIVE}$$

* Uncertain samples are recommended to be re-examined in two wells of the ELISA kit. If the results are again uncertain, a new sample should be selected and analyzed in 2-4 weeks. In case of repeated indeterminate results, such samples shall be considered negative.

11. PERFORMANCE CHARACTERISTICS

11.1. Analytical performance characteristics

Precision of measurement

Intra assay repeatability

The coefficient of variation (CV) for two sera with different levels of specific antibodies was evaluated in 32 replicates on one series of ELISA kits.

Sample No.	OD _{av}	IP _{av}	CV, %
G12	0,554	1,59	4,1
E5	0,999	2,87	6,5

Inter assay reproducibility

The coefficient of variation (CV) for three sera with different levels of specific antibodies was evaluated for 4 days in 4 sets of analysis, 8 replicates in each analysis.

Sample No.	OD _{av}	IP _{av}	CV, %
G12	0,561	1,56	4,7
E5	1,042	2,89	7,1

Analytical specificity

The test results are not affected by bilirubin at up to 0.21 mg/mL (361.8 µmol/L), haemoglobin at up to 10 mg/mL and triglycerides at up to 10 mg/mL (11.3 mmol/l) present in the sample.

11.2. Diagnostic characteristics

Studies of the characteristics of the method in comparison with a similar commercial ELISA kit were performed on a sample of characterized serum and a sample of donor blood serum. The relative sensitivity of «EQUI anti-Trichinella spiralis» ELISA kits was determined from 24 serum samples that were tested for antibodies to anti-*Trichinella spiralis* and characterized as positive in a commercial ELISA kit. According to the results of the analysis, the relative sensitivity of the sets «EQUI anti-Trichinella spiralis» is 96%. Studies of the relative specificity of ELISA kits were performed on a complex of 260 samples of donor blood serum. According to research, the relative specificity of ELISA kits «EQUI anti-Trichinella spiralis» was 98%. The coincidence percentage is 97.86%.

12. LIMITATIONS OF ASSAY

A positive result in the «EQUI anti-Trichinella spiralis» ELISA kit is an indication that the patient has antibodies of IgG and / or IgA specific to *Trichinella spiralis*. The presence of IgG antibodies in newborns is not an evidence of *Trichinella spiralis* invasion.

Uncertain results may indicate a history of *Trichinella spiralis* invasion. A negative result in the «EQUI anti-Trichinella spiralis» ELISA kit indicates the absence of antibodies specific for *Trichinella spiralis* in the test sample or the concentration of specific antibodies below the sensitivity limit of the analysis. Specific antibodies may not be detected at the beginning of the clinical manifestations of the invasion. In this case, it is recommended to re-obtain and test serum samples from patients with clinical signs of trichinosis in one or two weeks.

The final diagnosis cannot be established solely on the basis of serological test results. When making a diagnosis the results of a set of laboratory and instrumental studies, should be taken into account as well as clinical manifestations of the disease. Cross-reactions with antibodies to antigens of other helminths cannot be completely ruled out. To exclude a false-positive result, it is recommended to conduct a verification study of positive samples by immunoblotting.

13. DIFFICULTIES THAT CAN OCCUR DURING THE ASSAY PROCEDURE

Possible reasons	Solution
<i>High background in all wells</i>	
Contaminated washer	Clean the washer head and rinse according to the instructions for use
Poor quality or contaminated water	Use purified water with specific resistance $\geq 10 \text{ M}\Omega \cdot \text{cm}$

Use of poorly washed glassware	Use chemically clean utensils
Use of chlorinated disinfectants	Do not use chlorine disinfectants
Use of contaminated tips	Use new tips
Increased incubation times or change in the temperature conditions	Adhere to the incubation regime according to the instructions for use
<i>High background in a row of wells</i>	
Repeat application of TMB solution	TMB solution should be applied once
Contamination of the automatic pipette nozzle with conjugate solution	Clean the pipette and dial carefully liquid
Contamination of one of the washer's channel	Clean the flush channel, rinse washer
<i>Received OD of the positive control is below the border value</i>	
One of the reagents (conjugate solution or TMB solution) was not prepared in a correct way or was not added	Re-conduct ELISA, pay attention to the correctness of the introduction of these reagents
Reduced incubation times at any stage	Incubate according to instructions for use
<i>The colour density of the wells fails to meet the obtained optical density value</i>	
This may suggest that the optical beam has been displaced	Check the correct operation of the reader

14. TECHNICAL ASSISTANCE AND CUSTOMER SERVICE

In case of technical problems, you can obtain assistance by contacting the manufacturer.

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ASSAY PROCEDURE SCHEME

Keep all reagents for 30 min at temperature 18-25°C before use

Dispense 90 µl [DIL|SAMPLE] into the wells (purple)

Add to 10 µl of controls and samples into the wells:
A1 – [CONTROL|+], B1, C1, D1 – [CONTROL|-],
other wells – examined samples
(change of colour from purple to blue)

Cover strips with an adhesive film, incubate for **30 min at 37°C**

Rinse the wells 5 times with prepared 1:20 (1+19) washing solution TWEEN (300 µl per well)

Add 100 µl of [SOLN|CONJ] into all wells (green)

Cover strips with an adhesive film, incubate for **30 min at 37°C**

Rinse the wells 5 times with prepared 1:20 (1+19) washing solution TWEEN (300 µl per well)

Add 100 µl of [SOLN|TMB] into all wells

Incubate for **30 min** in the dark at **18-25°C**

Add 100 µl of [SOLN|STOP] into all wells (change of colour from blue to yellow)

Measure the optical density (OD) with an ELISA microplate reader at 450/620-695 nm

CALCULATION OF RESULTS

$$\bar{Nc} = (Nc1 + Nc2 + Nc3)/3;$$

$$CO = \bar{Nc} + 0,3;$$

$$IP_{\text{sample}} = OD_{\text{sample}}/CO$$

\bar{Nc} - the average value of OD 3-x [CONTROL|-]

CO - Cut off

IP_{sample} - sample positivity index

INTERPRETATION OF RESULTS


$IP_{\text{sample}} > 1,1$	POSITIVE
$0,9 \leq IP_{\text{sample}} \leq 1,1$	BORDERLINE
$IP_{\text{sample}} < 0,9$	NEGATIVE

anti-Lamblia

ELISA kit for the qualitative detection of antibodies
to *Giardia lamblia (intestinalis)*

Instructions for use

IVD

REF
EI-606

Σ 96
tests

CE

EQUI anti-Lamblia

ELISA kit for the qualitative detection of antibodies to
до *Giardia lamblia (intestinalis)*

1. INTENDED USE

The «EQUI anti-Lamblia» is ELISA kit intended to qualitatively detect antibodies to *Giardia lamblia (intestinalis)* in human serum or plasma by enzyme-linked immunosorbent assay (ELISA) to diagnose giardiasis. The testing procedure is designed for both manual arrangement with automatic pipettes and standard equipment, and for automated «open» immunoassay analysers.

Target group: children, pet owners, citizens of rural areas, summer house owners.

Usage: ELISA kit is used in clinical diagnostic laboratories and other institutions engaged in *in vitro* diagnostics.

2. CLINICAL SIGNIFICANCE

Giardiasis is considered one of the most common parasitic diseases of the small intestine in the world. This infection is a major cause of acute and chronic diarrhea, especially in children. The etiological agent of giardiasis is *Giardia lamblia*, which is also called *Giardia intestinalis* or *Giardia duodenalis*.

Giardia lamblia are unicellular flagellate protozoa that parasitize in the intestines of humans and some other mammals. During the life cycle of these parasites, two stages alternate: cysts, resistant to external conditions, and a vegetative form - trophozoites. Infection occurs when cysts enter the human gastrointestinal tract. After experiencing the effects of gastric acid, cysts in the duodenum turn into trophozoites, which parasitize in the upper parts of the small intestine. They absorb nutrients from the intestinal lumen, block parietal digestion and disrupt the motility of the intestine.

Humans get infected via fecal-oral routes through cyst-contaminated food, water, unwashed hands, and so on. *Giardia* can also be transmitted to humans from infected cats, dogs, and livestock. Giardiasis is especially common in regions with poor sanitation. In addition, human-to-human transmission is common in preschools.

In many cases, the invasion of *Giardia* occurs without clinical manifestations. In other cases, the first symptoms of giardiasis appear in 1-3 weeks after infection. They are most often manifested by spasms, bloating, nausea and diarrhea, which leads to dehydration and weight loss. The acute form of the disease can last up to two weeks and end in recovery without additional treatment or become chronic. Chronic giardiasis develops when the duration of the invasion is longer than 2 month and the exacerbation of clinical manifestations (diarrhea) is cyclical. *Giardia lamblia* parasitism can lead to malabsorption syndrome, which disrupts the absorption of carbohydrates and fats, as well as the metabolism of vitamins B12, A and C.

Immune response to invasion and non-immune factors are important to control the development of the disease and the severity of clinical manifestations. Both

humoral and cellular immunity play the part in the eradication of the pathogen, the role of which is still subjected to scientific research. In addition, partial resistance to re-infection is formed due to protective mechanisms of the body.

Typically, to diagnose giardiasis, the duodenal contents and feces are examined for trophozoites and cysts of giardiasis. In case of the chronic course of the disease, cysts get excreted periodically, and, considering this, the additional tests should be performed regularly for several weeks. Another method of diagnosing giardiasis is to detect *Giardia lamblia* antigens in the feces. However, serodiagnosis with the detection of specific antibodies to *Giardia* antigens is an important step in assessing the immune response of patients. Detection of specific IgM antibodies suggests an acute stage of giardiasis. However, the detection of specific IgG and IgA antibodies should be interpreted with caution: in some regions they persist for a long time after infection, while in others their level decreases after eradication of the pathogen.

3. ANALYSIS PRINCIPLE

The procedure of testing for *Giardia lamblia* specific antibodies in «EQUI anti-Lamblia» ELISA kit is based on «indirect» solid-phase ELISA with a two-stage incubation. Recombinant *Giardia lamblia* antigens are entrapped in the wells. During the first step of incubation of the test samples in the wells of the ELISA plate, *Giardia lamblia*-specific antibodies, if present in the samples, bind to the solid phase antigens. The wells are washed to remove unbound antibodies and have only specific antigen-antibody complexes left. Then, a conjugate of anti-species (anti-IgG and anti-IgA) monoclonal antibodies with horseradish peroxidase is added, which binds to solid-phase immune complexes. Unbound components are removed by washing. Antigen-antibody complexes are detected by adding a solution of chromogen 3,3',5,5'-tetramethylbenzidine (TMB) with hydrogen peroxide. After 30-minute incubation, the reaction is stopped by adding the stop solution. The optical density (OD) in the wells is determined using a spectrophotometer at 450/620-695 nm. The intensity of the yellow colour is proportional to the level of antibodies in the sample.

4. MATERIALS AND EQUIPMENT

4.1. Contents of the ELISA kit

STRIPS	1 x 96 wells	Microplate Each plate well is coated with <i>Giardia lamblia</i> purified antigens. The wells are detachable. After the first opening, store unused strips in the package at 2-8 °C for a maximum of 6 months
CONTROL +	1 x 0,35 ml	Positive control Conjugated specific monoclonal antibody solution with preservative (pink). Store at 2-8 °C
CONTROL -	1 x 1,2 ml	Negative control Negative human serum with a preservative (yellow). Store at 2-8 °C

DIL SAMPLE	1 x 11 ml	Serum dilution solution Buffer solution with a milk extract, a detergent and a preservative (purple). Store at 2-8 °C
SOLN CONJ	1 x 13 ml	Conjugate solution (ready to use) Buffer solution of monoclonal antibodies to human IgG and IgA, conjugated with horseradish peroxidase, with stabilizers and preservative (green). Store at 2-8 °C
SOLN TMB	1 x 13 ml	TMB solution (ready to use) TMB solution, H ₂ O ₂ , a stabilizer, a preservative (colourless). Store at 2-8 °C
TWEEN WASH 20x	1 x 50 ml	Washing solution TWEEN (20x concentrated) 20-fold phosphate buffer concentrate with Tween-20 (colourless). Dilute TWEEN detergent (20x) at 1:20 with distilled or deionized water (e. g., 5 mL of concentrate + 95 mL of water for 8 wells) before use. Store the diluted solution at 2-8 °C for a maximum of 7 days
SOLN STOP	1 x 13 ml	Stop Solution (ready to use) 0.5 mol H ₂ SO ₄ solution (colourless). Store at 2-8 °C

The ELISA kit also includes adhesive films (2 items), sample application plan (1 item), checklist, and instruction for use.

4.2. Optional reagents, materials and equipment

Automatic single and multichannel pipettes 10–1000 µL, tips, volumetric laboratory glassware (10–1,000 mL), deionized or distilled water, thermostat at 37 °C, automatic or semi-automatic plate washer, spectrophotometer (reader) for microplates at 450/620–695 nm, appropriate containers for potentially contaminated waste, timer, filter paper, disposable powder-free gloves, disinfectants.

5. PRECAUTIONS AND SAFETY

5.1. Precautions

Be sure to read the instructions for use carefully before the test. The validity of the test results depends on strict following of the test procedure.

- do not use the ELISA kit components after the expiry date;
- do not use for analysis or mix components of different batches, components of kits for different nosologies, or reagents from other manufacturers with the «EQUI anti-Lambliia» ELISA kit;
- do not freeze the ELISA kit or its contents;
- after using a reagent, close each vial with its cap;
- when washing, control filling and complete aspiration of solution from the wells;
- use a new pipette tip each time you add samples or reagents;
- prevent direct sunlight from reaching the reagents from the ELISA kit;
- SOLN|TMB solution must be colourless before use. Do not use the solution if its colour is blue or yellow. Avoid contact of SOLN|TMB with metals or metal ions. Use only clean glassware thoroughly rinsed with distilled water;

- do not use reagents with colour not in line with para. 4.1;
- under no circumstances should the same glassware be used for SOLN|CONJ and SOLN|TMB;
- do not evaluate the test results visually (without a reader);
- any optional equipment that is in direct contact with biological material or kit components should be considered contaminated and requires cleaning and decontamination;
- the ELISA kit includes materials for 96 tests. Dispose of the used components as well as any remaining unused components.

5.2. Safety requirements

- all reagents in the ELISA kit are for laboratory professional use for *in vitro* diagnosis only and may only be used by qualified personnel;
- conduct the tests in disposable powder-free gloves and goggles only;
- do not eat, drink, smoke, or apply make-up in the test room;
- do not mouth-pipette the solutions;
- controls from the «EQUI anti-Lambliа» ELISA kit have been tested and found to be for anti-HIV1/2, anti-HCV and anti-*Treponema pallidum* antibodies and HBsAg negative; however, controls and test samples should be handled as potentially hazardous infectious materials;
- some of the kit components contain low concentrations of harmful substances and can damage skin or mucosa. In case of contact of SOLN|TMB, SOLN|STOP and SOLN|CONJ with mucous membranes or skin, immediately wash the affected area with plenty of water;
- in case of spillage of acid-free solutions, e. g. sera, treat the surface with a disinfectant solution and then wipe dry with filter paper. Otherwise first neutralize acid with sodium bicarbonate solution and then wipe the surface dry as described above.

5.3. Waste inactivation and disposal

- the liquid waste must be inactivated, for example, with hydrogen peroxide solution at the final concentration of 6% for 3 hours at room temperature, or with sodium hypochlorite at the final concentration of 5% for 30 minutes, or with other approved disinfectants;
- the solid waste must be inactivated by autoclaving at a temperature not less than 132°C;
- do not autoclave the solutions that contain sodium azide or sodium hypochlorite;
- disposal of inactivated waste must be conducted due to national laws and regulations.

6. STORAGE AND STABILITY

ELISA kit is stable up to the expiry date stated on the label when stored at 2-8°C. The kit should be transported at 2-8°C. Single transportation at a

temperature up to 23°C for two days is possible.

7. SAMPLE COLLECTION, TRANSPORTATION AND STORAGE GUIDELINES

Collect blood from the vein into the sterile test tube. Test tube must be marked with patient ID and date of sample collecting. Blood before serum separation can be stored at 2-8 °C for 24 hours, avoiding freezing.

Serum or plasma can be stored at 2-8 °C for maximum 3 days. Frozen serum can be stored for longer periods of time at -20 °C or -70 °C. Thaw frozen samples and keep them at room temperature for 30 minutes before use. After thawing, the stir samples to achieve homogeneity. Avoid repeated freezing-thawing cycles for test samples. If serum (or plasma) is turbid, remove insoluble inclusions by centrifugation at 3000 rpm for 10-15 minutes. Do not use serum samples with hyperlipidemia, hemolysis, and bacterial growth.

Transport serum samples in insulated containers. To do that, put closed labelled tubes in a plastic bag, tightly seal it and place in the centre of an insulated container. Put the frozen cold packs on the bottom, along the side walls of the insulated container and on top of the serum samples.

8. REAGENT PREPARATION

NOTE! It is very important to keep all ELISA kit components for at least 30 min at room temperature 18-25 °C before the assay!

8.1. Microplate preparation

To prevent water condensation in the wells, keep the **STRIPS** for 30 minutes at a room temperature before opening. Open the vacuum pack, detach the appropriate number of wells, and carefully pack the remaining wells with a desiccant and store tightly zip-locked at 2-8 °C. Storing the packed plate this way ensures its stability for 6 months.

8.2. Washing solution preparation

To prepare detergent, dilute **TWEEN|WASH|20x** at 1:20 (1+19) with distilled or deionized water and stir. E. g., 5 mL of concentrate + 95 mL of water, which is enough for 8 wells. If there are crystals present in the detergent concentrate, heat the vial at 37 °C until the crystals dissolve completely (15–20 minutes). Store the diluted solution at 2-8 °C for a maximum of 7 days.

9. ASSAY PROCEDURE

9.1. Prepare the necessary number of wells (four wells for controls and a necessary number of wells for test samples) and insert them into the ELISA plate frame.

Be sure to add control wells in every test run.

9.2. Fill in the sample application plan.

9.3. Prepare the detergent as per para. 8.2.

9.4. Add 80 µL of **DIL|SAMPLE** into each plate well.

9.5. Add 20 µL of controls and test samples into the wells:

CONTROL|+ – into well A1,

CONTROL|- – into wells B1, C1 and D1,

and test samples into the remaining wells.

At the time of adding, the solution changes its colour from brown to blue. Pipette the mix in the wells carefully to avoid foaming.

- 9.6. Cover the strips up with adhesive film and incubate for **30 minutes at 37 °C**.
- 9.7. Remove and discard the adhesive film and wash all wells 5 times with automatic washer or 8-channel pipette as follows:
- aspirate the content of all wells into a liquid waste container;
 - add a minimum of 300 µl of diluted washing solution to each well, soak each well for 30 seconds;
 - aspirate the content of all wells again. The residual volume after every aspiration should be less than 5 µl;
 - repeat the washing step 4 more times;
 - after the final aspiration, eliminate extra moisture by tapping the plate against a piece of filter paper.
- 9.8. Add 100 µL of [SOLN|CONJ] into each well. Cover the strips with a new piece of adhesive film and incubate for **30 minutes at 37 °C**.
- 9.9. Following incubation, remove the film carefully and wash the wells five times as described in para. 9.7.
- 9.10. Add 100 µL of [SOLN|TMB] into the wells; do not touch the bottom and the walls of the plate wells.
- 9.11. Incubate the strips for **30 minutes** in a dark place at a room temperature of 18-25 °C. Do not use adhesive film at this stage.
- 9.12. Add 100 µL of [SOLN|STOP] into each strip well to stop the enzymatic reaction; adhere to the same sequence of actions as when adding [SOLN|TMB]. At the time of adding, the solution colour changes from blue to yellow, and clear solution slightly changes its shade.
- 9.13. Measure the optical density (OD) of the wells at 450/620-695 nm wavelength using an ELISA microplate reader within 5 minutes after stopping the reaction. Pay attention to the cleanness of the plate bottom and the absence of bubbles in the wells before reading.

Measurement at the single wavelength of 450 nm is possible, in that case, it is needed to leave one well for blank (only [SOLN|TMB] and [SOLN|STOP] must be added in blank well).

10. CALCULATION AND INTERPRETATION OF RESULTS

10.1. Calculation of results

Calculate the average OD for the negative control (\bar{Nc}), Cut off (CO) and a sample positivity index (IP_{sample}).

$$\bar{Nc} = (Nc1 + Nc2 + Nc3)/3; \quad CO = \bar{Nc} + 0,25$$

$$IP_{\text{sample}} = OD_{\text{sample}}/CO, \text{ where } OD_{\text{sample}} \text{ is the OD sample.}$$

10.2. Quality control (assay validation)

The test results are considered valid if they meet the following requirements:

CONTROL + OD \geq 1,0

CONTROL - OD \leq 0,150

CONTROL - $\bar{N}c \times 0,5 \leq N_{cn} \leq \bar{N}c \times 2,0$

where N_{cn} is the OD for each N_c run

If any of the OD values for the negative control is beyond the above interval, it should be discarded, and N_c is calculated based on the remaining OD values for the negative control. If several OD values for the negative control fail to meet the above requirements, the test is considered invalid and requires a new run.

10.3. Interpretation of results

$IP_{\text{sample}} > 1,1$ POSITIVE
 $0,9 \leq IP_{\text{sample}} \leq 1,1$ BORDERLINE*
 $IP_{\text{sample}} < 0,9$ NEGATIVE

* Uncertain samples are recommended to be re-examined in two wells of the ELISA kit. If the results are again uncertain, a new sample should be selected and analyzed in 2-4 weeks. In case of repeated indeterminate results, such samples shall be considered negative.

11. PERFORMANCE CHARACTERISTICS

11.1. Analytical performance characteristic

Precision of measurement

Intra assay repeatability

The coefficient of variation (CV) for two sera with different levels of specific antibodies was evaluated in 32 replicates on one series of ELISA kits.

Sample No.	OD _{av}	IP _{av}	CV, %
14L	0,679	2,47	6,5
16L	0,490	1,79	6,6

Inter assay reproducibility

The coefficient of variation (CV) for three sera with different levels of specific antibodies was evaluated for 3 days in 3 sets of analysis, 8 replicates in each analysis.

Sample No.	OD _{av}	IP _{av}	CV, %
14L	0,670	2,39	5,55
16L	0,463	1,65	7,06

Analytical specificity

The test results are not affected by bilirubin at up to 0.21 mg/mL (361.8 μ mol/L), haemoglobin at up to 10 mg/mL and triglycerides at up to 10 mg/mL (11.3 mmol/l) present in the sample.

11.2. Diagnostic characteristics

Studies of the characteristics of the method in comparison with a similar commercial ELISA kit were performed on a sample of characterized sera, the target group of children and a group of donors. The relative sensitivity of «EQUI anti-Lambli» ELISA kits was determined from a group of 23 serum samples that were tested for antibodies to *Giardia lamblia* and characterized as positive in a commercial ELISA kit. All sera were also determined to be positive in «EQUI anti-Lambli» kits, so the relative sensitivity equals 100%. For 148 serum samples of children that were tested and characterized in commercial analogues, the relative specificity of «EQUI anti-Lambli» ELISA kits was 92.86%, the percentage of coincidence - 93.24%. According to a similar principle, for 238 serum samples of donor blood, the relative specificity was 97% and the percentage of coincidence was 96.64%.

12. LIMITATIONS OF ASSAY

The final diagnosis cannot be made solely on the basis of serological test results, since clinical manifestations of the disease and laboratory data (such as the detection of cysts in faecal samples or trophozoites in duodenal contents; the results of detection of *Giardia lamblia* antigen in faeces) should be taken into account as well.

Additionally, cross-reactions with antibodies to antigens of other parasites cannot be completely ruled out.

Giardia lamblia-specific antibodies may not be detected in case of children with persistent and prolonged giardiasis.

It should be noted that IgG antibodies to *Giardia lamblia* can be detected via ELISA for a long time, even after successful treatment.

13. DIFFICULTIES THAT CAN OCCUR DURING THE ASSAY PROCEDURE

Possible reasons	Solution
High background in all wells	
Contaminated washer	Clean the washer head and rinse according to the instructions for use
Poor quality or contaminated water	Use purified water with specific resistance $\geq 10 \text{ M}\Omega \cdot \text{cm}$
Use of poorly washed glassware	Use chemically clean utensils
Use of chlorinated disinfectants	Do not use chlorine disinfectants
Use of contaminated tips	Use new tips
Increased incubation times or change in the temperature conditions	Adhere to the incubation regime according to the instructions for use
High background in a row of wells	

Repeat application of TMB solution	TMB solution should be applied once
Contamination of the automatic pipette nozzle with conjugate solution	Clean the pipette and dial carefully liquid
Contamination of one of the washer's channel	Clean the flush channel, rinse washer
<i>Received OD of the positive control is below the border value</i>	
One of the reagents (conjugate solution or TMB solution) was not prepared in a correct way or was not added	Re-conduct ELISA, pay attention to the correctness of the introduction of these reagents
Reduced incubation times at any stage	Incubate according to instructions for use
<i>The colour density of the wells fails to meet the obtained optical density value</i>	
This may suggest that the optical beam has been displaced	Check the correct operation of the reader

14. TECHNICAL ASSISTANCE AND CUSTOMER SERVICE

In case of technical problems, you can obtain assistance by contacting the manufacturer.

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	Manufacturer
	Authorized Representative in the European Community
	In vitro diagnostic medical device
	Catalogue number
	Date of manufacture
	Use by date
	Batch code
	Temperature limit
	Contains sufficient for <n> tests
	Caution
	Non-Sterile
	Consult instructions for use
	Keep away from sunlight
	Keep dry
	Compliance with EU safety requirements

Edition 7, 18.02.2022

For questions and suggestions regarding the ELISA kit contact:

 Obelis s.a.
 Bd Général Wahis 53
 1030 Brussels
 Belgium
 Tel: +(32)2 732-59-54
 Fax: +(32)2 732-60-03
 mail@obelis.net

 Ekvitestlab LLC
 Velyka Vasylkivska St. 114, Kyiv, Ukraine, 03150
 Tel: 0(800)31-89-87, +38 (044)334-89-87,
 e-mail: info@equitest.com.ua, www.equitest.com.ua

ASSAY PROCEDURE SCHEME

Keep all reagents for 30 min at temperature 18-25°C before use

Dispense 80 µl [DIL|SAMPLE] into the wells
(purple)

Add to 20 µl of controls and samples into the wells:
A1 – [CONTROL|+], B1, C1, D1 – [CONTROL|-],
other wells – examined samples
(change of colour from purple to blue)

Cover strips with an adhesive film, incubate for **30 min at 37°C**

Rinse the wells 5 times with prepared 1:20 (1+19) washing solution TWEEN
(300 µl per well)

Add 100 µl of [SOLN|CONJ] into all wells
(green)

Cover strips with an adhesive film, incubate for **30 min at 37°C**

Rinse the wells 5 times with prepared 1:20 (1+19) washing solution TWEEN
(300 µl per well)

Add 100 µl of [SOLN|TMB] into all wells

Incubate for **30 min** in the dark at **18-25°C**

Add 100 µl of [SOLN|STOP] into all wells
(change of colour from blue to yellow)

Measure the optical density (OD) with an ELISA microplate reader at
450/620-695 nm

CALCULATION OF RESULTS

$$\bar{Nc} = (Nc1 + Nc2 + Nc3)/3;$$

$$CO = \bar{Nc} + 0,25;$$

$$IP_{\text{sample}} = OD_{\text{sample}}/CO$$

\bar{Nc} - the average value of OD 3-x [CONTROL|-]

CO - Cut off

IP_{sample} - sample positivity index

INTERPRETATION OF RESULTS

$IP_{\text{sample}} > 1,1$	POSITIVE
$0,9 \leq IP_{\text{sample}} \leq 1,1$	BORDERLINE
$IP_{\text{sample}} < 0,9$	NEGATIVE

Fecal Occult Blood Rapid Test Cassette (Feces)

INTENDED USE

Fecal Occult Blood Rapid Test Cassette (Feces) is a rapid chromatographic immunoassay for the qualitative detection of human occult blood in feces by professional laboratories or physician's offices. It is useful to detect bleeding caused by a number of gastrointestinal disorders, e.g., diverticulitis, colitis, polyps, and colorectal cancer.

Fecal Occult Blood Rapid Test Cassette (Feces) is recommended for use in 1) routine physical examinations, 2) hospital monitoring for bleeding in patients, and 3) screening for colorectal cancer or gastrointestinal bleeding from any source.

INTRODUCTION

Most of diseases can cause hidden blood in the stool. In the early stages, gastrointestinal problems such as colon cancer, ulcers, polyps, colitis, diverticulitis, and fissures may not show any visible symptoms, only occult blood. Traditional guaiac-based method lacks sensitivity and specificity, and has diet-restriction prior to the testing.

Fecal Occult Blood Rapid Test Cassette (Feces) is a rapid test to qualitatively detect low levels of fecal occult blood in feces. The test uses double antibody-sandwich assay to selectively detect as low as 50 ng/mL of hemoglobin or 6 µg hemoglobin/g feces. In addition, unlike the guaiac assays, the accuracy of the test is not affected by the diet of the patients.

PRINCIPLE

Fecal Occult Blood Rapid Test Cassette (Feces) is a lateral flow chromatographic immunoassay based on the principle of the double antibody-sandwich technique. The membrane is pre-coated with anti-hemoglobin antibodies on the test line region of the device. During testing, the specimen reacts with the colloidal gold coated with anti-hemoglobin antibodies. The mixture migrates upward on the membrane chromatographically by capillary action to react with anti-hemoglobin antibodies on the membrane and generate a colored line. The presence of this colored line in the test region indicates a positive result, while its absence indicates a negative result. To serve as a procedural control, a colored line will always appear in the control line region indicating that proper volume of specimen has been added and membrane wicking has occurred.

MATERIALS PROVIDED

- 20 Test cassettes
- 20 Specimen collection tubes with buffer
- 1 Package insert

MATERIALS REQUIRED BUT NOT PROVIDED

1. Specimen collection containers
2. Clock or timer

STORAGE AND STABILITY

All reagents are ready to use as supplied. Store unused test device unopened at 2°C-30°C. If stored at 2°C-8°C, ensure that the test device is brought to room temperature before opening. The test is not stable out of the expiration date printed on the sealed pouch. Do not freeze the kit or expose the kit over 30°C.

PRECAUTIONS

1. For professional *in vitro* diagnostic use only.
2. This package insert must be read completely before performing the test. Failure to follow the insert gives inaccurate test results.
3. Do not use it if the tube/pouch is damaged or broken.
4. Test is for single use only. Do not re-use under any circumstances.
5. **Do not use specimen with visible blood for the testing.**
6. Handle all specimens as if they contain infectious agents. Observe established standard procedure for proper disposal of specimens.
7. Specimen extraction buffer contains Sodium Azide (0.1%). Avoid contact with skin or eyes. Do not ingest.
8. Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are assayed.
9. Humidity and temperature can adversely affect results.
10. Do not perform the test in a room with strong air flow, i.e. electric fan or strong air conditioning.

PATIENT PREPARATION

1. A specimen should not be collected from a patient with following conditions that may interfere with the test results:

- Menstrual bleeding
 - Bleeding hemorrhoids
 - Constipating bleeding
 - Urinary bleeding.
2. Dietary restrictions are not necessary.
 3. Alcohol and certain medications such as aspirin, indomethacin, phenylbutazone, reserpine, cortocosteroids, and nonsteroidal anti-inflammatory drugs may cause gastrointestinal irritation and subsequent bleeding, thus gives positive reactions. On the advice of the physician, such substances should be discontinued at least 48 hours prior to testing.

SPECIMEN COLLECTION AND PREPARATION

Consider any materials of human origin as infectious and handle them using standard biosafety procedures.


1. Collect a random sample of feces in a clean, dry receptacle.
2. Unscrew the top of the collection tube and remove the applicator stick.
3. Randomly pierce the fecal specimen in at least five (5) different sites.
4. Remove excess sample off the shaft and outer grooves. Be sure sample remains on inside grooves.
5. Replace the stick in the tube and tighten securely.
6. Shake the specimen collection bottle so that there is proper homogenisation of feces in buffer solution.

Note: Specimens prepared in the specimen collection tube may be stored at room temperature (15-30°C) for 3 days maximum, at 2-8°C for 7 days maximum or at -20°C for 3 months maximum if not tested within 1 hour after preparation.

TEST PROCEDURE

Allow the test cassette, specimen, and/or controls to reach room temperature (15-30°C) prior to testing.

1. Remove the test cassette from the foil pouch and use it as soon as possible. Best results will be obtained if the assay is performed within one hour.
2. Place the test cassette on a clean, flat surface.
3. Shake the specimen collection tube several times.
4. Hold the specimen collection tube upright and then unscrew and open the upper cap.
5. Squeeze 3 drops (~90 µL) of the sample solution in the sample well of the cassette and start the timer.
6. Wait for the colored line(s) to appear. Read results in 5 minutes. Do not interpret the result after 5 minutes.

INTERPRETATION OF RESULTS

(Please refer to the illustration above)

Positive: Two lines appear. One colored line should be in the control line region (C) and another apparent colored line should be in the test line region (T).

Negative: One colored line appears in the control line region (C). No line appears in the test line region (T).

Invalid: Control line fails to appear. The test should be repeated using a new cassette. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

NOTE:

1. The intensity of color in the test region (T) may vary depending on the concentration of analytes present in the specimen. Therefore, any shade of color in the test region should be considered positive. Note that this is a qualitative test only, and

Fecal Occult Blood Rapid Test Cassette (Feces)

cannot determine the concentration of analytes in the specimen.

2. Insufficient specimen volume, incorrect operating procedure or expired tests are the most likely reasons for control band failure.

QUALITY CONTROL

An internal procedural control is included in the test. A colored line appearing in the control line region (C) is an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique. Control standards are not supplied with this kit; however it is recommended that positive and negative controls be tested as a good laboratory practice to confirm the test procedure and to verify proper test performance.

LIMITATIONS

1. This test kit is to be used for the qualitative detection of human hemoglobin in fecal samples. A positive result suggests the presence of human hemoglobin in fecal samples. In addition to intestinal bleeding the presence of blood in stools may have other causes such as hemorrhoids, blood in urine etc.
2. Not all colorectal bleedings are due to precancerous or cancerous polyps. The information obtained by this test should be used in conjunction with other clinical findings and testing methods, such as colonoscopy gathered by the physician.
3. Negative results do not exclude bleeding since some polyps and colorectal region cancers can bleed intermittently or not at all. Additionally, blood may not be uniformly distributed in fecal samples. Colorectal polyps at an early stage may not bleed.
4. Urine and excessive dilution of sample with water from toilet bowl may cause erroneous test results. The use of a receptacle is recommended.
5. Feces specimens should not collect during the menstrual period and not three day before or afterwards, at bleeding due to constipation, bleeding haemorrhoids, or at taking rectally administered medication. It could cause false positive results.
6. This test may be less sensitive for detecting upper g.i. Bleeding because blood degrades as it passes through the g.i. Track.
7. The Fecal Occult Blood Rapid Test Cassette (Feces) is to aid in diagnosis and is not intended to replace other diagnostic procedures such as G.I. fibroscope, endoscopy, colonoscopy, or X-ray analysis. Test results should not be deemed conclusive with respect to the presence or absence of gastrointestinal bleeding or pathology. A positive result should be followed up with additional diagnostic procedures to determine the exact cause and source for the occult blood in the feces.

PERFORMANCE CHARACTERISTICS

1. Sensitivity: 99.6%

Fecal Occult Blood Rapid Test Cassette (Feces) can detect the levels of human occult blood as low as 50 ng/mL hemoglobin or 6 µg hemoglobin/g feces.

2. Prozone Effect:

It is observed that this FOB test can detect 2 mg/mL hemoglobin.

3. Specificity: 99.9%


Fecal Occult Blood Rapid Test Cassette (Feces) is specific to human hemoglobin. Specimen containing the following substances at the standard concentration was tested on both positive and negative controls and showed no effects on test results at standard concentration.

Substances	Concentrations (Diluted with the extraction buffer)
Beef hemoglobin	2 mg/mL
Chicken hemoglobin	0.5 mg/mL
Pig hemoglobin	0.5 mg/mL
Goat hemoglobin	0.5 mg/mL
Horse hemoglobin	20 mg/mL
Rabbit hemoglobin	0.06 mg/mL


REFERENCES

1. Simon J.B. Occult Blood Screening for Colorectal Carcinoma: A Critical Review, Gastroenterology, Vol. 1985;88:820.
2. Blebae J. and Napherson RA. False-Positive Guaiac Testing With Iodine, Arch Pathol Lab Med, 1985;109:437-40.

INDEX OF SYMBOLS

	Consult instructions for use		Tests per kit		Authorized Representative
	For <i>in vitro</i> diagnostic use only		Use by		Do not reuse
	Store between 2~30°C		Lot Number		Catalog#

 Zhejiang Orient Gene Biotech Co.,Ltd
Address: 3787#, East Yangguang Avenue, Dipu Street,
Anji 313300, Huzhou, Zhejiang, China
Tel: +86-572-5226111 Fax: +86-572-5226222
Website: www.orientgene.com

 Shanghai International Holding Corp. GmbH (Europe)
Add: Eiffestrasse 80, 20537 Hamburg, Germany

 GEFOB-602b

H. pylori Ag Rapid Test Cassette (Feces)

INTENDED USE

H. pylori Ag Rapid Test Cassette (Feces) is a sandwich lateral flow chromatographic immunoassay for the qualitative detection of H. Pylori antigen in feces. It is for professional *in vitro* diagnostic use only.

INTRODUCTION

H. Pylori is associated with a variety of gastrointestinal diseases included non-ulcer dyspepsia, duodenal and gastric ulcer and active, chronic gastritis.^{1,2} The prevalence of H. pylori infection could exceed 90% in patients with signs and symptoms of gastrointestinal diseases. Recent studies indicate an association of H. Pylori infection with stomach cancer.³ H. Pylori colonizing in the gastrointestinal system elicits specific antibody responses^{4,5,6} which aids in the diagnosis of H. Pylori infection and in monitoring the prognosis of the treatment of H. Pylori related diseases. Antibiotics in combination with bismuth compounds have been shown to be effective in treating active H. Pylori infection. Successful eradication of H. pylori is associated with clinical improvement in patients with gastrointestinal diseases providing a further evidence.⁷

PRINCIPLE

H. pylori Ag Rapid Test Cassette (Feces) is a lateral flow chromatographic immunoassay based on the principle of the double antibody–sandwich technique. The test cassette consists of: 1) a burgundy colored conjugate pad containing H. Pylori antibodies conjugated with color particles (H. Pylori conjugates). 2) a nitrocellulose membrane strip containing a test band (T band) and a control band (C band). The T band is pre-coated with non-conjugated H. Pylori antibodies.

When an adequate volume of test specimen is dispensed into the sample well of the cassette, the specimen migrates by capillary action across the cassette. The antigen of H. Pylori if present in the specimen will bind to the H. Pylori antibodies conjugates. The immunocomplex is then captured on the membrane by the pre-coated H. Pylori antibodies, forming a burgundy colored T band, indicating a H. Pylori antigen positive test result. To serve as a procedural control, a colored line will always appear in the control line region indicating that proper volume of specimen has been added and membrane wicking has occurred. Otherwise, the test result is invalid and the specimen must be retested with another device.

PRODUCT CONTENTS

H. pylori Ag Rapid Test Cassette (Feces) containing anti- H. pylori antibodies particles and anti-H. pylori antibodies coated on the membrane.

MATERIALS SUPPLIED

- 20 Sealed pouches each containing a test cassette and a desiccant
- 20 Specimen collection tubes with extraction buffer, 2.0 mL
- 1 Package insert

MATERIAL REQUIRED BUT NOT PROVIDED

- 1. Clock or timer
- 2. Specimen collection containers.

STORAGE AND STABILITY

All reagents are ready to use as supplied. Store unused test device unopened at 2°C-30°C. If stored at 2°C-8°C, ensure that the test device is brought to room temperature before opening. The test is not stable out off the expiration date printed on the sealed pouch. Do not freeze the kit or expose the kit over 30°C.

WARNINGS AND PRECAUTIONS

- 1. For professional *in vitro* diagnostic use only.
- 2. Do not use it if the tube/pouch is damaged or broken.
- 3. Test is for single use only. Do not re- use under any circumstances.
- 4. Handle all specimens as if they contain infectious agents. Observe established standard procedure for proper disposal of specimens
- 5. Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are assay.
- 6. Humidity and temperature can adversely affect results

SPECIMEN COLLECTION

Collect sufficient quantity of feces (1-2 mL or 1-2 g) in a clean, dry specimen collection container to obtain maximum antigens (if present). Best results will be obtained if the assay is performed within 6 hours after collection. Specimen collected may be stored for 3 days at 2-8°C if not tested within 6 hours. For long term storage, specimens should be kept below -20°C.


To process fecal specimens:

• For Solid Specimens:

Unscrew the cap of the specimen collection tube, then randomly stab the specimen collection applicator into the fecal specimen in at least 3 different sites to collect approximately 50 mg of feces (equivalent to 1/4 of a pea). Do not scoop the fecal specimen.

• For Liquid Specimens:


Hold the dropper vertically, aspirate fecal specimens, and then transfer 2 drops (approximately 80 µL) into the specimen collection tube containing the dilution buffer. Screw on and tighten the cap onto the specimen collection tube, then shake the specimen collection tube vigorously to mix the specimen and the dilution buffer. Leave the tube alone for 2 minutes.

TEST PROCEDURE

1. Remove the test device from its foil pouch by tearing along the notch and use it as soon as possible.
2. Specimen collection. See also specimen collection.
3. Holding the sample collection device upright, carefully break off the tip of collection device.
4. Squeeze 2 drops (~80 µL) of the sample solution in the sample well of the cassette, as in the illustration.
5. Read the test results in 10 minutes. It is important that the background is clear before the result is read. Do not read results after 10 minutes. To avoid confusion, discard the test device after interpreting the result.

INTERPRETATION OF RESULTS

H. pylori Ag Rapid Test Cassette (Feces)

Positive: Two lines appear. One colored line should be in the control line region (C) and another apparent colored line should be in the test line region (T).

Negative: One colored line appears in the control line region(C). No line appears in the test line region (T).

Invalid: Control line fails to appear.

QUALITY CONTROL

A procedural control is included in the test. A colored line appearing in the control line region (C) is an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.

Control standards are not supplied with this kit; however, it is recommended that positive and negative controls be tested as a good laboratory practice to confirm the test procedure and to verify proper test performance.

LIMITATIONS

1. The Assay Procedure and the Assay Result Interpretation must be followed closely when testing the presence of H. Pylori antigen in feces from individual subjects. Failure to follow the procedure may give inaccurate results.

2. H. pylori Ag Rapid Test Cassette (Feces) is limited to the qualitative detection of H. Pylori antigen in feces. The intensity of the test band does not have linear correlation with the antigen titer in the specimen.

3. A negative result for an individual subject indicates absence of detectable H. Pylori antigen. However, a negative test result does not preclude the possibility of exposure to or infection with H. Pylori.

4. A negative result can occur if the quantity of the H. Pylori antigen present in the specimen is below the detection limits of the assay, or the antigen that are detected are not present during the stage of disease in which a sample is collected.

5. The results obtained with this test should only be interpreted in conjunction with other diagnostic procedures and clinical findings.

PERFORMANCE CHARACTERISTICS

A study was performed with 165 patient feces samples including both symptomatic gastrointestinal disorders and samples from non-symptomatic patients and 100 normal feces samples. Comparison for all subjects with H. pylori Ag Rapid Test Cassette (Feces) and reference ELISA kit is showed in the following table:

Method		EIA		Total Results
H.P Test Cassette	Results	Positive	Negative	
	Positive	163	0	163
	Negative	2	100	102
Total Results		165	100	265

Relative sensitivity: 98.8%

Relative specificity: 100%

Accuracy:98.9%

REFERENCE

1. Marshall,B.J.et.al. Pyloric Campylobacter infection and gastroduodenal disease. Med. J. Australia.149:439-44, 1985.

2. Marshall,B.J.et.al. Prospective double-blind trial of duodenal ulcer relapse after eradication of Campylobacter pylori. Lancet. Dec.1437-42,1988.

3. Megraud,F.et.al. Seroepidemiology of Campylobacter pylori infection in virious populations J.Clin.Microbiology. 27:1870-3,1989.

4. Soll,A.H. Pathogenesis of peptic ulcer and implications for therapy. New England J. Med.322:909-916,1990.

5. Parsonnet,J.et.al. Helicobacter pylori infection and the risk of gastric carcinoma. New England J.Med. 325:1127-31,1991.

6. Ansong,R. et.al. Evaluation of techniques for isolation, subcultivation and preservation of Helicobacter pylori. J.Clin.Micro. 29:51-53,1991.

7. Pronovost,A.P.et.al. Evaluation of a new immunodiagnostic assay for Helicobacter pylori antibody detection: Correlation with histopathological and microbiological results. J.Clin.Microbiol.32:46-50,1994.

INDEX OF SYMBOLS

	Consult instructions for use		Tests per kit		Authorized Representative
	For <i>in vitro</i> diagnostic use only		Use by		Do not reuse
	Store between 2~30°C		Lot Number		Catalog#

Zhejiang Orient Gene Biotech Co.,Ltd
Address: 3787#, East Yangguang Avenue, Dipu Street,
Anji 313300, Huzhou, Zhejiang, China.
TEL: +86-572-5226111 FAX: +86-572-5226222
Website: www.orientgene.com

Shanghai International Holding Corp. GmbH (Europe)
Add: Eiffestrasse 80, 20537 Hamburg, Germany

GCHP-602a

Revision Date: 2022-03-08
B20435-03

Malaria P.f./P.v. Ag Rapid Test Cassette (Whole Blood)

INTENDED USE

The Malaria P.f./P.v. Ag Rapid Test Cassette (Whole Blood) is a rapid lateral flow chromatographic immunoassay for the simultaneous detection and differentiation of Malaria P.falciparum specific histidine rich protein-2 (Pf HRP-II) and Malaria P.vivax specific lactate dehydrogenase (Pv-LDH) in human blood specimen as an aid in the diagnosis of Malaria infection. It is for *In-Vitro* Diagnostic use only.

INTRODUCTION

Malaria is a serious, sometimes fatal, parasitic disease characterized by fever, chills, and anemia and is caused by a parasite that is transmitted from one human to another by the bite of infected Anopheles mosquitoes. There are four kinds of malaria that can infect humans: Plasmodium falciparum, P. vivax, P. ovale, and P. malariae. In humans, the parasites (called sporozoites) migrate to the liver where they mature and release another form, the merozoites. The disease now occurs in more than 90 countries worldwide, and it is estimated that there are over 500 million clinical cases and 2.7 million malaria-caused deaths per year. At the present, malaria is diagnosed by looking for the parasites in a drop of blood. Blood will be put onto a microscope slide and stained so that the parasites will be visible under a microscope.

PRINCIPLE

The Malaria P.f./P.v. Ag Rapid Test Cassette (Whole Blood) contains a membrane, which is precoated with mouse monoclonal antibodies specific to HRP-II of P. falciparum on test line Pf region and with mouse monoclonal antibodies specific to lactate dehydrogenase of P.vivax species on test line Pv region respectively. Conjugate pad is dispensed with monoclonal antibodies conjugated to colloidal gold, which are specific to P.falciparum histidine rich protein-2 (Pf HRP-II) and specific to the lactate dehydrogenase of P.vivax.

During the assay, an adequate volume of the blood specimen is dispensed into the sample well (S) of the test cassette, a lysis buffer is added to the buffer well (B). The buffer contains a detergent that lyses the red blood cells and releases various antigens, which migrate by capillary action across the strip held in the cassette. Pv-LDH if presents in the specimen will bind to the Pv-LDH-gold conjugates. The immunocomplex is then captured on the membrane by the pre-coated anti-Pv-LDH antibody, forming a burgundy colored Pv band, indicating a Pv positive test result.

Alternatively, pHRP-II if presents in the specimen will bind to the pHRP-II-gold conjugates. The immunocomplex is then captured on the membrane by the pre-coated anti-pHRP-II antibodies, forming a burgundy colored Pf band, indicating a Pf positive test result.

Absence of any T bands suggests a negative result. The test contains an internal control (C band) which should exhibit a burgundy colored band of the immunocomplex of goat anti- mouse IgG I mouse IgG (anti-Pv-LDH and anti-pHRP-II)-gold conjugates regardless of the color development on any of the T bands. Otherwise, the test result is invalid and the specimen must be retested with another device.

MATERIALS SUPPLIED

25 Sealed pouches each containing a test cassette, a dropper and a desiccant
1 Buffer, 7.0 mL
1 Package insert

MATERIAL REQUIRED BUT NOT PROVIDED

1. Clock or timer
2. Collection by venipuncture: collection tube (containing EDTA, citrate or heparin)
3. Collection using a lancet: sterile lancet

STORAGE AND STABILITY

All reagents are ready to use as supplied. Store unused test device unopened, preferably at 2°C-30°C. Do not expose the kit over 30°C. Do not freeze the kit. Ensure that the test device is brought to room temperature before opening. The test device is stable through the expiration date printed on the sealed pouch if it is stored at 2°C-30°C.

WARNINGS AND PRECAUTIONS

1. For professional *in vitro* diagnostic use only. Do not use after expiration date.
2. The instruction must be followed exactly to get accurate results. Failure to follow the insert gives inaccurate test results.
3. Do not eat, drink or smoke in the area where the specimens or kits are handled.
4. Handle all specimens as if they contain infectious agents. Observe established precautions against microbiological hazards throughout testing and follow the standard procedures for proper disposal of specimens.

5. Hemolized blood may be used for the testing, but do not take precipitants.
6. Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are being tested.
7. Humidity and temperature can adversely affect results.
8. Do not perform the test in a room with strong air flow, ie. an electric fan or strong airconditioning.

SPECIMEN COLLECTION

Collection by venipuncture:

- 1) Collect whole blood into a collection tube (containing EDTA, citrate or heparin) by venipuncture.
- 2) If specimens are not immediately tested, they should be refrigerated at 2-8°C. For storage periods greater than three days, freezing is recommended. They should be brought to room temperature prior to use. Using the specimen after long-term storage of more than three days can cause non-specific reaction.
- 3) When stored at 2-8°C, the whole blood sample should be used within three days.

Collection using a lancet:

- 1) Clean the area to be lanced with an alcohol swab.
- 2) Squeeze the end of the fingertip and pierce with a sterile lancet.
- 3) Wipe away the first drop of blood with sterile gauze or cotton.
- 4) Using the dropper provided, while gently squeezing the tube, immerse the open end in the blood drop and then gently release the pressure to draw blood into the dropper.

TEST PROCEDURE

Allow the test device, specimen, buffer, and/or controls to equilibrate to room temperature (15-30°C) prior to testing.


1. Remove the test cassette from the foil pouch and use it as soon as possible. Best results will be obtained if the assay is performed within one hour.
2. Place the test cassette on a clean and level surface. Be sure to label the device with specimen's ID number.
3. With a 5 µL mini plastic dropper provided, draw whole blood specimen to exceed the specimen line as showed in the following image and then transfer drawn whole blood into the sample well (S). Then add 3 drops (about 120 µL) of Lysis Buffer to the buffer well (B) immediately.

Note: Practice a few times prior to testing if you are not familiar with the mini dropper. For better precision, transfer specimen by pipette capable to deliver 5 µL of volume.

4. Set up timer.

If preferred, after 5 minutes of adding specimen and buffer, you may add one more drop of Lysis Buffer to help the background become clearer.

5. Results can be read in 20 to 30 minutes. It may take more than 20 minutes to have the background become clearer. Don't read results after 30 minutes. To avoid confusion, discard the test cassette after interpreting the result.

Malaria P.f./P.v. Ag Rapid Test Cassette (Whole Blood)

INTERPRETATION OF RESULTS

(Please refer to the illustration above)

POSITIVE:

P.f Positive: One line appears in the control region, and one line appears in P.f. line region.

P.v Positive: One line appears in the control region and one line appears in Pv line region.

P.f and P.v Positive: One line appears in the control region, one line appears in Pv line region and one line appears in P.f. line region.

NEGATIVE: Only one colored line appears in the control region.

INVALID: Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test device. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

QUALITY CONTROL

Internal procedural controls are included in the test. A colored line appearing in the control region (C) is an internal procedural control. It confirms sufficient specimen volume and correct procedural technique. Control standards are not supplied with this kit; however, it is recommended that positive and negative controls be tested as a good laboratory practice to confirm the test procedure and to verify proper test performance.

LIMITATIONS

- The Malaria P.f./P.v. Ag Rapid Test Cassette (Whole Blood) is for *in vitro* diagnostic use only. This test should be used for the detection of P.f and P.v antigens in whole blood specimens only. Neither the quantitative value nor the rate of increase in P.f and P.v concentration can be determined by this qualitative test.
- The Malaria P.f./P.v. Ag Rapid Test Cassette (Whole Blood) will only indicate the presence of antigens of P.f and / or P.v in the specimen and should not be used as the sole criterion for the diagnosis of malaria infection.
- As known relevant interference, haemolytic samples, rheumatoid factors-contained samples and lipaemic, icteric samples can lead to impair the test results.
- The test is limited to the detection of antigen to Malaria Plasmodium sp. Although the test is very accurate in detecting HRP-II specific to P.f or pLDH specific to P.v, a low incidence of false results can occur. Other clinically available tests are required if questionable results are obtained.
- If the test result is negative and clinical symptoms persist, additional testing using other clinical methods is recommended. A negative result does not at any time preclude the possibility of malaria infection.
- As with all diagnostic tests, a definitive clinical diagnosis should not be based on the results of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated.

PERFORMANCE CHARACTERISTICS

1. Clinical Performance for P.f Ag test:

A total of 352 samples from susceptible subjects were tested by the Malaria P.f./P.v. Ag Rapid Test Cassette (Whole Blood) and by thick blood smear test.

Method		Smear Test		Total Results
Malaria Pf/Pv Ag Rapid Test	Results	Positive	Negative	
	Positive	50	4	54
	Negative	0	298	298
Total Results		50	302	352

Relative Sensitivity: 100%

Relative Specificity: 98.7%

Overall Agreement: 98.9%

2. Clinical Performance for P.v Ag test:

A total of 289 samples from susceptible subjects were tested by the Malaria P.f./P.v. Ag Rapid Test Cassette (Whole Blood) and by thick blood smear test.

Method		Smear Test		Total Results
Malaria Pf/Pv Ag Rapid Test	Results	Positive	Negative	
	Positive	63	3	66
	Negative	0	223	223
Total Results		63	226	289

Relative Sensitivity: 100%

Relative Specificity: 98.7%

Overall Agreement: 99.0%

3. Precision: Within-run and between-run have been determined by the testing 10 replicates of four specimens: a negative, a low positive, a medium positive and a strong positive. All values were correctly identified 100% of the time.

4. Interference: To evaluate the interference of Malaria P.f./P.v. Ag Rapid Test Cassette (Whole Blood) with known relevant interfering specimens, the haemolytic samples, rheumatoid factors-contained samples and lipaemic, icteric samples were investigated. In these studies, those specimens did not interfere with the Malaria P.f./P.v. Ag Rapid Test Cassette (Whole Blood).

REFERENCE

- Leonard K. Basco, Frederique Marquet, Michael M. Makler, and Jacques Le Bras.: Plasmodium falciparum and Plasmodium vivax: Lactate Dehydrogenase Activity and its Application for *in vitro* Drug Susceptibility Assay. Experimental Parasitology 80, 260-271 (1995)
- David L. Vander Jagt, Lucy A. Hunsaker and John E. Heidrich : Partial Purification and Characterization of Lactate Dehydrogenase from Plasmodium falciparum. Molecular and Biochemical Parasitology, 4 (1981) 255-264
- David J. Bzik, Barbara A. Fox and Kenneth Gonyer : Expression of Plasmodium falciparum lactate dehydrogenase in Escherichia coli Molecular and Biochemical Parasitology, 59(1993) 155-166
- Histidine-Rich Protein II: a Novel Approach to Malaria Drug Sensitivity Testing ANTIMICROBIAL AGENTS AND CHEMOTHERAPY, June 2002, p. 1658©1664 Vol. 46, No. 6

INDEX OF SYMBOLS

	Consult instructions for use		Tests per kit		Authorized Representative
	For <i>in vitro</i> diagnostic use only		Use by		Do not reuse
	Store between 2~30°C		Lot Number		Catalog#

Zhejiang Orient Gene Biotech Co., Ltd
Address: 3787#, East Yangguang Avenue, Dipu Street,
Anji 313300, Huzhou, Zhejiang, China
Tel: +86-572-5226111 Fax: +86-572-5226222
Website: www.orientgene.com

QARAD BV
Cipalstraat 3, 2440 Geel BELGIUM

GCMAL(pf/pv)-402a

Vitrotest SARS-CoV-2 IgM

ELISA test kit for the qualitative determination of IgM class antibodies to coronavirus SARS-CoV-2

REF ELO34-96

Σ 96

IVD

1. INTENDED USE

The test kit Vitrotest SARS-CoV-2 IgM is an enzyme linked immunosorbent assay (ELISA) for the qualitative determination of IgM class antibodies to N (nucleocapsid) and S (spike) antigens of SARS-CoV-2 coronavirus in human serum or plasma.

The test kit might be applied for the ELISA using both automatic pipettes and standard equipment as well as open system automated ELISA analyzers.

2. CLINICAL VALUE

COVID-19 is an infectious disease caused by a new SARS-CoV-2 coronavirus which had not previously been detected in humans.

The viral infection leads to the development of a respiratory flu-like disease with symptoms such as cough and fever. In more severe cases pneumonia can develop. The average incubation period of the COVID-19 is 6.5 days, but it can range from 3 to 21 days.

Coronavirus SARS-CoV-2 is an RNA-containing virus with a characteristic envelope with spikes in the form of a "corona". The main structural proteins of the virus include the nucleocapsid protein and transmembrane S (spike) protein with a receptor-binding domain (RBD), which binds to human cell receptors ACE2, causing infection of the mucosal epithelial cells. Both proteins are highly immunogenic antigens for humans.

Specific antibodies (IgM, IgG, IgA) to the SARS-CoV-2 N- and S-proteins appear on the 7-11th day from the moment the virus enters / contacts the body of an infected person. Anti-SARS-CoV-2 IgM and IgA antibodies can be detected as early as 4 days after the first symptoms of the disease. N- and S-specific IgM in not requiring intensive care patients may peak in the second week after COVID-19 symptoms appear. According to the literature, the simultaneous determination of N- and S-specific IgM antibodies makes it possible to identify up to 75 % of patients infected with SARS-CoV-2 in the first week of clinical manifestations.

At the same time, the level of specific IgG continues to grow and can be detected in almost 100 % of patients (except for those with immunosuppression) by the 20th day from the onset of COVID-19 symptoms.

3. PRINCIPLE OF THE TEST

Determination of IgM antibodies specific to SARS-CoV-2 in the test kit Vitrotest SARS-CoV-2 IgM is based on "IgM-capture" solid phase ELISA in a two-step incubation procedure. Microwells are coated with monoclonal antibodies to human IgM. During the first incubation step, IgM antibodies will be bound to the solid phase precoated monoclonal antibodies. The wells are washed to remove unbound components. Recombinant antigens of SARS-CoV-2 (N- and S- antigens), which is conjugated to horseradish peroxidase (HRP), is added next and binds to the specific IgM captured on the solid phase. Unbound components are removed by washing. Antigen-antibody complexes are revealed by addition of chromogen solution containing 3,3',5,5'- tetramethylbenzidine (TMB) and hydrogen peroxide. After 30 minutes the reaction has been stopped, the absorbance values are read using a spectrophotometer at 450/620-695 nm. The colour intensity is proportional to the amount of the specific antibodies present in the sample.

4. MATERIALS AND EQUIPMENT

4.1. Composition of the test kit

ELISA STRIPS	1x96 wells	Microplate (12 strips x 8 wells) Each well is coated with monoclonal antibodies to human IgM. The wells can be separated.
CONTROL +	1x0.5 ml	Positive control Solution of specific monoclonal immunoglobulins with preservative (pink).
CONTROL -	1x0.5 ml	Negative control Buffer solution with detergent and preservative (yellow).
CONTROL CUT-OFF	1x0.5 ml	Cut-off control Solution of specific monoclonal immunoglobulins with preservative (orange).

SAMPLE DILUENT	1x12 ml	Sample diluent Buffer solution with detergent and preservative (violet).
CONJUGATE SOLUTION	1x12 ml	Conjugate solution Buffer solution of recombinant antigens of SARS-CoV-2 conjugated to HRP with stabilizers and preservative (green), ready to use.
TMB SOLUTION	1x12 ml	TMB solution TMB, H ₂ O ₂ , stabilizers, preservative (colourless), ready to use.
WASH TWEEN 20X	1x50 ml	Washing solution Tw20 (20x concentrate) 20X concentrated of PBS buffer with Tween-20 and NaCl (colourless)
STOP SOLUTION	1x12 ml	Stop Solution 0,5 mol/l H ₂ SO ₄ (colourless), ready to use.

Adhesive films (2), sera identification plan (1), instruction for use and certificate of analysis.

4.2. Material required but not provided

- variable volume automatic pipettes (10µl–1000 µl) and disposable pipette tips;
- plate reader (single wavelength 450 nm or dual wavelength 450/620–695 nm);
- volumetric laboratory glassware (10–1000 ml);
- distilled/DI water;
- incubator thermostatically controlled at 37 °C;
- automatic/semiautomatic plate washer;
- appropriate waste containers for potentially contaminated materials;
- timer;
- absorbent paper;
- disposable gloves;
- disinfectants;
- protective clothes.

5. PRECAUTIONS AND SAFETY

5.1. Precautions

The ELISA assays are time and temperature sensitive. Strictly follow the test procedure and do not modify it.

- do not use expired reagents;
- do not use for analyses and do not mix reagents from different lots or from test kits of different nosology as well as other manufacturer's reagents with Vitrotest kits;

*Note: it is possible to use **WASH TWEEN 20X**, **TMB SOLUTION** and **STOP SOLUTION** from other Vitrotest ELISA kits.*

- close reagents after use only with appropriate caps;
- control the filling and full aspiration of the solution in the wells;
- use a new tip for each sample and reagent;
- avoid exposure of kit reagents to direct sunlight;
- **TMB SOLUTION** must be colourless before use. If **TMB SOLUTION** is blue or yellow it cannot be used. Avoid any contact of **TMB SOLUTION** with metals or metal ions. Use glassware thoroughly washed and rinsed with distilled/DI water;
- never use the same glassware for **CONJUGATE SOLUTION** and **TMB SOLUTION**;

The manufacturer is not responsible or liable for any incorrect results and/or incidents taking place as a result of any violation of the instruction. The manufacturer is not responsible for visual readings of samples (without using a plate reader).

5.2. Safety

- all components of test kit are intended for *in vitro* diagnostic use only;
- all materials of human or animal origin should be regarded and handled as potentially infectious;
- the ELISA is only designed for qualified personnel;
- disposable gloves and safety glasses must be worn at all times while performing analysis;
- never eat, drink, smoke or apply cosmetics in the assay laboratory;
- never pipette solutions by mouth;
- controls do not contain of human origin components;
- avoid contact with **STOP SOLUTION** containing 0.5 mol/l H₂SO₄. It may cause skin irritation and burns;
- some components of the test kit contain low concentrations of harmful compounds and could cause irritation of the skin and the mucosa. In the case of contact of **TMB SOLUTION**, **STOP SOLUTION** or **CONJUGATE SOLUTION** with skin or mucosa, the place of contact should be immediately rinsed with large amounts of water;

- in case of spilling of solutions that do not contain acid, e.g. sera, rinse the surface with disinfectant, then dry it with absorbent paper. In other case acid first must be neutralized by sodium bicarbonate and then wiped out as described above;
- for information on hazardous substances included in the kit please refer to Safety Data Sheets. Safety Data Sheets for this product are available upon request.

5.3. Waste treatment

Patient specimens, controls, and incubated microplate strips should be treated as infectious waste, residues of chemicals and preparations are generally considered as hazardous waste. The disposal of this kind of waste is regulated through national and regional laws and regulations. Contact your local authorities or waste management companies which will give advice on how to dispose hazardous waste.

6. STORAGE AND STABILITY

Reagents are stable until stated expiration date on the label when stored refrigerated (2-8 °C). Do not freeze. The kit should be shipped at 2-8 °C. Single transportation at the temperature up to 23 °C for two days is acceptable.

After the first opening of the packaging, the components of the ELISA kits are stable within 3 months, except for those specified in p. 8 of this Instruction.

7. SPECIMEN COLLECTION

The fresh serum or plasma (EDTA, lithium-heparin, sodium citrate) samples can be stored for 3 days at 2-8 °C, or frozen for longer periods at -20 – -70 °C. Frozen samples must be thawed and kept at room temperature for at least 30 min before use. Do not use preheated samples. Mix thawed samples thoroughly to homogeneity. Avoid repeated freezing/thawing. Samples containing aggregates must be clarified by centrifugation (3000 rpm for 10-15 min). Do not use hyperlipemic, hyperhaemolysed or contaminated by microorganisms serum specimens. The presence of bilirubin up to concentration of 0.21 mg/ml (361.8 µmol/l), haemoglobin up to concentration of 10 mg/ml and triglycerides up to concentration of 10 mg/ml (11.3 mmol/l) are allowed.

8. REAGENT PREPARATION

It is very important to keep all test components for at least 30 min at room temperature (18-25 °C) before the assay!

8.1. [ELISA STRIPS] preparation

Before opening the bag with [ELISA STRIPS], keep it at room temperature for 30 min to avoid water condensation inside the wells. Open the vacuum bag and take out the necessary number of the wells. Once opened the bag with the remaining strips must be **resealed with zip-lock** immediately and kept refrigerated at 2-8 °C for no more than 3 months.

8.2. Washing solution preparation

Check the [WASH TWEEN]20X for the presence of salt crystals. If crystals have formed, re-solubilise by warming at 37°C, until crystals dissolve (15-20 min). Dilute the [WASH TWEEN]20X 1:20 (1+19) with distilled/DI water before use. For example, 4 ml concentrate + 76 ml water is sufficient for 8 wells. Once diluted it is stable at 2-8 °C for 7 days.

9. ASSAY PROCEDURE

- 9.1. Take out from the protective bag the support frame and the necessary number of the wells [ELISA STRIPS] (the number of specimens + 4 for controls). Place the wells into the frame. Wells with the controls must be included in every test.
- 9.2. Complete the sera identification plan.
- 9.3. Prepare washing solution (see 8.2).
- 9.4. Dispense 90 µl of [SAMPLE DILUENT] into each well [ELISA STRIPS].
- 9.5. Dispense 10 µl of controls and patient samples into the wells in the following order: A1 – [CONTROL +], B1, C1 – [CONTROL CUT-OFF], and D1 – [CONTROL -], other wells – patient samples. Mix gently to avoid foaming. The colour of the sample diluent changes from violet to blue.
- 9.6. Cover strips with an adhesive film and incubate for 30 min at 37 °C.
- 9.7. Remove and discard the adhesive film and wash all wells 5 times with automatic washer or 8-channel pipette as follows:
 - aspirate the contents of all wells into a liquid waste container and add immediately a minimum of 300 µl of diluted washing solution to each well;
 - soak each well for 30 s between each wash cycle;
 - aspirate again. The residual volume must be lower than 5 µl;
 - repeat the washing step 4 times;
 - after the final washing cycle, turn down the plate onto an absorbent paper and tap it to remove any residual buffer.

- 9.8. Dispense 100 µl of [CONJUGATE SOLUTION] per well. Cover strips with a new adhesive film, incubate for 30 min at 37 °C.
- 9.9. Remove and discard the adhesive film and wash all wells five times as described above (see 9.7).
- 9.10. Dispense 100 µl [TMB SOLUTION] into all wells. Do not touch the walls and bottoms of the wells to avoid contamination.
- 9.11. Incubate the strips for 30 min at room temperature (18-25 °C) in the dark. Do not use adhesive film in this step.
- 9.12. Add 100 µl [STOP SOLUTION] to each well in the same order and at the same rate as for [TMB SOLUTION].
- 9.13. Read the optical density (OD) of the wells at 450/620-695 nm using a microplate reader within 5 min after adding the [STOP SOLUTION]. Pay attention to the cleanness of the plate bottom and absence of bubbles in the wells before reading.

Measurement in the single-wave procedure at 450 nm is possible. Reserve blank well to adjust spectrophotometer in such analysis. Only [TMB SOLUTION] and [STOP SOLUTION] must be added in blank well.

10. CALCULATION AND INTERPRETATION OF RESULTS

10.1. Validation of the test

The test run may be considered valid provided the following criteria are met:

[CONTROL +]	OD ≥ 1.2
[CONTROL CUT-OFF]	OD in a range 0.25-0.65
[CONTROL -]	OD ≤ 0.150

If one of the control cut-off absorbances does not match the above criteria, this value should be discarded and a mean value should be calculated using the remaining cut-off value. If both control cut-off absorbance do not meet the criteria, the test is invalid and must be re-tested.

10.2. Calculation of results

The cut-off (CO) is the mean optical density (OD) of the wells containing [CONTROL CUT-OFF]:

$$CO = (OD_{\text{CONTROL CUT-OFF } 1} + OD_{\text{CONTROL CUT-OFF } 2})/2;$$

The sample result is reported as a Ratio:

$$\text{Ratio}_{\text{sample}} = OD_{\text{sample}}/CO, \quad OD_{\text{sample}} - \text{optical density of the well containing sample}$$

10.3. Interpretation of results

Ratio _{sample} > 1.1	POSITIVE
0.9 ≤ Ratio _{sample} ≤ 1.1	DOUBTFUL*
Ratio _{sample} < 0.9	NEGATIVE

* If the result is doubtful, repeat the test. If it remains doubtful, collect a new serum sample.

11. PERFORMANCE CHARACTERISTICS

11.1. Specificity and sensitivity

To assess the diagnostic characteristics of the Vitrotest SARS-CoV-2 IgM test kit we used Anti-SAR-SCoV-2 Verification Panel for Serology Assays (manufactured by NIBSC, UK), which contains 23 blood plasma samples from COVID-19 convalescents and 14 negative blood plasma samples. The results obtained in the Vitrotest SARS-CoV-2 IgM test kit completely coincide with the panel passport data.

In a comparative studies of the sensitivity of the Vitrotest SARS-CoV-2 IgM test kit on 85 samples, which were determined as positive in a commercial analogue recommended by the FDA, the percentage of agreement of the obtained results was 96.5 %.

To determine the sensitivity of the Vitrotest SARS-CoV-2 IgM test kit in the early stages of SARS-CoV-2 infection, 1117 blood sera obtained 5-15 days from the onset of clinical manifestations of COVID-19 were tested. Specific IgM antibodies were detected in 972 samples.

The specificity of the Vitrotest SARS-CoV-2 IgM test kit on 706 human blood serum samples obtained during the first half of 2019 (before the start of the COVID-19 pandemic) was 99.4 %.

11.2. Accuracy

Intra assay repeatability

Coefficient of variation (CV) was calculated by measuring 2 samples with various specific antibody levels in 32-replicate determinations using 1 lot of the test kit.

Serum No.	OD	Ratio	CV, %
336	0.755	3.45	2.7
345	2.408	11.02	1.6

Inter assay reproducibility

Coefficient of variation (CV) was calculated by measuring 2 samples with various specific antibody levels in 4 ELISA performances during 4 days, in 8-replicate determinations.

Serum No.	OD	Ratio	CV, %
336	0.758	3.46	2.4
345	2.409	10.99	1.2

12. LIMITATIONS OF THE PROCEDURE

If the test sample was obtained in the first days after infection IgM antibodies may not be detected. Therefore, a negative result of the test for IgM antibodies to SARS-CoV-2 does not exclude the infection of the patient with the virus. In the presence of clinical manifestations of the disease it is recommended to repeat the test at least in 1-2 weeks, for example when using test kits Vitrotest SARS-CoV-2 IgM, Vitrotest SARS-CoV-2 IgG and Vitrotest SARS-CoV-2 Total Ab.

Negative test results in immunosuppressed individuals should also be interpreted with caution.

The final diagnosis cannot be established solely on the basis of serological test results. The diagnosis should take into consideration clinical history, symptomatology, as well as the results of other laboratory tests (including PCR).

13. TROUBLESHOOTING

<i>Possible causes</i>	<i>Solutions</i>
<i>High background in all wells</i>	
Contaminated washer	Clean the washer head, then rinse it with 30 % ethanol and distilled water
Low quality water or contaminated water	Use distilled/DI with resistivity ≥ 10 M Ω -cm.
Using contaminated glassware	Use clean glassware
Using chlorine based disinfectants	Use disinfectants without chlorine
Using contaminated tips	Use new tips
Increased time of incubation or temperature regimen was changed	Follow incubation regimen according to instruction for use
<i>High background in a few wells</i>	
TMB solution was added more than once	Add TMB solution once
Pipette shaft was contaminated with conjugate solution	Clean the pipette; pipette the liquids carefully
One the channels of the washer was contaminated	Clean the washer channel, clean the washer
<i>OD of the positive control below normal</i>	
Conjugate solution/tmb solution was prepared improperly or not added	Run ELISA repeatedly, prepared conjugate solution / TMB solution properly
Reduced incubation time in one of the stages	Follow incubation regimen according to the instruction for use
<i>Visual colour intensity of the wells does not correspond to optical density</i>	
The optical beam or another component of the reader is misaligned or malfunctioning	Test the absorbance reader's performance

REFERENCE

1. Baoqing Sun , Ying Feng , Xiaoneng Mo et.al. Kinetics of SARS-CoV-2 specific IgM and IgG responses in COVID-19 patients.// *Emerging Microbes & Infections*, - 2020 - V.9 – p.940-948.
2. Juanjuan Zhao Jr., Quan Yuan et. al. Antibody responses to SARS-CoV-2 in patients of novel coronavirus disease 2019.// *Clinical Infectious Diseases.*, - 2020 Mar. 20 doi: 10.1093/cid/ciaa344.
3. Laboratory testing for coronavirus disease (COVID-19) in suspected human cases // WHO. Interim guidance 19 March 2020. WHO reference number: WHO/COVID-19/laboratory/2020.5.
4. Marco Cascella ; Michael Rajnik et.al. Features, Evaluation and Treatment Coronavirus (COVID-19). // NCBI Bookshelf. StatPearls Publishing; 2020 – P.16.
5. Patrick C. Y. Woo, Susanna K. P. Lau. et.al. Differential Sensitivities of Severe Acute Respiratory Syndrome (SARS) Coronavirus Spike Polypeptide Enzyme-Linked Immunosorbent Assay (ELISA) and SARS Coronavirus Nucleocapsid Protein ELISA for Serodiagnosis of SARS Coronavirus Pneumonia. // *J. Clin. Microbiol.*, - 2005 – V. 43 N.7 - p. 3054–3058.
6. Quan-Xin Long, Bai-Zhong Liu et. al. Antibody responses to SARS-CoV-2 in patients with COVID-19.// *Nature Medicine.*, - 2020 April 29 doi: 10.1038/s41591-020-0897-1.
7. Quan-xin Long, Hai-jun Deng et.al. Antibody responses to SARS-CoV-2 in COVID-19 patients: the perspective application of serological tests in clinical practice. // medRxiv preprint doi: <https://doi.org/10.1101/2020.03.18.20038018>.
8. Shu-Yuan Xiao, Yingjie Wu, Huan Liu. Evolving status of the 2019 novel corona virus infection: Proposal of conventional serologic assays for disease diagnosis and infection monitoring. // *J Med Virol.*, - 2020 - 92(5) – p.464-467. doi: 10.1002/jmv.25702. Epub 2020 Feb 17.
9. Wu, L.-P., Wang, N.-C. et.al. (2007). Duration of Antibody Responses after Severe Acute Respiratory Syndrome. // *Emerging Infectious Diseases*, - 2007 - 13(10) - p.1562-1564.

Catalogue number

Consult instructions for use

In vitro diagnostic medical device

Manufacturer

Caution

Contains sufficient for <n> tests

Temperature limit

Batch code

Use-by date

Date of manufacture

Keep away from sunlight

Signifies European conformity (CE) mark

Inst_SARS-CoV-2_IgM_EL034-96_V01_ENG
Edition 1st, 16.12.2021

Vitrotest Europe Sp. z o.o.
Krakowska str., 139-155, 50-428, Wrocław, Poland
tel.: +48 882 950 379,
e-mail: info@vitrotest.pl, www.vitrotest.pl

Vitrotest SARS-CoV-2 IgM

ASSAY PROCEDURE

Keep all reagents and specimens for at least 30 min at 18-25 °C before use

Dispense 90 µl of [SAMPLE DILUENT] into the wells of [ELISA STRIPS].
(violet colour)

Dispense 10 µl of controls and samples into the wells in the following order:
A1 – [CONTROL +],
B1, C1 – [CONTROL CUT-OFF],
D1 – [CONTROL -],
E1 and other wells – patient samples
(colour changes from violet to blue)

Cover wells with an adhesive film, incubate for 30 min at 37 °C

Rinse the wells 5 times with diluted 1:20 (1+19) washing solution Tween-20 (300 µl per well)

Add 100 µl of [CONJUGATE SOLUTION] into the wells
(green colour)

Cover wells with an adhesive film, incubate for 30 min at 37 °C

Rinse the wells 5 times with diluted 1:20 (1+19) washing solution Tween-20 (300 µl per well)

Add 100 µl of [TMB SOLUTION] into the wells

Incubate for 30 min in the dark at 18-25 °C

Add 100 µl of [STOP SOLUTION]
(colour changes from blue to yellow)

Determine the optical density (OD) at 450/620-695nm

CALCULATION

$$CO = (OD_{\text{CONTROL CUT-OFF 1}} + OD_{\text{CONTROL CUT-OFF 2}}) / 2;$$

$$\text{Ratio}_{\text{sample}} = OD_{\text{sample}} / CO$$

INTERPRETATION

$\text{Ratio}_{\text{sample}} > 1.1$	POSITIVE
$0.9 \leq \text{Ratio}_{\text{sample}} \leq 1.1$	DOUBTFUL
$\text{Ratio}_{\text{sample}} < 0.9$	NEGATIVE

Vitrotest SARS-CoV-2 IgG QuantiSpike

ELISA test kit for the quantitative determination of IgG class antibodies to coronavirus SARS-CoV-2 Spike protein

REF EL040-96

IVD

Σ
96

1. INTENDED USE

The test kit Vitrotest SARS-CoV-2 IgG QuantiSpike is an enzyme linked immunosorbent assay (ELISA) for the quantitative determination of IgG class antibodies to SARS-CoV-2 Spike protein synthesized in humans due to the disease or vaccination in serum or plasma.

The test kit might be applied for the ELISA using both automatic pipettes and standard equipment as well as open system automated ELISA analyzers.

2. CLINICAL VALUE

COVID-19 is an infectious disease caused by a new SARS-CoV-2 coronavirus which had not previously been detected in humans.

The viral infection leads to the development of a respiratory flu-like disease with symptoms such as cough and fever. In more severe cases pneumonia can develop. The average incubation period of the COVID-19 is 6.5 days, but it can range from 3 to 21 days.

SARS-CoV-2 is an RNA-virus with a specific envelope with spikes in the form of a "corona". The main structural proteins of the virus include envelope protein (E), membrane protein (M), spike (S) glycoprotein, and nucleocapsid (N) protein. S protein on the surface of the SARS-CoV-2 virion mediates the receptor recognition and cell membrane fusion with ACE2 molecules, which are mainly expressed on type II pneumocytes, colon and kidney epithelial cells. It contains three fragments, namely the ectodomain, the transmembrane domain and the short intracellular segment. The ectodomain consists of a receptor-binding subunit S1 containing the RBD domain and a fusion subunit (S2). During viral infection, S1 C-terminal domain binds to the extracellular peptidase (PD) domain of ACE2 to ensure that the virus attaches to the surface of the target cell. The S1 N-terminal domain binds to glycans causing the cleavage of S protein between S1 and S2 fragments by cellular proteases, which, in turn, initiates the fusion of viral and cell membranes by the S2 subunit.

Although most viral proteins are able to induce the production of specific antibodies after SARS-CoV-2 infection, and antibodies to N- and S-protein are widely used in the serological diagnosis of COVID-19, antibodies targeting viral S-protein are more noteworthy because they can block SARS-CoV-2 entry into the host cells. And since most vaccines induce antibodies to the spike protein the determination of IgG specific to this antigen also makes it possible to assess the presence of protective antibodies after the disease or vaccination against COVID-19.

3. PRINCIPLE OF THE TEST

Determination of IgG antibodies to S-protein of SARS-CoV-2 in the test kit Vitrotest SARS-CoV-2 IgG QuantiSpike is based on a solid phase, indirect ELISA in a two-step incubation procedure. Microwells are coated with the recombinant antigen, SARS-CoV-2 S-protein analogue. During the first incubation step, the specific antibodies to SARS-CoV-2 S-protein, if present in the sample, will be bound to the solid phase precoated antigens. The wells are washed to remove unbound antibodies. A secondary antibody (anti-IgG), which is conjugated to horseradish peroxidase (HRP), is added next and binds to the immune complexes on the solid phase. Unbound components are removed by washing. Antigen-antibody complexes are revealed by addition of chromogen solution containing 3,3',5,5'-tetramethylbenzidine (TMB) and hydrogen peroxide. After 15 min the reaction has been stopped, the absorbance values are read using a spectrophotometer at 450/620-695 nm. The colour intensity is proportional to the amount of the antibodies present in the sample.

Internal calibrators of the Vitrotest SARS-CoV-2 IgG QuantiSpike test kit are standardized according to the First WHO International Standard for anti-SARS-CoV-2 immunoglobulin (human) code: 20/136 (NIBSC, UK), which contains 1000 binding antibody units (BAU) per ml.

4. MATERIALS AND EQUIPMENT

4.1. Composition of the test kit

ELISA STRIPS	1x96 wells	Microplate (12 strips x 8 wells) Each well is coated with the recombinant antigens, SARS-CoV-2 S-protein analogues. The wells can be separated.
--------------	------------	--

PREDILUTION PLATE	1x96 wells	Microplate for preliminary dilution of sera
CAL 0	1x0.3 ml	Calibrator 0 Buffer solution with detergent and preservative (yellow).
CAL 25	1x0.3 ml	Calibrator 25 Solution of specific monoclonal immunoglobulins to S-protein (25 BAU/ml) with stabilizers and preservative (green).
CAL 50	1x0.3 ml	Calibrator 50 Solution of specific monoclonal immunoglobulins to S-protein (50 BAU/ml) with stabilizers and preservative (orange).
CAL 100	1x0.3 ml	Calibrator 100 Solution of specific monoclonal immunoglobulins to S-protein (100 BAU/ml) with stabilizers and preservative (pink).
CAL 200	1x0.3 ml	Calibrator 200 Solution of specific monoclonal immunoglobulins to S-protein (200 BAU/ml) with stabilizers and preservative (violet).
CONTROL +	1x0.5 ml	Positive control Solution of specific monoclonal immunoglobulins with a known concentration to S-protein with stabilizers and preservative (red).
SAMPLE PREDILUENT	1x50 ml	Sample prediluent Buffer solution with detergent and preservative (brown-green).
SAMPLE DILUENT	1x12 ml	Sample diluent Buffer solution with detergent and preservative (yellow).
CONJUGATE SOLUTION	1x12 ml	Conjugate solution Buffer solution of monoclonal antibodies to human IgG conjugated to HRP with stabilizers and preservative (violet), ready to use.
TMB SOLUTION	1x12 ml	TMB solution TMB, H ₂ O ₂ , stabilizers, preservative (colourless), ready to use.
WASH TWEEN 20X	1x50 ml	Washing solution Tw20 (20x concentrate) 20X concentrated of phosphate buffer with Tween-20 and NaCl (colourless).
STOP SOLUTION	1x12 ml	Stop Solution 0.5 mol/l H ₂ SO ₄ (colourless), ready to use.

Adhesive films (2), sera identification plan (1), blank calibration curve (1), instruction for use and certificate of analysis.

4.2. Material required but not provided

- Variable volume automatic pipettes (10 µl–1000 µl) and disposable pipette tips;
- plate reader (single wavelength 450 nm or dual wavelength 450/620–695 nm);
- volumetric laboratory glassware (10–1000 ml);
- distilled/DI water;
- incubator thermostatically controlled at 37 °C;
- automatic/semiautomatic plate washer;
- appropriate waste containers for potentially contaminated materials;
- timer;
- absorbent paper;
- disposable gloves;
- disinfectants;
- protective clothes.

5. PRECAUTIONS AND SAFETY

5.1. Precautions

The ELISA assays are time and temperature sensitive. Strictly follow the test procedure and do not modify it.

- do not use expired reagents;
- do not use for analyses and do not mix reagents from different lots or from test kits of different nosology as well as other manufacturer's reagents with Vitrotest kits;

Note: it is possible to use [WASH TWEEN|20X], [TMB SOLUTION], [STOP SOLUTION] and [SAMPLE PREDILUENT] from other lots.

- close reagents after use only with appropriate caps;
- control the filling and full aspiration of the solution in the wells;
- use a new tip for each sample and reagent;
- avoid exposure of kit reagents to direct sunlight;
- [TMB SOLUTION] must be colourless before use. If [TMB SOLUTION] is blue or yellow it cannot be used. Avoid any contact of [TMB SOLUTION] with metals or metal ions. Use glassware thoroughly washed and rinsed with distilled/DI water;
- never use the same glassware for [CONJUGATE SOLUTION] and [TMB SOLUTION].

The manufacturer is not responsible or liable for any incorrect results and/or incidents taking place as a result of any violation of the instruction. The manufacturer is not responsible for visual readings of samples (without using a plate reader).

5.2. Safety

- all components of test kit are intended for *in vitro* diagnostic use only;
- all materials of human or animal origin should be regarded and handled as potentially infectious;
- the ELISA is only designed for qualified personnel;
- disposable gloves and safety glasses must be worn at all times while performing analysis;
- never eat, drink, smoke or apply cosmetics in the assay laboratory;
- never pipette solutions by mouth;
- controls do not contain of human origin components;
- avoid contact with [STOP SOLUTION] containing 0.5 mol/l H₂SO₄. It may cause skin irritation and burns;
- some components of the test kit contain low concentrations of harmful compounds and could cause irritation of the skin and the mucosa. In the case of contact of [TMB SOLUTION], [STOP SOLUTION] or [CONJUGATE SOLUTION] with skin or mucosa, the place of contact should be immediately rinsed with large amounts of water;
- in case of spilling of solutions that do not contain acid, e.g. sera, rinse the surface with disinfectant, then dry it with absorbent paper. In other case acid first must be neutralized by sodium bicarbonate and then wiped out as described above;
- for information on hazardous substances included in the kit please refer to Safety Data Sheets. Safety Data Sheets for this product are available upon request.

5.3. Waste treatment

Patient specimens, calibrators, control and incubated microplate strips should be treated as infectious waste, residues of chemicals and preparations are generally considered as hazardous waste. The disposal of this kind of waste is regulated through national and regional laws and regulations. Contact your local authorities or waste management companies which will give advice on how to dispose hazardous waste.

6. STORAGE AND STABILITY

Reagents are stable until stated expiration date on the label when stored refrigerated (2-8 °C). Do not freeze. The kit should be shipped at 2-8 °C. Single transportation at the temperature up to 23 °C for two days is acceptable.

After the first opening of the packaging, the components of the ELISA kits are stable within 3 months, except for those specified in p. 8 of this Instruction.

7. SPECIMEN COLLECTION

The fresh serum or plasma (EDTA, lithium-heparin) samples can be stored for 3 days at 2-8 °C, or frozen for longer periods at -20 – -70 °C. Frozen samples must be thawed and kept at room temperature for at least 30 minutes before use. Do not use preheated samples. Mix thawed samples thoroughly to homogeneity. Avoid repeated freezing/thawing. Samples containing aggregates must

be clarified by centrifugation (3000 rpm for 10-15 min). Do not use hyperlipemic, hyperhaemolysed or contaminated by microorganisms serum specimens. The presence of bilirubin up to concentration of 0.21 mg/ml (361.8 $\mu\text{mol/l}$), haemoglobin up to concentration of 10 mg/ml and triglycerides up to concentration of 10 mg/ml (11.3 mmol/l) are allowed.

8. REAGENT PREPARATION

It is very important to keep all test components for at least 30 min at room temperature (18-25 °C) before the assay!

8.1. [ELISA STRIPS] preparation

Before opening the bag with the [ELISA STRIPS], keep it at room temperature for 30 min to avoid water condensation inside the wells. Open the vacuum bag and take out the necessary number of wells. Once opened, the bag with the remaining strips and desiccant must be *resealed with the zip-lock* immediately and kept refrigerated at 2-8 °C for no more than 3 months.

8.2. Washing solution preparation

Check the [WASH TWEEN 20X] for the presence of salt crystals. If crystals have formed, re-solubilise them by warming the vial at 37 °C, until crystals have been fully dissolved (15-20 min). Dilute the [WASH TWEEN 20X] 1:20 (1+19) with distilled/DI water before use. For example, 4ml concentrate + 76 ml water is sufficient for 8 wells. Once diluted it is stable at 2-8 °C for 7 days.

8.3. Predilution of samples, calibrators and positive control

Predilute patient samples, calibrators and positive control 1:10 with [SAMPLE PREDILUENT] immediately before test. Dispense 90 μl of [SAMPLE PREDILUENT] in the wells of [PREDILUTION PLATE], add 10 μl of samples, calibrators and positive control. Gently mix the content in the wells. After addition of the sample colour of the sample prediluent changes from brown-green to blue.

The procedure for dilution of samples, controls and calibrators should be carried out immediately before analysis.

9. ASSAY PROCEDURE

- 9.1. Take out from the protective bag the support frame and the necessary number of the wells [ELISA STRIPS] (the number of specimens, 1 well for the positive control and 5 wells for calibrators). Place the wells into the frame. Wells with calibrators and positive control must be included in every test.
- 9.2. Complete the sera identification plan.
- 9.3. Prepare washing solution (see 8.2.).
- 9.4. Predilute patient samples, calibrators and positive control (see 8.3).
- 9.5. Dispense 90 μl of [SAMPLE DILUENT] in the wells of [ELISA STRIPS].
- 9.6. Add 10 μl of prediluted (1:10) calibrators, positive control and patient samples to the wells in the following order: A1- [CAL 200], B1 - [CAL 100], C1 - [CAL 50], D1 - [CAL 25], E1 - [CAL 0] and F1 - [CONTROL +] respectively; other wells – patient samples. The final dilution in the wells is 1:100. Pipette gently to avoid foaming. The colour of the sample diluent changes from yellow to green.

Given the technical features of the equipment used for analysis, the order of dispense calibrators can be reversed: A1 - [CAL 0], B1 - [CAL 25], C1 - [CAL 50], D1 - [CAL 100], E1 - [CAL 200].

- 9.7. Cover strips with an adhesive film and incubate for 30 min at 37 °C.
- 9.8. Remove and discard the adhesive film and wash all wells 5 times with automatic washer or 8-channel pipette as follows:
 - aspirate the contents of all wells into a liquid waste container and add immediately a minimum of 300 μl of diluted washing solution to each well;
 - soak each well for 30 s between each wash cycle;
 - aspirate again. The residual volume must be lower than 5 μl ;
 - repeat the washing step 4 times;
 - after the final washing cycle, turn down the plate onto an absorbent paper and tap it to remove any residual buffer.
- 9.9. Dispense 100 μl [CONJUGATE SOLUTION] per well. Cover strips with a new adhesive film, incubate for 30 min at 37 °C.
- 9.10. Remove and discard the adhesive film and wash all wells five times as described above (see 9.8).
- 9.11. Dispense 100 μl [TMB SOLUTION] into all wells. Do not touch the walls and bottoms of the wells to avoid contamination.

- 9.12. Incubate the strips for 15 min at room temperature (18-25 °C) in the dark. Do not use adhesive film in this step.
- 9.13. Add 100 μ [STOP SOLUTION] to each well in the same order and at the same rate as for [TMB SOLUTION].
- 9.14. Read the optical density (OD) of the wells at 450/620-695 nm using a microplate reader within 5 min after adding the [STOP SOLUTION]. Pay attention to the cleanness of the plate bottom and absence of bubbles in the wells before reading.

Measurement in the single-wave procedure at 450 nm is possible. Reserve blank well to adjust spectrophotometer in such analysis. Only [TMB SOLUTION] and [STOP SOLUTION] must be added in blank well.

10. CALCULATION AND INTERPRETATION OF RESULTS


10.1. Validation of the test

The test run may be considered valid provided the following criteria are met:

[CAL 0]	OD \leq 0.100
[CAL 25]	OD \geq 0.120
[CAL 200]	OD \geq 1.500
[CONTROL +]	Within the concentration range on the tube label and on the certificate of analysis

10.2. Calculation of results

To determine the concentration of specific IgG antibodies in BAU/ml build a calibration curve. Plot the absorbances (OD) for 5 calibrators [CAL 0], [CAL 25], [CAL 50], [CAL 100] and [CAL 200] on Y axis and their corresponding concentrations in BAU/ml (0, 25, 50, 100 and 200 BAU/ml) respectively, on X axis. Using the absorbance values for each sample and positive control determine the corresponding concentration (BAU/ml) from the calibration curve.

Note that the calibration curve is intended for use only as an example, not for calculation of your results.

In case of OD of patient samples is higher of [CAL 200] the test result will be ">200 BAU/ml". Such samples can be retested at a dilution of 1:800. For diluted samples (1:800) multiply calculated results by 8.

$$\text{final concentration} = \text{concentration by the calibration curve} \times 8$$

If the optical density of the samples at a dilution of 1:800 is still above the value [CAL 200] it is recommended to retest such samples at a dilution of 1:1000 and 1:4000. In this case, the concentration of specific antibodies should be multiplied by the dilution rate of 10 and 40, respectively.

It is recommended to use computer software to read and calculate the obtained results.

10.3 Interpretation of results

Ig G concentration	Interpretation
> 25 BAU/ml	POSITIVE
20 - 25 BAU/ml	DOUBTFUL*
< 20 BAU/ml	NEGATIVE

* If the result is doubtful, repeat the test. If it remains doubtful, collect a new serum sample.

11. PERFORMANCE CHARACTERISTICS

11.1. Specificity and sensitivity

The specificity of the test kit Vitrotest SARS-CoV-2 IgG QuantiSpike using 352 samples of human sera obtained in the first half of 2019 (before the COVID-19 pandemic) was 100 %.

The sensitivity of the test kit Vitrotest SARS-CoV-2 IgG QuantiSpike, using 49 serum samples of COVID-19 convalescents (obtained 2-10 months after the disease), was 100 %. In addition, when testing 18 serum samples of vaccinated individuals all samples contained IgG antibodies to spike protein at a concentration of more than 1000 BAU/ml.

Diagnostic characteristics of the test kit Vitrotest SARS-CoV-2 IgG QuantiSpike were also evaluated on the verification panel Anti-SARSCoV-2 Verification Panel for Serology Assays code: 20/B770 (manufactured by NIBSC, UK), which consists of 23 characterized plasma samples of COVID-19 convalescents containing antibodies to SARS-CoV-2, and 14 characterized blood plasma samples not containing antibodies to SARS-CoV-2. The sensitivity and specificity of the test kit on this panel was 100%.

In the study of samples of the First WHO International Reference Panel for anti-SARS-CoV-2 immunoglobulin code: 20/268 (manufactured by NIBSC, UK) in the test kit Vitrotest SARS-CoV-2 IgG QuantiSpike obtained results matched the panel passport data.

11.2. Accuracy

Intra assay repeatability

Coefficient of variation (CV) was calculated by measuring 2 samples with various specific antibody levels in 32-replicate determinations using 1 lot of the test kit.

Serum No.	OD	Concentration, BAU/ml	CV, %
783	0.671	35.9	5.5
977	2.222	92.3	7.5

Inter assay reproducibility

Coefficient of variation (CV) was calculated by measuring 2 samples with various specific antibody levels in 4 ELISA performances during 4 days, in 8-replicate determinations.

Serum No.	OD	Concentration, BAU/ml	CV, %
783	0.643	34.8	5.3
977	2.143	93.4	6.9

11.3 Analytical sensitivity

The limit of determination (LOD), the lowest concentration of the analyte in the sample, which is detected with the declared probability for the test kit Vitrotest SARS-CoV-2 IgG QuantiSpike is 3.5 BAU/ml.

11.4. Linearity range

The linearity range of the test kit Vitrotest SARS-CoV-2 IgG QuantiSpike is within 10-164 BAU/ml.

11.5. Compliance of the test-kit calibrators with the International Standard

Vitrotest SARS-CoV-2 IgG QuantiSpike calibrators comply with the First WHO International Standard for anti-SARS-CoV-2 immunoglobulin (human) code: 20/136 (NIBSC, UK). The coefficient of determination (R²) is 0.99.

12. LIMITATIONS OF THE PROCEDURE

If the test sample was obtained in the first days after the infection IgG antibodies may not be detected. Therefore a negative result does not exclude the SARS-CoV-2 infection. In the presence of clinical manifestations of the disease it is recommended to repeat the test in 1-2 weeks.

Negative test results in immunosuppressed individuals should also be interpreted with caution.

The diagnosis should take into consideration clinical history, symptomatology, as well as the results of other laboratory tests.

13. TROUBLESHOOTING

<i>Possible causes</i>	<i>Solutions</i>
<i>High background in all wells</i>	
Contaminated washer	Clean the washer head, then rinse it with 30 % ethanol and distilled water
Low quality water or contaminated water	Use distilled/DI with resistivity $\geq 10 \text{ M}\Omega\text{-cm}$.
Using contaminated glassware	Use clean glassware
Using chlorine based disinfectants	Use disinfectants without chlorine
Using contaminated tips	Use new tips
Increased time of incubation or temperature regimen was changed	Follow incubation regimen according to instruction for use
<i>High background in a few wells</i>	
TMB solution was added more than once	Add TMB solution once
Pipette shaft was contaminated with conjugate solution	Clean the pipette; pipette the liquids carefully
One the channels of the washer was contaminated	Clean the washer channel, clean the washer
<i>OD of the positive control below normal</i>	
Conjugate solution/tmb solution was prepared improperly or not added	Run ELISA repeatedly, prepared conjugate solution / TMB solution properly
Reduced incubation time in one of the stages	Follow incubation regimen according to the instruction for use
<i>Visual colour intensity of the wells does not correspond to optical density</i>	
The optical beam or another component of the reader is misaligned or malfunctioning	Test the absorbance reader's performance

REFERENCE

1. Juanjuan Zhao Jr., Quan Yuan et. al. Antibody responses to SARS-CoV-2 in patients of novel coronavirus disease 2019.// *Clinical Infectious Diseases.*, - 2020 Mar. 20 doi: 10.1093/cid/ciaa344.
2. Laboratory testing for coronavirus disease (COVID-19) in suspected human cases // WHO. Interim guidance 19 March 2020. WHO reference number: WHO/COVID-19/laboratory/2020.5.
3. Marco Cascella ; Michael Rajnik et.al. Features, Evaluation and Treatment Coronavirus (COVID-19). // NCBI Bookshelf. StatPearls Publishing; 2020 – P16.
4. Patrick C. Y. Woo, Susanna K. P. Lau. et.al. Differential Sensitivities of Severe Acute Respiratory Syn- drome (SARS) Coronavirus Spike Polypeptide Enzyme-Linked Immunosorbent Assay (ELISA) and SARS Coronavirus Nucleocapsid Protein ELISA for Serodiagnosis of SARS Coronavirus Pneumonia. // *J. Clin. Microbiol.*, - 2005 – V. 43 N.7 - p. 3054–3058.
5. Quan-Xin Long, Bai-Zhong Liu et. al. Antibody responses to SARS-CoV-2 in patients with COVID-19.// *Nature Medicine.*, - 2020 April 29 doi: 10.1038/s41591-020-0897-1.
6. Quan-xin Long, Hai-jun Deng et.al. Antibody responses to SARS-CoV-2 in COVID-19 patients: the perspective application of serological tests in clinical practice. // medRxiv preprint doi: <https://doi.org/10.1101/2020.03.18.20038018>.
7. Shu-Yuan Xiao, Yingjie Wu, Huan Liu. Evolving status of the 2019 novel corona virus infection: Proposal of conventional serologic assays for disease diagnosis and infection monitoring. // *J Med Virol.*, - 2020 - 92(5) – p.464-467. doi: 10.1002/jmv.25702. Epub 2020 Feb 17.
8. Wu, L.-P., Wang, N.-C. et.al. (2007). Duration of Antibody Responses after Severe Acute Respiratory Syndrome. // *Emerging Infectious Diseases.*, - 2007 - 13(10) - p.1562-1564.
9. Bao Y., Ling Y., Chen Y. et. al. Dynamic anti-spike protein antibody profiles in COVID-19 patients.// *International Journal of Infectious Diseases.*- 103 (2021) - p.540–548.
10. Brochet E., Demey B. et. al. Anti-spike, Anti-nucleocapsid and Neutralizing Antibodies in SARS-CoV-2 Inpatients and Asymptomatic Individuals.// *Front. Microbiol.*, 19 October 2020/<https://doi.org/10.3389/fmicb.2020.584251>

Catalogue number

Consult instructions for use

In vitro diagnostic medical device

Manufacturer

Caution

Contains sufficient for <n> tests

Temperature limit

Batch code

Use-by date

Date of manufacture

Keep away from sunlight

Signifies European conformity (CE) mark

Inst_SARS-CoV-2-IgG_QuantiSpike_IgG_EL040-96_V01_ENG
Edition 1st, 16.12.2021

Vitrotest Europe Sp. z o.o.
Krakowska str., 139-155, 50-428, Wrocław, Poland
tel.: +48 882 950 379,
e-mail: info@vitrotest.pl, www.vitrotest.pl

Vitrotest SARS-CoV-2 IgG QuantiSpike

ASSAY PROCEDURE

Keep all reagents and specimens for at least 30 min at 18-25 °C before use

Dispense 90 µl of [SAMPLE PREDILUENT] (brown-green colour) into the wells of [PREDILUTION PLATE] and add 10 µl of calibrators, positive control and the samples
(colour changes from brown-green to blue)

Dispense 90 µl of [SAMPLE DILUENT] (yellow colour) into the wells of [ELISA STRIPS] and add 10 µl of prediluted calibrators, positive control and samples into the wells of [ELISA STRIPS] in the following order:
A1 - [CAL 200], B1 - [CAL 100], C1 - [CAL 50], D1 - [CAL 25], E1 - [CAL 0] (or vice versa: A1 - [CAL 0], B1 - [CAL 25], C1 - [CAL 50], D1 - [CAL 100], E1 - [CAL 200]),
F1 - [CONTROL +], G1 and other wells – patient samples
(colour changes from yellow to green)

Cover wells with an adhesive film, incubate for 30 min at 37 °C

Rinse the wells 5 times with diluted 1:20 (1+19) washing solution Tween-20 (300 µl per well)

Add 100 µl of [CONJUGATE SOLUTION] to each well
(violet colour)

Cover wells with an adhesive film, incubate for 30 min at 37 °C

Rinse the wells 5 times with diluted 1:20 (1+19) washing solution Tween-20 (300 µl per well)

Add 100 µl of [TMB SOLUTION] to each well

Incubate the plate for 15 min in the dark at 18-25 °C

Stop the reaction by adding 100 µl of [STOP SOLUTION]
(colour changes from blue to yellow)

Determine the optical density (OD) at 450/620-695 nm

Build a calibration curve, determine the concentration of IgG specific antibodies to SARS-CoV-2 S-protein in the samples and interpret the test results according to the table:

IgG concentration	Interpretation
> 25 BAU/ml	POSITIVE
20-25 BAU/ml	DOUBTFUL
< 20 BAU/ml	NEGATIVE

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«HSV 2 IgG-ИФА»

**A SOLID-PHASE ENZYME IMMUNOASSAY
FOR THE QUALITATIVE DETERMINATION
OF IgG ANTIBODIES TO HERPES SIMPLEX
VIRUSES 2 (HSV 2) IN HUMAN SERUM OR PLASMA**

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Authorized Representative in EU:
Polmed.de
Steinacker 20, D-73773
Aichwald, Germany
e-mail: info@polmed.de

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xemahelp

xemahelp@gmail.com

Instruction for use
A solid-phase enzyme immunoassay kit
for the quantitative determination of
total IgE in human serum or plasma

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
K200

CONTENT

1. INTENDED USE	2
2. GENERAL INFORMATION	2
3. TEST PRINCIPLE	2
4. KIT COMPONENTS	3
5. EQUIPMENT AND MATERIAL REQUIRED BUT NOT PROVIDED	4
6. WARNING AND PRECAUTIONS	4
7. SPECIMEN COLLECTION, TRANSPORTATION AND STORAGE OF SAMPLES	5
8. TRANSPORTATION AND STORAGE TERMS OF KIT, WASTE DISPOSAL	5
9. REAGENTS PREPARATION	6
10. ASSAY PROCEDURE	7
11. TEST VALIDITY	8
12. EXPECTED VALUES	8
13. PERFORMANCE CHARACTERISTICS	9
14. REFERENCES	10
SAMPLES IDENTIFICATION PLAN	11

Instruction for use
A solid-phase enzyme immunoassay kit
for the quantitative determination of
total IgE in human serum or plasma
Total IgE EIA

1. INTENDED USE

The Total IgE EIA kit is an enzyme immunoassay, intended for the quantitative determination of total IgE concentration in human serum or plasma.

The field of application is clinical laboratory diagnostics.

2. GENERAL INFORMATION

Total immunoglobulin E (IgE) serum level is widely reported as the laboratory marker of atopic diseases such as atopic asthma, atopic dermatitis, and pollenosis. An atopic (IgE-dependent) mechanism can also underlie gastroenterocolitis, urticaria, other forms of vasculitis (including systemic), cholecystitis, vulvovaginitis, and cystitis. Part of the drug allergy (mainly to penicillin and protein drugs) also develops according to the IgE-dependent mechanism. In all of the conditions listed above, the production of high titers of specific IgE antibodies can lead to an increase in the level of total IgE in the serum. A particularly high level of total IgE is characteristic of atopic dermatitis. In addition to atopic diseases, total serum IgE is significantly increased in parasitic infestations and mycoses (especially systemic), rarely in systemic autoimmune diseases and immunodeficiency states (especially in hyper-IgE syndrome), as well as in mastocytosis (mast cell tumor) and extremely rare IgE-myeloma. A decrease in the level of total IgE in serum (below 15 IU/ml in adults) is a rare and little-studied phenomenon described in hypogammaglobulinemia, some autoimmune diseases, ulcerative colitis, and primary biliary cirrhosis.

3. TEST PRINCIPLE

The determination of the total IgE is based on the two-site sandwich enzyme immunoassay principle. On the inner surface of the microplate wells are immobilized specific murine monoclonal antibodies to human IgE. Second antibodies – rabbit polyclonal antibodies to IgE conjugated to the horseradish peroxidase is used as enzyme conjugate. The analysis procedure includes three stages of incubation:

- during the first stage the total IgE from the specimen is captured by the monoclonal antibodies coated onto the microwell surface;
- during the second stage horseradish peroxidase-conjugated with rabbit polyclonal antibodies bind to free epitopes of immobilized total IgE, fixed in the formed at the previous stage complexes;
- during the third stage, the complexes formed due to the reaction with the chromogen 3,3',5,5'-tetramethylbenzidine are visualized.

After stopping the reaction with a stop solution, the intensity of the color of the microwells is measured. The optical density in the microwell is directly related to the quantity of the measured total IgE in the serum specimen (plasma).

The concentration is determined according to the calibration graph of the dependence of the optical density on the content of total IgE in the calibration samples.

4. KIT COMPONENTS

Code of component	Symbol	Name	Volume	Qty, pcs.	Description
P200Z	SORB MTP	Microplate	-	1	96-well polystyrene strip microplate coated with murine monoclonal antibodies to total IgE; ready to use
C200Z	CAL 1	Calibrator C1	0.8 mL	1	Solution based on phosphate buffer, free of total IgE, with preservative, ready to use (yellow liquid)
C200Z	CAL 1-5	Calibrators	0.8 mL	5	Solutions based on phosphate buffer, containing 50; 200; 500 and 1000 IU/mL of total IgE, ready to use (red liquids)
Q200Z	CONTROL	Control serum	0.8 mL	1	Solution based on human serum, containing of known total IgE content, with preservative, ready to use (colourless liquid)
T200Z	CONJ HRP	Conjugate Solution	14 mL	1	Solution of rabbit polyclonal antibodies to human total IgE conjugated to the horseradish peroxidase; ready to use (red liquid)
S011Z	DIL	EIA Buffer	14 mL	1	Buffer solution with detergent and preservative, ready to use (blue liquid)
R055Z	SUBS TMB	Substrate Solution	14 mL	1	Tetramethylbenzidine (TMB) substrate solution; ready to use (colourless liquid)
S008Z	BUF WASH 26X	26x Concentrate Washing Solution	30 mL	1	Buffer solution with detergent, 26x concentrate (colourless liquid)
R050Z	STOP	Stop Solution	14 mL	1	5.0% solution of sulphuric acid; ready to use (colourless liquid)
The kit also includes instruction for use, quality control data sheet and plate sealing tape (3 pcs.)					

5. EQUIPMENT AND MATERIAL REQUIRED BUT NOT PROVIDED

- microplate photometer with 450 nm wavelength or 450\620-680 nm;
- dry thermostat for 37 °C±2 °C;
- automatic plate washer (optional);
- micropipettes with variable volume, range volume 5-1000 µL;
- graduated cylinder of 1000 mL capacity;
- distilled or deionized water;
- timer;
- vortex mixer;
- disposable gloves;
- absorbent paper.

6. WARNING AND PRECAUTIONS

In order to prevent incorrect results, strictly follow the recommended order and duration of the analysis procedure.

6.1. The kit is for *in vitro* diagnostic use only. For professional laboratory use.

6.2. Follow the rules mentioned below during the kit using:

- do not use kit beyond expire date;
- do not use the kit if its packaging is damaged;
- in order to avoid contamination, use new tips to pipette samples and reagents;
- use only verified equipment;
- close each vial with its own cap, after using the reagent;
- do not use components of other kits or reagents of other manufacturers;
- do not let wells dry after completing the rinsing step; immediately proceed to the next stage;
- avoid bubbles when adding reagents.

ATTENTION! The TMB substrate solution is light sensitive. Avoid prolonged exposure of the component to light.

6.3. Some kit components, such as stop solution, substrate solution, and washing solution, may cause toxic or irritant effects. If they get on the skin or mucosa, the affected area should be washed with plenty of running water.

6.4. All human products, including patient samples, should be considered potentially infectious. Handling and disposal should be in accordance with the procedures defined by an appropriate national biohazard safety guidelines or regulations.

6.5. The Calibrators and Control Serum included in the kit are negative for antibodies to HIV 1,2, hepatitis C virus and HBsAg, but the reagents should be considered as potentially infectious material and handled carefully.

6.6. Specimens must not contain any azide compounds, as they inhibit activity of peroxidase.

6.7. Wear protective gloves, protective clothing, eye protection, face protection.

6.8. Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.

6.9. Safety Data Sheet for this product is available upon request directly from XEMA LLC.

6.10. Serious incidents related to the kit must be reported to the manufacturer, Authorized Representative, and to the Competent Authority of the EU member state(s) where the incident has occurred.

7. SPECIMEN COLLECTION, TRANSPORTATION AND STORAGE OF SAMPLES

7.1. Blood sampling should be carried out from the cubital vein with a disposable needle using a vacuum blood sampling system. Serum or plasma specimens should be clearly labeled and identified. Serum must be separated from the clot as early as possible to avoid hemolysis of red blood cells. If there are any visible particles in the sample, they should be removed by centrifugation at 3000-5000 rpm for 20 minutes at room temperature or by filtration.

Don't use samples with high lipidemia, hemolysis as they may give false test results.

7.2. Specimen should be stored at +2...+8°C up to 3 days. Specimen held for a longer time, should be placed in a freezer at -15°C or below; do not refreeze/thaw samples.

7.3. For the transportation of samples, it is recommended to use triple packaging. The primary package is the labeled tube containing the sample. Secondary packaging is a polyethylene bag that is hermetically closed with a zip-lock. The outer packaging is a heat-insulating container, while the secondary packaging is placed in the outer packaging for transportation in the center of the thermal container. Frozen refrigerants are placed on the bottom, along the side walls of the thermal container, and cover the samples with them.

8. TRANSPORTATION AND STORAGE TERMS OF KIT, WASTE DISPOSAL

Information about the singularity storage conditions, transportation of the kit, and disposal of waste should be taken into account by all persons who participate in these processes.

8.1. Transportation

The Total IgE EIA kit should be transported in the manufacturer's packaging at +2...+8°C. Single transportation at the temperature up to 25°C for 5 days is acceptable.

8.2. Storage

The Total IgE EIA kit should be stored in the manufacturer's packaging at +2...+8°C. Do not freeze.

The kit contains reagents sufficient for 96 determinations including Calibrators and Control Serum.

Once opened test-kit is stable for 2 months when stored properly as intended by manufacturer at 2-8°C.

In case of partial use of the kit, the components should be stored in the following way:

- strips that remain unused must be carefully sealed with the plate sealing tape and stored at +2...+8°C within 2 months;
- EIA Buffer, Substrate Solution, Stop Solution, and Washing Solution concentrate after opening the vial, can be stored tightly closed at +2...+8°C until the kit's shelf life;
- Conjugate Solution, Calibrators and Control Serum after opening the vial, can be stored tightly closed at +2...+8°C within 2 months;
- diluted Washing Solution can be stored at room temperature (+18...+25°C) for up to 5 days or at +2...+8°C for up to 14 days.

Kits that were stored in violation of the storage condition cannot be used.

8.3. Disposal

Expired kit components, used reagents and materials, as well as residual samples must be inactivated and disposed of in accordance with legal requirements.

9. REAGENTS PREPARATION

9.1. All reagents (including microstrips) and test samples should be allowed to reach room temperature (+18...+25 °C) for at least 30 minutes before use.

9.2. Microplate preparation

Open the package with the microplate and install the required number of strips into the frame. Unused strips must be sealed with plate sealing tape to prevent moisture from affecting the plate's holes and placed back in the bag.

9.3. Washing solution preparation

Add the contents of the 30 mL washing solution concentrate vial to 750 mL of distilled or deionized water and mix thoroughly. In case of partial use of the kit, take the necessary amount of washing solution concentrate and dilute it 26 times with distilled or deionized water.

The spending of the components in case of partial use of the kit is given in the table:

Quantity of strips	1	2	3	4	5	6	7	8	9	10	11	12
Volume of the washing solution concentrate, mL	2.5	5	7.5	10	12.5	15	17.5	20	22.5	25	27.5	30
Volume of water, mL	62.5	125	187.5	250	312.5	375	437.5	500	562.5	625	687.5	750

10. ASSAY PROCEDURE

- 10.1 Put the desired number of strips into the frame based on the number of test samples in 2 replicates and 12 wells for Calibrators and Control Serum (2 wells for each Calibrator (CAL 1-5) and 2 wells for Control Serum (Q)).
- 10.2 Dispense **50 µL of EIA Buffer** to all wells.
- 10.3 Dispense **50 µL of Calibrators and Control Serum as well as 50 µL of test serum/plasma samples (SAMP)** to the wells of the microplate according to the scheme below. The introduction of Calibrators, Control Serum and test samples should be carried out within 5 minutes to ensure equal incubation time for the first and last samples.

NOTE: during performing several independent series of tests, Calibrators, and Control Serum should be used each time.

Scheme of introduction of samples

	1	2	3	4	5	6	7	8	9	10	11	12
A	CAL1	CAL1	SAMP3	SAMP3	SAMP11	SAMP11						
B	CAL2	CAL2	SAMP4	SAMP4	SAMP12	SAMP12						
C	CAL3	CAL3	SAMP5	SAMP5	SAMP13	SAMP13						
D	CAL4	CAL4	SAMP6	SAMP6	SAMP14	SAMP14						
E	CAL5	CAL5	SAMP7	SAMP7	SAMP15	SAMP15						
F	Q	Q	SAMP8	SAMP8								
G	SAMP1	SAMP1	SAMP9	SAMP9								
H	SAMP2	SAMP2	SAMP10	SAMP10								

- 10.4 Carefully mix the contents of the microplate in a circular motion on a horizontal surface, cover strips with a plate sealing tape and incubate for **30 minutes at +37 °C**.
- 10.5 At the end of the incubation period, remove and discard the plate cover. Aspirate and wash each well 3 times using an automatic washer or an 8-channel dispenser. For each washing, add 300 µL of Washing Solution (see 9.3) to all wells, then remove the liquid by aspiration or decantation. The residual volume of the Washing Solution after each aspiration or decantation should be no more than 5µL. After washing, carefully remove the remaining liquid from the wells on the absorbent paper. For the automatic washer/analyzer, the Washing Solution volume can be increased to 350 µL.
- 10.6 Add **100 µL of Conjugate Solution** to all wells.
- 10.7 Cover strips with a plate sealing tape and incubate for **30 minutes at +37 °C**.
- 10.8 At the end of the incubation period, aspirate and wash each well 5 times as described in 10.5.
- 10.9 Add **100 µL of Substrate Solution** to all wells. The introduction of the Substrate Solution into the wells must be carried out within 2-3 minutes. Incubate the microplate in the dark **at room temperature (+18...+25°C) for 20 minutes**.
- 10.10 Add **100 µL of Stop Solution** to all wells in the same order as the Substrate Solution. After adding the Stop Solution, the contents of the wells turn yellow.
- 10.11 Read the optical density (OD) of the wells at 450nm and reference light filters 620–680 nm using a microplate photometer within 5 minutes of adding the stop solution. Set photometer blank on CAL1.
- 10.12 Plot a calibration curve in linear coordinates: (x) is the concentration of total IgE in the Calibrators IU/mL, (y) – OD versus concentration of total IgE (OD 450 nm / 620–680 nm). Manual or computerized data reduction is applicable at this stage. Point-by-point or linear data reduction is recommended due to non-linear shape of curve.
- 10.13 Determine the corresponding concentration of total IgE in tested samples from the calibration curve.

11. TEST VALIDITY

The test run shall be considered valid if the OD of CAL1 is above 0.15, and the values of the Control Serum fall into the required range (see Quality control Data Sheet).

12. EXPECTED VALUES

12.1. Therapeutical consequences should not be based on the results of IVD methods alone – all available clinical and laboratory findings should be used by a physician to elaborate therapeutically measures. Each laboratory should establish its own normal range for total IgE. Based on data obtained by XEMA LLC, the following normal range is recommended (see below).

NOTE: values of total IgE concentrations in the tested samples that are below the LoD (3 IU/mL) and also exceed the value of the upper calibrator (1000 IU/mL) should be provided in the following form : «the total IgE concentration of tested sample X is «lower than 3 IU/mL» or «higher than 1000 IU/mL».

12.2. The calibrators concentration values of the Total IgE EIA kit are expressed in IU/mL. To calculate concentrations in ng/mL, the received concentration value in IU/mL shall be multiplied by 2.4.

$$1 \text{ IU/mL} = 2.4 \text{ ng/mL.}$$

Sex, age	Units, IU/mL		Units alternative, ng/mL	
	Lower limit	Upper limit	Lower limit	Upper limit
< 6 months	-	12	-	28.8
6-12 months	-	30	-	72.0
1-3 yrs	-	45	-	108.0
4-6 yrs	-	70	-	168.0
7-9 yrs	-	90	-	216.0
10-15 yrs	-	120	-	288.0
>15 yrs	-	130	-	312.0

13. PERFORMANCE CHARACTERISTICS

13.1. Analytical performance characteristics

13.1.1 Precision of Measurement

Repeatability (Intra assay repeatability) was determined by evaluation the coefficient of variation (CV) for 2 different samples during 1 day in 24 replicates on one series of ELISA kit.

Sample	Concentration, IU/mL	CV, %
1	10.6	4.33
2	116.2	5.47

Reproducibility (Inter assay reproducibility) was determined by evaluating the coefficients of variation for 2 samples during 5 days in 8-replicate determinations.

Sample	Concentration, IU/mL	CV, %
1	12.5	8.36
2	113.4	1.47

Reproducibility between lots was investigated by testing samples for one day on three lots. Each sample was run in 8 replicates.

Sample	Concentration1, IU/mL	Concentration2, IU/mL	Concentration3, IU/mL	CV, %
1	12.7	13.3	12.3	3.66
2	115.5	117.8	115.1	1.25

13.1.2 Trueness

The trueness of measurement is the degree of closeness of the average value obtained from a large number of measurement results to the true value. The bias of the measurement result (bias of measurements) is the difference between the mathematical expectation of the measurement result and the true value of the measurand. The bias was calculated for each sample and it was determined that it corresponds to the specified limits of ± 10%.

13.1.3 Linearity

Linearity was determined using sera samples with known total IgE concentration (low and high) and mixing them with each other and buffer solution in different proportions. According to the measurements, linear range of kit is 50-1000 IU/mL \pm 10%.

13.1.4 Analytical sensitivity

Limit of detection (LoD) – the lowest total IgE concentration in the serum or plasma sample that is detected by the Total IgE EIA kit is no lower than 3 IU/mL.

Limit of quantification (LoQ) – the lowest concentration of the analyte in the sample that is determined quantitatively with the declared trueness for Total IgE EIA kit is 50IU/mL.

13.1.5 Analytical specificity

For the analysis result is not affected by the presence in the sample of bilirubin in a concentration of up to 0.21 mg/mL and hemoglobin in a concentration of up to 10 mg/mL.

The cross-reactivity of total IgE with other analytes is shown in the table:

Analyte	Concentration, IU/mL	Cross-reactivity, %
IgA	1000	Not detected
IgM	1000	Not detected
IgG	1000	Not detected

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LOT _____ DATE _____

	Manufacturer
	<i>In vitro</i> diagnostic medical device
	Catalogue number
 YYYY-MM	Use-by date
	Batch code
	Temperature limit
	Contains sufficient for <n> tests
	Caution
	Consult instructions for use
	Conformity Marking with technical regulations in Ukraine
	Authorized representative in the European Community/European Union
	CE Conformity Marking

**For any issues related to operation of the kit and technical support,
please contact by telefon number**

+38 044 294-69-78

or write to:

qa@xema.com.ua

XEMA LLC
Akademika Yefremova St. 23
03179, Kyiv, Ukraine
tel.:+38 044 422-62-16
tel.:+38 044 294-69-78
E-mail: qa@xema.com.ua
www.xema.in.ua

Instruction for use
A solid-phase enzyme immunoassay kit
for the quantitative determination of
triiodothyronine in human serum or plasma

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ASSAY PROCEDURE

Preparation of working conjugate solution

Conjugate Dilution Buffer + Conjugate Concentrate = Working conjugate solution

1 : 1

450 μ L 450 μ L

Dispensing of Calibrators, Control Serum and test samples

test samples 25 μ L

CAL1 CAL2 CONTROL

SORB MTP

Working Conjugate Solution

100 μ L

Incubation 1

+37 $^{\circ}$ C

60'

continuous shaking

OD measuring, calculation of results

450 nm

Stop Solution

100 μ L

Incubation 2

+18-25 $^{\circ}$ C

15'

Substrate Solution

100 μ L

SUBS TMB

Washing 5 times

CONTENT

1. INTENDED USE	2
2. GENERAL INFORMATION	2
3. TEST PRINCIPLE	2
4. KIT COMPONENTS	3
5. EQUIPMENT AND MATERIAL REQUIRED BUT NOT PROVIDED	4
6. WARNING AND PRECAUTIONS	4
7. SPECIMEN COLLECTION, TRANSPORTATION AND STORAGE OF SAMPLES	5
8. TRANSPORTATION AND STORAGE TERMS OF KIT, WASTE DISPOSAL	5
9. REAGENTS PREPARATION	6
10. ASSAY PROCEDURE	6
11. TEST VALIDITY	7
12. EXPECTED VALUES	8
13. PERFORMANCE CHARACTERISTICS	8
14. REFERENCES	9
SAMPLES IDENTIFICATION PLAN	10

Instruction for use
A solid-phase enzyme immunoassay kit
for the quantitative determination of
triiodothyronine in human serum or plasma
T3 EIA

1. INTENDED USE

The T3 EIA kit is an enzyme immunoassay, intended for the quantitative determination of triiodothyronine in human serum or plasma.

The field of application is clinical laboratory diagnostics.

2. GENERAL INFORMATION

Triiodothyronine (T3) is a hormone with a molecular weight of 651 Da, 58% of which is iodine. Thyroid hormones thyroxin (T4) and 3,5,3'-triiodothyronine (T3) exert regulatory influences on growth, differentiation, cellular metabolism and development of skeletal and organ systems. T4 and T3 in blood are found both in free and bound form – mostly, they are bound to thyroxin binding globulin (TBG). Only free forms of T3 and T4 exert hormonal activity also their percentage is very low – 0.3% for T3 and 0.03% for T4.

The concentration of T3 is much less than that of T4 but its metabolic activity is about 3 times greater. About 80% of T3 is produced in peripheral tissues by deiodination of T4, and only 20% is secreted by thyroid gland. That is why in hypothyroid patients T3 level may for a long time remain on the lower limit of the normal range, because its loss may be compensated by enhanced conversion of T4 into T3.

Determination of T3 level is most useful in T3-hyperthyroidism because 5-10% of such patients do not show significant changes in T4 level while concentration of T3 is highly elevated. Elevated T3 levels are seen in early thyroid hypofunction, after intake of estrogens, oral contraceptives, heroin, methadone, during pregnancy.

Decreased concentrations of T3 are found in initial stage of hyperthyroidism, acute and subacute thyroiditis, after intake of androgens, dexamethasone, salicylates. Decreased concentrations of T3 are found in initial stage of hyperthyroidism, acute and subacute thyroiditis, after intake of androgens, dexamethasone, salicylates.

3. TEST PRINCIPLE

The determination of triiodothyronine is based on the competition principle of the enzyme immunoassay. On the inner surface of the microplate wells are immobilized specific rabbit polyclonal to T3 antibodies. T3 conjugated to the horseradish peroxidase is used as enzyme conjugate. The analysis procedure includes two stages of incubation:

- during the first stage T3 from the specimen competes with the conjugated T3 for coating antibodies. As a result, a complex bounded to the solid phase and containing peroxidase is formed.

- during the second stage, the complexes formed due to the reaction with the chromogen 3,3',5,5'-tetramethylbenzidine are visualized.

After stopping the reaction with a stop solution, the intensity of the color of the microwells is measured. The optical density in the microwell is inversely related to the quantity of the measured T3 in the serum specimen (plasma).

The concentration is determined according to the calibration graph of the dependence of the optical density on the content of T3 in the calibration samples.

4. KIT COMPONENTS

Code of component	Symbol	Name	Volume	Qty, pcs.	Description
P211Z	SORB MTP	Microplate	-	1	96-well polystyrene strip microplate coated with rabbit polyclonal antibodies to T3, ready to use;
C211Z	CAL 1	Calibrator C1	0.5 mL	1	Solution based on tris buffer (pH 7.2-7.4), free of T3, with preservative, ready to use (yellow liquid)
C211Z	CAL 2-5	Calibrators	0.5 mL	4	Solutions based on tris buffer (pH 7.2-7.4), containing 0,75; 1,5; 7,5 and 15 nmol/L of T3, with preservative, ready to use (blue liquids)
Q211Z	CONTROL	Control serum	0.5 mL	1	Solution based on human plasma, containing of known T3 content, with preservative, ready to use (colourless liquid)
T211XZ	CONJ 2X	Conjugate Concentrate	7 mL	1	Solution of T3 conjugated to the horseradish peroxidase; 2x concentrate (purple liquid)
ST211Z	DIL CONJ	Conjugate Dilution Buffer	7 mL	1	Buffer solution with detergent ready to use (red liquid)
R055Z	SUBS TMB	Substrate Solution	14 ml (MЛ)	1	Tetramethylbenzidine (TMB) substrate solution; ready to use (colourless liquid)
S008Z	BUF WASH 26X	26x Concentrate Washing Solution	22 ml (MЛ)	1	Buffer solution with detergent, 26x concentrate (colourless liquid)
R050Z	STOP	Stop Solution	14 ml (MЛ)	1	5,0% solution of sulphuric acid; ready to use (colourless liquid)
The kit also includes instruction for use, quality control data sheet and plate sealing tape (2 pcs.)					

5. EQUIPMENT AND MATERIAL REQUIRED BUT NOT PROVIDED

- microplate photometer with 450 nm wavelength;
- shaker maintaining a speed of 500 rpm for +37 °C±2°C;
- automatic plate washer (optional);
- micropipettes with variable volume, range volume 5-1000 µL;
- graduated cylinder of 1000 mL capacity;
- distilled or deionized water;
- timer;
- vortex mixer;
- disposable gloves;
- absorbent paper.

6. WARNING AND PRECAUTIONS

In order to prevent incorrect results, strictly follow the recommended order and duration of the analysis procedure.

6.1. The kit is for *in vitro* diagnostic use only. For professional laboratory use.

6.2. Follow the rules mentioned below during the kit using:

- do not use kit beyond expire date;
- do not use the kit if its packaging is damaged;
- in order to avoid contamination, use new tips to pipette samples and reagents;
- use only verified equipment;
- close each vial with its own cap, after using the reagent;
- do not use components of other kits or reagents of other manufacturers;
- do not let wells dry after completing the rinsing step; immediately proceed to the next stage;
- avoid bubbles when adding reagents.

ATTENTION! The TMB substrate solution is light sensitive. Avoid prolonged exposure of the component to light.

6.3. Some kit components, such as stop solution, substrate solution, and washing solution, may cause toxic or irritant effects. If they get on the skin or mucosa, the affected area should be washed with plenty of running water.

6.4. All human products, including patient samples, should be considered potentially infectious. Handling and disposal should be in accordance with the procedures defined by an appropriate national biohazard safety guidelines or regulations.

6.5. The Calibrators and Control Serum included in the kit are negative for antibodies to HIV 1,2, hepatitis C virus and HBsAg, but the reagents should be considered as potentially infectious material and handled carefully.

6.6. Specimens must not contain any azide compounds, as they inhibit activity of peroxidase.

6.7. Wear protective gloves, protective clothing, eye protection, face protection.

6.8. Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.

6.9. Safety Data Sheet for this product is available upon request directly from XEMA LLC.

6.10. Serious incidents related to the kit must be reported to the manufacturer, Authorized Representative, and to the Competent Authority of the EU member state(s) where the incident has occurred.

7. SPECIMEN COLLECTION, TRANSPORTATION AND STORAGE OF SAMPLES

7.1. Blood sampling should be carried out from the cubital vein with a disposable needle using a vacuum blood sampling system. Serum or plasma specimens should be clearly labeled and identified. Serum must be separated from the clot as early as possible to avoid hemolysis of red blood cells. If there are any visible particles in the sample, they should be removed by centrifugation at 3000-5000 rpm for 20 minutes at room temperature or by filtration.

Don't use samples with high lipidemia, hemolysis as they may give false test results.

7.2. Specimen should be stored at +2...+8°C up to 3 days. Specimen held for a longer time, should be placed in a freezer at -15°C or below; do not refreeze/thaw samples.

7.3. For the transportation of samples, it is recommended to use triple packaging. The primary package is the labeled tube containing the sample. Secondary packaging is a polyethylene bag that is hermetically closed with a zip-lock. The outer packaging is a heat-insulating container, while the secondary packaging is placed in the outer packaging for transportation in the center of the thermal container. Frozen refrigerants are placed on the bottom, along the side walls of the thermal container, and cover the samples with them.

8. TRANSPORTATION AND STORAGE TERMS OF KIT, WASTE DISPOSAL

Information about the singularity storage conditions, transportation of the kit, and disposal of waste should be taken into account by all persons who participate in these processes.

8.1. Transportation

The T3 EIA kit should be transported in the manufacturer's packaging at +2...+8°C. Single transportation at the temperature up to 25°C for 5 days is acceptable.

8.2. Storage

The T3 EIA kit should be stored in the manufacturer's packaging at +2...+8°C. Do not freeze.

The kit contains reagents sufficient for 96 determinations including Calibrators and Control Serum.

Once opened test-kit is stable for 2 months when stored properly as intended by manufacturer at 2-8°C.

In case of partial use of the kit, the components should be stored in the following way:

- strips that remain unused must be carefully sealed with the plate sealing tape and stored at +2...+8°C within 2 months;
- Substrate Solution, Stop Solution, and Washing Solution concentrate after opening the vial, can be stored tightly closed at +2...+8°C until the kit's shelf life;
- Conjugate Concentrate, Conjugate Dilution Buffer, Calibrators and Control Serum after opening the vial, can be stored tightly closed at +2...+8°C within 2 months;
- diluted Washing Solution can be stored at room temperature (+18...+25°C) for up to 5 days or at +2...+8°C for up to 14 days.

Kits that were stored in violation of the storage condition cannot be used.

8.3. Disposal

Expired kit components, used reagents and materials, as well as residual samples must be inactivated and disposed of in accordance with legal requirements.

9. REAGENTS PREPARATION

9.1. All reagents (including microstrips) and test samples should be allowed to reach room temperature (+18...+25 °C) for at least 30 minutes before use.

9.2. Microplate preparation

Open the package with the microplate and install the required number of strips into the frame. Unused strips must be sealed with plate sealing tape to prevent moisture from affecting the plate's holes and placed back in the bag.

9.3. Washing Solution preparation

Add the contents of the 22 mL Washing Solution concentrate vial to 550 mL of distilled or deionized water and mix thoroughly. In case of partial use of the kit, take the necessary amount of Washing Solution concentrate and dilute it 26 times with distilled or deionized water.

9.4. Working conjugate solution preparation

Prepare a working conjugate solution by 2 dilutions of Conjugate Concentrate in Conjugate Dilution Buffer (eg, 450 µL of concentrate + 450 µL of Conjugate Dilution Buffer). In the case of partial use of the kit, take the necessary amount of Conjugate Concentrate and dilute it 2 times with Conjugate Dilution Buffer, since the working conjugate solution in a diluted form is not stored for a long time.

The spending of the components in case of partial use of the kit is given in the table:

Quantity of strips	1	2	3	4	5	6	7	8	9	10	11	12
Volume of the Washing Solution concentrate, mL	1.8	3.6	5.4	7.2	9	10.8	12.6	14.4	16.2	18	19.8	22
Volume of water, mL	45	90	135	180	225	270	315	360	405	450	495	550
Volume of Conjugate Concentrate, mL	0.45	0.9	1.35	1.8	2.25	2.7	3.15	3.6	4.05	4.5	4.95	5.4
Volume of Conjugate Dilution Buffer, mL	0.45	0.9	1.35	1.8	2.25	2.7	3.15	3.6	4.05	4.5	4.95	5.4

10. ASSAY PROCEDURE

- Put the desired number of strips into the frame based on the number of test samples in 2 replicates and 12 wells for Calibrators and Control Serum (2 wells for each calibrator (CAL 1-5) and 2 wells for control serum (Q)).
- Prepare Working conjugate solution as described in 9.4.
- Dispense **25 µL of Calibrators and Control Serum as well as 25 µL of test serum/plasma samples (SAMP)** to the wells of the microplate according to the scheme below. The introduction of Calibrators, Control Serum and test samples should be carried out within 5 minutes to ensure equal incubation time for the first and last samples.

Note: during performing several independent series of tests, Calibrators, and Control Sample should be used each time.

Scheme of introduction of samples

	1	2	3	4	5	6	7	8	9	10	11	12
A	CAL1	CAL1	SAMP3	SAMP3	SAMP11	SAMP11						
B	CAL2	CAL2	SAMP4	SAMP4	SAMP12	SAMP12						
C	CAL3	CAL3	SAMP5	SAMP5								
D	CAL4	CAL4	SAMP6	SAMP6								
E	CAL5	CAL5	SAMP7	SAMP7								
F	Q	Q	SAMP8	SAMP8								
G	SAMP1	SAMP1	SAMP9	SAMP9								
H	SAMP2	SAMP2	SAMP10	SAMP10								

- 10.4 Dispense **100 µL of Working conjugate solution** to all wells.
- 10.5 Carefully mix the contents of the microplate in a circular motion on a horizontal surface, cover strips with a plate sealing tape and incubate for **60 minutes at +37°C with continuous shaking 500 rpm**.
- 10.6 At the end of the incubation period, remove and discard the plate cover. Aspirate and wash each well 5 times using an automatic washer or an 8-channel dispenser. For each washing, add 300 µL of Washing Solution (see 9.3) to all wells, then remove the liquid by aspiration or decantation. The residual volume of the Washing Solution after each aspiration or decantation should be no more than 5µL. After washing, carefully remove the remaining liquid from the wells on the absorbent paper. For the automatic washer/analyzer, the Washing Solution volume can be increased to 350 µL.
- 10.7 Add **100 µL of Substrate Solution** to all wells. The introduction of the substrate solution into the wells must be carried out within 2-3 minutes. Incubate the microplate in the dark **at room temperature (+18...+25°C) for 15 minutes**.
- 10.8 Add **100 µL of Stop Solution** to all wells in the same order as the substrate solution. After adding the Stop Solution, the contents of the wells turn yellow.
- 10.9 Read the optical density (OD) of the wells at 450nm using a microplate photometer within 5 minutes of adding the Stop Solution.
- 10.10 Plot a calibration curve in semi-logarithmic coordinates: (x) is the decimal logarithm of the T3 concentration in the calibrators nmol/L, (y) – OD versus T3 concentration (OD 450 nm). Manual or computerized data reduction is applicable at this stage. Point-by-point or linear data reduction is recommended due to non-linear shape of curve. Adjust the concentration of CAL1 to an infinitesimally small value, for example, 0.001 nmol/L.
- 10.11 Determine the corresponding concentration of T3 in tested samples from the calibration curve.

11. TEST VALIDITY

The test run shall be considered valid if the OD of CAL1 is above 1.2, and the values of the Control Serum fall into the required range (see Quality control Data Sheet).

12. EXPECTED VALUES

Therapeutical consequences should not be based on results of IVD methods alone – all available clinical and laboratory findings should be used by a physician to elaborate therapeutically measures. Each laboratory should establish its own normal range for T3. Based on data obtained by XEMA, the following normal range is recommended (see below). NOTE: the patients that have received murine monoclonal antibodies for radioimaging or immunotherapy develop high titered anti-mouse antibodies (HAMA). The presence of these antibodies may cause false results in the present assay. Sera from HAMA positive patients should be treated with depleting adsorbents before assaying.

NOTE: values of T3 concentrations in the tested samples that are below the LoD (0.2 nmol/L) and also exceed the value of the upper calibrator (15 nmol/L) should be provided in the following form: «the T3 concentration of tested sample X is «lower than 0.2 nmol/L» or «higher than 15 nmol/L».

The concentration values of the T3 EIA kit calibrators are expressed in nmol/L. To convert the concentration in ng/mL it is necessary to multiply by 0.65 the obtained concentration value in nmol/L.

$$1 \text{ nmol/L} = 0.65 \text{ ng/mL}$$

Sex, age	Units, nmol/L		Units alternative, ng/mL	
	Lower limit	Upper limit	Lower limit	Upper limit
Healthy donors	1.2	3.2	0.8	2.1

13. PERFORMANCE CHARACTERISTICS

13.1. Analytical performance characteristics

13.1.1 Precision of Measurement

Repeatability (Intra assay repeatability) was determined by evaluation the coefficient of variation (CV) for 2 different samples during 1 day in 24 replicates on one series of ELISA kit.

Sample	Concentration, nmol/L	CV, %
1	2.32	9.16
2	1.45	9.66

Reproducibility (Inter assay reproducibility) was determined by evaluating the coefficients of variation for 2 samples during 5 days in 8-replicate determinations.

Sample	Concentration, nmol/L	CV, %
1	1.38	9.89
2	1.75	8.41

Reproducibility between lots was investigated by testing samples for one day on three lots. Each sample was run in 8 replicates.

Sample	Concentration1, nmol/L	Concentration2, nmol/L	Concentration3, nmol/L	CV, %
1	2.12	2.02	2.27	13.9
2	1.56	1.44	1.81	15.6

13.1.2 Trueness

The trueness of measurement is the degree of closeness of the average value obtained from a large number of measurement results to the true value. The bias of the measurement result (bias of measurements) is the difference between the mathematical expectation of the measurement result and the true value of the measurand. The bias was calculated for each sample and it was determined that it corresponds to the specified limits of $\pm 10\%$.

13.1.3 Linearity

Linearity was determined using sera samples with known T3 concentration (low and high) and mixing them with each other and buffer solution in different proportions. According to the measurements, linear range of kit is 0.75 –15 nmol/L $\pm 10\%$.

13.1.4 Analytical sensitivity

Limit of detection (LoD) – the lowest T3 concentration in the serum or plasma sample that is detected by the T3 EIA kit is no lower than 0.2 nmol/L.

Limit of quantification (LoQ) – the lowest concentration of the analyte in the sample that is determined quantitatively with the declared trueness for T3 EIA kit is 0.55 nmol/L.

3.1.5 Analytical specificity

For the analysis result is not affected by the presence in the sample of bilirubin in a concentration of up to 0.21 mg/mL and hemoglobin in a concentration of up to 10 mg/mL.

The cross-reactivity of T3 with other analytes is shown in the table:

Analyte	Cross-reactivity, %
L-Thyroxin	0.01
D-Thyroxin	0.04

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LOT _____ DATE _____

	Manufacturer
	<i>In vitro</i> diagnostic medical device
	Catalogue number
 YYYY-MM	Use-by date
	Batch code
	Temperature limit
	Contains sufficient for <n> tests
	Caution
	Consult instructions for use
	Conformity Marking with technical regulations in Ukraine
	Authorized representative in the European Community/European Union
	CE Conformity Marking

**For any issues related to operation of the kit and technical support,
please contact by telefon number**

+38 044 294-69-78

or write to:

qa@xema.com.ua

XEMA LLC
Akademika Yefremova St. 23
03179, Kyiv, Ukraine
tel.:+38 044 422-62-16
tel.:+38 044 294-69-78
E-mail: qa@xema.com.ua
www.xema.com.ua

Instruction for use
A solid-phase enzyme immunoassay kit
for the quantitative determination of
thyroxin in human serum or plasma

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ASSAY PROCEDURE

Dispensing of Calibrators, Control Serum and test samples

test samples
CAL
CONTROL
25 µL
SORB MTP

Conjugate Solution

CONJ HRP
100 µL

Incubation 1

+37 °C
60'

Washing 5 times

Washing 5 times

OD measuring, calculation of results

450 nm

Stop Solution

STOP
100 µL

Incubation 2

+18-25 °C
15'

Substrate Solution

SUBS TMB
100 µL

CONTENT

1. INTENDED USE	2
2. GENERAL INFORMATION	2
3. TEST PRINCIPLE	3
4. KIT COMPONENTS	4
5. EQUIPMENT AND MATERIAL REQUIRED BUT NOT PROVIDED	5
6. WARNING AND PRECAUTIONS	5
7. SPECIMEN COLLECTION, TRANSPORTATION AND STORAGE OF SAMPLES	6
8. TRANSPORTATION AND STORAGE TERMS OF KIT, WASTE DISPOSAL	6
9. REAGENTS PREPARATION	7
10. ASSAY PROCEDURE	7
11. TEST VALIDITY	9
12. EXPECTED VALUES	9
13. PERFORMANCE CHARACTERISTIC	9
14. REFERENCES	10
SAMPLES IDENTIFICATION PLAN	11

Instruction for use
A solid-phase enzyme immunoassay kit
for the quantitative determination of
thyroxine in human serum or plasma
T4 EIA

1. INTENDED USE

The T4 EIA kit is an enzyme immunoassay, intended for the quantitative determination of thyroxine in human serum or plasma.

The field of application is clinical laboratory diagnostics.

2. GENERAL INFORMATION

Thyroxine (T4) and triiodothyronine (T3) are hormones that are produced by the thyroid gland and circulate in the blood both free and bound - mainly with thyroxine-binding globulin (TBG). Only free T3 and T4 are characterized by Hormonal activity, but their share is very small: 0.03% of the total content for T4 and 0.3% - for T3. Concentration of T4 in serum blood is the most accepted indicator of thyroid gland function, which allows you to clearly distinguish between hyper-, hypo- and euthyroidism.

Increase of total T4 concentration is observed with hyperthyroidism, with pituitary tumors, with conditions with elevated TSH levels (pregnancy, acute or chronic active hepatitis, estrogen-secreting tumors or estrogen intake, genetically conditional increase), while taking oral contraceptives, heroin, methadone, thyroid drugs, TSH, thyroliberin.

Decrease of total T4 concentration is observed in hypothyroidism, panhypopituitarism, states of low levels of TSH (acromegaly, nephrotic syndrome, hypoproteinemia, chronic liver disease, androgen-secreting tumors, or androgens, genetically determined decrease), hemolysis, exercise, when taking amino salicylic and acetylsalicylic acids, glucocorticoids, sulfonamides, cholestyramine, reserpine, potassium iodide, triiodothyronine.

3. TEST PRINCIPLE

Determination of the thyroxine is based on competition principle of the enzyme immunoassay. Microwells plate is coated with specific murine monoclonal to thyroxine antibodies. Thyroxine conjugated to the horseradish peroxidase is used as enzyme conjugate. The analysis procedure includes two stages of incubation:

- during the first stage thyroxine from the specimen competes with the conjugated thyroxine for coating antibodies. As a result, a complex bounded to the solid phase and containing peroxidase is formed.
- during the second stage, the complexes formed due the reaction with the chromogen 3,3',5,5'-tetramethylbenzidine are visualized.

After stopping the reaction with a stop solution, the intensity of the color of the microwells is measured. Optical density in the microwell is inversely related to the quantity of the measured thyroxine in the specimen of the serum (plasma).

The concentration is determined according to the calibration graph of the dependence of the optical density on the content of thyroxine in the calibration samples.

4. KIT COMPONENTS

Code of component	Symbol	Name	Volume	Qty, pcs.	Description
P212Z	SORB MTP	Microplate	-	1	96-well polystyrene strip microplate coated with murine monoclonal antibodies to T4; ready to use
C212Z	CAL 1	Calibrator C1	0.5 mL	1	Solution based on human plasma, free of thyroxin, with preservative, ready to use (yellow liquid)
C212Z	CAL 2-5	Calibrators	0.5 mL	4	Solutions based on human plasma, containing 32; 64; 160 and 320 nmol/L of thyroxin, with preservative, ready to use (red liquids)
Q212Z	CONTROL	Control Serum	0.5 mL	1	Solution based on human plasma, containing of known thyroxin content, with preservative, ready to use (colourless liquid)
T212XZ	CONJ	Conjugate Solution	14 mL	1	Solution of thyroxin conjugated to the horseradish peroxidase; ready to use (red liquid)
R055Z	SUBS TMB	Substrate Solution	14 mL	1	Tetramethylbenzidine (TMB) substrate solution; ready to use (colourless liquid)
S008Z	BUF WASH 26X	26x Concentrate Washing Solution	22 mL	1	Buffer solution with detergent, 26x concentrate (colourless liquid)
R050Z	STOP	Stop Solution	14 mL	1	5.0% solution of sulphuric acid; ready to use (colourless liquid)
The kit also includes instruction for use, quality control data sheet and plate sealing tape (2 pcs.)					

5. EQUIPMENT AND MATERIAL REQUIRED BUT NOT PROVIDED

- microplate photometer with 450 nm wavelength;
- dry thermostat for $+37^{\circ}\text{C}\pm 2^{\circ}\text{C}$;
- automatic plate washer (optional);
- micropipettes with variable volume, range volume 5-1000 μL ;
- graduated cylinder of 1000 mL capacity;
- distilled or deionized water;
- timer;
- vortex mixer;
- disposable gloves;
- absorbent paper.

6. WARNING AND PRECAUTIONS

In order to prevent incorrect results, strictly follow the recommended order and duration of the analysis procedure.

6.1. The kit is for *in vitro* diagnostic use only. For professional laboratory use.

6.2. Follow the rules mentioned below during the kit using:

- do not use kit beyond expire date;
- do not use the kit if its packaging is damaged;
- in order to avoid contamination, use new tips to pipette samples and reagents;
- use only verified equipment;
- close each vial with its own cap, after using the reagent;
- do not use components of other kits or reagents of other manufacturers;
- do not let wells dry after completing the rinsing step; immediately proceed to the next stage;
- avoid bubbles when adding reagents.

ATTENTION! The TMB substrate solution is light sensitive. Avoid prolonged exposure of the component to light.

6.3. Some kit components, such as stop solution, substrate solution, and washing solution, may cause toxic or irritant effects. If they get on the skin or mucosa, the affected area should be washed with plenty of running water.

6.4. All human products, including patient samples, should be considered potentially infectious. Handling and disposal should be in accordance with the procedures defined by an appropriate national biohazard safety guidelines or regulations.

6.5. The Calibrators and Control Serum included in the kit are negative for antibodies to HIV 1,2, hepatitis C virus and HBsAg, but the reagents should be considered as potentially infectious material and handled carefully.

6.6. Specimens must not contain any azide compounds, as they inhibit activity of peroxidase.

6.7. Wear protective gloves, protective clothing, eye protection, face protection.

6.8. Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.

6.9. Safety Data Sheet for this product is available upon request directly from XEMA LLC.

6.10. Serious incidents related to the kit must be reported to the manufacturer, Authorized Representative, and to the Competent Authority of the EU member state(s) where the incident has occurred.

7. SPECIMEN COLLECTION, TRANSPORTATION AND STORAGE OF SAMPLES

7.1. Blood sampling should be carried out from the cubital vein with a disposable needle using a vacuum blood sampling system. Serum or plasma specimens should be clearly labeled and identified. Serum must be separated from the clot as early as possible to avoid hemolysis of red blood cells. If there are any visible particles in the sample, they should be removed by centrifugation at 3000-5000 rpm for 20 minutes at room temperature or by filtration.

Don't use samples with high lipidemia, hemolysis as they may give false test results.

7.2. Specimen should be stored at +2...+8°C up to 3 days. Specimen held for a longer time, should be placed in a freezer at -15°C or below, do not refreeze/thaw samples.

7.3. For the transportation of samples, it is recommended to use triple packaging. The primary package is the labeled tube containing the sample. Secondary packaging is a polyethylene bag that is hermetically closed with a zip-lock. The outer packaging is a heat-insulating container, while the secondary packaging is placed in the outer packaging for transportation in the center of the thermal container. Frozen refrigerants are placed on the bottom, along the side walls of the thermal container, and cover the samples with them.

8. TRANSPORTATION AND STORAGE TERMS OF KIT, WASTE DISPOSAL

Information about the singularity storage conditions, transportation of the kit, and disposal of waste should be taken into account by all persons who participate in these processes.

8.1. Transportation

The T4 EIA kit should be transported in the manufacturer's packaging at +2...+8°C. Single transportation at the temperature up to 25°C for 5 days is acceptable.

8.2. Storage

The T4 EIA kit should be stored in the manufacturer's packaging at +2...+8°C. Do not freeze.

The kit contains reagents sufficient for 96 determinations including Calibrators and Control Serum.

Once opened test-kit is stable for 2 months when stored properly as intended by manufacturer at 2-8°C.

In case of partial use of the kit, the components should be stored in the following way:

- strips that remain unused must be carefully sealed with the plate sealing tape and stored at +2...+8°C within 2 months;
- Substrate Solution, Stop Solution, and Washing Solution concentrate after opening the vial, can be stored tightly closed at +2...+8°C until the kit's shelf life;
- Conjugate Solution, Calibrators and Control Serum after opening the vial, can be stored tightly closed at +2...+8°C within 2 months;
- *NOTE: Single freezing of Calibrators and Control Serum in aliquots is allowed*
- diluted washing solution can be stored at room temperature (+18...+25°C) for up to 5 days or at +2...+8°C for up to 14 days.

Kits that were stored in violation of the storage condition cannot be used.

8.3. Disposal

Expired kit components, used reagents and materials, as well as residual samples must be inactivated and disposed of in accordance with legal requirements.

9. REAGENTS PREPARATION

9.1. All reagents (including microstrips) and test samples should be allowed to reach room temperature (+18...+25 °C) for at least 30 minutes before use.

9.2. Microplate preparation

Open the package with the microplate and install the required number of strips into the frame. Unused strips must be sealed with plate sealing tape to prevent moisture from affecting the plate's holes and placed back in the bag.

9.3. Washing Solution preparation

Add the contents of the 22 mL Washing Solution concentrate vial to 550 mL of distilled or deionized water and mix thoroughly. In case of partial use of the kit, take the necessary amount of washing solution concentrate and dilute it 26 times with distilled or deionized water.

The spending of the components in case of partial use of the kit is given in the table:

Quantity of strips	1	2	3	4	5	6	7	8	9	10	11	12
Volume of the Washing Solution concentrate, mL	1.8	3.6	5.4	7.2	9	10.8	12.6	14.4	16.2	18	19.8	22
Volume of water, mL	45	90	135	180	225	270	315	360	405	450	495	550

10. ASSAY PROCEDURE

- 10.1 Put the desired number of strips into the frame based on the number of test samples in 2 replicates and 12 wells for Calibrators and Control Serum (2 wells for each calibrator (CAL 1-5) and 2 wells for control serum (Q)).
- 10.2 Dispense **25 µL of Calibrators and Control Serum as well as 25 µL of test serum/plasma samples (SAMP)** to the wells of the microplate according to the scheme below. The introduction of Calibrators, Control Serum and test samples should be carried out within 5 minutes to ensure equal incubation time for the first and last samples.

Note: during performing several independent series of tests, Calibrators, and Control Sample should be used each time.

Scheme of introduction of samples

	1	2	3	4	5	6	7	8	9	10	11	12
A	CAL1	CAL1	SAMP3	SAMP3	SAMP11	SAMP11						
B	CAL2	CAL2	SAMP4	SAMP4	SAMP12	SAMP12						
C	CAL3	CAL3	SAMP5	SAMP5								
D	CAL4	CAL4	SAMP6	SAMP6								
E	CAL5	CAL5	SAMP7	SAMP7								
F	Q	Q	SAMP8	SAMP8								
G	SAMP1	SAMP1	SAMP9	SAMP9								
H	SAMP2	SAMP2	SAMP10	SAMP10								

- 10.3 Add **100 µL of the Conjugate Solution** to all wells.
- 10.4 Carefully mix the contents of the microplate in a circular motion on a horizontal surface, cover strips with a plate sealing tape and incubate for **60 minutes at +37°C**.
- 10.5 At the end of the incubation period, remove and discard the plate cover. Aspirate and wash each well 5 times using an automatic washer or an 8-channel dispenser. For each washing, add 300 µL of Washing Solution (see 9.3) to all wells, then remove the liquid by aspiration or decantation. The residual volume of the Washing Solution after each aspiration or decantation should be no more than 5µL. After washing, carefully remove the remaining liquid from the wells on the absorbent paper. For the automatic washer/analyzer, the Washing Solution volume can be increased to 350 µL.
- 10.6 Add **100 µL of Substrate Solution** to all wells. The introduction of the substrate solution into the wells must be carried out within 2-3 minutes. Incubate the microplate in the dark **at room temperature (+18...+25°C) for 15 minutes**.
- 10.7 Add **100 µL of Stop Solution** to all wells in the same order as the substrate solution. After adding the Stop Solution, the contents of the wells turn yellow.
- 10.8 Read the optical density (OD) of the wells at 450nm using a microplate photometer within 5 minutes of adding the Stop Solution.
- 10.9 Plot a calibration curve in semi-logarithmic coordinates: (x) is the decimal logarithm of the T4 concentration in the calibrators nmol/L, (y) – OD versus T4 concentration (OD 450 nm). Manual or computerized data reduction is applicable at this stage. Point-by-point or linear data reduction is recommended due to non-linear shape of curve. Adjust the concentration of CAL1 to an infinitesimally small value, for example, 0.001 nmol/L.
- 10.10 Determine the corresponding concentration of T4 in tested samples from the calibration curve.

11. TEST VALIDITY

The test run shall be considered valid if the OD of CAL1 is above 1.2, and the values of the Control Serum fall into the required range (see Quality control Data Sheet).

12. EXPECTED VALUES

12.1. Therapeutical consequences should not be based on results of IVD methods alone – all available clinical and laboratory findings should be used by a physician to elaborate therapeutically measures. Each laboratory should establish its own normal range for T4. Based on data obtained by XEMA, the following normal range is recommended (see below). NOTE: the patients that have received murine monoclonal antibodies for radioimaging or immunotherapy develop high titered anti-mouse antibodies (HAMA). The presence of these antibodies may cause false results in the present assay. Sera from HAMA positive patients should be treated with depleting adsorbents before assaying

NOTE: values of T4 concentrations in the tested samples that are below the LoD (3.0 nmol/L) and also exceed the value of the upper calibrator (320 nmol/L) should be provided in the following form : «the T4 concentration of tested sample X is «lower than 3.0 nmol/L» or «higher than 320 nmol/L».

12.2. The calibrators concentration values of the T4 EIA kit are expressed in nmol/L. To calculate concentrations in µg/dl, the received concentration value in nmol/L shall be multiplied by 0.0775.

1 nmol/L = 0.0775 µg/dl

Sex, age	Units, nmol/L		Units alternative, µg/dl	
	Lower limit	Upper limit	Lower limit	Upper limit
Healthy donors	60	160	4.7	12.4
Males				
>61 yrs	60	129	4.7	10.0
Females				
>61 yrs	70	135	5.4	10.5
Children				
1-5 yrs	90	190	7.0	14.7
6-10 yrs	83	170	6.4	13.2
>10 yrs	60	160	4.7	12.4

13. PERFORMANCE CHARACTERISTICS

13.1. Analytical performance characteristics

3.1.1 Precision of Measurement

Repeatability (Intra assay repeatability) was determined by evaluation the coefficient of variation (CV) for 2 different samples during 1 day in 24 replicates on one series of ELISA kit.

Sample	Concentration, nmol/L	CV, %
1	17.5	4.36
2	110.7	3.67

Reproducibility (Inter assay reproducibility) was determined by evaluating the coefficients of variation for 2 samples during 5 days in 8-replicate determinations.

Sample	Concentration, nmol/L	CV, %
1	16.4	1.17
2	111.1	5.43

Reproducibility between lots was investigated by testing samples for one day on three lots. Each sample was run in 8 replicates.

Sample	Concentration1, nmol/L	Concentration2, nmol/L	Concentration3, nmol/L	CV, %
1	14.59	13.67	15.39	5.92
2	116.23	114.53	120.13	2.45

13.1.2 Trueness

The trueness of measurement is the degree of closeness of the average value obtained from a large number of measurement results to the true value. The bias of the measurement result (bias of measurements) is the difference between the mathematical expectation of the measurement result and the true value of the measurand. The bias was calculated for each sample and it was determined that it corresponds to the specified limits of $\pm 10\%$.

13.1.3 Linearity

Linearity was determined using sera samples with known T4 concentration (low and high) and mixing them with each other and buffer solution in different proportions. According to the measurements, linear range of kit is 0.75–15 nmol/L $\pm 10\%$.

13.1.4 Analytical sensitivity

Limit of detection (LoD) – the lowest T4 concentration in the serum or plasma sample that is detected by the T4 EIA kit is no lower than 3 nmol/L.

Limit of quantification (LoQ) – the lowest concentration of the analyte in the sample that is determined quantitatively with the declared trueness for T4 EIA kit is 32 nmol/L.

13.1.5 Analytical specificity

For the analysis result is not affected by the presence in the sample of bilirubin in a concentration of up to 0.21 mg/mL and hemoglobin in a concentration of up to 10 mg/mL.

The cross-reactivity of T4 with other analytes is shown in the table:

Analyte	Cross-reactivity, %
T3	0.5
D-Thyroxin	30

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LOT _____ DATE _____

	Manufacturer
	<i>In vitro</i> diagnostic medical device
	Catalogue number
 YYYY-MM	Use-by date
	Batch code
	Temperature limit
	Contains sufficient for <n> tests
	Caution
	Consult instructions for use
	Conformity Marking with technical regulations in Ukraine
	Authorized representative in the European Community/European Union
	CE Conformity Marking

**For any issues related to operation of the kit and technical support,
please contact by telefon number**

+38 044 294-69-78

or write to:

qa@xema.com.ua

XEMA LLC
Akademika Yefremova St. 23
03179, Kyiv, Ukraine
tel.:+38 044 422-62-16
tel.:+38 044 294-69-78
E-mail: qa@xema.com.ua
www.xema.in.ua

Instruction for use
A solid-phase enzyme immunoassay kit
for the quantitative determination of
free triiodothyronine in human serum or plasma

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







ASSAY PROCEDURE

CONTENT

1. INTENDED USE	2
2. GENERAL INFORMATION	2
3. TEST PRINCIPLE	2
4. KIT COMPONENTS	3
5. EQUIPMENT AND MATERIAL REQUIRED BUT NOT PROVIDED	4
6. WARNING AND PRECAUTIONS	4
7. SPECIMEN COLLECTION, TRANSPORTATION AND STORAGE OF SAMPLES	5
8. TRANSPORTATION AND STORAGE TERMS OF KIT, WASTE DISPOSAL	5
9. REAGENTS PREPARATION	6
10. ASSAY PROCEDURE	6
11. TEST VALIDITY	8
12. EXPECTED VALUES	8
13. PERFORMANCE CHARACTERISTIC	8
14. REFERENCES	10
SAMPLES IDENTIFICATION PLAN	11

Instruction for use
A solid-phase enzyme immunoassay kit
for the quantitative determination of
free triiodothyronine in human serum or plasma
fT3 EIA

1. INTENDED USE

The fT3 EIA kit is an enzyme immunoassay, intended for the quantitative determination of free triiodothyronine in human serum or plasma.

The field of application is clinical laboratory diagnostics.

2. GENERAL INFORMATION

Thyroxine (T4) and triiodothyronine (T3) are hormones that are produced by the thyroid gland and circulate in the blood both free and bound - mainly with thyroxine-binding globulin (TBG). Only free T3 and T4 are characterized by Hormonal activity, but their share is very small: 0.03% of the total content for T4 and 0.3% - for T3.

The concentration of T3 is much less than that of T4 but its metabolic activity is about 3 times greater. About 80% of T3 is produced in peripheral tissues by deiodination of T4, and only a small amount of it is secreted by thyroid gland. That is why in hypothyroid patients T3 level may for a long time remain on the lower limit of the normal range, because its loss may be compensated by enhanced conversion of T4 into T3.

The determination of total and free T3 concentration is carried out at the initial stage of hyperthyroidism, in case of recurrence of hyperthyroidism, in the differential diagnosis of hyperthyroidism, in case of a symptomatic increase of the T3 level, in case of acute hyperthyroidism after suppressive therapy with L-thyroxine.

3. TEST PRINCIPLE

Determination of the fT3 is based on competition principle of the enzyme immunoassay. Microwells plate is coated with specific rabbit polyclonal to T3 antibodies. fT3 conjugated to the horseradish peroxidase is used as enzyme conjugate. The analysis procedure includes two stages of incubation:

- during the first stage fT3 from the specimen competes with the conjugated fT3 for coating antibodies. As a result, a complex bounded to the solid phase and containing peroxidase is formed.

- during the second stage, the complexes formed due the reaction with the chromogen 3,3',5,5'-tetramethylbenzidine are visualized.

After stopping the reaction with a stop solution, the intensity of the color of the microwells is measured. Optical density in the microwell is inversely related to the quantity of the measured fT3 in the specimen of the serum (plasma).

The concentration is determined according to the calibration graph of the dependence of the optical density on the content of fT3 in the calibration samples.

4. KIT COMPONENTS

Code of component	Symbol	Name	Volume	Qty, pcs.	Description
P213Z	SORB MTP	Microplate	-	1	96-well polystyrene strip microplate coated with rabbit polyclonal antibodies to T3; ready to use
C213Z	CAL 1	Calibrator C1	0.5 mL	1	Solution based on human plasma, free of FT3, with preservative, ready to use (yellow liquid)
C213Z	CAL 2-6	Calibrators	0.5 mL	5	Solutions based on human plasma, containin 2,5; 5; 10; 20 and 40 pmol/L of FT3, with preservative, ready to use (blue liquids)
Q213Z	CONTROL	Control Serum	0.5 mL	1	Solution based on human plasma, containing of known FT3 content, with preservative, ready to use (colourless liquid)
T213Z	CONJ HRP	Conjugate Solution	14 mL	1	Solution of FT3 conjugated to the horseradish peroxidase; ready to use (blue liquid)
R055Z	SUBS TMB	Substrate Solution	14 mL	1	Tetramethylbenzidine (TMB) substrate solution; ready to use (colourless liquid)
S008Z	BUF WASH 26X	26x Concentrate Washing Solution	22 mL	1	Buffer solution with detergent, 26x concentrate (colourless liquid)
R050Z	STOP	Stop Solution	14 mL	1	5.0% solution of sulphuric acid; ready to use (colourless liquid)
The kit also includes instruction for use, quality control data sheet and plate sealing tape (2 pcs.)					

5. EQUIPMENT AND MATERIAL REQUIRED BUT NOT PROVIDED

- microplate photometer with 450 nm wavelength;
- dry thermostat for $+37^{\circ}\text{C}\pm 2^{\circ}\text{C}$;
- automatic plate washer (optional);
- micropipettes with variable volume, range volume 5-1000 μL ;
- graduated cylinder of 1000 mL capacity;
- distilled or deionized water;
- timer;
- vortex mixer;
- disposable gloves;
- absorbent paper.

6. WARNING AND PRECAUTIONS

In order to prevent incorrect results, strictly follow the recommended order and duration of the analysis procedure.

6.1. The kit is for *in vitro* diagnostic use only. For professional laboratory use.

6.2. Follow the rules mentioned below during the kit using:

- do not use kit beyond expire date;
- do not use the kit if its packaging is damaged;
- in order to avoid contamination, use new tips to pipette samples and reagents;
- use only verified equipment;
- close each vial with its own cap, after using the reagent;
- do not use components of other kits or reagents of other manufacturers;
- do not let wells dry after completing the rinsing step; immediately proceed to the next stage;
- avoid bubbles when adding reagents.

ATTENTION! The TMB substrate solution is light sensitive. Avoid prolonged exposure of the component to light.

6.3. Some kit components, such as stop solution, substrate solution, and washing solution, may cause toxic or irritant effects. If they get on the skin or mucosa, the affected area should be washed with plenty of running water.

6.4. All human products, including patient samples, should be considered potentially infectious. Handling and disposal should be in accordance with the procedures defined by an appropriate national biohazard safety guidelines or regulations.

6.5. The Calibrators and Control Serum included in the kit are negative for antibodies to HIV 1,2, hepatitis C virus and HBsAg, but the reagents should be considered as potentially infectious material and handled carefully.

6.6. Specimens must not contain any azide compounds, as they inhibit activity of peroxidase.

6.7. Wear protective gloves, protective clothing, eye protection, face protection.

6.8. Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.

6.9. Safety Data Sheet for this product is available upon request directly from XEMA LLC.

6.10. Serious incidents related to the kit must be reported to the manufacturer, Authorized Representative, and to the Competent Authority of the EU member state(s) where the incident has occurred.

7. SPECIMEN COLLECTION, TRANSPORTATION AND STORAGE OF SAMPLES

7.1. Blood sampling should be carried out from the cubital vein with a disposable needle using a vacuum blood sampling system. Serum or plasma specimens should be clearly labeled and identified. Serum must be separated from the clot as early as possible to avoid hemolysis of red blood cells. If there are any visible particles in the sample, they should be removed by centrifugation at 3000-5000 rpm for 20 minutes at room temperature or by filtration.

Don't use samples with high lipidemia, hemolysis as they may give false test results.

7.2. Specimen should be stored at +2...+8°C up to 3 days. Specimen held for a longer time, should be placed in a freezer at -15°C or below, do not refreeze/thaw samples.

7.3. For the transportation of samples, it is recommended to use triple packaging. The primary package is the labeled tube containing the sample. Secondary packaging is a polyethylene bag that is hermetically closed with a zip-lock. The outer packaging is a heat-insulating container, while the secondary packaging is placed in the outer packaging for transportation in the center of the thermal container. Frozen refrigerants are placed on the bottom, along the side walls of the thermal container, and cover the samples with them.

8. TRANSPORTATION AND STORAGE TERMS OF KIT, WASTE DISPOSAL

Information about the singularity storage conditions, transportation of the kit, and disposal of waste should be taken into account by all persons who participate in these processes.

8.1. Transportation

The FT3 EIA kit should be transported in the manufacturer's packaging at +2...+8°C. Single transportation at the temperature up to 25°C for 5 days is acceptable.

8.2. Storage

The FT3 EIA kit should be stored in the manufacturer's packaging at +2...+8°C. Do not freeze.

The kit contains reagents sufficient for 96 determinations including Calibrators and Control Serum.

Once opened test-kit is stable for 2 months when stored properly as intended by manufacturer at 2-8°C.

In case of partial use of the kit, the components should be stored in the following way:

- strips that remain unused must be carefully sealed with the plate sealing tape and stored at +2...+8°C within 2 months;
- Substrate Solution, Stop Solution, and Washing Solution concentrate after opening the vial, can be stored tightly closed at +2...+8°C until the kit's shelf life;
- Conjugate Solution, Calibrators and Control Serum after opening the vial, can be stored tightly closed at +2...+8°C within 2 months;
- diluted washing solution can be stored at room temperature (+18...+25°C) for up to 5 days or at +2...+8°C for up to 14 days.

Kits that were stored in violation of the storage condition cannot be used.

8.3. Disposal

Expired kit components, used reagents and materials, as well as residual samples must be inactivated and disposed of in accordance with legal requirements.

9. REAGENTS PREPARATION

All reagents (including microstrips) and test samples should be allowed to reach room temperature (+18...+25 °C) for at least 30 minutes before use.

9.2. Microplate preparation

Open the package with the microplate and install the required number of strips into the frame. Unused strips must be sealed with plate sealing tape to prevent moisture from affecting the plate's holes and placed back in the bag.

9.3. Washing Solution preparation

Add the contents of the 22 mL Washing Solution concentrate vial to 550 mL of distilled or deionized water and mix thoroughly. In case of partial use of the kit, take the necessary amount of washing solution concentrate and dilute it 26 times with distilled or deionized water.

The spending of the components in case of partial use of the kit is given in the table:

Quantity of strips	1	2	3	4	5	6	7	8	9	10	11	12
Volume of the Washing Solution concentrate, mL	1.8	3.6	5.4	7.2	9	10.8	12.6	14.4	16.2	18	19.8	22
Volume of water, mL	45	90	135	180	225	270	315	360	405	450	495	550

10. ASSAY PROCEDURE

- 10.1 Put the desired number of strips into the frame based on the number of test samples in 2 replicates and 14 wells for Calibrators and Control Serum (2 wells for each calibrator (CAL 1-6) and 2 wells for control serum (Q)).
- 10.2 Dispense **25 µL of Calibrators and Control Serum as well as 25 µL of test serum/plasma samples (SAMP)** to the wells of the microplate according to the scheme below. The introduction of Calibrators, Control Serum and test samples should be carried out within 5 minutes to ensure equal incubation time for the first and last samples.

Note: during performing several independent series of tests, Calibrators, and Control Sample should be used each time.

Scheme of introduction of samples

	1	2	3	4	5	6	7	8	9	10	11	12
A	CAL1	CAL1	SAMP2	SAMP2	SAMP10	SAMP10						
B	CAL2	CAL2	SAMP3	SAMP3	SAMP11	SAMP11						
C	CAL3	CAL3	SAMP4	SAMP4	SAMP12	SAMP12						
D	CAL4	CAL4	SAMP5	SAMP5								
E	CAL5	CAL5	SAMP6	SAMP6								
F	CAL6	CAL6	SAMP7	SAMP7								
G	Q	Q	SAMP8	SAMP8								
H	SAMP1	SAMP1	SAMP9	SAMP9								

- 10.3 Add **100 µL of the Conjugate Solution** to all wells.
- 10.4 Carefully mix the contents of the microplate in a circular motion on a horizontal surface, cover strips with a plate sealing tape and incubate for **60 minutes at +37°C**.
- 10.5 At the end of the incubation period, remove and discard the plate cover. Aspirate and wash each well 5 times using an automatic washer or an 8-channel dispenser. For each washing, add 300 µL of Washing Solution (see 9.3) to all wells, then remove the liquid by aspiration or decantation. The residual volume of the Washing Solution after each aspiration or decantation should be no more than 5µL. After washing, carefully remove the remaining liquid from the wells on the absorbent paper. For the automatic washer/analyzer, the Washing Solution volume can be increased to 350 µL.
- 10.6 Add **100 µL of Substrate Solution** to all wells. The introduction of the substrate solution into the wells must be carried out within 2-3 minutes. Incubate the microplate in the dark **at room temperature (+18...+25°C) for 15 minutes**.
- 10.7 Add **100 µL of Stop Solution** to all wells in the same order as the substrate solution. After adding the Stop Solution, the contents of the wells turn yellow.
- 10.8 Read the optical density (OD) of the wells at 450nm using a microplate photometer within 5 minutes of adding the Stop Solution.
- 10.9 Plot a calibration curve in semi-logarithmic coordinates: (x) is the decimal logarithm of the FT3 concentration in the calibrators pmol/L, (y) – OD versus FT3 concentration (OD 450 nm). Manual or computerized data reduction is applicable at this stage. Point-by-point or linear data reduction is recommended due to non-linear shape of curve. Adjust the concentration of CAL1 to an infinitesimally small value, for example, 0.001 pmol/L.
- 10.10 Determine the corresponding concentration of FT3 in tested samples from the calibration curve.

11. TEST VALIDITY

The test run shall be considered valid if the OD of CAL1 is above 1.2, and the values of the Control Serum fall into the required range (see Quality control Data Sheet).

12. EXPECTED VALUES

Therapeutical consequences should not be based on results of IVD methods alone – all available clinical and laboratory findings should be used by a physician to elaborate therapeutically measures. Each laboratory should establish its own normal range for FT3. Based on data obtained by XEMA, the following normal range is recommended (see below). NOTE: the patients that have received murine monoclonal antibodies for radioimaging or immunotherapy develop high titered anti-mouse antibodies (HAMA). The presence of these antibodies may cause false results in the present assay. Sera from HAMA positive patients should be treated with depleting adsorbents before assaying.

NOTE: values of FT3 concentrations in the tested samples that are below the LoD (0.5 pmol/L) and also exceed the value of the upper calibrator (40 pmol/L) should be provided in the following form: «the FT3 concentration of tested sample X is «lower than 0.5 pmol/L» or «higher than 40pmol/L».

Sex, age	Units, pmol/L	
	Lower limit	Upper limit
Healthy donors	2.5	5.8

13. PERFORMANCE CHARACTERISTICS

13.1. Analytical performance characteristics

3.1.1 Precision of Measurement

Repeatability (Intra assay repeatability) was determined by evaluation the coefficient of variation (CV) for 2 different samples during 1 day in 24 replicates on one series of ELISA kit.

Sample	Concentration, pmol/L	CV, %
1	4.32	7.44
2	6.87	5.14

Reproducibility (Inter assay reproducibility) was determined by evaluating the coefficients of variation for 2 samples during 5 days in 8-replicate determinations..

Sample	Concentration, pmol/L	CV, %
1	2,34	7,12
2	3,83	6,41

Reproducibility between lots was investigated by testing samples for one day on three lots. Each sample was run in 8 replicates.

Sample	Concentration1, pmol/L	Concentration2, pmol/L	Concentration3, pmol/L	CV, %
1	5.17	5.42	4.78	6.54
2	3.61	3.78	3.45	9.6

13.1.2 Trueness

The trueness of measurement is the degree of closeness of the average value obtained from a large number of measurement results to the true value. The bias of the measurement result (bias of measurements) is the difference between the mathematical expectation of the measurement result and the true value of the measurand. The bias was calculated for each sample and it was determined that it corresponds to the specified limits of $\pm 10\%$.

13.1.3 Linearity

Linearity was determined using sera samples with known fT3 concentration (low and high) and mixing them with each other and buffer solution in different proportions. According to the measurements, linear range of kit is 2.5–40 pmol/L $\pm 10\%$.

13.1.4 Analytical sensitivity

Limit of detection (LoD) – the lowest fT3 concentration in the serum or plasma sample that is detected by the fT3 EIA kit is no lower than 2.0 pmol/L.

Limit of quantification (LoQ) – the lowest concentration of the analyte in the sample that is determined quantitatively with the declared trueness for fT3 EIA kit is 2.25 pmol/L.

13.1.5 Analytical specificity

For the analysis result is not affected by the presence in the sample of bilirubin in a concentration of up to 0.21 mg/mL and hemoglobin in a concentration of up to 10 mg/mL.

The cross-reactivity of fT3 with other analytes is shown in the table:

Analyte	Cross-reactivity, %
L-Thyroxin	0,01
D-Thyroxin	0,04

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LOT _____ DATE _____

	Manufacturer
	<i>In vitro</i> diagnostic medical device
	Catalogue number
 YYYY-MM	Use-by date
	Batch code
	Temperature limit
	Contains sufficient for <n> tests
	Caution
	Consult instructions for use
	Conformity Marking with technical regulations in Ukraine
	Authorized representative in the European Community/European Union
	CE Conformity Marking

**For any issues related to operation of the kit and technical support,
please contact by telefon number**

+38 044 294-69-78

or write to:

qa@xema.com.ua

XEMA LLC
Akademika Yefremova St. 23
03179, Kyiv, Ukraine
tel.:+38 044 422-62-16
tel.:+38 044 294-69-78
E-mail: qa@xema.com.ua
www.xema.in.ua

Instruction for use
A solid-phase enzyme immunoassay kit
for the quantitative determination of
free thyroxin in human serum or plasma

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ASSAY PROCEDURE

Dispensing of Calibrators, Control Serum and test samples

test samples
CAL
CONTROL
SORB MTP
25 µl

Conjugate Solution

CONJ HRP
100 µl

Incubation 1

+37 °C
60'

Washing 5 times

Washing 5 times

OD measuring, calculation of results

450 nm

Stop Solution

STOP
100 µl

Incubation 2

+18-25 °C
15'

Substrate Solution

SUBS TMB
100 µl

CONTENT

1. INTENDED USE	2
2. GENERAL INFORMATION	2
3. TEST PRINCIPLE	2
4. KIT COMPONENTS	3
5. EQUIPMENT AND MATERIAL REQUIRED BUT NOT PROVIDED	4
6. WARNING AND PRECAUTIONS	4
7. SPECIMEN COLLECTION, TRANSPORTATION AND STORAGE OF SAMPLES	5
8. TRANSPORTATION AND STORAGE TERMS OF KIT, WASTE DISPOSAL	5
9. REAGENTS PREPARATION	6
10. ASSAY PROCEDURE	6
11. TEST VALIDITY	7
12. EXPECTED VALUES	7
13. PERFORMANCE CHARACTERISTICS	8
14. REFERENCES	9
SAMPLES IDENTIFICATION PLAN	10

Instruction for use
A solid-phase enzyme immunoassay kit
for the quantitative determination of
free thyroxin in human serum or plasma
fT4 EIA

1. INTENDED USE

The fT4 EIA kit is an enzyme immunoassay, intended for the quantitative determination of free thyroxin in human serum or plasma.

The field of application is clinical laboratory diagnostics.

2. GENERAL INFORMATION

Thyroid hormones thyroxin (T4) and 3,5,3'-triiodothyronine (T3) exert regulatory influences on growth, differentiation, cellular metabolism and development of skeletal and organ systems. T4 and T3 in blood are found both in free and bound form – mostly, they are bound to thyroxin binding globulin (TBG). Only free forms of T3 and T4 exert hormonal activity also their percentage is very low – 0.3% for T3 and 0.03% for T4.

The concentration of T4 is generally accepted as an index of thyroid function which provide enough information to differentiate between hyper-, hypo- and euthyroidism.

Elevation of total T4 is found in hyperthyroidism, in patients with tumours of pituitary gland, in subjects with elevated TBG level (pregnancy, acute or chronic active hepatitis, estrogen-secreting tumours or estrogen intake, hereditary elevation of TBG), in patients taking oral contraceptives, heroin, methadone, thyroid preparations, TSH, thyroliberin.

Low total T4 is found in hypothyroidism, in patients with panhypopituitarism, in subjects with low TBG level (acromegaly, nephritic syndrome, hypoproteinemia, chronic liver diseases, androgen-secreting tumours, hereditary reduction), in patients taking aminosalicilic and acetylsalicilic acids, cholestyramine, reserpine, potassium iodide, triiodothyronine.

3. TEST PRINCIPLE

Determination of free thyroxin is based on competition principle of the enzyme immunoassay. Microwells plate is coated with specific murine monoclonal antibodies to T4. fT4 conjugated to the horseradish peroxidase is used as enzyme conjugate. The analysis procedure includes two stages of incubation:

- during the first stage fT4 from the specimen competes with the conjugated fT4 for coating antibodies. As a result, a complex bounded to the solid phase and containing peroxidase is formed.

- during the second stage, the complexes formed due the reaction with the chromogen 3,3',5,5'-tetramethylbenzidine are visualized.

After stopping the reaction with a stop solution, the intensity of the color of the microwells is measured. Optical density in the microwell is inversely related to the quantity of the measured fT4 in the specimen of the serum (plasma).

The concentration is determined according to the calibration graph of the dependence of the optical density on the content of fT4 in the calibration samples.

4. KIT COMPONENTS

Code of component	Symbol	Name	Volume	Qty, pcs.	Description
P214Z	SORB MTP	Microplate	-	1	96-well polystyrene strip microplate coated with murine monoclonal antibodies to T4; ready to use
C214Z	CAL 1	Calibrator C1	0.5 mL	1	Solution based on human plasma, free of FT4, with preservative, ready to use (yellow liquid)
C214Z	CAL 2-6	Calibrators	0.5 ml	5	Solutions based on human plasma, containing 5; 10; 25, 50 and 100 pmol/L of FT4, with preservative, ready to use (red liquids)
Q214Z	CONTROL	Control Serum	0.5 ml	1	Solution based on human plasma, containing of known FT4 content, with preservative, ready to use (colourless liquid)
T214Z	CONJ HRP	Conjugate Solution	12 ml	1	Solution of FT4 conjugated to the horseradish peroxidase; ready to use (red liquid)
R055Z	SUBS TMB	Substrate Solution	12 ml	1	Tetramethylbenzidine (TMB) substrate solution; ready to use (colourless liquid)
S008Z	BUF WASH 26X	26x Concentrate Washing Solution	22 ml	1	Buffer solution with detergent, 26x concentrate (colourless liquid)
R050Z	STOP	Stop Solution	12 ml	1	5.0% solution of sulphuric acid; ready to use (colourless liquid)
The kit also includes instruction for use, quality control data sheet and plate sealing tape (1 pcs.)					

5. EQUIPMENT AND MATERIAL REQUIRED BUT NOT PROVIDED

- microplate photometer with 450 nm wavelength;
- dry thermostat for $+37^{\circ}\text{C}\pm 1^{\circ}\text{C}$;
- automatic plate washer (optional);
- micropipettes with variable volume, range volume 5-1000 μL ;
- graduated cylinder of 1000 mL capacity;
- distilled or deionized water;
- timer;
- vortex mixer;
- disposable gloves;
- absorbent paper.

6. WARNING AND PRECAUTIONS

In order to prevent incorrect results, strictly follow the recommended order and duration of the analysis procedure.

6.1. The kit is for *in vitro* diagnostic use only. For professional laboratory use.

6.2. Follow the rules mentioned below during the kit using:

- do not use kit beyond expire date;
- do not use the kit if its packaging is damaged;
- in order to avoid contamination, use new tips to pipette samples and reagents;
- use only verified equipment;
- close each vial with its own cap, after using the reagent;
- do not use components of other kits or reagents of other manufacturers;
- do not let wells dry after completing the rinsing step; immediately proceed to the next stage;
- avoid bubbles when adding reagents.

ATTENTION! The TMB substrate solution is light sensitive. Avoid prolonged exposure of the component to light.

6.3. Some kit components, such as stop solution, substrate solution, and washing solution, may cause toxic or irritant effects. If they get on the skin or mucosa, the affected area should be washed with plenty of running water.

6.4. All human products, including patient samples, should be considered potentially infectious. Handling and disposal should be in accordance with the procedures defined by an appropriate national biohazard safety guidelines or regulations.

6.5. The Calibrators and Control Serum included in the kit are negative for antibodies to HIV 1,2, hepatitis C virus and HBsAg, but the reagents should be considered as potentially infectious material and handled carefully.

6.6. Specimens must not contain any azide compounds, as they inhibit activity of peroxidase.

6.7. Wear protective gloves, protective clothing, eye protection, face protection.

6.8. Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.

6.9. Safety Data Sheet for this product is available upon request directly from XEMA LLC.

6.10. Serious incidents related to the kit must be reported to the manufacturer, Authorized Representative, and to the Competent Authority of the EU member state(s) where the incident has occurred.

7. SPECIMEN COLLECTION, TRANSPORTATION AND STORAGE OF SAMPLES

7.1. Blood sampling should be carried out from the cubital vein with a disposable needle using a vacuum blood sampling system. Serum or plasma specimens should be clearly labeled and identified. Serum must be separated from the clot as early as possible to avoid hemolysis of red blood cells. If there are any visible particles in the sample, they should be removed by centrifugation at 3000-5000 rpm for 20 minutes at room temperature or by filtration.

Don't use samples with high lipidemia, hemolysis as they may give false test results.

7.2. Specimen should be stored at +2...+8°C up to 3 days. Specimen held for a longer time, should be placed in a freezer at -15°C or below; do not refreeze/thaw samples.

7.3. For the transportation of samples, it is recommended to use triple packaging. The primary package is the labeled tube containing the sample. Secondary packaging is a polyethylene bag that is hermetically closed with a zip-lock. The outer packaging is a heat-insulating container, while the secondary packaging is placed in the outer packaging for transportation in the center of the thermal container. Frozen refrigerants are placed on the bottom, along the side walls of the thermal container, and cover the samples with them.

8. TRANSPORTATION AND STORAGE TERMS OF KIT, WASTE DISPOSAL

Information about the singularity storage conditions, transportation of the kit, and disposal of waste should be taken into account by all persons who participate in these processes.

8.1. Transportation

The ft4 EIA kit should be transported in the manufacturer's packaging at +2...+8°C. Single transportation at the temperature up to 25°C for 5 days is acceptable.

8.2. Storage

The ft4 EIA kit should be stored in the manufacturer's packaging at +2...+8°C. Do not freeze.

The kit contains reagents sufficient for 96 determinations including Calibrators and Control Serum.

Once opened test-kit is stable for 2 months when stored properly as intended by manufacturer at 2-8°C.

In case of partial use of the kit, the components should be stored in the following way:

- the remaining strips should be immediately resealed in the bag along with the silica gel, closed with the zip-lock, and stored at +2...+8°C within 2 months
- Substrate Solution, Stop Solution, and Washing Solution concentrate after opening the vial, can be stored tightly closed at +2...+8°C until the kit's shelf life;
- Conjugate Solution, Calibrators and Control Serum after opening the vial, can be stored tightly closed at +2...+8°C within 2 months;
- *NOTE: Single freezing of Calibrators and Control Serum in aliquots is allowed.*
- diluted washing solution can be stored at room temperature (+18...+25°C) for up to 5 days or at +2...+8°C for up to 14 days.

Kits that were stored in violation of the storage condition cannot be used.

8.3. Disposal

Expired kit components, used reagents and materials, as well as residual samples must be inactivated and disposed of in accordance with legal requirements.

9. REAGENTS PREPARATION

9.1. All reagents (including microstrips) and test samples should be allowed to reach room temperature (+18...+25 °C) for at least 30 minutes before use.

9.2. Microplate preparation

Open the package with the microplate and install the required number of strips into the frame. The remaining strips should be immediately resealed in the bag along with the silica gel and closed with the zip-lock to prevent moisture from affecting the plate's strips.

9.3. Washing Solution preparation

Add the contents of the 22 mL Washing Solution concentrate vial to 550 mL of distilled or deionized water and mix thoroughly. In case of partial use of the kit, take the necessary amount of washing solution concentrate and dilute it 26 times with distilled or deionized water.

The spending of the components in case of partial use of the kit is given in the table:

Quantity of strips	1	2	3	4	5	6	7	8	9	10	11	12
Volume of the Washing Solution concentrate, mL	1.8	3.6	5.4	7.2	9	10.8	12.6	14.4	16.2	18	19.8	22
Volume of water, mL	45	90	135	180	225	270	315	360	405	450	495	550

10. ASSAY PROCEDURE

- 10.1 Put the desired number of strips into the frame based on the number of test samples in 2 replicates and 14 wells for Calibrators and Control Serum (2 wells for each calibrator (CAL 1-6) and 2 wells for control serum (Q)).
- 10.2 Dispense **25 µL of Calibrators and Control Serum as well as 25 µL of test serum/plasma samples (SAMP)** to the wells of the microplate according to the scheme below. The introduction of Calibrators, Control Serum and test samples should be carried out within 5 minutes to ensure equal incubation time for the first and last samples.

Note: during performing several independent series of tests, Calibrators, and Control Sample should be used each time.

Scheme of introduction of samples

	1	2	3	4	5	6	7	8	9	10	11	12
A	CAL1	CAL1	SAMP2	SAMP2	SAMP10	SAMP10						
B	CAL2	CAL2	SAMP3	SAMP3	SAMP11	SAMP11						
C	CAL3	CAL3	SAMP4	SAMP4	SAMP12	SAMP12						
D	CAL4	CAL4	SAMP5	SAMP5								
E	CAL5	CAL5	SAMP6	SAMP6								
F	CAL6	CAL6	SAMP7	SAMP7								
G	Q	Q	SAMP8	SAMP8								
H	SAMP1	SAMP1	SAMP9	SAMP9								

- 10.3 Add **100 µL of the Conjugate Solution** to all wells.
- 10.4 Carefully mix the contents of the microplate in a circular motion on a horizontal surface, cover strips with a plate sealing tape and incubate for **60 minutes at +37°C**.
- 10.5 At the end of the incubation period, remove and discard the plate cover. Aspirate and wash each well 5 times using an automatic washer or an 8-channel dispenser. For each washing, add 300 µL of Washing Solution (see 9.3) to all wells, then remove the liquid by aspiration or decantation. The residual volume of the Washing Solution after each aspiration or decantation should be no more than 5µL. After washing, carefully remove the remaining liquid from the wells on the absorbent paper. For the automatic washer/analyzer, the Washing Solution volume can be increased to 350 µL
- 10.6 Add **100 µL of Substrate Solution** to all wells. The introduction of the substrate solution into the wells must be carried out within 2-3 minutes. Incubate the microplate in the dark **at room temperature (+18...+25°C) for 15 minutes**.
- 10.7 Add **100 µL of Stop Solution** to all wells in the same order as the substrate solution. After adding the Stop Solution, the contents of the wells turn yellow.
- 10.8 Read the optical density (OD) of the wells at 450nm using a microplate photometer within 5 minutes of adding the Stop Solution.
- 10.9 Plot a calibration curve in semi-logarithmic coordinates: (x) is the decimal logarithm of the fT4 concentration in the calibrators pmol/L, (y) – OD versus fT4 concentration (OD 450 nm). Manual or computerized data reduction is applicable at this stage. Point-by-point or linear data reduction is recommended due to non-linear shape of curve. Adjust the concentration of CAL1 to an infinitesimally small value, for example, 0.001 pmol/L.
- 10.10 Determine the corresponding concentration of fT4 in tested samples from the calibration curve.

11. TEST VALIDITY

The test run shall be considered valid if the OD of CAL1 is above 1.2, and the values of the Control Serum fall into the required range (see Quality control Data Sheet).

12. EXPECTED VALUES

Therapeutical consequences should not be based on results of IVD methods alone – all available clinical and laboratory findings should be used by a physician to elaborate therapeutically measures. Each laboratory should establish its own normal range for fT4. Based on data obtained by XEMA, the following normal range is recommended (see below). NOTE: the patients that have received murine monoclonal antibodies for radioimaging or immunotherapy develop high titered anti-mouse antibodies (HAMA). The presence of these antibodies may cause false results in the present assay. Sera from HAMA positive patients should be treated with depleting adsorbents before assaying.

NOTE: values of fT4 concentrations in the tested samples that are below the LoD (0.75 pmol/L) and also exceed the value of the upper calibrator (100 pmol/L) should be provided in the following form : «the fT4 concentration of tested sample X is «lower than 0.75 pmol/L» or «higher than 100 pmol/L».

Sex, age	Units, pmol/L	
	Lower limit	Upper limit
Healthy donors		
< 60 yrs	10	25
> 60 yrs	10	21
Pregnancy week		
1st trimester	9	26
2nd trimester	6	21
3rd trimester	6	21

13. PERFORMANCE CHARACTERISTICS

13.1. Analytical performance characteristics

13.1.1 Precision of Measurement

Repeatability (Intra assay repeatability) was determined by evaluation the coefficient of variation (CV) for 2 different samples during 1 day in 24 replicates on one series of ELISA kit.

Sample	Concentration, pmol/L	CV, %
1	54.4	5.83
2	85.23	3.67

Reproducibility (Inter assay reproducibility) was determined by evaluating the coefficients of variation for 2 samples during 5 days in 8-replicate determinations.

Sample	Concentration, pmol/L	CV, %
1	54.36	1.15
2	85.73	3.23

Reproducibility between lots was investigated by testing samples for one day on three lots. Each sample was run in 8 replicates.

Sample	Concentration1, pmol/L	Concentration2, pmol/L	Concentration3, pmol/L	CV, %
1	54.59	52.67	60.39	7.19
2	85.23	87.53	85.13	1.58

13.1.2 Trueness

The trueness of measurement is the degree of closeness of the average value obtained from a large number of measurement results to the true value. The bias of the measurement result (bias of measurements) is the difference between the mathematical expectation of the measurement result and the true value of the measurand. The bias was calculated for each sample and it was determined that it corresponds to the specified limits of $\pm 10\%$.

13.1.3 Linearity

Linearity was determined using sera samples with known fT4 concentration (low and high) and mixing them with each other and buffer solution in different proportions. According to the measurements, linear range of kit is 5-100 pmol/L $\pm 10\%$.

13.1.4 Analytical sensitivity

Limit of detection (LoD) – the lowest fT4 concentration in the serum or plasma sample that is detected by the fT4 EIA kit is no lower than 0.75 pmol/L.

Limit of quantification (LoQ) – the lowest concentration of the analyte in the sample that is determined quantitatively with the declared trueness for fT4 EIA kit is 5 pmol/L.

13.1.5 Analytical specificity

For the analysis result is not affected by the presence in the sample of bilirubin in a concentration of up to 0.21 mg/mL and hemoglobin in a concentration of up to 10 mg/mL.

The cross-reactivity of fT4 with other analytes is shown in the table:

Analyte	Cross-reactivity, %
L-Thyroxin	100
D-Thyroxin	94
3,3',5'-Triiodo-L-Thyronine (Reverse T3)	86
3,3',5'-Triiodo-L-Thyronine (T3)	3.3
3,3',5'-Triiodo-D-Thyronine	1.8
3,3',5'-Triiodothyropropionic acid	0.6

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LOT _____ DATE _____

	Manufacturer
	<i>In vitro</i> diagnostic medical device
	Catalogue number
 YYYY-MM	Use-by date
	Batch code
	Temperature limit
	Contains sufficient for <n> tests
	Caution
	Consult instructions for use
	Conformity Marking with technical regulations in Ukraine
	Authorized representative in the European Community/European Union
	CE Conformity Marking

**For any issues related to operation of the kit and technical support,
please contact by telefon number**

+38 044 294-69-78

or write to:

qa@xema.com.ua

XEMA LLC
Akademika Yefremova St. 23
03179, Kyiv, Ukraine
tel.:+38 044 422-62-16
tel.:+38 044 294-69-78
E-mail: qa@xema.com.ua
www.xema.in.ua

Instruction for use
A solid-phase enzyme immunoassay kit
for the quantitative determination of
antigen CYFRA 21-1 in human serum or plasma

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
ASSAY PROCEDURE

CONTENT

1. INTENDED USE	2
2. GENERAL INFORMATION	2
3. TEST PRINCIPLE	3
4. KIT COMPONENTS	4
5. EQUIPMENT AND MATERIAL REQUIRED BUT NOT PROVIDED	5
6. WARNING AND PRECAUTIONS	5
7. SPECIMEN COLLECTION, TRANSPORTATION AND STORAGE OF SAMPLES	6
8. TRANSPORTATION AND STORAGE TERMS OF KIT, WASTE DISPOSAL	6
9. REAGENTS PREPARATION	7
10. ASSAY PROCEDURE	7
11. TEST VALIDITY	9
12. EXPECTED VALUES	9
13. PERFORMANCE CHARACTERISTICS	10
14. REFERENCES	11
SAMPLES IDENTIFICATION PLAN	12

Instruction for use
A solid-phase enzyme immunoassay kit
for the quantitative determination of
antigen CYFRA 21-1 in human serum or plasma
CYFRA 21-1 EIA

1. INTENDED USE

The CYFRA 21-1 EIA kit is an enzyme immunoassay, intended for the quantitative determination of antigen CYFRA 21-1 in human serum or plasma.

The field of application is clinical laboratory diagnostics.

2. GENERAL INFORMATION

The CYFRA 21-1 antigen is a fragment of cytokeratin 19, which is formed as a result of proteolysis and, unlike the main cytokeratin structure, is able to change into a soluble form and enter the systemic bloodstream.

The precursor molecule of the CYFRA 21-1 antigen - cytokeratin 19 - is expressed in all normal tissues, but a particularly high level of expression is observed in lung or bladder wall tumor cells.

Increased content of CYFRA 21-1 is observed in the blood of patients with lung tumors (mainly squamous cell carcinoma, less often adenocarcinoma and other histological forms) and bladder tumors. Determination of the level of the CYFRA 21-1 antigen is useful for monitoring the effectiveness of treatment and monitoring the course of these tumors; however, the results of the measurement of the CYFRA 21-1 antigen should always be interpreted in conjunction with the results of other research methods and clinical data.

3. TEST PRINCIPLE

The determination of the CYFRA 21-1 is based on the two-site sandwich enzyme immunoassay principle. On the inner surface of the microplate wells are immobilized specific murine monoclonal antibodies to soluble cytokeratin 8/19 (CYFRA 21-1). Second antibodies – murine monoclonal antibodies to human CYFRA 21-1 conjugated to the horseradish peroxidase is used as enzyme conjugate. The analysis procedure includes two stages of incubation:

- during the first stage CYFRA 21-1 from the specimen is captured by the antibodies coated onto the microwell surface, as well as horseradish peroxidase-conjugated monoclonal antibodies bind to free epitopes of immobilized CYFRA 21-1;

- during the second stage, the complexes formed due to the reaction with the chromogen 3,3',5,5'-tetramethylbenzidine are visualized.

After stopping the reaction with a stop solution, the intensity of the color of the microwells is measured. The optical density in the microwell is directly related to the quantity of the measured CYFRA 21-1 in the serum specimen (plasma). The concentration is determined according to the calibration graph of the dependence of the optical density on the content of CYFRA 21-1 in the calibration samples.

4. KIT COMPONENTS

Code of component	Symbol	Name	Volume	Qty, pcs.	Description
P236Z	SORB MTP	Microplate	-	1	96-well polystyrene strip microplate coated with murine monoclonal antibodies to CYFRA 21-1; ready to use
C236Z	CAL 1	Calibrator C1	2 mL	1	Solution based on phosphate buffer (pH 7.2-7.4), free of CYFRA 21-1, with preservative, ready to use (colourless liquid)
C236Z	CAL 2-5	Calibrators	0.8 mL	4	Solution based on phosphate buffer (pH 7.2-7.4), containing 3; 10; 25 and 50 ng/mL of CYFRA 21-1, with preservative, ready to use (red liquids)
Q236Z	CONTROL	Control Serum	0.8 mL	1	Solution based on human serum, containing of known CYFRA 21-1 content, with preservative, ready to use (colourless liquid)
T236Z	CONJ HRP	Conjugate Solution	6 mL	1	Solution of murine monoclonal antibodies to CYFRA 21-1 conjugated to the horseradish peroxidase; ready to use (red liquid)
R055Z	SUBS TMB	Substrate Solution	14 mL	1	Tetramethylbenzidine (TMB) substrate solution; ready to use (colourless liquid)
S008Z	BUF WASH 26X	26x Concentrate Washing Solution	22 mL	1	Buffer solution with detergent, 26x concentrate (colourless liquid)
R050Z	STOP	Stop Solution	14 mL	1	5.0% solution of sulphuric acid; ready to use (colourless liquid)
The kit also includes instruction for use, quality control data sheet and plate sealing tape (2 pcs.)					

5. EQUIPMENT AND MATERIAL REQUIRED BUT NOT PROVIDED

- microplate photometer with 450 nm wavelength or 450\620-680 nm;
- dry thermostat for 37 °C±2 °C;
- automatic plate washer (optional);
- micropipettes with variable volume, range volume 5-1000 µL;
- graduated cylinder of 1000 mL capacity;
- distilled or deionized water;
- timer;
- vortex mixer;
- disposable gloves;
- absorbent paper.

6. WARNING AND PRECAUTIONS

In order to prevent incorrect results, strictly follow the recommended order and duration of the analysis procedure.

6.1. The kit is for *in vitro* diagnostic use only. For professional laboratory use.

6.2. Follow the rules mentioned below during the kit using:

- do not use kit beyond expire date;
- do not use the kit if its packaging is damaged;
- in order to avoid contamination, use new tips to pipette samples and reagents;
- use only verified equipment;
- close each vial with its own cap, after using the reagent;
- do not use components of other kits or reagents of other manufacturers;
- do not let wells dry after completing the rinsing step; immediately proceed to the next stage;
- avoid bubbles when adding reagents.

ATTENTION! The TMB substrate solution is light sensitive. Avoid prolonged exposure of the component to light.

6.3. Some kit components, such as stop solution, substrate solution, and washing solution, may cause toxic or irritant effects. If they get on the skin or mucosa, the affected area should be washed with plenty of running water.

6.4. All human products, including patient samples, should be considered potentially infectious. Handling and disposal should be in accordance with the procedures defined by an appropriate national biohazard safety guidelines or regulations.

6.5. The Calibrators and Control Serum included in the kit are negative for antibodies to HIV 1,2, hepatitis C virus and HBsAg, but the reagents should be considered as potentially infectious material and handled carefully.

6.6. Specimens must not contain any azide compounds, as they inhibit activity of peroxidase.

6.7. Wear protective gloves, protective clothing, eye protection, face protection.

6.8. Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.

6.9. Safety Data Sheet for this product is available upon request directly from XEMA LLC.

6.10. Serious incidents related to the kit must be reported to the manufacturer, Authorized Representative, and to the Competent Authority of the EU member state(s) where the incident has occurred.

7. SPECIMEN COLLECTION, TRANSPORTATION AND STORAGE OF SAMPLES

7.1. Blood sampling should be carried out from the cubital vein with a disposable needle using a vacuum blood sampling system. Serum or plasma specimens should be clearly labeled and identified. Serum must be separated from the clot as early as possible to avoid hemolysis of red blood cells. If there are any visible particles in the sample, they should be removed by centrifugation at 3000-5000 rpm for 20 minutes at room temperature or by filtration.

Don't use samples with high lipidemia, hemolysis as they may give false test results.

7.2. Specimen should be stored at +2...+8°C up to 3 days. Specimen held for a longer time, should be placed in a freezer at -15°C or below, do not refreeze/thaw samples.

7.3. For the transportation of samples, it is recommended to use triple packaging. The primary package is the labeled tube containing the sample. Secondary packaging is a polyethylene bag that is hermetically closed with a zip-lock. The outer packaging is a heat-insulating container, while the secondary packaging is placed in the outer packaging for transportation in the center of the thermal container. Frozen refrigerants are placed on the bottom, along the side walls of the thermal container, and cover the samples with them.

8. TRANSPORTATION AND STORAGE TERMS OF KIT, WASTE DISPOSAL

Information about the singularity storage conditions, transportation of the kit, and disposal of waste should be taken into account by all persons who participate in these processes.

8.1. Transportation

The CYFRA 21-1 EIA kit should be transported in the manufacturer's packaging at +2...+8°C. Single transportation at the temperature up to 25°C for 5 days is acceptable.

8.2. Storage

The CYFRA 21-1 EIA kit should be stored in the manufacturer's packaging at +2...+8°C. Do not freeze.

The kit contains reagents sufficient for 96 determinations including Calibrators and Control Serum.

Once opened test-kit is stable for 2 months when stored properly as intended by manufacturer at 2-8°C.

In case of partial use of the kit, the components should be stored in the following way:

- strips that remain unused must be carefully sealed with the plate sealing tape and stored at +2...+8°C within 2 months;
- Substrate Solution, Stop Solution, and Washing Solution concentrate after opening the vial, can be stored tightly closed at +2...+8°C until the kit's shelf life;
- Conjugate Solution, Calibrators and Control Serum after opening the vial, can be stored tightly closed at +2...+8°C within 2 months;
- diluted Washing Solution can be stored at room temperature (+18...+25°C) for up to 5 days or at +2...+8°C for up to 14 days.

Kits that were stored in violation of the storage condition cannot be used.

8.3. Disposal

Expired kit components, used reagents and materials, as well as residual samples must be inactivated and disposed of in accordance with legal requirements.

9. REAGENTS PREPARATION

9.1. All reagents (including microstrips) and test samples should be allowed to reach room temperature (+18...+25 °C) for at least 30 minutes before use.

9.2. Microplate preparation

Open the package with the microplate and install the required number of strips into the frame. Unused strips must be sealed with plate sealing tape to prevent moisture from affecting the plate's holes and placed back in the bag.

9.3. Washing Solution preparation

Add the contents of the 22 mL Washing Solution concentrate vial to 550 mL of distilled or deionized water and mix thoroughly. In case of partial use of the kit, take the necessary amount of Washing Solution concentrate and dilute it 26 times with distilled or deionized water.

The spending of the components in case of partial use of the kit is given in the table:

Quantity of strips	1	2	3	4	5	6	7	8	9	10	11	12
Volume of the Washing Solution concentrate, mL	1.8	3.6	5.4	7.2	9	10.8	12.6	14.4	16.2	18	19.8	22
Volume of water, mL	45	90	135	180	225	270	315	360	405	450	495	550

9.4. Samples preparation

If suggested analyte concentration in the sample exceeds the 50 ng/mL, additionally dilute this sample accordingly, using (Calibrator C1). Use of other buffers or reagents for sample dilution may lead to incorrect measurement.

NOTE: in order to obtain reliable results, we recommend to use several successive dilutions of the blood serum (plasma) sample

Do not dilute Control Serum and Calibrators!

10. ASSAY PROCEDURE

10.1 Put the desired number of strips into the frame based on the number of test samples in 2 replicates and 12 wells for Calibrators and Control Serum (2 wells for each Calibrator (CAL 1-5) and 2 wells for Control Serum (Q)).

10.2 If necessary, dilute the test samples as described in 9.4.

10.3 Dispense **50 µL of Calibrators and Control Serum as well as 50 µL of test serum/plasma samples (SAMP)** to the wells of the microplate according to the scheme below. The introduction of Calibrators, Control Serum and test samples should be carried out within 5 minutes to ensure equal incubation time for the first and last samples.

NOTE: during performing several independent series of tests, Calibrators, and Control Serum should be used each time.

10.4 Dispense **50 µL of Conjugate Solution** to all wells.

Scheme of introduction of samples

	1	2	3	4	5	6	7	8	9	10	11	12
A	CAL1	CAL1	SAMP3	SAMP3	SAMP11	SAMP11						
B	CAL2	CAL2	SAMP4	SAMP4	SAMP12	SAMP12						
C	CAL3	CAL3	SAMP5	SAMP5								
D	CAL4	CAL4	SAMP6	SAMP6								
E	CAL5	CAL5	SAMP7	SAMP7								
F	Q	Q	SAMP8	SAMP8								
G	SAMP1	SAMP1	SAMP9	SAMP9								
H	SAMP2	SAMP2	SAMP10	SAMP10								

- 10.5 Carefully mix the contents of the microplate in a circular motion on a horizontal surface, cover strips with a plate sealing tape and incubate for **60 minutes at +37°C**.
- 10.6 At the end of the incubation period, remove and discard the plate cover. Aspirate and wash each well 5 times using an automatic washer or an 8-channel dispenser. For each washing, add 300 µL of Washing Solution (see 9.3) to all wells, then remove the liquid by aspiration or decantation. The residual volume of the Washing Solution after each aspiration or decantation should be no more than 5µL. After washing, carefully remove the remaining liquid from the wells on the absorbent paper. For the automatic washer/analyzer, the wash solution volume can be increased to 350 µL.
- 10.7 Add **100 µL of Substrate Solution** to all wells. The introduction of the Substrate Solution into the wells must be carried out within 2-3 minutes. Incubate the microplate in the dark **at room temperature (+18...+25°C) for 15 minutes**.
- 10.8 Add **100 µL of Stop Solution** to all wells in the same order as the Substrate Solution. After adding the Stop Solution, the contents of the wells turn yellow.
- 10.9 Read the optical density (OD) of the wells at 450nm and reference light filters 620–680 nm using a microplate photometer within 5 minutes of adding the Stop Solution. Set photometer blank on CAL1.
- 10.10 Plot a calibration curve in linear coordinates: (x) is the CYFRA 21-1 concentration in the calibrators ng/mL, (y) – OD versus CYFRA 21-1 concentration (OD 450 nm / 620–680 nm). Manual or computerized data reduction is applicable at this stage. Point-by-point or linear data reduction is recommended due to non-linear shape of curve.
- 10.11 Determine the corresponding concentration of CYFRA 21-1 in tested samples from the calibration curve. In the case of preliminary dilution of the test sample (see 9.4), the obtained result should be multiplied by the dilution factor.

11. TEST VALIDITY

The test run shall be considered valid if the OD of CAL1 is above 0.15, and the values of the Control Serum fall into the required range (see Quality control Data Sheet).

12. EXPECTED VALUES

Therapeutical consequences should not be based on results of IVD methods alone – all available clinical and laboratory findings should be used by a physician to elaborate therapeutically measures. Each laboratory should establish its own normal range for CYFRA 21-1. Based on data obtained by XEMA, the following normal range is recommended (see below). NOTE: the patients that have received murine monoclonal antibodies for radioimaging or immunotherapy develop high titered antimouse antibodies (HAMA). The presence of these antibodies may cause false results in the present assay. Sera from HAMA positive patients should be treated with depleting adsorbents before assaying.

NOTE: values of CYFRA 21-1 concentrations in the tested samples that are below the LoD (0.5 ng/mL) and also exceed the value of the upper Calibrator (50 ng/mL) should be provided in the following form : «the CYFRA 21-1 concentration of tested sample X is «lower than 0.5 ng/mL» or «higher than 50 ng/mL».

Sex, age	Units, ng/mL	
	Lower limit	Upper limit
Healthy donors	-	3.0

13. PERFORMANCE CHARACTERISTICS

13.1. Analytical performance characteristics

13.1.1 Precision of Measurement

Repeatability (Intra assay repeatability) was determined by evaluation the coefficient of variation (CV) for 2 different samples during 1 day in 24 replicates on one series of ELISA kit.

Sample	Concentration, ng/mL	CV, %
1	12.3	6.2
2	25	3.3

Reproducibility (Inter assay reproducibility) was determined by evaluating the coefficients of variation for 2 samples during 5 days in 8-replicate determinations.

Sample	Concentration, ng/mL	CV, %
1	12.27	4.3
2	25.89	5.2

Reproducibility between lots was investigated by testing samples for one day on three lots. Each sample was run in 8 replicates.

Sample	Concentration1, ng/mL	Concentration2, ng/mL	Concentration3, ng/mL	CV, %
1	12.32	12.02	12.81	5.2
2	25.02	25.6	26.0	2.9

13.1.2 Trueness

The trueness of measurement is the degree of closeness of the average value obtained from a large number of measurement results to the true value. The bias of the measurement result (bias of measurements) is the difference between the mathematical expectation of the measurement result and the true value of the measurand. The bias was calculated for each sample and it was determined whether it corresponds to the specified limits of $\pm 10\%$.

13.1.3 Linearity

Linearity was determined using sera samples with known CYFRA 21-1 concentration (low and high) and mixing them with each other and buffer solution in different proportions. According to the measurements, linear range of kit is 3-25 ng/mL $\pm 10\%$.

13.1.4 Analytical sensitivity

Limit of detection (LoD) – the lowest CYFRA 21-1 concentration in the serum or plasma sample that is detected by the CYFRA 21-1 EIA kit is no lower than 0.5 ng/mL.

Limit of quantification (LoQ) – the lowest concentration of the analyte in the sample that is determined quantitatively with the declared trueness for CYFRA 21-1 EIA kit is 3 ng/mL.

13.1.5 Hook Effect

Hook effect is absent for all samples up to reasonably foreseen concentrations 50 ng/mL.

13.1.6 Analytical specificity

For the analysis result is not affected by the presence in the sample of bilirubin in a concentration of up to 0.21 mg/mL and hemoglobin in a concentration of up to 10 mg/mL.

The cross-reactivity of CYFRA 21-1 with other analytes is shown in the table:

Analyte	Cross-reactivity, %
CA 15-3	<0.1
CA 125	<0.1
CA 19-9	<0.1
AFP	<0.1
PSA	<0.1

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LOT _____ DATE _____

	Manufacturer
	<i>In vitro</i> diagnostic medical device
	Catalogue number
 YYYY-MM	Use-by date
	Batch code
	Temperature limit
	Contains sufficient for <n> tests
	Caution
	Consult instructions for use
	Conformity Marking with technical regulations in Ukraine
	Authorized representative in the European Community/European Union
	CE Conformity Marking

**For any issues related to operation of the kit and technical support,
please contact by telefon number**

+38 044 294-69-78

or write to:

qa@xema.com.ua

XEMA LLC
Akademika Yefremova St. 23
03179, Kyiv, Ukraine
tel.:+38 050 422-62-16
tel.:+38 044 294-69-78
E-mail: qa@xema.com.ua
www.xema.in.ua

Instruction for use
A solid-phase enzyme immunoassay kit
for the quantitative determination of
carbohydrate antigen 242
in human serum or plasma

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







K243

CONTENT

1. INTENDED USE	2
2. GENERAL INFORMATION	2
3. TEST PRINCIPLE	3
4. KIT COMPONENTS	4
5. EQUIPMENT AND MATERIAL REQUIRED BUT NOT PROVIDED	5
6. WARNING AND PRECAUTIONS	5
7. SPECIMEN COLLECTION, TRANSPORTATION AND STORAGE OF SAMPLES	6
8. TRANSPORTATION AND STORAGE TERMS OF KIT, WASTE DISPOSAL	6
9. REAGENTS PREPARATION	7
10. ASSAY PROCEDURE	7
11. TEST VALIDITY	9
12. EXPECTED VALUES	9
13. PERFORMANCE CHARACTERISTICS	10
14. REFERENCES	10
SAMPLES IDENTIFICATION PLAN	11

Instruction for use
A solid-phase enzyme immunoassay kit
for the quantitative determination of
carbohydrate antigen 242 in human serum or plasma
CA 242 EIA

1. INTENDED USE

The CA 242 EIA kit is an enzyme immunoassay, intended for the quantitative determination of carbohydrate antigen 242 in human serum or plasma.

The field of application is clinical laboratory diagnostics.

2. GENERAL INFORMATION

The carbohydrate antigen CA 242 is one of the most advanced markers of gastrointestinal cancer. CA 242 is found on cells of colonic mucosa as well as on apical part of cells lining pancreatic ducts.

CA 242 is one of the most important markers used in oncology. For differential diagnostics between pancreatic cancer (PC) and chronic pancreatitis, diagnostic specificity of CA 242 is 1.4 fold higher than that of CA 19-9. In patients with PC, a positive prognostic value of CA 242 determination is higher than that of CA 19-9 at any stage of the disease.

3. TEST PRINCIPLE

The determination of the CA 242 is based on the two-site sandwich enzyme immunoassay principle. On the inner surface of the microplate wells are immobilized specific murine monoclonal antibodies to CA 242/CA 19-9. Second antibodies – murine monoclonal antibodies to human CA 242 conjugated to the horseradish peroxidase is used as enzyme conjugate. The analysis procedure includes two stages of incubation:

- during the first stage CA 242 from the specimen is captured by the antibodies coated onto the microwell surface, as well as horseradish peroxidase-conjugated monoclonal antibodies bind to free epitopes of immobilized CA 242;

- during the second stage, the complexes formed due to the reaction with the chromogen 3,3',5,5'-tetramethylbenzidine are visualized.

After stopping the reaction with a stop solution, the intensity of the color of the microwells is measured. The optical density in the microwell is directly related to the quantity of the measured CA 242 in the serum specimen (plasma). The concentration is determined according to the calibration graph of the dependence of the optical density on the content of CA 242 in the calibration samples.

4. KIT COMPONENTS

Code of component	Symbol	Name	Volume	Qty, pcs.	Description
P223Z	SORB MTP	Microplate	-	1	96-well polystyrene strip microplate coated with murine monoclonal antibodies to CA 242/CA 19-9; ready to use
C243Z	CAL 1	Calibrator C1	2 mL	1	Solution based on tris buffer (pH 7.2-7.4), free of CA 242, with preservative, ready to use (colourless liquid)
C243Z	CAL 2-5	Calibrators	0.5 mL	4	Solutions based on tris buffer (pH 7.2-7.4), containing 15; 50; 100 and 200 U/mL of CA 242, with preservative, ready to use (blue liquids)
Q243Z	CONTROL	Control Serum	0.5 mL	1	Solution based on human serum, containing of known CA 242 content, with preservative, ready to use (colourless liquid)
T243Z	CONJ HRP	Conjugate Solution	14 mL	1	Solution of murine monoclonal antibodies to human CA 242 conjugated to the horseradish peroxidase; ready to use (red liquid)
R055Z	SUBS TMB	Substrate Solution	14 mL	1	Tetramethylbenzidine (TMB) substrate solution; ready to use (colourless liquid)
S008Z	BUF WASH 26X	26x Concentrate Washing Solution	22 mL	1	Buffer solution with detergent, 26x concentrate (colourless liquid)
R050Z	STOP	Stop Solution	14 mL	1	5.0% solution of sulphuric acid; ready to use (colourless liquid)

The kit also includes instruction for use, quality control data sheet and plate sealing tape (2 pcs.)

5. EQUIPMENT AND MATERIAL REQUIRED BUT NOT PROVIDED

- microplate photometer with 450 nm wavelength or 450\620-680 nm;
- dry thermostat for +37°C±2°C;
- automatic plate washer (optional);
- micropipettes with variable volume, range volume 5-1000 µL;
- graduated cylinder of 1000 mL capacity;
- distilled or deionized water;
- timer;
- vortex mixer;
- disposable gloves;
- absorbent paper.

6. WARNING AND PRECAUTIONS

In order to prevent incorrect results, strictly follow the recommended order and duration of the analysis procedure.

6.1. The kit is for *in vitro* diagnostic use only. For professional laboratory use.

6.2. Follow the rules mentioned below during the kit using:

- do not use kit beyond expiry date;
- do not use the kit if its packaging is damaged;
- in order to avoid contamination, use new tips to pipette samples and reagents;
- use only verified equipment;
- close each vial with its own cap, after using the reagent;
- do not use components of other kits or reagents of other manufacturers;
- do not let wells dry after completing the rinsing step; immediately proceed to the next stage;
- avoid bubbles when adding reagents.

ATTENTION! The TMB substrate solution is light sensitive. Avoid prolonged exposure of the component to light.

6.3. Some kit components, such as stop solution, substrate solution, and washing solution, may cause toxic or irritant effects. If they get on the skin or mucosa, the affected area should be washed with plenty of running water.

6.4. All human products, including patient samples, should be considered potentially infectious. Handling and disposal should be in accordance with the procedures defined by an appropriate national biohazard safety guidelines or regulations.

6.5. The Calibrators and Control Serum included in the kit are negative for antibodies to HIV 1,2, hepatitis C virus and HBsAg, but the reagents should be considered as potentially infectious material and handled carefully.

6.6. Specimens must not contain any azide compounds, as they inhibit activity of peroxidase.

6.7. Wear protective gloves, protective clothing, eye protection, face protection.

6.8. Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.

6.9. Safety Data Sheet for this product is available upon request directly from XEMA LLC.

6.10. Serious incidents related to the kit must be reported to the manufacturer, Authorized Representative, and to the Competent Authority of the EU member state(s) where the incident has occurred.

7. SPECIMEN COLLECTION, TRANSPORTATION AND STORAGE OF SAMPLES

7.1. Blood sampling should be carried out from the cubital vein with a disposable needle using a vacuum blood sampling system. Serum or plasma specimens should be clearly labeled and identified. Serum must be separated from the clot as early as possible to avoid hemolysis of red blood cells. If there are any visible particles in the sample, they should be removed by centrifugation at 3000-5000 rpm for 20 minutes at room temperature or by filtration.

Don't use samples with high lipidemia, hemolysis as they may give false test results.

7.2. Specimen should be stored at +2...+8°C up to 3 days. Specimen held for a longer time, should be placed in a freezer at -15°C or below, do not refreeze/thaw samples.

7.3. For the transportation of samples, it is recommended to use triple packaging. The primary package is the labeled tube containing the sample. Secondary packaging is a polyethylene bag that is hermetically closed with a zip-lock. The outer packaging is a heat-insulating container, while the secondary packaging is placed in the outer packaging for transportation in the center of the thermal container. Frozen refrigerants are placed on the bottom, along the side walls of the thermal container, and cover the samples with them.

8. TRANSPORTATION AND STORAGE TERMS OF KIT, WASTE DISPOSAL

Information about the singularity storage conditions, transportation of the kit, and disposal of waste should be taken into account by all persons who participate in these processes.

8.1. Transportation

The CA 242 EIA kit should be transported in the manufacturer's packaging at +2...+8°C. Single transportation at the temperature up to 25°C for 5 days is acceptable.

8.2. Storage

The CA 242 EIA kit should be stored in the manufacturer's packaging at +2...+8°C. Do not freeze.

The kit contains reagents sufficient for 96 determinations including Calibrators and Control Serum.

Once opened test-kit is stable for 2 months when stored properly as intended by manufacturer at 2-8°C.

In case of partial use of the kit, the components should be stored in the following way:

- strips that remain unused must be carefully sealed with the plate sealing tape and stored at +2...+8°C within 2 months;
- Substrate Solution, Stop Solution, and Washing Solution concentrate after opening the vial, can be stored tightly closed at +2...+8°C until the kit's shelf life;
- Conjugate Solution, Calibrators and Control Serum after opening the vial, can be stored tightly closed at +2...+8°C within 2 months;
- diluted Washing Solution can be stored at room temperature (+18...+25°C) for up to 5 days or at +2...+8°C for up to 14 days.

NOTE: Single freezing of Calibrators and Control Serum in aliquots is allowed.

Kits that were stored in violation of the storage condition cannot be used.

8.3. Disposal

Expired kit components, used reagents and materials, as well as residual samples must be inactivated and disposed of in accordance with legal requirements.

9. REAGENTS PREPARATION

9.1. All reagents (including microstrips) and test samples should be allowed to reach room temperature (+18...+25 °C) for at least 30 minutes before use.

9.2. Microplate preparation

Open the package with the microplate and install the required number of strips into the frame. Unused strips must be sealed with plate sealing tape to prevent moisture from affecting the plate's holes and placed back in the bag.

9.3. Washing Solution preparation

Add the contents of the 22 mL Washing Solution concentrate vial to 550 mL of distilled or deionized water and mix thoroughly. In case of partial use of the kit, take the necessary amount of Washing Solution concentrate and dilute it 26 times with distilled or deionized water.

The spending of the components in case of partial use of the kit is given in the table:

Quantity of strips	1	2	3	4	5	6	7	8	9	10	11	12
Volume of the Washing Solution concentrate, mL	1.8	3.6	5.4	7.2	9	10.8	12.6	14.4	16.2	18	19.8	22
Volume of water, mL	45	90	135	180	225	270	315	360	405	450	495	550

9.4. Samples preparation

If suggested analyte concentration in the sample exceeds the 200 U/mL, additionally dilute this sample accordingly, using (Calibrator C1). Use of other buffers or reagents for sample dilution may lead to incorrect measurement.

NOTE: in order to obtain reliable results, we recommend to use several successive dilutions of the blood serum (plasma) sample

10. ASSAY PROCEDURE

- 10.1 Put the desired number of strips into the frame based on the number of test samples in 2 replicates and 12 wells for Calibrators and Control Serum (2 wells for each Calibrator (CAL 1-5) and 2 wells for Control Serum (Q)).
- 10.2 If necessary, dilute the test samples as described in 9.4.
- 10.3 Dispense **25 µL of Calibrators and Control Serum as well as 25 µL of test serum/plasma samples (SAMP)** to the wells of the microplate according to the scheme below. The introduction of Calibrators, Control Serum and test samples should be carried out within 5 minutes to ensure equal incubation time for the first and last samples.

NOTE: during performing several independent series of tests, Calibrators, and Control Serum should be used each time.

Scheme of introduction of samples

	1	2	3	4	5	6	7	8	9	10	11	12
A	CAL1	CAL1	SAMP3	SAMP3	SAMP11	SAMP11						
B	CAL2	CAL2	SAMP4	SAMP4	SAMP12	SAMP12						
C	CAL3	CAL3	SAMP5	SAMP5								
D	CAL4	CAL4	SAMP6	SAMP6								
E	CAL5	CAL5	SAMP7	SAMP7								
F	Q	Q	SAMP8	SAMP8								
G	SAMP1	SAMP1	SAMP9	SAMP9								
H	SAMP2	SAMP2	SAMP10	SAMP10								

- 10.4 Dispense **100 µL of Conjugate Solution** to all wells.
- 10.5 Carefully mix the contents of the microplate in a circular motion on a horizontal surface, cover strips with a plate sealing tape and incubate for **60 minutes at +37°C**.
- 10.6 At the end of the incubation period, remove and discard the plate cover. Aspirate and wash each well 5 times using an automatic washer or an 8-channel dispenser. For each washing, add 300 µL of Washing Solution (see 9.3) to all wells, then remove the liquid by aspiration or decantation. The residual volume of the Washing Solution after each aspiration or decantation should be no more than 5µL. After washing, carefully remove the remaining liquid from the wells on the absorbent paper. For the automatic washer/analyzer, the wash solution volume can be increased to 350 µL.
- 10.9 Add **100 µL of Substrate Solution** to all wells. The introduction of the Substrate Solution into the wells must be carried out within 2-3 minutes. Incubate the microplate in the dark **at room temperature (+18...+25°C) for 15 minutes**.
- 10.10 Add **100 µL of Stop Solution** to all wells in the same order as the Substrate Solution. After adding the Stop Solution, the contents of the wells turn yellow.
- 10.11 Read the optical density (OD) of the wells at 450nm and reference light filters 620–680 nm using a microplate photometer within 5 minutes of adding the Stop Solution. Set photometer blank on CAL1.
- 10.13 Plot a calibration curve in linear coordinates: (x) is the CA 242 concentration in the calibrators U/mL, (y) – OD versus CA 242 concentration (OD 450 nm / 620–680 nm). Manual or computerized data reduction is applicable at this stage. Point-by-point or linear data reduction is recommended due to non-linear shape of curve.
- 10.14 Determine the corresponding concentration of CA 242 in tested samples from the calibration curve. In the case of preliminary dilution of the test sample (see 9.4), the obtained result should be multiplied by the dilution factor.

11. TEST VALIDITY

The test run shall be considered valid if the OD of CAL1 is above 0.15, and the values of the Control Serum fall into the required range (see Quality control Data Sheet).

12. EXPECTED VALUES

12.1. Therapeutical consequences should not be based on results of IVD methods alone - all available clinical and laboratory findings should be used by a physician to elaborate therapeutically measures. Each laboratory should establish its own normal range for CA 242. Based on data obtained by XEMA LLC, the following normal range is recommended (see below).

NOTE: values of CA 242 concentrations in the tested samples that are below the LoD (0.5 U/mL) and also exceed the value of the upper Calibrator (200 U/mL) should be provided in the following form : «the CA 242 concentration of tested sample X is «lower than 0.5 U/mL» or «higher than 200 U/mL».

Sex, age	Units, U/mL	
	Lower limit	Upper limit
Males	-	20
Females	-	20

13. PERFORMANCE CHARACTERISTICS

13.1. Analytical performance characteristics

13.1.1 Precision of Measurement

Repeatability (Intra assay repeatability) was determined by evaluation the coefficient of variation (CV) for 2 different samples during 1 day in 24 replicates on one series of ELISA kit.

Sample	Concentration, U/mL	CV, %
1	10.12	3.2
2	53.64	2.8

Reproducibility (Inter assay reproducibility) was determined by evaluating the coefficients of variation for 2 samples during 5 days in 8-replicate determinations.

Sample	Concentration, U/mL	CV, %
1	10.27	7.0
2	53.87	6.1

Reproducibility between lots was investigated by testing samples for one day on three lots. Each sample was run in 8 replicates.

Sample	Concentration1, U/mL	Concentration2, U/mL	Concentration3, U/mL	CV, %
1	10.32	10.02	10.81	3.8
2	53.71	53.56	54.32	0.6

13.1.2 Trueness

The trueness of measurement is the degree of closeness of the average value obtained from a large number of measurement results to the true value. The bias of the measurement result (bias of measurements) is the difference between the mathematical expectation of the measurement result and the true value of the mezhurand. The bias was calculated for each sample and it was determined that it corresponds to the specified limits of $\pm 10\%$.

13.1.3 Linearity

Linearity was determined using sera samples with known CA 242 concentration (low and high) and mixing them with each other and buffer solution in different proportions. According to the measurements, linear range of kit is 15-100 U/mL \pm 10%.

13.1.4 Analytical sensitivity

Limit of detection (LoD) – the lowest CA 242 concentration in the serum or plasma sample that is detected by the CA 242 EIA kit is no lower than 0.5 U/mL.

Limit of quantification (LoQ) – the lowest concentration of the analyte in the sample that is determined quantitatively with the declared trueness for CA 242 EIA kit is 15 U/mL.

13.1.5 Hook Effect

Hook effect is absent for all samples up to reasonably foreseen concentrations 200 U/mL.

13.1.6 Analytical specificity

For the analysis result is not affected by the presence in the sample of bilirubin in a concentration of up to 0.21 mg/mL and hemoglobin in a concentration of up to 10 mg/mL.

The cross-reactivity of CA 242 with other analytes is shown in the table:

Analyte	Cross-reactivity, %
CEA	<0.1
CA 15-3	<0.1

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LOT _____ DATE _____

	Manufacturer
	<i>In vitro</i> diagnostic medical device
	Catalogue number
 YYYY-MM	Use-by date
	Batch code
	Temperature limit
	Contains sufficient for <n> tests
	Caution
	Consult instructions for use
	Conformity Marking with technical regulations in Ukraine
	Authorized representative in the European Community/European Union
	CE Conformity Marking

**For any issues related to operation of the kit and technical support,
please contact by telefon number**

+38 044 294-69-78

or write to:

qa@xema.com.ua

XEMA LLC
Akademika Yefremova St. 23
03179, Kyiv, Ukraine
tel.:+38 044 422-62-16
tel.:+38 044 294-69-78
E-mail: qa@xema.com.ua
www.xema.com.ua

Instruction for use
A solid-phase enzyme immunoassay kit
for the quantitative determination of
carbohydrate antigen 72-4
in human serum or plasma

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








ASSAY PROCEDURE

During performing several independent series of tests, Calibrators, and Control Serum should be used each time.

CONTENT

1. INTENDED USE	2
2. GENERAL INFORMATION	2
3. TEST PRINCIPLE	2
4. KIT COMPONENTS	3
5. EQUIPMENT AND MATERIAL REQUIRED BUT NOT PROVIDED	4
6. WARNING AND PRECAUTIONS	4
7. SPECIMEN COLLECTION, TRANSPORTATION AND STORAGE OF SAMPLES	5
8. TRANSPORTATION AND STORAGE TERMS OF KIT, WASTE DISPOSAL	5
9. REAGENTS PREPARATION	6
10. ASSAY PROCEDURE	6
11. TEST VALIDITY	7
12. EXPECTED VALUES	8
13. PERFORMANCE CHARACTERISTICS	8
14. REFERENCES	9
SAMPLES IDENTIFICATION PLAN	10

Instruction for use
A solid-phase enzyme immunoassay kit
for the quantitative determination of
carbohydrate antigen 72-4
in human serum or plasma
CA 72-4 EIA

1. INTENDED USE

The CA 72-4 EIA kit is an enzyme immunoassay, intended for the quantitative determination of carbohydrate antigen 72-4 in human serum or plasma.

Determination of CA 72-4 antigen concentration in serum (plasma) is used as an auxiliary method of early diagnosis, monitoring the effectiveness of therapy in malignant tumors of glandular tissue, such as gastric carcinoma, colon or ovarian cancer, for all population groups.

The field of application is clinical laboratory diagnostics.

2. GENERAL INFORMATION

CA 72-4, or a carbohydrate antigen 72-4, is a high MM (230-1000 kD) antigen (epitope) associated to gastric and ovarian cancer as well as some other malignancies and not expressed in noticeable quantities in tissues of healthy adult individuals.

Quantitative determination of CA 72-4 in serum or plasma is helpful (particularly, in combination with CA 19-9 – see XEMA LLC, Cat.# K223) for monitoring of gastric cancer and its therapy, while combined determination of CA 72-4 and CA 125 (see XEMA LLC, Cat.# K222) is used for monitoring of ovarian cancer.

Elevated levels of CA 72-4 are often seen in adenocarcinomas of the gastro-intestinal tract, ovaries (mucinous type) and lungs. Besides, raised CA 72-4 is sometimes also seen in patients with benign pathology (chronic inflammation, cysts, fibrosis). That is why, results of CA 72-4 determination should always be interpreted in conjunction with other laboratory and clinical data.

Functional purpose. Determination of the concentration of CA 72-4 antigen in serum (plasma) is used as an auxiliary method for early diagnosis, monitoring the effectiveness of therapy for malignant glandular tumours, such as gastric carcinoma, colon or ovarian cancer, for all population groups.

3. TEST PRINCIPLE

The determination of carbohydrate antigen 72-4 is based on the two-site sandwich enzyme immunoassay principle. On the inner surface of the microplate wells are immobilized specific murine monoclonal antibodies to human CA 72-4. Second antibodies – murine monoclonal antibodies to human CA 72-4 conjugated to the horseradish peroxidase is used as enzyme conjugate. The analysis procedure includes two stages of incubation:

- during the first stage CA 72-4 from the specimen is captured by the antibodies coated onto the microwell surface, as well as horseradish peroxidase-conjugated monoclonal antibodies bind to free epitopes of immobilized CA 72-4;

- during the second stage, the complexes formed due to the reaction with the chromogen 3,3',5,5'-tetramethylbenzidine are visualized.

After stopping the reaction with a stop solution, the intensity of the color of the microwells is measured. The optical density in the microwell is directly related to the quantity of the measured CA 72-4 in the serum specimen (plasma). The concentration is determined according to the calibration graph of the dependence of the optical density on the content of CA 72-4 in the calibration samples.

4. KIT COMPONENTS

Code of component	Symbol	Name	Volume	Qty, pcs.	Description
P244Z	SORB MTP	Microplate	-	1	96-well polystyrene strip microplate coated with murine monoclonal antibodies to human CA 72-4, ready to use
C244Z	CAL 1	Calibrator C1	0.5 mL	1	Solution based on human serum, free of CA 72-4, with preservative, ready to use (yellow liquid)
C244Z	CAL 2-5	Calibrators	0.5 mL	4	Solutions based on human serum, containing 5; 15; 50 and 200 U/mL of CA 72-4, with preservative, ready to use (blue liquids)
Q244Z	CONTROL	Control Serum	0.5 mL	1	Solution based on human serum, containing of known CA 72-4 content, with preservative, ready to use (yellow liquid)
T244XZ	CONJ 11X	Conjugate Concentrate	1.2 mL	1	Solution of murine monoclonal antibodies to human CA 72-4 conjugated to the horseradish peroxidase, 11x concentrate (green liquid)
ST244Z	DIL CONJ	Conjugate Dilution Buffer	12 mL	1	Buffer solution with detergent, ready to use (blue liquid)
R055Z	SUBS TMB	Substrate Solution	12 mL	1	Tetramethylbenzidine (TMB) substrate solution, ready to use (colourless liquid)
S008Z	BUF WASH 26X	26x Concentrate Washing Solution	22 mL	1	Buffer solution with detergent, 26x concentrate (colourless liquid)
R050Z	STOP	Stop Solution	12 mL	1	5.0% solution of sulphuric acid, ready to use (colourless liquid)
The kit also includes instruction for use, quality control data sheet and plate sealing tape (1 pcs.)					

5. EQUIPMENT AND MATERIAL REQUIRED BUT NOT PROVIDED

- microplate photometer with 450\620-680 nm wavelength;
- dry thermostat for 37°C±1°C;
- automatic plate washer (optional);
- micropipettes with variable volume, range volume 5-1000 µL;
- graduated cylinder of 1000 mL capacity;
- distilled or deionized water;
- timer;
- vortex mixer;
- disposable gloves;
- absorbent paper.

6. WARNING AND PRECAUTIONS

In order to prevent incorrect results, strictly follow the recommended order and duration of the analysis procedure.

6.1. The kit is for *in vitro* diagnostic use only. For professional laboratory use.

6.2. Follow the rules mentioned below during the kit using:

- do not use kit beyond expire date;
- do not use the kit if its packaging is damaged;
- in order to avoid contamination, use new tips to pipette samples and reagents;
- use only verified equipment;
- close each vial with its own cap, after using the reagent;
- do not use components of other kits or reagents of other manufacturers;
- do not let wells dry after completing the rinsing step; immediately proceed to the next stage;
- avoid bubbles when adding reagents.

ATTENTION! The TMB substrate solution is light sensitive. Avoid prolonged exposure of the component to light.

6.3. Some kit components, such as stop solution, substrate solution, and washing solution, may cause toxic or irritant effects. If they get on the skin or mucosa, the affected area should be washed with plenty of running water.

6.4. All human products, including patient samples, should be considered potentially infectious. Handling and disposal should be in accordance with the procedures defined by an appropriate national biohazard safety guidelines or regulations.

6.5. The Calibrators and Control Serum included in the kit are negative for antibodies to HIV 1,2, hepatitis C virus and HBsAg, but the reagents should be considered as potentially infectious material and handled carefully.

6.6. Specimens must not contain any azide compounds, as they inhibit activity of peroxidase.

6.7. Wear protective gloves, protective clothing, eye protection, face protection.

6.8. Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.

6.9. Safety Data Sheet for this product is available upon request directly from XEMA LLC.

6.10. Serious incidents related to the kit must be reported to the manufacturer, Authorized Representative, and to the Competent Authority of the EU member state(s) where the incident has occurred.

7. SPECIMEN COLLECTION, TRANSPORTATION AND STORAGE OF SAMPLES

7.1. Blood sampling should be carried out from the cubital vein with a disposable needle using a vacuum blood sampling system. Serum or plasma specimens should be clearly labeled and identified. Serum must be separated from the clot as early as possible to avoid hemolysis of red blood cells. If there are any visible particles in the sample, they should be removed by centrifugation at 3000-5000 rpm for 20 minutes at room temperature or by filtration.

Don't use samples with high lipidemia, hemolysis as they may give false test results.

7.2. Specimen should be stored at +2...+8°C up to 3 days. Specimen held for a longer time, should be placed in a freezer at -15°C or below; do not refreeze/thaw samples.

7.3. For the transportation of samples, it is recommended to use triple packaging. The primary package is the labeled tube containing the sample. Secondary packaging is a polyethylene bag that is hermetically closed with a zip-lock. The outer packaging is a heat-insulating container, while the secondary packaging is placed in the outer packaging for transportation in the center of the thermal container. Frozen refrigerants are placed on the bottom, along the side walls of the thermal container, and cover the samples with them.

8. TRANSPORTATION AND STORAGE TERMS OF KIT, WASTE DISPOSAL

Information about the singularity storage conditions, transportation of the kit, and disposal of waste should be taken into account by all persons who participate in these processes.

8.1. Transportation

The CA 72-4 EIA kit should be transported in the manufacturer's packaging at +2...+8°C. Single transportation at the temperature up to 25°C for 5 days is acceptable.

8.2. Storage

The CA 72-4 EIA kit should be stored in the manufacturer's packaging at +2...+8°C. Do not freeze.

The kit contains reagents sufficient for 96 determinations including Calibrators and Control Serum.

Once opened test-kit is stable for 2 months when stored properly as intended by manufacturer at 2-8°C.

In case of partial use of the kit, the components should be stored in the following way:

- the remaining strips should be immediately resealed in the bag along with the silica gel, closed with the zip-lock, and stored at +2...+8°C within 2 months;
- Substrate Solution, Stop Solution, and Washing Solution concentrate after opening the vial, can be stored tightly closed at +2...+8°C until the kit's shelf life;
- Conjugate Concentrate, Conjugate Dilution Buffer, Calibrators and Control Serum after opening the vial, can be stored tightly closed at +2...+8°C within 2 months;
- diluted Washing Solution can be stored at room temperature (+18...+25°C) for up to 5 days or at +2...+8°C for up to 14 days.

Kits that were stored in violation of the storage condition cannot be used.

8.3. Disposal

Expired kit components, used reagents and materials, as well as residual samples should be inactivated and disposed of in accordance with legal requirements.

9. REAGENTS PREPARATION

9.1. All reagents (including microstrips) and test samples should be allowed to reach room temperature (+18...+25 °C) for at least 30 minutes before use.

9.2. Microplate preparation

Open the package with the microplate and install the required number of strips into the frame. The remaining strips should be immediately resealed in the bag along with the silica gel and closed with the zip-lock to prevent moisture from affecting the plate's strips.

9.3. Washing Solution preparation

Add the contents of the 22 mL Washing Solution concentrate vial to 550 mL of distilled or deionized water and mix thoroughly. In case of partial use of the kit, take the necessary amount of Washing Solution concentrate and dilute it 26 times with distilled or deionized water.

9.4. Working conjugate solution preparation

Prepare in a different container a working conjugate solution by 11 dilutions of Conjugate Concentrate in Conjugate Dilution Buffer (eg, 90 µL of concentrate + 900 µL of Conjugate Dilution Buffer). In the case of partial use of the kit, take the necessary amount of Conjugate Concentrate and dilute it 11 times with Conjugate Dilution Buffer, since the working conjugate solution in a diluted form is not stored for a long time.

The spending of the components in case of partial use of the kit is given in the table:

Quantity of strips	1	2	3	4	5	6	7	8	9	10	11	12
Volume of the Washing Solution concentrate, mL	1.8	3.6	5.4	7.2	9	10.8	12.6	14.4	16.2	18	19.8	22
Volume of water, mL	45	90	135	180	225	270	315	360	405	450	495	550
Volume of Conjugate Concentrate, mL	0.09	0.18	0.27	0.36	0.45	0.54	0.63	0.72	0.81	0.9	0.99	1.08
Volume of Conjugate Dilution Buffer, mL	0.9	1.8	2.7	3.6	4.5	5.4	6.3	7.2	8.1	9	9.9	10.8

10. ASSAY PROCEDURE

- Put the desired number of strips into the frame based on the number of test samples in 2 replicates and 12 wells for Calibrators and Control Serum (2 wells for each Calibrator (CAL 1-5) and 2 wells for Control Serum (Q)).
- Dispense **20 µL of Calibrators and Control Serum as well as 20 µL of test serum/plasma samples (SAMP)** to the wells of the microplate according to the scheme below. The introduction of Calibrators, Control Serum and test samples should be carried out within 5 minutes to ensure equal incubation time for the first and last samples.

NOTE: during performing several independent series of tests, Calibrators, and Control Serum should be used each time.

Scheme of introduction of samples

	1	2	3	4	5	6	7	8	9	10	11	12
A	CAL1	CAL1	SAMP3	SAMP3	SAMP11	SAMP11						
B	CAL2	CAL2	SAMP4	SAMP4	SAMP12	SAMP12						
C	CAL3	CAL3	SAMP5	SAMP5								
D	CAL4	CAL4	SAMP6	SAMP6								
E	CAL5	CAL5	SAMP7	SAMP7								
F	Q	Q	SAMP8	SAMP8								
G	SAMP1	SAMP1	SAMP9	SAMP9								
H	SAMP2	SAMP2	SAMP10	SAMP10								

- 10.3. Dispense **100 µL of Working conjugate solution** to all wells (see 9.4).
- 10.4. Carefully mix the contents of the microplate in a circular motion on a horizontal surface, cover strips with a plate sealing tape and incubate for **60 minutes at +37°C**.
- 10.5. At the end of the incubation period, remove and discard the plate cover. Aspirate and wash each well **5 times** using an automatic washer or an 8-channel dispenser. For each washing, add 300 µL of Washing Solution (see 9.3) to all wells, then remove the liquid by aspiration or decantation. The residual volume of the Washing Solution after each aspiration or decantation should be no more than 5µL. After washing, carefully remove the remaining liquid from the wells on the absorbent paper. For the automatic washer/analyzer, the Washing Solution volume can be increased to 350 µL.
- 10.6. Add **100 µL of Substrate Solution** to all wells. The introduction of the Substrate Solution into the wells must be carried out within 2-3 minutes. Incubate the microplate in the dark **at room temperature (+18...+25°C) for 30 minutes**.
The incubation time can be varied depending on the intensity of the blue colour development.
- 10.7. Add **100 µL of Stop Solution** to all wells in the same order as the Substrate Solution. After adding the Stop Solution, the contents of the wells turn yellow.
- 10.8. Read the optical density (OD) of the wells at 450 nm and reference light filters 620–680 nm using a microplate photometer within 5 minutes of adding the stop solution. Set photometer blank on CAL1.
- 10.9. Plot a calibration curve in linear coordinates: (x) is the concentration of CA 72-4 in the Calibrators U/mL, (y) – OD versus concentration of CA 72-4 (OD 450 nm / 620–680 nm). Manual or computerized data reduction is applicable at this stage. For the algorithm calculation (approximation) of the calibration curve, using the interval (segment-linear, point-to-point) method is recommended.
- 10.10. Determine the corresponding concentration of CA 72-4 in tested samples from the calibration curve.

11. TEST VALIDITY

The test run shall be considered valid if the OD of CAL1 is above 0.15, and the values of the Control Serum fall into the required range (see Quality control Data Sheet).

12. EXPECTED VALUES

Therapeutical consequences should not be based on results of IVD methods alone – all available clinical and laboratory findings should be used by a physician to elaborate therapeutically measures. Each laboratory should establish its own normal range for CA 74-2. Based on data obtained by XEMA, the following normal range is recommended (see below). NOTE: the patients that have received murine monoclonal antibodies for radioimaging or immunotherapy develop high titered antimouse antibodies (HAMA). The presence of these antibodies may cause false results in the present assay. Sera from HAMA positive patients should be treated with depleting adsorbents before assaying.

NOTE: values of CA 74-2 concentrations in the tested samples that are below the LoD (0.3 U/mL) and also exceed the value of the upper calibrator (200 U/mL) should be provided in the following form : «the CA 74-2 concentration of tested sample X is «lower than 0.3 U/mL» or «higher than 200 U/mL».

Sex, age	Units, U/mL	
	Lower limit	Upper limit
Healthy donors	-	6.0

13. PERFORMANCE CHARACTERISTICS

13.1. Analytical performance characteristics

13.1.1. Precision of Measurement

Reproducibility. The coefficient of variation of determining the content of CA 72-4 in the same sample of blood serum (plasma) using the kit CA 72-4 EIA does not exceed 10%.

13.1.2. Trueness

The trueness of measurement is the degree of closeness of the average value obtained from a large number of measurement results to the true value. The bias of the measurement result (bias of measurements) is the difference between the mathematical expectation of the measurement result and the true value of the measurand. The bias was calculated for each sample and it was determined whether it corresponds to the specified limits of $\pm 10\%$.

13.1.3. Linearity

Linearity was determined using sera samples with known CA 72-4 concentration (low and high) and mixing them with each other and buffer solution in different proportions. According to the measurements, linear range of kit is 5–200 U/mL $\pm 10\%$.

13.1.4. Analytical sensitivity

Limit of detection (LoD) – the lowest CA 72-4 concentration in the serum or plasma sample that is detected by the CA 72-4 EIA kit is no lower than 0.3 U/mL.

Limit of quantification (LoQ) – the lowest concentration of the analyte in the sample that is determined quantitatively with the declared trueness for CA 72-4 EIA kit is 5 U/mL.

13.1.5. Hook Effect

Hook effect is absent for all samples up to reasonably foreseen concentrations 20000 U/mL.

13.1.6. Analytical specificity

For the analysis result is not affected by the presence in the sample of bilirubin in a concentration of up to 0.21 mg/mL and hemoglobin in a concentration of up to 10 mg/mL.

The cross-reactivity of CA 74-2 with other analytes is shown in the table:

Analyte	Cross-reactivity, %
CEA	<0.1
CA 125	<0.1
CA 19-9	<0.1

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LOT _____

DATE _____

	Manufacturer
	<i>In vitro</i> diagnostic medical device
	Catalogue number
 YYYY-MM	Use-by date
	Batch code
	Temperature limit
	Contains sufficient for <n> tests
	Caution
	Consult instructions for use
	Conformity Marking with technical regulations in Ukraine

**For any issues related to operation of the kit and technical support,
please contact by telefon number**

+38 044 294-69-78

or write to:

qa@xema.com.ua

XEMA LLC
Akademika Yefremova St. 23
03179, Kyiv, Ukraine
tel.:+38 044 422-62-16
tel.:+38 044 294-69-78
E-mail: qa@xema.com.ua
www.xema.com.ua

**INSTRUCTION FOR USE
A SOLID-PHASE ENZYME IMMUNOASSAY KIT
FOR THE QUANTITATIVE DETERMINATION
OF C-PEPTIDE IN HUMAN SERUM OR PLASMA**

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

ASSAY PROCEDURE


Dispensing of Calibrators, Control Serum and test samples


Conjugate Solution


Incubation 1


Washing 5 times


OD measuring, calculation of results


Stop Solution

Incubation 2

Substrate Solution

K267C

CONTENT

1. INTENDED USE	2
2. GENERAL INFORMATION	2
3. TEST PRINCIPLE	2
4. KIT COMPONENTS	3
5. EQUIPMENT AND MATERIAL REQUIRED BUT NOT PROVIDED	4
6. WARNING AND PRECAUTIONS	4
7. SPECIMEN COLLECTION, TRANSPORTATION AND STORAGE OF SAMPLES	5
8. TRANSPORTATION AND STORAGE TERMS OF KIT, WASTE DISPOSAL	5
9. REAGENTS PREPARATION	6
10. ASSAY PROCEDURE	6
11. TEST VALIDITY	7
12. EXPECTED VALUES	8
13. PERFORMANCE CHARACTERISTICS	8
14. REFERENCES	9
SAMPLES IDENTIFICATION PLAN	10

**INSTRUCTION FOR USE
A SOLID-PHASE ENZYME IMMUNOASSAY KIT
FOR THE QUANTITATIVE DETERMINATION OF C-PEPTIDE
IN HUMAN SERUM OR PLASMA**

C-peptide EIA

1. INTENDED USE

A solid-phase enzyme immunoassay for the quantitative determination of C-peptide in blood serum or plasma.

The field of application is clinical laboratory diagnostics.

2. GENERAL INFORMATION

C-peptide is a component of the endocrine secretion of the pancreas. This protein is necessary for insulin synthesis in pancreatic cells, a multi-step process during which inactive proinsulin is broken down to release active insulin. C-peptide and insulin are secreted in equimolar amounts, so determining the level of C-peptide allows to estimate insulin secretion. It should be noted that although the number of C-peptide and insulin molecules formed during secretion into the bloodstream is the same, the molar concentration of C-peptide in the blood is approximately 5 times higher than the molar concentration of insulin, which is most likely due to the different rates of excretion of these substances from the bloodstream.

C-peptide determination has a number of advantages over insulin determination: the half-life of C-peptide in the bloodstream is longer than that of insulin, so the level of C-peptide is a more stable indicator than the concentration of insulin. In immunoassays, C-peptide does not cross-react with insulin, which makes it possible to measure C-peptide to assess insulin secretion even when taking exogenous insulin and in the presence of autoantibodies to insulin, which is important in the examination of patients with diabetes. The level of C-peptide changes in accordance with fluctuations in insulin levels. The ratio of these indicators may change in case of liver and kidney disease, as insulin is metabolised mainly by the liver, while C-peptide is metabolised and excreted by the kidneys. In this regard, the determination of this indicator can be useful for the correct interpretation of changes in blood insulin levels in case of liver function disorders.

3. TEST PRINCIPLE

The determination of C-peptide is based on the two-site sandwich enzyme immunoassay principle. On the inner surface of the microplate wells are immobilized specific murine monoclonal antibodies to human C-peptide. Second antibodies – murine monoclonal antibodies to human C-peptide conjugated to the horseradish peroxidase is used as enzyme conjugate. The analysis procedure includes two stages of incubation:

- during the first stage C-peptide from the specimen is captured by the antibodies coated onto the microwell surface, as well as horseradish peroxidase-conjugated monoclonal antibodies bind to free epitopes of immobilized C-peptide;
- during the second stage, the complexes formed due to the reaction with the chromogen 3,3',5,5'-tetramethylbenzidine are visualized.

After stopping the reaction with a stop solution, the intensity of the color of the microwells is measured. The optical density in the microwell is directly related to the quantity of the measured C-peptide in the serum specimen (plasma). The concentration is determined according to the calibration graph of the dependence of the optical density on the content of C-peptide in the calibration samples.

4. KIT COMPONENTS

Code of component	Symbol	Name	Volume	Qty, pcs.	Description
P267CZ	SORB MTP	Microplate	-	1	96-well polystyrene strip microplate coated with murine monoclonal antibodies to human C-peptide, ready to use
C267CZ	CAL 1	Calibrator C1	0.5 mL	1	Solution based on tris buffer (pH 7.2-7.4), free of human C-peptide, lyophilized
C267CZ	CAL 2-6	Calibrators	0.5 mL	5	Solutions based on tris buffer (pH 7.2-7.4), containing 0.3; 0.7; 3; 7 and 30 ng/mL of human C-peptide, contains blue dye, lyophilized Note: Concentrations of C-peptide in Calibrators may differ from the specified values, the exact values are indicated on the component labels
Q267CZ	CONTROL	Control Serum	0.5 mL	1	Solution based on human serum, containing of known C-peptide content, with preservative, lyophilized
T267CZ	CONJ HRP	Conjugate Solution	12 mL	1	Solution of murine monoclonal antibodies to human C-peptide conjugated to the horseradish peroxidase, ready to use (blue liquid)
R055Z	SUBS TMB	Substrate Solution	12 mL	1	Tetramethylbenzidine (TMB) substrate solution, ready to use (colourless liquid)
S008Z	BUF WASH 26X	26x Concentrate Washing Solution	22 mL	1	Buffer solution with detergent, 26x concentrate (colourless liquid)
R050Z	STOP	Stop Solution	12 mL	1	5.0% solution of sulphuric acid, ready to use (colourless liquid)
The kit also includes instruction for use, quality control data sheet and plate sealing tape (2 pcs.)					

5. EQUIPMENT AND MATERIAL REQUIRED BUT NOT PROVIDED

- microplate photometer with 450 nm and 620-680 nm wavelength;
- automatic plate washer (optional);
- micropipettes with variable volume, range volume 5-1000 µL;
- graduated cylinder of 1000 mL capacity;
- distilled or deionized water;
- timer;
- disposable gloves;
- absorbent paper.

6. WARNING AND PRECAUTIONS

In order to prevent incorrect results, strictly follow the recommended order and duration of the analysis procedure.

6.1. The kit is for *in vitro* diagnostic use only. For professional laboratory use.

6.2. Follow the rules mentioned below during the kit using:

- do not use kit beyond expire date;
- do not use the kit if its packaging is damaged;
- in order to avoid contamination, use new tips to pipette samples and reagents;
- use only verified equipment;
- close each vial with its own cap, after using the reagent;
- do not use components of other kits or reagents of other manufacturers;
- do not let wells dry after completing the rinsing step, immediately proceed to the next stage;
- avoid bubbles when adding reagents.

ATTENTION! The TMB substrate solution is light sensitive. Avoid prolonged exposure of the component to light.

6.3. Some kit components, such as stop solution, substrate solution, and washing solution, may cause toxic or irritant effects. If they get on the skin or mucosa, the affected area should be washed with plenty of running water.

6.4. All human products, including patient samples, should be considered potentially infectious. Handling and disposal should be in accordance with the procedures defined by an appropriate national biohazard safety guidelines or regulations.

6.5. The Calibrators and Control Serum included in the kit are negative for antibodies to HIV 1,2, hepatitis C virus and HBsAg, but the reagents should be considered as potentially infectious material and handled carefully.

6.6. Specimens must not contain any azide compounds, as they inhibit activity of peroxidase.

6.7. Wear protective gloves, protective clothing, eye protection, face protection.

6.8. Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.

6.9. Safety Data Sheet for this product is available upon request directly from XEMA LLC.

6.10. Serious incidents related to the kit must be reported to the manufacturer, Authorized Representative, and to the Competent Authority of the EU member state(s) where the incident has occurred.

7. SPECIMEN COLLECTION, TRANSPORTATION AND STORAGE OF SAMPLES

7.1. Blood sampling should be carried out from the cubital vein with a disposable needle using a vacuum blood sampling system. Serum or plasma specimens should be clearly labeled and identified. Serum must be separated from the clot as early as possible to avoid hemolysis of red blood cells. If there are any visible particles in the sample, they should be removed by centrifugation at 3000-5000 rpm for 20 minutes at room temperature or by filtration.

Don't use samples with high lipidemia, hemolysis as they may give false test results.

7.2. Specimen should be stored at +2...+8°C up to 3 days. Specimen held for a longer time, should be placed in a freezer at -15°C or below, do not refreeze/thaw samples.

7.3. For the transportation of samples, it is recommended to use triple packaging. The primary package is the labeled tube containing the sample. Secondary packaging is a polyethylene bag that is hermetically closed with a zip-lock. The outer packaging is a heat-insulating container, while the secondary packaging is placed in the outer packaging for transportation in the center of the thermal container. Frozen refrigerants are placed on the bottom, along the side walls of the thermal container, and cover the samples with them.

8. TRANSPORTATION AND STORAGE TERMS OF KIT, WASTE DISPOSAL

Information about the singularity storage conditions, transportation of the kit, and disposal of waste should be taken into account by all persons who participate in these processes.

8.1. Transportation

The C-peptide EIA kit should be transported in the manufacturer's packaging at +2...+8°C. Single transportation at the temperature up to 25°C for 5 days is acceptable.

8.2. Storage

The C-peptide EIA kit should be stored in the manufacturer's packaging at +2...+8°C. Do not freeze.

The kit contains reagents sufficient for 96 determinations including Calibrators and Control Serum.

Once opened test-kit is stable for 2 months when stored properly as intended by manufacturer at 2-8°C.

In case of partial use of the kit, the components should be stored in the following way:

- strips that remain unused must be carefully sealed with the plate sealing tape and stored at +2...+8°C within 2 months;
- Substrate Solution, Stop Solution and Washing Solution concentrate after opening the vial, can be stored tightly closed at +2...+8°C until the kit's shelf life;
- Conjugate Solution after opening the vial, can be stored tightly closed at +2...+8°C within 2 months;
- Calibrators and Control Serum after dissolving should be stored frozen in aliquots below **-15°C**.

NOTE: only one freezing/thawing cycle of Calibrators and Control Serum is allowed.

Kits that were stored in violation of the storage condition cannot be used.

8.3. Disposal

Expired kit components, used reagents and materials, as well as residual samples must be inactivated and disposed of in accordance with legal requirements.

9. REAGENTS PREPARATION

9.1. All reagents (including microstrips) and test samples should be allowed to reach room temperature (+18...+25 °C) for at least 30 minutes before use.

9.2. Microplate preparation

Open the package with the microplate and install the required number of strips into the frame. Unused strips must be sealed with plate sealing tape to prevent moisture from affecting the plate's holes and placed back in the bag.

9.3. Washing Solution preparation

Add the contents of the 22 mL Washing Solution concentrate vial to 550 mL of distilled or deionized water and mix thoroughly. In case of partial use of the kit, take the necessary amount of Washing Solution concentrate and dilute it 26 times with distilled or deionized water.

The spending of the components in case of partial use of the kit is given in the table:

Quantity of strips	1	2	3	4	5	6	7	8	9	10	11	12
Volume of the Washing Solution concentrate, mL	1.8	3.6	5.4	7.2	9	10.8	12.6	14.4	16.2	18	19.8	22
Volume of water, mL	45	90	135	180	225	270	315	360	405	450	495	550

9.4. Calibrators and Control Serum preparation

Before first use of the kit dissolve the Calibrators and Control Serum: add 0.5 mL deionized water to each vial and mix thoroughly avoiding foam formation. Liquid Calibrators and Control Serum should be assayed within **72 hours**. For next assays Liquid Calibrators and Control Serum should be aliquoted and stored frozen below -15°C **immediately**.

10. ASSAY PROCEDURE

- 10.1. Put the desired number of strips into the frame based on the number of test samples in 2 replicates and 14 wells for Calibrators and Control Serum (2 wells for each Calibrator (CAL 1-6) and 2 wells for Control Serum (Q)).
- 10.2. Prepare Calibrators and Control Serum as described in 9.4.
- 10.3. Dispense **25 µL of Calibrators and Control Serum as well as 25 µL of test serum/plasma samples** (SAMP) to the wells of the microplate according to the scheme below. The introduction of Calibrators, Control Serum and test samples should be carried out within 5 minutes to ensure equal incubation time for the first and last samples.

NOTE: during performing several independent series of tests, Calibrators and Control Serum should be used each time.

Scheme of introduction of samples

	1	2	3	4	5	6	7	8	9	10	11	12
A	CAL1	CAL1	SAMP2	SAMP2	SAMP10	SAMP10						
B	CAL2	CAL2	SAMP3	SAMP3	SAMP11	SAMP11						
C	CAL3	CAL3	SAMP4	SAMP4	SAMP12	SAMP12						
D	CAL4	CAL4	SAMP5	SAMP5								
E	CAL5	CAL5	SAMP6	SAMP6								
F	CAL6	CAL6	SAMP7	SAMP7								
G	Q	Q	SAMP8	SAMP8								
H	SAMP1	SAMP1	SAMP9	SAMP9								

- 10.4. Dispense **100 µL of Conjugate Solution** to all wells.
- 10.5. Carefully mix the contents of the microplate in a circular motion on a horizontal surface, cover strips with a plate sealing tape and incubate for **60 minutes at room temperature (+18...+25°C)**.
- 10.6. At the end of the incubation period, remove and discard the plate cover. Aspirate and wash each well **5 times** using an automatic washer or an 8-channel dispenser. For each washing, add 300 µL of Washing Solution (see 9.3) to all wells, then remove the liquid by aspiration or decantation. The residual volume of the Washing Solution after each aspiration or decantation should be no more than 5 µL. After washing, carefully remove the remaining liquid from the wells on the absorbent paper. For the automatic washer/analyzer, the wash solution volume can be increased to 350 µL.
- 10.7. Add **100 µL of Substrate Solution** to all wells. The introduction of the Substrate Solution into the wells must be carried out within 2-3 minutes. Incubate the microplate in the dark **at room temperature (+18...+25°C) for 15 minutes**.
- 10.8. Add **100 µL of Stop Solution** to all wells in the same order as the Substrate Solution. After adding the Stop Solution, the contents of the wells turn yellow.
- 10.9. Read the optical density (OD) of the wells at 450 nm and reference light filters 620–680 nm using a microplate photometer within 5 minutes of adding the Stop Solution.
- 10.10. Plot a calibration curve in linear coordinates: (x) is the C-peptide concentration in the calibrators ng/mL, (y) – OD versus C-peptide concentration (OD 450 nm / 620–680 nm). Manual or computerized data reduction is applicable at this stage. For the algorithm calculation (approximation) of the calibration curve, using the interval (segment-linear, point-to-point) method is recommended.
- 10.11. Determine the corresponding concentration of C-peptide in tested samples from the calibration curve.

11. TEST VALIDITY

The test run shall be considered valid if the OD of CAL1 is above 0.15, and the values of the Control Serum fall into the required range (see Quality control Data Sheet).

12. EXPECTED VALUES

Therapeutical consequences should not be based on the results of IVD methods alone – all available clinical and laboratory findings should be used by a physician to elaborate therapeutically measures. Each laboratory should establish its own normal range for C-peptide. Based on data obtained by XEMA LLC, the following normal range is recommended (see below).

NOTE: values of C-peptide concentrations in the tested samples that are below the LoD (0.015 ng/mL) and also exceed the value of the upper Calibrator (30 ng/mL) should be provided in the following form : «the C-peptide concentration of tested sample X is «lower than 0.015 ng/mL» or «higher than 30* ng/mL».*

* - concentration of the C-peptide of the upper calibration sample may slightly differ from the specified value, the exact value is indicated on the component label.

The calibrators concentration values of the C-peptide EIA kit are expressed in ng/mL. To calculate concentrations in pmol/L, the received concentration value in ng/mL shall be multiplied by 331.

$$1 \text{ ng/mL} = 331 \text{ pmol/L}$$

Sex, age	Units, ng/mL		Units alternative, pmol/L	
	Lower limit	Upper limit	Lower limit	Upper limit
Healthy donors	0.9	5.0	298	1655

13. PERFORMANCE CHARACTERISTICS

13.1. Analytical performance characteristics

13.1.1. Precision of Measurement

Repeatability (Intra assay repeatability) was determined by evaluation the coefficient of variation (CV) for 2 different samples during 1 day in 24 replicates on one series of ELISA kit.

Sample	Concentration, ng/mL	CV, %
1	5.12	3.2
2	3.32	2.8

Reproducibility (Inter assay reproducibility) was determined by evaluating the coefficients of variation for 2 samples during 5 days in 8-replicate determinations.

Sample	Concentration, ng/mL	CV, %
1	5.27	4.0
2	3.87	2.4

Reproducibility between lots was investigated by testing samples for one day on three lots. Each sample was run in 8 replicates.

Sample	Concentration 1, ng/mL	Concentration 2, ng/mL	Concentration 3, ng/mL	CV, %
1	5.32	5.02	5.81	7.2
2	3.71	3.56	3.32	6.6

13.1.2. Trueness

The trueness of measurement is the degree of closeness of the average value obtained from a large number of measurement results to the true value. The bias of the measurement result (bias of measurements) is the difference between the mathematical expectation of the measurement result and the true value of the measurand. The bias was calculated for each sample and it was determined that it corresponds to the specified limits of $\pm 10\%$.

13.1.3. Linearity

Linearity was determined using sera samples with known C-peptide concentration (low and high) and mixing them with each other and buffer solution in different proportions. According to the measurements, linear range of kit is 5-20 ng/mL \pm 10%.

13.1.4. Analytical sensitivity

Limit of detection (LoD) – the lowest C-peptide concentration in the serum or plasma sample that is detected by the C-peptide EIA kit is no lower than 0.015 ng/mL.

Limit of quantification (LoQ) – the lowest concentration of the analyte in the sample that is determined quantitatively with the declared trueness for C-peptide EIA kit is 0.2 ng/mL.

13.1.5. Hook Effect

Hook effect is absent for all samples up to reasonably foreseen concentrations 40 ng/mL.

13.1.6. Analytical specificity

For the analysis result is not affected by the presence in the sample of bilirubin in a concentration of up to 0.21 mg/mL and hemoglobin in a concentration of up to 10 mg/mL.

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**A SOLID-PHASE ENZYME IMMUNOASSAY
FOR THE QUANTITATIVE DETERMINATION
OF TOTAL IgG IN HUMAN BIOLOGICAL FLUIDS**

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Instruction for use

A SOLID-PHASE ENZYME IMMUNOASSAY FOR THE QUANTITATIVE DETERMINATION OF TOTAL IgG IN HUMAN BIOLOGICAL FLUIDS

1. INTENDED USE

A solid-phase enzyme immunoassay for the quantitative determination of total IgG in biological fluids.

This kit is designed for measurement of total IgG in biological fluids. For possibility of use with other sample types, please, refer to Application Notes (on request). The kit contains reagents sufficient for 96 determinations and allows to analyze 42 unknown samples in duplicates.

2. SUMMARY AND EXPLANATION

Immunoglobulin G (IgG) is the main part of serum γ - globulin fraction. IgG is secreted during secondary immune response and plays a key role in humoral immunity. Decrease of serum IgG concentration below 5 g/l is a marker of severe life-threatening immunodeficiency. Determination of serum IgG concentration and IgG/IgA/IgM ratios can be used for monitoring of humoral immune status. Marked elevation of serum IgG may be observed in chronic inflammation, autoimmune diseases and myeloma.

3. PRINCIPLE OF THE TEST

This test is based on two-site sandwich enzyme immunoassay principle. Tested specimen is placed into the microwells coated by specific murine monoclonal to human total IgG-antibodies. Antigen from the specimen is captured by the antibodies coated onto the microwell surface. Unbound material is removed by washing procedure. Second antibodies – murine monoclonal to human total IgG, labelled with peroxidase enzyme, are then added into the microwells. After subsequent washing procedure, the remaining enzymatic activity bound to the microwell surface is detected and quantified by addition of chromogen-substrate mixture, stop solution and photometry at 450 nm. Optical density in the microwell is directly related to the quantity of the measured analyte in the specimen.

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5.1. Contents of the Kit

5. KIT COMPONENTS

Symbol	Description	Qty	Units	Colour code	Stability of opened/diluted components
1	SORB MTP total IgG EIA strips, 8x12 wells	1	pcs		until exp.date
2	CAL 1-5 Calibrator set, 1 ml each. The set contains 5 calibrators: 0; 1; 5; 10; 25 g/l	5	pcs	blue (Cl – colourless)	2 months
3	CONTROL Control serum (1 ml)	1	pcs	colourless	2 months
4	CONJ HRP Conjugate, 14 ml	1	pcs	red	until exp.date
5	DIL SPE EIA sample buffer 100 ml	1	pcs	blue	until exp.date
6	SUBS TMB Substrate solution, 14 ml	1	pcs	colourless	until exp.date
7	BUF WASH 26X Washing solution concentrate 26X, 22 ml	1	pcs	colourless	Concentrate – until exp.date Diluted washing solution – 45 days at 2-8 °C or 15 days at RT
8	STOP Stop solution, 14 ml	1	pcs	colourless	until exp.date
9	N003 Plate sealing tape	2	pcs		N/A
10	K2711 Instruction total IgG EIA	1	pcs		N/A
11	K271Q QC data sheet total IgG EIA	1	pcs		N/A

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- Distilled or deionized water;
- Automatic or semiautomatic multichannel micropipettes, 100–250 µl, is useful but not essential;
- Calibrated micropipettes with variable volume, range volume 10–250 µl;
- Dry thermostat for 37 °C ±0.1 °C
- Calibrated microplate photometer with 450 nm wavelength and OD measuring range 0–3.0

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7.4. Assay procedure

1	Put the desired number of microstrips into the frame; allocate 12 wells for the calibrators CAL 1-5 and control samples CONTROL and two wells for each unknown sample. DO NOT REMOVE ADHESIVE SEALING TAPE FROM UNUSED STRIPS.
2	Dilute samples using buffer DIL SPE (EIA sample buffer) 5000 fold. See table M for dilution modes and factors for different types of analyzed material. Do not dilute control sample and calibrators.
3	If suggested analyte concentration in the sample exceeds the highest calibrator, additionally dilute this sample accordingly, using DIL SPE (EIA sample buffer). Use of other buffers or reagents for sample dilution may lead to incorrect measurement.
4	Pipet 100 µl of calibrators CAL 1-5 and control samples CONTROL into allocated wells. For testing of blood serum or plasma pipet 100 µl of the unknown diluted sample (DILUTION 2) into the allocated wells. See table M for the volumes of other materials. Pipetting should be made within 3 minutes, to ensure an uniform incubation time for all samples. Carefully mix the contents of the wells by short horizontal rotating of the plate for 5-7 seconds and cover the wells by plate adhesive tape (included into the kit).
5	Incubate 30 minutes at +37 °C.
6	Prepare washing solution by 26X dilution of washing solution concentrate BUF WASH 26X by distilled water. Minimal quantity of washing solution should be 250 µl per well. Wash strips 3 times.
7	Dispense 100 µl of CONJ HRP into the wells. Cover the wells by plate adhesive tape.
8	Incubate 30 minutes at +37 °C.
9	Wash the strips 5 times.
10	Dispense 100 µl of SUBS TMB into the wells.
11	Incubate 10-20 minutes at +18...+25 °C.
12	Dispense 100 µl of STOP into the wells.
13	Measure OD (optical density) at 450 nm.
14	Set photometer blank on first calibrator.
15	Apply point-by-point method for data reduction. Use Calculation factor listed in table M to calculate analyte concentration in different material types.

7.5. Sample processing

Material type	Notes on material collection, storage and handling	Sample dilution example	EIA sample buffer into the well, µl	Sample into the well, µl	Calculation factor
blood serum or plasma	Grossly hemolytic, lipemic, or turbid samples should be avoided and should be treated by centrifugation before testing.	10 µl of sample + 990 µl of diluent = DILUTION 1. 10 µl of DILUTION 1 + 490 µl of diluent = DILUTION 2	0	100	1
saliva	Grossly hemolytic, lipemic, or turbid samples should be avoided and should be treated by centrifugation before testing.		90	10	0.002
urine	Grossly hemolytic, lipemic, or turbid samples should be avoided and should be treated by centrifugation before testing.		50	50	0.0004
cerebrospinal fluid	Grossly hemolytic, lipemic, or turbid samples should be avoided and should be treated by centrifugation before testing.	10 µl of sample + 500 µl of diluent	0	100	0.01

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
9.1.2. Plot a calibration curve on graph paper: OD versus total IgG concentration.

9.1.3. Determine the corresponding concentration of total IgG in unknown samples from the calibration curve. Manual or computerized data reduction is applicable on this stage. Point-by-point or linear data reduction is

recommended due to non-linear shape of curve.

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10. EXPECTED VALUES

Therapeutical consequences should not be based on results of IVD methods alone – all available clinical and laboratory findings should be used by a physician to elaborate therapeutically measures. Each laboratory should establish its own normal range for total IgG. Based on data obtained by XEMA, the following normal range is recommended (see below). NOTE: the patients that have received murine monoclonal antibodies for radioimaging or immunotherapy develop high titered anti-mouse antibodies (HAMA). The presence of these antibodies may cause false results in the present assay. Sera from HAMA positive patients should be treated with depleting adsorbents before assaying.

Sex, age	Units, g/l	
	Lower limit	Upper limit
newborn	7.0	15
1-3 month	2.7	8.0
4-6 month	1.8	8.5
7-12 month	3.5	12
1-6 yrs	6.5	18
7-11 yrs	8.5	15
> 11 yrs	9.0	20

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11.3. Linearity. Linearity was checked by assaying dilution series of 5 samples with different total IgG concentrations. Linearity percentages obtained ranged within 90 to 110%.

11.4. Recovery. Recovery was estimated by assaying 5 mixed samples with known total IgG concentrations. The recovery percentages ranged from 90 to 110%.

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Instruction for use

A SOLID-PHASE ENZYME IMMUNOASSAY FOR THE QUANTITATIVE DETERMINATION OF TOTAL IGA IN HUMAN BIOLOGICAL FLUIDS

1. INTENDED USE

A solid-phase enzyme immunoassay for the quantitative determination of total IgA in biological fluids.

This kit is designed for measurement of total IgA in biological fluids. For possibility of use with other sample types, please, refer to Application Notes (on request). The kit contains reagents sufficient for 96 determinations and allows to analyze 42 unknown samples in duplicates.

2. SUMMARY AND EXPLANATION

Immunoglobulin A (IgA) is a main factor of mucosal immune response to bacteria and viruses. Selective IgA deficiency is one of the most frequent hereditary disorders causing chronic infections inflammation in gastrointestinal, urinary and respiratory systems. Determination of IgA concentration in serum and other biological fluids can be used as screening for selective IgA deficiency and other immunodeficiency syndromes. Marked elevation of serum IgA is observed in some autoimmune diseases and IgA myeloma.

3. PRINCIPLE OF THE TEST

This test is based on two-site sandwich enzyme immunoassay principle. Tested specimen is placed into the microwells coated by specific murine monoclonal to human total IgA-antibodies. Antigen from the specimen is captured by the antibodies coated onto the microwell surface. Unbound material is removed by washing procedure. Second antibodies – murine monoclonal to human total IgA, labelled with peroxidase enzyme, are then added into the microwells. After subsequent washing procedure, the remaining enzymatic activity bound to the microwell surface is detected and quantified by addition of chromogen-substrate mixture, stop solution and photometry at 450 nm. Optical density in the microwell is directly related to the quantity of the measured analyte in the specimen.

4. WARNINGS AND PRECAUTIONS

4.1. For professional use only.

4.2. This kit is intended for in vitro diagnostic use only.

4.3. INFECTION HAZARD: There is no available test methods that can absolutely assure that Hepatitis B and C viruses, HIV-1/2, or other infectious agents are not present in the reagents of this kit. All human products, including patient samples, should be considered potentially infectious. Handling and disposal should be in accordance with the procedures defined by an appropriate national biohazard safety guidelines or regulations.

4.4. Avoid contact with stop solution containing 5.0% H_2SO_4 . It may cause skin irritation and burns.

4.5. Wear disposable latex gloves when handling specimens and reagents. Microbial contamination of reagents may give false results.

4.6. Do not use the kit beyond the expiration date.

4.7. All indicated volumes have to be performed according to the protocol. Optimal test results are only obtained when using calibrated pipettes and microplate readers.

4.8. Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.

4.9. Chemicals and prepared or used reagents have to be treated as hazardous waste according to the national biohazard safety guidelines or regulations.

4.10. Do not mix reagents from different lots.

4.11. Replace caps on reagents immediately. Do not swap caps.

4.12. Do not pipette reagents by mouth.

4.13. Specimens must not contain any AZIDE compounds – they inhibit activity of peroxidase.

4.14. Material Safety Data Sheet for this product is available upon request directly from XEMA Co., Ltd.

4.15. The Material Safety Data Sheet fit the requirements of EU Guideline 91/155 EC.

5.1. Contents of the Kit

5. KIT COMPONENTS

Symbol	Description	Qty	Units	Colour code	Stability of opened/diluted components
1	SORB MTP total IgA EIA strips, 8x12 wells	1	pcs		until exp.date
2	CAL 1-5 Calibrator set, 1 ml each. The set contains 5 calibrators: 0; 0.1; 0.5; 2; 5 g/l	5	pcs	blue (C1 – colourless)	2 months
3	CONTROL Control serum (1 ml)	1	pcs	colourless	2 months
4	CONJ HRP Conjugate, 14 ml	1	pcs	blue	until exp.date
5	DIL SPE EIA buffer 100 ml	1	pcs	blue	until exp.date
6	SUBS TMB Substrate solution, 14 ml	1	pcs	colourless	until exp.date
7	BUF WASH 26X Washing solution concentrate 26X, 22 ml	1	pcs	colourless	Concentrate – until exp.date Diluted washing solution - 45 days at 2-8 °C or 15 days at RT
8	STOP Stop solution, 14 ml	1	pcs	colourless	until exp.date
9	N003 Plate sealing tape	2	pcs		N/A
10	K275I Instruction total IgA EIA	1	pcs		N/A
11	K275Q QC data sheet total IgA EIA	1	pcs		N/A

5.2. Equipment and material required but not provided

- Distilled or deionized water;
- Automatic or semiautomatic multichannel micropipettes, 100–250 μl , is useful but not essential;
- Calibrated micropipettes with variable volume, range volume 5–250 μl ;
- Dry thermostat for +37 °C \pm 0.1 °C
- Calibrated microplate photometer with 450 nm wavelength and OD measuring range 0–3.0

5.3. Storage and stability of the Kit

Store the whole kit at +2...+8 °C upon receipt until the expiration date.

After opening the pouch keep unused microtiter wells TIGHTLY SEALED BY ADHESIVE TAPE (INCLUDED) to minimize exposure to moisture.

6. SPECIMEN COLLECTION AND STORAGE

This kit is intended for use with serum or plasma (ACD- or heparinized). Grossly hemolytic, lipemic, or turbid samples should be avoided.

Specimens may be stored for up to 48 hours at +2...+8 °C before testing.

7. TEST PROCEDURE**7.1. Reagent Preparation**

- All reagents (including unsealed microstrips) should be allowed to reach room temperature (+18...+25 °C) before use.
- All reagents should be mixed by gentle inversion or vortexing prior to use. Avoid foam formation.
- It is recommended to spin down shortly the tubes with calibrators on low speed centrifuge.
- Prepare washing solution from the concentrate BUF WASH 26X by 26 dilutions in distilled water.

7.2. Procedural Note:

It is recommended that pipetting of all calibrators and samples should be completed within 3 minutes.

7.3. Assay flowchart

See the example of calibration graphic in Quality Control data sheet.

7.4. Assay procedure

1	Put the desired number of microstrips into the frame; allocate 12 wells for the calibrators CAL 1-5 and control samples CONTROL and two wells for each unknown sample. DO NOT REMOVE ADHESIVE SEALING TAPE FROM UNUSED STRIPS.
2	Dilute samples using buffer DIL SPE (EIA buffer) 5000 fold. See table M for dilution modes and factors for different types of analyzed material. Do not dilute control sample and calibrators.
3	If suggested analyte concentration in the sample exceeds the highest calibrator, additionally dilute this sample accordingly, using DIL SPE (EIA buffer). Use of other buffers or reagents for sample dilution may lead to incorrect measurement.
4	Pipet 100 µl of calibrators CAL 1-5 and control samples CONTROL into allocated wells. For testing of blood serum or plasma pipet 100 µl of the diluted unknown sample into the allocated wells. See table M for the volumes of other materials. Pipetting should be made within 3 minutes, to ensure an uniform incubation time for all samples. Carefully mix the contents of the wells by short horizontal rotating of the plate for 5-7 seconds and cover the wells by plate adhesive tape (included into the kit).
5	Incubate 30 minutes at +37 °C.
6	Prepare washing solution by 26X dilution of washing solution concentrate BUF WASH 26X by distilled water. Minimal quantity of washing solution should be 250 µl per well. Wash strips 3 times.
7	Dispense 100 µl of CONJ HRP into the wells.
8	Incubate 30 minutes at +37 °C.
9	Wash the strips 5 times.
10	Dispense 100 µl of SUBS TMB into the wells.
11	Incubate 10-20 minutes at +18...+25 °C.
12	Dispense 100 µl of STOP into the wells.
13	Measure OD (optical density) at 450 nm.
14	Set photometer blank on first calibrator.
15	Apply point-by-point method for data reduction. Use Calculation factor listed in table M to calculate analyte concentration in different material types.

7.5. Sample processing

Material type	Notes on material collection, storage and handling	Sample dilution example	EIA buffer into the well, µl	Sample into the well, µl	Calculation factor
blood serum or plasma	Grossly hemolytic, lipemic, or turbid samples should be avoided and should be treated by centrifugation before testing.	10 µl of sample + 990 µl of diluent = DILUTION 1. 10 µl of DILUTION 1 + 490 µl of diluent = DILUTION 2	0	100	1
saliva	Grossly hemolytic, lipemic, or turbid samples should be avoided and should be treated by centrifugation before testing.	5 µl of sample + 500 µl of diluent	90	10	0.2
urine	Grossly hemolytic, lipemic, or turbid samples should be avoided and should be treated by centrifugation before testing.		80	20	0.001
cerebrospinal fluid	Grossly hemolytic, lipemic, or turbid samples should be avoided and should be treated by centrifugation before testing.		50	50	0.0004

8. QUALITY CONTROL

It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day to day validity of results.

The test must be performed exactly as per the manufacturer's instructions for use. Moreover the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable federal, state, and local standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications.

9. CALCULATION OF RESULTS


9.1. Calculate the mean absorbance values (OD450) for each pair of calibrators and samples.

9.2. Plot a calibration curve on graph paper: OD versus total IgA concentration.

9.3. Determine the corresponding concentration of total IgA in unknown samples from the calibration curve. Manual or computerized data reduction is applicable on this stage. Point-by-point or linear data reduction is recommended due to non-linear shape of curve.

9.4. Below is presented a typical example of a standard curve with the XEMA Co. Not for calculations!

Calibrators	Value	Absorbance Units (450 nm)
CAL 1	0 g/l	0.08
CAL 2	0.1 g/l	0.40
CAL 3	0.5 g/l	1.32
CAL 4	2 g/l	2.45
CAL 5	5 g/l	2.92

10. EXPECTED VALUES

Therapeutical consequences should not be based on results of IVD methods alone – all available clinical and laboratory findings should be used by a physician to elaborate therapeutically measures. Each laboratory should establish its own normal range for total IgA. Based on data obtained by XEMA, the following normal range is recommended (see below). NOTE: the patients that have received murine monoclonal antibodies for radioimaging or immunotherapy develop high titered anti-mouse antibodies (HAMA). The presence of these antibodies may cause false results in the present assay. Sera from HAMA positive patients should be treated with depleting adsorbents before assaying.

Sex, age	Units, g/l	
	Lower limit	Upper limit
Healthy donors	0.9	5.0
>61 yr	1.0	6.5
newborn	-	0.05
1-3 month	0.06	0.6
4-6 month	0.1	1.0
7-12 month	0.35	1.7
1-6 yrs	0.8	2.2
7-11 yrs	0.9	2.6

11. PERFORMANCE CHARACTERISTICS

11.1. Analytical specificity / Cross reactivity

Analyte	Cross-reactivity, % wt/wt
IgG	<0.1
IgM	<0.1
IgE	<0.1

11.2. Analytical sensitivity. Sensitivity of the assay was assessed as being 0.06 g/l.

11.3. Linearity. Linearity was checked by assaying dilution series of 5 samples with different total IgA concentrations. Linearity percentages obtained ranged within 90 to 110%.

11.4. Recovery. Recovery was estimated by assaying 5 mixed samples with known total IgA concentrations. The recovery percentages ranged from 90 to 110%.

12. LITERATURE

1. Heiddis B. Valdimarsdottir and Arthur A. Stone – Psychosocial Factors and Secretory Immunoglobulin A. Critical Reviews in Oral Biology & Medicine, Jan 1997; 8: 461 – 474.

2. Amir H Abdul Latiff and Michael A Kerr – The clinical significance of immunoglobulin A deficiency. Ann Clin Biochem, Mar 2007; 44: 131 – 139.

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Instruction for use

A SOLID-PHASE ENZYME IMMUNOASSAY FOR THE QUANTITATIVE DETERMINATION OF TOTAL IgM IN HUMAN BIOLOGICAL FLUIDS

1. INTENDED USE

A solid-phase enzyme immunoassay for the quantitative determination of total IgM in biological fluids.

This kit is designed for measurement of total IgM in biological fluids. For possibility of use with other sample types, please, refer to Application Notes (on request). The kit contains reagents sufficient for 96 determinations and allows to analyze 42 unknown samples in duplicates.

2. SUMMARY AND EXPLANATION

Immunoglobulin M (IgM) is secreted during primary immune response and exists in monomeric and pentameric forms. Elevated serum IgM is observed in chronic inflammation, macroglobulinemia and IgM myeloma. Decreased IgM level may occur in some immunodeficiency syndromes.

3. PRINCIPLE OF THE TEST

This test is based on two-site sandwich enzyme immunoassay principle. Tested specimen is placed into the microwells coated by specific murine monoclonal to human total IgM-antibodies. Antigen from the specimen is captured by the antibodies coated onto the microwell surface. Unbound material is removed by washing procedure. Second antibodies – murine monoclonal to human total IgM, labelled with peroxidase enzyme, are then added into the microwells. After subsequent washing procedure, the remaining enzymatic activity bound to the microwell surface is detected and quantified by addition of chromogen-substrate mixture, stop solution and photometry at 450 nm. Optical density in the microwell is directly related to the quantity of the measured analyte in the specimen.

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4.4. Avoid contact with stop solution containing 5.0% H₂SO₄. It may cause skin irritation and burns.

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5.1. Contents of the Kit

5. KIT COMPONENTS

Symbol	Description	Qty	Units	Colour code	Stability of opened/diluted components
1	total IgM EIA strips, 8x12 wells	1	pcs		until exp.date
2	total IgM EIA strips, 8x12 wells Calibrator set, 1 ml each. The set contains 5 calibrators: 0; 0.5; 2; 5; 10 g/l	5	pcs	purple (C1 – colourless)	2 months
3	Control serum (1 ml)	1	pcs	colourless	2 months
4	Conjugate, 14 ml	1	pcs	purple	until exp. date
5	EIA buffer 100 ml	1	pcs	blue	until exp. date
6	Substrate solution, 14 ml	1	pcs	colourless	until exp.date
7	Washing solution concentrate 26X, 22 ml	1	pcs	colourless	Concentrate – until exp.date Diluted washing solution – 45 days at 2-8 °C or 15 days at RT
8	Stop solution, 14 ml	1	pcs	colourless	until exp. date
9	Plate sealing tape	2	pcs		N/A
10	Instruction total IgM EIA	1	pcs		N/A
11	QC data sheet total IgM EIA	1	pcs		N/A

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- Distilled or deionized water;
- Automatic or semiautomatic multichannel micropipettes, 100–250 µl, is useful but not essential;
- Calibrated micropipettes with variable volume, range volume 5–250 µl;
- Dry thermostat for 37 °C ±0.1 °C
- Calibrated microplate photometer with 450 nm wavelength and OD measuring range 0–3.0

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7. TEST PROCEDURE**7.1.** Reagent Preparation

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7.4. Assay procedure

1	Put the desired number of microstrips into the frame; allocate 12 wells for the calibrators CAL 1–5 and control samples CONTROL and two wells for each unknown sample. DO NOT REMOVE ADHESIVE SEALING TAPE FROM UNUSED STRIPS.
2	Dilute samples using buffer DIL SPE (EIA sample buffer) 5000 fold. See table M for dilution modes and factors for different types of analyzed material. Do not dilute control sample and calibrators.
3	If suggested analyte concentration in the sample exceeds the highest calibrator, additionally dilute this sample accordingly, using DIL SPE (EIA sample buffer). Use of other buffers or reagents for sample dilution may lead to incorrect measurement.
4	Pipet 100 µl of calibrators CAL 1–5 and control samples CONTROL into allocated wells. For testing of blood serum or plasma pipet 100 µl of the unknown diluted sample (DILUTION 2) into the allocated wells. See table M for the volumes of other materials. Pipetting should be made within 3 minutes, to ensure an uniform incubation time for all samples. Carefully mix the contents of the wells by short horizontal rotating of the plate for 5–7 seconds and cover the wells by plate adhesive tape (included into the kit).
5	Incubate 30 minutes at 37 °C.
6	Prepare washing solution by 26X dilution of washing solution concentrate BUF WASH 26X by distilled water. Minimal quantity of washing solution should be 250 µl per well. Wash strips 3 times.
7	Dispense 100 µl of CONJ HRP into the wells.
8	Incubate 30 minutes at 37 °C.
9	Wash the strips 5 times.
10	Dispense 100 µl of SUBS TMB into the wells.
11	Incubate 10–20 minutes at +18...+25 °C.
12	Dispense 100 µl of STOP into the wells.
13	Measure OD (optical density) at 450 nm.
14	Set photometer blank on first calibrator.
15	Apply point-by-point method for data reduction. Use Calculation factor listed in table M to calculate analyte concentration in different material types.

7.5. Sample processing

Material type	Notes on material collection, storage and handling	Sample dilution example	EIA sample buffer into the well, μl	Sample into the well, μl	Calculation factor
blood serum or plasma	Grossly hemolytic, lipemic, or turbid samples should be avoided and should be treated by centrifugation before testing.	10 μl of sample + 990 μl of diluent = DILUTION 1. 10 μl of DILUTION 1 + 490 μl of diluent = DILUTION 2	0	100	1
saliva	Grossly hemolytic, lipemic, or turbid samples should be avoided and should be treated by centrifugation before testing.		90	10	0.002
urine	Grossly hemolytic, lipemic, or turbid samples should be avoided and should be treated by centrifugation before testing.		50	50	0.0004
cerebrospinal fluid	Grossly hemolytic, lipemic, or turbid samples should be avoided and should be treated by centrifugation before testing.		80	20	0.001

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
9.2. Plot a calibration curve on graph paper: OD versus total IgM concentration.

9.3. Determine the corresponding concentration of total IgM in unknown samples from the calibration curve. Manual or computerized data reduction is applicable on this stage. Point-by-point or linear data reduction is recommended due to non-linear shape of curve.

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10. EXPECTED VALUES

Therapeutical consequences should not be based on results of IVD methods alone – all available clinical and laboratory findings should be used by a physician to elaborate therapeutically measures. Each laboratory should establish its own normal range for total IgM. Based on data obtained by XEMA, the following normal range is recommended (see below). NOTE: the patients that have received murine monoclonal antibodies for radioimaging or immunotherapy develop high titered anti-mouse antibodies (HAMA). The presence of these antibodies may cause false results in the present assay. Sera from HAMA positive patients should be treated with depleting adsorbents before assaying.

Sex, age	Units, g/l	
	Lower limit	Upper limit
Healthy donors	0.7	3.7
newborn	0.1	0.35
1-3 month	0.12	0.9
4-6 month	0.25	1.2
7-12 month	0.35	1.0
1-6 yrs	0.55	2.2
7-11 yrs	0.65	1.7

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11.3. Linearity. Linearity was checked by assaying dilution series of 5 samples with different total IgM concentrations. Linearity percentages obtained ranged within 90 to 110%.

11.4. Recovery. Recovery was estimated by assaying 5 mixed samples with known total IgM concentrations. The recovery percentages ranged from 90 to 110%.

12. LITERATURE

Erik J. Wiersma, Cathy Collins, Shafie Fazel, and Marc J. Shulman Structural and Functional Analysis of J Chain-Deficient IgM J. Immunol., Jun 1998; 160: 5979 – 5989.

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