STORZ

KARL STORZ—ENDOSKOPE

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40 7120 50 – Shaver-Handstück DRILLCUT-X® II

40 7120 35 – Shaver-Handstück DRILLCUT-X® II-35

INSTRUCTION MANUAL

40 7120 50 – Shaver Handpiece DRILLCUT-X® II

40 7120 35 – Shaver Handpiece DRILLCUT-X® II-35

MANUAL DE INSTRUCCIONES

40 7120 50 – Pieza manual del shaver DRILLCUT-X® II

40 7120 35 – Pieza manual del shaver DRILLCUT-X® II-35

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Thank you for your expression of confidence in the KARL STORZ brand name. Like all of our other products, this medical device is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern, high-quality piece of equipment from KARL STORZ.

This instruction manual is intended to aid you in the proper installation, connection, and operation of this medical device.

- Failure to observe and follow this instruction manual and the instruction manuals of medical devices used in combination can result in injury to patients, users and third parties as well as damage to the medical device.
- Please read all relevant instruction manuals and always follow the instructions given. Check the functioning of the medical devices used in combination.
- Take note of all the instructions on the device.
- Keep this instruction manual available for ready reference in a conspicuous location near the device.

KARL STORZ reserves the right to make engineering modifications, in particular in the context of ongoing development and improvement of our devices.

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It is recommended that the suitability of the products for the intended procedure be checked prior to use.

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1	Normal use	3
1.1	Intended use	3
1.2	Indications	3
1.3	Contraindications	3
1.4	Complications and other side effects	4
1.5	Target groups	4
1.6	User qualification	4
1.7	Patient population	4
1.8	Intended application areas	4
1.9	Intended conditions of use	4
1.10	User profile of physician and medical support staff	5
2	Safety	6
3	Description	8
3.1	Device overview	8
3.2	Symbols employed	11
3.2.1	Safety instructions	11
3.2.2	Pictograms/packaging symbols	12
4	Initial operation	13
4.1	Unpacking	13
4.2	Scope of supply	13
4.3	Required accessories	13
4.4	Possible combinations	14
4.5	Preparing the device for use	15
4.5.1	Assembling the handle (optional)	15
4.5.2	Connecting the suction and irrigation tubes	16
4.5.3	Assembling the reusable shaver attachments	16
4.5.4	Assembling the single-use shaver attachments	18
4.5.5	Assembling the single-use sinus shaver	19
5	Operation	20
5.1	Basic information on operation	20
5.2	Disassembly	21
5.3	Removing the handle from the handpiece (optional)	21
6	Reprocessing	22
6.1	Overview table	22
6.2	General safety instructions	23
6.3	Accessories	24

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6.4	Preparation for cleaning and disinfection .	24
6.5	Manual precleaning.....	24
6.6	Manual cleaning.....	25
6.7	Manual disinfection.....	25
6.8	Machine cleaning and disinfection.....	25
6.9	Connecting.....	26
6.10	Assembly, inspection and care.....	26
6.11	Packaging systems.....	27
6.12	Sterilization.....	27
6.13	Limits of reprocessing.....	27
7	Service, repair and disposal.....	28
7.1	Warranty.....	28
7.2	Maintenance.....	28
7.3	Repairing device parts.....	28
7.4	Repair program.....	29
7.5	Infection prevention.....	29
7.6	Disposal.....	29
8	Technical description.....	30
8.1	Technical data.....	30
8.2	Operating/storage conditions.....	31
8.3	Directive compliance for DRILLCUT-X® II (40 7120 50) and DRILLCUT® X II-35 (40 7120 35).....	31
9	Troubleshooting.....	32
10	Spare parts and accessories.....	33
10.1	Recommended accessories.....	33
11	Subsidiaries.....	36

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1 Normal use

1.1 Intended use DRILLCUT X® II and DRILLCUT X® II-35 shaver handpiece

The shaver handpiece serves to accommodate and drive rotating and oscillating instruments by means of mechanical transmission and conversion of speed and torque, as well as the rotational direction for shaver and drill applications. The shaver handpiece is non-invasive and is designed for transient use in surgical invasive interventions.

The use of the handpiece in any way other than indicated above is not permitted for safety reasons.

Universal handle 40 7120 90

The handle can be fixed on a corresponding handpiece to allow better handling of the instrument.

1.2 Indications

The shaver handpiece is employed in the following applications:

- Minimally invasive interventions in ENT surgery as well as in plastic and neuroendoscopic procedures for controlled excision, morcellation, removal, ablation, and transection of bone and tissue.

1.3 Contraindications

There are no known contraindications relating directly to the device. The DRILLCUT X® II and DRILLCUT X® II-35 shaver handpiece must not be used if, in the opinion of the responsible physician, the patient's health is endangered through its use, e.g., due to the patient's general condition or if the endoscopic method as such is contraindicated.

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1.4 Complications and other side effects

- Bleeding
- Infections
- Injury to surrounding tissue and bones
- Burns
- Transection of or injury to nerves
- Undesirable ablation of tissue

1.5 Target groups

This instruction manual is intended for physicians, medical support staff and sterilization assistants.

1.6 User qualification

Only physicians and medical support staff with a relevant specialist qualification and who have received training on the device may use the DRILLCUT X® II and DRILLCUT X® II-35 handpiece. Only technically suitable persons may provide training.

1.7 Patient population

Gender: irrelevant
Age: irrelevant
Weight: irrelevant
State of health: suitable for endoscopic interventions, in the physician's judgement

1.8 Intended application areas

Can be used universally for ENT and NEURO.

1.9 Intended conditions of use

The shaver handpiece may only be used in hospitals and doctors' offices in suitable ambient conditions. For ambient conditions see the technical data.

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**1.10 User profile of physician and
medical support staff**

- Recognized medical skills in the field of application (specialist physician, qualified medical staff).
- Adequate powers of comprehension to rationally assess the surgical situation in hand.
- Adequate language skills in at least one of the languages used in the instruction manual.
- Be thoroughly trained in the operation and use of the device.
- Knowledge of the contents of the instruction manual.
- No physical impairments that could diminish perception of activation and alarm signals (visual and acoustic).

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2 Safety

Risks of infection

The instruments are not sterile when delivered. The use of non-sterile instruments poses a risk of infection for patients, users and third parties.

- *Inspect reusable instruments for visible contamination. Reprocess contaminated instruments.*
- *Reprocess reusable instruments before and after each use. When doing so, use validated procedures.*

With sterile single-use instruments, check the expiration date and ensure that the packaging is undamaged.

- *If the expiration date has passed or if the packaging is damaged, dispose of the instruments.*

Combination of system components

The DRILLCUT X® II/DRILLCUT X® II-35 shaver handpiece can only be used as a system in combination with other instruments and devices from KARL STORZ.

- *Suitable system components are listed in Section 4.4 'Possible combinations'.*
- *Combination with non-suitable instruments and devices can result in uncontrolled conduct and injury to the patient, user or third parties.*

Improper use

Incorrect application of medical instruments poses a risk of injury for patients.

- *Users of medical instruments must have an appropriate medical qualification and be acquainted with the application.*

Risks due to damaged parts

Incorrectly assembled and damaged instruments can lead to injuries to the patient.

- *Care must be taken not to leave missing or broken-off components inside the patient.*
- *Check the following properties of instruments*

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- before and after every use to ensure that they:*
- *are free from damage,*
 - *are in full working order,*
 - *have no unintentional rough surfaces,*
 - *have no sharp corners, burred edges,*
 - *have no projecting parts,*
 - *are complete.*

Risks from incorrect maintenance and installation

Maintenance and installation errors endanger the safety of physicians, medical support staff and patients. The instrument works unreliably or only with limited output.



- *Only staff authorized by KARL STORZ may perform the following work:*
 - *Assembly,*
 - *Expansion,*
 - *Readjustments,*
 - *Modifications,*
 - *Repairs.*
- *The electrical installation of the location in which the medical device is connected and operated must conform to the applicable statutes and standards*

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




3.2 Symbols employed

3.2.1 Safety instructions

HAZARD symbol: This is a hazard symbol. It warns of the risk of injury. To prevent injuries, follow all the measures marked with the hazard symbol.

Signal words:

WARNING! designates a possible imminent risk. If this is not avoided, it **could** lead to death or serious injuries.

CAUTION! designates a possible imminent risk. If this is not avoided, it could lead to minor injuries.

ATTENTION! designates a possibly harmful situation. If it is not avoided, the medical devcies could be damaged.

NOTE! indicates special information on the operation of the device or important information.

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3.2.2 Pictograms/packaging symbols

	Manufacturer
	Consult instructions for use
	Catalogue number
	Serial number
	Number of products in the product packaging
	Non-sterile
	Keep away from sunlight
	Keep dry
	Applied part type BF
	This device has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE).
	Federal (USA) law restricts this device to sale by or on the order of a physician

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4 Initial operation

4.1 Unpacking

Carefully remove the shaver handpiece and the accessories from the packaging.

Check for missing items and for evidence of shipping damage.

If necessary: notify the manufacturer or the supplier.

TIP: If possible, retain the original packing materials for later use; these can be used if the device has to be transported, for example.

4.2 Scope of supply

- DRILLCUT X® II shaver handpiece
Model 407120 50 with integrated connecting cable
- DRILLCUT X® II-35 shaver handpiece
Model 407120 35 with integrated connecting cable

4.3 Required accessories

ATTENTION! Non-compatible accessories can lead to malfunctioning or damage to the DRILLCUT X® II/DRILLCUT X® II-35 handpiece.

- Only use the compatible accessories designated by KARL STORZ.

DRILLCUT X® II/DRILLCUT X® II-35 shaver handpiece requires additional accessories for normal use:

- Motor control unit
- Shaver attachments
- Sinus shaver

For compatible accessories, see Section 4.4 Possible combinations.

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
* For a detailed list of compatible shaver attachments and sinus shavers, please refer to our catalogs ENDOSCOPES AND INSTRUMENTS FOR ENT and NEUROSURGERY .

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4.5 Preparing the device for use

4.5.1 Assembling the handle (optional)

ATTENTION! Damage to the clamping jaw ⑪ in the handle.


- Do not turn the rotary knob counterclockwise while the handle is disassembled.

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NOTE! Depending on the application and desired holding position, the handpiece can either be used without the handle or with the optional handle.


NOTE! The handle can be aligned to either point to the front [Fig. 2] or to the back [Fig. 3].

NOTE: Before assembling, ensure that the clamping jaw ⑪ does not project beyond the guide channel.

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1. Insert the handle ⑧ with the guide in the guide grooves ③ of the DRILLCUT X® II/DRILLCUT X® II-35 shaver handpiece and push into the desired position.

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2. Fix the handle ⑧ into place: Turn the rotary knob ⑨ clockwise as far as it will go.
3. Check that the handpiece is securely in place before use.

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4.5.2 Connecting the suction and irrigation tubes

Prerequisite: The suction and irrigation tubes must be connected to the respective suction/irrigation pump. Observe the instruction manual for the device used.

1. Connect the irrigation tube to the irrigation port ④ of the shaver handpiece [Fig. 1].

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
2. Connect the suction tube to the suction port ⑤ of the shaver handpiece [Fig. 2].

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4.5.3 Assembling the reusable shaver attachments


ATTENTION! If sheaths with different sheath codes are assembled, this can result in excessive heat generation at the tip of the shaver attachments as well as excessive wear.

- Ensure that the outer and inner sheaths have the same sheath code [Fig. 3].

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ATTENTION! Damage to the shaver blade due to incorrect alignment of the inner and outer sheaths.

- Make sure that the cutting window of the inner sheath is aligned with the orientation pin of the set wheel.

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
2. Pull back and hold the locking device ② of the DRILLCUT X® II//DRILLCUT X® II-35 shaver handpiece.

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3. Insert the shaver attachment into the handpiece as far as it will go.

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4. Verriegelung ② loslassen.


4. Release the locking device ②.

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☛ The shaver attachment is locked in place.

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5. Check that it is firmly in place by gently pulling axially on the set wheel ⑮.

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6. With **straight** shaver attachments, check the rotation function by turning the set wheel ⑮. It must not be possible to rotate **curved** shaver attachments when locked.

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4.5.4 Assembling the single-use shaver attachments

! WARNING!

Risk of infection from non-sterile instruments:

- Check the packaging for damage and to ensure the expiration date has not passed. If the expiration date has passed or if the packaging is damaged, dispose of the instruments.
- Do not remove the instrument from the packaging until immediately before use.

1. Insert the inner sheath into the outer sheath from behind.

2. Pull back and hold the locking device ② of the DRILLCUT X® II/DRILLCUT X® II-35 shaver handpiece.

3. Insert the shaver insert into the handpiece as far as it will go.
4. Release the locking device ②.
- ☛ The shaver attachment is locked in place.
5. Check that it is firmly in place by gently pulling axially on the set wheel ⑮.
6. Check the rotating function by rotating the set wheel ⑮.

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ATTENTION: The use of an unsuitable operating mode may damage the sinus drill.

- Only operate the DRILLCUT X® II-35 in 35k operating mode with the 35k sinus drills with red connector specified for this purpose.
- Only operate the DRILLCUT X® II-35 with other (12k) sinus drills in 12k operating mode.

1. Insert the inner sheath into the outer sheath from behind.
2. Pull back and hold the locking device ② of the DRILLCUT X® II/DRILLCUT X® II-35 shaver handpiece.

3. Insert the sinus shaver insert into the handpiece as far as it will go.
4. Release the locking device ②.
- ☛ The sinus shaver attachment is locked in place.
5. Check that it is firmly in place by gently pulling axially on the set wheel ⑮.
6. Check the rotating function by rotating the set wheel ⑮.

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5.1 Basic information on operation

! WARNUNG!

If an incorrectly attached handle becomes detached, the shaver handpiece can move in an uncontrolled fashion. This can result in injury to the patient.

- Ensure that the handle is correctly located in place in the guide grooves ③ and has no play.

! WARNUNG!

The use of applied parts outside of the field of vision can cause injury to tissue and damage to accessories.

- Only activate applied parts under visual control.

! CAUTION!

Rotating instruments can generate heat. This can cause the shaver handpiece to heat up and result in burns to the patient and user.

- Always work with suction and irrigation.

Check system combination before use

NOTE: Check the rinsing and suction function in vitro before using the system combination. To do this, connect the DRILLCUT-X® II shaver handpiece with suction drill/shaver attachment to the control unit and activate the suction/rinsing function [Fig. 1].
(Refer to the instruction manual for the control unit used, 96056025D, 96600526D or 960260003D.)

Rinse shaver/drill attachments after use

NOTE: Rinse the suction drill/shaver attachment with the DRILLCUT X® II handpiece (with the cutting opening open) immediately after use by activating the suction function, preferably with distilled water in order to prevent the attachment from becoming blocked [Fig. 1].

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5.2 Disassembly

1. Pull back and hold the locking device ② of the shaver handpiece.
2. Remove the shaver attachment or sinus shaver from the instrument holder on the DRILLCUT-X® II/ DRILLCUT-X® II-35 and release the locking device again [Fig. 1].

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3. Remove the suction tube from the suction port ⑤ on the DRILLCUT-X® II/ DRILLCUT-X® II-35 [Fig. 2].
4. Remove the irrigation tube from the irrigation port ④ of the DRILLCUT-X® II/ DRILLCUT-X® II-35 [Fig. 3].

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5. Pull the inner sheath of the shaver attachment or the sinus shaver out of the outer sheath toward proximal [Fig. 4].

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5.3 Removing the handle from the handpiece (optional)

1. Loosen the rotary knob ⑨ on the handle ⑧: Turn the rotary knob counterclockwise [Fig. 5].
2. Pull the handle axially from the handpiece via the guide grooves ③.

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* with use of the irrigation adaptor 41250 RA

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- This reprocessing step is part of the validated process
 - This reprocessing step may be carried out additionally in order to aid cleaning
 - X On the grounds of potential damage, this reprocessing step must NEVER be performed
- No symbol:** This reprocessing step has neither been validated in respect of efficacy nor in respect of material compatibility by KARL STORZ.

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6.2 General safety instructions

WARNING:

Risk of infection: When carrying out any work on contaminated medical devices, the guidelines of the Employers' Liability Insurance Association and equivalent organizations striving to ensure personal safety as well as other country-specific laws must be observed.

WARNING!

Incorrectly reprocessed instruments expose patients, users and third parties to a risk of infection as well as the danger of an instrument malfunction. The product-accompanying documents must be observed.

- An overview of the procedures validated by KARL STORZ can be found in the manual 'Cleaning, disinfection, sterilization and care of KARL STORZ instruments'.
- This manual can be downloaded or requested by visiting www.karlstorz.com.

WARNING!

Reprocessing of disposable articles can jeopardize the functionality and safety of the medical device and thus the patient and user and lead to infections.

- Do not reprocess disposable products. Dispose of disposable products in accordance with the applicable regulations.

ATTENTION: Using chemicals not approved by KARL STORZ may damage medical devices.

- Only use chemicals approved by KARL STORZ for reprocessing. A complete list can be found on the Internet at www.karlstorz.com.

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ATTENTION: Prolonged immersion and incorrect concentration of the cleaning/disinfection solutions may result in damage.

- When preparing and using solutions, follow the chemical manufacturer's specifications, paying close attention to proper concentration, exposure time and service life.
- The microbiological range of action of the chemicals used must be taken into account.

6.3 Accessories

Accessories required for carrying out reprocessing:

- Brushes: 27652, 27651 UB
- Adaptor: 41250 RA

6.4 Preparation for cleaning and disinfection

Heavy soiling, corrosive solutions and pharmaceuticals must be removed from the medical device immediately after use. To this end, preclean the medical device, by wiping it down and rinsing it, for example. As a general rule, KARL STORZ recommends manual precleaning under cold running water.

6.5 Manual precleaning Brushing the surfaces

Any visible contamination or heavy soiling must be removed by cleaning the surfaces under cold running water with the aid of a brush or sponge [Fig. 1].

Brushing the lumina

The working channels, lumina and hollow spaces must be precleaned under cold running water with the aid of a suitable brush [Fig. 2].

Rinsing out with a water jet pistol

Pull back the locking device ② and insert the cleaning adaptor 41250 RA ⑫ into the instrument holder of the DRILLCUT-X® II/DRILLCUT-X® II-35 [Fig. 3]. Release the locking device ② again [Fig. 4].

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The working channel, lumina and hollow spaces must be rinsed using a water jet pistol. Use a suitable irrigation attachment for this. For rinsing, use pulses of at least 3 x 5 seconds from a jet of water with a pressure of at least 2.5 bar. The cleaning gun with accessories (Art. no. 27660) is well-suited for this purpose [Fig. 5].

6.6 Manual cleaning

The medical device must be completely immersed in a cleaning solution. Through disassembling and opening the instrument, and filling the lumina in a targeted manner, it must be ensured that even surfaces with restricted access are covered and no air bubbles are present. After the necessary exposure time as stated by the manufacturer, clean the instrument mechanically with the aid of brushes or a sponge. Finally, it must be rinsed with cold water to ensure neutralization.

6.7 Manual disinfection

The medical device must be completely immersed in a disinfectant solution. Through disassembling and opening the instrument, and filling the lumina in a targeted manner, it must be ensured that even surfaces with restricted access are covered and no air bubbles are present. At the end of the necessary exposure time as stated by the manufacturer, the medical device must be rinsed several times in order to remove all chemical residues. To this end, water of the best possible quality must be used, taking into account the national regulations. Finally, all of the surfaces, joints, openings, channels and lumina are dried completely with (preferably sterile) compressed air in accordance with the national regulations. The cleaning gun with accessories (Art. no. 27660) is well-suited for this purpose.

6.8 Machine cleaning and disinfection

The following methods for machine decontamination have been validated and approved subject to compliance with the process parameters described in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments':

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• **Machine cleaning/thermal disinfection**

Thermal disinfection is preferred. The relevant national requirements and the A₀ value must be taken into account when using this method.

The selection of a suitable slide-in tray or instrument holder, which should ensure that the medical device is thoroughly rinsed out or through, must take place in consultation with the manufacturer of the device.

NOTE: If necessary, the instrument must be dried off afterwards by hand.

6.9 Connecting

In order to ensure effective machine cleaning and disinfection, the instrument must be connected to the washer and disinfector. Connect the DRILLCUT X® II/DRILLCUT-X® II-35 up to the washer and disinfector via the cleaning adapter 12 (Art. no. 41250 RA) (Fig. 1).

6.10 Assembly, inspection and care

The cleaned and disinfected medical device must be visually inspected for cleanliness, completeness, damage and dryness:

- If residues or contamination are still present, the medical device must be manually cleaned and subjected to a full cleaning and disinfection process once more.
- Damaged or corroded medical devices must be withdrawn from use.
- Joints, threads and glide surfaces must be lubricated locally with instrument oil (Art. no. 27656 B) [Fig. 2].
Treat the handpiece with a care oil (universal spray, Art. no. 280053 B) by spraying it into the instrument holder from the front [Fig. 3].
- Afterwards, a functional check must be carried out.

NOTE: The oil used for this purpose must be suitable for the subsequent sterilization procedure (silicone-free and paraffin- or white oil-based).

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NOTE: Use items from the catalog 'Care, Sterilization and Storage Techniques' during care procedures.

6.11 Packaging systems

Only standardized and approved packaging materials and systems may be used (EN 868 Parts 2-10, EN ISO 11607 Parts 1 + 2, DIN 58953).1

6.12 Sterilization

The procedures as well as the process-relevant parameters for the individually validated methods are described in detail in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments'. The method must be selected taking into account the respective applicable national requirements and in consultation with the device manufacturers.

The following sterilization procedures have been validated and approved by KARL STORZ for this medical device:

- **Steam sterilization using the fractionated prevacuum procedure**

The medical device must be sterilized in its fully assembled state using the fractionated prevacuum procedure (DIN EN ISO 17665-1) at 134°C–137°C with a minimum exposure time of 4 minutes and a maximum of 18 minutes. This procedure is only suitable for thermostable instruments.

NOTE: Greased components must be disassembled for sterilization, so that steam penetration is ensured.

6.13 Limits of reprocessing

The end of the product's service life and correct functioning are determined by mechanical stress and chemical influences within the scope of reprocessing and application.

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7 Service, repair and disposal

7.1 Warranty

The warranty conditions provided can be found in the Standard Conditions of Business of KARL STORZ.

The medical device must always be sent to your local KARL STORZ subsidiary (see the section entitled 'Subsidiaries'), even during the warranty period.

Opening of the equipment or performance of any repairs or modifications to the device by unauthorized persons shall relieve KARL STORZ of any liability for its performance. Any such opening, repair, or modification performed during the warranty period shall invalidate any warranty.

7.2 Maintenance

No further repair measures are required beyond reprocessing.

However, the legal requirements as per the Medical Devices Operator Ordinance (MPBetreibV) for devices to be used in combination apply (see Section 4.4). The relevant warnings in the instruction manuals must be observed.

7.3 Repairing device parts

Only personnel authorized by the manufacturer may repair device parts. The operations described in this manual are exempt from this rule.

The authorized personnel may only repair device parts which the manufacturer designates as repairable.

Upon request, the manufacturer will provide authorized personnel with information on repairs, e.g., circuit diagrams, spare parts lists, setting instructions.

The fact of having such documents at one's disposal does not correspond to authorization.

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7.4 Repair program

For this medical device, individual repair is necessary. To bridge the repair period, you will receive a device on loan, which must then be returned as soon as you receive the repaired medical device.

In the case of repairs, please contact your local KARL STORZ subsidiary (see section ‘Subsidiaries’) or authorized dealer.

7.5 Infection prevention

To prevent infections, contaminated devices and instruments must not be returned.

Decontaminate instruments and devices before returning. KARL STORZ returns contaminated instruments/devices to the sender.

7.6 Disposal

This device has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE).

At the end of its useful operating life, dispose of the device as electronic scrap.

Please ask KARL STORZ, a KARL STORZ subsidiary (see the section ‘Subsidiaries’), or your authorized dealer about the collection point in your area.

Within the scope of application of this directive, KARL STORZ is responsible for the proper disposal of this device.

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8 Technical description

8.1 Technical data

(sequential in combination with the KARL STORZ control
unit **40 7014 20/40 7016 20-1/40 7017 20-1**)

Power consumption	150 – 200 W
Operating temperature	10 – 30°C
Operating mode	Intermittent duty
Duty cycle	50%; t = 1 min
Weight with cable	371 g
Max. dimensions (W × H × D)	140 × 34 × 25 mm
Average life cycle	200 reprocessing cycles
Shaver speed range (DRILLCUT-X® II)	
Minimum speed	300 rpm
Maximum speed	10,000 rpm
Start-up maximum speed	4000 rpm
Rotation mode	Oscillation
Shaver speed range (DRILLCUT-X® II-35)	
Minimum speed	900 rpm
Maximum speed	10,000 rpm
Start-up maximum speed	4000 rpm
Rotation mode	Oscillation
Sinus drill speed range 12k (DRILLCUT-X® II)	
Minimum speed	300 rpm
Maximum speed	12,000 rpm
Start-up maximum speed	12,000 rpm
Rotation mode	Rotation

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**Sinus drill speed range 35k
(DRILLCUT-X® II-35)**

Minimum speed	1.500 U/min
Maximum speed	35.000 rpm
Start-up maximum speed	35.000 rpm
Rotation mode	Rotation

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8.2 Operating/storage conditions

	Temperature	Rel. humidity
Storage	-20°C – +60°C	10% – 90%
Operation	+10°C – +30°C	30% – 70%
Air pressure	800 hPa – 1080 hPa	

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**8.3 Directive compliance for
DRILLCUT-X® II (40 7120 50) and
DRILLCUT® X II-35 (40 7120 35)**

**In accordance with the Medical Device
Directive (MDD):**

This medical device bears a CE mark in accordance with the Medical Devices Directive (MDD) 93/42/EEC.

A code number after the CE mark indicates the responsible notified body.

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(Not included in the scope of supply)

Shaver blades:

straight, 35° curved, 65° curved non-serrated, single-serrated, double-serrated, reusable or disposable (sterile), pack of 5.

Sinus Drill, single-use:

15° curved, 40° curved, 70° curved

For article numbers see catalogs
ENDOSCOPES AND INSTRUMENTS FOR ENT
and NEUROSURGERY

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KARL STORZ Endoscopy Canada Ltd.
7171 Millcreek Drive, Mississauga
Ontario, L5N 3R3, Canada
Phone: +1 905 816-4500, Fax: +1 905 858-4599
Toll free phone: 1-800-268-4880 (Canada only)
Toll free fax: 1-800-482-4198 (Canada only)
E-Mail: info.canada@karlstorz.com

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100, Fax: +1 424 218-8526
Toll free phone: 800 421-0837 (USA only)
Toll free fax: 800 321-1304 (USA only)
E-Mail: info@kseal.com

KARL STORZ Veterinary Endoscopy-America, Inc.
1 South Los Carneros Road
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorzet.com

KARL STORZ Endoscopia Latino-America, Inc.
815 N. W. 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@ksela.com

KARL STORZ Endoscopia México S.A. de C.V.
Edificio Atlantic, Oficina 3G
Calle D e/ 1ra y 3ra
10400 Vedado, Havana, Cuba
Phone: +537 836 95 06, Fax: +537 836 97 76
E-Mail: kstorzcuba@gmail.com

KARL STORZ Endoscopia México S.A. de C.V.
Lago Constanza No 326,
Col. Granada Del. Miguel Hidalgo
C.P. 11520 México D.F.
Phone: +52 (55) 1101 1520
E-Mail: mx-info@karlstorz.com

KARL STORZ Marketing América Do Sul Ltda.
Rua Joaquim Floriano, nº. 413, 20º andar – Itaim Bibi,
CEP-04534-011 São Paulo, Brasil
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680
E-Mail: info@karlstorz.com.br

KARL STORZ Endoscopia Argentina S.A.
Zufriategui 627 6° Piso, B1638 CAA - Vicente Lopez
Provincia de Buenos Aires, Argentina
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773
E-Mail: info@karlstorz.com.ar

KARL STORZ Endoskopi Norge AS
Østensjøveien 15B
0661 Oslo, Norway
Phone: +47 6380 5600, Fax: +47 6380 5601
post@karlstorz.no

KARL STORZ Endoskop Sverige AB
Storsåtragränd 14
12739 Skärholmen, Sweden
14108 Kungens Kurva, Sweden
Phone: +46 8 505 648 00
E-Mail: kundservice@karlstorz.se

KARL STORZ Endoscopy Suomi OY
Taivaltie 5
01610 Vantaa, Finland
Phone: +358 (0)96824774, Fax: +358 (0)96824775
E-Mail: asiakaspalvelu@karlstorz.fi

KARL STORZ SE & Co. KG
Representation Office
Kestucio g. 59
08124 Vilnius, Lithuania
Phone: +370 68567000
E-Mail: LT-LV-info@karlstorz.com

KARL STORZ Endoskopi Danmark A/S
Skovlytoften 33
2840 Holte, Denmark
Phone: +45 45162600, Fax: +45 45162609
E-Mail: marketing@karlstorz.dk

KARL STORZ Endoscopy (UK) Ltd.
415 Perth Avenue, Slough
Berkshire, SL1 4TQ, United Kingdom
Phone: +44 1753 503500, Fax: +44 1753 578124
E-Mail: info-uk@karlstorz.com

KARL STORZ Endoscopie Nederland B.V.
Displayweg 2
3821 BT Amersfoort, Netherlands
Phone: +31 (0)33 4545890
E-Mail: info-nl@karlstorz.com

KARL STORZ Endoscopy Belgium N.V.
Phone: +31 (0)33 4545890
E-Mail: info-be@karlstorz.com

KARL STORZ Endoscopie France S. A. S.
12, rue Georges Guynemer, Quartier de l'Europe
78280 Guyancourt, France
Phone: +33 1 30484200, Fax: +33 1 30484201
E-Mail: marketing-fr@karlstorz.com

KARL STORZ Endoskop Austria GmbH
Landstraßer Hauptstr. 148/1/G1
1030 Wien, Austria
Phone: +43 1 71 56 0470, Fax: +43 1 71 56 0479
E-Mail: storz-austria@karlstorz.at

KARL STORZ Endoscopia Ibérica S. A.
Parque Empresarial San Fernando
Edificio Munich – Planta Baja
28830 Madrid, Spain
Phone: +34 91 6771051, Fax: +34 91 6772981
E-Mail: info-es@karlstorz.com

KARL STORZ Endoscopia Italia S. r. l.
Via dell'Artigianato, 3
37135 Verona, Italy
Phone: +39 045 8222000, Fax: +39 045 8222001
E-Mail: info-ita@karlstorz.com

KARL STORZ Adria Eos d.o.o.
Capraška 6
10000 Zagreb, Croatia
Phone: +385 1 6406 070, Fax: +385 1 6406 077
E-Mail: info@karlstorz.hr

KARL STORZ Endoskopija d.o.o.
Cesta v Gorice 34b
1000 Ljubljana, Slovenia
Phone: +386 1 620 5880, Fax: + 386 1 620 5882
E-Mail: pisarna@karlstorz.si

KARL STORZ Endoskop Polska Marketing Sp. z o.o.
ul. Bojkowska 47
44-100 Gliwice, Poland
Phone: +48 32 706 13 00, Fax: +48 32 706 13 07
E-Mail: info-pl@karlstorz.com

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KARL STORZ ENDOSKOPE
East Mediterranean and Gulf S.A.L.
Block M, 3rd Floor
Beirut Souks, Weygand Street
2012 3301 Beirut, Lebanon
Phone: +961 1 999390, Fax +961 1 999391
E-Mail: info@karlstorz-emg.com

KARL STORZ Endoscopy (South Africa) (Pty) Ltd.
P.O. 6061, Roggebaai 8012
Cape Town, South Africa
Phone: +27 21 417 2600, Fax: +27 21 421 5103
E-Mail: info@karlstorz.co.za

TOO KARL STORZ Endoscopy Kasachstan
Saryarka, 6, BC "Arman", off. 910
010000 Astana, Republic of Kazakhstan
Phone: +7 7172 552-549, 552-788, Fax: -444
E-Mail: info@karlstorz.kz

KARL STORZ ENDOSKOPE
East Mediterranean & Gulf (branch)
Building West Side 7A – Unit 7WA – 3008
Dubai Airport Free Zone, P.O. Box 54983
Dubai - United Arab Emirates
Phone: +971 (0)4 2958887, Fax: +971 (0)4 3205282
Service Hotline: +971 (0)4 3415882
E-Mail: info-gne@karlstorz-emg.com

KARL STORZ Endoscopy India Private Limited
11th Floor, Dr. Gopal Das Bhawan
28, Barakhamba Road
New Delhi 110001, India
Phone: +91 11 4374 3000, Fax: +91 11 4374 3010
E-Mail: corporate@karlstorz.in

KARL STORZ SE & Co. KG
Resident Representative Office
14th Floor, MPlaza Saigon
39 Le Duan, District 1
Ho Chi Minh City, Vietnam
Phone: +84 8 3823 8000, Fax: +84 8 3823 8039
E-Mail: infovietnam@karlstorz.com

KARL STORZ Endoscopy Australia Pty. Ltd.
15 Orion Road Lane Cove NSW 2066
P O Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (0)2 9420 0695
E-Mail: info@karlstorz.au

KARL STORZ Endoscopy China Ltd., Hong Kong
Room 2503-05, 25F AXA Tower, Landmark East,
No. 100 How Ming Street, Kwun Tong, Kowloon,
Hong Kong, People's Republic of China
Phone: +852 28 65 2411, Fax: +852 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd.
Beijing Branch,
Room 1805-1807, Building B, 18F Beijing IFC
No. 8, Jianguomenwai Street, Chaoyang District,
Beijing, 100022, People's Republic of China
Phone: +86 10 5638188, Fax: +86 10 5638199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shanghai Branch
36F, Tower 1 Grand Gateway
No. 1 Hongguiao Road, Xuhui District,
Shanghai, 200030, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch
Room 805, 8F Jin Jiang International Building
No. 1 West Linjiang Road, Wuhou District,
Chengdu, 6100414, People's Republic of China
Phone: +86 28 86587977, Fax: +86 28 86587975
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shenyang Branch,
Room 2001-2005, 20F N-MEDIA International Center,
No. 167 Youth Avenue, Shenhe District,
Shenyang, 110014, People's Republic of China
Phone: +86 24 23181118, Fax: +86 24 23181119
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Guangzhou Branch
Room 02B, 35F Teem Tower,
No. 208 Tianhe Road, Tianhe District,
Guangzhou, 510620, People's Republic of China
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd.
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Singapore Sales Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Korea Co. Ltd.
I-Park Bundang1, 102dong, 8F
239 Jeongjail-Ro, Bundang-Gu
Seongnam, 13556, Korea
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299
E-Mail: infokorea@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.,
12F, No. 192, Sec. 2, Chung Hsin Rd.,
Sindian District, New Taipei City, Taiwan
Phone: +886 933 014 160, Fax: +886 2 8672 6399
E-Mail: mingwang@karlstorz.tw

KARL STORZ SE & Co. KG
Representative Office Philippines
1901 Picadilly Star Bldg., 4th Avenue, BGC
Taguig City 1636, Philippines
Phone: +63 2 317 45 00, Fax: +63 2 317 45 11
E-Mail: philippines@karlstorz.com

KARL STORZ Endoscopy Japan K. K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 3 6380-8633
E-Mail: info@karlstorz.co.jp

www.karlstorz.com

KARL STORZ SE & Co. KG

Dr.-Karl-Storz-Straße 34
78532 Tuttlingen
Postfach 230
78503 Tuttlingen
Germany

Telefon: +49 (0)7461 708-0
Telefax: +49 (0)7461 708-105
E-Mail: info@karlstorz.com
Web: www.karlstorz.com