Medtronic

Abre™

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
Figure 1. Abre™ venous self-expanding stent

Device description

The Abre™ venous self-expanding stent system (Abre™ system) is a vascular self-expanding nitinol stent system. The Abre™ stent is intended for permanent implant and comes premounted on a 9 Fr (3.1 mm) delivery system inserted over a 0.89 mm (0.035 in) wire.

- The Abre™ self-expanding stent, as shown in Figure 1, is cut from a nickel titanium alloy (nitinol) tube in an open lattice design (1), with integral nitinol markers at the trailing end (2) and the leading end (3) of the stent. Upon deployment, the Abre stent exerts an outward force to establish patency.
- The Abre™ delivery system, as shown in Figure 1, has a triaxial shaft design, comprised of an inner shaft assembly (4), a retractable sheath (silver, 5), an isolation sheath (blue, 6), and an ergonomic handle (7). The inner shaft assembly terminates with a flexible catheter tip (9) and originates at the luer hub (10). The deployment handle (7) has a removable locking pin (11), thumbwheel (12), and luer hub (10).
- The radiopaque markers, at the trailing end and the leading end of the stent, guide stent positioning at the target lesion before Abre stent deployment. There is a single luer hub on the deployment handle. A thumbwheel on the deployment handle rotates to pull back the retractable sheath (silver). A locking pin prevents the stent from being deployed before use, and must be removed to actuate the thumbwheel. The Abre stent is fully deployed when the radiopaque marker (8) in the retractable sheath (silver) reaches beyond the integral nitinol markers on the trailing end of the stent. A uniquely designed isolation sheath (blue, attached to the deployment handle) improves control and accuracy of stent delivery.

Indications for use

The Abre™ venous self-expanding stent system is intended for use in the iliofemoral veins for the treatment of symptomatic venous outflow obstruction.

Contraindications

- Do not use the Abre system with patients in whom anticoagulant or antiplatelet therapy is contraindicated.
- Do not use the Abre system with patients with known hypersensitivity to nickel titanium (nitinol).

Warnings

- The device is provided sterile, for single use only. Do not reprocess or resterilize. Reprocessing or resterilizing could increase the risk of patient infection and the risk of compromised device performance.
- If unusually high resistance is encountered when advancing the Abre delivery system over the guidewire, assess the cause of the resistance before proceeding.
- If high resistance is felt when initially rotating the thumbwheel, do not force deployment. Carefully withdraw the system and do not use it.

Precautions

- Use of the Abre stent in vessel beds other than the iliofemoral vein region has not been studied.
- Abre stent system implant procedures should only be performed by physicians that have experience with interventional techniques in the vascular system.
- Appropriate stent size selection is crucial. Use Table 1 for guidance in selecting the appropriate stent size. Proper size selection reduces stent migration and ensures appropriate stent apposition to the vessel wall.
- Carefully inspect the sterile package and device before use to verify that no damage occurred during shipment. Do not use

the stent system if it is damaged or compromised.

- Carefully remove the system from the tray and tube without kinking the system.
- Do not use the stent system if the Abre stent is partially deployed before starting the procedure.
- Gaining access too close to the intended treatment site can lead to difficulty with stent deployment.
- Always use an introducer sheath during the implant procedure to protect both the vessel and puncture site.
- Ensure that the locking pin remains in locked position until the target site is reached.
- If the retractable sheath (silver) is exposed beyond the hemostatic valve, friction between the retractable sheath and the hemostatic valve can impair deployment. To avoid friction between the retractable sheath (silver) and the hemostatic valve, pull the introducer sheath back so that the hemostatic valve covers the end of the isolation sheath.
- Do not reposition the Abre stent after establishing apposition against the vessel wall.
- The Abre delivery system is not designed for recapturing the stent.
- The Abre stent is not designed to be lengthened or shortened from its nominal length. Excessive stent lengthening or shortening can increase the risk of stent fracture.
- If resistance is felt when withdrawing the system, do not force withdrawal.
- Use care when crossing a deployed Abre stent with any adjunct device to avoid stent dislodgement or damage to the adjunct device.

Potential adverse events

- Access failure
- Access site infection
- Allergic reaction to contrast medium or procedure medications
- Allergic reaction to nitinol or other device materials
- Arrhythmia
- AV fistula
- Bleeding
- Bruising
- Death
- Device breakage
- Device maldeployment
- Edema
- Fever
- Hematoma
- Hypotension, nausea, or other vasovagal response
- Infection
- Myocardial infarction
- Pain
- Pseudoaneurysm
- Pulmonary embolism
- Renal insufficiency/renal failure (new or worsening)
- Sepsis
- Stent fracture
- Stent malapposition
- Stent malposition
- Stent migration
- Stroke/paradoxical embolism/transient ischemic attack/intracerebral hemorrhage
- Tissue necrosis
- Transfusion reaction following blood transfusion for treatment

of major bleeding

- Vessel damage, including perforation or rupture
- Venous occlusion/thrombosis, within or outside of stented segment

Stent size selection

Considering the estimated anatomic vessel diameter, use Table 1 to select the Abre™ stent diameter size. Choose a stent length that extends beyond both ends of the target lesion, with at least 1 cm on each side of the lesion to reduce the risk of restenosis.

Table 1. Sizing guide

Stent diameter (mm)	Estimated anatomic vessel diameter (mm)	Stent length (mm)
10	7.5–9.5	40, 60, 80, 100, 120, 150
12	9.5–11.5	60, 80, 100, 120, 150
14	11.5–13.5	60, 80, 100, 120, 150
16	13.5–15.5	60, 80, 100, 120, 150
18	15.5–17.5	60, 80, 100, 120, 150
20	17.5–19.0	60, 80, 100, 120, 150

Caution: Appropriate stent size selection is crucial. Use Table 1 for guidance in selecting the appropriate stent size. Proper size selection reduces stent migration and ensures appropriate stent apposition to the vessel wall.

How supplied

Warning: The device is provided sterile, for single use only. Do not reprocess or resterilize. Reprocessing or resterilizing could increase the risk of patient infection and the risk of compromised device performance.

Directions for use

Required items for implant procedure

- 5–10 cc syringe filled with saline
- 0.89 mm (0.035 in) guidewire
- 9 Fr minimum hemostatic introducer sheath
- Balloon catheter

Preparation

1. Open the shelf box and remove the contents. Carefully inspect the pouch, looking for damage to the sterile barrier; then carefully peel open the pouch and extract the device by the handle.
Caution: Carefully inspect the sterile package and device before use to verify that no damage occurred during shipment. Do not use the stent system if it is damaged or compromised.
Caution: Carefully remove the delivery system from the tray and tube without kinking the system.
2. Attach a syringe filled with saline to the luer lock injection hub. Inject the saline through the guidewire lumen until it comes out of the catheter tip.
3. Examine the catheter to ensure that the Abre stent is fully contained within the retractable sheath (silver).
Caution: Do not use the stent system if the Abre stent is partially deployed before starting the procedure.

Stent deployment

1. Insert the introducer sheath and the guidewire.
 - a. Gain access at the appropriate site, using a compatible introducer sheath with a hemostatic valve.
Caution: Gaining access too close to the intended treatment site can lead to difficulty with stent deployment.
 - b. Insert a guidewire of appropriate length across the target lesion via the introducer sheath.
Caution: Always use an introducer sheath during the implant procedure to protect both the vessel and puncture site.
2. Predilate the lesion with a balloon catheter using standard techniques. Remove the balloon from the patient while maintaining access with the guidewire.
3. Introduce the Abre™ stent delivery system over the guidewire through the hemostatic valve and introducer sheath.

Warning: If unusually high resistance is encountered, when advancing the Abre™ delivery system over the guidewire, assess the cause of the resistance before proceeding.

Caution: Ensure that the locking pin remains in locked position until the target site is reached.

4. Advance the Abre delivery system until the leading edge of the stent is beyond the target lesion.
5. Ensure that the isolation sheath (blue) remains inside the introducer sheath valve.
Caution: If the retractable sheath (silver) is exposed beyond the hemostatic valve, friction between the retractable sheath and the hemostatic valve can impair deployment. To avoid friction between the retractable sheath (silver) and the hemostatic valve, pull the introducer sheath back so that the hemostatic valve covers the end of the isolation sheath.
6. Remove the locking pin by holding the device handle with one hand and gently pulling the locking pin with the other hand.
7. While holding the device with one hand, use the other hand to gently support and position the isolation sheath at the hemostatic valve (see Figure 2). Do not constrict the retractable sheath during deployment.


Figure 2. Hand placement for stent deployment

8. Initiate Abre stent deployment by slowly rotating the thumbwheel in the direction indicated by the arrow on the handle.
9. After the leading end of the stent has emerged from the retractable sheath, but before achieving vessel apposition, reposition the Abre stent as needed.

Warning: If high resistance is felt when initially rotating the thumbwheel, do not force deployment. Carefully withdraw the system and do not use it.

Caution: Do not reposition the Abre stent after establishing apposition against the vessel wall.

Caution: The Abre delivery system is not designed for recapturing the stent.

10. Continue to rotate the thumbwheel until the stent is fully deployed.

Caution: The Abre stent is not designed to be lengthened or shortened from its nominal length. Excessive stent lengthening or shortening can increase the risk of stent fracture.

Note: If placement of sequential stents is necessary, ensure that there is sufficient overlap.

Post stent deployment

1. Withdraw the entire delivery system from the patient.
Caution: If resistance is felt when withdrawing the system, do not force withdrawal.
2. Perform post-deployment balloon dilation as needed, using an appropriately sized balloon catheter with conventional dilation techniques.
Caution: Use care when crossing a deployed Abre stent with any adjunct device to avoid dislodging the stent or damaging the adjunct device.
3. Confirm that the Abre stent is fully expanded; then remove the balloon catheter from the patient.
Note: Always visualize the stent by using imaging techniques to verify full vessel wall apposition.
4. Remove the guidewire and introducer sheath from the patient and discard the delivery system, guidewire, and introducer sheath.

MRI safety information

MR Conditional

Nonclinical testing demonstrated that the Abre stent in single and overlapped conditions is MR Conditional for stents up to 150 mm. A patient with this device can be scanned safely, immediately after stent placement, under the following conditions:

- Static magnetic field of 1.5 Tesla or 3.0 Tesla
- Maximum spatial gradient magnetic field of 4000 Gauss/cm or

less (40 T/m)

- Maximum MR system-reported, whole-body-averaged specific absorption rate (SAR) of 2.0 W/kg (normal operating mode)

MR image quality may be compromised if the area of interest is in the exact location or close to the position of the Abre stent.

MRI-related temperature rise

Under the scan conditions as defined in the previous section, "MR Conditional," the Abre stent is expected to produce a maximum temperature rise less than or equal to 5.2 °C after 15 minutes of continuous scanning (per pulse sequence). The effect of temperature rise in the MRI environment for stents with fractured struts is not known.


It is recommended that patients register conditions under which the implant can be scanned safely with the MedicAlert Foundation (www.medicalert.org) or equivalent organization.

Artifact information

In nonclinical testing, the maximum artifact size as seen on the gradient echo pulse sequence at 3.0 Tesla extends approximately 5 mm, relative to the size and shape of the Abre™ stent. The lumen of the stent can be visualized using the T1-weighted spin echo pulse sequence and the T1-weighted gradient echo pulse sequence at 3.0 Tesla.

Warranty

This product is subject to standard Medtronic warranty terms. For standard terms and conditions of supply for Australia and New Zealand, go to www.medtronic.com.

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© 2018 Medtronic. All rights reserved.
Medtronic and Medtronic logo are trademarks of Medtronic. All other brands are trademarks of a Medtronic company.
Medtronic Australasia Pty Ltd
Medtronic, Inc.