

PRECISE RESULTS

MINI TORIC allows precise correction of astigmatism.

- The Mini Toric Calculator software, developed on the basis of a vectorial calculation algorithm, ensures elevated **reliability of results**.

Pre-operative planning

Refractive surprise analysis

Surgeon

Patient

Pre-op Parameters

Implant Parameters (IOL plane)

Print

Clear all

Flat

Step

Step Meridian

Keratometry (anterior)

42.00

43.50

@

110

D

mm

Corneal Refractive Index

1.3375

Use total corneal power (e.g., CorT total)

Axial Length

24.00

mm

Phaco Incision Flattening

0.20

D

Phaco Incision Meridian

90

°

Phaco Target

1.35

@

113

IOL Type

MINI TORIC Ready

Details ...

Sph Equiv

20.00

Cylinder

2.00

Planned +ve Axis

113

IOL Power

Sphere

18.00

SRK/T A Constant: 119

Effective lens position: 5.27 mm

Set lens constant

Order IOL ...

OD

OS

ASSORT

Phaco Incision

Phaco Flattening

Phaco Target

IOL Power

Expected Refraction

MONOFOCAL IOL FOR THE SURGICAL TREATMENT OF ASTIGMATISM

SAFETY

MINI TORIC offers surgeons very high intra and post-operative safety.

High sterility

Lower risk of IOL damaging

PRELOADED INJECTION SYSTEM FOR MINI INCISION

Reduction of the induced astigmatism

Reduced surgery time

INSTRUCTION FOR USE:

- Insert the cartridge, containing the IOL, into the injector.
- Apply viscoelastic solution first through the front hole (until it is filled) and then a little through the rear hole.
- Remove the clip of safety block by lifting it using your index and thumb fingers.
- Close the cartridge sides together until the "click-lock" mechanism engages. Push forward the plunger softly and ensure that the silicone cushion correctly enters into the cartridge loading chamber.

TECHNICAL SPECIFICATIONS

Material	Copolymer	Cylindrical diopter powers (D) on IOL optics	Cylindrical diopter powers (D) on cornea
Positioning	Bag		
Total diameter (mm)	10.75		
Optics diameter (mm)	6		
Vaulting	5°	1,5	1
Optics shape	Biconvex aspheric (front) toric aspheric (back)	2	1,3
Posterior edge	Double	2,5	1,7
Estimated A Constant	118.6	3	2
Estimated A.C.D. (mm)	5.32	3,5	2,3
Spherical dioptric range (D)	OD to +30D (OD to +10D incr. 1D; +10D to +30D incr. 0.5D)	4	2,7
		4,5	3

Bibliography 1 Mukherjee R. et al., Micron 43 (9) 2012, 937-94
2 Anastasi E. et al., Poster presentato a XXXIV European Society of Cataract and Refractive Surgery Copenhagen, Denmark 2016, FP 24606
3 Chang F., J Cataract Refract Surg. 2003 May 29 (5): 935-40
4 Ruhswurm I. et al., J Cataract Refract Surg. 2000 Jul 26 (7): 1022-1027

medicalinformation@sifigroup.com
www.sifigroup.com

MONOFOCAL IOL FOR THE SURGICAL TREATMENT OF ASTIGMATISM

THE BEST COMBINATION BETWEEN
SAFETY AND STABILITY

ROTATIONAL STABILITY

QUALITY OF VISION


PRECISE RESULTS

SAFETY

www.minitoriccalculator.com

ROTATIONAL STABILITY


MINI TORIC has a design which ensures the IOL **rotational stability** over time.




- Micro-roughness and the 120° contact of the 4 haptics support the anchorage to the capsular bag (**gripping ability**) reducing the risk of rotation.

- 5° vaulting loops and **constant thickness of the optics edge** synergically contribute for a better adherence of the IOL to the posterior capsule.

- The symmetrical design of the loops facilitates both the **self-centering** (evidenced by the technological platform of the MINI IOLs) and **360° rotation in either direction** simplifying the IOL alignment on the target axis.




C-loop IOLs
can be rotated clockwise only




IOLs with 4 fenestrated haptics,
like MINI TORIC,
can be rotated clockwise
and counterclockwise easily
and in extreme safety

A trial on **MINI TORIC** carried out on 12 patients, 2 months after the implant, showed:

- excellent rotational stability over time²**



Postoperative absolute rotation <5° in 92% of **MINI TORIC** implanted



2° MINI TORIC
mean rotation is not significant^{3,4}

- centration stability** in the capsular bag due to the 4 haptics design²




Slit lamp examination

PROVED ROTATIONAL STABILITY


QUALITY OF VISION²

MINI TORIC is designed to ensure **high quality of vision**.




- Anti-scattering optics** edge avoids the secondary light rays divert on the retina creating effects of visual confusion.

SCATTERING PHENOMENA



ANTI-SCATTERING OPTICS



Optical simulations performed by Zemax software