№ 06/22


Date: 19.01.2022 p.

STATEMENT

We, Vitrotest Bioreagent LLC, having a registered office at M. Boychuka street 18b/56, Kyiv 01103 Ukraine, assign SRL SANMEDICO having a registered office at A. Corobceanu street 7A, apt. 9, Chişinău MD-2012, Moldova, as authorized representative.

We declare that the company mentioned above is authorized to register, notify, renew or modify the registration of medical devices on the territory of the Republic of Moldova.

Signature: _____ Director Ihor Nikolaienko Ph.D

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DECLARATION OF CONFORMITY №UA-TK039

Manufacturer: Vitrotest Bioreagent LLC
State registration № 42149820

Legal address: M.Boychuka 18b, of.56, Kyiv 01103 Ukraine
Manufacturer's address: Kurortnaya 11, Kyiv, 04075, Ukraine

Description of the product:

Name	Catalog Number
ELISA test-kit for the qualitative and semiquantitative determination of IgG antibodies to SARS-CoV-2 «Vitrotest SARS-CoV-2 IgG»	TK039

Classification:

According to medical devices technical regulation for *in vitro* diagnostics, approved by Cabinet of Ministers decree from 02.10.2013 №754

Is not a part of A and B lists, is not a device for self-testing, not for performance assessment.

Conformity assessment procedure:

Annex 3 of medical devices technical regulation for *in vitro* diagnostics, approved by Cabinet of Ministers decree from 02.10.2013 №754

Vitrotest Bioreagent declares the execution of all demands regarding the device, that was mentioned above, according to medical devices technical regulation for *in vitro* diagnostics, approved by Cabinet of Ministers decree from 02.10.2013 №754, and the requirements of further regulations:

ДСТУ EN ISO 13485:2018
ДСТУ EN ISO 14971:2015
ДСТУ EN 13641:2015
ДСТУ EN ISO 15223-1:2018 (EN ISO 15223-1:2016, IDT; ISO 15223-1:2016, Corrected version 2017-03, IDT)
ДСТУ EN ISO 23640:2015 (EN ISO 23640:2015, IDT; ISO 23640:2011, IDT)
ДСТУ EN 13612:2015
ДСТУ EN ISO 18113-1:2018 (EN ISO 18113-1:2011, IDT; ISO 18113-1:2009, IDT)
ДСТУ EN ISO 18113-2:2018 (EN ISO 18113-2:2011, IDT; ISO 18113-2:2009, IDT)
ДСТУ EN 980:2007

The declaration is made under sole responsibility of the manufacturer.

Date of issue: 17.04.2020

Validity of declaration till: 17.04.2025

Director

(signature)

Ihor Nikolaienko, Ph.D.
(ПІП)

Edition 1 from 17.04.2020

DECLARATION OF CONFORMITY №UA-TK033

Manufacturer: Vitrotest Bioreagent LLC
State registration № 42149820

Legal address: M.Boychuka 18b, of.56, Kyiv 01103 Ukraine
Manufacturer's address: Kurortnaya 11, Kyiv, 04075, Ukraine

Description of the product:

Name	Catalog Number
ELISA test-kit for the qualitative and semiquantitative determination of IgA antibodies to SARS-CoV-2 «Vitrotest SARS-CoV-2 IgA»	TK033
ELISA test-kit for the qualitative determination of IgM antibodies to SARS-CoV-2 «Vitrotest SARS-CoV-2 IgM»	TK034, TK042

Classification:

According to medical devices technical regulation for *in vitro* diagnostics, approved by Cabinet of Ministers decree from 02.10.2013 №754

Is not a part of A and B lists, is not a device for self-testing, not for performance assessment.

Conformity assessment procedure:

Annex 3 of medical devices technical regulation for *in vitro* diagnostics, approved by Cabinet of Ministers decree from 02.10.2013 №754

Vitrotest Bioreagent declares the execution of all demands regarding the device, that was mentioned above, according to medical devices technical regulation for *in vitro* diagnostics, approved by Cabinet of Ministers decree from 02.10.2013 №754, and the requirements of further regulations:

ДСТУ EN ISO 13485:2018
ДСТУ EN ISO 14971:2015
ДСТУ EN 13641:2015
ДСТУ EN ISO 15223-1:2018 (EN ISO 15223-1:2016, IDT; ISO 15223-1:2016, Corrected version 2017-03, IDT)
ДСТУ EN ISO 23640:2015 (EN ISO 23640:2015, IDT; ISO 23640:2011, IDT)
ДСТУ EN 13612:2015
ДСТУ EN ISO 18113-1:2018 (EN ISO 18113-1:2011, IDT; ISO 18113-1:2009, IDT)
ДСТУ EN ISO 18113-2:2018 (EN ISO 18113-2:2011, IDT; ISO 18113-2:2009, IDT)
ДСТУ EN 980:2007

The declaration is made under sole responsibility of the manufacturer.

Date of issue: 20.11.2020

Validity of declaration till: 20.11.2025

Director

Ihor Nikolaienko, Ph.D.
(п/п)

Vitrotest® SARS-CoV-2 IgG

ELISA test-kit for the qualitative and semiquantitative determination of IgG antibodies to SARS-CoV-2

TK039

96 tests

1. INTENDED USE

The test-kit «Vitrotest SARS-CoV-2 IgG» is an enzyme linked immunosorbent assay (ELISA) for the qualitative and semiquantitative determination of IgG antibodies to SARS-CoV-2 coronavirus in human serum or plasma. The test-kit might be applied for the ELISA using both automatic pipettes and standard equipment as well as open system automated ELISA analyzers.

2. CLINICAL VALUE

COVID-19 is an infectious disease caused by a new SARS-CoV-2 coronavirus which had not previously been detected in humans.

The viral infection leads to the development of a respiratory flu-like disease with symptoms such as cough and fever. In more severe cases pneumonia can develop. The average incubation period of the COVID-19 is 6,5 days, but it can range from 3 to 21 days.

Coronavirus SARS-CoV-2 is an RNA-containing virus with a characteristic envelope with spikes in the form of a "corona". The main structural proteins of the virus include the nucleocapsid protein and transmembrane S (spike) protein with a receptor-binding domain (RBD), which binds to human cell receptors ACE2, causing infection of the mucosal epithelial cells. Both proteins are highly immunogenic antigens for humans.

Specific antibodies (IgM, IgG, IgA) against the virus appear 7-11th days from the moment of infection and/or contact of the body with the virus. Anti-SARS-CoV-2 IgM antibodies can be detected beginning from day 4 of the first symptoms of the disease. Specific antibodies are detected in less than 40% of patients with COVID-19 during the first week of the onset of symptoms. The levels of these antibodies increase rapidly and are detected in almost all patients by the end of the second week of clinical manifestations. On the 20th day, from the onset of the symptoms of the COVID-19, specific IgG are detected in almost 100% of patients (except for individuals with immunosuppression).

According to some scientific data, there is a clear correlation between the severity of the disease and the level of specific IgG in the blood in the second week of clinical manifestations.

3. PRINCIPLE OF THE TEST

Determination of IgG antibodies specific to SARS-CoV-2 in the test-kit «Vitrotest SARS-CoV-2 IgG» is based on a solid phase, indirect ELISA in a two-step incubation procedure. Microwells are coated with the nucleocapsid recombinant antigens of SARS-CoV-2. During the first incubation step, the specific antibodies to SARS-CoV-2, if present in the sample, are bound to the solid phase precoated antigens. The wells are washed to remove unbound antibodies. A secondary antibodies (anti-IgG), which are conjugated to horseradish peroxidase (HRP), are added next and bind to the immune complexes on the solid phase. Unbound components are removed by washing. Antigen-antibody complexes are revealed by addition of chromogen solution containing 3,3',5,5'- tetramethylbenzidine (TMB) and hydrogen peroxide. After 30 minutes the reaction has been stopped, the absorbance values are read using a spectrophotometer at 450/620-695 nm. The colour intensity is proportional to the amount of the specific antibodies present in the sample.

4. MATERIALS AND EQUIPMENT

4.1. Composition of the test-kit

ELISA STRIPS	1x96 wells	Microplate (12 strips x 8 wells) Each well is coated with the recombinant antigens SARS-CoV-2. The wells can be separated.
PREDILUTION PLATE	1x96 wells	Microplate for preliminary dilution of sera
CONTROL +	1x0.3 ml	Positive control Solution of conjugated specific antibodies with preservative (pink).
CONTROL -	1x0.5 ml	Negative control Albumin solution with preservative (yellow).
SAMPLE PREDILUENT	1x20 ml	Sample prediluent Buffer solution with detergent and preservative (brown-green).
SAMPLE DILUENT	1x12 ml	Sample diluent Buffer solution with detergent and preservative (yellow).

CONJUGATE SOLUTION	1x12 ml	Conjugate solution (ready to use) Buffer solution of monoclonal antibodies to human IgG conjugated to HRP with stabilizers and preservative (violet).
TMB SOLUTION	1x12 ml	TMB solution (ready to use) TMB, H ₂ O ₂ , stabilizers and preservative (colourless).
WASH TWEEN 20X	1x50 ml	Washing solution Tw20 (20x concentrated) 20X concentrate of PBS buffer with Tween-20 and NaCl (colourless)
STOP SOLUTION	1x12 ml	Stop Solution (ready to use) 0.5 M H ₂ SO ₄ (colourless).

Adhesive films (2), sera identification plan (1) and instruction for use.

4.2. Material required but not provided

- Variable volume automatic pipettes (10µl–1000µl) and disposable pipette tips;
- plate reader (single wavelength 450 nm or dual wavelength 450/620–695 nm);
- volumetric laboratory glassware (10–1000ml);
- distilled/DI water;
- incubator thermostatically controlled at 37°C;
- automatic/semiautomatic plate washer;
- appropriate waste containers for potentially contaminated materials;
- timer;
- absorbent paper;
- disposable gloves;
- disinfectants;
- protective clothes.

5. PRECAUTIONS AND SAFETY

5.1. Precautions

The ELISA assays are time and temperature sensitive. Strictly follow the test procedure and do not modify it.

- do not use expired reagents;
- do not use for analyses and do not mix reagents from different lots or from test-kits of different nosology as well as other manufacturer's reagents with Vitrotest® kits;

*Note: it is possible to use **WASH TWEEN 20X**, **TMB SOLUTION**, **SAMPLE PREDILUENT** and **STOP SOLUTION** from other Vitrotest® ELISA kits.*

- close reagents after use only with appropriate caps;
- control the filling and full aspiration of the solution in the wells;
- use a new tip for each sample and reagent;
- avoid exposure of kit reagents to direct sunlight;
- **TMB SOLUTION** must be colourless before use. If **TMB SOLUTION** is blue or yellow it cannot be used. Avoid any contact of **TMB SOLUTION** with metals or metal ions. Use glassware thoroughly washed and rinsed with distilled/DI water;
- never use the same glassware for **CONJUGATE SOLUTION** and **TMB SOLUTION**.

The manufacturer is not responsible or liable for any incorrect results and/or incidents taking place as a result of any violation of the instruction. The manufacturer is not responsible for visual readings of samples (without using a plate reader).

5.2. Safety

- all reagents included in the kit are intended for in vitro diagnostic use only;
- the test-kit is designed for use by qualified personnel only;
- protective clothes, disposable gloves and safety glasses must be worn at all times while performing analysis;
- never eat, drink, smoke or apply cosmetics in the assay laboratory;
- never pipette solutions by mouth;
- positive and negative controls do not contain components of human origin;
- the liquid waste must be inactivated, for example, with hydrogen peroxide solution at the final concentration of 6% for 3 hours at room temperature, or with sodium hypochlorite at the final concentration of 5% for 30 minutes, or with other approved disinfectants;
- the solid waste must be inactivated by autoclaving at 121°C for 1 hour;
- dispose of inactivated waste in accordance of national laws and regulations;
- do not autoclave the solutions that contain sodium azide or sodium hypochlorite;

- some components of the test-kit contain low concentrations of harmful compounds and could cause irritation of the skin and the mucosa. In the case of contact of TMB SOLUTION STOP SOLUTION or CONJUGATE SOLUTION with skin or mucosa, the place of contact should be immediately rinsed with large amounts of water;
- in case of spilling of solutions that do not contain acid, e.g. sera, rinse the surface with disinfectant, then dry it with absorbent paper. In other case acid first must be neutralized by sodium bicarbonate and then wiped out as described above.

6. STORAGE AND STABILITY

Reagents are stable until stated expiration date on the label when stored refrigerated (2-8°C). Do not freeze. The kit should be shipped at 2-8°C. Single transportation at the temperature up to 23°C for two days is acceptable.

7. SPECIMEN COLLECTION

The fresh serum or plasma samples can be stored for 3 days at 2-8 °C, or frozen for longer periods at -20 – -70°C. Frozen samples must be thawed and kept at room temperature for at least 30 minutes before use. Do not use preheated samples. Mix thawed samples thoroughly to homogeneity. Avoid repeated freezing/thawing. Samples containing aggregates must be clarified by centrifugation (3000 rpm for 10-15 min). Do not use hyperlipemic, hyperhaemolysed or contaminated by microorganisms serum specimens. The presence of bilirubin up to concentration of 0.21 mg/ml (361.8 µM/l), haemoglobin up to concentration of 10 mg/ml and triglycerides up to concentration of 10 mg/ml (11.3 mM/l) are allowed.

8. REAGENT PREPARATION

It is very important to keep all test components for at least 30 min at room temperature (18-25 °C) before the assay!

8.1. ELISA STRIPS preparation

Before opening the bag with ELISA STRIPS, keep it at room temperature for 30 minutes to avoid water condensation inside the wells. Open the vacuum bag and take out the necessary number of the wells. Once opened the bag with the remaining strips must be resealed with zip-lock immediately and kept refrigerated at 2-8°C for no more than 3 months.

8.2. Washing solution preparation

Check the WASH TWEEN|20X for the presence of salt crystals. If crystals have formed, re-solubilise by warming at 37°C, until crystals dissolve (15-20min). Dilute the WASH TWEEN|20X 1:20 (1+19) with distilled/DI water before use. For example, 4 ml concentrate + 76 ml water is sufficient for 8 wells. Once diluted it is stable at 2-8°C for 1 week.

8.3. Predilution of samples and controls

Predilute patient samples and controls 1:10 with SAMPLE PREDILUENT immediately before test. Dispense 90 µl of SAMPLE PREDILUENT in the wells of PREDILUTION PLATE, add 10 µl of samples and controls. Gently mix the content in the wells. After addition of the sample colour of the sample prediluent changes from brown-green to blue.

After dilution and transfer of the samples, the used PREDILUTION PLATE wells must be disinfected by soaking in a disinfectant solution or by autoclaving.

The procedure for dilution of samples and controls should be carried out immediately before analysis.

9. ASSAY PROCEDURE

1. Prepare the necessary number (the number of specimens and controls) of the wells ELISA STRIPS (see 8.1). Place the wells into the frame. Wells with the controls must be included in every test.
2. Complete the sera identification plan.
3. Prepare washing solution (see 8.2).
4. Predilute patient samples and controls (see 8.3).
5. Dispense 90 µl of SAMPLE DILUENT into each well.
6. Add 10 µl of prediluted controls and patient samples to the wells in the following order: A1 – CONTROL|+, B1, C1 and D1 – CONTROL|-, other wells – patient samples. The final dilution in the wells is 1:100. Pipette gently to avoid foaming. The colour of the sample diluent changes from yellow to green.
7. Cover strips with an adhesive film and incubate for 30 min at 37°C.
8. Remove and discard the adhesive film and wash all wells 5 times with automatic washer or 8-channel pipette as follows:
 - aspirate the contents of all wells into a liquid waste container and add immediately a minimum of 300 µl of diluted Washing solution to each well;
 - soak each well for 30 seconds between each wash cycle;
 - aspirate again. The residual volume must be lower than 5µl;
 - repeat the washing step 4 times;
 - after the final washing cycle, turn down the plate onto an absorbent paper and tap it to remove any residual buffer.

- 9.9. Dispense 100 µl of [CONJUGATE SOLUTION] per well. Cover strips with a new adhesive film, incubate for 30 min at 37°C.
- 9.10. Remove and discard the adhesive film and wash all wells five times as described above (see 9.8).
- 9.11. Dispense 100 µl [TMB SOLUTION] into all wells. Do not touch the walls and bottoms of the wells to avoid contamination.
- 9.12. Incubate the strips for 30 minutes at room temperature (18-25°C) in the dark. Do not use adhesive film in this step.
- 9.13. Add 100 µl [STOP SOLUTION] to each well in the same order and at the same rate as for [TMB SOLUTION].
- 9.14. Read the optical density (OD) of the wells at 450/620-695 nm using a microplate reader within 5 minutes after adding the [STOP SOLUTION]. Pay attention to the cleanness of the plate bottom and absence of bubbles in the wells before reading.

Measurement in the single-wave procedure at 450 nm is possible. Reserve blank well to adjust spectrophotometer in such analysis. Only [TMB SOLUTION] and [STOP SOLUTION] must be added in blank well.

10. CALCULATION AND INTERPRETATION OF RESULTS

10.1. Calculation of results

Calculate the mean absorbance value for 3 negative controls (Nc), Cut off value (CO) and Sample Index of Positivity (IP_{sample}),

$$Nc = (Nc1 + Nc2 + Nc3)/3; \quad CO = Nc + 0.3; \quad IP_{sample} = OD_{sample}/CO$$

10.2. Validation of the test

The test run may be considered valid provided the following criteria are met:

[CONTROL +]	OD ≥ 1.200
[CONTROL -]	OD ≤ 0.150
[CONTROL -]	$Nc \times 0.5 \leq Ncn \leq Nc \times 2.0$

If one of the negative control absorbances does not match the above criteria, this value should be discarded and a mean value should be calculated using the other two values. If more than one negative control absorbance does not meet the criteria, the test is invalid and must be re-tested.

10.3. Interpretation of results

IP _{sample} > 1.1	POSITIVE
$0.9 \leq IP_{sample} \leq 1.1$	DOUBTFUL*
IP _{sample} < 0.9	NEGATIVE

* If the result is doubtful, repeat the test. If it remains doubtful, collect a new serum sample.

Using IP_{sample} allows performing semiquantitative comparative analysis of the level of specific antibodies in dynamics in paired sera obtained from patients with an interval of 2-4 weeks.

11. PERFORMANCE CHARACTERISTICS

11.1. Specificity and sensitivity

The specificity of the test-kit «Vitrotest SARS-CoV-2 IgG» was evaluated using 304 samples of human serum obtained during 1-6 months in 2019 (before the COVID-19 pandemic). IgG to SARS-CoV-2 was not detected in all tested samples, the specificity was 100%.

For sensitivity determinations 27 serum samples of patients who have recovered from COVID-19 infection were evaluated. IgG antibodies to SARS-CoV-2 were detected in all tested sera in the «Vitrotest SARS-CoV-2 IgG» test-kit.

11.2. Accuracy

Intra assay repeatability

Coefficient of variation (CV) was calculated by measuring 2 samples with various specific antibody levels in 32-replicate determinations using 1 lot of the test-kit.

Serum No.	OD _{mean}	IP _{mean}	CV, %
27	0.966	2.99	5.4
30	2.684	8.31	4.6

Inter assay reproducibility

Coefficient of variation (CV) was calculated by measuring 2 samples with various specific antibody levels in 4 ELISA performances during 4 days, in 8-replicate determinations.

Serum No.	OD _{mean}	IP _{mean}	CV, %
27	0.937	2.89	5.4
30	2.642	8.14	4.0

12. LIMITATIONS OF THE PROCEDURE

If the test sample was obtained in the first days after infection IgG antibodies may not be detected. Therefore, a negative result of the test for IgG antibodies to SARS-CoV-2 does not exclude the infection of the patient with the virus. In the presence of clinical manifestations of the disease it is recommended to repeat the test at least in two weeks.

Negative test results in immunosuppressed individuals should also be interpreted with caution.

The final diagnosis cannot be established solely on the basis of serological test results. The diagnosis should take into consideration clinical history, symptomatology, as well as the results of other laboratory tests (including PCR).

13. TROUBLESHOOTING

Possible causes	Solutions
<i>High background in all wells</i>	
Contaminated washer	Clean the washer head, then rinse it with 30% ethanol and distilled water
Low quality water or contaminated water	Use distilled/DI with resistivity ≥ 10 M Ω -cm.
Using contaminated glassware	Use clean glassware
Using chlorine based disinfectants	Use disinfectants without chlorine
Using contaminated tips	Use new tips
Increased time of incubation or temperature regimen was changed	Follow incubation regimen according to instruction for use
<i>High background in a few wells</i>	
TMB solution was added more than once	Add TMB solution once
Pipette shaft was contaminated with conjugate solution	Clean the pipette; pipette the liquids carefully
One the channels of the washer was contaminated	Clean the washer channel, clean the washer
<i>OD of the positive control below normal</i>	
Conjugate solution/tmb solution was prepared improperly or not added	Run ELISA repeatedly, prepare conjugate solution / TMB solution properly
Reduced incubation time in one of the stages	Follow incubation regimen according to the instruction for use
<i>Visual colour intensity of the wells does not correspond to optical density</i>	
The optical beam or another component of the reader is misaligned or malfunctioning	Test the absorbance reader's performance

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1. Juanjuan Zhao Jr., Quan Yuan et. al. Antibody responses to SARS-CoV-2 in patients of novel coronavirus disease 2019. // Clinical Infectious Diseases., - 2020 Mar. 20 doi: 10.1093/cid/ciaa344.
2. Laboratory testing for coronavirus disease (COVID-19) in suspected human cases // WHO. Interim guidance 19 March 2020. WHO reference number: WHO/COVID-19/laboratory/2020.5.
3. Marco Cascella ; Michael Rajnik et.al. Features, Evaluation and Treatment Coronavirus (COVID-19). // NCBI Bookshelf. StatPearls Publishing; 2020 – P.16.
4. Patrick C. Y. Woo, Susanna K. P. Lau. et.al. Differential Sensitivities of Severe Acute Respiratory Syndrome (SARS) Coronavirus Spike Polypeptide Enzyme-Linked Immunosorbent Assay (ELISA) and SARS Coronavirus Nucleocapsid Protein ELISA for Serodiagnosis of SARS Coronavirus Pneumonia. // J. Clin. Microbiol., - 2005 – V. 43 N.7 - p. 3054–3058.
5. Quan-Xin Long, Bai-Zhong Luet et. al. Antibody responses to SARS-CoV-2 in patients with COVID-19. // Nature Medicine., - 2020 April 29 doi: 10.1038/s41591-020-0897-1.
6. Quan-xin Long, Hai-jun Deng et.al. Antibody responses to SARS-CoV-2 in COVID-19 patients: the perspective application of serological tests in clinical practice. // medRxiv preprint doi: <https://doi.org/10.1101/2020.03.18.20038018>.
7. Shu-Yuan Xiao, Yingjie Wu, Huan Liu. Evolving status of the 2019 novel coronavirus infection: Proposal of conventional serologic assays for disease diagnosis and infection monitoring. // J Med Virol., - 2020 - 92(5) – p.464-467. doi: 10.1002/jmv.25702. Epub 2020 Feb 17.
8. Wu, L.-P., Wang, N.-C. et.al. (2007). Duration of Antibody Responses after Severe Acute Respiratory Syndrome. // Emerging Infectious Diseases., - 2007 - 13(10) - p.1562-1564.

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For questions and suggestions regarding the kit, contact the manufacturer:

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Vitrotest® SARS-CoV-2 IgG

ASSAY PROCEDURE

Keep all reagents and specimens for at least 30 min at 18-25°C before use

Dispense 90 µl of **SAMPLE PREDILUENT** (*brown-green colour*) into the wells of **PREDILUTION PLATE** and add 10 µl of controls and the samples (*colour changes from brown-green to blue*)

Dispense 90µl of **SAMPLE DILUENT** (*yellow colour*) into the wells of **ELISA STRIPS** and 10µl prediluted samples and controls:

A1 – **CONTROL +**,

B1, C1, D1 – **CONTROL -**,

E1 and other wells – patient samples

(*colour changes from yellow to green*)

Cover wells with an adhesive film, incubate for 30 min at 37°C

Rinse the wells 5 times with diluted 1:20 (1+19) washing solution Tween-20 (300µl per well)

Add 100µl of **CONJUGATE SOLUTION** into the wells (*green colour*)

Cover wells with an adhesive film, incubate for 30 min at 37°C

Rinse the wells 5 times with diluted 1:20 (1+19) washing solution Tween-20 (300µl per well)

Add 100µl of **TMB SOLUTION** into the wells

Cover wells with an adhesive film, incubate for 30 min in the dark at 18-25°C

Add 100µl of **STOP SOLUTION** (*colour changes from blue to yellow*)

Determine the optical density (OD) at 450/620-695nm

CALCULATION

$$Nc = (Nc1 + Nc2 + Nc3)/3;$$

$$CO = Nc + 0.3;$$

$$IP_{\text{sample}} = OD_{\text{sample}} / CO;$$

$$Nc - OD_{\text{mean}} \text{ for 3 } \mathbf{CONTROL -}$$

CO - Cut off, IP- Index of Positivity

INTERPRETATION

$IP_{\text{sample}} > 1.1$	POSITIVE
$0.9 \leq IP_{\text{sample}} \leq 1.1$	DOUBTFUL
$IP_{\text{sample}} < 0.9$	NEGATIVE