**SYPHILIS SEROLOGY KIT**  
**DIRECTIONS FOR USE**

**RPR CARBON KIT: For Detection Of Syphilis.**

**SUMMARY**

At one time, syphilis was a major medical disease with a host of different manifestations transmitted primarily through sexual contact. The advent of penicillin in 1943 changed this. The etiologic agent of syphilis is *Treponema pallidum*, a spiral bacterium (spirochete). The spirochete causes some damage to the heart and the liver, releasing some tissue fragments. The patient's immune system produces antibodies, called reagins, against these fragments. There are two different techniques for the detection of syphilis. TPHA tests, which detect antibodies to *Treponema pallidum*, and non-treponemal serologic tests, which detect Reagin in infected people.

**INTENDED PURPOSE**

The reagent is a test reagent intended to be used to qualitatively and semi-quantitatively determine the presence or absence of Reagin (antibodies against Syphilis) in the serum or plasma of patients when tested in accordance with the recommended techniques stated in this IFU.

**PRINCIPLE**

When used by the recommended techniques, the reagent will agglutinate (clump) in the presence of reagin. No agglutination usually indicates the absence of reagin (see **Limitations**).

**KIT DESCRIPTION**

Lorne RPR Carbon Kit is a non-treponemal serologic test for the detection of syphilis. The RPR Carbon Antigen contains micro particulate carbon, which aids in the microscopic reading of results. The reagents do not contain or consist of CMR substances, or endocrine disrupting substances or that could result in sensitisation or an allergic reaction by the user. All the reagents are supplied at optimum dilution for use with all recommended techniques without the need for further dilution or addition. For lot reference number and expiry date see **Vial Labels**.

**STORAGE**

Do not freeze. Reagent vials should be stored at 2 - 8°C on receipt. Prolonged storage at temperatures outside this range may result in accelerated loss of reagent reactivity.

**SPECIMEN COLLECTION**

Specimens should be drawn with or without anticoagulant using an aseptic phlebotomy technique. If testing is delayed specimens can be stored at 2-8°C for 7 days or for up to 3 months at or below -20°C. Specimens must be free from bacterial contamination, fibrin, haemolysis and lipaemia.

**PRECAUTIONS**

1. The kit is for *in vitro* diagnostic use only.
2. Do not use kit past expiration date (see **Vial and Box Labels**).
3. Protective clothing should be worn when handling the reagents, such as disposable gloves and a laboratory coat.
4. The reagents in this kit have been processed to reduce the bio-burden, but are not supplied sterile. Once a vial has been opened the contents should remain viable up until the expiry date.
5. No known tests can guarantee products derived from human or animal sources are free from infectious agents. Care must be taken in the use and disposal of each vial and its contents.
6. RPR Positive Control: H319 - Causes serious eye irritation. Follow the precautionary statement given in the SDS.

**DISPOSAL OF KIT REAGENT AND DEALING WITH SPILLAGES**

For information on disposal of kit reagent and decontamination of a spillage site see **Material Safety Data Sheets**, available on request.

**CONTROLS AND ADVICE**

1. It is recommended that the RPR Positive and Negative Controls are tested in parallel with each batch of tests. Tests must be considered invalid if controls do not show expected results.
2. Shake all the reagents well before use to ensure homogeneity.
3. Do not interchange components between different kits.
4. All the reagents must be allowed to reach 18-25°C before use.
5. The circles on the agglutination cards should never be touched with fingers, as this may invalidate the test results.
6. Use of kit and interpretation of results must be carried out by properly trained and qualified personnel in accordance with the requirements of country where reagents are in use.
7. The user must determine suitability of the kit for use in other techniques.

**KIT COMPONENTS PROVIDED**

- 1) RPR Carbon Reagent (White cap, 1x3 mL (150 tests) or 2x5 mL (500 tests)): Carbon particles coated with a lipid complex (cardiolipin, lecithin and cholesterol) in phosphate buffer 20 mmol/L, pH 7.0 containing a preservative.
- 2) RPR Positive Control (Red cap, 1 mL): Artificial serum with reagin titer  $\geq 1/4$ .
- 3) RPR Negative Control (Blue cap, 1 mL): Animal serum containing a preservative
- 4) Dispensing bottle (Green cap, 1 x 2 ml).
- 5) Dispensing Needle (x1).
- 6) Disposable agglutination slides.
- 7) Plastic stirrers.

**MATERIALS AND EQUIPMENT REQUIRED BUT NOT SUPPLIED**

- a) Pipette capable of accurately delivering 50  $\mu$ l
- b) Mechanical rotating table capable of rotating at 80-100 rpm.
- c) 9 g/L saline solution.

**QUALITATIVE TECHNIQUE**

1. Allow the reagents and samples to reach room temperature. The sensitivity of the test may be reduced at low temperatures.
2. Place 50  $\mu$ L of the sample and one drop of each Positive and Negative Controls into separate circles on the slide test.
3. Swirl the RPR Carbon Reagent gently before using. Invert the dropper assembly and press gently to remove air bubbles from the micropipette.
4. Place the micropipette in a vertical position and perpendicular to the slide, and add one drop (20  $\mu$ L) of this reagent next to the samples to be tested.
5. Mix the drops with a stirrer, spreading them over the entire surface of the circle. Use different stirrers for each sample
6. Place the slide on a mechanical rotating table at 80-100 r.p.m. for 8 min. False positive results could appear if the test is read after more than 8 minutes.

**INTERPRETATION OF QUALITATIVE RESULTS**

1. **Reactive:** Visible agglutination (medium to large clumps) constitutes a positive result and within the accepted limitations of the test procedure, indicates the presence of reagin.
2. **Weak-Reactive:** Weak agglutination (small clumps) around the periphery of the test area constitutes a weak positive result and within the accepted limitations of the test procedure, indicates the presence of reagin.
3. **Negative:** No agglutination constitutes a negative result and within the accepted limitations of the test procedure, indicates the absence of reagin.

**SEMI QUANTITATIVE TECHNIQUE**

1. The semi-quantitative test can be performed in the same way as the quantitative technique using dilutions of the serum in 9 g/L saline solution.
2. Make doubling dilutions of specimen as follows:

Dilution	Serum	Saline
1/2	100 µl undiluted serum	100 µl
1/4	100 µl 1/2 diluted serum	100 µl
1/8	100 µl 1/4 diluted serum	100 µl
1/16	100 µl 1/8 diluted serum	100 µl

**Lorne Laboratories Limited**  
Unit 1 Cutbush Park Industrial Estate  
Danehill  
Lower Earley  
Berkshire, RG6 4UT  
United Kingdom  
Tel: +44 (0) 118 921 2264  
Fax: +44 (0) 118 986 4518  
E-mail: info@lornelabs.com

- Test the specimen dilutions in the same way as for the quantitative technique above.
- Read the test and note the last positive dilution series.

#### STABILITY OF THE REACTIONS

Slide tests should be interpreted straight after the 8-minute rotating period to avoid the possibility that a negative result may be incorrectly interpreted as positive due to drying of the reagent.

#### LIMITATIONS

- RPR carbon test is non-specific for syphilis. All Reactive samples should be retested with treponemic methods such as TPHA and FTA-Abs to confirm the results.
- A Non Reactive result by itself does not exclude a diagnosis of syphilis. Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.
- False positive results have been reported in diseases such as infectious mononucleosis, viral pneumonia, toxoplasmosis, pregnancy and autoimmune diseases.
- Bilirubin ( $\leq 20$  mg/dL), hemoglobin ( $\leq 10$  g/L) and lipids ( $\leq 10$  g/L), do not interfere. Rheumatoid factors ( $\geq 300$  IU/mL), interfere. Other substances may interfere<sup>5</sup>.
- False positive or negative results may also occur due to:
  - Not expelling air from end of needle
  - Not maintaining dispensing bottle and needle in a vertical position when dispensing the antigen.
  - When transferring the specimen from the collecting tube some of the specimen being drawn up in to the teat
  - Contamination of test materials
  - Improper storage of test materials or omission of reagents
  - Deviation from the recommended techniques

#### SPECIFIC PERFORMANCE CHARACTERISTICS

- The kit has been characterised by all the procedures mentioned in the **Recommended Techniques**.
- Prior to release, each lot of Lorne RPR Syphilis Kit is tested by the **Recommended Techniques** to ensure suitable reactivity.
- The reagent sensitivity is calibrated against the WHO 1<sup>st</sup> International Standard for human syphilitic plasma (NIBSC reference number 05/132).
- Prozone effect:** No prozone effect was detected up to titers  $\geq 1/128$ .
- Diagnostic sensitivity:** 100%
- Diagnostic specificity:** 100 %.

#### DISCLAIMER

- The user is responsible for the performance of the kit by any method other than those mentioned in the **Recommended Techniques**.
- Any deviations should be validated prior to use using established laboratory procedures.

#### BIBLIOGRAPHY

- David S.Jacobs et al. Laboratory Test Handbook, 3<sup>rd</sup> edition, Lexi-Comp Inc, 1994.

#### AVAILABLE KIT SIZES

Kit Size	Catalogue Number
150 Tests Per Kit	044150A
500 Tests Per Kit	044500A

Advena Ltd. Tower Business Centre, 2<sup>nd</sup> Flr.,  
Tower Street, Swatar, BKR 4013, Malta

# Certificate

Standard **ISO 9001:2015**

Certificate Registr. No. 01 100 1810008

Certificate Holder: **MACHEREY-NAGEL GmbH & Co. KG**  
Neumann-Neander-Str. 6-8  
52355 Düren  
Germany

including the locations according to annex

Scope: Design and development, production and distribution  
of products for filtration, rapid tests, water analysis,  
chromatography and bioanalysis

Proof has been furnished by means of an audit that the  
requirements of ISO 9001:2015 are met.

Validity: The certificate is valid from 2018-08-29 until 2020-05-28.

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Standard **ISO 9001:2015**

Certificate Registr. No. 01 100 1810008

No.	Location	Scope
/01	MACHEREY-NAGEL GmbH & Co. KG Neumann-Neander-Str. 6-8 52355 Düren Germany	Design and development, production and distribution of products for filtration, rapid tests, water analysis, chromatography and bioanalysis
/02	MACHEREY-NAGEL GmbH & Co. KG Valenciener Str. 11 52355 Düren Germany	Design and development, production, distribution, service and administration
/03	MACHEREY-NAGEL GmbH & Co. KG Papiermühle 50 52349 Düren Germany	Waste disposal
/04	MACHEREY-NAGEL GmbH & Co. KG Bahnstr. 120 52355 Düren Germany	Storage
/05	MACHEREY-NAGEL GmbH & Co. KG Monschauer Str. 64 52355 Düren Germany	Production

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003285

## EC-Declaration of Conformity

According to Directive 98/79/EC on in-vitro-diagnostic devices, Annex III


Manufacturer: nal von minden GmbH, Carl-Zeiss Str.12, 47445 Moers  
 Classification: Other Products

We herewith declare on our sole responsibility that all batches of below mentioned In-vitro-diagnostic devices are conform with the Essential Requirements Annex I of the directive 98/79/EC of the European Parliament and of the Council of 27 October 1998 on in vitro diagnostic medical devices. The products are suitable for the intended application (only professional users).

Relevant standards and guidelines are applied.

1010002N-20	NADAL® Cholera O1/O139 Test
1010002N-20_SSL	NADAL® Cholera O1/O139 Test
1120003N-20	NADAL® Candida albicans Test
1130002N-30	NADAL® HSV 1 IgG/IgM Test
1130002N-30_SSL	NADAL® HSV 1 IgG/IgM Test
1130003N-30	NADAL® HSV 2 IgG/IgM Test
1130003N-30_SSL	NADAL® HSV 2 IgG/IgM Test
1200001	NADAL® Lactoferrin Test
1201004N-10	NADAL® Ferritin Test
1201004N-10_SSL	NADAL® Ferritin Test
1212001	NADAL® Calprotectin Test
1212001_SSL	NADAL® Calprotectin Test
1212002	NADAL® Claprotectin/Lactoferrin Test
1222001	NADAL® Enterovirus Test
1222001_SSL	NADAL® Enterovirus Test
1232001	NADAL® Campylobacter Test
1242001	NADAL® Salmonella spp. Test
1242001_SSL	NADAL® Salmonella spp. Test
1242002	NADAL® Salmonella typhi Test


1242002_SSL	NADAL® Salmonella typhi Test
1252001	NADAL® Listeria Test
1252001_SSL	NADAL® Listeria Test
1262001	NADAL® Shigella Test
1262001_SSL	NADAL® Shigella Test
1262002	NADAL® Shigella dysenteriae Test
1062002_SSL	NADAL® Shigella dysenteriae Test

161001	NADAL® hLH Ovulation Test
161001_SSL	NADAL® hLH Ovulation Test
162001	NADAL® hLH Ovulation Test
162001_SSL	NADAL® hLH Ovulation Test
164001	NADAL® hLH Ovulation Test
165001	NADAL® hLH Ovulation Test
165001_SSL	NADAL® hLH Ovulation Test
165003	NADAL® hLH Ovulation Test
165003_SSL	NADAL® hLH Ovulation Test
166001	NADAL® hLH Ovulation Test
172003N-10	NADAL® FSH Menopause Test
194002	NADAL® pH-Test
1941333	NADAL® Infectious Diseases Puffer
201001	NADAL® Syphilis Test
201001_SSL	NADAL® Syphilis Test
202001	NADAL® Syphilis Test
202001_SSL	NADAL® Syphilis Test
203001	NADAL® Syphilis Test
203002	NADAL® Syphilis Test
203002_SSL	NADAL® Syphilis Test
2090001	NADAL® Entamoeba Test
2090001_SSL	NADAL® Entamoeba Test
2201001N-10	NADAL® IgE Test
2210001N-20	NADAL® Ebola Test
221001A	NADAL® Strep A Test
221005	NADAL® Strep A Reagenz 1
221006	NADAL® Strep A Reagenz 2
221050N-50	NADAL® Strep A plus Test
221050N-50_SSL	NADAL® Strep A plus Test

222001A	NADAL® Strep A Test
222007	NADAL® Strep A plus Test
222008	NADAL® Strep A plus Test
222008_SSL	NADAL® Strep A plus Test
222011	NADAL® Strep A plus Test
222049NBUL-20	NADAL® Strep A Scan Test
232001	NADAL® Strep B Test
232001_SSL	NADAL® Strep B Test
232005	NADAL® Strep B Reagenz 1
232006	NADAL® Strep B Reagenz 2
241005N-10	NADAL® Influenza A+B Test
241006N-25	NADAL® Influenza A+B Test
242001	NADAL® Influenza A+B Test
242006N-10	NADAL® Influenza A/B Test
252001	NADAL® Mononucleosis Test
252003	NADAL® Mononucleosis Test
252003N-20	NADAL® Mononucleosis Test
252017N-05	NADAL® Mononucleosis Test
262001	NADAL® H.Pylori Ab Test
262001_SSL	NADAL® H.Pylori Ab Test
262002	NADAL® H.Pylori Ag Test
262004NBUL-10	NADAL® H. pylori Ag Scan Test
272001	NADAL® FOB Test
272001_SSL	NADAL® FOB Test
272007	NADAL® FOB & Hb/Hp Extraktionsröhrchen
272011N-25	NADAL® Hb/Hp Complex plus Test
272015	NADAL® Hb/Hp Complex Test
272015_SSL	NADAL® Hb/Hp Complex Test
272016	NADAL® FOB plus Test
272031RU	NADAL® Hb/Hp Complex Test
272031RU-01	NADAL® Hb/Hp Complex Test
272031RU_SSL	NADAL® Hb/Hp Complex Test
272035	NADAL® FOB60 plus Test
272035_SSL	NADAL® FOB60 plus Test


272037	NADAL® FOB & Hb/Hp patient set
272040	NADAL® FOB II Test
272041	NADAL® FOB60 Test
272041_SSL	NADAL® FOB60 Test
272042	NADAL® FOB75 Test
272043	NADAL® FOB75 II Tet
282000	NADAL® Troponin I Test
282001	NADAL® Troponin I Test
282001_SSL	NADAL® Troponin I Test
282003	NADAL® Cardiac Combo Test
282015	NADAL® Troponin I Test
292001N-05	NADAL® MyoglobinTest
302001	NADAL® CK-MB Test
302001_SSL	NADAL® CK-MB Test
311003N-10	NADAL® CrP Test
311003N-20	NADAL® CrP Test
311004	NADAL® CRP plus Test
311006	NADAL® CRP plus Test
312001	NADAL® CrP Test
312002	NADAL® High Sensitivity CrP Test
312017	NADAL® CRP Quant RFID Chip
312017_SSL	NADAL® CRP Quant RFID Chip
312021NBUL-20	NADAL® CRP Quant Test
312021NBUL-40	NADAL® CRP Quant Test
322003N-30	NADAL® Tuberkulose IgG/IgM Test
322003N-30_SSL	NADAL® Tuberkulose IgG/IgM Test
331001	NADAL® Microalbumin Test
331004N-50	NADAL® Mikroalbumin Test
351006	NADAL® D-Dimer Test
351006_SSL	NADAL® D-Dimer Test
351006N-25	NADAL® D-Dimer Test
351007	NADAL® D-Dimer cassette
431001N-03	NADAL® PROM Amniotic fluid Test
431001N-03_SSL	NADAL® PROM Amniotic fluid Test

431001N-10	NADAL® PROM Amniotic fluid Test
431001N-10_SSL	NADAL® PROM Amniotic fluid Test
431001N-20	NADAL® PROM Amniotic fluid Test
431006N-01	NADAL® PROM Amniotic fluid Test
431006N-03	NADAL® PROM Amniotic fluid Test
431006N-10	NADAL® PROM Amniotic fluid Test
431006N-20	NADAL® PROM Amniotic fluid Test
432000N-03	NADAL® fFN Test
432000N-10	NADAL® fFN Test
432000N-20	NADAL® fFN Test
472003N-10	NADAL® Malaria 4 Species Test
472003N-25	NADAL® Malaria 4 Species Test
472003N-25_SSL	NADAL® Malaria 4 Species Test
472008	NADAL® Malaria 4 species Test
472009	NADAL® Malaria Pf/Pv Ab Test
472030N-10	NADAL® Malaria Pf/Pan Ag 4 Species Test
472030N-25	NADAL® Malaria Pf/Pan Ag 4 Species Test
472036N-25	NADAL® Malaria Pf Ag 4 Species Test
472036N-25_SSL	NADAL® Malaria Pf Ag 4 Species Test
481008	NADAL® Adenovirus Respiratory Test
481013	NADAL® Adenovirus positive control
481015	NADAL® Rota-Adenovirus Test
481015_SSL	NADAL® Rota-Adenovirus Test
481015N-20	NADAL® Rota-Adenovirus Test
481016	NADAL® Adenovirus Test
481016_SSL	NADAL® Adenovirus Test
481017	NADAL® Rotavirus Test
481049NBUL-10	NADAL® Rota-Adenovirus Scan

	Test
491000N-10	NADAL® RSV Test
491000N-10_SSL	NADAL® RSV plus Test
491003N-25	NADAL® RSV Test
491003N-25_SSL	NADAL® RSV Test
491005	NADAL® RSV Test
491008NBUL-20	NADAL® RSV Scan Test
491009	NADAL® RSV-Adenovirus Respiratory Test
495001	NADAL® MRSA Latextest
495001_SSL	NADAL® MRSA Latextest
501006	NADAL® E.coli O157 Test
501006_SSL	NADAL® E.coli O157 Test
501012	NADAL® EHEC Verotoxin 1/2 Test
501012_SSL	NADAL® EHEC Verotoxin 1/2 Test
511006	NADAL® Cryptosporidium Test
511006_SSL	NADAL® Cryptosporidium Test
521001	NADAL® Giardia Test
521001_SSL	NADAL® Giardia Test
521009	NADAL® Giardia Test
521009_SSL	NADAL® Giardia Test
521010	NADAL® Crypto/Giardia Test
521010_SSL	NADAL® Crypto/Giardia Test
532001_SSL	NADAL® Dengue IgG/IgM Test
532001N-25	NADAL® Dengue IgG/IgM Test
532002N-25	NADAL® Dengue Ag Test
532003N-25	NADAL® Dengue Ag+IgG/IgM Test
532003N-25_SSL	NADAL® Dengue Ag+IgG/IgM Test
532004N-25	NADAL® Dengue IgG/IgM Test
532004N-25_SSL	NADAL® Dengue IgG/IgM Test
532011N-25	NADAL® Dengue IgG/IgM Test
532012N-25	NADAL® Dengue NS1 Ag Test
532016N-25	NADAL® Dengue NS1 Ag+IgG/IgM Test
532016N-25_SSL	NADAL® Dengue NS1 Ag+IgG/IgM Test

25_SSL	Test
542001N-25	NADAL® Tetanus Test
552005	NADAL® Legionella Test
552005_SSL	NADAL® Legionella Test
552006	NADAL® Legionella Test
552006_SSL	NADAL® Legionella Test
552020	NADAL® Legionella Test
552020_SSL	NADAL® Legionella Test
562003N-10	NADAL® BCA/Hb Combo Test
562003RU	NADAL® BCA/Hb Combo Test
572004N-10	NADAL® S. pneumoniae Test
572005	NADAL® Legionella/S. pneumoniae Test
580005N-25	NADAL® TSH Test
582003	NADAL® C. difficile Toxin A/B Test
582003_SSL	NADAL® C. difficile Toxin A/B Test
582004	NADAL® C. difficile GDH Ag Test
582008	NADAL® C. difficile Toxin A/B Test
582009	NADAL® C. perfringens Ag Test
582009_SSL	NADAL® C. perfringens Ag Test
582016N-10	NADAL® C. difficile Toxin A/B+GDH Test
600002N-30	NADAL® Rheumatoid Factors Test
600002N-30_SSL	NADAL® Rheumatoid Factors Test
611005N-10	NADAL® Gonorrhea Test
612004N-25	NADAL® CCA (Bilharzia) Test
622001N-30	NADAL® HAV IgM Test
622040N-30	NADAL® HEV IgM Test
622070N-30	NADAL® HAV IgG/IgM Test
622070N-30_SSL	NADAL® HAV IgG/IgM Test
652001N-30	NADAL® Chagas IgG Test
652001N-30_SSL	NADAL® Chagas IgG Test
662001N-30	NADAL® Leishmania Test
672001N-30	NADAL® Filariasis IgG/IgM Test
672001N-30_SSL	NADAL® Filariasis IgG/IgM Test


682002N-20	NADAL® Chikungunya IgM Test
692001N-30	NADAL® Typhoid IgG/IgM Test
692001N-30_SSL	NADAL® Typhoid IgG/IgM Test
712001	NADAL® AFP Test
712001_SSL	NADAL® AFP Test
722003	NADAL® CEA Test
722003_SSL	NADAL® CEA Test
790001	NADAL® Celiac Disease tTG Test
790002	NADAL® Celiac Disease tTG/Gliadine Test
795002	NADAL® ASO Latextest
795003	NADAL® ASO Latextest
795005	NADAL® CRP Latextest
795006	NADAL® CRP Latextest
795008	NADAL® Rheumatoid Factors Latextest
795008_SSL	NADAL® Rheumatoid Factors Latextest
795009	NADAL® Rheumatoid Factors Latextest
795010	NADAL® RPR Carbon Latextest
795011	NADAL® RPR Carbon Latextest
795015	NADAL® VDRL Latextest
795016	NADAL® VDRL Latextest
795016_SSL	NADAL® VDRL Latextest
795017	NADAL® Waaler Rose Latextest
795018	NADAL® Waaler Rose Latextest
795018_SSL	NADAL® Waaler Rose Latextest
795024	NADAL® IM Latextest
795024_SSL	NADAL® IM Latextest
795027	NADAL® TPHA Test

795027_SSL	NADAL® TPHA Test
795028	NADAL® TPHA Test
795028_SSL	NADAL® TPHA Test
795030	NADAL® Rose Bengale Latextest
840003N-10	NADAL® Trichomonas vaginalis Test
850003	NADAL® Astrovirus Test
850003_SSL	NADAL® Astrovirus Test
860001	NADAL® Bence Jones Protein Test
860001_SSL	NADAL® Bence Jones Protein Test
920001	NADAL® Norovirus Test
920001_SSL	NADAL® Norovirus Test
920002	NADAL® Norovirus GI/GII Test
920002_SSL	NADAL® Norovirus GI/GII Test
920002N-01	NADAL® Norovirus GI/GII Test
920002N-01_SSL	NADAL® Norovirus GI/GII Test

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### • test tubes

- **untreated test tubes**
  - w/o closing stopper
- **treated test tubes**
  - for platelet count from venous or capillary blood
  - for reticulocyte count
  - for haematological analyses
  - for preparing the blood-citrate mixture for the Hepato Quick
  - for preparing the blood-citrate mixture for coagulation physiological analyses
  - for serum collection
  - for plasma collection
  - for stabilising the enzyme of acid phosphatase
  - for blood sugar determination
  - for determination of the catecholamine

### • untreated reaction vessels

- w/o closing stoppers

### • closing stoppers

for test tubes and reaction vessels

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meets the provisions of the directive 98/79/EC and its transpositions in national laws which apply to it.

This declaration is valid for the batches released by KABE LABORTECHNIK GmbH.

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# Certificate

The Certification Body of  
TÜV Rheinland LGA Products GmbH

hereby certifies that the organization

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has established and applies a quality management system for medical devices  
for the following scope:

**Design and development, production and distribution of  
in vitro diagnostic devices and consumption materials  
for sample withdrawal, preparation and storage  
as well as single-use medical devices**

Proof has been furnished that the requirements specified in

**EN ISO 13485:2016**

are fulfilled. The quality management system is subject to yearly surveillance.

Effective Date: 2018-10-16  
Certificate Registration No.: SX 60133221 0001  
An audit was performed. Report No.: 21234760 009  
This Certificate is valid until: 2021-10-15

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**"VACUTEST KIMA" vacuum blood and biological liquids  
collection tubes in plastic.**

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**Hereby we declare**

under our sole responsibility that the above mentioned devices meet the applicable provisions of the Directive 98/79/EC as amended on "In Vitro Diagnostic Medical Devices".

All the supporting documents, as required by Annex III, in order to prove conformity to the Essential Requirements as listed in Annex I, are retained under the premises of the Manufacturer

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*IQNet, the association of the world's first class certification bodies, is the largest provider of management System Certification in the world. IQNet is composed of more than 30 bodies and counts over 150 subsidiaries all over the globe.*

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*Trading of the products of the Group: diagnostic kits, culture media for microbiology,  
plastic disposable labware, test tubes with predetermined vacuum and sterile needles.*

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*For timely and updated information about any changes in the certification status referred to in this certificate,  
please contact the number +39 02 725341 or email address info@icim.it.*

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**006376**

## EC Certificate

Directive 98/79/EC Annex IV, excluding Sections 4 and 6  
Full Quality Assurance System  
In Vitro Diagnostic Medical Devices

Registration No.: HL 60119814 0001

Report No.: 21265422 001

**Manufacturer:** Macherey-Nagel GmbH & Co. KG  
Neumann-Neander-Str. 6-8  
52355 Düren  
Deutschland

**Products:** Products for self-testing  
(see attachment for products and sites included)  
Replaces Certificate, Registration No.: HL 60076687 0001

**Expiry Date:** 2022-05-28

The Notified Body hereby declares that the requirements of Annex IV, excluding section 4 and 6 of the directive 98/79/EC have been met for the listed products. The above named manufacturer has established and applies a quality assurance system, which is subject to periodic surveillance, defined by Annex IV, section 5 of the aforementioned directive. For placing on the market of List A devices covered by this certificate an EC design-examination certificate according to Annex IV, section 4 and a verification of manufactured products according to section 6 is required.

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TÜV Rheinland LGA Products GmbH is a Notified Body according to Directive 98/79/EC concerning in vitro diagnostic medical devices with the identification number 0197.

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**Products for self-testing:**

- Single and multi-parameter disposable test strips for urine analysis
- indicator test strips and papers for measurement of pH in urine

**Additional site for warehousing and logistics:**

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