en

Iron2

04T98

H06712R02

B4T980

Iron2

FOR USE WITH

Alinity c

Read Highlighted Changes: Revised March 2022.

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Iron2

## INTENDED USE

The Iron2 assay is used for the direct colorimetric determination of iron without deproteinization in human serum or plasma on the Alinity c system.

The Iron2 assay is to be used as an aid in the diagnosis and treatment of diseases such as iron deficiency anemia, hemochromatosis (a disease associated with widespread deposit in the tissues of two iron-containing pigments, hemosiderin and hemofuscin, and characterized by pigmentation of the skin), and chronic renal disease.

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Iron exists in biological fluids as a component of hemoglobin and myoglobin and is bound in serum and plasma to transferrin, which acts as a carrier protein. Increased iron concentrations are seen in hemolytic anemias, hemochromatosis, and acute liver disease. Decreased iron concentrations are seen in iron deficiency and anemia of chronic disease, such as in chronic renal disease.<sup>1</sup> Major causes of iron deficiency include gastrointestinal and menstrual bleeding. For the assessment of the body's iron status, the measurement of transferrin and ferritin can provide more accurate information.<sup>2</sup>

## PRINCIPLES OF THE PROCEDURE

The Iron2 assay is an automated clinical chemistry assay.

At an acidic pH, iron is released from transferrin to which it is bound, and then quantitatively reduced to a ferrous state. The iron forms with ferene-S (3-(2-pyridyl)-5,6-bis-[2-(5-furylsulfonic acid)]-1,2,4-triazine), a stable colored complex of which the color intensity is proportional to the amount of iron in the sample.

Particular reaction conditions and a specific masking agent almost entirely eliminate the interference from copper.

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## Warnings and Precautions

- IVD**
- For *In Vitro* Diagnostic Use
- Rx ONLY**

### Safety Precautions

**CAUTION:** This product requires the handling of human specimens. It is recommended that all human-sourced materials and all consumables contaminated with potentially infectious materials be considered potentially infectious and handled in accordance with the OSHA Standard on Bloodborne Pathogens. Biosafety Level 2 or other appropriate regional, national, and institutional biosafety practices should be used for materials that contain, are suspected of containing, or are contaminated with infectious agents.<sup>3-6</sup>

|                                                                                                                                                                          |                                                                                   |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|
| The following warnings and precautions apply to: <b>R1</b>                                                                                                               |                                                                                   |
|   |                                                                                   |
| <b>WARNING</b>                                                                                                                                                           | Contains guanidine hydrochloride, acetic acid, thiourea and methylisothiazolones. |
| H302                                                                                                                                                                     | Harmful if swallowed.                                                             |
| H332                                                                                                                                                                     | Harmful if inhaled.                                                               |
| H315                                                                                                                                                                     | Causes skin irritation.                                                           |
| H319                                                                                                                                                                     | Causes serious eye irritation.                                                    |
| H317                                                                                                                                                                     | May cause an allergic skin reaction.                                              |
| H351                                                                                                                                                                     | Suspected of causing cancer.                                                      |
| H361                                                                                                                                                                     | Suspected of damaging fertility or the unborn child.                              |
| H402*                                                                                                                                                                    | Harmful to aquatic life.                                                          |
| H412                                                                                                                                                                     | Harmful to aquatic life with long lasting effects.                                |
| <b>Prevention</b>                                                                                                                                                        |                                                                                   |
| P202                                                                                                                                                                     | Do not handle until all safety precautions have been read and understood.         |
| P261                                                                                                                                                                     | Avoid breathing mist / vapors / spray.                                            |
| P264                                                                                                                                                                     | Wash hands thoroughly after handling.                                             |
| P271                                                                                                                                                                     | Use only outdoors or in a well-ventilated area.                                   |
| P273                                                                                                                                                                     | Avoid release to the environment.                                                 |
| P272                                                                                                                                                                     | Contaminated work clothing should not be allowed out of the workplace.            |
| P280                                                                                                                                                                     | Wear protective gloves / protective clothing / eye protection.                    |

| Response       |                                                                                                                                  |
|----------------|----------------------------------------------------------------------------------------------------------------------------------|
| P301+P330+P312 | IF SWALLOWED: Rinse mouth. Call a POISON CENTER or doctor / physician if you feel unwell.                                        |
| P302+P352      | IF ON SKIN: Wash with plenty of water.                                                                                           |
| P305+P351+P338 | IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. |
| P304+P340      | IF INHALED: Remove person to fresh air and keep comfortable for breathing.                                                       |
| P308+P313      | IF exposed or concerned: Get medical advice / attention.                                                                         |
| P333+P313      | If skin irritation or rash occurs: Get medical advice / attention.                                                               |
| P337+P313      | If eye irritation persists: Get medical advice / attention.                                                                      |
| P362+P364      | Take off contaminated clothing and wash it before reuse.                                                                         |
| Disposal       |                                                                                                                                  |
| P501           | Dispose of contents / container in accordance with local regulations.                                                            |

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| The following warnings and precautions apply to: <b>R2</b>                        |                                                                        |
|-----------------------------------------------------------------------------------|------------------------------------------------------------------------|
|  |                                                                        |
| <b>WARNING</b>                                                                    | Contains methylisothiazolones.                                         |
| H317                                                                              | May cause an allergic skin reaction.                                   |
| H402*                                                                             | Harmful to aquatic life.                                               |
| H412                                                                              | Harmful to aquatic life with long lasting effects.                     |
| Prevention                                                                        |                                                                        |
| P261                                                                              | Avoid breathing mist / vapors / spray.                                 |
| P272                                                                              | Contaminated work clothing should not be allowed out of the workplace. |
| P280                                                                              | Wear protective gloves / protective clothing / eye protection.         |
| P273                                                                              | Avoid release to the environment.                                      |
| Response                                                                          |                                                                        |
| P302+P352                                                                         | IF ON SKIN: Wash with plenty of water.                                 |
| P333+P313                                                                         | If skin irritation or rash occurs: Get medical advice / attention.     |
| P362+P364                                                                         | Take off contaminated clothing and wash it before reuse.               |
| Disposal                                                                          |                                                                        |
| P501                                                                              | Dispose of contents / container in accordance with local regulations.  |

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#### Reagent Handling

- Upon receipt, place reagent cartridges in an upright position for 1 hour before use to allow bubbles that may have formed to dissipate.

- If a reagent cartridge is dropped, place in an upright position for 1 hour before use to allow bubbles that may have formed to dissipate.
- Reagents are susceptible to the formation of foam and bubbles. Bubbles may interfere with the detection of the reagent level in the cartridge and cause insufficient reagent aspiration that may adversely affect results.

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#### Reagent Storage

|                 | Storage Temperature | Maximum Storage Time  | Additional Storage Instructions                                                                                                                                         |
|-----------------|---------------------|-----------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Unopened</b> | 2 to 8°C            | Until expiration date | Store in upright position.                                                                                                                                              |
| <b>Onboard</b>  | System Temperature  | 30 days               |                                                                                                                                                                         |
| <b>Opened</b>   | 2 to 8°C            | Until expiration date | Store in upright position. Do not reuse original reagent caps or replacement caps due to the risk of contamination and the potential to compromise reagent performance. |

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#### INSTRUMENT PROCEDURE

The Iron2 assay file must be installed on the Alinity c system prior to performing the assay.

The Alinity ci-series system software version 3.2.0 or higher must be installed on the Alinity c system prior to performing the assay.

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Other specimen types and collection tube types have not been verified with this assay.

| Specimen Types | Collection Tubes          |
|----------------|---------------------------|
| Serum          | Serum                     |
|                | Serum separator           |
| Plasma         | Lithium heparin           |
|                | Lithium heparin separator |
|                | Sodium heparin            |

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### Specimen Conditions

- Do not use:
  - heat-inactivated specimens
  - pooled specimens
  - grossly hemolyzed specimens
  - specimens with obvious microbial contamination
  - specimens with fungal growth
- For accurate results, serum and plasma specimens should be free of fibrin, red blood cells, and other particulate matter. Serum specimens from patients receiving anticoagulant or thrombolytic therapy may contain fibrin due to incomplete clot formation.
- To prevent cross contamination, use of disposable pipettes or pipette tips is recommended.

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NOTE: If fibrin, red blood cells, or other particulate matter are observed, mix by low speed vortex or by inverting 10 times prior to recentrifugation.

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- Frozen specimens must be completely thawed before mixing.
- Mix thawed specimens thoroughly by low speed vortex or by inverting 10 times.
- Visually inspect the specimens. If layering or stratification is observed, mix until specimens are visibly homogeneous.
- If specimens are not mixed thoroughly, inconsistent results may be obtained.
- Recentrifuge specimens.

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### Specimen Storage

| Specimen Type | Temperature                   | Maximum Storage Time   |
|---------------|-------------------------------|------------------------|
| Serum/Plasma  | Room temperature (20 to 25°C) | 10 hours <sup>7</sup>  |
|               | 2 to 8°C                      | 7 days <sup>8</sup>    |
|               | -20°C                         | 12 months <sup>9</sup> |

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For additional information on sample handling and processing, refer to CLSI GP44-A4.<sup>10</sup> The storage information provided here is based on references.

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Stored specimens must be inspected for particulates. If present, mix with a low speed vortex or by inversion and centrifuge the specimen to remove particulates prior to testing.

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- Iron2 assay file
- 04V6201 Consolidated Chemistry Calibrator
- Controls containing iron
- Saline (0.85% to 0.90% NaCl) for specimen dilution

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- If using primary or aliquot tubes, refer to the Alinity ci-series Operations Manual, Section 4 to ensure sufficient specimen is present.
- Minimum sample cup volume is calculated by the system and printed on the Order List report. To minimize the effects of evaporation, verify adequate sample cup volume is present prior to running the test.
- Minimum sample volume requirements:
  - Sample volume for single test: 20.0 µL.

NOTE: This amount does not include the dead volume plus the additional over-aspiration volume. For total sample volume requirements, refer to the Alinity ci-series Operations Manual, Section 4.
- Refer to the Consolidated Chemistry Calibrator package insert **[REF]** 04V6201 and/or commercially available control material package insert for preparation and usage.
- For general operating procedures, refer to the Alinity ci-series Operations Manual, Section 5.
- For optimal performance, it is important to perform routine maintenance as described in the Alinity ci-series Operations Manual, Section 9. Perform maintenance more frequently when required by laboratory procedures.

## Sample Dilution Procedures

Sample dilutions have not been evaluated for the Iron2 assay. Samples with an iron value exceeding 1143 µg/dL (204.6 µmol/L) are flagged with the code "> 1143 µg/dL" ("> 204.6 µmol/L"). The standard dilution factor for the Iron2 assay is 1:1.45.

For details on configuring automated dilutions, refer to the Alinity ci-series Operations Manual, Section 2.

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Calibration is stable for approximately 15 days (360 hours) but is required with each change in reagent lot. Verify calibration with at least 2 levels of controls according to the established quality control requirements for your laboratory. If control results fall outside acceptable ranges, recalibration may be necessary.

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- At least 2 levels of controls (low and high) are to be run every 24 hours.
- If more frequent control monitoring is required, follow the established quality control procedures for your laboratory.
- If quality control results do not meet the acceptance criteria defined by your laboratory, sample results may be suspect. Follow the established quality control procedures for your laboratory. Recalibration may be necessary. For troubleshooting information, refer to the Alinity ci-series Operations Manual, Section 10.
- Review quality control results and acceptance criteria following a change of reagent or calibrator lot.

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## Quality Control Guidance

Refer to "Basic QC Practices" by James O. Westgard, Ph.D. for guidance on laboratory quality control practices.<sup>11</sup>

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## RESULTS

### Calculation

The Iron2 assay utilizes the Linear data reduction method to generate a calibration and results.

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## Reportable Interval

Based on representative data for the limit of quantitation (LoQ) and the limit of detection (LoD), the ranges over which results can be reported are provided below according to the definitions from CLSI EP34, 1st ed.<sup>12</sup>

|                                                  | µg/dL    | µmol/L      |
|--------------------------------------------------|----------|-------------|
| Analytical Measuring Interval (AMI) <sup>a</sup> | 7 - 1143 | 1.3 - 204.6 |
| Reportable Interval <sup>b</sup>                 | 4 - 1143 | 0.7 - 204.6 |

<sup>a</sup> AMI: The AMI extends from the LoQ to the upper limit of quantitation (ULoQ). This is determined by the range of values in µg/dL (µmol/L) that demonstrated acceptable performance for linearity, imprecision, and bias.

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NOTE: The default Low Linearity value of the assay file corresponds to the lower limit of the reportable interval of 4 µg/dL (0.7 µmol/L). To flag values using the lower limit of the analytical measuring interval of 7 µg/dL (1.3 µmol/L), the operator must edit the Low Linearity value, adjusted by the standard dilution factor.

For detailed information on editing the result settings of assay parameters, refer to the Alinity ci-series Operations Manual, Section 2.

## LIMITATIONS OF THE PROCEDURE

- Results should be used in conjunction with other data; e.g., symptoms, results of other tests, and clinical impressions.
- Falsely elevated iron results may be observed at the low end of the analytical measuring interval in samples with triglyceride concentrations above 200 mg/dL.
- Falsely elevated iron results may be observed at the low end of the analytical measuring interval in samples with unconjugated bilirubin concentrations above 25 mg/dL.
- Iron dextran treatment can result in elevated total iron results.
- Use of the Iron2 assay for patients undergoing treatment with deferoxamine or other iron chelating compounds is not recommended.
- Transiently elevated iron levels can be observed post ingestion of supplements/vitamins that contain iron.<sup>13</sup>
- Rifampicin levels above 5 mg/L may produce artificially low results with the Iron2 assay.
- Substances that demonstrated interference with the Iron2 assay are listed in the SPECIFIC PERFORMANCE CHARACTERISTICS, Analytical Specificity, Interference section of this package insert.
- Potential interference has not been evaluated for substances other than those described in the SPECIFIC PERFORMANCE CHARACTERISTICS, Analytical Specificity, Interference section of this package insert.

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### Reference Range

| Age                           | Range (µg/dL) | Range (µmol/L) |
|-------------------------------|---------------|----------------|
| <b>Pediatric<sup>14</sup></b> |               |                |
| 0 to < 14 years               | 16 - 128      | 2.9 - 22.9*    |
| 14 to < 19 years (Female)     | 20 - 162      | 3.6 - 29.0*    |
| 14 to < 19 years (Male)       | 31 - 168      | 5.5 - 30.1*    |
| <b>Adult<sup>15</sup></b>     |               |                |
| Female                        | 50 - 170      | 9.0 - 30.4     |
| Male                          | 65 - 175      | 11.6 - 31.3    |

\* Alternate result units were calculated by Abbott and are not included in the citation provided.

Abbott has not evaluated reference ranges in the pediatric population.

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### Precision

#### Within-Laboratory Precision

A study was performed based on guidance from CLSI EP05-A3.<sup>16</sup> Testing was conducted using 3 lots of the Iron2 reagents, 3 lots of the Consolidated Chemistry Calibrator, 1 lot of commercially available controls, and 1 instrument. Two controls and 3 human serum panels were tested in a minimum of 2 replicates twice per day on 20 days on 3 reagent lot/calibrator lot combinations, where a unique reagent lot and a unique calibrator lot are paired. The performance from a representative combination is shown in the following table.

| Sample          | n  | Mean<br>(µg/dL) | Within-Run<br>(Repeatability) |     | Within-Laboratory <sup>a</sup> |                              |
|-----------------|----|-----------------|-------------------------------|-----|--------------------------------|------------------------------|
|                 |    |                 | SD                            | %CV | SD<br>(Range <sup>b</sup> )    | %CV<br>(Range <sup>b</sup> ) |
| Control Level 1 | 80 | 245             | 0.7                           | 0.3 | 1.4<br>(1.4 - 4.0)             | 0.6<br>(0.6 - 1.7)           |
| Control Level 2 | 80 | 69              | 0.5                           | 0.7 | 0.8<br>(0.7 - 1.2)             | 1.1<br>(1.0 - 1.8)           |
| Panel A         | 80 | 14              | 0.7                           | 4.6 | 0.8<br>(0.6 - 0.9)             | 5.6<br>(3.8 - 6.1)           |
| Panel B         | 80 | 375             | 0.7                           | 0.2 | 1.6<br>(1.6 - 4.2)             | 0.4<br>(0.4 - 1.2)           |
| Panel C         | 80 | 975             | 2.4                           | 0.2 | 2.8<br>(2.8 - 8.8)             | 0.3<br>(0.3 - 0.9)           |

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#### Reproducibility

A study was performed based on guidance from CLSI EP05-A3.<sup>16</sup> Testing was conducted using 1 lot of the Iron2 reagents, 1 lot of the Consolidated Chemistry Calibrator, 1 lot of commercially available controls, and 3 instruments. Each instrument was operated by a different technician, and each technician prepared an individual sample set. Two controls and 3 human serum panels were tested in a minimum of 3 replicates at 2 separate times per day on 5 different days.

| Sample          | n  | Mean<br>(µg/dL) | Repeatability |     | Within-Laboratory <sup>a</sup> |     | Reproducibility <sup>b</sup> |     |
|-----------------|----|-----------------|---------------|-----|--------------------------------|-----|------------------------------|-----|
|                 |    |                 | SD            | %CV | SD                             | %CV | SD                           | %CV |
| Control Level 1 | 90 | 245             | 0.7           | 0.3 | 0.9                            | 0.4 | 1.0                          | 0.4 |
| Control Level 2 | 90 | 69              | 0.4           | 0.5 | 0.4                            | 0.6 | 0.9                          | 1.3 |
| Panel A         | 90 | 14              | 0.4           | 2.8 | 0.4                            | 2.9 | 0.6                          | 4.3 |
| Panel B         | 90 | 374             | 0.9           | 0.2 | 1.0                            | 0.3 | 1.8                          | 0.5 |
| Panel C         | 90 | 970             | 2.4           | 0.2 | 2.7                            | 0.3 | 7.2                          | 0.7 |

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### Accuracy

A study was performed to estimate the bias of the Iron2 assay relative to standard reference material NIST SRM 3126. Testing was conducted using 3 concentrations of the standard across 3 lots of the Iron2 reagents, 2 lots of the Consolidated Chemistry Calibrator, and 1 instrument. The bias ranged from -3.7% to 2.4% across concentrations of the standard, calibrator, and reagent lots.

### Lower Limits of Measurement

A study was performed based on guidance from CLSI EP17-A2.<sup>17</sup> Testing was conducted using 3 lots of the Iron2 reagents on each of 2 instruments over a minimum of 3 days. The limit of blank (LoB), limit of detection (LoD), and limit of quantitation (LoQ) values are summarized below. These representative data support the lower limit of the analytical measuring interval.

|                  | µg/dL | µmol/L |
|------------------|-------|--------|
| LoB <sup>a</sup> | 1     | 0.2    |
| LoD <sup>b</sup> | 4     | 0.7    |
| LoQ <sup>c</sup> | 7     | 1.3    |

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<sup>b</sup> The LoD presented in the table is in alignment with the Iron2 assay on the ARCHITECT c System. The observed LoD on the Alinity c system was 3 µg/dL (0.5 µmol/L) and represents the lowest concentration at which the analyte can be detected with 95% probability based on n ≥ 60 replicates of low-analyte level samples.

<sup>c</sup> The LoQ presented in the table is in alignment with the LoQ for the Iron2 on the ARCHITECT c System. The observed LoQ on the Alinity c system was 5 µg/dL (0.9 µmol/L). This LoQ is defined as the lowest concentration at which a maximum allowable precision of 20 %CV was met and was determined from n ≥ 60 replicates of low-analyte level samples.

## Linearity

A study was performed based on guidance from CLSI EP06-A.<sup>18</sup>

This assay is linear across the analytical measuring interval of 7 to 1143 µg/dL (1.3 to 204.6 µmol/L).

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#### Potentially Interfering Endogenous Substances

A study was performed based on guidance from CLSI EP07, 3rd ed.<sup>19</sup> Each substance was tested at 2 levels of the analyte (approximately 50 µg/dL and 180 µg/dL).

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**Interference beyond ± 10% (based on 95% Confidence Interval [CI])** was observed at the concentrations shown below for the following substances.

| Interference beyond ± 10% (based on 95% Confidence Interval [CI]) |                   |                 |               |                 |                         |
|-------------------------------------------------------------------|-------------------|-----------------|---------------|-----------------|-------------------------|
| Potentially Interfering Substance                                 | Interferent Level |                 | Analyte Level |                 | % Interference (95% CI) |
|                                                                   | Default Units     | Alternate Units | Default Units | Alternate Units |                         |
| Bilirubin, unconjugated                                           | 30 mg/dL          | 513 µmol/L      | 50 µg/dL      | 8.95 µmol/L     | 12%<br>(11%, 13%)       |
| Triglycerides                                                     | 250 mg/dL         | 2.82 mmol/L     | 50 µg/dL      | 8.95 µmol/L     | 11%<br>(10%, 12%)       |

#### Potentially Interfering Exogenous Substances

A study was performed based on guidance from CLSI EP07, 3rd ed.<sup>19</sup> Each substance was tested at 2 levels of the analyte (approximately 50 µg/dL and 180 µg/dL).

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### No Significant Interference (Interference within ± 10%)

| Potentially Interfering Substance   | Interferent Level |                 |
|-------------------------------------|-------------------|-----------------|
|                                     | Default Units     | Alternate Units |
| Mebendazole                         | 2 mg/L            | 6.78 µmol/L     |
| Methylidopa                         | 25 mg/L           | 118 µmol/L      |
| Metronidazole                       | 130 mg/L          | 759 µmol/L      |
| Phenylbutazone                      | 330 mg/L          | 1069 µmol/L     |
| Rifampicin                          | 5 mg/L            | 6.10 µmol/L     |
| Sodium heparin                      | 4 U/mL            | N/A             |
| Stanozolol                          | 60 mg/L           | 183 µmol/L      |
| Theophylline (1,3-dimethylxanthine) | 60 mg/L           | 333 µmol/L      |

N/A = Not applicable

**Interference beyond ± 10% (based on 95% Confidence Interval [CI])** was observed at the concentrations shown below for the following substances.

| Interference beyond ± 10% (based on 95% Confidence Interval [CI]) |                   |                 |               |                 |                         |
|-------------------------------------------------------------------|-------------------|-----------------|---------------|-----------------|-------------------------|
| Potentially Interfering Substance                                 | Interferent Level |                 | Analyte Level |                 | % Interference (95% CI) |
|                                                                   | Default Units     | Alternate Units | Default Units | Alternate Units |                         |
| Deferoxamine                                                      | 0.05 mg/dL        | 0.89 µmol/L     | 50 µg/dL      | 8.95 µmol/L     | -11%<br>(-11%, -10%)    |
| Iron dextran                                                      | 17 mg/L           | 111 µmol/L      | 50 µg/dL      | 8.95 µmol/L     | 110%<br>(108%, 111%)    |
| Rifampicin                                                        | 10 mg/L           | 12.2 µmol/L     | 50 µg/dL      | 8.95 µmol/L     | -14%<br>(-15%, -13%)    |

Interferences from medication or endogenous substances may affect results.<sup>20</sup>

## Method Comparison

A study was performed based on guidance from CLSI EP09-A3,<sup>21</sup> using the Passing-Bablok regression method.

| Iron2 on the Alinity c system vs Iron on the ARCHITECT c System |     |                   |                         |              |       |                           |
|-----------------------------------------------------------------|-----|-------------------|-------------------------|--------------|-------|---------------------------|
|                                                                 | n   | Units             | Correlation Coefficient | Intercept    | Slope | Concentration Range       |
| Serum                                                           | 113 | µg/dL<br>(µmol/L) | 1.00                    | -1<br>(-0.2) | 1.05  | 16 - 879<br>(2.9 - 157.2) |

| Iron2 on the Alinity c system vs Iron2 on the ARCHITECT c System |     |                   |                         |            |       |                           |
|------------------------------------------------------------------|-----|-------------------|-------------------------|------------|-------|---------------------------|
|                                                                  | n   | Units             | Correlation Coefficient | Intercept  | Slope | Concentration Range       |
| Serum                                                            | 115 | µg/dL<br>(µmol/L) | 1.00                    | 1<br>(0.1) | 1.00  | 16 - 933<br>(2.9 - 166.9) |

## BIBLIOGRAPHY

- Chen TK, Knically DH, Grams ME. Chronic kidney disease diagnosis and management: a review. *JAMA* 2019;322(13):1294-1304
- Jacobs DS, DeMott WR, Oxley DK. *Laboratory Test Handbook*, 5th ed. Hudson, OH: Lexi-Comp; 2001:203-204.
- US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
- US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 6th ed. Washington, DC: US Government Printing Office; June 2020.
- World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
- Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.
- Henriksen LO, Faber NR, Moller MF, et al. Stability of 35 biochemical and immunological routine tests after 10 hours storage and transport of human whole blood at 21°C. *Scand J Clin Lab Invest* 2014;74(7):603-610.
- Klat RL, Stein CDS, Santos MD, et al. Effects of freeze-thaw cycles and assessment of short-term storage stability on serum iron, ferritin, and transferrin. *Clin Lab* 2021;67(2).
- Chen H, Sternberg MR, Schleicher RL, et al. Long-term stability of 18 nutritional biomarkers stored at -20 °C and 5 °C for up to 12 months. *J App Lab Med* 2018;3(1):100-108.
- Clinical and Laboratory Standards Institute (CLSI). *Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests; Approved Guideline—Fourth Edition*. CLSI Document GP44-A4. Wayne, PA: CLSI; 2010.
- Westgard JO. *Basic QC Practices; Training in Statistical Quality Control for Medical Laboratories*. 4th ed. Westgard QC, Inc.; 2016.
- Clinical and Laboratory Standards Institute (CLSI). *Establishing and Verifying an Extended Measuring Interval Through Specimen Dilution and Spiking*. 1st ed. CLSI Guideline EP34. Wayne, PA: CLSI; 2018.
- Burtis CA, Bruns DE, editors. *Tietz Fundamentals of Clinical Chemistry and Molecular Diagnostics*. 7th ed. St. Louis, MO: Saunders Elsevier; 2015.
- Colantonio DA, Kyriakopoulou L, Chan MK, et al. Closing the gaps in pediatric laboratory reference intervals: a CALIPER database of 40 biochemical markers in a healthy and multiethnic population of children. *Clin Chem* 2012;58:854-868.
- Wu AHB, editor. *Tietz Clinical Guide to Laboratory Tests*. 4th ed. St. Louis, MO: Elsevier Saunders; 2006.
- Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Precision of Quantitative Measurement Procedures: Approved Guideline—Third Edition*. CLSI Document EP05-A3. Wayne, PA: CLSI; 2014.
- Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition*. CLSI Document EP17-A2. Wayne, PA: CLSI; 2012.
- Clinical and Laboratory Standards Institute (CLSI). *Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline*. CLSI Document EP06-A. Wayne, PA: CLSI; 2003.
- Clinical and Laboratory Standards Institute (CLSI). *Interference Testing in Clinical Chemistry*. 3rd ed. CLSI Guideline EP07. Wayne, PA: CLSI; 2018.
- Young DS. Laboratory test listings. In: *Effects of Drugs on Clinical Laboratory Tests*. 5th ed. AACC Press; 2000:chap 3.
- Clinical and Laboratory Standards Institute (CLSI). *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition*. CLSI Document EP09-A3. Wayne, PA: CLSI; 2013.

## Key to Symbols

| ISO 15223 Symbols                |                                                                                         |
|----------------------------------|-----------------------------------------------------------------------------------------|
|                                  | Consult instructions for use                                                            |
|                                  | Manufacturer                                                                            |
|                                  | Sufficient for                                                                          |
|                                  | Temperature limitation                                                                  |
|                                  | Use by/Expiration date                                                                  |
| <b>IVD</b>                       | In Vitro Diagnostic Medical Device                                                      |
| <b>LOT</b>                       | Lot Number                                                                              |
| <b>REF</b>                       | List Number                                                                             |
| <b>SN</b>                        | Serial number                                                                           |
| Other Symbols                    |                                                                                         |
| <b>DISTRIBUTED IN THE USA BY</b> | Distributed in the USA by                                                               |
| <b>FOR USE WITH</b>              | Identifies products to be used together                                                 |
| <b>INFORMATION FOR USA ONLY</b>  | Information needed for United States of America only                                    |
| <b>PRODUCT OF IRELAND</b>        | Product of Ireland                                                                      |
| <b>R1</b>                        | Reagent 1                                                                               |
| <b>R2</b>                        | Reagent 2                                                                               |
| <b>Rx ONLY</b>                   | For use by or on the order of a physician only (applicable to USA classification only). |

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## BIBLIOGRAFIE

- US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
- US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
- World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
- Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.
- Burtis CA, Ashwood ER, Bruns DE, editors. *Tietz Textbook of Clinical Chemistry and Molecular Diagnostics*, 4th ed. St. Louis, MO: Elsevier Saunders; 2006:1912.
- Taylor EC, Sethi B. Stability of 27 biochemistry analytes in storage at a range of temperatures after centrifugation. *Br J Biomed Sci* 2011;68:147-157.
- Wu W, Yang D, Tiselius HG, et al. Collection and storage of urine specimens for measurement of urolithiasis risk factors. *Urology* 2015;85:299-303.
- Remer T, Montenegro-Bethancourt G, Shi L. Long-term urine biobanking: storage stability of clinical chemical parameters under moderate freezing conditions without use of preservatives. *Clin Biochem* 2014;47:307-311.

- Clinical and Laboratory Standards Institute (CLSI). *Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests; Approved Guideline—Fourth Edition*. CLSI Document GP44-A4. Wayne, PA: CLSI; 2010.
- Westgard JO. *Basic QC Practices*. 3rd ed. Madison, WI: Westgard Quality Corporation; 2010.
- Wu AHB. *Tietz Clinical Guide to Laboratory Tests*, 4th ed. Philadelphia, PA: WB Saunders; 2006:706-708.
- Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition*. CLSI Document EP05-A2. Wayne, PA: CLSI; 2004.
- Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition*. CLSI Document EP17-A2. Wayne, PA: CLSI; 2012.
- Clinical and Laboratory Standards Institute (CLSI). *Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline*. CLSI Document EP06-A. Wayne, PA: CLSI; 2003.
- Clinical and Laboratory Standards Institute (CLSI). *Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition*. CLSI Document EP07-A2. Wayne, PA: CLSI; 2005.
- Young DS. Laboratory test listings. In: *Effects of Drugs on Clinical Laboratory Tests*. 5th ed. AACC Press; 2000:chap 3.
- Clinical and Laboratory Standards Institute (CLSI). *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition*. CLSI Document EP09-A3. Wayne, PA: CLSI; 2013.

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## BIBLIOGRAFIE

- Burtis CA, Ashwood ER, Bruns DE, editors. *Tietz Textbook of Clinical Chemistry and Molecular Diagnostics*, 5th ed. St. Louis, MO: Elsevier Saunders; 2012:807-810, 812-813.
- Jacobs DS, DeMott WR, Oxley DK, et al. *Laboratory Test Handbook*. Hudson, OH: Lexi-Comp; 2001:144-147, 258-260, 275-279.
- US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
- US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
- World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
- Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.
- Kaplan LA, Pesce AJ, Kazmierczak SC. *Clinical Chemistry – Theory, Analysis, and Correlation*, 3rd ed. St Louis, MO: Mosby, 1996:461–463.
- McPherson RA, Bidkorphe EK, Castellani WJ, et al. *Analysis of Body Fluids in Clinical Chemistry; Approved Guideline (C49-A)*. Wayne, PA: The National Committee for Clinical Laboratory Standards, 2007.
- Guder WG, da Fonseca-Wollheim F, Heil W, et al. *The Quality of Diagnostic Samples*. Darmstadt, Germany: GIT Verlag; 2001:22,40,44,52.
- US Pharmacopeial Convention, Inc. General notices. In: *US Pharmacopeia National Formulary*. 1995 ed (USP 23/NF18). Rockville, MD: The US Pharmacopeial Convention, Inc; 1994:11.
- Pesce AJ, Kaplan LA, editors. *Methods in Clinical Chemistry*. St Louis, MO: CV Mosby, 1987:75.
- Westgard JO. *Basic QC Practices*. 3rd ed. Madison, WI: Westgard Quality Corporation; 2010.
- Burtis CA, Ashwood ER, Bruns DE, editors. *Tietz Textbook of Clinical Chemistry and Molecular Diagnostics*, 5th ed. St. Louis, MO: Elsevier Saunders; 2012:2139,2164,2168.
- Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition*. CLSI Document EP05-A2. Wayne, PA: CLSI; 2004.
- Clinical and Laboratory Standards Institute (CLSI). *Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline*. CLSI Document EP06-A. Wayne, PA: CLSI; 2003.
- Clinical and Laboratory Standards Institute (CLSI). *Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition*. CLSI Document EP07-A2. Wayne, PA: CLSI; 2005.
- Young DS. *Effects of Drugs on Clinical Laboratory Tests*, 5th ed. Washington, DC: AACC Press, 2000:3-176, 3-644, 3-714.

18. Clinical and Laboratory Standards Institute (CLSI). *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition*. CLSI Document EP09-A3. Wayne, PA: CLSI; 2013.

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For laboratory professional use only.

## NAME

Phosphorus2 (also referred to as Phos2)

## INTENDED USE

The Phosphorus2 assay is used for the quantitation of phosphorus in human serum, plasma, or urine on the Alinity c system.

Measurements of phosphorus (inorganic) are used in the diagnosis and treatment of various disorders, including parathyroid gland and kidney diseases, and vitamin D imbalance.

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Phosphorus is a major component of bone mineral, phospholipids in cell membranes, and nucleic acids. Phosphorus acts as a major pH buffer in serum and urine.<sup>1</sup> Phosphorus is essential for normal muscle contractility, neurologic function, electrolyte transport, and oxygen carrying by hemoglobin (2,3 diphosphoglycerate).<sup>2</sup>

Serum phosphate is filtered at the glomerulus and is reabsorbed primarily in the proximal tubule. Parathyroid hormone (PTH) increases renal phosphate excretion by inhibiting the sodium-phosphate cotransporter in the proximal tubule, whereas vitamin D enhances intestinal phosphate absorption.<sup>1</sup>

Hypophosphatemia may be caused by decreased intake, impaired intestinal absorption, redistribution into cells or bones, and renal losses (as in Fanconi syndrome).<sup>1</sup> Clinical manifestation of severe hypophosphatemia include encephalopathy, dilated cardiomyopathy, generalized muscle weakness leading to respiratory failure, destruction of muscles (rhabdomyolysis), and hemolysis.<sup>1</sup> Hemolysis of red blood cells may be seen with serum phosphorus levels < 0.5 mg/dL.<sup>3</sup>

Hyperphosphatemia is caused by excessive phosphate intake, increased intestinal absorption, redistribution from intracellular stores, or impaired renal excretion. Acute levels can increase the risk for precipitation of calcium phosphate in the kidney and soft tissues.<sup>1</sup> In chronic hyperphosphatemia (renal insufficiency) with a phosphate level > 6.5 mg/dL, mortality is higher due to the increased risk for the development of coronary and other vascular calcification, leading to increased systolic blood pressure and left ventricular hypertrophy.<sup>1</sup> A rapid increase in serum phosphorus can result in hypocalcemia, which can cause tetany, hypotension, seizures, and cardiac arrhythmias.<sup>1</sup>

A urine phosphorus-creatinine ratio and urine fractional excretion of phosphorus may provide useful information about therapies altering intestinal absorption or urine phosphorus handling.<sup>4</sup>

## PRINCIPLES OF THE PROCEDURE

The Phosphorus2 assay is an automated clinical chemistry assay. Inorganic phosphate reacts with ammonium molybdate to form a heteropolyacid complex.<sup>5</sup> The absorbance at 340 nm is directly proportional to the inorganic phosphorus level in the sample.

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### Kit Contents

Phosphorus2 Reagent Kit 04U03

NOTE: Some kit sizes may not be available. Please contact your local distributor.

Volumes (mL) listed in the following table indicate the volume per cartridge.

| REF                          | 04U0320 | 04U0330 |
|------------------------------|---------|---------|
| Tests per cartridge          | 280     | 700     |
| Number of cartridges per kit | 4       | 4       |
| Tests per kit                | 1120    | 2800    |
| R1                           | 20.7 mL | 44.4 mL |
| R2                           | 20.0 mL | 44.0 mL |

R2 Active ingredient: ammonium molybdate tetrahydrate 2.90 g/L.  
Preservative: ProClin 300.

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### Safety Precautions

**CAUTION:** This product requires the handling of human specimens. It is recommended that all human-sourced materials and all consumables contaminated with potentially infectious materials be considered potentially infectious and handled in accordance with the OSHA Standard on Bloodborne Pathogens. Biosafety Level 2 or other appropriate regional, national, and institutional biosafety practices should be used for materials that contain, are suspected of containing, or are contaminated with infectious agents.<sup>6-9</sup>

|                                                                                     |                                                                |
|-------------------------------------------------------------------------------------|----------------------------------------------------------------|
| The following warnings and precautions apply to: <b>R1</b>                          |                                                                |
|  |                                                                |
| <b>DANGER</b>                                                                       | Contains sulfuric acid.                                        |
| H314                                                                                | Causes severe skin burns and eye damage.                       |
| H290                                                                                | May be corrosive to metals.                                    |
| <b>Prevention</b>                                                                   |                                                                |
| P234                                                                                | Keep only in original container.                               |
| P260                                                                                | Do not breathe mist / vapors / spray.                          |
| P264                                                                                | Wash hands thoroughly after handling.                          |
| P280                                                                                | Wear protective gloves / protective clothing / eye protection. |

| Response       |                                                                                                                                  |
|----------------|----------------------------------------------------------------------------------------------------------------------------------|
| P301+P330+P331 | IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.                                                                               |
| P305+P351+P338 | IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. |
| P303+P361+P353 | IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water / shower.                            |
| P310           | Immediately call a POISON CENTER or doctor / physician.                                                                          |
| P390           | Absorb spillage to prevent material damage.                                                                                      |
| Disposal       |                                                                                                                                  |
| P501           | Dispose of contents / container in accordance with local regulations.                                                            |

The following warnings and precautions apply to: **R2**

| WARNING    |                                                                        |
|------------|------------------------------------------------------------------------|
| H317       | Contains methylisothiazolones.                                         |
| H402*      | May cause an allergic skin reaction.                                   |
| H402*      | Harmful to aquatic life.                                               |
| H412       | Harmful to aquatic life with long lasting effects.                     |
| Prevention |                                                                        |
| P261       | Avoid breathing mist / vapors / spray.                                 |
| P272       | Contaminated work clothing should not be allowed out of the workplace. |
| P273       | Avoid release to the environment.                                      |
| P280       | Wear protective gloves / protective clothing / eye protection.         |
| Response   |                                                                        |
| P302+P352  | IF ON SKIN: Wash with plenty of water.                                 |
| P333+P313  | If skin irritation or rash occurs: Get medical advice / attention.     |
| P362+P364  | Take off contaminated clothing and wash it before reuse.               |
| Disposal   |                                                                        |
| P501       | Dispose of contents / container in accordance with local regulations.  |

\* Not applicable where regulation EC 1272/2008 (CLP) has been implemented.

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#### Reagent Handling

- Upon receipt, place reagent cartridges in an upright position for 1 hour before use to allow bubbles that may have formed to dissipate.
- If a reagent cartridge is dropped, place in an upright position for 2 hours before use to allow bubbles that may have formed to dissipate.
- Reagents are susceptible to the formation of foam and bubbles. Bubbles may interfere with the detection of the reagent level in the cartridge and cause insufficient reagent aspiration that may adversely affect results.

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#### Reagent Storage

|          | Storage Temperature | Maximum Storage Time  | Additional Storage Instructions                                                                                                                                         |
|----------|---------------------|-----------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unopened | 2 to 8°C            | Until expiration date | Store in upright position.                                                                                                                                              |
| Onboard  | System Temperature  | 28 days               |                                                                                                                                                                         |
| Opened   | 2 to 8°C            | Until expiration date | Store in upright position. Do not reuse original reagent caps or replacement caps due to the risk of contamination and the potential to compromise reagent performance. |

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#### Indications of Reagent Deterioration

Deterioration of the reagents may be indicated when a calibration error occurs or a control value is out of the specified range. Associated test results are invalid, and samples must be retested. Assay recalibration may be necessary.

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#### INSTRUMENT PROCEDURE

The Phosphorus2 assay file must be installed on the Alinity c system prior to performing the assay.

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#### SPECIMEN COLLECTION AND PREPARATION FOR ANALYSIS

##### Specimen Types

The specimen types listed below were verified for use with this assay.

Other specimen types, collection tube types, and anticoagulants have not been verified with this assay.

| Specimen Type | Collection Vessel               | Special Conditions |
|---------------|---------------------------------|--------------------|
| Serum         | Serum tubes                     |                    |
|               | Serum separator tubes           |                    |
| Plasma        | Dipotassium EDTA tubes          |                    |
|               | Lithium heparin tubes           |                    |
|               | Lithium heparin separator tubes |                    |
|               | Sodium heparin tubes            |                    |

| Specimen Type                                                                                            | Collection Vessel                | Special Conditions                                                                                                            |
|----------------------------------------------------------------------------------------------------------|----------------------------------|-------------------------------------------------------------------------------------------------------------------------------|
| Urine (random specimens or timed specimens collected over intervals shorter than 24 hours) <sup>10</sup> | Clean glass or plastic container |                                                                                                                               |
| Urine (24 hour) <sup>10</sup>                                                                            | Clean glass or plastic container | Urine specimens should be collected in 6 mol/L HCl, 20 to 30 mL, to avoid precipitation of phosphate complexes. <sup>10</sup> |

- Liquid anticoagulants may have a dilution effect resulting in lower concentration values for individual specimens.

The instrument does not provide the capability to verify specimen types. It is the responsibility of the operator to verify that the correct specimen types are used in the assay.

### Specimen Conditions

- Do not use:
  - heat-inactivated specimens
  - pooled specimens
  - grossly hemolyzed specimens
  - specimens with obvious microbial contamination
  - specimens with fungal growth
- For accurate results, serum and plasma specimens should be free of fibrin, red blood cells, and other particulate matter. Serum specimens from patients receiving anticoagulant or thrombolytic therapy may contain fibrin due to incomplete clot formation. For accurate results, plasma specimens should be free of platelets and other particulate matter. Ensure centrifugation is adequate to remove platelets.
- To prevent cross contamination, use of disposable pipettes or pipette tips is recommended.

### Preparation for Analysis

- Follow the tube manufacturer's processing instructions for collection tubes. Gravity separation is not sufficient for specimen preparation.
- Specimens should be free of bubbles. Remove bubbles with an applicator stick before analysis. Use a new applicator stick for each specimen to prevent cross contamination.

To ensure consistency in results, recentrifuge specimens prior to testing if

- they contain fibrin, red blood cells, or other particulate matter.

NOTE: If fibrin, red blood cells, or other particulate matter are observed, mix by low-speed vortex or by inverting 10 times prior to recentrifugation.

Prepare frozen specimens as follows:

- Frozen specimens must be completely thawed before mixing.
- Mix thawed specimens thoroughly by low-speed vortex or by inverting 10 times.
- Visually inspect the specimens. If layering or stratification is observed, mix until specimens are visibly homogeneous.
- If specimens are not mixed thoroughly, inconsistent results may be obtained.
- Recentrifuge specimens.

Recentrifugation of Specimens

- Transfer specimens to a centrifuge tube and centrifuge.
- Transfer clarified specimen to a sample cup or secondary tube for testing. For centrifuged specimens with a lipid layer, transfer only the clarified specimen and not the lipemic material.

### Specimen Storage

| Specimen Type    | Temperature                   | Collection Vessel                                                                                                 | Maximum Storage Time  | Special Instructions                                    |
|------------------|-------------------------------|-------------------------------------------------------------------------------------------------------------------|-----------------------|---------------------------------------------------------|
| Serum/<br>Plasma | Room temperature (20 to 25°C) | Serum tubes                                                                                                       | 24 hours              |                                                         |
|                  |                               | Serum separator tubes                                                                                             |                       |                                                         |
|                  |                               | Lithium heparin plasma tubes<br>Lithium heparin separator tubes<br>Sodium heparin tubes<br>Dipotassium EDTA tubes |                       |                                                         |
|                  | 2 to 8°C                      | Serum tubes                                                                                                       | 3 days                |                                                         |
|                  |                               | Serum separator tubes                                                                                             |                       |                                                         |
|                  |                               | Lithium heparin plasma tubes<br>Lithium heparin separator tubes<br>Sodium heparin tubes<br>Dipotassium EDTA tubes |                       |                                                         |
|                  | -20°C                         | Serum tubes                                                                                                       | 30 days <sup>11</sup> | Avoid multiple freeze/thaw cycles. <sup>11</sup>        |
|                  |                               | Serum separator tubes                                                                                             |                       |                                                         |
|                  |                               | Lithium heparin plasma tubes<br>Lithium heparin separator tubes<br>Sodium heparin tubes<br>Dipotassium EDTA tubes |                       |                                                         |
| Urine            | Room temperature (20 to 25°C) | Clean glass or plastic container                                                                                  | 4 days                | Acidify to pH < 5                                       |
|                  | 2 to 8°C                      | Clean glass or plastic container                                                                                  | 7 days                | Acidify to pH < 5                                       |
|                  | -20°C                         | Clean glass or plastic container                                                                                  | 1 month               | Acidify to pH < 5<br>Avoid multiple freeze/thaw cycles. |

It is the responsibility of the individual laboratory to determine specific specimen stability criteria for their laboratory per their laboratory workflow.

For additional information on sample handling and processing, refer to CLSI GP44-A4.<sup>12</sup> The storage information provided here is based on references or data maintained by the manufacturer.

Each laboratory may establish a range around -20°C from either the freezer manufacturer's specifications or your laboratory standard operating procedure(s) for specimen storage.

Stored specimens must be inspected for particulates. If present, mix with a low-speed vortex or by inversion and centrifuge the specimen to remove particulates prior to testing.

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- Phosphorus2 assay file
- 04V6201 Consolidated Chemistry Calibrator
- Controls containing phosphorus
- Saline (0.85% to 0.90% NaCl) for specimen dilution

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- If using primary or aliquot tubes, refer to the Alinity ci-series Operations Manual, Section 4 to ensure sufficient specimen is present.
- Minimum sample cup volume is calculated by the system and printed on the Order List report. To minimize the effects of evaporation, verify adequate sample cup volume is present prior to running the test.
- Minimum sample volume requirements:
  - Sample volume for single test: 3.6 µL (serum/plasma); 16.0 µL (urine).

NOTE: This amount does not include the dead volume plus the additional over-aspiration volume. For total sample volume requirements, refer to the Alinity ci-series Operations Manual, Section 4.
- Refer to the Consolidated Chemistry Calibrator package insert <sup>REF</sup> 04V6201 and/or commercially available control material package insert for preparation and usage.
- For general operating procedures, refer to the Alinity ci-series Operations Manual, Section 5.
- For optimal performance, it is important to perform routine maintenance as described in the Alinity ci-series Operations Manual, Section 9. Perform maintenance more frequently when required by laboratory procedures.

### Sample Dilution Procedures

#### Serum/Plasma

The standard dilution factor, applied automatically by the system software to all results for the Phosphorus2 assay, is 1:1.73. Samples with a phosphorus value exceeding 25.3 mg/dL (8.17 mmol/L) are flagged with the code "> 25.3 mg/dL" (> 8.17 mmol/L). Automated or manual sample dilutions have not been evaluated for the Phosphorus2 assay.

#### Urine

The standard dilution factor, applied automatically by the system software to all results for the Phosphorus2 assay, is 1:17.33. Samples with a phosphorus value exceeding 186.2 mg/dL (60.14 mmol/L) are flagged with the code "> 186.2 mg/dL" (> 60.14 mmol/L). Automated or manual sample dilutions have not been evaluated for the Phosphorus2 assay.

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Calibration is stable for approximately 14 days (336 hours) but is required with each change in reagent cartridge. Verify calibration with at least 2 levels of controls according to the established quality control requirements for your laboratory. If control results fall outside acceptable ranges, recalibration may be necessary.

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- At least 2 levels of controls (low and high) are to be run every 24 hours.
- Run both levels of control with each cartridge change.
- If more frequent control monitoring is required, follow the established quality control procedures for your laboratory.
- If quality control results do not meet the acceptance criteria defined by your laboratory, sample results may be suspect. Follow the established quality control procedures for your laboratory. Recalibration may be necessary. For troubleshooting information, refer to the Alinity ci-series Operations Manual, Section 10.
- Review quality control results and acceptance criteria following a change of reagent or calibrator lot.

Controls should be used according to the guidelines and recommendations of the control manufacturer. Concentration ranges provided in the control package insert should be used only for guidance.

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### Quality Control Guidance

Refer to "Basic QC Practices" by James O. Westgard, Ph.D. for guidance on laboratory quality control practices.<sup>13</sup>

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## RESULTS

### Calculation

The Phosphorus2 (Phos2) assay utilizes the Linear data reduction method to generate a calibration and results for both the serum/plasma and urine applications.

Urine sample quantification (Phos2-U) is achieved using the calibration generated with the Phosphorus2 (Phos2) assay parameter file.

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### Reportable Interval

Based on representative data for the limit of quantitation (LoQ) and the limit of detection (LoD), the ranges over which results can be reported are provided below according to the definitions from CLSI EP34, 1st ed.<sup>14</sup>

#### Serum

|                                                  | mg/dL      | mmol/L      |
|--------------------------------------------------|------------|-------------|
| Analytical Measuring Interval (AMI) <sup>a</sup> | 0.5 - 25.3 | 0.16 - 8.17 |
| Reportable Interval <sup>b</sup>                 | 0.2 - 25.3 | 0.06 - 8.17 |

<sup>a</sup> AMI: The AMI extends from the LoQ to the upper limit of quantitation (ULoQ). This is determined by the range of values in mg/dL (mmol/L) that demonstrated acceptable performance for linearity, imprecision, and bias.

<sup>b</sup> The reportable interval extends from the LoD to the upper limit of the AMI.

NOTE: The default Low Linearity value of the assay file corresponds to the lower limit of the AMI. Samples with a phosphorus value below 0.5 mg/dL (0.16 mmol/L) are reported as "< 0.5 mg/dL" (< 0.16 mmol/L).

## Urine

|                                                  | mg/dL       | mmol/L       |
|--------------------------------------------------|-------------|--------------|
| Analytical Measuring Interval (AMI) <sup>a</sup> | 2.4 - 186.2 | 0.78 - 60.14 |
| Reportable Interval <sup>b</sup>                 | 1.0 - 186.2 | 0.32 - 60.14 |

<sup>a</sup> AMI: The AMI extends from the LoQ to the upper limit of quantitation (ULoQ). This is determined by the range of values in mg/dL (mmol/L) that demonstrated acceptable performance for linearity, imprecision, and bias.

<sup>b</sup> The reportable interval extends from the LoD to the upper limit of the AMI.

NOTE: The default Low Linearity value of the assay file corresponds to the lower limit of the AMI. Samples with a phosphorus value below 2.4 mg/dL (0.78 mmol/L) are reported as "< 2.4 mg/dL" ("< 0.78 mmol/L").

## LIMITATIONS OF THE PROCEDURE

- Results should be used in conjunction with other data; e.g., symptoms, results of other tests, and clinical impressions.
- Substances that demonstrated interference with the Phosphorus2 assay are listed in the SPECIFIC PERFORMANCE CHARACTERISTICS, Analytical Specificity, Interference section of this package insert.
- Potential interference has not been evaluated for substances other than those described in the SPECIFIC PERFORMANCE CHARACTERISTICS, Analytical Specificity, Interference section of this package insert.
- In very rare cases, gammopathy may cause unreliable results.<sup>15</sup>

## EXPECTED VALUES

It is recommended that each laboratory determine its own reference range based upon its particular locale and population characteristics.

### Reference Range (Serum)

| Category              | Range (mg/dL) | Range (mmol/L) |
|-----------------------|---------------|----------------|
| Children <sup>3</sup> | 4.0-7.0       | 1.29-2.26      |
| Adults* <sup>3</sup>  | 2.5-4.5       | 0.81-1.45      |

\* Adult plasma reference intervals are approximately 0.2 mg/dL (0.06 mmol/L) lower than for serum.<sup>16</sup>

NOTE: The system will flag results based only on the serum reference range in the assay file.

### Reference Range (Urine)<sup>17</sup>

| Category               | Range         | Range              |
|------------------------|---------------|--------------------|
| 24 Hour (adults)       | 0.4-1.3 g/day | 12.9-42.0 mmol/day |
| Random (adult males)   | 5-189 mg/dL   | 1.6-61 mmol/L      |
| Random (adult females) | 7-148 mg/dL   | 2.3-48 mmol/L      |

### 24-Hour Urinary Excretion

To convert results from mg/dL to g/day (24-hour urinary excretion)

$$24\text{-hour excretion} = [(V \times c) \div 100\,000] \text{ g/day}$$

Where:

V = 24-hour urine volume (mL)

c = analyte concentration (mg/dL)

To convert results from mmol/L to mmol/day (24-hour urinary excretion)

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Representative performance data are provided in this section.

Results obtained in individual laboratories may vary.

The Alinity c system and the ARCHITECT c System utilize the same reagents and sample/reagent ratios.

Unless otherwise specified, all studies were performed on the Alinity c system.

## Precision

### Within-Laboratory Precision

#### Serum

A study was performed based on guidance from CLSI EP05-A3.<sup>18</sup> Testing was conducted using 3 lots of the Phosphorus2 reagents, 3 lots of the Consolidated Chemistry Calibrator, 1 lot of commercially available controls, and 1 instrument. Two controls and 3 human serum panels were tested in a minimum of 2 replicates, twice per day on 20 days on 3 reagent lot/calibrator lot combinations, where a unique reagent lot and a unique calibrator lot are paired. The performance from a representative combination is shown in the following table.

| Sample          | n  | Mean (mg/dL) | Within-Run (Repeatability) |     | Within-Laboratory <sup>a</sup> |                           |
|-----------------|----|--------------|----------------------------|-----|--------------------------------|---------------------------|
|                 |    |              | SD                         | %CV | SD (Range <sup>b</sup> )       | %CV (Range <sup>b</sup> ) |
| Control Level 1 | 80 | 3.4          | 0.03                       | 0.8 | 0.03 (0.03-0.05)               | 1.0 (1.0-1.4)             |
| Control Level 2 | 80 | 7.4          | 0.03                       | 0.3 | 0.06 (0.06-0.06)               | 0.8 (0.8-0.8)             |
| Panel A         | 80 | 1.1          | 0.02                       | 2.0 | 0.03 (0.02-0.05)               | 2.5 (1.8-4.2)             |
| Panel B         | 80 | 8.9          | 0.04                       | 0.5 | 0.06 (0.06-0.06)               | 0.7 (0.7-0.7)             |
| Panel C         | 80 | 22.8         | 0.13                       | 0.6 | 0.16 (0.14-0.19)               | 0.7 (0.6-0.8)             |

<sup>a</sup> Includes within-run, between-run, and between-day variability.

<sup>b</sup> Minimum and maximum SD or %CV across the 3 reagent lot/calibrator lot combinations.

| Sample          | n  | Mean (mmol/L) | Within-Run (Repeatability) |     | Within-Laboratory <sup>a</sup> |                           |
|-----------------|----|---------------|----------------------------|-----|--------------------------------|---------------------------|
|                 |    |               | SD                         | %CV | SD (Range <sup>b</sup> )       | %CV (Range <sup>b</sup> ) |
| Control Level 1 | 80 | 1.10          | 0.006                      | 0.5 | 0.009 (0.009-0.009)            | 0.8 (0.8-0.8)             |
| Control Level 2 | 80 | 2.38          | 0.008                      | 0.3 | 0.019 (0.015-0.019)            | 0.8 (0.6-0.8)             |
| Panel A         | 80 | 0.36          | 0.005                      | 1.5 | 0.008 (0.007-0.008)            | 2.4 (2.1-2.4)             |
| Panel B         | 80 | 2.88          | 0.009                      | 0.3 | 0.021 (0.017-0.021)            | 0.7 (0.6-0.7)             |
| Panel C         | 80 | 7.37          | 0.044                      | 0.6 | 0.053 (0.046-0.061)            | 0.7 (0.6-0.8)             |

<sup>a</sup> Includes within-run, between-run, and between-day variability.

<sup>b</sup> Minimum and maximum SD or %CV across the 3 reagent lot/calibrator lot combinations.

#### Urine

A study was performed based on guidance from CLSI EP05-A3.<sup>18</sup> Testing was conducted using 3 lots of the Phosphorus2 reagents, 3 lots of the Consolidated Chemistry Calibrator, 1 lot of commercially available controls, and 1 instrument. Two controls and 3 urine panels were tested in a minimum of 2 replicates, twice per day on 20 days on 3 reagent lot/calibrator lot combinations, where a unique reagent lot and a unique calibrator lot are paired. The performance from a representative combination is shown in the following table.

| Sample          | n  | Mean (mg/dL) | Within-Run (Repeatability) |     | Within-Laboratory <sup>a</sup> |                           |
|-----------------|----|--------------|----------------------------|-----|--------------------------------|---------------------------|
|                 |    |              | SD                         | %CV | SD (Range <sup>b</sup> )       | %CV (Range <sup>b</sup> ) |
| Control Level 1 | 80 | 24.8         | 0.12                       | 0.5 | 0.24 (0.24-0.26)               | 0.9 (0.9-1.0)             |
| Control Level 2 | 80 | 48.1         | 0.26                       | 0.5 | 0.31 (0.26-0.41)               | 0.7 (0.5-0.8)             |

| Sample  | n  | Mean<br>(mg/dL) | Within-Run<br>(Repeatability) |     | Within-Laboratory <sup>a</sup> |                              |
|---------|----|-----------------|-------------------------------|-----|--------------------------------|------------------------------|
|         |    |                 | SD                            | %CV | SD<br>(Range <sup>b</sup> )    | %CV<br>(Range <sup>b</sup> ) |
| Panel A | 80 | 5.2             | 0.08                          | 1.5 | 0.21<br>(0.18-0.21)            | 4.1<br>(3.4-4.1)             |
| Panel B | 80 | 92.2            | 0.63                          | 0.7 | 1.09<br>(1.09-2.15)            | 1.2<br>(1.2-2.3)             |
| Panel C | 80 | 161.3           | 0.93                          | 0.6 | 1.11<br>(1.11-1.78)            | 0.7<br>(0.7-1.1)             |

<sup>a</sup> Includes within-run, between-run, and between-day variability.

<sup>b</sup> Minimum and maximum SD or %CV across the 3 reagent lot/calibrator lot combinations.

| Sample          | n  | Mean<br>(mmol/L) | Within-Run<br>(Repeatability) |     | Within-Laboratory <sup>a</sup> |                              |
|-----------------|----|------------------|-------------------------------|-----|--------------------------------|------------------------------|
|                 |    |                  | SD                            | %CV | SD<br>(Range <sup>b</sup> )    | %CV<br>(Range <sup>b</sup> ) |
| Control Level 1 | 80 | 8.02             | 0.037                         | 0.5 | 0.076<br>(0.076-0.082)         | 0.9<br>(0.9-1.0)             |
| Control Level 2 | 80 | 15.52            | 0.083                         | 0.5 | 0.099<br>(0.085-0.132)         | 0.6<br>(0.5-0.8)             |
| Panel A         | 80 | 1.66             | 0.026                         | 1.6 | 0.066<br>(0.057-0.066)         | 4.0<br>(3.4-4.0)             |
| Panel B         | 80 | 29.78            | 0.202                         | 0.7 | 0.351<br>(0.351-0.695)         | 1.2<br>(1.2-2.3)             |
| Panel C         | 80 | 52.10            | 0.301                         | 0.6 | 0.360<br>(0.360-0.574)         | 0.7<br>(0.7-1.1)             |

<sup>a</sup> Includes within-run, between-run, and between-day variability.

<sup>b</sup> Minimum and maximum SD or %CV across the 3 reagent lot/calibrator lot combinations.

#### Reproducibility

##### Serum

A study was performed based on guidance from CLSI EP05-A3.<sup>18</sup> Testing was conducted using 1 lot of the Phosphorus2 reagents, 1 lot of the Consolidated Chemistry Calibrator, 1 lot of commercially available controls, and 3 instruments. Each instrument was operated by a different technician, and each technician prepared an individual sample set. Two controls and 3 human serum panels were tested in a minimum of 3 replicates at 2 separate times per day on 5 different days.

| Sample          | n  | Mean<br>(mg/dL) | Repeatability |     | Within-Laboratory <sup>a</sup> |     | Reproducibility <sup>b</sup> |     |
|-----------------|----|-----------------|---------------|-----|--------------------------------|-----|------------------------------|-----|
|                 |    |                 | SD            | %CV | SD                             | %CV | SD                           | %CV |
| Control Level 1 | 90 | 3.4             | 0.02          | 0.5 | 0.02                           | 0.7 | 0.02                         | 0.7 |
| Control Level 2 | 90 | 7.4             | 0.03          | 0.4 | 0.04                           | 0.5 | 0.04                         | 0.5 |
| Panel A         | 90 | 1.1             | 0.02          | 1.4 | 0.02                           | 1.6 | 0.02                         | 1.7 |
| Panel B         | 90 | 8.9             | 0.03          | 0.4 | 0.04                           | 0.4 | 0.05                         | 0.5 |
| Panel C         | 90 | 22.8            | 0.06          | 0.3 | 0.15                           | 0.6 | 0.15                         | 0.7 |

<sup>a</sup> Includes repeatability (within-run), between-run, and between-day variability.

<sup>b</sup> Includes repeatability (within-run), between-run, between-day, and between-instrument variability.

| Sample          | n  | Mean<br>(mmol/L) | Repeatability |     | Within-Laboratory <sup>a</sup> |     | Reproducibility <sup>b</sup> |     |
|-----------------|----|------------------|---------------|-----|--------------------------------|-----|------------------------------|-----|
|                 |    |                  | SD            | %CV | SD                             | %CV | SD                           | %CV |
| Control Level 1 | 90 | 1.11             | 0.006         | 0.5 | 0.008                          | 0.7 | 0.008                        | 0.7 |
| Control Level 2 | 90 | 2.39             | 0.010         | 0.4 | 0.012                          | 0.5 | 0.013                        | 0.5 |

| Sample  | n  | Mean<br>(mmol/L) | Repeatability |     | Within-Laboratory <sup>a</sup> |     | Reproducibility <sup>b</sup> |     |
|---------|----|------------------|---------------|-----|--------------------------------|-----|------------------------------|-----|
|         |    |                  | SD            | %CV | SD                             | %CV | SD                           | %CV |
| Panel A | 90 | 0.35             | 0.006         | 1.6 | 0.006                          | 1.7 | 0.006                        | 1.8 |
| Panel B | 90 | 2.88             | 0.010         | 0.4 | 0.011                          | 0.4 | 0.016                        | 0.5 |
| Panel C | 90 | 7.36             | 0.020         | 0.3 | 0.047                          | 0.6 | 0.048                        | 0.7 |

<sup>a</sup> Includes repeatability (within-run), between-run, and between-day variability.

<sup>b</sup> Includes repeatability (within-run), between-run, between-day, and between-instrument variability.

##### Urine

A study was performed based on guidance from CLSI EP05-A3.<sup>18</sup> Testing was conducted using 1 lot of the Phosphorus2 reagents, 1 lot of the Consolidated Chemistry Calibrator, 1 lot of commercially available controls, and 3 instruments. Each instrument was operated by a different technician, and each technician prepared an individual sample set. Two controls and 3 urine panels were tested in a minimum of 3 replicates at 2 separate times per day on 5 different days.

| Sample          | n  | Mean<br>(mg/dL) | Repeatability |     | Within-Laboratory <sup>a</sup> |     | Reproducibility <sup>b</sup> |     |
|-----------------|----|-----------------|---------------|-----|--------------------------------|-----|------------------------------|-----|
|                 |    |                 | SD            | %CV | SD                             | %CV | SD                           | %CV |
| Control Level 1 | 90 | 25.6            | 0.25          | 1.0 | 0.26                           | 1.0 | 0.39                         | 1.5 |
| Control Level 2 | 90 | 48.7            | 0.35          | 0.7 | 0.37                           | 0.8 | 0.59                         | 1.2 |
| Panel A         | 90 | 4.8             | 0.16          | 3.2 | 0.18                           | 3.8 | 0.25                         | 5.2 |
| Panel B         | 90 | 91.4            | 0.49          | 0.5 | 0.53                           | 0.6 | 1.38                         | 1.5 |
| Panel C         | 90 | 160.0           | 0.74          | 0.5 | 0.91                           | 0.6 | 2.51                         | 1.6 |

<sup>a</sup> Includes repeatability (within-run), between-run, and between-day variability.

<sup>b</sup> Includes repeatability (within-run), between-run, between-day, and between-instrument variability.

| Sample          | n  | Mean<br>(mmol/L) | Repeatability |     | Within-Laboratory <sup>a</sup> |     | Reproducibility <sup>b</sup> |     |
|-----------------|----|------------------|---------------|-----|--------------------------------|-----|------------------------------|-----|
|                 |    |                  | SD            | %CV | SD                             | %CV | SD                           | %CV |
| Control Level 1 | 90 | 8.26             | 0.080         | 1.0 | 0.085                          | 1.0 | 0.125                        | 1.5 |
| Control Level 2 | 90 | 15.72            | 0.113         | 0.7 | 0.121                          | 0.8 | 0.189                        | 1.2 |
| Panel A         | 90 | 1.56             | 0.050         | 3.2 | 0.059                          | 3.8 | 0.080                        | 5.1 |
| Panel B         | 90 | 29.52            | 0.158         | 0.5 | 0.174                          | 0.6 | 0.446                        | 1.5 |
| Panel C         | 90 | 51.67            | 0.237         | 0.5 | 0.294                          | 0.6 | 0.809                        | 1.6 |

<sup>a</sup> Includes repeatability (within-run), between-run, and between-day variability.

<sup>b</sup> Includes repeatability (within-run), between-run, between-day, and between-instrument variability.

##### Accuracy

A study was performed to estimate the bias of the Phosphorus2 assay relative to standard reference material (NIST SRM 2186I Standard). Testing was conducted using 3 concentrations of the standard across 1 lot of the Phosphorus2 reagents, 1 lot of the Consolidated Chemistry Calibrator, and 1 instrument. The bias ranged from -0.2% to 0.0% for serum and -5.1% to -0.1% for urine.

##### Lower Limits of Measurement

A study was performed based on guidance from CLSI EP17-A2.<sup>19</sup> Testing was conducted using 3 lots of the Phosphorus2 reagents on each of 2 instruments over a minimum of 3 days. The maximum observed limit of blank (LoB), limit of detection (LoD), and limit of quantitation (LoQ) values are summarized below.

**Serum**

|                  | mg/dL | mmol/L |
|------------------|-------|--------|
| LoB <sup>a</sup> | 0.1   | 0.03   |
| LoD <sup>b</sup> | 0.2   | 0.06   |
| LoQ <sup>c</sup> | 0.5   | 0.16   |

<sup>a</sup> The LoB represents the 95th percentile from n ≥ 60 replicates of zero-analyte samples.

<sup>b</sup> The LoD represents the lowest concentration at which the analyte can be detected with 95% probability based on n ≥ 60 replicates of low-analyte level samples.

<sup>c</sup> The LoQ is defined as the lowest concentration at which a maximum allowable precision of 20 %CV was met and was determined from n ≥ 60 replicates of low-analyte level samples.

**Urine**

|                  | mg/dL | mmol/L |
|------------------|-------|--------|
| LoB <sup>a</sup> | 0.5   | 0.16   |
| LoD <sup>b</sup> | 1.0   | 0.32   |
| LoQ <sup>c</sup> | 2.4   | 0.78   |

<sup>a</sup> The LoB represents the 95th percentile from n ≥ 60 replicates of zero-analyte samples.

<sup>b</sup> The LoD represents the lowest concentration at which the analyte can be detected with 95% probability based on n ≥ 60 replicates of low-analyte level samples.

<sup>c</sup> The LoQ presented in the table is in alignment with the LoQ for the Phosphorus2 assay on the ARCHITECT c System. The observed LoQ on the Alinity c system was 2.2 mg/dL (0.71 mmol/L). This LoQ is defined as the lowest concentration at which a maximum allowable precision of 20 %CV was met and was determined from n ≥ 60 replicates of low-analyte level samples.

**Linearity**

A study was performed based on guidance from CLSI EP06-A.<sup>20</sup>

This assay is linear across the analytical measuring interval of 0.5 to 25.3 mg/dL (0.16 to 8.17 mmol/L) for serum, and 2.4 to 186.2 mg/dL (0.78 to 60.14 mmol/L) for urine.

**Analytical Specificity****Interference**

These studies were performed on the ARCHITECT c System.

**Serum****Potentially Interfering Endogenous Substances**

A study was performed based on guidance from CLSI EP07, 3rd ed.<sup>21</sup> Each substance was tested at 2 levels of the analyte (approximately 3.5 mg/dL and 5.0 mg/dL).

**No significant interference (interference within ± 10%)** was observed at the following concentrations.

| No Significant Interference (Interference within ± 10%) |                   |                 |
|---------------------------------------------------------|-------------------|-----------------|
| Potentially Interfering Substance                       | Interferent Level |                 |
|                                                         | Default Units     | Alternate Units |
| Bilirubin - conjugated                                  | 30 mg/dL          | 356 µmol/L      |
| Bilirubin - unconjugated                                | 60 mg/dL          | 1026 µmol/L     |
| Hemoglobin                                              | 125 mg/dL         | 1.25 g/L        |
| Total protein                                           | 15 g/dL           | 150 g/L         |
| Triglycerides                                           | 750 mg/dL         | 8.47 mmol/L     |

**Interference beyond ± 10% (based on 95% Confidence Interval**

**[CI])** was observed at the concentrations shown below for the following substances.

| Interference beyond ± 10% (based on 95% Confidence Interval [CI]) |                   |                 |               |                 |                         |
|-------------------------------------------------------------------|-------------------|-----------------|---------------|-----------------|-------------------------|
| Potentially Interfering Substance                                 | Interferent Level |                 | Analyte Level |                 | % Interference (95% CI) |
|                                                                   | Default Units     | Alternate Units | Default Units | Alternate Units |                         |
| Bilirubin - conjugated                                            | 40 mg/dL          | 474 µmol/L      | 3.5 mg/dL     | 1.13 mmol/L     | -12% (-13%, -11%)       |
| Hemoglobin                                                        | 250 mg/dL         | 2.50 g/L        | 3.5 mg/dL     | 1.13 mmol/L     | 14% (14%, 14%)          |
| Triglycerides                                                     | 1500 mg/dL        | 16.9 mmol/L     | 3.5 mg/dL     | 1.13 mmol/L     | 15% (14%, 16%)          |
| Triglycerides                                                     | 1500 mg/dL        | 16.9 mmol/L     | 5.0 mg/dL     | 1.62 mmol/L     | 10% (9%, 11%)           |

**Potentially Interfering Exogenous Substances**

A study was performed based on guidance from CLSI EP07, 3rd ed.<sup>21</sup> Each substance was tested at 2 levels of the analyte (approximately 3.5 mg/dL and 5.0 mg/dL).

**No significant interference (interference within ± 10%)** was observed at the following concentrations.

| No Significant Interference (Interference within ± 10%) |                   |                 |
|---------------------------------------------------------|-------------------|-----------------|
| Potentially Interfering Substance                       | Interferent Level |                 |
|                                                         | Default Units     | Alternate Units |
| Acetaminophen                                           | 160 mg/L          | 1059 µmol/L     |
| Acetylcysteine                                          | 150 mg/L          | 920 µmol/L      |
| Acetylsalicylic acid                                    | 30 mg/L           | 167 µmol/L      |
| Aminosalicylic acid (p-aminosalicylic acid)             | 47 mg/dL          | 3074 µmol/L     |
| Amiodarone                                              | 5 mg/dL           | 77.5 µmol/L     |
| Ampicillin-Na                                           | 80 mg/L           | 215 µmol/L      |
| Ascorbic acid                                           | 60 mg/L           | 341 µmol/L      |
| Biotin                                                  | 4250 ng/mL        | 17.4 µmol/L     |
| Ca-dobesilate                                           | 60 mg/L           | 143 µmol/L      |
| Cefotaxime                                              | 53 mg/dL          | 1166 µmol/L     |
| Cefoxitin                                               | 6600 mg/L         | 15.4 mmol/L     |
| Cyclosporine                                            | 2 mg/L            | 1.66 µmol/L     |
| Desacetylcefotaxime                                     | 6 mg/dL           | 145 µmol/L      |
| Doxycycline                                             | 20 mg/L           | 45.0 µmol/L     |
| Ibuprofen                                               | 220 mg/L          | 1067 µmol/L     |
| Levodopa                                                | 8 mg/L            | 40.6 µmol/L     |
| Mannitol                                                | 4 g/L             | 22.0 mmol/L     |
| Methicillin (sodium)                                    | 8 mg/dL           | 198 µmol/L      |
| Methotrexate                                            | 140 mg/dL         | 3080 µmol/L     |
| Methylidopa                                             | 25 mg/L           | 118 µmol/L      |
| Metronidazole                                           | 130 mg/L          | 759 µmol/L      |
| Naproxen                                                | 36 mg/dL          | 1562 µmol/L     |
| Nitroglycerin                                           | 0.015 mg/L        | 0.0660 µmol/L   |
| Omeprazole                                              | 0.9 mg/dL         | 26.1 µmol/L     |
| Phenylbutazone                                          | 330 mg/L          | 1069 µmol/L     |
| Primidone                                               | 6 mg/dL           | 275 µmol/L      |
| Promethazine                                            | 0.03 mg/dL        | 1.06 µmol/L     |
| Ranitidine                                              | 2 mg/dL           | 63.6 µmol/L     |
| Rifampicin                                              | 50 mg/L           | 61.0 µmol/L     |
| Sodium heparin                                          | 4 U/mL            | NA              |
| Theophylline (1,3-dimethylxanthine)                     | 60 mg/L           | 333 µmol/L      |

NA = Not applicable

**Interference beyond  $\pm 10\%$  (based on 95% Confidence Interval [CI])** was observed at the concentrations shown below for the following substances.

| Potentially Interfering Substance | Interference beyond $\pm 10\%$ (based on 95% Confidence Interval [CI]) |                 |               |                 |                         |
|-----------------------------------|------------------------------------------------------------------------|-----------------|---------------|-----------------|-------------------------|
|                                   | Interferent Level                                                      |                 | Analyte Level |                 | % Interference (95% CI) |
|                                   | Default Units                                                          | Alternate Units | Default Units | Alternate Units |                         |
| Mannitol                          | 5 g/L                                                                  | 27.5 mmol/L     | 3.5 mg/dL     | 1.13 mmol/L     | -11%                    |
|                                   |                                                                        |                 |               |                 | (-11%, -10%)            |
| Mannitol                          | 10 g/L                                                                 | 54.9 mmol/L     | 5.0 mg/dL     | 1.62 mmol/L     | -10%                    |
|                                   |                                                                        |                 |               |                 | (-11%, -10%)            |

## Urine

### Potentially Interfering Endogenous Substances

A study was performed based on guidance from CLSI EP07, 3rd ed.<sup>21</sup> Each substance was tested at 2 levels of the analyte (approximately 30 mg/dL and 90 mg/dL).

**No significant interference (interference within  $\pm 10\%$ )** was observed at the following concentrations.

| Potentially Interfering Substance | No Significant Interference (Interference within $\pm 10\%$ ) |                  |
|-----------------------------------|---------------------------------------------------------------|------------------|
|                                   | Interferent Level                                             |                  |
|                                   | Default Units                                                 | Alternate Units  |
| Acetic acid (8.5N)                | 6.25 mL/dL                                                    | 531 mmol/L       |
| Ascorbate                         | 220 mg/dL                                                     | 12.6 mmol/L      |
| Boric acid                        | 250 mg/dL                                                     | 40.5 mmol/L      |
| Glucose                           | 1000 mg/dL                                                    | 55.5 mmol/L      |
| Hydrochloric acid (6N)            | 2.5 mL/dL                                                     | 150 mmol/L       |
| Nitric acid (6N)                  | 5.0 mL/dL                                                     | 300 mmol/L       |
| Protein                           | 50 mg/dL                                                      | 0.500 g/L        |
| Sodium carbonate                  | 1.25 g/dL                                                     | 118 mmol/L       |
| Sodium fluoride                   | 400 mg/dL                                                     | 95.2 mmol/L      |
| Sodium oxalate                    | 60 mg/dL                                                      | 4476 $\mu$ mol/L |

### Potentially Interfering Exogenous Substances

A study was performed based on guidance from CLSI EP07, 3rd ed.<sup>21</sup> Each substance was tested at 2 levels of the analyte (approximately 30 mg/dL and 90 mg/dL).

**No significant interference (interference within  $\pm 10\%$ )** was observed at the following concentrations.

| Potentially Interfering Substance | No Significant Interference (Interference within $\pm 10\%$ ) |                  |
|-----------------------------------|---------------------------------------------------------------|------------------|
|                                   | Interferent Level                                             |                  |
|                                   | Default Units                                                 | Alternate Units  |
| Acetaminophen                     | 16 mg/dL                                                      | 1059 $\mu$ mol/L |
| Acetylcysteine                    | 15 mg/dL                                                      | 920 $\mu$ mol/L  |
| Biotin                            | 4250 ng/mL                                                    | 17.4 $\mu$ mol/L |
| Ibuprofen                         | 22 mg/dL                                                      | 1067 $\mu$ mol/L |
| Mannitol                          | 18 g/L                                                        | 98.8 mmol/L      |

Interferences from medication or endogenous substances may affect results.<sup>22</sup>

## Method Comparison

A study was performed based on guidance from CLSI EP09c, 3rd ed.<sup>23</sup> using the Passing-Bablok regression method.

| Phosphorus2 on the Alinity c system vs. Phosphorus on the ARCHITECT c System |     |                   |                         |                 |       |                           |
|------------------------------------------------------------------------------|-----|-------------------|-------------------------|-----------------|-------|---------------------------|
|                                                                              | n   | Units             | Correlation Coefficient | Intercept       | Slope | Concentration Range       |
| Serum                                                                        | 108 | mg/dL<br>(mmol/L) | 1.00                    | 0.0             | 1.05  | 1.0-22.8<br>(0.34-7.36)   |
| Urine                                                                        | 110 | mg/dL<br>(mmol/L) | 1.00                    | -0.6<br>(-0.21) | 1.05  | 6.2-167.9<br>(2.00-54.22) |

| Phosphorus2 on the Alinity c system vs. Phosphorus2 on the ARCHITECT c System |     |                   |                         |                 |                |                           |
|-------------------------------------------------------------------------------|-----|-------------------|-------------------------|-----------------|----------------|---------------------------|
|                                                                               | n   | Units             | Correlation Coefficient | Intercept       | Slope          | Concentration Range       |
| Serum                                                                         | 108 | mg/dL<br>(mmol/L) | 1.00                    | 0.0<br>(-0.02)  | 1.00<br>(1.01) | 1.1-23.7<br>(0.36-7.67)   |
| Urine                                                                         | 111 | mg/dL<br>(mmol/L) | 1.00                    | -0.3<br>(-0.08) | 1.01           | 4.4-174.1<br>(1.42-56.23) |

## BIBLIOGRAPHY

- Goldman L, Schafer AL, editors. *Goldman-Cecil Medicine*. 25th ed. Elsevier/Saunders; 2016.
- McPherson RA, Pincus MR, editors. *Henry's Clinical Diagnosis and Management by Laboratory Methods*. 23rd ed. St. Louis, Missouri: Elsevier; 2017.
- Burtis CA, Bruns DE, editors. *Tietz Fundamentals of Clinical Chemistry and Molecular Diagnostics*. 7th ed. St. Louis, MO: Saunders Elsevier; 2015.
- Dominguez JR, Kestenbaum B, Chonchol M, et al. Relationships between serum and urine phosphorus with all-cause and cardiovascular mortality: the osteoporotic fractures in men (MrOS) study. *Am J Kidney Dis* 2013;61(4):555-563.
- Burtis CA, Ashwood ER, editors. *Tietz Fundamentals of Clinical Chemistry*. 5th ed. WB Saunders; 2001:chap 38.
- US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
- US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 6th ed. Washington, DC: US Government Printing Office; June 2020.
- World Health Organization. *Laboratory Biosafety Manual*. 4th ed. Geneva: World Health Organization; 2020.
- Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.
- Rifai N, Horvath AR, Wittwer CT editors. *Tietz Textbook of Clinical Chemistry and Molecular Diagnostics*. 6th ed. Saunders Elsevier; 2018:1444.
- Kachhawa K, Kachhawa P, Varma M, et al. Study of the stability of various biochemical analytes in samples stored at different predefined storage conditions at an accredited laboratory of India. *J Lab Physicians* 2017;9(1):11-15.
- Clinical and Laboratory Standards Institute (CLSI). *Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests; Approved Guideline—Fourth Edition*. CLSI Document GP44-A4. Wayne, PA: CLSI; 2010.
- Westgard JO. *Basic QC Practices; Training in Statistical Quality Control for Medical Laboratories*. 4th ed. Westgard QC, Inc.; 2016.
- Clinical and Laboratory Standards Institute (CLSI). *Establishing and Verifying an Extended Measuring Interval Through Specimen Dilution and Spiking*. 1st ed. CLSI Guideline EP34. Wayne, PA: CLSI; 2018.
- Valentina Molinaris, Mario G. Bianchetti, Gregorio P. Milani, et al. Interferences in the measurement of circulating phosphate: a literature review. *Clin Chem Lab Med* 2020; 58(12): 1971–1977.
- Rifai N, Horvath AR, Wittwer CT, editors. *Tietz Textbook of Clinical Chemistry and Molecular Diagnostics*. 6th ed. St. Louis, MO: Elsevier; 2018.
- Wu AHB, editor. *Tietz Clinical Guide to Laboratory Tests*. 4th ed. St. Louis, MO: Elsevier Saunders; 2006.
- Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Precision of Quantitative Measurement Procedures: Approved Guideline—Third Edition*. CLSI Document EP05-A3. Wayne, PA: CLSI; 2014.
- Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition*. CLSI Document EP17-A2. Wayne, PA: CLSI; 2012.
- Clinical and Laboratory Standards Institute (CLSI). *Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline*. CLSI Document EP06-A. Wayne, PA: CLSI; 2003.
- Clinical and Laboratory Standards Institute (CLSI). *Interference Testing in Clinical Chemistry*. 3rd ed. CLSI Guideline EP07. Wayne, PA: CLSI; 2018.
- Young DS. Laboratory test listings. In: *Effects of Drugs on Clinical Laboratory Tests*. 5th ed. AACC Press; 2000:chap 3.
- Clinical and Laboratory Standards Institute (CLSI). *Measurement Procedure Comparison and Bias Estimation Using Patient Samples*. 3rd ed. CLSI Guideline EP09c. Wayne, PA: CLSI; 2018.

## ■ Key to Symbols

### ISO 15223 Symbols

|                                                                                   |                                           |
|-----------------------------------------------------------------------------------|-------------------------------------------|
|  | Consult instructions for use              |
|  | Manufacturer                              |
|  | Sufficient for                            |
|  | Temperature limitation                    |
|  | Use by/Expiration date                    |
| <b>IVD</b>                                                                        | <i>In Vitro</i> Diagnostic Medical Device |
| <b>LOT</b>                                                                        | Lot Number                                |
| <b>REF</b>                                                                        | List Number                               |
| <b>SN</b>                                                                         | Serial number                             |

### Other Symbols

|                                  |                                                                                         |
|----------------------------------|-----------------------------------------------------------------------------------------|
| <b>DISTRIBUTED IN THE USA BY</b> | Distributed in the USA by                                                               |
| <b>FOR USE WITH</b>              | Identifies products to be used together                                                 |
| <b>INFORMATION FOR USA ONLY</b>  | Information needed for United States of America only                                    |
| <b>PRODUCT OF IRELAND</b>        | Product of Ireland                                                                      |
| <b>R1</b>                        | Reagent 1                                                                               |
| <b>R2</b>                        | Reagent 2                                                                               |
| <b>Rx ONLY</b>                   | For use by or on the order of a physician only (applicable to USA classification only). |

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For customers in the European Union: if, in the course of using this device, you have reason to believe that a serious incident has occurred, report it to the manufacturer and to your national authority.

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## BIBLIOGRAFIE

1. US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
2. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
3. World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
4. Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.

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Revised March 2022.

REF 08P4320

Instructions must be carefully followed. Reliability of assay results cannot be guaranteed if there are any deviations from these instructions.

**WARNING:** The Alinity c Hemoglobin A1c assay has significant interference with the fetal hemoglobin (HbF). Hemoglobin A1c results are invalid for patients with abnormal amounts of HbF, including those with known Hereditary Persistence of Fetal Hemoglobin. For more information regarding the specific concentrations of HbF that were found to interfere with the Alinity c Hemoglobin A1c assay, refer to the SPECIFIC PERFORMANCE CHARACTERISTICS, Analytical Specificity section of this package insert.

## NAME

Alinity c Hemoglobin A1c Reagent Kit (also referred to as HbA1c)

## INTENDED USE

The Alinity c Hemoglobin A1c assay is used in clinical laboratories for the quantitative *in vitro* measurement of percent hemoglobin A1c (NGSP) or HbA1c fraction mmol/mol (IFCC) in human whole blood and hemolysate on the Alinity c analyzer.

Hemoglobin A1c measurements are used as an aid in the diagnosis of diabetes mellitus, as an aid to identify patients who may be at risk for developing diabetes mellitus, and for the monitoring of long-term blood glucose control in individuals with diabetes mellitus.

## SUMMARY AND EXPLANATION OF THE TEST

HbA1c is the fraction of hemoglobin A that is first reversibly, then irreversibly glycosylated at one or both *N*-terminal valines of the  $\beta$ -chain.<sup>1</sup> The longer red blood cells are in circulation and the higher the ambient glucose levels, the higher the concentration of HbA1c. HbA1c reflects the average blood glucose level during the preceding 2 to 3 months. The HbA1c assay is useful as an aid in the:

- diagnosis of diabetes mellitus,
- identification of patients at risk for developing diabetes, and
- monitoring of patients with diabetes mellitus.<sup>2-6</sup>

For monitoring diabetic patients, it is recommended that glycemic goals are individualized following current professional society recommendations.<sup>7</sup> As recommended by the American Diabetes Association (ADA), patients in the range of 5.7 to 6.4 %HbA1c (39 to 46 mmol/mol) would be in the category of increased risk for diabetes and results  $\geq 6.5\%$  (48 mmol/mol) may aid in the diagnosis of diabetes.<sup>7</sup>

Several studies, including the Diabetes Control and Complications Trial (DCCT), have shown that long-term control of diabetes can prevent complications such as cardiovascular disease, retinopathy, nephropathy, and neuropathy. Measurement of HbA1c can be invaluable in the monitoring of glycemic control of diabetic patients.<sup>8-10</sup>

This method is certified to the National Glycohemoglobin Standardization Program (NGSP), standardized to International Federation of Clinical Chemistry and Laboratory Medicine (IFCC), and traceable to DCCT.

## PRINCIPLES OF THE PROCEDURE

The Alinity c Hemoglobin A1c assay consists of two separate concentration measurements: glycosylated hemoglobin (HbA1c) and total hemoglobin (THb). The two concentrations are used to determine the percent HbA1c (NGSP units) or the hemoglobin fraction in mmol/mol (IFCC units).

The individual concentration values of HbA1c and THb generated by the Alinity c Hemoglobin A1c assay are used only for calculating the percent hemoglobin A1c or HbA1c fraction, and must not be used individually for diagnostic purposes.

The anticoagulated whole blood specimen is lysed automatically on the system for the Whole Blood application or may be lysed manually using the Hemoglobin A1c Diluent (A1cDIL) for the Hemolysate application.

### Glycosylated Hemoglobin (HbA1c)

The Alinity c Hemoglobin A1c assay utilizes an enzymatic method that specifically measures *N*-terminal fructosyl dipeptides of the  $\beta$ -chain of HbA1c.

- In the pretreatment process, the erythrocytes are lysed and the hemoglobin is transformed to methemoglobin by reaction with sodium nitrite.
- With the addition of Reagent 1 (R1) to the sample, the glycosylated *N*-terminal dipeptide (fructosyl-VH) of the  $\beta$ -chain of hemoglobin is cleaved by the action of protease. The hemoglobin is transformed to stable methemoglobin azide by the action of sodium azide and the concentration of the hemoglobin is determined by measuring absorbance.
- Addition of Reagent 2 (R2) starts a reaction and fructosyl peptide oxidase (FPOX) is allowed to react with fructosyl-VH. The HbA1c concentration is measured by determining the resultant hydrogen peroxide.

### Total Hemoglobin (THb)

The hemoglobin is oxidized to stable methemoglobin azide by the action of sodium nitrite and sodium azide and the concentration of the hemoglobin is determined by measuring absorbance (sample + R1).

### Hemoglobin A1c Calculations<sup>11</sup>

The final result is expressed as %HbA1c (NGSP) or mmol/mol HbA1c (IFCC) and is automatically calculated by the system from the HbA1c/THb ratio as follows:

#### mmol/mol HbA1c IFCC:

$$\text{HbA1c (mmol/mol)} = (\text{HbA1c/THb}) \times 1000$$

#### %HbA1c DCCT/NGSP:

$$\text{HbA1c (\%)} = \text{IFCC} \times 0.09148 + 2.152$$

Methodology: Enzymatic

For additional information on system and assay technology, refer to the Alinity ci-series Operations Manual, Section 3.

## REAGENTS

### Kit Contents

Alinity c Hemoglobin A1c Reagent Kit 08P43

NOTE: This product is composed of 3 components, which are packaged as a 2 cartridge reagent set. Both cartridges are required to perform the assay.

Volumes (mL) listed in the table below indicate the volume per cartridge set.

|                                                                                                                                                                                                                                                                                      |                |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------|
| <b>REF</b>                                                                                                                                                                                                                                                                           | <b>08P4320</b> |
| Tests per cartridge set                                                                                                                                                                                                                                                              | 260            |
| Number of cartridge sets per kit                                                                                                                                                                                                                                                     | 5              |
| Tests per kit                                                                                                                                                                                                                                                                        | 1300           |
| <b>R1</b>                                                                                                                                                                                                                                                                            | 46.0 mL        |
| <b>R2</b>                                                                                                                                                                                                                                                                            | 18.3 mL        |
| <b>A1cDIL</b>                                                                                                                                                                                                                                                                        | 67.6 mL        |
| <b>R1</b> Active ingredients: 10-(carboxymethylaminocarbonyl)-3,7-bis(dimethylamino) phenothiazine sodium salt (0.000817%), protease (bacterial) (< 1 MU/dL). Inactive ingredient: sodium azide (< 0.1%) as a stabilizer. Preservatives: sodium azide (< 0.1%), ProClin 300 (0.05%). |                |
| <b>R2</b> Active ingredients: peroxidase (horseradish) (5 to 15 kU/dL), fructosyl-peptide-oxidase ( <i>E. coli</i> , recombinant) (300 to 900 U/dL). Preservative: ofloxacin (0.001%).                                                                                               |                |
| <b>A1cDIL</b> Active ingredient: sodium nitrite (> 0.05 to < 0.3%). Preservative: ProClin 300 (0.01%).                                                                                                                                                                               |                |

### Warnings and Precautions

- **IVD**
- For *In Vitro* Diagnostic Use
- **Rx ONLY**

### Safety Precautions

**CAUTION:** This product requires the handling of human specimens. It is recommended that all human-sourced materials be considered potentially infectious and handled in accordance with the OSHA Standard on Bloodborne Pathogens. Biosafety Level 2 or other appropriate biosafety practices should be used for materials that contain or are suspected of containing infectious agents.<sup>12-15</sup>

|                                                                                     |                                                                                                                     |
|-------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------|
| The following warnings and precautions apply to: <b>R1</b>                          |                                                                                                                     |
|  |                                                                                                                     |
| <b>DANGER</b>                                                                       | Contains morpholinoethanesulfonic acid, monohydrate*, N,N-dimethylformamide, methylisothiazolones and sodium azide. |
| H317                                                                                | May cause an allergic skin reaction.                                                                                |
| H360                                                                                | May damage fertility or the unborn child.                                                                           |
| H316*                                                                               | Causes mild skin irritation.                                                                                        |
| EUH032                                                                              | Contact with acids liberates very toxic gas.                                                                        |
| <b>Prevention</b>                                                                   |                                                                                                                     |
| P201                                                                                | Obtain special instructions before use.                                                                             |
| P261                                                                                | Avoid breathing mist / vapors / spray.                                                                              |
| P280                                                                                | Wear protective gloves / protective clothing / eye protection.                                                      |
| P272                                                                                | Contaminated work clothing should not be allowed out of the workplace.                                              |
| <b>Response</b>                                                                     |                                                                                                                     |
| P302+P352                                                                           | IF ON SKIN: Wash with plenty of water.                                                                              |
| P333+P313                                                                           | If skin irritation or rash occurs: Get medical advice / attention.                                                  |
| P362+P364                                                                           | Take off contaminated clothing and wash it before reuse.                                                            |
| P308+P313                                                                           | IF exposed or concerned: Get medical advice / attention.                                                            |
| <b>Disposal</b>                                                                     |                                                                                                                     |
| P501                                                                                | Dispose of contents / container in accordance with local regulations.                                               |

\* Not applicable where regulation EC 1272/2008 (CLP) or OSHA Hazard Communication 29 CFR 1910.1200 (HCS) 2012 have been implemented.

|                                                            |                                                            |
|------------------------------------------------------------|------------------------------------------------------------|
| The following warnings and precautions apply to: <b>R2</b> |                                                            |
| H316*                                                      | Causes mild skin irritation.                               |
| P332+P313*                                                 | If skin irritation occurs: Get medical advice / attention. |

\* Not applicable where regulation EC 1272/2008 (CLP) or OSHA Hazard Communication 29 CFR 1910.1200 (HCS) 2012 have been implemented.

|                                                                                   |                                                                        |
|-----------------------------------------------------------------------------------|------------------------------------------------------------------------|
| The following warnings and precautions apply to: <b>A1cDIL</b>                    |                                                                        |
|  |                                                                        |
| <b>WARNING</b>                                                                    | Contains sodium nitrite*, maleic acid and methylisothiazolones.        |
| H317                                                                              | May cause an allergic skin reaction.                                   |
| H402*                                                                             | Harmful to aquatic life.                                               |
| <b>Prevention</b>                                                                 |                                                                        |
| P261                                                                              | Avoid breathing mist / vapors / spray.                                 |
| P272                                                                              | Contaminated work clothing should not be allowed out of the workplace. |
| P273*                                                                             | Avoid release to the environment.                                      |
| P280                                                                              | Wear protective gloves / protective clothing / eye protection.         |
| <b>Response</b>                                                                   |                                                                        |
| P302+P352                                                                         | IF ON SKIN: Wash with plenty of water.                                 |
| P333+P313                                                                         | If skin irritation or rash occurs: Get medical advice / attention.     |
| P362+P364                                                                         | Take off contaminated clothing and wash it before reuse.               |
| <b>Disposal</b>                                                                   |                                                                        |
| P501                                                                              | Dispose of contents / container in accordance with local regulations.  |

\* Not applicable where regulation EC 1272/2008 (CLP) or OSHA Hazard Communication 29 CFR 1910.1200 (HCS) 2012 have been implemented.

Follow local chemical disposal regulations based on your location along with recommendations and content in the Safety Data Sheet to determine the safe disposal of this product.

For the most current hazard information, see the product Safety Data Sheet.

Safety Data Sheets are available at [www.corelaboratory.abbott](http://www.corelaboratory.abbott) or contact your local representative.

For a detailed discussion of safety precautions during system operation, refer to the Alinity ci-series Operations Manual, Section 8.

### Reagent Handling

- Reagents are shipped refrigerated or on cold packs.
- Upon receipt, place reagent cartridges in an upright position for 8 hours before use to allow bubbles that may have formed to dissipate.
- If a reagent cartridge is dropped, place in an upright position for 8 hours before use to allow bubbles that may have formed to dissipate.
- Reagents are susceptible to the formation of foam and bubbles. Bubbles may interfere with the detection of the reagent level in the cartridge and cause insufficient reagent aspiration that may adversely affect results.

For a detailed discussion of reagent handling precautions during system operation, refer to the Alinity ci-series Operations Manual, Section 7.

## Reagent Storage

|                 | Storage Temperature | Maximum Storage Time  | Additional Storage Instructions                                                                                                                                         |
|-----------------|---------------------|-----------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Unopened</b> | 2 to 8°C            | Until expiration date | Store in upright position.                                                                                                                                              |
| <b>Onboard</b>  | System Temperature  | 50 days               |                                                                                                                                                                         |
| <b>Opened</b>   | 2 to 8°C            | Until expiration date | Store in upright position. Do not reuse original reagent caps or replacement caps due to the risk of contamination and the potential to compromise reagent performance. |

Reagents may be stored on or off the system. If removed from the system, store reagents with new replacement caps in an upright position at 2 to 8°C. For reagents stored off the system, it is recommended that they be stored in their original trays or boxes to ensure they remain upright.

For information on unloading reagents, refer to the Alinity ci-series Operations Manual, Section 5.

**A1cDIL** can be replaced independently of **R1** and **R2**.

## Indications of Reagent Deterioration

Deterioration of the reagents may be indicated when a calibration error occurs or a control value is out of the specified range.

Associated test results are invalid, and samples must be retested. Assay recalibration may be necessary.

For troubleshooting information, refer to the Alinity ci-series Operations Manual, Section 10.

## INSTRUMENT PROCEDURE

The Alinity c Hemoglobin A1c assay files must be installed on the Alinity c analyzer prior to performing the assay.

For the Whole Blood application, the following assay files are required:

- THbWB
- HbA1cWB
- A1cWB (for results in mmol/mol (IFCC) units) or %A1cWB (for results in %HbA1c (NGSP) units)

For the Hemolysate application, the following assay files are required:

- THbH
- HbA1cH
- A1cH (for results in mmol/mol (IFCC) units) or %A1cH (for results in %HbA1c (NGSP) units)

The Alinity ci-series software version 2.0 or higher must be installed prior to performing the assay.

For detailed information on assay file installation and viewing and editing assay parameters, refer to the Alinity ci-series Operations Manual, Section 2.

For information on printing assay parameters, refer to the Alinity ci-series Operations Manual, Section 5.

For a detailed description of system procedures, refer to the Alinity ci-series Operations Manual.

## SPECIMEN COLLECTION AND PREPARATION FOR ANALYSIS

### Specimen Types

The specimen type used for this assay is whole blood only.

The specimen type listed below was verified for use with this assay.

Other specimen types and anticoagulants have not been verified with this assay.

| Specimen Type | Collection Vessel                                                                                                                                                    | Special Conditions                                                                                                                                                                        |
|---------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Whole Blood   | Collection tubes<br>Acceptable anticoagulants are:<br>Dipotassium EDTA<br>Lithium heparin<br>Sodium heparin<br>Sodium fluoride/dipotassium EDTA<br>Tripotassium EDTA | Do not overfill specimen collection tubes. Whole blood samples greater than 78 mm in height from the bottom of tube will result in an instrument error and results will not be generated. |

- Use only whole blood specimens collected by standard venipuncture techniques into plastic tubes.
- The instrument does not provide the capability to verify specimen types. It is the responsibility of the operator to verify that the correct specimen types are used in the assay.

### Specimen Conditions

- Analyze fresh specimens if possible.
- Do not overfill specimen collection tubes. Whole blood samples greater than 78 mm in height from the bottom of tube will result in an instrument error and results will not be generated. Refer to the Alinity ci-series Operations Manual, Section 10.
- To prevent cross contamination, use of disposable pipettes or pipette tips is recommended.

### Preparation for Analysis

- **Do not centrifuge samples.**
- Follow the tube manufacturer's collection instructions for specimen collection tubes.
- For testing whole blood samples less than 600 µL, use 12 x 75 mm polypropylene conical bottom tubes.
- Visually inspect the specimens. If fibrin clots or particulate matter is observed, remove with a clean applicator stick.
- Mix all specimens thoroughly by low speed vortexing or gently inverting 10 times prior to loading onto the Alinity c analyzer.

Prepare frozen whole blood specimens as follows:

- Frozen specimens must be completely thawed before mixing.
- Mix thawed specimens thoroughly by low speed vortex or by inverting 10 times.
- Visually inspect the specimens. If layering or stratification is observed, mix until specimens are visibly homogeneous.
- If specimens are not mixed thoroughly, inconsistent results may be obtained.

Do not use Alinity c Sample Cups for whole blood samples. Refer to the Assay Procedure section of this package insert for further information.

### Specimen Storage

Analyze fresh specimens if possible.

| Specimen Type | Temperature      | Maximum Storage Time | Special Instructions               |
|---------------|------------------|----------------------|------------------------------------|
| Whole Blood   | Room temperature | 8 hours              |                                    |
|               | 2 to 8°C         | 7 days               |                                    |
| Hemolysate    | Room temperature | 4 hours              | Do not freeze hemolyzed specimens. |
|               | 2 to 8°C         | 24 hours             |                                    |

**Whole Blood:** If testing will be delayed more than 7 days, store at -70°C or colder.

**CAUTION: Whole blood specimens that require freezing must be stored at -70°C or colder.**

Avoid more than one freeze/thaw cycle.

Stored specimens must be inspected for particulates. If present, mix with a low speed vortex or by inversion.

## Specimen Shipping

Package and label specimens in compliance with applicable state, federal, and international regulations covering the transport of clinical specimens and infectious substances.

Do not exceed the storage limitations listed above.

## PROCEDURE

### Materials Provided

08P43 Alinity c Hemoglobin A1c Reagent Kit

### Materials Required but not Provided

- Alinity c Hemoglobin A1c assay files
- 08P4301 Alinity c Hemoglobin A1c Calibrators

Whole Blood application:

- Commercially available whole blood controls
- 12 x 75 mm polypropylene conical bottom tubes

Hemolysate application:

- 08P4310 Alinity c Hemoglobin A1c Controls or other commercially available controls
- Calibrated adjustable pipette capable of measuring 222 µL
- Calibrated micropipette capable of measuring 10 µL
- Vortex (optional)

For information on materials required for operation of the instrument, refer to the Alinity ci-series Operations Manual, Section 1.

For information on materials required for maintenance procedures, refer to the Alinity ci-series Operations Manual, Section 9.

### Assay Procedure

For a detailed description of how to run an assay, refer to the Alinity ci-series Operations Manual, Section 5.

Refer to the following table.

**NOTE:** No more than one replicate can be sampled from a sample cup or tube.

For information on ordering patient samples and controls, refer to the Alinity ci-series Operations Manual, Section 5.

|                                                                                           | Sample Type                              |                                         |
|-------------------------------------------------------------------------------------------|------------------------------------------|-----------------------------------------|
|                                                                                           | Whole Blood                              | Hemolysate                              |
| 1. Order the application:                                                                 | Whole Blood                              | Hemolysate                              |
| 2. Prepare samples according to the instructions from this section of the package insert: | Assay Procedure, Whole Blood Preparation | Assay Procedure, Hemolysate Preparation |
| 3. Load, then run:                                                                        | Whole Blood samples                      | Prepared, Hemolyzed samples             |
| 4. Report results in percent (NGSP) or mmol/mol (IFCC) hemoglobin A1c.                    |                                          |                                         |

### Whole Blood Preparation

- WARNING: Do not use Alinity c sample cups.**
- WARNING: Do not overfill specimen collection tubes. Whole blood samples greater than 78 mm in height from the bottom of tube will result in an instrument error and results will not be generated. Refer to Section 10 of the Alinity ci-series Operations Manual.**
- Select the appropriate sample vessel using the following table:

| If sample is:                     | Then:                                                                                                                                          |
|-----------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------|
| < 200 µL                          | Do not use the Whole Blood application. Follow the instructions in the Assay Procedure, Hemolysate Preparation section of this package insert. |
| 200 µL to 600 µL                  | Use 12 x 75 mm polypropylene conical bottom tubes only.                                                                                        |
| > 600 µL                          | Use suitable tubes or polypropylene conical bottom tubes.                                                                                      |
| > 78 mm in height in a large tube | Pipette 600 µL of sample into a suitable tube or polypropylene conical bottom tube.                                                            |

### Hemolysate Preparation

- The minimum sample volume requirement is 150 µL for Alinity c sample cups.
  - NOTE:** This volume supersedes the calculated minimum sample volume displayed on the User Interface.

Prepare the hemolysate samples as follows:

- Using a calibrated pipette, dispense 222 µL **A1cDIL** into a tube or sample cup.
- Using a calibrated micropipette, withdraw 10 µL of the well-mixed whole blood patient specimen.
- Wipe excess blood from the exterior of the pipette to ensure accurate transfer of the sample.
- Insert the pipette into the tube or sample cup containing the **A1cDIL** allowing the tip of the pipette to just make contact with the surface of the **A1cDIL** and dispense the 10 µL sample (1:23.2 dilution).
- Withdraw and dispense twice to rinse the pipette, always keeping the tip of the pipette in contact with the fluid in the tube.
- Mix hemolysate thoroughly by low speed vortexing or by gently inverting 10 times. Avoid foaming.
- Allow the hemolysate to stand for a minimum of 1 minute at room temperature prior to testing.
- If the hemolysate is prepared in a tube, transfer to a sample cup and place the cup on the instrument.

**NOTE:** The number of tests per kit is based on the 222 µL of **A1cDIL** and 10 µL of specimen volumes stated in steps 1 and 2 above. However, alternate volumes may be used for the 1:23.2 dilution, such as 555 µL of **A1cDIL** and 25 µL of specimen.

- If using primary or aliquot tubes, refer to the Alinity ci-series Operations Manual, Section 4 to ensure sufficient specimen is present.
- To minimize the effects of evaporation, verify adequate sample cup volume is present prior to running the test.
- Do not use Alinity c Sample Cups for whole blood samples.
- Minimum sample volume requirements:
  - Sample volume for single test: 200 µL (Whole Blood), 150 µL (Hemolysate).

**NOTE:** This amount does not include the dead volume plus the additional over-aspiration volume. For total sample volume requirements, refer to the Alinity ci-series Operations Manual, Section 4.
- Refer to the Alinity c Hemoglobin A1c Calibrators package insert and the Alinity c Hemoglobin A1c Controls package insert or commercially available control material package insert for preparation and usage.
- For general operating procedures, refer to the Alinity ci-series Operations Manual, Section 5.
- For optimal performance, it is important to perform routine maintenance as described in the Alinity ci-series Operations Manual, Section 9. Perform maintenance more frequently when required by laboratory procedures.

### Sample Dilution Procedures

Samples cannot be diluted for the Alinity c Hemoglobin A1c assay.

## Calibration

For instructions on performing a calibration, refer to the Alinity ci-series Operations Manual, Section 5.

Calibration is stable for approximately 50 days (1200 hours), but is required with each change in reagent lot. Verify calibration with 2 levels of controls according to the established quality control requirements for your laboratory. If control results fall outside acceptable ranges, recalibration may be necessary.

This assay may require recalibration after maintenance to critical parts or subsystems or after service procedures have been performed.

Both the Whole Blood and Hemolysate applications use the Alinity c Hemoglobin A1c Calibrators which are supplied separately.

Alinity c Hemoglobin A1c Calibrators are traceable to NGSP and IFCC reference methods.

## Quality Control Procedures

As appropriate, refer to your laboratory standard operating procedure(s) and/or quality assurance plan for additional quality control requirements and potential corrective actions.

- Two levels of controls are to be run every 24 hours. The Alinity c Hemoglobin A1c assay uses:
  - Commercially available whole blood controls for the Whole Blood application.
  - Alinity c Hemoglobin A1c Controls or other commercially available controls for the Hemolysate application.  
Refer to the Alinity c Hemoglobin A1c Controls Value Sheet for NGSP and IFCC ranges of the Alinity c Hemoglobin A1c Controls.
- If more frequent control monitoring is required, follow the established quality control procedures for your laboratory.
- If quality control results do not meet the acceptance criteria defined by your laboratory, sample results may be suspect. Follow the established quality control procedures for your laboratory. Recalibration may be necessary. For troubleshooting information, refer to the Alinity ci-series Operations Manual, Section 10.
- Review quality control results and acceptance criteria following a change of reagent or calibrator lot.

Commercial controls should be used according to the guidelines and recommendations of the control manufacturer. Concentration ranges provided in the control package insert should be used only for guidance.

For any control material in use, the laboratory should ensure that the matrix of the control material is suitable for use in the assay per the assay package insert.

## Quality Control Guidance

Refer to "Basic QC Practices" by James O Westgard, Ph.D. for guidance on laboratory quality control practices.<sup>16</sup>

## Verification of Assay Claims

For protocols to verify package insert claims, refer to Verification of Assay Claims in the Alinity ci-series Operations Manual.

## RESULTS

- The individual concentrations of the Whole Blood and Hemolysate applications are measured by the analyzer.  
**NOTE: HbA1c or THb concentrations must not be used individually for clinical purposes.**
- Refer to the LIMITATIONS OF THE PROCEDURE section of this package insert.

## Conventional Units (NGSP)

The percent HbA1c (%HbA1c) is automatically calculated by the system per the calculation provided in the Principles of the Procedure section.

## SI Units (IFCC)

The hemoglobin A1c fraction (mmol/mol HbA1c) is automatically calculated by the system per the calculation provided in the Principles of the Procedure section.

## Calculation

The Alinity c Hemoglobin A1c assay utilizes the Linear data reduction method to generate a calibration and results.

## Flags

Some results may contain information in the Flags field. For a description of the flags that may appear in this field, refer to the Alinity ci-series Operations Manual, Section 5.

## Measuring Interval

Measuring interval is defined as the range of values in %HbA1c (mmol/mol) which meets the limits of acceptable performance for linearity, imprecision, and bias.

The measuring interval of the Alinity c Hemoglobin A1c assay is 4.0 to 14.0 %HbA1c (20.22 to 129.51 mmol/mol HbA1c).

## LIMITATIONS OF THE PROCEDURE

- Do not use Alinity c Sample Cups for whole blood samples. Refer to the Assay Procedure section of this package insert for further information.
- This assay must be performed by qualified laboratory personnel, under appropriate laboratory conditions, solely for the intended use of the assay.
- Do not centrifuge samples.
- Do not freeze specimens that have been hemolyzed with the **A1cDIL**.
- Whole blood specimens that require freezing must be stored at -70°C or colder.
- Do not overfill specimen collection tubes. Whole blood samples greater than 78 mm in height from the bottom of tube will result in an instrument error and results will not be generated. Refer to Section 10 of the Alinity ci-series Operations Manual.
- Use specimen collection tubes, or for sample volumes < 600 µL use the 12 x 75 mm polypropylene conical bottom tubes as recommended in the PROCEDURE, Materials Required but not Provided section of this package insert.
- WARNING:** The Alinity c Hemoglobin A1c assay should not be used to diagnose diabetes during pregnancy. Hemoglobin A1c reflects the average blood glucose levels over the preceding 3 months (i.e., the average life span of a red blood cell) and therefore may be falsely low during pregnancy or any other condition associated with recent onset of hyperglycemia and/or decreased red blood cell survival.<sup>17-20</sup>
- Blood transfusions may impact the HbA1c concentration in the patient sample.
- The Alinity c Hemoglobin A1c assay should not be used to diagnose or monitor diabetes in patients with the following conditions:<sup>17-20</sup>
  - hemoglobinopathies except as demonstrated to produce acceptable performance (e.g., sickle cell trait - refer to the SPECIFIC PERFORMANCE CHARACTERISTICS section of this package insert)
  - abnormal red blood cell turnover (e.g., anemias from hemolysis and iron deficiency)
  - malignancies, and severe chronic hepatic and renal disease
- In cases of rapidly evolving Type 1 diabetes, the increase of HbA1c values might be delayed compared to the acute increase in glucose concentrations. In these conditions, diabetes mellitus must be diagnosed based on plasma glucose concentrations and/or the typical clinical symptoms.
- This test should not replace glucose testing for patients with Type 1 diabetes, pediatric patients, or pregnant women.
- The Alinity c Hemoglobin A1c assay is susceptible to interference effects from conjugated bilirubin at > 15.0 mg/dL (180 µmol/L) and unconjugated bilirubin at > 10.0 mg/dL (171 µmol/L).
- The observed bias for samples containing HbC, HbD, HbE, HbS and HbA2 may be impacted by the method used to determine the reference Hemoglobin A1c concentration.

- The Alinity c Hemoglobin A1c assay is susceptible to interference effects from HbF at > 5%. Glycated HbF is not detected by the Alinity c Hemoglobin A1c assay as it does not contain the  $\beta$ -chain that characterizes HbA1c. However, HbF is measured in the total hemoglobin constituent assay and as a consequence, specimens containing high amounts of HbF (> 5%) may result in lower than expected mmol/mol HbA1c values (IFCC) and %HbA1c values (NGSP).
- Refer to the SPECIMEN COLLECTION AND PREPARATION FOR ANALYSIS section of this package insert for specimen limitations.

## EXPECTED VALUES

It is recommended that each laboratory determine its own reference range based upon its particular locale and population characteristics. For monitoring diabetic patients, it is recommended that glycemic goals are individualized following current professional society recommendations.<sup>7</sup> The American Diabetes Association (ADA) recommendations<sup>7</sup> are summarized in the following table.

| HbA1c Value                | Glycemic Goal                 |
|----------------------------|-------------------------------|
| < 8 %HbA1c (64 mmol/mol)   | Less stringent                |
| < 7 %HbA1c (53 mmol/mol)   | General (non-pregnant adults) |
| < 6.5 %HbA1c (48 mmol/mol) | More stringent                |

HbA1c values above 6.5 %HbA1c (48 mmol/mol) are an indication of hyperglycemia during the preceding 2 to 3 months or longer. According to the recommendations of the ADA, HbA1c values above 6.5 %HbA1c (48 mmol/mol) are suitable for the diagnosis of diabetes mellitus. Patients with HbA1c values in the range of 5.7 to 6.4 %HbA1c (39 to 46 mmol/mol) may be at a risk of developing diabetes.<sup>3, 21</sup>

## SPECIFIC PERFORMANCE CHARACTERISTICS

Representative performance data are provided in this section. Results obtained in individual laboratories may vary.

The Alinity c analyzer and the ARCHITECT c System utilize the same reagents and sample/reagent ratios.

Unless otherwise specified, all studies were performed on the Alinity c analyzer.

### Precision

#### Within-Laboratory Precision

A study was performed based on guidance from CLSI EP05-A2.<sup>22</sup> Testing was conducted using 3 lots of the Alinity c Hemoglobin A1c Reagent Kit, 3 lots of the Alinity c Hemoglobin A1c Calibrators, and 1 lot of commercially available controls and 3 instruments. Three controls and 3 human whole blood panels were assayed in a minimum of 2 replicates at 2 separate times per day on 20 different days.

#### Whole Blood (%HbA1c)

| Sample                      | Instrument | n   | Mean<br>(%HbA1c<br>(NGSP)) | Within-Run<br>(Repeatability) |     | Within-<br>Laboratory<br>(Total) <sup>a</sup> |     | Precision<br>with<br>Additional<br>Component of<br>Between-Lot <sup>b</sup> |     |
|-----------------------------|------------|-----|----------------------------|-------------------------------|-----|-----------------------------------------------|-----|-----------------------------------------------------------------------------|-----|
|                             |            |     |                            | SD                            | %CV | SD                                            | %CV | SD                                                                          | %CV |
| Control 1                   | 1          | 240 | 5.1                        | 0.01                          | 0.1 | 0.05                                          | 1.0 | 0.10                                                                        | 1.9 |
|                             | 2          | 240 | 5.1                        | 0.04                          | 0.7 | 0.06                                          | 1.2 | 0.10                                                                        | 1.9 |
|                             | 3          | 240 | 5.1                        | 0.02                          | 0.5 | 0.06                                          | 1.1 | 0.11                                                                        | 2.1 |
| Control 2                   | 1          | 240 | 9.2                        | 0.03                          | 0.3 | 0.07                                          | 0.7 | 0.12                                                                        | 1.3 |
|                             | 2          | 240 | 9.3                        | 0.03                          | 0.4 | 0.07                                          | 0.7 | 0.13                                                                        | 1.5 |
|                             | 3          | 240 | 9.3                        | 0.03                          | 0.4 | 0.06                                          | 0.7 | 0.11                                                                        | 1.2 |
| Control 3                   | 1          | 240 | 13.5                       | 0.03                          | 0.2 | 0.06                                          | 0.5 | 0.15                                                                        | 1.1 |
|                             | 2          | 240 | 13.6                       | 0.04                          | 0.3 | 0.08                                          | 0.6 | 0.16                                                                        | 1.2 |
|                             | 3          | 240 | 13.6                       | 0.04                          | 0.3 | 0.06                                          | 0.5 | 0.12                                                                        | 0.9 |
| Panel<br>near 4.0<br>%HbA1c | 1          | 240 | 4.4                        | 0.03                          | 0.7 | 0.06                                          | 1.3 | 0.07                                                                        | 1.7 |
|                             | 2          | 240 | 4.4                        | 0.03                          | 0.7 | 0.05                                          | 1.2 | 0.07                                                                        | 1.5 |
|                             | 3          | 240 | 4.4                        | 0.03                          | 0.7 | 0.06                                          | 1.3 | 0.08                                                                        | 1.8 |
| Panel<br>near 6.5<br>%HbA1c | 1          | 240 | 6.2                        | 0.02                          | 0.3 | 0.03                                          | 0.5 | 0.07                                                                        | 1.2 |
|                             | 2          | 240 | 6.2                        | 0.02                          | 0.4 | 0.04                                          | 0.6 | 0.07                                                                        | 1.2 |
|                             | 3          | 240 | 6.2                        | 0.03                          | 0.4 | 0.04                                          | 0.6 | 0.08                                                                        | 1.3 |

| Sample                      | Instrument | n   | Mean<br>(%HbA1c<br>(NGSP)) | Within-Run<br>(Repeatability) |     | Within-<br>Laboratory<br>(Total) <sup>a</sup> |     | Precision<br>with<br>Additional<br>Component of<br>Between-Lot <sup>b</sup> |     |
|-----------------------------|------------|-----|----------------------------|-------------------------------|-----|-----------------------------------------------|-----|-----------------------------------------------------------------------------|-----|
|                             |            |     |                            | SD                            | %CV | SD                                            | %CV | SD                                                                          | %CV |
| Panel<br>near 8.0<br>%HbA1c | 1          | 240 | 7.9                        | 0.03                          | 0.3 | 0.05                                          | 0.6 | 0.09                                                                        | 1.2 |
|                             | 2          | 240 | 7.9                        | 0.04                          | 0.5 | 0.06                                          | 0.7 | 0.09                                                                        | 1.2 |
|                             | 3          | 240 | 7.9                        | 0.04                          | 0.5 | 0.05                                          | 0.7 | 0.11                                                                        | 1.3 |

<sup>a</sup> Within-Laboratory variability contains within-run, between-run, and between-day variance components.

<sup>b</sup> Contains within-run, within-day, between-day, and between-lot variance components.

#### Whole Blood (mmol/mol)

| Sample                         | Instrument | n   | Mean<br>(mmol/mol<br>(IFCC)) | Within-Run<br>(Repeatability) |     | Within-<br>Laboratory<br>(Total) <sup>a</sup> |     | Precision<br>with<br>Additional<br>Component of<br>Between-Lot <sup>b</sup> |     |
|--------------------------------|------------|-----|------------------------------|-------------------------------|-----|-----------------------------------------------|-----|-----------------------------------------------------------------------------|-----|
|                                |            |     |                              | SD                            | %CV | SD                                            | %CV | SD                                                                          | %CV |
| Control 1                      | 1          | 240 | 32.15                        | 0.127                         | 0.4 | 0.534                                         | 1.7 | 1.009                                                                       | 3.1 |
|                                | 2          | 240 | 32.40                        | 0.195                         | 0.6 | 0.546                                         | 1.7 | 0.959                                                                       | 3.0 |
|                                | 3          | 240 | 32.38                        | 0.169                         | 0.5 | 0.537                                         | 1.7 | 1.089                                                                       | 3.4 |
| Control 2                      | 1          | 240 | 77.17                        | 0.153                         | 0.2 | 0.627                                         | 0.8 | 1.332                                                                       | 1.7 |
|                                | 2          | 240 | 77.64                        | 0.241                         | 0.3 | 0.639                                         | 0.8 | 1.344                                                                       | 1.7 |
|                                | 3          | 240 | 77.80                        | 0.182                         | 0.2 | 0.582                                         | 0.7 | 1.204                                                                       | 1.5 |
| Control 3                      | 1          | 240 | 124.49                       | 0.197                         | 0.2 | 0.640                                         | 0.5 | 1.641                                                                       | 1.3 |
|                                | 2          | 240 | 125.22                       | 0.373                         | 0.3 | 0.734                                         | 0.6 | 1.693                                                                       | 1.4 |
|                                | 3          | 240 | 125.52                       | 0.297                         | 0.2 | 0.605                                         | 0.5 | 1.220                                                                       | 1.0 |
| Panel<br>near 20.0<br>mmol/mol | 1          | 240 | 24.81                        | 0.195                         | 0.8 | 0.457                                         | 1.8 | 0.745                                                                       | 3.0 |
|                                | 2          | 240 | 24.90                        | 0.146                         | 0.6 | 0.428                                         | 1.7 | 0.682                                                                       | 2.7 |
|                                | 3          | 240 | 24.93                        | 0.107                         | 0.4 | 0.364                                         | 1.5 | 0.781                                                                       | 3.1 |
| Panel<br>near 47.5<br>mmol/mol | 1          | 240 | 43.82                        | 0.152                         | 0.3 | 0.400                                         | 0.9 | 0.811                                                                       | 1.9 |
|                                | 2          | 240 | 43.93                        | 0.214                         | 0.5 | 0.422                                         | 1.0 | 0.804                                                                       | 1.8 |
|                                | 3          | 240 | 43.97                        | 0.192                         | 0.4 | 0.382                                         | 0.9 | 0.891                                                                       | 2.0 |
| Panel<br>near 63.9<br>mmol/mol | 1          | 240 | 62.94                        | 0.145                         | 0.2 | 0.396                                         | 0.6 | 0.886                                                                       | 1.4 |
|                                | 2          | 240 | 63.07                        | 0.238                         | 0.4 | 0.456                                         | 0.7 | 0.905                                                                       | 1.4 |
|                                | 3          | 240 | 63.13                        | 0.175                         | 0.3 | 0.385                                         | 0.6 | 0.959                                                                       | 1.5 |

<sup>a</sup> Within-Laboratory variability contains within-run, between-run, and between-day variance components.

<sup>b</sup> Contains within-run, within-day, between-day, and between-lot variance components.

### Lower Limits of Measurement

A study was performed based on guidance from CLSI EP17-A2.<sup>23</sup> Testing was conducted using 3 lots of the Alinity c Hemoglobin A1c Reagent Kit on each of 2 instruments over a minimum of 3 days. The maximum observed Limit of Blank (LoB) and Limit of Detection (LoD) are summarized below.

#### Whole Blood

|                  | Glycated Hemoglobin (HbA1c)<br>μmol/L | Total Hemoglobin (THb)<br>μmol/L |
|------------------|---------------------------------------|----------------------------------|
| LoB <sup>a</sup> | 43.8095                               | 0.0000                           |
| LoD <sup>b</sup> | 48.1568                               | 6.7745                           |

<sup>a</sup> The LoB represents the 95th percentile from  $n \geq 60$  replicates of zero-analyte samples.

<sup>b</sup> The LoD represents the lowest concentration at which the analyte can be detected with 95% probability based on  $n \geq 60$  replicates of low-analyte level samples.

## Hemolysate

|                  | Glycated Hemoglobin (HbA1c)<br>μmol/L | Total Hemoglobin (THb)<br>μmol/L |
|------------------|---------------------------------------|----------------------------------|
| LoB <sup>a</sup> | 39.0741                               | 12.2147                          |
| LoD <sup>b</sup> | 41.0713                               | 16.3614                          |

<sup>a</sup> The LoB represents the 95th percentile from  $n \geq 60$  replicates of zero-analyte samples.

<sup>b</sup> The LoD represents the lowest concentration at which the analyte can be detected with 95% probability based on  $n \geq 60$  replicates of low-analyte level samples.

## Linearity

A study was performed based on guidance from CLSI EP06-A.<sup>24</sup>

This assay is linear across the measuring interval of 4.0 to 14.0 %HbA1c (20.22 to 129.51 mmol/mol HbA1c).

## Analytical Specificity

### Hemoglobin Derivatives

This study was performed on the ARCHITECT c System.

A specificity study was performed based on guidance from CLSI EP07-A2.<sup>25</sup> Specificity was assessed by comparing test samples containing the potential interferents listed below to reference samples. No interference was observed for:

- Acetylated Hemoglobin with  $\geq 50$  mg/dL of acetylsalicylic acid (ASA) (aspirin)
- Carbamylated Hemoglobin with  $\geq 10$  mmol/L of Cyanate
- Labile Hemoglobin with  $\geq 1000$  mg/dL of Glucose

## IFCC

The Alinity c Hemoglobin A1c assay had a difference within  $\pm 7\%$  for samples with concentrations  $\geq 38.78$  mmol/mol HbA1c.

## NGSP

The Alinity c Hemoglobin A1c assay had a difference within  $\pm 5\%$  for samples with concentrations  $\geq 5.7$  %HbA1c.

## Hemoglobin Variants

A study was performed based on guidance from CLSI documents EP07-A2.<sup>25</sup>

Specificity was assessed by comparing samples containing abnormal hemoglobins to a comparator assay. Heterozygous hemoglobin variants (HbS, HbC, HbD, HbE, HbA2) do not interfere with the Alinity c Hemoglobin A1c assay.

## IFCC

The Alinity c Hemoglobin A1c assay had a difference within  $\pm 7\%$ .

## NGSP

The Alinity c Hemoglobin A1c assay had a difference within  $\pm 5\%$ .

| Heterozygous Hemoglobin Variant | Relative % Difference from Reference Concentration                                    |                    |                                               |                    |
|---------------------------------|---------------------------------------------------------------------------------------|--------------------|-----------------------------------------------|--------------------|
|                                 | ~ 6.0 %HbA1c (5.5 to 6.5 %HbA1c)                                                      |                    | ~9.0 %HbA1c (8.0 to 10.0 %HbA1c) <sup>a</sup> |                    |
|                                 | Relative % Difference                                                                 | Range <sup>b</sup> | Relative % Difference                         | Range <sup>b</sup> |
| HbC                             | -4.2                                                                                  | -8.9, -1.6         | -2.5                                          | -6.2, 0.0          |
| HbD                             | 1.9                                                                                   | -1.7, 4.8          | 2.3                                           | 1.1, 3.5           |
| HbE                             | 1.4                                                                                   | -3.2, 6.6          | -0.9                                          | -4.1, 2.5          |
| HbS                             | -0.3                                                                                  | -5.5, 1.8          | -1.9                                          | -4.9, 1.1          |
| HbA2                            | -0.8                                                                                  | -4.9, 3.6          | -0.6                                          | -3.7, 5.4          |
| HbF                             | Difference exceeds -5% when the amount of HbF in the sample exceeds ~5%. <sup>c</sup> |                    |                                               |                    |
|                                 | See the Heterozygous Hemoglobin Variant HbF interference in the following table.      |                    |                                               |                    |

<sup>a</sup> The HbA2 results at ~9.0 %HbA1c consisted of samples between 7.1 to 14.0 %HbA1c.

<sup>b</sup> The range is defined as the minimum and maximum relative % difference at each concentration level (~6.0 and ~9.0 %HbA1c).

<sup>c</sup> Refer to the LIMITATIONS OF THE PROCEDURE section of this package insert for further information.

**NOTE:** The presence of multiple variants in a sample may impact the % difference.

Heterozygous Hemoglobin Variant HbF interference is summarized in the following table.

A negative % difference with HbF is proportional in magnitude to the % HbF present in the sample. For example, when the amount of HbF in the sample was 20.4%, the % difference was -21.3% on the Alinity c analyzer.

| Heterozygous Hemoglobin Variant | % HbF | Individual % Interference from Comparator |
|---------------------------------|-------|-------------------------------------------|
| HbF                             | 3.2   | -3.1                                      |
|                                 | 4.6   | -4.9                                      |
|                                 | 6.2   | -4.6                                      |
|                                 | 8.6   | -7.4                                      |
|                                 | 11    | -9.8                                      |

Refer to the LIMITATIONS OF THE PROCEDURE section of this package insert for further information regarding the specific concentrations of HbF that were found to interfere with the Alinity c Hemoglobin A1c assay.

## Interference

This study was performed on the ARCHITECT c8000 and c4000 Systems.

### Potentially Interfering Substances

A study was performed based on guidance from CLSI EP07-A2.<sup>25</sup> Interference effects were assessed by comparing test samples containing the potential interferents listed below to reference samples.

## IFCC

The Alinity c Hemoglobin A1c assay had a difference within  $\pm 7\%$  for samples with concentrations  $\geq 38.78$  mmol/mol HbA1c.

## NGSP

The Alinity c Hemoglobin A1c assay had a difference within  $\pm 5\%$  for samples with concentrations  $\geq 5.7$  %HbA1c. The data in %HbA1c (NGSP) showing the highest level of interference observed on the ARCHITECT c8000 and c4000 Systems is presented in the table.

| Potentially Interfering Substance     | Interferent Concentration |               | % Interference <sup>a</sup> |                   |
|---------------------------------------|---------------------------|---------------|-----------------------------|-------------------|
|                                       | Conventional Units        | SI Units      | 6.0 to 7.0 %HbA1c           | $\geq 8.0$ %HbA1c |
| Ascorbic Acid                         | 3.0 mg/dL                 | 0.15 mmol/L   | 0.0                         | 0.0               |
| Bilirubin (Conjugated) <sup>b</sup>   | 15.0 mg/dL                | 180 μmol/L    | -3.2                        | -3.3              |
| Bilirubin (Unconjugated) <sup>c</sup> | 10.0 mg/dL                | 171 μmol/L    | -3.0                        | -2.7              |
| Glucose                               | 1000 mg/dL                | 55.5 mmol/L   | 0.0                         | 0.0               |
| Rheumatoid Factor                     | 200 IU/mL                 | 200 IU/mL     | 0.0                         | 0.0               |
| Triglycerides                         | 3000 mg/dL                | 33.9 mmol/L   | -1.6                        | -4.5              |
| Total Protein                         | 22 g/dL <sup>d</sup>      | 220 g/L       | 0.0                         | -1.1              |
| Urea                                  | 667 mg/dL                 | 111.06 mmol/L | 0.0                         | 0.0               |
| Vitamin E                             | 8.6 mg/dL                 | 200 μmol/L    | 1.6                         | 0.0               |

<sup>a</sup>

$$\% \text{ Interference} = \frac{\text{Test Result} - \text{Control Result}}{\text{Control Result}} \times 100$$

<sup>b</sup> Samples containing conjugated bilirubin at  $> 15.0$  mg/dL (180 μmol/L) demonstrated interference. Refer to the LIMITATIONS OF THE PROCEDURE section of this package insert for further information.

<sup>c</sup> Samples containing unconjugated bilirubin at  $> 10.0$  mg/dL (171 μmol/L) demonstrated interference. Refer to the LIMITATIONS OF THE PROCEDURE section of this package insert for further information.

<sup>d</sup> The total protein concentration of 22 g/dL includes serum protein as well as hemoglobin.

| Potentially Interfering Substance | Interferent Concentration |                        | % Interference <sup>a</sup> |              |
|-----------------------------------|---------------------------|------------------------|-----------------------------|--------------|
|                                   | Conventional Units        | SI Units               | 6.0 to 7.0 %HbA1c           | ≥ 8.0 %HbA1c |
| Acarbose                          | 50 mg/dL                  | 0.77 mmol/L            | 0.0                         | 0.0          |
| Acetaminophen                     | 200 mg/L                  | 1324.5 µmol/L          | -0.59                       | -0.80        |
| N-acetyl-4-benzoquinone imine     | 20 mg/L                   | 134.2 µmol/L           | -1.5                        | -0.8         |
| N-acetyl-L-cysteine               | 1600 mg/L                 | 9816 µmol/L            | -3.4                        | -2.7         |
| Acetylsalicylate                  | 50.8 mg/dL                | 2.82 mmol/L            | 0.0                         | 0.0          |
| Atorvastatin                      | 0.06 mg/dL                | 600 µg Eq/L            | 1.6                         | 0.0          |
| Captopril                         | 0.5 mg/dL                 | 23 µmol/L              | -1.5                        | -1.1         |
| Chlorpropamide                    | 74.7 mg/dL                | 2.7 mmol/L             | 0.0                         | -1.1         |
| Cyanate                           | 50 mg/dL                  | 6.16 mmol/L            | 0.0                         | 1.1          |
| Dipyrrone                         | 100 mg/L                  | 300.3 µmol/L           | 0.2                         | 0.1          |
| Furosemide                        | 6.0 mg/dL                 | 181 µmol/L             | 0.0                         | 1.1          |
| Gemfibrozil                       | 7.5 mg/dL                 | 300 µmol/L             | 0.0                         | 0.0          |
| Ibuprofen                         | 50 mg/dL                  | 2425 µmol/L            | 0.0                         | 1.1          |
| Insulin                           | 450 micro units per mL    | 450 micro units per mL | 0.8                         | 0.0          |
| Losartan                          | 5 mg/dL                   | 0.11 mmol/L            | 0.0                         | 0.0          |
| Metformin                         | 5.1 mg/dL                 | 310 µmol/L             | 0.0                         | 0.0          |
| Nicotinic Acid                    | 61 mg/dL                  | 4.95 mmol/L            | -1.5                        | -0.5         |
| Propranolol                       | 0.2 mg/dL                 | 7.71 µmol/L            | 0.0                         | -0.5         |
| Repaglinide                       | 0.006 mg/dL               | 132.57 nmol/L          | 0.8                         | 0.0          |
| 4-acetamido antipyrine            | 40 mg/L                   | 163.3 µmol/L           | -0.4                        | -0.4         |
| 4-aminoantipyrine                 | 40 mg/L                   | 197.0 µmol/L           | -0.4                        | -0.2         |
| 4-formylamino antipyrine          | 40 mg/L                   | 173.2 µmol/L           | -0.1                        | -0.6         |
| 4-methylamino antipyrine          | 40 mg/L                   | 184.3 µmol/L           | -0.3                        | -0.4         |

a

$$\% \text{ Interference} = \frac{\text{Test Result} - \text{Control Result}}{\text{Control Result}} \times 100$$

Refer to the LIMITATIONS OF THE PROCEDURE section of this package insert for further information.

### Method Comparison

A study was performed based on guidance from CLSI EP09-A3 using the Passing-Bablok regression method.<sup>26</sup>

|                                                      |             | Units           | n   | Correlation Coefficient | Intercept | Slope | Concentration Range |
|------------------------------------------------------|-------------|-----------------|-----|-------------------------|-----------|-------|---------------------|
| Alinity c Hemoglobin A1c vs NGSP Secondary           | Whole Blood | %HbA1c (NGSP)   | 134 | 0.99                    | -0.10     | 1.00  | 4.4 - 14.0          |
|                                                      | Hemolysate  | %HbA1c (NGSP)   | 134 | 0.99                    | -0.10     | 1.00  | 4.4 - 14.0          |
| Alinity c Hemoglobin A1c vs ARCHITECT Hemoglobin A1c | Whole Blood | %HbA1c (NGSP)   | 130 | 1.00                    | 0.10      | 1.00  | 4.5 - 13.1          |
|                                                      |             | mmol/mol (IFCC) | 130 | 1.00                    | 0.10      | 1.01  | 25.72 - 119.99      |
|                                                      | Hemolysate  | %HbA1c (NGSP)   | 133 | 1.00                    | 0.00      | 1.00  | 4.7 - 13.1          |
|                                                      |             | mmol/mol (IFCC) | 133 | 1.00                    | -0.39     | 1.02  | 28.22 - 119.44      |

### Predicted Bias

#### IFCC

The Alinity c Hemoglobin A1c assay is designed to have a predicted bias of ≤ 5% at 42.06, 47.53, and 53.00 mmol/mol HbA1c using Deming regression.

#### NGSP


The Alinity c Hemoglobin A1c assay is designed to have a predicted bias of ≤ 3% at 6.0, 6.5, and 7.0 %HbA1c using Deming regression. The predicted bias in %HbA1c (NGSP) ranged from -2.0% to -1.5%.

### Allowable Total Difference (ATD) Zone

The Alinity c Hemoglobin A1c assay is designed to have > 95% of observations in the ATD zone and the low limit of the two-sided 95% Confidence Interval (CI) > 89.5%. The percentage of observations in the ATD zone was 99.3% (133/134) and the lower limit of the two-sided 95% CI was 95.9%. The ATD zone plots are presented below.

#### ATD Zone- Alinity c analyzer

##### Hemoglobin A1c in Whole Blood, mmol/mol (IFCC)


x = NGSP Secondary Reference Laboratory (mmol/mol)

y = Alinity c Hemoglobin A1c (mmol/mol)

#### ATD Zone- Alinity c analyzer

##### Hemoglobin A1c in Whole Blood, %HbA1c (NGSP)

x = NGSP Secondary Reference Laboratory (%HbA1c)

y = Alinity c Hemoglobin A1c (%HbA1c)

## BIBLIOGRAPHY

- Bry L, Chen PC, Sacks DB. Effects of hemoglobin variants and chemically modified derivatives on assays for glycohemoglobin. *Clin Chem* 2001;47(2):153-163.
- Sacks DB. The diagnosis of diabetes is changing: How implementation of hemoglobin A1c will impact clinical laboratories. *Clin Chem* 2009;55(9):1612-1614.
- American Diabetes Association Workgroup Report: International Expert Committee report on the role of the A1c assay in the diagnosis of diabetes. *Diabetes Care* 2009;32(7):1327-1334.
- Fonseca V, Inzucchi SE, Ferrannini E. Redefining the diagnosis of diabetes using glycated hemoglobin. *Diabetes Care* 2009;32(7):1344-1345.
- American Diabetes Association Executive Summary: Standards of medical care in diabetes - 2010. *Diabetes Care* 2010;33(Suppl 1):S4-S10.
- Sacks DB, Bruns DE, Goldstein DE, et al. Guidelines and recommendations for laboratory analysis in the diagnosis and management of diabetes mellitus. *Clin Chem* 2002;48(3):436-472.
- American Diabetes Association. Position Statement: Standards of Medical Care in Diabetes - 2012. In: *Diabetes Care* 2012;35(Suppl 1):S11-S63.
- Rochman H. Hemoglobin A1c and diabetes mellitus. *Ann Clin Lab Sci* 1980;10(2):111-115.
- The Diabetes Control and Complications Trial Research Group. The effect of intensive treatment of diabetes on the development and progression of long-term complications in insulin-dependent diabetes mellitus. *N Engl J Med* 1993;329(14):977-986.
- Larsen ML, Hørder M, Mogensen EF. Effect of long-term monitoring of glycosylated hemoglobin levels in insulin-dependent diabetes mellitus. *N Engl J Med* 1990;323(15):1021-1025.
- Geistanger A, Arends S, Berding C, et al. Statistical methods for monitoring the relationship between the IFCC reference measurement procedure for hemoglobin A1c and the designated comparison methods in the United States, Japan, and Sweden. *Clin Chem* 2008;54(8):1379-1385.
- US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
- US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
- World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
- Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.
- Westgard JO. *Basic QC Practices*. 3rd ed. Madison, WI: Westgard Quality Corporation; 2010.
- Goldstein DE, Wiedmeyer HM, England JD, et al. Recent advances in glycosylated hemoglobin measurements. *Crit Rev Clin Lab Sci* 1984;21(3):187-228.
- Peacock I. Glycosylated haemoglobin: measurement and clinical use. *J Clin Pathol* 1984;37:841-851.
- Horton BF, Huisman TH. Studies on the heterogeneity of haemoglobin, VII. Minor haemoglobin components in haematological diseases. *Br J Haematol* 1965;11(3):296-304.
- Lind T, Cheyne GA. Effect of normal pregnancy upon the glycosylated haemoglobins. *Br J Obstet Gynaecol* 1979;86(3):210-213.
- American Diabetes Association. Position Statement: Diagnosis and classification of diabetes mellitus. *Diabetes Care* 2010;33(Suppl 1):S62-S69.
- Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition*. CLSI Document EP05-A2. Wayne, PA: CLSI; 2004.
- Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition*. CLSI Document EP17-A2. Wayne, PA: CLSI; 2012.
- Clinical and Laboratory Standards Institute (CLSI). *Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline*. CLSI Document EP06-A. Wayne, PA: CLSI; 2003.
- Clinical and Laboratory Standards Institute (CLSI). *Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition*. CLSI Document EP07-A2. Wayne, PA: CLSI; 2005.

- Clinical and Laboratory Standards Institute (CLSI). *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition*. CLSI Document EP09-A3. Wayne, PA: CLSI; 2013.

Note for number formatting:

- A space is used as thousands separator (example: 10 000 specimens).
- A period is used to separate the integer part from the fractional part of a number written in decimal form (example: 3.12%).

## Key to Symbols

| ISO 15223 Symbols                                                                 |                                                                                         |
|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------|
|  | Consult instructions for use                                                            |
|  | Manufacturer                                                                            |
|  | Sufficient for                                                                          |
|  | Temperature limitation                                                                  |
|  | Use by/Expiration date                                                                  |
| <b>IVD</b>                                                                        | <i>In Vitro</i> Diagnostic Medical Device                                               |
| <b>LOT</b>                                                                        | Lot Number                                                                              |
| <b>REF</b>                                                                        | List Number                                                                             |
| <b>SN</b>                                                                         | Serial number                                                                           |
| Other Symbols                                                                     |                                                                                         |
| <b>A1cDIL</b>                                                                     | Hemoglobin A1c Diluent                                                                  |
| <b>CONTAINS: AZIDE</b>                                                            | Contains Sodium Azide. Contact with acids liberates very toxic gas.                     |
| <b>DISTRIBUTED IN THE USA BY</b>                                                  | Distributed in the USA by                                                               |
| <b>INFORMATION FOR USA ONLY</b>                                                   | Information needed for United States of America only                                    |
| <b>PRODUCED FOR ABBOTT BY</b>                                                     | Produced for Abbott by                                                                  |
| <b>PRODUCT OF CANADA</b>                                                          | Product of Canada                                                                       |
| <b>R1</b>                                                                         | Reagent 1                                                                               |
| <b>R2</b>                                                                         | Reagent 2                                                                               |
| <b>Rx ONLY</b>                                                                    | For use by or on the order of a physician only (applicable to USA classification only). |

Alinity, ARCHITECT and related brand marks are trademarks of Abbott. Other trademarks are the property of their respective owners.

Abbott GmbH  
Max-Planck-Ring 2  
65205 Wiesbaden  
Germany  
+49-6122-580

**PRODUCED FOR ABBOTT BY**

Sekisui Diagnostics P.E.I. Inc.  
70 Watts Avenue  
Charlottetown  
Prince Edward Island  
C1E 2B9 Canada

**DISTRIBUTED IN THE USA BY**

Abbott Laboratories  
Abbott Park, IL 60064 USA

**Customer Service: Contact your local representative or find country-specific contact information on [www.corelaboratory.abbott](http://www.corelaboratory.abbott)**

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## BIBLIOGRAFIE

- Jacobs DS, DeMott WR, Oxley DK, et al. *Laboratory Test Handbook*. Hudson, OH: Lexi-Comp; 2001:131–133.
- Bendz H, Sjödin I, Toss G, et al. Hyperparathyroidism and long-term lithium therapy—a cross-sectional study and the effect of lithium withdrawal. *J Intern Med* 1996;240:357–365.
- US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
- US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
- World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
- Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.
- Burtis CA, Ashwood ER, editors. *Tietz Fundamentals of Clinical Chemistry*, 5th ed. Philadelphia, PA: WB Saunders; 2001:800.
- Guder WG, Narayanan S, Wisser H, et al. List of Analytes—Preamplified Variables. Annex In: *Samples: From the Patient to the Laboratory*. Darmstadt, Germany: GIT Verlag; 1996:Annex 10–1, 36–7.
- US Pharmacopeial Convention, Inc. General notices. In: *US Pharmacopeia National Formulary*, 1995 ed (USP 23/NF 18). Rockville, MD: The US Pharmacopeial Convention, Inc; 1994:11.
- Guder WG, da Fonseca-Wollheim F, Heil W, et al. *The Quality of Diagnostic Samples*. Darmstadt, Germany: GIT Verlag; 2001:14–15.
- Westgard JO. *Basic QC Practices*. 3rd ed. Madison, WI: Westgard Quality Corporation; 2010.
- Burtis CA, Ashwood ER, editors. *Tietz Textbook of Clinical Chemistry*, 2nd ed. Philadelphia, PA: WB Saunders; 1994:2180.
- National Committee for Clinical Laboratory Standards (NCCLS). *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition*. NCCLS Document EP5-A2. Wayne, PA: NCCLS; 2004.
- Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition*. CLSI Document EP17-A2. Wayne, PA: CLSI; 2012.
- Clinical and Laboratory Standards Institute (CLSI). *Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline*. CLSI Document EP06-A. Wayne, PA: CLSI; 2003.
- Young DS. *Effects of Drugs on Clinical Laboratory Tests*, 4th ed. Washington, DC: AACC Press; 1995:3-115–3-124.
- Clinical and Laboratory Standards Institute (CLSI). *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition*. CLSI Document EP09-A3. Wayne, PA: CLSI; 2013.

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## REFERINȚE

1. Friedman RB, Young DS. *Effects of Disease on Clinical Laboratory Tests*, 2nd ed. Washington, DC: AACC Press, 1989.
2. Thomas L. Immunochemical techniques. In: Thomas L, editor. *Clinical Laboratory Diagnostics: Use and Assessment of Clinical Laboratory Results*, 1st ed. Frankfurt/Main, Germany: TH-Books Verlagsgesellschaft; 1998:1427–41.
3. Simó JM, Joven J, Clivillé X, et al. Automated latex agglutination immunoassay of serum ferritin with a centrifugal analyzer. *Clin Chem* 40;1994:625–9.
4. US Department of Labor, Occupational Safety and Health Administration. 29 CFR Part 1910.1030. Bloodborne pathogens.
5. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*, 5th ed. Washington, DC: US Government Printing Office, December 2009.
6. World Health Organization. *Laboratory Biosafety Manual*, 3rd ed. Geneva: World Health Organization, 2004.
7. Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.
8. Westgard JO, Barry PL. *Cost-effective Quality Control: Managing the Quality and Productivity of Analytical Processes*. Washington, DC: AACC Press, 1986.
9. Young DS. *Effects of Drugs on Clinical Laboratory Tests*, 5th ed. Washington, DC: AACC Press, 2000.

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| Configure assay parameters — General                 |                                        |                                       |                               |
|------------------------------------------------------|----------------------------------------|---------------------------------------|-------------------------------|
| <input checked="" type="radio"/> General             | <input type="radio"/> Calibration      | <input type="radio"/> SmartWash       | <input type="radio"/> Results |
| Assay: <b>FER</b>                                    | Type: <b>Photometric</b>               | Version: †                            |                               |
| Number: <b>2906</b>                                  |                                        |                                       |                               |
| Run controls for onboard reagents by:†† <b>Lot</b>   |                                        |                                       |                               |
| <input checked="" type="radio"/> Reaction definition | <input type="radio"/> Reagent / Sample | <input type="radio"/> Validity checks |                               |
| Reaction mode: <b>End up</b>                         |                                        |                                       |                               |
| Primary                                              |                                        | Secondary                             | Read times                    |
| Wavelength: <b>572</b> / <b>None</b>                 |                                        | Main: <b>32 – 33</b>                  |                               |
| Last required read: <b>33</b>                        |                                        |                                       |                               |
| Absorbance range: <b>— —</b>                         |                                        | Color correction: <b>— —</b>          |                               |
| Sample blank type: <b>Self</b>                       |                                        | Blank: <b>19 – 20</b>                 |                               |

| O Reaction definition                |               | ● Reagent / Sample           |            | O Validity checks |                 |
|--------------------------------------|---------------|------------------------------|------------|-------------------|-----------------|
| Reagent: <b>FER0B</b>                |               | Reagent volume: <b>160</b>   |            | R1                | R2              |
| Diluent: <b>Saline</b>               |               | Water volume: <b>—</b>       |            |                   | <b>50</b>       |
| Diluent dispense mode: <b>Type 0</b> |               | Dispense mode: <b>Type 0</b> |            | <b>Type 0</b>     | <b>Type 2</b>   |
| Dilution name                        | Sample        | Diluted sample               | Diluent    | Water             | Dilution factor |
| <b>STANDARD</b>                      | <b>: 25.0</b> | <b>—</b>                     | <b>—</b>   | <b>=</b>          | <b>1:1.00</b>   |
| <b>Dil 1</b>                         | <b>: 20.0</b> | <b>15.0</b>                  | <b>100</b> | <b>=</b>          | <b>1:9.57</b>   |
|                                      |               |                              |            |                   |                 |
|                                      |               |                              |            |                   |                 |

| O Reaction definition                  |  | O Reagent / Sample                         |  | ● Validity checks |   |
|----------------------------------------|--|--------------------------------------------|--|-------------------|---|
| Reaction check: <b>Rate Ratio</b>      |  |                                            |  | A                 | B |
|                                        |  | Read time: <b>24 – 26</b>                  |  | <b>18 – 20</b>    |   |
|                                        |  | Calculation limits: <b>0.5000 – 9.0000</b> |  |                   |   |
|                                        |  | Minimum: <b>0.0750</b>                     |  |                   |   |
| Maximum absorbance variation: <b>—</b> |  |                                            |  |                   |   |

| Configure assay parameters — Calibration     |                                              |                                 |                                       |
|----------------------------------------------|----------------------------------------------|---------------------------------|---------------------------------------|
| <input type="radio"/> General                | <input checked="" type="radio"/> Calibration | <input type="radio"/> SmartWash | <input type="radio"/> Results         |
| Assay: <b>FER</b>                            | Calibration method: <b>Spline</b>            |                                 |                                       |
| <input checked="" type="radio"/> Calibrators | <input type="radio"/> Volumes                | <input type="radio"/> Intervals | <input type="radio"/> Validity checks |
| Calibrator set: <b>FER</b>                   |                                              |                                 |                                       |
| Blank: <b>Water</b>                          |                                              | Concentration: <b>0.0†</b>      |                                       |
| Cal 1: <b>FER1</b>                           |                                              | ††                              |                                       |
| Cal 2: <b>FER2</b>                           |                                              | ††                              |                                       |
| Cal 3: <b>FER3</b>                           |                                              | ††                              |                                       |
| Cal 4: <b>FER4</b>                           |                                              | ††                              |                                       |

| O Calibrators          |  | ● Volumes        |  | O Intervals |                | O Validity checks |          |
|------------------------|--|------------------|--|-------------|----------------|-------------------|----------|
| Calibrator: <b>FER</b> |  | Calibrator level |  | Sample      | Diluted sample | Diluent           | Water    |
| Blank: <b>Water</b>    |  | <b>25.0</b>      |  | <b>—</b>    | <b>—</b>       | <b>—</b>          | <b>—</b> |
| Cal 1: <b>FER1</b>     |  | <b>25.0</b>      |  | <b>—</b>    | <b>—</b>       | <b>—</b>          | <b>—</b> |
| Cal 2: <b>FER2</b>     |  | <b>25.0</b>      |  | <b>—</b>    | <b>—</b>       | <b>—</b>          | <b>—</b> |
| Cal 3: <b>FER3</b>     |  | <b>25.0</b>      |  | <b>—</b>    | <b>—</b>       | <b>—</b>          | <b>—</b> |
| Cal 4: <b>FER4</b>     |  | <b>25.0</b>      |  | <b>—</b>    | <b>—</b>       | <b>—</b>          | <b>—</b> |

| O Calibrators                     |  | O Volumes |  | ● Intervals |  | O Validity checks |  |
|-----------------------------------|--|-----------|--|-------------|--|-------------------|--|
| Calibration intervals:            |  |           |  |             |  |                   |  |
| Full interval: <b>720</b> (hours) |  |           |  |             |  |                   |  |
| Calibration type:                 |  |           |  |             |  |                   |  |
| Adjust type: <b>None</b>          |  |           |  |             |  |                   |  |

| O Calibrators                    |  | O Volumes    |  | O Intervals    |  | ● Validity checks |  |
|----------------------------------|--|--------------|--|----------------|--|-------------------|--|
| Blank absorbance range:          |  | <b>—</b>     |  | <b>—</b>       |  | <b>—</b>          |  |
| Span:                            |  | <b>Blank</b> |  | <b>— Blank</b> |  |                   |  |
| Span absorbance range:           |  | <b>—</b>     |  | <b>—</b>       |  | <b>—</b>          |  |
| Expected cal factor:             |  | <b>0.00</b>  |  |                |  |                   |  |
| Expected cal factor tolerance %: |  | <b>0</b>     |  |                |  |                   |  |

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| Configure assay parameters — Results |                                   |                                 |                                          |                                      |
|--------------------------------------|-----------------------------------|---------------------------------|------------------------------------------|--------------------------------------|
| <input type="radio"/> General        | <input type="radio"/> Calibration | <input type="radio"/> SmartWash | <input checked="" type="radio"/> Results | <input type="radio"/> Interpretation |
| Assay: <b>FER</b>                    |                                   |                                 | Assay number: <b>2906</b>                |                                      |
| Dilution default range:              |                                   |                                 | Result units: <b>ng/mL</b>               |                                      |
| Low-Linearity: <b>10.0</b>           |                                   |                                 |                                          |                                      |
| High-Linearity: <b>500.0</b>         |                                   |                                 |                                          |                                      |
| Gender and age specific ranges:*     |                                   |                                 |                                          |                                      |
| GENDER                               | AGE (UNITS)                       | NORMAL                          | EXTREME                                  |                                      |

| Configure result units            |               |
|-----------------------------------|---------------|
| Assay: <b>FER</b>                 |               |
| Version: <b>†</b>                 |               |
| Result units: <b>ng/mL</b>        |               |
| Decimal places: <b>1</b>          | [Range 0 – 4] |
| Correlation factor: <b>1.0000</b> |               |
| Intercept: <b>0.0000</b>          |               |

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| Configure assay parameters — Results |                                   |                                 |                                          |                                      |
|--------------------------------------|-----------------------------------|---------------------------------|------------------------------------------|--------------------------------------|
| <input type="radio"/> General        | <input type="radio"/> Calibration | <input type="radio"/> SmartWash | <input checked="" type="radio"/> Results | <input type="radio"/> Interpretation |
| Assay: <b>FER</b>                    |                                   |                                 | Assay number: <b>2906</b>                |                                      |
| Dilution default range:              |                                   |                                 | Result units: <b>µg/L</b>                |                                      |
| Low-Linearity: <b>10.0</b>           |                                   |                                 |                                          |                                      |
| High-Linearity: <b>500.0</b>         |                                   |                                 |                                          |                                      |
| Gender and age specific ranges:*     |                                   |                                 |                                          |                                      |
| GENDER                               | AGE (UNITS)                       | NORMAL                          | EXTREME                                  |                                      |

| Configure result units            |               |
|-----------------------------------|---------------|
| Assay: <b>FER</b>                 |               |
| Version: <b>†</b>                 |               |
| Result units: <b>µg/L</b>         |               |
| Decimal places: <b>1</b>          | [Range 0 – 4] |
| Correlation factor: <b>1.0000</b> |               |
| Intercept: <b>0.0000</b>          |               |

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# QUANTIA FERRITIN STANDARD

## INTENDED USE

The Quantia Ferritin Standard is intended for use in establishing the calibration for the Quantia Ferritin assay by turbidimetry. This standard is intended to be used together with the Quantia Ferritin reagents. Carefully read the instructions with the corresponding reagents before use.

## PRINCIPLE

When the standards, which contains ferritin, are mixed with the Quantia Ferritin reagent, a clear agglutination occurs which can be measured by turbidimetry.

## CONTENTS

[REF] 6K49-03 Quantia Ferritin Standard, [CAL 1-4] 1 x 1 mL

This product is of human origin and is standardized to the WHO International Standard. The concentrations in ng/mL are indicated on the vial labels. Contains sodium azide (< 0.1%).

**NOTE:** Substituting components from other manufacturers may affect the results of the assay.

## PRECAUTIONS

- [IVD]
- For *In Vitro* Diagnostic Use.
- Do not use components beyond the expiration date.
- Do not mix materials from different kit lot numbers.
- Do not mix fresh reagents with in-use reagents.
-  **CAUTION:** This product contains human sourced and/or potentially infectious components. Refer to the **CONTENTS** section of this package insert. No known test method can offer complete assurance that products derived from human sources or inactivated microorganisms will not transmit infection. Therefore, all human sourced materials should be considered potentially infectious. It is recommended that these standards and human specimens be handled in accordance with the OSHA Standard on Bloodborne Pathogens.<sup>1</sup> Biosafety Level 2<sup>2</sup> or other appropriate biosafety practices<sup>3,4</sup> should be used for materials that contain or are suspected of containing infectious agents. Donor units of human sourced material have been tested and found to be nonreactive for HBsAg, anti-HIV-1/HIV-2, and anti-HCV.
- The following warnings and precautions apply to Quantia Ferritin Standard:
  - Contains sodium azide.
  - EUH032 Contact with acids liberates very toxic gas.
  - P501 Dispose of contents / container in accordance with local regulations.
- Safety Data Sheets are available at [www.abbottiagnostics.com](http://www.abbottiagnostics.com) or contact your local representative.
- For a detailed discussion of safety precautions during system operation, refer to the **ARCHITECT System Operations Manual, Section 8**.

en

Quantia Ferritin Standard

[REF] 6K49-03

G10813/R01

S6K4Q0

[FOR USE WITH]

ARCHITECT

## STORAGE

Unopened standards are stable until the expiration date shown on the label if stored at 2 to 8°C.

All Quantia products contain preservatives, but are still susceptible to contamination. Handle with normal precautions.

## PREPARATION

The standards are ready-to-use.

## QUALITY CONTROL

After calibration, the control values should be within the established range. Otherwise, corrective measures should be taken by the user.

## REFERENCES

1. US Department of Labor, Occupational Safety and Health Administration. 29 CFR Part 1910.1030. Bloodborne pathogens.
2. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*, 5th ed. Washington, DC: US Government Printing Office, December 2009.
3. World Health Organization. *Laboratory Biosafety Manual*, 3rd ed. Geneva: World Health Organization, 2004.
4. Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.

| Key to Symbols |                                                                     |
|----------------|---------------------------------------------------------------------|
|                | Calibrator 1                                                        |
|                | Calibrators 1 through 4                                             |
|                | Calibrator 2                                                        |
|                | Calibrator 3                                                        |
|                | Calibrator 4                                                        |
|                | Contains sodium azide. Contact with acids liberates very toxic gas. |
|                | Distributed by                                                      |
|                | Identifies products to be used together                             |
|                | <i>In Vitro</i> Diagnostic Medical Device                           |
|                | Batch code/Lot number                                               |
|                | Product of Spain                                                    |
|                | Catalog number/List number                                          |
|                | Caution                                                             |
|                | Consult instructions for use                                        |
|                | Manufacturer                                                        |
|                | Temperature limitation                                              |
|                | Use by/Expiration date                                              |

# QUANTIA FERRITIN STANDARD

6K49-03

Biokit, S.A.  
Av. Can Montcau 7  
08186 Lliçà d'Amunt  
Barcelona, Spain

ABBOTT  
65205 Wiesbaden, Germany

September 2017

# QUANTIA FERRITIN/MYOGLOBIN/IgE CONTROL

October 2017

## INTENDED USE

The Quantia Ferritin/Myoglobin/IgE Control is intended for use in monitoring the quality control of results obtained with the Quantia (Ferritin, Myoglobin, and IgE) assays by turbidimetry. This control is intended to be used together with the Quantia (Ferritin, Myoglobin, and IgE) reagents. Carefully read the instructions with the corresponding reagents before use.

## PRINCIPLE

When the controls, which contain ferritin, myoglobin, and IgE are mixed with the corresponding Quantia reagents, a clear agglutination occurs which can be measured by turbidimetry.

## CONTENTS

**[REF]** 6K56-02 Quantia Ferritin/Myoglobin/IgE Control:

|                     |            |               |          |
|---------------------|------------|---------------|----------|
| <b>[CONTROL I]</b>  | low level  | (lyophilized) | 1 x 3 mL |
| <b>[CONTROL II]</b> | high level | (lyophilized) | 1 x 3 mL |

This product is of human origin. Refer to the control data sheet for assayed values. Contains sodium azide (0.2% lyophilized, <0.01% reconstituted).

**NOTE:** Substituting components from other manufacturers may affect the results of the assay.

## PRECAUTIONS

- **[IVD]**
- For *In Vitro* Diagnostic Use.
- Do not use components beyond the expiration date.
- Do not mix materials from different kit lot numbers.
- Do not mix fresh reagents with in-use reagents.
-  **CAUTION:** This product contains human sourced and/or potentially infectious components. Refer to the **CONTENTS** section of this package insert. No known test method can offer complete assurance that products derived from human sources or inactivated microorganisms will not transmit infection. Therefore, all human sourced materials should be considered potentially infectious. It is recommended that these controls and human specimens be handled in accordance with the OSHA Standard on Bloodborne Pathogens.<sup>1</sup> Biosafety Level 2<sup>2</sup> or other appropriate biosafety practices<sup>3,4</sup> should be used for materials that contain or are suspected of containing infectious agents.  
Donor units of human-sourced materials used in the controls have been tested and found to be nonreactive for HBsAg, anti-HIV-1/HIV-2, and anti-HCV.
- The following warnings and precautions apply to **[CONTROL I]** and **[CONTROL II]**:  
Contains sodium azide.  
EUH032 Contact with acids liberates very toxic gas.  
P501 Dispose of contents / container in accordance with local regulations.
- Safety Data Sheets are available at [www.abbottdiagnostics.com](http://www.abbottdiagnostics.com) or contact your local representative.
- For a detailed discussion of safety precautions during system operation, refer to the **ARCHITECT System Operations Manual, Section 8**.

en

Quantia Ferritin/Myoglobin/IgE Control

**[REF]** 6K56-02

**G10834/R01**  
**C6K5T0**

**FOR USE WITH**

**ARCHITECT**

## STORAGE

Quantia Ferritin/Myoglobin/IgE Control is stable until the expiration date shown on the label if stored at 2 to 8°C.

All Quantia products contain preservatives, but are still susceptible to contamination. Handle with normal precautions.

Once reconstituted, the controls can be used for 7 days if stored at 2 to 8°C.

## PREPARATION

The controls should be reconstituted with 3 mL of distilled water. Allow reconstituted material to stand for 5 minutes. Prior to use, gently swirl the vial.

## QUALITY CONTROL

The control values should be within the established range. Otherwise, corrective measures should be taken by the user.

## REFERENCES

1. US Department of Labor, Occupational Safety and Health Administration. 29 CFR Part 1910.1030. Bloodborne pathogens.
2. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*, 5th ed. Washington, DC: US Government Printing Office, December 2009.
3. World Health Organization. *Laboratory Biosafety Manual*, 3rd ed. Geneva: World Health Organization, 2004.
4. Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.

### Key to Symbols

|                                                                                     |                                                                     |
|-------------------------------------------------------------------------------------|---------------------------------------------------------------------|
| <b>CONTAINS: AZIDE</b>                                                              | Contains sodium azide. Contact with acids liberates very toxic gas. |
| <b>CONTROL</b>                                                                      | Control                                                             |
| <b>CONTROL I</b>                                                                    | Control I                                                           |
| <b>CONTROL II</b>                                                                   | Control II                                                          |
| <b>DISTRIBUTED BY</b>                                                               | Distributed by                                                      |
| <b>FOR USE WITH</b>                                                                 | Identifies products to be used together                             |
| <b>IVD</b>                                                                          | <i>In Vitro</i> Diagnostic Medical Device                           |
| <b>LOT</b>                                                                          | Batch code/Lot number                                               |
| <b>MEAN</b>                                                                         | Mean Value                                                          |
| <b>PRODUCT OF JAPAN</b>                                                             | Product of Japan                                                    |
| <b>RANGE</b>                                                                        | Range                                                               |
| <b>REAGENT</b>                                                                      | Reagent                                                             |
| <b>REF</b>                                                                          | Catalog number/List number                                          |
| <b>SI UNITS</b>                                                                     | SI Units                                                            |
|  | Caution                                                             |
|  | Consult instructions for use                                        |
|  | Manufacturer                                                        |
|  | Temperature limitation                                              |
|  | Use by/Expiration date                                              |

## QUANTIA FERRITIN/MYOGLOBIN/IgE CONTROL

**REF 6K56-02**

**CE**  Biokit, S.A.  
Av. Can Montcau 7  
08186 Lliçà d'Amunt  
Barcelona, Spain

**DISTRIBUTED BY**

ABBOTT  
65205 Wiesbaden, Germany

 **Abbott**

October 2017

# Multichem® S Plus

Replace up to four competitor chemistry QC products with Multichem S Plus

Corresponding Competitor Products

Multichem  
S Plus

## Features

- Human-based matrix proprietary formulation.
- Mimics performance of patient samples.
- Targeted at key clinical decision points.
- Liquid for ease of use.

## Specifications

- 36 month closed vial stability at -20°C to -80°C.
- 10 day open vial stability at 2°C to 8°C.

## Instruments

Multichem S Plus, independent clinical chemistry QC, is suitable for use on all major platforms.

| Order Information                             |              |           |
|-----------------------------------------------|--------------|-----------|
| Multichem S Plus                              | Config.      | Part Code |
| Assayed Single Level (Level 1)                | 15 x 10mL    | CH101PLA  |
| Assayed Single Level (Level 2)                | 15 x 10mL    | CH102PLA  |
| Assayed Single Level (Level 3)                | 15 x 10mL    | CH103PLA  |
| Unassayed Single Level (Level 1)              | 15 x 10mL    | CH101CRP  |
| Unassayed Single Level (Level 2)              | 15 x 10mL    | CH102CRP  |
| Unassayed Single Level (Level 3)              | 15 x 10mL    | CH103CRP  |
| Mini kit<br>(contains one vial of each level) | 3 x 1 x 10mL | CH100MCRP |

## Analyte List

### Chemistry

Albumin  
Bilirubin, Direct  
Bilirubin, Total  
Calcium  
Carbon Dioxide (Bicarbonate)  
Chloride  
Creatinine  
Glucose  
Iron  
Lactic acid  
Magnesium  
Phosphorous  
Potassium  
Protein, Total  
Sodium  
Total Iron Binding Capacity\*  
Unsaturated Iron Binding Capacity (UIBC)  
Urea  
Uric Acid

### Esoterics

ACE\*  
Bile Acids  
Bili, Indirect\*  
Caffeine\*  
Calcium, Ionized\*  
Copper\*  
Cortisol  
Ethanol  
Fructosamine\*  
NT-proBNP\*  
Osmolality\*  
Protein Electrophoresis\*  
Triiodothyronine\*  
Thyroxine, (Total T4)  
Troponin T\*  
Zinc\*

### Immunoproteins

Alpha-1 Aminoglycoprotein  
Alpha-1 Antitrypsin  
Alpha 2-macroglobulin\*  
ASO\*  
ADNase B (Anti-Streptococcal DNase B)\*  
Antithrombin III\*  
Apolipoprotein A1  
Apolipoprotein B  
Beta-2 Microglobulin  
C1 Inhibitor\*  
CH50 (Total hemolytic Complement)\*  
Cystatin C\*  
Complement C3  
Complement C4  
Ceruloplasmin  
C-Reactive Protein  
Ferritin\*  
Haptoglobin  
Hemopexin\*  
Immunoglobulin A  
Immunoglobulin G  
Immunoglobulin M  
IgE\*  
IgG1, Subclass\*  
IgG2, Subclass\*  
IgG3, Subclass\*  
IgG4, Subclass\*  
Kappa Light Chain\*  
Lambda Light Chain\*  
Lipoprotein (a)\*  
Prealbumin  
Properdin Factor B\*  
Retinol Binding Protein\*  
Rheumatoid Factor\*  
Transferrin  
sTfR (Soluble Transferrin Receptor)\*

### Enzymes

Acid Phosphatase  
Alanine Aminotransferase (ALT)  
Alkaline Phosphatase (ALP)  
Amylase (Pancreatic)  
Amylase (Total)  
Aspartate Aminotransferase (AST)  
Alpha-hydroxybutyrate Dehydrogenase  
(α-HBDH)\*  
Beta-hydroxybutyrate Dehydrogenase\*  
Cholinesterase  
Creatine Phosphokinase  
CK-MB\*  
Gamma Glutamyltransferase  
Lactate Dehydrogenase (L to P)  
Lipase  
Prostatic Acid Phosphatase\*

### Lipids

Cholesterol, HDL  
Cholesterol, LDL  
Cholesterol, Total  
Phospholipids\*  
Triglycerides

### Therapeutic Drugs

Acetaminophen  
Amikacin  
Carbamazepine  
Digoxin  
Gentamicin  
Lithium  
Phenobarbital  
Phenytoin  
Salicylic Acid  
Theophylline  
Tobramycin  
Valproic Acid  
Vancomycin

\*Please refer to lot specific package inserts for stability and performance claims.

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- 04R1001 Alinity ci-series Calibrator/Control Replacement Caps

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## BIBLIOGRAFIE

- US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
- US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
- World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.

- Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.

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©2021, 2022 Abbott Laboratories