

# Elecsys Tacrolimus

cobas®

REF



SYSTEM

07251254190

07251254501

300

cobas e 402  
cobas e 801

## English

For use in the USA only

## System information

Short name	ACN (application code number)
TCL	10022

## Intended use

Immunoassay for the in vitro quantitative determination of tacrolimus in EDTA human whole blood. The assay is used as an aid in the management of liver and kidney transplant patients receiving tacrolimus therapy.

The electrochemiluminescence immunoassay "ECLIA" is intended for use on **cobas e** immunoassay analyzers.

## Summary

Tacrolimus (also referred to as FK506) is a macrolide antibiotic identified as a product of the actinobacterium *Streptomyces tsukubaensis* in Japan in 1984.<sup>1,2,3</sup> Studies demonstrated that tacrolimus is 10-100 times more active than cyclosporine at inhibiting several immune responses.<sup>4</sup>

The main mechanism through which tacrolimus exerts its immunosuppressive effect is believed to be via the inhibition of T cell activation and proliferation. Intracellular tacrolimus binds an immunophilin called FK506-binding protein (FKBP-12) and these complexes then inhibit the enzymatic activity of calcineurin.<sup>5</sup> The inhibition of calcineurin restricts the dephosphorylation and nuclear translocation of nuclear factor of activated T cells (NFAT), which regulates transcription of several cytokines, including IL-2, IL-4, TNF- $\alpha$ , and interferon- $\gamma$ , and therefore limits lymphocyte activation and proliferation.<sup>6,7,8,9,10</sup>

Tacrolimus is highly lipophilic and absorption is incomplete and variable. Following absorption, tacrolimus is highly bound to proteins and erythrocytes, with 99 % of the drug within the plasma being bound to albumin or  $\alpha$ -1-glycoprotein.<sup>11</sup>

The bioavailability and metabolism of tacrolimus are predominantly influenced by the activity of the cytochrome P450 isozymes CYP3A4 and CYP3A5, as well as the efflux pump p-glycoprotein, which show significant inter- and intra-individual variability in expression and function.<sup>12,13,14</sup>

Tacrolimus displays a high degree of inter- and intra-patient variability, as well as potentially severe side effects from doses that are either too low or too high. Inadequate tacrolimus concentrations might result in rejection of the transplanted organ. High levels may lead to severe adverse effects. Principle adverse effects associated with tacrolimus include nephrotoxicity, neurotoxicity, gastrointestinal disturbances, diabetogenesis, hypertension and malignant complications.<sup>15,16</sup>

The application of therapeutic drug monitoring (TDM) and concentration-controlled dosing in order to maintain each patient's drug exposure within a narrow therapeutic window is part of standard clinical practice for many years and is a major support to patient management.<sup>16,17</sup> Trough concentration (C<sub>0</sub>) monitoring is still widely used as a guide to individualizing tacrolimus dose requirements, even though some controversies remain about the relationship between C<sub>0</sub> and clinical outcome. The area under the concentration-time curve (AUC<sub>0-12</sub>) is generally considered the best marker of exposure but is expensive and impractical. To assess the efficacy of alternative strategies to C<sub>0</sub>, multicenter prospective trials would be needed.<sup>16</sup>

## Test principle

Manual precipitation:

Before testing with the Elecsys Tacrolimus assay, samples, calibrators and controls must be **pretreated** with Elecsys ISD Sample Pretreatment.

The reagent lyses the cells, extracts tacrolimus, and precipitates most of the blood proteins. The **pretreated** samples are centrifuged, and an aliquot of the resulting supernatant containing tacrolimus is then assayed using the Elecsys Tacrolimus assay.

Competition principle. Total duration of assay: 18 minutes.

- 1st incubation: 21  $\mu$ L of pretreated sample is incubated with a tacrolimus-specific biotinylated antibody and a ruthenium complex<sup>a)</sup> labeled tacrolimus-derivate. Depending on the analyte concentration in the sample and the formation of the respective immune complex, the labeled antibody binding site is occupied in part with sample analyte and in part with ruthenylated hapten.
- 2nd incubation: After addition of streptavidin-coated microparticles, the complex becomes bound to the solid phase via interaction of biotin and streptavidin.
- The reaction mixture is aspirated into the measuring cell where the microparticles are magnetically captured onto the surface of the electrode. Unbound substances are then removed with ProCell II M. Application of a voltage to the electrode then induces chemiluminescent emission which is measured by a photomultiplier.
- Results are determined via a calibration curve which is instrument-specifically generated by 2-point calibration and a master curve provided via the **cobas link**.

a) Tris(2,2'-bipyridyl)ruthenium(II)-complex (Ru(bpy)<sub>3</sub><sup>2+</sup>)

## Reagents - working solutions

The **cobas e** pack is labeled as TCL.

- M Streptavidin-coated microparticles, 1 bottle, 12.4 mL:  
Streptavidin-coated microparticles 0.72 mg/mL; preservative.
- R1 Anti-Tacrolimus-S-Ab-biotin, 1 bottle, 21 mL:  
Biotinylated monoclonal anti-tacrolimus-antibody (sheep) 15  $\mu$ g/L;  
phosphate buffer 100 mmol/L, pH 7.8; preservative.
- R2 Tacrolimus~Ru(bpy)<sub>3</sub><sup>2+</sup>, 1 bottle, 14.8 mL:  
Tacrolimus-derivative labeled with ruthenium complex 4  $\mu$ g/L; citrate  
buffer 10 mmol/L, pH 3.3; preservative.

## Precautions and warnings

For in vitro diagnostic use for health care professionals. Exercise the normal precautions required for handling all laboratory reagents.

Infectious or microbial waste:

Warning: handle waste as potentially biohazardous material. Dispose of waste according to accepted laboratory instructions and procedures.

Environmental hazards:

Apply all relevant local disposal regulations to determine the safe disposal.

Safety data sheet available for professional user on request.

For USA: Caution: Federal law restricts this device to sale by or on the order of a physician.

This kit contains components classified as follows in accordance with the Regulation (EC) No. 1272/2008:



## Warning

- H317 May cause an allergic skin reaction.

## Prevention:

- P261 Avoid breathing mist or vapours.
- P272 Contaminated work clothing should not be allowed out of the workplace.
- P280 Wear protective gloves.

## Response:

- P333 + P313 If skin irritation or rash occurs: Get medical advice/attention.

# Elecsys Tacrolimus



P362 + P364 Take off contaminated clothing and wash it before reuse.

## Disposal:

P501 Dispose of contents/container to an approved waste disposal plant.

Product safety labeling follows EU GHS guidance.

Contact phone: 1-800-428-2336

Avoid foam formation in all reagents and sample types (specimens, calibrators and controls).

## Reagent handling

The reagents in the kit have been assembled into a ready-for-use unit that cannot be separated.

All information required for correct operation is available via the **cobas** link.

## Storage and stability

Store at 2-8 °C.

Do not freeze.

Store the **cobas e** pack **upright** in order to ensure complete availability of the microparticles during automatic mixing prior to use.

Stability:	
unopened at 2-8 °C	up to the stated expiration date
on the analyzers	16 weeks

## Specimen collection and preparation

Only the specimens listed below were tested and found acceptable.

K<sub>2</sub>-EDTA and K<sub>3</sub>-EDTA whole blood.

Specimens collected in EDTA tubes may be stored for up to 5 days at 15-25 °C or 7 days at 2-8 °C prior to being tested. If testing will be delayed by more than 7 days, store frozen at -20 °C (± 5 °C) or lower for up to 3 months. Freeze only once. Specimens must be mixed thoroughly after thawing to ensure consistency of the results.

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Mix thawed specimens thoroughly by hand or on a roller mixer or rocker. Visually inspect the specimens. If layering or stratification is observed, continue mixing until specimens are visibly homogeneous.

Do not use heat-inactivated samples.

Do not use samples and controls stabilized with azide.

Ensure the samples and calibrators are at 20-25 °C prior to measurement.

Use caution when handling patient specimens to prevent cross contamination. Use of disposable pipettes or pipette tips is recommended.

**Pretreated samples can be stored in closed tubes for up to 4 hours at 20-25 °C.**

**Due to evaporation effects, pretreated samples should be analyzed/measured within 30 minutes after opening the vials and loading the samples on the analyzer. Avoid delays between loading and measurement to ensure the 30 minute stability of pretreated samples.**

**A re-run requires repeating of the manual pretreatment procedure.**

Sample stability claims were established by experimental data by the manufacturer or based on reference literature and only for the temperatures/time frames as stated in the method sheet. It is the responsibility of the individual laboratory to use all available references and/or its own studies to determine specific stability criteria for its laboratory.

## Materials provided

See "Reagents – working solutions" section for reagents.

## Materials required (but not provided)

- [REF] 05889073190, ISD Sample Pretreatment, 1 x 30 mL
- [REF] 05889065190, Tacrolimus CalSet, for 6 x 1.0 mL

- [REF] 05889081160, PreciControl ISD, for 3 x 3.0 mL
- [REF] 11776576322, CalSet Vials, 2 x 56 empty snap-cap bottles
- [REF] 07299001190, Diluent Universal, 36 mL sample diluent
- General laboratory equipment
- Precision pipettes (use only positive displacement pipettes for ISD Sample Pretreatment reagent handling)
- Microcentrifuge tubes (2.0 mL capacity)
- Microcentrifuge (at least 10000 g)
- Vortex mixer
- Roller mixer or rocker
- cobas e** analyzer

Additional materials for the **cobas e** 402 and **cobas e** 801 analyzers:

- [REF] 06908799190, ProCell II M, 2 x 2 L system solution
- [REF] 04880293190, CleanCell M, 2 x 2 L measuring cell cleaning solution
- [REF] 07485409001, Reservoir Cup, 8 cups to supply ProCell II M and CleanCell M
- [REF] 06908853190, PreClean II M, 2 x 2 L wash solution
- [REF] 05694302001, Assay Tip/Assay Cup tray, 6 magazines x 6 magazine stacks x 105 assay tips and 105 assay cups, 3 wasteliners
- [REF] 07485425001, Liquid Flow Cleaning Cup, 2 adaptor cups to supply ISE Cleaning Solution/Elecsys SysClean for Liquid Flow Cleaning Detection Unit
- [REF] 07485433001, PreWash Liquid Flow Cleaning Cup, 1 adaptor cup to supply ISE Cleaning Solution/Elecsys SysClean for Liquid Flow Cleaning PreWash Unit
- [REF] 11298500160, ISE Cleaning Solution/Elecsys SysClean, 5 x 100 mL system cleaning solution

## Manual specimen pretreatment

Follow the steps listed below to pretreat calibrators, controls and/or specimens. **The technical notes hereafter are an essential part of the instructions and must be read thoroughly before completing each step.** Follow steps 1 through 7 to pretreat calibrators, controls and/or specimens.

Steps	Technical notes
1. Equilibrate all reagents, calibrators, controls and specimens to 20-25 °C. Mix all calibrators, controls and specimens gently but thoroughly just before use.	Do not vortex. The liquids may be mixed by hand or on a roller mixer or rocker.  The calibrators and controls are a whole-blood hemolysate and may be slightly different in appearance from whole-blood samples.
2. Label one microcentrifuge tube for each calibrator, control and/or specimen to be pretreated.	none
3. Using a precision pipette, transfer 300 µL of each calibrator, control and/or specimen to the appropriately labeled micro-centrifuge tube.	Use a fresh pipette tip for each calibrator, control and/or specimen.
4. Using a positive displacement pipette, add 300 µL of ISD Sample Pretreatment reagent to each microcentrifuge tube. Immediately cap each tube and immediately proceed to step 5.	Note: ISD Sample Pretreatment is highly volatile. Keep tightly closed when not in use to prevent evaporation.

# Elecsys Tacrolimus

Steps	Technical notes
5. Vortex each microcentrifuge tube for at least 10 seconds. Failure to perform this step may result in a supernatant that appears red. See step 6, technical note.	Note: Failure to vortex each tube immediately after addition of the ISD Sample Pretreatment reagent will lead to erroneous assay results. Sample and reagent mixture should be completely homogeneous immediately after vortexing. Visual inspection is required.
6. Centrifuge the samples for at least 4 minutes in a micro-centrifuge ( $\geq 10000$ g).	The centrifuged samples should have well-defined pellets and clear supernatant. The supernatant should not appear cloudy or red. If the supernatant is red, discard and replace it with a newly extracted sample.
7. Transfer each supernatant directly into an appropriate vial and immediately cap each vial. The samples are ready to be assayed.	Pretreated samples can be stored in closed tubes for up to 4 hours at 20-25 °C. <b>Please note:</b> <b>Due to evaporation effects, pretreated samples should be analyzed/measured within 30 minutes after opening the vials and loading the samples on the system. Avoid delays between loading and measurement to ensure the 30 minutes stability of pretreated samples.</b> <b>This is supported by running the tacrolimus samples in batch mode:</b> <b>Based on average system sample processing time, no more than 35 tacrolimus samples may be loaded per calibrated measuring cell onto the analyzers at the same time.</b>

## Assay

For optimum performance of the assay follow the directions given in this document for the analyzer concerned. Refer to the appropriate operator's manual for analyzer-specific assay instructions.

Resuspension of the microparticles takes place automatically prior to use.

Place the cooled (stored at 2-8 °C) **cobas e** pack on the reagent manager. Avoid foam formation. The system automatically regulates the temperature of the reagents and the opening/closing of the **cobas e** pack.

## Calibration

**Traceability:** This method has been standardized against reference standards traceable to tacrolimus reference material (USP = United States Pharmacopeia) by weight.

The predefined master curve is adapted to the analyzer using the relevant CalSet.

**Tacrolimus CalSet must be pretreated freshly before calibration.**

**Calibration frequency:** Calibration must be performed once per reagent lot using fresh reagent (i.e. not more than 24 hours since the **cobas e** pack was registered on the analyzer).

Calibration interval may be extended based on acceptable verification of calibration by the laboratory.

Renewed calibration is recommended as follows:

- after 12 weeks when using the same reagent lot
- after 28 days when using the same **cobas e** pack on the analyzer

- as required: e.g. quality control findings outside the defined limits

## Quality control

For quality control, use PreciControl ISD.

**PreciControl ISD must be pretreated freshly before measurement.**

In addition, other suitable control material can be used.

Controls for the various concentration ranges should be run individually at least once every 24 hours when the test is in use, once per **cobas e** pack, and following each calibration.

The control intervals and limits should be adapted to each laboratory's individual requirements. Values obtained should fall within the defined limits. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

If necessary, repeat the measurement of the samples concerned.

Follow the applicable government regulations and local guidelines for quality control.

## Calculation

The analyzer automatically calculates the analyte concentration of each sample (either in ng/mL, nmol/L, µg/L).

Conversion factors:  $\text{ng/mL} \times 1.0 = \mu\text{g/L}$

$\text{ng/mL} \times 1.2438 = \text{nmol/L}$

## Limitations - interference

The effect of the following endogenous substances, pharmaceutical compounds and clinical conditions on assay performance was tested. Interferences were tested up to the listed concentrations and no impact to results was observed.

Criterion: Recovery within  $\pm 0.3$  ng/mL (concentration range  $> 0.75$ -3 ng/mL) or within  $\pm 10$  % (concentration range  $> 3$ -30 ng/mL) of initial value.

## Endogenous substances

Compound	Concentration tested
Albumin	$\leq 12.0$ g/dL
Bilirubin	$\leq 1026$ µmol/L or $\leq 60.0$ mg/dL
Cholesterol	$\leq 500$ mg/dL
HASAb)	$\leq 10.0$ µg/mL
Hematocrit	25-60 %
IgG	$\leq 12.0$ g/dL
Intralipid	$\leq 1500$ mg/dL
Rheumatoid factors	up to 500 IU/mL
Uric acid	$\leq 20.0$ mg/dL

b) Human anti-sheep antibody

## Biotin interference

% Bias for samples containing various concentrations of biotin					
Samples (ng/mL)	Biotin concentration (ng/mL)				
	100	200	250	300	350
2.24	+1.2 %	+7.1 %	+11.4 %	+15.4 %	+17.6 %
18.2	+2.6 %	+4.5 %	+7.9 %	+10.2 %	+12.5 %

% Bias for samples containing various concentrations of biotin				
Samples (ng/mL)	Biotin concentration (ng/mL)			
	400	600	1000	1200
2.24	+20.4 %	+57.9 %	+122.6 %	+147.9 %
18.2	+12.1 %	+26.7 %	+64.0 %	+79.1 %

Specimens with biotin concentrations up to 100 ng/mL demonstrated  $\leq 10$  % bias in results. Biotin concentrations greater than 100 ng/mL led to higher positive bias Tacrolimus results. Some studies have shown that serum concentrations of biotin can reach up to 355 ng/mL within the first

# Elecsys Tacrolimus

hour after biotin ingestion for subjects consuming supplements of 20 mg biotin per day<sup>18</sup> and up to 1160 ng/mL for subjects after a single dose of 300 mg biotin.<sup>19</sup>

Do not test samples from patients who are taking biotin.

## Pharmaceutical compounds:

In vitro tests were performed on 16 commonly used pharmaceutical compounds. Criterion: Recovery within  $\pm 10\%$  of initial value.

25 special drugs were additionally tested. The compounds tested and shown in the table below were found not to interfere with the assay except Itraconazole (INN international non-proprietary name) recovery was found to interfere at 10  $\mu\text{g/mL}$  (116 %).

Drug	Concentration tested
Acetaminophen	200 mg/L
N-Acetylcysteine	1660 mg/L
Acetylsalicylic acid	1000 mg/L
Acyclovir	3.2 $\mu\text{g/mL}$
Amphotericin B	5.8 $\mu\text{g/mL}$
Ampicillin-Na	1000 mg/L
Ascorbic acid	300 mg/L
Cefoxitin	2500 mg/L
Ciprofloxacin	7.4 $\mu\text{g/mL}$
Cyclosporine	5000 ng/mL
Doxycycline	50 mg/L
Erythromycin	20 mg/dL
Everolimus	60 ng/mL
Fluconazole	30 $\mu\text{g/mL}$
Flucytosine	40 $\mu\text{g/mL}$
Gancyclovir	1000 $\mu\text{g/mL}$
Gentamicin	12 mg/dL
Heparin	5000 IU/L
Ibuprofen	500 mg/L
K <sub>2</sub> -EDTA	6 mg/mL
K <sub>3</sub> -EDTA	6 mg/mL
Kanamycin	100 $\mu\text{g/mL}$
Ketoconazole	50 $\mu\text{g/mL}$
Levodopa	20 mg/L
Lidocaine	6 mg/dL
Methyldopa+1.5	20 mg/L
Metronidazole	200 mg/L
MPA (mycophenolic acid) glucuronide	1800 $\mu\text{g/mL}$
Mycophenolic acid	500 $\mu\text{g/mL}$
Nitrofurantoin	6 $\mu\text{g/mL}$
Phenobarbital	15 mg/dL
Phenylbutazone	400 mg/L
Rifampicin	6 mg/dL
Sirolimus	60 ng/mL
Spectinomycin	100 $\mu\text{g/mL}$
Sulfomethoxazole	200 $\mu\text{g/mL}$
Theophylline	100 mg/L
Tobramycin	2 mg/dL
Trimethoprim	40 $\mu\text{g/mL}$

Drug	Concentration tested
Vancomycin	6 mg/dL

In rare cases, interference due to extremely high titers of antibodies to analyte-specific antibodies, streptavidin or ruthenium can occur. These effects are minimized by suitable test design.

For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

## Limits and ranges

### Measuring range

0.75-30 ng/mL (defined by the Limit of Quantitation and the maximum of the master curve). Values below the Limit of Quantitation are reported as < 0.75 ng/mL. Values above the measuring range are reported as > 30 ng/mL.

### Lower limits of measurement

*Limit of Blank, Limit of Detection and Limit of Quantitation*

Limit of Blank = 0.3 ng/mL

Limit of Detection = 0.5 ng/mL

Limit of Quantitation = 0.75 ng/mL with a total allowable error of  $\leq 25\%$ .

The Limit of Blank, Limit of Detection and Limit of Quantitation were determined in accordance with the CLSI (Clinical and Laboratory Standards Institute) EP17-A2 requirements.

The Limit of Blank is the 95<sup>th</sup> percentile value from  $n \geq 60$  measurements of analyte-free samples over several independent series. The Limit of Blank corresponds to the concentration below which analyte-free samples are found with a probability of 95 %.

The Limit of Detection is determined based on the Limit of Blank and the standard deviation of low concentration samples. The Limit of Detection corresponds to the lowest analyte concentration which can be detected (value above the Limit of Blank with a probability of 95 %).

The Limit of Quantitation is defined as the lowest amount of analyte in a sample that can be accurately quantitated with a total allowable relative error of  $\leq 25\%$ .

### Dilution

Specimens with a tacrolimus concentration of > 30 ng/mL will be flagged as "> 30 ng/mL" and may be diluted with the Manual Dilution Procedure.

Manual dilutions should be performed as follows:

- The suggested dilution for the Elecsys Tacrolimus assay is 1:3.
- Specimen must be diluted manually before pretreatment.
- Add 100  $\mu\text{L}$  of the patient specimen to 200  $\mu\text{L}$  of Diluent Universal, mix thoroughly and then proceed with the Manual Pretreatment Procedure in the procedure section.

The pre-diluted and pretreated sample must be handled like any other routine sample. Hence, customer will be given the "diluted" result which must be manually multiplied with the dilution factor. For detailed information, refer to the system Operator manual.

### Expected values

No firm therapeutic range exists for tacrolimus in whole blood. The complexity of the clinical state, individual differences in sensitivity to immunosuppressive and nephrotoxic effects of tacrolimus, coadministration of other immunosuppressants, type of transplant, time post-transplant, and a number of other factors contribute to different requirements for optimal blood levels of tacrolimus. Individual tacrolimus values cannot be used as the sole indicator for making changes in the treatment regimen. Each patient should be thoroughly evaluated clinically before treatment adjustments are made, and each assay user must establish his or her ranges based on clinical experience.

These ranges will vary according to the commercial in vitro diagnostic test used. Ranges must be established for each commercial test used.

### Specific performance data

Representative performance data on the analyzers are given below. Results obtained in individual laboratories may differ.

### Precision

Precision was determined using Elecsys reagents, samples and controls in a protocol (EP05-A3) of the CLSI (Clinical and Laboratory Standards

# Elecsys Tacrolimus

Institute): 2 runs per day in duplicate each for 21 days (n = 84). The following results were obtained:

cobas e 402 and cobas e 801 analyzers					
		Repeatability		Intermediate precision	
Sample	Mean ng/mL	SD ng/mL	CV %	SD ng/mL	CV %
HSP <sup>c)</sup> 1	1.27	0.056	4.4	0.069	5.4
HSP 2	3.68	0.114	3.1	0.123	3.3
HSP 3	8.03	0.251	3.1	0.303	3.8
HSP 4	14.4	0.329	2.3	0.389	2.7
HSP 5	29.3	0.783	2.7	1.05	3.6
PreciControl ISD 1	1.80	0.065	3.6	0.079	4.4
PreciControl ISD 2	8.91	0.193	2.2	0.254	2.9
PreciControl ISD 3	16.2	0.337	2.1	0.401	2.5

c) HSP = Human Sample Pool

## Method comparison

A comparison of the Elecsys Tacrolimus assay, [REF] 07251254190 (cobas e 801 analyzer; y) with the Elecsys Tacrolimus assay, [REF] 05889057190 (cobas e 601 analyzer; x) gave the following correlations (ng/mL):

Number of samples measured: 159

Passing/Bablok<sup>20</sup> Linear regression  
 $y = 0.902x - 0.263$   $y = 0.979x - 0.796$   
 $r = 0.921$   $r = 0.989$

The sample concentrations were between 0.885 - 29.9 ng/mL.

The Elecsys Tacrolimus assay is designed to have a correlation coefficient of  $\geq 0.9$  for specimens between 0.75-30 ng/mL when compared to an automated immunoassay.


A comparison of the Elecsys Tacrolimus assay (y) with the Abbott ARCHITECT Tacrolimus assay (x) using human whole blood EDTA specimens from liver and kidney transplant patients receiving tacrolimus therapy was performed. Regression analysis was performed using the Weighted Deming<sup>21</sup> regression method. The data from this study are summarized in the following table:

Elecsys Tacrolimus vs Abbott ARCHITECT Tacrolimus				
Number of samples		Slope (95% CI <sup>d)</sup>	Intercept (95% CI)	Correlation Coefficient (r)
Combined	553	0.99 (0.98, 1.01)	0.06 (-0.07, 0.18)	0.99
Kidney	346	1.01 (0.98, 1.03)	0.04 (-0.17, 0.24)	0.99
Liver	207	0.97 (0.95, 0.99)	0.13 (-0.01, 0.27)	0.99

d) CI = Confidence Interval

A bias analysis of the Elecsys Tacrolimus vs. Abbott ARCHITECT Tacrolimus assay was performed on the same 553 human whole blood EDTA samples in the range of 0.75 to 30 ng/mL. The following representative data in Bland-Altman Plot is provided to aid in understanding the differences between these two assays. The average ng/mL difference bias exhibited by Elecsys Tacrolimus vs. Abbott ARCHITECT Tacrolimus assay in this study was -0.01 ng/mL. The 95 % confidence interval of the ng/mL difference bias is -0.07 ng/mL to 0.06 ng/mL. The summary of the results is below.

Elecsys Tacrolimus ng/mL difference bias from Abbott ARCHITECT Tacrolimus assay




The Elecsys Tacrolimus assay is designed to have a correlation coefficient of  $\geq 0.9$  for specimens between 0.75-30 ng/mL when compared to the LC-MS/MS method.

A comparison of the Elecsys Tacrolimus assay (y) with an LC-MS/MS method (x) using clinical samples was performed. Regression analysis was performed using the Weighted Deming<sup>21</sup> regression method. The data from this study are summarized in the following table:

Elecsys Tacrolimus vs LC-MS/MS				
Number of samples		Slope (95% CI <sup>d)</sup>	Intercept (95% CI)	Correlation Coefficient (r)
Combined	554	0.92 (0.90, 0.95)	-0.01 (-0.16, 0.14)	0.96
Kidney	344	0.93 (0.90, 0.96)	0.04 (-0.19, 0.27)	0.97
Liver	210	0.91 (0.87, 0.95)	-0.05 (-0.25, 0.16)	0.95

A bias analysis of the Elecsys Tacrolimus vs. LC-MS/MS was performed on the same 554 human whole blood EDTA samples in the range of 0.75 to 30 ng/mL. The following representative data in Bland-Altman Plot is provided to aid in understanding the differences between these two methods. The average ng/mL difference bias exhibited by Elecsys Tacrolimus vs. LC-MS/MS method in this study was -0.81 ng/mL. The 95 % confidence interval of the ng/mL difference bias is -0.94 ng/mL to -0.69 ng/mL. The summary of the results is below.

Elecsys Tacrolimus ng/mL difference bias from LC-MS/MS





# Elecsys Tacrolimus

## Analytical specificity

A study was performed with the Elecsys Tacrolimus assay based on guidance from the CLSI document EP7-A2.

Metabolite	Maximum concentration of metabolite added ng/mL	Maximum cross-reactivity %
M I	50	1
M II	50	74
M III	50	3
M IV	50	1
M V	6.6	52
M VI	50	1
M VII	50	0
M VIII	50	4

## References

- Kino T, Hatanaka H, Hashimoto M, et al. FK506, a novel immunosuppressant isolated from a Streptomyces. I. Fermentation, isolation, and physico-chemical and biological characteristics. J Antibiot (Tokyo) 1987;40(9):1249-1255.
- Kino T, Hatanaka H, Miyata S, et al. FK-506, a novel immunosuppressant isolated from a Streptomyces. II. Immunosuppressive effect of FK-506 in vitro. J Antibiot (Tokyo) 1987;40(9):1256-1265.
- Wallemacq PE, Reding R. FK506 (tacrolimus), a novel immunosuppressant in organ transplantation: clinical, biomedical, and analytical aspects. Clin Chem 1993;39(11 Pt 1):2219-2228.
- Scott LJ, McKeage K, Keam SJ, et al. Tacrolimus: a further update of its use in the management of organ transplantation. Drugs 2003;63(12):1247-1297.
- Siekierka JJ, Hung SH, Poe M, et al. A cytosolic binding protein for the immunosuppressant FK506 has peptidyl-prolyl isomerase activity but is distinct from cyclophilin. Nature 1989;341(6244):755-757.
- Flanagan WM, Corthesy B, Bram RJ, et al. Nuclear association of a T-cell transcription factor blocked by FK-506 and cyclosporin A. Nature 1991;352:803-807.
- Jain J, McCaffrey PG, Miner Z, et al. The T-cell transcription factor NFATp is a substrate for calcineurin and interacts with Fos and Jun. Nature 1993;365:352-355.
- Shaw KT, Ho AM, Raghavan A, et al. Immunosuppressive drugs prevent a rapid dephosphorylation of transcription factor NFAT1 in stimulated immune cells. Proc Natl Acad Sci USA 1995;92:11205-11209.
- Clipstine NA, Crabtree GR. Identification of calcineurin as a key signalling enzyme in T-lymphocyte activation. Nature 1992;357:695-697.
- O'Keefe SJ, Tamura J, Kincaid RL, et al. FK-506 and CsA-sensitive activation of the interleukin-2 promoter by calcineurin. Nature 1992;357:692-694.
- Astellas Pharma. Prograf Package Insert. Available at: [http://www.accessdata.fda.gov/drugsatfda\\_docs/label/2012/050708s038lbl.pdf](http://www.accessdata.fda.gov/drugsatfda_docs/label/2012/050708s038lbl.pdf) [Last accessed 01/03/2012].
- Benet LZ, Cummins CL, Wu CY. Unmasking the dynamic interplay between intestinal efflux transporters and metabolic enzymes. Int J Pharm 2004;277:3-9.
- de Jonge H, Naesens M, Kuypers DR. New insights into the pharmacokinetics and pharmacodynamics of the calcineurin inhibitors and mycophenolic acid: possible consequences for therapeutic drug monitoring in solid organ transplantation. Ther Drug Monit 2009;31:416-435.
- Anglicheau D, Legendre C, Beaune P, et al. Cytochrome P450 3A polymorphisms and immunosuppressive drugs: an update. Pharmacogenomics 2007;8(7):835-849.
- Staatz CE, Tett SE. Clinical pharmacokinetics and pharmacodynamics of tacrolimus in solid organ transplant recipients. Clin Pharmacokinet 2004;43:623-653.
- Wallemacq P, Armstrong VW, Brunet M, et al. Opportunities to optimize tacrolimus therapy in solid organ transplantation: report of the European consensus conference. Ther Drug Monit 2009;31:139-152.
- Schiff J, Cole E, Cantarovich M. Therapeutic monitoring of calcineurin inhibitors for the nephrologist. Clin J Am Soc Nephrol 2007;2:374-384.
- Piketty ML, Prie D, Sedel F, et al. High-dose biotin therapy leading to false biochemical endocrine profiles: validation of a simple method to overcome biotin interference. Clin Chem Lab Med. 2017 May 1;55(6):817-825.
- Grimsey P, Frey N, Bendig G, et al. Population pharmacokinetics of exogenous biotin and the relationship between biotin serum levels and in vitro immunoassay interference. International Journal of Pharmacokinetics. 2017 Sept 14;2(4): 247-256.
- Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem 1988 Nov;26(11):783-790.
- Linnet K. (1990) Estimation of the Linear Relationship between the Measurements of two Methods with Proportional Errors. Statistics in Medicine, Vol. 9, 1463-1473.

For further information, please refer to the appropriate operator's manual for the analyzer concerned, the respective application sheets and the Method Sheets of all necessary components (if available in your country).

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

## Symbols

Roche Diagnostics uses the following symbols and signs in addition to those listed in the ISO 15223-1 standard (for USA: see [dialog. Roche.com](http://dialog. Roche.com) for definition of symbols used):

	Contents of kit
	Analyzers/Instruments on which reagents can be used
	Reagent
	Calibrator
	Volume for reconstitution
	Global Trade Item Number

## FOR US CUSTOMERS ONLY: LIMITED WARRANTY

Roche Diagnostics warrants that this product will meet the specifications stated in the labeling when used in accordance with such labeling and will be free from defects in material and workmanship until the expiration date printed on the label. THIS LIMITED WARRANTY IS IN LIEU OF ANY OTHER WARRANTY, EXPRESS OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR PARTICULAR PURPOSE. IN NO EVENT SHALL ROCHE DIAGNOSTICS BE LIABLE FOR INCIDENTAL, INDIRECT, SPECIAL OR CONSEQUENTIAL DAMAGES.

COBAS, COBAS E, ELECSYS and PRECICONTROL are trademarks of Roche. INTRALIPID is a trademark of Fresenius Kabi AB.

All other product names and trademarks are the property of their respective owners.

Additions, deletions or changes are indicated by a change bar in the margin.

© 2022, Roche Diagnostics



Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim  
[www.Roche.com](http://www.Roche.com)

+800 5505 6606



Distribution in USA by:  
Roche Diagnostics, Indianapolis, IN  
US Customer Technical Support 1-800-428-2336