SmartScore 4.0

Advanced imaging software that detects, quantifies and scores cardiac calcium plaque burden. Instantly.


Cardiovascular disease remains one of the most common health issues in the world today. As with many conditions, early detection and patient risk assessment are vital to preventing or minimizing long-term negative effects. But many cases of this disease aren't diagnosed until the patient presents with symptoms. Conventional procedures to assess risk can be expensive, time consuming, and uncomfortable. An accessible, patient-friendly risk assessment method could help clinicians devise a regimen for their patients that might lessen the chance for serious cardiac events.

Overview

SmartScore 4.0 is designed to identify the presence of regional and global coronary artery calcification from a CT scan, then measure and score the results. Scores can be calculated using a standard Agatston/Janowitz (AJ) method. When correlated with a patient's personal information, the score can yield an estimation of a patient's risk for coronary artery disease.

Highlights

- Non-invasive alternative to conventional assessment procedures.
- Score can be correlated to age group cohort to determine patient's risk per population
- Provides information on coronary artery wall calcium plaque buildup.
- Automatically detects calcium and highlights it in green.
- Streamlines workflow by networking patient demographics from CT scanner directly to SmartScore program.
- Free Hand Trace lets you outline specific ROIs.

CONFIDENTIAL

SECRET DE AFACERI

Visit us:

<https://www.gehealthcare.com/products/advanced-visualization/all-applications/smartscore-40>

Features

- AJ 130 scoring method uses conventional Agatston/Ianowitz technique with a threshold of 130 HU which is adjusted to the appropriate image slice thickness.
- Volume scoring calculates volumes (mm^3) of calcified plaque above the 130 HU threshold.
- Mass scoring calculates mass (mg) of calcified plaque above the 130 HU threshold.
- Individual and aggregate scores are computed for each artery type.
- Provides two methods of calcium scoring.
- Customizable parameters can be predefined based on your site's preferences.
- User interface streamlines your workflow and reduces reading time.
- Report options allow you to custom tailor your reports and distribute them in a variety of formats.

System Requirements

- EKG monitor with recording device and x-ray translucent lead cable.

- Gantry hardware upgrade kit for those scanner systems already in operation.

Compatible Platforms:

- AW: Z800 (and later)
- AW Server 3.2 Ext. 3.x with LT ML 350p G8** (and later)
- Recommended monitor resolution is up to dual 2MP (1600 x 1200) or a single 3MP (1536 x 2048)

Indications for Use

SmartScore is a non-invasive software option that can be used to evaluate calcified plaques in the coronary arteries, which may be a risk factor for coronary artery disease. SmartScore may be used to monitor the progression/regression of calcium in coronary arteries over time, which may aid in the prognosis of cardiac disease.

Regulatory Compliance

This product complies with the European Council Directive 93/42/EEC Medical Device Directive as amended by European Council Directive 2007/47/EC.

This product or its features may not be available in some other countries or regions. Please contact your sales associate.


Rx Only

CONFIDENTIAL
SECRET DE AFACERI

GE imagination at work

© 2020 General Electric Company.
All rights reserved. Data subject to change.
GE and GE Monogram are trademarks of General Electric Company.
* Trademark of General Electric Company.
DOC0793786

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
GE Healthcare

CardIQ Xpress 2.0

Bringing integration and automation to your CT Cardiac post processing needs.

Cardiac disease is one of the leading health concerns worldwide. Successful treatment of the many conditions that cause and perpetuate heart disease requires that physician's approach specific cardiac problems with as much information as today's technology can provide. To supply your referring physicians with this information you need software that gives insight into cardiac anatomy and tools for thorough evaluation into the extent of disease.

Overview

CardIQ Xpress 2.0 is an integrated post processing imaging analysis application dedicated to cardiovascular imaging on GE Advantage Workstation (AW) and AW Server Enterprise System (AWS). The CardIQ Xpress 2.0 software option can be used to display, reformat, and analyze 2D or 3D cardiac CT images for qualitative or quantitative assessment of heart anatomy and coronary artery vessels from a single or multiple cardiac phase image data sets. Cardiac motion is a very real challenge that can occur at any heart rate. CardIQ Xpress 2.0 is designed to work with SnapShot Freeze* images to automatically process and display images generated with reduced motion blur artifact.

Highlights

- Automatically segment coronary tree across phases.
- Automatically tracks and labels coronary arteries.
- Improved centerline editing tools for faster edits.
- Right mouse wheel menu for quick access to renaming, deleting, and editing centerlines.
- Plaque ID tool assists in visualizing and quantifying plaque burden
- Relative perfusion highlights and quantifies hypo-dense areas of myocardium

CONFIDENTIAL
SECRET DE AFACERI

Visit us:

<https://www.gehealthcare.com/products/advanced-visualization/all-applications/cardiq-xpress-20-reveal>

Features

- Pre-processing automatically recognizes cardiac datasets and performs all segmentations as data arrives on the system saving processing time
- Auto Coronary analysis automatically segments the coronary tree, tracks and labels the coronary arteries
- Three pre-defined orientations provide quick access to planes for best visualization of the coronaries in 2D
- Quick AVA allows access to vessel tracking at any time
- Relative perfusion color codes and quantifies percent of hypo-dense areas of myocardium with four selectable color maps and fusion overlays with the coronary tree
- Plaque ID provides volume measurements for four distinct Hounsfield ranges to aid with identification and visualization of coronary plaque
- Automatically display SnapShot Freeze¹ processed images for reduced motion blur
- Robust, automatic calculation of ejection fraction and stroke volumes from the 3D endocardium volumes
- Stenosis measurements, IVUS views and preset volume rendering models assist in communication of specific findings
- Measure ES and ED for ejection fraction & volume with automatic extraction of the left ventricle

- Create multiphase beating hearts
- Select oblique reformat views in the standard cath angles for easy analysis of the coronary vessels
- Display 4D valve views with a single click

System Requirements

- AW Workstation AW 4.7 Ext 14 or higher
- AW Server 3.2 Ext. 2.0 or higher

Regulatory Compliance

This product complies with Regulation (EU) 2017/745 of the European Parliament and of the Council on medical devices (MDR). This product or its features may not be available in some other countries or regions. Please contact your sales associate.

Rx Only

Indications for Use

CardIQ Xpress 2.0 is intended to provide an optimized non-invasive application to analyze cardiovascular anatomy and pathology and aid in determining treatment paths from a set of Computed Tomography (CT) Angiographic images. CardIQ Xpress 2.0 is a CT, non-invasive, image analysis software package, which aids in diagnosing of cardiovascular disease to include, coronary artery disease, functional parameters of the heart,

heart structures and follow-up for stent placement, bypasses and plaque imaging. CardIQ Xpress 2.0 offers unique tools such as automatic tracking, which will pre-process the CT data into multiple viewing ports to allow for an expedited read time improving workflow. With CardIQ Xpress 2.0, the user can color code the myocardial tissue to show hypo/hyper-dense areas in the myocardial tissue of the heart. With the IVUS-like view the user can color code the HU units of the plaque to better visualize the difference between calcified and non-calcified plaque in the wall of the vessel and the lumen to determine the amount of atherosclerosis. The user can see the different valve planes along with a variety of new layouts to align the heart. The IVUS-like view is created by applying GE's Volume Rendering on a cross-section perpendicular to the detected centerline. This view merely displays a cross section as in IVUS imaging and color codes like IVUS images. No new or additional diagnostic information is added. CardIQ Xpress 2.0 is for use on the Advantage Workstation (AW) platform, CT Scanner, PAC or Centricity stations, which can be used in the analysis of 2D or 3D CT angiography images/data derived from DICOM 3.0 CT scans.

CONFIDENTIAL
SECRET DE AFACERI

ion at work

© 2021 General Electric Company.
All rights reserved. Data subject to change.
GE and GE Monogram are trademarks of General Electric Company.
* Trademark of General Electric Company.
DOC1036477

Features Detail

Pre-Processing

- Automated processing of single or multiphase data
- Customizable menu to choose which protocols are pre-processed
- Pre-processed protocols include: Auto Coronary Analysis, Labeling of vessels, VR heart, Angiographic View, and Auto EF

Auto Launch

- Multiple exams can be loaded into the auto launch to ready them for review (AW only)
- Color-code display of ready-to-read exam listing (AW only)
- The ability to switch between exams without having to quit out of the application and reload


Automated Coronary Vessel Analysis (CVA)

- CVA allows the users to track, extract, visualize and measure coronary arteries from either single or multiple cardiac phase images
- Auto launch case selection and display of automatically tracked and labeled coronary vessels
- Simple, corrections to vessel branch tracking
- Curved, oblique, longitudinal, and cross-sectional reformatted views are automatically generated in the software for quick review

- Current-state tracking points within vessel analysis may be saved for future review and/or manipulation
- Color-coded plaque analysis with volume measurements for evaluating the change in plaque size over time
- Generate and save rotational movies from curved planar reformation, best section and cross-sectional (lumen) views

- Quantitative or qualitative coronary vessel assessment on user-selected vessel segments.
- Vessel measurements including:
 - Distance and volume
 - Cross-sectional area and mean diameter
 - Single or Dual reference point comparison
 - Relative percent stenosis
- Single-screen filming capabilities with multi-views within the screen to show the entire picture of the vessel
- Automatically display of the coronary vessel tree
- Using pre-set protocols, 2D or 3D coronary vessel tree models can be generated and displayed in an automatic fashion to give users a qualitative overview of coronary vessel structures

- New right mouse wheel menu for quick access to renaming, deleting and editing centerlines
- Improved centerline editing tools for faster edits
- The vessel tree models can be in the form of either 3D volume rendered or Maximum Intensity Projection (MIP)

2D Review

- Three pre-defined orientation protocols for easy review of the coronary vessels and chambers
- Dual reformat review allows automatic review of the coronaries from an axial image while linking to oblique views in longitudinal and cross-sectional planes
- Reformat review predefines workflow steps to automatically present thick, multi-planar views of cardiac anatomy
- One-click access to quick vessel analysis

GE imagination at work


CONFIDENTIAL
SECRET DE AFACERI

Quick/Dynamic AVA

- Single or dual click vessel analysis from any protocol; 3D, reformat or vessel tree
- Real-time tracking of the center line with instant visualization of a vessel for quick inspection
- Ability to extend tracking proximally or distally for full view of the vessel
- Ability to select a point above and below potential lesions for automated vessel analysis

Vessel Data Base

- A vessel name data base for quick access to vessel labels
- Anatomically relevant listings
- Ability to add new vessel names into the data base
- Automated processing of single or multiphase data

Relative Perfusion

- Hypo-dense areas of myocardium can be high-lighted with color to aid in the visualization of ischemic heart disease
- Pre-defined layouts for density enhancement
- A hybrid display view to show the relationship of the vessel tree with the perfusion defect
- Quantification of a hypo-perfused area as it relates to myocardial defects
- Four selectable color maps to help display the hypo-dense areas
- Exportable statistics via one touch


IVUS Views

- Interactive volume rendered images to better define the calcium, lumen and non-calcified plaque related to lumen narrowing or wall abnormalities
- Can be applied to any 2D MIP image to include; cross-section, MPR, best L-section images

- Display cross-sections perpendicular to the vessel centerline to create an IVUS-like image. No new diagnostic information is added

PlaqID

- Customizable color mapping to HU ranges for easy identification of plaque in axial, curved reformat or MPVR images
- Four distinct colors to aid in differentiating vessel lumen, non-calcified plaque and calcified plaque
- Easy modification of the color and transparency of plaque
- Volume, area to track size of plaque over time
- Smooth transition of color from one density to another
- Plaque volume and area automatically captured in summary table.

VR Heart

- One-touch automatic segmentation algorithm designed to extract the cardiac anatomy from within the chest image
- Optimize algorithms for the easy review of heart and bypass grafts
- Multiple optimized VR curve settings to enhance different structures within the heart.
- Automatic segmentation of both single or multi-phase cardiac image data sets
- 4D beating heart images with the ability to page through the phases, rotate the image with real-time functionality

Multi-Phase Image Review

- Multiphase images can be reviewed in any protocol
- 10 phase images show the heart throughout the complete cycle
- Quickly edit phase images to keep only the phase or phases needed for analysis of the coronary vessels once multi-phase review is complete

GE imagination at work

CONFIDENTIAL
SECRET DE AFACERI

Transparency View

- 3D vessel tree overlaid on the heart chambers with the ability to adjust opacity of chambers of heart, myocardium vessels tree and bones

3D Angiographic View

- One-touch automatic visualization of the coronary arteries in an x-ray angiographic view
- Single click cardiac cath views to orientate the angiographic view into RAO/LAO, or CRANIAL/CAUDAL orientation
- Added or removed vessels are easily archived with the Auto Select feature
- Toggle between gray scale or inverse video views

One Touch Cath Views

- Default cardiac cath views are provided with CardIQ Xpress 2.0 software to provide 3D VR heart, 3D vessel tree or reformatted image in the orientation of a standard cardiac catheterization
- Single click on the angulation area to change between one-touch views
- Interactive RAO/LAO & CRANIAL/CAUDAL orientations allow for manually positioning the view
- Customizable orientations can be easily generated and saved for future use

3D Ejection Fraction

- Automatic extraction of the left ventricle across all phases and the automatic detection of systole and diastole delivers calculated ejection fraction and stroke volumes
- Robust, automatic calculation of ejection fraction and stroke volumes from the 3D endocardium volumes

Layouts

- Configurable layouts
- Dual monitor support
- Configurable default layout for use at AVA launch
- Easily capture the screen layout for identical visualization on PACs or film

4D Movie Mode

- Movie mode allows the user to view and save multiple views of the heart and display as a movie sequence
- Load multiple phases from systole to diastole to create a beating heart movie
- Save movies as a DICOM image set or export in JPEG/MPEG format

SnapShot Freeze (SSF)

- Automatically display motion blur-reduced SnapShot Freeze processed images
- Reprocess SnapShot Freeze images for artifacts caused by motion blurring after manual vessel editing
- Motion blurring correction requires the use images produced by a CT using the optional SnapShot Freeze feature

Protocols

- CardIQ Xpress 2.0 is supplied with a set of pre-defined protocols that are easily adaptable to customized protocols
- The protocols have comprehensive instructions and tools which minimize the need to memorize the procedures or to refer continuously to user documentation.

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

SECRET DE AFACERI

CardIQ Function Xpress

Real time cardiac review the instant you're ready to read.

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SECRET DE AFACERI

Contents

1. Features.....	4
2. System Requirements.....	4
3. Image Requirements.....	4
4. Indications for Use.....	4
5. Regulatory Compliance.....	4
6. Features Detail.....	4
6.01 Features.....	4
6.02 Pre-processing.....	5
6.03 Auto Launch (AW only).....	5
6.04 Automatic Chamber Detection.....	5
6.05 Automatic Ejection Fraction.....	5
6.06 Display Table.....	5
6.07 Myocardial Analysis.....	6
6.08 Simplified Workflow.....	6

CONFIDENTIAL
SECRET DE AFACERI

Cardiac disease is one of the leading health concerns worldwide. Successful treatment of the many conditions that cause and perpetuate heart disease requires that physicians approach specific cardiac problems with as much information as today's technology can provide. To supply your referring physicians with this information you need software that gives you insight into cardiac anatomy and function supported by quantifiable data.


Overview

CardIQ Function Xpress post-processing software helps you evaluate cardiac function and aids in the diagnosis of cardiovascular disease with a high degree of confidence. Providing accurate and reproducible quantification of left and right ventricular volumes, ejection fractions, and myocardial mass, CardIQ Function Xpress is optimized to assess cardiac function using multiphase, multi-slice cardiac CT images.

The CardIQ Function Xpress option provides an easy-to-use and time-effective way for you to perform cardiovascular functional analysis.

Highlights

- Automatically pre-processes and loads exams (AW only).
- Automatically detects all chambers in all phases.
- Automatically performs ejection fraction and chamber volume analysis with 91% reliability on LV and RV automatic ejection fraction.
- Provides single click myocardial analysis with bull's eye images.
- Performs simultaneous case review in real time.
- Load multiphase data with no load time.

CONFIDENTIAL
SECRET DE AFACERI

Visit us:

<https://www.gehealthcare.com/products/advanced-visualization/all-applications/cardiq-function-xpress>

Features

- Performs behind-the-scenes function data process for real time ejection fraction review, volume analysis, and myocardial analysis.
- Enables you to extract, render, and display 3D volumetric models of the endocardium for ejection fraction calculation.
- Automatic detection of the endocardial and epicardial walls for wall motion, wall thickness, wall thickening and myocardial mass analysis.
- Automatic volume analysis of all chambers.
- Automatically calculates left atrium volume while excluding the pulmonary vein.
- Activate visual wall motion with short axis images in basal, mid, and distal orientations along with a two-chamber long axis view in a single click.
- Perform myocardial analysis with wall motion, wall thickness, and mass calculations.

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Image Requirements

CardIQ Function Xpress option can be applied to gated CT images. These images can be acquired on multi-slice CT scanners. The GE HealthCare family product line must have the CardIQ SnapShot option to acquire gated images.

Indications for Use

CardIQ Function Xpress is intended to provide an optimized noninvasive application to analyze cardiovascular anatomy and pathology and aid in determining treatment paths from a set of Computed Tomography (CT) Angiographic images. CardIQ Function Xpress in conjunction with CT cardiac images to automatically calculate and display various left ventricular and right ventricular functional parameters as ejection fraction, end systolic and end diastolic volumes, stroke volumes, wall motion, wall thickening, cardiac output, myocardial mass, systemic and pulmonary vascular resistance. Volume measurement of each chamber of the heart is also available. With CardIQ Function Xpress atrium volumes may be used to determine volume assessment of atrial disease to include but not limited to atrial fibrillation. CardIQ Function Xpress is a CT, non-invasive image analysis software package, which aids in the assessment of cardiac function and in determination of cardiovascular disease diagnosis and management. CardIQ Function Xpress is for use on the Advantage Workstation (AW) platform, CT Scanner, PAC or Centricity stations, which can be used in the analysis of 2D or 3D CT angiography images/data derived from DICOM 3.0 CT scans.

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

Features Detail

Features

CardIQ Function Xpress allows you to accurately and reproducibly quantify myocardial mass, left and right ventricular volumes and ejection fractions. The software automatically detects left ventricular endocardial and epicardial contours so that you can more easily assess left ventricular (LV) and right ventricular (RV) functional parameters. The program is optimized to perform assessment of cardiac function using GE HealthCare multiphase, multi-slice cardiac CT images.

LV Segmentation

RV Segmentation

CONFIDENTIAL
SECRET DE AFACERI

Features Detail *(cont.)*

Pre-processing

- The program automatically processes multi-phase data for functional analysis.
- You can select configurations to process only the chambers you require.
- For more streamlined reading, ejection fraction, volume analysis, and myocardial analysis are completed before you view the case.

Auto Launch (AW only)

- Load multiple cases into auto launch so they are ready to review when you are.
- Switch between exams without having to quit the application and reload.
- Exams ready for review display in green.

Automatic Chamber Detection

- The program loads single or multi-phase data sets and selects the heart's chambers for volume analysis.
- Analyze individual chambers with automatic contouring of the heart's four chambers.
- The system automatically detects the left/right ventricles and atriums and contours the endocardium for ejection fraction calculation.
- Works with single or multi-phase datasets


Automatic Ejection Fraction

- The system features a protocol that, when chosen, automatically detects the end diastolic and end systolic phases and calculates an accurate ejection fraction of the left ventricle, right ventricle, and left atrium.
- Ejection fraction measurements are available for the right and left ventricles.
- Stroke volume, end diastolic and end systolic volumes are automatically calculated and shown in a graphical representation for assessment of heart failure.
- With dual screens, you can display quantification assessments on the left monitor and a qualitative overview of the full chamber from apex to base on the right monitor.
- Customize layout designs to accommodate your viewing preferences:
 - Base to apex
 - Basal, mid, and distal along with a two-chamber long axis.
- You can change the selection of end diastolic and end systolic phase.
- While contour detection is automatic, scalpel and auto select tools are available with which you can edit the contours.

Display Table

A display table shows volume and provides a graphical update of end diastolic and end systolic


- Click on available tabs to select different chambers for analysis.
- Left and right ventricle measurements include:
 - Stroke volume
 - Ejection fraction
 - Cardiac output
 - Myocardial mass
 - End systolic volume
 - End diastolic volume
 - Pulmonary vascular resistance
 - Systemic vascular resistance
 - Body surface area
 - Heart rate
- Left atrium measurements include:
 - Min and max volume
 - Ejection fraction
 - Left atrium index
 - Body surface
 - Heart rate
- Right atrium measurements include:
 - Min and max volume
 - Heart rate
- You can film, save, and export the volume graphs.

CONFIDENTIAL
SECRET DE AFACERI

Features Detail *(cont.)*

Myocardial Analysis

- CardIQ Function Xpress automatically detects the endocardial and epicardial walls to analyze wall motion, wall thickness, wall thickening, and myocardial mass analysis.
- You can edit the contours with a freehand trace and click to spline tools.
- The program displays multiple bull's eye graphs to graphically represent these cardiac characteristics:
 - Wall motion
 - Wall thickness
 - Myocardial mass
 - Cardiac output
- Choose from four bull's eye graph overlays to show 17, 20, and 49 segment model, or a coronary artery territory model.
- Manual scoring is available for wall motion and thickening in 17 and 20 segments.
- Cord display is available for visual assessment of wall abnormalities.

Simplified Workflow

CardIQ Function Xpress provides you with the tools you need to accurately analyze, evaluate, and report on cardiac anatomy and function. Its ease of use will streamline your workflow and its accuracy will give you greater diagnostic confidence.

CardIQ Function Xpress simplifies your workflow with:

- Automatic ejection fraction protocol
- Streamlined multiphase beating images for easy visualization of cardiac wall motion
- Single click myocardial analysis with bull's eye images

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
GE HealthCare

CardEP

Streamlined post-processing for enhanced electrophysiology procedures.

CONFIDENTIAL

SECRET DE AFACERI

Contents

1. Features..... 4

2. Image Requirements..... 4

3. System Requirements..... 4

4. Indications for Use..... 4

5. Regulatory Compliance..... 4

6. Features Detail..... 5

6.01 Operation..... 5

6.02 3D VR Atrium..... 5

6.03 3D VR Heart..... 5

6.04 Image Reformation..... 5

6.05 Coronary Sinus and Pulmonary Vein Analysis..... 6

6.06 Filming..... 6

6.07 One Touch EP Views..... 6

6.08 Pulmonary Vein (PV) and Left Atrial Appendage Reformation..... 6

6.09 Phase Registration..... 6

6.10 EP Navigator Views..... 6

6.11 Batch Movie Mode..... 6

6.12 Protocols..... 7

6.13 Summary..... 7

CONFIDENTIAL

SECRET DE AFACERI


CT cardiac exams provide critical information to practitioners of cardiac therapy. Precise information pertaining to the left atrium's complex anatomy, the pulmonary veins, the coronary sinus, or the cardiac veins has a major impact on the efficacy of subsequent cardiac therapy: It can speed procedures and facilitate treatment. But analyzing and reporting the wealth of information CT cardiac studies provides can be time consuming. A program that automates many routine functions and gives you tools to easily quantify and qualify many aspects of cardiac function would streamline your workflow and give you greater diagnostic confidence.

Overview

CardEP is an integrated post processing image analysis software for the application of cardiovascular and electrophysiology imaging on the AW Workstation. With it, you can process display, reformat, and analyze 2D or 3D cardiac images for qualitative or quantitative assessment of heart anatomy and pulmonary veins from single or multi-phase cardiac image data sets.

Highlights

- Visualize the anatomical layout of the left atrium and pulmonary vessels.
- Visualize the origins of the pulmonary vessels.
- Automatic heart segmentation.
- Helps you visualize and quantify cardiac venous pathways.
- Offers you a variety of 2D, 3D, or reformatted protocols with which to perform image analysis.

CONFIDENTIAL

SECRET DE AFACERI

Visit us:

<https://www.gehealthcare.com/products/advanced-visualization/all-applications/cardep>

Features

CardEP offers a variety of features that give you image analysis flexibility:

- Volume render 3D heart and chamber views
- Obtain segmented 3D views of the left atrium and pulmonary veins
- View hollow 3D VR heart models
- Measure the pulmonary veins and left atrial appendage with predefined protocols
- Render and display 2D/3D views of the left chambers in single or multiple phases of the cardiac cycle with one-touch automatic extraction
- Reformat standard axial CT images of single or multiple cardiac phases into short or long axis and save the image series for further analysis.
- Automatically track, extract, and display pulmonary veins with the pulmonary vessel analysis tool
- Generate pulmonary vein measurements
- Automatically track, extract, and display coronary sinus and cardiac veins with the coronary sinus analysis tool
- Perform one-touch standard 2D reformation and 3D volume rendering, including hollow view, of cardiac anatomy with optimized AW Volume Viewer and 3D rendering tools
- Register images from different cardiac phases into a unique data set that you can save as a 3D object or use for further analysis
- Use predefined navigator views of the pulmonary veins to get a fly-through perspective of the vessels.

Image Requirements

CardEP accepts standard cardiac gated CT image sets acquired on CT scanners equipped with DICOM 3.0 standards. Images must meet the same image requirements as those for the basic Volume Viewer application.

System Requirements

Refer to AW Applications platform compatibility document.

Indications for Use

CardEP is a post-processing software option for the Advantage Workstation (AW) Platform. This product can be used for the analysis of CT angiographic images for the assessment of the heart to include the atria, pulmonary veins and coronary sinus. It provides quantitative analysis tools which include a number of display, measurement and model export capabilities. This product can be used to aid trained physicians in the visualization and assessment of cardiac anatomy.

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Features Detail

Operation


- On the menu page, multiple image processing tools for both single and multiple cardiac phase studies are provided. Access each pre-set protocol from the AW Volume Viewer application or from the application pull-down menu page on the browser.

- CardEP images present you with a detailed display of patient heart rate, cardiac scanning parameters, and each scan's temporal resolution.

3D VR Atrium

- Pre-set protocols enable you to semi-automatically generate and display 2D or 3D left atrial/ventricular views for a qualitative overview of the heart's chambers.
- The atrial/ventricular model can be either 3D volume rendered or Maximum Intensity Projection (MIP).
- With the endocardial view feature, you can view inside the atrium to analyze the relationship of the pulmonary vein ostia from the chamber's endocardial surface.
- With a one-touch protocol you can quickly change from a 3D model to an endoscopic view of the atrium and pulmonary veins.
- Anatomical landmark selection lets you deposit points on the 3D model, embed them, and then export the model for future registration use with a mapping system.
- Save predefined landmarks and models and export them to CD.
- Save current model state to the AW browser and reload for future analysis.
- Optimized VR and curve setting provide clearer contrast enhanced cardiac display.
- Choose from single or multi cardiac phase sets.
- Add/remove structure (Autoselect), scalping and brush threshold tools are all available for quick optimization of a 3D model.


3D VR Heart

- Automatically extracts cardiac anatomy from within the chest cavity.
- Optimized VR and curve setting provides clearer contrast enhanced cardiac display.
- Automatic segmentation for single or multi-phase cardiac data.
- The program supports multiphase thresholding, filtering, and scalping techniques.
- Create movies of rotating and beating hearts with easy export capabilities.
- Export enhanced 3D models, including anatomical markers.

Image Reformation

- With batch processing tools, you can reformat either single or multi-phase cardiac images in various non-axial orientations, e.g. short axis, long axis, chamber views. Save reformatted images in separate image series for further analysis (e.g. from PACs).


- Offers predefined protocols to generate automatically short and long axis views.
- Multi-phase image capability lets you quickly and efficiently assess phases or function.
- View the cross-sectional diameter of the pulmonary vein and left atrial appendage with specifically designed protocols.
- Measurement tools integrated with protocols enable you to quickly ascertain length and cross-sectional diameter measurements.

CONFIDENTIAL
SECRET DE AFACERI

Features Detail *(cont.)*

Coronary Sinus and Pulmonary Vein Analysis

- Load single or multi-phase cardiac data sets in the vessel analysis tool.

- With the semi-automated pulmonary vein/coronary sinus analysis tool (CSA), based on the Advanced Vessel Analysis (AVA) technique, you can track, extract, visualize, and measure the veins and sinus (physical dimension and stenosis sizing) from either single or multiple cardiac phase image data sets.
- Generate curved, oblique, longitudinal, and cross-sectional reformatted views of individual vessels.
- Generate 3D models of the atrium and pulmonary veins during your vessel analysis procedure.
- Save current state tracking points within your vessel analysis for future image review or manipulation.
- Generate and save rotational movies from curved planar reformation and cross-sectional (lumen) views.
- Generate and save movies from Best L Section reformatted images.
- Perform quantitative and qualitative coronary sinus assessments on segment you choose:
 - Distance measurement
 - Cross-sectional area
 - Single or dual reference point comparison
 - Relative percentage stenosis
 - Volume measurement
 - Mean diameter
- Generate beating, rotating movies for visualization and communicating findings.

Filming

CardEP offers you filming flexibility. Batch filming capability is provided, and any image can be queued for filming.

One Touch EP Views

- Default views provide you with visualization of the atria and pulmonary veins in correct anatomical orientation. You are provided the capability to create your own presets to save and apply to future analyses.
- Use the interactive RAO/LAO and CRA/CAU annotation feature to manually orient the view.

Pulmonary Vein (PV) and Left Atrial Appendage Reformation

- Predefined protocols allow you to generate reformation of the pulmonary veins and left atrial appendage in the long axis.
- Generate cross-sectional views of the veins or appendage to determine the ostium's diameter.

Phase Registration

- Suppress artifacts that may occur when the heart rate fluctuates during a scan with a semi-automated multi-phase tool that aligns cardiac chambers, veins, or sinus structures.
- You can automatically or manually register phases among images to improve image quality.
- Save the registered volume data sets for future review or manipulation.

EP Navigator Views

- With one-touch, the navigator displays endoscopic views of the left atrium.
- View the pulmonary vein and atrial appendage ostia from either the left or right side using a predefined orientation layout preset.
- Optimize interior atrium views with automated threshold settings.
- Global virtual endoscopic view provides you ostial displays of the right and left pulmonary veins and the left atrial appendage.

Batch Movie Mode

- Batch movie mode lets you define first and last images for movie creation.
- Create rotating and beating movies
- Save automatically generated movies to the viewer as DICOM image sets or send them to data export in JPEG/MPEG/AVI formats.
- Use the movie mode with MIP, 3D VR, and reformatted images.
- Generate and save movies from Best L Section Reformatted images, curved planar reformatted, cross-sectional, and lumen views.

CONFIDENTIAL
SECRET DE AFACERI

Features Detail *(cont.)*

Protocols

CardEP contains predefined protocols that you can edit or modify to adapt them to your specific customized protocols and save under a new name.

The CardEP protocols comprise comprehensive instructions and the tools you need to perform analyses. Tool icons are displayed in the protocol panels for fast and easy access.

Summary

CardEP gives you the tools and flexibility to quickly analyze cardiac studies. The information and statistics you can generate with this program may help referring physicians with their assessment of therapy options for their patients.

CONFIDENTIAL

SECRET DE AFACERI

About GE HealthCare

GE HealthCare is a leading global medical technology, pharmaceutical diagnostics, and digital solutions innovator, dedicated to providing integrated solutions, services, and data analytics to make hospitals more efficient, clinicians more effective, therapies more precise, and patients healthier and happier. Serving patients and providers for more than 100 years, GE HealthCare is advancing personalized, connected, and compassionate care, while simplifying the patient's journey across the care pathway. Together our Imaging, Ultrasound, Patient Care Solutions, and Pharmaceutical Diagnostics businesses help improve patient care from prevention and screening, to diagnosis, treatment, therapy, and monitoring. We are an \$18 billion business with 51,000 employees working to create a world where healthcare has no limits.

Follow us on [Facebook](#), [LinkedIn](#), [Twitter](#), [Instagram](#) and [Insights](#) for the latest news, or visit our website [gehealthcare.com](https://www.gehealthcare.com) for more information.

CONFIDENTIAL
SECRET DE AFACERI

Products mentioned in the material may be subject to government regulations and may not be available in all countries. Shipment and effective sale can only occur after approval from the regulator. Please check with local GE HealthCare representative for details.

©2024 GE HealthCare. GE is a trademark of General Electric Company used under trademark license.

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