# TERUMOBCT

## **40851, 40852, 40853, 40854, 40855, 40856**

Terumo - Platelet Additive Solution + (T-PAS+)

Instructions for Use

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## Description

---

The solution is presented as a single use container with a 200 mL, 220 mL, 250 mL, 280 mL, 300 mL, or 500 mL volume of T-PAS+ in a clear polyolefin bag covered by a laminated polypropylene overwrap. The solution has a pH of 7.2.

Each 1,000 mL contains:

- Natrii Citras dihydricus Ph Eur - 3.18 g
- Natrii acetat trihydricus Ph Eur - 4.42 g
- Natrii dihydrogenophosphas dihydricus Ph Eur - 1.05 g
- Dinatrii phosphas dodecahydricus Ph Eur - 7.69 g
- Kalii chloridum Ph Eur - 0.37 g
- Magnesii chloridum hexahydricum Ph Eur - 0.30 g
- Natrii chloridum Ph Eur - 4.05 g
- Aqua ad Injectibilia Ph Eur

## Storage Recommendations

---

Store between 0 °C and 25 °C.

The expiration date of the solution is printed on the bag and the product should not be used after this date.

## Indications

---

The Terumo - Platelet Additive Solution + (T-PAS+) is a platelet additive solution intended to partially replace plasma in the preparation and storage of a buffy coat-derived platelet concentrate or apheresis platelet units.

Please see the device operator's manual for a complete listing of warnings, cautions, and undesirable effects.

## Warnings

---

1. **Do Not Reuse/Not for Reuse:** Terumo BCT, Inc. products bearing the "Do Not Reuse" symbol are intended for single use only and are not intended to be reused or re-sterilized in any manner. Terumo BCT cannot ensure the functionality or sterility of the product if it is reused or re-sterilized.

Reuse of a single-use product could result in:

- Product performance issues due to a loss of product integrity, including but not limited to the following:
  - Fluid leaks
  - Parts that are warped or deformed
  - Plastics that are brittle and discolored
  - Filters that have reduced filtration capabilities
- Viral infections such as hepatitis or human immunodeficiency virus (HIV)
- Bacterial infections
- Cross-contamination

Any of these risks could result in serious injury or death. These risks are shared by product users, donors, patients, and recipients of end products of the device.

## Cautions

---

1. This product is sterilized with steam.
2. When you break frangible connectors, bend them in both directions to ensure that you break them completely. Failure to do so may result in restricted flow.

## Adverse Reactions

---

T-PAS+ is added to platelet concentrates after the collection is complete. It is not for direct intravenous infusion. T-PAS+ is not expected to cause adverse events other than those normally associated with platelet transfusion.

## Procedures for Use

---

- Not for direct intravenous infusion.
- Do not use unless the solution is clear and the individual package is intact.
- Do not remove the overwrap from the bag until immediately before use.
- Discard unused portion.
- Discard the empty bag that contained the solution. That bag is not for containing blood or blood components.
- Do not add any drug into the bag.

## Symbols and Certifications

---

| Symbol                                                                              | Definition                                                                                                                  |
|-------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------|
|   | Indicates the temperature limitations for product storage when accompanied by temperature(s).                               |
|  | Indicates that the product is sterilized using steam or dry heat.                                                           |
|  | Indicates that the product should not be reused.                                                                            |
|  | Indicates the product was manufactured in accordance with Annex II of the European Council Directive 93/42/EEC, as amended. |
|  | Indicates the product lot number when accompanied by a number.                                                              |
|  | Indicates the expiration date of the product when accompanied by a specific date.                                           |
|  | Indicates that the user should consult the instructions for use.                                                            |

| Symbol                                                                            | Definition                                                         |
|-----------------------------------------------------------------------------------|--------------------------------------------------------------------|
|  | Indicates the medical device manufacturer.                         |
|  | Indicates the Authorized representative in the European Community. |
|  | Indicates a medical device that is non-pyrogenic.                  |

## Return of Used Product

---

If for any reason this product must be returned to Terumo BCT, Inc., a returned goods authorization (an RGA number) is required from Terumo BCT prior to shipping.

Instructions for cleaning and materials, including appropriate shipping containers, proper labeling, and an RGA number, may be obtained from the Terumo BCT Quality Assurance Department.

**IT IS THE RESPONSIBILITY OF THE HEALTH CARE INSTITUTION TO ADEQUATELY PREPARE AND IDENTIFY THE PRODUCT FOR RETURN SHIPMENT.**

**Please contact your local representative for information regarding returned goods and product complaints.**

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©2019 Terumo BCT, Inc.