
# Mini-Spike®

---

---

**B | BRAUN**



<b>en</b>	Removal and injection spikes for single- and multi-dose vials.	<b>4</b>	<b>no</b>	Aspirasjons- og injeksjonspigger for enkelt- og flerdosebeholdere.	<b>28</b>
<b>de</b>	Entnahme- und Zuspitzspike für Einzel- und Mehrdosenfläschchen.	<b>5</b>	<b>pl</b>	Szpunkulce do pobierania i wstrzykiwania do fiolek jedno- i wielodawkowych.	<b>29</b>
<b>bg</b>	Шипове за изтегляне и инжектиране за едно- и многодозови флакони.	<b>6</b>	<b>pt</b>	Remoção e spikes de injeção para frascos de dose única ou multi-dose.	<b>31</b>
<b>cs</b>	Odběrové nebo injekční hroty pro jednodávkové a vícedávkové injekční lahvičky.	<b>8</b>	<b>pt</b> <b>BR</b>	Spikes de remoção e de injeção para frascos de dose única e multidoses.	<b>32</b>
<b>da</b>	Fjernelses- og injektionsspyd til hætteglas med en eller flere doser.	<b>9</b>	<b>ro</b>	Vârfuri pentru îndepărtare și injectare destinate flacoanelor unidoză și multidoză.	<b>34</b>
<b>el</b>	Ακίδες αφαίρεσης και έγχυσης για φιαλίδια μίας ή πολλών δόσεων.	<b>11</b>	<b>ru</b>	Наконечники для аспирации и инъекций, предназначенные для однодозных и мультidosных флаконов.	<b>35</b>
<b>es</b>	Punzones de extracción y de inyección para viales de dosis única o múltiple.	<b>12</b>	<b>sk</b>	Aspiračné a injekčné hroty pre jednodávkové a viacdávkové liekovky.	<b>37</b>
<b>et</b>	Ühe- ja mitmeannuseliste viaalide eemaldus- ja süsteteravikud.	<b>14</b>	<b>sl</b>	Odstranjevanje in konice za injiciranje za enodmerne in večodmerne viale.	<b>38</b>
<b>fi</b>	Poisto- ja injektointipiikit yksi- ja moniannosampulleille.	<b>15</b>	<b>sr</b>	Šiljci za ubrizgavanje ili izvlačenje tečnosti iz jednodoznih i višedoznih bočica.	<b>39</b>
<b>fr</b>	Perforateurs de retrait et d'injection pour des flacons mono- et multi-doses.	<b>16</b>	<b>sv</b>	Uttags- och injiceringspetsar för en- och flerdosflaskor.	<b>41</b>
<b>hu</b>	Tűskék felszíváshoz/befecskendezéshez egy- és multidózisos ampullákhoz.	<b>18</b>	<b>tr</b>	Tek ve çok dozlu flakonlar için sıvı çekme ve enjeksiyon spike'i.	<b>42</b>
<b>id</b>	Taji atau paku pelepasan dan injeksi untuk vial satu dosis dan multi-dosis.	<b>19</b>	<b>uk</b>	Аспіраційні канюлі для одноразової та багаторазової аспірації та ін'єкції у флакони.	<b>43</b>
<b>it</b>	Dispositivo di prelievo e iniezione per fiale monodose e multidose.	<b>21</b>	<b>vi</b>	Kim hút và bơm cho các loại ống thuốc một liều hoặc đa liều.	<b>45</b>
<b>kk</b>	Бір және көп дозалы құтыларға арналған алу және енгізу инелері.	<b>22</b>	<b>zh</b>	用于单剂量和多剂量小药瓶的抽吸和注入穿刺器。	<b>47</b>
<b>lt</b>	Nuėmimo ir injekcijos smaigai vienos dozės ar daugiadoziams flakonams.	<b>24</b>			
<b>lv</b>	Izņemšanas vai injekcijas adatas vienas un vairāku devu pudelītēm.	<b>25</b>			
<b>nl</b>	Verwijderings- en injectiespikes voor vials die een enkele of meerdere doses bevatten.	<b>26</b>			

**en** Directions  
**de** Anwendungshinweise  
**bg** Указания за употреба  
**cs** Pokyny pro použití  
**da** Betjeningsvejledning  
**el** Υποδείξεις εφαρμογής  
**es** Aplicaciones  
**et** Kasutusjuhendid

**fi** Käyttöohje  
**fr** Conseils d'utilisation  
**hu** Használati útmutatás  
**id** Petunjuk  
**it** Modalità d'uso  
**kk** Бағыттар  
**lt** Įspėjimai  
**lv** Norādījumi par lietošanu

**nl** Gebruiksaanwijzingen  
**no** Bruk  
**pl** Przepisy stosowania  
**pt** Instruções de uso  
**pt BR** Instruções de uso  
**ro** Instrucţiuni de folosire  
**ru** Инструкция по применению  
**sk** Aplikáčné pokyny

**sl** Napotki za uporabo  
**sr** Uputstvo  
**sv** Användningsanvisningar  
**tr** Kullanım talimatları  
**uk** Вказівки  
**vi** Hướng dẫn  
**zh** 说明

## **en** Instructions for use

### **Device description**

#### **Material used:**

Housing:	ABS/SAN
Closing Cap:	PP
Protective Cap and Valve Opener:	PE
Filter:	Acrylic copolymer / Nylon support
Valve Disk:	silicone

Mini-Spike® can be used with syringes which are in compliance with ISO 80369-7.

Mini-Spike® can be used with single- and multi-dose vials and fluid containers of different volumes.

Removal and injection spikes for single- and multi-dose vials.

### **Sterility**

Sterilized by Gamma radiation.

Do not use if package is damaged.

Visually inspect the product packaging to ensure that the sterile barrier system is intact.

### **Intended purpose**

Mini-Spike® is a vented dispensing pin for injecting/withdrawing fluids from single and multi-dose containers (i.e. glass vials, semi rigid and glass IV container).

Do not re-use.

### **Indication**

- Withdrawal and injection from/in multi-dose vials
- Withdrawal and injection from/in multi-dose vials with additional filtration of fluids (Mini-Spike® Filter and some Mini-Spike® Chemo variants only)
- Preparation of ready-for-use injections (according to the SmPC of the drug)
- Preparation of cytostatic injections (according to the SmPC of the drug)

#### **A: Venting filter 0.45 µm**

Protects the medication from contamination by bacteria.

#### **B: Venting filter 0.2 µm**

Protects the medication against contamination by bacteria and the user from contamination by aerosols.  
Specially suited for use with toxic substances.

#### **C: Micro-Tip**

Reduces the residual volume.

Specially suited to small volume, expensive medications.

#### **D: Liquid filter/particle filter 5 µm**

Protects against contamination by particles.

Specially suited for use with medications in powder form.

#### **E: Valve**

The valve also closes upside down containers and is drip proof.

### **Patient population**

Mini-Spike® is not directly used on patients as it is only used for drug admixture. No gender or age related limitations.

### **Intended User**

Mini-Spike® should only be used by healthcare professionals (e.g., professional nurses, licensed and/or certified physicians, healthcare workers, midwives, paramedics, pharmacists) who have been adequately trained within the educational background in this technique.

Upon medical assessment and adequate instruction, patients and/or caregivers may also be allowed to take over defined handling steps according to national guidelines. Healthcare workers are responsible for ensuring that patients and/or caregivers are instructed in the correct use of Mini-Spike® and that the instruction of caregiver and/or patient is documented by the healthcare worker according to local documentation guidelines.

### **Contraindications**

Withdrawal of cellular containing blood samples.

Reconstitution of solid pharmaceuticals with additional solvent. There are special transfer devices available for this indication.

### **Residual Risks/Side effects**

No side effects known.

The following complications may occur during the procedure of liquid transfer, reconstitution or dilution:

- Microbiological or particulate contamination as a result of poor aseptic technique up to infection
- Medication errors

### **Precautions**

Do not re-sterilize.

Always refer to the SmPCs of the drugs / solutions for admixture procedure.

### **Warning**

Do not re-use. Re-use of single-use devices creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability. Contamination and/or limited functionality of the device may lead to injury, illness, or death of the patient.

Every dispensing pin is intended to be used with a single container and is not made for re-use.

Do not use if packaging or product is damaged or contaminated. Use of damaged or contaminated device might lead to injury, illness, or death of the patient.

The duration of the application is applied in accordance with the used pharmaceutical and its SmPCs as well as in accordance with the national guidelines and/or hospital protocols.

If product is used with cytotoxic or antineoplastic drugs, please mind safe handling instructions and guidelines / institutional protocols.

Do not remove the dispensing pin/spike from the stopper / vial. The injection pin may cause injuries with careless handling.

Solution from a filled syringe into the container must only be injected upright.

Vials with inserted Mini-Spike® are to be stored in a vertical position (upright).

In case of repeated withdrawal of fluids from a multi-dose vial: disinfect the female luer connection / swabable valve as per institutional protocol.

For variants with Liquid filter / Particle filter: When using the following pharmaceuticals, the absorption is larger than 15% due to the liquid filter.

The resulting loss has to be taken into account:

- Corticotropin 30 %
- Cosyntropin 35 %
- Cyancobalamin 70 %
- Ergonovine Maleate 45 %
- Glucagon 60 %
- Insulin 90 %
- Vasopressin 95 %

Every withdrawal device is intended to be used with a single container and is not made for reuse, because even by applying the device according to its intended use there might occur fissuring or stress cracking due to mechanical and/or chemical influences.

Substances causing stress cracks:

Various inorganic and organic surface-active substances, solvents (e. g. macrogol/polyethylene glycol and/or benzyl alcohol) in medication or active substances (e.g. carboplatin, daunorubicin, fluorouracil, cisplatin, paclitaxel) can lead to stress and tension cracks in plastics. With these aggressive substances we recommend use of the spike for a single aspiration only.

If hygiene rules are disregarded, patient-cross-contamination is potentially possible, in spite of the use of the product.

### Operating Instructions

- 1) Remove protective cap from the medication container and disinfect the stopper.
- 2) Remove protective cap from the spike.
- 3) Insert spike.
- 4) Open snap cap.
- 5) Connect syringe and aspirate / inject required volume.
- 6) **Important:** Solution from a filled syringe into the container (Fig. 6) must only be injected right-side up. Risk of excess pressure and impaired drawing.

### Duration of use

Please refer to the SmPC of the medication as well as valid hygiene and safety regulations.

### Disposal

Disposal acc. to local guidelines and/or clinical protocols.

### Storage and handling conditions

Keep away from sunlight.

Keep dry.

### Notice to the user

If, during the use of this product or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorised representative and to your national authority.

If further Instructions for Use are required, they can be requested from the manufacturer or can be obtained from the B. Braun homepage: <https://eifu.bbraun.com/>.

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**AU** B. Braun Australia Pty. Ltd.  
Level 5, 7 – 9 Irvine Place  
Bella Vista NSW 2153  
Australia  
For sterilization type see primary  
packaging

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**CA** Distributed by: **B. Braun Medical Inc.**  
Bethlehem, PA USA 18018-3524

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**ID** Distributed by:  
PT. B. Braun Medical Indonesia  
Jakarta, Indonesia

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**US** Rx only  
Pictorials are for reference only  
Distributed by: **B. Braun Medical Inc.**  
Bethlehem, PA USA 18018-3524

**CH REP**  
B. Braun Medical AG  
Seesatz 17  
6204 Sempach  
Switzerland

Manufacturing site:  
B. Braun Melsungen AG  
Carl-Braun-Str. 1  
34212 Melsungen  
Germany

Made in Germany

B. Braun Melsungen AG  
34209 Melsungen  
Germany  
www.bbraun.com

0323 12410454