# CardIQ Xpress 2.0 Reveal

Bringing integration and automation to your CT Cardiac post processing needs.

Cardiac disease is one of the leading health concerns worldwide. Successful treatment of the many conditions that cause and perpetuate heart disease requires that physicians approach specific cardiac problems with as much information as today's technology can provide. To supply your referring physicians with this information you need software that gives you insight into cardiac anatomy and tools for thorough evaluation into the extent of disease.

## Overview

CardIQ Xpress 2.0 Reveal is an integrated post processing imaging analysis application dedicated to cardiovascular imaging on GE Advantage Workstation (AW) and AW Server (AWS). The CardIQ Xpress 2.0 Reveal software option can be used to display, reformat and analyze 2D or 3D cardiac CT images for qualitative or quantitative assessment of heart anatomy and coronary artery vessels from a single or multiple cardiac phase image data sets. Cardiac motion is a very real challenge that can occur at any heart rate. CardIQ Xpress 2.0 Reveal is designed to work with SnapShot Freeze\* images to automatically process and display images generated with reduced motion blur artifact.

## Highlights

- Automatically segment coronary tree across phases.
- Automatically tracks and labels coronary arteries.
- Improved centerline editing tools for faster edits.
- New right mouse wheel menu for quick access to renaming, deleting and editing centerlines.
- Plaque ID tool assists in visualizing and quantifying plaque burden
- Relative perfusion highlights and quantifies hypo-dense areas of myocardium

## Features

- Pre-processing automatically recognizes cardiac datasets and performs all segmentations as data arrives on the system saving processing time.
- Auto Coronary analysis automatically segments the coronary tree, tracks the centerline and labels the coronary arteries.
- For readers who prefer a 2D read, three pre-defined orientations provide quick access to planes for best visualization of the coronaries.
- Quick AVA allows access to vessel tracking at any time.
- Relative perfusion color codes and quantifies percent of hypo-dense areas of myocardium with four selectable color maps and fusion overlays with the coronary tree.
- Plaque ID provides volume measurements for four distinct Hounsfield ranges to aid with identification and visualization of coronary plaque in axial, curved or MPVR views.
- Automatically display SnapShot Freeze processed images for reduced motion blur<sup>1</sup>.
- Robust, automatic calculation of ejection fraction and stroke volumes from the 3D endocardium volumes.
- Portfolio of tools such as stenosis measurements, IVUS views and preset volume rendering models assist in communication of specific findings back to referring physicians.
- Measure ES and ED for ejection fraction & volume with automatic extraction of the left ventricle.

- Create multiphase beating heart images
- Select oblique reformat views in the standard cath angles for easy analysis of the coronary vessels
- Display 4D aortic valve and mitral valve views with a single click

## Indications for Use

CardIQ Xpress 2.0 is intended to provide an optimized non-invasive application to analyze cardiovascular anatomy and pathology and aid in determining treatment paths from a set of Computed Tomography (CT) Angiographic images. CardIQ Xpress 2.0 is a CT, non-invasive, image analysis software package, which aids in diagnosing of cardiovascular disease to include, coronary artery disease, functional parameters of the heart, heart structures and follow-up for stent placement, bypasses and plaque imaging. CardIQ Xpress 2.0 offers unique tools such as automatic tracking, which will pre-process the CT data into multiple viewing ports to allow for an expedited read time improving workflow. With CardIQ Xpress 2.0, the user can color code the myocardial tissue to show hypo/hyper-dense areas in the myocardial tissue of the heart. With the IVUS-like view the user can color code the HU units of the plaque to better visualize the difference between calcified and non-calcified plaque in the wall of the vessel and the lumen to determine the amount of atherosclerosis. The user can see the different valve

planes along with a variety of new layouts to align the heart. The IVUS-like view is created by applying GE's Volume Rendering on a cross-section perpendicular to the detected centerline. This view merely displays a cross section as in IVUS imaging and color codes like IVUS images. No new or additional diagnostic information is added. CardIQ Xpress 2.0 is for use on the Advantage Workstation (AW) platform, CT Scanner, PAC or Centricity stations, which can be used in the analysis of 2D or 3D CT angiography images/data derived from DICOM 3.0 CT scans.

## System Requirements

CardIQ Xpress 2.0 Reveal option can be installed on GE's AW workstation or AW Server with VolumeShare 5, VolumeShare 7 or later software.

### Compatible hardware VolumeShare 5:

- 8400, 8600 workstation 16GB RAM
- Z800 24GB RAM
- AW Server 2.0

### Compatible hardware VolumeShare 7:

- Z800, Z820 (and later) 24GB of RAM
- AW Server 3.1 and above (recommended monitor resolution is up to dual 2MP (1600 x 1200) or a single 3MP (1536 x 2048))
- Centricity™ Universal Viewer

## Regulatory Compliance

This product complies with the following requirements:  
European Council Directive 93/42/EEC concerning medical devices.

<sup>1</sup> SnapShot Freeze option is not compatible with Centricity Universal Viewer.

GE imagination at work

**CONFIDENTIAL**  
**SECRET DE AFACERI**

## Feature Details

### Pre-Processing

- Automated processing of single or multiphase data
- Customizable menu to choose which protocols are pre-processed
- Pre-processed protocols include: Auto Coronary Analysis, Labeling of vessels, VR heart, Angiographic View, and Auto EF

### Auto Launch

- Multiple exams can be loaded into the auto launch to ready them for review (AW only)
- Color-code display of ready-to-read exam listing (AW only)
- The ability to switch between exams without having to quit out of the application and reload


### Automated Coronary Vessel Analysis (CVA)

- CVA allows the users to track, extract, visualize and measure coronary arteries from either single or multiple cardiac phase images
- Auto launch case selection and display of automatically tracked and labeled coronary vessels
- Simple, corrections to vessel branch tracking

- Curved, oblique, longitudinal and cross-sectional reformatted views are automatically generated in the software for quick review
- Current-state tracking points within vessel analysis may be saved for future review and/or manipulation
- Color-coded plaque analysis with volume measurements for evaluating the change in plaque size over time
- Generate and save rotational movies from curved planar reformation, best section and cross-sectional (lumen) views
- Quantitative or qualitative coronary vessel assessment on user-selected vessel segments.
- Vessel measurements including:
  - Distance and volume
  - Cross-sectional area and mean diameter

- Single or Dual reference point comparison
- Relative percent stenosis
- Single-screen filming capabilities with multi-views within the screen to show the entire picture of the vessel
- Automatically display of the coronary vessel tree
- Using pre-set protocols, 2D or 3D coronary vessel tree models can be generate and displayed in an automatic fashion to give users a qualitative overview of coronary vessel structures

- New right mouse wheel menu for quick access to renaming, deleting and editing centerlines
- Improved centerline editing tools for faster edits
- The vessel tree models can be in the form of either 3D volume rendered or Maximum Intensity Projection (MIP)

CONFIDENTIAL

SECRET DE AFACERI

GE imagination at work

## 2D Review


- Three pre-defined orientation protocols for easy review of the coronary vessels and chambers
- Dual reformat review allows automatic review of the coronaries from an axial image while linking to oblique views in longitudinal and cross-sectional planes
- Reformat review predefines workflow steps to automatically present thick, multi-planer views of cardiac anatomy
- One-click access to quick vessel analysis

### Quick/Dynamic AVA

- Single or dual click vessel analysis from any protocol; 3D, reformat or vessel tree
- Real-time tracking of the center line with instant visualization of a vessel for quick inspection
- Ability to extend tracking proximally or distally for full view of the vessel
- Ability to select a point above and below potential lesions for automated vessel analysis

### Vessel Data Base

- A vessel name data base for quick access to vessel labels
- Anatomically relevant listings
- Ability to add new vessel names into the data base
- Automated processing of single or multiphase data

### Relative Perfusion

- Hypo-dense areas of myocardium can be high-lighted with color to aid in the visualization of ischemic heart disease
- Pre-defined layouts for density enhancement
- A hybrid display view to show the relationship of the vessel tree with the perfusion defect
- Quantification of a hypo-perfused area as it relates to myocardial defects
- Four selectable color maps to help display the hypo-dense areas
- Exportable statistics via one touch

### IVUS Views

- Interactive volume rendered images to better define the calcium, lumen and non-calcified plaque related to lumen narrowing or wall abnormalities

- Can be applied to any 2D MIP image to include; cross-section, MPR, best L-section images
- Display cross-sections perpendicular to the vessel centerline to create an IVUS-like image. No new diagnostic information is added

### PlaqID

- Customizable color mapping to HU ranges for easy identification of plaque in axial, curved reformat or MPVR images
- Four distinct colors to aid in differentiating vessel lumen, non-calcified plaque and calcified plaque
- Easy modification of the color and transparency of plaque
- Volume, area to track size of plaque over time
- Smooth transition of color from one density to another
- Plaque volume and area automatically captured in summary table.

GE imagination at work

SECRET DE AFACERI  
CONFIDENTIAL

### VR Heart

- One-touch automatic segmentation algorithm designed to extract the cardiac anatomy from within the chest image
- Optimize algorithms for the easy review of heart and bypass grafts
- Multiple optimized VR curve settings to enhance different structures within the heart.
- Automatic segmentation of both single or multi-phase cardiac image data sets
- 4D beating heart images with the ability to page through the phases, rotate the image with real-time functionality

### Multi-Phase Image Review

- Multiphase images can be reviewed in any protocol
- 10 phase images show the heart throughout the complete cycle
- Quickly edit phase images to keep only the phase or phases needed for analysis of the coronary vessels once multi-phase review is complete

### Transparency View

- 3D vessel tree overlaid on the heart chambers with the ability to adjust opacity of chambers of heart, myocardium vessels tree and bones

### 3D Angiographic View

- One-touch automatic visualization of the coronary arteries in an x-ray angiographic view
- Single click cardiac cath views to orientate the angiographic view into RAO/LAO, or CRANIAL/CAUDAL orientation
- Added or removed vessels are easily archived with the Auto Select feature
- Toggle between gray scale or inverse video views

### One Touch Cath Views

- Default cardiac cath views are provided with CardIQ Xpress 2.0 Reveal software to provide 3D VR heart, 3D vessel tree or reformatted image in the orientation of a standard cardiac catheterization
- Single click on the angulation area to change between one-touch views

- Interactive RAO/LAO & CRANIAL/CAUDAL orientations allow for manually positioning the view
- Customizable orientations can be easily generated and saved for future use

### 3D Ejection Fraction

- Automatic extraction of the left ventricle across all phases and the automatic detection of systole and diastole delivers calculated ejection fraction and stroke volumes
- Robust, automatic calculation of ejection fraction and stroke volumes from the 3D endocardium volumes

### Layouts

- Configurable layouts
- Dual monitor support
- Configurable default layout for use at AVA launch
- Easily capture the screen layout for identical visualization on PACs or film

### 4D Movie Mode

- Movie mode allows the user to view and save multiple views of the heart and display as a movie sequence

CONFIDENTIAL

SECRET DE AFACERI

GE imagination at work

- Load multiple phases from systole to diastole to create a beating heart movie
- Save movies as a DICOM image set or export in JPEG/MPEG format

#### SnapShot Freeze (SSF)

- Automatically display motion blur-reduced SnapShot Freeze processed images
- Reprocess SnapShot Freeze images for artifacts caused by motion blurring after manual vessel editing
- **Motion blurring correction** requires the use images produced by a CT using the optional SnapShot Freeze feature

#### Protocols

- CardIQ Xpress 2.0 Reveal is supplied with a set of pre-defined protocols that are easily adaptable to customized protocols
- The protocols have comprehensive instructions and tools which minimize the need to memorize the procedures or to refer continuously to user documentation.

#### Indications for Use

CardIQ Xpress 2.0 is intended to provide an optimized non-invasive application to analyze cardiovascular anatomy and pathology and aid in determining.

treatment paths from a set of Computed Tomography (CT) Angiographic images CardIQ Xpress 2.0 is a CT, non-invasive, image analysis software package, which aids in diagnosing of cardiovascular disease to include, coronary artery disease, functional parameters of the heart, heart structures and follow-up for stent placement, bypasses and plaque imaging.

CardIQ Xpress 2.0 offers unique tools such as automatic tracking, which will pre-process the CT data into multiple viewing ports to allow for an expedited read time improving workflow. With CardIQ Xpress 2.0, the user can color code the myocardial tissue to show hypo/hyper-dense areas in the myocardial tissue of the heart. With the IVUS-like view the user can color code the HU units of the plaque to better visualize the difference between calcified and non-calcified plaque in the wall of the vessel and the lumen to determine the amount of atherosclerosis.

The user can see the different valve planes along with a variety of new layouts to align the heart. The IVUS-like view is created by applying GE's Volume Rendering on a cross-section perpendicular to the detected centerline. This view merely displays a cross section as in IVUS imaging and

color codes like IVUS images. No new or additional diagnostic information is added. CardIQ Xpress 2.0 is for use on the Advantage Workstation (AW) platform, CT Scanner, PAC or Centricity stations, which can be used in the analysis of **2D or 3D CT angiography images**/data derived from DICOM 3.0 CT scans.

#### Regulatory Compliance

This product complies with the following requirements:

- European Council Directive 93/42/EEC concerning medical devices.

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
SECRET DE AFACERI

# CardIQ Function Xpress

Real time cardiac review the instant you're ready to read.


Cardiac disease is one of the leading health concerns world wide. Successful treatment of the many conditions that cause and perpetuate heart disease requires that physicians approach specific cardiac problems with as much information as today's technology can provide. To supply your referring physicians with this information you need software that gives you insight into cardiac anatomy and function supported by quantifiable data.

## Overview

Cardiac IQ Function Xpress post-processing software helps you evaluate cardiac function and diagnose cardiovascular disease with a high degree of confidence. Providing accurate and reproducible quantification of left and right ventricular volumes, ejection fractions, and myocardial mass, Cardiac IQ Function Xpress is optimized to assess cardiac function using multiphase, multi-slice cardiac CT images. The CardIQ Function Xpress option provides an easy-to-use and time-effective way for you to perform cardiovascular functional analysis.

## What's new

- Automatically pre-processes and loads exams.
- Automatically detects all chambers in all phases.
- Automatically performs ejection fraction and chamber volume analysis with 91% reliability on LV and RV automatic ejection fraction.
- Provides single click myocardial analysis with bull's eye images.
- Performs simultaneous case review in real time.
- Lets you load multiphase data with no load time.
- Lets you review multiple cases in real time.
- Lets you generate reports with imported analytics.

SECRET DE AFACERI

CONFIDENTIAL

Visit us:

[www.gehealthcare.com/aw/applications/cardiq-function-xpress/](http://www.gehealthcare.com/aw/applications/cardiq-function-xpress/)

### Features

- Apply the CardIQ Function Xpress application to standard axial CT images acquired on a GE scanner using the cardiac CT SnapShot imaging acquisition option.
- Automatically selects each heart chamber for individual volume analysis.
- Performs behind-the-scene function data process for real time ejection fraction review, volume analysis, and myocardial analysis.
- Enables you to extract, render, and display 3D volumetric models of the endocardium for ejection fraction calculation.
- Automatic detection of the endocardial and epicardial walls for wall motion, wall thickness, wall thickening and myocardial mass analysis. Lets you perform volume analysis of all chambers.
- You can automatically calculate left atrium volume while excluding the pulmonary vein.
- With a single click you can activate visual wall motion with short axis images in basal, mid, and distal orientations along with a two-chamber long axis view.
- You can perform myocardial analysis with wall motion, wall thickness, and mass calculations.
- The application's flexible reporting tool lets you include graphical representations.
- Display on screen a table of key functional parameters for instant visualization.

### Intended Use

CardIQ Function Xpress is intended to provide an optimized non-invasive application to analyze cardiovascular anatomy and pathology and aid in determining treatment paths from a set of Computed Tomography (CT) Angiographic images.

### Image Requirements

CardIQ Function Xpress accepts standard CT image sets acquired on qualified GE CT scanners using the appropriate cardiac imaging software. These images must meet the same image requirements as those for the basic Volume Share 2 application.

### System Requirements

- AW Volume Share Workstation with Volume Share 2 running 4 8400 hardware or higher (running 64 bit).
- Auto launch and preprocessing are available only on 8400 workstations with 16 GB RAM.
- 2 monitor configuration.
- Color landscape monitor.

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GE imagination at work

**CONFIDENTIAL**  
**SECRET DE AFACERI**


© 2012 General Electric Company.  
All rights reserved. Data subject to change.  
GE and GE Monogram are trademarks of General Electric Company.  
\* Trademark of General Electric Company.


## Overview

CardIQ Function Xpress allows you to accurately and reproducibly quantify myocardial mass, and left and right ventricular volumes and ejection fractions. The software automatically detects left ventricular endocardial and epicardial contours so that you can more easily assess left ventricular (LV) and right ventricular (RV) functional parameters. The program is optimized to perform assessment of cardiac function using GE multiphase, multi-slice cardiac CT images.

The CardIQ Function Xpress menu page consists of multiple image processing tools for single or multiphase cardiac studies. These protocols can be accessed either from the AW Volume Viewer application or directly from the application menu page on the browser.

CardIQ Function Xpress benefits to cardiac patients include:

- Minimally invasive alternative to help diagnose heart function.
- Data is acquired during a coronary CTA exam, so no additional dose or contrast is required.
- Offers a diagnostic option to patient who may not be able to undergo other non-invasive exams.

## Feature Details

CardIQ Function Xpress uses single or multiphase image data sets to measure, graphically plot, and tabulate a number of ventricular function and wall motion parameters.

### Pre-processing

- The program automatically processes multiphase data for functional analysis.
- You can select configurations to process only the chambers you require.
- For more streamlined reading, ejection fraction, volume analysis, and myocardial analysis are completed before you view the case.

### Auto Launch

- Load multiple cases into auto launch so they are ready to review when you are.
- Switch between exams without having to quit the application and reload.
- Exams ready for review display in green.

### Automatic Chamber Detection

- The program automatically loads single or multiphase data sets and selects the heart's chambers for volume analysis.
- Analyze individual chambers by automatically contouring any of the heart's four chambers.
- The system automatically detects the left and right ventricles and the left atrium and contours the endocardium for ejection fraction calculation.
- Choose single or multi-phase chamber detection.

### Automatic Ejection Fraction

- The system features a protocol that, when chosen, automatically detects the end diastolic and end systolic phases and calculates an accurate ejection fraction of the left ventricle, right ventricle, and left atrium.

- Ejection fraction measurements are available for

SECRET DE AFACERI  
CONFIDENTIAL


GE imagination at work

the right and left ventricles.

- Stroke volume, end diastolic and end systolic volumes are automatically calculated and shown in a graphical representation for assessment of heart failure.
- You can easily send ejection fractions and other quantification results to your report.
- With dual screens, you can display quantification assessments on the left monitor and a qualitative overview of the full chamber from apex to base on the right monitor.
- Customize layout designs to accommodate your viewing preferences:
  - Base to apex
  - Basal, mid, and distal along with a two-chamber long axis.
- You can change the selection of end diastolic and end systolic phase.
- While contour detection is automatic, scalpel and auto select tools are available with which you can edit the contours.

### Display Table

- A display table shows volume and provides a graphical update of end diastolic and end systolic locations.

- Click on available tabs to select different chambers for analysis.
- Left and right ventricle measurements include:
  - Stroke volume
  - Ejection fraction
  - Cardiac output
  - Myocardial mass
  - End systolic volume
  - End diastolic volume
  - Pulmonary vascular resistance
  - Systemic vascular resistance
  - Body surface area
  - Heart rate

- Left atrium measurements include:
  - Min and max volume
  - Ejection fraction
  - Left atrium index
  - Body surface
  - Heart rate
- Right atrium measurements include:
  - Min and max volume
  - Heart rate
- You can film, save, and export the volume graphs.

### Myocardial Analysis

- CardIQ Function Xpress automatically detects the endocardial and epicardial walls to analyze wall motion, wall thickness, wall thickening, and myocardial mass analysis.
- You can edit the contours with a freehand trace and click to spline tools.
- The program displays multiple bull's eye graphs to graphically represent these cardiac characteristics:
  - Wall motion
  - Wall thickness
  - Myocardial mass
  - Cardiac output
- Choose from four bull's eye graph overlays to show 17, 20, and 49 segment model, or a coronary artery territory model.
- Manual scoring is available for wall motion and thickening in 17 and 20 segment.
- Cord display is available for visual assessment of wall abnormalities.

### Reporting

- The application provides you with automatic report generation capabilities.
- You can import hospital or clinic logo and address to make your reports site specific.
- You can capture 3D volume rendering endocardium images with ejection fraction calculations displayed for inclusion in your report.
- You can export the volume graph of all heart chambers to your report.
- The program automatically fills these report fields:
  - Ejection fraction
  - Stroke volume
  - End diastolic volume
  - End systolic volume
  - Pulmonary vascular resistance
  - Systemic vascular resistance
  - Myocardial mass
  - Cardiac output
- You can create predefined conclusion templates for various types of findings.
  - Report templates are fully customizable, so you

SECRET DE AFACER  
CONFIDENTIAL

GE imagination at work

CardIQ Function Xpress - 2

can include logos, diagrams, signatures, patient history, and messages.

- Generate reports in PDF, HTML, or plain text.
- You can combine CardIQ Function Xpress reports with CardIQ Xpress 2.0 Elite reports.

### Simplified Workflow

CardIQ Function Xpress simplifies your workflow with:

- Automatic ejection fraction protocol
- Streamlined multiphase beating images for easy visualization of cardiac wall motion.
- Single click myocardial analysis with bull's eye images.
- Simplified reporting.

### Summary

CardIQ Function Xpress provides you with the tools you need to accurately analyze, evaluate, and report on cardiac anatomy and function. Its ease of use will streamline your workflow and its accuracy will give you greater diagnostic confidence.

### Indication for Use

CardIQ Function Xpress is used in conjunction with CT cardiac images to automatically calculate and display various left ventricular and right ventricular functional parameters such as ejection fraction, end systolic and end diastolic volumes, stroke volumes, wall motion, wall thickening, cardiac output, myocardial mass, and systemic and pulmonary vascular resistance. Volume measurement of each chamber of the heart is also available. With CardIQ Function Xpress atrium volumes may be used to determine volume assessment of atrial disease to including, but not limited to, atrial fibrillation. CardIQ Function Xpress is a CT, non-invasive image analysis software package, which aids in the assessment of cardiac function and in determination of cardiovascular disease diagnosis and management. CardIQ Function Xpress is available on the AW Workstation, CT Scanner, PAC, or Centricity stations, which can be used in the analysis of 2D or 3D CT angiography images/data derived from DICOM 3.0 CT scans.

CardIQ Function Xpress is an extension of Volume Viewer 4 and is dedicated to cardiovascular imaging. It provides an optimized non-invasive application to analyze cardiovascular anatomy and pathology and aid in determining treatment paths from a set of CT Angiographic images.

**SECRET DE AFACERI**  
**CONFIDENTIAL**


GE imagination at work

© 2012 General Electric Company.  
All rights reserved. Data subject to change.  
GE and GE Monogram are trademarks of General Electric Company.  
\* Trademark of General Electric Company.

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\*Marcă comercială a General Electric Company

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
**SECRET DE AFACERI**  
**CONFIDENTIAL**

# VessellIQ™ Xpress with AutoBone™ Xpress

Fast and Easy to Use Vessel Analysis with Automated Real-Time Tracking, Auto-Labeling and Simplified User Interface.

Peripheral vascular disease (PVD) is a condition similar to coronary artery disease and carotid artery disease. In PVD, fatty deposits build up in the inner linings of the artery walls. These blockages restrict blood circulation, mainly in arteries leading to the kidneys, stomach, arms, legs and feet.

There are different clinical techniques used to diagnose PVD. Computed Tomography Angiography (CTA) can potentially provide sufficient information for therapeutic decision-making and reduced imaging costs for the management of PVD. One of the advantages of computed tomography (CT) is that it shows both vascular & nonvascular anatomy and also provides 3D information for best view. Also, CT may be more useful for patients with limited or no peripheral access or with stents, who may not be candidates for magnetic resonance angiography (MRA). Complementary to the ability to perform CTA is to perform clinical assessment by a dedicated post-processing tool.

## Overview

VessellIQ™ Xpress with Autobone Xpress is a post-processing software package for the Advantage Workstation (AW) platform, AW Server, CT scanners and PACS reading stations. It is a tool to analyze CT angiographic data including stenosis analysis, thrombus, pre and post stent planning procedures, and directional vessel tortuosity visualization. Autobone Xpress provides zero-click bone segmentation for head, neck and other anatomy. The automation and ease of use streamlines workflow.

## Highlights

- Automated detection of aorta and iliacs with auto-labeling of vessels
- Automated real-time fast tracking for all vasculature
- New contextual wheel menu including easy access to advanced options
- New centerline edition capabilities
- Improved predefined measurements and edit capabilities
- New reporting interface which includes measurements in summary table
- New user interface

SECRET DE AFACER  
CONFIDENTIAL

Visit us:

<https://www.gehealthcare.com/products/advanced-visualization/all-applications/autobone-vessellq-xpress>

## Features

- Anatomy based protocols for guided workflow in the analysis for your peripheral vascular CTA acquisitions.
- Automatic bone removal for efficient review of vasculature.
- Automated abdominal aorta CTA protocol with vessel tracking and vessel labeling.
- Fast tracking detects entire length of vessel and displays the centerline based on a single start point.
- Advanced editing tools bringing new flexibility for faster editing based on targeted anatomy.
- Thrombus detection with measurements and comparison over time.
- Dynamic AVA provides access to vessel tracking at any time, from any protocol, on any vessel.
- Advanced measurement tools for stenosis detection and color identification of HU value for quantification of different densities.
- Create a queue of cases that are pre-loaded and pre-processed in the background while you read another exam.
- Measurements are organized in the summary table allowing for easy viewing and exporting.
- Fully customizable workflows with adjustable layouts, personalized parameter and settings, custom review steps.

- Save State lets you save and restore the state of the processed images at any stage

## System Requirements

- AW Workstation AW 4.7 ext 14 or higher
- AW Server 3.2 Ext. 3.2 or higher

## Indications for Use

VessellQ Xpress is an optional, non-invasive, optimized, post-processing application intended to provide images and tools to analyze vascular anatomy and pathology, aiding physicians in diagnosis and determination of treatment paths, from a set of Computed Tomography (CT) Angiographic images. VessellQ Xpress is an option for the Advantage Workstation (AW) platform, CT Scanner, and/or PACS, which **can be used in the analysis of 2D and 3D CT Angiography images/data** for the purpose of cardiovascular and vascular disease assessment. This software-only device is designed to support physician assessment for a wide variety of clinical uses such as stenosis analysis, pre/post stent planning, pre/post valve replacement planning, and directional vessel tortuosity visualization. VessellQ Xpress' automatic visualization tools provide users the capability to segment bony structures for accurate identification of the vessels. Additional tools enable analysis of the vascular anatomy including the aorta, valves, and branching vessels for: anatomical sizing;

density and volume analysis of segmented vasculature and calcified / non-calcified plaque; and measurements of abnormalities. The TAVI Analysis option for VessellQ Xpress is a planning tool used for Trans Aortic Valve Implantation (TAVI) procedures. It automatically segments the aorta and displays the aortic valve in multiple views for measurements of anatomic structures commonly needed for aortic annulus replacement planning. TAVI Analysis provides guided workflow and semi-automated tools to aid in evaluation of appropriate access pathways for interventional procedure planning.

AutoBone Xpress option is a software package that is intended to facilitate segmentation of bony structures and calcifications from abdominal and extremity CT Angiography data.

## Regulatory Compliance

This product complies with Regulation (EU) 2017/745 of the European Parliament and of the Council on medical devices (MDR). This product or its features may not be available in some other countries or regions. Please contact your sales associate.

**Rx Only**

**SECRET DE AFACERI**

**CONFIDENTIAL**

GE imagination at work

© 2021 General Electric Company.  
All rights reserved. Data subject to change.  
GE, GE Monogram, imagination at work, AutoBone Xpress, VessellQ Xpress are  
trademarks of General Electric Company.  
All other third-party trademarks are the property of their respective owners  
Doc1471990

## Features Detail

### Anatomy-based protocols

- The program offers you predefined protocols for fast, easy review of all vascular structures:
  - Intra-cerebral
  - Carotids
  - Thoracic Aorta
  - Abdominal Aorta
  - Runoff
  - Spinal
  - Extremity
- Protocols may be customized to meet your requirements and saved for quick access.

### AutoBone™ Xpress

- Anatomical driven protocols with automatic bone removal.
- Improve vessel visualization by removing obstructive bony detail


- Automatic calcium removal for vessel lumen visualization

### Auto Aorta


- With a single click, the program automatically removes bones and performs centerline tracking of the aorta and iliacs and labels vessels.

- Simply confirm the automated tracking or modify it by adding additional branch points, deleting a single point, or deleting multiple points from any view.

### Fast Tracking

- Semi-automatic detection of complete vessel from a single starting point

- Advanced editing with variable interpolation tool for quick centerline edits adjusted for length of vessel.

- Bulls eye adjustment capability on x-section views to perfect tracking of the smallest most tortuous vessels.
- Volume measurements of the thrombus allow you to track change over time. Additional measurements include :
  - Minimum diameter
  - Maximum diameter
  - Mean diameter
  - Cross-sectional area

### Vessel Label Data Base

- Each protocol is associated with a database of vessel names
- Customize vessel label database to meet your requirements for every anatomical region.
- A right mouse click provides access to the list of vessel names per anatomy.


### Dynamic AVA

- Real-time tracking of the vessel center line and instant visualization of an unfolded view for quick vessel inspection.
- Immediate access to vessel analysis tools.
- Edit contours with a selection of drawing tools with variable interpolation for fast editing.


CONFIDENTIAL

SECRET DE AFACERI

## Advanced Vessel Analysis Tools

- Semi-automated vessel analysis lets you track, extract, visualize, and size vessels.
- Track a vessel and its branches for quick review in curved, oblique longitudinal or cross-sectional views.
- Multi-point deposition capabilities for improved center line tracking of challenging cases.
- When you edit a vessel contour, the program automatically interpolates the edit to neighboring slices.
- Generate and save rotating movies
  - Curved planar reformation and cross-sectional (lumen) views
  - Best L-section reformatted images
- Select vessel segments for quantitative or qualitative analysis.


- Available measurements include:
  - Distance
  - Cross-sectional area

- Single or dual reference point comparison
- Relative percent stenosis / aneurysm
- Color Identification
- 3D Angle
- Thrombus Extraction
- Volume
- Mean, min, max diameter
- 
- Move freely between review steps and protocols without leaving the exam.

### Stenosis Detection & ColorID

- With the customizable one- or two-point deposit technique, you can automatically quantify degree of stenosis and segment length.
- By assigning customizable color mapping to CT number ranges within an image, you can identify plaque on an axial image, curved reformat, or MPVR image. Up to ten levels of color coding are available. Colors transition smoothly from one density to another.
- Make volume and area measurements to track plaque size over time.

- Custom sizeable 3D ROI lets you apply color over all anatomy.
- Generate multiple ROIs on one vessel or multiple vessels.
- Automatically include each defined plaque, its volume and area into the report.

### Quick AVA

- Dynamic AVA, provides access to vessel analysis from any protocol.
- Launch the Quick AVA feature from the main tool bar.
- Perform vessel analysis from any 3D or reformatted image with one- or two-point clicks.
- If you detect a lesion, you can analyze it by dropping a point above and below the section. The vessel analysis tool launches automatically for fast review in curved reformat, cross-sectional, lumen, and MRP view.
- Extend tracking proximally or distally for extended views at any time.

### Summary Table & Report

- Measurements are automatically captured in summary table for easy review.

- Simply click on any measurement in the summary table and you will be taken to that measurement location.

CONFIDENTIAL  
SECRET DE AFACERI

- Prepare report with a single click or generate an advanced report with customized options
- Send the report to database, film or PDF.

### Save State & Protocol

- Save the current state of tracking and measurements for easy reloading with a single mouse click.  
Save state captures:
  - 3D images
  - Vessel tracking
  - Measurements
  - Vessel contours.
- Transfer the Save State to the PACS or other AW workstations for future review.
- Save custom protocols including anatomical labels and measurements.

### Summary

VessellIQ Xpress provides user- friendly tools an intuitive interface to help streamline your vessel analysis workflow and help you generate more comprehensive reports for referring physicians.

### Multiple Study Comparison

- With VessellIQ Xpress, compare a previous study with the current study. You also have the flexibility to launch comparison within Quick AVA.

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**SECRET DE AFACERI**  
CONFIDENTIAL

## SmartScore 4.0

Advanced imaging software that detects, quantifies and scores cardiac calcium plaque burden. Instantly.


Cardiovascular disease remains one of the most common health issues in the world today. As with many conditions, early detection and patient risk assessment are vital to preventing or minimizing long-term negative effects. But many cases of this disease aren't diagnosed until the patient presents with symptoms. Conventional procedures to assess risk can be expensive, time consuming, and uncomfortable. An accessible, patient-friendly risk assessment method could help clinicians devise a regimen for their patients that might lessen the chance for serious cardiac events.

### Overview

SmartScore 4.0 is designed to identify the presence of regional and global coronary artery calcification from a CT scan, then measure and score the results. Scores can be calculated using a standard Agatston/Janowitz (AJ) method. When correlated with a patient's personal information, the score can yield an estimation of a patient's risk for coronary artery disease.

### What's new

- Non-invasive alternative to conventional assessment procedures.
- Score can be correlated to age group cohort to determine patient's risk per population
- Provides information on coronary artery wall calcium plaque buildup.
- Automatically detects calcium and highlights it in green.
- Streamlines workflow by networking patient demographics from CT scanner directly to SmartScore program.
- Free Hand Trace lets you outline specific ROIs.

SECRET DE AFACERI

CONFIDENTIAL

Visit us:

[www.gehealthcare.com/ow/applications/smartscore-4/](http://www.gehealthcare.com/ow/applications/smartscore-4/)

### Features

- AJ 130 scoring method uses conventional Agatston/Janowitz technique with a threshold of 130 HU which is adjusted to the appropriate image slice thickness.
- Volume scoring calculates volumes ( $\text{mm}^3$ ) of calcified plaque above the 130 HU threshold.
- Mass scoring calculates mass (mg) of calcified plaque above the 130 HU threshold.
- Individual and aggregate scores are computed for each artery type.
- Provides two methods of calcium scoring.
- Customizable parameters can be predefined based on your site's preferences.
- User interface streamlines your workflow and reduces reading time.
- Report options allow you to custom tailor your reports and distribute them in a variety of formats.

### System Requirements

- EKG monitor with recording device and x-ray translucent lead cable.
- Gantry hardware upgrade kit for those scanner systems already in operation.

### Minimum platform release:

- Software for Advantage Workstation 4.2P or higher.
- AW Server 2.

### Recommended Options

#### Postscript Printers:

- Codonics: NP-1660M
- Kodak: 3600 DMI
- Codonics 1660M, 1660MD or Horizon
- Lexmark Optra 1650N, 1855N, SC1275N, C710N, C72N, T612, or T614
- Seiko 1720D
- Kodak Dmi3600
- Quantum GL2101HD with film/thick paper
- Quantum GL2101HD with plain paper (see PI-008)
- Tally T8106
- HP LaserJet
- Xerox Phaser

### Intended Use

SmartScore is a non-invasive software option that can be used to evaluate calcified plaques in the coronary arteries, which may be a risk factor for coronary artery disease. SmartScore may be used to monitor the progression/regression of calcium in coronary arteries over time, which may aid in the prognosis of cardiac disease.

### Regulatory Compliance

This product complies with the European CE Marking regulation for Medical Devices Directive: Directive 93/42/EEC, dated 14 June 1993.

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
Vizitați-ne:  
[www.gehealthcare.com/aww/application/gsi-viewer/](http://www.gehealthcare.com/aww/application/gsi-viewer/)

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
*afodor*

## CT Perfusion 4D

**Comprehensive perfusion analysis available for neurology, oncology and cardiology.**

Easy-to-use, protocol-driven workflow for all organs.

Perfusion 4D has been expanded beyond stroke and oncology to now include dynamic myocardial perfusion. For any organ where blood perfusion is of interest Perfusion 4D is a complete package to quantitatively provide the necessary information to guide treatment decisions.

The protocol-driven design leads the user step-by-step through the process, reducing keystrokes and improving repeatability. Get the information you want quickly and reliably.

### Overview

GE's Perfusion 4D is a fast, easy-to-use automated software for analyzing CT Perfusion images related to stroke, tumor angiogenesis and dynamic myocardial perfusion. Its simple user interface and automated perfusion post-processing make it easy for you to diagnose quickly and accurately – and make treatment decisions more confidently.

GE leverages its deconvolution expertise, with its innovative delay correction algorithm. Perfusion 4D takes perfusion assessment to the next level, with the addition of CT Dynamic Myocardial Perfusion. Providing the classic color overlays with quantitative measurements to evaluate myocardial blood perfusion.

### What's new

- **Dynamic myocardial perfusion maps**
- Improved IQ of functional maps in presence of noise.
- Improved processing times
- Improved ROI mirroring for brain protocols
- Streamlined workflow for tissue classification.
- Permits injection rates of 4cc/sec.
- Incorporates the GE delay-corrected deconvolution method.
- Dynamic registration

**SECRET DE AFACERI**  
**CONFIDENTIAL**

Visit us:  
[http://www3.gehealthcare.com/en/products/categories/advanced\\_visualization/applications/ct\\_perfusion\\_4d\\_multi-organs](http://www3.gehealthcare.com/en/products/categories/advanced_visualization/applications/ct_perfusion_4d_multi-organs)

## Features

- Visualize all information in true volumetric form with the ability to employ all the volumetric-based image analysis tools in the AW Volume Viewer.
- Whole heart perfusion with semi-automatic guided workflow to aid in assessment of myocardial ischemia
- Smart Maps using the intelligent 4D noise suppression algorithm improves functional map image quality in the presence of noise.
- Streamlined tissue classification and visualization incorporates thresholding of selected functional maps.
- Dynamic non-rigid registration protocols for dynamic acquisitions within the body.
- Productivity enhancements include:
  - Ability to select automated vessel detection or manual vessel detection for Brain Perfusion
  - Faster processing times for functional map computations
  - Volumetric visualization of functional maps
  - Access from anywhere using AW Server.
  - Time stamps on functional maps and Regions of interest (ROIs).
  - Interactive large vessels exclusion.
- Many basic features, including the ability to:
  - Process double phase acquisition data.
  - Create 4D Regions of interest (ROIs).
- Display Regions of interest (ROI) statistical information.
- Display averaged time intensity information of the Regions of interest (ROI).
- Simultaneously review any number of functional maps.
- Customize protocols
- Save parameters.
- Includes these CT post-processing protocols:
  - Brain stroke
  - Brain tumor
  - Dynamic Myocardial Perfusion
  - Body tumor
  - Liver perfusion
  - Pancreas perfusion
  - Prostate perfusion
  - Kidney perfusion
  - Soft tissue perfusion
  - Spleen perfusion
  - Bone perfusion
  - Liver Dynamic Registration
  - Cardiac Dynamic Registration
- Functional maps:
  - Regional cerebral blood volume
  - Regional cerebral blood flow
  - Regional mean transit time
  - Contrast arrival delay
  - Transit time to peak of impulse residue function
  - Capillary permeability surface area for lengthened acquisition protocols
  - Mean Slope of Increase
  - Base image
  - Average image

## System Requirements

Advantage Workstation Volume Share

4.7 or higher, AW Server 3.2 or higher, Centricity™ Universal Viewer

- Z820, Z440 (and later) with 24GB or higher of RAM for the AW
- Recommended monitor resolution is up to dual 2MP (1600 x 1200) or a single 3MP (1536 x 2048) for the AWS

## Intended Use

CT Perfusion 4D is an image analysis software package that allows the user to produce dynamic image data and to generate information with regard to changes in image intensity over time. It supports the analysis of CT Perfusion images, obtained by cine imaging (in the head and body) after the intravenous injection of contrast, in calculation of the various perfusion-related parameters (i.e. regional blood flow, regional blood volume, mean transit time and capillary permeability).

The results are displayed in a user-friendly graphic format as parametric images. This software will aid in the assessment of the extent and type of perfusion, blood volume and capillary permeability changes, which may be related to stroke or tumor angiogenesis and the treatment thereof.

## Regulatory Compliance

This product complies with the European Council Directive 93/42/EEC Medical Device Directive as amended by European Council Directive 2007/47/EC.

GE imagination at work

SECRET DE AFACERI

CONFIDENTIAL

© 2017 General Electric Company.  
All rights reserved. Data subject to change.  
GE and GE Monogram are trademarks of General Electric Company.  
\* Trademark of General Electric Company  
DOC1996217

## Overview

CT Perfusion 4D – is a complete image analysis software package that allows the evaluation of dynamic CT data following an injection of a compact bolus of contrast material, and generating information regarding changes in image intensity over time.

It provides a quick and reliable assessment of the type and extent of perfusion disturbances by providing qualitative and quantitative information on various perfusion related parameters (i.e. regional blood flow, regional blood volume, mean transit time, capillary permeability and tissue classification), which may be related to acute stroke, tumor angiogenesis or myocardial ischemia and the treatment thereof.

The results are displayed in a user-friendly interactive format as functional images.

CT Perfusion 4D has been expanded to include dynamic myocardial perfusion. This protocol provides a complete workflow to guide the user through image preparation and generation of perfusion maps of the myocardium to aid in the assessment of myocardial perfusion abnormalities. Additionally, non-rigid registration protocols have been included within CT Perfusion 4D to assist in organ registration over time which is necessary in dynamic imaging. The dynamic registration protocol can be applied to any dynamic acquisition to minimize mis-registration between passes.

CT Perfusion 4D incorporates *smart smoothing*, a smoothing algorithm, which improves the image quality of the functional maps in the presence of noise.

Also included is a streamlined workflow for Tissue Classification for neuro perfusion and its visualization, incorporating interactive thresholding of selected functional maps. Tissue Classification may aid the clinician in determining the status of the tissue based on blood volume and one of blood flow, mean transit time, or transit time to IRF peak (Tmax) maps, where the first 6 hours after the onset of symptoms are critical in identifying the occurrence of stroke and the follow-up treatment.


## Inject at 4 cc / second

The deconvolution method employed is not dependent on the assumptions about the underlying hemodynamics of the circulation at the tissue level, and also permits slow injection rates of 4 cc/s, which makes the method more acceptable in patient studies.

The software uses the arterial input function to correct for the dispersion and re-circulation of contrast agent, and venous output function to normalize the calculated functional maps.

## Productivity

- **Tissue Classification:** Provides an interactive thresholding algorithm that aids clinicians in determining the status of the tissue based on Blood Volume and one of Blood Flow, Mean Transit Time, or Tmax maps.
- **Productivity enhancements:**
  - Robust automated vessel detection for artery and vein curves definition
  - Faster functional maps, calculated in a fraction of previous times
  - Volumetric visualization of functional maps
  - Access to CT Perfusion 4D from anywhere using AW Server
  - Interactive large vessels exclusion from computed functional maps.
  - Time stamps on functional maps and ROIs for better data management and tracking.

Blood Flow, Blood Volume,  
MTT and Tmax maps

## GE Delay-Corrected Deconvolution Method

Advantages of using the delay - corrected deconvolution algorithm include the fact that no assumptions need to be made regarding the absence of venous outflow, or the presence of the draining vein or the uniformity of the characteristics of perfusion in the brain. These assumptions may well lead to larger variability in the outcome of the perfusion measurements. The de-convolution method computes the required results using linear optimization techniques to ensure that the global minimum is obtained in contrast to non-linear optimization techniques. The delay corrected component of the algorithm will adjust the functional maps values based on the actual arrival of contrast to the tissue in question rather than overestimate the mean transit time by including the transit time of the contrast from the selected artery to the tissue.

CONFIDENTIAL

SECRET DE AFACERI

## Basic Features

- Process double phase acquisition data (data with minimal temporal sampling and data with longer temporal sampling in one acquisition to cover entire functional perfusion mechanism).
- Create 4D ROIs (multi-locations ROI with time-varying position).
- Display ROI statistical information (Area, Min., Max, Std Dev.).
- Display input ROI from several locations on the same graph.
- Display the averaged time-intensity information from an ROI (VOI).
- Simultaneously review any number of functional maps.
- Customize any viewport to series, graph or functional view.
- View series or functional data in axial, sagittal, coronal or oblique view.
- View series or functional data in 3D or VR mode
- Customize protocols (advanced settings, color ramp and/or W/L used for each generated functional map) and to save these user preferences.
- Save the state of a protocol including parameters, ROIs, functional maps and pre/post-processing results for further reloading in CT Perfusion.
- Save all images as processed images in the database.
- Save all images and graph data in Screensave and to print them to film or color paper printer.
- Save ROIs statistical information using Screensave.
- Print to film or color paper printer.
- Save parameters, used to generate functional maps, as DICOM objects and to restore them including ROIs with the corresponding maps in CT Perfusion.

## Protocols


CT Perfusion 4D has fully integrated protocols that assist the user in obtaining results efficiently for neurology, oncology and cardiology core areas.

CT Perfusion 4D package protocols are:

**CT Brain Stroke Protocol.** The CT brain stroke protocol will open a loading panel containing information about the selected perfusion series in a tabulated format. It is possible to remove unwanted time points from the loading sequence. Once the exam is loaded the SW will perform the following algorithmic steps:

- **Image Registration** to minimize the effects of patient movement.
- **Bone Removal** The cranium is automatically removed from the dataset.
- **Arterial Input.** Used as the input for the deconvolution algorithm. Artery is automatically selected and can be easily adjusted if needed.
- **Venous output** for partial volume correction. Vein is automatically selected and can be easily adjusted if needed.
- **Last pre-enhancement / First post-enhancement image.** These parameters are selected to ascertain correct input for functional maps calculations.
- **Final Settings and Compute.** The software will compute the functional maps with the saved settings, but can be made to stop and recompute the functional maps if changes to the algorithm settings were implemented. The changes may include any of the following:
  - Algorithm Resolution, Smart Smoothing size.
  - Maximum Blood Flow, Brain Density, Hematocrit Ratio.


The software will use deconvolution to generate all functional maps, *Blood Flow, Blood Volume, Mean Transit Time, and Transit Time to IRF Peak (Tmax)*. The Brain Stroke protocol has multiple review steps to guide the user through the functional information. Review settings: After loading the dynamic exam to process a review step is initiated. The general purpose of this step is to confirm the automatic selection of the input vessels for the brain stroke or tumor protocols, or to conduct a general review of the time resolved perfusion data and place ROI's to select vessels.


CONFIDENTIAL

SECRET DE AFACERI

- **Set Symmetry:** A review step in the brain stroke or tumor protocol that allows users to set a symmetry axis needed for mirroring the ROIs.


- **Review Functional Maps:** A review step enabling the viewing and assessment of functional maps such as Blood Flow, Blood Volume, Mean Transit Time, and Tmax. ROI mirroring functionality makes it easier to compare hemispheres.

- **Tissue Classification:** A review step specific to the brain stroke protocol, which enables the visualization of regions that are created as a result of segmentation of the tissue using the interactive thresholding of the Blood Volume and one of Blood flow, Mean Transit Time or Tmax maps.

The Brain Stroke protocol is typically used for the assessment of cerebral hemodynamic disturbances e.g. stroke and traumatic brain injury.

**CT Brain Tumor Protocol.** Essentially identical to the brain stroke protocol but without the Tissue Classification review step.

**CT Body Tumor Protocol.** The Body tumor protocol produces functional maps automatically once the user selects an artery using a simple and intuitive workflow.

The Body Tumor protocol also leverages most of the advanced algorithm settings that are part of the Brain Stroke protocol. The protocol has the following algorithms

- **Noise Threshold.** Thresholds can be customized to exclude air and bone from the calculations.
- **Artery Input.** System allows for an intuitive selection of the artery.
- **Last pre-enhancement / Last post-enhancement image / Last Second Phase Image** (last image of the phase 2 data). These can be selected to ascertain correct input into the parameter calculations.
- **Final Settings and Compute.** The user can compute immediately or customize the following advanced settings:
  - Algorithm Resolution, Smart Smoothing.
  - Maximum Blood Flow, Brain Density, Hematocrit Ratio.

The Compute command then engages the exclusive de-convolution algorithm to produce the parametric images.

**CT Pancreas Perfusion Protocol.** Identical to the Body Tumor Protocol, but optimized for pancreas tumors.

**CT Prostate Perfusion Protocol.** Identical to the Body Tumor Protocol, but optimized for prostate tumors.

**CT Kidney Perfusion Protocol.** Identical to the Body Tumor Protocol, but optimized for kidney tumors.

**CT Soft Tissue Perfusion Protocol.** Identical to the Body Tumor Protocol, with an additional step to define the vein output for partial volume correction (automatic or manual venous selection).

**CT Spleen Perfusion Protocol.** Identical to the Body Tumor Protocol, but optimized for spleen tumors.


**CT Bone Perfusion Protocol.** Identical to the Body Tumor Protocol, but optimized for bone tumors.

**SECRET DE AFACERI**


**CONFIDENTIAL**

**CT Non-Rigid Registration.** Fully automated protocols used for dynamic body acquisitions to register data for use in perfusion analysis. Provides simultaneous review of original data with the registered results for full confidence in the newly created series. Saved registered results are to be loaded into the appropriate body or cardiac perfusion protocol for functional analysis. Two registration protocols are provided. CT Cardiac Dynamic Registration and CT Liver Dynamic Registration.

**CT Liver Perfusion Protocol.** Identical to the Body Tumor Protocol, with an additional step to define the portal vein input (automatic or manual venous selection).

Processing times have improved exponentially, now with an average processing time of one to three minutes depending on slice thickness. **Hepatic Arterial Fraction maps have been improved and can be an efficient tool for tumor detection. A new map, hepatic arterial blood flow (HABF), has been added to assist in tumor assessment.**


#### Dynamic CT Myocardial Perfusion Protocol.


Complete guided workflow for assessing dynamic myocardial perfusion. The dynamic myocardial perfusion protocol is laid out in a series of three easy to follow steps. Start by selecting the series that has been generated from the non-rigid registration protocol and then follow the guided workflow:

- 1) **Cardiac Reformat** – images are automatically orientated into short and long axis views,

- 2) **Segmentation** – define the valve plane and the apex of the heart. The software will then calculate the endocardial and epicardial borders and provide a 17 segmented map of the myocardium along with a graphed view for each segment

- 3) **Functional Maps** – the third and final step produces the quantitative data with the color overlay and bullseye maps for blood flow, blood volume, mean slope of increase and mean transit time.

Quantitative data can be exported as a .csv file for a complete segmental and regional quantitative evaluation.

#### Functional Maps Available

The key perfusion parameters that CT Perfusion 4D package offers (measured in their appropriate physiological units) are:

1. regional Blood Volume (BV; ml/100g).  
Characterizes the amount of blood present in a local region; this is affected by the size of the blood vessels and the number of open capillaries.
2. Regional Blood Flow (BF; ml/min/100 g).  
It is a measure of the volume of blood that flows through a local tissue region within a period of time.
3. Regional Mean Transit Time (MTT; sec).  
Mean Transit time is the time that blood requires to traverse the tissue circulation from arterial inlet(s) to the venous outlet(s). Mathematically:  
 $rMTT = rCBV / rCBF$  (Cerebral), or

**CONFIDENTIAL**

**SECRET DE AFACERI**

$$rMTT = rTBV / rTBF \text{ (Tissue)}$$

#### 4. Contrast Arrival Delay (IRF T0)

It is the delay of arrival of contrast to any given location from its arrival at the artery ROI.

#### 5. Transit Time to IRF Peak (Tmax; sec)

It is the transit time to the peak of the Impulse Residue Function (IRF).

Mathematically:

$$T_{max} = IRF\ T0 + MTT/2$$

#### 6. Capillary Permeability Surface Area Product (PS; ml/min/100 g)

It is the unidirectional rate of transfer of blood solutes (contrast media) through the capillary endothelium into the interstitial space. PS is given in the same units as blood flow, ml/min/100g. Permeability Surface maps are intended to be used with lengthened acquisition times, typically 60 seconds are greater.

#### 7. Base Image corresponds to an average of all pre-enhancement images. It is used in baseline subtraction, and may also be as a reference image.

#### 8. Average Image corresponds to the as the average of all images at a given slice location.

#### 5. Overlay of a functional map onto a reference Image. Creation of a fused display of the colored parametric image and the selected CT image (gray scale or another colored parametric image).

#### 6. Display of functional data.

- Display of the statistical information for all ROIs in Summary Table.

#### 7. Saving and printing of processed data.

- All images, graph data, and functional data may be saved and printed to film or color paper printer.
- Customized protocol may be saved as a user preference for next sessions.
- If the session needs to be interrupted, the protocol state may be saved for further processing.

### Summary of Operations

#### 1. Selection and loading of image data set.

- Selection of dynamic\* data set or registered dynamic data.
- Loading of the selected series, including all slice locations and time points.

#### 2. Definition of cursor and ROI size and location.

- Consistent ROI definition process as in the AW viewer.
- Interrogation of time-intensity information via either a cursor or a region-of-interest (ROI).
- Automatic or semi-automatic selection of artery and vein (where applicable).

#### 3. Display and plotting of time-intensity information.

- Display and plotting of time-intensity information based on the algorithm chosen.

#### 4. Parametric image creation.

- Creation of functional maps; Image intensity values in the functional map represent the result of a calculation process, which generates a single-valued parameter at each pixel location in the original data set.
- Display of parametric images in color or gray scale.

SECRET DE AFACERI

CONFIDENTIAL

GE imagination at work

© 2017 General Electric Company.  
All rights reserved. Data subject to change.  
GE and GE Monogram are trademarks of General Electric Company.  
\* Trademark of General Electric Company  
Doc1996217.

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SECRET DE AFACERI

Visit us:

[http://www3.gehealthcare.com/en/products/categories/advanced\\_visualization/applications/ct\\_perfusion\\_4d\\_multi-organs](http://www3.gehealthcare.com/en/products/categories/advanced_visualization/applications/ct_perfusion_4d_multi-organs)

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SECRET DE AFACERI

CONFIDENTIAL

© 2012 General Electric Company  
All rights reserved. Design subject to change.  
GE and GE Monogram are trademarks of General Electric Company.  
© 2012 General Electric Company. All rights reserved. Data subject to change.  
GE and GE Monogram are trademarks of General Electric Company.  
\* Trademark of General Electric Company.

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
GE imagination at work

# Stroke VCAR

Part of a comprehensive acute Neuro workup—Semi-automated segmentation and assessment of Hematoma and Aneurysms in CT images.

## Clinical and Technical Background

Stroke and trauma to the head may result in intra cerebral or subdural hematomas. Patients under suspicion of such internal bleeding are triaged using CT imaging. It is important to determine both the presence as well as the amount of such bleeding to properly triage a patient for treatment or follow-up. Bleeding in the brain may also be caused by the rupture of aneurysms. It is important to detect the presence of the aneurysm and assess its potential hazard.

## Overview

**Stroke VCAR (Volume Computer Assisted Reading)** gives you a complete reading workflow solution for a comprehensive and robust analysis of hematoma and aneurysms. Hematoma assessment is done using semi-automated segmentation tools married with innovative interactive editing capability in the form of SmartMesh. Aneurysm assessment is done through an innovative user guided aneurysm segmentation and visualization. The program lets you generate a clear, concise clinical summary for sharing with referring physicians.

## Highlights

- The program allows you to quickly segment the hematomas in a non-contrast CT.
- Introducing SmartMesh: an intuitive volumetric editing tool.
- Intelligent user guided aneurysm segmentation.
- Fully integrated with Spectral CT when used with GSI datasets.

**SECRET DE AFACET!**  
**CONFIDE**

Visit us:

[http://www.gehealthcare.com/geo/products/StrokeVicar/advanced\\_visualization/combined\\_tomography\\_images\\_software](http://www.gehealthcare.com/geo/products/StrokeVicar/advanced_visualization/combined_tomography_images_software)

## Features

- Semi-automated segmentation and sizing of hematoma.
- Introducing SmartMesh—an interactive volumetric editing tool.
- Comprehensive reporting of segmented hematomas.
- Fully integrated with Spectral CT when using GSI data.
- Track hematoma changes both visually and quantitatively with longitudinal exam workflow.
- Simple user driven aneurysm detection workflow.
- Easy and intuitive workflow for aneurysm segmentation. Initial segmentation done with four guided clicks.
- A selection of application specific tools ensures ease of use.
- Create layouts to match your reading style.

## System Requirements

### Minimum platform release:

- AW VolumeShare 7 and AW Server 3.2 or later.
- Color monitors
- Single or Dual Display monitors for AW Server in landscape or portrait orientations.

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## Indications for Use

Stroke VCAR is a CT image analysis software package that allows for the analysis and visualization of Brain CT data derived from DICOM 3.0 compliant CT scans. Stroke VCAR is designed for the purpose of segmenting and assessing hemorrhages in the brain using automated tools on non contrast CT exams. Additionally Stroke VCAR provides a set of workflow tools for the segmentation and visualization of aneurysms in the brain from contrast enhanced CT exams. It is intended for use by clinicians to process, review, archive, print and distribute CT studies.

This software will assist the user by providing initial 3D segmentation, measurements and visualization of hemorrhages and aneurysm in the brain. The user has the ability to adjust, review and has to confirm the final segmentation.

GE imagination at work

**SECRET DE AFACER!**  
**CONFIDENTIAL**


© 2016 General Electric Company.  
All rights reserved. Data subject to change.  
GE and GE Monogram are trademarks of General Electric Company.  
\* Trademark of General Electric Company.

## Features Detail

### Hematoma


#### Segmentation

The Stroke VCAR hematoma tool uses an intelligent user guided algorithm for semi-automated segmentation of all types of hematomas in the brain. These segmentations can be converted into volume of interest to measure and assess their density.

#### SmartMesh editing

In cases when the semi-automated segmentation does not fully meet a user's expectation, SmartMesh allows for innovative editing capability to quickly complete the segmentation. A novel feature of the editing tool is to provide users with immediate visual feedback on what is included or excluded from the edition.

#### Reporting

Simple and intuitive reporting enables users to confidently manage and share hematoma findings in a consistent and standardized manner.

#### Follow up

Load repeat scans in Stroke VCAR Hematoma protocol and use its longitudinal analysis capability to track changes, both visually and quantitatively, due to progression or therapy response.


SECRET DE AFACER!  
CONFIDENTIAL


## Aneurysm

### Detection Workflow

A simple and intuitive workflow that helps users with a fast and easy assessment for the detection and treatment planning of aneurysms.


Define visualization planes to better see how the aneurysm is connected to the feeding vessel. Once the segmentation is complete, the user is able to interactively review 3D images of the segmented aneurysm and measurements.

### Segmentation

A simple guided workflow for the segmentation of an aneurysm. Fast initial segmentation is done with four guided clicks.

### Reporting

Simple and intuitive reporting enables users to share images and measurements of aneurysms in a consistent and standardized manner.

### Summary

Stroke VCAR provides workflow benefits using semi-automated hematoma and aneurysm segmentation tools enabling a standardized approach to your facilities review & analysis of these exams.

**CONFIDENTIAL**  
**SECRET DE AFACER!**

GE imagination at work

© 2016 General Electric Company.  
All rights reserved. Data subject to change.  
GE and GE Monogram are trademarks of General Electric Company.  
\* Trademark of General Electric Company.

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
# Bone VCAR

Automated spine labeling with curved reformat generation to improve ease of reading and reporting

## Background

Approximately six percent of the total CT exam procedure volume is dedicated spine imaging. Add to that the number of other routine procedures that include the spine and are referenced within an imaging report you will find a significant number of imaging series generated every day for every scanner where identification and review of the vertebral bodies and their discs are necessary. Bone VCAR will help simplify the reading experience and improve reporting efficiency by identifying and displaying the anatomy in an automated yet flexible implementation.

## Overview

Bone VCAR provides an automated spine labeling application based on a deep learning model. It can identify and label segments of the spine or the entire spine in a matter of seconds. Additionally, let Bone VCAR generate the curved reformats of the spine along with the perpendicular oblique views to easily see the true cross section of the vertebral bodies and disc spaces. Utilize these time saving features to assist in your reading experience. Bone VCAR can also improve the reporting efficiency for identifying and remembering key landmarks that are included in the report dictation.

## Highlights

- Automated spine labeling
- Automated curved reformats with oblique views parallel to the bodies and disc space
- Accessible for any exam type; trauma, oncology, dedicated spine, general imaging
- Accessible with all workflows

**SECRET DE AFACERI**  
CONFIDENTIAL

## Features

- Automated vertebrae labeling in 5 seconds or less for volumes of 300mm or less \*
- Deep learning application trained on global datasets acquired with a broad range of acquisition parameters achieving >90% labeling accuracy
- Automated generation of a 3D trace to generate curved reformat views
- Automated generation of oblique views perpendicular to vertebral bodies and disc spaces
- Compatible with axial source images as well as previously generated reformatted series
- Manually insert, edit or delete vertebral body labels as necessary
- Compatible with data acquisitions for dedicated spine work, trauma, oncology and routine imaging
- Load multiple series and vertebrae labeling will propagate to the different volumes
- Easily review multiple reconstructions with simultaneous display of different kernels

- Access from the Volume Viewer tool bar to easily apply during your standard read at any time
- Compatible with all DICOM CT data that meets image input requirements as described in the user guide

## System Requirements

### Minimum platform release:

- AW Workstations AW 4.7 Ext. 12 with VV14.0 Ext 6 or higher \*Timing specifications with Z440 platform configuration and normal anatomy
- AW Server 3.2 or higher
- Color Monitors
- Single or Dual Display Monitors for AW Server in landscape or portrait orientations

(recommended monitor resolution is up to dual 2MP (1600 x 1200) or a single 3MP (1536 x 2048))

## Regulatory Compliance

This product complies with the European Council Directive 93/42/EEC Medical Device Directive as amended by European Council Directive 2007/47/EC. This product or its features may not be available in some other countries or regions. Please contact your sales associate.

**Rx Only**

## Indications for Use

Bone VCAR is a post processing application for use in the analysis of CT images. The software is intended to support clinicians in the review of images that include the spine by providing tools to label the spine and optimize the display of anatomy within the CT image.

Bone VCAR is designed to support the clinician in visualization of the spine, by providing initial identification of vertebrae to assist in report dictation.

The software also assists the user by providing optimized display settings for easier identification of anatomy to facilitate fast image review and reporting of findings. Bone VCAR may be used for multiple care areas and is not specific to any disease state. It can be utilized during the review of various types of exams including trauma, oncology, and routine body.

CONFIDENTIAL

SECRET DE AFACERI

GE imagination at work

© 2019 General Electric Company.  
All rights reserved. Data subject to change.  
GE and GE Monogram are trademarks of General Electric Company.  
\* Trademark of General Electric Company.  
DOC2185114

## Features Detail

Bone VCAR is designed to support the clinician in easy visualization and identification of vertebral bodies. A refined layout including reformats through the vertebral bodies and disc spaces enables fast review of the anatomy to aid in the reading and reporting experience.

### Integration into VV


- Bone VCAR is integrated into the core toolset of Volume Viewer which makes it available to use in any reading review style desired.

- Compatible with reformatted series, Volume Rendered displays as well as GSI monochromatic datasets
- View labeled anatomy on two different reconstructions at the same time

### Automatic labeling

- Bone VCAR, when activated, will automatically label all the vertebral bodies within the volume for slice thickness & spacing of 2.5 or less

- Deep learning design, built with diverse database of scanner types, variable reconstruction parameters, scan ranges, contrast and non-contrast exams
- Easily edit, when needed, with auto propagation above and below the edited location
- Manually insert vertebrae and add label names to account for anatomical variations
- Load multiple series from the same acquisition and labels will propagate across the volumes

### Curved Reformats


- Curved reformats are automatically created and displayed in a review layout to allow visualization of standard and reformatted views simultaneously

- Oblique view perpendicular to the centerline to show easily the true cross section of the vertebral body and the disc space


- Easily edit the centerline to ensure accurate placement in the most challenging cases

SECRET DE AFACERI

CONFIDENTIAL

- Streamline the reading workflow by loading multiple series for simultaneous display with labels and curved reformat and sagittal views in a single layout

- Create batch reformats with the vertebrae labels and export to PACs to provide easy identification of vertebra bodies for expanded clinical team

### Summary

Bone VCAR is designed to provide an efficient and reliable method for labeling vertebral bodies and generating true sagittal, coronal and axial views via a 3D centerline detection through the center of the vertebral bodies. It is available for use with any DICOM CT datasets that meet the input requirements as defined in the user guide.


CONFIDENTIAL

SECRET DE AFACER!

GE imagination at work

© 2019 General Electric Company.  
All rights reserved. Data subject to change.  
GE and GE Monogram are trademarks of General Electric Company.  
\* Trademark of General Electric Company.  
DOC2185114

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
## Colon VCAR EC

A comprehensive reading workflow solution for colonic lesion detection.


Colorectal cancer is a major health concern for both men and women. Recent studies conclude that early detection of colonic polyps can reduce the mortality rate of colon cancer. You need a protocol that lets you quickly and accurately detect and report results of CT colon scans; a program that optimizes your workflow and improves your reading productivity.

### Overview

Colon VCAR EC gives you a comprehensive reading workflow solution for detecting colonic lesions with flexibility and efficient performance. With it you can visualize anatomy that would otherwise be obscured by tagged fluid. Plus digital contrast agent filter highlights potentially suspicious regions.

### What's new

- The program allows you to read and problem solve in 2D, 3D, or 360 degree dissection views.
- Prone/Supine polyp linking provides quick navigation between datasets.
- Digital Contrast Agent (DCA) provides automatic visualization of shapes characteristic of polyps.
- Electronic Cleansing (EC) subtracts tagged stool and fluid, making it easier to identify lesions.

**SECRET DE AFACERI**

**CONFIDENTIAL**


Visit us:

[www.gehealthcare.com/ow/applications/](http://www.gehealthcare.com/ow/applications/)

## Features

- Prone and supine synchronized reviewing for problem solving.
- You synchronize prone and supine navigation based on polyp bookmark location and distance from center line.
- Provides you with quick, easy to use polyp and center line measurement tools for complete reporting.
- You can leverage any combination of synchronized 2D, 3D, and 360 views.
- Dissection views make primary reading and problem solving more efficient.

- You can vary navigation speed.
- Lets you track entire colon from rectum to cecum for thorough, precise exams.
- You are provided with unrestricted viewing angles of all colon surfaces.
- You can apply color to bookmarked lesions for better visualization.
- The program allows you to assess luminal, mural and extracolonic anatomy and pelvic structures.

- A selection of application-specific tools ensures ease of use and thorough exams.
- Your workflow benefits from fast, interactive electronic patient reporting.

## System Requirements

- AW Workstation
- AW Server

## Indications for Use

Colon VCAR EC is a CT, non-invasive, image analysis software package that allows the visualization of 2D, 3D and dissected medical image data of the colon derived from DICOM 3.0 compliant CT scans. ColonVCAR is designed to aid the physician in evaluating the lumen and internal wall of the colon to confirm the presence or absence of colonic lesions (e.g. polyps). It provides functionality for 2D/3D rendering, bookmarking of suspected lesions, synchronized viewing of the 2D, 3D and 360 dissection views for data sets acquired in any position, and an object oriented endoluminal display. In comparison to Colonoscopy, this tool has an advantage of depth penetration due to its 3D presentation capability.

It is intended for use by clinicians to process, render, review, archive, print and distribute colon image studies.

The Colon VCAR DCA (Digital Contrast Agent) module is an automated highlight feature for the visual identification of spherical structures within the colon and is intended to be used as concurrent reading device. Digital Contrast Agent is a 3D filter that produces images that highlight spherical anatomical regions, such as polyps, and/or stool. Colon VCAR uses color to display these highlighted spheres. Images are made available to the physician to aid in characterization of potential polyps and thus, the patient management care decision process.

## Regulatory Compliance

This product complies with the European CE marking regulation following Medical Devices Directive: Directive 93/42/EEC.

This product has not been approved for sale in the United States by the United States Food and Drug Administration (FDA).

GE imagination at work

**CONFIDENTIAL**  
**SECRET DE AFACERI**

© 2015 General Electric Company.  
All rights reserved. Data subject to change.  
GE and GE Monogram are trademarks of General Electric Company.  
\* Trademark of General Electric Company.

## Features Detail

### Electronic Cleansing


Electronic Cleansing, accessible in both 2D and 3D review, subtracts tagged stool and fluid from the generated data making it easier to identify lesions that may be visible in only one view.

Tagged material is identified based on Hounsfield units. Processing tasks such as Auto Segmentation, 3D flythrough, and Dissection also take advantage of tagged material subtraction.

Auto Segmentation shows the tagged material while 3D flythrough enables quick anatomy visualization without the tagged material. With this method, the tagged material is removed yet shown in a different color at an adjustable opacity level.

Similarly, the tagged material can be removed from the dissection view as seen in the images below.

### User Interface

The Colon VCAR EC user interface has been optimized to facilitate fast, accurate examination of the colon's full extent. It provides seamless integration between automatic and manual interactions for a level of flexibility designed to reduce reading times.

### Simultaneous, Indexed Viewing

The program can display both prone and supine

views, which can be manually synchronized for simultaneous and indexed viewing. During review, you can "Lock to Target" with a single button click to quickly problem solve during the primary read. And you are able to bookmark locations of suspected polyps.

### 360-degree Dissection

This powerful tool provides a 360-degree view of the entire internal large bowel simultaneously displayed in supine and prone positions. This ability speeds viewing and analysis.

This feature lets you leverage all available views: 2D, 3D, and 360-degree dissection during the primary read. For more traditional reading methods, a 2D, 3D, or Dissection view can be laid out and reserved for primary reading with reference to other views for problem solving.

### Joystick Navigation

Colon VCAR EC lets you adjust the speed and direction of a flythrough simply by using the mouse. Flythroughs can also be set on autopilot for hands-off reviewing. Additionally, a single mouse click lets you lock onto a region of interest and fly around the target anomaly.

### Computer Aided Reading

Colon VCAR EC's Contrast Agent (DCA), a 3D filter, highlights spherical anatomy such as polyps in color. The resulting images facilitate the characterization of potential polyps.

### Virtual Biopsy

Problem solving while 3D reading is facilitated by the Virtual biopsy tool, which allows you to view areas in the navigator in color-coded or black/white scales. Different densities, determined by Hounsfield units, can be displayed using the color scale.

SECRET DE AFACERI


CONFIDENTIAL

GE imagination at work

Colon VCAR EC -1

### Auto-Segmentation


Colon VCAR EC quickly performs the time-intensive task of colon segmentation for both prone and supine data sets. All you do is select the protocol.

The program automatically defines a luminal centerline in both prone and supine positions. This can be applied to a high percentage of cases, depending on the level of insufflation and collapse.

Semi-automatic and manual segmentation modes are also available. Semi-automatic segmentation provides complementary, and often necessary, assistance in cases when the colon is collapsed and operator intervention is required. This mode is typically used to connect colon segments around a collapse to form a single centerline. Manual segmentation is typically used to extend the centerline to the end of the cecum for completeness. Manual editing is also available to refine automatic tracking if needed.

### Dissection (centerline) Contouring Verification

Using a color mapping technique, Dissection Contouring Verification outlines colon anatomy included in the centerline trace. You can review verification in any orthogonal plane.

### Automatic High Resolution Flythrough

Real time 512-matrix resolution automatic flythrough enhances navigational image quality, and provides operational flexibility, functionality, and simplicity. Anatomy remains sharp with minimal pixilation. Steering remains centered as you navigate through the colon. You can specify your navigation speed, and 3D navigation is correlated with 2D and 360-degree Dissection views for rapid synchronized referencing in primary or secondary reading.

512 matrix resolution

### Save State

This feature gives you the ability to save any centerline edits, bookmarks, and colored polyps for later reloading. It is also convenient for quick recall during consultation or exam review.

### Prone/Supine Polyp Linking


After bookmarking a polyp in prone and supine datasets, the program provides you with the capability to match and label them for easy tracking. This gives you the ability to read registered prone/supine datasets for a one pass read.

### Synchronized Prone/Supine 3D Dissection View Navigation

Once prone/supine polyp linkage is made, the program co-registers and links the datasets and polyps, letting you visualize opposing anatomy for a synchronized review.

### Small Bowel Extraction

This feature quickly segments the small bowel, giving you an unobstructed view of the large intestine.

### Analysis Mode

The program's Analysis Mode provides a variety of tools to help streamline your workflow:

- Lock Mode lets you target and possible lesion and stay locked on it at all times.
- Synchronized display can be customized to show two oblique images, axial, sagittal, coronal, volume rendered and/or 3D views within the colon in a Lock mode. Leveraging multiple synchronized planes and views in a problem-solving mode may enhance your reading productivity.

SECRET DE AFACERI


CONFIDENTIAL

GE imagination at work

## Color Centerline Tracking

With Colon VCAR EC you can display a real-time 3D model of the colon with a defined centerline to quickly orient the location of the area you are viewing to the overall colon. This synchronization is always linked and available for quick reference.

## Patient Report\*

Customizable, interactive reporting capability gives greater structure and flexibility in generating your reports. You can choose how many images to display per sheet, place your facility's logo, apply physician signatures, list patient demographics, and detail your comments and impressions. Export your report to CD, HTTP, or print it to paper.

The Quick Report panel gives you instant access to report input fields and provides access to automatic bookmark annotation on the colon map.

## Measurement Tools


With Colon VCAR EC, you can measure polyp distance from distal rectum to polyp location in millimeters, perform 2D measurements, and obtain statistics of ROI in Hounsfield units.

## Archive

You can define the number of images per polyp you want to archive to either a networked location or to magnetic optical disk.

## Customizable Review and Analysis Layouts

Create layouts that match your reading style and save them as protocols for future recall. Viewports can be shaped into a variety of square or rectangular configurations.

Five of the most popular protocol layouts are pre-built for quick loading. Choose from several dissection configurations:

- 360-dissection gives you convenience and productivity benefits.
- 120-dissection lets you read around distortions if necessary.

- Use the 2D read configuration if the exam has poor insufflation and multiple collapses. 2D can also be used for extracolonic review.

- Use the 3-series layout to load a prone, supine, and decubitus series in parallel to obtain a correlated view.

- You can also customize new layouts to suit your own reading style.
- If you work on a one-monitor Advantage Workstation, a single monitor display is also available.

## Movie Generator

You can generate movie clips of the area of interest simply by defining the start and end locations with a mouse click. A cross-sectional movie will be generated. You can also generate lock-to-target movies. Movies can be saved as mpeg, avi, HTTP and saved to CD.

SECRET DE AFACERI

GE imagination at work

Colon VCAR EC -3

### Entire Colon Movie

Generate a fly-through movie of the entire colon with one button click. Export the result as an mpeg or as screen captures to a PACS workstation for review or consultation.

### Additional Features

- Automatic center line tracking
- SmartCursor™ for easier navigation
- Synchronized reformatted views
- Multiple correlations of anatomy with CT colonography
- Exploration of suspicious areas in primary 2D reformatted sections
- Primary 3D reformatted sections. Obliques, Navigator, Quick Virtual Dissection, and axial views are all correlated to the same anatomy of interest.
- Black/white or color display

### Summary

Colon VCAR EC helps make colon evaluation a less traumatic procedure for patients. Preparation for the exam is less severe; the procedure is less time intensive and less expensive than optical colonoscopy. Moreover, the tools and features this program puts into your hands will streamline your workflow and potentially enable greater diagnostic confidence.

### Indications for Use

Colon VCAR EC is a CT, non-invasive, image analysis software package that allows the visualization of 2D, 3D and dissected medical image data of the colon derived from DICOM 3.0 compliant CT scans. ColonVCAR is designed to aid the physician in evaluating the lumen and internal wall of the colon to confirm the presence or absence of colonic lesions (e.g. polyps). It provides functionality for 2D/3D rendering, bookmarking of suspected lesions, synchronized viewing of the 2D, 3D and 360 dissection views for data sets acquired in any position, and an object oriented endoluminal display. In comparison to Colonoscopy, this tool has an advantage of depth penetration due to its 3D presentation capability. It is intended for use by clinicians to process, render, review, archive, print and distribute colon image studies.

The Colon VCAR DCA (Digital Contrast Agent) module is an automated highlight feature for the visual identification of spherical structures within the colon and is intended to be used as concurrent reading device. Digital Contrast Agent is a 3D filter that produces images that highlight spherical anatomical regions, such as polyps, and/or stool. Colon VCAR uses color to display these highlighted spheres. Images are made available to the physician to aid in characterization of potential polyps and thus, the patient management care decision process.

\*Report feature is only available on AW Workstation


GE imagination at work

**CONFIDENTIAL**  
**SECRET DE AFACERI**

© 2015 General Electric Company.  
All rights reserved. Data subject to change.  
GE and GE Monogram are trademarks of General Electric Company.  
\* Trademark of General Electric Company.

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
**SECRET DE AFACERI**  
**CONFIDENTIAL**

## Lung VCAR

For lung nodule analysis - automatic visualization, measurement, reporting and follow-up

### Clinical and Technical Background

CT scanning is utilized as the exam of choice to visualize and analyze complex lung pathology. The detection of pulmonary nodules and assessment of their evolution with CT are of major importance in chest imaging.

### Overview

Lung VCAR brings efficient CT pulmonary nodule assessment and diagnosis. The innovative Digital Contrast Agent (DCA) feature automatically visualizes lung nodules to help you confirm the presence or absence of suspicious lesions from 2 to 12 mm in size. Lung VCAR allows automated follow-up for lesion matching by the registration of two or more datasets, automatic lesion classification.

### Highlights

- Synchronized 2D, Digital Contrast Agent (DCA) and segmentation analysis.
- Automatic nodule visualization.
- Automatic nodule analysis (volume, doubling time, % growth).
- Automatic follow-up
- Reporting workflow

SECRET DE AFACERI  
CONFIDENTIAL

Visit us:

<https://www.gehealthcare.com/products/advanced-visualization/all-applications/lung-vcar>

## Features

### Review

- Ability to synchronize multiple images for nodule comparison.
- Ability to review single or multiple exams and compare axial, sagittal, oblique, coronal, and volume-rendered images.
- Automatically propagates previous exam bookmarks to current exam.
- Automatically segments both right and left lung to reduce visual distraction.
- Digital Contrast Agent (DCA) automatically highlights spherical shapes to enhance visualization of suspicious nodules.

### Analysis

- Performs automatic segmentation of all nodule types.
- Provides automatic nodule analysis, including % growth, doubling time
- Automatic Nodule Contour to verify pixels within the volume

### Indications for Use

Lung VCAR/AdvantageALA is intended to provide an optimized non-invasive application to measure abnormalities in the lung (for example, nodules, lesions, etc.) from a set of computed tomography (CT) images.

The software is designed to support the physician in confirming the presence or absence of physician identified lung lesions (e.g. nodules). The software allows measurement of volume over time using a consistent standardized measurement protocol, thus providing an estimation of the volume doubling time. The Lung VCAR/AdvantageALA software allows analysis and displays statistics for nodule characterization all the different nodule types. Lung VCAR optional Digital Contrast Agent (DCA) module is an automated highlight feature for the visual identification of possible lesions. Digital Contrast Agent (DCA) is a 3D filter that produces images that highlight spherical (S) or cylindrical (C) anatomical regions, such as nodules, cysts, scars, and vessels. Images are made available to the physician to aid in characterization of suspicious nodules and thus, the patient management care decision process. Lung VCAR/AdvantageALA provides the physician with additional information, meant to complement diagnosis based on classical techniques.

### System Requirements

Lung VCAR/AdvantageALA can be installed on AW Workstations AW 4.7 Ext. 14 or later and AW Server 3.2 Ext 3.2 or higher.

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SECRET DE AFACERI

GE imagination at work

© 2021 General Electric Company.  
All rights reserved. Data subject to change.  
GE and GE Monogram are trademarks of General Electric Company.  
\* Trademark of General Electric Company.  
DOC0751472

## Features Detail

### Scanning

Lung VCAR sensitivity and segmentation performance is a function of slice thickness and dose. Lung VCAR was designed to be compatible with slice thicknesses ranging from 0.625 – 1.25 mm with a dose range of >40.

### Review Correlated Workflow

The Lung VCAR user interface is designed to organize routine chest radiology reading and make it more efficient. Viewports and displayed anatomy are synchronized and correlated, facilitating reading by making review and problem solving more transparent.


Lung VCAR gives you the ability to:

- Synchronize multiple images for nodule comparison.
- Review single or multiple exams and compare axial, sagittal, coronal, oblique, and volume rendered images.
- Automatically propagate and display previous exam bookmarks on the current exam when two exams are loaded.

Lung VCAR provides you with tools to help improve your workflow and analysis precision.

### Lung Segmentation

Lung VCAR automatically segments both the right and left lungs. Segmentation reduces visual distractions associated with irrelevant anatomy during lung nodule evaluation. Total lung volume is captured in the patient report.

**Automatic Follow-up** synchronizes and displays current and prior patient exams and their corresponding bookmarks for quick and easy comparison review.

**Digital Contrast Agent (DCA)** automatically highlights spherical shapes to enhance visualization of suspicious nodules.

**DCA Color Mapping** highlights nodules in your choice of red, or for a colorblind-friendly display, orange.

**DCA Object Size (0-10mm)** gives you full control over visualization sensitivity to a desired nodule size and the number of false positives displayed.

**DCA Shape Choice (Sphere / Cylindrical)** lets you select from two different filters to help visualize specific shapes and assist in problem solving:

- Cylindrical for vessels
- Spheres for nodules

**DCA Toggle** enables you to toggle the DCA-highlighted nodules on and off.

**X-Reference / Correlation Bar** provides you with a quick reference for localizing a nodule's global location when reading 2D axial images. The x-reference bar is synchronized with both review and analysis layouts for immediate screen refresh when deposited at any location or on any anatomy.

CONFIDENTIAL

SECRET DE AFACEDU

## Bookmark Management

Lung VCAR provides you with custom-designed (select dot or small box bookmark format) bookmark management capability linked to review, paging and Digital Contrast Agent for a more organized nodule tracking and analysis workflow.

**Interface** lets you dynamically and easily view, page and deposit bookmarks within an exam. It provides:

- On/Off toggle capability for DCA-highlighted nodules
- Smooth, contiguous paging through large data sets
- Thick slab MIP paging
- Bookmark paging (next / prior and add / delete)

**Automatic Bookmark Propagation** lets you automatically propagate bookmarks from previous to current exam or current to previous exam.

**Deferred bookmarking** lets you queue bookmarked nodules for later segmentation (operator initiated segmentation).


**Immediate bookmarking** automatically segments the nodule and displays the analysis review layouts immediately after you place a bookmark

## Layouts

**Layout Preferences** let you create and save a layout to match your reading style. Layout elements such as viewport, image type in each viewport, slice thickness, window level/width, magnification, rendering mode, can be easily customized.


## Analysis

Analysis mode allows you to maximize use of a second monitor by displaying a segmented nodule of interest to an orientation useful for quick problem solving and analysis.

**Automatic segmentation of all nodule types** provides automated segmentation and sizing of all nodules types. This feature can:

- Display 3D volume rendered nodule view.
- On follow-up exams, display doubling time and percentage of nodule growth.
- Compare nodules from multiple exams

**Shutter view** provides you with a quick Volume Rendered framed view of a point of interest in relationship to surrounding structures for fast visualization and problem solving.

CONFIDENTIAL  
SECRET DE AFACER

**Manual or automatic nodule characterization** can automatically characterize and label nodule type and apply the necessary segmentation algorithm for accurate isolation. You can at any time manually override this feature and apply a different segmentation algorithm and labeling for any nodule.

**Automatic Nodule Analysis** provides the following nodule information:

- % growth
- Doubling time
- Volumes (solid component, non-solid component, total volume)
- Maximum distance in 3D (long axis)
- 2D measurements
- Hounsfield units (Min, Max, Avg)
- Scan interval

**Paging enhancements** provide a special scrolling tool that lets you page through large data sets while the image displays seamlessly.

**Automatic nodule contour** capability lets you verify pixels within the volume. It displays segmentation nodule contour on the 2D image for quick segmentation validation.

## Patient Report

DICOM SR allows a reporting workflow to present and archive images and measurements.

Report feature only available on AW Workstation.

## Summary

With automated nodule visualization and localization, increased detection sensitivity, accurate and repeatable segmentations analysis, Lung VCAR could streamline your reading workflow, enhance your efficiency, and increase your diagnostic confidence.

CONFIDENTIAL


SECRET DE AFACER!

GE imagination at work

© 2021 General Electric Company  
All rights reserved. Data subject to change.  
GE and GE Monogram are trademarks of General Electric Company.  
\* Trademark of General Electric Company  
DOC0751472

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GE, the GE Monogram, Centricity și Snapshot Freeze sunt trademarks of General Electric Company