ORGENTEC Diagnostika GmbH


Carl-Zeiss-Straße 49-51

55129 Mainz - Germany

Phone: +49 (0) 61 31 / 92 58-0

Fax: +49 (0) 61 31 / 92 58-58

Internet: www.orgentec.com

ORG 512 Anti-Sci-70

INTENDED PURPOSE

Anti-Sci-70 is an ELISA test system for the quantitative measurement of IgG class autoantibodies against Sci-70 in human serum or plasma. This product is intended for professional in vitro diagnostic use only.

Antibodies against Sci-70 (DNA topoisomerase I) are an accepted marker for progressive systemic sclerosis. They contribute to the differential diagnosis of sclerosis. Evaluation of a test result should always take into account all clinical and laboratory diagnostic findings.

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


Ready to use

PRINCIPLE OF THE TEST

Highly purified Sci-70 is bound to microwells.

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CONTENTS OF THE KIT

ORG 512		96	Sufficient for 96 determinations
MICROPLATE	1		One divisible microplate consisting of 12 modules of 8 wells each. Ready to use. Color code on module
CALIBRATOR A	1x 1.5 ml		Calibrator A 0 U/ml, containing serum/buffer matrix (PBS, BSA, detergent, NaN3 0.09%), yellow. Ready to use.
CALIBRATOR B	1x 1.5 ml		Calibrator B 12.5 U/ml, containing Scl-70 antibodies in a serum/buffer matrix (PBS, BSA, detergent, NaN3 0.09%), yellow. Ready to use.
CALIBRATOR C	1x 1.5 ml		Calibrator C 25 U/ml, containing Scl-70 antibodies in a serum/buffer matrix (PBS, BSA, detergent, NaN3 0.09%), yellow. Ready to use.
CALIBRATOR D	1x 1.5 ml		Calibrator D 50 U/ml, containing Scl-70 antibodies in a serum/buffer matrix (PBS, BSA, detergent, NaN3 0.09%), yellow. Ready to use.
CALIBRATOR E	1x 1.5 ml		Calibrator E 100 U/ml, containing Scl-70 antibodies in a serum/buffer matrix (PBS, BSA, NaN3 0.09%), yellow. Ready to use.
CALIBRATOR F	1x 1.5 ml		Calibrator F 200 U/ml, containing Scl-70 antibodies in a serum/buffer matrix (PBS, BSA, detergent, NaN3 0.09%), yellow. Ready to use.
CONTROL +	1x 1.5 ml		Control positive, containing Scl-70 antibodies in a serum/buffer matrix (PBS, BSA, detergent, NaN3 0.09%), yellow. Ready to use. The concentration is specified on the certificate of analysis.
CONTROL -	1x 1.5 ml		Control negative, containing Scl-70 antibodies in a serum/buffer matrix (PBS, BSA, detergent, NaN3 0.09%), yellow. Ready to use. The concentration is specified on the certificate of analysis.
DILUENT	20 ml		Sample Buffer P, containing PBS, BSA, detergent, preservative sodium azide 0.09%, yellow, concentrate (5 x).
CONJUGATE	15 ml		Enzyme Conjugate containing anti-human IgG antibodies, HRP labelled; PBS, BSA, detergent, preservative PROCLIN 0.05%, light red. Ready to use.
TMB	15 ml		TMB Substrate; containing 3,3', 5,5'- Tetramethylbenzidin, colorless. Ready to use.
STOP	15 ml		Stop solution; contains acid. Ready to use.
WASH	20 ml		Wash Buffer, containing Tris, detergent, preservative sodium azide 0.09%; 50 x conc.
	1		Certificate of Analysis

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CALCULATION OF RESULTS

For quantitative results plot the optical density of each calibrator versus the calibrator concentration to create a calibration curve. The concentration of patient samples may then be estimated from the calibration curve by interpolation.
Using data reduction software a 4-Parameter-Fit with lin-log coordinates for optical density and concentration is the data reduction method of choice.

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Reproducibility

Intra-assay precision: Coefficient of variation (CV) was calculated for each of three samples from the results of 24 determinations in a single run. Results for precision-within-assay are shown in the table below.
Inter-assay precision: Coefficient of variation (CV) was calculated for each of three samples from the results of 6 determinations in 5 different runs. Results for run-to-run precision are shown in the table below.

Intra-Assay		
Sample	Mean U/ml	CV %
1	45.7	4.0
2	90.4	3.2
3	184.1	3.4

Inter-Assay		
Sample	Mean U/ml	CV %
1	41.1	2.8
2	89.9	2.8
3	157.4	2.3

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		Clinical Diagnosis		
		POS	NEG	
ORG 512	POS	19	1	25 100 125
	NEG	6	99	
Sensitivity:		76.0	%	
Specificity:		99.0	%	
Overall agreement:		94.4	%	

LIMITATIONS OF THE PROCEDURE

This assay is a diagnostic aid. A definite clinical diagnosis should not be based on the results of a single test, but

should be made by the physician after all clinical and laboratory findings have been evaluated concerning the entire clinical picture of the patient. Also every decision for therapy should be taken individually.

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REFERENCES

- Alba P, Bento L, Cuadrado MJ, Karim Y, Tungekar MF, Abbs I et al. Anti-dsDNA, anti-Sm antibodies, and the lupus anticoagulant: significant factors associated with lupus nephritis. *Ann Rheum Dis* 2003; 62(6):556-560.
- Antico A, Platzgummer S, Bassetti D, Bizzaro N, Tozzoli R, Villalta D. Diagnosing systemic lupus erythematosus: new-generation immunoassays for measurement of anti-dsDNA antibodies are an effective alternative to the Farr technique and the Crithidia luciliae immunofluorescence test. *Lupus* 2010; 19(8):906-912.
- Brouwer R, Hengstman GJ, Vree EW, Ehrfeld H, Bozic B, Ghirardello A et al. Autoantibody profiles in the sera of European patients with myositis. *Ann Rheum Dis* 2001; 60(2):116-123.
- Castro C, Gourley M. Diagnostic testing and interpretation of tests for autoimmunity. *J Allergy Clin Immunol* 2010; 125(2 Suppl 2):S238-S247.
- Defendenti C, Atzeni F, Spina MF, Grosso S, Cereda A, Guercilena G et al. Clinical and laboratory aspects of Ro/SSA-52 autoantibodies. *Autoimmun Rev* 2011; 10(3):150-154.
- Eriksson C, Kokkonen H, Johansson M, Hallmans G, Wadell G, Rantapaa-Dahlqvist S. Autoantibodies predate the onset of Systemic Lupus Erythematosus in northern Sweden. *Arthritis Research & Therapy* 2011; 13(1):R30.
- Haugbro K, Nossent JC, Winkler T, Figenschau Y, Rekvig OP. Anti-dsDNA antibodies and disease classification in antinuclear antibody positive patients: the role of analytical diversity. *Ann Rheum Dis JID* - 0372355 2004; 63(4):386-394.
- Ippolito A, Wallace DJ, Gladman D, Fortin PR, Urowitz M, Werth V et al. Autoantibodies in systemic lupus erythematosus: comparison of historical and current assessment of seropositivity. *Lupus* 2011; 20(3):250-255.
- Isenberg DA, Manson JJ, Ehrenstein MR, Rahman A. Fifty years of anti-ds DNA antibodies: are we approaching journey's end? *Rheumatology (Oxford)* 2007; 46(7):1052-1056.
- Kattah NH, Kattah MG, Utz PJ. The U1-snRNP complex: structural properties relating to autoimmune pathogenesis in rheumatic diseases. *Immunol Rev* 2010; 233(1):126-145.
- Kumar Y, Bhatia A, Minz RW. Antinuclear antibodies and their detection methods in diagnosis of connective tissue diseases: a journey revisited. *Diagn Pathol* 2009; 4:1.
- Meroni PL, Schur PH. ANA screening: an old test with new recommendations. *Ann Rheum Dis* 2010; 69:1420-1422.
- Petri M, Magder L. Classification criteria for systemic lupus erythematosus: a review. *Lupus* 2004; 13(11):829-837.
- Poole BD, Schneider RI, Guthridge JM, Vette CA, Reichlin M, Harley JB et al. Early targets of nuclear RNP humoral autoimmunity in human systemic lupus erythematosus. *Arthritis Rheum* 2009; 60(3):848-859.
- Putova I, Dostal C, Becvar R. Prevalence of antinucleosome antibodies by enzyme-linked immunosorbent assays in patients with systemic lupus erythematosus and other autoimmune systemic diseases. *Ann N Y Acad Sci* 2007; 1109:275-286.
- Reveille JD. Predictive value of autoantibodies for activity of systemic lupus erythematosus. *Lupus JID* - 9204265 2004; 13(5):290-297.
- Simon JA, Cabiedes J, Ortiz E, Alcocer-Varela J, Sanchez-Guerrero J. Anti-nucleosome antibodies in patients with systemic lupus erythematosus of recent onset. Potential utility as a diagnostic tool and disease activity marker. *Rheumatology (Oxford)* 2004; 43(2):220-224.
- Sinclair D, Saas M, Williams D, Hart M, Goswami R. Can an ELISA replace immunofluorescence for the detection of anti-nuclear antibodies?--The routine use of anti-nuclear antibody screening ELISAs. *Clin Lab* 2007; 53(3-4):183-191.
- Tozzoli R, Bizzaro N, Tonutti E, Villalta D, Bassetti D, Manoni F et al. Guidelines for the laboratory use of autoantibody tests in the diagnosis and monitoring of autoimmune rheumatic diseases. *Am J Clin Pathol* 2002; 117(2):316-324.
- Maidhof W., Hilius O. Lupus: an overview of the disease and management options. *P T* 2012; 37(4):240-9.
- Hahn BH, McMahon MA, Wilkinson A, Wallace WD, Daikh DI, Fitzgerald JD et al. American College of Rheumatology guidelines for screening, treatment, and management of lupus nephritis. *Arthritis Care Res (Hoboken)* 2012; 64(6):797-808.

competent authority of the EU Member State in which the user and/or the patient is established .

Change Control

Former version: *ORG 512_IFU_EN_QM113139_2013-12-16_1.2* Reason for revision: *Introduction electronic IFU on homepage*

Notice to the user (European Union):

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the

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ORG 509 Anti-SS-B

INTENDED PURPOSE

Anti-SS-B is an ELISA test system for the quantitative measurement of IgG class autoantibodies against SS-B (La) in human serum or plasma. This product is intended for professional in vitro diagnostic use only.

Antibodies against SS-B are used for the differential diagnosis of systemic inflammatory autoimmune diseases. Autoantibodies against the SS-B protein are usually found together with anti-SS-A in cases of Sjogren's syndrome. Evaluation of a test result should always take into account all clinical and laboratory diagnostic findings.

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


Ready to use

PRINCIPLE OF THE TEST

Highly purified SS-B is bound to microwells.

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CONTENTS OF THE KIT

ORG 509		96	Sufficient for 96 determinations
MICROPLATE	1		One divisible microplate consisting of 12 modules of 8 wells each. Ready to use. Color code on module
CALIBRATOR A	1x 1.5 ml		Calibrator A 0 U/ml, containing serum/buffer matrix (PBS, BSA, detergent, NaN3 0.09%), yellow. Ready to use.
CALIBRATOR B	1x 1.5 ml		Calibrator B 12.5 U/ml, containing SS-B antibodies in a serum/buffer matrix (PBS, BSA, detergent, NaN3 0.09%), yellow. Ready to use.
CALIBRATOR C	1x 1.5 ml		Calibrator C 25 U/ml, containing SS-B antibodies in a serum/buffer matrix (PBS, BSA, detergent, NaN3 0.09%), yellow. Ready to use.
CALIBRATOR D	1x 1.5 ml		Calibrator D 50 U/ml, containing SS-B antibodies in a serum/buffer matrix (PBS, BSA, detergent, NaN3 0.09%), yellow. Ready to use.
CALIBRATOR E	1x 1.5 ml		Calibrator E 100 U/ml, containing SS-B antibodies in a serum/buffer matrix (PBS, BSA, NaN3 0.09%), yellow. Ready to use.
CALIBRATOR F	1x 1.5 ml		Calibrator F 200 U/ml, containing SS-B antibodies in a serum/buffer matrix (PBS, BSA, detergent, NaN3 0.09%), yellow. Ready to use.
CONTROL +	1x 1.5 ml		Control positive, containing SS-B antibodies in a serum/buffer matrix (PBS, BSA, detergent, NaN3 0.09%), yellow. Ready to use. The concentration is specified on the certificate of analysis.
CONTROL -	1x 1.5 ml		Control negative, containing SS-B antibodies in a serum/buffer matrix (PBS, BSA, detergent, NaN3 0.09%), yellow. Ready to use. The concentration is specified on the certificate of analysis.
DILUENT	20 ml		Sample Buffer P, containing PBS, BSA, detergent, preservative sodium azide 0.09%, yellow, concentrate (5 x).
CONJUGATE	15 ml		Enzyme Conjugate containing anti-human IgG antibodies, HRP labelled; PBS, BSA, detergent, preservative PROCLIN 0.05%, light red. Ready to use.
TMB	15 ml		TMB Substrate; containing 3,3', 5,5'- Tetramethylbenzidin, colorless. Ready to use.
STOP	15 ml		Stop solution; contains acid. Ready to use.
WASH	20 ml		Wash Buffer, containing Tris, detergent, preservative sodium azide 0.09%; 50 x conc.
	1		Certificate of Analysis

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CALCULATION OF RESULTS

For quantitative results plot the optical density of each calibrator versus the calibrator concentration to create a calibration curve. The concentration of patient samples may then be estimated from the calibration curve by interpolation.
Using data reduction software a 4-Parameter-Fit with lin-log coordinates for optical density and concentration is the data reduction method of choice.

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Reproducibility

Intra-assay precision: Coefficient of variation (CV) was calculated for each of three samples from the results of 24 determinations in a single run. Results for precision-within-assay are shown in the table below.
Inter-assay precision: Coefficient of variation (CV) was calculated for each of three samples from the results of 6 determinations in 5 different runs. Results for run-to-run precision are shown in the table below.

Intra-Assay		
Sample	Mean U/ml	CV %
1	28.8	5.6
2	67.1	5.8
3	143.2	5.2

Inter-Assay		
Sample	Mean U/ml	CV %
1	24.5	11.0
2	70.5	6.9
3	157.6	4.1

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		Clinical Diagnosis		
		POS	NEG	
ORG 509	POS	43	5	70 120 190
	NEG	27	115	
Sensitivity:		61.4	%	
Specificity:		95.8	%	
Overall agreement:		83.2	%	

LIMITATIONS OF THE PROCEDURE

This assay is a diagnostic aid. A definite clinical diagnosis should not be based on the results of a single test, but

should be made by the physician after all clinical and laboratory findings have been evaluated concerning the entire clinical picture of the patient. Also every decision for therapy should be taken individually.

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REFERENCES

1. Alba P, Bento L, Cuadrado MJ, Karim Y, Tungekar MF, Abbs I et al. Anti-dsDNA, anti-Sm antibodies, and the lupus anticoagulant: significant factors associated with lupus nephritis. *Ann Rheum Dis* 2003; 62(6):556-560.
2. Antico A, Platzgummer S, Bassetti D, Bizzaro N, Tozzoli R, Villalta D. Diagnosing systemic lupus erythematosus: new-generation immunoassays for measurement of anti-dsDNA antibodies are an effective alternative to the Farr technique and the Crithidia luciliae immunofluorescence test. *Lupus* 2010; 19(8):906-912.
3. Brouwer R, Hengstman GJ, Vree EW, Ehrfeld H, Bozic B, Ghirardello A et al. Autoantibody profiles in the sera of European patients with myositis. *Ann Rheum Dis* 2001; 60(2):116-123.
4. Castro C, Gourley M. Diagnostic testing and interpretation of tests for autoimmunity. *J Allergy Clin Immunol* 2010; 125(2 Suppl 2):S238-S247.
5. Defendenti C, Atzeni F, Spina MF, Grosso S, Cereda A, Guercilena G et al. Clinical and laboratory aspects of Ro/SSA-52 autoantibodies. *Autoimmun Rev* 2011; 10(3):150-154.
6. Eriksson C, Kokkonen H, Johansson M, Hallmans G, Wadell G, Rantapaa-Dahlqvist S. Autoantibodies predate the onset of Systemic Lupus Erythematosus in northern Sweden. *Arthritis Research & Therapy* 2011; 13(1):R30.
7. Haugbro K, Nossent JC, Winkler T, Figenschau Y, Rekvig OP. Anti-dsDNA antibodies and disease classification in antinuclear antibody positive patients: the role of analytical diversity. *Ann Rheum Dis JID* - 0372355 2004; 63(4):386-394.
8. Ippolito A, Wallace DJ, Gladman D, Fortin PR, Urowitz M, Werth V et al. Autoantibodies in systemic lupus erythematosus: comparison of historical and current assessment of seropositivity. *Lupus* 2011; 20(3):250-255.
9. Isenberg DA, Manson JJ, Ehrenstein MR, Rahman A. Fifty years of anti-ds DNA antibodies: are we approaching journey's end? *Rheumatology (Oxford)* 2007; 46(7):1052-1056.
10. Kattah NH, Kattah MG, Utz PJ. The U1-snRNP complex: structural properties relating to autoimmune pathogenesis in rheumatic diseases. *Immunol Rev* 2010; 233(1):126-145.
11. Kumar Y, Bhatia A, Minz RW. Antinuclear antibodies and their detection methods in diagnosis of connective tissue diseases: a journey revisited. *Diagn Pathol* 2009; 4:1.
12. Meroni PL, Schur PH. ANA screening: an old test with new recommendations. *Ann Rheum Dis* 2010; 69:1420-1422.
13. Petri M, Magder L. Classification criteria for systemic lupus erythematosus: a review. *Lupus* 2004; 13(11):829-837.
14. Poole BD, Schneider RI, Guthridge JM, Vette CA, Reichlin M, Harley JB et al. Early targets of nuclear RNP humoral autoimmunity in human systemic lupus erythematosus. *Arthritis Rheum* 2009; 60(3):848-859.
15. Putova I, Dostal C, Becvar R. Prevalence of antinucleosome antibodies by enzyme-linked immunosorbent assays in patients with systemic lupus erythematosus and other autoimmune systemic diseases. *Ann N Y Acad Sci* 2007; 1109:275-286.
16. Reveille JD. Predictive value of autoantibodies for activity of systemic lupus erythematosus. *Lupus JID* - 9204265 2004; 13(5):290-297.
17. Simon JA, Cabiedes J, Ortiz E, Alcocer-Varela J, Sanchez-Guerrero J. Anti-nucleosome antibodies in patients with systemic lupus erythematosus of recent onset. Potential utility as a diagnostic tool and disease activity marker. *Rheumatology (Oxford)* 2004; 43(2):220-224.
18. Sinclair D, Saas M, Williams D, Hart M, Goswami R. Can an ELISA replace immunofluorescence for the detection of anti-nuclear antibodies?--The routine use of anti-nuclear antibody screening ELISAs. *Clin Lab* 2007; 53(3-4):183-191.
19. Tozzoli R, Bizzaro N, Tonutti E, Villalta D, Bassetti D, Manoni F et al. Guidelines for the laboratory use of autoantibody tests in the diagnosis and monitoring of autoimmune rheumatic diseases. *Am J Clin Pathol* 2002; 117(2):316-324.
20. Maidhof W., Hilius O. Lupus: an overview of the disease and management options. *P T* 2012; 37(4):240-9.
21. Hahn BH, McMahon MA, Wilkinson A, Wallace WD, Daikh DI, Fitzgerald JD et al. American College of Rheumatology guidelines for screening, treatment, and management of lupus nephritis. *Arthritis Care Res (Hoboken)* 2012; 64(6):797-808.

Notice to the user (European Union):

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the

competent authority of the EU Member State in which the user and/or the patient is established .

Change Control

Former version: *ORG 509_IFU_EN_QM113136_2013-12-16_1.2*

Reason for revision: *Introduction electronic IFU on homepage*

- ① Pipet **100 µl** calibrator, control or patient sample
→ Incubate for **30 minutes** at room temperature
→ Discard the contents of the wells and wash 3 times with **300 µl** wash solution
- ② Pipet **100 µl** enzyme conjugate
→ Incubate for **15 minutes** at room temperature
→ Discard the contents of the wells and wash 3 times with **300 µl** wash solution
- ③ Pipet **100 µl** substrate solution
→ Incubate for **15 minutes** at room temperature
- ④ Add **100 µl** stop solution
→ Leave untouched for **5 minutes**
→ Read at **450 nm**

ORGENTEC Diagnostika GmbH


Carl-Zeiss-Straße 49-51

55129 Mainz - Germany

Phone: +49 (0) 61 31 / 92 58-0

Fax: +49 (0) 61 31 / 92 58-58

Internet: www.orgentec.com

ORG 601 Anti-CCP hs® (high sensitive)

INTENDED PURPOSE

Anti-CCP hs® (high sensitive) is an ELISA test system for the quantitative measurement of IgG class autoantibodies against cyclic citrullinated peptides (CCP) in human serum or plasma. This product is intended for professional in vitro diagnostic use only.

Measurement of anti-CCP antibodies may aid in the diagnosis of rheumatoid arthritis (RA), where anti-CCP antibody levels represent one parameter of a multi-criterion diagnostic process, encompassing both clinical and laboratory-based assessments.

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



Ready to use

PRINCIPLE OF THE TEST

Highly purified cyclic citrullinated vimentin peptides (CCP) is bound to microwells.

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CONTENTS OF THE KIT

ORG 601	 96	Sufficient for 96 determinations
MICROPLATE	1	One divisible microplate consisting of 12 modules of 8 wells each. Ready to use. Product code on module: CCP
CALIBRATOR A	1x 1.5 ml	Calibrator A 0 U/ml, containing serum/buffer matrix (PBS, BSA, detergent, NaN3 0.09%), yellow. Ready to use.
CALIBRATOR B	1x 1.5 ml	Calibrator B 20 U/ml, containing CCP antibodies in a serum/buffer matrix (PBS, BSA, detergent, NaN3 0.09%), yellow. Ready to use.
CALIBRATOR C	1x 1.5 ml	Calibrator C 40 U/ml, containing CCP antibodies in a serum/buffer matrix (PBS, BSA, detergent, NaN3 0.09%), yellow. Ready to use.
CALIBRATOR D	1x 1.5 ml	Calibrator D 100 U/ml, containing CCP antibodies in a serum/buffer matrix (PBS, BSA, detergent, NaN3 0.09%), yellow. Ready to use.
CALIBRATOR E	1x 1.5 ml	Calibrator E 300 U/ml, containing CCP antibodies in a serum/buffer matrix (PBS, BSA, NaN3 0.09%), yellow. Ready to use.
CALIBRATOR F	1x 1.5 ml	Calibrator F 1000 U/ml, containing CCP antibodies in a serum/buffer matrix (PBS, BSA, detergent, NaN3 0.09%), yellow. Ready to use.
CONTROL +	1x 1.5 ml	Control positive, containing CCP antibodies in a serum/buffer matrix (PBS, BSA, detergent, NaN3 0.09%), yellow. Ready to use. The concentration is specified on the certificate of analysis.
CONTROL -	1x 1.5 ml	Control negative, containing CCP antibodies in a serum/buffer matrix (PBS, BSA, detergent, NaN3 0.09%), yellow. Ready to use. The concentration is specified on the certificate of analysis.
DILUENT	20 ml	Sample Buffer P, containing PBS, BSA, detergent, preservative sodium azide 0.09%, yellow, concentrate (5 x).
CONJUGATE	15 ml	Enzyme Conjugate containing anti-human IgG antibodies, HRP labelled; PBS, BSA, detergent, preservative PROCLIN 0.05%, light red. Ready to use.
TMB	15 ml	TMB Substrate; containing 3,3', 5,5'- Tetramethylbenzidin, colorless. Ready to use.
STOP	15 ml	Stop solution; contains acid. Ready to use.
WASH	20 ml	Wash Buffer, containing Tris, detergent, preservative sodium azide 0.09%; 50 x conc.
 i	1	Certificate of Analysis

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1. Pipette **100 µl** of calibrators, controls and prediluted patient samples into the wells.
Incubate for **30 minutes** at room temperature (20-28 °C).
Discard the contents of the microwells and **wash 3 times** with **300 µl** of wash solution.
2. Dispense **100 µl** of enzyme conjugate into each well.
Incubate for **15 minutes** at room temperature.
Discard the contents of the microwells and **wash 3 times** with **300 µl** of wash solution.
3. Dispense **100 µl** of TMB substrate solution into each well.
Incubate for **15 minutes** at room temperature
4. **Add 100 µl** of stop solution to each well of the modules
Incubate for **5 minutes** at room temperature.
Read the optical density at 450 nm (reference 600-690nm) and calculate the results.
The developed colour is stable for at least 30 minutes. Read during this time.

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PERFORMANCE CHARACTERISTICS

Calibration

This assay system is calibrated in relative arbitrary units. It is calibrated against an external anti-CCP Assay, since no international reference sera for RA diagnostic are available so far.

Measuring range

The calculation range of this ELISA assay is 0 - 1000 U/ml

Expected values

In a normal range study with samples from healthy blood donors the following ranges have been established with this ELISA assay: Cut-off 20 U/ml

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		Clinical Diagnosis		
		POS	NEG	
ORG 601	POS	237	8	
	NEG	22	169	
		259	177	436
Sensitivity: 91.5 %				
Specificity: 95.5 %				
Overall agreement: 93.1 %				

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REFERENCES

- Bang H, Egerer K, Gauliard A, et al. Mutation and citrullination modifies vimentin to a novel autoantigen for rheumatoid arthritis. *Arthritis Rheum.* 56(8):2503, 2007.
- Aletaha D, et al. Rheumatoid arthritis classification criteria: An American College of Rheumatology/European League Against Rheumatism collaborative initiative. *Arthritis & Rheumatism* 62 (9):2569-2581, 2010.
- Prujin GJ, Wiik A, van Venrooij WJ. The use of citrullinated peptides and proteins for the diagnosis of rheumatoid arthritis. *Arthritis Res Ther* 12 (1):203, 2010.
- Snir O, Widhe M, Hermansson M, von Spee C, Lindberg J, Hensen S, Lundberg K, Engstrom A, Venables PJ, Toes RE, Holmdahl R, Klareskog L, Malmstrom V. Antibodies to several citrullinated antigens are enriched in the joints of rheumatoid arthritis patients. *Arthritis Rheum* 62 (1):44-52, 2009.
- Van Steendam K, Tillemann K, Deforce D. The relevance of citrullinated vimentin in the production of antibodies against citrullinated proteins and the pathogenesis of rheumatoid arthritis. *Rheumatology (Oxford)*, 2011.
- Van Steendam K, Tillemann K, De Ceuleneer M, De Keyser F, Elewaut D, Deforce D. Citrullinated vimentin as an important antigen in immune complexes from synovial fluid of rheumatoid arthritis patients with antibodies against citrullinated proteins. *Arthritis Res Ther* 12 (4):R132, 2010.
- Syversen SW, Goll GL, van der Heijde D, Landewe R, Lie BA, Odegard S, Uhlig T, Gaarder PI, Kvien TK. Prediction of radiographic progression in rheumatoid arthritis and the role of antibodies against mutated citrullinated vimentin: results from a ten-year prospective study. *Ann.Rheum.Dis.* 69 (2):345-351, 2009.
- Innala L, Kokkonen H, Eriksson C, Jidell E, Berglin E, Dahlqvist SR. Antibodies against mutated citrullinated vimentin are a better predictor of disease activity at 24 months in early rheumatoid arthritis than antibodies against cyclic citrullinated peptides. *J Rheumatol.* 35(6):1002, 2008.
- Mathsson L, Mullazehi M, Wick MC, et al. Antibodies against citrullinated vimentin in rheumatoid arthritis: higher sensitivity and extended prognostic value concerning future radiographic progression as compared with antibodies against cyclic citrullinated peptides. *Arthritis Rheum.* 58(1):36, 2008.
- Soós L, Szekanecz Z, Szabó Z, et al. Clinical evaluation of anti-mutated citrullinated vimentin by ELISA in rheumatoid arthritis. *J Rheumatol.* 34(8):1658, 2007.
- Bizzaro N, Tonutti E, Tozzoli R, Villalta D. Analytical and diagnostic characteristics of 11 2nd- and 3rd-generation immunoenzymatic methods for the detection of antibodies to citrullinated proteins. *Clin Chem.* 53 (8):1527, 2007.
- Coenen D, Verschuere P, Westhovens R, Bossuyt X. Technical and diagnostic performance of 6 assays for the measurement of citrullinated protein/peptide antibodies in the diagnosis of rheumatoid arthritis. *Clin Chem.* 53(3):498, 2007.
- Dejaco C, Klotz W, Larcher H, Duftner C, Schirmer M, Herold M. Diagnostic value of antibodies against a modified citrullinated vimentin in rheumatoid arthritis. *Arthritis Res Ther.* 8(4):R119, 2006.
- Keskin G, Inal A, Keskin D, et al. Diagnostic utility of anti-cyclic citrullinated peptide and anti-modified citrullinated vimentin antibodies in rheumatoid arthritis. *Protein Pept Lett.* 15(3):314 2008.
- Poulsom H, Charles PJ. Antibodies to citrullinated vimentin are a specific and sensitive marker for the diagnosis of rheumatoid arthritis. *Clin Rev Allergy Immunol.* 34(1):4 2008.
- Song JS, Park GB, Park AJ. Comparison of anti-mutated citrullinated vimentin with anti-cyclic citrullinated peptide and rheumatoid factors for the diagnostic value of rheumatoid arthritis [in Korean]. *J Korean Rheum Assoc.* 14(3):235 2007.

- Ursum J, Nielsen MM, van Schaardenburg D, et al. Antibodies to mutated citrullinated vimentin and disease activity score in early arthritis: a cohort study. *Arthritis Res Ther.* 10(1):R12, 2008.
- Besada E, Nikolaisen C, Nossent H. Diagnostic value of antibodies against mutated citrullinated vimentin for rheumatoid arthritis. *Clin Exp Rheumatol.* 29(1):85 2011.

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
Change Control

Former version: ORG 601_IFU_EN_QM113201_2016-04-18_2

Reason for revision: *Introduction electronic IFU on homepage*

ORGENTEC Diagnostika GmbH

Carl-Zeiss-Straße 49-51

55129 Mainz - Germany

Phone: +49 (0) 61 31 / 92 58-0

Fax: +49 (0) 61 31 / 92 58-58

Internet: www.orgentec.com

508_3

ORG 508 Anti-SS-A

INTENDED PURPOSE

Anti-SS-A is an ELISA test system for the quantitative measurement of IgG class autoantibodies against SS-A (52 and 60 kDa) in human serum or plasma. This product is intended for professional in vitro diagnostic use only.

This test is useful for the differential diagnosis and monitoring of systemic rheumatic inflammatory autoimmune diseases. Autoantibodies against the two antigens SS-A 52 and SS-A 60 are predominantly found in cases of Sjogren's syndrome. Evaluation of a test result should always take into account all clinical and laboratory diagnostic findings.

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Ready to use

PRINCIPLE OF THE TEST

Highly purified SS-A (52 and 60 kDa) is bound to microwells.

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CONTENTS OF THE KIT

ORG 508		96	Sufficient for 96 determinations
MICROPLATE		1	One divisible microplate consisting of 12 modules of 8 wells each. Ready to use. Color code on module
CALIBRATOR A		1x 1.5 ml	Calibrator A 0 U/ml, containing serum/buffer matrix (PBS, BSA, detergent, NaN3 0.09%), yellow. Ready to use.
CALIBRATOR B		1x 1.5 ml	Calibrator B 12.5 U/ml, containing SS-A antibodies in a serum/buffer matrix (PBS, BSA, detergent, NaN3 0.09%), yellow. Ready to use.
CALIBRATOR C		1x 1.5 ml	Calibrator C 25 U/ml, containing SS-A antibodies in a serum/buffer matrix (PBS, BSA, detergent, NaN3 0.09%), yellow. Ready to use.
CALIBRATOR D		1x 1.5 ml	Calibrator D 50 U/ml, containing SS-A antibodies in a serum/buffer matrix (PBS, BSA, detergent, NaN3 0.09%), yellow. Ready to use.
CALIBRATOR E		1x 1.5 ml	Calibrator E 100 U/ml, containing SS-A antibodies in a serum/buffer matrix (PBS, BSA, NaN3 0.09%), yellow. Ready to use.
CALIBRATOR F		1x 1.5 ml	Calibrator F 200 U/ml, containing SS-A antibodies in a serum/buffer matrix (PBS, BSA, detergent, NaN3 0.09%), yellow. Ready to use.
CONTROL +		1x 1.5 ml	Control positive, containing SS-A antibodies in a serum/buffer matrix (PBS, BSA, detergent, NaN3 0.09%), yellow. Ready to use. The concentration is specified on the certificate of analysis.
CONTROL -		1x 1.5 ml	Control negative, containing SS-A antibodies in a serum/buffer matrix (PBS, BSA, detergent, NaN3 0.09%), yellow. Ready to use. The concentration is specified on the certificate of analysis.
DILUENT		20 ml	Sample Buffer P, containing PBS, BSA, detergent, preservative sodium azide 0.09%, yellow, concentrate (5 x).
CONJUGATE		15 ml	Enzyme Conjugate containing anti-human IgG antibodies, HRP labelled; PBS, BSA, detergent, preservative PROCLIN 0.05%, light red. Ready to use.
TMB		15 ml	TMB Substrate; containing 3,3', 5,5'- Tetramethylbenzidin, colorless. Ready to use.
STOP		15 ml	Stop solution; contains acid. Ready to use.
WASH		20 ml	Wash Buffer, containing Tris, detergent, preservative sodium azide 0.09%; 50 x conc.
		1	Certificate of Analysis

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CALCULATION OF RESULTS

For quantitative results plot the optical density of each calibrator versus the calibrator concentration to create a calibration curve. The concentration of patient samples may then be estimated from the calibration curve by interpolation.
Using data reduction software a 4-Parameter-Fit with lin-log coordinates for optical density and concentration is the data reduction method of choice.

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Reproducibility

Intra-assay precision: Coefficient of variation (CV) was calculated for each of three samples from the results of 24 determinations in a single run. Results for precision-within-assay are shown in the table below.
Inter-assay precision: Coefficient of variation (CV) was calculated for each of three samples from the results of 6 determinations in 5 different runs. Results for run-to-run precision are shown in the table below.

Intra-Assay		
Sample	Mean U/ml	CV %
1	32.2	2.7
2	73.2	2.6
3	134.0	3.6

Inter-Assay		
Sample	Mean U/ml	CV %
1	33.8	6.4
2	71.3	6.2
3	133.1	1.1

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Study results

<u>Study population</u>	<u>n</u>	<u>n Pos</u>	<u>%</u>
Sjogren's syndrome	70	51	72.9
Normal human sera	100	7	7.0

		Clinical Diagnosis		
		POS	NEG	
ORG 508	POS	51	7	70 100 170
	NEG	19	93	
Sensitivity:		72.9	%	
Specificity:		93.0	%	
Overall agreement:		84.7	%	

LIMITATIONS OF THE PROCEDURE

This assay is a diagnostic aid. A definite clinical diagnosis should not be based on the results of a single test, but should be made by the physician after all clinical and laboratory findings have been evaluated concerning the entire

clinical picture of the patient. Also every decision for therapy should be taken individually.

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Change Control

Former version: *ORG 508_IFU_EN_QM113135_2013-12-16_1.2* Reason for revision: *Introduction electronic IFU on homepage*

REFERENCES

1. Alba P, Bento L, Cuadrado MJ, Karim Y, Tungekar MF, Abbs I et al. Anti-dsDNA, anti-Sm antibodies, and the lupus anticoagulant: significant factors associated with lupus nephritis. *Ann Rheum Dis* 2003; 62(6):556-560.
2. Antico A, Platzgummer S, Bassetti D, Bizzaro N, Tozzoli R, Villalta D. Diagnosing systemic lupus erythematosus: new-generation immunoassays for measurement of anti-dsDNA antibodies are an effective alternative to the Farr technique and the Crithidia luciliae immunofluorescence test. *Lupus* 2010; 19(8):906-912.
3. Brouwer R, Hengstman GJ, Vree EW, Ehrfeld H, Bozic B, Ghirardello A et al. Autoantibody profiles in the sera of European patients with myositis. *Ann Rheum Dis* 2001; 60(2):116-123.
4. Castro C, Gourley M. Diagnostic testing and interpretation of tests for autoimmunity. *J Allergy Clin Immunol* 2010; 125(2 Suppl 2):S238-S247.
5. Defendenti C, Atzeni F, Spina MF, Grosso S, Cereda A, Guercilena G et al. Clinical and laboratory aspects of Ro/SSA-52 autoantibodies. *Autoimmun Rev* 2011; 10(3):150-154.
6. Eriksson C, Kokkonen H, Johansson M, Hallmans G, Wadell G, Rantapaa-Dahlqvist S. Autoantibodies predate the onset of Systemic Lupus Erythematosus in northern Sweden. *Arthritis Research & Therapy* 2011; 13(1):R30.
7. Haugbro K, Nossent JC, Winkler T, Figenschau Y, Rekvig OP. Anti-dsDNA antibodies and disease classification in antinuclear antibody positive patients: the role of analytical diversity. *Ann Rheum Dis JID* - 0372355 2004; 63(4):386-394.
8. Ippolito A, Wallace DJ, Gladman D, Fortin PR, Urowitz M, Werth V et al. Autoantibodies in systemic lupus erythematosus: comparison of historical and current assessment of seropositivity. *Lupus* 2011; 20(3):250-255.
9. Isenberg DA, Manson JJ, Ehrenstein MR, Rahman A. Fifty years of anti-ds DNA antibodies: are we approaching journey's end? *Rheumatology (Oxford)* 2007; 46(7):1052-1056.
10. Kattah NH, Kattah MG, Utz PJ. The U1-snRNP complex: structural properties relating to autoimmune pathogenesis in rheumatic diseases. *Immunol Rev* 2010; 233(1):126-145.
11. Kumar Y, Bhatia A, Minz RW. Antinuclear antibodies and their detection methods in diagnosis of connective tissue diseases: a journey revisited. *Diagn Pathol* 2009; 4:1.
12. Meroni PL, Schur PH. ANA screening: an old test with new recommendations. *Ann Rheum Dis* 2010; 69:1420-1422.
13. Petri M, Magder L. Classification criteria for systemic lupus erythematosus: a review. *Lupus* 2004; 13(11):829-837.
14. Poole BD, Schneider RI, Guthridge JM, Vette CA, Reichlin M, Harley JB et al. Early targets of nuclear RNP humoral autoimmunity in human systemic lupus erythematosus. *Arthritis Rheum* 2009; 60(3):848-859.
15. Putova I, Dostal C, Becvar R. Prevalence of antinucleosome antibodies by enzyme-linked immunosorbent assays in patients with systemic lupus erythematosus and other autoimmune systemic diseases. *Ann N Y Acad Sci* 2007; 1109:275-286.
16. Reveille JD. Predictive value of autoantibodies for activity of systemic lupus erythematosus. *Lupus JID* - 9204265 2004; 13(5):290-297.
17. Simon JA, Cabiedes J, Ortiz E, Alcocer-Varela J, Sanchez-Guerrero J. Anti-nucleosome antibodies in patients with systemic lupus erythematosus of recent onset. Potential utility as a diagnostic tool and disease activity marker. *Rheumatology (Oxford)* 2004; 43(2):220-224.
18. Sinclair D, Saas M, Williams D, Hart M, Goswami R. Can an ELISA replace immunofluorescence for the detection of anti-nuclear antibodies?--The routine use of anti-nuclear antibody screening ELISAs. *Clin Lab* 2007; 53(3-4):183-191.
19. Tozzoli R, Bizzaro N, Tonutti E, Villalta D, Bassetti D, Manoni F et al. Guidelines for the laboratory use of autoantibody tests in the diagnosis and monitoring of autoimmune rheumatic diseases. *Am J Clin Pathol* 2002; 117(2):316-324.
20. Maidhof W., Hilius O. Lupus: an overview of the disease and management options. *P T* 2012; 37(4):240-9.
21. Hahn BH, McMahon MA, Wilkinson A, Wallace WD, Daikh DI, Fitzgerald JD et al. American College of Rheumatology guidelines for screening, treatment, and management of lupus nephritis. *Arthritis Care Res (Hoboken)* 2012; 64(6):797-808.

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- ① Pipet **100 µl** calibrator, control or patient sample
→ Incubate for **30 minutes** at room temperature
→ Discard the contents of the wells and wash 3 times with **300 µl** wash solution
- ② Pipet **100 µl** enzyme conjugate
→ Incubate for **15 minutes** at room temperature
→ Discard the contents of the wells and wash 3 times with **300 µl** wash solution
- ③ Pipet **100 µl** substrate solution
→ Incubate for **15 minutes** at room temperature
- ④ Add **100 µl** stop solution
→ Leave untouched for **5 minutes**
→ Read at **450 nm**

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- antiphospholipid syndrome? *J Autoimmun.* 2000 Sep;15(2):163-72.
7. Cervera R, Piette JC, Font J, Khamashta MA, Shoenfeld Y, *et al*; Euro-Phospholipid Project Group. Antiphospholipid syndrome: clinical and immunologic manifestations and patterns of disease expression in a cohort of 1,000 patients. *Arthritis Rheum.* 2002 Apr;46(4):1019-27.
 8. Linnemann B. Antiphospholipid syndrome - an update. *Vasa.* 2018 Oct;47(6):451-464.
 9. Calvo-Alén J, Toloza SM, Fernández M, Bastian HM, *et al*.; LUMINA Study Group. Systemic lupus erythematosus in a multiethnic US cohort (LUMINA). XXV. Smoking, older age, disease activity, lupus anticoagulant, and glucocorticoid dose as risk factors for the occurrence of venous thrombosis in lupus patients. *Arthritis Rheum.* 2005 Jul;52(7):2060-8.
 10. Janardhan V, Wolf PA, Kase CS, Massaro JM, *et al*. Anticardiolipin antibodies and risk of ischemic stroke and transient ischemic attack: the Framingham cohort and offspring study. *Stroke.* 2004 Mar;35(3):736-41.
 11. Harris EN, Pierangeli SS. Revisiting the anticardiolipin test and its standardization. *Lupus.* 2002;11(5):269-75.
 12. McClain MT, Arbuckle MR, Heinlen LD, Dennis GJ, *et al*. The prevalence, onset, and clinical significance of antiphospholipid antibodies prior to diagnosis of systemic lupus erythematosus. *Arthritis Rheum.* 2004 Apr;50(4):1226-32.
 13. Singh AK. Immunopathogenesis of the antiphospholipid antibody syndrome: an update. *Curr Opin Nephrol Hypertens.* 2001 May;10(3):355-8.
 14. de Laat B, de Groot PG. Autoantibodies directed against domain I of beta2-glycoprotein I. *Curr Rheumatol Rep.* 2011 Feb;13(1):70-6.
 15. Suh-Lailam BB, Cromar A, Davis KW, Tebo AE. APhL antibody ELISA as an alternative to anticardiolipin test for the diagnosis of antiphospholipid syndrome. *Int J Clin Exp Pathol.* 2012;5(3):210-5.
 16. Ofer-Shiber S, Molad Y. Frequency of vascular and pregnancy morbidity in patients with low vs. moderate-to-high titers of antiphospholipid antibodies. *Blood Coagul Fibrinolysis.* 2015 Apr;26(3):261-6.
 17. Meroni PL, Borghi MO, Raschi E, Tedesco F. Pathogenesis of antiphospholipid syndrome: understanding the antibodies. *Nat Rev Rheumatol.* 2011 Jun;7(6):330-9.
 18. de Laat B, Pengo V, Pabinger I, Musial J, *et al*. The association between circulating antibodies against domain I of beta2-glycoprotein I and thrombosis: an international multicenter study. *J Thromb Haemost.* 2009 Nov;7(11):1767-73.
 19. Sanmarco M, Roll P, Gayet S, Oksman F, *et al*. J. Combined search for anti-beta2-glycoprotein I and anticardiolipin antibodies in antiphospholipid syndrome: contribution to diagnosis. *J Lab Clin Med.* 2004 Sep;144(3):141-7.
 20. Otomo K, Atsumi T, Amengual O, Fujieda Y, *et al*. Efficacy of the antiphospholipid score for the diagnosis of antiphospholipid syndrome and its predictive value for thrombotic events. *Arthritis Rheum.* 2012 Feb;64(2):504-12.
 21. Basic QC Practices On-line Course; <http://www.Westgard.com>.
 22. Boscatto, LM. and Stuart, MC., 'Heterophilic antibodies: a problem for all immunoassays'. *Clin Chem*, 34, 1988, pp 27-33.

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Anti Cardiolipin Screen

for routine analysis

Quantitative determination of IgG or IgM autoantibodies against Cardiolipin in human serum or plasma

IVD

LOT

See external label

2°C 8°C

 $\Sigma = 96$ tests

REF DKO144

1. INTENDED PURPOSE

For *In Vitro* Diagnostic Use
For Laboratory Professional Use

Anti Cardiolipin Screen is a manual *in vitro* diagnostic device intended for the quantitative determination of both IgG and IgM class antibodies directed against cardiolipin in human serum or plasma. Results are to be used in conjunction with other clinical and laboratory data as an aid in the diagnosis of antiphospholipid syndrome (APS).

2. CLINICAL SIGNIFICANCE

Cardiolipin is a negatively charged phospholipid which is typically located in the inner mitochondrial membrane¹. Autoantibodies directed against cardiolipin are part of a group known as antiphospholipid antibodies which includes anti- β 2 glycoprotein 1 autoantibodies. Measurement of anti-cardiolipin autoantibodies is considered to be one of the most important markers to support diagnosis of antiphospholipid syndrome (APS)^{2,3}.

APS is a systemic autoimmune disorder characterised by a combination of arterial and/or venous thromboses, pregnancy complications, such as recurrent foetal loss, together with elevated levels of antiphospholipid antibodies⁴. APS was first described in patients with systemic lupus erythematosus (SLE), though it has subsequently been established that SLE may be independent of an underlying disease⁵.

APS can occur alone – Primary APS, or in association with other conditions, such as SLE – Secondary APS⁶ (6). However, it has been demonstrated that anti-cardiolipin antibodies (aCLs) can be detected in patients with SLE that do not develop Secondary APS⁷⁻⁹. However, thromboembolic events are the most common clinical manifestation of APS

Anti-cardiolipin antibodies can recognise both cardiolipin and portions of the phospholipid-protein complex β 2 glycoprotein 1-cardiolipin^{10,11}.

Studies have indicated an association of anti-cardiolipin IgG and IgM antibodies^{2,7,11-16}, with thrombotic events whereas others suggest these to be linked with IgG isotype, but not IgM^{6,10,17}.

aCL IgM antibodies have been shown to occur in infections such as chronic hepatitis C, leprosy, syphilis, but they are not directly involved in thrombotic events¹⁷.

It is likely that the presence of antiphospholipid antibodies, including aCL IgG and IgM constitutes the single most recognisable risk factor in cases of recurrent pregnancy loss and late placenta-mediated obstetric complications^{4-6,11,14-15,18,19}. Where patients can present with only adverse pregnancy outcomes with isolated vascular events or with both obstetric and thrombotic manifestations⁶. It has been suggested that anti- β 2-glycoprotein I – cardiolipin antibodies are able to recognise the antigen on placental tissues, inhibiting the growth and differentiation of trophoblasts which may eventually cause defective placentation²⁰.

3. PRINCIPLE OF THE METHOD

The Anti Cardiolipin Screen allows the determination of autoantibodies directed against the Cardiolipin- β 2-glycoprotein complex through two different calibration curves and enzyme conjugates (one specific for IgG test, one specific for IgM test) and one microplate. The principle of the method and assay procedure are the same for both assessments. Use reagents for IgG or reagents for IgM depending on the isotype which is under investigation.

The Anti Cardiolipin Screen is a two-step sandwich enzyme immunometric assay (ELISA) where patient samples, calibrators or controls are incubated on microtitre plates coated with the antigenic cardiolipin- β 2 glycoprotein complex. During the incubation, antibodies present in the test sample bind to the immobilised antigen complex. After the incubation, the bound/free separation is performed by a simple solid phase washing.

A subsequent incubation occurs with anti-human IgM or IgG conjugated with horseradish peroxidase (HRP), which binds to the immobilised antibodies. A further wash step is performed to remove excess conjugate. Then, a chromogenic substrate solution containing TMB is dispensed into the wells which reacts with the conjugated HRP and a blue colour develops that changes into yellow when the Stop Solution (H_2SO_4) is added. The colour intensity is directly proportional to the Anti Cardiolipin IgM or IgG concentration (depending on the conjugate used) in the sample.

Anti-cardiolipin antibody concentration in the sample is calculated through a calibration curve.

4. REAGENTS, MATERIALS AND INSTRUMENTATION

4.1. Reagents and materials supplied in the kit

For determination of IgG class antibodies

1. Anti Cardiolipin IgG Calibrators

(5 vials, 1.2 mL each)

Phosphate buffer 0.1M, $\text{NaN}_3 < 0.1\%$, human serum

CAL0 **REF** DCE002/11306-0

CAL1 **REF** DCE002/11307-0

CAL2 **REF** DCE002/11308-0

CAL3 **REF** DCE002/11309-0

CAL4 **REF** DCE002/11310-0

2. Controls (2 vials, 1.2 mL each, ready to use)

Phosphate buffer 0.1M, $\text{NaN}_3 < 0.1\%$, human serum

Negative Control **REF** DCE045/11301-0

Positive Control **REF** DCE045/11302-0

3. Conjugate IgG (1 vial, 15 mL)

Anti h-IgG conjugate with horseradish peroxidase (HRP),

BSA 0.1%, ProClin $< 0.0015\%$ **REF** DCE002/11302-0

For determination of IgM class antibodies

1. Calibrators (5 vials, 1.2 mL each)

Phosphate buffer 0.1M, $\text{NaN}_3 < 0.1\%$, human serum

CAL0 **REF** DCE002/11206-0

CAL1 **REF** DCE002/11207-0

CAL2 **REF** DCE002/11208-0

CAL3 **REF** DCE002/11209-0

CAL4 **REF** DCE002/11210-0

2. Controls (2 vials, 1.2 mL each, ready to use)

Phosphate buffer 0.1M, $\text{NaN}_3 < 0.1\%$, human serum

Negative Control **REF** DCE045/11201-0

Positive Control **REF** DCE045/11202-0

3. Conjugate (1 vial, 15 mL)

Anti h-IgM conjugate with horseradish peroxidase (HRP),

BSA 0.1%, ProClin $< 0.0015\%$ **REF** DCE002/11202-0

Common reagents

4. Sample Diluent (1 vial, 100 mL)

Phosphate buffer 0.1 M $\text{NaN}_3 < 0.1\%$

REF DCE053-0

5. Coated Microplate (1 breakable microplate)

Microplate coated with antigenic Cardiolipin- β 2-Glycoprotein complex

REF DCE002/14403-0

6. TMB Substrate (1 vial, 15 mL)

H_2O_2 -TMB (0.26 g/L) (avoid any skin contact)

REF DCE004-0

7. Stop Solution (1 vial, 15 mL)

Sulphuric acid 0.15M (avoid any skin contact)

REF DCE005-0

8. 10X Conc. Wash Solution (1 vial, 50 mL)

Phosphate buffer 0.2M, pH 7.4 **REF** DCE054-0

4.2. Materials required but not provided

Distilled water

4.3. Auxiliary materials and instrumentation

Automatic dispenser

Precision Pipetting Devices

Microplate reader (450 nm, 620-630 nm)

5. WARNINGS

- This kit is intended for *in vitro* use by professional persons only. Not for internal or external use in Humans or Animals.
- Use appropriate personal protective equipment while working with the reagents provided.
- Follow Good Laboratory Practice (GLP) for handling blood products.
- All human source material used in the preparation of the reagents has been tested and found negative for antibody to HIV 1&2, HbsAg, and HCV. No test method however can offer complete assurance that HIV, HBV, HCV or other infectious agents are absent. Therefore, the Calibrators and the Controls should be handled in the same manner as potentially infectious material.
- Material of animal origin used in the preparation of the kit has been obtained from animals certified as healthy and the bovine protein has been obtained from countries not infected by BSE, but these materials should be handled as potentially infectious.
- Some reagents contain small amounts of Sodium Azide (NaN_3) or ProClin™ 300 as preservative. Avoid contact with skin or mucosa.
- Sodium Azide may be toxic if ingested or absorbed through the skin or eyes; moreover, it may react with lead or copper plumbing to form potentially explosive metal azides. If you use a sink to remove the reagents, wash through large with amounts of water to prevent azide build-up.
- The TMB Substrate contains an irritant, which harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous, corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/ H_2O_2 to direct sunlight, metals or oxidants. Do not freeze the solution.

6. PRECAUTIONS

- Please adhere strictly to the sequence of pipetting steps provided in this protocol. The performance data represented here were obtained using specific reagents listed in this Instruction For Use.
- All reagents should be stored refrigerated at 2-8°C in their original container. Any exceptions are clearly indicated.
- Allow all kit components and specimens to reach room temperature (22-28°C) and mix well prior to use.
- Do not interchange kit components from different lots. The expiry date printed on box and vials labels must be observed. Do not use any kit component beyond their expiry date.
- **WARNING: the conjugate reagent is designed to ensure maximum dose sensitivity and may be contaminated by external agents if not used properly;** therefore, it is recommended to use disposable consumables (tips, bottles, trays, etc.). For divided doses, take the exact amount of conjugate needed and do not re-introduce any waste product into the original bottle. In addition, **for doses dispensed with the aid of automatic and semi-automatic devices,** before using the conjugate, it is advisable to clean the fluid handling system, ensuring that the

procedures of washing, deproteinisation and decontamination are effective in avoiding contamination of the conjugate; this procedure is highly recommended when the kit is processed using analysers which are not equipped with disposable tips.

For this purpose, DiaMetra supplies a separate decontamination reagent for cleaning needles.

- If you use automated equipment, the user has the responsibility to make sure that the kit has been appropriately tested.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background. To improve the performance of the kit on automatic systems is recommended to increase the number of washes.
- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate
- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.
- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.
- Maximum precision is required for reconstitution and dispensation of the reagents.
- Samples microbiologically contaminated, highly lipemic or haemolysed should not be used in the assay.
- Plate readers measure vertically. Do not touch the bottom of the wells.

7. REAGENT STORAGE AND STABILITY

Store the kit at 2 – 8°C in the dark.

- The kit is stable at 2 – 8°C until the expiry date stated on the external kit label.
- Once opened, the kit is stable at 2 – 8°C for 6 months.
- The diluted wash solution is stable for 30 days at 2-8°C.

Important note: open the bag containing the Coated Microplate only when it is at room temperature and close it immediately after use.

8. SAMPLE COLLECTION AND STORAGE

The assay should be performed using serum (standard sampling tubes or tubes containing serum separating gel) or plasma (lithium heparin, sodium heparin, or potassium EDTA) samples.

Sample Storage	Duration
2 – 8 °C	96 hours
Freeze/thaw cycles	3 cycles

9. PROCEDURE

9.1. Preparation of Calibrators and Controls

The calibrators and controls are ready to use.

The calibrators and approximately the following concentrations:

	C ₀	C ₁	C ₂	C ₃	C ₄
AU/mL	0	5	10	20	80

9.2. Preparation of the Wash Solution

Dilute the content of the vial “10X Conc. Wash Solution” with distilled water to a final volume of 500 mL prior to use. For smaller volumes respect the 1:10 dilution ratio.

It is possible to observe the presence of crystals within the concentrated wash solution; in this case mix at room temperature until the complete dissolution of crystals. For greater accuracy, dilute the whole bottle of concentrated wash solution to 500 mL, taking care also to transfer crystals completely by rinsing of the bottle, then mix until crystals are completely dissolved.

9.3. Preparation of Samples

All serum and plasma samples must be diluted 1:100 with sample diluent.

e.g., 10 µL of sample should be diluted with 990 µL of sample diluent.

9.4. Procedure

- **Allow all reagents to reach room temperature (22-28°C) for at least 30 minutes.** At the end of the assay, immediately store the reagents at 2-8°C: avoiding long exposure to room temperature.
- Unused coated microwell strips should be released securely in the foil pouch containing desiccant and stored at 2-8°C.
- To avoid potential microbial and/or chemical contamination, unused reagents should never be transferred into the original vials.
- As it is necessary to perform the determination in duplicate to improve accuracy of the test results, prepare two wells for each point of the calibration curve (C₀-C₄), two for each control, two for each sample, one for blank.

The following procedure is the same for both class IgG and IgM antibodies assay.

Reagents	Calibrator	Sample/ Controls	Blank
----------	------------	---------------------	-------

Use reagents for IgG or reagents for IgM depending on the isotype which is under investigation

Calibrator C ₀ -C ₄ (IgG or IgM)	100 µL		
Controls (IgG or IgM)		100 µL	
Diluted Sample		100 µL	

Incubate for 60 minutes at room temperature (22-28°C). Remove the content from each well; wash the wells 3 times with 300 µL diluted wash solution.

Important note: during each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel.

Automatic washer: if you use automated equipment, wash the wells at least 5 times.

Conjugate (IgG or IgM)	100 µL	100 µL	
---------------------------	--------	--------	--

Incubate for 60 minutes at room temperature (22-28°C). Remove the content from each well; wash the wells 3 times with 300 µL diluted wash solution.

Washing: follow the same indications of the previous point.

TMB Substrate	100 µL	100 µL	100 µL
---------------	--------	--------	--------

Incubate for 15 minutes in the dark at room temperature (22-28°C).

Stop Solution	100 µL	100 µL	100 µL
---------------	--------	--------	--------

Shake the microplate gently. Read the absorbance (E) at 450 nm against a reference wavelength of 620-630 nm or against Blank within 5 minutes.

10. QUALITY CONTROL

Good Laboratory Practice (GLP) requires the use of quality control specimens in each series of assays in order to check the performance of the assay. Controls should be treated as unknown samples, and the results analysed with appropriate statistical methods.

The kit controls provided in the kit should be tested as unknowns and are intended to assist in assessing the validity of results obtained with each assay plate.

The mean concentration of each control level is documented in the QC report included with each kit. These mean concentration levels are determined over several assays which are run in duplicate in multiple locations across each plate. This test is only valid if the optical density at 450 nm for the controls as well as for the calibrators (C₀-C₄) fall within with the respective range indicated on the Quality Control Certificate enclosed in each test kit.

DiaMetra recommends the users to maintain graphic records of the control values generated with each assay run, including the running means, SDs and %CVs. This

information will facilitate the controls trending analysis relating to the performance of current and historical control lots relative to the supplied Quality Control data. The trending will assist in the identification of assays which give control values significantly different from their average range.

When interpreting control data, users should note that this product was designed and developed as a manual product. The range stated on the QC certificate should be appropriate for assays that are performed manually and with strict adherence to the Assay Procedure described above. It is recognised by Quality Control professionals, that as a result of differences in conditions and practices, there will always be variability in the mean values and precision of control measurements between different laboratories²¹.

11. CALCULATION OF RESULTS

A variety of data reduction software packages are available, which may be employed to generate the mean calibration curve and to calculate the mean concentrations of unknown samples and controls. A 4-parameter logistic (4PL) curve fit with lin-log coordinates, **including Calibrator 0 is required**. A smoothed spline fit including Calibrator 0 can be used. Other curve fitting algorithms are not recommended.

Alternatively, a calibration curve may be prepared on semi-log graph paper by plotting mean absorbance on the Y-axis against concentration of analyte on the X-axis. Calibrator 0 should be included in the calibration curve. Read the mean absorbance value of each unknown sample off the curve.

12. MEASURING RANGE

12.1. For assessment of IgG class antibodies

The assay measuring range (AMR) is 2 – 80 AU/mL.

Any value that reads below 2 AU/mL should be reported as "< 2 AU/mL". Any value that reads above 80 AU/mL should be reported as "> 80 AU/mL".

12.2. For assessment of IgM class antibodies

The assay measuring range (AMR) is 2.08 – 80 AU/mL.

Any value that reads below 2.09 AU/mL should be reported as "< 2.08 AU/mL". Any value that reads above 80 AU/mL should be reported as "> 80 AU/mL".

13. METROLOGY AND TRACEABILITY

13.1. For assessment of IgG class antibodies

The calibrators of this kit are traceable to the Centres for Disease Control (CDC) Human IgG Anti-Cardiolipin Monoclonal Antibody HCAL - Catalogue [IS2717].

13.2. For assessment of IgM class antibodies

The calibrators of this kit are traceable to the Centres for Disease Control (CDC) Human IgM Anti-Cardiolipin Monoclonal Antibody EY2C9 [IS2718].

14. INTERPRETATION OF RESULTS

Concentration	Interpretation
< 8 AU/mL	The sample should be considered negative
8 – 10 AU/mL	The sample should be graded equivocal and repeat testing / sampling should be performed according to internal practices
>10 AU/mL	The sample should be considered positive

Determination of a range of expected values for a “normal” population of a given method is dependent on many factors, such as specificity and sensitivity of the method used and type of population under investigation. Therefore, each laboratory should consider the range given by the Manufacturer as a general indication and produce their own range of expected values based on the indigenous population.

Positive results should be verified concerning the entire clinical status of the patient, with the decision for therapy being taken on an individual basis. It is recommended that each laboratory establishes its own normal and pathological ranges of Anti-Cardiolipin antibody values.

15. PERFORMANCE CHARACTERISTICS

Representative performance data are shown. Results obtained at individual laboratories may vary.

15.1. For assessment of IgG class antibodies

15.1.1. Detection Capability

The limit of blank (LoB), limit of detection (LoD) and limit of quantitation (LoQ) were determined with guidance from CLSI EP17-A, “Protocols for Determination of Limits of Detection and Limits of Quantitation” using 6 blanks and 6 low level samples.

Sensitivity	Concentration
Limit of Blank (LoB)	0.59 AU/mL
Limit of Detection (LoD)	1.25 AU/mL
Limit of Quantitation (LoQ)	2.00 AU/mL

15.1.2. Trueness

Trueness of the Anti Cardiolipin Screen for assessment of IgG class antibodies has been demonstrated through performance of a recovery test using the CDC Human IgG Anti-Cardiolipin Monoclonal Antibody HCAL - Catalogue [IS2717].

15.1.3. Diagnostic Sensitivity and Specificity

The sensitivity and specificity were determined with guidance from CLSI EP-24 “Assessment of the Diagnostic Accuracy of Laboratory Tests Using Receiver Operating Characteristic Curves” using 50 negative and 51 positive samples run on two reagent lots.

		DKO144 - IgG		Total
		Positive	Negative	
True state	Positive	47	4	51
	Negative	0	57	57
Total		47	61	108

Diagnostic sensitivity: 92%

Diagnostic specificity: 100%

15.1.4. Precision

Precision of the Anti Cardiolipin Screen for determination of IgG class antibodies was determined by performing a complex precision study.

Repeatability: A total of 6 serum samples were assayed in 5 replicates, once a day for 5 days by 3 operators.

Data from one representative lot is shown below:

Sample	n	Mean Conc. (AU/mL)	Within run (Repeatability)	
			SD	CV
1	75	6.95	0.38	5.5%
2	75	11.03	0.51	4.6%
3	75	20.12	0.94	4.7%
4	75	30.26	1.86	6.1%
5	75	50.11	2.58	5.1%
6	75	71.71	2.53	3.5%

Reproducibility: A total of 6 serum samples were assayed in 5 replicates, once a day for 5 days by 3 operators.

Results for the combined data from two lots is shown below:

Sample	n	Mean Conc. (AU/mL)	Within Laboratory (Reproducibility)	
			SD	CV%
1	150	6.98	0.49	7.0%
2	150	11.17	0.84	7.5%
3	150	20.16	1.87	9.3%
4	150	30.38	3.06	10.1%
5	150	50.76	5.02	9.9%
6	150	72.30	3.93	5.4%

15.1.5. Linearity

Linearity was evaluated based on CLSI EP-06, “Evaluation of the Linearity of Quantitative Measurement Procedures”. For anti-cardiolipin IgG concentration by Anti Cardiolipin Screen, the measurement procedure shows linearity for the interval from 0.84 to 83.68 ng/mL within the allowable deviation of linearity (ADL) of $\pm 15\%$.

15.2. Analytical Specificity

The following substances do not interfere with a bias of $> \pm 15\%$ in the Anti Cardiolipin Screen assay when assessing IgG class antibodies when the concentrations are below the stated threshold presented in the following table.

Potentially Interfering Reagent	Threshold Concentration
Bilirubin, conjugated	15 mg/dL
Bilirubin, unconjugated	15 mg/dL
Haemoglobin	200 mg/dL
Total Protein	10 g/dL
Triglyceride	500 mg/dL

15.2.1. Serum-plasma study

The Anti Cardiolipin Screen matrix comparison study for assessment of IgG class antibodies was performed to evaluate the difference across tube types (serum separator tubes (SST), lithium heparin plasma, sodium heparin plasma and K2 EDTA plasma) versus the control samples (red top serum, without additive) following CLSI EP9-A3 guidelines. A total of 22 samples (18 native, 4 spiked) to cover the range were evaluated. Linear regression analysis was performed on the comparative data:

Sample type	Slope [95% CI]	Intercept (ng/mL) [95% CI]	Correlation coefficient (r)
SST	0.96 [0.92 – 0.98]	0.64 [-0.38 to 1.66]	1.00
Lithium Heparin	0.92 [0.88 to 0.96]	0.87 [-0.24 to 1.99]	1.00
Sodium Heparin	0.94 [0.89 to 0.98]	0.66 [-0.75 to 2.06]	0.99
EDTA	0.95 [0.92 to 0.99]	0.54 [-0.69 to 1.77]	1.00

15.3. For assessment of IgM class antibodies

15.3.1. Detection Capability

The limit of blank (LoB), limit of detection (LoD) and limit of quantitation (LoQ) were determined with guidance from CLSI EP17-A, "Protocols for Determination of Limits of Detection and Limits of Quantitation" using 6 blanks and 6 low level samples.

Sensitivity	Concentration
Limit of Blank (LoB)	0.76 AU/mL
Limit of Detection (LoD)	1.45 AU/mL
Limit of Quantitation (LoQ)	2.08 AU/mL

15.3.2. Trueness

Trueness of the Anti Cardiolipin Screen for assessment of IgM class antibodies demonstrated through performance of a recovery test using the CDC Human IgM Anti-Cardiolipin Monoclonal Antibody EY2C9 [IS2718].

15.3.3. Diagnostic Sensitivity and Specificity

The sensitivity and specificity were determined with guidance from CLSI EP-24 "Assessment of the Diagnostic Accuracy of Laboratory Tests Using Receiver Operating Characteristic Curves" using 73 negative and 62 positive samples run on two reagent lots.

		DKO144 IgM		Total
		Positive	Negative	
True state	Positive	50	12	62
	Negative	0	73	73
Total		50	85	135

Diagnostic sensitivity: 80%

Diagnostic specificity: 100%

15.3.4. Precision

Precision of the Anti Cardiolipin Screen for determination of IgM class antibodies determined by performing a complex precision study.

Repeatability: A total of 6 serum samples were assayed in 5 replicates, once a day for 5 days by 3 operators. Data from one representative lot is shown below:

Sample	n	Mean Conc. (AU/mL)	Within run (Repeatability)	
			SD	CV
1	75	7.74	0.40	5.2%
2	75	12.46	0.51	4.1%
3	75	21.06	0.99	4.7%
4	75	32.07	1.46	4.6%
5	75	55.15	1.19	2.2%
6	75	75.45	2.33	3.1%

Reproducibility: A total of 6 serum samples were assayed in 5 replicates, once a day for 5 days by 3 operators. Results for the combined data from two lots is shown below:

Sample	n	Mean Conc. (AU/mL)	Within Laboratory (Reproducibility)	
			SD	CV%
1	150	7.65	0.46	6.0%
2	150	12.33	0.74	6.0%
3	150	21.03	1.60	7.6%
4	150	31.89	1.98	6.2%
5	150	54.28	2.70	5.0%
6	150	75.46	2.63	3.5%

15.3.5. Linearity

Linearity was evaluated based on CLSI EP-06, "Evaluation of the Linearity of Quantitative Measurement Procedures". For anti-cardiolipin IgM concentration by Anti Cardiolipin Screen, the measurement procedure shows linearity for the interval from 0.82 to 86.88 AU/mL within the allowable deviation of linearity (ADL) of $\pm 15\%$.

15.3.6. Analytical Specificity

The following substances do not interfere with a bias of $> \pm 15\%$ in the Anti Cardiolipin Screen assay when assessing IgM class antibodies when the concentrations are below the stated threshold presented in the following table.

Potentially Interfering Reagent	Threshold Concentration
Bilirubin, conjugated	15 mg/dL
Bilirubin, unconjugated	15 mg/dL
Haemoglobin	200 mg/dL
Total Protein	10 g/dL
Triglyceride	500 mg/dL

15.3.7. Serum-plasma study

The Anti Cardiolipin Screen matrix comparison study for assessment of IgM class antibodies was performed to evaluate the difference across tube types (serum separator tubes (SST), lithium heparin plasma, sodium heparin plasma and K2 EDTA plasma) versus the control samples (red top serum, without additive) following CLSI EP9-A3 guidelines. A total of 20 samples (16 native, 4 spiked) to cover the assay range were evaluated. Linear regression analysis was performed on the comparative data:

Sample type	Slope [95% CI]	Intercept (ng/mL) [95% CI]	Correlation coefficient (r)
SST	1.02 [0.94 to 1.09]	0.19 [-1.34 to 1.72]	0.99
Lithium Heparin	0.93 [0.82 to 1.04]	0.73 [-1.49 to 2.95]	0.97
Sodium Heparin	0.93 [0.84 to 1.01]	0.64 [-1.00 to 2.28]	0.98

EDTA

0.96 [0.86 to 1.05]	0.65 [-1.19 to 2.49]	0.98
------------------------	-------------------------	------

16. LIMITATIONS OF USE

- As in the case of any diagnostic procedure, results must be interpreted in conjunction with the patient's clinical presentation and other information available to the physician.
- The performance characteristics of this assay have not been established in a paediatric population.
- Heterophilic antibodies in human serum can react with reagent immunoglobulins, interfering with *in vitro* immunoassays²². Patients routinely exposed to animals or to animal serum products can be prone to this interference and anomalous values may be observed.
- The presence of immune complexes or other immunoglobulin aggregates in the patient sample may cause an increased level of non-specific binding and produce false positives in this assay.

17. WASTE MANAGEMENT

Reagents must be disposed of in accordance with local regulations.

18. BIBLIOGRAPHY

1. Dusse LM, Silva FD, Freitas LG, Rios DR, *et al.* Antiphospholipid syndrome: a clinical and laboratorial challenge. *Rev Assoc Med Bras* (1992). 2014 Mar-Apr;60(2):181-6.
2. Roggenbuck D, Somma V, Schierack P, Borghi MO, Meroni PL. Autoantibody profiling in APS. *Lupus*. 2014 Oct;23(12):1262-4.
3. Banzato A, Pengo V. Clinical relevance of β_2 -glycoprotein-I plasma levels in antiphospholipid syndrome (APS). *Curr Rheumatol Rep*. 2014 Jun;16(6):424.
4. Amiral J, Peyrafitte M, Dunois C, Vissac AM, Seghatchian J. Anti-phospholipid syndrome: Current opinion on mechanisms involved, laboratory characterization and diagnostic aspects. *Transfus Apher Sci*. 2017 Aug;56(4):612-625.
5. Miyakis S, Lockshin MD, Atsumi T, Branch, *et al.* International consensus statement on an update of the classification criteria for definite antiphospholipid syndrome (APS). *J Thromb Haemost*. 2006 Feb;4(2):295-306.
6. Carreras LO, Forastiero RR, Martinuzzo ME. Which are the best biological markers of the antiphospholipid syndrome? *J Autoimmun*. 2000 Sep;15(2):163-72.
7. Cervera R, Piette JC, Font J, Khamashta MA, Shoenfeld Y, *et al.*; Euro-Phospholipid Project Group. Antiphospholipid syndrome: clinical and immunologic manifestations and patterns of disease expression in a cohort of 1,000 patients. *Arthritis Rheum*. 2002 Apr;46(4):1019-27.
8. Linnemann B. Antiphospholipid syndrome - an update. *Vasa*. 2018 Oct;47(6):451-464.
9. Calvo-Alén J, Toloza SM, Fernández M, Bastian HM, *et al.*; LUMINA Study Group. Systemic lupus erythematosus in a multiethnic US cohort (LUMINA). XXV. Smoking, older age, disease

- activity, lupus anticoagulant, and glucocorticoid dose as risk factors for the occurrence of venous thrombosis in lupus patients. *Arthritis Rheum.* 2005 Jul;52(7):2060-8.
10. Janardhan V, Wolf PA, Kase CS, Massaro JM, *et al.* Anticardiolipin antibodies and risk of ischemic stroke and transient ischemic attack: the Framingham cohort and offspring study. *Stroke.* 2004 Mar;35(3):736-41.
 11. Harris EN, Pierangeli SS. Revisiting the anticardiolipin test and its standardization. *Lupus.* 2002;11(5):269-75.
 12. McClain MT, Arbuckle MR, Heinlen LD, Dennis GJ, *et al.* The prevalence, onset, and clinical significance of antiphospholipid antibodies prior to diagnosis of systemic lupus erythematosus. *Arthritis Rheum.* 2004 Apr;50(4):1226-32.
 13. Singh AK. Immunopathogenesis of the antiphospholipid antibody syndrome: an update. *Curr Opin Nephrol Hypertens.* 2001 May;10(3):355-8.
 14. de Laat B, de Groot PG. Autoantibodies directed against domain I of beta2-glycoprotein I. *Curr Rheumatol Rep.* 2011 Feb;13(1):70-6.
 15. Suh-Lailam BB, Cromar A, Davis KW, Tebo AE. APhL antibody ELISA as an alternative to anticardiolipin test for the diagnosis of antiphospholipid syndrome. *Int J Clin Exp Pathol.* 2012;5(3):210-5.
 16. Ofer-Shiber S, Molad Y. Frequency of vascular and pregnancy morbidity in patients with low vs. moderate-to-high titers of antiphospholipid antibodies. *Blood Coagul Fibrinolysis.* 2015 Apr;26(3):261-6.
 17. Meroni PL, Borghi MO, Raschi E, Tedesco F. Pathogenesis of antiphospholipid syndrome: understanding the antibodies. *Nat Rev Rheumatol.* 2011 Jun;7(6):330-9.
 18. de Laat B, Pengo V, Pabinger I, Musial J, *et al.* The association between circulating antibodies against domain I of beta2-glycoprotein I and thrombosis: an international multicenter study. *J Thromb Haemost.* 2009 Nov;7(11):1767-73.
 19. Sanmarco M, Roll P, Gayet S, Oksman F, *et al.* J. Combined search for anti-beta2-glycoprotein I and anticardiolipin antibodies in antiphospholipid syndrome: contribution to diagnosis. *J Lab Clin Med.* 2004 Sep;144(3):141-7.
 20. Otomo K, Atsumi T, Amengual O, Fujieda Y, *et al.* Efficacy of the antiphospholipid score for the diagnosis of antiphospholipid syndrome and its predictive value for thrombotic events. *Arthritis Rheum.* 2012 Feb;64(2):504-12.
 21. Basic QC Practices On-line Course; <http://www.Westgard.com>.
 22. Boscato, LM. and Stuart, MC., 'Heterophilic antibodies: a problem for all immunoassays'. *Clin Chem*, 34, 1988, pp 27–33

19. REVISION IDENTIFIER

Additions or changes to the IFU are indicated by grey highlighting.

20. PRODUCT COMPLAINTS AND TECHNICAL SUPPORT

For a patient/user/third party in the European Union and in countries with similar regulatory regime (Regulation 2017/746/EU on IVD Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorised representative and to your national regulatory authority.

The manufacturer can be contacted through their customer service or technical support team. The contact details can be found below and on the company website: www.diametra.com.

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1. Dusse LM, Silva FD, Freitas LG, Rios DR, *et al.* Antiphospholipid syndrome: a clinical and laboratorial challenge. *Rev Assoc Med Bras* (1992). 2014 Mar-Apr;60(2):181-6.
2. Roggenbuck D, Somma V, Schierack P, Borghi MO, Meroni PL. Autoantibody profiling in APS. *Lupus*. 2014 Oct;23(12):1262-4.
3. Banzato A, Pengo V. Clinical relevance of β_2 -glycoprotein-I plasma levels in antiphospholipid syndrome (APS). *Curr Rheumatol Rep*. 2014 Jun;16(6):424.
4. Amiral J, Peyrafitte M, Dunois C, Vissac AM, Seghatchian J. Anti-phospholipid syndrome: Current opinion on mechanisms involved, laboratory characterization and diagnostic aspects. *Transfus Apher Sci*. 2017 Aug;56(4):612-625.
5. Miyakis S, Lockshin MD, Atsumi T, Branch, *et al.* International consensus statement on an update of the classification criteria for definite antiphospholipid syndrome (APS). *J Thromb Haemost*. 2006 Feb;4(2):295-306.
6. Carreras LO, Forastiero RR, Martinuzzo ME. Which are the best biological markers of the antiphospholipid syndrome? *J Autoimmun*. 2000 Sep;15(2):163-72.
7. Cervera R, Piette JC, Font J, Khamashta MA, Shoenfeld Y, *et al.*; Euro-Phospholipid Project Group. Antiphospholipid syndrome: clinical and immunologic manifestations and patterns of disease expression in a cohort of 1,000 patients. *Arthritis Rheum*. 2002 Apr;46(4):1019-27.
8. Linnemann B. Antiphospholipid syndrome - an update. *Vasa*. 2018 Oct;47(6):451-464.
9. Calvo-Alén J, Toloza SM, Fernández M, Bastian HM, *et al.*; LUMINA Study Group. Systemic lupus erythematosus in a multiethnic US cohort (LUMINA). XXV. Smoking, older age, disease activity, lupus anticoagulant, and glucocorticoid dose as risk factors for the occurrence of venous thrombosis in lupus patients. *Arthritis Rheum*. 2005 Jul;52(7):2060-8.
10. Janardhan V, Wolf PA, Kase CS, Massaro JM, *et al.* Anticardiolipin antibodies and risk of ischemic stroke and transient ischemic attack: the Framingham cohort and offspring study. *Stroke*. 2004 Mar;35(3):736-41.
11. Harris EN, Pierangeli SS. Revisiting the anticardiolipin test and its standardization. *Lupus*. 2002;11(5):269-75.
12. McClain MT, Arbuckle MR, Heinlen LD, Dennis GJ, *et al.* The prevalence, onset, and clinical significance of antiphospholipid antibodies prior to diagnosis of systemic lupus erythematosus. *Arthritis Rheum*. 2004 Apr;50(4):1226-32.
13. Singh AK. Immunopathogenesis of the antiphospholipid antibody syndrome: an update. *Curr Opin Nephrol Hypertens*. 2001 May;10(3):355-8.

14. de Laat B, de Groot PG. Autoantibodies directed against domain I of beta2-glycoprotein I. *Curr Rheumatol Rep*. 2011 Feb;13(1):70-6.
15. Suh-Lailam BB, Cromar A, Davis KW, Tebo AE. APhL antibody ELISA as an alternative to anticardiolipin test for the diagnosis of antiphospholipid syndrome. *Int J Clin Exp Pathol*. 2012;5(3):210-5.
16. Ofer-Shiber S, Molad Y. Frequency of vascular and pregnancy morbidity in patients with low vs. moderate-to-high titers of antiphospholipid antibodies. *Blood Coagul Fibrinolysis*. 2015 Apr;26(3):261-6.
17. Meroni PL, Borghi MO, Raschi E, Tedesco F. Pathogenesis of antiphospholipid syndrome: understanding the antibodies. *Nat Rev Rheumatol*. 2011 Jun;7(6):330-9.
18. de Laat B, Pengo V, Pabinger I, Musial J, *et al.* The association between circulating antibodies against domain I of beta2-glycoprotein I and thrombosis: an international multicenter study. *J Thromb Haemost*. 2009 Nov;7(11):1767-73.
19. Sanmarco M, Roll P, Gayet S, Oksman F, *et al.* J. Combined search for anti-beta2-glycoprotein I and anticardiolipin antibodies in antiphospholipid syndrome: contribution to diagnosis. *J Lab Clin Med*. 2004 Sep;144(3):141-7.
20. Otomo K, Atsumi T, Amengual O, Fujieda Y, *et al.* Efficacy of the antiphospholipid score for the diagnosis of antiphospholipid syndrome and its predictive value for thrombotic events. *Arthritis Rheum*. 2012 Feb;64(2):504-12.
21. Basic QC Practices On-line Course; <http://www.Westgard.com>.
22. Boscato, LM. and Stuart, MC., 'Heterophilic antibodies: a problem for all immunoassays'. *Clin Chem*, 34, 1988, pp 27-33.

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Anti Phospholipid Screen

for routine analysis

Quantitative determination of auto-antibodies against phospholipids in human serum or plasma

IVD

LOT

See external label

2°C 8°C

$\Sigma = 96$ test

REF DKO114

INTENDED USE

Anti Phospholipid Screen is an indirect solid phase immunoassay kit for the quantitative measurement of IgG and IgM class auto-antibodies directed against β -glycoprotein mediated anionic phospholipids in human serum or plasma, including Cardiolipin, Phosphatidyl serine, Phosphatidyl inositol, Phosphatidic acid, Phosphatidyl choline, Lyso-phosphatidyl choline, Phosphatidyl ethanolamine. The assay is intended for in vitro diagnostic use only as an aid in the diagnosis of increased risk of thrombosis in patients with Systemic Lupus Erythematosus (SLE) or similar disorders.

Anti Phospholipid Screen kit is intended for laboratory use only.

1. CLINICAL SIGNIFICANCE

The first study on the anti-phospholipid antibodies began in 1906 when Wasserman introduced a serological test for syphilis. In 1942 it was found the active component that is a phospholipid indicated by the name of Cardiolipin. In the 50's it was observed that a large number of people appeared to be positive for syphilis tests but did not show any evidence of disease. At the beginning the phenomenon was classified as a series of false positive syphilis test, then a more accurate analysis revealed, for this group of patients, a high prevalence of autoimmune disorders including Systemic Lupus Erythematosus (SLE) and Sjögrens syndrome.

The term lupus anticoagulant (LA), used for the first time in 1972, derives from experimental observations in which it was observed an increased risk of thrombosis, paradoxically, with the presence of some anticoagulants factors; the term LA is not totally correct, in fact the disease is present more frequently in patients without lupus and it is associated with thrombosis rather than to abnormal bleeding.

Some years later the role of a cofactor has been investigated, the β 2-glycoprotein I (apolipoprotein H) also said β 2GPI, and its interactions with anionic phospholipids in human serum / plasma. This cofactor is a β -globulin with a molecular weight of 50 kDa that has the concentration of 200 μ g / mL in plasma. The β 2GPI is involved in the regulation of blood coagulation, inhibiting the intrinsic way. β 2GPI in vivo is associated with negatively charged substances such as anionic phospholipids, heparin and lipoproteins. The region that binds phospholipids is in its fifth domain.

The acronym "aPL" (anti-phospholipid antibodies) indicates improperly antibodies directed against phospholipids negatively charged like Cardiolipin (CL), Phosphatidyl serine (PS) Phosphatidyl inositol (PI) and phosphatidic acid (PA); more correctly the term anti-phospholipid antibodies indicate those antibodies directed against the complex

between β 2GPI and anionic phospholipids that can bind to the fifth domain of β 2GPI. Among these, the Cardiolipin is the most commonly used phospholipid as an antigen for determining the aPL with ELISA method. Diagnostic laboratories measure the antibodies directed against the complex between β 2GPI and negatively charged phospholipids, as Phosphatidyl serine (PS) Phosphatidyl inositol (PI) and phosphatidic acid (PA). Some researchers suggest the use of PS instead of Cardiolipin in ELISA assays, for a more precise diagnosis. However, these antibodies against phospholipids are less commonly used, even if their use may increase the clinical sensitivity of patients samples with suspected Anti-phospholipid Syndrome (APS), but it can't replace the determination of autoantibodies anti-Cardiolipin.

2. PRINCIPLE

Anti Phospholipid Screen test is based on the binding of antibodies in human serum directed against the antigenic complex between anionic phospholipids (Cardiolipin, Phosphatidyl serine, Phosphatidyl inositol, Phosphatidic acid, Phosphatidyl choline, Lyso-phosphatidyl choline, Phosphatidyl ethanolamine) and β 2-Glycoprotein; these complexes are coated on the microplate. Any antibody of IgG class or IgM class in calibrators, controls or prediluted patient samples binds to its respective antigen.

After 60 minutes incubation, the microplate is washed with wash buffer for removing non-reactive serum components. An anti-human IgG conjugate solution (Conjugate IgG, reactive 4) or an anti-human IgM conjugate solution (Conjugate IgM, reactive 5) recognize IgG class or IgM class antibodies, respectively, bound to the immobilized antigens.

After a 30 minutes incubation any excess enzyme conjugate which is not specifically bound is washed away with the wash buffer.

A chromogenic substrate solution containing TMB is then dispensed into the wells. After a 15 minutes incubation the color development is stopped by adding the stop solution. The solutions color change into yellow. The amount of color is directly proportional to the concentration of IgG or IgM antibodies present in the original sample.

The concentration of IgG or IgM antibodies in the original sample is calculated through a calibration curve.

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Microplate reader (450 nm, 620-630 nm).

Notes

Store all reagents between 2-8°C in the dark.
Open the bag of reagent 6 (Coated Microplate) only when it is at room temperature and close it immediately after use; once opened, it is stable until expiry date of the kit.

4. WARNINGS

- This kit is intended for in vitro use by professional persons only. Not for internal or external use in Humans or Animals.
- Use appropriate personal protective equipment while working with the reagents provided.
- Follow Good Laboratory Practice (GLP) for handling blood products.
- Material of animal origin used in the preparation of the kit has been obtained from animals certified as healthy and the bovine protein has been obtained from countries not infected by BSE, but these materials should be handled as potentially infectious.
- All human source material used in the preparation of the reagents has been tested and found negative for antibody to HIV 1&2, HbsAg, and HCV. No test method however can offer complete assurance that HIV, HBV, HCV or other infectious agents are absent. Therefore, Calibrators and Controls should be handled in the same

manner as potentially infectious material.

- Some reagents contain small amounts of Sodium Azide (NaN_3) or Proclin 300^R as preservatives. Avoid the contact with skin or mucosa.
- Sodium Azide may be toxic if ingested or absorbed through the skin or eyes; moreover it may react with lead or copper plumbing to form potentially explosive metal azides. If you use a sink to remove the reagents, allow scroll through large amounts of water to prevent azide build-up.
- The TMB Substrate contains an irritant, which may be harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous and corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/ H_2O_2 to directed sunlight, metals or oxidants. Do not freeze the solution.

5. PRECAUTIONS

- Please adhere strictly to the sequence of pipetting steps provided in this protocol. The performance data represented here were obtained using specific reagents listed in this Instruction For Use.
- All reagents should be stored refrigerated at 2-8°C in their original container. Any exceptions are clearly indicated. The reagents are stable until the expiry date when stored and handled as indicated.
- Allow all kit components and specimens to reach room temperature (22-28°C) and mix well prior to use.
- Do not interchange kit components from different lots. The expiry date printed on box and vials labels must be observed. Do not use any kit component beyond their expiry date.
- **WARNING: the conjugate reagent is designed to ensure maximum dose sensitivity and may be contaminated by external agents if not used properly;** therefore, it is recommended to use disposable consumables (tips, bottles, trays, etc.). For divided doses, take the exact amount of conjugate needed and do not re-introduce any waste product into the original bottle. In addition, **for doses dispensed with the aid of automatic and semi-automatic devices,** before using the conjugate, it is advisable to clean the fluid handling system, ensuring that the procedures of washing, deproteinization and decontamination are effective in avoiding contamination of the conjugate; **this procedure is highly recommended when the kit is processed using analyzers which are not equipped with disposable tips.** For this purpose, Dia.Metra supplies a separate decontamination reagent for cleaning needles.
- If you use automated equipment, the user has the responsibility to make sure that the kit has been appropriately tested.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background. To improve the performance of the kit on automatic systems is recommended to increase the number of washes.
- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate
- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop

Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.

- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.
- Maximum precision is required for reconstitution and dispensation of the reagents.
- Samples microbiologically contaminated, highly lipemic or haemolysed should not be used in the assay.
- Plate readers measure vertically. Do not touch the bottom of the wells.

6. PROCEDURE

6.1. Preparation of Calibrators (C₀...C₄)

The Calibrators are ready to use and are mixed, so they have both IgG and IgM antibodies. The Calibrators have the following concentrations:

	C ₀	C ₁	C ₂	C ₃	C ₄
AU/mL	0	5	10	20	80

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6.2. Sample Preparation

Either human serum or plasma samples can be used for the test execution. Test samples should be clear. Contamination by lipemia should be avoided, although it does not interfere with this assay. Specimens may be refrigerated at 2-8°C for up to five days or stored at -20°C up to six months. Avoid repetitive freezing and thawing of serum samples. This may result in variable loss of autoantibody activity. Testing of heat-inactivated sera is not recommended.

All serum and plasma samples have to be diluted 1:100 with sample diluent; for example 10 μ L of sample may be diluted with 990 μ L of sample diluent.

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In concentrated wash solution it is possible to observe the presence of crystals. In this case mix at room temperature until complete dissolution of crystals is observed. For greater accuracy dilute the whole bottle of concentrated wash solution to 500 mL taking care also to transfer the crystals completely, then mix until the crystals are completely dissolved.

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Reagents	Calibrator	Sample/Controls	Blank
Calibrator C ₀ -C ₄	100 μ L		
Controls		100 μ L	
Diluted Sample		100 μ L	
Incubate 60 minutes at room temperature (22-28°C). Remove the content from each well, wash the wells 3 times with 300 μ L of diluted wash solution. Important note: during each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel. Automatic washer: if you use automated equipment, wash the wells at least 5 times.			
Conjugate (IgG or IgM)	100 μ L	100 μ L	
Incubate 30 minutes at room temperature (22-28°C). Remove the content from each well, wash the wells 3 times with 300 μ L of diluted wash solution. Washing: follow the same indications of the previous point.			
TMB Substrate	100 μ L	100 μ L	100 μ L
Incubate 15 minutes in the dark at room temperature (22-28°C).			
Stop Solution	100 μ L	100 μ L	100 μ L
Shake the microplate gently. Read the absorbance (E) at 450 nm against a reference wavelength of 620-630 nm or against Blank within 5 minutes.			

7. QUALITY CONTROL

- The anti-phospholipids Positive Control should be run with every batch of samples to ensure that all reagents and procedures perform properly.
- Because Positive Control is prediluted, it does not control for procedural methods associated with dilution of specimens.
- Additional suitable control sera may be prepared by aliquoting pooled human serum specimens and storing at < -20°C.
- In order for the test results to be considered valid, all of the criteria listed below must be met. If any of these are not met, the test should be considered invalid and the assay repeated:
 - The Positive Control are intended to monitor for substantial reagent failure and they will not ensure precision at the assay cut-off.
 - This test is only valid if the optical density at 450 nm for Positive Control as well as for the Calibrator (C₀-C₄) complies with the respective range indicated on the Quality Control Certificate enclosed to each test kit.

8. RESULTS

For Anti Phospholipid Screen kit a 4-Parameter-Fit with lin-log coordinates for optical density and concentration is the data reduction method of choice. Smoothed-Spline Approximation and log-log coordinates are also suitable. However we recommend using a Lin-Log Plot.

First calculate the averaged optical densities for each calibrator well. Use lin-log graph paper and plot the averaged optical density of each calibrator versus the concentration. Draw the best fitting curve approximating the path of all calibrator points. The calibrator points may

also be connected with straight line segments. The concentration of unknowns may then be estimated from the calibration curve by interpolation.

9. REFERENCE VALUES

In a normal range study with serum samples from healthy blood donors the following ranges have been established with the Anti Phospholipid Screen test:

	IgG (GPL AU/mL)	IgM (MPL AU/mL)
Normal	< 10	< 10
Elevated	≥ 10	≥ 10

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Positive results should be verified concerning the entire clinical status of the patient. Also every decision for therapy should be taken individually. It is recommended that each laboratory establishes its own normal and pathological ranges of seric Ab-Anti-Phospholipid.

10. LIMITATIONS OF PROCEDURE

The presence of immune complexes or other immunoglobulin aggregates in the patient sample may cause an increased level of non-specific binding and produce false positives in this assay

11. PERFORMANCE AND CHARACTERISTICS

11.1. Precision and reproducibility

Precision and reproducibility are evaluated by eight reply of two positive samples by two different runs with two different lots. Dispensing and washing operations were performed manually by an operator.

The results in terms of standard deviation and coefficient of variation were below:

Sample	IgG			
	1		2	
	SD	CV%	SD	CV%
Intra-assay	1.03	5.9	1.31	7.4
Inter-assay	0.26	9.2	5.25	11.7
Sample	IgM			
	1		2	
	SD	CV%	SD	CV%
Intra-assay	0.61	7.6	1.97	5.9
Inter-assay	0.15	7.1	2.98	6.6

11.2. Sensitivity

The clinical sensitivity of Anti Phospholipids Screen IgG assay is 92,3%.

The clinical sensitivity of Anti Phospholipids Screen IgM assay is 68,8%.

11.3. Specificity

The clinical specificity of Anti Phospholipids Screen IgG assay is 84,6%.

The clinical specificity of Anti Phospholipids Screen IgM assay is > 99,9%.

11.4. Detection Limit

The lowest concentration that can be distinguished from Calibrator zero is 0.03 AU/mL for IgG and 0.16 AU/mL for IgM.

12. WASTE MANAGEMENT

Reagents must be disposed off in accordance with local regulations..

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ERROR POSSIBLE CAUSES / SUGGESTIONS**No colorimetric reaction**

- no conjugate pipetted reaction after addition
- contamination of conjugates and/or of substrate
- errors in performing the assay procedure (e.g. accidental pipetting of reagents in a wrong sequence or from the wrong vial, etc.)

Too low reaction (too low ODs)

- incorrect conjugate (e.g. not from original kit)
- incubation time too short, incubation temperature too low

Too high reaction (too high ODs)

- incorrect conjugate (e.g. not from original kit)
- incubation time too long, incubation temperature too high
- water quality for wash buffer insufficient (low grade of deionization)
- insufficient washing (conjugates not properly removed)

Unexplainable outliers

- contamination of pipettes, tips or containers
- insufficient washing (conjugates not properly removed) too high within-run
- reagents and/or strips not pre-warmed to CV% Room Temperature prior to use
- plate washer is not washing correctly (suggestion: clean washer head)
- too high between-run - incubation conditions not constant (time, CV % temperature)
- controls and samples not dispensed at the same time (with the same intervals) (check pipetting order)
- person-related variation

DiaMetra	Packaging Information Sheet	Mod. PIS000-1
-----------------	------------------------------------	----------------------

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BIBLIOGRAFIA

1. Jennette, J.C. and Falk, R.J. Antineutrophil Cytoplasmic Autoantibodies and Associated Diseases: a Review. Am. J. Kidney Dis. 1990, Vol. XV, No. 6: 517 - 529.
2. Gross, W.L. et al. Antineutrophil Cytoplasmic Autoantibody-Associated Diseases: A Rheumatologist's Perspective. Am. J. Kidney Dis. 1991, Vol. XVIII, No. 2: 175 - 179.
3. Wieslander, J. How are Antineutrophil Cytoplasmic Autoantibodies Detected? Am. J. Kidney Dis. 1991, Vol. XVIII, No. 2: 154 - 158.
4. Lesavre, P. Antineutrophil cytoplasmic antibodies antigen specificity. Am. J. Kidney Dis. 1991, Vol. XVIII, No. 2: 159 - 163.
5. Hagen, E.C. et al. Antineutrophil cytoplasmic autoantibodies: a review of the antigens involved, the assays, and the clinical and possible pathogenic consequences. Blood 1993, Vol. 81: 1996 - 2000.
6. Gross, W.L. et al. Immunodiagnostische und immunopathogenetische Bedeutung von Anti-Neutrophilen-Cytoplasma-Antikörpern. Deutsche Medizinische Wochenschrift 1993, Vol. 118: 191 - 199.

Ed. 09/2018

DCM092-5

Dia.Metra S.r.l.

Via Pozzuolo 14,
06038 SPELLO (PG) Italy
Tel. +39-0742-24851
Fax +39-0742-316197
E-mail: info@diametra.com

DCM092-5
Ed. 09/2018

Anti MPO (p-ANCA)

for routine analysis

Quantitative determination of IgG class antibodies against myeloperoxidase (MPO) in human serum or plasma

IVD

LOT

See external label

$\Sigma = 96$ tests

REF DKO092

INTENDED USE

Anti MPO (p-ANCA) kit is an indirect solid phase enzyme immunoassay (ELISA) for the quantitative measurement of IgG class autoantibodies against myeloperoxidase (MPO) in human serum or plasma. The assay is intended for in vitro diagnostic use only as an aid in the diagnosis of certain autoimmune vasculitides such as Wegener's granulomatosis.

Anti MPO (p-ANCA) kit is intended for laboratory use only.

1. CLINICAL SIGNIFICANCE

Anti-neutrophilic-cytoplasm antibodies (ANCA) represent a group of autoantibodies directed towards the cytoplasmic components of the neutrophilic granulocytes and monocytes. The classical methods for the determination of ANCA are the immunofluorescent methods. With these indirect immunofluorescence techniques two main patterns are recognized, a cytoplasmic (c-ANCA) and a perinuclear (p-ANCA) type.

The main antigen for the c-ANCA is the proteinase 3 (PR3), which is a serine proteinase of the present in primary granules. Antibodies of p-ANCA positive sera are mainly directed to myeloperoxidase (MPO). Antibodies to other antigens e.g. lactoferrin, elastase, cathepsin-G and also lysozyme often result in a similar p-ANCA pattern. Beside different untypical variants of p-ANCA IF patterns granulocyte specific antinuclear antibodies (GS-ANA) is indistinguishable from p-ANCA. This makes a clear interpretation and classification of the IF patterns difficult.

Therefore every positive IF-ANCA findings esp. p-ANCA should be differentiated by ELISA techniques using purified antigens. A survey of documented clinical indications of specific ANCA is given in the table below. PR3-ANCA and MPO-ANCA are reliable serologic markers in the diagnostics of vasculitides.

PR3-ANCA is the classical autoantigen in Wegener's granulomatosis with a clinical specificity of more than 95%. c-ANCA is documented to be present in different diseases. Anti-MPO antibodies are highly specific for idiopathic and vasculitis associated crescentic glomerulonephritis and also for classic polyarteritis nodosa, Churg-Strauss syndrome and the polyangitis overlap syndrome without renal involvement. 7-10 With respect to sensitivity, either MPO or PR-3 antibodies were found in 77 to 100% of patients with idiopathic and vasculitis associated crescentic glomerulonephritis. In WG, anti-MPO antibodies were detected only occasionally and generally in patients negative for PR-3 antibodies.

The MPO and PR-3 specific ELISA methods can provide an important confirmatory result for two of the more important of the identified antigens. ELISA is also useful for interpreting "difficult" samples by IFA such as those which exhibit several antibodies simultaneously or those with high background fluorescence.

2. PRINCIPLE

Anti MPO (p-ANCA) test is based on the binding of the antibodies in the sample to the human neutrophil myeloperoxidase coated on the microplates. In the first step the antibodies in calibrators, controls or prediluted patient samples bind to the inner surface of the wells. After a 60 minutes incubation the microplate is washed with a wash buffer to remove the non-reactive serum components. In the second step an anti-human-IgG horseradish peroxidase conjugated solution recognizes the IgG class antibodies bound to the immobilized antigens. After a 30 minutes incubation any excess of enzyme conjugate which is not specifically bound is washed away with the wash buffer. Then a chromogenic substrate solution containing TMB is dispensed into the wells. After 15 minutes of incubation the color development is stopped by adding the stop solution. The solutions color change into yellow. The amount of color is directly proportional to the concentration of IgG antibodies present in the original sample. The concentration of IgG antibodies in the sample is calculated through a calibration curve.

3. REAGENTS, MATERIALS AND INSTRUMENTATION

3.1. Reagents and materials supplied in the kit

1. Calibrators (5 vials, 1.2 mL each)

Phosphate buffer 0,1M, NaN₃ < 0,1%, human serum

CAL0

REF DCE002/09206-0

CAL1

REF DCE002/09207-0

CAL2

REF DCE002/09208-0

CAL3

REF DCE002/09209-0

CAL4

REF DCE002/09210-0

2. Controls (2 vials, 1.2 mL each, ready for use)

Phosphate buffer 0,1M, NaN₃ < 0,1%, human serum

Negative control

REF DCE045/09201-0

Positive control

REF DCE045/09202-0

3. Sample Diluent (1 vial, 100 mL)

Phosphate buffer 0,1M, NaN₃ < 0,1%

REF DCE053-0

4. Conjugate (1 vial, 15 mL)

Anti h-IgG conjugated with horseradish peroxidase (HRP),

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Microplate coated with human neutrophil myeloperoxidase

REF DCE002/09203-0

6. TMB Substrate (1 vial, 15 mL)

H₂O₂-TMB 0.26 g/L (avoid any skin contact)

REF DCE004-0

7. Stop Solution (1 vial, 15 mL)

Sulphuric acid 0.15M (avoid any skin contact)

REF DCE005-0

8. 10X Conc. Wash Solution (1 vial, 50 mL)
Phosphate buffer 0,2M, pH 7.4 REF DCE054-0

3.2. Reagents necessary not supplied

Distilled water.

3.3. Auxiliary materials and instrumentation

Automatic dispenser.

Microplate reader (450 nm, 620-630 nm).

Notes

Store all reagents between 2-8°C in the dark.

Open the bag of reagent 5 (Coated Microplate) only when it is at room temperature and close it immediately after use; once opened, it is stable until expiry date of the kit.

4. WARNINGS

- This kit is intended for in vitro use by professional persons only. Not for internal or external use in Humans or Animals.
- Use appropriate personal protective equipment while working with the reagents provided.
- Follow Good Laboratory Practice (GLP) for handling blood products.
- All human source material used in the preparation of the reagents has been tested and found negative for antibody to HIV 1&2, HbsAg, and HCV. No test method however can offer complete assurance that HIV, HBV, HCV or other infectious agents are absent. Therefore, Calibrators and Controls should be handled in the same manner as potentially infectious material.
- Material of animal origin used in the preparation of the kit has been obtained from animals certified as healthy and the bovine protein has been obtained from countries not infected by BSE, but these materials should be handled as potentially infectious.
- Some reagents contain small amounts of Sodium Azide (NaN_3) or Proclin 300^R as preservatives. Avoid the contact with skin or mucosa.
- Sodium Azide may be toxic if ingested or absorbed through the skin or eyes; moreover it may react with lead or copper plumbing to form potentially explosive metal azides. If you use a sink to remove the reagents, allow scroll through large amounts of water to prevent azide build-up.
- The TMB Substrate contains an irritant, which may be harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous and corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/ H_2O_2 to directed sunlight, metals or oxidants. Do not freeze the solution.

5. PRECAUTIONS

- Please adhere strictly to the sequence of pipetting steps provided in this protocol. The performance data represented here were obtained using specific reagents listed in this Instruction For Use.
- All reagents should be stored refrigerated at 2-8°C in their original container. Any exceptions are clearly indicated. The reagents are stable until the expiry date when stored and handled as indicated.
- Allow all kit components and specimens to reach room temperature (22-28°C) and mix well prior to use.
- Do not interchange kit components from different lots. The expiry date printed on box and vials labels must be observed. Do not use any kit component beyond their expiry date.

- **WARNING: the conjugate reagent is designed to ensure maximum dose sensitivity and may be contaminated by external agents if not used properly;** therefore, it is recommended to use disposable consumables (tips, bottles, trays, etc.). For divided doses, take the exact amount of conjugate needed and do not re-introduce any waste product into the original bottle. In addition, **for doses dispensed with the aid of automatic and semi-automatic devices,** before using the conjugate, it is advisable to clean the fluid handling system, ensuring that the procedures of washing, deproteinization and decontamination are effective in avoiding contamination of the conjugate; **this procedure is highly recommended when the kit is processed using analyzers which are not equipped with disposable tips.**

For this purpose, Dia.Metra supplies a separate decontamination reagent for cleaning needles.

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6. PROCEDURE

6.1. Preparation of the Calibrators (C₀...C₄)

Since no international reference preparation for Anti-myeloperoxidase antibodies is available, the assay system is calibrated in relative arbitrary units. The Calibrators are to use and have the following concentration:

	C ₀	C ₁	C ₂	C ₃	C ₄
AU/mL	0	10	20	40	160

Once opened, the Calibrators are stable 6 months at 2-8°C.

6.2. Preparation of the Sample

Either human serum or plasma samples can be used for the test execution. Test samples should be clear. Contamination by lipemia is best avoided, but does not interfere with this assay.

Specimens may be refrigerated at 2-8°C for up to five days or stored at -20°C up to six months. Avoid repetitive freezing and thawing of serum samples. This may result in variable loss of autoantibody activity.

Testing of heat-inactivated sample is not recommended.
All serum and plasma samples must be prediluted 1:100 with sample diluents; for example 10 L of sample should be diluted with 990 L of sample diluent.
 The Controls are ready to use.

6.3. Preparation of the Wash Solution

Dilute the content of each vial of the "10X Conc. Wash Solution" with distilled water to a final volume of 500 mL prior to use. For smaller volumes respect the 1:10 dilution ratio. The diluted wash solution is stable for 30 days at 2-8°C.

In concentrated wash solution is possible to observe the presence of crystals; in this case mix at room temperature until the complete dissolution of crystals; for greater accuracy, dilute the whole bottle of concentrated wash solution to 500 mL, taking care to transfer completely the crystals, then mix until crystals are completely dissolved.

6.4. Procedure

- **Allow all reagents to reach room temperature (22-28°C) for at least 30 minutes.** At the end of the assay store immediately the reagents at 2-8°C: avoid long exposure to room temperature.
- Unused coated microwell strips should be released securely in the foil pouch containing desiccant and stored at 2-8°C.
- To avoid potential microbial and/or chemical contamination, unused reagents should never be transferred into the original vials.
- As it is necessary to perform the determination in duplicate in order to improve accuracy of the test results, prepare two wells for each point of the calibration curve (C₀-C₄), two for each Control, two for each sample, one for Blank.

Reagent	Calibrator	Sample/ Controls	Blank
Calibrator C ₀ -C ₄	100 µL		
Controls		100 µL	
Diluted Sample		100 µL	
Incubate 60 minutes at room temperature (22-28°C). Remove the contents from each well, wash the wells 3 times with 300 µL of diluted wash solution. Important note: during each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel. Automatic washer: if you use automated equipment, wash the wells at least 5 times.			
Conjugate	100 µL	100 µL	
Incubate 30 minutes at room temperature (22-28°C). Remove the contents from each well, wash the wells 3 times with 300 µL of diluted wash solution. Washing: follow the same indications of the previous point.			
TMB Substrate	100 µL	100 µL	100 µL
Incubate 15 minutes in the dark at room temperature (22-28°C).			

Stop Solution	100 µL	100 µL	100 µL
Shake the microplate gently. Read the absorbance (E) at 450 nm against a reference wavelength of 620-630 nm or against Blank within 5 minutes.			

7. QUALITY CONTROL

- The MPO IgG Positive and Negative Control should be run with every batch of samples to ensure that all reagents and procedures perform properly.
- Because Positive and the Negative Control are prediluted, they do not control for procedural methods associated with dilution of specimens.
- Additional suitable control sera may be prepared by aliquoting pooled human serum specimens and storing at < -20°C.
- In order for the test results to be considered valid, all of the criteria listed below must be met. If any of these are not met, the test should be considered invalid and the assay repeated:
 - The absorbance of the prediluted MPO IgG Positive must be greater than the absorbance of the prediluted Negative Control.
 - The Negative and Positive Control are intended to monitor for substantial reagent failure and they will not ensure precision at the assay cutoff.
 - This test is only valid if the optical density at 450 nm for Positive Control (1) and Negative Control (2) as well as for the Calibrator S0-S5 complies with the respective range indicated on the Quality Control Certificate enclosed to each test kit: If any of these criteria is not met, the results are invalid and the test should be repeated.

8. RESULTS

For Anti MPO (p-ANCA) test a 4-Parameter-Fit with lin-log coordinates for optical density and concentration is the data reduction method of choice. Smoothed-Spline Approximation and log-log coordinates are also suitable. Recommended Lin-Log Plot

First calculate the averaged optical densities for each calibrator well. Use lin-log graph paper and plot the averaged optical density of each calibrator versus the concentration. Draw the best fitting curve approximating the path of all calibrator points. The calibrator points may also be connected with straight line segments. The concentration of unknowns may then be estimated from the calibration curve by interpolation.

9. REFERENCE VALUES

In a normal range study with serum samples from healthy blood donors the following ranges have been established with the Anti MPO (p-ANCA) tests:

	Anti MPO (p-ANCA)
Cut-off	20 AU/mL

Please pay attention to the fact that the determination of a range of expected values for a "normal" population in a given method is dependent on many factors, such as specificity and sensitivity of the method used and type of population under investigation. Therefore each laboratory should consider the range given by the Manufacturer as a general indication and produce their own range of expected values based on the indigenous population where the laboratory works.

Positive results should be verified concerning the entire clinical status of the patient. Also every decision for therapy should be taken individually. It is recommended that each

laboratory establishes its own normal and pathological ranges of serum Anti-MPO.

10. LIMITATIONS OF PROCEDURE

The presence of immune complexes or other immunoglobulin aggregates in the patient sample may cause an increased level of non-specific binding and produce false positives in this assay.

11. PERFORMANCE AND CHARACTERISTICS

11.1. Specificity

Comparison test against a commercial reference kit, performed on 32 sera (8 of them positive sera and 24 negative sera) showed a 100% specificity.

11.2. Sensitivity

Comparison test against a commercial reference kit, performed on 32 sera (8 of them positive sera and 24 negative sera) showed a 100% sensitivity.

11.3. Detection Limit

The lowest concentration of anti MPO antibodies that can be distinguished from the Calibrator 0 is about 0.73 AU/mL with a confidence limit of 98%.

11.4. Precision and reproducibility

11.4.1. Intra-Assay

Within run variation was determined by replicate the measurements of three different sera with values in the range of the calibration curve. The within assay variability is $\leq 6.1\%$.

11.4.2. Inter-Assay

Between run variation was determined by replicate the measurements of two different control sera with different lots of kits and/or different mix of lots of reagents. The between assay variability is $\leq 12.3\%$.

12. WASTE MANAGEMENT

Reagents must be disposed off in accordance with local regulations.

BIBLIOGRAPHY

1. Jennette, J.C. and Falk, R.J. Antineutrophil Cytoplasmic Autoantibodies and Associated Diseases: a Review. Am. J. Kidney Dis. 1990, Vol. XV, No. 6: 517 - 529.
2. Gross, W.L. et al. Antineutrophil Cytoplasmic Autoantibody-Associated Diseases: A Rheumatologist's Perspective. Am. J. Kidney Dis. 1991, Vol. XVIII, No. 2: 175 - 179.
3. Wieslander, J. How are Antineutrophil Cytoplasmic Autoantibodies Detected? Am. J. Kidney Dis. 1991, Vol. XVIII, No. 2: 154 - 158.
4. Lesavre, P. Antineutrophil cytoplasmic antibodies antigen specificity. Am. J. Kidney Dis. 1991, Vol. XVIII, No. 2: 159 - 163.
5. Hagen, E.C. et al. Antineutrophil cytoplasmic autoantibodies: a review of the antigens involved, the assays, and the clinical and possible pathogenic consequences. Blood 1993, Vol. 81: 1996 - 2000.
6. Gross, W.L. et al. Immunodiagnostische und immunopathogenetische Bedeutung von Anti-Neutrophilen-Cytoplasma-Antikörpern. Deutsche Medizinische Wochenschrift 1993, Vol. 118: 191 - 199.

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Anti-LC-1 ELISA (IgG)

Test instruction


ORDER NO.	ANTIBODIES AGAINST	IG CLASS	SUBSTRATE	FORMAT
EA 1307-9601 G	LC-1	IgG	Ag-coated microplate wells	96 x 01 (96)

Indications: The ELISA test kit provides a semiquantitative in vitro assay for human autoantibodies of the IgG class against liver cytosol antigen type 1 (LC-1) in serum or plasma for the diagnosis of increase in transaminases for unclear reasons and suspected autoimmune hepatitis.

Application: The determination of autoantibodies against LC-1 is another important supplementing serological parameter in the diagnosis of autoimmune liver diseases. A reliable diagnosis of AIH is indispensable since untreated AIH rapidly turns into liver cirrhosis.

Principle of the test: The test kit contains microplate strips, each with 8 break-off reagent wells coated with LC-1. In the first reaction step, diluted patient samples are incubated in the wells. If samples are positive, specific IgG antibodies (also IgA and IgM) bind to the antigens. To detect the bound antibodies, a second incubation is carried out using an enzyme-labelled anti-human IgG (enzyme conjugate) catalysing a colour reaction.

Contents of the test kit:

Component	Colour	Format	Symbol
1. Microplate wells coated with antigens 12 microplate strips each containing 8 individual break-off wells in a frame, ready for use	---	12 x 8	
2. Calibrator (IgG, human), ready for use	dark red	1 x 2.0 ml	
3. Positive control (IgG, human), ready for use	blue	1 x 2.0 ml	
4. Negative control (IgG, human), ready for use	green	1 x 2.0 ml	
5. Enzyme conjugate peroxidase-labelled anti-human IgG (rabbit), ready for use	green	1 x 12 ml	
6. Sample buffer ready for use	light blue	1 x 100 ml	
7. Wash buffer 10x concentrate	colourless	1 x 100 ml	
8. Chromogen/substrate solution TMB/H ₂ O ₂ , ready for use	colourless	1 x 12 ml	
9. Stop solution 0.5 M sulphuric acid, ready for use	colourless	1 x 12 ml	
10. Test instruction	---	1 booklet	
11. Quality control certificate	---	1 protocol	
 Lot description			 Storage temperature
 In vitro diagnostic medical device			 Unopened usable until

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- **Coated wells:** Ready for use. Tear open the resealable protective wrapping of the microplate at the recesses above the grip seam. Do not open until the microplate has reached room temperature to prevent the individual strips from moistening. Immediately replace the remaining wells of a partly used microplate in the protective wrapping and tightly seal with the integrated grip seam (Do not remove the desiccant bag).
Once the protective wrapping has been opened for the first time, the wells coated with antigens can be stored in a dry place and at a temperature between +2°C and +8°C for 4 months.
- **Calibrator and controls:** Ready for use. The reagents must be mixed thoroughly before use.
- **Enzyme conjugate:** Ready for use. The enzyme conjugate must be mixed thoroughly before use.
- **Sample buffer:** Ready for use.
- **Wash buffer:** The wash buffer is a 10x concentrate. If crystallisation occurs in the concentrated buffer, warm it to +37°C and mix well before diluting. The quantity required should be removed from the bottle using a clean pipette and diluted with deionised or distilled water (1 part reagent plus 9 parts distilled water).
For example: For 1 microplate strip, 5 ml concentrate plus 45 ml water.
The working strength wash buffer is stable for 4 weeks when stored at +2°C to +8°C and handled properly.
- **Chromogen/substrate solution:** Ready for use. Close the bottle immediately after use, as the contents are sensitive to light . The chromogen/substrate solution must be clear on use. Do not use the solution if it is blue coloured.
- **Stop solution:** Ready for use.

Warning: The calibrator and controls of human origin have tested negative for HBsAg, anti-HCV, anti-HIV-1 and anti-HIV-2. Nonetheless, all materials should be treated as being a potential infection hazard and should be handled with care. Some of the reagents contain the agent sodium azide in a non-declarable concentration. Avoid skin contact.

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Sample dilution: Patient samples are diluted **1:101** in sample buffer. For example: dilute 10 µl sample in 1.0 ml sample buffer and mix well by vortexing (sample pipettes are not suitable for mixing).

NOTE: The calibrator and controls are prediluted and ready for use, do not dilute them.

Incubation

Sample incubation: (1st step)

Transfer 100 µl of the calibrator, positive or negative control or diluted patient samples into the individual microplate wells according to the pipetting protocol. Incubate for **30 minutes** at room temperature (+18°C to +25°C).

Washing:

Manual: Empty the wells and subsequently wash 3 times using 300 µl of working strength wash buffer for each wash.

Automatic: Wash the reagent wells 3 times with 450 µl of working strength wash buffer (program setting: e.g. TECAN Columbus Washer "Overflow Mode").

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Conjugate incubation: (2nd step)

Pipette 100 µl of enzyme conjugate (peroxidase-labelled anti-human IgG) into each of the microplate wells. Incubate for **30 minutes** at room temperature (+18°C to +25°C).

Washing:

Empty the wells. Wash as described above.

Substrate incubation: (3rd step)

Pipette 100 µl of chromogen/substrate solution into each of the microplate wells. Incubate for **15 minutes** at room temperature (+18°C to +25°C) protect from direct sunlight.

Stopping:

Pipette 100 µl of stop solution into each of the microplate wells in the same order and at the same speed as the chromogen/substrate solution was introduced.

Measurement:

Photometric measurement of the colour intensity should be made at a wavelength of 450 nm and a reference wavelength between 620 nm and 650 nm **within 30 minutes of adding the stop solution**. Prior to measuring, slightly shake the microplate to ensure a homogeneous distribution of the solution.

Pipetting protocol

	1	2	3	4	5	6	7	8	9	10	11	12
A	C	P 6	P 14	P 22								
B	pos.	P 7	P 15	P 23								
C	neg.	P 8	P 16	P 24								
D	P 1	P 9	P 17									
E	P 2	P 10	P 18									
F	P 3	P 11	P 19									
G	P 4	P 12	P 20									
H	P 5	P 13	P 21									

The above pipetting protocol is an example of the semiquantitative determination of antibodies in 24 patient samples (P 1 to P 24).

Calibrator (C), positive (pos.) and negative (neg.) control as well as the patient samples have been incubated in one well each. The reliability of the ELISA test can be improved by duplicate determinations of each sample.

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Calculation of results

Semiquantitative: Results can be evaluated semiquantitatively by calculating a ratio of the extinction value of the control or patient sample over the extinction value of calibrator. Use the following formula to calculate the ratio:

$$\frac{\text{Extinction of the control or patient sample}}{\text{Extinction of calibrator}} = \text{Ratio}$$

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Ratio <1.0: **negative**
Ratio ≥1.0: **positive**

For duplicate determinations the mean of the two values should be taken. If the two values deviate substantially from one another, EUROIMMUN recommends retesting the samples.

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Test characteristics

Calibration: As no international reference serum exists for antibodies against LC-1, the results are given as a ratio. This represents a relative measurement for the concentration of antibodies in serum or plasma.

For every group of tests performed, the extinction values of the calibrator and the ratio of the positive and negative control sera must lie within the limits stated for the relevant test kit lot. A quality control certificate containing these reference values is included. If the values specified for the control sera are not achieved, the test results may be inaccurate and the test should be repeated.

The binding activity of the antibodies and the activity of the enzyme used are temperature-dependent. It is therefore recommended using a thermostat in all three incubation steps. The higher the room temperature (+18°C to +25°C) during the incubation steps, the greater the extinction values will be. Corresponding variations apply also to the incubation times. However, the calibrator is subject to the same influences, so that such variations will be largely compensated in the calculation of the result.

Antigen: The microplate wells were coated with recombinant LC-1. The corresponding human cDNA was expressed in insect cells using a baculovirus vector. The specific target antigen of anti-LC-1 antibodies was identified in 1999 as the enzyme formiminotransferase cyclodeaminase (Lapierre et al.).

Detection limit: The lower detection limit is defined as the mean value of an analyte-free sample plus three times the standard deviation and is the smallest detectable antibody titer. The lower detection limit of the Anti-LC-1 ELISA (IgG) is ratio 0.05.

Cross reactivity: This ELISA specifically detects autoantibodies of class IgG against LC-1. When investigating patient sera for autoantibodies against SLA (n = 2) and LKM (n = 8) no cross reactions were found.

Interference: Haemolytic, lipaemic and icteric samples showed no influence on the result up to a concentration of 10 mg/ml for haemoglobin, 20 mg/ml for triglycerides and 0.4 mg/ml for bilirubin in this ELISA.

Reproducibility: The reproducibility of the test was investigated by determining the intra- and inter-assay coefficients of variation (CV) using 3 samples. The intra-assay CVs are based on 20 determinations and the inter-assay CVs on 4 determinations performed in 6 different test runs.

<i>Intra-assay variation, n = 20</i>		
Sample	Mean value (Ratio)	CV (%)
1	3.1	4.9
2	4.7	2.7
3	8.0	2.4

<i>Inter-assay variation, n = 4 x 6</i>		
Sample	Mean value (Ratio)	CV (%)
1	2.5	10.8
2	3.8	10.2
3	6.7	9.4

Prevalence and specificity: Sera of 93 patients with autoimmune hepatitis, 183 patients with hepatitis A, hepatitis B, toxic liver diseases, steatohepatitis or primary biliary cirrhosis (PBC) and 200 healthy blood donors were examined with the EUROIMMUN Anti-LC-1 ELISA (IgG). The prevalence of antibodies against LC-1 in autoimmune hepatitis was 5.4% with a specificity of 100%.

Reference range: The levels of the anti-LC-1 antibodies (IgG) were analysed with this EUROIMMUN ELISA in a panel of 200 healthy blood donors. With a cut-off ratio of 1.0, all blood donors were anti-LC-1 negative.

Literature references

1. Martini E, Abauf N, Cavalli F, Durand V, Johanet C, Homberg JC. **Antibody to liver cytosol (anti-LC-1) in patients with autoimmune chronic active hepatitis type 2.** Hepatology 8: 1662-1666 (1988).
2. Lapierre P, Hajoui O, Homberg JC, Alvarez F. **Formiminotransferase cyclodeaminase is an organ-specific autoantigen recognized by sera of patients with autoimmune hepatitis.** Gastroenterology 116: 643-649 (1999).

Anti-SLA/LP ELISA (IgG)


Test instruction

ORDER NO.	ANTIBODIES AGAINST	IG-CLASS	SUBSTRATE	FORMAT
EA 1302-9601 G	SLA/LP (soluble liver antigen/ liver-pancreas antigen)	IgG	Ag-coated microplate wells	96 x 01 (96)

Indications: Increase in transaminases for unclear reasons, suspected autoimmune hepatitis.

Principles of the test: The ELISA test kit provides a semiquantitative or quantitative in vitro assay for human autoantibodies of the IgG class against SLA/LP (soluble liver antigen/liver-pancreas antigen) in serum or plasma. The test kit contains microtiter strips each with 8 break-off reagent wells coated with SLA/LP. In the first reaction step, diluted patient samples are incubated in the wells. In the case of positive samples, specific IgG antibodies (also IgA and IgM) will bind to the antigens. To detect the bound antibodies, a second incubation is carried out using an enzyme-labelled anti-human IgG (enzyme conjugate) which is capable of promoting a colour reaction.

Contents of the test kit:

Component	Colour	Format	Symbol
1. Microplate wells, coated with antigens: 12 microplate strips each containing 8 individual break-off wells in a frame, ready for use	---	12 x 8	
2. Calibrator 1 200 RU/ml (IgG, human), ready for use	dark red	1 x 2.0 ml	
3. Calibrator 2 20 RU/ml (IgG, human), ready for use	red	1 x 2.0 ml	
4. Calibrator 3 2 RU/ml (IgG, human), ready for use	light red	1 x 2.0 ml	
5. Positive control (IgG, human), ready for use	blue	1 x 2.0 ml	
6. Negative control (IgG, human), ready for use	green	1 x 2.0 ml	
7. Enzyme conjugate peroxidase-labelled anti-human IgG (rabbit), ready for use	green	1 x 12 ml	
8. Sample buffer ready for use	light blue	1 x 100 ml	
9. Wash buffer 10x concentrate	colourless	1 x 100 ml	
10. Chromogen/substrate solution TMB/H ₂ O ₂ , ready for use	colourless	1 x 12 ml	
11. Stop solution 0.5 M sulphuric acid, ready for use	colourless	1 x 12 ml	
12. Test instruction	---	1 booklet	
13. Protocol with target values	---	1 protocol	
 Lot			 storage temperature
 In vitro determination			 unopened usable until

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Waste disposal: Patient samples, calibrators, controls and incubated microplate strips should be handled as infectious waste. All reagents should be disposed of according to official regulations.

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Note: All reagents must be brought to room temperature (+18°C to +25°C) around 30 minutes before use. After first use, the reagents are stable until the indicated expiry date if stored at +2°C to +8°C and protected from contamination, unless stated otherwise below.

- **Coated wells:** Ready for use. Tear open the resealable protective wrapping of the microplate at the recesses above the grip seam. Do not open until the microplate has reached room temperature to prevent the individual strips from moistening. Immediately replace the remaining wells of a partly used microplate in the protective wrapping and tightly seal with the integrated grip seam (Do not remove the desiccant bag).
Once the protective wrapping has been opened for the first time, the wells coated with antigens can be stored in a dry place and at a temperature between +2°C and +8°C for a minimum of 4 months.
- **Calibrators and controls:** Ready for use. The reagents must be mixed thoroughly before use.
- **Enzyme conjugate:** Ready for use. The enzyme conjugate must be mixed thoroughly before use.
- **Sample buffer:** Ready for use.
- **Wash buffer:** The wash buffer is a 10x concentrate. If crystallization occurs in the concentrated buffer, warm it to 37°C and mix well before diluting. The quantity required should be removed from the bottle using a clean pipette and diluted with deionized or distilled water (1 part reagent plus 9 parts distilled water).
For example, for 1 microplate: 5 ml concentrate plus 45 ml water.
The ready-to-use diluted wash buffer is stable for 4 weeks when stored at +2°C to +8°C and handled properly.
- **Chromogen/substrate solution:** Ready for use. Close the bottle immediately after use, as the contents are sensitive to light. The chromogen/substrate solution must be clear on use. Do not use the solution if it is blue coloured.
- **Stop solution:** Ready for use.

Warning: Calibrators and controls used have tested negative for HBsAg, anti-HCV, anti-HIV-1 and anti-HIV-2 using enzyme immunoassays and indirect immunofluorescence methods. Nonetheless, all materials should be treated as being a potential infection hazard and should be handled with care. Some of the reagents contain the toxic agent sodium azide. Avoid contact with the skin.

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Sample dilution: Patient samples are diluted **1:101** sample buffer. For example: dilute 10 µl serum in 1.0 ml sample buffer and mix well by vortexing (sample pipettes are not suitable for mixing).

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Sample incubation:
(1. step)

Transfer 100 µl calibrators, positive and negative controls or diluted patient samples into the individual microplate wells according to the pipetting protocol. Incubate for **30 minutes** at room temperature (+18°C to +25°C).

Washing:

Manual: Empty the wells and subsequently wash 3 times using 300 µl of working strength wash buffer for each wash.

Automatic: Wash reagent wells 3 times with 400 µl working strength wash buffer (program setting: e.g. TECAN Columbus Washer "Overflow Modus").

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Note: Residual liquid (> 10 µl) remaining in the reagent wells after washing can interfere with the substrate and lead to falsely low extinction values. Insufficient washing (e.g., less than 3 wash cycles, too small wash buffer volumes, or too short reaction times) can lead to falsely high extinction values.

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Conjugate incubation:
(2. step)

Pipette 100 µl of enzyme conjugate (peroxidase-labelled anti-human IgG) into each of the microplate wells. Incubate for **30 minutes** at room temperature (+18°C bis 25°C).

Washing:

Empty the wells. Wash as described above.

Substrate incubation:
(3. step)

Pipette 100 µl of chromogen/substrate solution into each of the microplate wells. Incubate for **15 minutes** at room temperature (+18°C bis 25°C) protect from direct sunlight.

Stopping the reaction:

Pipette 100 µl of stop solution into each of the microplate wells in the same order and at the same speed as the chromogen/substrate solution was introduced.

Measurement:

Photometric measurement of the colour intensity should be made at a wavelength of 450 nm and a reference wavelength of between 620 nm and 650 nm **within 30 minutes of adding the stop solution**. Prior to measuring, slightly shake the micro-plate to ensure a homogeneous distribution of the solution.

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The wells can be broken off individually from the strips. This makes it possible to adjust the number of test substrates used to the number of samples to be examined and minimizes reagent wastage.

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Calculation of results


Semiquantitative: Results can be evaluated semiquantitatively by calculating a ratio of the extinction value of the control or patient sample over the extinction value of calibrator 2. Calculate the ratio according the following formula:

$$\frac{\text{Extinction of the control or patient sample}}{\text{Extinction of calibrator 2}} = \text{Ratio}$$

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If the extinction of a serum sample lies above the value of calibrator 1 (200 RU/ml). The result should be given as ">200 RU/ml". It is recommended that the sample be re-tested at a dilution of 1:400. The result in IU/ml read from the calibration curve for this sample must then be multiplied by a factor of 4.

The upper limit of the normal range (**cut-off**) recommended by EUROIMMUN is 20 relative units (RU) /ml. EUROIMMUN recommends interpreting results as follows:

<20 RU/ml:	negative
≥20 RU/ml:	positive

For duplicate determinations the mean of the two values should be taken. If the two values deviate substantially from one another the sample should be retested.

For diagnosis, the clinical symptoms of the patient should always be taken into account alongside the serological results.

Test characteristics

Calibration: As no international reference serum exists for antibodies against SLA/LP, the calibration is performed in relative units (RU).


For every group of tests performed, the extinction values of the calibrators and the relative units and/or ratios determined for the positive and negative controls must lie within the limits stated for the relevant test kit lot. A protocol containing these target values is included. If the values specified for the controls are not achieved, the test results may be inaccurate and the test should be repeated.

The activity of the enzyme used is temperature-dependent and the extinction values may vary if a thermostat is not used. The higher the room temperature (+18°C to +25°C) during substrate incubation, the greater will be the extinction values. Corresponding variations apply also to the incubation times. However, the calibrators are subject to the same influences, with the result that such variations will be largely compensated in the calculation of the result.

Antigen: The microplate wells were coated with recombinant SLA/LP. The corresponding human cDNA was expressed in *E. coli*. Identification of SLA/LP at the DNA level succeeded in 1998 by cloning the target antigen (A. W. Lohse, University Hospital of Mainz). SLA/LP is probably a cytoplasmic molecule with a molecular weight of 50 kDa which is involved in the regulation of the protein biosynthesis (an UGA-suppressor-tRNA associated protein). Previous descriptions of SLA as being the liver cytochromes 8 and 18 or glutathione-S-transferase were apparently wrong.

Linearity: The linearity of the test was investigated by assaying serial dilutions of patient samples with high antibody concentrations. The chart below shows the typical linearity of samples on the basis of 3 patient sera. The Anti-SLA/LP ELISA (IgG) is linear in the measurement range 2 - 200 RU/ml.

Detection limit: The detection limit is defined as a value of three times the standard deviation of an analyte-free sample and is the smallest detectable antibody titer. The detection limit of the Anti-SLA/LP ELISA (IgG) is approximately 1 RU/ml.

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Interference: Haemolytic, lipaemic and icteric samples showed no influence at the result up to a concentration of 10 mg/ml for hemoglobin, 20 mg/ml for triglycerides and 0.4 mg/ml for bilirubin in this ELISA.

Reproducibility: The reproducibility of the test was investigated by determining the intra- and inter-assay coefficients of variation (CV) using 3 sera with values at different points on the calibration curve. The intra-assay CVs are based on 20 determinations and the inter-assay CVs on 4 determinations performed on 6 different days.

<i>Intra-Assay Variation, n = 20</i>		
Serum	Mean value (RU/ml)	CV (%)
1	41	2.6
2	110	2.3
3	154	3.7

<i>Inter-Assay Variation, n = 4 x 6</i>		
Serum	Mean value (RU/ml)	CV (%)
1	44	4.4
2	118	3.6
3	169	3.8

Prevalence and specificity: Sera from 454 patients with autoimmune hepatitis, 165 patients with other liver diseases and 200 healthy blood donors were examined with the Anti-SLA/LP ELISA. The prevalence anti-SLA/LP antibodies in non-Japanese AIH patients was in the range between 15% and 19%. The test showed a specificity of 100%.

Patient group (n=819)	Origin of sera	n	Anti-SLA/LP-positive
Autoimmune hepatitis	Prof. Lohse, Univ. Mainz	108	21 (19%)
	Brasilian panel	154	25 (16%)
	Japanese panel	43	2 (5%)
	American panel	149	23 (15%)
AIH (anti-LKM-1 positive)	Dr. Gruber, Univ. München	18	0
Primary biliary cirrhosis	Dr. Gruber, Univ. München PD Dr. Wick, Klinikum Großhadern	30	0
Hepatitis-B infection	Labor Dr. Stöcker, Lübeck	40	0
Hepatitis-C infection	Prof. Lohse, Univ. Mainz Labor Dr. Stöcker, Lübeck	39	0
Steatohepatitis	Prof. Lohse, Univ. Mainz	24	0
Toxic liver damage	Prof. Lohse, Univ. Mainz	14	0
Blood donors	Med. Univ. Lübeck	200	0

Reference range: The levels of the anti-SLA/LP antibodies (IgG) were analyzed with this EUROIMMUN ELISA in a panel of 200 healthy blood donors. With a cut-off of 20 RU/ml, all blood donors were anti-SLA/LP negative.

Clinical significance

The determination of autoantibodies against soluble liver antigen/liver-pancreas antigen (SLA/LP) is a new and important component in the diagnostics of autoimmune diseases of the liver.

Autoimmune diseases of the liver include

- autoimmune hepatitis (AIH),
- primary biliary liver cirrhosis (PBC) and
- primary sclerosing cholangitis (PSC).

Mainly women are affected by **autoimmune hepatitis (AIH)**, earlier designations: lupoid hepatitis, chronically active hepatitis). The disease manifests itself through increased levels of bilirubin, liver enzymes and immunoglobulins, through typical histological changes (liver biopsies show necrosis of the parenchymal cells with lymphocyte and plasma cell infiltrates) and through the occurrence of various autoantibodies. The disease can occur from infancy until the old age, however, it affects most frequently young adults.

The incidence of AIH in western Europe is 1.9 cases per 100,000 inhabitants per year. Untreated, autoimmune hepatitis soon develops into liver cirrhosis. However, if immunosuppressive low-dosed therapy is commenced in good time and continued consistently throughout the patient's life, the patient has a normal life expectancy. For differential diagnosis, a current infection with hepatitis viruses must be ruled out by investigation of the appropriate serological parameters.

Circulating autoantibodies have gained great significance for the diagnosis of AIH. They occur in the majority of patients, but their role in the pathogenesis is questionable. Also, there is no clear correlation between the activity or the prognosis of the disease and the antibody titre. Besides antibodies against SLA/LP the following autoantibodies are associated with AIH: antibodies against cell nuclei (**ANA**), **nDNA**, smooth muscle (**SMA**, with the most important target antigen being **F-actin**), liver-kidney micro-somes (**LKM-1**, target antigen: cytochrome P450 IID6) and granulocytes (**pANCA**). The **autoantibodies against SLA/LP** that can today be measured by EUROIMMUN enzyme immunoassays probably have the highest diagnostic accuracy of all antibodies involved in autoimmune hepatitis. Anti-SLA/LP occur in AIH either singly or together with other autoantibodies. Although their prevalence is only between 10% and 30%, the predictive value is almost 100%: Every positive result essentially provides evidence of autoimmune hepatitis (provided that the corresponding clinical symptoms are present).

Investigations of antibodies against SLA/LP have up until now only been possible in a few special laboratories. It was earlier suspected that the SLA antigen is identical with the liver cytokeratins 8 and 18 or with the enzyme glutathione-S-transferase. It has only recently been found that this is apparently not the case. The newly discovered target antigen SLA/LP has presently only been identified reliably at DNA level. It is probably a cytoplasmic molecule which is involved in the regulation of the protein biosynthesis (an UGA-suppressor-tRNA associated protein). The great diagnostic value of autoantibodies against SLA/LP has now been clearly proved.

Autoantibodies against cell nuclei (ANA) and against smooth muscle (SMA) are frequent in AIH, but also occur in 10% to 20% of patients with chronic viral hepatitis and in other diseases. Autoantibodies against LKM-1 can only be demonstrated in about 1% of adult AIH patients, but they are more frequent in children. Antibodies against LKM-1 can also be found 1% to 2% of patients with a positive hepatitis C serology. As opposed to all other autoantibodies, antibodies against SLA/LP are highly specific for AIH and have not been described in viral hepatitis.

Some authors classify autoimmune hepatitis in accordance with the autoantibody status: Subtype I (ANA, SMA), Subtype II (antibodies against LKM-1) and Subtype III (antibodies against SLA/LP). This classification has probably neither a clinical nor a therapeutic and prognostic significance.

The diagnosis of AIH is based on the clinical picture, the biochemical and serological test results and the histological evidence of an inflammation reaction. It requires the exclusion of other causes of chronic hepatitis, such as viruses, alcohol and drugs, and a delimitation to the other autoimmune liver diseases (primary biliary liver cirrhosis and primary sclerosing cholangitis). The diagnosis AIH is appropriate when 4 of the following 5 main criteria are fulfilled. If all 5 criteria are present, the diagnosis is definite. The diagnosis is confirmed if the patient responds to immunosuppressive therapy.

Main criteria for the diagnosis of autoimmune hepatitis:

1. Histological evidence of hepatitis
2. Detection of autoantibodies (ANA, SMA, LKM, SLA/LP)
3. Hypergammaglobulinaemia
4. Negative HBV and HCV serology
5. Detection of the HLA antigens B8, DR3 or DR4

The serological determination of **autoantibodies against SLA/LP** provides an exact delimitation to viral hepatitis in many patients with AIH. The autoimmune serology of hepatitis has thus been enriched by a further parameter whose significance must be rated more highly than that of antibodies against cell nuclei, DNA, smooth muscle and liver-kidney microsomes (**ANA, nDNA, ASMA, LKM**). Testing for these antibodies has considerable consequences for the hepatology clinic: the consequence of incorrect treatment of AIH with interferon can be just as fatal as the immunosuppressive therapy of the virus infection.

Primary biliary cirrhosis (PBC) is a liver cirrhosis which takes the form of a non-suppurative, destructive inflammation of the biliferous ducts. It occurs in woman ten times more frequently than in men. In the foreground of the clinical picture is cholestasis. In Europe, 13 out of every 100,000 persons develop this disease each year. As a cirrhosis does not occur in every case, and mostly not until the late stage, this disease is better described as “**chronic non-suppurative destructive cholangitis**”. Pathognomonic is the serological detection of autoantibodies against the mitochondrial fraction M2 (**AMA M2**) and against **Nuclear Dots** (acid protein SP100). Approximately 10% to 20% of PBC patients develop a secondary autoimmune hepatitis (also called **overlap syndrome**). In these cases, as in AIH, autoantibodies can frequently be detected. Antibodies against SLA/LP indicate a secondary AIH (overlap syndrome) which requires an immunosuppressive therapy.

The incidence of **primary sclerosing cholangitis (PSC)** is quoted as being 4 cases in every 100,000 inhabitants per year. The clinical picture is also characterised by cholestasis. Diagnostic pointers are the relevant laboratory findings, the histology and ERCP (endoscopic retrograde cholangiopancreatography). Men are mainly affected, and in half of the patients a concomitant Colitis ulcerosa is present (inversely, the prevalence of PSC in cases of Colitis ulcerosa is 4%). Most patients with PSC display autoantibodies against granulocytes (**pANCA**). Occasionally, pANCA can also occur in AIH and PBC, and their value for the differential diagnosis is limited in this respect - however, pANCA are an indication

of an autoimmune disease of the liver and can possibly be used as a means of delimitation to infectious forms of hepatitis. Even though overlaps between PSC and AIH have been reported, antibodies against SLA/LP have not been found in such patients.

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ANA Screen

for routine analysis

Qualitative determination of IgG class antibodies against nuclear antigens in human serum or plasma

IVD

LOT

See external label

2°C 8°C

$\Sigma = 96$ tests

REF DKO099

1. INTENDED PURPOSE

For *In Vitro* Diagnostic Use

For Laboratory Professional Use

The ANA Screen assay is a manual *in vitro* diagnostic device intended for the qualitative determination of IgG antibodies directed against Sm (Smith), RNP/Sm, Scl-70, SS-A (Ro) (52kDa e 60kDa), SS-B (La), Jo1, U1-SmRNP, CENP-B, dsDNA and Histones in human serum or plasma. Results are to be used in conjunction with other clinical and laboratory data as an aid in the diagnosis of systemic rheumatic autoimmune diseases.

2. CLINICAL SIGNIFICANCE

Anti-nuclear autoantibodies (ANAs) are a heterogeneous group of autoantibodies that recognise nuclear macromolecules and their complexes; antibodies to certain cytoplasmic proteins might also fall into this category¹⁻⁶. These antibodies commonly occur in the sera of patients with systemic autoimmune rheumatic diseases (SARDs) and can bind to DNA, RNA and proteins, as well as complexes of nucleic acids with proteins that are essential for different intracellular functions such as replication and transcription.

ANAs' targets can be divided into two groups in the context of systemic rheumatic autoimmune diseases: antibodies that recognise DNA, histones and nucleosomes^{1,2,5,6}; and those that bind to complexes of RNA with RNA-binding proteins (RBPs) and small nuclear ribonucleoproteins (snRNPs)¹.

Autoantibodies directed to nuclear antigens are serological hallmarks of SARDs and their detection is used by clinicians to support the diagnosis of disorders such as SLE, Sjögren syndrome, systemic sclerosis, polymyositis and mixed connective tissue disease. It is well documented that these are complex disorders and characterised by the presence of a wide spectrum of autoantibodies. Furthermore, it has widely been shown that ANAs are also frequently found in the sera of patients with non-rheumatic diseases and healthy subjects that do not develop a rheumatic disease.

Due to the complexity of SARDs, ANA tests should always be considered together with the evaluation of the patient's history, clinical manifestations and other laboratory tests to diagnose these conditions.

3. PRINCIPLE OF THE METHOD

The ANA Screen test is an indirect enzyme immunometric assay (ELISA) based on the binding of present antibodies with Sm, RNP/Sm, SS-A (Ro), SS-B (La), Scl-70, Jo1, U1-SmRNP, CENP-B, dsDNA and Histones antigens coated on the microplates. Any antibodies present in calibrators, controls or prediluted patient samples bind to the inner surface of the wells. After a 30 minutes incubation the microplate is washed with wash buffer for removing non-reactive serum components.

An anti-human-IgG horseradish peroxidase conjugate solution recognises IgG class antibodies bound to the immobilised antigens. After a 30 minutes incubation any excess enzyme conjugate, which is not specifically bound is washed away with wash buffer.

A chromogenic substrate solution containing TMB is dispensed into the wells. After 15 minutes of incubation the colour development is stopped by adding the stop solution. The solution turns yellow at this point. The level of colour is directly proportional to the concentration of IgG antibodies present in the original sample.

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- Calibrator (1 vial, 1.2 mL)
Phosphate buffer 0.1M, NaN₃ < 0.1%, human serum
CAL1 **REF DCE002/09906-0**
- Controls (2 vials, 1.2 mL each, ready to use)
Phosphate buffer 0.1M, NaN₃ < 0.1%, human serum
Negative Control **REF DCE045/09901-0**
Positive Control **REF DCE045/09902-0**
- Sample Diluent (1 vial, 100 mL)
Phosphate buffer 0.1M, NaN₃ < 0.1%
REF DCE053-0
- Conjugate (1 vial, 15 mL)
Anti h-IgG conjugated with horseradish peroxidase (HRP),
BSA 0.1%, ProClin™ < 0.0015%
REF DCE002/09902-0
- Coated Microplate (1 breakable microplate)
Microplate coated with ANA antigens
REF DCE002/09903-0
- TMB Substrate (1 vial, 15 mL)
H₂O₂-TMB (0,26 g/L) (avoid any skin contact)
REF DCE004-0
- Stop Solution (1 vial, 15 mL)
Sulphuric acid 0.15M (avoid any skin contact)
REF DCE005-0
- 10X Conc. Wash Solution (1 vial, 50 mL)
Phosphate buffer 0.2M pH 7.4 **REF DCE054-0**

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Precision Pipetting Devices Microplate reader (450 nm, 620-630 nm)

5. WARNINGS

- This kit is intended for *in vitro* use by professional persons only. Not for internal or external use in Humans or Animals.
- Use appropriate personal protective equipment while working with the reagents provided.
- Follow Good Laboratory Practice (GLP) for handling blood products.
- Material of animal origin used in the preparation of the kit has been obtained from animals certified as healthy and the bovine protein has been obtained from countries not infected by BSE, but these materials should be handled as potentially infectious.
- Some reagents contain small amounts of Sodium Azide (NaN_3) or ProClin™ 300 as preservative. Avoid the contact with skin or mucosa.
- All human source material used in the preparation of the reagents has been tested and found negative for antibody to HIV 1&2, HbsAg, and HCV. No test method however can offer complete assurance that HIV, HBV, HCV or other infectious agents are absent. Therefore, the Calibrators and the Controls should be handled in the same manner as potentially infectious material.
- Sodium Azide may be toxic if ingested or absorbed through the skin or eyes; moreover, it may react with lead or copper plumbing to form potentially explosive metal azides. If you use a sink to remove the reagents, wash through large amounts of water to prevent azide build-up.
- The TMB Substrate contains an irritant, which may be harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous and corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/ H_2O_2 to directed sunlight, metals or oxidants. Do not freeze the solution.

6. PRECAUTIONS

- Please adhere strictly to the sequence of pipetting steps provided in this protocol. The performance data represented here were obtained using specific reagents listed in this Instruction For Use.
- All reagents should be stored refrigerated at 2-8°C in their original container. Any exceptions are clearly indicated.
- Allow all kit components and specimens to reach room temperature (22-28°C) and mix well prior to use.
- Do not interchange kit components from different lots. The expiry date printed on box and vials labels must be observed. Do not use any kit component beyond their expiry date.
- If you use automated equipment, the user has the responsibility to make sure that the kit has been appropriately tested.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background. To improve the performance of the kit

on automatic systems is recommended to increase the number of washes.

- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate.
- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.
- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.
- Maximum precision is required for reconstitution and dispensation of the reagents.
- Samples microbiologically contaminated, highly lipemic or haemolysed should not be used in the assay.
- Plate readers measure vertically. Do not touch the bottom of the wells.
- **WARNING: the conjugate reagent is designed to ensure maximum dose sensitivity and may be contaminated by external agents if not used properly;** therefore, it is recommended to use disposable consumables (tips, bottles, trays, etc.). For divided doses, take the exact amount of conjugate needed and do not re-introduce any waste product into the original bottle. In addition, **for doses dispensed with the aid of automatic and semi-automatic devices,** before using the conjugate, it is advisable to clean the fluid handling system, ensuring that the procedures of washing, deproteinization and decontamination are effective in avoiding contamination of the conjugate. This **procedure is highly recommended when the kit is processed using analysers which are not equipped with disposable tips.**
For this purpose, DiaMetra supplies a separate decontamination reagent for cleaning needles.

7. REAGENT STORAGE AND STABILITY

Store the kit at 2÷8°C in the dark.

- The kit is stable at 2-8°C until the expiry date stated in the external kit label.
- Once opened, the kit is stable at 2-8°C for 6 months.
- The diluted wash solution is stable for 30 days at 2-8°C.

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8. SAMPLE COLLECTION AND STORAGE

The assay should be performed using serum (standard sampling tubes or tubes containing serum separating gel) or plasma (lithium heparin, sodium heparin or potassium EDTA) samples.

Sample Storage	Duration
2-8°C	96 hours
Freeze/thaw cycles	3 cycles

9. PROCEDURE

9.1. Preparation of the Calibrators (C₁)

Calibrators are ready to use; the concentration of the Calibrator is printed on the label.

9.2. Preparation of the Wash Solution

Dilute the contents of each vial of the buffered wash solution concentrate (10X) with distilled water to a final volume of 500 mL prior to use. For smaller volumes respect the 1:10 dilution ratio. The diluted wash solution is stable for 30 days at 2-8°C.

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9.3. Preparation of Samples

For determination of ANA antibodies, human serum or plasma are the preferred sample matrixes.

All serum and plasma samples must be prediluted with sample diluent 1:100; for example, 10 µL of sample may be diluted with 990 µL of sample diluent.

Fasting samples are not necessary and no special preparations are required. Collect blood by venepuncture into vacutainers and separate serum (after clot formation) or plasma from the cells by centrifugation.

Neither Bilirubin nor Haemolysis have significant effect on the procedure.

The Controls are ready to use.

9.4. Procedure

- **Allow all reagents to reach room temperature (22-28°C) for at least 30 minutes.** At the end of the assay store immediately the reagents at 2-8°C: avoid long exposure to room temperature.
- Unused coated microwell strips should be released securely in the foil pouch containing desiccant and stored at 2-8°C.
- To avoid potential microbial and/or chemical contamination, unused reagents should never be transferred into the original vials.
- As it is necessary to perform the determination in duplicate in order to improve accuracy of the test results, prepare two wells for each point of the calibration curve (C₁), two for each Control, two for each sample, one for Blank.

Reagent	Calibrator	Sample/ Controls	Blank
Calibrator C ₁	100 µL		
Controls		100 µL	
Diluted Sample		100 µL	

Incubate for 30 minutes at room temperature (22-28°C). Remove the content from each well, wash the wells 3 times with 300 µL of diluted wash solution.

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Incubate for 30 minutes at room temperature (22-28°C). Remove the content from each well, wash the wells 3 times with 300 µL of diluted wash solution.

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Shake the microplate gently.

Read the absorbance (E) at 450 nm against a reference wavelength of 620-630 nm or against Blank within 5 minutes.

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DiaMetra recommends the users to maintain graphic records of the control values generated with each assay run, including the running means, SDs and %CVs. This information will facilitate the controls trending analysis relating to the performance of current and historical control lots relative to the supplied Quality Control data. The trending will assist in the identification of assays which give control values significantly different from their average range.

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11. CALCULATION OF RESULTS

Determine the mean absorbance for each duplicate sample.

To obtain sample value: divide the mean absorbance of the sample by the absorbance of the Calibrator, then multiply by the concentration of the Calibrator, as shown below:

$$\text{Sample Conc.} = \frac{\text{OD sample}}{\text{OD C1}} \times \text{Conc. C1}$$

Reactivity is not connected in a linear way to the amount of antibodies present. Although an increase or a decrease in the concentration of antibodies results in increased or decreased responsiveness, the change is not in proportion (e.g., doubling the concentration of antibodies does not lead to a doubling of responsiveness).

To achieve greater accuracy in the determination of antibodies it is recommended to test serial dilutions of the sample. The final dilution that is positive in the test is the antibody concentration of the patient.

12. METROLOGY AND TRACEABILITY

The ANA Screen assay is traceable to the Antinuclear Antibodies (ANA) Reference Standards by Autoantibody Standardisation Committee (ASC) - Center for Disease Control (CDC) and WHO anti-dsDNA 15/174.

13. INTERPRETATION OF RESULTS

ANA Screen (AU/ mL)	Interpretation
< 25	The sample should be considered negative
≥ 25	The sample should be considered positive

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It is recommended that each laboratory establishes its own normal and pathological ranges of serum anti-ANA.

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14.1. Diagnostic Sensitivity and Specificity

The sensitivity and specificity were determined with guidance from CLSI EP-24 "Assessment of the Diagnostic Accuracy of Laboratory Tests Using Receiver Operating Characteristic Curves" using 50 negative and 53 positive samples run on a single reagent lot.

		DKO099		Total
		Positive	Negative	
True state	Positive	53	0	53
	Negative	0	50	50
Total		53	50	103

Diagnostic sensitivity: 100%

Diagnostic specificity: 100%

14.2. Precision

Precision of the ANA Screen kit was determined by performing a complex precision study.

Repeatability: A total of 3 serum samples were assayed using 2 lots of reagents in 10 replicates per sample, by 3 operators, once a day for 5 days. Data from one representative lot is shown below:

Repeatability results from one representative lot:

Sample	n	+ve results	-ve results	% +ve results	% -ve results
Cut off	150	82	68	54.7	45.3
Negative	150	1	149	0.7	99.3
Positive	150	150	0	100.0	0.0

Reproducibility results for the combined data from two lots is shown below:

Sample	n	+ve results	-ve results	% +ve results	% -ve results
Cut off	300	163	137	54.3	45.7
Negative	300	3	297	1.0	99.0
Positive	300	300	0	100.0	0.0

14.3. Serum-plasma study

The ANA Screen matrix comparison study was performed to evaluate the difference across tube types (serum separator tubes (SST), lithium heparin plasma, sodium heparin plasma and K2 EDTA plasma) versus the control samples (red top serum, without additive) following CLSI GP34-A "Validation and Verification of Tubes for Venous and Capillary Blood. Specimen Collection". A total of 20 samples (16 native, 4 contrived) were evaluated.

Sample Type	No. Samples	Results vs Serum	
		False pos	False neg
SST	20	0	0
K2 EDTA	20	0	0
Lithium heparin	20	0	0
Sodium heparin	20	0	0

14.4. Interferents

The following substances do not interfere in the ANA Screen assay when the concentrations presented in the following table are below the stated threshold.

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15. LIMITATIONS OF USE

- As in the case of any diagnostic procedure, results must be interpreted in conjunction with the patient's clinical presentation and other information available to the physician.
- The performance characteristics of this assay have not been established in a paediatric population.
- Heterophilic antibodies in human serum can react with reagent immunoglobulins, interfering with *in vitro* immunoassays⁸. Patients routinely exposed to animals or to animal serum products can be prone to this interference and anomalous values may be observed.

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19. PRODUCT COMPLAINTS AND TECHNICAL SUPPORT

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BIBLIOGRAFIA

1. Condemi, Jhon J. The autoimmune disease. The Journal of the American Medical Association 1987, Vo 258 n 20 2920 - 2929
2. Reichlin, M., Van Venrooij, W. J.: Autoantibodies to UNRP particles: relationship to clinical diagnosis and nephritis. Clin. exp. Immunol. 1991; 83:286-90.
3. Tan EM, et al.: The 1982 Revised Criteria for the classification of Systemic Lupus Erythematosus. Arthritis and Rheumatism 25: 1271-1277, 1982.
4. McCarty GA. Autoantibodies and their relation to rheumatic diseases. Medical Clinics of North America 70: 237-261, 1986.
5. Hardin JA. The lupus autoantigens and the pathogenesis of systemic lupus erythematosus. Arth Rheum 29: 457-460, 1986.
6. The first international Calibrator for antibodies to double stranded DNA. Annals of the Rheumatic Disease 198; Vo 47: 740 - 746
7. Homburger HA, et al: detection of Antinuclear antibodies, comparative evaluation of enzyme immunoassay and indirect immunofluorescence methods. Arch Pathol Lab Med 122:993-999 (1998)
8. Smeenk, R. et al. Avidity of Antibodies to dsDNA. Comparison of IFT on Crithidia Lucilaie, FARR assay and PEG assay. The journal of Immunology 1982 Vo 128 n.1. 73 -78

Ed. 09/2018

DCM095-7

Dia.Metra S.r.l.

Via Pozzuolo 14,
06038 SPELLO (PG) Italy
Tel. +39-0742-24851
Fax +39-0742-316197
E-mail: info@diametra.com

DCM095-7
Ed. 09/2018

Anti dsDNA IgG

for routine analysis


Quantitative determination of IgG class antibodies against dsDNA in human serum or plasma

IVD

LOT

See external label

$\Sigma = 96$ tests

REF DKO095

INTENDED USE

Anti dsDNA IgG kit is an indirect solid phase enzyme immunometric assay (ELISA) designed for the quantitative measurement of IgG class antibodies directed against dsDNA in human serum or plasma. Anti dsDNA IgG is intended for laboratory use only.

1. CLINICAL SIGNIFICANCE

Anti dsDNA IgG test is used for initial diagnosis of Systemic Lupus Erythematosus (SLE) and for diagnosis of SLE different diseases. Besides the determination of high titers of antinuclear antibody (ANA), the determination of autoantibodies against dsDNA is one of the ACR criteria (American College of Rheumatology) for the diagnosis of Systemic Lupus Erythematosus (SLE). The determination of the concentration of antibodies can be used to monitor treatment success and predict possible attacks of the disease (SLE)

2. PRINCIPLE

Anti dsDNA IgG test is based on the binding of serum or plasma IgG antibodies on dsDNA coated on the microplates. In the first step the antibodies in calibrators, controls or prediluted patient samples bind into the inner surface of the wells. After an incubation the microplate is washed with a wash buffer for removing non-reactive serum components. In the second step an anti human IgG horseradish peroxidase conjugated solution recognizes the IgG class antibodies bound to the immobilized dsDNA antigens. After an incubation any excess of enzyme conjugate, which is not specifically bound, is washed away with the wash buffer.

Then a chromogenic substrate solution containing TMB is dispensed into the wells. After 15 minutes of incubation the color development is stopped by adding the stop solution. The solution color changes into yellow. The amount of color is directly proportional to the concentration of the Anti dsDNA IgG antibodies present in the original sample.

The concentration of the anti dsDNA IgG antibodies in the sample are calculated through a calibration curve.

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Microplates reader (450 nm, 620-630 nm)

Notes

Store all reagents between $2\pm 8^\circ\text{C}$ in the dark.

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5. PRECAUTIONS

- Please adhere strictly to the sequence of pipetting steps provided in this protocol. The performance data represented here were obtained using specific reagents listed in this Instruction For Use.
- All reagents should be stored refrigerated at 2-8°C in their original container. Any exceptions are clearly indicated. The reagents are stable until the expiry date when stored and handled as indicated.
- Allow all kit components and specimens to reach room temperature (22-28°C) and mix well prior to use.
- Do not interchange kit components from different lots. The expiry date printed on box and vials labels must be observed. Do not use any kit component beyond their expiry date.
- **WARNING: the conjugate reagent is designed to ensure maximum dose sensitivity and may be contaminated by external agents if not used properly;** therefore, it is recommended to use disposable consumables (tips, bottles, trays, etc.).

For divided doses, take the exact amount of conjugate needed and do not re-introduce any waste product into the original bottle. In addition, **for doses dispensed with the aid of automatic and semi-automatic devices,** before using the conjugate, it is advisable to clean the fluid handling system, ensuring that the procedures of washing, deproteinization and decontamination are effective in avoiding contamination of the conjugate; **this procedure is highly recommended when the kit is processed using analyzers which are not equipped with disposable tips.**

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6. PROCEDURE

6.1. Preparation of the Calibrators ($\text{C}_0 \dots \text{C}_4$)

Calibration curve is ready to use and is calibrated against International Standard WHO Wo/80. The Calibrators have approximately the following concentration:

	C_0	C_1	C_2	C_3	C_4
IU/mL	0	12.5	25	50	200

NB: Calibrators concentration is lot-dependent; exact concentrations are indicated on Certificate of Analysis and labels.

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6.2. Preparation of the Sample

For determination of anti dsDNA antibodies, human serum or plasma are the preferred sample matrixes.

All serum and plasma samples have to be prediluted with sample diluent 1:100; for example 10 µL of sample may be diluted with 990 µL of sample diluent.

The patients need not to be fasting, and no special preparations are necessary. Collect blood by venipuncture into vacutainers and separate serum (after clot formation) or plasma from the cells by centrifugation.

Samples may be stored refrigerated at 2-8°C for at least 5 days. For longer storage of up to six months samples should be stored frozen at -20°C. To avoid repeated thawing and freezing the samples should be aliquoted.

Neither Bilirubin nor Hemolysis have significant effect on the procedure.

The Controls are ready to use.

6.3. Preparation of the Wash Solution

Dilute the contents of each vial of the buffered wash solution concentrate (10x) with distilled water to a final volume of 500 mL prior to use. For smaller volumes respect the 1:10 dilution ratio. The diluted wash solution is stable for 30 days at 2-8°C.

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6.4. Procedure

- **Allow all reagents to reach room temperature (22-28°C) for at least 30 minutes.** At the end of the assay, store immediately the reagents at 2-8°C: avoid long exposure to room temperature.
- Unused coated microwell strips should be released securely in the foil pouch containing desiccant and stored at 2-8°C.
- To avoid potential microbial and/or chemical contamination, unused reagents should never be transferred into the original vials.
- As it is necessary to perform the determination in duplicate in order to improve accuracy of the test results, prepare two wells for each point of the calibration curve (C₀-C₄), two for each Control, two for each sample, one for Blank.

Reagent	Calibrator	Sample/ Controls	Blank
Calibrator C ₀ -C ₄	100 µL		
Controls		100 µL	
Diluted Sample		100 µL	
Incubate 45 minutes at 37°C. Remove the content from each well; wash the wells 3 times with 300 µL of diluted wash solution. Important note: during each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel. Automatic washer: if you use automated equipment, wash the wells at least 5 times.			
Conjugate	100 µL	100 µL	
Incubate 45 minutes at 37°C. Remove the content from each well, wash the wells 3 times with 300 µL of diluted wash solution (if you use automated equipment, wash the wells at least 5 times). Washing: follow the same indications of the previous point.			
TMB Substrate	100 µL	100 µL	100 µL
Incubate 15 minutes in the dark at room temperature (22-28°C).			
Stop Solution	100 µL	100 µL	100 µL
Shake the microplate gently. Read the absorbance (E) at 450 nm against a reference wavelength of 620-630 nm or against Blank within 5 minutes.			

7. RESULTS

7.1. Calibration curve

For the test Anti dsDNA IgG a 4-Parameter-Fit with lin-log coordinates for optical density and concentration is the data reduction method of choice. Smoothed-Spline Approximation and log-log coordinates are also suitable. However we recommend using a Lin-Log curve.

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8. REFERENCE VALUES

Following ranges have been established with the Anti dsDNA IgG tests:

	Anti dsDNA IgG (IU/mL)
Normal	< 25
Elevated	> 25

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Positive results should be verified concerning the entire clinical status of the patient. Also every decision for therapy should be taken individually.

It is recommended that each laboratory establishes its own normal and pathological ranges of serum Anti dsDNA.

9. PERFORMANCE AND CHARACTERISTICS

9.1. Specificity

Comparison test against a commercial reference kit, performed on 100 sera (50 of them positive sera and 50 negative sera) showed a specificity > 99%.

9.2. Sensitivity

Comparison test against a commercial reference kit, performed on 100 sera (50 of them positive sera and 50 negative sera) showed a sensibility > 99%..

9.3. Detection limit

The lowest concentration of Anti dsDNA IgG that can be distinguished from zero Calibrator is 0.135 IU/mL.

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BIBLIOGRAPHY

1. Condemni, Jhon J. The autoimmune disease. The Journal of the American Medical Association 1987, Vo 258 n 20 2920 - 2929
2. Reichlin, M., Van Venrooij, W. J.: Autoantibodies to UNRP particles: relationship to clinical diagnosis and nephritis. Clin. exp. Immunol. 1991; 83:286-90.
3. Tan EM, et al.: The 1982 Revised Criteria for the classification of Systemic Lupus Erythematosus. Arthritis and Rheumatism 25: 1271-1277, 1982.
4. McCarty GA. Autoantibodies and their relation to rheumatic diseases. Medical Clinics of North America 70: 237-261, 1986.
5. Hardin JA. The lupus autoantigens and the pathogenesis of systemic lupus erythematosus. Arth Rheum 29: 457-460, 1986.

6. The first international Calibrator for antibodies to double stranded DNA. Annals of the Rheumatic Disease 198; Vo 47: 740 - 746
7. Homburger HA, et al: detection of Antinuclear antibodies, comparative evaluation of enzyme immunoassay and indirect immunofluorescence methods. Arch Pathol Lab Med 122:993-999 (1998)
8. Smeenk, R. et al. Avidity of Antibodies to dsDNA. Comparison of IFT on Crithidia Lucilaie, FARR assay and PEG assay. The journal of Immunology 1982 Vo 128 n.1. 73 -78

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1. Condem, Jhon J. The autoimmune disease. The Journal of the American Medical Association 1987, Vo 258 n 20 2920 - 2929
2. Reichlin, M., Van Venrooij, W. J.: Autoantibodies to UNRP particles: relationship to clinical diagnosis and nephritis. Clin. exp. Immunol. 1991; 83:286-90.
3. Tan EM, et al.: The 1982 Revised Criteria for the classification of Systemic Lupus Erythematosus. Arthritis and Rheumatism 25: 1271-1277, 1982.
4. McCarty GA. Autoantibodies and their relation to rheumatic diseases. Medical Clinics of North America 70: 237-261, 1986.
5. Hardin JA. The lupus autoantigens and the pathogenesis of systemic lupus erythematosus. Arth Rheum 29: 457-460, 1986.
6. The first international standard for antibody to double stranded DNA. Annals of the Rheumatic Disease 198; Vo 47: 740 - 746
7. Homburger HA, et al: detection of Antinuclear antibodies, comparative evaluation of enzyme immunoassay and indirect immunofluorescence methods. Arch Pathol Lab Med 122:993-999 (1998)
8. Smeenk, R. et al. Avidity of Antibodies to dsDNA. Comparison of IFT on Crithidia Lucilaie, FARR assay and PEG assay. The journal of Immunology 1982 Vo 128 n.1. 73 -78

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5. Jennette JC, Falk RJ, Bacon PA, Basu N, *et al.* 2012 revised International Chapel Hill Consensus Conference Nomenclature of Vasculitides. *Arthritis Rheum.* 2013 Jan; 65(1):1-11
6. Alberici F, Martorana D, Vaglio A. Genetic aspects of anti-neutrophil cytoplasmic antibody-associated vasculitis. *Nephrol Dial Transplant.* 2015 Apr;30 Suppl 1:i37-45.
7. Basic QC Practices On-line Course; <http://www.Westgard.com>.
8. Boscatto, LM. and Stuart, MC., 'Heterophilic antibodies: a problem for all immunoassays'. *Clin Chem*, 34, 1988, pp 27-33

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Anti PR-3 (c-ANCA)

for routine analysis

Semi quantitative determination of IgG class antibodies against proteinase 3 (PR-3) in human serum or plasma

IVD

LOT

See external label

2°C 8°C

$\Sigma = 96$ tests

REF DKO091

1. INTENDED PURPOSE

For *In Vitro* Diagnostic Use

For Laboratory Professional Use

Anti PR-3 (c-ANCA) assay is a manual *in vitro* diagnostic device intended for the quantitative determination of IgG class autoantibodies against proteinase 3 (PR-3) in human serum or plasma. Results are to be used in conjunction with other clinical and laboratory data as an aid in the diagnosis of certain autoimmune vasculitides such as Wegener's granulomatosis.

2. CLINICAL SIGNIFICANCE

Proteinase 3 (PR-3 or PR3) is a small serine protease located in the primary granules and secretory vesicles in the cytoplasm of neutrophils. PR3 is stored in the cytoplasmic granules in resting neutrophils, whilst the expression of the protein on cell surface increases during the activation of neutrophils and during apoptosis¹. PR3 is also found in neutrophil extracellular traps (NETs), that consist of chromatin released by neutrophils as defence against pathogens.

Autoantibodies against the PR3 are frequently detected in a specific class of autoimmune vasculitides known as Anti-neutrophil cytoplasmatic autoantibodies (ANCA) associated vasculitis (AAV). These vasculitides are characterised by the presence of autoantibodies against cytoplasmatic components of neutrophils that may have, cytoplasmatic staining (c-ANCA) or a perinuclear staining (p-ANCA)¹⁻³.

The mechanism of autoimmunity for ANCA antigens is not still being evaluated. Activated neutrophils can undergo apoptosis creating NETs, which are involved in the induction of inflammation. Persistence of NETs during the infection can lead to an autoimmune response^{3,4}. Increased expression of MPO and PR3 in neutrophils and monocytes can influence the pathogenic autoimmune response leading to ANCA-neutrophils activation³.

ANCA-associated vasculitides is classified as small vessel vasculitis (SVV), which predominantly affect the small vessels, intraparenchymal arteries, arterioles, capillaries and venules⁵. ANCA-associated vasculitides are necrotising vasculitis characterised by a lack of immune deposits in

vessel walls causing conditions such as Microscopic polyangiitis (MPA) and Granulomatosis with polyangiitis (GPA, formerly known as Wegner's syndrome).

The frequency of serotype (MPO or PR3), clinical manifestation and geographic distribution of ANCA-associated vasculitides varies, suggesting a role of the environment and the genetic background of patients affected by this category of pathologies^{1,3,6}.

3. PRINCIPLE OF THE METHOD

The Anti PR-3 (c-ANCA) assay is an indirect enzyme immunometric assay (ELISA) where anti- PR3 (antibody) in the sample binds to human neutrophil proteinase 3 coated into the microplate (solid phase). In the first step the antibodies in calibrators, controls or prediluted patient samples bind to the inner surface of the wells. After a 60-minute incubation, the microplate is washed to remove non-reactive serum components. In the second step an anti-human-IgG horseradish peroxidase (HRP) conjugated solution binds to the IgG class antibodies bound to the immobilised antigens. After a further incubation step, any excess enzyme conjugate which is not specifically bound is washed away.

A chromogenic substrate solution containing TMB is then dispensed into the wells. After a further incubation step, the enzyme HRP in the bound fraction reacts with the Substrate (H_2O_2) and the TMB Substrate and develops a blue colour that changes into yellow when the Stop Solution (H_2SO_4) is added. The colour intensity is directly proportional to the concentration of IgG antibodies in the sample.

The concentration of anti- PR3 IgG antibodies in the sample is calculated through a calibration curve.

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5. WARNINGS

- This kit is intended for *in vitro* use by professional persons only. Not for internal or external use in Humans or Animals.
- Use appropriate personal protective equipment while working with the reagents provided.
- Follow Good Laboratory Practice (GLP) for handling blood products.
- All human source material used in the preparation of the reagents has been tested and found negative for antibody to HIV 1&2, HbsAg, and HCV. No test method however can offer complete assurance that HIV, HBV, HCV or other infectious agents are absent. Therefore, the Calibrators and the Controls should be handled in the same manner as potentially infectious material.
- Material of animal origin used in the preparation of the kit has been obtained from animals certified as healthy and the bovine protein has been obtained from countries not infected by BSE, but these materials should be handled as potentially infectious.
- Some reagents contain small amounts of Sodium Azide (NaN_3) or ProClin™ 300 as preservative. Avoid contact with skin or mucosa.
- ProClin™ 300: Classification under CLP
Skin Sens. 1

Hazard statements

H317: May cause an allergic skin reaction

Precautionary statements

P261: Avoid breathing dust / fumes / gas / mist / vapours / spray.

P272: Contaminated work clothing should not be allowed out of the workplace.

P280: Wear protective gloves / protective clothing / eye protection / face protection.

P302/352 IF ON SKIN: Wash with plenty of water / soap and water.

P321: Specific treatment (see instructions on the label)

P333/313: If skin irritation or rash occurs: get medical attention.

- Sodium Azide may be toxic if ingested or absorbed through the skin or eyes; moreover, it may react with lead or copper plumbing to form potentially explosive metal azides. If you use a sink to remove the reagents, wash through large with amounts of water to prevent azide build-up.
- The TMB Substrate contains an irritant, which harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous, corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/ H_2O_2 to direct sunlight, metals or oxidants. Do not freeze the solution.

6. PRECAUTIONS

- Please adhere strictly to the sequence of pipetting steps provided in this protocol. The performance data represented here were obtained using specific reagents listed in this Instruction For Use.
- All reagents should be stored refrigerated at 2-8°C in their original container. Any exceptions are clearly indicated.
- Allow all kit components and specimens to reach room temperature (22-28°C) and mix well prior to use.
- Do not interchange kit components from different lots. The expiry date printed on box and vials labels must be observed. Do not use any kit component beyond their expiry date.
- If you use automated equipment, the user has the responsibility to make sure that the kit has been appropriately tested.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background. To improve the performance of the kit on automatic systems is recommended to increase the number of washes.
- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate
- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.
- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.
- Maximum precision is required for reconstitution and dispensation of the reagents.

- Samples microbiologically contaminated, highly lipemic or haemolysed should not be used in the assay.
- Plate readers measure vertically. Do not touch the bottom of the wells.
- **WARNING: the conjugate reagent is designed to ensure maximum dose sensitivity and may be contaminated by external agents if not used properly;** therefore, it is recommended to use disposable consumables (tips, bottles, trays, etc.). For divided doses, take the exact amount of conjugate needed and do not re-introduce any waste product into the original bottle. In addition, **for doses dispensed with the aid of automatic and semi-automatic devices,** before using the conjugate, it is advisable to clean the fluid handling system, ensuring that the procedures of washing, deproteinisation and decontamination are effective in avoiding contamination of the conjugate; **this procedure is highly recommended when the kit is processed using analysers which are not equipped with disposable tips.**
For this purpose, Diametra supplies a separate decontamination reagent for cleaning needles.

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Since no international reference preparation for Anti-proteinase 3 antibodies is available, the assay system is calibrated in relative arbitrary units. The Calibrators are ready to use and have the following concentration:

	C ₀	C ₁	C ₂	C ₃	C ₄
AU/mL	0	10	20	40	160

The controls are ready to use; the concentration of the control is printed on the label.

9.2. Preparation of the Wash Solution

Dilute the content of the vial "10X Conc. Wash Solution" with distilled water to a final volume of 500 mL prior to use. For smaller volumes respect the 1:10 dilution ratio.

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9.3. Preparation of Samples

Test samples should be clear. Contamination by lipemia is best avoided but does not interfere with this assay.

Testing of heat-inactivated serum or plasma samples is not recommended.

All serum or plasma samples must be prediluted 1:100 with sample diluent; for example, 10 µL of sample should be diluted with 990 µL of sample diluent.

9.4. Procedure

- **Allow all reagents to reach room temperature (22-28°C) for at least 30 minutes.** At the end of the assay, immediately store the reagents at 2-8°C: avoiding long exposure to room temperature.
- Unused coated microwell strips should be released securely in the foil pouch containing desiccant and stored at 2-8°C.
- To avoid potential microbial and/or chemical contamination, unused reagents should never be transferred into the original vials.
- As it is necessary to perform the determination in duplicate in order to improve accuracy of the test results, prepare two wells for each point of the calibration curve (C₀-C₄), two for each Control, two for each sample, one for Blank.

Reagent	Calibrator	Sample/ Controls	Blank
Calibrator C ₀ -C ₄	100 µL		
Controls		100 µL	
Diluted Sample		100 µL	
Incubate 60 minutes at room temperature (22-28°C). Remove the content from each well, wash the wells 3 times with 300 µL of diluted wash solution. Important note: during each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel. Automatic washer: if you use automated equipment, wash the wells at least 5 times.			
Conjugate	100 µL	100 µL	
Incubate 30 minutes at room temperature (22-28°C). Remove the contents from each well, wash the wells 3 times with 300 µL of diluted wash solution. Washing: follow the same indications of the previous point.			
TMB Substrate	100 µL	100 µL	100 µL
Incubate 15 minutes in the dark at room temperature (22-28°C).			
Stop Solution	100 µL	100 µL	100 µL
Shake the microplate gently Read the absorbance (E) at 450 nm against a reference wavelength of 620-630 nm or against Blank within 5 minutes.			

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The mean concentration of each control level is documented in the QC report included with each kit. These mean concentration levels are determined over several assays which are run in duplicate in multiple locations across each plate.

DiaMetra recommends the users to maintain graphic records of the control values generated with each assay run, including the running means, SDs and %CVs. This information will facilitate the controls trending analysis relating to the performance of current and historical control

lots relative to the supplied Quality Control data. The trending will assist in the identification of assays which give control values significantly different from their average range.

When interpreting control data, users should note that this product was designed and developed as a manual product. The range stated on the QC certificate should be appropriate for assays that are performed manually and with strict adherence to the Assay Procedure described above. It is recognised by Quality Control professionals, that as a result of differences in conditions and practices, there will always be variability in the mean values and precision of control measurements between different laboratories⁷.

In order for the test results to be considered valid, all of the criteria listed below must be met. If any of these are not met, the test should be considered invalid and the assay repeated:

- The absorbance of the prediluted PR-3 IgG Positive Control must be greater than the absorbance of the prediluted Negative Control.
- The negative and positive controls are intended to monitor for substantial reagent failure and they will not ensure precision at the assay cut-off.
- This test is only valid if the optical density at 450 nm for Negative Control and Positive Control as well as for the Calibrator C₀-C₄ complies with the respective range indicated on the Quality Control Certificate enclosed to each test kit: if any of these criteria is not met, the results are invalid and the test should be repeated.

11. CALCULATION OF RESULTS

A variety of data reduction software packages are available, which may be employed to generate the mean calibration curve and to calculate the mean concentrations of unknown samples and controls. A 4-parameter logistic (4PL) curve fit, **including Calibrator 0 is required**. A smoothed spline and log-log fit including Calibrator 0 can be used. Other curve fitting algorithms are not recommended.

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12. MEASURING RANGE

The assay measuring range (AMR) is 3 – 160 AU/mL. Any value that reads below 3 AU/mL should be reported as "< 3 AU/mL". Any value that reads above 160 AU/mL should be reported as "> 160 AU/mL".

13. METROLOGY AND TRACEABILITY

The calibrators of this kit are traceable to the Centres for Disease Control and Prevention (CDC) Human Reference Serum #16 for Anti-Proteinase 3 Ab Catalogue #IS2721.

14. EXPECTED VALUES

The following ranges were determined using the Anti PR-3 (c-ANCA) assay and are provided for information only. The 90 % reference interval for apparently healthy adults were

calculated by a non-parametric method following guidance from CLSI C28-A3 "Defining, Establishing and Verifying Reference Intervals in the Clinical Laboratory".

No. of subjects	73
Mean (AU/mL)	8.71
SD (AU/mL)	3.21
Median (AU/mL)	8.12
Reference Interval (AU/mL)	4.19 – 14.40

The above ranges should be considered as guidelines only; it is recommended that each laboratory establish its own expected range based upon its own patient population.

15. INTERPRETATION OF RESULTS

Concentration	Interpretation
< 16 AU/mL	The sample should be considered negative
16 – 20 AU/mL	The sample should be graded equivocal and repeat testing / sampling should be performed according to internal practices
≥ 20 AU/mL	The sample should be considered positive

The above index and interpretation should be considered as guidelines only. Positive results should be verified concerning the entire clinical status of the patient. Also, every decision for therapy should be taken on an individual patient basis.

It is recommended that each laboratory establishes its own normal and pathological ranges of anti- PR3.

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Sensitivity	Concentration
Limit of Blank (LoB)	0.6 AU/mL
Limit of Detection (LoD)	1.6 AU/mL
Limit of Quantitation (LoQ)	3.2 AU/mL

16.2. Trueness

Trueness has been demonstrated through a recovery test of the Anti PR-3 (c-ANCA) assay with the human Anti-Proteinase 3 reference serum #16 from the Centres for Disease Control and Prevention (CDC).

16.3. Precision

Precision of the Anti PR-3 (c-ANCA) assay was determined by performing a complex precision study.

Repeatability: A total of 6 samples were assayed in 5 replicates, once a day for 5 days by 3 operators. Data from one representative lot is shown below:

Sample	n	Mean Conc. (AU/mL)	Within run (Repeatability)	
			SD	CV%
1	75	16.51	0.96	5.8%
2	75	23.94	0.93	3.9%
3	75	40.03	1.23	3.1%
4	75	59.68	3.40	5.7%
5	75	86.96	2.80	3.2%
6	75	136.70	6.24	4.6%

Reproducibility: A total of 6 samples were assayed in 5 replicates, once a day for 5 days by 3 operators. Results for the combined data from two lots is shown below:

Sample	n	Mean Conc. (AU/mL)	Within Laboratory (Reproducibility)	
			SD	CV%
1	150	16.57	1.49	9.0%
2	150	23.95	1.53	6.4%
3	150	40.59	2.67	6.6%
4	150	60.11	4.15	6.9%
5	150	86.69	6.12	7.1%
6	150	137.19	9.63	7.0%

16.4. Linearity

Linearity was evaluated based on CLSI EP-06, "Evaluation of the Linearity of Quantitative Measurement Procedures". For anti- PR3 concentration by Anti PR-3 (c-ANCA) assay, the measurement procedure shows linearity for the interval from 1.73 to 162.62 AU/mL within the allowable deviation of linearity (ADL) of ±15 %.

16.5. Diagnostic Sensitivity and Specificity

The sensitivity and specificity were determined with guidance from CLSI EP-24 "Assessment of the Diagnostic Accuracy of Laboratory Tests Using Receiver Operating Characteristic Curves" using 73 negative and 50 positive samples run on a single reagent lot.

		DKO091		Total
		Positive	Negative	
True state	Positive	50	0	50
	Negative	2	71	73
Total		52	71	123

Diagnostic sensitivity: 100%

Diagnostic specificity: 97.3%

16.6. Interferents

The following substances do not interfere with a bias of $\pm 15\%$ in the Anti PR-3 (c-ANCA) assay when the concentrations are below the stated threshold presented in the following table.

Potentially Interfering Reagent	Threshold Concentration
Bilirubin, conjugated	15 mg/dL
Bilirubin, unconjugated	15 mg/dL
Haemoglobin	200 mg/dL
Total Protein	10 g/dL
Triglycerides	500 mg/dL

16.7. Serum-plasma study

The Anti PR-3 (c-ANCA) assay matrix comparison study was performed to evaluate the difference across tube types (serum separator tubes (SST), lithium heparin plasma, sodium heparin plasma and K2 EDTA plasma) versus the control samples (red top serum, without additive) following CLSI EP9-A3 guidelines. A total of 20 samples (16 native, 4 spiked or diluted) to cover the range of 2.80 to 74.40 AU/mL were evaluated. Linear regression analysis was performed on the comparative data:

Sample type	Slope [95% CI]	Intercept (AU/mL) [95% CI]	Correlation coefficient (r)
SST	1.0 [0.97 to 1.08]	0.14 [-1.17 to 1.45]	0.99
Lithium Heparin	1.1 [1.06 to 1.21]	-0.46 [-2.09 to 1.17]	0.99
Sodium Heparin	1.1 [1.05 to 1.14]	-0.91 [-1.88 to 0.06]	1.00
EDTA	1.1 [1.0 to 1.12]	0.00 [-1.79 to 1.78]	0.99

17. LIMITATIONS OF USE

- As in the case of any diagnostic procedure, results must be interpreted in conjunction with the patient's clinical presentation and other information available to the physician.
- The performance characteristics of this assay have not been established in a paediatric population.
- Heterophilic antibodies in human serum can react with reagent immunoglobulins, interfering with *in vitro* immunoassays¹¹. Patients routinely exposed to animals or to animal serum products can be prone to this interference and anomalous values may be observed.

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19. BIBLIOGRAPHY

1. Cornec D, Cornec-Le Gall E, Fervenza FC, Specks U. ANCA-associated vasculitis - clinical utility of using ANCA specificity to classify patients. *Nat Rev Rheumatol*. 2016 Oct;12(10):570-9.
2. Rasmussen N, Wiik A, Jayne DR. A historical essay on detection of anti-neutrophil cytoplasmic antibodies. *Nephrol Dial Transplant*. 2015 Apr;30 Suppl 1:i8-13.
3. Jennette JC, Nachman PH. ANCA Glomerulonephritis and Vasculitis. *Clin J Am Soc Nephrol*. 2017 Oct 6;12(10):1680-1691.
4. Grau RG. Drug-Induced Vasculitis: New Insights and a Changing Lineup of Suspects. *Curr Rheumatol Rep*. 2015 Dec;17(12):71.
5. Jennette JC, Falk RJ, Bacon PA, Basu N, et al. 2012 revised International Chapel Hill Consensus Conference Nomenclature of Vasculitides. *Arthritis Rheum*. 2013 Jan; 65(1):1-11
6. Alberici F, Martorana D, Vaglio A. Genetic aspects of anti-neutrophil cytoplasmic antibody-associated vasculitis. *Nephrol Dial Transplant*. 2015 Apr;30 Suppl 1:i37-45.
7. Basic QC Practices On-line Course; <http://www.Westgard.com>.
8. Boscatto, LM. and Stuart, MC., 'Heterophilic antibodies: a problem for all immunoassays'. *Clin Chem*, 34, 1988, pp 27-33

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21. PRODUCT COMPLAINTS AND TECHNICAL SUPPORT

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
ORG 511 Anti-RNP/Sm

INTENDED PURPOSE

Anti-RNP/Sm is an ELISA test system for the quantitative measurement of IgG class autoantibodies against RNP/Sm in human serum or plasma. This product is intended for professional in vitro diagnostic use only.

Antibodies against the RNP/Sm complex are useful in the diagnosis of mixed connective tissue disorder (MCTD, Sharp syndrome) and related autoimmune diseases. Antibodies against the 70 kDa protein of this complex are a very specific marker for Sharp syndrome. The Sm proteins are recognised by antibodies that may occur in cases of mixed connective tissue disorder and systemic lupus erythematosus. Evaluation of a test result should always take into account all clinical and laboratory diagnostic findings.

SYMBOLS USED ON LABELS

In vitro diagnostic medical device

Manufacturer

Catalogue number

Sufficient for 96 determinations

Batch code

Use by

Temperature limitation

Keep away from sunlight

Do not reuse

Date of manufacture

CE marked according to 98/79/EC


Consult instructions for use

Electronic Instruction For Use: version

Microplate

Calibrator

Calibrator


Calibrator

Calibrator

Calibrator

Calibrator

Control positive

Control negative

Sample Buffer P

Enzyme Conjugate

TMB Substrate

Stop solution

Wash Buffer

Ready to use

PRINCIPLE OF THE TEST

Highly purified RNP/Sm is bound to microwells.

The determination is based on an indirect enzyme linked immune reaction with the following steps:

Specific antibodies in the patient sample bind to the antigen coated on the surface of the reaction wells. After incubation, a washing step removes unbound and unspecifically bound serum or plasma components. Subsequently added enzyme conjugate binds to the immobilized antibody-antigen-complexes. After incubation, a second washing step removes unbound enzyme conjugate. After addition of substrate solution the bound enzyme conjugate hydrolyses the substrate forming a blue coloured product. Addition of an acid stops the reaction generating a yellow end-product. The intensity of the yellow color correlates with the concentration of the antibody-antigen-complex and can be measured photometrically at 450 nm.

WARNINGS AND PRECAUTIONS

- All reagents of this kit are intended for professional in vitro diagnostic use only.
- Components containing human serum were tested and found negative for HBsAg, HCV, HIV1 and HIV2 by FDA approved methods. No test can guarantee the absence of HBsAg, HCV, HIV1 or HIV2, and so all human serum based reagents in this kit must be handled as though capable of transmitting infection.
- Bovine serum albumin (BSA) used in components has been tested for BSE and found negative.
- Avoid contact with the substrate TMB (3,3',5,5'-Tetramethyl-benzidine).
- Stop solution contains acid, classification is non-hazardous. Avoid contact with skin.
- Control, sample buffer and wash buffer contain sodium azide 0.09% as preservative. This concentration is classified as non-hazardous.
- Enzyme conjugate contains ProClin 300 0.05% as preservative. This concentration is classified as non-hazardous.

During handling of all reagents, controls and serum samples observe the existing regulations for laboratory safety regulations and good laboratory practice:

- First aid measures: In case of skin contact, immediately wash thoroughly with water and soap. Remove contaminated clothing and shoes and wash before reuse. If system fluid comes into contact with skin, wash thoroughly with water. After contact with the eyes carefully rinse the opened eye with running water for at least 10 minutes. Get medical attention if necessary.
 - Personal precautions, protective equipment and emergency procedures:
Observe laboratory safety regulations. Avoid contact with skin and eyes. Do not swallow. Do not pipette by mouth. Do not eat, drink, smoke or apply makeup in areas where specimens or kit reagents are handled. When spilled, absorb with an inert material and put the spilled material in an appropriate waste disposal.
 - Exposure controls / personal protection: Wear protective gloves of nitril rubber or natural latex. Wear protective glasses. Used according to intended use no dangerous reactions known.
 - Conditions to avoid: Since substrate solution is light-sensitive. Store in the dark.
 - For disposal of laboratory waste the national or regional legislation has to be observed.
- Observe the guidelines for performing quality control in medical laboratories by assaying control sera.

CONTENTS OF THE KIT

ORG 511		96	Sufficient for 96 determinations
MICROPLATE	1		One divisible microplate consisting of 12 modules of 8 wells each. Ready to use. Color code on module
CALIBRATOR A	1x 1.5 ml		Calibrator A 0 U/ml, containing serum/buffer matrix (PBS, BSA, detergent, NaN3 0.09%), yellow. Ready to use.
CALIBRATOR B	1x 1.5 ml		Calibrator B 12.5 U/ml, containing RNP/Sm antibodies in a serum/buffer matrix (PBS, BSA, detergent, NaN3 0.09%), yellow. Ready to use.
CALIBRATOR C	1x 1.5 ml		Calibrator C 25 U/ml, containing RNP/Sm antibodies in a serum/buffer matrix (PBS, BSA, detergent, NaN3 0.09%), yellow. Ready to use.
CALIBRATOR D	1x 1.5 ml		Calibrator D 50 U/ml, containing RNP/Sm antibodies in a serum/buffer matrix (PBS, BSA, detergent, NaN3 0.09%), yellow. Ready to use.
CALIBRATOR E	1x 1.5 ml		Calibrator E 100 U/ml, containing RNP/Sm antibodies in a serum/buffer matrix (PBS, BSA, NaN3 0.09%), yellow. Ready to use.
CALIBRATOR F	1x 1.5 ml		Calibrator F 200 U/ml, containing RNP/Sm antibodies in a serum/buffer matrix (PBS, BSA, detergent, NaN3 0.09%), yellow. Ready to use.
CONTROL +	1x 1.5 ml		Control positive, containing RNP/Sm antibodies in a serum/buffer matrix (PBS, BSA, detergent, NaN3 0.09%), yellow. Ready to use. The concentration is specified on the certificate of analysis.
CONTROL -	1x 1.5 ml		Control negative, containing RNP/Sm antibodies in a serum/buffer matrix (PBS, BSA, detergent, NaN3 0.09%), yellow. Ready to use. The concentration is specified on the certificate of analysis.
DILUENT	20 ml		Sample Buffer P, containing PBS, BSA, detergent, preservative sodium azide 0.09%, yellow, concentrate (5 x).
CONJUGATE	15 ml		Enzyme Conjugate containing anti-human IgG antibodies, HRP labelled; PBS, BSA, detergent, preservative PROCLIN 0.05%, light red. Ready to use.
TMB	15 ml		TMB Substrate; containing 3,3', 5,5'- Tetramethylbenzidin, colorless. Ready to use.
STOP	15 ml		Stop solution; contains acid. Ready to use.
WASH	20 ml		Wash Buffer, containing Tris, detergent, preservative sodium azide 0.09%; 50 x conc.
	1		Certificate of Analysis

MATERIALS REQUIRED

- Microplate reader capable of endpoint measurements at 450 nm; optional: reference filter at 620 nm
- Data reduction software
- Multi-channel dispenser or repeatable pipette for 100 µl
- Vortex mixer
- Pipettes for 10 µl, 100 µl and 1000 µl
- Laboratory timing device
- Distilled or deionised water
- Measuring cylinder for 1000 ml and 100 ml
- Plastic container for storage of the wash solution

This ELISA assay is suitable for use on open automated ELISA processors. Each assay has to be validated on the respective automated system. Detailed information is provided upon request.

SPECIMEN COLLECTION, STORAGE AND HANDLING

- Collect whole blood specimens using acceptable medical techniques to avoid hemolysis.
- Allow blood to clot and separate the serum or plasma by centrifugation.
- Test serum should be clear and non-hemolyzed. Contamination by hemolysis or lipemia should be avoided, but does not interfere with this assay.
- Specimens may be refrigerated at 2-8°C for up to five days or stored at -20°C up to six months.
- Avoid repetitive freezing and thawing of serum or plasma samples. This may result in variable loss of antibody activity.
- Testing of heat-inactivated sera is not recommended.

STORAGE AND STABILITY

- Store test kit at 2-8°C in the dark.
- Do not expose reagents to heat, sun, or strong light during storage and usage.
- Store microplate sealed and dessicated in the clip bag provided.
- Shelf life of the unopened test kit is 18 months from day of production.
Unopened reagents are stable until expiration of the kit. See labels for individual batch.
- Diluted Wash Buffer and Sample Buffer are stable for at least 30 days when stored at 2-8°C.
We recommend consumption on the same day.

PROCEDURAL NOTES

- Do not use kit components beyond their expiration dates.
- Do not interchange kit components from different lots and products.
- All materials must be at room temperature (20-28°C) prior to use.
- Prepare all reagents and samples. Once started, perform the test without interruption.
- Double determinations may be done. By this means pipetting errors may become obvious.
- Perform the assay steps only in the order indicated.
- Always use fresh sample dilutions.
- Pipette all reagents and samples into the bottom of the wells.
- To avoid carryover or contamination, change the pipette tip between samples and different kit controls.
- Wash microwells thoroughly and remove the last droplets of wash buffer.
- All incubation steps must be accurately timed.
- Do not re-use microplate wells.

PREPARATION OF REAGENTS

WASH
Dilute the contents of one vial of the buffered wash solution concentrate (50x) with distilled or deionised water to a final volume of 1000 ml prior to use.

DILUENT
Sample Buffer P: Prior to use dilute the contents (20 ml) of one vial of sample buffer 5x concentrate with distilled or deionised water to a final volume of 100 ml.

Preparation of samples

Dilute patient samples 1:100 before the assay: Put 990 µl of prediluted sample buffer in a polystyrene tube and add 10 µl of sample. Mix well. Note: Calibrators / Controls are ready to use and need not be diluted.

TEST PROCEDURE

Prepare enough microplate modules for all calibrators / controls and patient samples.

- Pipette **100 µl** of calibrators, controls and prediluted patient samples into the wells.
Incubate for **30 minutes** at room temperature (20-28 °C).
Discard the contents of the microwells and **wash 3 times** with **300 µl** of wash solution.
- Dispense **100 µl** of enzyme conjugate into each well.
Incubate for **15 minutes** at room temperature.
Discard the contents of the microwells and **wash 3 times** with **300 µl** of wash solution.
- Dispense **100 µl** of TMB substrate solution into each well.
Incubate for **15 minutes** at room temperature
- Add 100 µl** of stop solution to each well of the modules
Incubate for **5 minutes** at room temperature.
Read the optical density at 450 nm (reference 600-690nm) and calculate the results.
The developed colour is stable for at least 30 minutes. Read during this time.

Example for a pipetting scheme:

	1	2	3	4	5	6	7	8	9	10	11	12
A	A	P1										
B	B	P2										
C	C	P3										
D	D											
E	E											
F	F											
G	C+											
H	C-											

P1, ... patient sample A-F calibrators C+, C- controls

VALIDATION

Test results are valid if the optical densities at 450 nm for calibrators / controls and the results for controls comply with the reference ranges indicated on the Certificate of Analysis enclosed in each test kit.
If these quality control criteria are not met the assay run is invalid and should be repeated.

CALCULATION OF RESULTS

For quantitative results plot the optical density of each calibrator versus the calibrator concentration to create a calibration curve. The concentration of patient samples may then be estimated from the calibration curve by interpolation.

Using data reduction software a 4-Parameter-Fit with lin-log coordinates for optical density and concentration is the data reduction method of choice.

PERFORMANCE CHARACTERISTICS

Calibration

The assay system is calibrated against the internationally recognized reference sera from CDC, Atlanta USA.

Measuring range

The calculation range of this ELISA assay is 0 - 200 U/ml

Expected values

In a normal range study with samples from healthy blood donors the following ranges have been established with this ELISA assay: Cut-off 25 U/ml

Interpretation of results

Negative: < 15 U/ml
Borderline: 15 - 25 U/ml
Positive: > 25 U/ml

Linearity

Patient samples containing high levels of specific antibody were serially diluted in sample buffer to demonstrate the dynamic range of the assay and the upper / lower end of linearity. Activity for each dilution was calculated from the calibration curve using a 4-Parameter-Fit with lin-log coordinates.

Sample	Dilution	Observed U/ml	Expected U/ml	O/E [%]
1	1:100	161.4	161.4	100
.	1:200	78.0	80.7	97
.	1:400	39.7	40.4	98
.	1:800	20.1	20.2	100
2	1:100	167.2	167.2	100
.	1:200	83.7	83.6	100
.	1:400	41.5	41.8	99
.	1:800	20.8	20.9	100

Limit of detection

Functional sensitivity was determined to be: 1 U/ml

Reproducibility

Intra-assay precision: Coefficient of variation (CV) was calculated for each of three samples from the results of 24 determinations in a single run. Results for precision-within-assay are shown in the table below.

Inter-assay precision: Coefficient of variation (CV) was calculated for each of three samples from the results of 6 determinations in 5 different runs. Results for run-to-run precision are shown in the table below.

Intra-Assay		
Sample	Mean U/ml	CV %
1	65.6	4.1
2	101.9	5.9
3	182.0	1.8

Inter-Assay		
Sample	Mean U/ml	CV %
1	33.3	4.2
2	109.0	3.1
3	176.8	2.9

Interfering substances

No interference has been observed with haemolytic (up to 1000 mg/dl) or lipemic (up to 3 g/dl triglycerides) sera or plasma, or bilirubin (up to 40 mg/dl) containing sera or plasma. Nor have any interfering effects been observed with the use of anticoagulants (Citrates, EDTA, Heparine). However for practical reasons it is recommended that grossly hemolyzed or lipemic samples should be avoided.

Study results

<u>Study population</u>	<u>n</u>	<u>n Pos</u>	<u>%</u>
SLE	70	37	52.9
MCTD	30	29	96.7
Rheumatoid arthritis	20	3	15.0
Normal human sera	100	2	2.0

		Clinical Diagnosis		
		POS	NEG	
ORG 511	POS	66	5	
	NEG	34	115	
		100	120	220
Sensitivity:		66.0	%	
Specificity:		95.8	%	
Overall agreement:		82.3	%	

LIMITATIONS OF THE PROCEDURE

This assay is a diagnostic aid. A definite clinical diagnosis should not be based on the results of a single test, but should be made by the physician after all clinical and laboratory findings have been evaluated concerning the entire clinical picture of the patient. Also every decision for therapy should be taken individually.

The above pathological and normal reference ranges for antibodies in patient samples should be regarded as recommendations only. Each laboratory should establish its own ranges according to ISO 15189 or other applicable laboratory guidelines.

REFERENCES

1. Alba P, Bento L, Cuadrado MJ, Karim Y, Tungekar MF, Abbs I et al. Anti-dsDNA, anti-Sm antibodies, and the lupus anticoagulant: significant factors associated with lupus nephritis. *Ann Rheum Dis* 2003; 62(6):556-560.
2. Antico A, Platzgummer S, Bassetti D, Bizzaro N, Tozzoli R, Villalta D. Diagnosing systemic lupus erythematosus: new-generation immunoassays for measurement of anti-dsDNA antibodies are an effective alternative to the Farr technique and the Crithidia luciliae immunofluorescence test. *Lupus* 2010; 19(8):906-912.
3. Brouwer R, Hengstman GJ, Vree EW, Ehrfeld H, Bozic B, Ghirardello A et al. Autoantibody profiles in the sera of European patients with myositis. *Ann Rheum Dis* 2001; 60(2):116-123.
4. Castro C, Gourley M. Diagnostic testing and interpretation of tests for autoimmunity. *J Allergy Clin Immunol* 2010; 125(2 Suppl 2):S238-S247.
5. Defendenti C, Atzeni F, Spina MF, Grosso S, Cereda A, Guercilena G et al. Clinical and laboratory aspects of Ro/SSA-52 autoantibodies. *Autoimmun Rev* 2011; 10(3):150-154.
6. Eriksson C, Kokkonen H, Johansson M, Hallmans G, Wadell G, Rantapaa-Dahlqvist S. Autoantibodies predate the onset of Systemic Lupus Erythematosus in northern Sweden. *Arthritis Research & Therapy* 2011; 13(1):R30.
7. Haugbro K, Nossent JC, Winkler T, Figenschau Y, Rekvig OP. Anti-dsDNA antibodies and disease classification in antinuclear antibody positive patients: the role of analytical diversity. *Ann Rheum Dis JID* - 0372355 2004; 63(4):386-394.
8. Ippolito A, Wallace DJ, Gladman D, Fortin PR, Urowitz M, Werth V et al. Autoantibodies in systemic lupus erythematosus: comparison of historical and current assessment of seropositivity. *Lupus* 2011; 20(3):250-255.
9. Isenberg DA, Manson JJ, Ehrenstein MR, Rahman A. Fifty years of anti-ds DNA antibodies: are we approaching journey's end? *Rheumatology (Oxford)* 2007; 46(7):1052-1056.
10. Kattah NH, Kattah MG, Utz PJ. The U1-snRNP complex: structural properties relating to autoimmune pathogenesis in rheumatic diseases. *Immunol Rev* 2010; 233(1):126-145.
11. Kumar Y, Bhatia A, Minz RW. Antinuclear antibodies and their detection methods in diagnosis of connective tissue diseases: a journey revisited. *Diagn Pathol* 2009; 4:1.
12. Meroni PL, Schur PH. ANA screening: an old test with new recommendations. *Ann Rheum Dis* 2010; 69:1420-1422.
13. Petri M, Magder L. Classification criteria for systemic lupus erythematosus: a review. *Lupus* 2004; 13(11):829-837.
14. Poole BD, Schneider RI, Guthridge JM, Velte CA, Reichlin M, Harley JB et al. Early targets of nuclear RNP humoral autoimmunity in human systemic lupus erythematosus. *Arthritis Rheum* 2009; 60(3):848-859.
15. Putova I, Dostal C, Becvar R. Prevalence of antinucleosome antibodies by enzyme-linked immunosorbent assays in patients with systemic lupus erythematosus and other autoimmune systemic diseases. *Ann N Y Acad Sci* 2007; 1109:275-286.
16. Reveille JD. Predictive value of autoantibodies for activity of systemic lupus erythematosus. *Lupus JID* - 9204265 2004; 13(5):290-297.
17. Simon JA, Cabiedes J, Ortiz E, Alcocer-Varela J, Sanchez-Guerrero J. Anti-nucleosome antibodies in patients with systemic lupus erythematosus of recent onset. Potential utility as a diagnostic tool and disease activity marker. *Rheumatology (Oxford)* 2004; 43(2):220-224.
18. Sinclair D, Saas M, Williams D, Hart M, Goswami R. Can an ELISA replace immunofluorescence for the detection of anti-nuclear antibodies?--The routine use of anti-nuclear antibody screening ELISAs. *Clin Lab* 2007; 53(3-4):183-191.
19. Tozzoli R, Bizzaro N, Tonutti E, Villalta D, Bassetti D, Manoni F et al. Guidelines for the laboratory use of autoantibody tests in the diagnosis and monitoring of autoimmune rheumatic diseases. *Am J Clin Pathol* 2002; 117(2):316-324.
20. Maidhof W., Hiliar O. Lupus: an overview of the disease and management options. *P T* 2012; 37(4):240-9.
21. Hahn BH, McMahon MA, Wilkinson A, Wallace WD, Daikh DI, Fitzgerald JD et al. American College of Rheumatology guidelines for screening, treatment, and management of lupus nephritis. *Arthritis Care Res (Hoboken)* 2012; 64(6):797-808.

Notice to the user (European Union):

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the EU Member State in which the user and/or the patient is established .

Change Control

Former version: *ORG 511_IFU_EN_QM113138_2013-12-16_1.2* Reason for revision: *Introduction electronic IFU on homepage*

- ① Pipet **100 µl** calibrator, control or patient sample
 - Incubate for **30 minutes** at room temperature
 - Discard the contents of the wells and wash 3 times with **300 µl** wash solution
- ② Pipet **100 µl** enzyme conjugate
 - Incubate for **15 minutes** at room temperature
 - Discard the contents of the wells and wash 3 times with **300 µl** wash solution
- ③ Pipet **100 µl** substrate solution
 - Incubate for **15 minutes** at room temperature
- ④ Add **100 µl** stop solution
 - Leave untouched for **5 minutes**
 - Read at **450 nm**

ORGENTEC Diagnostika GmbH

Carl-Zeiss-Straße 49-51

55129 Mainz - Germany

Phone: +49 (0) 61 31 / 92 58-0

Fax: +49 (0) 61 31 / 92 58-58

Internet: www.orgentec.com

633_3

ORG 633 Anti-Centromere B

INTENDED PURPOSE

Anti-Centromere B is an ELISA test system for the quantitative measurement of IgG class autoantibodies against centromere B in human serum or plasma. This product is intended for professional in vitro diagnostic use only.

The test is used as an aid in the differential diagnosis of inflammatory autoimmune diseases, e.g. CREST syndrome. Evaluation of a test result should always take into account all clinical and laboratory diagnostic findings.

SYMBOLS USED ON LABELS

In vitro diagnostic medical device

Manufacturer

Catalogue number

Sufficient for 96 determinations

Batch code

Use by

Temperature limitation

Keep away from sunlight

Do not reuse

Date of manufacture

CE marked according to 98/79/EC


Consult instructions for use

Electronic Instruction For Use: version

Microplate


Calibrator


Calibrator

Calibrator

Calibrator

Calibrator

Control positive

Control negative

Sample Buffer P

Enzyme Conjugate

TMB Substrate

Stop solution

Wash Buffer

Ready to use

PRINCIPLE OF THE TEST

Recombinant centromere protein B is bound to microwells.

The determination is based on an indirect enzyme linked immune reaction with the following steps:

Specific antibodies in the patient sample bind to the antigen coated on the surface of the reaction wells. After incubation, a washing step removes unbound and unspecifically bound serum or plasma components. Subsequently added enzyme conjugate binds to the immobilized antibody-antigen-complexes. After incubation, a second washing step removes unbound enzyme conjugate. After addition of substrate solution the bound enzyme conjugate hydrolyses the substrate forming a blue coloured product. Addition of an acid stops the reaction generating a yellow end-product. The intensity of the yellow color correlates with the concentration of the antibody-antigen-complex and can be measured photometrically at 450 nm.

WARNINGS AND PRECAUTIONS

- All reagents of this kit are intended for professional in vitro diagnostic use only.
- Components containing human serum were tested and found negative for HBsAg, HCV, HIV1 and HIV2 by FDA approved methods. No test can guarantee the absence of HBsAg, HCV, HIV1 or HIV2, and so all human serum based reagents in this kit must be handled as though capable of transmitting infection.
- Bovine serum albumin (BSA) used in components has been tested for BSE and found negative.
- Avoid contact with the substrate TMB (3,3',5,5'-Tetramethyl-benzidine).
- Stop solution contains acid, classification is non-hazardous. Avoid contact with skin.
- Control, sample buffer and wash buffer contain sodium azide 0.09% as preservative. This concentration is classified as non-hazardous.
- Enzyme conjugate contains ProClin 300 0.05% as preservative. This concentration is classified as non-hazardous.

During handling of all reagents, controls and serum samples observe the existing regulations for laboratory safety regulations and good laboratory practice:

- First aid measures: In case of skin contact, immediately wash thoroughly with water and soap. Remove contaminated clothing and shoes and wash before reuse. If system fluid comes into contact with skin, wash thoroughly with water. After contact with the eyes carefully rinse the opened eye with running water for at least 10 minutes. Get medical attention if necessary.
 - Personal precautions, protective equipment and emergency procedures:
Observe laboratory safety regulations. Avoid contact with skin and eyes. Do not swallow. Do not pipette by mouth. Do not eat, drink, smoke or apply makeup in areas where specimens or kit reagents are handled. When spilled, absorb with an inert material and put the spilled material in an appropriate waste disposal.
 - Exposure controls / personal protection: Wear protective gloves of nitril rubber or natural latex. Wear protective glasses. Used according to intended use no dangerous reactions known.
 - Conditions to avoid: Since substrate solution is light-sensitive. Store in the dark.
 - For disposal of laboratory waste the national or regional legislation has to be observed.
- Observe the guidelines for performing quality control in medical laboratories by assaying control sera.

CONTENTS OF THE KIT

ORG 633	▽ 96	Sufficient for 96 determinations
MICROPLATE	1	One divisible microplate consisting of 12 modules of 8 wells each. Ready to use. Product code on module: CEN
CALIBRATOR A	1x 1.5 ml	Calibrator A 0 U/ml, containing serum/buffer matrix (PBS, BSA, detergent, NaN3 0.09%), yellow. Ready to use.
CALIBRATOR B	1x 1.5 ml	Calibrator B 10 U/ml, containing centromere B antibodies in a serum/buffer matrix (PBS, BSA, detergent, NaN3 0.09%), yellow. Ready to use.
CALIBRATOR C	1x 1.5 ml	Calibrator C 30 U/ml, containing centromere B antibodies in a serum/buffer matrix (PBS, BSA, detergent, NaN3 0.09%), yellow. Ready to use.
CALIBRATOR D	1x 1.5 ml	Calibrator D 100 U/ml, containing centromere B antibodies in a serum/buffer matrix (PBS, BSA, detergent, NaN3 0.09%), yellow. Ready to use.
CALIBRATOR E	1x 1.5 ml	Calibrator E 300 U/ml, containing centromere B antibodies in a serum/buffer matrix (PBS, BSA, NaN3 0.09%), yellow. Ready to use.
CONTROL +	1x 1.5 ml	Control positive, containing centromere B antibodies in a serum/buffer matrix (PBS, BSA, detergent, NaN3 0.09%), yellow. Ready to use. The concentration is specified on the certificate of analysis.
CONTROL -	1x 1.5 ml	Control negative, containing centromere B antibodies in a serum/buffer matrix (PBS, BSA, detergent, NaN3 0.09%), yellow. Ready to use. The concentration is specified on the certificate of analysis.
DILUENT	20 ml	Sample Buffer P, containing PBS, BSA, detergent, preservative sodium azide 0.09%, yellow, concentrate (5 x).
CONJUGATE	15 ml	Enzyme Conjugate containing anti-human IgG antibodies, HRP labelled; PBS, BSA, detergent, preservative PROCLIN 0.05%, light red. Ready to use.
TMB	15 ml	TMB Substrate; containing 3,3', 5,5'- Tetramethylbenzidin, colorless. Ready to use.
STOP	15 ml	Stop solution; contains acid. Ready to use.
WASH	20 ml	Wash Buffer, containing Tris, detergent, preservative sodium azide 0.09%; 50 x conc.
i	1	Certificate of Analysis

MATERIALS REQUIRED

- Microplate reader capable of endpoint measurements at 450 nm; optional: reference filter at 620 nm
- Data reduction software
- Multi-channel dispenser or repeatable pipette for 100 µl
- Vortex mixer
- Pipettes for 10 µl, 100 µl and 1000 µl
- Laboratory timing device
- Distilled or deionised water
- Measuring cylinder for 1000 ml and 100 ml
- Plastic container for storage of the wash solution

This ELISA assay is suitable for use on open automated ELISA processors. Each assay has to be validated on the respective automated system. Detailed information is provided upon request.

SPECIMEN COLLECTION, STORAGE AND HANDLING

- Collect whole blood specimens using acceptable medical techniques to avoid hemolysis.
- Allow blood to clot and separate the serum or plasma by centrifugation.
- Test serum should be clear and non-hemolyzed. Contamination by hemolysis or lipemia should be avoided, but does not interfere with this assay.
- Specimens may be refrigerated at 2-8°C for up to five days or stored at -20°C up to six months.
- Avoid repetitive freezing and thawing of serum or plasma samples. This may result in variable loss of antibody activity.
- Testing of heat-inactivated sera is not recommended.

STORAGE AND STABILITY

- Store test kit at 2-8°C in the dark.
- Do not expose reagents to heat, sun, or strong light during storage and usage.

- Store microplate sealed and dessicated in the clip bag provided.
- Shelf life of the unopened test kit is 18 months from day of production.
Unopened reagents are stable until expiration of the kit. See labels for individual batch.
- Diluted Wash Buffer and Sample Buffer are stable for at least 30 days when stored at 2-8°C.
We recommend consumption on the same day.

PROCEDURAL NOTES

- Do not use kit components beyond their expiration dates.
- Do not interchange kit components from different lots and products.
- All materials must be at room temperature (20-28°C) prior to use.
- Prepare all reagents and samples. Once started, perform the test without interruption.
- Double determinations may be done. By this means pipetting errors may become obvious.
- Perform the assay steps only in the order indicated.
- Always use fresh sample dilutions.
- Pipette all reagents and samples into the bottom of the wells.
- To avoid carryover or contamination, change the pipette tip between samples and different kit controls.
- Wash microwells thoroughly and remove the last droplets of wash buffer.
- All incubation steps must be accurately timed.
- Do not re-use microplate wells.

PREPARATION OF REAGENTS

WASH

Dilute the contents of one vial of the buffered wash solution concentrate (50x) with distilled or deionised water to a final volume of 1000 ml prior to use.

DILUENT

Sample Buffer P: Prior to use dilute the contents (20 ml) of one vial of sample buffer 5x concentrate with distilled or deionised water to a final volume of 100 ml.

Preparation of samples

Dilute patient samples 1:100 before the assay: Put 990 µl of prediluted sample buffer in a polystyrene tube and add 10 µl of sample. Mix well. Note: Calibrators / Controls are ready to use and need not be diluted.

TEST PROCEDURE

Prepare enough microplate modules for all calibrators / controls and patient samples.

1. Pipette **100 µl** of calibrators, controls and prediluted patient samples into the wells.
Incubate for **30 minutes** at room temperature (20-28 °C).
Discard the contents of the microwells and **wash 3 times** with **300 µl** of wash solution.
2. Dispense **100 µl** of enzyme conjugate into each well.
Incubate for **15 minutes** at room temperature.
Discard the contents of the microwells and **wash 3 times** with **300 µl** of wash solution.
3. Dispense **100 µl** of TMB substrate solution into each well.
Incubate for **15 minutes** at room temperature
4. **Add 100 µl** of stop solution to each well of the modules
Incubate for **5 minutes** at room temperature.

Read the optical density at 450 nm (reference 600-690nm) and calculate the results.
The developed colour is stable for at least 30 minutes. Read during this time.

Example for a pipetting scheme:

[illegible]

P1, ... patient sample A-E calibrators C+, C- controls

VALIDATION

Test results are valid if the optical densities at 450 nm for calibrators / controls and the results for controls comply with the reference ranges indicated on the Certificate of Analysis enclosed in each test kit.
If these quality control criteria are not met the assay run is invalid and should be repeated.

CALCULATION OF RESULTS

For quantitative results plot the optical density of each calibrator versus the calibrator concentration to create a calibration curve. The concentration of patient samples may then be estimated from the calibration curve by interpolation. Using data reduction software a 4-Parameter-Fit with lin-log coordinates for optical density and concentration is the data reduction method of choice.

PERFORMANCE CHARACTERISTICS

Calibration

This assay system is calibrated in relative arbitrary units, since no international reference preparation is available for this assay.

Measuring range

The calculation range of this ELISA assay is 0 - 300 U/ml

Expected values

In a normal range study with samples from healthy blood donors the following ranges have been established with this ELISA assay: Cut-off 10 U/ml

Interpretation of results

Negative: < 10 U/ml
Positive: ≥ 10 U/ml

Linearity

Patient samples containing high levels of specific antibody were serially diluted in sample buffer to demonstrate the dynamic range of the assay and the upper / lower end of linearity. Activity for each dilution was calculated from the calibration curve using a 4-Parameter-Fit with lin-log coordinates.

Sample	Dilution	Observed U/ml	Expected U/ml	O/E [%]
1	1:100	136.8	136.8	100
.	1:200	67.1	68.4	98
.	1:400	35.2	34.2	103
.	1:800	16.9	17.1	99
2	1:100	285.0	285.0	100
.	1:200	139.2	142.5	98
.	1:400	73.5	71.3	103
.	1:800	37.0	35.6	104

Limit of detection

Functional sensitivity was determined to be: 1 U/ml

Reproducibility

Intra-assay precision: Coefficient of variation (CV) was calculated for each of three samples from the results of 24 determinations in a single run. Results for precision-within-assay are shown in the table below.

Inter-assay precision: Coefficient of variation (CV) was calculated for each of three samples from the results of 6 determinations in 5 different runs. Results for run-to-run precision are shown in the table below.

Intra-Assay		
Sample	Mean U/ml	CV %
1	15.2	5.4
2	122.0	4.4
3	220.0	4.7

Inter-Assay		
Sample	Mean U/ml	CV %
1	16.4	5.4
2	125.6	5.0
3	225.4	4.2

Interfering substances

No interference has been observed with haemolytic (up to 1000 mg/dl) or lipemic (up to 3 g/dl triglycerides) sera or plasma, or bilirubin (up to 40 mg/dl) containing sera or plasma. Nor have any interfering effects been observed with the use of anticoagulants (Citrate, EDTA, Heparine). However for practical reasons it is recommended that grossly hemolyzed or lipemic samples should be avoided.

Study results

<u>Study population</u>	<u>n</u>	<u>n Pos</u>	<u>%</u>
CREST syndrome	32	31	96.9
Rheumatoid arthritis	20	1	5.0
Normal human sera	100	6	6.0

		Clinical Diagnosis	
		POS	NEG
ORG 633	POS	31	7
	NEG	1	113
		32	120
			152

Sensitivity: 96.9 %

Specificity: 94.2 %

Overall agreement: 94.7 %

LIMITATIONS OF THE PROCEDURE

This assay is a diagnostic aid. A definite clinical diagnosis should not be based on the results of a single test, but should be made by the physician after all clinical and laboratory findings have been evaluated concerning the entire clinical picture of the patient. Also every decision for therapy should be taken individually.

The above pathological and normal reference ranges for antibodies in patient samples should be regarded as recommendations only. Each laboratory should establish its own ranges according to ISO 15189 or other applicable laboratory guidelines.

REFERENCES

- 1 Alba P, Bento L, Cuadrado MJ, Karim Y, Tungekar MF, Abbs I et al. Anti-dsDNA, anti-Sm antibodies, and the lupus anticoagulant: significant factors associated with lupus nephritis. *Ann Rheum Dis* 2003; 62(6):556-560.
- 2 Antico A, Platzgummer S, Bassetti D, Bizzaro N, Tozzoli R, Villalta D. Diagnosing systemic lupus erythematosus: new-generation immunoassays for measurement of anti-dsDNA antibodies are an effective alternative to the Farr technique and the Crithidia luciliae immunofluorescence test. *Lupus* 2010; 19(8):906-912.
- 3 Brouwer R, Hengstman GJ, Vree EW, Ehrfeld H, Bozic B, Ghirardello A et al. Autoantibody profiles in the sera of European patients with myositis. *Ann Rheum Dis* 2001; 60(2):116-123.
- 4 Castro C, Gourley M. Diagnostic testing and interpretation of tests for autoimmunity. *J Allergy Clin Immunol* 2010; 125(2 Suppl 2):S238-S247.
- 5 Defendenti C, Atzeni F, Spina MF, Grosso S, Cereda A, Guercilena G et al. Clinical and laboratory aspects of Ro/SSA-52 autoantibodies. *Autoimmun Rev* 2011; 10(3):150-154.
- 6 Eriksson C, Kokkonen H, Johansson M, Hallmans G, Wadell G, Rantapaa-Dahlqvist S. Autoantibodies predate the onset of Systemic Lupus Erythematosus in northern Sweden. *Arthritis Research & Therapy* 2011; 13(1):R30.
- 7 Haugbro K, Nossent JC, Winkler T, Figenschau Y, Rekvig OP. Anti-dsDNA antibodies and disease classification in antinuclear antibody positive patients: the role of analytical diversity. *Ann Rheum Dis JID* - 0372355 2004; 63(4):386-394.
- 8 Ippolito A, Wallace DJ, Gladman D, Fortin PR, Urowitz M, Werth V et al. Autoantibodies in systemic lupus erythematosus: comparison of historical and current assessment of seropositivity. *Lupus* 2011; 20(3):250-255.
- 9 Isenberg DA, Manson JJ, Ehrenstein MR, Rahman A. Fifty years of anti-ds DNA antibodies: are we approaching journey's end? *Rheumatology (Oxford)* 2007; 46(7):1052-1056.
10. Kattah NH, Kattah MG, Utz PJ. The U1-snRNP complex: structural properties relating to autoimmune pathogenesis in rheumatic diseases. *Immunol Rev* 2010; 233(1):126-145.
- 11.Kumar Y, Bhatia A, Minz RW. Antinuclear antibodies and their detection methods in diagnosis of connective tissue diseases: a journey revisited. *Diagn Pathol* 2009; 4:1.
12. Meroni PL, Schur PH. ANA screening: an old test with new recommendations. *Ann Rheum Dis* 2010; 69:1420-1422.
13. Petri M, Magder L. Classification criteria for systemic lupus erythematosus: a review. *Lupus* 2004; 13(11):829-837.
14. Poole BD, Schneider RI, Guthridge JM, Velte CA, Reichlin M, Harley JB et al. Early targets of nuclear RNP humoral autoimmunity in human systemic lupus erythematosus. *Arthritis Rheum* 2009; 60(3):848-859.
- 15.Putova I, Dostal C, Becvar R. Prevalence of antinucleosome antibodies by enzyme-linked immunosorbent assays in patients with systemic lupus erythematosus and other autoimmune systemic diseases. *Ann N Y Acad Sci* 2007; 1109:275-286.
16. Reveille JD. Predictive value of autoantibodies for activity of systemic lupus erythematosus. *Lupus JID* - 9204265 2004; 13(5):290-297.
17. Simon JA, Cabiedes J, Ortiz E, Alcocer-Varela J, Sanchez-Guerrero J. Anti-nucleosome antibodies in patients with systemic lupus erythematosus of recent onset. Potential utility as a diagnostic tool and disease activity marker. *Rheumatology (Oxford)* 2004; 43(2):220-224.
18. Sinclair D, Saas M, Williams D, Hart M, Goswami R. Can an ELISA replace immunofluorescence for the detection of anti-nuclear antibodies?--The routine use of anti-nuclear antibody screening ELISAs. *Clin Lab* 2007; 53(3-4):183-191.
19. Tozzoli R, Bizzaro N, Tonutti E, Villalta D, Bassetti D, Manoni F et al. Guidelines for the laboratory use of autoantibody tests in the diagnosis and monitoring of autoimmune rheumatic diseases. *Am J Clin Pathol* 2002; 117(2):316-324.
20. Maidhof W., Hilius O. Lupus: an overview of the disease and management options. *P T* 2012; 37(4):240-9.

Notice to the user (European Union):

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the EU Member State in which the user and/or the patient is established .

Change Control

Former version: ORG 633_IFU_EN_QM113209_2016-08-16_1.3 Reason for revision: *Introduction electronic IFU on homepage*

- 1 Pipet **100 µl** calibrator, control or patient sample
→ Incubate for **30 minutes** at room temperature
→ Discard the contents of the wells and wash 3 times with **300 µl** wash solution
- 2 Pipet **100 µl** enzyme conjugate
→ Incubate for **15 minutes** at room temperature
→ Discard the contents of the wells and wash 3 times with **300 µl** wash solution
- 3 Pipet **100 µl** substrate solution
→ Incubate for **15 minutes** at room temperature
- 4 Add **100 µl** stop solution
→ Leave untouched for **5 minutes**
→ Read at **450 nm**

ORGENTEC Diagnostika GmbH


Carl-Zeiss-Straße 49-51

55129 Mainz - Germany

Phone: +49 (0) 61 31 / 92 58-0

Fax: +49 (0) 61 31 / 92 58-58

Internet: www.orgentec.com

ORG 516 AMA-M2

INTENDED PURPOSE

AMA-M2 is an ELISA test system for the quantitative measurement of IgG class autoantibodies against mitochondrial M2 subtype antigen in human serum or plasma. This product is intended for professional in vitro diagnostic use only.

The test is used as an aid in the differential diagnosis of primary biliary cirrhosis (PBC). In patients with other autoimmune diseases occurrence of AMA antibodies may be related to the development or association of PBC. Evaluation of a test result should always take into account all clinical and laboratory diagnostic findings.

SYMBOLS USED ON LABELS

In vitro diagnostic medical device

Manufacturer

Catalogue number

Sufficient for 96 determinations

Batch code

Use by

Temperature limitation

Keep away from sunlight

Do not reuse

Date of manufacture

CE marked according to 98/79/EC


Consult instructions for use

Electronic Instruction For Use: version


Microplate

Calibrator

Calibrator


Calibrator

Calibrator

Calibrator

Calibrator

Control positive

Control negative

Sample Buffer P

Enzyme Conjugate

TMB Substrate

Stop solution

Wash Buffer

Ready to use

PRINCIPLE OF THE TEST

Highly purified mitochondrial M2 subtype (PDC-E2, BCOADC-E2, OGDC-E2) antigen is bound to microwells.

The determination is based on an indirect enzyme linked immune reaction with the following steps:

Specific antibodies in the patient sample bind to the antigen coated on the surface of the reaction wells. After incubation, a washing step removes unbound and unspecifically bound serum or plasma components. Subsequently added enzyme conjugate binds to the immobilized antibody-antigen-complexes. After incubation, a second washing step removes unbound enzyme conjugate. After addition of substrate solution the bound enzyme conjugate hydrolyses the substrate forming a blue coloured product. Addition of an acid stops the reaction generating a yellow end-product. The intensity of the yellow color correlates with the concentration of the antibody-antigen-complex and can be measured photometrically at 450 nm.



WARNINGS AND PRECAUTIONS

- All reagents of this kit are intended for professional in vitro diagnostic use only.
- Components containing human serum were tested and found negative for HBsAg, HCV, HIV1 and HIV2 by FDA approved methods. No test can guarantee the absence of HBsAg, HCV, HIV1 or HIV2, and so all human serum based reagents in this kit must be handled as though capable of transmitting infection.
- Bovine serum albumin (BSA) used in components has been tested for BSE and found negative.
- Avoid contact with the substrate TMB (3,3',5,5'-Tetramethyl-benzidine).
- Stop solution contains acid, classification is non-hazardous. Avoid contact with skin.
- Control, sample buffer and wash buffer contain sodium azide 0.09% as preservative. This concentration is classified as non-hazardous.
- Enzyme conjugate contains ProClin 300 0.05% as preservative. This concentration is classified as non-hazardous.

During handling of all reagents, controls and serum samples observe the existing regulations for laboratory safety regulations and good laboratory practice:

- First aid measures: In case of skin contact, immediately wash thoroughly with water and soap. Remove contaminated clothing and shoes and wash before reuse. If system fluid comes into contact with skin, wash thoroughly with water. After contact with the eyes carefully rinse the opened eye with running water for at least 10 minutes. Get medical attention if necessary.
 - Personal precautions, protective equipment and emergency procedures:
Observe laboratory safety regulations. Avoid contact with skin and eyes. Do not swallow. Do not pipette by mouth. Do not eat, drink, smoke or apply makeup in areas where specimens or kit reagents are handled. When spilled, absorb with an inert material and put the spilled material in an appropriate waste disposal.
 - Exposure controls / personal protection: Wear protective gloves of nitril rubber or natural latex. Wear protective glasses. Used according to intended use no dangerous reactions known.
 - Conditions to avoid: Since substrate solution is light-sensitive. Store in the dark.
 - For disposal of laboratory waste the national or regional legislation has to be observed.
- Observe the guidelines for performing quality control in medical laboratories by assaying control sera.

CONTENTS OF THE KIT

ORG 516		96	Sufficient for 96 determinations
	1		One divisible microplate consisting of 12 modules of 8 wells each. Ready to use. Product code on module: AMA
	1x 1.5 ml	Calibrator A	0 IU/ml, containing serum/buffer matrix (PBS, BSA, detergent, NaN3 0.09%), yellow. Ready to use.
	1x 1.5 ml	Calibrator B	12.5 IU/ml, containing AMA-M2 antibodies in a serum/buffer matrix (PBS, BSA, detergent, NaN3 0.09%), yellow. Ready to use.
	1x 1.5 ml	Calibrator C	25 IU/ml, containing AMA-M2 antibodies in a serum/buffer matrix (PBS, BSA, detergent, NaN3 0.09%), yellow. Ready to use.
	1x 1.5 ml	Calibrator D	50 IU/ml, containing AMA-M2 antibodies in a serum/buffer matrix (PBS, BSA, detergent, NaN3 0.09%), yellow. Ready to use.
	1x 1.5 ml	Calibrator E	100 IU/ml, containing AMA-M2 antibodies in a serum/buffer matrix (PBS, BSA, NaN3 0.09%), yellow. Ready to use.
	1x 1.5 ml	Calibrator F	200 IU/ml, containing AMA-M2 antibodies in a serum/buffer matrix (PBS, BSA, detergent, NaN3 0.09%), yellow. Ready to use.
	1x 1.5 ml	Control positive	containing AMA-M2 antibodies in a serum/buffer matrix (PBS, BSA, detergent, NaN3 0.09%), yellow. Ready to use. The concentration is specified on the certificate of analysis.
	1x 1.5 ml	Control negative	containing AMA-M2 antibodies in a serum/buffer matrix (PBS, BSA, detergent, NaN3 0.09%), yellow. Ready to use. The concentration is specified on the certificate of analysis.
	20 ml	Sample Buffer P	containing PBS, BSA, detergent, preservative sodium azide 0.09%, yellow, concentrate (5 x).
	15 ml	Enzyme Conjugate	containing anti-human IgG antibodies, HRP labelled; PBS, BSA, detergent, preservative PROCLIN 0.05%, light red. Ready to use.
	15 ml	TMB Substrate	; containing 3,3', 5,5'- Tetramethylbenzidin, colorless. Ready to use.
	15 ml	Stop solution	; contains acid. Ready to use.
	20 ml	Wash Buffer	containing Tris, detergent, preservative sodium azide 0.09%; 50 x conc.
	1		Certificate of Analysis

MATERIALS REQUIRED

- Microplate reader capable of endpoint measurements at 450 nm; optional: reference filter at 620 nm
- Data reduction software
- Multi-channel dispenser or repeatable pipette for 100 µl
- Vortex mixer
- Pipettes for 10 µl, 100 µl and 1000 µl
- Laboratory timing device
- Distilled or deionised water
- Measuring cylinder for 1000 ml and 100 ml
- Plastic container for storage of the wash solution

This ELISA assay is suitable for use on open automated ELISA processors. Each assay has to be validated on the respective automated system. Detailed information is provided upon request.

SPECIMEN COLLECTION, STORAGE AND HANDLING

- Collect whole blood specimens using acceptable medical techniques to avoid hemolysis.
- Allow blood to clot and separate the serum or plasma by centrifugation.
- Test serum should be clear and non-hemolyzed. Contamination by hemolysis or lipemia should be avoided, but does not interfere with this assay.
- Specimens may be refrigerated at 2-8°C for up to five days or stored at -20°C up to six months.
- Avoid repetitive freezing and thawing of serum or plasma samples. This may result in variable loss of antibody activity.
- Testing of heat-inactivated sera is not recommended.

STORAGE AND STABILITY

- Store test kit at 2-8°C in the dark.
- Do not expose reagents to heat, sun, or strong light during storage and usage.
- Store microplate sealed and dessicated in the clip bag provided.
- Shelf life of the unopened test kit is 18 months from day of production.
Unopened reagents are stable until expiration of the kit. See labels for individual batch.
- Diluted Wash Buffer and Sample Buffer are stable for at least 30 days when stored at 2-8°C.
We recommend consumption on the same day.

PROCEDURAL NOTES

- Do not use kit components beyond their expiration dates.
- Do not interchange kit components from different lots and products.
- All materials must be at room temperature (20-28°C) prior to use.
- Prepare all reagents and samples. Once started, perform the test without interruption.
- Double determinations may be done. By this means pipetting errors may become obvious.
- Perform the assay steps only in the order indicated.
- Always use fresh sample dilutions.
- Pipette all reagents and samples into the bottom of the wells.
- To avoid carryover or contamination, change the pipette tip between samples and different kit controls.
- Wash microwells thoroughly and remove the last droplets of wash buffer.
- All incubation steps must be accurately timed.
- Do not re-use microplate wells.

PREPARATION OF REAGENTS

 Dilute the contents of one vial of the buffered wash solution concentrate (50x) with distilled or deionised water to a final volume of 1000 ml prior to use.

 Sample Buffer P: Prior to use dilute the contents (20 ml) of one vial of sample buffer 5x concentrate with distilled or deionised water to a final volume of 100 ml.

Preparation of samples

Dilute patient samples 1:100 before the assay: Put 990 µl of prediluted sample buffer in a polystyrene tube and add 10 µl of sample. Mix well. Note: Calibrators / Controls are ready to use and need not be diluted.

TEST PROCEDURE

Prepare enough microplate modules for all calibrators / controls and patient samples.

- Pipette **100 µl** of calibrators, controls and prediluted patient samples into the wells.
Incubate for **30 minutes** at room temperature (20-28 °C).
Discard the contents of the microwells and **wash 3 times** with **300 µl** of wash solution.
- Dispense **100 µl** of enzyme conjugate into each well.
Incubate for **15 minutes** at room temperature.
Discard the contents of the microwells and **wash 3 times** with **300 µl** of wash solution.
- Dispense **100 µl** of TMB substrate solution into each well.
Incubate for **15 minutes** at room temperature
- Add 100 µl** of stop solution to each well of the modules
Incubate for **5 minutes** at room temperature.
Read the optical density at 450 nm (reference 600-690nm) and calculate the results.
The developed colour is stable for at least 30 minutes. Read during this time.

Example for a pipetting scheme:

	1	2	3	4	5	6	7	8	9	10	11	12
A	A	P1										
B	B	P2										
C	C	P3										
D	D											
E	E											
F	F											
G	C+											
H	C-											

P1, ... patient sample A-F calibrators C+, C- controls

VALIDATION

Test results are valid if the optical densities at 450 nm for calibrators / controls and the results for controls comply with the reference ranges indicated on the Certificate of Analysis enclosed in each test kit.
If these quality control criteria are not met the assay run is invalid and should be repeated.

CALCULATION OF RESULTS

For quantitative results plot the optical density of each calibrator versus the calibrator concentration to create a calibration curve. The concentration of patient samples may then be estimated from the calibration curve by interpolation.

Using data reduction software a 4-Parameter-Fit with lin-log coordinates for optical density and concentration is the data reduction method of choice.

PERFORMANCE CHARACTERISTICS

Calibration

The assay system is calibrated against the international reference preparation WHO 67/183 for AMA-M2 as 100 IU/ml.

Measuring range

The calculation range of this ELISA assay is 0 - 200 IU/ml

Expected values

In a normal range study with samples from healthy blood donors the following ranges have been established with this ELISA assay: Cut-off 10 IU/ml

Interpretation of results

Negative: < 10 IU/ml
Positive: ≥ 10 IU/ml

Linearity

Samples containing high levels of specific antibody were serially diluted in sample buffer to demonstrate the dynamic range of the assay and the upper / lower end of linearity. Activity for each dilution was calculated from the calibration curve using a 4-Parameter-Fit with lin-log coordinates.

Sample	Dilution	Observed IU/ml	Expected IU/ml	O/E [%]
WHO	1:100	108.5	100.0	109
.	1:200	51.2	50.0	102
.	1:400	25.2	25.0	101
.	1:800	12.8	12.5	102
.	1:1600	6.1	6.3	98
.	1:3200	3.1	3.1	99
1	1:100	49.5	49.5	100
.	1:200	25.0	24.8	101
.	1:400	12.2	12.4	99
.	1:800	5.9	6.2	95

Limit of detection

Functional sensitivity was determined to be: 1 IU/ml

Reproducibility

Intra-assay precision: Coefficient of variation (CV) was calculated for each of three samples from the results of 24 determinations in a single run. Results for precision-within-assay are shown in the table below.

Inter-assay precision: Coefficient of variation (CV) was calculated for each of three samples from the results of 6 determinations in 5 different runs. Results for run-to-run precision are shown in the table below.

Intra-Assay		
Sample	Mean IU/ml	CV %
1	39.8	7.0
2	81.3	3.8
3	177.3	3.6

Inter-Assay		
Sample	Mean IU/ml	CV %
1	40.1	6.2
2	84.6	11.8
3	180.4	3.8

Interfering substances

No interference has been observed with haemolytic (up to 1000 mg/dl) or lipemic (up to 3 g/dl triglycerides) sera or plasma, or bilirubin (up to 40 mg/dl) containing sera or plasma. Nor have any interfering effects been observed with the use of anticoagulants (Citrate, EDTA, Heparine). However for practical reasons it is recommended that grossly hemolyzed or lipemic samples should be avoided.

Study results

Study population	n	n Pos	%
Primary biliary cirrhosis (PBC)	143	139	97.2
Rheumatoid Arthritis	60	1	1.7
Normal human sera	267	18	6.7

		Clinical Diagnosis	
		POS	NEG
ORG 516	POS	139	19
	NEG	4	308
		143	327

Sensitivity: 97.2 %

Specificity: 94.2 %

Overall agreement: 95.1 %

LIMITATIONS OF THE PROCEDURE

This assay is a diagnostic aid. A definite clinical diagnosis should not be based on the results of a single test, but should be made by the physician after all clinical and laboratory findings have been evaluated concerning the entire clinical picture of the patient. Also every decision for therapy should be taken individually.

The above pathological and normal reference ranges for antibodies in patient samples should be regarded as recommendations only. Each laboratory should establish its own ranges according to ISO 15189 or other applicable laboratory guidelines.

REFERENCES

1. Berg, P.A. and Klein, R. Diagnose der primär-biliären Zirrhose. IVD Nachrichten 1990; 1/1: 6 -7.
2. Berg, P.A. and Klein, R. Heterogeneity of anti-mitochondrial antibodies. Sem. Liver Dis. 1989; 9: 103 - 116.
3. Berg, P.A. and Klein, R. Immunology of primary biliary cirrhosis. Ballière's Clin.Gastroenterol. 1987; 1: 675 - 706.
4. Baum, H. and Palmer, C. The PBC specific antigen. Mol. Aspects Med. 1985; 8: 201 - 234.
5. Fussey, S.P.M., Guest, J.R., James, O.F W. et al. Identification and analysis of the major M2 autoantigens in primary biliary cirrhosis. PNAS, USA 1988; 85: 8654 - 8658.

Notice to the user (European Union):

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the EU Member State in which the user and/or the patient is established .

Change Control

Former version: *ORG 516_IFU_EN_QM113145_2013-12-16_2.1* Reason for revision: *Introduction electronic IFU on homepage*

- 1 Pipet **100 µl** calibrator, control or patient sample
→ Incubate for **30 minutes** at room temperature
→ Discard the contents of the wells and wash 3 times with **300 µl** wash solution
- 2 Pipet **100 µl** enzyme conjugate
→ Incubate for **15 minutes** at room temperature
→ Discard the contents of the wells and wash 3 times with **300 µl** wash solution
- 3 Pipet **100 µl** substrate solution
→ Incubate for **15 minutes** at room temperature
- 4 Add **100 µl** stop solution
→ Leave untouched for **5 minutes**
→ Read at **450 nm**

Anti-LKM-1 ELISA (IgG)

Test instruction




ORDER NO.	ANTIBODIES AGAINST	IG CLASS	SUBSTRATE	FORMAT
EA 1321-9601 G	LKM-1	IgG	Ag-coated microplate wells	96 x 01 (96)

Indications: The ELISA test kit provides a semiquantitative or quantitative in vitro assay for human anti-bodies of the IgG class against liver-kidney microsomes (LKM) in serum or plasma for the diagnosis of inexplainable increase in transaminases, suspected autoimmune hepatitis.

Application: According to the simplified diagnostic criteria by EM Hennes and colleagues (*International Autoimmune Hepatitis Group*) published in 2008, the detection of autoantibodies against LKM belongs to the routine investigations performed to diagnose autoimmune hepatitis. Antibodies against LKM-1 are mostly observed in children, but may be also present in adult patients with AIH. For delimitation from a virus hepatitis, the parallel determination of the other autoantibodies associated with AIH, such as ANA, pANCA, ASMA or antibodies against LC-1 and SLA/LP is recommended.

Principles of the test: The test kit contains microtiter strips each with 8 break-off reagent wells coated with LKM-1. In the first reaction step, diluted patient samples are incubated in the wells. In the case of positive samples, specific IgG antibodies (also IgA and IgM) will bind to the antigens. To detect the bound antibodies, a second incubation is carried out using an enzyme-labelled anti-human IgG (enzyme conjugate) catalysing a colour reaction.

Contents of the test kit:

Component	Colour	Format	Symbol
1. Microplate wells coated with antigens 12 microplate strips each containing 8 individual break-off wells in a frame, ready for use	---	12 x 8	
2. Calibrator 1 200 RU/ml (IgG, human), ready for use	dark red	1 x 2.0 ml	
3. Calibrator 2 20 RU/ml (IgG, human), ready for use	red	1 x 2.0 ml	
4. Calibrator 3 2 RU/ml (IgG, human), ready for use	light red	1 x 2.0 ml	
5. Positive control (IgG, human), ready for use	blue	1 x 2.0 ml	
6. Negative control (IgG, human), ready for use	green	1 x 2.0 ml	
7. Enzyme conjugate peroxidase-labelled anti-human IgG (rabbit), ready for use	green	1 x 12 ml	
8. Sample buffer ready for use	light blue	1 x 100 ml	
9. Wash buffer 10x concentrate	colourless	1 x 100 ml	
10. Chromogen/substrate solution TMB/H ₂ O ₂ , ready for use	colourless	1 x 12 ml	
11. Stop solution 0.5 M sulphuric acid, ready for use	colourless	1 x 12 ml	
12. Test instruction	---	1 booklet	
13. Quality control certificate	---	1 protocol	
 Lot description			 Storage temperature
 In vitro diagnostic medical device			 Unopened usable until

Preparation and stability of the reagents

Note: All reagents must be brought to room temperature (+18°C to +25°C) approx. 30 minutes before use. After first use, the reagents are stable until the indicated expiry date if stored at +2°C to +8°C and protected from contamination, unless stated otherwise below.

- **Coated wells:** Ready for use. Tear open the resealable protective wrapping of the microplate at the recesses above the grip seam. Do not open until the microplate has reached room temperature to prevent the individual strips from moistening. Immediately replace the remaining wells of a partly used microplate in the protective wrapping and tightly seal with the integrated grip seam (Do not remove the desiccant bag).
Once the protective wrapping has been opened for the first time, the wells coated with antigens can be stored in a dry place and at a temperature between +2°C and +8°C for 4 months.
- **Calibrators and controls:** Ready for use. The reagents must be mixed thoroughly before use.
- **Enzyme conjugate:** Ready for use. The enzyme conjugate must be mixed thoroughly before use.
- **Sample buffer:** Ready for use.
- **Wash buffer:** The wash buffer is a 10x concentrate. If crystallisation occurs in the concentrated buffer, warm it to 37°C and mix well before diluting. The quantity required should be removed from the bottle using a clean pipette and diluted with deionised or distilled water (1 part reagent plus 9 parts distilled water).
For example: For 1 microplate strip, 5 ml concentrate plus 45 ml water.
The working strength wash buffer is stable for 4 weeks when stored at +2°C to +8°C and handled properly.
- **Chromogen/substrate solution:** Ready for use. Close the bottle immediately after use, as the contents are sensitive to light ☀. The chromogen/substrate solution must be clear on use. Do not use the solution if it is blue coloured.
- **Stop solution:** Ready for use.

Storage and stability: The test kit has to be stored at a temperature between +2°C to +8°C. Do not freeze. Unopened, all test kit components are stable until the indicated expiry date.

Waste disposal: Patient samples, calibrators, controls and incubated microplate strips should be handled as infectious waste. All reagents must be disposed of in accordance with local disposal regulations.

Warning: The calibrators and controls of human origin have tested negative for HBsAg, anti-HCV, anti-HIV-1 and anti-HIV-2. Nonetheless, all materials should be treated as being a potential infection hazard and should be handled with care. Some of the reagents contain the agent sodium azide in a non-declarable concentration. Avoid skin contact.

Preparation and stability of the patient samples

Samples: Human serum or EDTA, heparin or citrate plasma.

Stability: Patient samples to be investigated can generally be stored at +2°C to +8°C for up to 14 days. Diluted samples should be incubated within one working day.

Sample dilution: Patient samples are diluted **1:101** in sample buffer. For example: dilute 10 µl of sample in 1.0 ml sample buffer and mix well by vortexing (sample pipettes are not suitable for mixing).

NOTE: Calibrators and controls are prediluted and ready for use, do not dilute them.

Incubation

For **semiquantitative analysis** incubate **calibrator 2** along with the positive and negative controls and patient samples. For **quantitative analysis** incubate **calibrators 1, 2 and 3** along with the positive and negative controls and patient samples.

(Partly) manual test performance

Sample incubation:
(1st step) Transfer 100 µl of the calibrators, positive and negative controls or diluted patient samples into the individual microplate wells according to the pipetting protocol. Incubate for **30 minutes** at room temperature (+18°C to +25°C).

Washing: Manual: Empty the wells and subsequently wash 3 times using 300 µl of working strength wash buffer for each wash.
Automatic: Wash reagent wells 3 times with 450 µl of working strength wash buffer (program setting: e.g. TECAN Columbus Washer "Overflow Modus").

Leave the wash buffer in each well for 30 to 60 seconds per washing cycle, then empty the wells. After washing (manual and automated tests), thoroughly dispose of all liquid from the microplate by tapping it on absorbent paper with the openings facing downwards to remove all residual wash buffer.

Note: Residual liquid (> 10 µl) remaining in the reagent wells after washing can interfere with the substrate and lead to false low extinction values. Insufficient washing (e.g., less than 3 wash cycles, too small wash buffer volumes, or too short residence times) can lead to false high extinction values.

Free positions on the microplate strip should be filled with blank wells of the same plate format as that of the parameter to be investigated.

Conjugate incubation:
(2nd step) Pipette 100 µl of enzyme conjugate (peroxidase-labelled anti-human IgG) into each of the microplate wells. Incubate for **30 minutes** at room temperature (+18°C to +25°C).

Washing: Empty the wells. Wash as described above.

Substrate incubation:
(3rd step) Pipette 100 µl of chromogen/substrate solution into each of the microplate wells. Incubate for **15 minutes** at room temperature (+18°C to +25°C), protect from direct sunlight.

Stopping the reaction: Pipette 100 µl of stop solution into each of the microplate wells in the same order and at the same speed as the chromogen/substrate solution was introduced.

Measurement: **Photometric measurement** of the colour intensity should be made at a **wavelength of 450 nm** and a reference wavelength between 620 nm and 650 nm **within 30 minutes of adding the stop solution**. Prior to measuring, slightly shake the microplate to ensure a homogeneous distribution of the solution.

Test performance using fully automated analysis devices

Sample dilution and test performance are carried out fully automatically using an analysis device. The incubation conditions programmed in the respective software authorised by EUROIMMUN may deviate slightly from the specifications given in the ELISA test instruction. However, these conditions were validated in respect of the combination of the EUROIMMUN Analyzer I, Analyzer I-2P or the DSX from Dynex and this EUROIMMUN ELISA. Validation documents are available on enquiry.

Automated test performance using other fully automated, open system analysis devices is possible. However, the combination should be validated by the user.

Pipetting protocol

	1	2	3	4	5	6	7	8	9	10	11	12
A	C 2	P 6	P 14	P 22			C 1	P 4	P 12	P 20		
B	pos.	P 7	P 15	P 23			C 2	P 5	P 13	P 21		
C	neg.	P 8	P 16	P 24			C 3	P 6	P 14	P 22		
D	P 1	P 9	P 17				pos.	P 7	P 15	P 23		
E	P 2	P 10	P 18				neg.	P 8	P 16	P 24		
F	P 3	P 11	P 19				P 1	P 9	P 17			
G	P 4	P 12	P 20				P 2	P 10	P 18			
H	P 5	P 13	P 21				P 3	P 11	P 19			

The pipetting protocol for microtiter strips 1-4 is an example for the **semiquantitative analysis** of 24 patient samples (P 1 to P 24).

The pipetting protocol for microtiter strips 7-10 is an example for the **quantitative analysis** of 24 patient samples (P 1 to P 24).

The calibrators (C 1 to C 3), the positive (pos.) and negative (neg.) controls, and the patient samples have each been incubated in one well. The reliability of the ELISA test can be improved by duplicate determinations for each sample.

The wells can be broken off individually from the strips. This makes it possible to adjust the number of test substrates used to the number of samples to be examined and minimises reagent wastage.

Both positive and negative controls serve as internal controls for the reliability of the test procedure. They should be assayed with each test run.

Calculation of results


Semiquantitative: Results can be evaluated semiquantitatively by calculating a ratio of the extinction value of the control or patient sample over the extinction value of calibrator 2. Calculate the ratio according to the following formula:

$$\frac{\text{Extinction of the control or patient sample}}{\text{Extinction of calibrator 2}} = \text{Ratio}$$

EUROIMMUN recommends interpreting results as follows:

Ratio <1.0:	negative
Ratio ≥1.0:	positive

Quantitative: The standard curve from which the concentration of antibodies in the serum samples can be taken is obtained by point-to-point plotting of the extinction values measured for the 3 calibration sera against the corresponding units (linear/linear). Use "point-to-point" plotting for calculation of the standard curve by computer. The following plot is an example of a typical calibration curve. Please do not use this curve for the determination of antibody concentrations in patient samples.

If the extinction for a patient sample lies above the value of calibrator 1 (200 RU/ml), the result should be reported as ">200 RU/ml". It is recommended that the sample be re-tested at a dilution of e.g. 1:400. The result in RU/ml read from the calibration curve for this sample must then be multiplied by a factor of 4.

The upper limit of the normal range (**cut-off**) recommended by EUROIMMUN is 20 relative units (RU)/ml. EUROIMMUN recommends interpreting results as follows:

<20 RU/ml:	negative
≥20 RU/ml:	positive

For duplicate determinations the mean of the two values should be taken. If the two values deviate substantially from one another, EUROIMMUN recommends to retest the samples.

For diagnosis, the clinical picture of the patient always needs to be taken into account along with the serological findings.

Test characteristics

Calibration: As no international reference serum exists for antibodies against LKM-1, the calibration is performed in relative units (RU).

For every group of tests performed, the extinction values of the calibrators and the relative units and/or ratio determined for the positive and negative controls must lie within the limits stated for the relevant test kit lot. A quality control certificate containing these reference values is included. If the values specified for the controls are not achieved, the test results may be inaccurate and the test should be repeated.

The binding activity of the antibodies and the activity of the enzyme used are temperature-dependent. It is therefore recommended using a thermostat in all three incubation steps. The higher the room temperature during the incubation steps, the greater will be the extinction values. Corresponding variations apply also to the incubation times. However, the calibrators are subject to the same influences, with the result that such variations will be largely compensated in the calculation of the result.

Antigen: The reagent wells are coated with recombinant cytochrome P450 IID6 which constitutes the specific target antigen for antibodies against LKM1.

Linearity: The linearity of the Anti-LKM-1 ELISA (IgG) was determined by assaying 4 serial dilutions of different patient samples. The coefficient of determination R^2 for all sera was > 0.95 . The Anti-LKM-1 ELISA (IgG) is linear at least in the tested concentration range (2 RU/ml to 194 RU/ml).

Detection limit: The lower detection limit is defined as the mean value of an analyte-free sample plus three times the standard deviation and is the smallest detectable antibody titer. The lower detection limit of the Anti-LKM-1 ELISA (IgG) is 1.4 RU/ml.

Cross reactivity: This ELISA showed no cross reactivity.

Interference: Haemolytic, lipaemic and icteric samples showed no influences on the result up to a concentration of 10 mg/ml for haemoglobin, 20 mg/ml for triglycerides and 0.4 mg/ml for bilirubin in this ELISA.

Reproducibility: The reproducibility of the test was investigated by determining the intra- and inter-assay coefficients of variation (CV) using 3 sera. The intra-assay CVs are based on 20 determinations and the inter-assay CVs on 4 determinations performed in 6 different test runs.

<i>Intra-Assay Variation, n = 20</i>		
Serum	Mean value (RU/ml)	CV (%)
1	55	3.0
2	96	2.7
3	158	2.3

<i>Inter-Assay Variation, n = 4 x 6</i>		
Serum	Mean value (RU/ml)	CV (%)
1	58	3.2
2	96	3.7
3	160	2.5

Specificity and sensitivity: 18 patient samples suffering from autoimmune hepatitis and 489 patient samples from a reference laboratory were investigated with the EUROIMMUN Anti-LKM-1 ELISA (IgG). The EUROIMMUN-IIFT (IgG) was used as a reference method. The ELISA has a specificity of 99.4% and a sensitivity of 100% with reference to the EUROIMMUN IIFT.

Serum panel (n = 507)		IIFT (rat liver/rat kidney)	
		positive	negative
Anti-LKM-1 ELISA	positive	27	3
	negative	0	477

A patient sample which reacted positive in ELISA and negative in IIFT belongs to a patient with characterised AIH.

Reference range: The levels of the anti-LKM-1 antibodies (IgG) were analysed with this EUROIMMUN ELISA in a panel of 200 healthy blood donors. With a cut-off of 20 RU/ml, 0.5% of the blood donors were anti-LKM-1 positive.

Clinical significance

In Western Europe the incidence of AIH is 1.9 cases per 100,000 inhabitants per year. Untreated, AIH soon develops into liver cirrhosis. However, if low-dose immunosuppressive therapy is started early enough and continued lifelong, patients have a normal life expectancy.

Circulating autoantibodies have come to play a significant role in the diagnosis of AIH. Antibodies against the following antigens are associated with AIH: soluble liver antigen/liver-pancreas antigen (SLA/LP), cell nuclei (ANA), nDNA, smooth muscles (SMA, the most important target antigen being F actin), liver-kidney microsomes (LKM-1, target antigen: cytochrome P450 IID6), liver cytosolic antigen type 1 (LC-1, target antigen: formiminotransferase cyclo-deaminase) and granulocytes (pANCA). Anti-mito-chondrial antibodies (AMA) are also investigated in this context to exclude the possibility of primary biliary cirrhosis (PBC). AIH is sometimes classified according to the antibody status, i.e., subtype I (ANA, SMA), subtype II (antibodies against LKM-1 and LC-1), or subtype III (antibodies against SLA/LP). However, this classification is probably neither clinically nor therapeutically or prognostically relevant, since 10 to 20% of patients with PBC develop secondary AIH (overlap). In these cases the same autoantibodies as in AIH are frequently detected.

Autoantibodies against LKM-1 (LKM-1, antigen: cytochrome P450 IID6) are present in 1% of adults with AIH. In children they are more common. Antibodies against LKM-1 are also found in 1 to 2% of patients with hepatitis C-positive serology.

The highest diagnostic accuracy currently available for AIH is probably provided by the various EUROIMMUN enzyme immunoassays that detect autoantibodies against SLA/LP. Although SLA/LP autoantibodies have a prevalence of only 10 to 30% in AIH patients, the predictive value is nearly 100%. Every positive anti-SLA/LP result essentially indicates AIH (provided the relevant clinical symptoms are also present).

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HP Ag

**Enzyme Immunoassay for the
qualitative/quantitative determination
of Helicobacter pylori Antigen
in human stools**

- for “in vitro” diagnostic use only -

DIA.PRO

**Diagnostic Bioprobes Srl
Via G. Carducci n° 27
20099 Sesto San Giovanni
(Milano) - Italy**

Phone +39 02 27007161

Fax +39 02 44386771

e-mail: info@diapro.it

REF HPAG.CE.96
96 Tests

HP Ag

A. INTENDED USE

Enzyme Immunoassay (ELISA) for the one-step qualitative/quantitative determination of *Helicobacter pylori* Antigen (HP Ag) in human stools. The kit may be used for the follow-up of HP-infected patients and their pharmacological treatment. For "in vitro" diagnostic use only.

B. INTRODUCTION

Helicobacter pylori (Hp) is a Gram negative bacterium, firstly isolated in gastric mucosa by Marshall and Warren in 1983.

This bacterium is widely diffused in men, without limitations of sex and age; it has been found that infections can be transmitted directly by contact with contaminated biological fluids (saliva, stool, body secretions) and also from contaminated food and beverages.

H.pylori, and in particular some pathogenic strains (CagA +), is the etiological agent responsible of most of active infections and lesions of the gastric mucosa in man.

H.pylori infection moreover acts as cofactor in the development of tumoral pathologies of the gastric apparatus and it is suspected to be associated to some inflammatory pathologies of the genital female apparatus, evolving toward neoplastic transformation.

At the present time, the identification of *Helicobacter pylori* is mostly made with invasive histochemical techniques, with the determination of its urease activity on a isotopic substrate (breath test and mass analysis), with time-consuming bacteriological culture systems and with expensive molecular biology techniques (PCR).

ELISA for HP Ag have been only recently introduced as a specific, fast, non invasive (analysis of stools) and cheaper method of detection.

C. PRINCIPLE OF THE TEST

Stools from patients are used as a source of sample for the determination of HP antigen.

Microplates are coated with a cocktail of affinity purified mouse monoclonal antibodies directed to the most specific *Helicobacter pylori* antigens.

In the 1st incubation, the solid phase is treated with the sample, previously extracted from stools, and simultaneously with a mixture of monoclonal antibodies to Hp, conjugated with peroxidase (HRP).

After washing out all the other components of the sample, in the 2nd incubation the bound enzyme specifically present on the solid phase generates an optical signal that is proportional to the amount of *H.pylori* antigens present in the sample.

D. COMPONENTS

Code HPAG.CE.96 contains reagents to perform 96 tests.

1. Microplate MICROPLATE

8x12 microwell strips coated with anti HP Ag specific affinity purified mouse monoclonal antibodies and sealed into a bag with desiccant. Allow the microplate to reach room temperature before opening; reseal unused strips in the bag with desiccant and store at 4°C.

2. Calibration Set: CAL ...

n° 1 set of 4 vials - Lyophilized calibrators. To be dissolved with EIA grade water. When dissolved, Calibrators have the following concentrations: **0 - 0.1 - 0.5 - 1.0 ug/ml HP Ag**

They contain fetal bovine serum, inactivated HP Ag, 10 mM phosphate buffer pH 7.4+/-0.1, 0.02% gentamicine sulphate and 0.045% ProClin 300 as preservatives.

Important Note: Calibrators when dissolved are not stable. Proceed as described in the proper section for storage.

3. Wash buffer concentrate WASHBUF 20X

1x60ml/bottle - 20x concentrated solution. Once diluted, the wash solution contains 10 mM phosphate buffer pH 7.0+/-0.2 and 0.05% Tween 20 and 0.045% ProClin 300

4. Enzyme Conjugate CONJ

1x16ml/vial - Ready to use component. It contains Horseradish Peroxidase (HRP) labeled mouse monoclonal antibodies to HP Ag, 10 mM Tris buffer pH 6.8+/-0.1, 2% BSA, 0.045% ProClin 300 and 0.02% gentamicine sulphate as preservatives. The Enzyme Conjugate is color coded red.

5. Chromogen/Substrate SUBS TMB

1x25ml/vial - It contains a 50 mM citrate-phosphate buffered solution at pH 3.5-3.8, 0.03% tetra-methyl-benzidine (TMB) and 0.02% hydrogen peroxide (H₂O₂).

Note: To be stored protected from light as sensitive to strong illumination.

6. Specimen Diluent: DILSPE

2x60ml/vial - Buffered solution for the extraction of HP Ag from the specimen and preparation of the sample. It contains 10 mM Tris-HCl buffer pH 7.4+/-0.1, 2% BSA, 0.045% ProClin 300 and 0.02% gentamicine sulphate as preservatives. The component is color coded blue.

7. Sulphuric Acid H₂SO₄ 0.3 M

1x15ml/vial - It contains 0.3 M H₂SO₄ solution. Attention: Irritant (H315, H319; P280, P302+P352, P332+P313, P305+P351+P338, P337+P313, P362+P363).

8. Plate sealing foils: n° 2

9. Package insert: n° 1

Upon request:

HP Ag Extraction kit n° 1

The kit contains all what is necessary to prepare n° 50 samples extracted from stools collected by patients.

E. MATERIALS REQUIRED BUT NOT PROVIDED

1. Calibrated variable volume Micropipettes ranging 1000 ul and 200 ul; disposable plastic tips.
2. EIA grade water (double distilled or deionised, charcoal treated to remove oxidizing chemicals used as disinfectants).
3. Timer with 60 minute range or higher.
4. Absorbent paper tissues.
5. Calibrated ELISA microplate thermostatic incubator (dry or wet), set at +37°C.
6. Calibrated ELISA microwell reader with 450nm (reading) and with 620-630nm (blanking) filters.
7. Calibrated ELISA microplate washer.
8. Vortex or similar mixing tools.
9. Disposable plastic micro-spoon stools collection container (available upon request from Dia.Pro s.r.l.)

F. WARNINGS AND PRECAUTIONS

1. The kit has to be used by skilled and properly trained technical personnel only, under the supervision of a medical doctor responsible of the laboratory.
2. All the personnel involved in performing the assay have to wear protective laboratory clothes, talc-free gloves and glasses. The use of any sharp (needles) or cutting (blades) devices should be avoided. All the personnel involved should be trained in biosafety procedures, as recommended by the Center for

Disease Control, Atlanta, U.S. and reported in the National Institute of Health's publication: "Biosafety in Microbiological and Biomedical Laboratories", ed. 1984.

3. The laboratory environment should be controlled so as to avoid contaminants such as dust or air-borne microbial agents, when opening kit vials and microplates and when performing the test. Protect the Chromogen/Substrate from strong light and avoid vibration of the bench surface where the test is undertaken.

4. Upon receipt, store the kit at 2..8°C into a temperature controlled refrigerator or cold room.

5. Do not interchange components between different lots of the kits. It is recommended that components between two kits of the same lot should not be interchanged.

6. Check that the liquid components of the kit are clear and do not contain visible heavy particles or aggregates. If not, advise the laboratory supervisor to initiate the necessary procedures for kit replacement.

7. Avoid cross-contamination between samples by using disposable tips and changing them after each sample. Do not reuse disposable tips.

8. Avoid cross-contamination between kit components by using disposable tips and changing them between the use of each one. Do not reuse disposable tips.

9. Do not use the kit after the expiration date stated on the external container and internal (vials) labels.

10. Treat all specimens as potentially infective, according to national regulations and laws concerning biological sample handling and wasting.

11. The use of disposable plastic-ware is recommended in the preparation of the liquid components or in transferring components into automated workstations, in order to avoid cross contamination.

12. Wastes produced during the use of the kit have to be discarded in compliance with national directives and laws concerning laboratory waste of chemical and biological substances. In particular, liquid waste generated from the washing procedure, from residuals of controls and from samples has to be treated as potentially infective material and inactivated before waste. Suggested procedures of inactivation are treatment with a 10% final concentration of household bleach for 16-18 hrs or heat inactivation by autoclave at 121°C for 20 min..

13. Accidental spills from samples and operations have to be adsorbed with paper tissues soaked with household bleach and then with water. Tissues should then be discarded in proper containers designated for laboratory/hospital waste.

14. The Sulphuric Acid is an irritant. In case of spills, wash the surface with plenty of water

15. Other waste materials generated from the use of the kit (example: tips used for samples and controls, tools for the extraction of the sample from specimens, used microplates, etc.) should be handled as potentially infective and disposed according to national directives and laws concerning laboratory wastes.

G. SPECIMEN: COLLECTION, PREPARATION AND WARNINGS

1) It is recommended to collect fresh stools in the morning with the plastic collector provided on request together with the kit. Alternatively a conical bottomed disposable tube, provided by the laboratory to the patient, may be used.

2) The patient submitted to the test should not be under antibiotic or anti-bacterial treatments as this pharmaceutical therapy is known to affect H.pylori up to a certain extent, depending on the antibiotic used, giving origin to false interpretation.

3) The patient has to be asked to collect the specimen avoiding any possible contact with urine or water using the plastic spoon present in the stool collector and taking just the amount of specimen necessary to fill up the cavity of the spoon.

4) The patient is asked to deliver the specimen the same day to the laboratory. From the time of collection, the specimen can be stored in the laboratory up to 24 hr at 2..8°C or kept frozen at -20°C for longer time.

5) Specimens, and then samples derived from them, have to be clearly identified with codes or names in order to avoid misinterpretation of results. Bar code labeling and electronic reading is recommended when the number of samples on testing is pretty high.

Important Note: Degradation of HP antigen heavily occurs in stools after 24 hrs generating false negative results, even if the specimen is stored at 2..8°C.

The next following operations are described and represented in figures in the Instructions for Use of the Stool Extraction Kit provided together with the kit.

Operate according to the following instructions:

1) Open the stool collection device and introduce the extraction brush deeply into the specimen. Rotate the brush 2-4 times in order to collect the right amount of biological material (about 0.2 gr).

2) Transfer the brush carefully into the test tube supplied in the kit and then add 1 ml Specimen Diluent. Keeping the brush inside the tube, mix vigorously on vortex for 1 min +/-10% in order to dissolve H.pylori into solution.

3) Discard the brush and insert the filtering piston, supplied with the kit, into the tube. Push gently the piston down into the tube in order to collect not more than 150-200 ul of the liquid phase of the suspension, volume enough to carry out the test.

Important Notes:

a) Be careful not to apply a too strong manual pressure on the piston. The piston could break the tube and spills could be generated. If this should happen, use a paper towel soaked with an hospital disinfectant to clean up the contaminated surfaces.

b) Avoid any addition of preservatives to samples, especially sodium azide as this chemical would affect the enzymatic activity of the conjugate, generating false negative results.

H. PREPARATION OF COMPONENTS AND WARNINGS

A study conducted on an opened kit has not pointed out any relevant loss of activity up to 6 re-uses of the device and up to 3 months.

Microplates:

Allow the microplate to reach room temperature (about 1 hr) before opening the container. Check that the desiccant has not turned dark green, indicating a defect in storing.

In this case, call Dia.Pro's customer service.

Unused strips have to be placed back inside the aluminum pouch, with the desiccant supplied, firmly zipped and stored at +2..8°C.

Important Note: After first opening, remaining strips are stable until the humidity indicator inside the desiccant bag turns from yellow to green.

Calibration Set:

Add the volume of ELISA grade water reported in the label to the lyophilized powder of each Calibrator. Let fully dissolve the content and then gently mix on vortex.

Important Note: When dissolved, Calibrators are not stable. Store Calibrators frozen in aliquots at -20°C, carefully labeled with the content of HP Ag present in each of them.

Wash buffer concentrate:

The concentrated solution has to be diluted 20x with ELISA grade water and mixed gently end-over-end before use. During

preparation avoid foaming as the presence of bubbles could impact on the efficiency of the washing cycles.

Important Note: Once diluted, the wash solution is stable for 1 week at +2..8° C.

Enzyme Conjugate:

Ready to use. Mix well on vortex before use.

Chromogen/Substrate:

Ready to use. Mix well on vortex before use.

Avoid contamination of the liquid with oxidizing chemicals, air-driven dust or microbes. Do not expose to strong light, oxidizing agents and metallic surfaces.

If this component has to be transferred use only plastic, and if possible, sterile disposable container.

Specimen Diluent:

Ready to use. Mix well on vortex before use.

Sulphuric Acid:

Ready to use. Mix well on vortex before use.

Attention: Irritant (H315, H319; P280, P302+P352, P332+P313, P305+P351+P338, P337+P313, P362+P363).

Legenda:

Warning H statements:

H315 – Causes skin irritation.

H319 – Causes serious eye irritation.

Precautionary P statements:

P280 – Wear protective gloves/protective clothing/eye protection/face protection.

P302 + P352 – IF ON SKIN: Wash with plenty of soap and water.

P332 + P313 – If skin irritation occurs: Get medical advice/attention.

P305 + P351 + P338 – IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P337 + P313 – If eye irritation persists: Get medical advice/attention.

P362 + P363 – Take off contaminated clothing and wash it before reuse.

I. INSTRUMENTS AND TOOLS USED IN COMBINATION WITH THE KIT

1. Micropipettes have to be calibrated to deliver the correct volume required by the assay and must be submitted to regular decontamination (70% ethanol, 10% solution of bleach, hospital grade disinfectants) of those parts that could accidentally come in contact with the sample or the components of the kit. They should also be regularly maintained in order to show a precision of 1% and a trueness of $\pm 2\%$.
2. The ELISA incubator has to be set at +37°C (tolerance of $\pm 0.5^\circ\text{C}$) and regularly checked to ensure the correct temperature is maintained. Both dry incubators and water baths are suitable for the incubations, provided that the instrument is validated for the incubation of ELISA tests.
3. The ELISA washer is extremely important to the overall performances of the assay. The washer must be carefully validated in advance, checked for the delivery of the right dispensation volume and regularly submitted to maintenance according to the manufacturer's instructions for use. In particular the washer, at the end of the daily workload, has to be extensively cleaned out of salts with deionized water. Before use, the washer has to be extensively primed with the diluted Washing Solution. The instrument weekly has to be submitted to decontamination according to its manual (NaOH 0.1 M decontamination suggested).

5 washing cycles (aspiration + dispensation of 350ul/well of washing solution + 20 sec soaking = 1 cycle) are sufficient to ensure the assay with the declared performances. If soaking is not possible add one more cycle of washing.

An incorrect washing cycle or salt-blocked needles are the major cause of false positive reactions. The washer has to be correctly optimized using the kit controls/calibrator and reference panels, before using the kit for routine laboratory tests.

Important Note: Due to the nature of the sample used and the possible presence of particles in the sample, be careful to control that the needles of the washer do not get blocked by the presence of stool bodies.

4. Incubation times have a tolerance of $\pm 5\%$.
5. The ELISA microplate reader has to be equipped with a reading filter of 450nm and with a second filter of 620-630nm, mandatory for blanking purposes. Its standard performances should be (a) bandwidth ≤ 10 nm; (b) absorbance range from 0 to ≥ 2.0 ; (c) linearity to ≥ 2.0 ; repeatability $\geq 1\%$. Blanking is carried out on the well identified in the section "Internal Quality Control". The optical system of the reader has to be calibrated regularly to ensure that the correct optical density is measured. It should be regularly maintained according to the manufacturer 's instructions.
6. When using an ELISA automated workstation, all critical steps (dispensation, incubation, washing, reading, shaking, data handling) have to be carefully set, calibrated, controlled and regularly serviced in order to match the values reported in the sections "Internal Quality Control". The assay protocol has to be installed in the operating system of the unit and validated as for the washer and the reader. In addition, the liquid handling part of the station (dispensation and washing) has to be validated and correctly set. Particular attention must be paid to avoid carry over by the needles used for dispensing samples and for washing. This must be studied and controlled to minimize the possibility of contamination of adjacent wells due to strongly reactive samples, leading to false positive results. The use of ELISA automated work stations is recommended when the number of samples to be tested exceed 20-30 units per run.

Important Note: Due to the nature of the sample used and the possible presence of particles in the sample, be careful to control that the needles of the workstation do not get blocked by the presence of stool bodies. We strongly suggest to use disposable sample tips in order to avoid any block or damage of fix probes.

7. Dia.Pro's customer service offers support to the user in the setting and checking of instruments used in combination with the kit, in order to assure full compliance with the requirements described. Support is also provided for the installation of new instruments to be used with the kit.
8. Upon request, Dia.Pro srl offers a sample preparation device able to produce a particle free sample showing excellent performances in the assay. Please inquire.

L. PRE ASSAY CONTROLS AND OPERATIONS

1. Prepare the sample from stools as described in section G and represented in the Instructions for Use of the HP Ag Extraction Kit.
2. Check the expiration date of the kit printed on the external label of the kit box. Do not use if expired.
3. Check that the liquid components are not contaminated by naked-eye visible particles or aggregates. Check that the Chromogen/Substrate is colorless or pale blue by aspirating a small volume of it with a sterile transparent plastic pipette. Check that no breakage occurred in transportation and no spillage of liquid is present inside the box. Check that the aluminum pouch, containing the microplate, is not punctured or damaged.

- Dilute all the content of the 20x concentrated Wash Solution as described above.
- Dissolve the Calibrator Set as described above.
- Allow all the other components to reach room temperature (about 1 hr) and then mix as described.
- Set the ELISA incubator at +37°C and prepare the ELISA washer by priming with the diluted washing solution, according to the manufacturers instructions. Set the right number of washing cycles reported in the specific-section.
- Check that the ELISA reader has been turned on at least 20 minutes before reading.
- If using an automated workstation, turn it on, check settings and be sure to use the right assay protocol.
- Check that the micropipettes are set to the required volume.
- Check that all the other equipment is available and ready to use.
- In case of problems, do not proceed further with the test and advise the supervisor.

M. ASSAY PROCEDURE

The assay has to be carried out according to what reported below, taking care to maintain the same incubation time for all the samples in testing.

Two procedures are available: a quantitative method able to provide a quantification of HP Ag in the specimen and a qualitative method.

A. Quantitative Assay

- Place the required number of strips in the plastic holder and carefully identify the wells for calibrators and samples. Leave A1+B1 wells empty for blanking purposes.
- Pipette 100 µl Calibrators in duplicate into the calibration wells (see the example of dispensation reported below).
- With the Pasteur pipette supplied aspirate the liquid filtered up into the inner chamber of the piston and dispense 2 drops (about 100 µl) of sample into the sample well. Check for the presence of samples in wells by naked eye (there is a marked color difference between empty and full wells) or by reading at 450/620nm. (samples show OD values higher than 0.100).
- Dispense then 100 µl Enzymatic Conjugate in all wells, except for A1+B1, used for blanking operations.

Important note: Be careful not to touch the inner surface of the well with the pipette tip when the conjugate is dispensed. Contamination might occur.

- Following addition of the conjugate, check that the color of the samples has turned from brown to pale reddish and incubate the microplate for **120 min at +37°C**.

Important notes: Strips have to be sealed with the adhesive sealing foil only when the test is performed manually. Do not cover strips when using ELISA workstations.

- When the first incubation is over, wash the microwells as previously described (section I.3)
- Pipette 200 µl Chromogen/Substrate into all the wells, A1+B1 included. Incubate the microplate protected from light at **room temperature (18-24°C) for 20 min**.

Important note: Do not expose to strong direct light as a high background might be generated.

- Pipette 100 µl Sulphuric Acid into all the wells to stop the enzymatic reaction, using the same pipetting sequence as in step 6.

- Measure the color intensity of the solution in each well, as described in section I.5 using a 450nm filter (reading) and a 620-630nm filter (background subtraction, mandatory), blanking the instrument on A1 or B1 or both.

An example of dispensation scheme is reported below:

		Microplate					
		1	2	3	4	5	6
A	BLK	CAL4					
B	BLK	CAL4					
C	CAL1	S1					
D	CAL1	S2					
E	CAL2	S3					
F	CAL2	S4					
G	CAL3	S5					
H	CAL3	S6					

Legenda: BLK = Blank CAL = Calibrator S = Sample

B. Qualitative Assay

- Place the required number of strips in the plastic holder and carefully identify the wells for calibrators and samples. Leave A1well empty for blanking purposes.
- Pipette 100 µl Calibrator 1 in duplicate, 100 µl Calibrator 2 in duplicate, 100 µl Calibrator 4 in single and then 100 µl samples. Check for the presence of samples in wells as reported before.
- Dispense 100 µl Enzymatic Conjugate in all wells, except for A1, used for blanking operations.

Important note: Be careful not to touch the inner surface of the well with the pipette tip when the conjugate is dispensed. Contamination might occur.

- Following addition of the conjugate, check that the color of the samples has turned from brown to pale reddish and then incubate the microplate for **120 min at +37°C**.

Important notes: Strips have to be sealed with the adhesive sealing foil only when the test is performed manually. Do not cover strips when using ELISA workstations.

- When the first incubation is over, wash the microwells as previously described (section I.3)
- Pipette 200 µl Chromogen/Substrate into all the wells, A1 included. Incubate the microplate protected from light at **room temperature (18-24°C) for 20 min**.

Important note: Do not expose to strong direct light as a high background might be generated.

- Pipette 100 µl Sulphuric Acid into all the wells to stop the enzymatic reaction, using the same pipetting sequence as in step 6.
- Measure the color intensity of the solution in each well, as described in section I.5 using a 450nm filter (reading) and a 620-630nm filter (background subtraction, mandatory), blanking the instrument on A1.

An example of dispensation scheme is reported below:

Microplate						
	1	2	3	4	5	6
A	BLK	S3				
B	CAL1	S4				
C	CAL1	S5				
D	CAL2	S6				
E	CAL2	S7				
F	CAL4	S8				
G	S1	S9				
H	S2	S10				

Legenda: BLK = Blank CAL = Calibrator S = Sample

Important notes:

1. Ensure that no fingerprints or dust are present on the external bottom of the microwell before reading. They could generate false positive results on reading.
2. Reading should ideally be performed immediately after the addition of the acid solution but definitely no longer than 20 minutes afterwards. Some self-oxidation of the chromogen can occur leading to a higher background.

N. ASSAY SCHEME

Operations	Procedure
Calibrators&samples	100 ul
Enzyme Conjugate	100 ul
1st incubation	120 min
Temperature	+37°C
Washing steps	n° 5 cycles with 20" of soaking OR n° 6 cycles without soaking
Chromogen/Substrate	200ul
2nd incubation	20 min
Temperature	room
Sulphuric Acid	100 ul
Reading OD	450nm / 620-630nm

O. INTERNAL QUALITY CONTROL

A check is performed on the controls/calibrator any time the kit is used in order to verify whether the expected OD450nm/620-630nm or S/Co values have been matched in the analysis. Ensure that:

Parameter	Requirements
Blank well	< 0.100 OD450nm value
CAL 0 ug/ml	< 0.200 mean OD450nm value after blanking
CAL 0.1 ug/ml	OD450nm > OD450nm CAL 0 ug/ml + 0.100
CAL 1 ug/ml	> 1.000 OD450nm value

If the results of the test match the requirements stated above, proceed to the next section. If they do not, do not proceed any further and perform the following checks:

Problem	Check
Blank well > 0.100 OD450nm	1. that the Chromogen/Substrate solution has not become contaminated during the assay
CAL 0 ug/ml > 0.200 OD450nm after blanking	1. that the washing procedure and the washer settings are as validated in the pre qualification study; 2. that the proper washing solution has been used and the washer has been primed with it before use; 3. that no mistake has been done in the assay procedure (dispensation of positive calibrators instead of the negative one); 4. that no contamination of the calibrator or of the wells where the calibrator was dispensed has occurred due to spills of positive samples or of the enzyme conjugate; 5. that micropipettes have not become contaminated with positive samples or with the enzyme conjugate; 6. that the washer needles are not blocked or partially obstructed.
CAL 0.1 ug/ml OD450nm < Cal 0 ug/ml + 0.100	1. that the procedure has been correctly performed; 2. that no mistake has occurred during its distribution (ex.: dispensation of negative calibrator instead) 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the calibrator has occurred.
CAL 1 ug/ml < 1.000 OD450nm	1. that the procedure has been correctly performed; 2. that no mistake has occurred during the distribution of the calibrator (dispensation of negative calibrator instead). 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the calibrator has occurred.

If any of the above problems have occurred, report the problem to the supervisor for further actions.

Important note:

The analysis must be done proceeding as the reading step described in the section M, points A8 and B8.

P. CALCULATION OF RESULTS

Quantitative Assay:

Calculate the mean OD450nm/620-630nm value of the calibrators. Then draw a calibration curve possibly using a 4 parameters fitting curve system. Then calculate on the curve the concentration of HP antigen in the sample.

Qualitative Assay:

The test results are calculated by means of a cut-off value determined from the OD450nm/ 620-630nm value of the CAL 0 ug/ml (CAL 0) and the OD450nm/620-630nm of the CAL 0.1 ug/ml (CAL 0.1) with the following formula:

$$\text{Cut-Off} = (\text{CAL 0} + \text{CAL 0.1}) / 2$$

Important note: When the calculation of results is performed by the operating system of an ELISA automated work station, ensure that the proper formulation is used to calculate the cut-off value and generate the correct interpretation of results.

Q. INTERPRETATION OF RESULTS

In the quantitative assay, samples showing a concentration of H.pylori antigen higher than 0.05 ug/ml are considered positive. For the qualitative assay, test results are interpreted as a ratio of the sample OD450nm/620-630nm (S) and the Cut-Off value (Co), mathematically S/Co, according to the following table:

S/Co	Interpretation
< 1.0	Negative
1.0 – 1.1	Equivocal
> 1.1	Positive

A negative result indicates that the patient is not infected by H.pylori.

Any patient showing an equivocal result should be retested on a second sample.

A positive result is indicative of HP infection and therefore the patient should be treated accordingly.

Important notes:

1. Interpretation of results should be done under the supervision of the laboratory supervisor to reduce the risk of judgment errors and misinterpretations.
2. Any positive result should be confirmed first by repeating the test and then, if still positive, by an alternative method before a diagnosis of HP infection is confirmed.
3. When test results are transmitted from the laboratory to another department, attention must be paid to avoid erroneous data transfer.
4. Diagnosis of HP infection has to be taken and released to the patient by a suitably qualified medical doctor. This should be done taking also into account other diagnostic evidences of infection.

An example of qualitative method is reported below. (data obtained proceeding as the the reading step described in the section M, point B8):


Note: The following data must not be used instead of real figures obtained by the user.

Cal 0 ug/ml: 0.040 – 0.060 OD450nm
Mean Value: 0.050 OD450nm
Lower than 0.200 – Accepted

Cal 0.1 ug/ml: 0.210-0.230 OD450nm
Mean Value: 0.220 OD450nm
Higher than Cal 0 ug/ml + 0.100 – Accepted
Cut-Off = (CAL 0 + CAL 0.1) / 2 = 0.135
Calibrator 1 ug/ml: 2.000 OD450nm
OD450nm higher than 1.000 – Accepted

Sample 1: 0.028 OD450nm
Sample 2: 1.690 OD450nm
Sample 1 S/Co < 1.0 = negative
Sample 2 S/Co > 1.1 = positive

An Example of Calibration curve is reported below:

R. PERFORMANCE CHARACTERISTICS

Evaluation of Performances was conducted by testing negative and positive samples in an external clinical site and in the own laboratories as well.

1. Limit of detection

The limit of detection of the assay was calculated by examining serial dilution of HP antigen in Sample Diluent.

Results show that the analytical sensitivity of the assay is better than **0.05 ug/ml** when the limit of dilution is considered mean OD450nm/620-630nm CAL 0 ug/ml + 5 SD.

2. Diagnostic performances:

The diagnostic performances of the kit were studied in some external trials against (a) the "breath test", considered by the medical literature the Gold Standard for HP Ag determination and (b) against a commercial CE-marked ELISA, as well.

In the studies were the breath test was used as reference stools were collected from patients the same day of the breath test.

In the study carried out with a commercial ELISA as reference the same samples were extracted with the specific devise supplied in combination with kits.

In the studies conducted in reference to the **breath test** the following mean values were obtained from two centers, one in Italy and one in Spain:

- Sensitivity: 98%
- Specificity: 96%

No confirmation by ELISA was possible on discrepant samples.

In the studies conducted in reference to a CE-marked commercial **ELISA**, similarly based on monoclonal antibodies to H.pylori, the following mean values were obtained:

- Sensitivity: 100%
- Specificity: 93%

The discrepant samples (7% "false positive") turned to be positive in the breath test, revealing a better overall performance of our kit respect to the commercial reference one, when compared on the gold standard method.

Due to the highest specificity shown by the monoclonal antibodies used in the kit no cross-reaction was observed with Campilobacter species.

3. Precision:

The variability shown in the tables did not result in sample misclassification. CV values ranging 4-8%, depending on OD450nm/620-630nm values were observed.

S. LIMITATIONS

False negative results were obtained from samples extracted from specimens stored for more than 1 day at 2..8°C. This finding is explain by the fact that the devise detects "live" H.pylori as confirmed by the results of excellent comparison to the breath test.

False positive results were mostly obtained from samples of liquid stools, difficult to handle and extract.

ASSAY GRAPHICAL SCHEME


Add 100 µl Calibrators, samples and conjugate to the plate and then incubate for 120 min at +37°C


Wash as described in the proper section

Add 200 µl Chromogen/Substrate and incubate 20 min at r.t.

Add 100 µl Sulphuric Acid

Read the plate at 450nm (reading) and at 620-630nm (blanking)

BIBLIOGRAPHY

Fujisawa T. et al. J Clin Lab Anal. 2001; 15(3):154-9
 Braden B. et al. Ann Med. 2001 Mar; 33(2):91-7
 Dominguez-Bella MG. Et al. FEMS Microbiol Lett. 2001 Apr 20;198(1):15-6
 Ding HJ. Et al. J Gastroenterol. 2001 Apr; 36(4):237-41
 Monteiro L. et al. J Microbiol Methods. 2001 Jun; 45(2):89-94
 Kim N. et al. Korean J Intern Med. 2000 Dec; 15(3):187-94
 Monteiro L. et al. Am J Gastroenterol. 2001 Feb; 96(2):353-8
 Shimada T. et al. Nippon Rinsho. 2001 Feb; 59(2):280-5
 Brown LM. Epidemiol Rev. 2000; 22(2):283-97
 Ameriso SF. Et al. Stroke 2001 Feb; 32(2):385-91
 Vaira D. et al. Aliment Pharmacol Ther. 2000 Oct; 14 suppl (3):13-22

All the IVD Products manufactured by the company are under the control of a certified Quality Management System in compliance with ISO 13485 rule. Each lot is submitted to a quality control and released into the market only if conforming with the EC technical specifications and acceptance criteria.

Manufacturer:

Dia.Pro Diagnostic Bioprobes Srl
 Via G. Carducci n° 27 – Sesto San Giovanni (MI) – Italy

Alpha-Fetoprotein (AFP) Test System

Product Code: 1925-300

1.0 INTRODUCTION

Intended Use: The Quantitative Determination of Alpha-Fetoprotein (AFP) Concentration in Human Serum by a Microplate Enzyme Immunoassay, Colorimetric

2.0 SUMMARY AND EXPLANATION OF THE TEST

Alpha-Fetoprotein (AFP) is a glycoprotein with a molecular weight of 70 kDa. AFP is normally produced during fetal development by the hepatocytes, yolk sac and, to a lesser extent, the gastrointestinal tract. Serum concentrations reach a peak level of up to 10 mg/ml at twelve weeks of gestation.¹ This peak level gradually decreases to less than 25 ng/ml after one year of postpartum. Thereafter, the levels reduce further to less than 10 ng/ml.

Elevated levels of AFP are found in patients with primary hepatoma and yolk sac-derived germ tumors. AFP is the most useful marker for the diagnosis and management of hepatocellular carcinoma.² AFP is also elevated in pregnant women. Presence of abnormally high AFP concentrations in pregnant women provides a risk marker for Down syndrome.³

In this method, AFP calibrator, patient specimen or control is first added to a streptavidin coated well. Biotinylated and enzyme labeled monoclonal antibodies (directed against distinct and different epitopes of AFP) are added and the reactants mixed. Reaction between the various AFP antibodies and native AFP forms a sandwich complex that binds with the streptavidin coated to the well. After the completion of the required incubation period, the enzyme-AFP antibody bound conjugate is separated from the unbound enzyme-AFP conjugate by aspiration or decantation. The activity of the enzyme present on the surface of the well is quantitated by reaction with a suitable substrate to produce color.

The employment of several serum references of known alpha-fetoprotein (AFP) levels permits the construction of a dose response curve of activity and concentration. From comparison to the dose response curve, an unknown specimen's activity can be correlated with AFP concentration.


3.0 PRINCIPLE

Immunoenzymometric assay (TYPE 3):

The essential reagents required for an immunoenzymometric assay include high affinity and specificity antibodies (enzyme and immobilized), with different and distinct epitope recognition, **in excess**, and native antigen. In this procedure, the immobilization takes place during the assay at the surface of a microplate well through the interaction of streptavidin coated on the well and exogenously added biotinylated monoclonal anti-AFP antibody.

Upon mixing monoclonal biotinylated antibody, the enzyme-labeled antibody and a serum containing the native antigen, reaction results between the native antigen and the antibodies, without competition or steric hindrance, to form a soluble

sandwich complex. The interaction is illustrated by the following equation:

Simultaneously, the complex is deposited to the well through the high affinity reaction of streptavidin and biotinylated antibody. This interaction is illustrated below:
 $\text{EnzAb} - \text{Ag}_{\text{AFP}} - \text{B}^{\text{tr}}\text{Ab}_{(\text{m})} + \text{Streptavidin}_{\text{C.W.}} \Rightarrow \text{Immobilized complex}$
 $\text{Streptavidin}_{\text{C.W.}}$ = Streptavidin immobilized on well
Immobilized complex = sandwich complex bound to the well

After equilibrium is attained, the antibody-bound fraction is separated from unbound antigen by decantation or aspiration. The enzyme activity in the antibody-bound fraction is directly proportional to the native antigen concentration. By utilizing several different serum references of known antigen values, a dose response curve can be generated from which the antigen concentration of an unknown can be ascertained.

4.0 REAGENTS

Materials Provided:

- AFP Calibrators – 1 ml/vial – Icons A-F**
Six (6) vials of references AFP antigen at levels of 0 (A), 5 (B), 25 (C), 50 (D), 250 (E) and 500 (F) ng/ml. Store at 2-8°C. A preservative has been added.
- Note:** The calibrators, human serum based, were calibrated using a reference preparation, which was assayed against the WHO 1st IRP # 72/225.
- B. AFP Enzyme Reagent – 13ml/vial – Icon**
One (1) vial containing enzyme labeled antibody, biotinylated monoclonal mouse IgG in buffer, dye, and preservative. Store at 2-8°C.
- C. Streptavidin Coated Microplate – 96 wells – Icon**
One 96-well microplate coated with streptavidin and packaged in an aluminum bag with a drying agent. Store at 2-8°C.
- D. Wash Solution Concentrate – 20ml/vial – Icon**
One (1) vial containing a surfactant in buffered saline. A preservative has been added. Store at 2-8°C.
- E. Substrate A – 7ml/vial – Icon S^A**
One (1) vial containing tetramethylbenzidine (TMB) in buffer. Store at 2-8°C.
- F. Substrate B – 7ml/vial – Icon S^B**
One (1) vial containing hydrogen peroxide (H₂O₂) in buffer. Store at 2-8°C.
- G. Stop Solution – 8ml/vial – Icon**
One (1) vial containing a strong acid (1N HCl). Store at 2-8°C.
- H. Product Instructions.**

Note 1: Do not use reagents beyond the kit expiration date.

Note 2: Opened reagents are stable for sixty (60) days when stored at 2-8°C. **Kit and component stability are identified on the label.**

Note 3: Above reagents are for a single 96-well microplate.

4.1 Required But Not Provided:

- Pipette(s) capable of delivering 0.025 & 0.050ml (25 & 50µl) volumes with a precision of better than 1.5%.
- Dispenser(s) for repetitive deliveries of 0.100 & 0.350ml (100 & 350µl) volumes with a precision of better than 1.5%.
- Microplate washers or a squeeze bottle (optional).
- Microplate Reader with 450nm and 620nm wavelength absorbance capability.
- Absorbent Paper for blotting the microplate wells.
- Plastic wrap or microplate cover for incubation steps.
- Vacuum aspirator (optional) for wash steps.
- Timer.
- Quality control materials

5.0 PRECAUTIONS

For In Vitro Diagnostic Use

Not for Internal or External Use in Humans or Animals

All products that contain human serum have been found to be non-reactive for Hepatitis B Surface Antigen, HIV 1&2 and HCV Antibodies by FDA licensed reagents. Since no known test can offer complete assurance that infectious agents are absent, all human serum products should be handled as potentially hazardous and capable of transmitting disease. Good laboratory procedures for handling blood products can be found in the Center for Disease Control / National Institute of Health, "Biosafety in Microbiological and Biomedical Laboratories," 2nd Edition, 1988, HHS Publication No. (CDC) 88-8395.

Safe Disposal of kit components must be according to local regulatory and statutory requirement.

6.0 SPECIMEN COLLECTION AND PREPARATION

The specimens shall be blood, serum in type and the usual precautions in the collection of venipuncture samples should be observed. For accurate comparison to established normal values, a fasting morning serum sample should be obtained. The blood should be collected in a plain redtop venipuncture tube without additives or anti-coagulants. Allow the blood to clot. Centrifuge the specimen to separate the serum from the cells.

In patients receiving therapy with high biotin doses (i.e. >5mg/day), no sample should be taken until at least 8 hours after the last biotin administration, preferably overnight to ensure fasting sample.

Samples may be refrigerated at 2-8°C for a maximum period of five (5) days. If the specimen(s) cannot be assayed within this time, the sample(s) may be stored at temperatures of -20 °C for up to 30 days. Avoid use of contaminated devices. Avoid repetitive freezing and thawing. When assayed in duplicate, 0.050ml (50µl) of the specimen is required.

7.0 QUALITY CONTROL

Each laboratory should assay controls at levels in the low, normal and elevated range for monitoring assay performance. These controls should be treated as unknowns and values determined in every test procedure performed. Quality control charts should be maintained to follow the performance of the supplied reagents. Pertinent statistical methods should be employed to ascertain trends. Significant deviation from established performance can indicate unnoticed change in experimental conditions or degradation of kit reagents. Fresh reagents should be used to determine the reason for the variations.

8.0 REAGENT PREPARATION

- Wash Buffer**
Dilute contents of wash concentrate to 1000ml with distilled or deionized water in a suitable storage container. Store diluted buffer at 2-30°C for up to 60 days.
- Working Substrate Solution** – Stable for one (1) year
Pour the contents of the amber vial labeled Solution 'A' into the clear vial labeled Solution 'B'. Place the yellow cap on the clear vial for easy identification. Mix and label accordingly. Store at 2 - 8°C.

Note 1: Do not use the working substrate if it looks blue.
Note 2: Do not use reagents that are contaminated or have bacteria growth.

9.0 TEST PROCEDURE

Before proceeding with the assay, bring all reagents, serum reference calibrators and controls to room temperature (20-27 °C).
****Test Procedure should be performed by a skilled individual or trained professional****

- Format the microplates' wells for each serum reference calibrator, control and patient specimen to be assayed in duplicate. **Replace any unused microwell strips back into the aluminum bag, seal and store at 2-8°C.**

- Pipette 0.025ml (25µl) of the appropriate serum reference calibrator, control or specimen into the assigned well.
- Add 0.100ml (100µl) of the AFP Enzyme Reagent to each well. **It is very important to dispense all reagents close to the bottom of the coated well.**
- Mix (See Note) the microplate for 20-30 seconds until homogenous.
- Swirl the microplate gently for 20-30 seconds to mix and cover.
- Incubate 60 minutes at room temperature.
- Discard the contents of the microplate by decantation or aspiration. If decanting, tap and blot the plate dry with absorbent paper.
- Add 0.350ml (350µl) of wash buffer (see Reagent Preparation Section), decant (tap and blot) or aspirate. Repeat two (2) additional times for a total of three (3) washes. **An automatic or manual plate washer can be used. Follow the manufacturer's instruction for proper usage. If a squeeze bottle is employed, fill each well by depressing the container (avoiding air bubbles) to dispense the wash. Decant the wash and repeat two (2) additional times.**
- Add 0.100ml (100µl) of working substrate solution to all wells (see Reagent Preparation Section). **Always add reagents in the same order to minimize reaction time differences between wells.**
DO NOT SHAKE THE PLATE AFTER SUBSTRATE ADDITION
- Incubate at room temperature for fifteen (15) minutes.
- Add 0.050ml (50µl) of stop solution to each well and mix gently for 15-20 seconds. **Always add reagents in the same order to minimize reaction time differences between wells.**
- Read the absorbance in each well at 450nm (using a reference wavelength of 620-630nm to minimize well imperfections) in a microplate reader. **The results should be read within thirty (30) minutes of adding the stop solution.**

Note: Cycle (start and stop) mixing (4 cycles) for 5-8 seconds/cycle is more efficient than one continuous (20-30 seconds) cycle to achieve homogeneity. A plate mixer can be used to perform the mixing cycle.

10.0 CALCULATION OF RESULTS


A dose response curve is used to ascertain the concentration of AFP in unknown specimens.

- Record the absorbance obtained from the printout of the microplate reader as outlined in Example 1.
- Plot the absorbance for each duplicate serum reference versus the corresponding AFP concentration in ng/ml on linear graph paper (do not average the duplicates of the serum references before plotting).
- Draw the best-fit curve through the plotted points.
- To determine the concentration of AFP for an unknown, locate the average absorbance of the duplicates for each unknown on the vertical axis of the graph, find the intersecting point on the curve, and read the concentration (in ng/ml) from the horizontal axis of the graph (the duplicates of the unknown may be averaged as indicated). In the following example, the average absorbance (0.420) intersects the dose response curve at 33.2 ng/ml AFP concentration (See Figure 1).

Note: Computer data reduction software designed for ELISA assays may also be used for the data reduction. **If such software is utilized, the validation of the software should be ascertained.**

EXAMPLE 1

Sample I.D.	Well Number	Abs (A)	Mean Abs (B)	Value (ng/ml)
Cal A	A1	0.012	0.011	0
	B1	0.011		
	C1	0.100		
Cal B	D1	0.097	0.098	5
	E1	0.336		
Cal C	F1	0.333	0.335	25
	G1	0.612		
Cal D	H1	0.577	0.594	50
	A2	2.005		
Cal E	B2	1.975	1.990	250
	C2	2.664		
Cal F	D2	2.680	2.672	500
	E2	0.427		
Patient	F2	0.413	0.420	33.2

*The data presented in Example 1 and Figure 1 is for illustration only and **should not** be used in lieu of a dose response curve prepared with each assay.

11.0 QC PARAMETERS

In order for the assay results to be considered valid the following criteria should be met:

1. The absorbance (OD) of calibrator F should be ≥ 1.3 .
2. The absorbance (OD) of calibrator A should be ≤ 0.035 .
3. Four out of six quality control pools should be within the established ranges.

12.0 RISK ANALYSIS

The MSDS and Risk Analysis Form for this product are available on request from Monobind Inc.

12.1 Assay Performance

1. It is important that the time of reaction in each well is held constant to achieve reproducible results.
2. Pipetting of samples should not extend beyond ten (10) minutes to avoid assay drift.
3. Highly lipemic, hemolyzed or grossly contaminated specimen(s) should not be used.
4. If more than one (1) plate is used, it is recommended to repeat the dose response curve.
5. The addition of substrate solution initiates a kinetic reaction, which is terminated by the addition of the stop solution. Therefore, the substrate and stop solution should be added in the same sequence to eliminate any time-deviation during reaction.
6. Plate readers measure vertically. Do not touch the bottom of the wells.
7. Failure to remove adhering solution adequately in the aspiration or decantation wash step(s) may result in poor replication and spurious results.
8. Use components from the same lot. No intermixing of reagents from different batches.
9. Patient specimens with AFP concentrations above 500 ng/ml may be diluted (for example 1/10 or higher) with normal male serum (AFP < 10 ng/ml) and re-assayed. The sample's concentration is obtained by multiplying the result by the dilution factor (x10).
10. Accurate and precise pipetting, as well as following the exact time and temperature requirements prescribed are essential. Any deviation from Monobind's IFU may yield inaccurate results.
11. All applicable national standards, regulations and laws, including, but not limited to, good laboratory procedures, must be strictly followed to ensure compliance and proper device usage.
12. It is important to calibrate all the equipment e.g. Pipettes, Readers, Washers and/or the automated instruments used with this device, and to perform routine preventative maintenance.
13. Risk Analysis- as required by CE Mark IVD Directive 98/79/EC - for this and other devices, made by Monobind, can be requested via email from Monobind@monobind.com.

12.2 Interpretation

1. **Measurements and interpretation of results must be performed by a skilled individual or trained professional.**
2. Laboratory results alone are only one aspect for determining patient care and should not be the sole basis for therapy, particularly if the results conflict with other determinants.

3. The reagents for the test system procedures have been formulated to eliminate maximal interference; however, potential interaction between rare serum specimens and test reagents can cause erroneous results. Heterophilic antibodies often cause these interactions and have been known to be problems for all kinds of immunoassays (Boscato LM, Stuart MC. "Heterophilic antibodies: a problem for all immunoassays" *Clin.Chem.* 1988;34:27-33). For diagnostic purposes, the results from this assay should be used in combination with clinical examination, patient history and all other clinical findings.
4. For valid test results, adequate controls and other parameters must be within the listed ranges and assay requirements.
5. If test kits are altered, such as by mixing parts of different kits, which could produce false test results, or if results are incorrectly interpreted, **Monobind shall have no liability.**
6. If computer controlled data reduction is used to interpret the results of the test, it is imperative that the predicted values for the calibrators fall within 10% of the assigned concentrations.
7. AFP has a low clinical sensitivity and specificity as a tumor marker. Clinically an elevated **AFP value alone is not of diagnostic value as a test for cancer** and should only be used in conjunction with other clinical manifestations (observations) and diagnostic parameters. AFP levels are known to be elevated in a number of benign diseases and conditions including pregnancy and non-malignant liver diseases such as hepatitis and cirrhosis.

13.0 EXPECTED RANGE OF VALUES

Approximately 97-98% of the normal healthy population has AFP levels less than 8.5ng/ml.⁴ In high-risk patients, AFP values between 100-350 ng/ml suggest hepatocellular carcinoma. Concentrations over 350 ng/ml usually indicate the disease.

TABLE 1
Expected Values for the AFP AccuBind® ELISA Test System

Male and Female	<8.5ng/ml (97-98%)
-----------------	--------------------

Values for AFP for a normal, healthy population and pregnant women, during gestation cycle, are given in Table 2. The values depicted below represent limited in house studies in concordance with published literature.^{8,9,10}

TABLE 2
Median Values during Gestation.

Gestation (Week)	AFP (ng/ml)
15	40.14
16	42.91
17	52.34
18	61.50
19	75.57
20	83.31
21	90.46

It is important to keep in mind that establishment of a range of values, which can be expected to be found by a given method for a population of "normal" persons, is dependent upon a multiplicity of factors: the specificity of the method, the population tested and the precision of the method in the hands of the analyst. For these reasons, each laboratory should depend upon the range of expected values established by the Manufacturer only until an in-house range can be determined by the analysts using the method with a population indigenous to the area in which the laboratory is located.

14.0 PERFORMANCE CHARACTERISTICS

14.1 Precision

The within and between assay precision of the AFP AccuBind® ELISA test system were determined by analyses on three different levels of control sera. The number, mean value, standard deviation and coefficient of variation for each of these control sera are presented in Table 3 and Table 4.

TABLE 3
Within Assay Precision (Values in ng/ml)

Sample	N	X	σ	C.V.
Level 1	24	14.71	0.67	4.6
Level 2	24	71.89	2.68	3.7

TABLE 4
Between Assay Precision* (Values in ng/ml)

Sample	N	X	σ	C.V.
Level 1	30	16.20	1.41	8.7
Level 2	30	88.26	7.47	8.5
Level 3	30	188.43	11.92	6.3

*As measured in thirty experiments in duplicate.

14.2 Sensitivity

The AFP AccuBind® ELISA Test System has a sensitivity of 0.01 ng. This is equivalent to a sample containing 0.44 ng/ml AFP concentration. The sensitivity (detection limit) was ascertained by determining the variability of the '0 ng/ml' calibrator and using the 2σ (95% certainty) statistic to calculate the minimum dose.

14.3 Accuracy

The AFP AccuBind® ELISA Test System was compared with a reference method. Biological specimens with concentrations ranging from 1.0 to 41 ng/ml were assayed. The total number of such specimens was 42. The least square regression equation and the correlation coefficient were computed for the AFP procedure in comparison with the reference method. The data obtained is displayed in Table 5.

TABLE 5

Method	Mean	Least Square Regression Analysis	Correlation Coefficient
This Method (Y)	5.27	$y = 0.746(x) + 1.0007$	0.973
Reference (X)	5.72		

Only slight amounts of bias between the AFP AccuBind® ELISA Test System and the reference method are indicated by the closeness of the mean values. The least square regression equation and correlation coefficient indicates excellent method agreement.

14.4 Specificity

No interference was detected with the performance of AFP AccuBind® ELISA Test System upon addition of massive amounts of the following substances to a human serum pool.

SUBSTANCE	Cross Reactivity	Concentration
Acetyl salicylic Acid	ND	100 µg/ml
Amethopterin	ND	100 µg/ml
Ascorbic Acid	ND	100 µg/ml
Atropine	ND	100 µg/ml
Caffeine	ND	100 µg/ml
CEA	ND	10 µg/ml
PSA	ND	1.0 µg/ml
CA-125	ND	10,000 U/ml
hCG	ND	1000 IU/ml
hLH	ND	10 IU/ml
hTSH	ND	100 mIU/ml
hPRL	ND	100 µg/ml

14.5 Linearity & Hook Effect:

Three different lot preparations of the AFP AccuBind® ELISA test system reagents were used to assess the linearity and hook effect. Massive concentrations of AFP (> 100,000 ng/ml) were used for linear dilutions in pooled human patient sera.

The test showed no hook effect up to concentrations of 10,000 ng/ml and a with a dose recovery of 86.1 to 113.6%.

15.0 REFERENCES

1. Wild D, *The Immunoassay Handbook*, Stockton Press, 445 (1994).
2. Henry JB, *Clinical Diagnosis and Management by Laboratory Methods*, WB Saunders Company, 1075 (1996).
3. Wild D, *The Immunoassay Handbook*, Stockton Press p400-02. (1994)
4. Li D, Mallory T, Satomura S, "AFP; a new generation of tumor marker for hepatocellular carcinoma", *Clin Chem Acta*, 313, 15-9 (2001).
5. Mizejewski GJ, 'Alfa-fetoprotein structure and function; relevant to isoforms, epitopes and conformational variants' *Exp Biol Med*, 226, 337-408 (2001).
6. Johnson OJ, Williams R, 'Cirrhosis and etiology of hepatocellular carcinoma', *J Hepatology*, 4, 140-147 (1987).

7. Javadpour N, 'The role of biologic tumor markers in testicular cancer', *Cancer*, 45, 1755-61 (1980).
8. Canick JA, Rish S, 'The accuracy of assigned risks in maternal serum screening', *Prenatal Diagnosis*; 18:413-415 (1998).
9. NIH State-of-the Science Conference Statement on Management of Menopause-Related Symptoms. NIH Consensus State Sci Statements. Mar 21-23; 22(1), 1-38 (2005).
10. Tietz NW, ED: *Clinical Guide to Laboratory Tests* 3rd Ed, Philadelphia, WA Saunders Co (1995).

Effective Date: 2021-Sep-23
MP1925

Rev. 8 DCO: 1509
Product Code: 1925-300

Size	96(A)	192(B)
A)	1ml set	1ml set
B)	1 (13ml)	2 (13ml)
C)	1 plate	2 plates
D)	1 (20ml)	1 (20ml)
E)	1 (7ml)	2 (7ml)
F)	1 (7ml)	2 (7ml)
G)	1 (8ml)	2 (8ml)

For Orders and Inquires, please contact


Monobind Inc.
100 North Pointe Drive
Lake Forest, CA 92630 USA

Tel: +1 949.951.2665 Mail: info@monobind.com
Fax: +1 949.951.3539 Fax: www.monobind.com

Please visit our website to learn more about our products and services.

Glossary of Symbols (EN 980/ISO 15223)

NovaLisa®

Chlamydia pneumoniae IgG

ELISA

CE 0483

Only for in-vitro diagnostic use

English	2
Deutsch	7
Français	12
Italiano	18
Español	24
Português	30
Bibliography / Literatur / Bibliographie / Bibliografia / Bibliografía/ Bibliografia	38
Abbreviations / Abkürzungen / Abréviations / Abbreviazioni / Abreviaciones / Abreviaturas	38
Symbols Key / Symbolschlüssel / Explication des Symboles / Legenda / Símbolos / Tabela de símbolos	39
Summary of Test Procedure / Kurzanleitung Testdurchführung / Résumé de la procedure de test / Schema della procedura / Resumen de la técnica / Resumo do Procedimento de Teste	40

Product Number: CHLG0510 (96 Determinations)

ENGLISH

1. INTRODUCTION

Chlamydiae are non-motile, Gram-negative and obligatory intracellular growing bacteria which form characteristic inclusions within the cytoplasm of parasitized cells. They are easily visible in the light microscope. Three different Chlamydia species pathogenic for humans are known: Chlamydia trachomatis, Chlamydia pneumoniae and Chlamydia psittaci, and one species only pathogenic for animals (C. pecorum). Chlamydia trachomatis is the most prevalent agent of sexually transmitted diseases worldwide (400-500 million cases) and the number of infections is constantly growing. Pregnant women infected with C. trachomatis may transmit these bacteria during childbirth, causing conjunctivitis or pneumonia in newborns. Untreated cases of chlamydial infection can lead to chronic salpingitis, possibly resulting in ectopic pregnancy or infertility. In males, C. trachomatis is a major cause of non-gonococcal urethritis. A severe problem in Chlamydia infections is the frequent asymptomatic insidious course which may result in the initiation of chronic diseases. In many instances primary infections are not recognized and only the sequelae caused by ascended, persisting agents are diagnosed.

Species	Disease	Symptoms (e.g.)	Transmission route
Chlamydia trachomatis	Trachoma Lymphogranuloma venereum (LGV) Urogenital infections	Bilateral conjunctivitis with foreign body sensation, tearing eyes, and discharge of serous secretions from the eye. Genital lesion; swelling of the inguinal and femoral lymph nodes; formation of suppurating buboes. In men: Inflammation of the urethra (urethritis), inflammation of the prostate (prostatitis) and the epididymis (epididymitis). Sterility (infertility). In women: Inflammation of the urethra (urethritis), Bartholin's glands, inflammation of the cervix (cervicitis), uterine mucous membrane (endometritis) and fallopian tubes (salpingitis) and vaginal discharge. Sterility (Infertility)	Sexually transmitted. Can be passed from an infected mother to her baby during childbirth.
C. pneumoniae	Respiratory diseases (e.g. pneumonia, bronchitis) discussed: endocarditis, coronary heart diseases	Fever, cough (usually productive)	Man to Man Aerogenic, droplet infection
C. psittaci	Ornithosis (Psittacosis)	hepatosplenomegaly; unproductive cough, fever and severe headache	Inhalation of feces from infected birds; contact with infected avian viscera

Infection or presence of pathogen may be identified by:

- Microscopy
- PCR
- Serology: Detection of antigens by ELISA
Detection of antibodies by ELISA

2. INTENDED USE

The Chlamydia pneumoniae IgG ELISA is intended for the qualitative determination of IgG class antibodies against Chlamydia pneumoniae in human serum or plasma (citrate, heparin).

3. PRINCIPLE OF THE ASSAY

The qualitative immunoenzymatic determination of specific antibodies is based on the ELISA (Enzyme-linked Immunosorbent Assay) technique.

Microtiterplates are coated with specific antigens to bind corresponding antibodies of the sample. After washing the wells to remove all unbound sample material a horseradish peroxidase (HRP) labelled conjugate is added. This conjugate binds to the captured antibodies. In a second washing step unbound conjugate is removed. The immune complex formed by the bound conjugate is visualized by adding Tetramethylbenzidine (TMB) substrate which gives a blue reaction product.

The intensity of this product is proportional to the amount of specific antibodies in the sample. Sulphuric acid is added to stop the reaction. This produces a yellow endpoint colour. Absorbance at 450/620 nm is read using an ELISA Microtiterplate reader.

4. MATERIALS

4.1. Reagents supplied

- **Microtiterplate:** 12 break-apart 8-well snap-off strips coated with *Chlamydia pneumoniae* antigens; in resealable aluminium foil.
- **IgG Sample Dilution Buffer:** 1 bottle containing 100 mL of phosphate buffer (10 mM) for sample dilution; pH 7.2 ± 0.2; coloured yellow; ready to use; white cap; ≤ 0.0015% (v/v) CMIT/ MIT (3:1).
- **Stop Solution:** 1 bottle containing 15 mL sulphuric acid, 0.2 mol/L; ready to use; red cap.
- **Washing Buffer (20x conc.):** 1 bottle containing 50 mL of a 20-fold concentrated phosphate buffer (0.2 M), pH 7.2 ± 0.2, for washing the wells; white cap.
- **Conjugate:** 1 bottle containing 20 mL of peroxidase labelled antibody to human IgG in phosphate buffer (10 mM); coloured blue; ready to use; black cap.
- **TMB Substrate Solution:** 1 bottle containing 15 mL 3,3',5,5'-tetramethylbenzidine (TMB), < 0.1 %; ready to use; yellow cap.
- **Positive Control:** 1 vial containing 2 mL control; coloured yellow; ready to use; red cap; ≤ 0.02% (v/v) MIT.
- **Cut-off Control:** 1 vial containing 3 mL control; coloured yellow; ready to use; green cap; ≤ 0.02% (v/v) MIT.
- **Negative Control:** 1 vial containing 2 mL control; coloured yellow; ready to use; blue cap; ≤ 0.0015% (v/v) CMIT/ MIT (3:1).

For hazard and precautionary statements see 12.1

For potential hazardous substances please check the safety data sheet.

4.2. Materials supplied

- 1 Cover foil
- 1 Instruction for use (IFU)
- 1 Plate layout

4.3. Materials and Equipment needed

- ELISA Microtiterplate reader, equipped for the measurement of absorbance at 450/620 nm
- Incubator 37 °C
- Manual or automatic equipment for rinsing Microtiterplates
- Pipettes to deliver volumes between 10 and 1000 µL
- Vortex tube mixer
- Distilled water
- Disposable tubes

5. STABILITY AND STORAGE

Store the kit at 2...8 °C. The opened reagents are stable up to the expiry date stated on the label when stored at 2...8 °C.

6. REAGENT PREPARATION

It is very important to bring all reagents and samples to room temperature (20...25 °C) and mix them before starting the test run!

6.1. Microtiterplate

The break-apart snap-off strips are coated with *Chlamydia pneumoniae* antigens. Immediately after removal of the strips, the remaining strips should be resealed in the aluminium foil along with the desiccant supplied and stored at 2...8 °C.

6.2. Washing Buffer (20x conc.)

Dilute Washing Buffer 1 + 19; e. g. 10 mL Washing Buffer + 190 mL distilled water. The diluted buffer is stable for 5 days at room temperature (20...25 °C). In case crystals appear in the concentrate, warm up the solution to 37°C e.g. in a water bath. Mix well before dilution.

6.3. TMB Substrate Solution

The reagent is ready to use and has to be stored at 2...8 °C, away from the light. The solution should be colourless or could have a slight blue tinge. If the substrate turns into blue, it may have become contaminated and should be thrown away.

7. SAMPLE COLLECTION AND PREPARATION

Use human serum or plasma (citrate, heparin) samples with this assay. If the assay is performed within 5 days after sample collection, the samples should be kept at 2...8 °C; otherwise they should be aliquoted and stored deep-frozen (-70...-20 °C). If samples are stored frozen, mix thawed samples well before testing. Avoid repeated freezing and thawing. Heat inactivation of samples is not recommended.

7.1. Sample Dilution

Before assaying, all samples should be diluted 1+100 with IgG Sample Dilution Buffer. Dispense 10 µL sample and 1 mL IgG Sample Dilution Buffer into tubes to obtain a 1+100 dilution and thoroughly mix with a Vortex.

8. ASSAY PROCEDURE

Please read the instruction for use carefully **before** performing the assay. Result reliability depends on strict adherence to the instruction for use as described. The following test procedure is only validated for manual procedure. If performing the test on ELISA automatic systems we recommend increasing the washing steps from three up to five and the volume of Washing Buffer from 300 µL to 350 µL to avoid washing effects. Pay attention to chapter 12. Prior to commencing the assay, the distribution and identification plan for all samples and standards/controls (duplicates recommended) should be carefully established on the plate layout supplied in the kit. Select the required number of microtiter strips or wells and insert them into the holder.

Perform all assay steps in the order given and without any delays.

A clean, disposable tip should be used for dispensing each standard/control and sample.

Adjust the incubator to 37 ± 1 °C.

1. Dispense 100 µL standards/controls and diluted samples into their respective wells. Leave well A1 for the Substrate Blank.
2. Cover wells with the foil supplied in the kit.
3. **Incubate for 1 hour $\pm$ 5 min at 37 ± 1 °C.**
4. When incubation has been completed, remove the foil, aspirate the content of the wells and wash each well three times with 300 µL of Washing Buffer. Avoid overflows from the reaction wells. The interval between washing and aspiration should be > 5 sec. At the end carefully remove remaining fluid by tapping strips on tissue paper prior to the next step!
Note: Washing is important! Insufficient washing results in poor precision and false results.
5. Dispense 100 µL Conjugate into all wells except for the Substrate Blank well A1.
6. **Incubate for 30 min at room temperature (20...25 °C).** Do not expose to direct sunlight.
7. Repeat step 4.
8. Dispense 100 µL TMB Substrate Solution into all wells.
9. **Incubate for exactly 15 min at room temperature (20...25 °C) in the dark.** A blue colour occurs due to an enzymatic reaction.
10. Dispense 100 µL Stop Solution into all wells in the same order and at the same rate as for the TMB Substrate Solution, thereby a colour change from blue to yellow occurs.
11. Measure the absorbance at 450/620 nm within 30 min after addition of the Stop Solution.

8.1. Measurement

Adjust the ELISA Microtiterplate reader **to zero** using the **Substrate Blank**.

If - due to technical reasons - the ELISA Microtiterplate reader cannot be adjusted to zero using the Substrate Blank, subtract its absorbance value from all other absorbance values measured in order to obtain reliable results!

Measure the absorbance of all wells at **450 nm** and record the absorbance values for each standard/control and sample in the plate layout.

Bichromatic measurement using a reference wavelength of 620 nm is recommended.

Where applicable calculate the mean absorbance values of all duplicates.

9. RESULTS

9.1. Run Validation Criteria

In order for an assay run to be considered valid, these Instructions for Use have to be strictly followed and the following criteria must be met:

- **Substrate Blank:** Absorbance value < **0.100**
- **Negative Control:** Absorbance value < **0.200** and < **Cut-off**
- **Cut-off Control:** Absorbance value **0.150 – 1.300**
- **Positive Control:** Absorbance value > **Cut-off**

If these criteria are not met, the test is not valid and must be repeated.

9.2. Calculation of Results

The Cut-off is the mean absorbance value of the Cut-off Control determinations.

Example: Absorbance value Cut-off Control 0.44 + absorbance value Cut-off control 0.42 = 0.86 / 2 = 0.43
Cut-off = 0.43

9.2.1. Results in Units [NTU]

Sample (mean) absorbance value x 10 = [NovaTec Units = NTU]
Cut-off

Example: $\frac{1.591 \times 10}{0.43} = 37$ NTU (Units)

9.3. Interpretation of Results

Cut-off	10 NTU	-
Positive	> 11 NTU	Antibodies against the pathogen are present. There has been a contact with the antigen (pathogen resp. vaccine).
Equivocal	9 – 11 NTU	Antibodies against the pathogen could not be detected clearly. It is recommended to repeat the test with a fresh sample in 2 to 4 weeks. If the result is equivocal again the sample is judged as negative .
Negative	< 9 NTU	The sample contains no antibodies against the pathogen. A previous contact with the antigen (pathogen resp. vaccine) is unlikely.
Diagnosis of an infectious disease should not be established on the basis of a single test result. A precise diagnosis should take into consideration clinical history, symptomatology as well as serological data. In immunocompromised patients and newborns serological data only have restricted value.		

9.3.1. Antibody Isotypes and State of Infection

Serology	Significance
IgM	Characteristic of the primary antibody response High IgM titer with low IgG titer: → suggests a current or very recent infection Rare: → persisting IgM
IgG	Characteristic of the secondary antibody response May persist for several years High IgG titer with low IgM titer: → may indicate a past infection
IgA	Produced in mucosal linings throughout the body (⇒ protective barrier) Usually produced early in the course of the infection

10. SPECIFIC PERFORMANCE CHARACTERISTICS

The results refer to the groups of samples investigated; these are not guaranteed specifications.

For further information about the specific performance characteristics please contact NovaTec Immundiagnostica GmbH.

10.1. Precision

<u>Intraassay</u>	<u>n</u>	<u>Mean (E)</u>	<u>CV (%)</u>
#1	24	0,376	11,46
#2	24	0,887	7,69
#3	24	1,105	8,01
<u>Interassay</u>	<u>n</u>	<u>Mean (NTU)</u>	<u>CV (%)</u>
#1	12	26,31	3,37
#2	12	25,00	5,09
#3	12	5,66	9,75

10.2. Diagnostic Specificity

The diagnostic specificity is defined as the probability of the assay of scoring negative in the absence of the specific analyte. It is 95.12% (95% confidence interval: 83.47% - 99.4%).

10.3. Diagnostic Sensitivity

The diagnostic sensitivity is defined as the probability of the assay of scoring positive in the presence of the specific analyte. It is 91.67% (95% confidence interval: 81.61% - 97.24%).

10.4. Interferences

Interferences with hemolytic, lipemic or icteric samples are not observed up to a concentration of 10 mg/mL hemoglobin, 5 mg/mL triglycerides and 0.5 mg/mL bilirubin.

10.5. Cross Reactivity

Investigation of a sample panel with antibody activities to potentially cross-reacting parameters did not reveal evidence of false-positive results due to cross-reactions.

A cross reaction with Chlamydia trachomatis cannot be excluded with sera containing antibodies to LPS and MOMP.

11. LIMITATIONS OF THE PROCEDURE

Bacterial contamination or repeated freeze-thaw cycles of the sample may affect the absorbance values.

12. PRECAUTIONS AND WARNINGS

- The test procedure, the information, the precautions and warnings in the instructions for use have to be strictly followed. The use of the testkits with analyzers and similar equipment has to be validated. Any change in design, composition and test procedure as well as for any use in combination with other products not approved by the manufacturer is not authorized; the user himself is responsible for such changes. The manufacturer is not liable for false results and incidents for these reasons. The manufacturer is not liable for any results by visual analysis of the patient samples.
- Only for in-vitro diagnostic use.
- All materials of human or animal origin should be regarded and handled as potentially infectious.
- All components of human origin used for the production of these reagents have been tested for anti-HIV antibodies, anti-HCV antibodies and HBsAg and have been found to be non-reactive.
- Do not interchange reagents or Microtiterplates of different production lots.
- No reagents of other manufacturers should be used along with reagents of this test kit.
- Do not use reagents after expiry date stated on the label.
- Use only clean pipette tips, dispensers, and lab ware.
- Do not interchange screw caps of reagent vials to avoid cross-contamination.
- Close reagent vials tightly immediately after use to avoid evaporation and microbial contamination.
- After first opening and subsequent storage check conjugate and standard/control vials for microbial contamination prior to further use.
- To avoid cross-contamination and falsely elevated results pipette patient samples and dispense reagents without splashing accurately into the wells.
- The ELISA is only designed for qualified personnel following the standards of good laboratory practice (GLP).
- For further internal quality control each laboratory should additionally use known samples.

12.1. Safety note for reagents containing hazardous substances

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Therefore, the following hazard and precautionary statements apply.

Warning

H317	May cause an allergic skin reaction.
P261	Avoid breathing spray
P280	Wear protective gloves/ protective clothing.
P302+P352	IF ON SKIN: Wash with plenty of soap and water.
P333+P313	If skin irritation or rash occurs: Get medical advice/ attention.
P362+P364	Take off contaminated and Wash it before reuse.

Further information can be found in the safety data sheet.

12.2. Disposal Considerations

Residues of chemicals and preparations are generally considered as hazardous waste. The disposal of this kind of waste is regulated through national and regional laws and regulations. Contact your local authorities or waste management companies which will give advice on how to dispose hazardous waste.

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Ballard, Ronald C.; Papp, John R. (2006): Chlamydia trachomatis Infections of the Genital Tract (Oculogenital Infections and Lymphogranuloma Venereum). In Richard L. Guerrant, David H. Walker, Peter F. Weller (Eds.): Tropical infectious diseases. Principles, pathogens & practice. 2nd ed. Philadelphia: Churchill Livingstone, pp. 526–534.

Dawson, Chandler R.; Schachter, Julius (2006): Trachoma. In Richard L. Guerrant, David H. Walker, Peter F. Weller (Eds.): Tropical infectious diseases. Principles, pathogens & practice. 2nd ed. Philadelphia: Churchill Livingstone, 519–525.

Hoyme, U. B.; Spitzbart, H. (1996): Past and current prevalence of Chlamydia trachomatis in women in Germany. In Angelika Stary (Ed.): Proceedings of the third meeting of the European society for chlamydia research. Vienna, Austria, September 11–14, 1996. Bologna: Esculapio, p. 391.

Marrie, Thomas J. (2006): Psittacosis. In Richard L. Guerrant, David H. Walker, Peter F. Weller (Eds.): Tropical infectious diseases. Principles, pathogens & practice. 2nd ed. Philadelphia: Churchill Livingstone, pp. 535–538.

Paavonen, Jorma (1996): Chlamydia trachomatis: A major cause of mucopurulent cervicitis and pelvi inflammatory disease in women. In Peter Elsner (Ed.): Sexually transmitted diseases. Advances in diagnosis and treatment ; 50 tables. Basel: Karger (Current problems in dermatology, 24), pp. 110–122.

Petersen, Eiko E.; Clad, Andreas (1995): Genitale Chlamydien-Infektionen. In *Deutsches Ärzteblatt* 92 (5), pp. 277–282.

Weström, L. V. (1996): Chlamydia and its effect on reproduction. In *Journal of the British Fertility Society* 1 (1), pp. 23–30.

Weström, L. V. (1996): Consequences of genital Chlamydia infections in women. In Angelika Stary (Ed.): Proceedings of the third meeting of the European society for chlamydia research. Vienna, Austria, September 11–14, 1996. Bologna: Esculapio, pp. 137–140.

ABBREVIATIONS / ABKÜRZUNGEN / ABRÉVIATIONS / ABBREVIAZIONI / ABREVIACIONES / ABREVIATURAS

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Chlamydia pneumoniae IgG ELISA

Test Preparation

Prepare reagents and samples as described.
Establish the distribution and identification plan for all samples and standards/controls on the plate layout supplied in the kit.
Select the required number of microtiter strips or wells and insert them into the holder.

Assay Procedure

	Substrate Blank (A1)	Negative Control	Cut-off Control	Positive Control	Sample (diluted 1+100)
Negative Control	-	100 µL	-	-	-
Cut-off Control	-	-	100 µL	-	-
Positive Control	-	-	-	100 µL	-
Sample (diluted 1+100)	-	-	-	-	100 µL
Cover wells with foil supplied in the kit Incubate for 1 h at 37 ± 1 °C Wash each well three times with 300 µL of Washing Buffer					
Conjugate	-	100 µL	100 µL	100 µL	100 µL
Incubate for 30 min at room temperature (20...25 °C) Do not expose to direct sunlight Wash each well three times with 300 µL of Washing Buffer					
TMB Substrate Solution	100 µL	100 µL	100 µL	100 µL	100 µL
Incubate for exactly 15 min at room temperature (20...25 °C) in the dark					
Stop Solution	100 µL	100 µL	100 µL	100 µL	100 µL
Photometric measurement at 450 nm (reference wavelength: 620 nm)					

NovaTec Immundiagnostica GmbH

Waldstraße 23 A6
63128 Dietzenbach, Germany

Tel.: +49 (0) 6074-48760
Email: info@NovaTec-ID.com
Internet: www.NovaTec-ID.com

Fax: +49 (0) 6074-487629

CHLG0510-2020-06-04-Ka-gr.Lot 163

NovaLisa®

Chlamydia pneumoniae IgM

ELISA

CE 0483

Only for in-vitro diagnostic use

English	2
Deutsch	7
Français	12
Italiano	17
Español	22
Português	27
Bibliography / Literatur / Bibliographie / Bibliografia / Bibliografía Bibliografia	34
Abbreviations / Abkürzungen / Abréviations / Abbreviazioni / Abreviaciones / Abreviaturas	34
Symbols Key / Symbolschlüssel / Explication des Symboles / Legenda / Símbolos / Tabela de símbolos	35
Summary of Test Procedure / Kurzanleitung Testdurchführung / Résumé de la procedure de test / Schema della procedura / Resumen de la técnica / Resumo do Procedimento de Teste	36

Product Number: CHLM0510 (96 Determinations)

ENGLISH

1. INTRODUCTION

Chlamydiae are non-motile, Gram negative and obligatory intracellular growing bacteria which form characteristic inclusions within the cytoplasm of parasitized cells. They are easily visible in the light microscope. Three different Chlamydia species pathogenic for humans are known: Chlamydia trachomatis, Chlamydia pneumoniae and Chlamydia psittaci, and one species only pathogenic for animals (C. pecorum). Chlamydia trachomatis is the most prevalent agent of sexually transmitted diseases worldwide (400-500 million cases) and the number of infections is constantly growing. Pregnant women infected with C. trachomatis may transmit these bacteria during childbirth, causing conjunctivitis or pneumonia in newborns. Untreated cases of chlamydial infection can lead to chronic salpingitis, possibly resulting in ectopic pregnancy or infertility. In males, C. trachomatis is a major cause of non-gonococcal urethritis. A severe problem in Chlamydia infections is the frequent asymptomatic insidious course which may result in the initiation of chronic diseases. In many instances primary infections are not recognized and only the sequelae caused by ascended, persisting agents are diagnosed.

Species	Disease	Symptoms (e.g.)	Transmission route
Chlamydia trachomatis	Trachoma Lymphogranuloma venereum (LGV) Urogenital infections	Bilateral conjunctivitis with foreign body sensation, tearing eyes, and discharge of serous secretions from the eye. Genital lesion; swelling of the inguinal and femoral lymph nodes; formation of suppurating buboes. In men: Inflammation of the urethra (urethritis), inflammation of the prostate (prostatitis) and the epididymis (epididymitis). Sterility (infertility). In women: Inflammation of the urethra (urethritis), Bartholin's glands, inflammation of the cervix (cervicitis), uterine mucous membrane (endometritis) and fallopian tubes (salpingitis) and vaginal discharge. Sterility (Infertility)	Sexually transmitted. Can be passed from an infected mother to her baby during childbirth.
C. pneumoniae	Respiratory diseases (e.g. pneumonia, bronchitis) discussed: endocarditis, coronary heart diseases	Fever, cough (usually productive)	Man to Man Aerogenic, droplet infection
C. psittaci	Ornithosis (Psittacosis)	Hepatosplenomegaly; Unproductive cough, fever and severe headache	Inhalation of feces from infected birds; contact with infected avian viscera

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- Microscopy
- PCR
- Serology: Detection of antigens e.g. by ELISA
Detection of antibodies e.g. by ELISA

2. INTENDED USE

The Chlamydia pneumoniae IgM ELISA is intended for the qualitative determination of IgM class antibodies against Chlamydia pneumoniae in human serum or plasma (citrate, heparin).

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4. MATERIALS

4.1. Reagents supplied

- **Microtiterplate:** 12 break-apart 8-well snap-off strips coated with *Chlamydia pneumoniae* antigens; in resealable aluminium foil.
- **IgM Sample Dilution Buffer:** 1 bottle containing 100 mL of phosphate buffer (10 mM) for sample dilution; pH 7.2 ± 0.2 ; anti-human IgG (RF Absorbent); coloured green; ready to use; white cap; $\leq 0.0015\%$ (v/v) CMIT/ MIT (3:1).
- **Stop Solution:** 1 bottle containing 15 mL sulphuric acid, 0.2 mol/L; ready to use; red cap.
- **Washing Buffer (20x conc.):** 1 bottle containing 50 mL of a 20-fold concentrated phosphate buffer (0.2 M), pH 7.2 ± 0.2 , for washing the wells; white cap.
- **Conjugate:** 1 bottle containing 20 mL of peroxidase labelled antibody to human IgM in phosphate buffer (10 mM); coloured red; ready to use; black cap.
- **TMB Substrate Solution:** 1 bottle containing 15 mL 3,3',5,5'-tetramethylbenzidine (TMB), $< 0.1\%$; ready to use; yellow cap.
- **Positive Control:** 1 vial containing 2 mL control; coloured yellow; ready to use; red cap; $\leq 0.02\%$ (v/v) MIT.
- **Cut-off Control:** 1 vial containing 3 mL control; coloured yellow; ready to use; green cap; $\leq 0.02\%$ (v/v) MIT.
- **Negative Control:** 1 vial containing 2 mL control; coloured yellow; ready to use; blue cap; $\leq 0.0015\%$ (v/v) CMIT/ MIT (3:1).

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4.3. Materials and Equipment needed

- ELISA Microtiterplate reader, equipped for the measurement of absorbance at 450/620 nm
- Incubator 37 °C
- Manual or automatic equipment for rinsing Microtiterplates
- Pipettes to deliver volumes between 10 and 1000 μL
- Vortex tube mixer
- Distilled water
- Disposable tubes

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7.1. Sample Dilution

Before assaying, all samples should be diluted 1+100 with IgM Sample Dilution Buffer. Dispense 10 μL sample and 1 mL IgM Sample Dilution Buffer into tubes to obtain a 1+100 dilution and thoroughly mix with a Vortex.

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10.2. Diagnostic Specificity

The diagnostic specificity is defined as the probability of the assay of scoring negative in the absence of the specific analyte.

It is 99.26% (95% confidence interval: 95.97% - 99.98%).

10.3. Diagnostic Sensitivity

The diagnostic sensitivity is defined as the probability of the assay of scoring positive in the presence of the specific analyte.

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10.5. Cross Reactivity

Investigation of a sample panel with antibody activities to potentially cross-reacting parameters did not reveal significant evidence of false-positive results due to cross-reactions.

A cross reaction with *Chlamydia trachomatis* cannot be excluded with sera containing antibodies to LPS and MOMP.

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