Brilliance™ Candida Agar

REF **CM1002B**

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Bibliografie

- Public Health England. 2015a. 'Investigation of Nasal Samples'. UK Standards for Microbiology Investigations B 5 (7.1). <https://www.gov.uk/government/publications/smi-b-5-investigation-of-nose-swabs>.
- Müller, F. M.C., A. H. Groll, and T. J. Walsh. 1999. 'Current Approaches to Diagnosis and Treatment of Fungal Infections in Children Infected with Human Immuno Deficiency Virus'. *European Journal of Pediatrics* 158 (3): 187–99. <https://doi.org/10.1007/s004310051047>.
- Public Health England. 2013. 'Sexually Transmitted Infections'. UK Standards for Microbiology Investigations S 6 (1.2). <https://www.gov.uk/government/publications/smi-s-6-sexually-transmitted-infections>.
- Public Health England. 2015b. 'Investigation of Throat Related Specimens'. UK Standards for Microbiology Investigations B 9 (9). <https://www.gov.uk/government/publications/smi-b-9-investigation-of-throat-swabs>.
- Public Health England. 2016. 'Investigation of Dermatological Specimens for Superficial Mycoses'. UK Standards for Microbiology Investigation B 39 (3.1). <https://www.gov.uk/government/publications/smi-b-39-investigation-of-dermatological-specimens-for-superficial-mycoses>.
- Public Health England. 2017a. 'Genital Tract and Associated Specimens'. UK Standards for

Microbiology Investigation B 28 (4.6).

<https://www.gov.uk/government/publications/smi-b-28-investigation-of-genital-tract-and-associated-specimens>.

- Public Health England. 2017b. 'Introduction to the Preliminary Identification of Medically Important Bacteria and Fungi from Culture'. UK Standards for Microbiology Investigations ID 1 (2). <https://www.gov.uk/government/publications/smi-id-1-preliminary-identification-of-medically-important-bacteria>.
- Public Health England. 2017c. 'Investigation of Intravascular Cannulae and Associated Specimens'. UK Standards for Microbiology Investigations B 20 (6.1). <https://www.gov.uk/government/publications/smi-b-20-investigation-of-intravascular-cannulae-and-associated-specimens>.
- Public Health England. 2018. 'Swabs from Skin and Superficial Soft Tissue Infections'. UK Standards for Microbiology Investigations B 11 (6.5). <https://www.gov.uk/government/publications/smi-b-11-investigation-of-skin-superficial-and-non-surgical-wound-swabs>.
- Public Health England. 2019a. 'Investigation of Bronchoalveolar Lavage, Sputum and Associated Specimens'. UK Standards for Microbiology Investigations B 57 (3.5). <https://www.gov.uk/government/publications/smi-b-57-investigation-of-bronchoalveolar-lavage-sputum-and-associated-specimens>.
- Public Health England. 2019b. 'Investigation of Urine'. UK Standards for Microbiology Investigation B 41 (8.7). <https://www.gov.uk/government/publications/smi-b-41-investigation-of-urine>.
- Sganga, Gabriele. 2009. 'Clinical Aspects of Invasive Candidiasis in the Surgical Patient'. *Drugs* 69: 29–32. <https://doi.org/10.2165/11315600-000000000-00000>.
- Chen, Sharon C.A., D. Marriott, E. G. Playford, Q. Nguyen, D. Ellis, W. Meyer, T. C. Sorrell, et al. 2009. 'Candidaemia with Uncommon Candida Species: Predisposing Factors, Outcome, Antifungal Susceptibility, and Implications for Management'. *Clinical Microbiology and Infection* 15 (7): 662–69. <https://doi.org/10.1111/j.1469-0691.2009.02821.x>.
- Farooqi, J. Q., K. Jabeen, N. Saeed, N. Iqbal, B. Malik, S. R. Lockhart, A. Zafar, M. E. Brandt, and R. Hasan. 2013. 'Invasive Candidiasis in Pakistan: Clinical Characteristics, Species Distribution and Antifungal Susceptibility'. *Journal of Medical Microbiology* 62 (2): 259–68. <https://doi.org/10.1099/jmm.0.048785-0>.
- Pongrácz, Júlia, Emese Juhász, Miklós Iván, and Katalin Kristóf. 2015. 'Significance of Yeasts in Bloodstream Infection: Epidemiology and Predisposing Factors of Candidaemia in Adult Patients at a University Hospital (2010-2014)'. *Acta Microbiologica et Immunologica Hungarica* 62 (3): 317–29. <https://doi.org/10.1556/030.62.2015.3.9>.
- Khadka, Sundar, Jeevan Bahadur Sherchand, Bharat Mani Pokhrel, Keshab Parajuli, Shyam Kumar Mishra, Sangita Sharma, Niranjan Shah, et al. 2017. 'Isolation, Speciation and Antifungal Susceptibility Testing of Candida Isolates from Various Clinical Specimens at a Tertiary Care Hospital, Nepal'. *BMC Research Notes* 10 (1): 1–5. <https://doi.org/10.1186/s13104-017-2547-3>.
- Xiang, Guo Shi, Ling Bing Zeng, Jie Yu Zhang, Ying Ying, Xi Ren Deng, Xue Fei Hu, Yan Yan, Ke Hua Yu, Wei Wen Zou, and Xiao Tian Huang. 2018. 'In Vitro Antidrug Susceptibility Testing of Candida Species Isolated from Aseptic Body Fluids'. *Jundishapur Journal of Microbiology* 11 (8). <https://doi.org/10.5812/ijm.55547>.

18. Xiao, Meng, Xin Fan, Xin Hou, Sharon C.A. Chen, He Wang, Fanrong Kong, Zi Yong Sun, Yun Zhuo Chu, and Ying Chun Xu. 2018. 'Clinical Characteristics of the First Cases of Invasive Candidiasis in China Due to Pan-Echinocandin-Resistant Candida Tropicalis and Candida Glabrata Isolates with Delineation of Their Resistance Mechanisms'. Infection and Drug Resistance 11: 155–61. <https://doi.org/10.2147/IDR.S152785>.

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Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, UK

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©2022 Thermo Fisher Scientific Inc. Toate drepturile rezervate.
Toate celelalte mărci comerciale aparțin Thermo Fisher Scientific Inc. și subsidiarelor acesteia.