



PRODUCT INFORMATION

**Thermo Scientific  
GeneJET Genomic  
DNA Purification Kit  
#K0721, #K0722**

Pub. No. MAN0012663

Rev. Date 12 October 2016 (Rev. B.00)



Read Storage information (p. 2) before first use!

[www.thermofisher.com](http://www.thermofisher.com)

**For Research Use Only. Not for use in diagnostic procedures.**

#\_  
Lot\_  
Exp. \_

### CERTIFICATE OF ANALYSIS

Thermo Scientific GeneJET Genomic DNA Purification Kit is qualified by isolating genomic DNA from 200  $\mu$ L of blood and 5 mg of mammalian tissue following described protocols. The purified genomic DNA has an  $A_{260/280}$  ratio of  $\geq 1.7$ . A single band of more than 30 kb is seen after agarose gel electrophoresis and ethidium bromide staining. Functional quality of genomic DNA is evaluated by PCR amplification of a single-copy gene and digestion with restriction enzymes.

**Quality authorized by:**



Jurgita Zilinskiene

<b>CONTENTS</b>	page
COMPONENTS OF THE KIT .....	2
STORAGE.....	2
DESCRIPTION .....	2
PRINCIPLE .....	2
IMPORTANT NOTES .....	3
ADDITIONAL MATERIALS AND EQUIPMENT REQUIRED.....	3
GENOMIC DNA PURIFICATION PROTOCOLS.....	4
A. Mammalian Tissue and Rodent Tail Genomic DNA Purification Protocol .....	4
B. Cultured Mammalian Cells Genomic DNA Purification Protocol .....	6
C. Mammalian Blood Genomic DNA Purification Protocol .....	7
D. Gram-Negative Bacteria Genomic DNA Purification Protocol .....	8
E. Gram-Positive Bacteria Genomic DNA Purification Protocol .....	9
F. Yeast Genomic DNA Purification Protocol.....	10
G. DNA Purification from Buccal Swabs .....	11
TROUBLESHOOTING .....	12

## COMPONENTS OF THE KIT

GeneJET Genomic DNA Purification Kit	#K0721 50 preps	#K0722 250 preps
Proteinase K Solution	1.2 mL	5 × 1.2 mL
RNase A Solution	1 mL	5 × 1 mL
Digestion Solution	11 mL	55 mL
Lysis Solution	24 mL	2 × 60 mL
Wash Buffer I (concentrated)	10 mL	40 mL
Wash Buffer II (concentrated)	10 mL	40 mL
Elution Buffer (10 mM Tris-Cl, pH 9.0, 0.1 mM EDTA)	30 mL	150 mL
GeneJET Genomic DNA Purification Columns pre-assembled with Collection Tubes	50	250
Collection Tubes	50	250

### STORAGE

Proteinase K and RNase A solutions are stable at room temperature as long as not opened. After being opened they should be stored at -20 °C. Other components of the kit should be stored at room temperature (15-25 °C).

**Note. Close the bag with GeneJET Genomic DNA Purification Columns tightly after each use!**

### DESCRIPTION

The GeneJET™ Genomic DNA Purification Kit is designed for rapid and efficient purification of high quality genomic DNA from various mammalian cell culture and tissue samples, whole blood, bacteria and yeast. The kit utilizes silica-based membrane technology in the form of a convenient spin column, eliminating the need for expensive resins, toxic phenol-chloroform extractions, or time-consuming alcohol precipitation. The standard procedure takes less than 20 minutes following cell lysis and yields purified DNA of more than 30 kb in size. Isolated DNA can be used directly in PCR, Southern blotting and enzymatic reactions. See Table 1 for typical genomic DNA yields from various sources.

### PRINCIPLE

Depending on the starting material, samples are digested with Proteinase K in either the supplied Digestion or Lysis Solution. RNA is removed by treating the samples with RNase A. The lysate is then mixed with ethanol and loaded on the purification column where the DNA binds to the silica membrane. Impurities are effectively removed by washing the column with the prepared wash buffers. Genomic DNA is then eluted under low ionic strength conditions with the Elution Buffer.

**Table 1.** Typical genomic DNA yields from various sources.

Source	Quantity	Yield, µg
Mammalian blood	200 µL	4-6
Mouse heart	10 mg	10-15
Mouse tail	0.5 cm	8-10
Rat liver	10 mg	10-20
Rat spleen	5 mg	20-30
Rat kidney	10 mg	25-30
Rabbit ear	20 mg	5-10
<i>Bacillus pumilis</i> cells	2×10 <sup>9</sup> cells	10-15
<i>Escherichia coli</i> cells	2×10 <sup>9</sup> cells	10-15
HeLa cells	2×10 <sup>6</sup> cells	15-20
Jurkat cells	5×10 <sup>6</sup> cells	25-30
<i>Saccharomyces cerevisiae</i> cells	1×10 <sup>8</sup> cells	3-5

## IMPORTANT NOTES

- To minimize DNA degradation, avoid repeated freeze/thaw cycles of the samples and perform extractions from fresh material or material that has been immediately frozen and stored at -20 °C or -70 °C.
- Add the indicated volume of ethanol (96-100%) to **Wash Buffer I** (concentrated) and **Wash Buffer II** (concentrated) prior to first use:

	#K0721 50 preps		#K0722 250 preps	
	Wash Buffer I	Wash Buffer II	Wash Buffer I	Wash Buffer II
Concentrated wash solution	10 mL	10 mL	40 mL	40 mL
<b>Ethanol (96-100%)</b>	<b>30 mL</b>	<b>30 mL</b>	<b>120 mL</b>	<b>120 mL</b>
Total volume:	40 mL	40 mL	160 mL	160 mL

After the ethanol has been added, mark the check box on the bottle's cap to indicate the completed step.

- Check the **Digestion Solution** and **Lysis Solution** for salt precipitation before each use. Re-dissolve any precipitate by warming the solution at 37 °C, then cool back down to 25 °C before use.
- Wear gloves when handling the **Lysis Solution** and **Wash Buffer I** as these reagents contain irritants.

## ADDITIONAL MATERIALS AND EQUIPMENT REQUIRED

- Pipets and pipet tips
- Vortex
- Ethanol (96-100%)
- 1.5 mL microcentrifuge tubes
- Microcentrifuge
- Thermomixer, shaking water bath or rocking platform capable of heating up to 56 °C
- Disposable gloves

### Buffers

For mammalian cell lysate preparation:

- PBS (137 mM NaCl, 2.7 mM KCl, 10 mM Na<sub>2</sub>HPO<sub>4</sub>, 2 mM KH<sub>2</sub>PO<sub>4</sub>, pH 7.4)
- TE buffer (10 mM Tris-HCl, pH 8.0, 1 mM EDTA)

For gram-positive bacteria lysate preparation

- Gram-positive bacteria lysis buffer (20 mM Tris-HCl, pH 8.0, 2 mM EDTA, 1.2% Triton X-100, add lysozyme to 20 mg/mL immediately before use)

For yeast lysate preparation:

- Yeast lysis buffer (5 mg/mL zymolyase 20T, 1 M sorbitol, 0.1 M EDTA)

## GENOMIC DNA PURIFICATION PROTOCOLS

Protocols for genomic DNA purification from mammalian tissue and rodent tail, cultured mammalian cells, mammalian blood, gram-negative, gram-positive bacteria, yeast and buccal swabs are described on p.4-11.

### A. Mammalian Tissue and Rodent Tail Genomic DNA Purification Protocol

Step	Procedure														
1	Grind up to 20 mg of mammalian tissue (use up to 10 mg of spleen tissue), 0.6 cm (rat) or 0.5 cm (mouse) tail clip in liquid nitrogen using a mortar and pestle. Alternatively, cut the tissue into small pieces or disrupt it using a homogenizer.														
2	Collect the material into a 1.5 mL microcentrifuge tube (not provided) and resuspend in 180 $\mu$ L of Digestion Solution. Add 20 $\mu$ L of Proteinase K Solution and mix thoroughly by vortexing or pipetting to obtain a uniform suspension.														
3	<p>Incubate the sample at 56 °C until the tissue is completely lysed and no particles remain. During incubation vortex the vial occasionally or use a shaking water bath, rocking platform or thermomixer.</p> <p>Suggested incubation times:</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th style="text-align: center;">Quantity</th> <th style="text-align: center;">Suggested incubation time</th> </tr> </thead> <tbody> <tr> <td>5 mg of tissue (except spleen)</td> <td>1 hour</td> </tr> <tr> <td>10 mg of tissue (except spleen)</td> <td>2 hours</td> </tr> <tr> <td>20 mg of tissue (except spleen)</td> <td>3 hours</td> </tr> <tr> <td>5 mg of spleen tissue</td> <td>2 hours</td> </tr> <tr> <td>10 mg of spleen tissue</td> <td>3 hours</td> </tr> <tr> <td>Mouse tail (0.5 cm), rat tail (0.6 cm)</td> <td>6 hours</td> </tr> </tbody> </table> <p><b>Note.</b> Lysis time varies on the type and amount of tissue processed. In some cases incubation time should be prolonged to 6-8 hours or overnight (for rodent tail) until complete lysis occurs.</p>	Quantity	Suggested incubation time	5 mg of tissue (except spleen)	1 hour	10 mg of tissue (except spleen)	2 hours	20 mg of tissue (except spleen)	3 hours	5 mg of spleen tissue	2 hours	10 mg of spleen tissue	3 hours	Mouse tail (0.5 cm), rat tail (0.6 cm)	6 hours
Quantity	Suggested incubation time														
5 mg of tissue (except spleen)	1 hour														
10 mg of tissue (except spleen)	2 hours														
20 mg of tissue (except spleen)	3 hours														
5 mg of spleen tissue	2 hours														
10 mg of spleen tissue	3 hours														
Mouse tail (0.5 cm), rat tail (0.6 cm)	6 hours														
4	Add 20 $\mu$ L of RNase A Solution, mix by vortexing then incubate for 10 min at room temperature.														
5	Add 200 $\mu$ L of Lysis Solution. Mix thoroughly by vortexing for 15 s until a homogeneous mixture is obtained.														
6	Add 400 $\mu$ L of 50% ethanol and mix by pipetting or vortexing.														
7	<p>Transfer the prepared lysate to a GeneJET Genomic DNA Purification Column inserted in a collection tube. Centrifuge the column for 1 min at 6000 <math>\times</math> g. Discard the collection tube containing the flow-through solution. Place the GeneJET Genomic DNA Purification Column into a new 2 mL collection tube (included).</p> <p><b>Note. Close the bag with GeneJET Genomic DNA Purification Columns tightly after each use!</b></p>														
8	Add 500 $\mu$ L of Wash Buffer I (with ethanol added). Centrifuge for 1 min at 8000 $\times$ g. Discard the flow-through and place the purification column back into the collection tube.														

Step	Procedure
9	<p>Add 500 <math>\mu</math>L of Wash Buffer II (with ethanol added) to the GeneJET Genomic DNA Purification Column. Centrifuge for 3 min at maximum speed (<math>\geq 12000 \times g</math>).</p> <p><i>Optional.</i> If residual solution is seen in the purification column, empty the collection tube and re-spin the column for 1 min. at maximum speed.</p> <p>Discard the collection tube containing the flow-through solution and transfer the GeneJET Genomic DNA Purification Column to a sterile 1.5 mL microcentrifuge tube (not included).</p>
10	<p>Add 200 <math>\mu</math>L of Elution Buffer to the center of the GeneJET Genomic DNA Purification Column membrane to elute genomic DNA. Incubate for 2 min at room temperature and centrifuge for 1 min at <math>8000 \times g</math>.</p> <p><b>Note</b></p> <ul style="list-style-type: none"> <li>• For maximum DNA yield, repeat the elution step with additional 200 <math>\mu</math>L of Elution Buffer.</li> <li>• If more concentrated DNA is required or DNA is isolated from a small amount of starting material (e.g., <math>&lt;5</math> mg of tissue) the volume of the Elution Buffer added to the column can be reduced to 50-100 <math>\mu</math>L. Please be aware that smaller volumes of Elution Buffer will result in smaller final quantity of eluted DNA.</li> </ul>
11	<p>Discard the purification column. Use the purified DNA immediately in downstream applications or store at <math>-20</math> <math>^{\circ}</math>C.</p>

## B. Cultured Mammalian Cells Genomic DNA Purification Protocol

Step	Procedure
1	<p>a) <u>Suspension cells</u> Collect up to <math>5 \times 10^6</math> cells in a centrifuge tube. Pellet cells by centrifugation for 5 min at <math>250 \times g</math>. Discard the supernatant. Rinse cells once with PBS to remove residual medium and repeat the centrifugation step. Discard the supernatant.</p> <p>b) <u>Adherent cells</u> Remove the growth medium from a culture plate containing up to <math>2 \times 10^6</math> cells. Rinse cells once with PBS to remove residual medium. Discard PBS. Detach the cells from the culture plate by scraping in an appropriate volume of PBS or by trypsinization. Transfer the cells to a microcentrifuge tube and pellet them by centrifugation for 5 minutes at <math>250 \times g</math>. Discard supernatant.</p>
2	Resuspend the cells collected in step 1a or 1b in 200 $\mu$ L of TE buffer or PBS. Add 200 $\mu$ L of Lysis Solution and 20 $\mu$ L of Proteinase K Solution to the cell pellet. Mix thoroughly by vortexing or pipetting to obtain a uniform suspension.
3	Incubate the sample at 56 °C while vortexing occasionally or use a shaking water bath, rocking platform or thermomixer until the cells are completely lysed (10 min).
4	Add 20 $\mu$ L of RNase A Solution, mix by vortexing and incubate the mixture for 10 min at room temperature.
5	Add 400 $\mu$ L of 50% ethanol and mix by pipetting or vortexing.
6	Transfer the prepared lysate to a GeneJET Genomic DNA Purification Column inserted in a collection tube. Centrifuge the column for 1 min at $6000 \times g$ . Discard the collection tube containing the flow-through solution. Place the GeneJET Genomic DNA Purification Column into a new 2 mL collection tube (included). <b>Note. Close the bag with GeneJET Genomic DNA Purification Columns tightly after each use!</b>
7	Add 500 $\mu$ L of Wash Buffer I (with ethanol added). Centrifuge for 1 min at $8000 \times g$ . Discard the flow-through and place the purification column back into the collection tube.
8	Add 500 $\mu$ L of Wash Buffer II (with ethanol added) to the GeneJET Genomic DNA Purification Column. Centrifuge for 3 min at maximum speed ( $\geq 12000 \times g$ ). <i>Optional.</i> If residual solution is seen in the purification column, empty the collection tube and re-spin the column for 1 min at maximum speed. Discard the collection tube containing the flow-through solution and transfer the GeneJET Genomic DNA Purification Column to a sterile 1.5 mL microcentrifuge tube (not included).
9	Add 200 $\mu$ L of Elution Buffer to the center of the GeneJET Genomic DNA Purification Column membrane to elute genomic DNA. Incubate for 2 min at room temperature and centrifuge for 1 min at $8000 \times g$ . <b>Note</b> <ul style="list-style-type: none"> <li>For maximum DNA yield, repeat the elution step with additional 200 <math>\mu</math>L of Elution Buffer.</li> <li>If more concentrated DNA is required or DNA is isolated from a small amount of starting material (e.g., <math>\leq 1 \times 10^6</math> of cultured mammalian cells) the volume of the Elution Buffer added to the column can be reduced to 50-100 <math>\mu</math>L. Please be aware that smaller volumes of Elution Buffer will result in smaller final quantity of eluted DNA.</li> </ul>
10	Discard the purification column. Use the purified DNA immediately in downstream applications or store at -20 °C.



### C. Mammalian Blood Genomic DNA Purification Protocol

Step	Procedure
1	Add 400 $\mu$ L of Lysis Solution and 20 $\mu$ L of Proteinase K Solution to 200 $\mu$ L of whole blood, mix thoroughly by vortexing or pipetting to obtain a uniform suspension.
2	Incubate the sample at 56 °C while vortexing occasionally or use a shaking water bath, rocking platform or thermomixer until the cells are completely lysed (10 min).
3	Add 200 $\mu$ L of ethanol (96-100%) and mix by pipetting or vortexing.
4	Transfer the prepared lysate to a GeneJET Genomic DNA Purification Column inserted in a collection tube. Centrifuge the column for 1 min at 6000 $\times$ g. Discard the collection tube containing the flow-through solution. Place the GeneJET Genomic DNA Purification Column into a new 2 mL collection tube (included). <b>Note. Close the bag with GeneJET Genomic DNA Purification Columns tightly after each use!</b>
5	Add 500 $\mu$ L of Wash Buffer I (with ethanol added). Centrifuge for 1 min at 8000 $\times$ g. Discard the flow-through and place the purification column back into the collection tube.
6	Add 500 $\mu$ L of Wash Buffer II (with ethanol added) to the GeneJET Genomic DNA Purification Column. Centrifuge for 3 min at maximum speed ( $\geq$ 12000 $\times$ g). <i>Optional.</i> If residual solution is seen in the purification column, empty the collection tube and re-spin the column for 1 min. at maximum speed. Discard the collection tube containing the flow-through solution and transfer the GeneJET Genomic DNA Purification Column to a sterile 1.5 mL microcentrifuge tube (not included).
7	Add 200 $\mu$ L of Elution Buffer to the center of the GeneJET Genomic DNA Purification Column membrane to elute genomic DNA. Incubate for 2 min at room temperature and centrifuge for 1 min at 8000 $\times$ g. <b>Note</b> <ul style="list-style-type: none"> <li>• For maximum DNA yield, repeat the elution step with additional 200 <math>\mu</math>L of Elution Buffer.</li> <li>• If more concentrated DNA is required or DNA is isolated from a small amount of starting material (e.g., 50 <math>\mu</math>L) the volume of the Elution Buffer added to the column can be reduced to 50-100 <math>\mu</math>L. Please be aware that smaller volumes of Elution Buffer will result in smaller final quantity of eluted DNA.</li> </ul>
8	Discard the purification column. Use the purified DNA immediately in downstream applications or store at -20 °C.

## D. Gram-Negative Bacteria Genomic DNA Purification Protocol

Step	Procedure
1	Harvest up to $2 \times 10^9$ bacterial cells in a 1.5 or 2 mL microcentrifuge tube by centrifugation for 10 min at $5000 \times g$ . Discard the supernatant.
2	Resuspend the pellet in 180 $\mu$ L of Digestion Solution. Add 20 $\mu$ L of Proteinase K Solution and mix thoroughly by vortexing or pipetting to obtain a uniform suspension.
3	Incubate the sample at 56 °C while vortexing occasionally or use a shaking water bath, rocking platform or thermomixer until the cells are completely lysed (~30 min).
4	Add 20 $\mu$ L of RNase A Solution, mix by vortexing and incubate the mixture for 10 min at room temperature.
5	Add 200 $\mu$ L of Lysis Solution to the sample. Mix thoroughly by vortexing for about 15 s until a homogeneous mixture is obtained.
6	Add 400 $\mu$ L of 50% ethanol and mix by pipetting or vortexing.
7	Transfer the prepared lysate to a GeneJET Genomic DNA Purification Column inserted in a collection tube. Centrifuge the column for 1 min at $6000 \times g$ . Discard the collection tube containing the flow-through solution. Place the GeneJET Genomic DNA Purification Column into a new 2 mL collection tube (included). <b>Note. Close the bag with GeneJET Genomic DNA Purification Columns tightly after each use!</b>
8	Add 500 $\mu$ L of Wash Buffer I (with ethanol added). Centrifuge for 1 min at $8000 \times g$ . Discard the flow-through and place the purification column back into the collection tube.
9	Add 500 $\mu$ L of Wash Buffer II (with ethanol added) to the GeneJET Genomic DNA Purification Column. Centrifuge for 3 min at maximum speed ( $\geq 12000 \times g$ ). <i>Optional.</i> If residual solution is seen in the purification column, empty the collection tube and re-spin the column for 1 min. at maximum speed. Discard the collection tube containing the flow-through solution and transfer the GeneJET Genomic DNA Purification Column to a sterile 1.5 mL microcentrifuge tube (not included).
10	Add 200 $\mu$ L of Elution Buffer to the center of the GeneJET Genomic DNA Purification Column membrane to elute genomic DNA. Incubate for 2 min at room temperature and centrifuge for 1 min at $8000 \times g$ . <b>Note</b> <ul style="list-style-type: none"> <li>• For maximum DNA yield, repeat the elution step with additional 200 <math>\mu</math>L of Elution Buffer.</li> <li>• If more concentrated DNA is required or DNA is isolated from a small amount of starting material the volume of the Elution Buffer added to the column can be reduced to 50-100 <math>\mu</math>L. Please be aware that smaller volumes of Elution Buffer will result in smaller final quantity of eluted DNA.</li> </ul>
11	Discard the purification column. Use the purified DNA immediately in downstream applications or store at -20 °C.

## E. Gram-Positive Bacteria Genomic DNA Purification Protocol

### Before starting

Prepare Gram-positive bacteria lysis buffer: 20 mM Tris-HCl, pH 8.0, 2 mM EDTA, 1.2% Triton X-100, add lysozyme to 20 mg/mL immediately before use.

Step	Procedure
1	Harvest up to $2 \times 10^9$ bacterial cells in a 1.5 or 2 mL microcentrifuge tube by centrifugation for 10 min at $5000 \times g$ . Discard the supernatant.
2	Resuspend the pellet in 180 $\mu$ L of Gram-positive bacteria lysis buffer. Incubate for 30 min at 37 °C.
3	Add 200 $\mu$ L of Lysis Solution and 20 $\mu$ L of Proteinase K. Mix thoroughly by vortexing or pipetting to obtain a uniform suspension.
4	Incubate the sample at 56 °C while vortexing occasionally or use a shaking water bath, rocking platform or thermomixer until the cells are completely lysed (~30 min).
5	Add 20 $\mu$ L of RNase A Solution, mix by vortexing and incubate the mixture for 10 min at room temperature.
6	Add 400 $\mu$ L of 50% ethanol and mix by pipetting or vortexing.
7	Transfer the prepared lysate to a GeneJET Genomic DNA Purification Column inserted in a collection tube. Centrifuge the column for 1 min at $6000 \times g$ . Discard the collection tube containing the flow-through solution. Place the GeneJET Genomic DNA Purification Column into a new 2 mL collection tube (included). <b>Note. Close the bag with GeneJET Genomic DNA Purification Columns tightly after each use!</b>
8	Add 500 $\mu$ L of Wash Buffer I (with ethanol added). Centrifuge for 1 min at $8000 \times g$ . Discard the flow-through and place the purification column back into the collection tube.
9	Add 500 $\mu$ L of Wash Buffer II (with ethanol added) to the GeneJET Genomic DNA Purification Column. Centrifuge for 3 min at maximum speed ( $\geq 12000 \times g$ ). <i>Optional.</i> If residual solution is seen in the purification column, empty the collection tube and re-spin the column for 1 min. at maximum speed. Discard the collection tube containing the flow-through solution and transfer the GeneJET Genomic DNA Purification Column to a sterile 1.5 mL microcentrifuge tube (not included).
10	Add 200 $\mu$ L of Elution Buffer to the center of the GeneJET Genomic DNA Purification Column membrane to elute genomic DNA. Incubate for 2 min at room temperature and centrifuge for 1 min at $8000 \times g$ . <b>Note</b> <ul style="list-style-type: none"><li>• For maximum DNA yield, repeat the elution step with additional 200 <math>\mu</math>L of Elution Buffer.</li><li>• If more concentrated DNA is required or DNA is isolated from a small amount of starting material the volume of the Elution Buffer added to the column can be reduced to 50-100 <math>\mu</math>L. Please be aware that smaller volumes of Elution Buffer will result in smaller final quantity of eluted DNA.</li></ul>
11	Discard the purification column. Use the purified DNA immediately in downstream applications or store at -20 °C.

## F. Yeast Genomic DNA Purification Protocol

### Before starting

Prepare Yeast lysis buffer: 5 mg/mL zymolyase 20T, 1 M sorbitol, 0.1 M EDTA.

Step	Procedure
1	Harvest up to $1 \times 10^8$ yeast cells in a 1.5 or 2 mL microcentrifuge tube by centrifugation for 5-10 s at maximum speed $\geq 12000 \times g$ . Discard the supernatant.
2	Resuspend the pellet in 500 $\mu$ L of Yeast lysis buffer. Incubate for 1 hour at 37 °C.
3	Centrifuge cells for 10 min at $3000 \times g$ . Discard the supernatant.
4	Resuspend the pellet in 180 $\mu$ L of Digestion Solution. Add 20 $\mu$ L of Proteinase K Solution and mix thoroughly by vortexing or pipetting to obtain a uniform suspension.
5	Incubate the sample at 56 °C while vortexing occasionally or use a shaking water bath, rocking platform or thermomixer until the cells are completely lysed (~45 min).
6	Add 20 $\mu$ L of RNase A Solution, mix by vortexing and incubate the mixture for 10 min at room temperature.
7	Add 200 $\mu$ L of Lysis Solution. Mix thoroughly by vortexing for 15 s until a homogeneous mixture is obtained.
8	Add 400 $\mu$ L of 50% ethanol and mix by pipetting or vortexing.
9	Transfer the prepared lysate to a GeneJET Genomic DNA Purification Column inserted in a collection tube. Centrifuge the column for 1 min at $6000 \times g$ . Discard the collection tube containing the flow-through solution. Place the GeneJET Genomic DNA Purification Column into a new 2 mL collection tube (included). <b>Note.</b> Close the bag with GeneJET Genomic DNA Purification Columns tightly after each use!
10	Add 500 $\mu$ L of Wash Buffer I (with ethanol added). Centrifuge for 1 min at $8000 \times g$ . Discard the flow-through and place the purification column back into the collection tube.
11	Add 500 $\mu$ L of Wash Buffer II (with ethanol added) to the GeneJET Genomic DNA Purification Column. Centrifuge for 3 min at maximum speed ( $\geq 12000 \times g$ ). <i>Optional.</i> If residual solution is seen in the purification column, empty the collection tube and re-spin the column for 1 min. at maximum speed. Discard the collection tube containing the flow-through solution and transfer the GeneJET Genomic DNA Purification Column to a sterile 1.5 mL microcentrifuge tube (not included).
12	Add 200 $\mu$ L of Elution Buffer to the center of the GeneJET Genomic DNA Purification Column membrane to elute genomic DNA. Incubate for 2 min at room temperature and centrifuge for 1 min at $8000 \times g$ . <b>Note</b> <ul style="list-style-type: none"><li>• For maximum DNA yield, repeat the elution step with additional 200 <math>\mu</math>L of Elution Buffer.</li><li>• If more concentrated DNA is required or DNA is isolated from a small amount of starting material the volume of the Elution Buffer added to the column can be reduced to 50-100 <math>\mu</math>L. Please be aware that smaller volumes of Elution Buffer will result in smaller final quantity of eluted DNA.</li></ul>
13	Discard the purification column. Use the purified DNA immediately in downstream applications or store at -20 °C.

### G. DNA Purification from Buccal Swabs

Step	Procedure
1	To collect a sample, scrape the swab 5-6 times against the inside cheek.
2	Swirl the swab for 30-60 s in 200 $\mu$ L of 1 $\times$ PBS.
3	Go to step 1 of the standard Mammalian Blood Genomic DNA Purification Protocol (p. 7).

## TROUBLESHOOTING

Problem	Possible cause and solution
<b>Low yield of purified DNA</b>	<p><b>Excess sample used during lysate preparation.</b> Reduce the amount of starting material. Do not use more tissue or cells than indicated in lysis protocols.</p> <p><b>Starting material was not completely digested.</b> Extend the Proteinase K digestion at 56 °C until complete lysis occurs and no particles remain.</p> <p><b>Ethanol was not added to the lysate.</b> Make sure that the ethanol was added to the lysate before applying the sample to the Purification Column.</p> <p><b>Ethanol was not mixed with the lysate.</b> After the addition of ethanol to the lysate mix the sample by vortexing or pipetting.</p> <p><b>Ethanol was not added to Wash Buffers.</b> Make sure that ethanol was added to Wash Buffer I and Wash Buffer II before use. Follow the instructions for Wash Buffer preparation on p.3.</p>
<b>Purified DNA is degraded</b>	<p><b>Sample was frozen and thawed repeatedly.</b> Avoid repeated freeze / thaw cycles of the samples. Use a new sample for DNA isolation. Perform extractions from fresh material when possible.</p> <p><b>Inappropriate sample storage conditions.</b> Store mammalian tissues at -70 °C and bacteria at -20 °C until use. Whole blood can be stored at 4 °C for no longer than 1-2 days. For long term storage blood samples should be aliquoted in 200 µL portions and stored at -20 °C.</p>
<b>RNA contamination</b>	<p><b>RNase A treatment was not carried out.</b> Carry out RNase A treatment step described in the purification procedure.</p>
<b>Column becomes clogged during purification</b>	<p><b>Excess sample was used during lysate preparation.</b> Reduce the amount of starting material. A maximum of <math>2 \times 10^9</math> of bacteria cells, <math>5 \times 10^6</math> of suspension cells and 20 mg of mammalian tissue is recommended for lysate preparation.</p> <p><b>Tissue was not completely digested.</b> Extend the Proteinase K digestion at 56 °C until complete lysis occurs and no particles remain.</p>
<b>Inhibition of downstream enzymatic reactions</b>	<p><b>Purified DNA contains residual ethanol.</b> If residual solution is seen in the purification column after washing the column with Wash Buffer II, empty the collection tube and re-spin the column for an additional 1 min. at maximum speed (<math>\geq 12000 \times g</math>).</p> <p><b>Purified DNA contains residual salt.</b> Use the correct order for the Washing Buffers. Always wash the purification column with Wash Buffer I first and then proceed to washing with Wash Buffer II.</p>

### **LIMITED USE LABEL LICENSE: Internal Research and Development Use Only.**

The purchase of this product conveys to the buyer the limited, non-exclusive, non-transferable right (without the right to resell, repackage, or further sublicense) to use this product for internal research and development purposes. No other license is granted to the buyer whether expressly, by implication, by estoppel or otherwise. In particular, the purchase of the product does not include or carry any right or license to use, develop, or otherwise exploit this product commercially and no rights are conveyed to the buyer to use the product or components of the product for purposes including but not limited to provision of services to a third party, generation of commercial databases or clinical diagnostics. This product is sold pursuant to authorization from Thermo Fisher Scientific and Thermo Fisher Scientific reserves all other rights. For information on purchasing a license for uses other than internal research and development purposes, please contact [outlicensing@lifetech.com](mailto:outlicensing@lifetech.com) or Out Licensing, Life Technologies Inc., 5781 Van Allen Way, Carlsbad, California 92008.

### **PRODUCT USE LIMITATION**

This product is developed, designed and sold exclusively *for research purposes and in vitro use only*. The product was not tested for use in diagnostics or for drug development, nor is it suitable for administration to humans or animals.

Please refer to [www.thermofisher.com](http://www.thermofisher.com) for Material Safety Data Sheet of the product.

© 2016 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

# TaqMan<sup>®</sup> Genotyping Master Mix

Catalog Numbers 4371353, 4371355, 4381656, 4371357, and 4381657

Pub. No. 4374656 Rev. D

**Note:** For safety and biohazard guidelines, see the “Safety” appendix in the *TaqMan<sup>®</sup> Genotyping Master Mix Protocol* (Pub. No. 4371131). Read the Safety Data Sheets (SDSs) and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves.

## Contents and storage

Cat. No.	Number of reactions	Contents	Storage
4371353	40	1 mL	2–8°C for up to one year
4371355	400	10 mL	
4381656	800	2 × 10 mL	
4371357	2,000	50 mL	
4381657	4,000	2 × 50 mL	



## Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale found on Life Technologies' website at [www.thermofisher.com/us/en/home/global/terms-and-conditions.html](http://www.thermofisher.com/us/en/home/global/terms-and-conditions.html). If you have any questions, please contact Life Technologies at [www.thermofisher.com/support](http://www.thermofisher.com/support).



**Manufacturer:** Thermo Fisher Scientific Baltics UAB | V.A. Graiciuno 8, LT-02241 | Vilnius, Lithuania

The information in this guide is subject to change without notice.

**DISCLAIMER:** TO THE EXTENT ALLOWED BY LAW, THERMO FISHER SCIENTIFIC INC. AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

**Revision history:** Pub. No. 4374656

Revision	Date	Description
D	04 September 2018	Updated manufacturing address, branding, licensing, trademarks, general style and format.
C	September 2011	Baseline for this revision.

**Important Licensing Information:** These products may be covered by one or more Limited Use Label Licenses. By use of these products, you accept the terms and conditions of all applicable Limited Use Label Licenses.

©2018 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. TaqMan is a registered trademark of Roche Molecular Systems, Inc., used under permission and license.


# TaqMan SNP Genotyping Assays

Applied Biosystems™ TaqMan™ SNP Genotyping Assays from Thermo Fisher Scientific provide a highly flexible technology for detection of polymorphisms within any genome. Applied Biosystems™ TaqMan™ Assays have the simplest workflow available and are the quickest way to generate genotyping data. Based on powerful Applied Biosystems™ TaqMan™ probe and primer chemistry and designs, and coupled to dependable Applied Biosystems™ instruments and software, these made-to-order assays produce high-confidence results. TaqMan Assays are ideal for genotyping applications, including screening, association, candidate region, candidate gene, and fine-mapping studies.

Content-rich marker selection tools simplify study design and help you select from a library of human and mouse assays. This library includes over 7 million genome-wide human assays (of which 3.5 million are HapMap SNP-based assays, 160,000 are validated assays, and over 950,000 are coding region assays) and 10,000 mouse assays. We also offer 2,700 inventoried drug metabolism genotyping assays. Additionally, with Applied Biosystems™ Custom TaqMan™ SNP Genotyping Assays you can confidentially submit target SNP sequences for any genome to create your own assays. Let TaqMan SNP Genotyping Assays accelerate the pace of your discovery by eliminating time-consuming experimental design and optimization.

## Powerful, proven chemistry

Whether your genotyping studies require targeted detection of essential SNPs, or the flexibility to choose SNPs for mapping, TaqMan SNP Genotyping Assays are the technology of choice. Proven TaqMan



**Figure 1. Allelic discrimination is achieved by the selective annealing of TaqMan MGB probes.**

probes, which incorporate minor groove binder (MGB) technology at the 3' end, deliver superior allelic discrimination. The MGB molecule binds to the minor groove of the DNA helix, improving hybridization-based assays by stabilizing the MGB probe–template complex. This increased binding stability permits the use of probes as short as 13 bases for improved mismatch discrimination and greater flexibility when designing assays for difficult or variable sequences. In addition to SNP detection, TaqMan probes can be designed to detect multiple nucleotide polymorphisms (MNPs) and insertion/deletions (indels).

Detection is achieved with proven 5' nuclease chemistry by means of exonuclease cleavage of an allele-specific 5' dye label, which generates the permanent assay signal (Figure 1). All MGB probes include a nonfluorescent

quencher (NFQ) that virtually eliminates the background fluorescence associated with traditional quenchers, and provides a greater signal-to-noise ratio for superior assay sensitivity.

## TaqMan SNP Genotyping Assays collection

TaqMan SNP Genotyping Assays are the world's largest collection of single-tube, ready-to-use SNP assays available. The TaqMan SNP Genotyping Assays library consists of two collections of human assays and one of mouse assays, and can be supplemented with assays designed using our Custom TaqMan SNP Genotyping Assays Service.

## Over 7 million human SNP genotyping assays

This assay group contains over 7 million genome-wide SNPs, providing unprecedented marker coverage. Included in this collection are 160,000 validated assays that have approximately 10 kb spacing across gene regions. These assays were subjected to an extensive minor allele frequency test in 2–4 ethnic populations (45 individual samples per ethnic group) and as a result, offer the highest success rate. Also included are over 600,000 assays for the detection of nonsynonymous SNPs in coding regions, including many putative functional SNPs. Visit [thermofisher.com/taqmansnp](http://thermofisher.com/taqmansnp) for more information.

## Over 10,000 mouse SNP genotyping assays

The Applied Biosystems™ Mouse TaqMan™ Predesigned SNP Genotyping Assays collection consists of over 10,000 assays, and can be supplemented with assays designed using our Custom TaqMan SNP Genotyping Assays Service.

## TaqMan Drug Metabolism Genotyping Assays

The collection of Applied Biosystems™ TaqMan™ Drug Metabolism Genotyping Assays includes 2,700 assays that target high-value polymorphisms in 221 drug metabolism genes. These assays have proven performance in four different ethnic population samples, consisting of 45 individuals each. To enable easy identification, these assays have been mapped to the common public allele nomenclature websites where possible. Visit [thermofisher.com/taqmandme](http://thermofisher.com/taqmandme) for more information.

All TaqMan SNP Genotyping Assays are generated using next-generation algorithms from the Thermo Fisher Scientific bioinformatics pipeline. For all predesigned assays, bioinformatics evaluation of target SNP sequences includes the masking of adjacent SNPs and ambiguous bases so that assay design and subsequent performance is not affected by the poor quality of the underlying sequence. Lastly,

the assay designs are aligned to the human genome using BLAST to ensure that each assay binds uniquely to the intended polymorphism. As the Custom TaqMan SNP Genotyping Assay Service is confidential and secure, you simply perform your own bioinformatics analysis prior to submitting your sequence for assay design.

## Custom assay service for any possible SNP

Custom TaqMan SNP Genotyping Assays can be developed for any SNP in any organism. This service can generate assays for the detection of SNPs, MNPs, indels of up to 6 bases, or QSY™-labeled probes for multiplexing SNP assays.

Custom TaqMan SNP Genotyping Assays provide you with a complete service that includes secure and confidential ordering, assay design and manufacturing, and quality-control testing for synthesis accuracy and formulation completeness. Additionally, custom human assays are subjected to a functional test on 20 unique DNA samples.

Use the free Applied Biosystems™ Custom TaqMan™ Assay Design Tool to input and submit your sequence for assay design. This easy-to-use online resource lets you quickly submit your sequence information and start the ordering process securely and confidentially. Access the Custom TaqMan Assay Design Tool at [thermofisher.com/snpcadt](http://thermofisher.com/snpcadt)


## Quality design and manufacturing

Probes and primers used in TaqMan SNP Genotyping Assays are designed using our rigorous bioinformatics pipeline. This proprietary group of algorithms has generated millions of TaqMan Assay designs by utilizing heuristic design rules deduced from both manufacturing and assay performance data. All assays are designed to perform under universal reaction conditions, as calculated probe and primer melting temperatures are consistent and include contributions from associated probe conjugates (i.e., dyes and MGB).

After manufacturing, assay components undergo extensive laboratory testing at our state-of-the-art manufacturing facility. Quality-control testing includes mass spectrometry for sequence verification and formulation assessments of probe and primer concentrations. Additionally, all human SNP genotyping assays are functionally tested to ensure allelic discrimination.

## Simple workflow for quick results

TaqMan SNP Genotyping Assays constitute the simplest SNP genotyping technology available. We deliver your ready-to-use SNP genotyping assay at ambient temperature in a convenient, single-tube format. The rest is easy. Just combine the assay with Applied Biosystems™ TaqMan™ Genotyping Master Mix or TaqMan™ Universal PCR Master Mix and your purified DNA sample (Figure 2). There is no need to optimize probe, primer, salt concentrations, or temperature because all assays use universal reagent concentrations and thermal cycling conditions.



**Figure 2. A simple workflow and reliable instruments combine to generate fast, high-confidence results.**

After generating an endpoint read using a thermal cycler or real-time PCR instrument, no transfers, washes, or additional reagents are required, and the plate remains sealed; just read the plate and analyze the genotypes. This reduces the chance of contamination, sample mix-up, and sample loss. The simplicity of the chemistry allows you to easily automate the reaction for massively parallel genotyping studies, readily increasing the number of assays, number of samples, or both. Additionally, the analysis software allows you to auto-call genotypes, minimizing manual intervention.

### Reliable real-time PCR platforms

A suite of superior Applied Biosystems instrument platforms is available for processing and analyzing TaqMan SNP Genotyping Assays (Table 1). These instruments, which meet all throughput needs and budgets, include the 7500, 7500 Fast, 7900HT Fast, ViiA™ 7, StepOne™, and StepOnePlus™ Real-Time PCR Systems, and the QuantStudio™ 3, 5, 6, 7, and 12K Flex Real-Time PCR Systems (Figure 3). Following PCR amplification, an endpoint read can be performed on any Applied Biosystems real-time PCR system. All of these dependable instruments offer the advanced multicolor detection capabilities required for highly accurate and reproducible allelic discrimination assays.

### Data analysis software


The sophisticated SDS software package provided with all Applied Biosystems real-time PCR systems facilitates experimental setup, data collection, and assay performance analysis. The SDS software uses an advanced multicomponent algorithm to calculate the distinct signal contribution of each allele of a marker from the fluorescence measurements of each sample well during the assay plate read. The multicomponent data collected from the plate read are stored as SDS files, ready for genotype determination by the SDS software or optional Applied Biosystems™ TaqMan™ Genotyper Software (Figure 4).

**Table 1. Applied Biosystems instrument capacities.**

Instrument	Capacity
7500/7500 Fast Real-Time PCR System	96-well block (standard or Fast)
7900HT Fast Real-Time PCR System	96- and 384-well blocks (standard or Fast)
ViiA 7 Real-Time PCR System	96-well (standard or Fast), 384-well, and TaqMan Array Card blocks
StepOne Real-Time PCR System	48-well block (standard or Fast)
StepOnePlus Real-Time PCR System	96-well block (standard or Fast)
QuantStudio 3 Real-Time PCR System	96-well block (standard or Fast)
QuantStudio 5 Real-Time PCR System	96-well (standard or Fast) and 384-well blocks
QuantStudio 6 Real-Time PCR System	96-well (standard or Fast) and 384-well blocks
QuantStudio 7 Real-Time PCR System	96-well (standard or Fast), 384-well, and TaqMan Array Card blocks
QuantStudio 12K Flex Real-Time PCR System	96-well (standard or Fast), 384-well, TaqMan Array Card, and OpenArray plate blocks



**Figure 3. The QuantStudio 5 Real-Time PCR System (left) and the QuantStudio 12K Flex Real-Time PCR System (right), which offers the highest throughput of all Applied Biosystems real-time PCR instruments.**



**Figure 4. TaqMan Genotyper Software automatically determines sample genotypes and displays data.**

TaqMan Genotyper Software is a great resource for fast and accurate genotype calling. It is a free SNP genotyping data analysis tool for use with TaqMan SNP Genotyping Assays performed in 48-, 96-, or 384-well microtiter plates or on Applied Biosystems™ TaqMan™ OpenArray™ Plates. It has a state-of-the-art genotype-calling algorithm, an intuitive user interface, and enhanced study-based analysis features. The software enables multi-plate data analysis for high-throughput workflows and improved accuracy in genotype calling; versatile export features and comprehensive quality-control features facilitate streamlining of the entire workflow. TaqMan Genotyper Software can be downloaded at [thermofisher.com/taqmangenotyper](http://thermofisher.com/taqmangenotyper)

### Simple ordering

Selecting and ordering TaqMan SNP Genotyping Assays is as simple as “point and click.” Use SNPbrowser Software to select the most informative SNPs for your genotyping studies. As you identify SNPs of interest, simply upload your selected TaqMan SNP Genotyping Assays to our online ordering tool.

The TaqMan Assay online ordering tool (Figure 5) enables you to search, select, and order from our catalog of over 7 million made-to-order predesigned TaqMan SNP Genotyping Assays. You can search for SNPs using any of several criteria: National Center for Biotechnology Information (NCBI) gene ID, NCBI SNP reference ID (rs#), or gene symbol. You can further refine your search by using SNP type (i.e., intragenic, 5’ or 3’ UTR, chromosome, etc.).

Our Custom TaqMan SNP Genotyping Assays supply you with SNPs that are not available from our predesigned assay collection, including those from any nonhuman organism. This service designs assays for all possible SNP, MNP, and indel targets but without the up-front bioinformatic preparation used for the predesigned made-to-order assays. Our complementary Custom TaqMan Assay Design Tool conveniently formats your target sequence for submission to our manufacturing facilities. To order custom assays, simply prepare your target sequence according to the Design and Ordering Guide, and upload your submission file at [thermofisher.com/snpcaadt](http://thermofisher.com/snpcaadt)

The screenshot displays the TaqMan Assay search and order tool interface. On the left, there are search filters for Species (Human), Gene (BRCA1), Assay Type (Functionally Tested), and SNP Type (Intragenic, Silent Mutation, Transition Substitution). The main search bar contains the query 'c\_60612706\_20'. The results show a table with columns for SNP ID, Gene, Location, SNP Type, Assay Type, and Made To Order | Cat #. The selected SNP is rs34844365, located on Chromosome 17 at position 41251811. The SNP Type is Silent Mutation, Transition Substitution, Intragenic. The Assay Type is Functionally Tested. The Made To Order | Cat # is 4351379. Below the table, there are buttons for 'View Details', 'View Assay on Map', 'List Price (USD)', 'Check availability', and 'Add To Cart'. The detailed view includes Product Details (SNP ID: rs34844365, Assay Type: Functionally Tested, NCBI dbSNP Submissions: 2, Location: Chr. 17:41251811 on NCBI Build 37, Set Membership, Context Sequence [VIC/FAM], Phenotype, Polymorphism, Allele Nomenclature, Literature Links) and Gene Details (Gene: BRCA1, Gene Name: breast cancer 1, early onset). A table of Transcript Accession, SNP Location, SNP Type, Codon Change, Amine Acid Change, and Protein ID is also shown.

Transcript Accession	SNP Location	SNP Type	Codon Change	Amine Acid Change	Protein ID
NM_007254.3	760	Silent Mutation	ACA.ACG	T176T	NP_009225.1
NM_007257.3	760	Silent Mutation	ACA.ACG	T129T	NP_009228.2
NM_007258.3	760	Silent Mutation	ACA.ACG	T176T	NP_009229.2
NM_007299.3	760	Silent Mutation	ACA.ACG	T176T	NP_009230.3
NM_007350.3	760	Silent Mutation	ACA.ACG	T176T	NP_009231.2

**Figure 5.** Our TaqMan Assay search and order tool makes online ordering easy. For convenient online ordering and multiple search options for all our genotyping assays, including keyword, batch, and location searches, visit [thermofisher.com/taqmansnp](http://thermofisher.com/taqmansnp)



## Ordering information

Size	Human assays (Cat. No.)	Nonhuman assays (Cat. No.)	Number of SNPs	No. of 5 $\mu$ L reactions (384-well plate)	No. of 25 $\mu$ L reactions (96-well plate)	Assay mix formulation	Assay type
<b>Predesigned TaqMan SNP Genotyping Assays for Human and Mouse</b>							
Small	4351379	4351384*	>4.5 million	1,500	300	40X	made-to-order
Medium	4351376	4351382*	>4.5 million	5,000	1,000	40X	made-to-order
Large	4351374	4351380*	>4.5 million	12,000	2,400	80X	made-to-order
<b>Custom TaqMan SNP Genotyping Assays</b>							
Small	4331349	4332077	$\infty$	1,500	300	40X	made-to-order
Medium	4332072	4332075	$\infty$	5,000	1,000	40X	made-to-order
Large	4332073	4332076	$\infty$	12,000	2,400	80X	made-to-order
<b>TaqMan Drug Metabolism Genotyping Assays</b>							
Small	4362691	NA	2,700	750	150	20X	inventoried

All assays are quality-control tested using a mass spectrometer to verify sequence and yield. All assays have a VIC™ dye-labeled probe, a FAM™ dye-labeled probe, and two target-specific primers. All assays, excluding Custom TaqMan SNP Genotyping Assays, undergo bioinformatics evaluation of target SNP sequences.

Functional testing against 20 unique genomic DNA samples is performed on all custom and predesigned made-to-order human TaqMan SNP Genotyping Assays. Validation testing against four populations (45 samples/population) was performed on all 160,000 validated TaqMan SNP Genotyping Assays, and all TaqMan Drug Metabolism Genotyping Assays.

\* Over 10,000 mouse assays available.

applied  
biosystems

Find out more at [thermofisher.com/taqmansnp](https://thermofisher.com/taqmansnp)

**ThermoFisher**  
SCIENTIFIC

**For Research Use Only. Not for use in diagnostic procedures.** © 2016 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. TaqMan is a trademark of Roche Molecular Systems, Inc., used under permission and license. **CO020271 0116**

## Safety Data Sheet

Safety Data Sheet according to Regulation (EC) No. 1907/2006 (REACH)  
Classification according to Regulation (EC) No. 1272/2008 [CLP]

### SECTION 1: Identification of the substance/mixture and of the company/undertaking

#### Product identifier

**Product code** 4351376  
**Product name** Dye Labeled Oligonucleotide Assays  
**Chemical Name** Not Applicable  
**REACH registration number** Formamide : 01-2119496064-35-XXXX

#### Relevant identified uses of the substance or mixture and uses advised against

**Relevant identified uses** For research use only  
**Use Description Code** SU22 - Professional uses: Public domain (administration, education, entertainment, services, craftsmen), PROC15 - Use as laboratory reagent, PC21 - Laboratory chemicals, SU24 - Scientific research and development  
**Uses advised against** Not for consumer use.

#### Details of the supplier of the safety data sheet

##### Manufacturer / Supplier

LIFE TECHNOLOGIES EUROPE BV  
KWARTSWEG 2  
2665 NN BLEISWIJK  
NETHERLANDS  
31-(0)180 392 400  
Email: MSDS@lifetech.com

Life Technologies Limited  
3 Fountain Drive  
Inchinnan Business Park  
Paisley  
PA4 9RF, UK  
+44 (0)141 814 6100

**24 hour Emergency Response for Hazardous Materials** Within the USA + Canada: 1-800-424-9300 and  
**[or Dangerous Goods] Incident. Spill, Leak, Fire,** 1-703-527-3887  
**Exposure, or Accident. Call CHEMTREC** Outside the USA + Canada: 1-703-741-5970

#### Country Specific Emergency Number (if available):

CHEMTREC Ireland (Dublin) +(353)-19014670 (Greeting Language: English and Irish)  
CHEMTREC UK (London) +(44)-870-8200418 (Greeting Language: English)

### SECTION 2: Hazards identification

**Revision date** 05-Mar-2020  
**Product code** 4351376

**Page** 1 / 10  
**Product name** Dye Labeled Oligonucleotide Assays



## Classification of the substance or mixture

Classification according to Regulation (EC) No. 1272/2008 [CLP]

### Physical hazards

Not Hazardous

### Health hazards

Carcinogenicity	Category 2
Reproductive Toxicity	Category 1B
Specific target organ toxicity - Repeated exposure	Category 2

### Environmental hazards

Not Hazardous

### Additional information

No information available

## Label elements

Labelling according to Regulation (EC) No 1272/2008 [CLP]

### Hazard pictograms



### Signal Word

Danger

### Hazard Statements

H351 - Suspected of causing cancer if swallowed

H360 - May damage fertility or the unborn child if swallowed

H373 - May cause damage to organs through prolonged or repeated exposure

### Precautionary Statements

#### Prevention

P201 - Obtain special instructions before use

P202 - Do not handle until all safety precautions have been read and understood

P260 - Do not breathe dust/fume/gas/mist/vapours/spray

P280 - Wear protective gloves/protective clothing/eye protection/face protection

#### Response

P308 + P313 - IF exposed or concerned: Get medical advice/attention

P314 - Get medical advice/attention if you feel unwell

#### Storage

Not Applicable

#### Disposal

P501 - Dispose of contents/ container to an approved waste disposal plant

Revision date 05-Mar-2020  
Product code 4351376

Page 2 / 10  
Product name Dye Labeled Oligonucleotide Assays

**Other hazards**

Not Applicable

## SECTION 3: Composition/information on ingredients

Component	CAS No	EINECS-No.	Weight-%	REACH registration number	Classification according to Regulation (EC) No. 1272/2008 [CLP]
Formamide 75-12-7 ( 1-5 )	75-12-7	200-842-0	1-5	01-2119496064-35-X XXX	Repr. 1B - H360 Carc. 2 - H351 STOT RE 2 - H373

## SECTION 4: First aid measures

**Description of first aid measures**

<b>Skin contact</b>	Wash off immediately with plenty of water for at least 15 minutes. Remove and wash contaminated clothing and gloves, including the inside, before re-use. Immediate medical attention is required.
<b>Eye contact</b>	Rinse immediately with plenty of water, also under the eyelids, for at least 15 minutes. Immediate medical attention is required.
<b>Ingestion</b>	Never give anything by mouth to an unconscious person. Do not induce vomiting without medical advice. If swallowed, rinse mouth with water (only if the person is conscious). Risk of serious damage to the lungs (by aspiration). Get medical attention if symptoms occur.
<b>Inhalation</b>	Remove to fresh air. If not breathing, give artificial respiration. If symptoms persist, call a doctor.
<b>Notes to Physician</b>	Treat symptomatically.

**Most important symptoms and effects, both acute and delayed**

H360 - May damage fertility or the unborn child if swallowed H351 - Suspected of causing cancer if swallowed H373 - May cause damage to organs through prolonged or repeated exposure

**Indication of any immediate medical attention and special treatment needed**

IF exposed or concerned: Get medical advice/attention. Get medical advice/attention if you feel unwell.

## SECTION 5: Firefighting measures

**Extinguishing media**

<b>Suitable extinguishing media</b>	Foam. Dry powder. Dry chemical. Carbon dioxide (CO <sub>2</sub> ). Water spray.
<b>Unsuitable extinguishing media</b>	No information available.

**Special hazards arising from the substance or mixture**

None known

**Protective equipment and precautions for firefighters**

Wear self-contained breathing apparatus and protective suit.

## SECTION 6: Accidental release measures

### **Personal precautions, protective equipment and emergency procedures**

Avoid exposure to vapour  
Avoid breathing vapours or mists  
Ensure adequate ventilation  
Avoid contact with skin, eyes or clothing  
Use personal protection equipment  
See section 8 for more information

### **Environmental precautions**

No special environmental precautions required.

### **Methods and material for containment and cleaning up**

Soak up with inert absorbent material (e.g. sand, silica gel, acid binder, universal binder, sawdust). Sweep up and shovel into suitable containers for disposal. Clean contaminated surface thoroughly.

### **Reference to other sections**

See section 8 for more information.

## SECTION 7: Handling and storage

### **Precautions for safe handling**

Always wear recommended Personal Protective Equipment. Wash hands before breaks and immediately after handling the product. Do not get in eyes, on skin, or on clothing. Avoid breathing vapours or mists. If during normal use the material presents a respiratory hazard, use adequate ventilation and/or wear appropriate respirator. See section 8 for more information.

### **Conditions for safe storage, including any incompatibilities**

Keep in properly labelled containers. Keep in a dry, cool and well-ventilated place. Store in accordance with local regulations. Keep away from combustible material.

### **Specific end use(s)**

For research use only.

## SECTION 8: Exposure controls/personal protection

### Control parameters

Chemical Name	EU OEL (TWA)	EU OEL (STEL)	EU Skin Notation
Formamide 75-12-7	None	None	None
Chemical Name	Austria	Belgium (TWA)	Czech Republic
Formamide 75-12-7	9 ppm 16 mg/m <sup>3</sup>	10 ppm 18 mg/m <sup>3</sup>	None
Chemical Name	Denmark (TWA)	Finland OEL (TWA)	France OEL (VME)
Formamide 75-12-7	10 ppm 18 mg/m <sup>3</sup>	10 ppm 19 mg/m <sup>3</sup>	20 ppm 30 mg/m <sup>3</sup>
Chemical Name	Germany OEL (TWA)	Ireland (TWA)	Italy OEL (TWA)
Formamide 75-12-7	None	10 ppm 18 mg/m <sup>3</sup>	None
Chemical Name	Lithuania OEL (TWA)	Netherlands OEL (MAC)	Norway
Formamide 75-12-7	10 ppm 20 mg/m <sup>3</sup>	None	10 ppm TWA 18 mg/m <sup>3</sup> TWA 15 ppm STEL 27 mg/m <sup>3</sup> STEL
Chemical Name	Poland	Portugal	Spain OEL (TWA)
Formamide 75-12-7	23 mg/m <sup>3</sup> TWA Skin Notation	10 ppm TWA skin - potential for cutaneous exposure	10 ppm 19 mg/m <sup>3</sup>
Chemical Name	Sweden - Occupational Exposure Limits - TLVs (LLVs)	Switzerland	United Kingdom
Formamide 75-12-7	10 ppm TLV NGV; 20 mg/m <sup>3</sup> TLV NGV	10 ppm TWA 18 mg/m <sup>3</sup> TWA	20 ppm TWA; 37 mg/m <sup>3</sup> TWA

**Engineering Measures**      Ensure adequate ventilation, especially in confined areas.

### Exposure controls

#### Personal protection equipment

**Respiratory protection**      In case of insufficient ventilation wear respirators and components tested and approved under appropriate government standards.

**Hand protection**      Glove material: Nitrile rubber with thickness (mm) :5 Break through time (hours) :>1  
Recommended glove type has not been tested for use with this product.  
Information is based on professional knowledge

**Eye protection**      Tight sealing safety goggles.

**Skin and Body Protection**      Wear laboratory coat for body protection.

**Hygiene Measures**      Handle in accordance with good industrial hygiene and safety practice.

**Environmental exposure controls**  
No special environmental precautions required.

## SECTION 9: Physical and chemical properties

### Information on basic physical and chemical properties

<b>Appearance</b>	liquid	
<b>Odour</b>	No data	
<b>Odour Threshold</b>	No data	
<b>Molecular Weight</b>	No data	
<b>pH</b>	8	
<b>Melting point / melting range</b>	°C No data	°F No data
<b>Boiling point / boiling range</b>	°C No data	°F No data
<b>Flash point</b>	°C No data	°F No data
<b>Autoignition Temperature</b>	°C No data	°F No data
<b>Decomposition temperature</b>	°C No data	°F No data
<b>Evaporation rate</b>	No data	
<b>Flammability (solid, gas)</b>	No data	
<b>Upper explosion limit</b>	No data	
<b>Lower explosion limit</b>	No data	
<b>Vapour Pressure</b>	No data	
<b>Vapour density</b>	No data	
<b>Relative density</b>	No data	
<b>Specific gravity</b>	No data	
<b>Solubility</b>	No data	
<b>Partition coefficient: n-octanol/water</b>	No data	
<b>Viscosity</b>	No data	
<b>Explosive properties</b>	No data	
<b>Oxidising properties</b>	No data	

### Other information

No data.

## SECTION 10: Stability and reactivity

<b>Reactivity</b>	None known.
<b>Chemical stability</b>	Stable under normal conditions.
<b>Possibility of hazardous reactions</b>	Hazardous reaction has not been reported.
<b>Conditions to avoid</b>	No information available.
<b>Incompatible materials</b>	No information available.
<b>Hazardous decomposition products</b>	No data available.

## SECTION 11: Toxicological information

### Information on toxicological effects

Chemical Name	Oral LD50	Dermal LD50	Inhalation LC50
Formamide	3200 mg/kg	13500 mg/kg	3900 ppm/6H

### Principal Routes of Exposure

<b>Skin corrosion/irritation</b>	Data are conclusive but insufficient for classification
<b>Serious eye damage/irritation</b>	Data are conclusive but insufficient for classification
<b>Respiratory or skin sensitisation</b>	Data are conclusive but insufficient for classification
<b>Specific target organ toxicity (STOT) – single exposure</b>	Data are conclusive but insufficient for classification
<b>Specific target organ toxicity (STOT) – repeated exposure</b>	Target organ(s) : Cardiovascular System Hematopoietic System
<b>Carcinogenicity</b>	Contains a known or suspected carcinogen
<b>Germ cell mutagenicity</b>	Data are conclusive but insufficient for classification
<b>Reproductive Toxicity</b>	May cause adverse reproductive effects - such as birth defect, miscarriages, or infertility
<b>Aspiration Hazard</b>	Data are conclusive but insufficient for classification

## SECTION 12: Ecological information

### Ecotoxicity

The environmental impact of this product has not been fully investigated.

Chemical Name	Toxicity to algae	Toxicity to daphnia and other aquatic invertebrates	Toxicity to fish	Microtox Data	log Pow
Formamide	Desmodesmus subspicatus EC50>500 mg/L (72 h) Desmodesmus subspicatus EC50>500 mg/L (96 h)	Daphnia magna EC50>500 mg/L (48 h)	No data available	No data available	logPow-0.82

**Persistence and degradability** No information available.

**Bioaccumulative potential** No information available.

### Results of PBT and vPvB assessment

No information available.

### Other adverse effects

No information available.

## SECTION 13: Disposal considerations

### Waste treatment methods

The generation of waste should be avoided or minimized wherever possible. Empty containers or liners may retain some product residues. This material and its container must be disposed of in accordance with approved disposal technique. Disposal of this product, its solutions or of any by-products, shall comply with the requirements of all applicable local, regional or national/federal regulations.

## SECTION 14: Transport information

### IATA / ADR / DOT-US / IMDG

Not regulated in the meaning of transport regulations

<b>UN number</b>	Not Applicable
<b>UN proper shipping name</b>	Not Applicable
<b>Transport hazard class(es)</b>	Not Applicable
<b>Packing group</b>	Not Applicable

### Environmental hazards

Not Applicable

### Special precautions for user

Not Applicable

### Transport in bulk according to Annex II of MARPOL and the IBC Code

Not Applicable.

## SECTION 15: Regulatory information

### Safety, health and environmental regulations/legislation specific for the substance or mixture

#### Substances of Very High Concern

Chemical Name	Weight-%	EU - REACH (1907/2006) - Article 59(1) - Candidate List of Substances for Eventual Inclusion in Annex XIV
Formamide	1-5	Reason for inclusion Toxic for reproduction, Article 57c

### Substance subject to authorisation per REACH Annex XIV

None

### Restricted substances under EC 1907/2006, Annex XVII

Chemical Name	Weight-%	EU - REACH (1907/2006) - Annex XVII - Restrictions on Certain Dangerous Substances
Formamide	1-5	Use restricted. See item 30.

### Substances listed under Annex I of Regulation (EC) No 689/2008

None.

Revision date 05-Mar-2020  
Product code 4351376

Page 8 / 10  
Product name Dye Labeled Oligonucleotide Assays

**Restricted substances under Annex V of Regulation (EC) No 689/2008**

None.

**Substances under Regulation (EC) No 850/2004 of the European Parliament and of the Council of 29 April 2004 on persistent organic pollutants and amending Directive 79/117/EEC**

None.

**German Water hazard classes (Wassergefährdungsklassen)**

Chemical Name	Weight-%	Water hazard class (WGK)
Formamide	1-5	hazard class 1 - slightly hazardous to water

**Other International Inventories**

Chemical Name	EINECS (European Union)	ELINCS (European List of Notified Chemical Substances)	ENCS (Japan)	PICCS (Philippines)
Formamide	Listed	-	Listed	Listed

Chemical Name	AICS (Australia)	South Korea (KECL)	Canada (DSL)	NDSL
Formamide	Listed	Listed	Listed	-

**Chemical safety assessment**

No Chemical safety assessment has been carried out.



## SECTION 16: Other information

**Reason for revision** Update according to Commission Regulation (EU) No 830/2015  
**Revision number** 6  
**Revision date** 05-Mar-2020

### References

- ECHA: <http://echa.europa.eu/>
- TOXNET: <http://toxnet.nlm.nih.gov/>
- eChemPortal: <http://www.echemportal.org/>
- LOLI database: <https://www.chemadvisor.com/loli-database>

### Classification and procedure used to derive the classification for mixtures according to Regulation (EC) 1272/2008 [CLP]:

Carcinogenicity	Category 2	Calculation method
Reproductive Toxicity	Category 1B	Calculation method
Specific target organ toxicity - Repeated exposure	Category 2	Calculation method

"The above information was acquired by diligent search and/or investigation and the recommendations are based on prudent application of professional judgment. The information shall not be taken as being all inclusive and is to be used only as a guide. All materials and mixtures may present unknown hazards and should be used with caution. Since the Company cannot control the actual methods, volumes, or conditions of use, the Company shall not be held liable for any damages or losses resulting from the handling or from contact with the product as described herein. THE INFORMATION IN THIS SDS DOES NOT CONSTITUTE A WARRANTY, EXPRESSED OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR PURPOSE"

## Safety Data Sheet

Safety Data Sheet according to Regulation (EC) No. 1907/2006 (REACH)  
Classification according to Regulation (EC) No. 1272/2008 [CLP]

### SECTION 1: Identification of the substance/mixture and of the company/undertaking

#### Product identifier

**Product code** 10977035  
**Product name** DNASE/RNASE-FREE DISTILLED WATER

**Chemical Name** Not Applicable  
**REACH registration number** No registration number is given yet for this substance / substances in this mixture since the annual import quantity is less than one tonnage per annum or the transition period for its registration according to Article 23 of REACH has not yet expired.

#### Relevant identified uses of the substance or mixture and uses advised against

**Relevant identified uses** For research use only  
**Use Description Code** SU22 - Professional uses: Public domain (administration, education, entertainment, services, craftsmen), PROC15 - Use as laboratory reagent, PC21 - Laboratory chemicals, SU24 - Scientific research and development  
**Uses advised against** Not for consumer use.

#### Details of the supplier of the safety data sheet

##### Manufacturer / Supplier

LIFE TECHNOLOGIES EUROPE BV  
KWARTSWEG 2  
2665 NN BLEISWIJK  
NETHERLANDS  
31-(0)180 392 400  
Email: MSDS@lifetech.com

Life Technologies Limited  
3 Fountain Drive  
Inchinnan Business Park  
Paisley  
PA4 9RF, UK  
+44 (0)141 814 6100

**24 hour Emergency Response for Hazardous Materials** Within the USA + Canada: 1-800-424-9300 and  
**[or Dangerous Goods] Incident. Spill, Leak, Fire,** 1-703-527-3887  
**Exposure, or Accident. Call CHEMTREC** Outside the USA + Canada: 1-703-741-5970

#### Country Specific Emergency Number (if available):

CHEMTREC Ireland (Dublin) +(353)-19014670 (Greeting Language: English and Irish)  
CHEMTREC UK (London) +(44)-870-8200418 (Greeting Language: English)

**Revision date** 28-Oct-2019  
**Product code** 10977035

**Page** 1 / 9  
**Product name** DNASE/RNASE-FREE DISTILLED WATER

## SECTION 2: Hazards identification

### Classification of the substance or mixture

**Classification according to Regulation (EC) No. 1272/2008 [CLP]**

**Physical hazards**

Not Hazardous

**Health hazards**

Not Hazardous

**Environmental hazards**

Not Hazardous

**Additional information**

Not Applicable

**Label elements**

**Labelling according to Regulation (EC) No 1272/2008 [CLP]**

**Hazard pictograms**

None

**Signal Word**

None

**Hazard Statements**

Not Applicable

**EU Specific Hazard Statements**

Not Applicable

**Precautionary Statements**

**Prevention**

Not Applicable

**Response**

Not Applicable

**Storage**

Not Applicable

**Disposal**

Not Applicable

**Other hazards**

This mixture does not contain any substances that are assessed to be a PBT or a vPvB

### SECTION 3: Composition/information on ingredients

The product contains no substances which at their given concentration, are considered to be hazardous to health. We recommend handling all chemicals with caution.

### SECTION 4: First aid measures

#### Description of first aid measures

<b>Skin contact</b>	Rinse skin with water. Immediate medical attention is not required.
<b>Eye contact</b>	Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
<b>Ingestion</b>	Not expected to present a significant ingestion hazard under anticipated conditions of normal use. If you feel unwell, seek medical advice.
<b>Inhalation</b>	Not expected to be an inhalation hazard under anticipated conditions of normal use of this material. Consult a physician if necessary.
<b>Notes to Physician</b>	Treat symptomatically.

#### Most important symptoms and effects, both acute and delayed

Not Applicable

#### Indication of any immediate medical attention and special treatment needed

None.

### SECTION 5: Firefighting measures

#### Extinguishing media

<b>Suitable extinguishing media</b>	Water spray. Carbon dioxide (CO <sub>2</sub> ). Foam. Dry chemical.
<b>Unsuitable extinguishing media</b>	No information available.

#### Special hazards arising from the substance or mixture

Not known

#### Protective equipment and precautions for firefighters

Standard procedure for chemical fires.

## SECTION 6: Accidental release measures

### **Personal precautions, protective equipment and emergency procedures**

Ensure adequate ventilation  
Always wear recommended Personal Protective Equipment  
Use personal protection equipment  
See section 8 for more information

### **Environmental precautions**

No special environmental precautions required.

### **Methods and material for containment and cleaning up**

Soak up with inert absorbent material.

### **Reference to other sections**

See section 8 for more information.

## SECTION 7: Handling and storage

### **Precautions for safe handling**

Use personal protective equipment as required. No special handling advices are necessary.

### **Conditions for safe storage, including any incompatibilities**

Keep in a dry, cool and well-ventilated place. Keep in properly labelled containers.

### **Specific end use(s)**

For research use only.

## SECTION 8: Exposure controls/personal protection

### Control parameters

**Exposure Limits** Contains no substances with occupational exposure limit values.

**Engineering Measures** Ensure adequate ventilation, especially in confined areas.

### Exposure controls

#### **Personal protection equipment**

**Respiratory protection** In case of insufficient ventilation wear respirators and components tested and approved under appropriate government standards.

**Hand protection** Wear suitable gloves Glove material: Compatible chemical-resistant gloves.

**Eye protection** Tight sealing safety goggles.

**Skin and Body Protection** Wear suitable protective clothing.

**Hygiene Measures** Handle in accordance with good industrial hygiene and safety practice.

#### **Environmental exposure controls**

No special environmental precautions required.

## SECTION 9: Physical and chemical properties

### Information on basic physical and chemical properties

<b>Appearance</b>	liquid	
<b>Odour</b>	No data	
<b>Odour Threshold</b>	No data	
<b>Molecular Weight</b>	No data	
<b>pH</b>	6-8	
<b>Melting point / melting range</b>	°C No data	°F No data
<b>Boiling point / boiling range</b>	°C No data	°F No data
<b>Flash point</b>	°C No data	°F No data
<b>Autoignition Temperature</b>	°C No data	°F No data
<b>Decomposition temperature</b>	°C No data	°F No data
<b>Evaporation rate</b>	No data	
<b>Flammability (solid, gas)</b>	No data	
<b>Upper explosion limit</b>	No data	
<b>Lower explosion limit</b>	No data	
<b>Vapour Pressure</b>	No data	
<b>Vapour density</b>	No data	
<b>Relative density</b>	No data	
<b>Specific gravity</b>	No data	
<b>Solubility</b>	No data	
<b>Partition coefficient: n-octanol/water</b>	No data	
<b>Viscosity</b>	No data	
<b>Explosive properties</b>	No data	
<b>Oxidising properties</b>	No data	

### Other information

No data.

## SECTION 10: Stability and reactivity

<b>Reactivity</b>	None known.
<b>Chemical stability</b>	Stable under normal conditions.
<b>Possibility of hazardous reactions</b>	Hazardous reaction has not been reported.
<b>Conditions to avoid</b>	No information available.
<b>Incompatible materials</b>	No dangerous reaction known under conditions of normal use.
<b>Hazardous decomposition products</b>	No data available.

## SECTION 11: Toxicological information

### Information on toxicological effects

There is no evidence available indicating acute toxicity.

### Principal Routes of Exposure

**Skin corrosion/irritation** Data are conclusive but insufficient for classification

**Serious eye damage/irritation** Data are conclusive but insufficient for classification

**Respiratory or skin sensitisation** Data are conclusive but insufficient for classification

**Specific target organ toxicity (STOT) – single exposure** Data are conclusive but insufficient for classification

**Specific target organ toxicity (STOT) – repeated exposure** Data are conclusive but insufficient for classification

**Carcinogenicity** Data are conclusive but insufficient for classification

**Germ cell mutagenicity** Data are conclusive but insufficient for classification

**Reproductive Toxicity** Data are conclusive but insufficient for classification

**Aspiration Hazard** Data are conclusive but insufficient for classification

## SECTION 12: Ecological information

### Ecotoxicity

Contains no substances known to be hazardous to the environment or not degradable in waste water treatment plants.

**Persistence and degradability** No information available.

**Bioaccumulative potential** No information available.

### Results of PBT and vPvB assessment

This mixture does not contain any substances that are assessed to be a PBT or a vPvB.

### Other adverse effects

No information available.



## SECTION 13: Disposal considerations

### Waste treatment methods

The generation of waste should be avoided or minimized wherever possible. Empty containers or liners may retain some product residues. This material and its container must be disposed of in accordance with approved disposal technique. Disposal of this product, its solutions or of any by-products, shall comply with the requirements of all applicable local, regional or national/federal regulations.

## SECTION 14: Transport information

### IATA / ADR / DOT-US / IMDG

Not regulated in the meaning of transport regulations

<b>UN number</b>	Not Applicable
<b>UN proper shipping name</b>	Not Applicable
<b>Transport hazard class(es)</b>	Not Applicable
<b>Packing group</b>	Not Applicable

### Environmental hazards

Not Applicable

### Special precautions for user

Not Applicable

### Transport in bulk according to Annex II of MARPOL and the IBC Code

Not Applicable.

## SECTION 15: Regulatory information

### Safety, health and environmental regulations/legislation specific for the substance or mixture

#### Substances of Very High Concern

None.

#### Substance subject to authorisation per REACH Annex XIV

None

#### Restricted substances under EC 1907/2006, Annex XVII

None.

#### Substances listed under Annex I of Regulation (EC) No 689/2008

None.

#### Restricted substances under Annex V of Regulation (EC) No 689/2008

None.

Revision date 28-Oct-2019  
Product code 10977035

Page 8 / 9  
Product name DNASE/RNASE-FREE DISTILLED WATER

**Substances under Regulation (EC) No 850/2004 of the European Parliament and of the Council of 29 April 2004 on persistent organic pollutants and amending Directive 79/117/EEC**

None.

**German Water hazard classes (Wassergefährdungsklassen)**

Not classified.

**Other International Inventories**

No information available

**Chemical safety assessment**

No Chemical safety assessment has been carried out.

SECTION 16: Other information

<b>Reason for revision</b>	Update according to Commission Regulation (EU) No 830/2015
<b>Revision number</b>	3
<b>Revision date</b>	28-Oct-2019

**References**

- ECHA: <http://echa.europa.eu/>
- TOXNET: <http://toxnet.nlm.nih.gov/>
- eChemPortal: <http://www.echemportal.org/>
- LOLI database: <https://www.chemadvisor.com/loli-database>

**Classification and procedure used to derive the classification for mixtures according to Regulation (EC) 1272/2008 [CLP]:**

Not classified

"The above information was acquired by diligent search and/or investigation and the recommendations are based on prudent application of professional judgment. The information shall not be taken as being all inclusive and is to be used only as a guide. All materials and mixtures may present unknown hazards and should be used with caution. Since the Company cannot control the actual methods, volumes, or conditions of use, the Company shall not be held liable for any damages or losses resulting from the handling or from contact with the product as described herein. THE INFORMATION IN THIS SDS DOES NOT CONSTITUTE A WARRANTY, EXPRESSED OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR PURPOSE"

## Safety Data Sheet

(In accordance with COMMISSION REGULATION (EU) No 830/2015)

### SECTION 1: Identification of the substance/mixture and of the company/undertaking

#### Product identifier

**Product code** C5192  
**Product name** Wash buffer 1 (conc.)

**Chemical Name** Not Applicable  
**REACH registration number** No registration number is given yet for this substance / substances in this mixture since the annual import quantity is less than one tonnage per annum or the transition period for its registration according to Article 23 of REACH has not yet expired.

#### Relevant identified uses of the substance or mixture and uses advised against

**Relevant identified uses** For research use only. Not for use in diagnostic procedures  
**Use Description Code** SU22 - Professional uses: Public domain (administration, education, entertainment, services, craftsmen), PROC15 - Use as laboratory reagent, PC21 - Laboratory chemicals, SU24 - Scientific research and development  
**Uses advised against** Not for consumer use.

#### Details of the supplier of the safety data sheet

##### Manufacturer/Supplier

Thermo Fisher Scientific Baltics UAB  
V.Graiciuno 8  
LT-02241 Vilnius  
Lithuania  
Tel.: +370 5 2602131  
Fax: +370 5 2602142

Life Technologies Limited  
3 Fountain Drive  
Inchinnan Business Park  
Paisley  
PA4 9RF, UK  
+44 (0)141 814 6100

LIFE TECHNOLOGIES EUROPE BV  
KWARTSWEG 2  
2665 NN BLEISWIJK  
NETHERLANDS  
31-(0)180 392 400  
Email: MSDS@lifetech.com

**24 hour Emergency Response for Hazardous Materials** Within the USA + Canada: 1-800-424-9300 and +1  
**[or Dangerous Goods] Incident. Spill, Leak, Fire,** 703-527-3887  
**Exposure, or Accident. Call CHEMTREC** Outside the USA + Canada: +1 703-741-5970

#### Country Specific Emergency Number (if available):

**Revision date** 29-Mar-2019  
**Product code** C5192

**Page** 1 / 9  
**Product name** Wash buffer 1 (conc.)

SECTION 2: Hazards identification

**Classification of the substance or mixture**

**Classification according to Regulation (EC) No 1272/2008 [CLP]**

**Physical hazards**

Not Hazardous

**Health hazards**

Acute oral toxicity	Category 4
Skin corrosion/irritation	Category 2
Serious eye damage/eye irritation	Category 2

**Environmental hazards**

Not Hazardous

**Additional information**

No information available

**Label elements**

**Labelling according to Regulation (EC) No 1272/2008 [CLP]**

**Hazard pictograms**



**Signal Word**

Warning

**Hazard Statements**

- H302 - Harmful if swallowed
- H315 - Causes skin irritation
- H319 - Causes serious eye irritation

**Precautionary Statements**

**Prevention**

- P264 - Wash hands thoroughly after handling
- P280 - Wear protective gloves/protective clothing/eye protection/face protection
- P270 - Do not eat, drink or smoke when using this product

**Response**

- P301 + P310 - IF SWALLOWED: Immediately call a POISON CENTER or doctor/physician
- P302 + P352 - IF ON SKIN: Wash with plenty of soap and water
- P305 + P351 + P338 - IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if

present and easy to do. Continue rinsing  
P330 - Rinse mouth  
P332 + P313 - If skin irritation occurs: Get medical advice/attention

#### Storage

Not Applicable

#### Disposal

P501 - Dispose of contents/ container to an approved waste disposal plant

#### Other hazards

Not Applicable

### SECTION 3: Composition/information on ingredients

We recommend handling all chemicals with caution.

Component	CAS No	EINECS-No.	Weight-%	REACH registration number	Classification according to Regulation (EC) No 1272/2008 [CLP]
Guanidine hydrochloride 50-01-1 ( 40 - 70 )	50-01-1	-	40 - 70	-	Acute Tox. 4 - H302 Acute Tox. 4 - H332 Skin Irrit. 2 - H315 Eye Irrit. 2 - H319

### SECTION 4: First aid measures

#### Description of first aid measures

##### Skin contact

Wash off immediately with plenty of water for at least 15 minutes. Remove and wash contaminated clothing and gloves, including the inside, before re-use. Immediate medical attention is required.

##### Eye contact

Rinse immediately with plenty of water, also under the eyelids, for at least 15 minutes. Immediate medical attention is required.

##### Ingestion

Never give anything by mouth to an unconscious person. Do not induce vomiting without medical advice. Get medical attention if symptoms occur.

##### Inhalation

Remove to fresh air. If not breathing, give artificial respiration. If symptoms persist, call a doctor.

##### Notes to Physician

Treat symptomatically.

#### Most important symptoms and effects, both acute and delayed

#### Indication of any immediate medical attention and special treatment needed

IF INHALED: Remove person to fresh air and keep comfortable for breathing. If skin irritation occurs: Get medical advice/ attention. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

## SECTION 5: Firefighting measures

### Extinguishing media

**Suitable extinguishing media**

Water spray. Carbon dioxide (CO<sub>2</sub>). Foam. Dry chemical.

**Unsuitable extinguishing media**

No information available.

### Special hazards arising from the substance or mixture

Not known.

### Advice for fire-fighters

Wear self-contained breathing apparatus and protective suit.

## SECTION 6: Accidental release measures

### Personal precautions, protective equipment and emergency procedures

Ensure adequate ventilation

Avoid contact with skin, eyes or clothing

Use personal protection equipment

See Section 8 for more detail.

### Environmental precautions

Avoid discharge into drains and waterways whenever possible.

### Methods and material for containment and cleaning up

Soak up with inert absorbent material.

### Reference to other sections

See section 8 for more information.

## SECTION 7: Handling and storage

### Precautions for safe handling

Always wear recommended Personal Protective Equipment. See Section 8 for more detail. Do not get in eyes, on skin, or on clothing. Do not ingest. If during normal use the material presents a respiratory hazard, use adequate ventilation and/or wear appropriate respirator.

### Conditions for safe storage, including any incompatibilities

Keep in a dry, cool and well-ventilated place. Keep in properly labelled containers. Store in accordance with local regulations.

### Specific end use(s)

For research use only. Not for use in diagnostic procedures.

## SECTION 8: Exposure controls/personal protection

### Control parameters

Chemical Name	EU OEL (TWA)	EU OEL (STEL)	EU Skin Notation
Guanidine hydrochloride 50-01-1	None	None	None

Chemical Name	Austria	Belgium (TWA)	Denmark (TWA)	Finland OEL (TWA)
Guanidine hydrochloride 50-01-1	None	None	None	None

Chemical Name	France OEL (VME)	Germany OEL (TWA)	Ireland (TWA)	Italy OEL (TWA)
Guanidine hydrochloride 50-01-1	None	None	None	None

Chemical Name	Sweden - Occupational Exposure Limits - TLVs (LLVs)	Netherlands OEL (MAC)	Spain OEL (TWA)	United Kingdom
Guanidine hydrochloride 50-01-1	None	None	None	None

Chemical Name	European Union	France OEL (VME)	Germany OEL (TWA)
Guanidine hydrochloride 50-01-1	None	None	None

Chemical Name	Italy OEL (TWA)	Portugal	Netherlands OEL (MAC)	Finland OEL (TWA)
Guanidine hydrochloride 50-01-1	None	None	None	None

Chemical Name	Austria	Denmark	Poland	Switzerland
Guanidine hydrochloride 50-01-1	None	None	None	None

Chemical Name	Ireland	Norway	Lithuania OEL (TWA)	Spain OEL (TWA)
Guanidine hydrochloride 50-01-1	None	None	None	None

**Engineering measures**      Ensure adequate ventilation, especially in confined areas.

### Exposure controls

#### Personal protection equipment

**Respiratory protection**      In case of insufficient ventilation wear respirators and components tested and approved under appropriate government standards.

**Hand protection**      Glove material: Nitrile rubber. with thickness (mm) :.5. Break through time. (hours) :.1.  
Recommended glove type has not been tested for use with this product.  
Information is based on professional knowledge.

**Eye protection**      Tight sealing safety goggles.

**Skin and Body Protection**      Wear laboratory coat for body protection.

**Hygiene Measures**      Handle in accordance with good industrial hygiene and safety practice.

**Environmental exposure controls**  
No special environmental precautions required.

## SECTION 9: Physical and chemical properties

### Information on basic physical and chemical properties

<b>Appearance</b>	Colourless, liquid	
<b>Odour</b>	odourless	
<b>Odour Threshold</b>	Mixture has not been tested	
<b>pH</b>	Mixture has not been tested	
<b>Melting point / melting range</b>	°C Mixture has not been tested	°F Mixture has not been tested
<b>Boiling point / boiling range</b>	°C Mixture has not been tested	°F Mixture has not been tested
<b>Flash point</b>	°C Mixture has not been tested	°F Mixture has not been tested
<b>Autoignition Temperature</b>	°C Mixture has not been tested	°F Mixture has not been tested
<b>Decomposition temperature</b>	°C Mixture has not been tested	°F Mixture has not been tested
<b>Evaporation rate</b>	No data	
<b>Flammability (solid, gas)</b>	Not Applicable	
<b>Upper explosion limit</b>	Mixture has not been tested	
<b>Lower explosion limit</b>	Mixture has not been tested	
<b>Vapour Pressure</b>	Mixture has not been tested	
<b>Relative density</b>	Mixture has not been tested	
<b>Specific gravity</b>	No data	
<b>Solubility</b>	Soluble in water	
<b>Partition coefficient: n-octanol/water</b>	No data	
<b>Viscosity</b>	Mixture has not been tested	
<b>Explosive properties</b>	Mixture has not been tested	
<b>Oxidising properties</b>	Mixture has not been tested	

### Other information

No data.

## SECTION 10: Stability and reactivity

<b>Reactivity</b>	None known.
<b>Chemical stability</b>	Stable under normal conditions.
<b>Possibility of hazardous reactions</b>	Hazardous reaction has not been reported.
<b>Conditions to avoid</b>	Contact with acids or bleach liberates toxic gases. DO NOT ADD acids or bleach to any liquid wastes containing this product.
<b>Incompatible materials</b>	No dangerous reaction known under conditions of normal use.
<b>Hazardous decomposition products</b>	Carbon oxides. Nitrogen oxides (NOx). Sulphur oxides.



## SECTION 11: Toxicological information

### Information on toxicological effects

Chemical Name	LD50 (oral, rat/mouse)	LD50 (dermal, rat/rabbit)	LC50 (inhalation, rat/mouse)
Guanidine hydrochloride	= 475 mg/kg (Rat)	Skin (rabbit) 500 mg/24h SEVERE	No data available

### Principal Routes of Exposure,

<b>Skin corrosion/irritation</b>	Skin irritation
<b>Serious eye damage/irritation</b>	Irritating to eyes
<b>Respiratory or skin sensitisation</b>	Conclusive but not sufficient for classification
<b>Specific target organ toxicity (STOT) – single exposure</b>	May cause respiratory irritation
<b>Specific target organ toxicity (STOT) – repeated exposure</b>	Conclusive but not sufficient for classification
<b>Carcinogenicity</b>	Conclusive but not sufficient for classification
<b>Germ cell mutagenicity</b>	Conclusive but not sufficient for classification
<b>Reproductive Toxicity</b>	Conclusive but not sufficient for classification
<b>Aspiration Hazard</b>	Conclusive but not sufficient for classification

## SECTION 12: Ecological information

### Toxicity

The environmental impact of this product has not been fully investigated.

Chemical Name	Freshwater Algae Data	Water Flea Data	Freshwater Fish Species Data	Microtox Data	log Pow
Guanidine hydrochloride	No data available	No data available	No data available	No data available	logPow<=-1.7

**Persistence and degradability** No information available.

**Bioaccumulative potential** No information available.

### Results of PBT and vPvB assessment

No information available.

### Other adverse effects

No information available.

## SECTION 13: Disposal considerations

### Waste treatment methods

The generation of waste should be avoided or minimized wherever possible. Empty containers or liners may retain some product residues. This material and its container must be disposed of in accordance with approved disposal technique. Disposal of this product, its solutions or of any by-products, shall comply with the requirements of all applicable local, regional or national/federal regulations.

## SECTION 14: Transport information

### IATA / ADR / DOT-US / IMDG

Not regulated in the meaning of transport regulations

<b>UN number</b>	Not Applicable
<b>UN proper shipping name</b>	Not Applicable
<b>Transport hazard class(es)</b>	Not Applicable
<b>Packing group</b>	Not Applicable

### Environmental hazards

Not Applicable

### Special precautions for user

Not Applicable

### Transport in bulk according to Annex II of MARPOL and the IBC Code

Not Applicable.

## SECTION 15: Regulatory information

### Safety, health and environmental regulations/legislation specific for the substance or mixture

#### Substances of Very High Concern

None.

#### EU REACH (1907/2006) - Annex XIV - Substances Subject to Authorization

None

#### Restricted substances under EC 1907/2006, Annex XVII

None.

#### Substances listed under Annex I of Regulation (EC) No 689/2008

None.

#### Restricted substances under Annex V of Regulation (EC) No 689/2008

None.

#### Substances under Regulation (EC) No 850/2004 of the European Parliament and of the Council of 29 April 2004 on persistent organic pollutants and amending Directive 79/117/EEC

Revision date 29-Mar-2019  
Product code C5192

Page 8 / 9  
Product name Wash buffer 1 (conc.)

None.

### German Water hazard classes (Wassergefährdungsklassen)

Chemical Name	Weight-%	Germany - Water Classification - Substances According to AwSV Classified By or Based on the VwVwS
Guanidine hydrochloride	40 - 70	hazard class 1 - slightly hazardous to water

### Other International Inventories

Chemical Name	EINECS (European Union)	ELINCS (European List of Notified Chemical Substances)	ENCS (Japan)	PICCS (Philippines)
Guanidine hydrochloride	Listed	-	Listed	Listed

Chemical Name	AICS (Australia)	South Korea (KECL)	Canada (DSL)	NDSL
Guanidine hydrochloride	Listed	Listed	Listed	-

### Chemical safety assessment

No Chemical safety assessment has been carried out.

### SECTION 16: Other information

<b>Reason for revision</b>	Update according to Commission Regulation (EU) No 830/2015
<b>Revision number</b>	1
<b>Revision date</b>	29-Mar-2019

### References

- ECHA: <http://echa.europa.eu/>
- TOXNET: <http://toxnet.nlm.nih.gov/>
- eChemPortal: <http://www.echemportal.org/>
- LOLI database: <https://www.chemadvisor.com/loli-database>

### Classification and procedure used to derive the classification for mixtures according to Regulation (EC) 1272/2008 [CLP]:

Acute oral toxicity	Category 4	Calculation method
Skin corrosion/irritation	Category 2	Calculation method
Serious eye damage/eye irritation	Category 2	Calculation method
Calculation method		Calculation method

"The above information was acquired by diligent search and/or investigation and the recommendations are based on prudent application of professional judgment. The information shall not be taken as being all inclusive and is to be used only as a guide. All materials and mixtures may present unknown hazards and should be used with caution. Since the Company cannot control the actual methods, volumes, or conditions of use, the Company shall not be held liable for any damages or losses resulting from the handling or from contact with the product as described herein. THE INFORMATION IN THIS SDS DOES NOT CONSTITUTE A WARRANTY, EXPRESSED OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR PURPOSE"

Revision date 29-Mar-2019  
Product code C5192

Page 9 / 9  
Product name Wash buffer 1 (conc.)

# Safety Data Sheet

(In accordance with COMMISSION REGULATION (EU) No 830/2015)

## SECTION 1: Identification of the substance/mixture and of the company/undertaking

### Product identifier

**Product code** 21875034  
**Product name** RPMI 1640 Medium (1X), liquid  
**Chemical Name** Not Applicable  
**REACH registration number** No registration number is given yet for this substance / substances in this mixture since the annual import quantity is less than one tonnage per annum or the transition period for its registration according to Article 23 of REACH has not yet expired.

### Relevant identified uses of the substance or mixture and uses advised against

**Relevant identified uses** For in vitro diagnostic use  
**Use Description Code** SU22 - Professional uses: Public domain (administration, education, entertainment, services, craftsmen), PROC15 - Use as laboratory reagent, PC21 - Laboratory chemicals, SU20 - Health services  
**Uses advised against** Not for consumer use.

### Details of the supplier of the safety data sheet

#### Manufacturer/Supplier

LIFE TECHNOLOGIES EUROPE BV  
KWARTSWEG 2  
2665 NN BLEISWIJK  
NETHERLANDS  
31-(0)180 392 400  
Email: MSDS@lifetech.com

Life Technologies Limited  
3 Fountain Drive  
Inchinnan Business Park  
Paisley  
PA4 9RF, UK  
+44 (0)141 814 6100

**24 hour Emergency Response for Hazardous Materials** Within the USA + Canada: 1-800-424-9300 and +1  
**[or Dangerous Goods] Incident. Spill, Leak, Fire,** 703-527-3887  
**Exposure, or Accident. Call CHEMTREC** Outside the USA + Canada: +1 703-741-5970

#### Country Specific Emergency Number (if available):

CHEMTREC Ireland (Dublin) +(353)-19014670 (Greeting Language: English and Irish)  
CHEMTREC UK (London) +(44)-870-8200418 (Greeting Language: English)

## SECTION 2: Hazards identification

### Classification of the substance or mixture

#### Classification according to Regulation (EC) No 1272/2008 [CLP]

##### **Physical hazards**

Not Hazardous

##### **Health hazards**

Not Hazardous

##### **Environmental hazards**

Not Hazardous

##### **Additional information**

Not Applicable

### Label elements

#### Labelling according to Regulation (EC) No 1272/2008 [CLP]

##### **Hazard pictograms**

None

##### **Signal Word**

None

##### **Hazard Statements**

Not Applicable

##### **EU Specific Hazard Statements**

Not Applicable

##### **Precautionary Statements**

##### **Prevention**

Not Applicable

##### **Response**

Not Applicable

##### **Storage**

Not Applicable

##### **Disposal**

Not Applicable

##### **Other hazards**

This mixture does not contain any substances that are assessed to be a PBT or a vPvB

### SECTION 3: Composition/information on ingredients

The product contains no substances which at their given concentration, are considered to be hazardous to health. We recommend handling all chemicals with caution.

### SECTION 4: First aid measures

#### Description of first aid measures

<b>Skin contact</b>	Rinse with plenty of water. Immediate medical attention is not required.
<b>Eye contact</b>	Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
<b>Ingestion</b>	Not expected to present a significant ingestion hazard under anticipated conditions of normal use. If you feel unwell, seek medical advice.
<b>Inhalation</b>	Not expected to be an inhalation hazard under anticipated conditions of normal use of this material. Consult a physician if necessary.
<b>Notes to Physician</b>	Treat symptomatically.

#### Most important symptoms and effects, both acute and delayed

Not Applicable

#### Indication of any immediate medical attention and special treatment needed

None.

### SECTION 5: Firefighting measures

#### Extinguishing media

<b>Suitable extinguishing media</b>	Water spray. Carbon dioxide (CO <sub>2</sub> ). Foam. Dry chemical.
<b>Unsuitable extinguishing media</b>	No information available.

#### Special hazards arising from the substance or mixture

Not known.

#### Advice for fire-fighters

Standard procedure for chemical fires.

## SECTION 6: Accidental release measures

### **Personal precautions, protective equipment and emergency procedures**

Ensure adequate ventilation  
Always wear recommended Personal Protective Equipment  
Use personal protection equipment  
See Section 8 for more detail.

### **Environmental precautions**

No special environmental precautions required.

### **Methods and material for containment and cleaning up**

Soak up with inert absorbent material.

### **Reference to other sections**

See section 8 for more information.

## SECTION 7: Handling and storage

### **Precautions for safe handling**

Use personal protective equipment as required. No special handling advices are necessary.

### **Conditions for safe storage, including any incompatibilities**

Keep in a dry, cool and well-ventilated place. Keep in properly labelled containers.

### **Specific end use(s)**

For in vitro diagnostic use.

## SECTION 8: Exposure controls/personal protection

### Control parameters

**Exposure Limits** Contains no substances with occupational exposure limit values.

**Engineering measures** Ensure adequate ventilation, especially in confined areas.

### Exposure controls

#### **Personal protection equipment**

**Respiratory protection** In case of insufficient ventilation wear respirators and components tested and approved under appropriate government standards.

**Hand protection** Wear suitable gloves. Glove material: Compatible chemical-resistant gloves.

**Eye protection** Tight sealing safety goggles.

**Skin and Body Protection** Wear suitable protective clothing.

**Hygiene Measures** Handle in accordance with good industrial hygiene and safety practice.

#### **Environmental exposure controls**

No special environmental precautions required.



## SECTION 9: Physical and chemical properties

### Information on basic physical and chemical properties

<b>Appearance</b>	Solution	
<b>Odour</b>	No data available	
<b>pH</b>	7.2	
<b>Melting point / melting range</b>	°C Mixture has not been tested	°F Mixture has not been tested
<b>Boiling point / boiling range</b>	°C Mixture has not been tested	°F Mixture has not been tested
<b>Flash point</b>	°C Mixture has not been tested	°F Mixture has not been tested
<b>Autoignition Temperature</b>	°C Mixture has not been tested	°F Mixture has not been tested
<b>Decomposition temperature</b>	°C Mixture has not been tested	°F Mixture has not been tested
<b>Evaporation rate</b>	No data	
<b>Flammability (solid, gas)</b>	No data available	
<b>Upper explosion limit</b>	Mixture has not been tested	
<b>Lower explosion limit</b>	Mixture has not been tested	
<b>Vapour Pressure</b>	Mixture has not been tested	
<b>Relative density</b>	Mixture has not been tested	
<b>Specific gravity</b>	No data	
<b>Solubility</b>	No data available	
<b>Partition coefficient: n-octanol/water</b>	No data available	
<b>Viscosity</b>	No data	
<b>Explosive properties</b>	Mixture has not been tested	
<b>Oxidising properties</b>	No data	

### Other information

No data available.

## SECTION 10: Stability and reactivity

<b>Reactivity</b>	None known.
<b>Chemical stability</b>	Stable under normal conditions.
<b>Possibility of hazardous reactions</b>	Hazardous reaction has not been reported.
<b>Conditions to avoid</b>	No information available.
<b>Incompatible materials</b>	No dangerous reaction known under conditions of normal use.
<b>Hazardous decomposition products</b>	No data available.

## SECTION 11: Toxicological information

### Information on toxicological effects

There is no evidence available indicating acute toxicity.

### Principal Routes of Exposure,

<b>Skin corrosion/irritation</b>	Conclusive but not sufficient for classification
<b>Serious eye damage/irritation</b>	Conclusive but not sufficient for classification
<b>Respiratory or skin sensitisation</b>	Conclusive but not sufficient for classification
<b>Specific target organ toxicity (STOT) – single exposure</b>	Conclusive but not sufficient for classification
<b>Specific target organ toxicity (STOT) – repeated exposure</b>	Conclusive but not sufficient for classification
<b>Carcinogenicity</b>	Conclusive but not sufficient for classification
<b>Germ cell mutagenicity</b>	Conclusive but not sufficient for classification
<b>Reproductive Toxicity</b>	Conclusive but not sufficient for classification
<b>Aspiration Hazard</b>	Conclusive but not sufficient for classification

## SECTION 12: Ecological information

### Toxicity

Contains no substances known to be hazardous to the environment or not degradable in waste water treatment plants.

**Persistence and degradability** No information available.

**Bioaccumulative potential** No information available.

### Results of PBT and vPvB assessment

This mixture does not contain any substances that are assessed to be a PBT or a vPvB.

### Other adverse effects

No information available.

## SECTION 13: Disposal considerations

### Waste treatment methods

The generation of waste should be avoided or minimized wherever possible. Empty containers or liners may retain some product residues. This material and its container must be disposed of in accordance with approved disposal technique. Disposal of this product, its solutions or of any by-products, shall comply with the requirements of all applicable local, regional or national/federal regulations.

## SECTION 14: Transport information

### IATA / ADR / DOT-US / IMDG

Not regulated in the meaning of transport regulations

<b>UN number</b>	Not Applicable
<b>UN proper shipping name</b>	Not Applicable
<b>Transport hazard class(es)</b>	Not Applicable
<b>Packing group</b>	Not Applicable

### Environmental hazards

Not Applicable

### Special precautions for user

Not Applicable

### Transport in bulk according to Annex II of MARPOL and the IBC Code

Not Applicable.

## SECTION 15: Regulatory information

### Safety, health and environmental regulations/legislation specific for the substance or mixture

#### Substances of Very High Concern

None.

#### EU REACH (1907/2006) - Annex XIV - Substances Subject to Authorization

None

#### Restricted substances under EC 1907/2006, Annex XVII

None.

#### Substances listed under Annex I of Regulation (EC) No 689/2008

None.

#### Restricted substances under Annex V of Regulation (EC) No 689/2008

None.

#### Substances under Regulation (EC) No 850/2004 of the European Parliament and of the Council of 29 April 2004 on persistent organic pollutants and amending Directive 79/117/EEC

None.

#### German Water hazard classes (Wassergefährdungsklassen)

Not classified.

### Other International Inventories

Revision date 28-Mar-2019  
Product code 21875034

Page 8 / 9  
Product name RPM1 1640 Medium (1X), liquid

No information available

### Chemical safety assessment

No Chemical safety assessment has been carried out.

## SECTION 16: Other information

<b>Reason for revision</b>	Update according to Commission Regulation (EU) No 830/2015
<b>Revision number</b>	3
<b>Revision date</b>	28-Mar-2019

### References

- ECHA: <http://echa.europa.eu/>
- TOXNET: <http://toxnet.nlm.nih.gov/>
- eChemPortal: <http://www.echemportal.org/>
- LOLI database: <https://www.chemadvisor.com/loli-database>

### Classification and procedure used to derive the classification for mixtures according to Regulation (EC) 1272/2008 [CLP]:

Not classified

"The above information was acquired by diligent search and/or investigation and the recommendations are based on prudent application of professional judgment. The information shall not be taken as being all inclusive and is to be used only as a guide. All materials and mixtures may present unknown hazards and should be used with caution. Since the Company cannot control the actual methods, volumes, or conditions of use, the Company shall not be held liable for any damages or losses resulting from the handling or from contact with the product as described herein. THE INFORMATION IN THIS SDS DOES NOT CONSTITUTE A WARRANTY, EXPRESSED OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR PURPOSE"

---


<b>Revision date</b>	28-Mar-2019
<b>Product code</b>	21875034

<b>Product name</b>	RPMI 1640 Medium (1X), liquid
---------------------	-------------------------------

# RPMI 1640 Medium

For various human clinical samples

Pub. No. MAN0018935 Rev. 1.0

 **WARNING!** Read the Safety Data Sheets (SDSs) and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves. Safety Data Sheets (SDSs) are available from [thermofisher.com/support](https://www.thermofisher.com/support).

## Intended Use

For *in vitro* diagnostic use

The isolation of human viruses from clinical samples using cell culture remains necessary because it is the only technique capable of providing a viable isolate that can be used for antiviral susceptibility testing. An additional advantage is that in contrast to most antigen and nucleic acid detection methods, viral culture allows detection of multiple viruses, not all of which may have been suspected at the time diagnostic culture was requested.

RPMI 1640 cell culture media products are for professional use. They are used in medical laboratories by personnel who have received specialized education and training with regard to procedures utilizing In Vitro Diagnostic products. IVD products of this type are not intended as sole determinant in a diagnostic situation. Test results are interpreted by a healthcare professional as part of the clinical management of a patient.

## Principle and explanation of procedure

RPMI is a commonly used cell culture media for diagnostic virology (1, 2). RPMI 1640 Medium was originally developed to culture human leukemic cells in suspension and as a monolayer. Roswell Park Memorial Institute (RPMI) 1640 Medium has since been found suitable for a variety of mammalian cells, including HeLa, Jurkat, MCF-7, PC12, PBMC, astrocytes, and carcinomas.

RPMI 1640 Medium is unique from other media because it contains the reducing agent glutathione and high concentrations of vitamins. RPMI 1640 Medium contains biotin, vitamin B12, and PABA, which are not found in Eagle's Minimal Essential Medium or Dulbecco's Modified Eagle Medium. In addition, the vitamins inositol and choline are present in very high concentrations. RPMI 1640 Medium contains no proteins, lipids, or growth factors. Therefore, RPMI 1640 Medium requires supplementation, commonly with 10% Fetal Bovine Serum (FBS). RPMI 1640 Medium uses a sodium bicarbonate buffer system (2.0 g/L), and therefore requires a 5–10% CO<sub>2</sub> environment to maintain physiological pH.

## Contents and storage

All quality control testing results are reported on lot-specific Certificate of Analysis available on our website: [thermofisher.com](https://www.thermofisher.com).


Product	Cat. No.	Storage	Shelf life <sup>[1]</sup>
RPMI 1640 (1X) [+] L-Glutamine [-] Phenol Red	11835030 <sup>[2]</sup> 11835055 <sup>[2]</sup> 11835063 <sup>[3]</sup>	2°C to 8°C Protect from light	12 months
RPMI Medium (1X) 1640 [+] L-Glutamine	11875085 <sup>[4]</sup> 11875093 <sup>[4]</sup> 11875101 <sup>[2]</sup> 11875119 <sup>[4]</sup> 11875127 <sup>[2]</sup> 11875135 <sup>[4]</sup> 11875168 <sup>[2]</sup> 11875176 <sup>[2]</sup>	2°C to 8°C Protect from light	12 months





For In Vitro Diagnostic Use.

Product	Cat. No.	Storage	Shelf life <sup>[1]</sup>
RPMI Medium (1X) 1640 [-] L-Glutamine	21870076 <sup>[4]</sup> 21870084 <sup>[4]</sup> 21870092 <sup>[4]</sup> 21870100 <sup>[4]</sup>	2°C to 8°C Protect from light	24 months
RPMI 1640 W/GLUT (1X) (CE)	21875034 <sup>[3]</sup> 21875042 <sup>[3]</sup> 21875059 <sup>[3]</sup>	2°C to 8°C Protect from light	12 months
RPMI 1640 (1X) [+] L-Glutamine [+] HEPES	22400071 <sup>[2]</sup> 22400089 <sup>[4]</sup> 22400097 <sup>[2]</sup> 22400105 <sup>[4]</sup> 22400121 <sup>[2]</sup> 22400197 <sup>[2]</sup>	2°C to 8°C Protect from light	12 months
RPMI 1640 W/O L-GLUTAMINE (CE)	31870017 <sup>[3]</sup> 31870025 <sup>[3]</sup>	2°C to 8°C Protect from light	12 months
RPMI 1640 W/O PHENOL RED	32404014 <sup>[3]</sup>	2°C to 8°C Protect from light	12 months
RPMI 1640 W/25MMHEPES W/OL-GLUT	42401018 <sup>[3]</sup> 42402016 <sup>[3]</sup>	2°C to 8°C Protect from light	12 months
RPMI MEDIUM 1640 W/HEPES (CE)	52400017 <sup>[3]</sup> 52400025 <sup>[3]</sup>	2°C to 8°C Protect from light	12 months
RPMI 1640 W/GLUTAMAX-I (1X)	61870010 <sup>[3]</sup> 61870036 <sup>[2]</sup> 61870127 <sup>[2]</sup> 61870143 <sup>[2]</sup> 61870150 <sup>[2]</sup>	2°C to 8°C Protect from light	12 months
RPMI 1640 W/HEPES W/GLUTAMAX-I	72400013 <sup>[3]</sup> 72400021 <sup>[3]</sup>	2°C to 8°C Protect from light	12 months
RPMI (1X) + GlutaMAX -I	72400047 <sup>[2]</sup> 72400120 <sup>[2]</sup> 72400146 <sup>[2]</sup> 72400153 <sup>[2]</sup>	2°C to 8°C Protect from light	12 months

<sup>[1]</sup> Shelf life is determined from Date of Manufacture. Do not use beyond the labelled expiration date.

<sup>[2]</sup>  Manufacturer: Life Technologies Corporation | 3175 Staley Road | Grand Island, NY 14072

<sup>[3]</sup>  Manufacturer: Life Technologies™ Ltd. | 3 Fountain Drive, Inchinnan Business Park | Paisley PA49RF, Scotland, United Kingdom | Tel: +44 (0)141 81416305

<sup>[4]</sup>  Dual manufactured.

## Precautions

Do not use the product if packaging, including bottles and vials, have been compromised and/or show evidence of microbial contamination, cloudy appearance, discoloration, drying, cracking, or other signs of deterioration.



**CAUTION!** Human samples are potentially biohazardous. Follow standard precautions for handling, storage and disposal.



**WARNING!** Do not use for injection or infusion! Please report any serious incidents in relation to the device to the manufacturer and the Competent Authority of the EU Member State in which the user and/or patient is established.

- Once opened, use RPMI 1640 Medium within 14 days for maximal growth performance.
- Avoid repeated warming/cooling and prolonged exposure to light.
- Do not use beyond labeled expiration date.
- All solutions that come into contact with clinical samples must be sterile. Always use proper aseptic techniques and work inside a laminar flow hood. Consult our **Gibco Cell Culture Basics** for aseptic handling.

## Test protocol

There is no single type of cell culture that can support the growth of all medically relevant viruses. As such, virology laboratories must maintain several different cell culture types. The choice of cell line used for a specific specimen is determined by the information communicated from the ordering physician to the laboratory and by knowledge of the specimens usually isolated from a given specimen type.

Ready to-use commercial cell culture media undergoes strict quality control to ensure sterility, but may become contaminated while handling. Follow the below guidelines for sterile handling to avoid contamination.

- Always wipe your hand and work area with 70% ethanol.
- Wipe the outside of the containers, flasks, plates, and dishes with 70% ethanol before placing them in the cell culture hood.
- Avoid pouring media and reagents directly from bottles or flasks.
- Use sterile pipette tips and pipettes to work with liquids, and use each pipette tip only once to avoid cross-contamination. Do not unwrap sterile pipettes until they are ready to be used. Keep pipettes and tips within the clean work area.
- Do not talk while performing sterile procedures and perform your cell culture as rapidly as possible to minimize contamination.

## Quality control

Standard evaluations for cell culture media are pH, osmolality, endotoxins and sterility testing for liquid products. All quality control testing results are reported on lot specific Certificate of Analysis available on our website: [thermofisher.com](http://thermofisher.com).

## Related products

Product	Source
Gentamicin 50 mg/mL	15750078
Gibco Amphotericin B	15290018
Penicillin Streptomycin 10,000 U/mL	15140122
PBS, pH 7.4	10010031
Phytohemagglutinin, M form (PHA-M)	10576015
FBS	16000044

## References

1. Winn, W. C., & Koneman, E. W. (2006). Koneman's color atlas and textbook of diagnostic microbiology (6th ed.). Philadelphia: Lippincott Williams & Wilkins.
2. WHO Guidelines on the Establishment of Virology Laboratories in Developing Countries, 2008.
3. Griffith, B P. "Principles of laboratory isolation and identification of the human immunodeficiency virus (HIV)" Yale journal of biology and medicine vol. 60,6 (1987): 575-87.
4. Krowicka, Halina et al. "Use of tissue culture cell lines to evaluate HIV antiviral resistance" AIDS research and human retroviruses vol. 24,7 (2008): 957-67.

## Labeling symbols

The symbols present on the product label are explained in the following table.

	MANUFACTURER		USE BY
	IN VITRO DIAGNOSTIC MEDICAL DEVICE		CONSULT INSTRUCTIONS FOR USE
	CATALOG NUMBER		CAUTION, CONSULT ACCOMPANYING DOCUMENTS
	BATCH CODE		UPPER AND LOWER LIMITS OF TEMPERATURE
	Sterilized using aseptic processing technique		PROTECT FROM LIGHT
	European Mark of Conformity		AUTHORISED REPRESENTATIVE IN THE EUROPEAN COMMUNITY

## Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and

Conditions of Sale at [www.thermofisher.com/us/en/home/global/terms-and-conditions.html](http://www.thermofisher.com/us/en/home/global/terms-and-conditions.html). If you have any questions, please contact Life Technologies at [www.thermofisher.com/support](http://www.thermofisher.com/support).



**Manufacturer:**  
Life Technologies Corporation |  
3175 Staley Road |  
Grand Island, NY 14072



European Regulatory Affairs  
Life Technologies Europe B.V.  
Kwartsweg 2, 2665 NN Bleiswijk  
The Netherlands  
Tel: +31 (0) 10 714 5000



**Manufacturer:**  
Life Technologies™ Ltd. |  
3 Fountain Drive, Inchinnan Business Park |  
Paisley PA49RF, Scotland, United Kingdom |  
Tel: +44 (0)141 81416305



**Manufacturer:**  
Dual manufactured products

The information in this guide is subject to change without notice.

**DISCLAIMER:** TO THE EXTENT ALLOWED BY LAW, THERMO FISHER SCIENTIFIC INC. AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

**Revision history:** Pub. No. MAN0018935

Revision	Date	Description
1.0	12 November 2019	New document

**Important Licensing Information:** These products may be covered by one or more Limited Use Label Licenses. By use of these products, you accept the terms and conditions of all applicable Limited Use Label Licenses.

©2019 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.