Enteroclon Anti-Salmonella Gr. B, C, D, E Enteroclon und Testserum Anti-Salmonella O Enteroclon Anti-Salmonella Vi

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for slide agglutination

Instructions for professional use

Intended use

The Enteroclons and test sera are intended for use either in identifying group-specific O-antigens, and thus for the assignment of an isolate to the corresponding group of the Kauffmann-White scheme, or in the detection of additional O-antigens, insofar as these are required to determine or check the serovar.

Additional serological differentiation must be performed using the Enteroclons and the test sera Anti-Salmonella H in accordance with the Kauffmann-White scheme.

Principle of the test

If the Salmonella strain possesses an antigen that is covered by the Enteroclon or test serum, this antigen becomes bound when mixed with the specific antibody. The antigen-antibody reaction results in clearly visible agglutination of the strain.

Composition

The Enteroclons are produced from cell culture supernatants of hybridoma cell lines which secrete antibodies against the respective Salmonella antigens.

The term Enteroclon is applied by the company to products that are prepared exclusively from monoclonal antibodies.

The test sera are immunized rabbit sera from which un-specific agglutinins have been removed by absorption.

Preservative: Sodium azide (NaN₃) 0.9 mg/ml

Available specificities Anti-Salmonella

Specificity	REF	Vol	Form
Gr. B	TR 1201	1 ml	liquid
	TR 5201	5 ml	liquid
Gr. C	TR 1202	1 ml	liquid
	TR 1203	1 ml	liquid
Gr. D	TR 5203	5 ml	liquid
	TR 1204	1 ml	liquid
O 2	TR 1301	1 ml	liquid
O 4	TR 1302	1 ml	liquid
	TR 5302	5 ml	liquid
O 5	TR 1303	1 ml	liquid
	TR 5303	5 ml	liquid
O 6 ¹⁾	TR 1304	1 ml	lyophilized
O 7	TR 1305	1 ml	liquid
O 8	TR 1306	1 ml	liquid
O 9	TR 1307	1 ml	liquid
	TR 5307	5 ml	liquid

Specificity	REF	Vol	Form
O 10	TR 1308	1 ml	liquid
O 11	TR 1323	1 ml	lyophilized
O 13	TR 1325	1 ml	lyophilized
O 14 ²⁾	TR 1309	1 ml	lyophilized
O 15	TR 1310	1 ml	liquid
O 16	TR 1328	1 ml	lyophilized
O 17	TR 1329	1 ml	lyophilized
O 18	TS 1330	1 ml	lyophilized
O 19	TR 1311	1 ml	liquid
O 20	TR 1312	1 ml	liquid
O 21	TR 1331	1 ml	lyophilized
O 22	TS 1332	1 ml	lyophilized
O 25	TR 1335	1 ml	lyophilized
O 27	TR 1313	1 ml	liquid
O 28	TR 1336	1 ml	lyophilized
O 30	TR 1339	1 ml	lyophilized
O 34 ³⁾	TR 1314	1 ml	liquid
O 35	TR 1341	1 ml	lyophilized
O 38	TR 1344	1 ml	lyophilized
O 39	TR 1345	1 ml	lyophilized
O 40	TR 1346	1 ml	lyophilized
O 41	TR 1347	1 ml	lyophilized
O 42	TR 1348	1 ml	lyophilized
O 43	TR 1349	1 ml	lyophilized
O 44	TR 1350	1 ml	lyophilized
O 45	TR 1351	1 ml	lyophilized
O 46	TR 1315	1 ml	liquid
O 47	TR 1353	1 ml	lyophilized
O 48	TR 1354	1 ml	lyophilized
O 50	TR 1355	1 ml	lyophilized
O 51	TR 1356	1 ml	lyophilized
O 52	TR 1357	1 ml	lyophilized
O 53	TR 1358	1 ml	lyophilized
O 54	TR 1359	1 ml	lyophilized
O 55 ⁴⁾	TR 1360	1 ml	lyophilized
O 56	TR 1361	1 ml	lyophilized
O 57	TR 1362	1 ml	lyophilized
O 58	TR 1363	1 ml	lyophilized
O 59	TR 1364	1 ml	lyophilized
O 60	TR 1365	1 ml	lyophilized
O 61	TR 1366	1 ml	lyophilized
O 62	TR 1367	1 ml	lyophilized
O 63	TR 1368	1 ml	lyophilized
O 65	TR 1369	1 ml	lyophilized
O 66	TR 1370	1 ml	lyophilized
O 67	TR 1371	1 ml	lyophilized
Vi	TR 1316	1 ml	liquid

¹⁾ To be used for testing serovars of the group C₂-C₃ for the presence or absence of the factor O 6.

²⁾ Identifies O 14 of the group H but not of the group C₁.

³⁾ May only be used for differentiation within group E, since it cross-reacts with O 12₂, O 28 and O 55.

⁴⁾ Reacts with some strains of the group O 59.

Shelf life and conditions for storage

Enteroclons that are frequently needed in the microbiological laboratory are available in liquid form. Enteroclons and test sera for specificities that are needed less often are lyophilized (freeze-dried).

If stored at 2...8 °C both before and after opening, liquid test reagents may be used up to the date given on the label. The reagents are ready to use. After use, the bottle must be properly closed.

Lyophilized Enteroclons / test sera must be dissolved in 1 ml distilled water before use. They must be properly closed with the enclosed pipette. If stored at 2...8 °C, they remain usable for at least 18 months. However, they must not be used after the date given on the label.

The test reagents/ test sera may sometimes show turbidity which is not microbially caused. Such turbidity does not impair effectivity and the test reagents/ test sera can be clarified by centrifugation or filtration.

The test reagents/ test sera must have their temperatures adjusted to room temperature (18...26 °C) before use.

Warnings and precautions

The biotechnological manufacture of the Enteroclons means that the risk of contamination by infectious agents can be virtually excluded. Because they contain animal materials (fetal calf serum, stabilizer), they should be treated as potentially infectious and handled accordingly.

The test sera are biological materials which should be treated as potentially infectious and handled accordingly.

As these materials contain sodium azide, contact with the skin and mucous membranes must be avoided. In case of contact, rinse with plenty of water.

Since the performance of the slide agglutination test involves working with native pathogenic materials, all necessary work protection procedures must be adhered to (risk of infection)!

Materials and equipment not supplied

Glass slides, stirring rods, physiological saline (NaCl 9 g/l), disposal containers for infectious material, Kauffmann-White scheme.

Test material und methodology

Transfer a small amount of bacterial mass from a 16-20 hour old subculture (e.g. nutrient agar or Kligler culture medium) onto a slide and mix well with a drop of test reagent/ test serum (ca. 25 µl), so that a homogenous, slightly milky suspension results. The slide should be placed on a dark surface.

The reaction is read with the naked eye by holding the slide in front of a light source against a black background and rocking it (tilting it back and forth).

To control for spontaneous agglutination, a negative control (NC) must be performed at the same time using physiological saline in place of the Enteroclon.


Evaluation

The test can only be evaluated if the negative control (NC) remains milky-opaque.


Positive: visible agglutination after the sample has been tilted back and forth less than 20 times. In a strongly positive reaction (PR), agglutination (coarsely or finely

flocculent) appears as soon as the bacterial mass is mixed in.

In a weakly positive result, agglutination only appears after the slide has been tilted back and forth 10-20 times.

Negative: if the suspension remains milky-opaque or the reaction begins to occur only after the slide has been tilted back and forth more than 20 times, the result is negative (NR).

Limitations of the method

The test reagents or test sera react solely with Salmonella strains which contain antigens of the specificity named in the declaration.

In exceptional cases, there is the possibility of cross-reactions with other genera of the Enterobacteriaceae family (e. g. with strains of *Citrobacter* spp., *Hafnia alvei*, *Proteus* spp. or *E. coli*) due to antigen identities or related antigens. For this reason, biochemical examination is required to establish membership of the genus Salmonella.

Explanation of the symbols used

IVD In Vitro Diagnostic Medical Device

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Information for professional use

Intended use

The Enteroclons are intended for use in the serological identification of Salmonella strains isolated from test material of varied origin with respect to the Kauffmann-White scheme, using slide agglutination.

Additional serological differentiation must be performed using the group-specific Enteroclons and the Enteroclons Anti-Salmonella O, Vi in accordance with the Kauffmann-White scheme.

Principle of the test

If the isolated strain possesses an antigen that is covered by the Enteroclon, this antigen becomes bound when mixed with the specific antibody. The antigen-antibody reaction results in clearly visible agglutination of the strain.

Composition

The Enteroclons Anti-Salmonella A-67, Anti-Salmonella I and Anti-Salmonella II are mixtures of monoclonal (mouse) antibodies.

They are produced from cell culture supernatants of hybridoma cell lines which secrete antibodies against the

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The Enteroclons Anti-Salmonella A-67 und Anti-Salmonella I are available in liquid form. Enteroclon Anti-Salmonella II is lyophilized (freeze-dried).

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The lyophilized Enteroclon Anti-Salmonella II must be dissolved in 1 ml or 5 ml distilled water in accordance with the declaration before use. It can be properly closed with the enclosed pipette. If stored at 2...8 °C, it remains usable for at least 2 years. However, it must not be used after the date given on the label.

The test reagents may sometimes show turbidity which is not microbially caused. Such turbidity does not impair effectivity and the test reagents can be clarified by centrifugation or filtration.

The test reagents must have their temperatures adjusted to room temperature (18...26 °C) before use.

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As these materials contain sodium azide, contact with the skin and mucous membranes must be avoided! In case of contact, rinse with plenty of water.

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Glass slides, stirring rods, physiological saline, disposal containers for infectious material, Kauffmann-White scheme.

Test material und methodology

Transfer a small amount of bacterial mass from a colony under test onto a slide and mix well with a drop of test reagent (ca. 25 µl), so that a homogenous, slightly milky suspension results. The slide should be placed on a dark surface.

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
In exceptional circumstances, selective media may adversely affect the agglutinability of the bacteria. This problem can be avoided by removing the bacteria from the nutrient or blood agar or from the Kligler culture medium.

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
Positive: visible agglutination after the sample has been tilted back and forth less than 20 times. In a strongly positive reaction (PR), agglutination (coarsely or finely flocculent) appears as soon as the bacterial mass is mixed in. In a weakly positive result, agglutination only appears after the slide has been tilted back and forth 10-20 times.

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Limitations of the method

The test reagents react specifically with individual Salmonella strains, i. e. they only agglutinate Salmonella strains belonging to the groups named in the declaration.

In exceptional cases, there is the possibility of cross-reactions with other genera of the *Enterobacteriaceae* family (e. g. with strains of *Citrobacter* spp., *Hafnia alvei*, *Proteus* spp. or *E. coli*) due to antigen identities or related antigens. For this reason, biochemical examination is required to establish membership of the genus Salmonella.

Explanation of the symbols used

- In Vitro Diagnostic Medical Device
- Batch code (Lot)
- Catalogue number
- Temperature limitation
- Use by YYYY-MM (MM = end of month)

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
Information for professional use

Intended use

The test reagents are intended for use in the serological identification of *Salmonella* strains isolated from test material of varied origin, using slide agglutination.

Additional serological differentiation must be performed using the group-specific test reagents (Enteroclons) and the test reagents (Enteroclons) Anti-Salmonella O, Vi in accordance with the Kauffmann-White scheme.

Principle of the test

If the isolated strain possesses an antigen that is covered by the test reagent, this *Salmonella* antigen becomes bound when mixed with the specific antibody. The antigen-antibody reaction results in clearly visible agglutination of the strain.

Composition

The test reagents are mixtures of monoclonal (mouse) antibodies. They are produced from cell culture supernatants of hybridoma cell lines which secrete antibodies against the respective group specific *Salmonella* O-antigens.

For the coverage of the O-Group K the test reagent Anti-Salmonella OMC contains test serum of immunized rabbit sera from which unspecific agglutinins have been removed by absorption.

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The test reagents are available in liquid form.

If stored at 2...8 °C before and after opening, they may be used up to the date given on the label. After use, the bottle must be properly closed.

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Warnings and precautions

The biotechnological manufacture of the test reagents means that the risk of contamination by infectious agents can be virtually excluded. Because they contain animal materials (fetal calf serum, stabilizer) and added test serum, they should be treated as potentially infectious and handled accordingly.

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Test material und methodology

The test material is streaked onto non-selective culture media (Nutrient Agar, Blood Agar) or Kligler Medium. Incubate for 16 – 20 hours at 35...37 °C.

Transfer a small amount of bacterial mass from a suspicious colony under test onto a slide and mix well with a drop of test reagent (ca. 25 µl), so that a homogenous, slightly milky suspension results. The slide should be placed on a dark surface.

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To control for spontaneous agglutination, a negative control (NC) must be performed at the same time using physiological saline in place of the test reagent.

Test procedure


First carry out the test with the O-Group Pool Test Reagents Anti-Salmonella OMA and Anti-Salmonella OMB. Approx. 98 % of *Salmonella* can be assigned to Anti-Salmonella OMA and Anti-Salmonella OMB.

If the strain does not agglutinate with these both test reagents, it is recommended to test this strain with Anti-Salmonella Vi, [REF] TR 1316. If this reaction is also negative the strain must then be agglutinated with the O-Group Pool Test Reagents Anti-Salmonella C, Anti-Salmonella D, Anti-Salmonella E, Anti-Salmonella F, Anti-Salmonella G.

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
Limitations of the method

The test reagents react specifically with individual *Salmonella* strains, i. e. they only agglutinate *Salmonella* strains belonging to the O-Group Pools named in the declaration.

In exceptional cases, there is the possibility of cross-reactions with other genera of the *Enterobacteriaceae* family (e. g. with strains of *Citrobacter* spp., *Hafnia alvei*, *Proteus* spp. or *E. coli*) due to antigen identities or related antigens. For this reason, biochemical examination is required to establish membership of the genus *Salmonella*.

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GmbH Berlin,