TDR Resin Anaerobic

【Package】

25 bottles / kit

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TDR Resin Anaerobic culture bottles are used with automated blood culture systems in qualitative procedures for recovery and detection of anaerobic and facultative anaerobic microorganisms from blood and other normally sterile body fluids. It is intended to be used for aiding to diagnose bacteremia.

【Principle】

If anaerobic and facultative anaerobic are present in the sample to be tested, the microorganisms are able to metabolize nutrients in the medium to produce carbon dioxide (CO₂); an increase in CO₂ content in the bottle will cause the sensor at the bottom of the bottle changing from blue-green gradually to yellow. The automatic microbial culture system determines whether or not the culture bottle is positive or whether or not the test sample contains viable microorganisms by identifying the speed or amount of color change.

【Composition】

Culture bottles contain sensor, complex medium, resin and an atmosphere of N₂ and CO₂. The complex medium consists of the following reactive components: combination of peptones/biological extracts (≥1.85%, W/V), carbon sources (≥0.14%, W/V), anticoagulants (≥0.05%, W/V), vitamins and amino acids (≥0.001%, W/V), reducing agent (≥0.05%, W/V), other trace elements and growth factor.

【Storage and Stability】

Stable for 12 months when stored at 2°C~30°C and protected from light, freeze is forbidden.

【Matched Instrument】

Review the appropriate Mindray automated blood culture systems User Manual before use.

【Sample Requirement】

- Obtain samples prior to initiating antibiotic therapy.
- Strictly observed the sterile procedure during sample collection to prevent contamination.
- Blood collection tubes containing anticoagulants are forbidden to be used for sample collection and culture.
- The recommended specimen volume is 5–10 ml.

【Test Procedures and Result Explanation】

1. Culture bottle preparation: label patient information on the culture bottle; open the plastic

bottle cap and carefully disinfect the stopper;

2. Blood collection procedure: collect blood by venipuncture and inoculate directly or collect blood using a syringe and inoculate;
 - (a) If venipuncture is used to directly collect blood and inoculate various types of culture bottles, first inoculate aerobic culture bottles and then anaerobic ones to avoid the oxygen remaining in the pipeline being transferred to the anaerobic culture bottles;
 - (b) If a syringe is used to collect blood and inoculate various types of culture bottles, first inoculate anaerobic culture bottles and then aerobic ones to avoid the oxygen remaining in the pipeline being transferred to the anaerobic bottles;
 - (c) The blood collection volume for each bottle is 5-10mL (maximum 10ml) and a volume close to 10ml is beneficial to improve the detection rate of microorganisms;
 - (d) An inoculated amount exceeding the upper limit should be avoided, and the scale line on the label of the culture bottle can be referenced for monitoring the sample volume in the culture bottle.
3. Culture: load the inoculated culture bottles into the automatic microbial culture system for culturing as soon as possible; if an unavoidable delay occurs, the inoculated bottles should be stored at room temperature before loading to the instrument and the time delay for loading should not exceed 24 hours.
4. Interpretation of results:
 - (a) Under the monitoring of the automatic microbial culture system, the bottles have been incubated for 5 days or detected as positive, all the positive bottles are smeared and sub-cultured.
 - (b) A negative smear for a positive bottle indicates a possible false positive and the bottle should be re-loaded into the instrument for incubation until growth occurs or the retest result is positive.
 - (c) The culture bottle which is initially determined as false positive and retested as positive should be smeared and sub-cultured.
 - (d) In some cases (where there was a high clinical suspicion of bloodstream infection with a negative result), the negative culture bottle should be checked by smear and subculture before discarding it.

【Limitations of Detection Method】

1. This medium is not suitable for the culture of obligate aerobes, fungi and mycobacterium tuberculosis.
2. *Streptococcus pneumoniae* is prone to autolysis if it is not transferred immediately after a positive is reported.
3. Rare cases that may occur: there are bacteria growing in the culture bottle, but it cannot be reported as positive due to an insufficient concentration of CO₂ produced.
4. In the case of antibiotics in the sample, there are occasional cases that bacteria grow in the culture bottle, but it cannot be reported as positive.
5. The patient sample which is confirmed as positive may contains microorganisms which can be reported as positive in the smear test but cannot grow on the conventional sub-culture mediums. If it is suspected to be this case, the sample sub-culture should be conducted on

a special medium. Similarly, positive samples may also contain microorganisms that cannot be observed by conventional smear staining methods, which may require a special staining and subculture medium for detection and culture.

6. If the specimen volume inoculated in the culture bottle is insufficient, some strains that may be susceptible to anticoagulants (sodium polyanethole sulfonate), such as *Peptostreptococcus anaerobius*, may not grow or produce a small amount of CO₂.
7. False positive results may be indicated for some negative bottles in the Gram staining. In this case, contamination from medium components, staining reagents or small amounts of inactive organisms contained in the slide should be excluded.

【Quality Control】

The quality control (QC) has been performed throughout the manufacturing period of the product, and the user should also do it according to local regulations.

QC strain: *Clostridium histolyticum* ATCC 19401, *Clostridium perfringens* bacteria ATCC 13124, *Bacteroides fragilis* ATCC 25285, *Common bacteroides* ATCC 8482, *Streptococcus pneumoniae* ATCC 6305, *Staphylococcus aureus* ATCC 25923, *Escherichia coli* ATCC 25922. Procedure: prepare a bacterial suspension with a concentration of 100 CFU/ml; inoculate 100 CFU into the culture bottle (add 5ml sterile defibrillated sheep blood); load the inoculated culture bottle into the instrument and the positive result should be reported within 72h.

【Performance Characteristic Of The Test】

Repeatability

Table 1 is the result of an internal inoculation study. Different batches of TDR Resin Anaerobic were used and 10 replicate tests were performed for each microorganism for each batch of culture bottles. After the culture bottles were reported as positive, one bottle of the replicate bottles was taken randomly for sub-culture and Gram staining microscopy. The detected strain was confirmed as the same as the inoculated one.

Table 1 Results of Repeatability Test

Microorganism	Strain ID (ATCC)	Inoculated Range (CFU/bottle)	Time to Detection		Positive rate
			Mean (H)	Range (H)	
<i>Escherichia coli</i>	25922	16~24	11.09	10.7~11.6	100%
<i>Bacteroides fragilis</i>	25285	33	49.07	38.5~73.3	100%
<i>Clostridium perfringens</i>	13124	19	12.97	12.2~13.9	100%
<i>Clostridium histolyticum</i>	19401	8	70.35	62.4~82.1	100%
<i>Bacteroides vulgatus</i>	8482	29	58.85	42.7~78.2	100%
<i>Staphylococcus aureus</i>	25923	19~29	15.83	14~17.4	100%
<i>Streptococcus pneumoniae</i>	6305	8	20.61	18.5~24.8	100%

Limit of Detection (LoD)

Table 2 is the result of an internal inoculation study. 20 replicate tests were confirmed for each microorganism, and at least 95% of the test results were achieved at LoD. In this inoculation study, one bottle of the replicate bottles was taken randomly for subculture and Gram staining microscopy. The detected strain was confirmed as the same as the inoculated one.

Table 2 Results of LoD

Microorganism	Strain ID (ATCC)	LoD (CFU/bottle)	Positive rate
<i>Escherichia coli</i>	25922	4	100%
<i>Bacteroides fragilis</i>	25285	9	100%
<i>Clostridium perfringens</i>	13124	6	100%
<i>Staphylococcus aureus</i>	25923	3	100%
<i>Enterobacter faecalis</i>	29212	5	100%
<i>Streptococcus pneumoniae</i>	6305	4	100%

Growth Performance

Table 3 demonstrated the results of the internal inoculation study with/without use of sterile defibrillated sheep blood. Different batches of culture bottles were used for detection with a target final concentration of inoculated microorganisms for each bottle of 30 CFU/bottle, while the actual range for inoculum concentration was 8-135 CFU/bottle. In this inoculation study, one bottle of the replicate bottles was taken randomly for subculture and Gram staining microscopy. The detected strain was confirmed as the same as the inoculated one.

Table 3 Results of Growth Performance Test

Microorganism	Strain ID (ATCC)	Range (CFU/bottle)	Without blood			Blood (10ml)		
			Detection rate	Time to Detection		Detection rate	Time to Detection	
				Mean (H)	Range (H)		Mean (H)	Range (H)
<i>Staphylococcus epidermidis</i>	ATCC 12228	8~25	100%	41.5	29.5~49.6	100%	15.2	14.5~15.6
<i>Listeria monocytogenes</i>	ATCC 15313	73	100%	26.9	20.9~37.5	100%	18.4	11.2~25.6
<i>Escherichia coli</i>	ATCC 25922	35~41	100%	11.0	10.8~11.4	100%	10.4	10.1~10.7
<i>Klebsiella pneumoniae</i>	ATCC 700603	34~72	100%	11.7	11.4~12.1	100%	10.5	10.4~10.5
<i>Streptococcus pneumoniae</i>	ATCC 6305	16~54	100%	21.1	18.9~24	100%	14.4	14.3~14.5
<i>Clostridium perfringens</i>	ATCC 13124	23~51	100%	13.3	13.1~13.4	100%	12.5	10.1~13.9
<i>Bacteroides fragilis</i>	ATCC	73~135	100%	40.5	36.6~42.6	100%	31.6	17.1~39.9

	25285							
<i>Enterobacter faecalis</i>	ATCC 29212	15	100%	22.5	20.6~24.4	100%	11.1	10.2~12.7
<i>Staphylococcus aureus</i>	ATCC 25923	8~24	100%	17.4	16.1~19.1	100%	13.2	12.4~14.9
<i>Proteus mirabilis</i>	ATCC 35659	32~35	100%	11.7	11.5~12.1	100%	10.9	10.5~11.7
<i>Streptococcus agalactiae</i>	ATCC 13813	30~35	100%	14.8	14.7~15.1	100%	9.4	8.7~9.8
<i>Enterobacter cloacea</i>	CICC 21539	18~29	100%	12.3	10.8~13.4	100%	11.1	10.1~11.6

Neutralization of Antimicrobials

The adsorption capacity of the resin in the culture medium for antimicrobial drugs varies depending on the drug concentration in the sample. An internal inoculation study demonstrated that antibacterial drugs were effectively neutralized by the medium in the TDR Resin Anaerobic (based on the 100% detection rate of susceptible strains at peak serum concentration). In this trial, sensitive strains were inoculated and antibacterial drugs were directly added at their peak serum concentration into culture bottles, using culture bottles without antibiotics as positive growth controls.

The following classes of antibacterial drugs can be effectively neutralized by the medium at their peak serum concentration: penicillins, aminoglycosides, glycopeptides, polymyxins, lipopeptides, oxazolidinones, quinolones.

Clinical Study Results (Blood)

The results of the blood culture in the TDR Resin Anaerobic and the subculture results of the medium after culture were compared. A clinical study conducted in Hunan Province compared the clinical performance of TDR Resin Anaerobic in blood sample detection. If the culture result was positive for a culture bottle and isolate growth was found in the subculture, then this bottle was confirmed as true positive. For a culture bottle with a negative result after 5-day incubation, if no isolate growth was seen in the subculture, then it was confirmed as true negative.

In this clinical study, totally 471 blood samples were collected and in the clinical culture test, 38 culture bottles were positive and 433 were negative. 37 bottles culture bottles reported positive were confirmed as true positive in the subculture, 1 false positives, accounting for 7.9% (37/471) and 0.2% (1/471) of the total bottles in the study, respectively. All the culture bottles reported negative obtained negative results in the subculture, accounting for 91.9% (433/471) of the total bottles in the study.

Table 4 summarized the comparison results of the blood sample detection using the TDR Resin Anaerobic and the subculture results.

Table 4 Clinical Comparison Results (Blood Samples)

TDR Resin Anaerobic	Subculture		Total
	Positive	Negative	

TDR	Positive	37	1	38
	Negative	0	433	433
Total		37	434	471

Table 5 summarized the results of the subculture to obtain isolates after the blood detection using the TDR Resin Anaerobic.

Table 5 Amount of Microorganisms Isolated (Blood Samples)

Serial	Species of microorganisms	Quantity
1	Obligate anaerobes	2
2	<i>Enterobacter aerogenes</i>	2
3	<i>Escherichia coli</i>	12
4	<i>Klebsiella pneumoniae</i>	6
5	<i>Enterobacter faecalis</i>	3
6	<i>Staphylococcus lentus</i>	1
7	<i>Staphylococcus aureus</i>	5
8	<i>Candida tropicalis</i>	2
9	<i>Staphylococcus hominis</i>	1
10	<i>Staphylococcus haemolyticus</i>	1
11	<i>Micrococcus luteus</i>	1
12	<i>Staphylococcus capitis</i>	1
Total		37

Clinical Study Results (Sterile Body Fluid)

The results of the sterile body fluid sample in the TDR Resin Anaerobic and the subculture results of the medium after culture were compared. A clinical study conducted in Hunan Province compared the clinical performance of TDR Resin Anaerobic in sterile body fluid sample detection. If the culture result was positive for a culture bottle and isolate growth was found in the subculture, then this bottle was confirmed as true positive. For a culture bottle with a negative result after 5-day incubation, if no isolate growth was seen in the subculture, then it was confirmed as true negative.

In this clinical study, totally 27 sterility body fluid samples were collected, and in the clinical culture test, 9 culture bottles were positive and 18 were negative. 8 bottles culture bottles reported positive were confirmed as true positive in the subculture, and 1 false positive was found, accounting for 29.6% (8/27) and 3.7% (1/27) of the total bottles in the study, respectively. All the culture bottles reported negative obtained confirmed negative results in the subculture, accounting for 66.7% (18/27) of the total bottles in the study.

Table 6 summarized the comparison results of the sterile body fluid sample detection using the TDR Resin Anaerobic and the subculture results.

Table 6 Clinical Comparison Results (Sterile Body Fluid Samples)

TDR Resin Anaerobic		Subculture		Total
		Positive	Negative	
TDR	Positive	8	1	9
	Negative	0	18	18
Total		8	19	27

Table 7 summarized the results of the subculture to obtain isolates after the sterile body fluid sample detection using the TDR Resin Anaerobic.

Table 7 Amount of Microorganisms Isolated (Sterility Body Fluid Samples)

Serial	Species of microorganisms	Quantity
1	<i>Staphylococcus epidermidis</i>	1
2	<i>Escherichia coli</i>	2
3	<i>Klebsiella pneumoniae</i>	1
4	<i>Staphylococcus aureus</i>	2
5	<i>Salmonella</i>	1
6	<i>Streptococcus constellatus</i>	1
Total		8

Table 8 summarized the distribution of the types of the sterile body fluid samples detected by TDR Resin Anaerobic.

Table 8 Distribution of Types of Sterile Body Fluid Samples

Serial	Sample type	Quantity	No. of positive samples
1	Puncture fluid	2	0
2	Bile	3	3
3	Catheter tip	1	1
4	Peritoneal fluid	11	2
5	Cerebrospinal fluid	7	0
6	Hydrothorax	2	1
7	Drainage fluid	1	1
Total		27	8

Percentage of Positive Cultures in Clinical Studies (Blood and Sterile Body Fluid)

In this clinical study, an overall positive culture percentage of 8.1% (38/471) with a 95% CI of 6%-11% was observed in all TDR Resin Anaerobic added with blood samples.

In this clinical study, the overall positive rate of all TDR Resin Anaerobic added with sterile body fluid was 33.3% (9/27), with a 95% CI of 16%-51%.

The expected percentage of positives varies based on various factors, such as subject

population, microorganism distribution, different trial sites, etc. The positive culture percentage provided was based on the clinical research data obtained in this clinical trial.

【Precautions】

- For IVD professional use only.
- Culture bottle should be utilized by trained laboratory personnel only.
- Before use, if the culture bottle is broken or leaking, the fluid in the culture bottle is turbid or the sensor at the bottle bottom becomes yellow, please don't use it and dispose as waste.
- Great care must be taken to prevent contamination of the patient sample during venipuncture and during inoculation into the culture bottle, since contamination could lead to a specimen being determined positive when a clinically relevant isolate is not actually present.
- If inoculated culture bottles have been delayed in their receipt into the laboratory or have been incubated prior to entry into blood culture system, they should be visually inspected for indications of microbial growth. If microbial growth is evident, treat the bottles as positive and do not place the bottle in the instrument.
- Microorganisms are often few in quantity and may appear intermittently in the blood stream; therefore, consecutive blood samples should be collected from each patient.
- The positive culture bottle may be overfilled or contains microorganisms with high gas production, leading to a high internal pressure of the positive culture bottle. Therefore, special care should be taken in subculture of positive bottles. Before Gram staining or disposal, brief venting should be done for positive culture bottles to release the gas produced during microbial metabolism.
- Attention should be paid to the inoculated amount, which is recommended to be no less than the specified minimum amount to reduce the risk of missed detection.
- The growth of some fastidious bacteria requires growth factors in the medium. For example, *Haemophilus influenzae* requires NAD and V factors and blood samples contain these growth factors. Therefore, if fastidious bacteria requiring special growth factors are suspected in non-blood samples or small blood samples, it is recommended to supplement the essential growth factors for growth support, for example, adding 10ml sterile blood such as sterile defibrillated sheep blood.
- All used culture bottles, materials of blood collection, and other waste should be disinfected per the lab requirements.
- Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the country in which the user and/or the patient is established.

【References】

1. Liu Xiguang. Modern Microbiology Diagnosis[M]. Beijing: People's Medical Publishing House, 2002.
2. Li Zhongxing, Zheng Jiaqi, Li Jiahong. Diagnostic Bacteriology[M]. Hong Kong: The Yellow River Culture Publishing House, 1992.
3. Chen Tianshou. Manufacture and Application of Microbiological Culture Media[M]. Beijing:

- China Agriculture Press, 1995.
- 4. Clinical and Laboratory Standards Institute, M22-A3, Quality Control for Commercially Prepared Microbiological Culture Media; Approved Standard, Vol.24 No.19.
- 5. Clinical and Laboratory Standards Institute, M47-A, Principles and Procedures for Blood Cultures; Approved Guideline. Vol. 27 No.17.

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Cat. No.	Description
105-006860-00	TDR Resin Anaerobic

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
Manufacturer: Hunan Mindray Medical Technology Co., Ltd.
 Manufacturing address: 3/F, West, Building C, Luvalley Science & Technology Innovation and Entrepreneurship Park, No.1698, Yuelu West Avenue, High-tech Development Zone, Changsha, 410221, P.R.China

【EC-Representative】

Shanghai International Holding Corp. GmbH (Europe)
 Eiffestraße 80, 20537 Hamburg, Germany

【Revision table】

Revision Date	Revision Version	Revision Description
2022-04	A	Initial version
<u>2022-06</u>	<u>B</u>	<u>Added 【Additional Materials Required】 and symbol “Keep away from sunlight”;</u> <u>Revised 【Intended Purpose】 and 【Composition】</u>

-  Manufacturer
-  Use by, YYYY-MM-DD
-  Authorized Representative in the European Community
-  In Vitro Diagnostic Medical Device
-  Temperature limitation

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1. Liu Xiguang. Modern Microbiology Diagnosis[M]. Beijing: People's Medical Publishing House, 2002.
2. Li Zhongxing, Zheng Jiaqi, Li Jiahong. Diagnostic Bacteriology[M]. Hong Kong: The Yellow River Culture Publishing House, 1992.
3. Chen Tianshou. Manufacture and Application of Microbiological Culture Media[M]. Beijing: China Agriculture Press, 1995.
4. Clinical and Laboratory Standards Institute, M22-A3, Quality Control for Commercially Prepared Microbiological Culture Media; Approved Standard, Vol.24 No.19.
5. Clinical and Laboratory Standards Institute, M47-A, Principles and Procedures for Blood Cultures; Approved Guideline. Vol. 27 No.17.

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Producent: Hunan Mindray Medical Technology Co., Ltd.