

## SABOURAUD CHLORAMPHENICOL GENTAMICIN AGAR

### AGAR / SELECTIVE MEDIUM FOR THE ISOLATION OF FUNGI

56594 - 63774

IVD

#### 1- INTENDED USE

Sabouraud Chloramphenicol Gentamicin Agar is recommended for the isolation of all species of yeasts and filamentous fungi from biological specimens presenting mixed fungal and bacterial flora.

#### 2- PRINCIPLE

The growth of fungi is promoted by the nutrients provided by peptones and glucose. Gentamicin and chloramphenicol inhibit the growth of most bacteria such as *Enterobacteriaceae*, *Pseudomonas*, and *Staphylococci*.

#### 3- HOW SUPPLIED

- Ready to use medium:
  - box of 20 Petri dishes (90 mm) (SCG) code 63774
  - 25 x 8 ml slant tubes code 56594

#### 4- THEORETICAL COMPOSITION (g/l of distilled water)

This agar is prepared according to the theoretical formula of medium C of the European Pharmacopoeia (1), with the addition of chloramphenicol and gentamicin.

|         |    |
|---------|----|
| Peptone | 10 |
| Glucose | 40 |
| Agar    | 12 |

#### 5- STORAGE

- Medium in tubes: at +2-8°C.
- Medium in Petri dishes: at +2-8°C.

The expiry date and batch number are indicated on the packaging.

#### 6- INSTRUCTIONS

##### Material:

- Material provided: Sabouraud Gentamicin Chloramphenicol agar.

##### Inoculation:

All biological specimens can be inoculated onto Sabouraud Gentamicin Chloramphenicol agar. Refer to current recommendations for storage of biological specimens (2).

##### Incubation:

The incubation temperature and incubation time vary according to the characteristics and the type of fungus examined. These parameters are determined by the user.

##### Temperature:

As a rule, the optimal temperature for the detection of fungi in biological specimens is 30-35°C.

The incubation of several culture media at different temperatures (25-27°C, 30-35°C) may be justified when testing for particular fungal species.

##### Duration:

The incubation time varies according to the nature of the specimen and the fungal species suspected.

Recommendations are generally as follows:

- For all specimens examined for the presence of yeasts (except *C. neoformans*): incubation for 24 to 72 hours with daily examination of the culture.
- For deep specimens: incubation for 1 to 2 weeks with daily examination of the culture.
- Culture of integuments looking for *Dermatophytes* requires an incubation for 3 to 4 weeks with twice-weekly examination.
- When *C. neoformans* or a dimorphic fungus is suspected in CSF or any other specimen: incubation for 1 month with daily examination of the culture.

##### Reading - Interpretation:

- After incubation, observe the growth of micro-organisms.
- Perform identification of the fungus (or fungi) isolated from clearly defined colonies: macroscopic and microscopic morphological examination, biochemical tests\*, complementary tests.

\* The identification of yeast colonies (matt white appearance, 1 to 2 mm in diameter) can be performed with the Auxacolor® 2 gallery (code 56513).

#### 7- PERFORMANCE/QUALITY CONTROL OF THE TEST

- Appearance of the ready to use medium: **amber-coloured** agar.
- The growth performances of Sabouraud Gentamicin Chloramphenicol agar are verified with the following strains:

| STRAINS                                  | INCUBATION CONDITIONS  | CULTURE RESULT |
|--|------------------------|----------------|
| <i>Candida albicans</i> ATCC 26790       | 24-48 hours at 30-35°C | Good growth    |
| <i>Candida tropicalis</i> ATCC 750       | 24-48 hours at 30-35°C | Good growth    |
| <i>Candida glabrata</i>                  | 24-48 hours at 30-35°C | Good growth    |
| <i>Cryptococcus neoformans</i>           | 24-48 hours at 30-35°C | Good growth    |
| <i>Staphylococcus aureus</i> ATCC 25923  | 24-48 hours at 35-38°C | No growth      |
| <i>Pseudomonas aeruginosa</i> ATCC 27853 | 24-48 hours at 35-38°C | No growth      |

#### 8- QUALITY CONTROL OF THE MANUFACTURER

All manufactured reagents are prepared according to our Quality System, starting from reception of raw material to the final commercialization of the product. Each lot is submitted to quality control assessments and is only released to the market, after conforming to pre-defined acceptance criteria. The records relating to production and control of each single lot are kept within Bio-Rad.

#### 9- LIMITS OF USE

- Antibacterial agents (chloramphenicol and gentamicin) can inhibit the growth of certain pathogenic moulds.

#### 10- REFERENCES

1. Pharmacopée Européenne, Addendum 4.2 de la IVe édition, 2002, p. 2841.
2. Basic Laboratory Procedures Clinical Bacteriology. World Health Organization. Geneva. 1991. 1<sup>st</sup> edition.



**Bio-Rad**

3, boulevard Raymond Poincaré  
92430 Marnes-la-Coquette France  
Tel. : +33 (0) 1 47 95 60 00  
Fax : +33 (0) 1 47 41 91 33



11/2003