Mannitol Salt Agar

Selective medium for isolation and enumeration of staphylococci from clinical samples and other materials, according to USP/EP/JP.

DESCRIPTION

Mannitol Salt Agar is a selective medium used for isolating pathogenic staphylococci from clinical samples, food and other materials of sanitary importance.

This medium is prepared according to recommendations of the harmonized USP/EP/JP method for the detection of *S. aureus* in non sterile pharmaceutical products.

TYPICAL FORMULA	(g/l)
Pancreatic Digest of Casein	5.0
Peptic Digest of Animal Tissue	5.0
Beef Extract	1.0
D-Mannitol	10.0
Sodium Chloride	75.0
Phenol Red	0.025
Agar	15.0
Final pH 7.4 ± 0.2 at 25°C	

METHOD PRINCIPLE

Pancreatic digest of casein, peptic digest of animal tissue and beef extract provide amino acids, nitrogen, carbon, vitamins and minerals for organisms growth. Mannitol is the fermentable carbohydrate. The high salt content of 7.5% inhibits most bacteria other than staphylococci. Phenol red is the pH indicator. Agar is the solidifying agent.

PREPARATION

<u>Dehydrated medium</u>	Suspend 111 g of the powder in 1 liter of distilled or deionized water. Mix well. Heat to boil for 1 minute shaking frequently until completely dissolved. Sterilize in autoclave at 121°C for 15 minutes.
<u>Medium in bottles</u>	Melt the content of the bottle in a water bath at 100°C (loosing the cap partially removed) until completely dissolved. Then screw the cap and check the homogeneity of the dissolved medium, if it is the case turning the bottle upside down. Cool at 45-50°C, mix well avoiding foam formation and aseptically distribute into Petri dishes.

TEST PROCEDURE

Inoculate plates by the direct streaking of the material to be examined over the agar surface. Incubate aerobically at 35 ± 2°C for 24-48 hours.

Harmonized USP/EP/JP method for microbiological examination of non sterile products recommends to inoculate the sample in Tryptic Soy Broth (ref. 24444). Subculture on a plate of Mannitol Salt Agar and incubate at 30-35°C for 18-72 hours.

INTERPRETING RESULTS

S. aureus cultivates with yellow or white colonies surrounded by a yellow zone. Confirm by identification tests*. Coagulase-negative Staphylococci form small colorless to red colonies with no color change to the medium

*Suspect colonies can be subcultured to a moderately selective medium such as Baird Parker RPF Agar (ref. 10521, 402210) for the determination of coagulase activity (ISO 6888-2).

APPEARANCE OF THE MEDIUM

Dehydrated medium: free-flowing, homogeneous, beige-pink.

Prepared medium: slightly opalescent, pinkish-red.

STORAGE

The powder is very hygroscopic, store the powder at 10-30°C, in a dry environment, in its original container tightly closed. Store bottles and prepared plates at 10-25°C away from light. Do not use the product beyond its expiry date on the label or if product shows any evidence of contamination or any sign of deterioration.

SHELF LIFE

Dehydrated medium: 4 years.
 Medium in bottles: 2 years.
 Ready-to-use plates: 6 months.

QUALITY CONTROL

Plates are inoculated with the microbial strains indicated in the QC table.
 Inoculum for productivity: 10-100 CFU
 Inoculum for selectivity: 10⁴-10⁶ CFU
 Incubation conditions: aerobically at 35 ± 2°C for 24-48 hours.
 *30-35°C for 18-72 h (USP/EP/JP Growth Promotion Testing).

QC Table.

Microorganism		Growth	Specification
<i>Staphylococcus aureus</i>	ATCC® 25923	Good	Yellow colonies with yellow zone
<i>Staphylococcus aureus</i> *	ATCC® 6538	Good	Yellow colonies with yellow zone
<i>Staphylococcus epidermidis</i>	ATCC® 12228	Good	Red colonies
<i>Escherichia coli</i>	ATCC® 25922	Inhibited	---
<i>Escherichia coli</i> *	ATCC® 8739	Inhibited	---

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BIBLIOGRAPHY

1. European Pharmacopoeia 6.5 (2009). 2.6.13 Microbiological examination of non-sterile products: Test for specified microorganisms.
2. United States Pharmacopoeia 32 NF 27 (2009). <62> Microbiological examination of non-sterile products: Test for specified microorganisms.
3. Japanese Pharmacopoeia 4.05 (2008). Microbiological examination of non-sterile products: Test for specified microorganisms.
4. ISO 6888-2:1999 + A1:2003. Microbiology of food and animal feeding stuffs – Horizontal method for the enumeration of coagulase-positive staphylococci (*Staphylococcus aureus* and other species) – Part 2: Technique using rabbit plasma fibrinogen agar medium.
5. Kloos, W.E., and T.L. Bannerman (1995) *Staphylococcus* and *Micrococcus*. In Manual of clinical microbiology, 6th ed.
6. Chapman, G.H. (1945) The significance of sodium chloride in studies of staphylococci. J. Bacteriol. 50:201-203.

PRESENTATION		Contents	Ref.
Mannitol Salt Agar	90 mm ready-to-use plates	20 plates	10030
Mannitol Salt Agar	90 mm ready-to-use plates	100 plates	10030*
Mannitol Salt Agar	Bottles	6 x 500 ml bottles	470080
Mannitol Salt Agar	Bottles	6 x 200 ml bottles	412290
Mannitol Salt Agar	Bottles	6 x 100 ml bottles	402290
Mannitol Salt Agar	Dehydrated medium	500 g of powder	610029
Mannitol Salt Agar	Dehydrated medium	100 g of powder	620029
Mannitol Salt Agar	Dehydrated medium	5 kg of powder	6100295

TABLE OF SYMBOLS

LOT Batch code	IVD <i>In vitro</i> Medical Diagnostic Device	 Manufacturer	 Use by	 Fragile, handle with care
REF Catalogue number	 Temperature limitation	 Contains sufficient for <n> tests	 Caution, consult Instruction For Use	 Do not reuse

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*30-35°C per 18-72 ore (USP/EP/JP Growth Promotion Testing).

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liofilchem@liofilchem.net

Sabouraud Dextrose Agar

Medium for the cultivation and enumeration of yeasts and moulds from different materials, according to EN ISO 11133 and USP/EP/JP.

DESCRIPTION

Sabouraud Dextrose Agar (SDA) is a non selective isolation medium used for the growth and maintenance of pathogenic and non-pathogenic fungi from clinical and nonclinical specimens. It is also used for recovery and total counting of yeasts and moulds in environmental monitoring.

This medium complies with EN ISO 11133 for microbiological examination of food, animal feed and water, where it is described as the main reference medium to carry out quantitative testing on culture media intended for fungi.

Its formula conforms to the recommendations of the harmonized method in the United States Pharmacopoeia (USP), European Pharmacopoeia (EP) and Japanese Pharmacopoeia (JP) for the microbiological examination of non sterile products. The medium is also available as gamma-irradiated triple bagged plates, particularly suitable for use in restricted areas like isolators and clean rooms.

TYPICAL FORMULA	(g/l)
Pancreatic Digest of Casein	5.0
Peptic Digest of Animal Tissue	5.0
Dextrose	40.0
Agar	15.0
Final pH 5.6 ± 0.2 at 25°C	

METHOD PRINCIPLE

Pancreatic digest of casein and peptic digest of animal tissue provide amino acids, nitrogen, carbon, vitamins and minerals for organisms growth. Dextrose is an energy source. Agar is the solidifying agent. The high concentration of dextrose and the acidic pH of the medium permit selectivity of fungi.

The medium can be supplemented with chloramphenicol to increase bacterial inhibition and recovery of dermatophytes.

PREPARATION

<u>Dehydrated medium</u>	Suspend 65 g of the powder in 1 liter of distilled or deionized water. Mix well. Heat to boil shaking frequently until completely dissolved. Sterilize in autoclave at 121°C for 15 minutes.
<u>Medium in bottles</u>	Melt the content of the bottle in a water bath at 100°C (loosing the cap partially removed) until completely dissolved. Then screw the cap and check the homogeneity of the dissolved medium, if it is the case turning the bottle upside down. Cool at 45-50°C, mix well avoiding foam formation and aseptically distribute into Petri dishes.

TEST PROCEDURE

For use in medical microbiology

Streak the specimen as soon as possible after it is received in the laboratory to obtain isolated colonies. Prepared tubed slants primarily are intended for use with pure cultures for maintenance or other purposes. Incubation conditions may vary according to the type of specimen and the microorganisms being tested for.

For use in food, animal feed and water testing

Refer to EN ISO 11133 for specific instructions.

For use in industrial microbiology

Control of non-sterile products

Refer to the procedure described in the harmonized chapters of the Pharmacopoeia.

Passive Air Monitoring

Take the lid off the settle plate and leave the medium exposed to the air for a period of time no longer than 4 hours (settling plates filled with 30 ml of medium may compensate for water loss during extended incubation periods). Plates can be placed according to the 1/1/1 scheme (for 1 h, about 1 above the floor, at least 1 m from the walls or any obstacle).

Surfaces and Personnel Hygiene Monitoring

Take a swab sample for irregular surfaces or use the sampling template 10x10 (ref. 96762) to sample a well defined area of the test surface. Inoculate a 90 mm plate by streaking the swab over the agar surface. Furthermore, the medium is suitable for personnel hygiene monitoring to detect microbial contamination of gloves or hands e.g. in a 5-finger-print.

Incubate the plates at 20-25°C for 5-7 days or at 30-35°C for 24-48 hours.

INTERPRETING RESULTS

Transfer of growth from slants to plated media may be required in order to obtain pure cultures of fungi. Examine for fungal colonies exhibiting typical microscopic and colonial morphology. Biochemical tests may be required for final identification.

The total combined yeasts/moulds count (TYMC) is considered to be equal to the number of CFU found per each plate. When an acceptable criterion for microbiological quality is prescribed it is interpreted as follows:

- 10¹ CFU: maximum acceptable count = 20;
- 10² CFU: maximum acceptable count = 200;
- 10³ CFU: maximum acceptable count = 2000, and so forth.

In procedures intended for environmental and personnel hygiene monitoring, observe daily for the formation of colonies.

APPEARANCE

Dehydrated medium: free-flowing, homogeneous, light beige.

Prepared medium: slightly opalescent, light amber.

STORAGE

The powder is very hygroscopic, store the powder at 10-30°C, in a dry environment, in its original container tightly closed. Store bottles, tubes and prepared plates at 10-25°C away from light. Do not use the product beyond its expiry date on the label or if product shows any evidence of contamination or any sign of deterioration.

SHELF LIFE

Dehydrated medium: 4 years.

Medium in bottles: 2 years.

Medium in tubes: 1 year.

Ready-to-use plates (90 and 60 mm): 6 months.

Contact plates (55 mm): 9 months

QUALITY CONTROL

The medium is inoculated with the microbial strains indicated in the QC table.

Inoculum for productivity: 50-100 CFU.

Incubation conditions: 32.5 ± 2.5°C for 24-48 h (*C. albicans*) and at 22.5 ± 2.5°C for up to 5 days (all listed organisms), under aerobic atmosphere.

QC Table.

Microorganism		Growth
<i>Candida albicans</i>	WDCM 00054	Good
<i>Aspergillus brasiliensis</i>	WDCM 00053	Good
<i>Saccharomyces cerevisiae</i>	WDCM 00058	Good

WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is intended for *in vitro* diagnostic use and must be used only by properly trained operators.

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BIBLIOGRAPHY

1. EN ISO 11133:2014+Amd1:2018. Microbiology of food, animal feed and water – Preparation, production, storage and performance testing of culture media.
2. European Pharmacopoeia 6.5 (2009) 2.6.13. Microbiological examination of non-sterile products: Test for specified microorganisms.
3. United States Pharmacopoeia 32 NF 27 (2009) <62> Microbiological examination of non-sterile products: Test for specified microorganisms.
4. Japanese Pharmacopoeia 4.05 (2008) Microbiological examination of non-sterile products: Test for specified microorganisms.
5. Sabouraud, R. (1892) Ann. Dermatol. Syphilol. 3:1061.

PRESENTATION	Category	Packaging	Ref.
Sabouraud Dextrose Agar	90 mm plates	20 plates	10035
Sabouraud Dextrose Agar	90 mm plates	100 plates	10035*
Sabouraud Dextrose Agar	90 mm plates (triple-wrapped and gamma-irradiated)	20 plates	10035S
Sabouraud Dextrose Agar	90 mm plates (triple-wrapped and gamma-irradiated, 30 ml filling volume)	20 plates	10114S f
Sabouraud Dextrose Agar	60 mm plates	20 plates	163402 f
Sabouraud Dextrose Agar	60 mm plates	450 plates	173402 f
Sabouraud Dextrose Agar	Tubes - Bottles	10 x 9 ml slant tubes	30093
Sabouraud Dextrose Agar	55 mm contact plates	20 plates	15327 f
Sabouraud Dextrose Agar	55 mm contact plates irradiated	20 plates	15327S f
Sabouraud Dextrose Agar	Tubes - Bottles	20 x 9 ml slant tubes	31093
Sabouraud Dextrose Agar	Tubes - Bottles	6 x 500 ml bottles	470040
Sabouraud Dextrose Agar	Tubes - Bottles	6 x 200 ml bottles	412280
Sabouraud Dextrose Agar	Tubes - Bottles	25 x 200 ml bottles	452280
Sabouraud Dextrose Agar	Tubes - Bottles	6 x 100 ml bottles	402280
Sabouraud Dextrose Agar	Dehydrated culture medium	500 g of powder	610103
Sabouraud Dextrose Agar	Dehydrated culture medium	100 g of powder	620103
Sabouraud Dextrose Agar	Dehydrated culture medium	5 kg of powder	6101035

f: Not CE Marked

TABLE OF SYMBOLS

LOT Batch code	IVD <i>In vitro</i> Diagnostic Medical Device	 Manufacturer	 Use by	 Fragile, handle with care
REF Catalogue number	 Temperature limitation	 Contains sufficient for <n> tests	 Caution, consult Instruction For Use	 Do not reuse

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INTERPRETING RESULTS

Transfer of growth from slants to plated media may be required in order to obtain pure cultures of fungi.

Examine for fungal colonies exhibiting typical microscopic and colonial morphology. Biochemical tests may be required for final identification.

The total combined yeasts/moulds count (TYMC) is considered to be equal to the number of CFU found per each plate.

When an acceptable criterion for microbiological quality is prescribed it is interpreted as follows:

- 10¹ CFU: maximum acceptable count = 20;
- 10² CFU: maximum acceptable count = 200;
- 10³ CFU: maximum acceptable count = 2000, and so forth.

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Nutrient Agar ISO 16266

Medium for cultivating non-fastidious organisms and confirming *Pseudomonas aeruginosa*, according to ISO 16266.

DESCRIPTION

Nutrient Agar ISO 16266 is a medium used for the cultivation of non-fastidious organisms from clinical specimens and environmental samples.

This medium is formulated according to ISO 16266 for the detection and enumeration of *Pseudomonas aeruginosa* in water by the membrane filtration technique.

TYPICAL FORMULA

	(g/l)
Peptone	5.0
Meat Extract	1.0
Yeast Extract	2.0
Sodium Chloride	5.0
Agar	15.0
Final pH 7.4 ± 0.2 at 25°C	

METHOD PRINCIPLE

Peptone and meat extract provide amino acids, nitrogen, carbon, vitamins and minerals for organisms growth. Yeast extract is a source of vitamins, particularly of B-group. Sodium chloride maintains the osmotic balance of the medium. Agar is the solidifying agent.

PREPARATION

<u>Dehydrated medium</u>	Suspend 28 g of the powder in 1 liter of distilled or deionized water. Mix well. Heat to boil shaking frequently until completely dissolved. Sterilize in autoclave at 121°C for 15 minutes.
<u>Medium in bottles</u>	Melt the content of the bottle in a water bath at 100°C (loosing the cap partially removed) until completely dissolved. Then screw the cap and check the homogeneity of the dissolved medium, if it is the case turning the bottle upside down. Cool at 45-50°C, mix well avoiding foam formation and aseptically distribute into Petri dishes.

TEST PROCEDURE

According to ISO 16266, transfer the membrane and presumptive *Pseudomonas aeruginosa* to a plate of Nutrient Agar ISO 16266. Incubate aerobically at 36 ± 2°C for 20-24 hours.

Alternatively, the medium can be inoculated by spread plating or direct streaking of the sample over the agar surface.

INTERPRETING RESULTS

Observe for colony growth. Confirm *P. aeruginosa* by performing the oxidase test (ref. 88029).

APPEARANCE

Dehydrated medium: free-flowing, homogeneous, beige.

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Inoculum for productivity: 50-100 CFU

Incubation conditions: aerobically at $36 \pm 1^\circ\text{C}$ for 20-24 hours.

QC Table.

Microorganism	Growth
<i>Pseudomonas aeruginosa</i> ATCC® 27853	Good
<i>Escherichia coli</i> ATCC® 25922	Good

WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is intended for *In vitro* diagnostic use and must be used only by properly trained operators.

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BIBLIOGRAPHY

- ISO 16266:2008. Water Quality – Detection and enumeration of *Pseudomonas aeruginosa* – Method by membrane filtration.
- Marshall, R.T. (1993) Standard methods for the microbiological examination of dairy products, 16th ed.

PRESENTATION	Contents	Ref.
Nutrient Agar ISO 16266 90 mm ready-to-use plates	20 plates	10044
Nutrient Agar ISO 16266 90 mm ready-to-use plates	100 plates	10044*
Nutrient Agar ISO 16266 Slant tubes	20 x 7 ml tubes	31083
Nutrient Agar ISO 16266 Slant tubes	10 x 7 ml tubes	30083
Nutrient Agar ISO 16266 Bottles	6 x 500 ml bottles	470060
Nutrient Agar ISO 16266 Bottles	6 x 200 ml bottles	412190
Nutrient Agar ISO 16266 Bottles	6 x 100 ml bottles	402190
Nutrient Agar ISO 16266 Dehydrated medium	100 g of powder	620036
Nutrient Agar ISO 16266 Dehydrated medium	500 g of powder	610036
Nutrient Agar ISO 16266 Dehydrated medium	5 kg of powder	6100365

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liofilchem@liofilchem.net

Selenite Broth

Liquid medium for selective enrichment of *Salmonella* spp, from clinical and nonclinical samples, according to APHA.

DESCRIPTION

Selenite Broth is an enrichment medium used for the selective isolation of *Salmonella* and some species of *Shigella*.

This medium is prepared according to the original formula described as Selenite F Broth by Leifson and recommended by the American Public Health Association for the examination of food.

TYPICAL FORMULA

	(g/l)
Enzymatic Digest of Casein	5.0
Lactose	4.0
Sodium Phosphate	10.0
Sodium Selenite	4.0

Final pH 7.0 ± 0.2 at 25°C

METHOD PRINCIPLE

Enzymatic digest of casein provides amino acids, nitrogen, carbon, vitamins and minerals for organisms growth. Lactose is the fermentable carbohydrate. Sodium phosphate is the buffer. Sodium selenite is the selective agent inhibiting many species of Gram-positive and Gram-negative bacteria including enterococci and coliforms.

PREPARATION

Dehydrated medium Suspend 23 g of the powder in 1 liter of distilled or deionized water. Mix well. Heat to boil shaking frequently until completely dissolved. Dispense into suitable containers (bottles or tubes). A depth of at least 5 cm is recommended, as salmonellae survive better at low oxygen tensions. DO NOT AUTOCLAVE.

TEST PROCEDURE

Inoculate the tube with 1-2 g of stool specimen or other solid material (approximately 10-15% by volume) and emulsify in the broth. For urines, the broth should be used at double concentration and inoculated with its own volume of the specimen. Incubate at 35 ± 2°C for 12-24 hours (coliforms may overgrow the pathogens if incubated for longer than 24 hours).

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Subculture to a selective and differential enteric plated medium, such as XLD Agar (ref. 10056), Hektoen Enteric Agar (ref. 10043) or MacConkey Agar (ref. 10029), streaking for isolation. Examine for typical colony morphology. Confirm with further biochemical tests.

APPEARANCE

Dehydrated medium: free-flowing, homogeneous, white to light beige.

Prepared medium: clear, very pale yellow.

STORAGE

The powder is very hygroscopic, store the powder at 10-30°C, in a dry environment, in its original container tightly closed. Store bottles and tubes at 2-8°C away from light. Do not use the product beyond its expiry date on the label or if product shows any evidence of contamination or any sign of deterioration.

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Inoculum for productivity: ≤ 100 CFU.

Inoculum for selectivity: $>10^3$ CFU.

Incubation conditions: aerobically at $35 \pm 2^\circ\text{C}$ for 18-24 hours.

QC Table.

Microorganism		Growth
<i>Salmonella</i> Typhimurium	ATCC® 14028	Good
<i>Shigella sonnei</i>	ATCC® 25931	Good
<i>Escherichia coli</i>	ATCC® 25922	Partially to completely inhibited

WARNING AND PRECAUTIONS

The product contains hazardous substances and is classified as dangerous. It is recommended to consult the safety data sheet for its correct use. The product is intended for *In vitro* diagnostic use and must be used only by properly trained operators.

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BIBLIOGRAPHY

1. Versalovic, J., K.C. Carroll, G. Funke, J.H. Jorgensen, M.L. Landry, and D.W. Warnock (2011) Manual of Clinical Microbiology. 10th ed. ASM Press, Washington, D.C.
2. Quality Control for Commercially Prepared Microbiological Media (2004) - 3rd ed. M22-A3. Clinical and Laboratory Standards Institute - CLSI (NCCLS), Wayne, PA.
3. Vanderzant, C., and D.F. Splittstoesser (eds.). Compendium of methods for the microbiological examination of foods, 3rd ed. American Public Health Association, Washington, D.C.
4. Leifson, E. (1939) New selenite selective enrichment medium for the isolation of typhoid and paratyphoid bacilli. Am. J. Hyg. 24:423-432.

PRESENTATION		Contents	Ref.
Selenite Broth	Tubes	20 x 10 ml tubes	24110
Selenite Broth	Tubes	20 x 5 ml tubes	24143
Selenite Broth	Bottles	6 x 100 ml bottles	402050
Selenite Broth	Bottles	6 x 200 ml bottles	412050
Selenite Broth (Double Concentration)	Bottles	6 x 200 ml bottles	432050
Selenite Broth	Bottles	6 x 500 ml bottles	470020
Selenite Broth	Bottles	6 x 1000 ml bottles	463130
Selenite Broth	Dehydrated medium	500 g of powder	610145
Selenite Broth	Dehydrated medium	100 g of powder	620145
Selenite Broth	Dehydrated medium	5 kg of powder	6101455

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COAGULASE TEST

Freeze-dried rabbit plasma for coagulase test.

TYPICAL FORMULA

Freeze-dried rabbit plasma with EDTA..... 4.0 ml

DESCRIPTION

COAGULASE TEST is constituted by rabbit freeze-dried plasma with EDTA, for the determination of enzyme coagulase produced by *Staphylococcus aureus*.

PRINCIPLE

The test shows the coagulant activity of potentially pathogenic staphylococci on plasma. This activity is due to two factors at least: free coagulase (extracellular enzyme) and bound coagulase or clumping factor (antigen of cellular wall).

The enzyme is present in the most of biotypes belonging to the species of *Staphylococcus aureus* and in biotypes belonging to the species *S. intermedius* and *S. hyicus*, pathogenic opportunist for animals, and is always absent in saprophyte and commensal species.

TECHNIQUE

1. Take a bottle of freeze-dried plasma from the kit and aseptically reconstitute with 4 ml of sterile physiological solution. Mix until completely dissolved, avoiding foam formation.
2. Prepare a culture broth in Brain Heart Infusion, taking one or more colonies from selective culture media for *Staphylococcus aureus* isolation, and incubate at 36+/-1 °C for 4-6 hours.
3. Mix 0.5 ml of Coagulase Test with 0.5 ml of culture broth and incubate at 36+/-1 °C for 1-2-4-8-24 hours.
4. Examine for the formation of coagulum, carefully inclining the tube to one side, without shaking. Do not incubate after 24 hours because cases of fibrinolysis may occur.

INTERPRETATION OF RESULTS

In the most or in the whole culture medium the presence of enzyme coagulase is shown by the formation of a well defined coagulum. The test determines both free and bound coagulase.

STORAGE

2-8 °C away from light, until the expiry date on the label or until signs of deterioration or contamination are evident.

WARNING and PRECAUTIONS

The product contains dangerous substances according to directives 1999/45/CE and 2001/60/CE or for which exist recognized exposure limits.

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COAGULASE TEST

PRESENTATION

Glass bottles containing 4 ml rabbit freeze-dried plasma with EDTA.

STORAGE

2-8 °C

PACKAGING

Code	Content	Packaging
88030	5 bottles x 40 test	• 5 bottles in box

USE

COAGULASE TEST is constituted by rabbit freeze-dried plasma with EDTA, for the determination of enzyme coagulase produced by *Staphylococcus aureus*.

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APPEARANCE of the MEDIUM

Compact freeze-dried, beige to light pink.

SHELF LIFE

2 years

QUALITY CONTROL

- Control of general characteristics, label and print
- Microbiological control
Incubation conditions: 12-24 h at 36 ± 1 °C, in aerobiosis

Microorganisms		Coagulation
<i>Escherichia coli</i>	ATCC 25922	-
<i>Staphylococcus aureus</i>	ATCC 25923	+

TABLE OF SYMBOLS

Symbol	Meanings
	Catalogue number
	In vitro Diagnostic Medical Device
	Manufacturer
	Temperature limitation
	Use by
	Batch code
	Consult accompanying documents

Sabouraud Dextrose Broth

Liquid medium for the cultivation of yeasts and moulds from different materials, according to USP/EP/JP.

DESCRIPTION

Sabouraud Dextrose Broth (SDB) is a liquid medium recommended for use in qualitative procedures for isolation of yeasts and moulds and for the culture or subculture of fungi from clinical and nonclinical specimens.

This medium conforms to the requirements of the harmonized method in the United States Pharmacopoeia (USP), European Pharmacopoeia (EP) and Japanese Pharmacopoeia (JP) for the microbiological examination of non sterile products.

TYPICAL FORMULA (g/l)

Pancreatic Digest of Casein	5.0
Peptic Digest of Animal Tissue	5.0
Dextrose	20.0
Final pH 5.6 ± 0.2 at 25°C	

METHOD PRINCIPLE

Pancreatic digest of casein and peptic digest of animal tissue provide amino acids, nitrogen, carbon, vitamins and minerals for organisms growth. Dextrose is an energy source. The high concentration of dextrose and the acidic pH of the medium permit selectivity of fungi.

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PREPARATION

Dehydrated medium Suspend 30 g of the powder in 1 liter of distilled or deionized water. Mix well. Heat to boil shaking frequently until completely dissolved. Dispense into appropriate containers. Sterilize in autoclave at 121°C for 15 minutes.

TEST PROCEDURE

For use in medical microbiology

Inoculate the specimen directly into the broth. Incubate aerobically at 25°C for 2-7 days (incubation conditions may vary according to the type of specimen and the microorganisms being tested for).

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To prepare the fungal test strains grow *C. albicans* or *A. brasiliensis* at 20-25°C for 48-72 hours or 5-7 days, respectively.

To test for *C. albicans*, inoculate the preparation of the product to be examined 1:100 in SDB and incubate at 30-35°C for 3-5 days. Subculture on a plate of Sabouraud Dextrose Agar (ref. 10035).

INTERPRETING RESULTS

Turbidity indicates microbial growth.

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Prepared medium: clear, light amber, may have a slight precipitate.

STORAGE

The powder is very hygroscopic, store the powder at 10-30°C, in a dry environment, in its original container tightly closed. Store tubes and bottles at 10-25°C away from light. Do not use the product beyond its expiry date on the label or if product shows any evidence of contamination or any sign of deterioration.

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Inoculum for productivity: ≤ 100 CFU.

Incubation conditions: $32.5 \pm 2.5^\circ\text{C}$ for 48-72 h (*C. albicans*) and at $22.5 \pm 2.5^\circ\text{C}$ for up to 5 days (all listed organisms), under aerobic atmosphere.

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BIBLIOGRAPHY

1. European Pharmacopoeia 6.5 (2009) 2.6.13. Microbiological examination of non-sterile products: Test for specified microorganisms.
2. United States Pharmacopoeia 32 NF 27 (2009) <62> Microbiological examination of non-sterile products: Test for specified microorganisms.
3. Japanese Pharmacopoeia 4.05 (2008) Microbiological examination of non-sterile products: Test for specified microorganisms.
4. Sabouraud, R. (1892) Ann. Dermatol. Syphilol. 3:1061.

PRESENTATION

		Contents	Ref.
Sabouraud Dextrose Broth	Tubes	20 x 10 ml tubes	24109
Sabouraud Dextrose Broth	Bottles	6 x 100 ml bottles	402040
Sabouraud Dextrose Broth	Bottles	25 x 100 ml bottles	452040
Sabouraud Dextrose Broth	Bottles	6 x 500 ml bottles	471070
Sabouraud Dextrose Broth	Dehydrated medium	500 g of powder	610104
Sabouraud Dextrose Broth	Dehydrated medium	100 g of powder	620104

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liofilchem@liofilchem.net

Bile Aesculin Azide Agar

Selective medium for detection and enumeration of enterococci in water and other materials, according to ISO 7899-2.

TYPICAL FORMULA	(g/l)
Tryptone	17.0
Peptone	3.0
Yeast Extract	5.0
Ox-bile	10.0
Sodium Chloride	5.0
Aesculin	1.0
Ferric Ammonium Citrate	0.5
Sodium Azide	0.15
Agar	15.0
Final pH 7.1 ± 0.1 at 25°C	

DESCRIPTION

Bile Aesculin Azide Agar is a selective medium used for isolating and enumerating enterococci from environmental samples of sanitary importance and clinical specimens.

This medium complies with ISO 7899-2 for rapid confirmation of typical colonies on the primary isolation Slanetz Bartley Agar.

PRINCIPLE

Tryptone and peptone provide amino acids, nitrogen, carbon, vitamins and minerals for organisms growth. Yeast extract is a source of vitamins, particularly of B-group. Ox-bile inhibits the growth of numerous accompanying bacteria. Sodium chloride maintains the osmotic balance of the medium. The glycoside aesculin is hydrolyzed from enterococci to aesculetin and glucose. The aesculetin reacts with iron ions forming a dark brown or black complex. Sodium azide suppress the growth of Gram-negative bacteria. Agar is the solidifying agent.

PREPARATION

Suspend 56.7 g of powder in 1 liter of deionized or distilled water. Bring to boil and shake until completely dissolved. Mix well. Sterilize in autoclave at 121°C for 15 minutes. Cool up to 45-50°C. Pour in Petri dishes.

TECHNIQUE

ISO 7899-2 recommends to filter the water sample through a filter membrane (0.45 µm pore diameter), transfer the membrane onto a Slanetz Bartley Agar plate (ref. 163462) and incubate aerobically at 36 ± 2°C for 40-48 h.

Confirm red-maroon-pink colonies by transferring the membrane and the colonies onto a plate of Aesculin Azide Bile Agar which has been preheated to 44°C. Incubate at 44 ± 0.5°C for 2 h.

Alternatively, sample can be inoculated by spread plating, pour plating or by direct streaking on the medium surface. Incubate at 35 ± 2°C for 18-24 h.

INTERPRETATION OF RESULTS

Enterococci typically produce colonies showing a tan-black color in the surrounding medium.

STORAGE

The powder is very hygroscopic, store the powder at 10-30°C, in a dry environment, in its original container tightly closed and use it before the expiry date on the label or until signs of deterioration or contamination are evident. Store prepared plates at 2-8°C away from light.

WARNING AND PRECAUTIONS

The product contains hazardous substances and is classified as dangerous. It is recommended to consult the safety data sheet for its correct use. The product is designed for *in vitro* diagnostic use only and must be used by properly trained operators.

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REFERENCES

- ISO 7899-2:2000. Water quality – Detection and enumeration of intestinal enterococci – Part 2: Membrane filtration method.
- Facklam R.R. and M. Moody (1970) Presumptive identification of group D streptococci: the bile-aesculin test. *App. Microbiol.* 20:245-250.
- Isenberg H.D. and D. Goldber (1970) Laboratory studies with a selective Enterococcus medium. *Appl. Microbiol.* 20:433-436
- Slanetz L.W. and C.H. Bartley (1957) Numbers of enterococci in water, sewage and faeces determined by the membrane filtration technique with an improved medium. *J. Bact.* 74:591-595.

PRODUCT SPECIFICATIONS

NAME

Bile Aesculin Azide Agar

PRESENTATION

Dehydrated medium

STORAGE

10-30°C

PACKAGING

Ref.	Content	Packaging
610001	500 g	500 g of powder in plastic bottle
620001	100 g	100 g of powder in plastic bottle
6100015	5 Kg	5 kg of powder in plastic bottle

pH OF THE MEDIUM

7.1 ± 0.1

USE

Bile Aesculin Azide Agar is a selective medium used for confirmation and enumeration of enterococci from water and other samples according to ISO 7899-2

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4 years

QUALITY CONTROL

- Control of general characteristics, label and print
- Microbiological control
Inoculum for productivity: 50-100 CFU
Inoculum for selectivity: 10^4 - 10^6 CFU
Incubation Conditions: 18-24 h at $35 \pm 2^\circ\text{C}$, in aerobiosis

Microorganism		Growth	Specification
<i>Enterococcus faecalis</i>	ATCC® 19433	Good	Blackening
<i>Enterococcus faecium</i>	ATCC® 19434	Good	Blackening
<i>Escherichia coli</i>	ATCC® 25922	Inhibited	---
<i>Streptococcus pyogenes</i>	ATCC® 19615	Inhibited	---

TABLE OF SYMBOLS

LOT	Batch code	IVD	In vitro Diagnostic Medical Device		Manufacturer		Use by		Fragile, handle with care
REF	Catalogue number		Temperature limitation		Contains sufficient for <n> tests		Caution, consult instructions for use		Do not reuse

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1. ISO 7899-2:2000. Water quality – Detection and enumeration of intestinal enterococci – Part 2: Membrane filtration method.
2. Facklam R.R. and M. Moody (1970) Presumptive identification of group D streptococci: the bile-aesculin test. App. Microbiol. 20:245-250.
3. Isenberg H.D. and D. Goldber (1970) Laboratory studies with a selective Enterococcus medium. Appl. Microbiol. 20:433-436
4. Slanetz L.W. and C.H. Bartley (1957) Numbers of enterococci in water, sewage and faeces determined by the membrane filtration technique with an improved medium. J. Bact. 74:591-595.

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COLUMBIA AGAR BASE

Medium for fastidious microorganisms isolation from clinical samples.

TYPICAL FORMULA (g/l)

Peptospecial	23.0
Starch	1.0
Sodium Chloride	5.0
Agar	14.0
Final pH = 7.3 ± 0.2 at 25 °C.	

DIRECTIONS

Suspend 43.0 g of powder in 1 liter of distilled or deionized water. Heat to boiling until completely dissolved. Sterilize in autoclave at 121 °C for 15 minutes. Cool to 45-50 °C and aseptically add 5% defibrinated sterile sheep blood. Mix well. Dispense in petri dishes.

Columbia Agar Base can be also enriched in various way:

- with 2 vials of CNA (Staf / Strep) supplement (colistin sulphate 5 mg/vial, nalidixic acid 8 mg/vial, code 81048), each one reconstituted with 5 ml of sterile distilled water; final medium will contain colistin sulphate 10 mg/l and nalidixic acid 16 mg/l.
- with 2 vials of *Gardnerella vaginalis* supplement (gentamicin 3 mg/vial, amphotericin B 1mg/vial, nalidixic acid 15 mg/vial, code 81040), each one reconstituted with 5 ml of a 1:1 solution of ethyl alcohol and sterile distilled water; final medium will contain gentamicin 6 mg/l, amphotericin B 2 mg/l and nalidixic acid 30 mg/l.

DESCRIPTION

COLUMBIA AGAR BASE, enriched with sterile sheep blood (5%), is suitable for isolation and growth of fastidious microorganisms such as streptococci, staphylococci, pneumococci and listeriae from clinical samples.

TECHNIQUE

Inoculate the medium with the specimen streaking by a sterile loop and incubate at 36 ± 1 °C for 18-48 hours aerobically, anaerobically or under conditions of increased CO₂ (5-10%), in accordance with established laboratory procedures.

Examine plates for growth and hemolytic reactions. Four types of hemolysis on blood agar media can be described:

1. α-hemolysis is the reduction of hemoglobin to methemoglobin in the medium surrounding the colony, causing a greenish discolorization of the medium.
2. β-hemolysis is the lysis of red blood cells, producing a clear zone surrounding the colony.
3. γ-hemolysis indicates no destruction of red blood cells and no change in the color of the medium.
4. δ-hemolysis indicates a partial lysis.

QUALITY CONTROL

Dehydrated medium

Appearance: free-flowing, homogeneous.

Color: beige.

Prepared medium

Appearance: opaque.

Color: cherry red.

Incubation conditions: 36 ± 1 °C for 18-48 hours at 5-10% CO₂.

Microorganism	ATCC	Growth	Characteristics
<i>Streptococcus pyogenes</i>	19615	good	β-hemolysis
<i>Streptococcus pneumoniae</i>	6303	good	α-hemolysis
<i>Staphylococcus aureus</i>	25923	good	β-hemolysis
<i>Gardnerella vaginalis</i>	14018	good	β-hemolysis

LIOFILCHEM s.r.l.

Via Scozia, Zona Ind.le - 64026, Roseto D.A. (TE) - ITALY

Phone +390858930745 Fax +390858930330

Website: www.liofilchem.net E-mail: liofilchem@liofilchem.net

PERFORMANCE AND LIMITATIONS

When this medium is enriched with 10% sterile sheep blood, heated at 80 °C for 10 minutes until a chocolate color is obtained, and an antibiotic mixture is added (vancomycin, colimycin, trimethoprim, amphoterycin B) it is suitable for the selective isolation of the pathogens neisseria. If used without the addition of blood, the medium is suitable for growing of *Brucella abortus*, *Yersinia pestis*, *Clostridium perfringens* and *enterobacteria*. Hemolytic reactions of some strains of Group D streptococci have been shown to be affected by differences in animal blood. Such strains are beta –hemolytic on horse and rabbit blood agar and alpha-hemolytic on sheep blood agar.

STORAGE

The powder is very hygroscopic: store the powder at 10-30 °C, in a dry environment, in its original container tightly closed and use it before the expiry date on the label or until signs of deterioration or contamination are evident.
Store prepared plates at 2-8 °C.

REFERENCES

1. Ellner, P.D., C.J. Stoessel., E. Drakeford, and F. Vasi (1966). A new culture medium for medical bacteriology. Am. J.Clin. Path. **45**, 502-504.
2. Isenberg, H.D. (ed.) (1992). Clinical microbiology procedures handbook, vol. 1 American Society for Microbiology, Washington, DC.

PRESENTATION

Product	REF	
COLUMBIA AGAR BASE (11.6 l)	610013	500 g
COLUMBIA AGAR BASE (2.3 l)	620013	100 g
COLUMBIA AGAR BASE (116.2 l)	6100135	5 Kg
SHEEP BLOOD DEFIBRINATED	83296	50 ml
CNA (Staf / Strep) supplement	81048	10 vials
Gardnerella vaginalis supplement	81040	10 vials

TABLE OF SYMBOLS

 LOT Batch code	 Caution, consult accompanying documents	 Manufacturer	 Contains sufficient for <n> tests	 IVD In Vitro Diagnostic Medical Device
 REF Catalogue number	 Fragile, handle with care	 Use by	 Temperature limitation	 Keep away from heat source

**LIOFILCHEM s.r.l.**

Via Scozia, Zona Ind.le - 64026, Roseto D.A. (TE) - ITALY

Phone +390858930745 Fax +390858930330

Website: www.liofilchem.net E-mail: liofilchem@liofilchem.net

ENDO AGAR

Medium for coliforms confirmatory test.

TYPICAL FORMULA (g/l)

Peptone	10.0
Lactose	10.0
Dipotassium Phosphate	3.5
Agar	15.0
Sodium Sulphite	2.5
Basic Fuchsin	0.5
Final pH = 7.5 ± 0.2 at 25 °C.	

DIRECTIONS

Suspend 41.5 g of powder in 1 liter of distilled or deionized water. Heat to boiling with frequent and careful overturnings until complete dissolution. Autoclave at 121 °C for 15 minutes. Evenly disperse the precipitate when dispensing. Use immediately.

DESCRIPTION

ENDO AGAR is used for confirming the presence of coliforms organisms.

TECHNIQUE

For the confirmation of presumptive tests with liquid media, subculture tubes showing gas, or acid and gas formation, onto an Endo Agar plate. Incubate at 36 ± 1 °C for 24 hours. Lactose fermenting coliforms (e.g. *E. coli*) give rise to deep red colonies which color the surrounding medium and possess a golden metallic sheen. Non-lactose fermenters form colorless translucent colonies, against the pink to colorless medium.

QUALITY CONTROL

Dehydrated medium

Appearance: free-flowing, homogeneous.

Color: medium purple.

Prepared medium

Appearance: opalescent with precipitates.

Color: pink.

Incubation conditions: 36 ± 1 °C for 24 ± 2 hours.

Microorganism	ATCC	Growth	Characteristics
<i>Staphylococcus aureus</i>	25923	markedly to completely inhibited	
<i>Escherichia coli</i> □ □	25922	good	red colonies w / green metallic sheen
<i>Salmonella typhimurium</i>	14028	good	colorless to pink colonies

PERFORMANCE AND LIMITATIONS

If the medium is to be used the same day it is rehydrated, it does not need to be autoclaved. Boil to dissolve completely before dispensing into plates.

STORAGE

The powder is very hygroscopic: store the powder at 10-30 °C, in a dry environment, in its original container tightly closed and use it before the expiry date on the label or until signs of deterioration or contamination are evident. The medium should be used the day it is prepared: if it is necessary store in the dark at 2-8 °C for no more than 3 days.

REFERENCES

- Endo, S. (1904). Über ein Verfahren zum Nachweis der Typhusbacillen. Centr. Bakt., Abt 1, Orig. 35:109-110.
- American Public Health Association. (1975). Standard methods for the examination of water and wastewater, 14th ed.

PRESENTATION

Product	REF	Σ
ENDO AGAR (12.0 l)	610020	500 g
ENDO AGAR (2.4 l)	620020	100 g

TABLE OF SYMBOLS

LOT Batch code	 Caution, consult accompanying documents	 Manufacturer	 Contains sufficient for <n> tests	 Keep away from heat source
REF Catalogue number	 Fragile, handle with care	 Use by	 Temperature limitation	

LIOFILCHEM s.r.l.

Via Scozia, Zona Ind.le - 64026, Roseto D.A. (TE) - ITALY

Phone +390858930745 Fax +390858930330

Website: www.liofilchem.net E-mail: liofilchem@liofilchem.net

KLIGLER IRON AGAR

Differential medium for enterobacteria identification.

TYPICAL FORMULA (g/l)

Proteose Peptone	20.0
Sodium Chloride	5.0
Yeast Extract	3.0
Meat Extract	3.0
Ferrous Sulfate	0.2
Sodium Thiosulphate	0.3
Lactose	10.0
Glucose	1.0
Phenol Red	0.024
Agar	11.0
Final pH = 7.4 ± 0.2 at 25 °C.	

DIRECTIONS

Suspend 53.5 g of powder in 1 liter of distilled or deionized water. Heat to boiling until completely dissolved. Dispense into final tubes. Sterilize in autoclave at 121°C for 15 minutes. Cool in a slanting position.

DESCRIPTION

KLIGLER IRON AGAR is a solid medium used to distinguish between *Enterobacteriaceae* on the basis of their ability to ferment lactose and / or glucose and to produce hydrogen sulphide.

TECHNIQUE

Inoculate by stabbing the butt and abundantly streaking the slope. Incubate at 36 ± 1°C for 18-24 hours and check the color of the medium both in the butt and at the slope. Also check for the presence of gas in the butt and the presence of the black precipitate (H₂S).

QUALITY CONTROL

Dehydrated medium

Appearance: free-flowing, homogeneous.

Color: pinkish beige.

Prepared medium

Appearance: slightly opalescent, slight precipitate.

Color: slightly orange-red.

Incubation conditions: 36 ± 1°C for 18-24 hours.

Microorganism	ATCC	Growth	Slant/butt	Gas	H ₂ S
<i>Citrobacter freundii</i>	8090	good	acid/acid	+	+
<i>Escherichia coli</i>	25922	good	acid/acid	+	-
<i>Proteus vulgaris</i>	6380	good	alkaline/acid	-	+

PERFORMANCE AND LIMITATIONS

A pure culture is essential when inoculating Kligler Iron Agar. If inoculated with a mixed culture, irregular observations may occur.

STORAGE

The powder is very hygroscopic: store the powder at 10-30 °C, in a dry environment, in its original container tightly closed and use it before the expiry date on the label or until signs of deterioration or contamination are evident.

Store prepared tubes at 2-8°C.

REFERENCES

1. MacFaddin, J.F. (1976). Biochemical tests for identification of medical bacteria.
2. Kligler, I.J. (1918). J. Exp. Med. **28**: 319-322.

LIOFILCHEM s.r.l.

Via Scozia, Zona Ind.le - 64026, Roseto D.A. (TE) - ITALY

Phone +390858930745 Fax +390858930330

Website: www.liofilchem.net E-mail: liofilchem@liofilchem.net

PRESENTATION

Product	REF	
KLIGLER IRON AGAR (9.3 l)	610023	500 g
KLIGLER IRON AGAR (1.8 l)	620023	100 g

TABLE OF SYMBOLS

LOT Batch code	 Caution, consult accompanying documents	 Manufacturer	 Contains sufficient for <n> tests	IVD <i>In Vitro</i> Diagnostic Medical Device
REF Catalogue number	 Fragile, handle with care	 Use by	 Temperature limitation	 Keep away from heat source

LIOFILCHEM s.r.l.

Via Scozia, Zona Ind.le - 64026, Roseto D.A. (TE) - ITALY

Phone +390858930745 Fax +390858930330

Website: www.liofilchem.net E-mail: liofilchem@liofilchem.net

MUELLER HINTON AGAR

Medium for susceptibility test (Kirby-Bauer method).

TYPICAL FORMULA (g/L)

Meat Extract	2.0
Casamino Acids, Technical	17.5
Starch	1.5
Agar	15.0
Final pH 7.3 ± 0.1	

DESCRIPTION

MUELLER HINTON AGAR is used for antimicrobial susceptibility testing of rapidly growing aerobic microorganisms by the disk diffusion technique.

PRINCIPLE

Casamino acids and meat extract are a source of amino acids, nitrogen, minerals, vitamins, carbon and other factors which increase the growth of microorganisms. Starch acts as a protective substance against toxic molecules which can be present in the medium. Hydrolysis of starch during sterilization supplies a little amount of glucose which represents a source of energy. Agar is the solidifying agent. Kirby-Bauer method is based on the diffusion, through the agar, of antimicrobial substances which soaks paper disks: microorganism growth shows an inhibition halo around the disk and the diameter of the halo is correlated to the Minimal Inhibiting Concentration (MIC).

PREPARATION

Suspend 36.0 g of powder in 1 litre of distilled or deionized water. Heat to boiling and shake until completely dissolved. Sterilize in autoclave at 121°C for 15 minutes. Dispense in final containers.

TECHNIQUE

Transfer 4-5 colonies in an appropriate broth.

Place it in a 37°C incubator until an opacity is obtained equivalent to the standard opacity of 0.5 on the MacFarland scale. Introduce a sterile swab into the inoculum and inoculate the agar passing 2 or 3 times onto the entire surface.

Press the disk containing the antimicrobial on the agar surface.

Incubate at 36±1°C for 18 hours, measure the inhibition zone with a compass and compare to the NCCLS recommended zone ranges.

INTERPRETATION OF RESULTS

Compare obtained values of inhibition halo diameter with the values reported on NCCLS M100(M2) document.

STORAGE

10-30°C away from light, until the expiry date on the label or until signs of deterioration or contamination are evident.

WARNING and PRECAUTIONS

The product is not classified as hazardous by current legislation and does not contain harmful substances in concentrations of ≥1%. The product is designed for *In vitro* diagnostic use and must be used only by properly trained operators.

DISPOSAL of WASTE

Disposal of waste must be carried out according to national and local regulations in force.

REFERENCES

1. Bauer et al. (1966). J. Clin. Pathol. 45:493-496.
2. Mueller, J.H., and Hinton. 1941. Proc. Soc. Exp. Biol. Med. 48: 330-333.
3. NCCLS. Performance standards for susceptibility testing; Twelve Informational Supplement. NCCLS Document M100-S12, January 2002.

Liofilchem s.r.l.

Via Scozia-Zona industriale - 64026 Roseto degli Abruzzi Tel. +39.085.8930745 - Fax +39.085.8930330
Web site: <http://www.liofilchem.net> E-mail: liofilchem@liofilchem.net

PRODUCT SPECIFICATIONS

NAME

MUELLER HINTON AGAR

PRESENTATION

Dehydrated culture medium

STORAGE

10-30°C

PACKAGING

Code	Content	Packaging
610033	500 g	500 g of powder in plastic bottle
620033	100 g	100 g of powder in plastic bottle
6100335	5 kg	5 kg of powder in plastic container

pH OF THE MEDIUM

7.3 ± 0.1

USE

MUELLER HINTON AGAR is used for antimicrobial susceptibility testing of rapidly growing aerobic microorganisms by the disk diffusion technique.

TECHNIQUE

Refer to technical sheet of the product.

APPEARANCE OF THE MEDIUM

Amber medium, slightly opalescent.

SHELF LIFE

4 years

QUALITY CONTROL

- Control of general characteristics, label and print
- Sterility control
7 days at 25 ± 1°C, in aerobiosis
7 days at 36 ± 1°C, in aerobiosis
- Microbiological control
Incubation conditions: 18-24 h at 36 ± 1°C

Microorganism	ATCC	Growth	Characteristics
<i>Enterococcus faecalis</i>	ATCC 29212	Good	White colonies
<i>Escherichia coli</i>	ATCC 25922	Good	Colorless colonies
<i>Proteus mirabilis</i>	ATCC 25933	Good	Colorless colonies
<i>Staphylococcus aureus</i>	ATCC 25923	Good	White colonies

TABLE of SYMBOLS

IVD In vitro Diagnostic Medical Device	LOT Batch code	 Manufacturer	 Contains sufficient for <n> tests
REF Catalogue number	 Temperature limitation	 Use by	 Caution, consult accompanying documents

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S.S. AGAR (MODIFIED)

Selective medium for the isolation of *Salmonella* spp. and *Shigella* spp.

TYPICAL FORMULA	(g/l)
Peptone	5.5
Meat Extract	5.0
Lactose	10.0
Sodium Tiosulfate	8.5
Yeast Extract	5.0
Sodium Citrate	1.0
Bile Salts N.3	1.5
Ferric Ammonium Citrate	1.5
Brilliant Green	0.33 mg
Neutral Red	0.025
Agar	14.0
Final pH 7.0 ± 0.2	

DESCRIPTION

S.S. AGAR (MODIFIED) is a highly selective medium for the isolation of *Salmonella* spp. and some species of *Shigella* from clinical specimens and food.

PRINCIPLE

Gram-positive microorganisms and coliforms are inhibited by selective components: brilliant green, bile salts n.3, sodium tiosulfate and citrate. The differentiations of microorganisms is obtained through the introduction of lactose in the medium. Lactose fermented bacteria cause acidification, thus formation of red colonies for the presence of neutral red. Not-fermented microorganisms form instead colourless colonies. Sodium tiosulfate in combination with iron acts as indicator for sulphur production causing the blackening of the colony center.

PREPARATION

Suspend 52.0 g of the powder in 1 litre of distilled or deionized water. Mix well. Heat to boil shaking frequently until dissolved completely. DO NOT AUTOCLAVE. Cool to 45-50°C. In aseptic conditions dispense in Petri dishes and let solidify the medium with the lids of the plates partially removed.

TECHNIQUE

Inoculate the plate streaking the sample onto the agar surface to isolate pure colonies from samples containing a mixed flora. Incubate at 36±1°C for 18-24 hours.

INTERPRETATION OF RESULTS

Salmonella spp. and other lactose not-fermented microorganisms can produce opaque, translucent or transparent colonies, with or without black center. *Shigella* colonies are colourless. The few lactose fermented microorganisms, that are able to growth on the medium, show reddish mucoid colonies.

STORAGE

10-30°C away from light, until the expiry date on the label. Eliminate if signs of deterioration or contamination are evident.

WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is designed for *in vitro* diagnostic use only and must be used by properly trained operators.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to the national and local regulations in force.

REFERENCES

1. Gray L.D. (1995). *Escherichia*, *Salmonella*, *Shigella* and *Yersinia*, p. 450-456. In *Manual of clinical microbiology*, 6th ed. American society of microbiology.
2. Leifson E. (1935). *J. Pathol. Bacteriol.* 40: 581.
3. Rose, H.M., and M.H. Kolodny (1942). *J. Lab. Clin. Med.* 27: 1081-1083.

LIOFILCHEM® S.r.l.

Via Scozia, Zona Ind.le - 64026, Roseto degli Abruzzi (TE) - ITALY
Tel +39 0858930745 Fax +39 0858930330 Website: www.liofilchem.net E-mail: liofilchem@liofilchem.net

PRODUCT SPECIFICATIONS

NAME

S.S. AGAR (MODIFIED)

PRESENTATION

Dehydrated medium

STORAGE

10-30°C

PACKAGING

Ref.	Content	Packaging
610042	500 g	500 g of powder in plastic bottle
620042	100 g	100 g of powder in plastic bottle
6100425	5 kg	5 kg of powder in plastic container

pH OF THE MEDIUM

7.0 ± 0.2

USE

S.S. AGAR (MODIFIED) is a highly selective medium for the isolation of *Salmonella* spp. and some species of *Shigella* from clinical specimens and foods

TECHNIQUE

Refer to technical sheet of the product

APPEARANCE OF THE MEDIUM

Dehydrated medium

Appearance: free-flowing, homogeneous

Colour: light-pink

Prepared medium

Appearance: opalescent

Colour: purple

SHELF LIFE

4 years

QUALITY CONTROL

- Control of general characteristics, label and print
- Microbiological control
Inoculum for productivity: 10-100 UFC/ml
Inoculum for selectivity: 10⁴-10⁵ UFC/ml
Inoculum for specificity: ≤10⁴ UFC/ml
Incubation Conditions: 18-24 h at 35 ± 2°C, in aerobiosis

Microorganism

Shigella flexneri
Salmonella typhimurium
Enterococcus faecalis
Staphylococcus aureus
Escherichia coli

ATCC® 12022
ATCC® 14028
ATCC® 29212
ATCC® 25923
ATCC® 25922

Growth

Good
Good
Inhibited
Inhibited
Partially Inhibited

Features

Colourless colonies
Colourless colonies with or without black center

Pink or red colonies

TABLE OF SYMBOLS

LOT	Batch code	IVD	In vitro Diagnostic Medical Device		Manufacturer		Use by
REF	Catalogue number		Temperature limitation		Contains sufficient for <n> tests		Caution, consult instructions for use

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TYPICAL FORMULA	(g/l)
Magnesium Sulfate	0.2
Ammonium Dihydrogen Phosphate	1.0
Dipotassium Phosphate	1.0
Sodium Citrate	2.0
Sodium Chloride	5.0
Brom Thymol Blue	0.08
Agar	15.0

Final pH = 6.8 ± 0.2 at 25 °C.

DIRECTIONS

Suspend 24.3 g of powder in 1 liter of distilled or deionized water. Heat to boiling until completely dissolved. Dispense into final tubes and sterilize in the autoclave at 121°C for 15 minutes. Allow the medium to solidify in a slant position.

DESCRIPTION

SIMMONS CITRATE AGAR is recommended for the differentiation and identification of *Enterobacteriaceae* on the basis of citrate utilization.

TECHNIQUE

Inoculate the medium with the specimen by stabbing the butt and streaking the slope. Incubate at 36 ± 1 °C for 24-48 hours. Organisms able to utilize ammonium dihydrogen phosphate and sodium citrate as the sole sources of nitrogen and carbon respectively will grow on this medium and produce an alkaline reaction as evidenced by a change in the color of the bromthymol blue indicator from green (neutral) to blue (alkaline).

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Color: yellow, may have green tinge.

Prepared medium

Appearance: slightly opalescent, may have a slight precipitate.

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Website: www.liofilchem.net E-mail: liofilchem@liofilchem.net

REFERENCES

1. Ewig W.H. and Edwards P.R. (1960). Bull. Bact. Nomen. And Taxon. **10**:1-12.
2. American Public Health Association (1981). Standard Methods for the Examination of Water and Wastewater, 15th ed. APHA Inc, Washington DC.
3. Matsen J.M., and Sherris J.C. (1969) Appl. Microbiol. **18**: 452-454.

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Website: www.liofilchem.net E-mail: liofilchem@liofilchem.net

SCHAEDLER BROTH

Basal medium for aerobic and anaerobic bacteria isolation.

TYPICAL FORMULA (g/l)

Tryptic Soy Broth	10.0
Tryptone	5.0
Yeast Extract	5.0
Glucose	5.0
Hemin	0.01
L-Cystine	0.4

Final pH = 7.6 ± 0.2 at 25 °C.

DIRECTIONS

Suspend 25.4 g of powder in 1 liter of distilled or deionized water. Heat to boiling until completely dissolved. Dispense into final containers. Sterilize in the autoclave at 121 °C for 15 minutes.

DESCRIPTION

SCHAEDLER BROTH is a nonselective medium used for fastidious aerobic and anaerobic bacteria, such as lactobacilli, streptococci, clostridia and *Bacteroides*. The medium is highly nutritious due to its content of peptones, glucose and yeast extract. Hemin supplies the X factor required by many fastidious microorganisms.

TECHNIQUE

Cool the medium at 25°C. If the medium is not used immediately after its preparation, it must be regenerated at 100°C for 20 minutes in order to restore conditions of anaerobiosis. Do not repeat this operation more than once. Inoculate the material to examine into the tube using a sterile bent glass rod. Incubate at 36±1°C for 18-48 hours in aero/anaerobic conditions required by the species being studied.

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Color: light tan.

Prepared medium

Appearance: clear to slightly opalescent, may have a slight precipitate.

Color: light to medium amber.

Incubation conditions: 36 ± 1 °C for 18-48 hours.

Microorganism	ATCC	Growth
<i>Clostridium perfringens</i>	13124	good
<i>Streptococcus pyogenes</i>	19615	good
<i>Escherichia coli</i>	25922	good
<i>Staphylococcus aureus</i>	25923	good

PERFORMANCE AND LIMITATIONS

Since the nutritional requirements of organisms are different, some strains may be encountered that fail to grow or grow poorly on this medium.

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Store prepared tubes at 2-8 °C.

REFERENCES

- MacFaddin, J.D. (1985). Media for isolation-cultivation-identification-maintenance of medical bacteria, p. 695-699, vol. 1. Williams & Wilkins, Baltimore, MD.
- Schaedler, R.W., R. Dubos, and R. Costello. (1965). The development of bacterial flora in the gastrointestinal tract of mice. J. Exp. Med. **122**: 59.

PRESENTATION

Product	REF	
SCHAEDLER BROTH (19.6 l)	610137	500 g
SCHAEDLER BROTH (3.9 l)	620137	100 g

TABLE OF SYMBOLS

 Batch code	 Caution, consult accompanying documents	 Manufacturer	 Contains sufficient for <n> tests	 In Vitro Diagnostic Medical Device
 Catalogue number	 Fragile, handle with care	 Use by	 Temperature limitation	 Keep away from heat source

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EC DECLARATION OF CONFORMITY

The company Liofilchem® S.r.l., registered office in Via Scozia, 64026 Roseto degli Abruzzi (TE) Italy, as a manufacturer of the *in vitro* medical-diagnostic device listed in the attached table, Revision 31.0 of 08.01.2016

hereby certifies under its own responsibility

1. that the above mentioned device complies with all the applicable provisions of Directive 98/79/EC (Annex III) and its relevant transposition into national law;
2. the above mentioned is not included in Annex II, List A and B of Directive 98/79/EC;
3. that the technical documentation referred to at Annex III of the Directive 98/79/EC is available for the national authorities in its facility and that this documentation shall be kept for 5 years after the last product has been manufactured;
4. that the manufacturing process follows suitable principles of quality assurance;
5. that, has implemented and keep up to date, a post-production surveillance system for monitoring the products;
6. that the device in question, was introduced into the market provided with CE mark.

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Rev. 31.0 del 08.01.2016

621021	HEART INFUSION BROTH
621022	MIDDLEBROOK 7H10 AGAR BASE
621210	WURTZ LACTOSE AGAR
621265	ISOSENSITEST AGAR
621367	BILE BACTERIOLOGICAL
621401	IRON SULPHITE AGAR
621402	CARY BLAIR TRANSPORT MEDIUM
621601	GLUCOSE
621618	CHROMATIC™ MH
621619	CHROMATIC™ CRE AGAR BASE
621701	PEPTONE BACTERIOLOGICAL
622202	STREPTOCOCCUS SELECTIVE AGAR
630026	LOWENSTEIN JENSEN MEDIUM w GLYCEROL 1 litre
71618	ENTEROSYSTEM 18R 20 Tests
71630	STAF SYSTEM 18 R 20 Tests
71670	COPRO SYSTEM 40 Tests
71675	COPRO SYSTEM Plus 20 Tests
71678	PATHOGENIC SYSTEM DOUBLE 40 Tests
71679	PATHOGENIC SYSTEM 20 Tests
71681	PATHOGENIC SYSTEM AST
71714	INTEGRAL SYSTEM ENTEROBATTERI 20 Tests
71718	INTEGRAL SYSTEM STAFILOCOCCI 20 Tests
71720	INTEGRAL SYSTEM STREPTOCOCCI 20 Tests
71724	INTEGRAL SYSTEM GARDNERELLA 20 TESTS
71822	INTEGRAL SYSTEM YEASTS Plus 20 Tests
72560	STREPTO SYSTEM 12 R 40 Tests
72592	MYCOPLASMA SYSTEM Plus 20 Tests
74156	A.F. GENITAL SYSTEM 20 Tests
74160	URIN SYSTEM Plus 20 Tests
74161	URIN SYSTEM Chrom 20 Tests
76010	Sensi Test gram-negative 20 Tests
76020	Sensi Test gram-positive 20 Tests
76031	SensiQuattro Gram-negative 20 Tests
76032	SensiQuattro Gram-positive 20 Tests
76033	SensiQuattro Candida EU 20 Tests
78618	ENTERO PLURI TEST 10 Tests
78619	ENTERO PLURI TEST 25 Tests
78620	OXI/FERM PLURI TEST 10 Tests
78621	OXI/FERM PLURI TEST 25 Tests
79010	Sensi Test gram-negative 4 Tests
79020	Sensi Test gram-positive 4 Tests
79031	SensiQuattro Gram-negative 4 Tests
79032	SensiQuattro Gram-positive 4 Tests
79033	SensiQuattro Candida EU 4 Tests
79156	A.F. GENITAL SYSTEM 4 Tests
79160	URIN SYSTEM Plus 4 Tests
79161	URIN SYSTEM Chrom 4 Tests
79560	STREPTO SYSTEM 12 R 8 Tests
79592	MYCOPLASMA SYSTEM Plus 4 Tests
79618	ENTEROSYSTEM 18R 4 Tests
79630	STAF SYSTEM 18 R 4 Tests
79670	COPRO SYSTEM 8 Tests
79675	COPRO SYSTEM Plus 4 Tests
79678	PATHOGENIC SYSTEM DOUBLE 8 Tests
79679	PATHOGENIC SYSTEM 4 Tests
79681	PATHOGENIC SYSTEM AST
79714	INTEGRAL SYSTEM ENTEROBATTERI 4 Tests
79718	INTEGRAL SYSTEM STAFILOCOCCI 4 Tests
79720	INTEGRAL SYSTEM STREPTOCOCCI 4 Tests
79724	INTEGRAL SYSTEM GARDNERELLA 4 Tests

79822	INTEGRAL SYSTEM YEASTS Plus 4 Tests
80009	IODINE MKTT SOLUTION 10 x 10 ml
80010	XLT 4 supplement 2 x 50 ml
80021	GLYCEROL supplement 4 x 50 ml
80022	POTASSIUM TELLURITE 1% suppl. 5 x 10 ml
80031	TWEEN 80 supplement 2 x 50 ml
80040	CHROMATIC™ SALMONELLA Supplement 2x50 ml
80047	MULLER KAUFFMANN 3X50 ml (Iodio/B.G.O.1%)
80053	VITAMIN K 1% supplement 5 x 5 ml
80056	LEGIONELLA growth supplement 10 vials
80057	H2O2 REAGENT 1 x 10 ml
80060	DECONTAM-KIT
80110	UREA 40% 6X100 ml
80219	EGG YOLK emulsion 4 x 50 ml
80252	ENTEROSYSTEM 18R REAGENT 100/200 Tests
80253	COPRO SYSTEM REAGENTS (antisera)
80257	LISTERIA SYSTEM 18R -REAG 100/200 Tests
80258	AF GENITAL SYSTEM REAGENT
80260	IDENTIF. SYSTEM-REAGENT 100/200 Tests
80271	KOVAC'S REAGENT 4x25 ml
80272	FERRIC CHLORIDE 10% 2x 25 ml
80273	NINHYDRIN 7% 10 ml
80275	MIF COLOR KIT 50 Tests
80276	ZIEHL-NEESEN 3 x 250 ml
80277	METHYLENE BLUE Solution 250 ml
80279	VASELINE OIL 4 x 50 ml
80280	V.P. TEST-Reagent 10x10ml
80281	V.P. TEST EP 10 x 10 mL
80282	Kit May-Grünwald Giemsa
80290	SAFRANIN SOLUTION 1000 ml
80291	POTASSIUM TELLURITE 3.5% suppl.5x10 ml
80292	UREA 40 % supplement 10 x 5 ml
80293	GRAM COLOR KIT 4 x 250 ml
80294	KIT COLOR ALBERT 2 x 250 ml
80295	DECOLOURIZING SOLUTION 1000 ml
80296	LUGOL PVP SOLUTION 1000 ML
80297	SAFRANIN SOLUTION 500 ml
80298	LUGOL PVP SOLUTION 250 ml
80299	CRYSTAL VIOLET SOLUTION 1000 ml
80300	TTC 1% supplement 5 x 10 ml
80350	ANTIBIOTIC TEST
80351	RAPID ANTIBIOTIC TEST 50 Tests
80380	KINYOUN COLOR KIT 2 x 250 ml
80390	FIXUR 1
80409	IODINE SOLUTION 10 x 10 ml
80410	XLT 4 SUPPLEMENT 4 x 50 ml
80422	POTASSIUM TELLURITE 1% Supplement 10 x 10 ml
80430	TTC 1% supplement 10 x 10 ml
80431	TWEEN 80 Supplement 4 x 50 ml
80453	VITAMIN K 1% SUPPLEMENT 10 x 5 ml
80491	POTASSIUM TELLURITE 3,5% Supplement 10 x 10
81001	AMPICILLIN supplement 10 vials
81002	LEGIONELLA (BMPA) supplement 10 vials
81003	BRUCELLA supplement 10 vials
81004	CAMPYLOBACTER Preston supplem 10 vials
81006	CN (Pseudomonas) supplement 10 vials
81007	CLOSTRIDIUM difficile suppleme 10 vials
81008	LEGIONELLA (GVPC) supplement 10 vials
81009	IODINE solution 5 x 10 ml
81011	CLOSTRIDIUM perfringens (T.S.C.) sup.10 v.

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91200	DISC DISPENSER 8 CARTRIDGES
91203	DISC DISPENSER 6 CARTRIDGES
92000	AMOX*/SULB 2/1 AXS 0.016-256* 30 MIC Tests
920000	AMOX*/SULB 2/1 AXS 0.016-256* 100 MIC Tests
92001	RIFAMPICIN RD 0.002-32 30 MIC Tests
920010	RIFAMPICIN RD 0.002-32 100 MIC Tests
920011	RIFAMPICIN RD 0.002-32 10 MIC Tests
92002	FUSIDIC ACID FU 0.016-256 30 MIC Tests
920020	FUSIDIC ACID FU 0.016-256 100 MIC Tests
920021	FUSIDIC ACID FU 0.016-256 10 MIC Tests
92003	AMPICILLIN AMP 0.016-256 30 MIC Tests
920030	AMPICILLIN AMP 0.016-256 100 MIC Tests
920031	AMPICILLIN AMP 0.016-256 10 MIC Tests
92004	POLYMYXIN B PB 0.064-1024 30 MIC Tests
920040	POLYMYXIN B PB 0.064-1024 100 MIC Tests
920041	POLYMYXIN B PB 0.064-1024 10 MIC Tests
92005	CEFPODOXIME PX 0.016-256 30 MIC Tests
920050	CEFPODOXIME PX 0.016-256 100 MIC Tests
920051	CEFPODOXIME PX 0.016-256 10 MIC Tests
92006	CEFOTAXIME CTX 0.016-256 30 MIC Tests
920060	CEFOTAXIME CTX 0.016-256 100 MIC Tests
920061	CEFOTAXIME CTX 0.016-256 10 MIC Tests
92007	CEFOTAXIME CTX 0.002-32 30 MIC Tests
920070	CEFOTAXIME CTX 0.002-32 100 MIC Tests
920071	CEFOTAXIME CTX 0.002-32 10 MIC Tests
92008	CEFPIROME CR 0.016-256 30 MIC Tests
920080	CEFPIROME CR 0.016-256 100 MIC Tests
920081	CEFPIROME CR 0.016-256 10 MIC Tests
92009	GENTAMICIN CN 0.016-256 30 MIC Tests
920090	GENTAMICIN CN 0.016-256 100 MIC Tests
920091	GENTAMICIN CN 0.016-256 10 MIC Tests
92010	GENTAMICIN CN 0.064-1024 30 MIC Tests
920100	GENTAMICIN CN 0.064-1024 100 MIC Tests
920101	GENTAMICIN CN 0.064-1024 10 MIC Tests
92011	GATIFLOXACIN GAT 0.002-32 30 MIC Tests
920110	GATIFLOXACIN GAT 0.002-32 100 MIC Tests
920111	GATIFLOXACIN GAT 0.002-32 10 MIC Tests
92012	TEICOPLANIN TEC 0.016-256 30 MIC Tests
920120	TEICOPLANIN TEC 0.016-256 100 MIC Tests
920121	TEICOPLANIN TEC 0.016-256 10 MIC Tests
92013	ENROFLOXACIN ENR 0.002-32 30 MIC Tests
920130	ENROFLOXACIN ENR 0.002-32 100 MIC Tests
920131	ENROFLOXACIN ENR 0.002-32 10 MIC Tests
92014	SPECTINOMYCIN SPC 0.064-1024 30 MIC Tests
920140	SPECTINOMYCIN SPC 0.064-1024 100 MIC Tests
920141	SPECTINOMYCIN SPC 0.064-1024 10 MIC Tests
92015	OXACILLIN OX 0.016-256 30 MIC Tests
920150	OXACILLIN OX 0.016-256 100 MIC Tests
920151	OXACILLIN OX 0.016-256 10 MIC Tests
92016	CEFTIZOXIME CZX 0.016-256 30 MIC Tests
920160	CEFTIZOXIME CZX 0.016-256 100 MIC Tests
920161	CEFTIZOXIME CZX 0.016-256 10 MIC Tests
92017	MECILLINAM MEC 0.002-32 30 MIC Tests
920170	MECILLINAM MEC 0.002-32 100 MIC Tests
920171	MECILLINAM MEC 0.002-32 10 MIC Tests
92018	AMIKACIN AK 0.016-256 30 MIC Tests
920180	AMIKACIN AK 0.016-256 100 MIC Tests
920181	AMIKACIN AK 0.016-256 10 MIC Tests
92019	BACITRACIN BA 0.016-256 30 MIC Tests
920190	BACITRACIN BA 0.016-256 100 MIC Tests

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96021	SALMONELLA PARATYPHI C TOTALE 20 ml
96022	BRUCELLA MELITENSIS 20 ml
96023	BRUCELLA SUIS 20 ml
96031	SALMONELLA TYPHI H SLIDE 5 ml
96032	SALMONELLA TYPHI O SLIDE 5 ml
96033	SALMONELLA TYPHI TOTALE 5 ml SLIDE
96034	SALMONELLA PARATYPHI AH SLIDE 5 ml
96035	SALMONELLA PARATYPHI AO 5 ml SLIDE
96036	SALMONELLA PARATYPHI A TOTALE 5ml SLIDE
96037	SALMONELLA PARATYPHI BH 5 ml SLIDE
96038	SALMONELLA PARATYPHI BO 5 ml SLIDE
96039	SALMONELLA PARATYPHI B TOTALE 5ml SLIDE
96040	SALMONELLA PARATYPHI CH 5 ml SLIDE
96041	SALMONELLA PARATYPHI CO 5 ml SLIDE
96042	SALMONELLA PARATYPHI C TOTALE 5 ml SLIDE
96043	BRUCELLA TOTALE SLIDE 5 ml SLIDE
96044	BRUCELLA ABORTUS 5 ml SLIDE
96045	BRUCELLA MELITENSIS SLIDE 5 ml
96046	BRUCELLA BENGAL ROSE SLIDE 5 ml
96047	PROTEUS OX2 5 ml SLIDE
96048	PROTEUS OX19 5 ml SLIDE
96049	PROTEUS OXK 5 ml SLIDE
96093	CONTROLLO NEGATIVO/NEGATIVE CONTROL 0.5ml
96096	POSITIVE CONTROL FOR SALMONELLA 0.5ml
96097	POSITIVE CONTROL FOR PROTEUS 0.5ml
96098	POSITIVE CONTROL FOR BRUCELLA 0.5ml
96142	Legionella Latex Kit
96143	CAMPYLOBACTER LATEX KIT
96144	CLOSTRIDIUM DIFFICILE LATEX KIT
96148	SHIGELLA ANTISERUM
96150	E. COLI O157 LATEX KIT
96151	SALMONELLA LATEX KIT
96153	STREPTO B LATEX KIT
96154	STREPTO A LATEX KIT
96155	BENCE JONES LATEX TEST
96316	Clostridium difficile GDH Card
96317	Clostridium Difficile Toxin A+B Card
96318	Giardia Card
96319	Listeria Monocytogenes Card
96320	Salmonella Ag Card
96321	O157 E.coli Card
96401	ONE STEP AMP DRMG SCREEN 20 CARDS
96404	ONE STEP COC DRMG SCREEN
96405	ONE STEP THC DRMG SCREEN
96406	ONE STEP M-AMP DRMG SCREEN 20 CARDS
96411	ONE STEP BRUPRENORPHINE DRMG SCREEN 20 CARDS
96415/20	FECAL OCCULT BLOOD CARD
96418	STREPTO A CARD 30 CARDS
96441	Gonorrhea Ag Card
96442	Gardnerella Vaginalis Card
96443	Trichomonas Vaginalis Card
96444	B.J. Free Kappa/Lambda Dipstick
96455	H.PYLORI CARD 20 CARD
96460	HCG URINE/SERUM CARD 50 CARD
96461	HCG URINE/SERUM CARD 100 CARD
96462	MICROALBUMIN CARD URINE 20 Cards
96465	AFP -ALFA FETO CARD 20 CARDS
96468	TUBERCOLOSI CARD 20 CARDS
96480	IgE TOTAL CARD
96485	CEA CARD 20 Cards

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Rev. 31.0 del 08.01.2016

96487	MYOGLOBIN
96488	TROPONIN 20 CARDS
96490	FERRITIN CARD
96495	SIFILIDE CARD 20 CARDS
96498	IM MONONUCLEOSIS INFECTION 20 CARDS
96590	URINE STRIP
96900	GIOTTO READER
96909	BIOMIC V3
96914	BIOMIC V3 AST
96915	BIOMIC V3 ID
96916	BIOMIC V3 CC
96919	AST Software
96931	ID Software
96932	CC Software
96933	Micropiastre 96 pozzetti Software

97800	ROTASTICK ONE STEP KIT 20 Tests
97801	RSV STICK ONE STEP 20 Tests
97802	ROTA/ADENO COMBI STICK ONE STEP 20 Tests
97803	H.PYLORI FECAL Ag ONE STEP 20 Tests
97805	STREP B STICK ONE STEP ASSAY 20 Tests
97807	ADENOSTICK ONE STEP ASSAY 20 Tests
9999	Blank Discs
99003	KPC&MBL disc kit (acc. to EUCAST)
99004	ESBL disc kit (acc. to EUCAST)
99005	ESBL disc kit (acc. to CLSI)
99006	ESBL (Chromos. Ind. AmpC) disc kit (acc. to EUCAST)
99007	KPC&MBL&OXA-48 disc kit (acc. to EUCAST)
99008	ESBL+AmpC screen disc kit
99009	AmpC disc kit

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