Chromatic™ VRE

ENGLISH

Chromogenic medium for the detection of Vancomycin Resistant Enterococci (VRE) directly from clinical specimens.

DESCRIPTION

Chromatic™ VRE is a chromogenic medium used for the detection of vancomycin resistant *Enterococcus faecium* and *Enterococcus faecalis* directly from clinical specimens in 24 hours.

Vancomycin Resistant Enterococci (VRE) were first isolated in Europe in the late 1980s and recently have been associated with nosocomial infections around the world.

Six types of vancomycin resistance have been characterized on both a phenotypic and genotypic basis in enterococci. Five of this types (VanA, B, C, D, E and G) correspond to acquired resistance, presumably by transfer from other bacteria in the intestinal tract; one type (VanC) is an intrinsic property of non-pathogenic enterococcal species, *E. gallinarum* and *E. casseliflavus*. The *vanA* and *vanB* operons are located on plasmids or in the chromosome, whereas the *vanC*, *vanD*, *vanE* and *vanG* operons have, thus far, been found only in the chromosome. VanA and VanB types of glycopeptide resistance have been associated with outbreaks of vancomycin resistant enterococci, mostly *E. faecium* and *E. faecalis*, and can be transferred to other organisms, including *S. aureus*. VanA-type strains display high levels of inducible resistance to both vancomycin and teicoplanin, whereas VanB-type has variable levels of inducible resistance to vancomycin only. The major concern of the emergence and dissemination of high-level resistance to vancomycin in enterococci is related to the risk for the spread of glycopeptide resistance to more pathogenic bacteria which in turn can lead to clinical isolates resistant to all antibiotics.

Chromatic™ VRE medium can be used to promptly detect VRE, differentiating between *E. faecalis* and *E. faecium*, reducing time and labour, increasing specificity and sensitivity and eliminating pre-enrichment and isolation of pure colonies.

TYPICAL FORMULA (g/l)

Peptone Mix	40.0
Sodium Chloride	5.0
Chromogenic and Selective Mix	6.5
Agar	15.0

Final pH 6.9 ± 0.2 at 25°C

METHOD PRINCIPLE

Peptone mix supply amino acids, nitrogen, carbon, minerals, vitamins and other nutrients which support the growth of microorganisms. Sodium chloride maintains the osmotic balance of the medium. The chromogenic and selective mix allows the identification of bacteria on the basis of the colony color and morphology, while inhibiting most of yeasts and bacteria with the exception of VRE. Agar is the solidifying agent

TEST PROCEDURE

Faecal screening swabs, faecal samples, isolated colony or microbial suspensions can be inoculated by direct streaking on the plate. Incubate aerobically at 37°C for 18-24 hours. Negative plates may be re-incubated for an additional 24 h.

INTERPRETING RESULTS

After incubation observe the color of the colonies and interpret the results as indicated in the ID table.

Identifications can be confirmed using Liofilchem® Strep-Check Kit (ref. 96016), taking a well-isolated colony directly from the Chromatic™ VRE plate.

ID Table.

Microorganism	Typical colony color
Vancomycin resistant <i>E. faecalis</i>	Purple-mauve
Vancomycin resistant <i>E. faecium</i>	Blue-green

See pictures in Appendix I.

APPEARANCE

Slightly opalescent, amber.

STORAGE

2-8°C away from light. Do not use the product beyond its expiry date on the label or if product shows any evidence of contamination or any sign of deterioration.

SHELF LIFE

4 months.

QUALITY CONTROL

Plates are inoculated with the microbial strains indicated in the QC table.

Inoculum for productivity: 10-100 CFU/ml.

Inoculum for selectivity: 10⁴-10⁵ CFU/ml.

Incubation conditions: aerobically at 35±2°C for 18-24 h.

QC Table.

Microorganism		Growth	Specification
<i>Enterococcus faecalis</i> (VRE)	ATCC® 51299	Good	Purple-mauve colonies
<i>Enterococcus faecium</i> (VRE)	DSM 13590	Good	Blue-green colonies
<i>Enterococcus faecalis</i>	ATCC® 29212	Inhibited	---
<i>Enterococcus faecium</i>	ATCC® 19434	Inhibited	---
<i>Escherichia coli</i>	ATCC® 25922	Inhibited	---

WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is intended for *In vitro* diagnostic use and must be used only by properly trained operators.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAPHY

1. Uttley AH, Collins CH, Naidoo J, George RC (1988) Vancomycin-resistant enterococci. *Lancet* 1 (8575-6):57-8.
2. Leclercq R, Derlot E, Duval J, Courvalin P (1988) Plasmid-mediated resistance to vancomycin and teicoplanin in *Enterococcus faecium*. *N. Engl. J. Med.* 319(3): 157-61.
3. Chavers LS, Moser SA, Benjamin WH, Banks SE, Steinhauer JR, Smith AM, Johnson CN, Funkhouser E, Chavers LP, Stamm AM, Waiters KB (2003) Vancomycin-resistant enterococci: 15 years and counting. *J. Hosp. Infect.* 53(3): 159-71.
4. Courvalin P (2006) Vancomycin Resistance in Gram-positive Cocci. *Clin. Infect. Dis.* 42:S25-34.

PRESENTATION

		Contents	Ref.
Chromatic™ VRE	90 mm ready-to-use plates	20 plates	11621

TABLE OF SYMBOLS

LOT Batch code	IVD <i>In vitro</i> Diagnostic Medical Device	 Manufacturer	 Use by	 Fragile, handle with care	 Keep away from sunlight
REF Catalogue number	 Temperature limitation	 Contains sufficient for <n> tests	 Caution, consult Instruction For Use	 Do not reuse	

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