TBX Agar (Tryptone Bile X-Glucuronide)

Catalog #	Description
3564035	TBX Agar (Tryptone Bile X-Glucuronide) , dehydrated, 500 g
3555309	TBX Agar (Tryptone Bile X-Glucuronide) , 100 ml x 6 bottles

For laboratory use only.

Intended Use

Selective chromogenic medium used for the direct enumeration (without colony confirmation) of glucuronidase positive *Escherichia coli* in products intended for human and animal consumption, primary production samples and environmental samples.

Principle

The principle of the medium relies on the detection of the β -D-glucuronidase activity of *E. coli*. The medium contains a complex chromogen linked to a sugar (5-bromo-4-chloro-3-indolyl- β -glucuronide [BCIG]) which is specific to β -D-glucuronidase. Absorbed by the target microorganism, the BCIG is hydrolyzed by this enzyme. The sugar is consumed by the bacteria and the chromogenic agent accumulates in this same cell, giving *E. coli* colonies a blue or blue-green color. With the combined action of high incubation temperature and the presence of bile salts, this medium is selective against Gram positive bacteria and interfering flora.

Theoretical Composition

Base Medium

Enzymatic casein digest	20 g
Bile salts #3	1.5 g
BCIG	0.075 g
Agar	12 g
Distilled water	1,000 ml
Final pH at 25°C = 7.2 $\pm$ 0.2	

Shelf Life and Storage

Store ready-to-use medium 2–8°C in a dark place. Store dehydrated medium at 15–25°C in carefully sealed bottles in a cool, dry place. Medium prepared by user from dehydrated product has a shelf life of 30 days when stored at 2–8°C in a dark place.

Required Materials Not Supplied

This list is not exhaustive.

Equipment

- All usual laboratory equipment
- Incubators or incubation room
- Scales
- Stirrer/homogenizer
- Vortexer
- Water-bath capable of maintaining $\pm$ 1°C
- Sterile forceps for handling membranes

Supplies

- Sterile Petri dishes ($\varnothing$ = 90 mm)
- Sterile disposable platinum loop or inoculating loops
- Sterile pipettes (1 ml, etc)
- Sterile membrane filters ($\varnothing$ = 85 mm, 0.45 μ m to 1.2 μ m)
- Sterile spreaders

Precautions

- Respect Good Laboratory Practice (EN ISO 7218). Appropriate protection, such as gloves and lab coats, should be worn when working with potentially infectious live bacteria
- Media that have come in contact with food, environmental or primary production samples should be considered contaminated and should be disposed of in accordance with local rules and regulations
- Some *E. coli* strains:
 - are β -D-glucuronidase negative (around 3–4 %) and do not produce a blue color (e.g. *E. coli* O157)
 - do not develop at the high temperature of 44°C (e.g. *E. coli* O157:H7)
- Some serovars of *Salmonella* and a few species of *Shigella* also possess the enzyme β -D-glucuronidase (<1.5%)
- If the presence of stressed *E. coli* is suspected when inoculating without a membrane, incubate 4 hr at $37 \pm 1^\circ\text{C}$, then 21 ± 3 hr at $44 \pm 1^\circ\text{C}$
- Incubation should not exceed 24 hr and 45°C
- For media which are ready-to-use or prepared in advance, avoid any prolonged overheating during autoclave sterilization. To maintain optimal quality, the medium must not undergo more than 1 cycle of heating and solidification after initial preparation
- The medium may present a frothy appearance after solidification in bottles. It nevertheless retains all its qualities when the appearance changes after melting and shaking
- As development of colonies at the bottom of the Petri dish may interfere with reading (due to staining of colonies), the time lapse between depositing the inoculum in Petri dishes and addition of the medium must not exceed 15 min
- Avoid trapping any air bubbles underneath the membrane when depositing it on the agar. If necessary, gently and carefully flatten the membrane using sterile forceps
- For SDS product safety information and certificate of analysis, visit bio-rad.com

Quality Control

Every product manufactured and marketed by Bio-Rad is subject to a quality assurance procedure at all stages, from reception of raw materials through to marketing of the finished products. Each batch of finished product undergoes quality control according to EN ISO 11133 and is marketed only if it satisfies the acceptability criteria. Documentation relative to the production and quality control of each batch is kept on file.

Protocol

Dehydrated Base Medium Preparation

- Shake the bottle before use
- Dissolve 33.6 g of powder in 1 L of sterile distilled water. Wait for 5 min, then mix until a homogeneous suspension is obtained
- Heat gently, swirling frequently, and bring to a boil until completely dissolved
- Dispense and sterilize at $121 \pm 3^\circ\text{C}$ for 15 min
- Cool to 44–47°C

Reconstitution ratio: 33.6 g/1 L (500 g of powder makes 14.8 L of medium base)

Sample Preparation and Enrichment

- Prepare and enrich sample according to the standard method applicable to the product concerned

Inoculation and Incubation

Membrane filtration (for products containing stressed bacteria; ISO 16649-1)

- Prepare Petri dishes with minerals modified glutamate medium (MMGA) and with TBX medium (approximately 15 ml/dish)
- Under aseptic conditions, place the membrane on the dry surface of the MMGA agar
- Transfer 1 ml of sample and/or 1 ml of its decimal dilutions to the center of the membrane
- Spread the inoculum uniformly over the entire surface of the membrane, avoiding distributing it beyond the membrane
- Leave the inoculated dishes in a horizontal position at room temperature for about 15 min, until the inoculum has been absorbed by the agar
- Invert the dishes and incubate at $37 \pm 1^\circ\text{C}$ for 4 hr $\pm$ 15 min
- After resuscitation of injured bacteria, using sterile tweezers, transfer the membranes from MMGA medium to Petri dishes containing TBX medium
- Incubate at $44 \pm 1^\circ\text{C}$ for 20–24 hr

Pour plate inoculation (stressed or non-stressed bacteria; ISO 16649-2)

- Using sterile pipettes, transfer 1 ml of sample and/or 1 ml of its decimal dilutions to sterile Petri dishes
- Quickly pour about 15 ml of tempered (44–47°C) TBX medium to the Petri dish. Homogenize by swirling and leave to solidify on a cool, level surface
- Incubate at $44 \pm 1^\circ\text{C}$ for 18–24 hr

Surface inoculation (detection and MNP method; ISO 16649-3)

- Transfer 1 ml of sample and/or its decimal dilutions to 10 ml of MMGA simple and/or double strength MMGA according to the standard protocol
- Incubate at $37 \pm 1^\circ\text{C}$ for 24 ± 2 hr
- Using a sterile loop, streak each positive tube onto TBX dishes
- Incubate at $44 \pm 1^\circ\text{C}$ for 18–24 hr

Reading and Interpretation

- After incubation, count any characteristic glucuronidase positive *E. coli* colonies, which are blue or blue-green. Other colonies are white/beige
- For calculation and expression of results, refer to standards ISO 16649-1, -2 and -3 or ISO 7218

References

ISO 7218:2007/AMD1:2013. Microbiology of food and animal feeding stuffs — General requirements and guidance for microbiological examinations

ISO 16649-1:2018. Microbiology of the food chain — Horizontal method for the enumeration of beta-glucuronidase-positive *Escherichia coli* — Part 1: Colony-count technique at 44 degrees C using membranes and 5-bromo-4-chloro-3-indolyl beta-D-glucuronide.

ISO 16649-2:2001. Microbiology of food and animal feeding stuffs — Horizontal method for the enumeration of beta-glucuronidase-positive *Escherichia coli* — Part 2: Colony-count technique at 44 degrees C using 5-bromo-4-chloro-3-indolyl beta-D-glucuronide.

ISO 16649-3:2015. Microbiology of the food chain — Horizontal method for the enumeration of beta-glucuronidase-positive *Escherichia coli* — Part 3: Detection and most probable number technique C using 5-bromo-4-chloro-3-indolyl beta-D-glucuronide.

Revision History

Release date	Document number	Change
July 2022	5096 Ver A	- Major change - New document design - Document number change — previous version: V4_12/08/11

BIO-RAD is a trademark of Bio-Rad Laboratories, Inc. All trademarks used herein are the property of their respective owner.

TBX Agar (Tryptone Bile X-Glucuronide)

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3564035 **TBX Agar (Tryptone Bile X-Glucuronide)**, base déshydratée, 500 g
3555309 **TBX Agar (Tryptone Bile X-Glucuronide)**, 100 ml x 6 flacons

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Riferimenti

ISO 6579-1:2017/AMD 1:2020. Microbiology of the food chain — Horizontal method for the detection, enumeration, and serotyping of *Salmonella* spp.

ISO 7218:2007/AMD1:2013. Microbiology of food and animal feeding stuffs - General requirements and guidance for microbiological examinations

ISO 19250:2010. Water quality — Detection of *Salmonella* spp.

NF U47-100 July 2007. Animal health analysis methods – Isolation and identification of *Salmonella* or search for particular serovar's in the animal production sector.

NF U47-101 November 2007. Animal health analysis methods – Isolation and identification of any salmonella serotype or of specified salmonella serotypes among birds.

NF U47-102 January 2008. Animal health analysis methods – Isolation and identification of any salmonella serotype or of specified salmonella serotypes among mammals.

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ISO 16649-1:2018. Microbiology of the food chain — Horizontal method for the enumeration of beta-glucuronidase-positive Escherichia coli — Part 1: Colony-count technique at 44 degrees C using membranes and 5-bromo-4-chloro-3-indolyl beta-D-glucuronide.

ISO 16649-2:2001. Microbiology of food and animal feeding stuffs — Horizontal method for the enumeration of beta-glucuronidase-positive Escherichia coli — Part 2: Colony-count technique at 44 degrees C using 5-bromo-4-chloro-3-indolyl beta-D-glucuronide.

ISO 16649-3:2015. Microbiology of the food chain — Horizontal method for the enumeration of beta-glucuronidase-positive Escherichia coli — Part 3: Detection and most probable number technique C using 5-bromo-4-chloro-3-indolyl beta-D-glucuronide.

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