Anti-LC-1 ELISA (IgG)

Test instruction



ORDER NO.	ANTIBODIES AGAINST	IG CLASS	SUBSTRATE	FORMAT
EA 1307-9601 G	LC-1	IgG	Ag-coated microplate wells	96 x 01 (96)

Indications: The ELISA test kit provides a semiquantitative in vitro assay for human autoantibodies of the IgG class against liver cytosol antigen type 1 (LC-1) in serum or plasma for the diagnosis of increase in transaminases for unclear reasons and suspected autoimmune hepatitis.

Application: The determination of autoantibodies against LC-1 is another important supplementing serological parameter in the diagnosis of autoimmune liver diseases. A reliable diagnosis of AIH is indispensable since untreated AIH rapidly turns into liver cirrhosis.

Principle of the test: The test kit contains microplate strips, each with 8 break-off reagent wells coated with LC-1. In the first reaction step, diluted patient samples are incubated in the wells. If samples are positive, specific IgG antibodies (also IgA and IgM) bind to the antigens. To detect the bound antibodies, a second incubation is carried out using an enzyme-labelled anti-human IgG (enzyme conjugate) catalysing a colour reaction.

Contents of the test kit:

Component	Colour	Format	Symbol
1. Microplate wells coated with antigens 12 microplate strips each containing 8 individual break-off wells in a frame, ready for use	---	12 x 8	
2. Calibrator (IgG, human), ready for use	dark red	1 x 2.0 ml	
3. Positive control (IgG, human), ready for use	blue	1 x 2.0 ml	
4. Negative control (IgG, human), ready for use	green	1 x 2.0 ml	
5. Enzyme conjugate peroxidase-labelled anti-human IgG (rabbit), ready for use	green	1 x 12 ml	
6. Sample buffer ready for use	light blue	1 x 100 ml	
7. Wash buffer 10x concentrate	colourless	1 x 100 ml	
8. Chromogen/substrate solution TMB/H ₂ O ₂ , ready for use	colourless	1 x 12 ml	
9. Stop solution 0.5 M sulphuric acid, ready for use	colourless	1 x 12 ml	
10. Test instruction	---	1 booklet	
11. Quality control certificate	---	1 protocol	
 Lot description			 Storage temperature
 In vitro diagnostic medical device			 Unopened usable until

Storage and stability: The test kit has to be stored at a temperature between +2°C to +8°C. Do not freeze. Unopened, all test kit components are stable until the indicated expiry date.

Waste disposal: Patient samples, calibrators, controls and incubated microplate strips should be handled as infectious waste. All reagents must be disposed of in accordance with local disposal regulations.

Preparation and stability of the reagents

Note: All reagents must be brought to room temperature (+18°C to +25°C) approx. 30 minutes before use. After first use, the reagents are stable until the indicated expiry date if stored at +2°C to +8°C and protected from contamination, unless stated otherwise below.

- **Coated wells:** Ready for use. Tear open the resealable protective wrapping of the microplate at the recesses above the grip seam. Do not open until the microplate has reached room temperature to prevent the individual strips from moistening. Immediately replace the remaining wells of a partly used microplate in the protective wrapping and tightly seal with the integrated grip seam (Do not remove the desiccant bag).
Once the protective wrapping has been opened for the first time, the wells coated with antigens can be stored in a dry place and at a temperature between +2°C and +8°C for 4 months.
- **Calibrator and controls:** Ready for use. The reagents must be mixed thoroughly before use.
- **Enzyme conjugate:** Ready for use. The enzyme conjugate must be mixed thoroughly before use.
- **Sample buffer:** Ready for use.
- **Wash buffer:** The wash buffer is a 10x concentrate. If crystallisation occurs in the concentrated buffer, warm it to +37°C and mix well before diluting. The quantity required should be removed from the bottle using a clean pipette and diluted with deionised or distilled water (1 part reagent plus 9 parts distilled water).
For example: For 1 microplate strip, 5 ml concentrate plus 45 ml water.
The working strength wash buffer is stable for 4 weeks when stored at +2°C to +8°C and handled properly.
- **Chromogen/substrate solution:** Ready for use. Close the bottle immediately after use, as the contents are sensitive to light ☀. The chromogen/substrate solution must be clear on use. Do not use the solution if it is blue coloured.
- **Stop solution:** Ready for use.

Warning: The calibrator and controls of human origin have tested negative for HBsAg, anti-HCV, anti-HIV-1 and anti-HIV-2. Nonetheless, all materials should be treated as being a potential infection hazard and should be handled with care. Some of the reagents contain the agent sodium azide in a non-declarable concentration. Avoid skin contact.

Preparation and stability of the patient samples

Samples: Human serum or EDTA, heparin or citrate plasma.

Stability: Patient samples to be investigated can generally be stored at +2°C to +8°C for up to 14 days. Diluted samples should be incubated within one working day.

Sample dilution: Patient samples are diluted **1:101** in sample buffer. For example: dilute 10 µl sample in 1.0 ml sample buffer and mix well by vortexing (sample pipettes are not suitable for mixing).

NOTE: The calibrator and controls are prediluted and ready for use, do not dilute them.

Incubation

Sample incubation: (1st step)

Transfer 100 µl of the calibrator, positive or negative control or diluted patient samples into the individual microplate wells according to the pipetting protocol. Incubate for **30 minutes** at room temperature (+18°C to +25°C).

Washing:

Manual: Empty the wells and subsequently wash 3 times using 300 µl of working strength wash buffer for each wash.

Automatic: Wash the reagent wells 3 times with 450 µl of working strength wash buffer (program setting: e.g. TECAN Columbus Washer "Overflow Mode").

Leave the wash buffer in each well for 30 to 60 seconds per washing cycle, then empty the wells. After washing (manual and automated tests), thoroughly dispose of all liquid from the microplate by tapping it on absorbent paper with the openings facing downwards to remove all residual wash buffer.

Note: Residual liquid (> 10 µl) remaining in the reagent wells after washing can interfere with the substrate and lead to false low extinction values. Insufficient washing (e.g., less than 3 wash cycles, too small wash buffer volumes, or too short residence times) can lead to false high extinction values.

Free positions on the microplate strip should be filled with blank wells of the same plate format as that of the parameter to be investigated.

Conjugate incubation: (2nd step)

Pipette 100 µl of enzyme conjugate (peroxidase-labelled anti-human IgG) into each of the microplate wells. Incubate for **30 minutes** at room temperature (+18°C to +25°C).

Washing:

Empty the wells. Wash as described above.

Substrate incubation: (3rd step)

Pipette 100 µl of chromogen/substrate solution into each of the microplate wells. Incubate for **15 minutes** at room temperature (+18°C to +25°C) protect from direct sunlight.

Stopping:

Pipette 100 µl of stop solution into each of the microplate wells in the same order and at the same speed as the chromogen/substrate solution was introduced.

Measurement:

Photometric measurement of the colour intensity should be made at a wavelength of 450 nm and a reference wavelength between 620 nm and 650 nm **within 30 minutes of adding the stop solution**. Prior to measuring, slightly shake the microplate to ensure a homogeneous distribution of the solution.

Pipetting protocol

	1	2	3	4	5	6	7	8	9	10	11	12
A	C	P 6	P 14	P 22								
B	pos.	P 7	P 15	P 23								
C	neg.	P 8	P 16	P 24								
D	P 1	P 9	P 17									
E	P 2	P 10	P 18									
F	P 3	P 11	P 19									
G	P 4	P 12	P 20									
H	P 5	P 13	P 21									

The above pipetting protocol is an example of the semiquantitative determination of antibodies in 24 patient samples (P 1 to P 24).

Calibrator (C), positive (pos.) and negative (neg.) control as well as the patient samples have been incubated in one well each. The reliability of the ELISA test can be improved by duplicate determinations of each sample.

The wells can be broken off individually from the strips. This makes it possible to adjust the number of test substrates used to the number of samples to be examined and minimises reagent wastage.

Both positive and negative controls serve as internal controls for the reliability of the test procedure. They should be assayed with each test run.

Calculation of results

Semiquantitative: Results can be evaluated semiquantitatively by calculating a ratio of the extinction value of the control or patient sample over the extinction value of calibrator. Use the following formula to calculate the ratio:

$$\frac{\text{Extinction of the control or patient sample}}{\text{Extinction of calibrator}} = \text{Ratio}$$

EUROIMMUN recommends interpreting results as follows:

Ratio <1.0: **negative**
Ratio ≥1.0: **positive**

For duplicate determinations the mean of the two values should be taken. If the two values deviate substantially from one another, EUROIMMUN recommends retesting the samples.

For diagnosis, the clinical picture of the patient always needs to be taken into account along with the serological findings.

Test characteristics

Calibration: As no international reference serum exists for antibodies against LC-1, the results are given as a ratio. This represents a relative measurement for the concentration of antibodies in serum or plasma.

For every group of tests performed, the extinction values of the calibrator and the ratio of the positive and negative control sera must lie within the limits stated for the relevant test kit lot. A quality control certificate containing these reference values is included. If the values specified for the control sera are not achieved, the test results may be inaccurate and the test should be repeated.

The binding activity of the antibodies and the activity of the enzyme used are temperature-dependent. It is therefore recommended using a thermostat in all three incubation steps. The higher the room temperature (+18°C to +25°C) during the incubation steps, the greater the extinction values will be. Corresponding variations apply also to the incubation times. However, the calibrator is subject to the same influences, so that such variations will be largely compensated in the calculation of the result.

Antigen: The microplate wells were coated with recombinant LC-1. The corresponding human cDNA was expressed in insect cells using a baculovirus vector. The specific target antigen of anti-LC-1 antibodies was identified in 1999 as the enzyme formiminotransferase cyclodeaminase (Lapierre et al.).

Detection limit: The lower detection limit is defined as the mean value of an analyte-free sample plus three times the standard deviation and is the smallest detectable antibody titer. The lower detection limit of the Anti-LC-1 ELISA (IgG) is ratio 0.05.

Cross reactivity: This ELISA specifically detects autoantibodies of class IgG against LC-1. When investigating patient sera for autoantibodies against SLA (n = 2) and LKM (n = 8) no cross reactions were found.

Interference: Haemolytic, lipaemic and icteric samples showed no influence on the result up to a concentration of 10 mg/ml for haemoglobin, 20 mg/ml for triglycerides and 0.4 mg/ml for bilirubin in this ELISA.

Reproducibility: The reproducibility of the test was investigated by determining the intra- and inter-assay coefficients of variation (CV) using 3 samples. The intra-assay CVs are based on 20 determinations and the inter-assay CVs on 4 determinations performed in 6 different test runs.

<i>Intra-assay variation, n = 20</i>		
Sample	Mean value (Ratio)	CV (%)
1	3.1	4.9
2	4.7	2.7
3	8.0	2.4

<i>Inter-assay variation, n = 4 x 6</i>		
Sample	Mean value (Ratio)	CV (%)
1	2.5	10.8
2	3.8	10.2
3	6.7	9.4

Prevalence and specificity: Sera of 93 patients with autoimmune hepatitis, 183 patients with hepatitis A, hepatitis B, toxic liver diseases, steatohepatitis or primary biliary cirrhosis (PBC) and 200 healthy blood donors were examined with the EUROIMMUN Anti-LC-1 ELISA (IgG). The prevalence of antibodies against LC-1 in autoimmune hepatitis was 5.4% with a specificity of 100%.

Reference range: The levels of the anti-LC-1 antibodies (IgG) were analysed with this EUROIMMUN ELISA in a panel of 200 healthy blood donors. With a cut-off ratio of 1.0, all blood donors were anti-LC-1 negative.

Literature references

1. Martini E, Abauf N, Cavalli F, Durand V, Johanet C, Homberg JC. **Antibody to liver cytosol (anti-LC-1) in patients with autoimmune chronic active hepatitis type 2.** Hepatology 8: 1662-1666 (1988).
2. Lapierre P, Hajoui O, Homberg JC, Alvarez F. **Formiminotransferase cyclodeaminase is an organ-specific autoantigen recognized by sera of patients with autoimmune hepatitis.** Gastroenterology 116: 643-649 (1999).



REF ORG 548 Anti-MCV

INTENDED PURPOSE

Anti-MCV is an ELISA test system for the quantitative measurement of IgG class autoantibodies against mutated citrullinated vimentin (MCV) in human serum or plasma. This product is intended for professional in vitro diagnostic use only.

Measurement of anti-MCV antibodies contributes to early diagnosis of rheumatoid arthritis (RA), where anti-MCV antibody levels represent one parameter of a multi-criterion diagnostic process, encompassing both clinical and laboratory-based assessments.

SYMBOLS USED ON LABELS

	In vitro diagnostic medical device		Microplate
	Manufacturer		Calibrator
	Catalogue number		Calibrator
	Sufficient for ... determinations		Calibrator
	Batch code		Calibrator
	Use by		Calibrator
	Temperature limitation		Calibrator
	Keep away from sunlight		Control positive
	Do not reuse		Control negative
	Date of manufacture		Sample Buffer P
	CE marked according to 98/79/EC		Enzyme Conjugate
	Consult electronic Instructions For Use		TMB Substrate
	Electronic Instruction For Use: version		Stop solution
	Applicable for U.S.A.: Prescription in vitro diagnostic product		Wash Buffer
	Applicable for U.S.A.: in vitro diagnostic product		Ready to use
			50 x concentrate

PRINCIPLE OF THE TEST

Mutated citrullinated vimentin (MCV) is bound to microwells.

The determination is based on an indirect enzyme linked immune reaction with the following steps:

Specific antibodies in the patient sample bind to the antigen coated on the surface of the reaction wells. After incubation, a washing step removes unbound and unspecifically bound serum or plasma components. Subsequently added enzyme conjugate binds to the immobilized antibody-antigen-complexes. After incubation, a second washing step removes unbound enzyme conjugate. After addition of substrate solution the bound enzyme conjugate hydrolyses the substrate forming a blue coloured product. Addition of an acid stops the reaction generating a yellow end-product. The intensity of the yellow color correlates with the concentration of the antibody-antigen-complex and can be measured photometrically at 450 nm.

WARNINGS AND PRECAUTIONS

- All reagents of this kit are intended for professional in vitro diagnostic use only.
- Components containing human serum were tested and found negative for HBsAg, HCV, HIV1 and HIV2 by FDA approved methods. No test can guarantee the absence of HBsAg, HCV, HIV1 or HIV2, and so all human serum based reagents in this kit must be handled as though capable of transmitting infection.
- Bovine serum albumin (BSA) used in components has been tested for BSE and found negative.
- Avoid contact with the substrate TMB (3,3',5,5'-Tetramethyl-benzidine).
- Stop solution contains acid, classification is non-hazardous. Avoid contact with skin.
- Control, sample buffer and wash buffer contain sodium azide 0.09% as preservative. This concentration is classified as non-hazardous.
- Enzyme conjugate contains ProClin 300 0.05% as preservative. This concentration is classified as non-hazardous.

During handling of all reagents, controls and serum samples observe the existing regulations for laboratory safety regulations and good laboratory practice:

- First aid measures: In case of skin contact, immediately wash thoroughly with water and soap. Remove contaminated clothing and shoes and wash before reuse. If system fluid comes into contact with skin, wash thoroughly with water. After contact with the eyes carefully rinse the opened eye with running water for at least 10 minutes. Get medical attention if necessary.

- Personal precautions, protective equipment and emergency procedures:

Observe laboratory safety regulations. Avoid contact with skin and eyes. Do not swallow. Do not pipette by mouth. Do not eat, drink, smoke or apply makeup in areas where specimens or kit reagents are handled. When spilled, absorb with an inert material and put the spilled material in an appropriate waste disposal.

- Exposure controls / personal protection: Wear protective gloves of nitril rubber or natural latex. Wear protective glasses. Used according to intended use no dangerous reactions known.
- Conditions to avoid: Since substrate solution is light-sensitive. Store in the dark.
- For disposal of laboratory waste the national or regional legislation has to be observed.

Observe the guidelines for performing quality control in medical laboratories by assaying control sera.

CONTENTS OF THE KIT

ORG 548	▽ 96	Sufficient for 96 determinations
MICROPLATE	1	One divisible microplate consisting of 12 modules of 8 wells each. Ready to use. Product code on module: MCV
CALIBRATOR A	1x 1.5 ml	Calibrator A 0 U/ml, containing serum/buffer matrix (PBS, BSA, detergent, NaN3 0.09%), yellow. Ready to use.
CALIBRATOR B	1x 1.5 ml	Calibrator B 20 U/ml, containing MCV antibodies in a serum/buffer matrix (PBS, BSA, detergent, NaN3 0.09%), yellow. Ready to use.
CALIBRATOR C	1x 1.5 ml	Calibrator C 40 U/ml, containing MCV antibodies in a serum/buffer matrix (PBS, BSA, detergent, NaN3 0.09%), yellow. Ready to use.
CALIBRATOR D	1x 1.5 ml	Calibrator D 100 U/ml, containing MCV antibodies in a serum/buffer matrix (PBS, BSA, detergent, NaN3 0.09%), yellow. Ready to use.
CALIBRATOR E	1x 1.5 ml	Calibrator E 300 U/ml, containing MCV antibodies in a serum/buffer matrix (PBS, BSA, NaN3 0.09%), yellow. Ready to use.
CALIBRATOR F	1x 1.5 ml	Calibrator F 1000 U/ml, containing MCV antibodies in a serum/buffer matrix (PBS, BSA, detergent, NaN3 0.09%), yellow. Ready to use.
CONTROL +	1x 1.5 ml	Control positive, containing MCV antibodies in a serum/buffer matrix (PBS, BSA, detergent, NaN3 0.09%), yellow. Ready to use. The concentration is specified on the certificate of analysis.
CONTROL -	1x 1.5 ml	Control negative, containing MCV antibodies in a serum/buffer matrix (PBS, BSA, detergent, NaN3 0.09%), yellow. Ready to use. The concentration is specified on the certificate of analysis.
DILUENT	20 ml	Sample Buffer P, containing PBS, BSA, detergent, preservative sodium azide 0.09%, yellow, concentrate (5 x).
CONJUGATE	15 ml	Enzyme Conjugate containing anti-human IgG antibodies, HRP labelled; PBS, BSA, detergent, preservative PROCLIN 0.05%, light red. Ready to use.
TMB	15 ml	TMB Substrate; containing 3,3', 5,5'- Tetramethylbenzidin, colorless. Ready to use.
STOP	15 ml	Stop solution; contains acid. Ready to use.
WASH	20 ml	Wash Buffer, containing Tris, detergent, preservative sodium azide 0.09%; 50 x conc.

MATERIALS REQUIRED

- Microplate reader capable of endpoint measurements at 450 nm; optional: reference filter at 620 nm
- Data reduction software
- Multi-channel dispenser or repeatable pipette for 100 µl
- Vortex mixer
- Pipettes for 10 µl, 100 µl and 1000 µl
- Laboratory timing device
- Distilled or deionised water
- Measuring cylinder for 1000 ml and 100 ml
- Plastic container for storage of the wash solution

This ELISA assay is suitable for use on open automated ELISA processors. Each assay has to be validated on the respective automated system. Detailed information is provided upon request.

SPECIMEN COLLECTION, STORAGE AND HANDLING

- Collect whole blood specimens using acceptable medical techniques to avoid hemolysis.
- Allow blood to clot and separate the serum or plasma by centrifugation.
- Test serum should be clear and non-hemolyzed. Contamination by hemolysis or lipemia should be avoided, but does not interfere with this assay.
- Specimens may be refrigerated at 2-8°C for up to five days or stored at -20°C up to six months.
- Avoid repetitive freezing and thawing of serum or plasma samples. This may result in variable loss of antibody activity.
- Testing of heat-inactivated sera is not recommended.

STORAGE AND STABILITY

- Store test kit at 2-8°C in the dark.
- Do not expose reagents to heat, sun, or strong light during storage and usage.
- Store microplate sealed and dessicated in the clip bag provided.
- Shelf life of the unopened test kit is 18 months from day of production.
Unopened reagents are stable until expiration of the kit. See labels for individual batch.
- Diluted Wash Buffer and Sample Buffer are stable for at least 30 days when stored at 2-8°C.
We recommend consumption on the same day.

PROCEDURAL NOTES

- Do not use kit components beyond their expiration dates.
- Do not interchange kit components from different lots and products.
- All materials must be at room temperature (20-28°C) prior to use.
- Prepare all reagents and samples. Once started, perform the test without interruption.
- Double determinations may be done. By this means pipetting errors may become obvious.
- Perform the assay steps only in the order indicated.
- Always use fresh sample dilutions.
- Pipette all reagents and samples into the bottom of the wells.
- To avoid carryover or contamination, change the pipette tip between samples and different kit controls.
- Wash microwells thoroughly and remove the last droplets of wash buffer.
- All incubation steps must be accurately timed.
- Do not re-use microplate wells.

PREPARATION OF REAGENTS

WASH

Dilute the contents of one vial of the buffered wash solution concentrate (50x) with distilled or deionised water to a final volume of 1000 ml prior to use.

DILUENT

Sample Buffer P: Prior to use dilute the contents (20 ml) of one vial of sample buffer 5x concentrate with distilled or deionised water to a final volume of 100 ml.

Preparation of samples

Dilute patient samples 1:100 before the assay: Put 990 µl of prediluted sample buffer in a polystyrene tube and add 10 µl of sample. Mix well. Note: Calibrators / Controls are ready to use and need not be diluted.

TEST PROCEDURE

Prepare enough microplate modules for all calibrators / controls and patient samples.

1. Pipette **100 µl** of calibrators, controls and prediluted patient samples into the wells.
Incubate for **30 minutes** at room temperature (20-28 °C).
Discard the contents of the microwells and **wash 3 times** with **300 µl** of wash solution.
2. Dispense **100 µl** of enzyme conjugate into each well.
Incubate for **15 minutes** at room temperature.
Discard the contents of the microwells and **wash 3 times** with **300 µl** of wash solution.
3. Dispense **100 µl** of TMB substrate solution into each well.
Incubate for **15 minutes** at room temperature
4. **Add 100 µl** of stop solution to each well of the modules
Incubate for **5 minutes** at room temperature.
Read the optical density at 450 nm (reference 600-690nm) and calculate the results.
The developed colour is stable for at least 30 minutes. Read during this time.

Example for a pipetting scheme:

	1	2	3	4	5	6	7	8	9	10	11	12
A	A	P1										
B	B	P2										
C	C	P3										
D	D											
E	E											
F	F											
G	C+											
H	C-											

P1, ... patient sample A-F calibrators C+, C- controls

VALIDATION

Test results are valid if the optical densities at 450 nm for calibrators / controls and the results for controls comply with the reference ranges indicated on the Certificate of Analysis enclosed in each test kit.
If these quality control criteria are not met the assay run is invalid and should be repeated.

CALCULATION OF RESULTS

For quantitative results plot the optical density of each calibrator versus the calibrator concentration to create a calibration curve. The concentration of patient samples may then be estimated from the calibration curve by interpolation.
Using data reduction software a 4-Parameter-Fit with lin-log coordinates for optical density and concentration is the data reduction method of choice.

PERFORMANCE CHARACTERISTICS

CALIBRATION

This assay system is calibrated in relative arbitrary units, since no international reference preparation is available for this assay.

Measuring range

The calculation range of this ELISA assay is 0 - 1000 U/ml

Expected values

In a normal range study with samples from healthy blood donors the following ranges have been established with this ELISA assay: Cut-off 20 U/ml

Interpretation of results

Negative: < 20 U/ml
Positive: ≥ 20 U/ml

Linearity

Patient samples containing high levels of specific antibody were serially diluted in sample buffer to demonstrate the dynamic range of the assay and the upper / lower end of linearity. Activity for each dilution was calculated from the calibration curve using a 4-Parameter-Fit with lin-log coordinates.

Sample	Dilution	Observed U/ml	Expected U/ml	O/E [%]
1	1:100	882.8	882.8	100
.	1:200	386.0	441.4	87
.	1:400	205.2	220.7	93
.	1:800	110.7	110.4	100
.	1:1600	52.2	55.2	95
.	1:3200	23.4	27.6	85
2	1:100	932.1	932.1	100
.	1:200	486.0	466.1	104
.	1:400	250.1	233.0	107
.	1:800	126.6	116.5	109
.	1:1600	61.7	58.3	106
.	1:3200	28.2	29.1	97
3	1:100	727.9	727.9	100
.	1:200	362.4	364.0	100
.	1:400	178.2	182.0	98
.	1:800	85.7	91.0	94
.	1:1600	47.1	45.5	104
.	1:3200	19.2	22.7	85

Limit of detection

Functional sensitivity was determined to be: 1 U/ml

Reproducibility

Intra-assay precision: Coefficient of variation (CV) was calculated for each of three samples from the results of 24 determinations in a single run. Results for precision-within-assay are shown in the table below.
Inter-assay precision: Coefficient of variation (CV) was calculated for each of three samples from the results of 6 determinations in 5 different runs. Results for run-to-run precision are shown in the table below.

Intra-Assay		
Sample	Mean U/ml	CV %
1	22.7	6.2
2	118.8	6.4
3	548.1	4.6

Inter-Assay		
Sample	Mean U/ml	CV %
1	20.2	5.3
2	111.0	9.2
3	451.6	7.7

Interfering substances

No interference has been observed with haemolytic (up to 1000 mg/dl) or lipemic (up to 3 g/dl triglycerides) sera or plasma, or bilirubin (up to 40 mg/dl) containing sera or plasma. Nor have any interfering effects been observed with the use of anticoagulants (Citrate, EDTA, Heparine). However for practical reasons it is recommended that grossly hemolyzed or lipemic samples should be avoided.

Study results

<u>Study population</u>	<u>n</u>	<u>n Pos</u>	<u>%</u>
Rheumatoid arthritis	490	398	81.2
Other diseases	522	14	2.7
Normal human sera	234	1	0.4

		Clinical Diagnosis		
		POS	NEG	
ORG 548	POS	398	15	
	NEG	92	741	
		490	756	1246
Sensitivity: 81.2 %				
Specificity: 98.0 %				
Overall agreement: 91.4 %				

LIMITATIONS OF THE PROCEDURE

This assay is a diagnostic aid. A definite clinical diagnosis should not be based on the results of a single test, but should be made by the physician after all clinical and laboratory findings have been evaluated concerning the entire clinical picture of the patient. Also every decision for therapy should be taken individually.

The above pathological and normal reference ranges for antibodies in patient samples should be regarded as recommendations only. Each laboratory should establish its own ranges according to ISO 15189 or other applicable laboratory guidelines.

REFERENCES

1. F.Bobbio-Pallavicini, C.Alpini, R.Caporali, S.Avalle, S.Bugatti, C.Montecuccio. Autoantibody profile in rheumatoid arthritis during long-term infliximab treatment. Arthritis Res Ther 2004, 6:R264-R272 (DOI 10.1186/ar1173)
2. E.R.Vossenaar, N.Deprés, E.Lapointe,A.van der Heijden, M.Lora,T.Senshu,W.J.van Venrooij, H.A.Ménard. Rheumatoid arthritis specific anti Sa antibodies target citrullinated vimentin. Arthritis Research & Therapie Vol.6 No.2
3. O.Vittecoq, S.Poplin, K.Krzanowska, F.Jouen-Beades, J.F.Ménard, A.Daragon, F.Tron, X.Loet. Rheumatoid Factor is the strongest predictor of radiological progression of rheumatoid arthritis in a three-year prospective study in community-recruited patients. Rheumatology 2003; 42:939-946
4. W.J. van Venrooij, J.M.Hazes, H.Visser. Anticitrullinated protein/peptide antibody and its role in the diagnosis and prognosis of early rheumatoid arthritis. The Netherlands Journal of Medicine.
5. M.Escalon, F.J.Lópes-Longo, C.M. González, I.Monteagudo, M.Rodríguez-Mahou, R.Grau, L.Carreno. Anti-Sa Sera from patients with Rheumatoid Arthritis contain at least 2 different subpopulations of Anti-Sa antibodies. The Journal of Rheumatology 2002; 29:10 2053-60
6. Ch.Vincent, L.Nogueira, M.Sebba, S.Chapuy-Regaud, M.Arnaud, O.Letourneur, D.Rolland, B.Rounie, A. Cantagrel, M.Jolivet, G.Serre. Detection of antibodies to determined recombinant tat flaggrin by Enzyme-Linked Immunosorbent Assay. Arthritis & Rheumatism Vo. 46, No.8, August 2002, pp. 2051 -58
7. G.Steiner, J.Smolén. Antibodies in rheumatoid arthritis and their clinical significance. Arthritis Res 2002,4 (suppl 2):S1-S5
8. R.Goldbach-Mansky, J.Lee, A.McCoy, J.Hoxworth, C.Yarboro, J.S.Smolén, G.Steiner, A.Rosen, C.Zhang, H.A. Ménard, Z.J.Zhou, T.Palosuo, W.J.Van Venrooij, R.L.Wilder, J.H.Klippel, H.R.Schumacher Jr., H.S.El-Gabalawy. Rheumatoid arthritis associated antibodies in patients with synovitis of recent onset. Arthritis Res 2000,2:236–243
9. H.Ménard, E.Lapointe, M.D.Rochdi, Z.J.Zhou. Insights into rheumatoid arthritis derived from the Sa immune system. Arthritis Research 2000,2:429-432
10. G.Hayem, P.Chazerain, B.Combe, A.Elias, T.Haim, P.Nicaise, K.Benali, J-F Eliau, M-F Kahn, J.Sany, O.Meyer. Anti-Sa antibody is an accurate diagnostic and prognostic Marker in adult Rheumatoid Arthritis.The Journal of Rheumatology 1999;26:7-13
11. N. Deprés, G.Boire, F.J. Lopez-Longo, H.A. Ménard. The Sa System : A novel antigen antibody system specific for rheumatoid arthritis. The Journal of Rheumatology 1994; 21 :1 027- 33
12. E.R.Vossenaar,T.R.D. Radstake,A. van der Heijden, M.A.M. van Mansum, C. Dieteren, D.-J. de Rooij, P. Barrera, A.J.W. Zendman,W.J. van Venrooij. Expression and activity of citrullinated peptidylarginine deiminase enzymes in monocytes and macrophages. Ann Rheum Dis 2004; 63:373-381
13. H. Burkhardt, B. Sehnert, R. Bockermann, A. Engström, J.R. Kalden, R. Holmdahl. Humoral immune response to citrullinated collagen type II determinants in early rheumatoid arthritis. Eur. J. Immunol. 2005.35: 1643-1652

14. Ch.Vincent, L. Nogueira, C. Clavel, M. Sebbag, G. Serre. Autoantibodies to citrullinated pro-teins: ACPA. Autoimmunity, February 2005; 38 (1): 17-24
15. M. Sebbag, S. Chapuy-Regaud, I. Auger, E. Petit-Teixeira, C. Clavel, L. Nogueira, Ch.Vincent, F. Cornélis,J. Roudier, G. Serre.Clinical and pathophysiological significance of the autoimmune response to citrullinated proteins in rheumatoid arthritis. Joint bone Spine 71 (2004) 493-502

References Concerning Classification of RA

16. Liao KP, Batra KL, Chibnik L, Schur PH, Costenbader KH. Anti-CCP revised criteria for the classification of rheumatoid arthritis. Ann Rheum Dis 2008.
17. Arnett FC, Edworthy SM, Bloch DA, McShane DJ, Fries JF, Cooper NS et al. The American Rheumatism Association 1987 revised criteria for the classification of rheumatoid arthritis. Arthritis Rheum 1988; 31(3):315 -24.Further Reading on Diagnostic Performance of Anti-MCV vs. Predicate Assay:
18. Gyetvai A, Szekanecz Z, Soos L, Szabo Z, Fekete A, Kapitány A et al. New classification of the shared epitope in rheumatoid arthritis: impact on the production of various anti-citrullinated protein antibodies. Rheumatology (Oxford) 2010; 49(1):25-33.
19. Puijn G, Wiik A, van Venrooij W. The use of citrullinated peptides and proteins for the diagnosis of rheumatoid arthritis. Arthritis Research & Therapy 2010; 12(1):203.
20. Aggarwal R, Liao K, Nair R, Ringold S, Costenbader KH. Anti-citrullinated peptide antibody assays and their role in the diagnosis of rheumatoid arthritis. Arthritis Rheum 2009; 61(11):1472-83.
21. Feitsma AL, van der Voort EI, Franken KL, El BH, Elferink BG, Drijfhout JW et al. Identification of citrullinated vimentin peptides as T cell epitopes in HLA-DR4-positive patients with rheumatoid arthritis. Arthritis Rheum 2009; 62(1):117-25.
22. Klareskog L, Catrina AI, Paget S. Rheumatoid arthritis. Lancet 2009; 373(9664):659-72.
23. Liu X, Jia R, Zhao J, Li Z. The Role of Anti-Mutated Citrullinated Vimentin Antibodies in the Diagnosis of Early Rheumatoid Arthritis. J Rheumatol 2009.
24. Luime JJ, Colm EM, Hazes JM, Lubberts E. Does anti-MCV has additional value as serological marker in the diagnostic and prognostic work-up of patients with rheumatoid arthritis? A systematic review. Ann Rheum Dis 2009.
25. Mjaavatten MD, Uhlig T, Haugen AJ, Nygaard H, Sidenval G, Helgetveit K et al. Positive anti-citrullinated protein antibody status and small joint arthritis are consistent predictors of chronic disease in patients with very early arthritis: results from the NOR-VEAC cohort. Arthritis Res Ther 2009; 11(5):R146.
26. Engelmann R, Brandt J, Eggert M, Karberg K, Krause A, Neeck G et al. The Anti-mutated Citrullinated Vimentin Response Classifies Patients with Rheumatoid Arthritis into Broad and Narrow Responders. J Rheumatol 2009.
27. Raza K, Mathsson L, Buckley CD, Filer A, Ronnelid J. Anti-modified citrullinated vimentin (MCV) antibodies in patients with very early synovitis. Ann Rheum Dis 2009.
28. Snir O, Widhe M, Hermansson M, von SC, Lindberg J, Hensen S et al. Antibodies to several citrullinated antigens are enriched in the joints of rheumatoid arthritis patients. Arthritis Rheum 2009; 62(1):44-52.
29. Syversen SW, Goll GL, van der Heijde D, Landewe R, Lie BA, Odegard S et al. Prediction of radiographic progression in rheumatoid arthritis and the role of antibodies against mutated citrullinated vimentin: results from a ten-year prospective study. Ann Rheum Dis 2009.
30. van der Linden MP, van der Woude D, Ioan-Facsinay A, Levarht EW, Stoeken-Rijsbergen G, Huizinga TW et al. Value of anti-modified citrullinated vimentin and third-generation anti-cyclic citrullinated peptide compared with second-generation anti-cyclic citrullinated peptide and rheumatoid factor in predicting disease outcome in undifferentiated arthritis and rheumatoid arthritis. Arthritis Rheum 2009; 60(8):2232-41.
31. Wagner E, Skoumal M, Bayer PM, Klaushofer K. Antibody against mutated citrullinated vimentin: a new sensitive marker in the diagnosis of rheumatoid arthritis. Rheumatol Int 2009.
32. Keskin G, Inal A, Keskin D, Pekel A, Baysal O, Dizer U et al. Diagnostic Utility of Anti-Cyclic Citrullinated Peptide and Anti-Modified Citrullinated Vimentin Antibodies in Rheumatoid Arthritis. Protein Pept Lett 2008; 15 (3):314-7.
33. Szekanecz Z, Lakos G. Rheumatoid arthritis diagnosis with antimutated citrullinated vimentin ELISA by Orgentec Diagnostika. Exp Opin Med Diagnostics 2008; 2(9):1083-90.
34. Szekanecz Z, Soos L, Szabo Z, Fekete A, Kapitány A, Vegvari A et al. Anti-Citrullinated Protein Antibodies in Rheumatoid Arthritis: As Good as it Gets? Clin Rev Allergy Immunol 2008; 34(1):26-31.

Notice to the user (European Union):

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the EU Member State in which the user and/or the patient is established .

Change Control

Former version: *ORG 548_IFU_EN_QM113183_2023-08-02_5*

Reason for revision: Definition of symbols used and symbols updated

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Anti-SLA/LP ELISA (IgG)


Test instruction

ORDER NO.	ANTIBODIES AGAINST	IG-CLASS	SUBSTRATE	FORMAT
EA 1302-9601 G	SLA/LP (soluble liver antigen/ liver-pancreas antigen)	IgG	Ag-coated microplate wells	96 x 01 (96)

Indications: Increase in transaminases for unclear reasons, suspected autoimmune hepatitis.

Principles of the test: The ELISA test kit provides a semiquantitative or quantitative in vitro assay for human autoantibodies of the IgG class against SLA/LP (soluble liver antigen/liver-pancreas antigen) in serum or plasma. The test kit contains microtiter strips each with 8 break-off reagent wells coated with SLA/LP. In the first reaction step, diluted patient samples are incubated in the wells. In the case of positive samples, specific IgG antibodies (also IgA and IgM) will bind to the antigens. To detect the bound antibodies, a second incubation is carried out using an enzyme-labelled anti-human IgG (enzyme conjugate) which is capable of promoting a colour reaction.

Contents of the test kit:

Component	Colour	Format	Symbol
1. Microplate wells, coated with antigens: 12 microplate strips each containing 8 individual break-off wells in a frame, ready for use	---	12 x 8	
2. Calibrator 1 200 RU/ml (IgG, human), ready for use	dark red	1 x 2.0 ml	
3. Calibrator 2 20 RU/ml (IgG, human), ready for use	red	1 x 2.0 ml	
4. Calibrator 3 2 RU/ml (IgG, human), ready for use	light red	1 x 2.0 ml	
5. Positive control (IgG, human), ready for use	blue	1 x 2.0 ml	
6. Negative control (IgG, human), ready for use	green	1 x 2.0 ml	
7. Enzyme conjugate peroxidase-labelled anti-human IgG (rabbit), ready for use	green	1 x 12 ml	
8. Sample buffer ready for use	light blue	1 x 100 ml	
9. Wash buffer 10x concentrate	colourless	1 x 100 ml	
10. Chromogen/substrate solution TMB/H ₂ O ₂ , ready for use	colourless	1 x 12 ml	
11. Stop solution 0.5 M sulphuric acid, ready for use	colourless	1 x 12 ml	
12. Test instruction	---	1 booklet	
13. Protocol with target values	---	1 protocol	
 Lot			 storage temperature
 In vitro determination			 unopened usable until

Storage and stability: The test kit has to be stored at a temperature between +2°C to +8°C. Do not freeze. Unopened, all test kit components are stable until the indicated expiry date.

Waste disposal: Patient samples, calibrators, controls and incubated microplate strips should be handled as infectious waste. All reagents should be disposed of according to official regulations.

Preparation and stability of the reagents

Note: All reagents must be brought to room temperature (+18°C to +25°C) around 30 minutes before use. After first use, the reagents are stable until the indicated expiry date if stored at +2°C to +8°C and protected from contamination, unless stated otherwise below.

- **Coated wells:** Ready for use. Tear open the resealable protective wrapping of the microplate at the recesses above the grip seam. Do not open until the microplate has reached room temperature to prevent the individual strips from moistening. Immediately replace the remaining wells of a partly used microplate in the protective wrapping and tightly seal with the integrated grip seam (Do not remove the desiccant bag).
Once the protective wrapping has been opened for the first time, the wells coated with antigens can be stored in a dry place and at a temperature between +2°C and +8°C for a minimum of 4 months.
- **Calibrators and controls:** Ready for use. The reagents must be mixed thoroughly before use.
- **Enzyme conjugate:** Ready for use. The enzyme conjugate must be mixed thoroughly before use.
- **Sample buffer:** Ready for use.
- **Wash buffer:** The wash buffer is a 10x concentrate. If crystallization occurs in the concentrated buffer, warm it to 37°C and mix well before diluting. The quantity required should be removed from the bottle using a clean pipette and diluted with deionized or distilled water (1 part reagent plus 9 parts distilled water).
For example, for 1 microplate: 5 ml concentrate plus 45 ml water.
The ready-to-use diluted wash buffer is stable for 4 weeks when stored at +2°C to +8°C and handled properly.
- **Chromogen/substrate solution:** Ready for use. Close the bottle immediately after use, as the contents are sensitive to light. The chromogen/substrate solution must be clear on use. Do not use the solution if it is blue coloured.
- **Stop solution:** Ready for use.

Warning: Calibrators and controls used have tested negative for HBsAg, anti-HCV, anti-HIV-1 and anti-HIV-2 using enzyme immunoassays and indirect immunofluorescence methods. Nonetheless, all materials should be treated as being a potential infection hazard and should be handled with care. Some of the reagents contain the toxic agent sodium azide. Avoid contact with the skin.

Preparation and stability of the patient samples

Sample material: Human serum or EDTA, heparin or citrate plasma.

Stability: Patient samples to be investigated can generally be stored at +2°C to +8°C for up to 14 days. Diluted samples should be incubated within one working day.

Sample dilution: Patient samples are diluted **1:101** sample buffer. For example: dilute 10 µl serum in 1.0 ml sample buffer and mix well by vortexing (sample pipettes are not suitable for mixing).

NOTE: Calibrators and controls are prediluted and ready for use, do not dilute them.

Incubation

For **semiquantitative analysis** incubate **calibrator 2** along with the positive and negative controls and patient samples. For **quantitative analysis** incubate **calibrators 1, 2 and 3** along with the positive and negative controls and patient samples.

Sample incubation:
(1. step)

Transfer 100 µl calibrators, positive and negative controls or diluted patient samples into the individual microplate wells according to the pipetting protocol. Incubate for **30 minutes** at room temperature (+18°C to +25°C).

Washing:

Manual: Empty the wells and subsequently wash 3 times using 300 µl of working strength wash buffer for each wash.

Automatic: Wash reagent wells 3 times with 400 µl working strength wash buffer (program setting: e.g. TECAN Columbus Washer "Overflow Modus").

Leave the wash buffer in each well for 30 to 60 seconds per washing cycle, then empty the wells. After washing (manual and automated tests), thoroughly dispose of all liquid from the microplate by tapping it on absorbent paper with the openings facing downwards to remove all residual wash buffer.

Note: Residual liquid (> 10 µl) remaining in the reagent wells after washing can interfere with the substrate and lead to falsely low extinction values. Insufficient washing (e.g., less than 3 wash cycles, too small wash buffer volumes, or too short reaction times) can lead to falsely high extinction values.

Free positions on the microplate strip should be filled with blank wells of the same plate format as that of the parameter to be investigated.

Conjugate incubation:
(2. step)

Pipette 100 µl of enzyme conjugate (peroxidase-labelled anti-human IgG) into each of the microplate wells. Incubate for **30 minutes** at room temperature (+18°C bis 25°C).

Washing:

Empty the wells. Wash as described above.

Substrate incubation:
(3. step)

Pipette 100 µl of chromogen/substrate solution into each of the microplate wells. Incubate for **15 minutes** at room temperature (+18°C bis 25°C) protect from direct sunlight.

Stopping the reaction:

Pipette 100 µl of stop solution into each of the microplate wells in the same order and at the same speed as the chromogen/substrate solution was introduced.

Measurement:

Photometric measurement of the colour intensity should be made at a wavelength of 450 nm and a reference wavelength of between 620 nm and 650 nm **within 30 minutes of adding the stop solution**. Prior to measuring, slightly shake the micro-plate to ensure a homogeneous distribution of the solution.

Pipetting protocol

	1	2	3	4	5	6	7	8	9	10	11	12
A	C 2	P 6	P 14	P 22			C 1	P 4	P 12	P 20		
B	pos.	P 7	P 15	P 23			C 2	P 5	P 13	P 21		
C	neg.	P 8	P 16	P 24			C 3	P 6	P 14	P 22		
D	P 1	P 9	P 17				pos.	P 7	P 15	P 23		
E	P 2	P 10	P 18				neg.	P 8	P 16	P 24		
F	P 3	P 11	P 19				P 1	P 9	P 17			
G	P 4	P 12	P 20				P 2	P 10	P 18			
H	P 5	P 13	P 21				P 3	P 11	P 19			

The pipetting protocol for microtiter strips 1-4 is an example for the **semiquantitative analysis** of 24 patient samples (P 1 to P 24).

The pipetting protocol for microtiter strips 7-10 is an example for the **quantitative analysis** of 24 patient samples (P 1 to P 24).

The calibrators (C 1 to C 3), the positive (pos.) and negative (neg.) controls, and the patient samples have each been incubated in one well. The reliability of the ELISA test can be improved by duplicate determinations for each sample.

The wells can be broken off individually from the strips. This makes it possible to adjust the number of test substrates used to the number of samples to be examined and minimizes reagent wastage.

Both positive and negative controls serve as internal controls for the reliability of the test procedure. They should be assayed with each test run.

Calculation of results


Semiquantitative: Results can be evaluated semiquantitatively by calculating a ratio of the extinction value of the control or patient sample over the extinction value of calibrator 2. Calculate the ratio according the following formula:

$$\frac{\text{Extinction of the control or patient sample}}{\text{Extinction of calibrator 2}} = \text{Ratio}$$

EUROIMMUN recommends interpreting results as follows:

Ratio <1.0:	negative
Ratio ≥1.0:	positive

Quantitative: The standard curve from which the concentration of antibodies in the serum samples can be taken is obtained by point-to-point plotting of the extinction values measured for the 3 calibration sera against the corresponding units (linear/linear). Use "point-to-point" plotting for calculation of the standard curve by computer. The following plot is an example of a typical calibration curve. Please do not use this curve for the determination of antibody concentrations in patient samples.

If the extinction of a serum sample lies above the value of calibrator 1 (200 RU/ml). The result should be given as ">200 RU/ml". It is recommended that the sample be re-tested at a dilution of 1:400. The result in IU/ml read from the calibration curve for this sample must then be multiplied by a factor of 4.

The upper limit of the normal range (**cut-off**) recommended by EUROIMMUN is 20 relative units (RU) /ml. EUROIMMUN recommends interpreting results as follows:

<20 RU/ml:	negative
≥20 RU/ml:	positive

For duplicate determinations the mean of the two values should be taken. If the two values deviate substantially from one another the sample should be retested.

For diagnosis, the clinical symptoms of the patient should always be taken into account alongside the serological results.

Test characteristics

Calibration: As no international reference serum exists for antibodies against SLA/LP, the calibration is performed in relative units (RU).


For every group of tests performed, the extinction values of the calibrators and the relative units and/or ratios determined for the positive and negative controls must lie within the limits stated for the relevant test kit lot. A protocol containing these target values is included. If the values specified for the controls are not achieved, the test results may be inaccurate and the test should be repeated.

The activity of the enzyme used is temperature-dependent and the extinction values may vary if a thermostat is not used. The higher the room temperature (+18°C to +25°C) during substrate incubation, the greater will be the extinction values. Corresponding variations apply also to the incubation times. However, the calibrators are subject to the same influences, with the result that such variations will be largely compensated in the calculation of the result.

Antigen: The microplate wells were coated with recombinant SLA/LP. The corresponding human cDNA was expressed in *E. coli*. Identification of SLA/LP at the DNA level succeeded in 1998 by cloning the target antigen (A. W. Lohse, University Hospital of Mainz). SLA/LP is probably a cytoplasmic molecule with a molecular weight of 50 kDa which is involved in the regulation of the protein biosynthesis (an UGA-suppressor-tRNA associated protein). Previous descriptions of SLA as being the liver cytochromes 8 and 18 or glutathione-S-transferase were apparently wrong.

Linearity: The linearity of the test was investigated by assaying serial dilutions of patient samples with high antibody concentrations. The chart below shows the typical linearity of samples on the basis of 3 patient sera. The Anti-SLA/LP ELISA (IgG) is linear in the measurement range 2 - 200 RU/ml.

Detection limit: The detection limit is defined as a value of three times the standard deviation of an analyte-free sample and is the smallest detectable antibody titer. The detection limit of the Anti-SLA/LP ELISA (IgG) is approximately 1 RU/ml.

Cross reactivity: This ELISA showed no cross reactivity.

Interference: Haemolytic, lipaemic and icteric samples showed no influence at the result up to a concentration of 10 mg/ml for hemoglobin, 20 mg/ml for triglycerides and 0.4 mg/ml for bilirubin in this ELISA.

Reproducibility: The reproducibility of the test was investigated by determining the intra- and inter-assay coefficients of variation (CV) using 3 sera with values at different points on the calibration curve. The intra-assay CVs are based on 20 determinations and the inter-assay CVs on 4 determinations performed on 6 different days.

<i>Intra-Assay Variation, n = 20</i>		
Serum	Mean value (RU/ml)	CV (%)
1	41	2.6
2	110	2.3
3	154	3.7

<i>Inter-Assay Variation, n = 4 x 6</i>		
Serum	Mean value (RU/ml)	CV (%)
1	44	4.4
2	118	3.6
3	169	3.8

Prevalence and specificity: Sera from 454 patients with autoimmune hepatitis, 165 patients with other liver diseases and 200 healthy blood donors were examined with the Anti-SLA/LP ELISA. The prevalence anti-SLA/LP antibodies in non-Japanese AIH patients was in the range between 15% and 19%. The test showed a specificity of 100%.

Patient group (n=819)	Origin of sera	n	Anti-SLA/LP-positive
Autoimmune hepatitis	Prof. Lohse, Univ. Mainz	108	21 (19%)
	Brasilian panel	154	25 (16%)
	Japanese panel	43	2 (5%)
	American panel	149	23 (15%)
AIH (anti-LKM-1 positive)	Dr. Gruber, Univ. München	18	0
Primary biliary cirrhosis	Dr. Gruber, Univ. München PD Dr. Wick, Klinikum Großhadern	30	0
Hepatitis-B infection	Labor Dr. Stöcker, Lübeck	40	0
Hepatitis-C infection	Prof. Lohse, Univ. Mainz Labor Dr. Stöcker, Lübeck	39	0
Steatohepatitis	Prof. Lohse, Univ. Mainz	24	0
Toxic liver damage	Prof. Lohse, Univ. Mainz	14	0
Blood donors	Med. Univ. Lübeck	200	0

Reference range: The levels of the anti-SLA/LP antibodies (IgG) were analyzed with this EUROIMMUN ELISA in a panel of 200 healthy blood donors. With a cut-off of 20 RU/ml, all blood donors were anti-SLA/LP negative.

Clinical significance

The determination of autoantibodies against soluble liver antigen/liver-pancreas antigen (SLA/LP) is a new and important component in the diagnostics of autoimmune diseases of the liver.

Autoimmune diseases of the liver include

- autoimmune hepatitis (AIH),
- primary biliary liver cirrhosis (PBC) and
- primary sclerosing cholangitis (PSC).

Mainly women are affected by **autoimmune hepatitis (AIH)**, earlier designations: lupoid hepatitis, chronically active hepatitis). The disease manifests itself through increased levels of bilirubin, liver enzymes and immunoglobulins, through typical histological changes (liver biopsies show necrosis of the parenchymal cells with lymphocyte and plasma cell infiltrates) and through the occurrence of various autoantibodies. The disease can occur from infancy until the old age, however, it affects most frequently young adults.

The incidence of AIH in western Europe is 1.9 cases per 100,000 inhabitants per year. Untreated, autoimmune hepatitis soon develops into liver cirrhosis. However, if immunosuppressive low-dosed therapy is commenced in good time and continued consistently throughout the patient's life, the patient has a normal life expectancy. For differential diagnosis, a current infection with hepatitis viruses must be ruled out by investigation of the appropriate serological parameters.

Circulating autoantibodies have gained great significance for the diagnosis of AIH. They occur in the majority of patients, but their role in the pathogenesis is questionable. Also, there is no clear correlation between the activity or the prognosis of the disease and the antibody titre. Besides antibodies against SLA/LP the following autoantibodies are associated with AIH: antibodies against cell nuclei (**ANA**), **nDNA**, smooth muscle (**SMA**, with the most important target antigen being **F-actin**), liver-kidney micro-somes (**LKM-1**, target antigen: cytochrome P450 IID6) and granulocytes (**pANCA**). The **autoantibodies against SLA/LP** that can today be measured by EUROIMMUN enzyme immunoassays probably have the highest diagnostic accuracy of all antibodies involved in autoimmune hepatitis. Anti-SLA/LP occur in AIH either singly or together with other autoantibodies. Although their prevalence is only between 10% and 30%, the predictive value is almost 100%: Every positive result essentially provides evidence of autoimmune hepatitis (provided that the corresponding clinical symptoms are present).

Investigations of antibodies against SLA/LP have up until now only been possible in a few special laboratories. It was earlier suspected that the SLA antigen is identical with the liver cytokeratins 8 and 18 or with the enzyme glutathione-S-transferase. It has only recently been found that this is apparently not the case. The newly discovered target antigen SLA/LP has presently only been identified reliably at DNA level. It is probably a cytoplasmic molecule which is involved in the regulation of the protein biosynthesis (an UGA-suppressor-tRNA associated protein). The great diagnostic value of autoantibodies against SLA/LP has now been clearly proved.

Autoantibodies against cell nuclei (ANA) and against smooth muscle (SMA) are frequent in AIH, but also occur in 10% to 20% of patients with chronic viral hepatitis and in other diseases. Autoantibodies against LKM-1 can only be demonstrated in about 1% of adult AIH patients, but they are more frequent in children. Antibodies against LKM-1 can also be found 1% to 2% of patients with a positive hepatitis C serology. As opposed to all other autoantibodies, antibodies against SLA/LP are highly specific for AIH and have not been described in viral hepatitis.

Some authors classify autoimmune hepatitis in accordance with the autoantibody status: Subtype I (ANA, SMA), Subtype II (antibodies against LKM-1) and Subtype III (antibodies against SLA/LP). This classification has probably neither a clinical nor a therapeutic and prognostic significance.

The diagnosis of AIH is based on the clinical picture, the biochemical and serological test results and the histological evidence of an inflammation reaction. It requires the exclusion of other causes of chronic hepatitis, such as viruses, alcohol and drugs, and a delimitation to the other autoimmune liver diseases (primary biliary liver cirrhosis and primary sclerosing cholangitis). The diagnosis AIH is appropriate when 4 of the following 5 main criteria are fulfilled. If all 5 criteria are present, the diagnosis is definite. The diagnosis is confirmed if the patient responds to immunosuppressive therapy.

Main criteria for the diagnosis of autoimmune hepatitis:

1. Histological evidence of hepatitis
2. Detection of autoantibodies (ANA, SMA, LKM, SLA/LP)
3. Hypergammaglobulinaemia
4. Negative HBV and HCV serology
5. Detection of the HLA antigens B8, DR3 or DR4

The serological determination of **autoantibodies against SLA/LP** provides an exact delimitation to viral hepatitis in many patients with AIH. The autoimmune serology of hepatitis has thus been enriched by a further parameter whose significance must be rated more highly than that of antibodies against cell nuclei, DNA, smooth muscle and liver-kidney microsomes (**ANA, nDNA, ASMA, LKM**). Testing for these antibodies has considerable consequences for the hepatology clinic: the consequence of incorrect treatment of AIH with interferon can be just as fatal as the immunosuppressive therapy of the virus infection.

Primary biliary cirrhosis (PBC) is a liver cirrhosis which takes the form of a non-suppurative, destructive inflammation of the biliferous ducts. It occurs in woman ten times more frequently than in men. In the foreground of the clinical picture is cholestasis. In Europe, 13 out of every 100,000 persons develop this disease each year. As a cirrhosis does not occur in every case, and mostly not until the late stage, this disease is better described as “**chronic non-suppurative destructive cholangiitis**”. Pathognomonic is the serological detection of autoantibodies against the mitochondrial fraction M2 (**AMA M2**) and against **Nuclear Dots** (acid protein SP100). Approximately 10% to 20% of PBC patients develop a secondary autoimmune hepatitis (also called **overlap syndrome**). In these cases, as in AIH, autoantibodies can frequently be detected. Antibodies against SLA/LP indicate a secondary AIH (overlap syndrome) which requires an immunosuppressive therapy.

The incidence of **primary sclerosing cholangitis (PSC)** is quoted as being 4 cases in every 100,000 inhabitants per year. The clinical picture is also characterised by cholestasis. Diagnostic pointers are the relevant laboratory findings, the histology and ERCP (endoscopic retrograde cholangiopancreatography). Men are mainly affected, and in half of the patients a concomitant Colitis ulcerosa is present (inversely, the prevalence of PSC in cases of Colitis ulcerosa is 4%). Most patients with PSC display autoantibodies against granulocytes (**pANCA**). Occasionally, pANCA can also occur in AIH and PBC, and their value for the differential diagnosis is limited in this respect - however, pANCA are an indication

of an autoimmune disease of the liver and can possibly be used as a means of delimitation to infectious forms of hepatitis. Even though overlaps between PSC and AIH have been reported, antibodies against SLA/LP have not been found in such patients.

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 ORG 510 Anti-Sm

INTENDED PURPOSE

Anti-Sm is an ELISA test system for the quantitative measurement of IgG class autoantibodies against Sm in human serum or plasma. This product is intended for professional in vitro diagnostic use only.

The detection of autoantibodies against Sm proteins is a component of the multi-parametric ACR criteria for the diagnosis of systemic lupus erythematosus (SLE). The detection of Sm antibodies serves as a prognostic marker for SLE, there is a relationship between the appearance of Sm antibodies and severe organ manifestations of the disease. Evaluation of a test result should always take into account all clinical and laboratory diagnostic findings.

SYMBOLS USED ON LABELS

	In vitro diagnostic medical device		Microplate
	Manufacturer		Calibrator
	Catalogue number		Calibrator
	Sufficient for ... determinations		Calibrator
	Batch code		Calibrator
	Use by		Calibrator
	Temperature limitation		Calibrator
	Keep away from sunlight		Control positive
	Do not reuse		Control negative
	Date of manufacture		Sample Buffer P
	CE marked according to 98/79/EC		Enzyme Conjugate
	Consult electronic Instructions For Use		TMB Substrate
	Electronic Instruction For Use: version		Stop solution
			Wash Buffer
			Ready to use
			50 x concentrate

PRINCIPLE OF THE TEST

Highly purified Sm is bound to microwells.

The determination is based on an indirect enzyme linked immune reaction with the following steps:
Specific antibodies in the patient sample bind to the antigen coated on the surface of the reaction wells. After incubation, a washing step removes unbound and unspecifically bound serum or plasma components. Subsequently added enzyme conjugate binds to the immobilized antibody-antigen-complexes. After incubation, a second washing step removes unbound enzyme conjugate. After addition of substrate solution the bound enzyme conjugate hydrolyses the substrate forming a blue coloured product. Addition of an acid stops the reaction generating a yellow end-product. The intensity of the yellow color correlates with the concentration of the antibody-antigen-complex and can be measured photometrically at 450 nm.

WARNINGS AND PRECAUTIONS

- All reagents of this kit are intended for professional in vitro diagnostic use only.
- Components containing human serum were tested and found negative for HBsAg, HCV, HIV1 and HIV2 by FDA approved methods. No test can guarantee the absence of HBsAg, HCV, HIV1 or HIV2, and so all human serum based reagents in this kit must be handled as though capable of transmitting infection.
- Bovine serum albumin (BSA) used in components has been tested for BSE and found negative.
- Avoid contact with the substrate TMB (3,3',5,5'-Tetramethyl-benzidine).
- Stop solution contains acid, classification is non-hazardous. Avoid contact with skin.
- Control, sample buffer and wash buffer contain sodium azide 0.09% as preservative. This concentration is classified as non-hazardous.
- Enzyme conjugate contains ProClin 300 0.05% as preservative. This concentration is classified as non-hazardous.

During handling of all reagents, controls and serum samples observe the existing regulations for laboratory safety regulations and good laboratory practice:

- First aid measures: In case of skin contact, immediately wash thoroughly with water and soap. Remove contaminated clothing and shoes and wash before reuse. If system fluid comes into contact with skin, wash thoroughly with water. After contact with the eyes carefully rinse the opened eye with running water for at least 10 minutes. Get medical attention if necessary.
- Personal precautions, protective equipment and emergency procedures:

Observe laboratory safety regulations. Avoid contact with skin and eyes. Do not swallow. Do not pipette by mouth. Do not eat, drink, smoke or apply makeup in areas where specimens or kit reagents are handled. When spilled, absorb with an inert material and put the spilled material in an appropriate waste disposal.

- Exposure controls / personal protection: Wear protective gloves of nitril rubber or natural latex. Wear protective glasses. Used according to intended use no dangerous reactions known.
- Conditions to avoid: Since substrate solution is light-sensitive. Store in the dark.
- For disposal of laboratory waste the national or regional legislation has to be observed.

Observe the guidelines for performing quality control in medical laboratories by assaying control sera.

CONTENTS OF THE KIT

ORG 510	▽ 96	Sufficient for 96 determinations
MICROPLATE	1	One divisible microplate consisting of 12 modules of 8 wells each. Ready to use. Color code on module
CALIBRATOR A	1x 1.5 ml	Calibrator A 0 U/ml, containing serum/buffer matrix (PBS, BSA, detergent, NaN3 0.09%), yellow. Ready to use.
CALIBRATOR B	1x 1.5 ml	Calibrator B 12.5 U/ml, containing Sm antibodies in a serum/buffer matrix (PBS, BSA, detergent, NaN3 0.09%), yellow. Ready to use.
CALIBRATOR C	1x 1.5 ml	Calibrator C 25 U/ml, containing Sm antibodies in a serum/buffer matrix (PBS, BSA, detergent, NaN3 0.09%), yellow. Ready to use.
CALIBRATOR D	1x 1.5 ml	Calibrator D 50 U/ml, containing Sm antibodies in a serum/buffer matrix (PBS, BSA, detergent, NaN3 0.09%), yellow. Ready to use.
CALIBRATOR E	1x 1.5 ml	Calibrator E 100 U/ml, containing Sm antibodies in a serum/buffer matrix (PBS, BSA, NaN3 0.09%), yellow. Ready to use.
CALIBRATOR F	1x 1.5 ml	Calibrator F 200 U/ml, containing Sm antibodies in a serum/buffer matrix (PBS, BSA, detergent, NaN3 0.09%), yellow. Ready to use.
CONTROL +	1x 1.5 ml	Control positive, containing Sm antibodies in a serum/buffer matrix (PBS, BSA, detergent, NaN3 0.09%), yellow. Ready to use. The concentration is specified on the certificate of analysis.
CONTROL -	1x 1.5 ml	Control negative, containing Sm antibodies in a serum/buffer matrix (PBS, BSA, detergent, NaN3 0.09%), yellow. Ready to use. The concentration is specified on the certificate of analysis.
DILUENT	20 ml	Sample Buffer P, containing PBS, BSA, detergent, preservative sodium azide 0.09%, yellow, concentrate (5 x).
CONJUGATE	15 ml	Enzyme Conjugate containing anti-human IgG antibodies, HRP labelled; PBS, BSA, detergent, preservative PROCLIN 0.05%, light red. Ready to use.
TMB	15 ml	TMB Substrate; containing 3,3', 5,5'- Tetramethylbenzidin, colorless. Ready to use.
STOP	15 ml	Stop solution; contains acid. Ready to use.
WASH	20 ml	Wash Buffer, containing Tris, detergent, preservative sodium azide 0.09%; 50 x conc.

MATERIALS REQUIRED

- Microplate reader capable of endpoint measurements at 450 nm; optional: reference filter at 620 nm
- Data reduction software
- Multi-channel dispenser or repeatable pipette for 100 µl
- Vortex mixer
- Pipettes for 10 µl, 100 µl and 1000 µl
- Laboratory timing device
- Distilled or deionised water
- Measuring cylinder for 1000 ml and 100 ml
- Plastic container for storage of the wash solution

This ELISA assay is suitable for use on open automated ELISA processors. Each assay has to be validated on the respective automated system. Detailed information is provided upon request.

SPECIMEN COLLECTION, STORAGE AND HANDLING

- Collect whole blood specimens using acceptable medical techniques to avoid hemolysis.
- Allow blood to clot and separate the serum or plasma by centrifugation.
- Test serum should be clear and non-hemolyzed. Contamination by hemolysis or lipemia should be avoided, but does not interfere with this assay.
- Specimens may be refrigerated at 2-8°C for up to five days or stored at -20°C up to six months.
- Avoid repetitive freezing and thawing of serum or plasma samples. This may result in variable loss of antibody activity.
- Testing of heat-inactivated sera is not recommended.

STORAGE AND STABILITY

- Store test kit at 2-8°C in the dark.
- Do not expose reagents to heat, sun, or strong light during storage and usage.
- Store microplate sealed and dessicated in the clip bag provided.
- Shelf life of the unopened test kit is 18 months from day of production.
Unopened reagents are stable until expiration of the kit. See labels for individual batch.
- Diluted Wash Buffer and Sample Buffer are stable for at least 30 days when stored at 2-8°C.
We recommend consumption on the same day.

PROCEDURAL NOTES

- Do not use kit components beyond their expiration dates.
- Do not interchange kit components from different lots and products.
- All materials must be at room temperature (20-28°C) prior to use.
- Prepare all reagents and samples. Once started, perform the test without interruption.
- Double determinations may be done. By this means pipetting errors may become obvious.
- Perform the assay steps only in the order indicated.
- Always use fresh sample dilutions.
- Pipette all reagents and samples into the bottom of the wells.
- To avoid carryover or contamination, change the pipette tip between samples and different kit controls.
- Wash microwells thoroughly and remove the last droplets of wash buffer.
- All incubation steps must be accurately timed.
- Do not re-use microplate wells.

PREPARATION OF REAGENTS

WASH
Dilute the contents of one vial of the buffered wash solution concentrate (50x) with distilled or deionised water to a final volume of 1000 ml prior to use.

DILUENT
Sample Buffer P: Prior to use dilute the contents (20 ml) of one vial of sample buffer 5x concentrate with distilled or deionised water to a final volume of 100 ml.

Preparation of samples

Dilute patient samples 1:100 before the assay: Put 990 µl of prediluted sample buffer in a polystyrene tube and add 10 µl of sample. Mix well. Note: Calibrators / Controls are ready to use and need not be diluted.

21. Hahn BH, McMahon MA, Wilkinson A, Wallace WD, Daikh DI, Fitzgerald JD et al. American College of Rheumatology guidelines for screening, treatment, and management of lupus nephritis. *Arthritis Care Res (Hoboken)* 2012; 64(6):797-808.

Notice to the user (European Union):

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the EU Member State in which the user and/or the patient is established .

Change Control

Former version: *ORG 510_IFU_EN_QM113137_2018-01-02_3*

Reason for revision: *Definition of symbols used and symbols updated*

LIMITATIONS OF THE PROCEDURE

This assay is a diagnostic aid. A definite clinical diagnosis should not be based on the results of a single test, but should be made by the physician after all clinical and laboratory findings have been evaluated concerning the entire clinical picture of the patient. Also every decision for therapy should be taken individually.

The above pathological and normal reference ranges for antibodies in patient samples should be regarded as recommendations only. Each laboratory should establish its own ranges according to ISO 15189 or other applicable laboratory guidelines.

REFERENCES

1. Alba P, Bento L, Cuadrado MJ, Karim Y, Tungekar MF, Abbs I et al. Anti-dsDNA, anti-Sm antibodies, and the lupus anticoagulant: significant factors associated with lupus nephritis. *Ann Rheum Dis* 2003; 62(6):556-560.
2. Antico A, Platzgummer S, Bassetti D, Bizzaro N, Tozzoli R, Villalta D. Diagnosing systemic lupus erythematosus: new-generation immunoassays for measurement of anti-dsDNA antibodies are an effective alternative to the Farr technique and the Crithidia luciliae immunofluorescence test. *Lupus* 2010; 19(8):906-912.
3. Brouwer R, Hengstman GJ, Vree EW, Ehrfeld H, Bozic B, Ghirardello A et al. Autoantibody profiles in the sera of European patients with myositis. *Ann Rheum Dis* 2001; 60(2):116-123.
4. Castro C, Gourley M. Diagnostic testing and interpretation of tests for autoimmunity. *J Allergy Clin Immunol* 2010; 125(2 Suppl 2):S238-S247.
5. Defendenti C, Atzeni F, Spina MF, Grosso S, Cereda A, Guercilena G et al. Clinical and laboratory aspects of Ro/SSA-52 autoantibodies. *Autoimmun Rev* 2011; 10(3):150-154.
6. Eriksson C, Kokkonen H, Johansson M, Hallmans G, Wadell G, Rantapaa-Dahlqvist S. Autoantibodies predate the onset of Systemic Lupus Erythematosus in northern Sweden. *Arthritis Research & Therapy* 2011; 13(1):R30.
7. Haugbro K, Nossent JC, Winkler T, Figenschau Y, Rekvig OP. Anti-dsDNA antibodies and disease classification in antinuclear antibody positive patients: the role of analytical diversity. *Ann Rheum Dis JID* - 0372355 2004; 63(4):386-394.
8. Ippolito A, Wallace DJ, Gladman D, Fortin PR, Urowitz M, Werth V et al. Autoantibodies in systemic lupus erythematosus: comparison of historical and current assessment of seropositivity. *Lupus* 2011; 20(3):250-255.
9. Isenberg DA, Manson JJ, Ehrenstein MR, Rahman A. Fifty years of anti-ds DNA antibodies: are we approaching journey's end? *Rheumatology (Oxford)* 2007; 46(7):1052-1056.
10. Kattah NH, Kattah MG, Utz PJ. The U1-snRNP complex: structural properties relating to autoimmune pathogenesis in rheumatic diseases. *Immunol Rev* 2010; 233(1):126-145.
11. Kumar Y, Bhatia A, Minz RW. Antinuclear antibodies and their detection methods in diagnosis of connective tissue diseases: a journey revisited. *Diagn Pathol* 2009; 4:1.
12. Meroni PL, Schur PH. ANA screening: an old test with new recommendations. *Ann Rheum Dis* 2010; 69:1420-1422.
13. Petri M, Magder L. Classification criteria for systemic lupus erythematosus: a review. *Lupus* 2004; 13(11):829-837.
14. Poole BD, Schneider RI, Guthridge JM, Velte CA, Reichlin M, Harley JB et al. Early targets of nuclear RNP humoral autoimmunity in human systemic lupus erythematosus. *Arthritis Rheum* 2009; 60(3):848-859.
15. Putova I, Dostal C, Becvar R. Prevalence of antinucleosome antibodies by enzyme-linked immunosorbent assays in patients with systemic lupus erythematosus and other autoimmune systemic diseases. *Ann N Y Acad Sci* 2007; 1109:275-286.
16. Reveille JD. Predictive value of autoantibodies for activity of systemic lupus erythematosus. *Lupus JID* - 9204265 2004; 13(5):290-297.
17. Simon JA, Cabiedes J, Ortiz E, Alcocer-Varela J, Sanchez-Guerrero J. Anti-nucleosome antibodies in patients with systemic lupus erythematosus of recent onset. Potential utility as a diagnostic tool and disease activity marker. *Rheumatology (Oxford)* 2004; 43(2):220-224.
18. Sinclair D, Saas M, Williams D, Hart M, Goswami R. Can an ELISA replace immunofluorescence for the detection of anti-nuclear antibodies?--The routine use of anti-nuclear antibody screening ELISAs. *Clin Lab* 2007; 53(3-4):183-191.
19. Tozzoli R, Bizzaro N, Tonutti E, Villalta D, Bassetti D, Manoni F et al. Guidelines for the laboratory use of autoantibody tests in the diagnosis and monitoring of autoimmune rheumatic diseases. *Am J Clin Pathol* 2002; 117(2):316-324.
20. Maidhof W., Hillas O. Lupus: an overview of the disease and management options. *P T* 2012; 37(4):240-9.

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en

Anti-Salmonella O-Group Pools

Test Reagents for slide agglutination

INFORMATION FOR PROFESSIONAL USE

CE

IVD

	O-Group-Pool	Vol.		O-Group-Pool	Vol.
Anti-Salmonella OMA	A, B, D, E, L	1 ml, 3 ml	Anti-Salmonella OME	X, Y, Z, 51, 52, 53	1 ml
Anti-Salmonella OMB	C, F, G, H	1 ml, 3 ml	Anti-Salmonella OMF	54-59	1 ml
Anti-Salmonella OMC	I, J, K, M, N, O, P	1 ml	Anti-Salmonella OMG	60-63, 65-67	1 ml
Anti-Salmonella OMD	Q, R, S, T, U, V, W	1 ml			

Intended use
The test reagents are intended for use in the serological identification of *Salmonella* strains isolated from test material of varied origin, using slide agglutination. Additional serological differentiation must be performed using the group-specific test reagents and the Anti-Salmonella O, Vi test reagents in accordance with the White-Kauffmann-Le-Minor scheme (Kauffmann-White scheme).

Principle of the test
If the isolated strain possesses an antigen that is covered by the test reagent, this *Salmonella* antigen becomes bound when mixed with the specific antibody. The antigen-antibody reaction results in clearly visible agglutination of the strain.

Composition
The test reagents are mixtures of monoclonal antibodies. They are produced from cell culture supernatants of hybridoma cell lines which secrete antibodies against the respective group specific *Salmonella* O-antigens. For the coverage of the O-Group K, the Anti-Salmonella OMC test reagent contains serum of immunized rabbits from which unspecific agglutinins have been removed by absorption. Preservative: Sodium azide (NaN₃) 0.9 mg/ml

Form in which product is supplied, shelf life and conditions for storage
The test reagents are available in liquid form ready to use. If stored at 2...8 °C before and after opening, they may be used up to the date given on the label. After use, the bottle must be properly closed. The test reagents may sometimes show turbidity which is not microbially caused. Such turbidity does not impair effectivity and the test reagents can be clarified by centrifugation or filtration. The test reagents must have their temperatures adjusted to room temperature (18...26 °C) before use.


Warnings and precautions
The biotechnological manufacture of the test reagents means that the risk of contamination by infectious agents can be virtually excluded. Because they contain animal materials (fetal calf serum, stabilizer) and because of the addition of rabbit serum, they should be treated as potentially infectious and handled accordingly. As these materials contain sodium azide, contact with the skin and mucous membranes must be avoided! In case of contact, rinse with plenty of water. Since the performance of the slide agglutination test involves working with native pathogenic materials, all necessary work protection procedures must be adhered to (risk of infection)!

Materials and equipment not supplied
Glass slides, stirring rods, physiological saline, disposal containers for infectious material, White-Kauffmann-Le-Minor scheme.


Test material and methodology
The test material is streaked onto non-selective culture media (Nutrient Agar, Blood Agar) or Kligler Medium. Incubate for 16 – 20 hours at 35 ...37 °C. Transfer a small amount of bacterial mass from a suspicious colony under test onto a slide and mix well with a drop of test reagent (ca. 25 µl), so that a homogenous, slightly milky suspension results. The slide should be placed on a dark surface. The reaction is read with the naked eye by holding the slide in front of a light source against a black background and rocking it (tilting it back and forth). To control for spontaneous agglutination, a negative control (NC) must be performed at the same time using physiological saline in place of the test reagent.

Test procedure
First carry out the test with the Anti-Salmonella OMA and Anti-Salmonella OMB test reagents. Approx. 98 % of *Salmonella* can be assigned to Anti-Salmonella OMA and Anti-Salmonella OMB. If the strain does not agglutinate with these both test reagents, it is recommended to test this strain with Anti-Salmonella Vi, [REF](#) TR 1316. If this reaction is also negative the strain must then be agglutinated with the Anti-Salmonella OMC, Anti-Salmonella OMD, Anti-Salmonella OME, Anti-Salmonella OMF, and Anti-Salmonella OMG test reagents.

Evaluation
The test can only be evaluated if the negative control (NC) remains milky-opaque.
Positive: visible agglutination after the sample has been tilted back and forth less than 20 times. In a strongly positive reaction (PR), agglutination (coarsely or finely flocculent) appears as soon as the bacterial mass is mixed in. In a weakly positive result, agglutination only appears after the slide has been tilted back and forth 10-20 times.

Negative: if the suspension remains milky-opaque or the reaction begins to occur only after the slide has been tilted back and forth more than 20 times, the result is negative (NR).

Quality assurance during the testing procedure
For the quality control of serological determination of *Salmonella* type-O antigens by slide agglutination test, the good expression of the strain's cell surface antigen is important. The use of strains from interlaboratory tests, or field strains of defined origin that have been fully characterised by an external laboratory, or the use of suitable commercially available *Salmonella* test antigens is therefore recommended for quality control.

Limitations of the method
The test reagents react specifically with individual *Salmonella* strains, i. e. they only agglutinate *Salmonella* strains belonging to the O-Group Pools named in the declaration. In exceptional cases, there is the possibility of cross-reactions with other genera of *Enterobacteriaceae* (e. g. with strains of *Citrobacter* spp., *Hafnia alvei*, *Proteus* spp. or *E. coli*) due to antigen identities or related antigens. For this reason, biochemical examination is required to establish membership of the genus *Salmonella*.

LOT	Batch code (Lot)		Use by YYYY-MM (MM = end of month)
REF	Catalogue number		Temperature limitation
IVD	In Vitro Diagnostic Medical Device		Consult instructions for use
TR	Test reagent		Slide agglutination

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Alkaline Saline Peptone Water ISO

Enrichment medium for detection of *Vibrio* spp, according to ISO 21872.

Instructions For Use

ENGLISH

DESCRIPTION

Alkaline Saline Peptone Water ISO (ASPW) is a liquid medium used for the selective enrichment of *Vibrio parahaemolyticus*, *Vibrio cholerae* and *Vibrio vulnificus* from food, animal feed and environmental samples. This medium complies to ISO 21872-1.

TYPICAL FORMULA*

	(g/litre)
Peptone	20.0
Sodium Chloride	20.0

Final pH 8.6 ± 0.2 at 25°C

*Adjusted and/or supplemented as required to meet performance criteria.

METHOD PRINCIPLE

Peptone provide carbon, nitrogen and amino acids for bacterial growth. The high concentration of sodium chloride along with the alkaline pH provide a favourable environment for growth of *Vibrio* spp while inhibiting other organisms.

PREPARATION

Dehydrated medium Suspend 40.0 g of the powder in 1 liter of distilled or deionized water. Mix well. Heat to boil and shake until completely dissolved. Dispense into final containers, in quantities required for the examination. Sterilize in autoclave at 121°C for 15 minutes.

TEST PROCEDURE

Following the procedure described by ISO 21872-1

1. Take test portions (25 g or 25 ml) and homogenize in a bottle containing 225 ml of ASPW.
2. Incubate at 41.5 ± 1°C and/or 37 ± 1°C for 6 ± 1 h ^(a).
3. Transfer 1 ml of the primary enrichment culture into a tube containing 10 ml of ASPW ^(b).
4. Incubate at 41.5 ± 1°C (*V. parahaemolyticus* and *V. cholerae*) and/or 37 ± 1°C (*V. vulnificus*) for 18 ± 1 h ^(c).
5. From the cultures obtained in the ASPW, inoculate with 1 µl sampling loop the surface of a TCBS agar plate (ref. 11195) to obtain well-isolated colonies (for more details see the IFU of TCBS agar).

Footnotes

- a** For primary enrichment, incubation temperature depends on the target species and state of product:
 - For detection of *V. parahaemolyticus* and *V. cholerae* in fresh product, incubate at 41.5 ± 1°C.
 - For detection of *V. vulnificus* in fresh product, incubate at 37 ± 1°C.
 - For detection of all *Vibrio* species in deep frozen, dried or salted product, incubate at 37 ± 1°C.
- b** Do NOT agitate the sample while taking the aliquot from the surface of the liquid medium.
- c** When choosing the incubation temperature for secondary enrichment consider the target species only.

INTERPRETING RESULTS

Growth in ASPW is indicated by turbidity compared to an uninoculated control.

For interpreting the test results on TCBS agar see the relevant IFU.

APPEARANCE OF THE MEDIUM

Dehydrated medium: free-flowing, homogeneous, beige.

Prepared medium: clear, light amber.

STORAGE

The powder is very hygroscopic, store the powder at 10-30°C, in a dry environment, in its original container tightly closed. Store tubes and bottles at 10-25°C away from light. Do not use the product beyond its expiry date on the label or if product shows any evidence of contamination or any sign of deterioration.

SHELF LIFE

Dehydrated medium: 4 years.

Tubes/Bottles: 2 years.

QUALITY CONTROL

To check the performance of the medium the following reference strains can be used.

Strain		Inoculum	Incubation	Growth
<i>Vibrio parahaemolyticus</i>	WDCM 00185	≤100 CFU	37 ± 1 °C / 18 ± 1 h	Good
<i>Vibrio furnissii</i>	WDCM 00186			Good
<i>Escherichia coli</i>	WDCM 00013	>10 ³ CFU		Inhibited

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The product is available in the various configurations listed below. There may be additional product ref. numbers as well. For an updated listing of available products, visit liofilchem.com

Product	Format	Packaging	Ref.
Alkaline Saline Peptone Water ISO (ASPW)	Tube	100 x 10 ml	26486
Alkaline Saline Peptone Water ISO (ASPW)	Bottle	6 x 225 ml	414130
Alkaline Saline Peptone Water ISO (ASPW)	Dehydrated medium	500 g of powder	610377

This document is available from the online Support Center:

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liofilchem.com/ifu-sds

BIBLIOGRAPHY / BIBLIOGRAFIA

1. EN ISO 11133:2014+Amd1:2018. Microbiology of food, animal feed and water – Preparation, production, storage and performance testing of culture media.
2. ISO 21872-1:2017. Microbiology of the food chain – Horizontal method for the determination of *Vibrio* spp. – Part 1: Detection of potentially enteropathogenic *Vibrio parahaemolyticus*, *Vibrio cholerae* and *Vibrio vulnificus*.
3. Shread P., Donovan T. J., and Lee J. V. (1991) Soc. Gen. Microbiol. Q. 8:184.
4. Cruickshank R. (1968) Medical Microbiology. 11th ed. Livingstone Ltd, London, UK.

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liofilchem@liofilchem.com

Antibiotic Disc

Antibiotic discs for susceptibility tests

Contents	Go to	Page
English	Click here	1
Italiano	Click here	2
Español	Click here	4
Deutsch	Click here	5
Français	Click here	7
Ελληνικά	Click here	8
Português	Click here	10
Polski	Click here	12
Product list	Click here	15
References	Click here	19
AST media	Click here	21

Antibiotic Disc

ENGLISH

Antibiotic discs for susceptibility tests

DESCRIPTION

Antibiotic Disc are paper discs with special features, that are impregnated with antibiotic and used for the susceptibility test according to the Kirby-Bauer antibiotic testing (KB testing or disk diffusion antibiotic sensitivity testing).

Antibiotic Disc are available in a large variety of configurations. Each configuration is available in packages of 50 and 250 tests.

CONTENTS OF THE PACKAGES

Discs in cartridge

The 50-test box contains 1 cartridge with 50 discs packed in a desiccant envelope.

The 250-test box contains 5 cartridges of 50 discs, each cartridge individually packed in a desiccant envelope.

Each package also contains a transparent resealable bag.

Discs in canister

The canister contains 250 discs and a desiccant tablet.

METHOD PRINCIPLE

The discs are applied to the surface of a culture medium inoculated with a pure colony suspension of the microorganism under examination. After incubation, the plates are examined, the inhibition halos around each disc are examined and compared with the standard inhibition haloes: in this way the microorganisms are defined as being susceptible, intermediate or resistant to the tested antimicrobial agents.

COMPOSITION

Liofilchem antimicrobial susceptibility test discs are made of high-quality paper in compliance with WHO and FDA specifications.

The discs are manufactured under the quality systems UNI EN ISO 9001 and EN ISO 13485, and to DIN specification for potency, i.e. the concentration of each antibiotic is within 90-125% of the concentration stated on the disc.

GATHERING AND KEEPING SAMPLES

The colonies that are to be subjected to the susceptibility test are taken up by culture media that have been previously swabbed with the sample under examination. In the case of mixed colonies the bacterial strains must be purified before they are swabbed on the plates for the susceptibility test.

TEST PROCEDURE

1. Allow discs to equilibrate to room temperature before opening the container (cartridge or canister) in order to minimize the condensation on the discs, which could affect long-term stability.
2. Make a suspension of the test organism to the density of a 0.5 McFarland turbidity standard.
3. Using a swab, inoculate a suitable agar plate medium by uniformly spreading the suspension over the entire agar surface.
4. Apply discs firmly to the surface of the inoculated agar plate.
5. Incubate plates in an inverted position at the appropriate temperature, atmosphere and time, following the methodology chosen (e.g. CLSI, EUCAST).
6. Return unused discs to the refrigerator/freezer as soon as the application of the discs has been completed (see STORAGE).

NOTE 1: The medium to be used depends on the organism under investigation and the methodology followed, and must be validated by the media manufacturer for antimicrobial susceptibility testing. A list of recommended agar media can be found at the end of this IFU.

NOTE 2: It is recommended to use the inoculum suspension within 15 minutes of preparation, apply discs within 15 minutes of inoculation and incubate plates within 15 minutes of disc application.

For more details, please refer to the current published standards.

EVALUATING THE RESULTS / QUALITY CONTROL

At the end of the incubation period, measure the inhibition halos and interpret according to the current reference standards:

[Antibiotic Disc Interpretative Criteria and Quality Control \(pdf file\)](#)

CLINICAL INTERPRETATION

The susceptibility test carried out *in vitro* cannot exactly reproduce *in vivo* conditions. Nevertheless, it shows the effect of the concentration of the antibiotic, which varies in the culture medium in relation to the growth of the microbial population.

The final choice of antibiotic to administer to the patient is the responsibility of the clinician who possesses all the information on the patient.

LIMITS

Diffusion susceptibility tests use an *in vitro* technique and cannot therefore reproduce the extremely complex *in vivo* conditions. Nevertheless, it is a useful and important tool that helps the clinician choose the correct therapy. Many variable factors influence the final result of the diffusion susceptibility test. The main ones are: the culture medium used, impregnation of the discs, inoculation of the medium, temperature, time and incubation atmosphere of the plates, pre-incubation and pre-diffusion conditions, depth of the medium, etc.

PRECAUTIONS

The Antibiotic Disc cannot be classified as being hazardous according to current legislation. Antibiotic Disc are disposable products. Antibiotic Disc are only for diagnostic *in vitro* use and are intended for professional use. They must be used in the laboratory by properly trained operators using approved aseptic and safety methods for pathogenic agents.

STORAGE

The unopened package of Antibiotic Disc can be stored in most cases at -20°C to $+8^{\circ}\text{C}$ till the expiry date. Some products have to be stored at -20°C as maximum storage temperature. The recommended temperature limits can be found both on the product envelop and on the box label. Leftover discs from an opened CARTRIDGE need to be stored at $2-8^{\circ}\text{C}$ for no more than 7 days. The cartridge containing unused discs should be returned into its desiccant envelope and then inserted into the resealable bag. Discs in a CANISTER can be used for up to 2 months from first opening and must be stored at the label storage temperature. Dispose of expired discs.

ELIMINATING USED MATERIAL

After use, Antibiotic Disc and the material that comes into contact with the sample must be decontaminated and disposed of in accordance with current laboratory techniques for the decontamination and disposal of potentially infected material.

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The following list of products might be out-of-date

View the complete range of Antimicrobial Discs in [cartridge](#) and [canister](#) on Liofilchem's website

Antibiotic discs in cartridge and canister			CLSI ¹	EUCAST ^{3,4}	Ref.*
Description		µg			
Amikacin	AK	30	✓	✓	9004
Amoxicillin	AML	2			9151
Amoxicillin	AML	10		✓	9133
Amoxicillin	AML	25			9179
Amoxicillin	AML	30			9005
Amoxicillin-clavulanic acid	AUG	3 (2/1)		✓	9191
Amoxicillin-clavulanic acid	AUG	7.5			9255
Amoxicillin-clavulanic acid	AUG	30 (20/10)	✓	✓	9048
Amoxicillin 2 + Clavulanic acid 0.1	AC	2.1 (2/0.1)			9273 ♦
Amoxicillin 2 + Clavulanic acid 0.5	AC	2.5 (2/0.5)			9274 ♦
Amoxicillin 5 + Clavulanic acid 0.1	AC	5.1 (5/0.1)			9275 ♦
Amoxicillin 5 + Clavulanic acid 0.5	AC	5.5 (5/0.5)			9276 ♦
Amoxicillin 5 + Clavulanic acid 1	AC	6 (5/1)			9277 ♦
Amoxicillin 10 + Clavulanic acid 0.1	AC	10.1 (10/0.1)			9278 ♦
Amoxicillin 10 + Clavulanic acid 0.5	AC	10.5 (10/0.5)			9279 ♦
Amoxicillin 10 + Clavulanic acid 1	AC	11 (10/1)			9280 ♦
Ampicillin	AMP	2		✓	9115
Ampicillin	AMP	10	✓	✓	9006
Ampicillin-sulbactam	AMS	20 (10/10)	✓	✓	9031
Ampliclox (Ampicillin + Cloxacillin)	ACL	30 (25/5)			9122
Azithromycin	AZM	15	✓	✓	9105
Azlocillin	AZL	75	✓		9007
Aztreonam	ATM	30	✓	✓	9008
Bacitracin	BA	10 units			9051
Carbenicillin	CAR	100	✓		9009
Cefaclor	CEC	30	✓	✓	9010
Cefadroxil	CDX	30		✓	9052
Cefamandole	MA	30	✓		9014
Cefazolin	KZ	30	✓	✓	9015
Cefepime	FEP	10			9220
Cefepime	FEP	30	✓	✓	9104
Cefepime + Clavulanic acid	FEL	40 (30/10)		✓ ¹⁰	9143
Cefiderocol	FDC	30	✓	✓	9266
Cefixime	CFM	5	✓	✓	9089
Cefoperazone	CFP	30			9016
Cefoperazone	CFP	75	✓		9108
Cefotaxime	CTX	5		✓	9152
Cefotaxime	CTX	30	✓		9017
Cefotaxime	CTX	75			9134
Cefotaxime + Clavulanic acid	CTL	40 (30/10)		✓ ¹⁰	9182
Cefotaxime + Clavulanic acid + Cloxacillin	CTLCL			✓ ¹⁰	9203
Cefotaxime + Cloxacillin	CTC	230 (30/200)		✓ ¹⁰	9224
Cefotetan	CTT	30	✓		9081
Cefoxitin	FOX	30	✓	✓	9018
Cefoxitin + Cloxacillin	FOC	230 (30/200)			9144
Cefpirome	CR	30			9185
Cefpodoxime	PX	10	✓	✓	9064
Cefpodoxime + Clavulanic acid	PXL	11 (10/1)			9190
Cefprozil	CPR	30	✓		9112
Cefsulodin	CSD	30			9053
Ceftaroline	CPT	5		✓	9195
Ceftaroline	CPT	30	✓		9198
Ceftazidime	CAZ	10		✓	9153
Ceftazidime	CAZ	30	✓		9019
Ceftazidime-avibactam	CZA	14 (10/4)		✓	9206

Antibiotic discs in cartridge and canister						
Description		µg		CLSI ¹	EUCAST ^{3,4}	Ref.*
Ceftazidime-avibactam	CZA	50	(30/20)	✓		9205
Ceftazidime + Clavulanic acid	CAL	20	(20/10)		✓ ⁹	9258 ♦
Ceftazidime + Clavulanic acid	CAL	40	(30/10)		✓ ¹⁰	9145
Ceftazidime + Clavulanic acid + Cloxacillin	CALC				✓ ¹⁰	9204
Ceftazidime + Cloxacillin	CAC				✓ ¹⁰	9225
Ceftibuten	CTB	30		✓	✓	9101
Ceftizoxime	CZX	30		✓		9054
Ceftobiprole	BPR	5		✓	✓	9242
Ceftolozane-tazobactam	C/T	40	(30/10)	✓	✓	9246
Ceftriaxone	CRO	30		✓	✓	9020
Cefuroxime	CXM	1				9232
Cefuroxime	CXM	5				9236 ♦
Cefuroxime	CXM	30		✓	✓	9021
Cephalexin	CL	30			✓	9011
Cephalothin	KF	30		✓		9013
Cephradine	CE	30				9055
Chloramphenicol	C	10				9128
Chloramphenicol	C	30		✓	✓	9022
Cinoxacin	CIN	100		✓		9057
Ciprofloxacin	CIP	5		✓	✓	9056
Clarithromycin	CLR	15		✓		9098
Clavulanic acid	CLA	1				9229 ♦
Clavulanic acid	CLA	2				9228 ♦
Clavulanic acid	CLA	5				9230 ♦
Clavulanic acid	CLA	10				9231 ♦
Clindamycin	CD	2		✓	✓	9047
Clindamycin	CD	10				9146
Cloxacillin	CX	5				9058
Colistin sulfate	CS	10		✓		9023 ♦
Colistin sulfate	CS	25				9184 ♦
Colistin sulfate	CS	30	units			9141 ♦
Daptomycin (includes Ca ²⁺)	DAP	30				9090
Dicloxacillin	DCX	1				9093
Dipicolinic acid	DP					9194
Doripenem	DOR	10		✓	✓	9154
Doxycycline	DXT	30		✓		9059
EDTA	ED					9087
Eravacycline	ERV	20		✓	✓	9238
Eravacycline	ERV	50				9237 ♦
Ertapenem	ETP	10		✓	✓	9061
Ertapenem + Cloxacillin	ET+CL					9199
Ertapenem + Phenylboronic acid	ET+BO					9202
Erythromycin	E	2				9180
Erythromycin	E	15		✓	✓	9024
Fosfomycin	FOS	50				9025
Fosfomycin (includes G6P 50)	FOS	100	(100/50)			9121
Fosfomycin (includes G6P 50)	FOS	200	(200/50)	✓	✓	9109
Fosfomycin 200 + G6P 200	FGP	200	(200/200)			9214 ♦
Furazolidone	FR	50				9099
Fusidic acid	FC	10		✓	✓	9049
Fusidic acid	FC	30				9111
Gatifloxacin	GAT	5		✓		9169
Gentamicin	CN	10		✓	✓	9026
Gentamicin	CN	30			✓	9125
Gentamicin	CN	120		✓		9124
Gentamicin	CN	500				9288 ♦
Imipenem	IMI	10		✓	✓	9079
Imipenem-relebactam	I/R	35	(10/25)	✓	✓	9253
Imipenem + Cloxacillin	IMI+CL					9086
Imipenem + EDTA	IM+ED					9183

Antibiotic discs in cartridge and canister						
Description		µg		CLSI ¹	EUCAST ^{3,4}	Ref.*
Imipenem + Phenylboronic acid	IMI+BO					9085
Kanamycin	K	30		✓		9027
Lefamulin	LMU	5			✓	9249 ♦
Lefamulin	LMU	20		✓		9250 ♦
Levofloxacin	LEV	5		✓	✓	9102
Levonadifloxacin	LND	10		✓		9267
Lincomycin	MY	2				9028
Lincomycin	MY	15				9116
Linezolid	LNZ	10			✓	9155
Linezolid	LNZ	30		✓		9136
Lomefloxacin	LOM	10		✓		9113
Mecillinam	MEC	10		✓	✓	9156
Meropenem	MRP	10		✓	✓	9068
Meropenem + Cloxacillin	MR+CL				✓ ¹⁰	9175
Meropenem + EDTA	MR+ED				✓ ¹⁰	9178
Meropenem + Phenylboronic acid	MR+BO				✓ ¹⁰	9176
Methicillin	MET	5				9029
Metronidazole	MTZ	5			✓	9076
Metronidazole	MTZ	50				9119
Mezlocillin	MEZ	75				9062
Minocycline	MN	30		✓	✓	9030
Moxifloxacin	MXF	5		✓	✓	9103
Mupirocin	MUP	5				9189
Mupirocin	MUP	200		✓	✓	9157
Nafcillin	NAF	1		✓		9174
Nalidixic acid	NA	30		✓	✓	9001
Neomycin	N	10			✓	9287 ♦
Netilmicin	NET	10			✓	9170
Netilmicin	NET	30		✓		9033
Nitrofurantoin	F	50				9181
Nitrofurantoin	F	100			✓	9158
Nitrofurantoin	F	300		✓		9034
Nitroxoline	NI	30			✓	9209
Norfloxacin	NOR	10		✓	✓	9035
Novobiocin	NO	5				9117
Novobiocin	NO	30				9063
Ofloxacin	OFX	5		✓	✓	9080
Omadacycline	OMC	30		✓		9252 ♦
Oritavancin	ORI	5				9201
Oxacillin	OX	1		✓	✓	9036
Oxacillin	OX	5				9135
Oxolinic acid	OA	2				9002
Oxytetracycline	OT	30				9065
Pefloxacin	PEF	5		✓	✓	9091
Penicillin G	P	1 unit			✓	9130
Penicillin G	P	2 units				9127
Penicillin G	P	10 units		✓		9037
Penicillin G 1 + Clavulanic acid 2	PC	2 (1 unit/2)				9227 ♦
Penicillin G 1 + Clavulanic acid 5	PC	5 (1 unit/5)				9240 ♦
Penicillin G 1 + Clavulanic acid 10	PC	10 (1 unit/10)				9241 ♦
Penicillin G 1 + Clavulanic acid 20	PC	20 (1 unit/20)				9254 ♦
Phenoxymethylpenicillin	PV	10				9171
Phenylboronic acid	BO					9193
Pipemidic acid	PI	20				9003
Piperacillin	PRL	30			✓	9159
Piperacillin	PRL	100		✓		9038
Piperacillin-tazobactam	TZP	36 (30/6)			✓	9160
Piperacillin-tazobactam	TZP	110 (100/10)		✓		9100
Polymyxin B	PB	100 units				9066
Polymyxin B	PB	300 units		✓		9120

Antibiotic discs in cartridge and canister					
Description		µg	CLSI ¹	EUCAST ^{3,4}	Ref.*
Quinupristin-dalfopristin	QDA	15	✓	✓	9161
Rifampicin	RD	5	✓	✓	9118
Rifampicin	RD	30			9039
Roxithromycin	RXT	15			9060
Sisomicin	SIS	30			9046
Sodium Fusidate	FC	30			9131
Spectinomycin	SPC	100	✓		9067
Streptomycin	S	10	✓		9040
Streptomycin	S	300	✓	✓	9162
Sulbactam	SU	20			9129
Sulfadiazine	SUZ	300			9150
Sulfafurazole	SF	300			9041
Sulfamethoxazole	SMX	50			9084
Sulfamethoxazole	SMX	100			9187
Sulfaprim	SXT	50			9132
Sulfonamide	S3	300	✓		9126
Tedizolid	TZD	2	✓	✓	9243
Tedizolid	TZD	20			9245
Teicoplanin	TEC	30	✓	✓	9050
Telithromycin	TEL	15	✓	✓	9172
Temocillin	TMO	30		✓	9186
Tetracycline	TE	30	✓	✓	9043
Ticarcillin	TC	75	✓	✓	9070
Ticarcillin-clavulanic acid	TTC	85 (75/10)	✓	✓	9096
Tigecycline	TGC	15	✓	✓	9147
Tobramycin	TOB	10	✓	✓	9044
Tobramycin	TOB	30			9163
Trimethoprim	TM	2.5			9083
Trimethoprim	TM	5	✓	✓	9110
Trimethoprim-sulfamethoxazole	SXT	2.5 (0.125/2.375)			9299 ♦
Trimethoprim-sulfamethoxazole	SXT	25 (1.25/23.75)	✓	✓	9042
Vancomycin	VA	5		✓	9164
Vancomycin	VA	30	✓		9045

●, to be stored at -20°C ♦, not CE marked

Storage Unless otherwise indicated, discs may be stored at any temperature between -20°C and +8°C.

CE Marking Unless otherwise indicated, all products intended for clinical applications are CE marked according to the European Directive 98/79/EC for In Vitro Diagnostic Medical Devices.

* **Packaging** 5 cartridges of 50 discs (5x50 discs = 250 discs).

Single cartridges of 50 discs are available: add /1 to the catalogue ref. no. to indicate the relevant item.

Example: ref. 9045/1 indicates Vancomycin 30 µg in one single cartridge of 50 discs.

250 discs in a CANISTER: add /2 to the catalogue ref. no. to indicate the relevant item.

Example: ref. 9045/2 indicates Vancomycin 30 µg in a canister of 250 discs.

Note that testing and reporting antimicrobial agents for which there are no interpretive criteria are the responsibility of the chief microbiologist and such decisions should be made with input from the infectious disease clinicians.

For full details on specific organism/agent combinations refer to current CLSI and EUCAST recommendations.

Also see the EUCAST guidance at http://www.eucast.org/clinical_breakpoints/when_there_are_no_breakpoints/

Antibiotic discs in cartridge and canister				
Veterinary				
Description		µg	CLSI ⁷	Ref.*
Aminosidine	AM	60		9301
Apramycin	AP	15	✓	9300
Ceftiofur	FUR	30	✓	9251
Enrofloxacin	ENR	5	✓	9233
Florfenicol	FFC	30	✓	9234
Flumequine	UB	30		9208
Lincomycin-spectinomycin	MSP	100		9302
Marbofloxacin	MAR	5	✓	9297
Neomycin	N	30	✓	9032
Spectinomycin	SPC	100	✓	9067
Spiramycin	SP	100		9088
Tiamulin	T	30	✓	9094
Tilmicosin	TIL	15	✓	9298
Tylosin	TY	30	✓	9082

* **Packaging** 5 cartridges of 50 discs (5x50 discs = 250 discs).
 Single cartridges of 50 discs are available: add /1 to the catalogue ref. no. to indicate the relevant item.
 Example: ref. 9082/1 indicates Tylosin 30 µg in one single cartridge of 50 discs.

250 discs in a CANISTER: add /2 to the catalogue ref. no. to indicate the relevant item.
 Example: ref. 9082/2 indicates Tylosin 30 µg in a canister of 250 discs.

REFERENCES

1. CLSI. Performance Standards for Antimicrobial Susceptibility Testing; 33rd ed. CLSI Supplement M100S. Wayne, PA: Clinical and Laboratory Standards Institute; 2023.
2. CLSI. Performance Standards for Antimicrobial Disk Susceptibility Tests; 13th ed. CLSI standard M02. Wayne, PA: Clinical and Laboratory Standards Institute; 2018.
3. The European Committee on Antimicrobial Susceptibility Testing. Breakpoint tables for interpretation of MICs and zone diameters. Version 13.2, 2023. <http://www.eucast.org>
4. The European Committee on Antimicrobial Susceptibility Testing. Routine and extended internal quality control for MIC determination and disk diffusion as recommended by EUCAST, Version 13.2, 2023. <http://www.eucast.org>
5. EUCAST Disk Diffusion Method for Antimicrobial Susceptibility Testing - Version 11.0, January 2023.
6. EUCAST Disk Diffusion Anaerobic Bacteria - Version 2.0, January 2023.
7. CLSI. Performance Standards for Antimicrobial Disk and Dilution Susceptibility Tests for Bacteria Isolated From Animals, 4th ed. CLSI supplement VET08. Wayne, PA: Clinical and Laboratory Standards Institute; 2018.
8. CLSI. Performance Standards for Antimicrobial Disk and Dilution Susceptibility Tests for Bacteria Isolated From Animals; 5th ed. CLSI standard VET01. Wayne, PA: Clinical and Laboratory Standards Institute; 2018.
9. EUCAST technical guidance on the use of the combination disk test (CDT) for confirmation of ESBL in Enterobacterales. New disk potencies for combination disks containing cefotaxime and ceftazidime without and with clavulanic acid (February 2019).
10. EUCAST guidelines for detection of resistance mechanisms and specific resistances of clinical and/or epidemiological importance. Version 2.0, July 2017.
11. DIN 58940-2 Medical microbiology - Susceptibility testing of microbial pathogens to antimicrobial agents - Part 2: Active substance carriers for the agar diffusion test; 2007-10.
12. FDA (1978) Codes of Fed.Rebs. 21.Part 460.
13. WHO (1977) Tech rep.Ser.n°610.

Agar media for AST of bacteria (*)

Description		CLSI	EUCAST	Packaging	REF
Mueller Hinton II Agar (cation-adjusted), for non-fastidious bacteria	MHA	✓	✓	20 plates 90 mm 10 plates 140 mm	10031 10231
Mueller Hinton II Agar with 5% Sheep Blood, for fastidious organisms, including <i>Pasteurella multocida</i> and <i>Mannheimia haemolytica</i>	MHA + 5% blood	✓		20 plates 90 mm 10 plates 140 mm	10131 11231
Mueller Hinton Fastidious Agar: MHA with 20 mg/l β-NAD and 5% Horse Blood, for <i>Streptococcus</i> spp., <i>Haemophilus</i> spp., and other fastidious organisms	MH-F		✓	20 plates 90 mm 10 plates 140 mm	10132 11132
Haemophilus Test Agar, for <i>Haemophilus</i> spp.	HTM	✓		20 plates 90 mm	10080
GC Agar Base with 1% defined growth supplement, for <i>Neisseria gonorrhoeae</i>		✓			
Mueller Hinton Chocolate Agar, for <i>Actinobacillus pleuropneumoniae</i> and <i>Histophilus somni</i>	Chocolate MHA	✓		20 plates 90 mm	10335
Fastidious Anaerobe Agar with defibrinated Horse Blood, for selected rapidly growing anaerobic bacteria (i.e. <i>Bacteroides</i> spp., <i>Prevotella</i> spp., <i>F. necrophorum</i> , <i>C. perfringens</i> and <i>C. acnes</i>)	FAA-HB		✓	20 plates 90 mm	10062

*Available as ready to use plated media in the Liofilchem Catalog.

NOTE: MHA, MH-F and FAA-HB into Petri dishes should have a level depth of 4 ± 0.5 mm.

Find out more on Disc Diffusion Antimicrobial Susceptibility Tests

Liofilchem® and the Liofilchem company logo are registered trademarks of LIOFILCHEM s.r.l.
© Copyright LIOFILCHEM 2023

F01511

LIOFILCHEM® s.r.l.

Via Scozia, 64026 Roseto degli Abruzzi (TE) Italy
Tel. +39 0858930745 Fax +39 0858930330 www.liofilchem.com

Bile Aesculin Agar

Selective medium for the preliminary identification of group D streptococci.

TYPICAL FORMULA (g/l)

Beef Extract	3.0
Meat Peptone	5.0
Esculin	1.0
Oxgall	40.0
Ferric Citrate	0.5
Agar	14.0
Final pH 6.4 ± 0.2 at 25°C	

DESCRIPTION

Bile Aesculin Agar is used for the selective isolation and differentiation of enterococci and streptococci (serology group D) from food and pharmaceutical products. This medium is not intended for use in the diagnosis of disease or other conditions in humans.

PRINCIPLE

Beef extract and meat peptone provide amino acids, nitrogen, carbon, vitamins and minerals. Oxgall is ox bile purified and dehydrated. It contains a mix of biliary salts and is used in media for enterobacteria, as selective agent, because inhibit Gram-positive bacteria other than Group D streptococci. Esculin is a glucoside which is hydrolysed by Group D streptococci to form aesculetin and dextrose. Aesculetin combines with ferric citrate in the medium to form a dark brown or black complex which is indicative of a positive result. Agar is the solidifying agent.

PREPARATION

Suspend 63.5 g of powder in 1 liter of deionized or distilled water. Bring to boil and shake until completely dissolved. Sterilize at 121°C for 15 minutes. Cool up to 45-50°C. Dispense in Petri dishes.

TECHNIQUE

Inoculate the medium by streaking directly the sample onto the agar surface.

Incubate at 36 ± 1°C for 18-24 hours in aerobic atmosphere.

INTERPRETATION OF RESULTS

Group D streptococci grow on this medium producing a dark brown color around the colonies.

STORAGE

The powder is very hygroscopic, store the powder at 10-30°C, in a dry environment, in its original container tightly closed and use it before the expiry date on the label or until signs of deterioration or contamination are evident. Store prepared plates at 2-8°C away from light.

WARNING AND PRECAUTIONS

For professional use only. Operators must be trained and have certain experience in the laboratory methods. Please read the instructions carefully before using this product. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this document.

Consult the Safety Data Sheet (SDS) for information regarding hazards and safe handling practices.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to the national and local regulations in force.

REFERENCES

- European Pharmacopoeia 9.0 (2017) - 2.6.12. Microbiological examination of non-sterile products: Microbial enumeration tests. 2.6.13. Microbiological examination of non-sterile products: Test for specified microorganisms.
- Swan, A. 1954. The use of bile-esculin medium and of Maxted's technique of Lancefield grouping in the identification of enterococci (group D streptococci). J. Clin. Pathol. 7:160.
- Facklam, R. R., and M. D. Moody. 1970. Presumptive identification of group D streptococci: the bileesculin test. Appl. Microbiol. 20:245.
- Rochaix, A. 1924. Milieux a leculine pour le diagnostid differentiel des bacteries du grojps strepto-enteropneumocoque. Compt. Rend. Soc. Biol. 90:771-772.
- Meyer, K., and H. Schönfeld. 1926. Über die Unterscheidung des Enterococcus vom Streptococcus viridans und die Beziehung beider zum Strptoccus lactis. Zentralb. Bakteriell Parasitenkd. Infektionskr. Hyg. Abt. I orig. 99:402-416.
- Schleifer, K. H., and R. Kilpper-Balz. 1987. Molecular and chemotaxonomic approaches to the classification of streptococci, enterococci and lactococci: a review. Syst. Appl. Microbiol. 10:1-19.
- Vanderzant, C., and D. F. Splittstoesser (eds.). Compendium of methods for the microbiological examination of foods, 4th ed. American Public Health Association, Washington, D.C.
- FDA-BAM. www.fda.gov/Food/ScienceResearch/LaboratoryMethods/BacteriologicalAnalytical manualBAM/ default.htm.
- Marshall, R. T. (ed.). 2004. Standard methods for the examination of dairy products, 17th ed. American Public Health Association, Washington, D.C.

LIOFILCHEM® S.r.l.

Via Scozia, 64026 Roseto degli Abruzzi (TE) ITALY
Tel +39 0858930745 Fax +39 0858930330 www.liofilchem.com liofilchem@liofilchem.com

PRODUCT SPECIFICATIONS

NAME

Bile Aesculin Agar

PRESENTATION

Dehydrated medium

STORAGE

10-30°C

PACKAGING

Ref.	Content	Packaging
620210	100 g	100 g of powder in plastic bottle
610210	500 g	500 g of powder in plastic bottle
6102105	5 kg	5 kg of powder in plastic bottle

pH OF THE MEDIUM

6.4 ± 0.2

USE

Bile Aesculin Agar is used for the selective isolation and differentiation of enterococci and streptococci (serology group D) from food and pharmaceutical products

TECHNIQUE

Refer to technical sheet of the product

APPEARANCE OF THE MEDIUM

Powder medium

Appearance: free-flowing, homogeneous

Colour: greenish, light to medium beige

Ready-to-use medium

Appearance: slightly opalescent

Colour: greenish to medium amber

SHELF LIFE

4 years

QUALITY CONTROL

- Control of general characteristics, label and print
- Microbiological control
Inoculum for productivity: 50-100 CFU
Inoculum for selectivity: 10⁴-10⁶ CFU
Incubation Conditions: 24 h/ 36 ± 1°C

Microorganism

Enterococcus faecalis ATCC 19433

Streptococcus pyogenes ATCC 19615

Growth

Good

Inhibited

Characteristics

Blackening of the medium

TABLE OF SYMBOLS

LOT	Batch code		Consult instructions for use		Manufacturer		Use by
REF	Catalogue number		Temperature limitation		Contains sufficient for <n> tests		Keep away from sunlight

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MacConkey Broth

Differential medium for detection of coliform bacteria.

DESCRIPTION

MacConkey Broth is a liquid medium used for the detection and enumeration of coliforms from water, and milk.

TYPICAL FORMULA*	(g/l)
Peptone	20.0
Lactose	10.0
Sodium Chloride	5.0
Bile Salts No.3	5.0
Bromocresol Purple	0.01
Final pH 7.4 ± 0.2 at 25°C	

*Adjusted and/or supplemented as required to meet performance specifications.

METHOD PRINCIPLE

Peptone provides amino acids, nitrogen, carbon, vitamins and minerals for organisms growth. Lactose is the fermentable carbohydrate. Sodium Chloride maintain the osmotic balance of the medium. Bile salts inhibit the growth of Gram-positive bacteria. Bromocresol purple is the pH indicator.

PREPARATION

Dehydrated medium

Suspend 40.0 g of the powder in 1 liter of distilled or deionized water. Mix well. Heat to boil shaking frequently until completely dissolved. Dispense into tubes containing Durham tubes. Autoclave at 121°C for 15 minutes.

Ensure there are no air bubbles in the Durham tube prior to use as it may lead to a misinterpretation (false-positive reaction in gas production) of results. If needed, invert the tube to allow the air being released from the Durham tube.

TEST PROCEDURE

For the presumptive coliform examination, inoculate a series of tubes with the test sample and its dilutions and incubate at 35-37°C for 18-24 hours.

The number of total coliforms is estimated by counting the number of tubes giving a positive reaction for acid (yellow colouration of the broth) and gas production and comparing the pattern of positive results with standard charts to find MPN (Most Probable Number) of total 100 ml sample.

For the differential test, subculture from each tube showing positive reaction into a fresh tube of MacConkey Broth. Incubate at 42-44°C for up to 48 h to detect *Escherichia coli*.

INTERPRETING RESULTS

Turbidity of MacConkey Broth indicates microbial growth.

Acid production due to lactose fermentation causes a color change of the medium to yellow. Gas is also produced, collecting in Durham tubes.

The ability to form gas in MacConkey Broth at 44°C is specific for *E. coli*, which is indicative of faecal pollution of the original sample.

Non-fermenting Gram-negative organisms produce good growth but without significant changes in color nor gas formation.

STORAGE

The powder is very hygroscopic, store the powder at 10-30°C, in a dry environment, in its original container tightly closed. Do not use the product beyond its expiry date on the label or if product shows any evidence of contamination or any sign of deterioration.

SHELF LIFE

Dehydrated medium: 4 years.

QUALITY CONTROL

Appearance of Dehydrated Medium: Free-flowing, homogeneous, light beige.

Appearance of Prepared Medium: Clear, purple.

Expected Cultural Response:

Control strain		Inoculum	Incubation	Growth	Acid	Gas
Salmonella Typhimurium	WDCM 00031 (ATCC® 14028; NCTC 12023)	≤100 CFU	18-24 h / 35-37°C	Good	Neg. (purple)	Neg.
Escherichia coli	WDCM 00012 (ATCC® 8739; NCTC 12923)	≤100 CFU	18-24 h / 42-44°C	Good	Pos. (yellow)	Pos.
Staphylococcus aureus	WDCM 00032 (ATCC® 6538; NCTC 10788)	1000 CFU	48 h / 35-37°C	Inhibition	—	—

Please refer to the actual batch related Certificate of Analysis (CoA).

WARNING AND PRECAUTIONS

For professional use only. Operators must be trained and have certain experience in the laboratory methods. Please read the instructions carefully before using this product. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this document.

Consult the Safety Data Sheet (SDS) for information regarding hazards and safe handling practices.

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BIBLIOGRAPHY

See the references at the end of this document.

TABLE OF SYMBOLS

See the table of symbols at the end of this document.

The product is available in the configurations listed below. There may be additional product ref. numbers as well. For an updated listing of available products, visit [liofilchem.com](https://www.liofilchem.com)

Product	Format	Packaging	Ref.
MacConkey Broth	Dehydrated medium	100 g	620171
		500 g	610171

Significant changes from previous version:

Document	Release Date	Change Summary
610171_IFU-0	2022-11-04	Document creation

This IFU document and the SDS are available from the online Support Center:

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liofilchem.com/ifu-sds

References / Riferimenti

1. European Pharmacopoeia (EP) 2.6.13. Microbiological examination of non-sterile products: Test for specified microorganisms.
2. United States Pharmacopoeia (USP) <62> Microbiological examination of non-sterile products: Test for specified microorganisms.
3. Japanese Pharmacopoeia (JP) 4.05 Microbiological examination of non-sterile products: Test for specified microorganisms.
4. Departments of the Environment, Health, Social Security and Public Health Laboratory Service (1982) The Bacteriological Examination of Drinking Water Supplies, Report No. 71, HMSO London.
5. World Health Organization (1963) International Standards for Drinking Water, 2nd ed., WHO, Geneva.
6. Murray, Baron, Jorgensen, Landry and Pfaller ed. (2007) Manual of clinical microbiology, 9th ed. American Society for Microbiology, Washington, D.C.
7. MacConkey, A. (1908) Bile salt media and their advantages in some bacteriological examinations. J. Hyg. 8: 322-334.
8. MacConkey, A. (1905) Lactose-fermenting bacteria in faeces. J. Hyg. 8: 333-379.

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CIC-C1q ELISA (IgG)

Test instruction

ORDER NO.	Circulating Immune Complexes (CIC-C1q)	SUBSTRATE	FORMAT
EA 1818-9601 G		C1q-coated microplate wells	96 x 01 (96)

Indication: The ELISA test kit provides a semiquantitative or quantitative in vitro assay against C1q-binding-circulating immune complexes (CIC) containing IgG in human serum or plasma for the diagnosis of rheumatic, autoimmune and infectious diseases, allergy, haematological and neoplastic disorders.

Application: Increased levels of circulating immune complexes can be detected in patients with different autoimmune, microbial and neoplastic diseases. Immune complexes are produced as a consequence of the body's immune response. Investigation of these complexes can therefore be useful for assessing the disease activity, for therapy control and for establishing a prognosis. Moreover, circulating immune complexes can become pathogenetically active through complement activation.

Principle of the test: The test kit contains microplate strips, each with 8 break-off reagent wells coated with complement fragment C1q. In the first reaction step, diluted patient samples are incubated in the wells. If samples are positive, circulating immune complexes bind to the C1q via the Fc fragment of the immunoglobulins. To detect the bound immune complexes, a second incubation is carried out using an enzyme-labelled anti-human IgG (enzyme conjugate) catalysing a colour reaction. The intensity of the colour formed is proportional to the concentration of circulating immune complexes against C1q.

Contents of the test kit:

Component	Colour	Format	Symbol
1. Microplate wells coated with C1q 12 microplate strips each containing 8 individual break-off wells in a frame, ready for use	---	12 x 8	
2. Calibrator 1 200 RU/ml (CIC-C1q, human), lyophilised	colourless	1 x 2.0 ml	
3. Calibrator 2 20 RU/ml (CIC-C1q, human), lyophilised	colourless	1 x 2.0 ml	
4. Calibrator 3 2 RU/ml (CIC-C1q, human), lyophilised	colourless	1 x 2.0 ml	
5. Positive control (CIC-C1q, human), lyophilised	colourless	1 x 2.0 ml	
6. Negative control (CIC-C1q, human), lyophilised	colourless	1 x 2.0 ml	
7. Enzyme conjugate peroxidase-labelled anti-human IgG (rabbit), ready for use	light green	1 x 12 ml	
8. Sample buffer ready for use	bluish	1 x 100 ml	
9. Wash buffer 10x concentrated	colourless	1 x 100 ml	
10. Chromogen/substrate solution TMB/H ₂ O ₂ , ready for use	colourless	1 x 12 ml	
11. Stop solution 0.5 M sulphuric acid, ready for use	colourless	1 x 12 ml	
12. Plastic foil	---	2 pieces	
13. Test instruction	---	1 booklet	
14. Quality control certificate	---	1 protocol	
Lot description			Storage temperature
In vitro diagnostic medical device			Unopened usable until

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- **Coated wells:** Ready for use. Tear open the resealable protective wrapping of the microplate at the recesses above the grip seam. Do not open until the microplate has reached room temperature to prevent the individual strips from moistening. Immediately replace the remaining wells of a partly used microplate in the protective wrapping and tightly seal with the integrated grip seam (Do not remove the desiccant bag).
Once the protective wrapping has been opened for the first time, the wells coated with C1q can be stored in a dry place and at a temperature between +2°C and +8°C for 4 months.
- **Calibrators and controls:** Lyophilised. Reconstitute the content of each vial with 2 ml deionised or distilled water. The reconstituted calibrators and controls must be mixed thoroughly before use.
The reconstituted calibrators and controls are stable at +18°C to +25°C for up to 4 hours. If the reconstituted calibrators and controls are to be stored for a longer period of time, they should be kept aliquotated at -20°C or lower until expiry date. Repeated freeze-thaw cycles of calibrators and controls should be avoided.
- **Enzyme conjugate:** Ready for use. The enzyme conjugate must be mixed thoroughly before use.
- **Sample buffer:** Ready for use.
- **Wash buffer:** The wash buffer is a 10x concentrate. If crystallisation occurs in the concentrated buffer, warm it to +37°C and mix well before diluting. The quantity required should be removed from the bottle using a clean pipette and diluted with deionised or distilled water (1 part reagent plus 9 parts distilled water).
For example: For 1 microplate strip, 5 ml concentrate plus 45 ml water.
The working strength wash buffer is stable for 4 weeks when stored at +2°C to +8°C and handled properly.
- **Chromogen/substrate solution:** Ready for use. Close the bottle immediately after use, as the contents are sensitive to light. The chromogen/substrate solution must be clear on use. Do not use the solution if it is blue coloured.
- **Stop solution:** Ready for use.

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Warning: The calibrators and controls of human origin have tested negative for HBsAg, anti-HCV, anti-HIV-1 and anti-HIV-2. Nonetheless, all materials should be treated as being a potential infection hazard and should be handled with care. Some of the reagents contain the agent sodium azide in a non-declarable concentration. Avoid skin contact.

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Sample dilution: Patient samples are diluted **1:51** in sample buffer. For example: dilute 10 µl sample in 0.5 ml sample buffer and mix well by vortexing (sample pipettes are not suitable for mixing).

NOTE: The calibrators and controls are prediluted and ready for use after reconstitution, do not dilute them.

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(Partly) manual test performance

Sample incubation: (1st step)

Transfer 100 µl of the calibrators, positive and negative controls or diluted patient samples into the individual microplate wells according to the pipetting protocol. Cover the wells and incubate for **30 minutes** at room temperature (+18°C to +25°C) on a **microplate shaker (600 rpm**, horizontal rotary shaker with a shaking amplitude of approx. 3 mm, e.g. microplate shaker TPM 4 from Sarstedt or Compact Digital Microplate Shaker from Thermo Scientific).*

Washing:

Manual: Remove the plastic foil. Empty the wells and subsequently wash 3 times using 350 µl of working strength wash buffer for each wash.

Automatic: Remove the plastic foil and wash the reagent wells 3 times with 450 µl of working strength wash buffer (program setting: e.g. TECAN Columbus Washer "Overflow Mode").

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Note: Residual liquid (> 10 µl) in the reagent wells after washing can interfere with the substrate and lead to false low extinction values. Insufficient washing (e.g. less than 3 wash cycles, too small wash buffer volumes, or too short residence times) can lead to false high extinction values.

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Conjugate incubation: (2nd step)

Pipette 100 µl of enzyme conjugate (peroxidase-labelled anti-human IgG) into each of the microplate wells. Incubate for **30 minutes** at room temperature (+18°C to +25°C) on a **microplate shaker (600 rpm**, horizontal rotary shaker with a shaking amplitude of approx. 3 mm, e.g. microplate shaker TPM 4 from Sarstedt or Compact Digital Microplate Shaker from Thermo Scientific).

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Substrate incubation: (3rd step)

Pipette 100 µl of chromogen/substrate solution into each of the microplate wells. Incubate for **15 minutes** at room temperature (+18°C to +25°C) (protect from direct sunlight).

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Measurement:

Photometric measurement of the colour intensity should be made at a **wavelength of 450 nm** and a reference wavelength between 620 nm and 650 nm **within 30 minutes of adding the stop solution**. Prior to measuring, slightly shake the microplate to ensure a homogeneous distribution of the solution.

* Please note that in case of deviating device specifications customers are required to carry out their own validation.

Test performance using fully automated analysis devices

Sample dilution and test performance are carried out fully automatically using an analysis device. The incubation conditions programmed in the respective software authorised by EUROIMMUN may deviate slightly from the specifications given in the ELISA test instruction. However, these conditions were validated in respect of the combination of the DSX from Dynex and this EUROIMMUN ELISA. Validation documents are available on enquiry.

Automated test performance using other fully automated, open-system analysis devices is possible. However, the combination should be validated by the user.

Pipetting protocol

	1	2	3	4	5	6	7	8	9	10	11	12
A	C 2	P 6	P 14	P 22			C 1	P 4	P 12	P 20		
B	pos.	P 7	P 15	P 23			C 2	P 5	P 13	P 21		
C	neg.	P 8	P 16	P 24			C 3	P 6	P 14	P 22		
D	P 1	P 9	P 17				pos.	P 7	P 15	P 23		
E	P 2	P 10	P 18				neg.	P 8	P 16	P 24		
F	P 3	P 11	P 19				P 1	P 9	P 17			
G	P 4	P 12	P 20				P 2	P 10	P 18			
H	P 5	P 13	P 21				P 3	P 11	P 19			

The pipetting protocol for microtiter strips 1 to 4 is an example for the **semiquantitative analysis** of 24 patient samples (P 1 to P 24).

The pipetting protocol for microtiter strips 7 to 10 is an example for the **quantitative analysis** of 24 patient samples (P 1 to P 24).

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Calculation of results


Semiquantitative: Results can be evaluated semiquantitatively by calculating a ratio of the extinction value of the control or patient sample over the extinction value of calibrator 2. Calculate the ratio according to the following formula:

$$\frac{\text{Extinction of the control or patient sample}}{\text{Extinction of calibrator 2}} = \text{Ratio}$$

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Quantitative: The calibration curve from which the concentration of antibodies in the patient samples can be taken is obtained by point-to-point plotting of the extinction values measured for the 3 calibrators against the corresponding units (linear/linear). Use "point-to-point" plotting for calculation of the calibration curve by computer. The following plot is an example of a typical calibration curve. Please do not use this curve for the determination of circulating immune complex concentrations in patient samples.

If the extinction for a patient sample lies above the value of calibrator 1 (200 RU/ml), the result should be reported as ">200 RU/ml". It is recommended that the sample be re-tested at a dilution of 1:201. The result in RU/ml read from the calibration curve for this sample must then be multiplied by a factor of 4.

The upper limit of the normal range (cut-off value) recommended by EUROIMMUN is 20 relative units per millilitre (RU/ml). This corresponds to 4 µg Eq/ml (µg equivalents of heat-aggregated human IgG per ml). EUROIMMUN recommends interpreting results as follows:

<20 RU/ml:	negative
≥20 RU/ml:	positive

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Test characteristics

Calibration: The calibration is performed in relative units (RU). Five RU/ml correspond to one µg Eq/ml (µg equivalents of heat-aggregated human IgG per ml).

For every group of tests performed, the extinction values of the calibrators and the relative units and/or ratio determined for the positive and negative controls must lie within the limits stated for the relevant test kit lot. A quality control certificate containing these reference values is included. If the values specified for the controls are not achieved, the test results may be inaccurate and the test should be repeated.

The binding activity of the antibodies and the activity of the enzyme used are temperature-dependent. It is therefore recommended using a thermostat in all three incubation steps. The higher the room temperature (+18°C to +25°C) during the incubation steps, the greater will be the extinction values. Corresponding variations apply also to the incubation times. However, the calibrators, controls and patient samples are subject to the same influences, with the result that such variations will be largely compensated in the calculation of the result.

Antigen: The microplate wells were coated with the complement factor C1q prepared from human serum. The circulating immune complexes bind to the C1q via the Fc fragment of the immunoglobulins. C1q is a glycoprotein composed of 18 polypeptide chains and consisting of three nonidentical subunits with M.W. of 29,000, 26,000 and 19,000 respectively. C1q is complexed with two C1r and two C1s molecules to form the first component of complement (C1). Activation of complement via the classical pathway is triggered by binding of C1q to immune complexes or to a variety of other activating substances. Subsequent to C1q binding, C1r and C1s are converted to proteolytic enzymes, which are responsible for continuation of activation via the classical pathway.

Linearity: The linearity of the CIC-C1q ELISA (IgG) was determined by assaying at least 4 serial dilutions of different patient samples. The CIC-C1q ELISA (IgG) is linear at least in the tested concentration range (8 RU/ml to 193 RU/ml).

Detection limit: The lower detection limit is defined as the mean value of an analyte-free sample plus three times the standard deviation and is the smallest detectable antibody titer. The lower detection limit of the CIC-C1q ELISA (IgG) is 0.1 RU/ml

Analytical specificity: The ELISA presented here specifically detects circulating immune complexes containing IgG with the ability to bind to human complement fragment C1q. Elevated concentrations of CIC-C1q are observed in up to 10% of blood donors without any clinical manifestations.

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Sensitivity and specificity: The CIC-C1q levels were determined in 55 serum samples from patients with systemic lupus erythematosus (SLE) or rheumatoid arthritis (RA) and from apparently healthy blood donors (origin: Europe, China) using the EUROIMMUN ELISA and another commercially available ELISA (Bühlmann Laboratories) as reference. Of the analysed sera, the same 22 were positive and the same 30 were negative in both assays. Two sera were found positive by the EUROIMMUN ELISA only, and one serum by the reference ELISA (Bühlmann Laboratories) only. The sensitivity amounted to 96%, with a specificity of 94%.

Reference range: The CIC-C1q concentrations were determined in serum samples from a total of 250 apparently healthy blood donors (origin: Germany) using this EUROIMMUN ELISA. CIC-C1q was detected at concentrations ranging from 1 to 296 RU/ml. The 90.0 percentile of the reference interval was calculated to be 15 RU/ml and the 95.0 percentile was found to be 31 RU/ml.

Clinical significance

Immune complexes are protein aggregates that form via the binding of antibodies to antigens. They are continually formed as a result of the immune response. Since they are rapidly degraded, circulating immune complexes (synonym: circulating Ab-Ag complexes, CIC) are not normally detected in serum or only in low concentrations. If CIC are present in disease in such high concentrations that they cannot be quickly eliminated from the body, e.g. if phagocyte capacity is overwhelmed, accumulation of the immune complexes occurs. This can influence the coagulation system and lead to increased thrombocyte and erythrocyte aggregation. Through activation of the complement system CIC can induce inflammatory reactions and cause damage to organs, as is the case, for example, in glomerulonephritis, Crohn's disease, ulcerative colitis, rheumatic diseases and viral or bacterial infections.

In the classical activation pathway the antigen-antibody complexes bind to the complement components C1, C4 and C2 and trigger the enzyme cascade. This is continued and finished with the components C3, C5, C6, C7, C8 and C9. The factors, which include C1q, act as indicators of immune complex diseases with complement use.

Circulating or extravascular immune complexes can become pathogenetically active. When immune complexes accumulate on or in tissues this can lead to inflammatory processes with corresponding clinical symptoms. The involvement of immune complexes can be detected in a range of immune complex-associated autoimmune diseases, such as systemic lupus erythematosus (CLE), rheumatoid arthritis (RA) and immune complex glomerulonephritis. Elevated CIC serum concentrations are also found in bacterial, viral and parasitic infections as post- and parainfectious immune complex diseases, as well as in allergies, chronic skin diseases and neurological diseases.

Although the detection of CIC is not specific for a particular disease these analyses can yield useful information about the immunopathology, course and prognosis of disease. The presence of CIC is an indicator of an overwhelmed immune defence or an autoimmune conflict, which should be clarified by targeted tests. In autoimmune, microbial or neoplastic diseases measurement of CIC is used as a marker for disease activity and also for assessing organ manifestation and for monitoring treatment. The test method of choice is ELISA, which allows quantitative in vitro determination of C1q-binding, IgG-containing CIC. Complement CIC are bound to immobilised purified C1q proteins. Persistently elevated CIC-C1q indicates chronic active disease. Normalisation is considered an indicator of successful treatment. It is often possible to correlate the CIC level with clinical symptoms and relapses. The WHO recommends always performing a second determination several weeks later in order to evaluate the persistence of the immune complexes in the circulation and thereby the pathological significance. The second determination also allows treatment to be monitored. With appropriate therapy a clear decrease in CIC should occur.

According to the WHO the following diseases are CIC-associated: rheumatic diseases (rheumatoid arthritis, systemic lupus erythematosus, Sjögren's syndrome, mixed connective tissue disease, periarteritis nodosa, Felty's syndrome, Bechterew's disease, Reiter's disease), viral infections (hepatitis B, cytomegaly, mononucleosis, subacute and sclerosing panencephalitis), bacterial infections (infectious endocarditis, disseminated gonorrhoea, syphilis, streptococcus and meningococcus infections), parasitosis (malaria, schistosomiasis, trypanosomiasis, toxoplasmosis), glomerulonephritis, neoplastic diseases and other diseases such as Crohn's disease, ulcerative colitis, idiopathic interstitial pneumonia, cystic fibrosis, multiple sclerosis, thrombotic thrombocytopenic purpura, Behcet's syndrome Churg-Strauss syndrome and chronic liver diseases.

Literature references

1. Aziz M, Akhtar S, Malik A. **Evaluation of cell-mediated immunity and circulating immune complexes as prognostic indicators in cancer patients.** Cancer detection and prevention 22.2 (1997): 87-99.
2. Braun A, Sis J, Max R, Mueller K, Fiehn C, Zeier M, Andrassy K. **Anti-chromatin and anti-C1q antibodies in systemic lupus erythematosus compared to other systemic autoimmune diseases.** Scand J Rheumatol 36 (2007) 291-298.
3. EUROIMMUN AG. Stöcker W, Schlumberger W, Krüger C. **Alle Beiträge zum Thema Autoimmundiagnostik.** In: Gressner A, Arndt T (Hrsg.) Lexikon der Medizinischen Laboratoriumsdiagnostik. 2. Auflage. Springer Medizin Verlag, Heidelberg (2012).
4. Flierman R, Daha MR. **Pathogenic role of anti-C1q autoantibodies in the development of lupus nephritis--a hypothesis.** Mol Immunol 44 (2007) 133-138.
5. Oelzner P, Deliyska B, Fünfstück R, Hein G, Herrmann D, Stein G. **Anti-C1q antibodies and antiendothelial cell antibodies in systemic lupus erythematosus - relationship with disease activity and renal involvement.** Clin Rheumatol 22 (2003) 271-278.
6. O'Flynn J, Flierman R, van der Pol P, Rops A, Satchell SC, Mathieson PW, van Kooten C, van der Vlag J, Berden JH, Daha MR. **Nucleosomes and C1q bound to glomerular endothelial cells serve as targets for autoantibodies and determine complement activation.** Mol Immunol 49 (2011) 75-83.

7. Saadoun, David, et al. **Anti-C1q antibodies in hepatitis C virus infection.** Clinical & Experimental Immunology 145.2 (2006): 308-312.
8. Shoenfeld Y, Szyper-Kravitz M, Witte T, Doria A, Tsutsumi A, Tatsuya A, Dayer JM, Roux-Lombard P, Fontao L, Kallenberg CG, Bijl M, Matthias T, Fraser A, Zandman-Goddard G, Blank M, Gilburd B, Meroni PL. **Autoantibodies against protective molecules--C1q, C-reactive protein, serum amyloid P, mannose-binding lectin, and apolipoprotein A1: prevalence in systemic lupus erythematosus.** Ann N Y Acad Sci 1108 (2007) 227-239.
9. Smykał-Jankowiak K, Niemir ZI, Polcyn-Adamczak M. **Do circulating antibodies against C1q reflect the activity of lupus nephritis?** Pol Arch Med Wewn 121 (2011) 287-295.
10. Trouw LA, Groeneveld TW, Seelen MA, Duijs JM, Bajema IM, Prins FA, Kishore U, Salant DJ, Verbeek JS, van Kooten C, Daha MR. **Anti-C1q autoantibodies deposit in glomeruli but are only pathogenic in combination with glomerular C1q-containing immune complexes.** J Clin Invest 114 (2004) 679-688.
11. Walport MJ. **Complement and systemic lupus erythematosus.** Arthritis Res 4 (2002) 279-293.
12. Zabaleta-Lanz ME, Muñoz LE, Tapanes FJ, Vargas-Arenas RE, Daboin I, Barrios Y, Pinto JA, Bianco NE. **Further description of early clinically silent lupus nephritis.** Lupus 15 (2006) 845-851.

EE Broth-Mossel

Liquid medium for the cultivation and selective enrichment of Enterobacteriaceae from different types of samples, according to USP/EP/JP.

DESCRIPTION

Enterobacteriaceae Enrichment Broth-Mossel is a selective medium used for the detection of bile-tolerant Gram-negative bacteria in food and other materials of sanitary importance.

This medium complies with the recommendations of the harmonized method in the United States Pharmacopoeia (USP), European Pharmacopoeia (EP) and Japanese Pharmacopoeia (JP) for the microbiological examination of nonsterile products.

TYPICAL FORMULA	(g/l)
Pancreatic Digest of Gelatin	10.0
Glucose Monohydrate	5.0
Dehydrated Ox Bile	20.0
Potassium Dihydrogen Phosphate	2.0
Disodium Hydrogen Phosphate, Anhydrous	6.4*
Brilliant Green	0.015
Final pH 7.2 ± 0.2 at 25°C	

* Equivalent to 8.0 g of Disodium Hydrogen Phosphate Dihydrate.

METHOD PRINCIPLE

Pancreatic digest of gelatin provides amino acids, nitrogen, carbon, vitamins and minerals for organisms growth. Glucose is the fermentable carbohydrate. Ox bile and brilliant green are selective agents effective against Gram-positive cocci. Potassium phosphate and sodium phosphate act as buffer.

PREPARATION

Dehydrated medium Suspend 43.4 g of the powder in 1 liter of distilled or deionized water. Mix well. Heat to boil shaking frequently until completely dissolved. DO NOT AUTOCLAVE.

TEST PROCEDURE

As in the Pharmacopoeia, prepare the sample using a 1 in 10 dilution of not less than 1 g of the product to be examined by choosing Tryptic Soy Broth (ref. 24513 or 452080) as diluent and incubate at 20-25°C for 2-5 hour to resuscitate bacteria.

For qualitative test (test for absence), transfer the volume of the pre-enrichment broth corresponding to 1 g of the product to be examined to EE Broth-Mossel.

For quantitative test, enumerate Enterobacteriaceae found per milliliter or per gram of test sample by using the Most Probable Number (MPN) technique. Use the volume of the pre-enrichment broth containing 0.1 g, 0.01 g and 0.001 g (or 0.1 ml, 0.01 ml and 0.001 ml) of the product to be examined to inoculate EE Broth-Mossel.

For both types of test incubate EE Broth-Mossel at 30-35°C for 24-48 h and continue analysis by subculturing on Violet Red Bile Glucose Agar (ref. 11184). Incubate plates aerobically at 30-35°C for 18-24 hours.

INTERPRETING RESULTS

Turbidity of EE Broth-Mossel indicates microbial growth; acid production causes a color change of the medium to yellow.

No growth of colonies on Violet Red Bile Glucose Agar is reported as absence of bile-tolerant Gram-negative bacteria.

Growth of colonies constitutes a positive result and the probable number of bacteria is determined from the table below.

MPN Table.

Results for each quantity of product			Probable number of bacteria per gram or per milliliter of product
0.1 g or 0.1 ml	0.01 g or 0.01 ml	0.001 g or 0.001 ml	
+	+	+	>10 ³
+	+	—	10 ³ - 10 ²
+	—	—	10 ² - 10
—	—	—	<10

APPEARANCE

Dehydrated medium: free-flowing, homogeneous, light beige to light green.

Prepared medium: clear, green.

STORAGE

The powder is very hygroscopic, store the powder at 10-30°C, in a dry environment, in its original container tightly closed. Store bottles and prepared plates at 10-25°C away from light. Do not use the product beyond its expiry date on the label or if product shows any evidence of contamination or any sign of deterioration.

SHELF LIFE

Dehydrated medium: 4 years.

Medium in tubes/bottles: 1 year.

QUALITY CONTROL

The medium is inoculated with the microbial strains indicated in the QC table.

Inoculum for productivity: ≤100 CFU.

Inoculum for selectivity: >100 CFU.

Incubation conditions: 18-24 h at 30-35°C (Pharmacopoeia growth promotion).

QC Table.

Microorganism		Specification
<i>Escherichia coli</i>	ATCC® 8739	Good growth
<i>Pseudomonas aeruginosa</i>	ATCC® 9027	Good growth
<i>Staphylococcus aureus</i>	ATCC® 6538	Inhibition

WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is intended for professional use only and must be used by properly trained operators.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAPHY

1. EN ISO 11133:2014. Microbiology of food, animal feed and water – Preparation, production, storage and performance testing of culture media.
2. European Pharmacopoeia 6.5 (2009) 2.6.13. Microbiological examination of non-sterile products: Test for specified microorganisms.
3. United States Pharmacopoeia 32 NF 27 (2009) <62> Microbiological examination of non-sterile products: Test for specified microorganisms.
4. Japanese Pharmacopoeia 4.05 (2008) Microbiological examination of non-sterile products: Test for specified microorganisms.
5. ISO 21528-1:2004. Microbiology of food and animal feeding stuffs – Horizontal method for the detection and enumeration of Enterobacteriaceae – Detection and enumeration by MPN technique with pre-enrichment.
6. ISO 21528-2:2004. Microbiology of food and animal feeding stuffs – Horizontal method for the detection and enumeration of Enterobacteriaceae – Colony count method.
7. Davidson, Roth, and Gambrel-Lenarz (2004) In Wehr and Frank (ed.) Standard methods for the microbiological examination of dairy products, 17th ed. American Public Health Association, Washington, D.C.
8. Kornacki and Johnson (2001) In Downes and Ito (ed.) Compendium of methods for the microbiological examination of foods, 4th ed. American Public Health Association, Washington D.C.
9. Mossel, Vissar, and Cornellisen (1963) J. Appl. Bacteriol. 26:444.

PRESENTATION		Contents	Ref.
EE Broth-Mossel	Tubes	20 x 10 ml tubes	24096
EE Broth-Mossel	Bottles (screw cap)	6 x 100 ml bottles	402480
EE Broth-Mossel	Bottles (flip-off cap)	25 x 100 ml bottles	453080
EE Broth-Mossel	Bottles (perforable cap)	6 x 100 ml bottles	495000
EE Broth-Mossel	Dehydrated medium	500 g of powder	610017
EE Broth-Mossel	Dehydrated medium	100 g of powder	620017

TABLE OF SYMBOLS

LOT Batch code	 Keep away from sunlight	 Manufacturer	 Use by	 Fragile, handle with care
REF Catalogue number	 Temperature limitation	 Contains sufficient for <n> tests	 Caution, consult Instruction For Use	 Do not reuse

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IVD

INTEGRAL SYSTEM ENTEROBATTERI

Colorimetric system for biochemical identification and susceptibility testing of enterobacteria

DESCRIPTION

INTEGRAL SYSTEM ENTEROBATTERI is a 24-well system containing desiccated biochemical and antibiotic substrates for the biochemical identification and susceptibility testing of enterobacteria.

The system is inoculated with the bacterial suspension of the micro-organism under examination and incubated at $36\text{ }^{\circ}\text{C} \pm 1^{\circ}\text{C}$ for 18-24 hours. The tests for identification and for the susceptibility testing are interpreted by assessing the change in colour of the various wells.

KIT CONTENTS

The kit contains:

- 20 systems of INTEGRAL SYSTEM ENTEROBATTERI
- 20 Vials of Suspension Medium (7.0 ml/vial)
- 20 Oxidase Disc
- 1 Instruction sheet

ITEMS NECESSARY BUT NOT INCLUDED IN THE KIT

- Vaseline oil for microbiology use (Vaselin oil 2 bottles of 50 ml)cod. 80278
- Physiological Solution (20 vials of 7.0 ml/vial).....cod. 20095
- Kovac's Reagent (2 bottles of 25 ml).....cod. 80270
- Sundry microbiology laboratory materials

CONFIGURATION

The configuration of the system is shown in Table 1.

Table 1

Well	BIOCHEMICAL IDENTIFICATION
1-LDC	Decarboxylation of lysine
2-ODC	Decarboxylation of ornithine
3-H ₂ S	Production of hydrogen sulphide
4-IND	Indole test
5-LAC	Fermentation of lactose
6-DUL	Fermentation of dulcitol
7-PA	Deamination of phenylalanine
8-UR	Hydrolysis of urea
9-CIT	Utilisation of citrate
10-OX	Oxidase test
Well	SUSCEPTIBILITY TESTING (*)
11-AK	Amikacin - 32 µg/ml
12-CN	Gentamycin - 8 µg/ml
13-TOB	Tobramycin - 8 µg/ml
14-TZP	Piperacillin+Tazobactam - 128/4 µg/ml
15-FOS	Fosfomicin - 200 µg/ml
16-CFP	Cefoperazone - 64 µg/ml
17-CTX	Cefotaxime - 64 µg/ml
18-CAZ	Ceftazidime - 32 µg/ml
19-AMS	Ampicillin + Sulbactam - 32/16 µg/ml
20-NA	Nalidixic acid - 32 µg/ml
21-CIP	Ciprofloxacin - 4 µg/ml
22-LEV	Levofloxacin - 8 µg/ml
23-SXT	Co-Trimoxazole - 8 µg/ml
24-C	Growth control for the susceptibility testing

(*): The concentration of each antibiotic is according to the updated NCCLS-January 2004, Vol 24 N°1.⁽¹⁾

PRINCIPLE OF THE METHOD

INTEGRAL SYSTEM ENTEROBATTERI makes possible simultaneous biochemical identification and assessment of the susceptibility testing of enterobacteria isolated from clinical samples.

- The identification is based on biochemical tests performed on culture media containing specific substrates in the wells from **1-LDC** to **10-OX**.
- The susceptibility testing is assessed on the basis of the growth or inhibition of the micro-organisms in the media containing the antibiotic and a growth indicator in the wells from **11-AK** to **23-SXT**. Well **24-C** does not contain antibiotics but only the culture medium and the indicator; it serves as a microbe growth control for assessment of the susceptibility testing.

COMPOSITION

Table n°2

<i>Well</i>	<i>Content</i>
1-LDC	Culture medium for showing lysine decarboxylase
2-ODC	Culture medium for showing ornithine decarboxylase
3-H₂S	Culture medium for showing production of hydrogen sulphide
4-IND	Culture medium for showing production of indole
5-LAC	Culture medium for showing fermentation of lactose
6-DUL	Culture medium for showing fermentation of dulcitol
7-PA	Culture medium for showing phenylalanine deamination
8-UR	Culture medium for showing hydrolysis of urea
9-CIT	Culture medium for showing utilisation of citrate
10-OX	Culture medium for showing enzyme cytochrome-oxidase
11-AK	Culture medium containing Amikacin - 32 µg/ml
12-CN	Culture medium containing Gentamycin - 8 µg/ml
13-TOB	Culture medium containing Tobramycin - 8 µg/ml
14-TZP	Culture medium containing Piperacillin+Tazobactam - 128/4 µg/ml
15-FOS	Culture medium containing Fosfomycin - 200 µg/ml
16-CFP	Culture medium containing Cefoperazone - 64 µg/ml
17-CTX	Culture medium containing Cefotaxime - 64 µg/ml
18-CAZ	Culture medium containing Ceftazidime - 32 µg/ml
19-AMS	Culture medium containing Ampicillin + Sulbactam - 32/16 µg/ml
20-NA	Culture medium containing Nalidixic acid - 32 µg/ml
21-CIP	Culture medium containing Ciprofloxacin - 4 µg/ml
22-LEV	Culture medium containing Levofloxacin - 8 µg/ml
23-SXT	Culture medium containing Co-Trimoxazole - 8 µg/ml
24-C	Medium without antibiotics

Suspension Medium (g/l): Yeast extract **5 g**; Beef peptone **3 g**; Glucose **2 g**; Distilled water **1000.0 ml**;
pH 6.8±0.2

COLLECTION AND STORAGE OF THE SAMPLES

The colonies to be subjected to the biochemical identification test and susceptibility testing assessment with INTEGRAL SYSTEM ENTEROBATTERI must be taken from a culture medium – selective or non selective – used for the isolation of enterobacteria.

TEST PROCEDURE

1. Take a system from the kit.
2. After checking that the colonies that have developed on the culture media presumably belong to the enterobacteria group, take one or more morphologically similar, well isolated, colonies from a solid culture medium and suspend them in 5 ml of Physiological Solution for microbiological use in such a way as to obtain a turbidity equivalent to 0.5 McFarland (**Bacterial suspension**).
3. Trasfer
 - 0.2 ml of **Bacterial suspension** into each of the first 10 wells and to wells **1-LDC**, **2-ODC**, **3-H₂S** and **8-UR** add 2 drops of Vaseline Oil for microbiological use (BIOCHEMICAL IDENTIFICATION).
 - 0.01 ml of **Bacterial suspension** to the vial of **Suspension Medium** contained into the kit and distribute 0.2 ml into the wells from **11-AK** to **24-C** (SUSCEPTIBILITY TESTING).
4. Cover the system with the lid provided and incubate at 36 °C ± 1°C for 18-24 hours.
5. After incubation. add:
 - 2 drops of Kovac's Reagent to well **4-IND** (INDOLE TEST).
 - 1 OXIDASE DISC to well **10-OX** (OXIDASE TEST).

INTERPRETATION OF THE RESULTS

BIOCHEMICAL IDENTIFICATION

Interpret the results of the first 10 wells using Table n° 3 and form the 4-digit code following the instructions given to the paragraph **FORMATION OF NUMERICAL CODE**. Then use the code and Table n°5 to identify the bacteria.

SUSCEPTIBILITY TESTING

Observe the changes of colour of the wells from **11-AK** to **23-SXT** and interpret the results using Table n°3. The control well (**24-C**) must be positive (yellow). If it is negative (blue or grey), it is necessary to check the vitality of the inoculum and whether it has been correctly prepared; then to repeat the tests using a new system.

Table n°3

Well	BIOCHEMICAL IDENTIFICATION	Well colour	
		Positive reaction	Negative reaction
1-LDC	Decarboxylation of lysine	Violet	Yellow-Brown
2-ODC	Decarboxylation of ornithine	Violet	Yellow-Brown
3 -H₂S	Production of hydrogen sulphide	Black	Yellow
4-IND	Indole test	Pink-Red ring	Yellow
5-LAC	Fermentation of lactose	Yellow	Blue-Green
6-DUL	Fermentation of dulcitol	Yellow	Blue-Green
7-PA	Deamination of phenylalanine	Black-Brown	Yellow
8-UR	Hydrolysis of urea	Red-Fuchsia	Yellow-Orange
9-CIT	Utilisation of citrate	Blue-Dark Green	Green
10-OX	Oxidase test	Blue-Purple (immediate reaction)	Colourless

SUSCEPTIBILITY TESTING		
WELL COLOUR	BACTERIAL GROWTH	INTERPRETATION
Blue	Inhibited	S = Sensitive
Grey	Intermediate	I = Intermediate sensitivity
Yellow	Good	R = Resistant

FORMATION OF THE NUMERICAL CODE

1) The 10 biochemical tests are divided into 3 groups each containing 3 tests and one group containing 1 test, each one is indicated with a positiveness value of 1,2,4.

- Value 1 : first test positive in each group
- (LDC, IND, PA,OX)
- Value 2 : second test positive in each group
- (ODC,LAC,UR)
- Value 4 : third test positive in each group
- (H₂S, DUL,CIT)
- Value zero: negative reaction in every group

2) Adding the number of positive reactions in each group, one obtains a 4-digit code which is used to identify the micro-organism from the table of numerical codes as in the example.

	Group I			Group II			Group III			Group IV
Well	LDC	ODC	H ₂ S	IND	LAC	DUL	PA	UR	CIT	OX
Positiveness code	1	2	4	1	2	4	1	2	4	1
Results	+	+	-	+	+	-	-	-	-	-
Codes addition	3			3			0			0
NUMERICAL CODE: 3300				IDENTIFICATION: <i>E.coli</i>						

TABLE OF BIOCHEMICAL REACTIONS

Table n°4

Micro-organisms	Group I			Group II			Group III			Group IV
	LDC	ODC	H ₂ S	IND	LAC	DUL	PA	UR	CIT	OX
<i>E.coli</i>	V	V	-	+	+	V	-	-	-	-
<i>E.coli inactive</i>	V	V	-	+	-	V	-	-	-	-
<i>Shigella spp.</i>	-	V	-	V	-	-	-	-	-	-
<i>Edwardsiella spp.</i>	+	+	+	+	-	-	-	-	-	-
<i>Citrobacter spp.</i>	-	V	+	V	V	V	-	V	V	-
<i>Salmonella spp.</i>	+	+	V	-	-	+	-	-	V	-
<i>Arizona spp.</i>	+	+	+	V	V	-	-	-	+	-
<i>K. pneumoniae</i>	+	-	-	-	+	V	-	V	V	-
<i>Klebsiella oxytoca</i>	+	-	-	+	+	V	-	+	V	-
<i>Enterobacter cloacae</i>	V	+	-	-	V	-	-	V	+	-
<i>E. aerogenes</i>	+	+	-	-	+	-	-	-	+	-
<i>Enterobacter hafnia</i>	+	+	-	-	V	-	-	-	+	-
<i>Serratia spp.</i>	+	+	-	-	-	-	-	V	+	-
<i>Proteus mirabilis</i>	-	+	+	-	-	-	+	+	V	-
<i>Proteus vulgaris</i>	-	-	V	+	-	-	+	+	V	-
<i>Morganella morganii</i>	-	+	-	+	-	-	+	+	-	-
<i>Providencia stuartii</i>	-	-	-	+	-	-	+	-	+	-
<i>Providencia rettgeri</i>	-	-	-	+	-	-	+	+	+	-
<i>Pseudomonas spp.</i>	+	+	-	-	-	-	-	-	+	+

+ : Positive Reaction

- : Negative Reaction

V: Variable Reaction

*Dulcitol: Some strains of *Salmonella* ferment dulcitol slowly; they may therefore be identified initially as *Arizona spp.* A confirmation test is thus required using a polyvalent antiserum for *Salmonella*.

TABLE OF NUMERICAL CODES

Table n°5

Code	Micro-organism	Code	Micro-organism	Code	Micro-organism	Code	Micro-organism
0000	<i>Shigella</i> spp.	2240	<i>Enterobacter cloacae</i>	4220	<i>Citrobacter</i> spp.	6160	<i>Citrobacter</i> spp.
0100	<i>E.coli</i> inactive 45%	2260	<i>Enterobacter cloacae</i>	4240	<i>Citrobacter</i> spp.	6200	<i>Citrobacter</i> spp.
	<i>Shigella</i> spp. 55%	2300	<i>E.coli</i>	4260	<i>Citrobacter</i> spp.	6220	<i>Citrobacter</i> spp.
0130	<i>Proteus vulgaris</i>	2500	<i>E.coli</i> inactive	4300	<i>Citrobacter</i> spp.	6240	<i>Citrobacter</i> spp.
0150	<i>Providencia stuartii</i>	2700	<i>E.coli</i>	4320	<i>Citrobacter</i> spp.	6260	<i>Citrobacter</i> spp.
0170	<i>Proteus vulgaris</i> 0.5%	3040	<i>Enterobacter hafniae</i> 60%	4340	<i>Citrobacter</i> spp.	6300	<i>Citrobacter</i> spp.
	<i>Providencia rettgeri</i> 99.5%		<i>Serratia</i> spp. 37%	4360	<i>Citrobacter freundii</i>	6320	<i>Citrobacter</i> spp.
0300	<i>E.coli</i>		<i>Enterobacter cloacae</i> 3%	4400	<i>Citrobacter</i> spp.	6340	<i>Citrobacter</i> spp.
0500	<i>E.coli</i> inactive	3041	<i>Pseudomonas</i> spp.	4420	<i>Citrobacter</i> spp.	6360	<i>Citrobacter</i> spp.
0700	<i>E.coli</i>	3060	<i>Serratia</i> spp. 98.9 %	4440	<i>Citrobacter</i> spp.	6400	<i>Citrobacter</i> spp.
1100	<i>E.coli</i> inactive		<i>Enterobacter cloacae</i> 1.1%	4460	<i>Citrobacter</i> spp.	6420	<i>Citrobacter</i> spp.
1200	<i>Klebsiella pneumoniae</i>	3100	<i>E.coli</i> inactive	4500	<i>Citrobacter</i> spp.	6440	<i>Citrobacter</i> spp.
1220	<i>Klebsiella pneumoniae</i>	3240	<i>Enterobacter cloacae</i> 7.7%	4520	<i>Citrobacter</i> spp.	6460	<i>Citrobacter</i> spp.
1240	<i>Klebsiella pneumoniae</i>		<i>Enterobacter aerogenes</i> 75%	4540	<i>Citrobacter</i> spp.	6500	<i>Citrobacter</i> spp.
1260	<i>Klebsiella pneumoniae</i>		<i>Enterobacter hafniae</i> 17.3%	4560	<i>Citrobacter</i> spp.	6520	<i>Citrobacter</i> spp.
1300	<i>E.coli</i>	3260	<i>Enterobacter cloacae</i>	4600	<i>Citrobacter</i> spp.	6540	<i>Citrobacter</i> spp.
1320	<i>Klebsiella oxytoca</i>	3300	<i>E.coli</i>	4620	<i>Citrobacter</i> spp.	6560	<i>Citrobacter</i> spp.
1360	<i>Klebsiella oxytoca</i>	3400	<i>Salmonella</i> spp.	4640	<i>Citrobacter</i> spp.	6600	<i>Citrobacter</i> spp.
1500	<i>E.coli</i> inactive	3440	<i>Salmonella</i> spp.	4660	<i>Citrobacter</i> spp.	6620	<i>Citrobacter</i> spp.
1600	<i>Klebsiella pneumoniae</i>	3500	<i>E.coli</i> inactive	4700	<i>Citrobacter</i> spp.	6640	<i>Citrobacter</i> spp.
1620	<i>Klebsiella pneumoniae</i>	3700	<i>E.coli</i>	4720	<i>Citrobacter</i> spp.	6660	<i>Citrobacter</i> spp.
1640	<i>Klebsiella pneumoniae</i>	4000	<i>Citrobacter</i> spp.	4740	<i>Citrobacter</i> spp.	6700	<i>Citrobacter</i> spp.
1660	<i>Klebsiella pneumoniae</i>	4020	<i>Citrobacter</i> spp.	4760	<i>Citrobacter</i> spp.	6720	<i>Citrobacter</i> spp.
1700	<i>E.coli</i>	4040	<i>Citrobacter</i> spp.	6000	<i>Citrobacter</i> spp.	6740	<i>Citrobacter</i> spp.
1720	<i>Klebsiella oxytoca</i>	4060	<i>Citrobacter</i> spp.	6020	<i>Citrobacter</i> spp.	6760	<i>Citrobacter</i> spp.
1760	<i>Klebsiella oxytoca</i>	4100	<i>Citrobacter</i> spp.	6030	<i>Proteus mirabilis</i>	7040	<i>Arizona</i> spp.
2000	<i>Shigella</i> spp.	4120	<i>Citrobacter</i> spp.	6040	<i>Citrobacter</i> spp.	7100	<i>Edwardsiella</i> spp.
2040	<i>Enterobacter cloacae</i>	4130	<i>Proteus vulgaris</i>	6060	<i>Citrobacter</i> spp.	7140	<i>Arizona</i> spp.
2060	<i>Enterobacter cloacae</i>	4140	<i>Citrobacter</i> spp.	6070	<i>Proteus mirabilis</i>	7240	<i>Arizona</i> spp.
2100	<i>E.coli</i> inactive 44.6%	4160	<i>Citrobacter</i> spp.	6100	<i>Citrobacter</i> spp.	7340	<i>Arizona</i> spp.
	<i>Shigella</i> spp. 55.4%	4170	<i>Proteus vulgaris</i>	6120	<i>Citrobacter</i> spp.	7400	<i>Salmonella</i> spp.
2130	<i>Morganella morganii</i>	4200	<i>Citrobacter</i> spp.	6140	<i>Citrobacter</i> spp.	7440	<i>Salmonella</i> spp.

QUALITY CONTROL

Each batch of INTEGRAL SYSTEM ENTEROBATTERI is subjected to quality control using the following ATCC reference strains.

<i>Enterobacter cloacae</i>	ATCC 13047	<i>Proteus mirabilis</i>	ATCC 25933
<i>Escherichia coli</i>	ATCC 25922	<i>Pseudomonas aeruginosa</i>	ATCC 27853
<i>Klebsiella pneumoniae</i>	ATCC 13883	<i>Salmonella typhimurium</i>	ATCC 14028

FACTORS THAT MAY INVALIDATE THE RESULTS

Poor standardisation of the inoculum; application of the method to micro-organisms not in the Enterobacteria group; mixed or contaminated cultures; use of expired systems or expired supplementary reagents; incorrect application of the technique.

LIMITS AND WARNINGS

The identification of *Salmonella spp.* and of *Shigella spp.* must be confirmed using appropriate antisera. In some cases, supplementary confirmation tests may be required for definitive identification of the micro-organisms.

PERFORMANCE

The bacterial identification results obtained with the INTEGRAL SYSTEM ENTEROBATTERI agree with those obtained using traditional methods as indicated by Piccolomini et al.⁽⁴⁾.

The antibiogram results obtained with the INTEGRAL SYSTEM ENTEROBATTERI agree with those obtained using the radial diffusion methods according to Bauer et al.⁽⁵⁾ recommended by the Food and Drug Administration (FDA) and by the National Committee for Clinical Laboratory Standards (NCCLS) U.S.A.^(1,8).

PRECAUTIONS

The product, INTEGRAL SYSTEM ENTEROBATTERI, cannot be classified as hazardous under current legislation nor does it contain harmful substances in concentrations $\geq 1\%$. It therefore does not require a Safety Data Sheet to be available. INTEGRAL SYSTEM ENTEROBATTERI is a disposable device to be used only for diagnostic use *in vitro*; it is intended for use in a professional environment and must be used in the laboratory by properly trained personnel, using approved asepsis and safety methods for handling pathogenic agents.

STORAGE

Store at 2-8°C in the original packaging. Keep away from sources of heat and avoid excessive changes in temperature. In such conditions, the product will remain valid until the expiry date indicated on the label. Do not use beyond that date. Eliminate without using if there are signs of deterioration.

DISPOSAL OF USED MATERIAL

After use, INTEGRAL SYSTEM ENTEROBATTERI and material that has come into contact with the sample must be decontaminated and disposed of in accordance with the techniques used in the laboratory for decontamination and disposal of potentially infected material.

PRESENTATION

Product	Code	Kit
INTEGRAL SYSTEM ENTEROBATTERI	71714	20 tests

TABLE OF SYMBOLS

SYMBOL	MEANING	SYMBOL	MEANING
	In Vitro Diagnostic Medical Device		Do not reuse
	Manufacturer		Contains sufficient for <n> tests
	Catalogue number		Fragile, handle with care
	Use by		Caution, consult accompanying documents
	Temperature limitation		Batch code

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- 1.
2. NCCLS – Performance Standards for antimicrobial Susceptibility Testing, eleventh Informational Supplement- January 2004, Vol 24 N°1.
3. Ewing H.E., 1986 Edwards and Ewings “Identification of Enterobacteriaceae”, 4th ed. Burgess Publishing co.Minneapolis
4. Kelly M.T., D.J. Brenne and J.J. Farmer,1985. “Enterobacteriaceae”, p.263-277 E.H.Lennette, A.Balows, W.J. Hauster,Jr.and H.J. Shadomy (ed.) Manual of Clinical Microbiology, 4th ed. American Society for Microbiology Washington, D.C.
5. Piccolomini R., Di Girolamo A., Catamo G., Cellini L., Allocati N. and Ravagnan G.’ Enterosistem 18-R: description and Comparative Evaluation with Conventional Methods for Identification of Members of the Family *Enterobacteriaceae*. 1991 *J.Clin.Microbiol.* **29**(10), 2300-2304.
6. Baker C. et al. Inoculum Standardization in Antimicrobial susceptibility testing:evaluation of overnight agar cultures and the rapid inoculum standardization system”. 1983 *J.Clin.Microbiol.* **17**, 450-457.
7. Lampe M/F et al. Relationship of early readings of minimal inhibitory concentrations to the results of overnight tests. 1975. *Ant. Agents Chernoth.* **8**, 429-433.
8. Bauer, A.W., W.M.M. Kirby, J.C. Sherris and M. Turk. 1966. Antimicrobial susceptibility testing by a standardized single disk method. *AM. J. Clin. Pathol.* 45:493-496.
9. Dati di archivio Liofilchem s.r.l. “Integral System Enterobatteri” (Giugno 2003).

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