Mycobacteria Growth Indicator Tube 7 mL Cu BACTEC MGIT 960 Supplement Kit

L000180JAA(01)
2013-07
Română

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REFERINȚE

- Bloom, B.R., and C.J.L. Murray. 1992. Tuberculosis: commentary on a reemergent killer. *Science* 257:1055-1064.
- Horsburg, C.R., Jr., 1991. *Mycobacterium avium* complex infection in the acquired immunodeficiency syndrome. *N. Engl. J. Med.* 324:1332-1338.
- Tenover, F.C., et al, 1993. The resurgence of tuberculosis: is your laboratory ready? *J. Clin. Microbiol.* 31:767-770.
- Kent, P.T., and G.P. Kubica. 1985. Public health mycobacteriology: a guide for the level III laboratory. USDHHS, Centers for Disease Control, Atlanta.
- Cohn, M.L., R.F. Waggoner and J.K. McClatchy. 1968. The 7H11 medium for the cultivation of mycobacteria. *Am. Rev. Respir. Dis.* 98:295-296.
- Youmans, G.P. 1979. Cultivation of mycobacteria, the morphology and metabolism of mycobacteria, p. 25-35. *Tuberculosis*. W.B. Saunders Co., Philadelphia.
- Clinical and Laboratory Standards Institute. 2005. Approved Guideline M29-A3. Protection of laboratory workers from occupationally acquired infections, 3rd ed. CLSI, Wayne, Pa.
- Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. *Infect. Control Hospital Epidemiol.* 17:53-80.
- U.S. Department of Health and Human Services. 2007. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 5th ed. U.S. Government Printing Office, Washington, D.C.
- Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). *Official Journal L262*, 17/10/2000, p. 0021-0045.
- Isenberg, Henry D. (ed.) 1992. *Clinical microbiology procedures handbook*. vol. 1. American Society for Microbiology, Washington, D.C.
- Clinical and Laboratory Standards Institute. 2004. Approved Standard M22-A3. Quality control for commercially prepared microbiological culture media, 3rd ed., CLSI, Wayne, Pa.

Manufacturer / Výrobce / Producent / Fabrikant / Tootja / Valmistaja / Fabricant / Hersteller / Κατασκευαστής / Gyártó / Ditta produttrice / Gamintojas / Producent / Fabricante / Vyroba / Tillverkare / Производител / Producător / Üretici / Proizvođač / Производитель / Атқарушы

Use by / Spotřebujte do / Anvendes før / Houdbaar tot / Kasutada enne / Viimeinkäyttöpäivä / A utiliser avant / Verwendbar bis / Ημερομηνία λήξης / Felhasználhatóság dátuma / Usare entro / Naudokite iki / Brukes før / Stosować do / Utilizar em / Použite do / Usar antes de / Använd före / Используйте до / A se utiliza până la / Son kullanna tarihi / Upotrebiti do / Использовать до / дейін пайдаланура / Upotrijebiti do / YYYY-MM-DD / YYYY-MM (MM = end of month) / RRRR-MM-DD / RRRR-MM (MM = konec měsíce) / ÁÁÁÁ-MM-DD / ÁÁÁÁ-MM (MM = slutning af måned) / JJJJ-MM-DD / JJJJ-MM (MM = einde maand) / AAAA-KK-PP / AAAA-KK (KK = kuu lõpp) / VVVV-KK-PP / VVVV-KK (kuukauden loppuun mennessä) / AAAA-MM-JJ / AAAA-MM (MM = fin du mois) / JJJJ-MM-TT / JJJJ-MM (MM = Monatsende) / EEEE-MM-HH / EEEE-MM (MM = τέλος του μήνα) / ÉÉÉÉ-HH-NN / ÉÉÉÉ-HH (HH = hónap utolsó napja) / AAAA-MM-GG / AAAA-MM (MM = fine mese) / MMMM-MM-DD / MMMM-MM (MM = mensesio pabaiga) / ÁÁÁÁ-MM-DD / ÁÁÁÁ-MM (MM = slutten av måneden) / RRRR-MM-DD / RRRR-MM (MM = koniec miesiąca) / AAAA-MM-DD / AAAA-MM (MM = fim do mês) / RRRR-MM-DD / RRRR-MM (MM = koniec mesiacu) / aaaa-mm-dd / aaaa-mm (mm = fin del mes) / ÁÁÁÁ-MM-DD / ÁÁÁÁ-MM (MM = slutet på månaden) / ГГГГ-ММ-ДД / ГГГГ-ММ (ММ = края на месеца) / АААА-ЛЛ-ЗЗ / АААА-ЛЛ (ЛЛ = sfârșitul lunii) / YYYY-AA-GG / YYYY-AA (AA = ayın sonu) / GGGG-MM-DD / GGGG-MM (MM = kraj meseca) / ГГГГ-ММ-ДД / ГГГГ-ММ (ММ = конец месяца) / ЖЖЖЖ-АА-КК / ЖЖЖЖ-АА (АА = айдың соңы) / GGGG-MM-DD / GGGG-MM (MM = kraj mjeseca)

Catalog number / Katalogové číslo / Katalognummer / Catalogusnummer / Kataloogi number / Tuotenumero / Numéro catalogue / Bestellnummer / Αριθμός καταλόγου / Katalogusszám / Numero di catalogo / Katalogo numeris / Numer katalogowy / Número do catálogo / Katalogové číslo / Número de catálogo / Каталоген номер / Număr de catalog / Katalog numarası / Kataloški broj / Номер по каталогу / Каталог нөмірі

Authorized Representative in the European Community / Autorizovaný zástupce pro Evropskou unii / Autoriseret repræsentant i EU / Erkend vertegenwoordiger in de Europese Unie / Volitatus esindaja Euroopa Nõukogus / Valtuutettu edustaja Euroopan yhteisössä / Représentant agréé pour la C.E.E. / Autorisierte EG-Vertretung / Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα / Hivatalos képviselő az Európai Unióban / Rappresentante autorizzato nella Comunità europea / Įgaliotasis atstovas Europos Bendrijoje / Autorisert representant i EU / Autoryzowane przedstawicielstwo w Unii Europejskiej / Representante autorizado na União Europeia / Autorizovaný zástupca v Európskom spoločenstve / Representante autorizado en la Comunidad Europea / Auktoriserad representant i EU / Оторизирован представитель в ЕС / Reprezentant autorizat în Uniunea Europeană / Avrupa Topluluğu Yetkili Temsilcisi / Ovlašćení predstavnik u Evropskoj zajednici / Уполномоченный представитель в Европейском сообществе / Европа қауымдастығындағы уәкілетті өкіл / Autorizuirani predstavnik u EU

In Vitro Diagnostic Medical Device / Lékařské zařízení určené pro diagnostiku in vitro / In vitro diagnostisk medicinsk anordning / Medisch hulpmiddel voor in vitro diagnose / In vitro diagnostika meditsiiniaparatuur / Lääkinnällinen in vitro -diagnostiikkalaitte / Dispositif médical de diagnostic in vitro / Medizinisches In-vitro-Diagnostikum / In vitro διαγνωστική ιατρική συσκευή / In vitro diagnosztikai orvosi eszköz / Dispositivo medico diagnostico in vitro. / In vitro diagnostikos prietaisas / In vitro diagnostisk medisinsk utstyr / Urządzenie medyczne do diagnostyki in vitro / Dispositivo médico para diagnóstico in vitro / Medicínska pomôcka na diagnostiku in vitro / Dispositivo médico de diagnóstico in vitro / Medicinsk anordning för in vitro-diagnostik / Медицински уред за диагностика ин витро / Aparatură medicală de diagnosticare in vitro / In Vitro Diagnostik Tibbi Cihaz / Medicinski uređaj za in vitro dijagnostiku / Медицинский прибор для диагностики ин витро / Жасанды жағдайда жүзрізетін медициналық диагностика аспабы / Medicínska pomagala za In Vitro Dijagnostiku

Temperature limitation / Teplotní omezení / Temperaturbegrænsning / Temperatuurlimiet / Temperatuuri piirang / Lämpötilarajoitus / Température limite / Zulässiger Temperaturenbereich / Οπο θερμότητας / Hőmérsékleti határ / Temperatura limite / Laikymo temperatūra / Temperaturbegrænsning / Ograniczenie temperatury / Limitação da temperatura / Ohraničenje toploty / Limitación de temperatura / Temperaturbegrænsning / Температури ограничения / Limitare de temperatură / Sıcaklık sınırlaması / Ograničenje temperature / Ограничение температуры / Температураны шектеу / Dozvoljena temperatura

Batch Code (Lot) / Kód (číslo) šarže / Batch kode (Lot) / Chargennummer (lot) / Partii kood / Eräkoodi (LOT) / Code de lot (Lot) / Chargencode (Chargenbezeichnung) / Κωδικός παρτίδας (Παρτίδα) / Tétel száma (Lot) / Codice del lotto (partita) / Partijos numeris (Lot) / Batch-kode (Serie) / Kod partii (seria) / Código do lote (Lote) / Kód série (šarža) / Código de lote (Lote) / Satskod (parti) / Код (Партида) / Număr lot (Lotul) / Parti Kodu (Lot) / Kod serije / Код партии (лот) / Топтама коды / Lot (kod)

Consult Instructions for Use / Prostudujte pokyny k použití / Læs brugsanvisningen / Raadpleeg gebruiksaanwijzing / Lugeda kasutusjuhendit / Tarkista käyttöohjeista / Consulter la notice d'emploi / Gebrauchsanweisung beachten / Συμβουλευτείτε τις οδηγίες χρήσης / Olvassa el a használati utasítást / Consultare le istruzioni per l'uso / Skaitykite naudojimo instrukcijas / Se i brugsanvisningen / Zobacz instrukcja użytkowania / Consulte as instruções de utilização / Pozri Pokyny na používanie / Consultar las instrucciones de uso / Se brugsanvisningen / Направете справка в инструкциите за употреба / Consultați instrucțiunile de utilizare / Kullanım Talimatları'na başvurun / Pogledajte uputstvo za upotrebu / См. руководство по эксплуатации / Пайдалану нұсқаулығымен танысып алыңыз / Koristi upute za upotrebu

Becton, Dickinson and Company
7 Loveton Circle
Sparks, MD 21152 USA
800-638-8663
www.bd.com/ds

Benex Limited
Pottery Road, Dun Laoghaire
Co. Dublin, Ireland

ATCC is a trademark of the American Type Culture Collection.
BD, BD Logo and all other trademarks are property of Becton, Dickinson and Company. © 2013 BD.

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Kit Components

Kit Product No.	Kit Product Description
245124	BACTEC MGIT 960 Growth Supplement Kit

Kit Component(s)	Kit Component(s) Description
800632JAA	MGIT PANTA, Lyophilized
WP245124JAA	BACTEC MGIT 960 Growth Supplement

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1. Identification

Product identifier

Product No.:	Product name:	Common name(s), synonym(s)
800632JAA	MGIT PANTA, Lyophilized	

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Hazardous Decomposition Products: Not known.

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Ingestion: No harmful effects expected in amounts likely to be ingested by accident.

Inhalation: Limited inhalation hazard at normal work temperatures.

Skin Contact: Negligible irritation to skin at ambient temperatures.

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CERCLA Hazardous Substance List (40 CFR 302.4):

<u>Chemical Identity</u>	<u>Reportable quantity</u>
Phosphoric acid, sodium salt (1:2)	5000 lbs.

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<u>Chemical Identity</u>	<u>Reportable quantity</u>
Phosphoric acid, sodium salt (1:2)	5000 lbs.

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<u>Chemical Identity</u>	<u>Reportable quantity</u>
Phosphoric acid, sodium salt (1:2)	Reportable quantity: 5000 lbs.

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Chemical Identity
Phosphoric acid, sodium salt (1:2)

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Chemical Identity
Phosphoric acid, sodium salt (1:2)

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Chemical Identity
Phosphoric acid, sodium salt (1:2)

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1. Identification

Product identifier

Product No.:	Product name:	Common name(s), synonym(s)
WP245124JAA	BACTEC MGIT 960 Growth Supplement	

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Other hazards which do not result in GHS classification: None.

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Mixtures

Chemical Identity	Common name and synonyms	CAS number	Content in percent (%)*
Ethanol		64-17-5	28.158%
9-Octadecenoic acid (9Z)-		112-80-1	1.224%

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6. Accidental release measures

Personal precautions, protective equipment and emergency procedures: No special precautionary health measures should be needed under anticipated conditions of use.

Methods and material for containment and cleaning up: No specific clean-up procedure noted.

Environmental Precautions: Avoid release to the environment.

7. Handling and storage

Precautions for safe handling: When using do not eat, drink or smoke. Read and follow manufacturer's recommendations. Use personal protective equipment as required.

Conditions for safe storage, including any incompatibilities: Store in a cool, dry place. Keep container tightly closed.

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Control Parameters

Occupational Exposure Limits

Chemical Identity	Type	Exposure Limit Values	Source
Ethanol	TWA	1,000 ppm 1,900 mg/m3	US. OSHA Table Z-1-A (29 CFR 1910.1000) (1989)
	TWA	1,000 ppm 1,900 mg/m3	US. Tennessee. OELs. Occupational Exposure Limits, Table Z1A (06 2008)
	AN ESL	1,000 ppb	US. Texas. Effects Screening Levels (Texas Commission on Environmental Quality) (12 2010)
	ST ESL	10,000 ppb	US. Texas. Effects Screening Levels (Texas Commission on Environmental Quality) (12 2010)
	AN ESL	1,880 µg/m3	US. Texas. Effects Screening Levels (Texas Commission on Environmental Quality) (12 2010)
	ST ESL	18,800 µg/m3	US. Texas. Effects Screening Levels (Texas Commission on Environmental Quality) (12 2010)
	TWA PEL	1,000 ppm 1,900 mg/m3	US. California Code of Regulations, Title 8, Section 5155. Airborne Contaminants (08 2010)
	STEL	1,000 ppm	US. ACGIH Threshold Limit Values (12 2010)

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	REL	1,000 ppm 1,900 mg/m3	US. NIOSH: Pocket Guide to Chemical Hazards (2005)
	PEL	1,000 ppm 1,900 mg/m3	US. OSHA Table Z-1 Limits for Air Contaminants (29 CFR 1910.1000) (02 2006)

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Relative density:	No data available.
Solubility(ies)	
Solubility in water:	Completely Soluble
Solubility (other):	Water.: No data available.
Partition coefficient (n-octanol/water):	No data available.
Auto-ignition temperature:	No data available.
Decomposition temperature:	No data available.
Viscosity:	Not determined.

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Ethanol soil - Very mobile liquid
9-Octadecenoic acid (9Z)- No data available.

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14. Transport information

DOTUN Number: Not regulated.
UN Proper Shipping Name: Not regulated.
Transport Hazard Class(es)
Class: Not regulated.
Label(s): Not regulated.
Packing Group: Not regulated.
Marine Pollutant: Not regulated.
Limited quantity Not regulated.
Excepted quantity Not regulated.

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1. Bloom, B.R., and C.J.L. Murray. 1992. Tuberculosis: commentary on a reemergent killer. *Science* 257:1055-1064.
2. Horsburg Jr., C.R. 1991. *Mycobacterium avium* complex infection in the acquired immunodeficiency syndrome. *N. Engl. J. Med.* 324:1332-1338.
3. Tenover, F.C., et al. 1993. The resurgence of tuberculosis: Is your laboratory ready? *J. Clin. Microbiol.* 31:767-770.
4. Cohn, M.L., R.F. Waggoner, and J.K. McClatchy. 1968. The 7H11 medium for the cultivation of mycobacteria. *Am. Rev. Resp. Dis.* 98:295-296.
5. Youmans, G.P. 1979. Cultivation of mycobacteria, the morphology and metabolism of mycobacteria, p. 25-35. *Tuberculosis*. W.B. Saunders Company, Philadelphia.
6. Kent, P.T., and G.P. Kubica. 1985. Public health mycobacteriology: A guide for the level III laboratory. USDHHS, Centers for Disease Control, Atlanta.
7. Bloodborne pathogens. Code of Federal Regulations, Title 29, Part 1910.1030, Federal Register 1991, 56:64175-64182.
8. Isenberg, Henry D. 1992. Clinical microbiology procedures handbook, vol. 1. American Society for Microbiology, Washington, D.C.

Manufacturer / Výrobce / Producent / Fabrikant / Tootja / Valmistaja / Fabricant / Hersteller / Κατασκευαστής / Gyártó / Ditta produttrice / Gamintojas / Producent / Fabricante / Výrobca / Tillverkare / Производител / Producător / Üretici / Proizvođač / Производитель

Use by / Spotřebujte do / Anvendes før / Houdbaar tot / Kasutada enne / Viimeinkäyttöpäivä / A utiliser avant / Verwendbar bis / Ημερομηνία λήξης / Felhasználhatóság dátuma / Usare entro / Naudokite iki / Brukes før / Stosować do / Utilizar em / Použite do / Usar antes de / Använd före / Используйте до / A se utiliza până la / Son kullanna tarihi / Uptrebiti do / Upotrebiti do / Использовать до

YYYY-MM-DD / YYYY-MM (MM = end of month) /
RRRR-MM-DD / RRRR-MM (MM = konec měsíce) /
AAAA-MM-DD / AAAA-MM (MM = slutning af måned) /
JJJJ-MM-DD / JJJJ-MM (MM = einde maand) /
AAAA-KK-PP / AAAA-KK (KK = kuu lõpp) /
VVVV-KK-PP / VVVV-KK (kuukauden loppuun mennessä) /
AAAA-MM-JJ / AAAA-MM (MM = fin du mois) /
JJJJ-MM-TT / JJJJ-MM (MM = Monatsende) /
EEEE-MM-HH / EEEE-MM (MM = τέλος του μήνα) /
EEEE-HH-NN / EEEE-HH (HH = hónap utolsó napja) /
AAAA-MM-GG / AAAA-MM (MM = fine mese) /
MMMM-MM-DD / MMMM-MM (MM = mensesio pabaiga) /
AAAA-MM-DD / AAAA-MM (MM = slutten av månaden) /
RRRR-MM-DD / RRRR-MM (MM = koniec miesiąca) /
AAAA-MM-DD / AAAA-MM (MM = fim do mês) /
RRRR-MM-DD / RRRR-MM (MM = koniec mesiacu) /
aaaa-mm-dd / aaaa-mm (mm = fin del mes) /
AAAA-MM-DD / AAAA-MM (MM = slutet på månaden) /
ГГГГ-ММ-ДД / ГГГГ-ММ (ММ = края на месеца) /
AAAA-LL-ZZ / AAAA-LL (LL = sfârșitul lunii) /
YYYY-AA-GG / YYYY-AA (AA = ayın sonu) /
GGGG-MM-DD / GGGG-MM (MM = kraj meseca) /
ГГГГ-ММ-ДД / ГГГГ-ММ (ММ = конец месяца)

Catalog number / Katalogové číslo / Katalognummer / Catalogusnummer / Kataloogi number / Tuotenumero / Numéro catalogue / Bestellnummer / Αριθμός καταλόγου / Katalógusszám / Numero di catalogo / Katalogo numeris / Numer katalogowy / Número do catálogo / Katalogové číslo / Número de catálogo / Каталоген номер / Număr de catalog / Katalog numarası / Kataloški broj / Номер по каталогу

Authorized Representative in the European Community / Autorizovaný zástupce pro Evropskou unii / Autoriseret repræsentant i EU / Erkend vertegenwoordiger in de Europese Unie / Volitatud esindaja Euroopa Nõukogus / Valtuutettu edustaja Euroopan yhteisössä / Représentant agréé pour la C.E.E. / Autorisierte EG-Vertretung / Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα / Hivatalos képviselő az Európai Unióban / Rappresentante autorizzato nella Comunità europea / Įgaliotasis atstovas Europos Bendrijoje / Autorisert representant i EU / Autoryzowane przedstawicielstwo w Unii Europejskiej / Representante autorizado na União Europeia / Autorizovaný zástupca v Európskom spoločenstve / Representante autorizado en la Comunidad Europea / Auktoriserad representant i EU / Оторизован представитель в EU / Représentant autorizat în Uniunea Europeană / Avrupa Topluluğu Yetkili Temsilcisi / Ovlašćeni predstavnik u Evropskoj zajednici / Уполномоченный представитель в Европейском сообществе

In Vitro Diagnostic Medical Device / Lékařské zařízení určené pro diagnostiku in vitro / In vitro diagnostisk medicinsk anordning / Medisch hulpmiddel voor in vitro diagnose / In vitro diagnostika meditsiiniaparatuur / Lääkinnällinen in vitro -diagnostiikkalaitte / Dispositif médical de diagnostic in vitro / Medizinisches In-vitro-Diagnostikum / In vitro διαγνωστική ιατρική συσκευή / In vitro diagnosztikai orvosi eszköz / Dispositivo medico diagnostico in vitro. / In vitro diagnostikos prietaisas / In vitro diagnostisk medisinsk utstyr / Urządzenie medyczne do diagnostyki in vitro / Dispositivo médico para diagnóstico in vitro / Medicínska pomôcka na diagnostiku in vitro / Dispositivo médico de diagnóstico in vitro / Medicínsk anordning för in vitro-diagnostik / Медицински уред за диагностика ин витро / Aparatură medicală de diagnosticare in vitro / In Vitro Diagnostik Tibbi Cihaz / Medicinski uređaj za in vitro dijagnostiku / Медицинский прибор для диагностики ин витро

Temperature limitation / Teplotní omezení / Temperaturbegrænsning / Temperatuurlimiet / Temperatuuri piirang / Lämpötilarajoitus / Température limite / Zulässiger Temperaturenbereich / Όριο θερμοκρασίας / Hőmérsékleti határ / Temperatura limite / Laikymo temperatūra / Temperaturbegrænsning / Ograniczenie temperatury / Limitação da temperatura / Ohraničenje teploty / Limitación de temperatura / Temperaturbegrænsning / Температурни ограничения / Limitare de temperatură / Sicaklık sınırlaması / Ograničenje temperature / Ограничение температуры

Batch Code (Lot) / Kód (číslo) šarže / Batch code (Lot) / Chargennummer (lot) / Partii kood / Eräkoodi (LOT) / Code de lot (Lot) / Chargencode (Chargenbezeichnung) / Κωδικός παρτίδας (Παρτίδα) / Tétel száma (Lot) / Codice del lotto (partita) / Partijos numeris (Lot) / Batch-kode (Serie) / Kod partii (seria) / Código do lote (Lote) / Kód série (šarža) / Código de lote (Lote) / Satskod (parti) / Код (Партида) / Număr lot (Lotul) / Parti Kodu (Lot) / Код партии (лот)

Consult Instructions for Use / Prostudujte pokyny k použití / Læs brugsanvisningen / Raadpleeg gebruiksaanwijzing / Lugeda kasutusjuhendit / Tarkista käyttöohjeista / Consulter la notice d'emploi / Gebrauchsanweisung beachten / Συμβουλευτείτε τις οδηγίες χρήσης / Olvassa el a használati utasítást / Consultare le istruzioni per l'uso / Skaitykite naudojimo instrukcijas / Se i brugsanvisningen / Zobacz instrukcja użytkowania / Consulte as instruções de utilização / Pozri Pokyny na používanie / Consultar las instrucciones de uso / Se brugsanvisningen / Направте справка в инструкциите за употреба / Consultați instrucțiunile de utilizare / Kullanım Talimatları'na başvurun / Pogledajte uputstvo za upotrebu / См. руководство по эксплуатации

Do not reuse / Nepoužívejte opakovaně / Må ikke genbruges / Niet opnieuw gebruiken / Mitte kasutada korduvalt / Ei saa käyttää uudelleen / Usage unique / Nicht wiederverwenden / Μην το ξαναχρησιμοποιείτε / Egyszer használatos / Non riutilizzare / Tik vienkartiniam naudojimui / Må ikke gjenbrukes / Nie stosować powtórnie / Não reutilizar / Nepouzívajte opakovane / No reusar / Får ej återanvändas / Не используйте отново / A nu se reutiliza / Tekrar kullanmayın / Ne upotrebljavajte ponovo / Не использовать повторно

Becton, Dickinson and Company
7 Loveton Circle
Sparks, Maryland 21152 USA
800-638-8663

BENEX Limited
Bay K 1a/d, Shannon Industrial Estate
Shannon, County Clare, Ireland
Tel: 353-61-47-29-20
Fax: 353-61-47-25-46

BACTEC™ MGIT™ 960 PZA Kit

For the Antimycobacterial Susceptibility Testing of *Mycobacterium tuberculosis*

I. INTENDED USE

The **BACTEC™ MGIT™ 960 PZA Kit** is used as a rapid qualitative procedure for susceptibility testing of *Mycobacterium tuberculosis*, from culture, to pyrazinamide (PZA). The **BACTEC MGIT 960 PZA Kit** is used with the **BACTEC MGIT 960 System**.

II. SUMMARY AND EXPLANATION

Antimycobacterial susceptibility testing is necessary for the proper treatment of patients with tuberculosis. The treatment of tuberculosis is commonly through a multiple drug regimen that includes the antimycobacterial drug pyrazinamide. It is important that the antimycobacterial drug prescribed shows appropriate activity against *Mycobacterium tuberculosis*, i.e., susceptibility of the isolate to the drug.

Multidrug resistant *Mycobacterium tuberculosis* (MDR-TB) has recently become a serious public health problem.¹ Resistance to any of the primary drugs, including pyrazinamide, makes the disease more difficult and expensive to treat. The rapid detection of these resistant isolates is critical to the effective patient management.

Two methods have been widely used for antimycobacterial susceptibility testing. The first method, known as the Method of Proportion,² uses Middlebrook and Cohn 7H10 Agar. It compares colony counts on drug-containing and drug-free media. The testing for pyrazinamide requires some modification from the general methods because the drug is active *in vitro* only at lower pH values.³ A modification to the method of proportion method was developed using a 7H10 agar medium at pH 5.5, with a drug concentration of 25-50 µg/mL.⁴ A limitation of the method is that at a pH of 5.5, many isolates of *M. tuberculosis* either fail to grow or grow poorly. Agar-based methods such as the agar proportion method have not proven to be satisfactory for PZA susceptibility testing because of failure of many isolates to grow when the agar has been acidified for the PZA test.

The second method, known as the **BACTEC 460TB** radiometric susceptibility method,⁵ is based on the production of radioactive ¹⁴C-labeled carbon dioxide by the growing mycobacteria, manifested by a Growth Index increase in the system. A modification to the **BACTEC 460TB** susceptibility method was developed using a modified 7H12 radiometric medium, **BACTEC PZA Test Medium**, with a reduced pH of 6.0.⁶ At this pH, PZA activity against mycobacteria can be determined without inhibiting the growth of most *M. tuberculosis* isolates. The **BACTEC 460TB PZA** susceptibility test uses a pyrazinamide drug concentration of 100 µg/mL. Susceptibility testing in the **BACTEC 460TB System** has proven to be satisfactory and is presently considered the reference method for PZA susceptibility testing. The National Committee for Clinical Laboratory Standards (NCCLS) recommends the **BACTEC 460TB** method for PZA susceptibility testing.²

Use of the **BACTEC MGIT 960 System** in combination with the **BACTEC MGIT 960 PZA kit** is a non-radiometric method of determining antimycobacterial susceptibility to PZA. The **BACTEC MGIT 960 PZA Kit** has been developed to allow susceptibility testing at a PZA concentration of 100 µg/mL. This concentration correlates with the concentration used in the **BACTEC 460TB System**.

III. PRINCIPLES OF THE PROCEDURE

BACTEC MGIT 960 PZA Medium is a tube containing a modified Middlebrook 7H9 Broth, which supports the growth and detection of mycobacteria at a reduced pH of 5.9. The **MGIT 960 PZA Medium** tube contains a fluorescent compound embedded in silicone on the bottom of a 16 x 100 mm round-bottom tube. The fluorescent compound is sensitive to the presence of oxygen dissolved in the broth. The initial concentration of dissolved oxygen quenches the emission from the compound, and little fluorescence can be detected. Later, actively growing and respiring microorganisms consume the oxygen, which allows the compound to fluoresce.

The **BACTEC MGIT 960 PZA Kit** is a 4 – 21 day qualitative test. The test is based on growth of the *M. tuberculosis* isolate in a drug-containing tube compared to a drug-free tube (Growth Control). The **BACTEC MGIT 960** instrument continually monitors tubes for increased fluorescence. Analysis of fluorescence in the drug-containing tube compared to the fluorescence of the Growth Control tube is used by the instrument to determine susceptibility results.

The **BACTEC MGIT 960** instrument automatically interprets these results and reports a susceptible or resistant result.

IV. REAGENTS

The **BACTEC MGIT PZA Medium** contains 110 µL of fluorescent indicator and 7 mL of PZA broth. The indicator contains Tris 4, 7-diphenyl-1, 10 phenanthroline ruthenium chloride pentahydrate in a silicone rubber base. The tubes are capped with a polypropylene cap. The pH is adjusted to 5.9.

Approximate Formula* Per L of Purified Water:

Modified Middlebrook 7H9 broth	5.9 g
Casein peptone	1.25 g

BACTEC MGIT 960 PZA Kit contains two lyophilized vials of pyrazinamide and six vials of PZA Supplement.

Approximate Formula* Per Vial Lyophilized drug:

Pyrazinamide	20,000 µg
--------------	-----------

BACTEC MGIT 960 PZA Supplement contains 15 mL of enrichment

Approximate Formula* Per L Purified Water:

Bovine albumin	50.0 g
Catalase	0.03 g
Dextrose	20.0 g
Oleic Acid	0.1 g
Polyoxyethylene stearate (POES)	1.1 g

**Adjusted and/or supplemented as required to meet performance criteria.*

Storage and reconstitution of reagents:

BACTEC MGIT 960 PZA Medium - On receipt, store at 2 – 25°C. DO NOT FREEZE. Broth should appear clear and colorless. Do not use if turbid. Minimize exposure to light. Tubes stored as labeled prior to use, may be inoculated up to the expiration date.

BACTEC MGIT 960 PZA Drug vials – On receipt, store the lyophilized drug vials at 2 – 8°C. Once reconstituted, the antibiotic solution may be frozen and stored at -20°C or colder up to six months, not to exceed the original expiration date. Once thawed, use immediately. Discard unused portions.

BACTEC MGIT 960 PZA Supplement – On receipt, store in dark at 2 – 8°C. Avoid freezing or overheating. Open and use prior to the expiration date. Minimize exposure to light.

Directions For Use:

Reconstitute each **BACTEC MGIT 960 PZA** lyophilized drug vial with **2.5 mL** of sterile distilled/deionized water to make a stock solution of 8000 µg/mL.

Warnings and Precautions: For *in vitro* Diagnostic Use.

POTENTIALLY INFECTIOUS TEST SPECIMEN: Pathogenic microorganisms, including hepatitis viruses and Human Immunodeficiency Virus, may be present in clinical specimens. “Standard Precautions”⁷⁻¹⁰ and institutional guidelines should be followed in handling all items contaminated with blood and other body fluids.

Working with *M. tuberculosis* growth in culture requires Biosafety Level (BSL) 3 practices, containment equipment and facilities.

Read and follow directions contained in all appropriate package inserts including the **BBL™ MGIT™** 7 mL Mycobacteria Growth Indicator Tube.

Prior to use, the user should examine the tubes and vials for evidence of contamination or damage. Discard any tubes or vials if they appear unsuitable. Dropped tubes should be examined carefully. If damage is seen, the tube should be discarded.

In the event of tube breakage: 1) Close the instrument drawers; 2) Turn off the instrument; 3) Vacate the area immediately; 4) Consult your facility/CDC guidelines. An inoculated, leaking or broken tube may produce an aerosol of mycobacteria; appropriate handling should be observed.

Autoclave all inoculated **MGIT** tubes prior to disposal.

V. SPECIMEN PREPARATION

All preparations detailed below must be from pure cultures of *M. tuberculosis*. The laboratory should confirm, by appropriate identification techniques, that the isolate to be tested is a pure culture of *M. tuberculosis*.

Preparation of the Isolate from Solid Media:

1. Add 4 mL of **BBL** Middlebrook 7H9 Broth (or **BBL MGIT** broth) to a 16.5 x 128 mm sterile tube with cap containing 8–10 glass beads.
2. Scrape with a sterile loop as many colonies as possible from growth no more than fourteen days old, trying not to remove any solid medium. Suspend the colonies in the Middlebrook 7H9 Broth.
3. Vortex the suspension for 2–3 minutes to break up the larger clumps. The suspension should exceed a 1.0 McFarland standard in turbidity.
4. Let the suspension sit for 20 min without disturbing.
5. Transfer the supernatant fluid to another 16.5 x 128 mm sterile tube with cap (avoid transferring any of the sediment) and let sit for another 15 min.
6. Transfer the supernatant fluid (it should be smooth, free of any clumps) to a third 16.5 x 128 mm sterile tube. NOTE: The organism suspension should be greater than a 0.5 McFarland standard at this step.
7. Adjust the suspension to a 0.5 McFarland standard by visual comparison to a 0.5 McFarland turbidity standard. Do not adjust below a 0.5 McFarland standard.
8. Dilute 1 mL of the adjusted suspension in 4 mL of sterile saline (1:5 dilution). Use this as the AST inoculum and proceed to "Inoculation Procedure for **BACTEC MGIT** 960 PZA Susceptibility Test."

Preparation from a Positive BACTEC MGIT Tube:

1. The first day of an instrument positive **BACTEC MGIT** tube is considered Day 0.
2. For the preparation of the test inoculum, a positive 7 mL **MGIT** tube should be used the day **after** it first becomes positive on the **BACTEC MGIT** 960 instrument (Day 1), up to and including the fifth day (Day 5) after instrument positivity. A tube which has been positive longer than five days should be subcultured to a fresh 7 mL **MGIT** tube containing **BACTEC MGIT** 960 Growth Supplement and tested on the **BACTEC MGIT** 960 instrument until positive, and used from one to five days following positivity.
3. If the tube is a Day 1 or Day 2 positive, no dilution is required. Use this as the AST inoculum and proceed to "Inoculation Procedure for **BACTEC MGIT** 960 PZA Susceptibility Test."
4. If the tube is a Day 3, Day 4, or Day 5 positive, then dilute 1 mL of the positive broth in 4 mL of sterile saline (1:5 dilution). Use this as the AST inoculum and proceed to "Inoculation Procedure for **BACTEC MGIT** 960 PZA Susceptibility Test."

VI. PROCEDURE

Materials Provided: **BACTEC MGIT** 960 PZA Kit containing two vials each lyophilized drug and six vials of PZA Supplement (approximately 50 tests per kit).

Materials Required But Not Provided: **BACTEC MGIT** 960 PZA Medium (25 tubes per carton), ancillary culture media, reagents, quality control organisms and laboratory equipment as required for this procedure.

Inoculation Procedure for BACTEC MGIT 960 PZA Susceptibility Test:

Important considerations when preparing the PZA AST Set are the proper reconstitution of the lyophilized drug, use of pure culture and the proper dilution of the organism for the Growth Control and PZA tube. It is important to add drug only to the corresponding **MGIT** tube labeled “PZA”. Use only the **MGIT** 960 PZA supplement supplied with the kit and **MGIT** 960 PZA Medium tubes when performing the PZA AST set.

1. Label two 7 mL **MGIT** 960 PZA Medium tubes for each test isolate. Label one as GC (Growth Control), one as PZA. Place the tubes in the correct sequence in the two tube AST set carrier (see **BACTEC MGIT** 960 User’s Manual, AST Instructions).
2. Aseptically add 0.8 mL of **BACTEC MGIT** 960 PZA Supplement to each tube.
3. Using a micropipet, aseptically pipet 100 µL of 8000 µg/mL **BACTEC MGIT** 960 PZA drug solution to the appropriately labeled **MGIT** PZA tube. No PZA drug solution should be added to the **MGIT** GC tube.

Drug	Concentration of Drug after Reconstitution*	Volume Added to MGIT Tubes for Test	Final Concentration in MGIT Tubes
MGIT PZA	8000 µg/mL	100 µL	100 µg/mL*

*PZA must be reconstituted using 2.5 mL sterile/deionized water to achieve the concentration indicated.

4. **Growth Control tube preparation and inoculation:** Aseptically pipet 0.5 mL of the AST inoculum (see “SPECIMEN PREPARATION”) into 4.5 mL of sterile saline to prepare the **1:10** Growth Control suspension. Mix the Growth Control suspension thoroughly. Inoculate 0.5 mL of the **1:10** Growth Control suspension into the **MGIT** tube labeled “GC.”

NOTE: It is important to use an appropriately prepared **1:10** dilution for the “GC” tube to ensure accurate AST results and avoid PZA AST set errors.

5. **Drug-containing tube inoculation:** Aseptically pipet 0.5 mL of the AST inoculum (see “SPECIMEN PREPARATION”) into the **MGIT** tube labeled “PZA”.
6. Tightly recap the tubes. Mix tubes thoroughly by gentle inversion three to four times.
7. Enter the PZA set into the **BACTEC MGIT** 960 instrument using the AST set entry feature (refer to the **BACTEC MGIT** 960 User’s Manual, AST Instructions). Ensure that the Growth Control tube is in the first left tube position. Select PZA as the drug in the 2 tube AST set carrier definition when performing the AST set entry.
8. Streak 0.1 mL of the organism suspension to a **Trypticase**[™] Soy Agar with 5% Sheep Blood (TSA II) plate. Enclose in a plastic bag. Incubate at 35 – 37°C.
9. Check the blood agar plate at 48 h for bacterial contamination. If the blood agar plate shows no growth, then allow PZA testing to proceed. If the blood agar plate shows growth, discard the PZA set (refer to the **BACTEC MGIT** 960 User’s Manual, AST Instructions) and repeat testing with a pure culture of *Mycobacterium tuberculosis*.

VII USER QUALITY CONTROL

Upon receipt of a new shipment or lot number of **BACTEC MGIT 960 PZA** Kit vials or **BACTEC MGIT PZA** Medium, it is recommended that the control organism shown below be tested. The control organism should be a pure culture and the culture should be prepared according to "SPECIMEN PREPARATION" instructions.

The quality control (QC) AST Set should be prepared according to the "Inoculation Procedure for **BACTEC MGIT 960 PZA** Susceptibility Test" instructions. Important considerations when preparing the QC AST Set are the proper reconstitution of the lyophilized drug, use of pure culture and the proper dilution of the QC organism for the Growth Control and PZA tubes. It is important to add drug only to the corresponding **MGIT** tube labeled "PZA".

The same control organism should be run as batch QC once each week when susceptibility testing is performed. Observation of the proper results, as shown below, within 4 - 20 days indicates that the **BACTEC MGIT 960 PZA** reagents are ready for use in testing patient isolates.

If the proper results are not observed, do not report patient results. Repeat QC and any patient isolates affected by the initial QC failure. If the repeat QC does not perform as expected, do not report patient results. Do not use the product until you have contacted Technical Services at (800) 638-8663 (United States Only).

Strain	GC	MGIT PZA
<i>M. tuberculosis</i> ATCC™ 27294	Positive	Susceptible

During the external evaluation of the **BACTEC MGIT 960 PZA** Kit the average time to result for the control organism was seven days with a range of four to eleven days. The most common causes of QC failures during the external evaluation were over-inoculated PZA Sets and contaminated QC cultures.

VIII. RESULTS

The **BACTEC MGIT 960** instrument will monitor AST sets until a susceptible or resistant determination is made. Once the set testing is completed, the results are reported by the **BACTEC MGIT 960** instrument (refer to the **BACTEC MGIT 960** User's Manual, AST Instructions). The **BACTEC MGIT 960** instrument will report an AST Set result as an Error ("X"), no susceptibility interpretation, when certain conditions occur that may affect the test results. Conditions that may result in an Error ("X") result are described in the AST Instructions,

Section 6 – Troubleshooting of the **BACTEC MGIT 960** User's Manual.

It is important to include the test method, drug name and concentration when reporting results. The Pulmonary and/or Infectious Disease specialist in TB control should be consulted concerning the appropriate therapeutic regimen and dosages.

Mono-resistance to pyrazinamide is uncommon, therefore in the event of unexpected resistant results, verify purity and identification of the isolate tested as *M. tuberculosis*. Guidelines for mycobacterial purity checks can be found in the NCCLS M24 standard.²

BACTEC MGIT 960 PZA result reporting

Drug (concentration)	MGIT 960 result	Recommended Report	Action
PZA (100 µg/mL)	Susceptible	Isolate tested with BACTEC MGIT 960 [PZA/100 µg/mL] and result is susceptible.	No action.
	Resistant	Isolate tested with BACTEC MGIT 960 [PZA/100 µg/mL] and result is resistant.	If isolate is mono-resistant to PZA, confirm that isolate tested is a pure culture of <i>Mycobacterium tuberculosis</i> .
	Error "X"	No report.	Repeat test.

IX. LIMITATIONS OF THE PROCEDURE

The **BACTEC MGIT 960** PZA susceptibility test does not interpret the degree of susceptibility of the isolate being tested. Results are reported as either susceptible or resistant.

The **BACTEC MGIT 960** PZA susceptibility test can only be performed using the **BACTEC MGIT 960** instrument. The PZA Sets cannot be read manually.

Use only pure cultures of *M. tuberculosis*. Cultures that are contaminated or that may contain multiple species of mycobacteria may give erroneous results and should not be tested. Direct testing from clinical specimens is not recommended.

Suspensions made from solid media must be allowed to settle for the prescribed times prior to standardization. Inoculum preparations made from solid media should be visually compared to a 0.5 McFarland turbidity standard; failure to do so may give inaccurate results or cause an AST Set error.

Failure to use the 1:5 dilution of the organism suspension, when indicated, to inoculate the drug containing tubes may give inaccurate results.

Failure to use a 1:10 dilution of the organism suspension for the inoculation of the Growth Control tube may give inaccurate results or cause an AST Set error.

Failure to reconstitute the PZA drug with the appropriate volume of sterile distilled / deionized water may give inaccurate results.

Thorough mixing of inoculated tubes is important. Failure to mix the tubes adequately may lead to false resistant results.

Failure to load the tubes of the AST Set into the AST Set Carrier in the proper sequence may give inaccurate results. Failure to select the appropriate set carrier drug definition may result in invalid or inaccurate results.

Failure to load the AST Set into the instrument correctly will result in an anonymous condition that must be resolved within eight hours. If condition is not resolved within eight hours, the AST Set must be discarded and set up again.

Failure to use the **BACTEC MGIT 960 PZA Supplement** in the PZA AST set may give inaccurate results. DO NOT add **BACTEC MGIT 960 SIRE Supplement** or **BACTEC MGIT 960 Growth Supplement** to the PZA AST set.

Failure to use **BACTEC MGIT PZA Medium** for the PZA AST set may give inaccurate results. DO NOT substitute **BBL MGIT 7 mL Mycobacteria Growth Indicator Tubes** for **BACTEC MGIT PZA Medium**.


X. EXPECTED VALUES

A total of 118 clinical isolates of *M. tuberculosis* were tested with the **BACTEC MGIT 960 PZA** susceptibility test at four geographically diverse sites. The testing included both fresh clinical and subcultured isolates from both liquid and solid culture sources. A total of 228 PZA susceptibility tests (liquid and solid) were performed.

During the external evaluation of the **BACTEC MGIT 960 PZA Kit**, there were nine PZA tests from clinical trial isolates that required repeat testing due to contamination (six isolates) or over-inoculation/procedural errors (three isolates).

The average time-to-result for the **BACTEC MGIT 960 PZA** susceptibility test is seven days with a range from four to seventeen days. The data are shown in Figure 1.

Figure 1: Distribution of **BACTEC MGIT 960** PZA AST Time to Result

XI. PERFORMANCE CHARACTERISTICS

Analytical Studies

Liquid and Solid Media AST Inoculum Ranges:

Liquid media - The recommended procedure for preparing a PZA Set from a positive **MGIT** 7 mL tube uses a direct inoculum on Day 1 and Day 2 post-positivity and a dilute (1:5) inoculum on Day 3 to Day 5 post-positivity. Internal studies show that inocula prepared from a Day 1 to Day 5 positive **MGIT** 7 mL tube range between 2.0×10^4 to 7.5×10^6 CFU/mL.

Solid media - The recommended procedure for preparing a PZA Set from growth on solid media (up to 14 days after first visible growth is seen) uses a 1:5 dilution of an organism suspension equivalent to a 0.5 McFarland Standard. Internal studies show that inocula prepared from solid medium culture range between 2.1×10^5 to 3.9×10^6 CFU/mL.

Lot Reproducibility:

Lot reproducibility was evaluated using twenty-five *M. tuberculosis* strains (including three ATCC™ strains). Each strain was tested in triplicate with the **BACTEC MGIT 960 PZA** susceptibility test. Each replicate represented a separate test condition differentiated by lot of PZA drug, PZA Supplement and PZA medium used (three lots each).

Observed results were compared to the expected results. The overall reproducibility for the **BACTEC MGIT 960 PZA** susceptibility test is 96.8%.

CDC Challenge Panel Testing:

The performance of the **BACTEC MGIT 960 PZA** susceptibility test was evaluated using a panel of challenge strains obtained from the Centers for Disease Control and Prevention (CDC), Atlanta, GA, USA. The panel consisted of nine strains of *M. tuberculosis* with known susceptibility patterns (using **BACTEC 460TB**). The panel was tested in triplicate with the **BACTEC MGIT 960 PZA** susceptibility test. The **BACTEC MGIT 960 PZA** results were compared to the CDC expected results. The overall agreement with CDC expected results for the **BACTEC MGIT 960 PZA** test is 98.7%.

Clinical Evaluation

The **BACTEC MGIT 960 PZA** susceptibility test was evaluated at four geographically diverse clinical sites composed of regional reference centers and university hospital-based laboratories, including two ex-US sites. The **BACTEC MGIT 960 PZA** susceptibility test was compared to the **BACTEC 460TB PZA** susceptibility test method.

Reproducibility Testing:

The reproducibility of the **BACTEC MGIT 960 PZA** susceptibility test was evaluated at the clinical sites using a panel of five qualified strains. The **BACTEC MGIT 960 PZA** test results were compared to the expected results. The overall reproducibility for the **BACTEC MGIT 960 PZA** susceptibility test is 94%.

CDC Challenge Panel Testing:

The performance of the **BACTEC MGIT 960 PZA** susceptibility test was evaluated at each of the four clinical sites using a panel of challenge strains obtained from the Centers for Disease Control and Prevention (CDC), Atlanta, GA, USA. The panel consisted of nine strains of *M. tuberculosis* with known susceptibility patterns (using **BACTEC 460TB**). Of the thirty-six PZA results collected with the **BACTEC MGIT 960 PZA** susceptibility test, thirty-three agreed with the CDC expected results. The calculated percent agreement to the CDC expected results for the **BACTEC MGIT 960 PZA** susceptibility test is 91.7%.

Clinical Isolate Testing:

A total of 118 clinical isolates of *M. tuberculosis* were tested with the **BACTEC MGIT 960** PZA susceptibility test and the **BACTEC 460TB** PZA susceptibility test. This included testing of both fresh clinical and subcultured isolates from both liquid and solid culture sources. This generated a total of 228 test results.

Table 1 presents the results from clinical isolate testing for PZA drug at 100 µg/mL from liquid source cultures, from solid source cultures and both source cultures combined.

Table 1: Clinical Isolate Results – BACTEC MGIT 960 PZA test Compared to BACTEC 460TB test

Source	# Tests	BACTEC 460TB System		BACTEC MGIT 960 System			
		Expected PZA Results		Susceptible Results		Resistant Results	
		S	R	# agree	Category agreement % (95% CI)	# agree	Category agreement % (95% CI)
LIQUID	112	89	23	88	98.9% (93.9-100)	22	95.7% (78.1-99.9)
SOLID	113*	90	23	88	97.8% (92.2-99.7)	20	87.0% (66.4-97.2)
ALL	225*	179	46	176	98.3% (95.2-99.7)	42	91.3% (79.2-97.6)

*Three **BACTEC 460TB** borderline results are not included in this table.

All isolates with discordant **BACTEC MGIT 960** PZA results were tested using the **BACTEC 460TB** PZA susceptibility test at two independent sites. Discordant results were those strains where the **BACTEC MGIT 960** PZA result differed from the **BACTEC 460TB** PZA result. Borderline results are not included in the performance calculations for the **BACTEC MGIT 960** PZA Kit.

Of the four discordant PZA susceptible (S-**MGIT 960**, R-**BACTEC 460TB**) isolates tested, one had susceptible results from both independent sites and the other three had resistant results from both independent sites. Of the three discordant PZA resistant (R-**MGIT 960**, S-**BACTEC 460TB**) isolates tested, all isolates had susceptible results from both independent sites.

Two of the three **BACTEC 460TB** borderline PZA results (S-**MGIT 960**, B-**BACTEC 460TB**) had susceptible results from both independent sites. One of the three **BACTEC 460TB** borderline PZA results (R-**MGIT 960**, B-**BACTEC 460TB**) had one independent site determine a susceptible result. The other independent site determined a borderline result.

XII. AVAILABILITY

Cat. No. Description

245128 **BACTEC™ MGIT™ 960** PZA Kit, carton of 2 lyophilized drug vials and 6 PZA Supplements.

245115 **BACTEC™ MGIT™ PZA** Medium, carton of 25 tubes.

XIII. REFERENCES

1. Barenfanger, J. 1993. Making your lab safe against multi-drug resistant *Mycobacterium tuberculosis*. Clin. Microbiol. Newsl. 15: 76-80.
2. National Committee for Clinical Laboratory Standards. 2001. Approved standard: M24-A. Susceptibility testing of mycobacteria, *Nocardia*, and other aerobic actinomycetes . NCCLS, Wayne, Pa.
3. Butler, W.R. and Kilburn. 1982. Improved method for testing susceptibility of *Mycobacterium tuberculosis* to pyrazinamide. J.Clin.Microbiol. 16:1106-1109.
4. Heifets, L.B. and Iseman, M.D. 1985. Radiometric method for testing susceptibility of *Mycobacterium tuberculosis* to pyrazinamide in 7H12 broth. J.Clin.Microbiol. 21:200-204.
5. BD Diagnostic Systems. **BACTEC**™ 460TB System Product and Procedure Manual.
6. Salfinger, M. et al. 1989. Rapid radiometric method for pyrazinamide susceptibility testing of *Mycobacterium tuberculosis*. Res. Microbiol 1989. 140:301-309.
7. National Committee for Clinical Laboratory Standards. 2002. Approved Guideline M29-A2. Protection of laboratory workers from occupationally acquired infections. NCCLS, Wayne, Pa.
8. Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. Infect. Control Hospital Epidemiol. 17:53-80.
9. U.S. Department of Health and Human Services. 1999. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC) 4th ed. U.S. Government Printing Office, Washington, D.C.
10. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). Official Journal L262, 17/10/2000, p. 0021-0045.

XIV. TECHNICAL INFORMATION

In the United States, telephone Technical Services, toll free (800) 638-8663.

Approved by:

Supervisor:_____

Date:_____

Director:_____

Date:_____

Effective Date:_____

Reviewed by

Date:

PI Rev 07/02

Rev 07/02

Becton, Dickinson and Company
BD, Franklin Lakes, NJ 07417
USA
www.bd.com

Kit Components

Kit Product No.	Kit Product Description
245128	Kit Drug Pza Mgit

Kit Component(s)	Kit Component(s) Description
8019891	MGIT PZA
WP245128JAA	SUPPLEMENT PZA

IMDG

Special precautions for user: Not regulated.

IATA

Special precautions for user: Not regulated.

Please note: If a listed component does not have a corresponding document included, this means that the product is not hazardous and does not require a Safety Data Sheet.

Becton, Dickinson and
Company
BD, Franklin Lakes, NJ
07417 USA
www.bd.com

SAFETY DATA SHEET

1. Identification

Product identifier

Product No.:	Product name:	Common name(s), synonym(s)
8019891	MGIT PZA	

Other means of identification

SDS number: 088100156805

Recommended use and restriction on use

Recommended use: Laboratory Chemicals

Restrictions on use: None known.

Manufacturer/Importer/Supplier/Distributor Information

Manufacturer

Company Name: BD Diagnostic Systems
Address: 7 Loveton Circle
21152 Sparks, MD USA
Telephone: 1 410 771 0100 or 1 800 638 8663
Fax:
Contact Person: Tech Services

Emergency telephone number: ChemTrec 1 800 424 9300

2. Hazard(s) identification

Hazard Classification

Not classified

Label Elements

Hazard Symbol: No symbol
Signal Word: No signal word.
Hazard Statement: Not applicable
Precautionary Statements: Not applicable

Other hazards which do not
result in GHS classification: None.

3. Composition/information on ingredients

Becton, Dickinson and Company
BD, Franklin Lakes, NJ
07417 USA
www.bd.com

Mixtures

Composition Comments: The components are not hazardous or are below required disclosure limits.

4. First-aid measures

General information: Get medical attention if symptoms occur.

Ingestion: Get medical attention if symptoms occur.

Inhalation: Provide fresh air, warmth and rest, preferably in comfortable upright sitting position.

Skin Contact: Wash contact areas with soap and water. Remove contaminated clothing. Launder contaminated clothing before reuse.

Eye contact: Flush thoroughly with water. If irritation occurs, get medical assistance.

Most important symptoms/effects, acute and delayed

Symptoms: No data available.

Indication of immediate medical attention and special treatment needed

Treatment: No data available.

5. Fire-fighting measures

General Fire Hazards: Extinguish all ignition sources. Avoid sparks, flames, heat and smoking. Ventilate. Use water spray to keep fire-exposed containers cool.

Suitable (and unsuitable) extinguishing media

Suitable extinguishing media: Water spray, fog, CO₂, dry chemical, or alcohol resistant foam.

Unsuitable extinguishing media: None known.

Specific hazards arising from the chemical: None known.

Special protective equipment and precautions for firefighters

Special fire fighting procedures: No unusual fire or explosion hazards noted.

Special protective equipment for fire-fighters: Firefighters must use standard protective equipment including flame retardant coat, helmet with face shield, gloves, rubber boots, and in enclosed spaces, SCBA.

Becton, Dickinson and Company
BD, Franklin Lakes, NJ
07417 USA
www.bd.com

6. Accidental release measures

Personal precautions, protective equipment and emergency procedures:	No special precautionary health measures should be needed under anticipated conditions of use.
Methods and material for containment and cleaning up:	No specific clean-up procedure noted.
Environmental Precautions:	Avoid release to the environment.

7. Handling and storage

Precautions for safe handling:	When using do not eat, drink or smoke. Read and follow manufacturer's recommendations. Use personal protective equipment as required.
Conditions for safe storage, including any incompatibilities:	Store in a cool, dry place. Keep container tightly closed.

8. Exposure controls/personal protection

Control Parameters

Occupational Exposure Limits

None of the components have assigned exposure limits.

Appropriate Engineering Controls

No special requirements under ordinary conditions of use and with adequate ventilation.

Individual protection measures, such as personal protective equipment

General information:	Always observe good personal hygiene measures, such as washing after handling the material and before eating, drinking, and/or smoking. Routinely wash work clothing to remove contaminants. Discard contaminated footwear that cannot be cleaned.
Eye/face protection:	Wear safety glasses with side shields (or goggles).
Skin Protection	
Hand Protection:	Chemical resistant gloves Suitable gloves can be recommended by the glove supplier. Wash hands after contact.
Other:	Wear a lab coat or similar protective clothing.
Respiratory Protection:	Respiratory protection not required.
Hygiene measures:	Observe good industrial hygiene practices.

Becton, Dickinson and
Company
BD, Franklin Lakes, NJ
07417 USA
www.bd.com

9. Physical and chemical properties

Appearance

Physical state:	solid
Form:	solid
Color:	According to product specification.
Odor:	Characteristic
Odor threshold:	No data available.
pH:	No data available.
Melting point/freezing point:	No data available.
Initial boiling point and boiling range:	No data available.
Flash Point:	No data available.
Evaporation rate:	No data available.
Flammability (solid, gas):	No data available.
Upper/lower limit on flammability or explosive limits	
Flammability limit - upper (%):	No data available.
Flammability limit - lower (%):	No data available.
Explosive limit - upper (%):	No data available.
Explosive limit - lower (%):	No data available.
Vapor pressure:	No data available.
Vapor density:	No data available.
Relative density:	No data available.
Solubility(ies)	
Solubility in water:	No data available.
Solubility (other):	No data available.
Partition coefficient (n-octanol/water):	No data available.
Auto-ignition temperature:	No data available.
Decomposition temperature:	No data available.
Viscosity:	No data available.

10. Stability and reactivity

Reactivity:	Stable under normal temperature conditions and recommended use.
Chemical Stability:	Material is stable under normal conditions.
Possibility of hazardous reactions:	Not known.
Conditions to avoid:	Avoid exposure to high temperatures or direct sunlight.
Incompatible Materials:	Strong oxidizers.

Becton, Dickinson and Company
BD, Franklin Lakes, NJ
07417 USA
www.bd.com

Hazardous Decomposition Products:	Not known.
--	------------

11. Toxicological information

General information:	No data on possible toxicity effects have been found.
-----------------------------	---

Information on likely routes of exposure

Ingestion:	No harmful effects expected in amounts likely to be ingested by accident.
-------------------	---

Inhalation:	Limited inhalation hazard at normal work temperatures.
--------------------	--

Skin Contact:	Negligible irritation to skin at ambient temperatures.
----------------------	--

Eye contact:	Do not get in eyes.
---------------------	---------------------

Symptoms related to the physical, chemical and toxicological characteristics

Ingestion:	No data available.
-------------------	--------------------

Inhalation:	No data available.
--------------------	--------------------

Skin Contact:	No data available.
----------------------	--------------------

Eye contact:	No data available.
---------------------	--------------------

Information on toxicological effects

Acute toxicity (list all possible routes of exposure)

Oral

Product:	No data available.
-----------------	--------------------

Dermal

Product:	No data available.
-----------------	--------------------

Inhalation

Product:	No data available.
-----------------	--------------------

Repeated dose toxicity

Product:	No data available.
-----------------	--------------------

Skin Corrosion/Irritation

Product:	No data available.
-----------------	--------------------

Serious Eye Damage/Eye Irritation

Product:	No data available.
-----------------	--------------------

Becton, Dickinson and Company
BD, Franklin Lakes, NJ
07417 USA
www.bd.com

Respiratory or Skin Sensitization

Product: No data available.

Carcinogenicity

Product: No data available.

IARC Monographs on the Evaluation of Carcinogenic Risks to Humans:

No carcinogenic components identified

US. National Toxicology Program (NTP) Report on Carcinogens:

No carcinogenic components identified

US. OSHA Specifically Regulated Substances (29 CFR 1910.1001-1050):

No carcinogenic components identified

Germ Cell Mutagenicity

In vitro

Product: No data available.

In vivo

Product: No data available.

Reproductive toxicity

Product: No data available.

Specific Target Organ Toxicity - Single Exposure

Product: No data available.

Specific Target Organ Toxicity - Repeated Exposure

Product: No data available.

Aspiration Hazard

Product: No data available.

Other effects: None known.

12. Ecological information

Ecotoxicity:

Acute hazards to the aquatic environment:

Fish

Product: No negative effects on the aquatic environment are known.

Becton, Dickinson and Company
BD, Franklin Lakes, NJ
07417 USA
www.bd.com

Aquatic Invertebrates

Product: No negative effects on the aquatic environment are known.

Chronic hazards to the aquatic environment:

Fish

Product: No negative effects on the aquatic environment are known.

Aquatic Invertebrates

Product: No negative effects on the aquatic environment are known.

Toxicity to Aquatic Plants

Product: No negative effects on the aquatic environment are known.

Persistence and Degradability

Biodegradation

Product: Expected to be readily biodegradable.

BOD/COD Ratio

Product: No data available.

Bioaccumulative potential

Bioconcentration Factor (BCF)

Product: No data available.

Partition Coefficient n-octanol / water (log K_{ow})

Product: No data available.

Mobility in soil: No data available.

Known or predicted distribution to environmental compartments

Other adverse effects: The product is not expected to be hazardous to the environment.

13. Disposal considerations

General information: Dispose of waste and residues in accordance with local authority requirements.

Disposal instructions: No specific disposal method required.

Contaminated Packaging: No data available.

Becton, Dickinson and Company
BD, Franklin Lakes, NJ
07417 USA
www.bd.com

14. Transport information

DOTUN Number:	Not regulated.
UN Proper Shipping Name:	Not regulated.
Transport Hazard Class(es)	
Class:	Not regulated.
Label(s):	Not regulated.
Packing Group:	Not regulated.
Marine Pollutant:	Not regulated.
Limited quantity	Not regulated.
Excepted quantity	Not regulated.
Special precautions for user:	Not regulated.

IMDG

UN Number:	Not regulated.
UN Proper Shipping Name:	Not regulated.
Transport Hazard Class(es)	
Class:	Not regulated.
Subsidiary risk:	Not regulated.
EmS No.:	Not regulated.
Packing Group:	Not regulated.
Environmental Hazards	
Marine Pollutant:	Not regulated.
Special precautions for user:	Not regulated.

IATA

UN Number:	Not regulated.
Proper Shipping Name:	Not regulated.
Transport Hazard Class(es):	
Class:	Not regulated.
Subsidiary risk:	Not regulated.
Packing Group:	Not regulated.
Environmental Hazards	
Marine pollutant:	Not regulated.
Special precautions for user:	Not regulated.

15. Regulatory information

US Federal Regulations

TSCA Section 12(b) Export Notification (40 CFR 707, Subpt. D)

None present or none present in regulated quantities.

US. OSHA Specifically Regulated Substances (29 CFR 1910.1001-1050)

None present or none present in regulated quantities.

**Becton, Dickinson and
Company**
BD, Franklin Lakes, NJ
07417 USA
www.bd.com

CERCLA Hazardous Substance List (40 CFR 302.4):

None present or none present in regulated quantities.

Superfund Amendments and Reauthorization Act of 1986 (SARA)

Hazard categories

Not classified

Not classified

SARA 302 Extremely Hazardous Substance

None present or none present in regulated quantities.

SARA 304 Emergency Release Notification

None present or none present in regulated quantities.

SARA 311/312 Hazardous Chemical

None present or none present in regulated quantities.

SARA 313 (TRI Reporting)

None present or none present in regulated quantities.

Clean Water Act Section 311 Hazardous Substances (40 CFR 117.3)

None present or none present in regulated quantities.

Clean Air Act (CAA) Section 112(r) Accidental Release Prevention (40 CFR 68.130):

None present or none present in regulated quantities.

US State Regulations

US. California Proposition 65

No ingredient regulated by CA Prop 65 present.

US. New Jersey Worker and Community Right-to-Know Act

No ingredient regulated by NJ Right-to-Know Law present.

US. Massachusetts RTK - Substance List

No ingredient regulated by MA Right-to-Know Law present.

US. Pennsylvania RTK - Hazardous Substances

No ingredient regulated by PA Right-to-Know Law present.

US. Rhode Island RTK

No ingredient regulated by RI Right-to-Know Law present.

16. Other information, including date of preparation or last revision
--

Issue Date: 07/17/2018

Version #: 1.0

**Becton, Dickinson and
Company**
BD, Franklin Lakes, NJ
07417 USA
www.bd.com

Revision Information: No data available.

Further Information: No data available.

Disclaimer: Disclaimer:
The information contained herein has been obtained from various sources and is believed to be correct as of the date issued. However, neither BD nor any of its subsidiaries assumes any liabilities whatsoever for the accuracy or completeness of the information contained herein. Final determination of suitability for a particular use of any material is the sole responsibility of the user. All materials may present unknown hazards and should be used with caution. Although certain hazards are described herein, we cannot guarantee that these are the only hazards that exist. BD provides SDS in electronic form so the information may be more easily accessed. Due to the possibility of errors during transmission, BD makes no representations as to the completeness or accuracy of the information.

Becton, Dickinson and
Company
BD, Franklin Lakes, NJ
07417 USA
www.bd.com

SAFETY DATA SHEET

1. Identification

Product identifier

Product No.:	Product name:	Common name(s), synonym(s)
WP245128JAA	SUPPLEMENT PZA	

Other means of identification

SDS number: 088100177573

Recommended use and restriction on use

Recommended use: Laboratory Chemicals

Restrictions on use: None known.

Manufacturer/Importer/Supplier/Distributor Information

Manufacturer

Company Name: BD Diagnostic Systems
Address: 7 Loveton Circle
21152 Sparks, MD USA
Telephone: 1 410 771 0100 or 1 800 638 8663
Fax:
Contact Person: Tech Services

Emergency telephone number: ChemTrec 1 800 424 9300

2. Hazard(s) identification

Hazard Classification

Not classified

Label Elements

Hazard Symbol: No symbol
Signal Word: No signal word.
Hazard Statement: Not applicable
Precautionary Statements: Not applicable

Other hazards which do not
result in GHS classification: None.

3. Composition/information on ingredients

Becton, Dickinson and
Company
BD, Franklin Lakes, NJ
07417 USA
www.bd.com

Mixtures

Chemical Identity	Common name and synonyms	CAS number	Content in percent (%)*
Ethanol		64-17-5	0.9177%

* All concentrations are percent by weight unless ingredient is a gas. Gas concentrations are in percent by volume.

4. First-aid measures

General information:	Get medical attention if symptoms occur.
Ingestion:	Get medical attention if symptoms occur.
Inhalation:	Provide fresh air, warmth and rest, preferably in comfortable upright sitting position.
Skin Contact:	Wash contact areas with soap and water. Remove contaminated clothing. Launder contaminated clothing before reuse.
Eye contact:	Flush thoroughly with water. If irritation occurs, get medical assistance.

Most important symptoms/effects, acute and delayed

Symptoms: No data available.

Indication of immediate medical attention and special treatment needed

Treatment: No data available.

5. Fire-fighting measures

General Fire Hazards: Extinguish all ignition sources. Avoid sparks, flames, heat and smoking. Ventilate. Use water spray to keep fire-exposed containers cool.

Suitable (and unsuitable) extinguishing media

Suitable extinguishing media: Water spray, fog, CO₂, dry chemical, or alcohol resistant foam.

Unsuitable extinguishing media: None known.

Specific hazards arising from the chemical: None known.

Special protective equipment and precautions for firefighters

Special fire fighting procedures: No unusual fire or explosion hazards noted.

Becton, Dickinson and Company
BD, Franklin Lakes, NJ
07417 USA
www.bd.com

Special protective equipment for fire-fighters:

Firefighters must use standard protective equipment including flame retardant coat, helmet with face shield, gloves, rubber boots, and in enclosed spaces, SCBA.

6. Accidental release measures

Personal precautions, protective equipment and emergency procedures:

No special precautionary health measures should be needed under anticipated conditions of use.

Methods and material for containment and cleaning up:

No specific clean-up procedure noted.

Environmental Precautions:

Avoid release to the environment.

7. Handling and storage

Precautions for safe handling:

When using do not eat, drink or smoke. Read and follow manufacturer's recommendations. Use personal protective equipment as required.

Conditions for safe storage, including any incompatibilities:

Store in a cool, dry place. Keep container tightly closed.

8. Exposure controls/personal protection

Control Parameters

Occupational Exposure Limits

Chemical Identity	Type	Exposure Limit Values	Source
Ethanol	TWA	1,000 ppm 1,900 mg/m3	US. OSHA Table Z-1-A (29 CFR 1910.1000) (1989)
	TWA	1,000 ppm 1,900 mg/m3	US. Tennessee. OELs. Occupational Exposure Limits, Table Z1A (06 2008)
	AN ESL	1,000 ppb	US. Texas. Effects Screening Levels (Texas Commission on Environmental Quality) (12 2010)
	ST ESL	10,000 ppb	US. Texas. Effects Screening Levels (Texas Commission on Environmental Quality) (12 2010)
	AN ESL	1,880 µg/m3	US. Texas. Effects Screening Levels (Texas Commission on Environmental Quality) (12 2010)
	ST ESL	18,800 µg/m3	US. Texas. Effects Screening Levels (Texas Commission on Environmental Quality) (12 2010)
	TWA PEL	1,000 ppm 1,900 mg/m3	US. California Code of Regulations, Title 8, Section 5155. Airborne Contaminants (08 2010)
	STEL	1,000 ppm	US. ACGIH Threshold Limit Values (12 2010)
	REL	1,000 ppm 1,900 mg/m3	US. NIOSH: Pocket Guide to Chemical Hazards (2005)

Becton, Dickinson and
Company
BD, Franklin Lakes, NJ
07417 USA
www.bd.com

	PEL	1,000 ppm 1,900 mg/m3	US. OSHA Table Z-1 Limits for Air Contaminants (29 CFR 1910.1000) (02 2006)
--	-----	-----------------------	---

Appropriate Engineering Controls

No special requirements under ordinary conditions of use and with adequate ventilation.

Individual protection measures, such as personal protective equipment

General information:	Always observe good personal hygiene measures, such as washing after handling the material and before eating, drinking, and/or smoking. Routinely wash work clothing to remove contaminants. Discard contaminated footwear that cannot be cleaned.
Eye/face protection:	Wear safety glasses with side shields (or goggles).
Skin Protection	
Hand Protection:	Chemical resistant gloves Suitable gloves can be recommended by the glove supplier. Wash hands after contact.
Other:	Wear a lab coat or similar protective clothing.
Respiratory Protection:	Respiratory protection not required.
Hygiene measures:	Observe good industrial hygiene practices.

9. Physical and chemical properties

Appearance

Physical state:	solid
Form:	Solid or Flake
Color:	According to product specification.
Odor:	Characteristic
Odor threshold:	No data available.
pH:	No data available.
Melting point/freezing point:	No data available.
Initial boiling point and boiling range:	No data available.
Flash Point:	Not applicable
Evaporation rate:	No data available.
Flammability (solid, gas):	No data available.
Upper/lower limit on flammability or explosive limits	
Flammability limit - upper (%):	No data available.
Flammability limit - lower (%):	No data available.
Explosive limit - upper (%):	No data available.
Explosive limit - lower (%):	No data available.
Vapor pressure:	No data available.
Vapor density:	No data available.
Relative density:	No data available.

Becton, Dickinson and Company
BD, Franklin Lakes, NJ
07417 USA
www.bd.com

Solubility(ies)

Solubility in water:	Completely Soluble
Solubility (other):	Water.: No data available.
Partition coefficient (n-octanol/water):	No data available.
Auto-ignition temperature:	No data available.
Decomposition temperature:	No data available.
Viscosity:	Not determined.

10. Stability and reactivity

Reactivity:	Stable under normal temperature conditions and recommended use.
Chemical Stability:	Material is stable under normal conditions.
Possibility of hazardous reactions:	Not known.
Conditions to avoid:	Avoid exposure to high temperatures or direct sunlight.
Incompatible Materials:	Strong oxidizers.
Hazardous Decomposition Products:	Not known.

11. Toxicological information

General information: No data on possible toxicity effects have been found.

Information on likely routes of exposure

Ingestion:	No harmful effects expected in amounts likely to be ingested by accident.
Inhalation:	Limited inhalation hazard at normal work temperatures.
Skin Contact:	Negligible irritation to skin at ambient temperatures.
Eye contact:	Do not get in eyes.

Symptoms related to the physical, chemical and toxicological characteristics

Ingestion:	No data available.
Inhalation:	No data available.
Skin Contact:	No data available.
Eye contact:	No data available.

Becton, Dickinson and Company
BD, Franklin Lakes, NJ
07417 USA
www.bd.com

Information on toxicological effects

Acute toxicity (list all possible routes of exposure)

Oral

Product: No data available.

Dermal

Product: No data available.

Inhalation

Product: No data available.

Repeated dose toxicity

Product: No data available.

Specified substance(s):

Ethanol

Based on available data, the classification criteria are not met.
LOAEL (Rat(Female, Male), Inhalation, 7,318 - 7,496 h): 1.3 mg/l Inhalation
Read-across from supporting substance (structural analogue or surrogate),
Weight of Evidence study
NOAEL (Guinea pig, Inhalation, 10.5 Weeks): 3,000 ppm(m) Inhalation
Experimental result, Supporting study
LOAEL (Rat(Male), Inhalation, 1 - 6 Weeks): 13.3 mg/l Inhalation Read-
across from supporting substance (structural analogue or surrogate),
Supporting study
LOAEL (Monkey, Inhalation, 5 - 20 d): 3.99 mg/l Inhalation Read-across
from supporting substance (structural analogue or surrogate), Supporting
study

Skin Corrosion/Irritation

Product: No data available.

Specified substance(s):

Ethanol

in vivo (Rabbit): Not irritant Experimental result, Key study

Serious Eye Damage/Eye Irritation

Product: No data available.

Specified substance(s):

Ethanol

in vivo (Rabbit, 24 - 72 hrs): Not irritating EU

Respiratory or Skin Sensitization

Product: No data available.

Becton, Dickinson and Company
BD, Franklin Lakes, NJ
07417 USA
www.bd.com

Specified substance(s):

Ethanol

Based on available data, the classification criteria are not met.
Skin sensitization:, in vivo (Guinea pig): Non sensitising

Carcinogenicity

Product:

No data available.

Specified substance(s):

Ethanol

Based on available data, the classification criteria are not met.

IARC Monographs on the Evaluation of Carcinogenic Risks to Humans:

No carcinogenic components identified

US. National Toxicology Program (NTP) Report on Carcinogens:

No carcinogenic components identified

US. OSHA Specifically Regulated Substances (29 CFR 1910.1001-1050):

No carcinogenic components identified

Germ Cell Mutagenicity

In vitro

Product:

No data available.

Specified substance(s):

Ethanol

Based on available data, the classification criteria are not met.

In vivo

Product:

No data available.

Specified substance(s):

Ethanol

Based on available data, the classification criteria are not met.

Reproductive toxicity

Product:

No data available.

Specified substance(s):

Ethanol

Based on available data, the classification criteria are not met.

Specific Target Organ Toxicity - Single Exposure

Product:

No data available.

Specified substance(s):

Ethanol

Based on available data, the classification criteria are not met.

Specific Target Organ Toxicity - Repeated Exposure

Product:

No data available.

Specified substance(s):

Ethanol

Based on available data, the classification criteria are not met.

Aspiration Hazard

Product:

No data available.

**Becton, Dickinson and
Company**
BD, Franklin Lakes, NJ
07417 USA
www.bd.com

Other effects: None known.

12. Ecological information

Ecotoxicity:

Acute hazards to the aquatic environment:

Fish

Product: No negative effects on the aquatic environment are known.

Aquatic Invertebrates

Product: No negative effects on the aquatic environment are known.

Chronic hazards to the aquatic environment:

Fish

Product: No negative effects on the aquatic environment are known.

Aquatic Invertebrates

Product: No negative effects on the aquatic environment are known.

Toxicity to Aquatic Plants

Product: No negative effects on the aquatic environment are known.

Persistence and Degradability

Biodegradation

Product: Expected to be readily biodegradable.

BOD/COD Ratio

Product: No data available.

Bioaccumulative potential

Bioconcentration Factor (BCF)

Product: No data available.

Specified substance(s):

Becton, Dickinson and Company
BD, Franklin Lakes, NJ
07417 USA
www.bd.com

Ethanol

Potential to bioaccumulate is low.
Cyprinus carpio, Bioconcentration Factor (BCF): 4.5 Aquatic sediment Read-across from supporting substance (structural analogue or surrogate), Supporting study
Cyprinus carpio, Bioconcentration Factor (BCF): 3 Aquatic sediment Read-across from supporting substance (structural analogue or surrogate), Supporting study
Leuciscus idus, Bioconcentration Factor (BCF): 0.2 Aquatic sediment Read-across from supporting substance (structural analogue or surrogate), Not specified
Cyprinus carpio, Bioconcentration Factor (BCF): 1 Aquatic sediment Read-across from supporting substance (structural analogue or surrogate), Supporting study

Partition Coefficient n-octanol / water (log Kow)

Product: Log Kow: No data available.

Mobility in soil: No data available.

Known or predicted distribution to environmental compartments

Ethanol soil - Very mobile liquid

Other adverse effects: The product is not expected to be hazardous to the environment.

13. Disposal considerations

General information: Dispose of waste and residues in accordance with local authority requirements.

Disposal instructions: No specific disposal method required.

Contaminated Packaging: No data available.

14. Transport information

DOTUN Number:	Not regulated.
UN Proper Shipping Name:	Not regulated.
Transport Hazard Class(es)	
Class:	Not regulated.
Label(s):	Not regulated.
Packing Group:	Not regulated.
Marine Pollutant:	Not regulated.
Limited quantity	Not regulated.
Excepted quantity	Not regulated.
Special precautions for user:	Not regulated.

Becton, Dickinson and Company
BD, Franklin Lakes, NJ
07417 USA
www.bd.com

IMDG

UN Number:	Not regulated.
UN Proper Shipping Name:	Not regulated.
Transport Hazard Class(es)	
Class:	Not regulated.
Subsidiary risk:	Not regulated.
EmS No.:	Not regulated.
Packing Group:	Not regulated.
Environmental Hazards	
Marine Pollutant:	Not regulated.
Special precautions for user:	Not regulated.

IATA

UN Number:	Not regulated.
Proper Shipping Name:	Not regulated.
Transport Hazard Class(es):	
Class:	Not regulated.
Subsidiary risk:	Not regulated.
Packing Group:	Not regulated.
Environmental Hazards	
Marine pollutant:	Not regulated.
Special precautions for user:	Not regulated.

15. Regulatory information

US Federal Regulations

TSCA Section 12(b) Export Notification (40 CFR 707, Subpt. D)

None present or none present in regulated quantities.

US. OSHA Specifically Regulated Substances (29 CFR 1910.1001-1050)

None present or none present in regulated quantities.

CERCLA Hazardous Substance List (40 CFR 302.4):

<u>Chemical Identity</u>	<u>Reportable quantity</u>
Ethanol	100 lbs.

Superfund Amendments and Reauthorization Act of 1986 (SARA)

Hazard categories

Not classified
Not classified

SARA 302 Extremely Hazardous Substance

None present or none present in regulated quantities.

Becton, Dickinson and Company
BD, Franklin Lakes, NJ
07417 USA
www.bd.com

SARA 304 Emergency Release Notification

<u>Chemical Identity</u>	<u>Reportable quantity</u>
Ethanol	100 lbs.

SARA 311/312 Hazardous Chemical

<u>Chemical Identity</u>	<u>Threshold Planning Quantity</u>
Ethanol	10000 lbs

SARA 313 (TRI Reporting)

None present or none present in regulated quantities.

Clean Water Act Section 311 Hazardous Substances (40 CFR 117.3)

None present or none present in regulated quantities.

Clean Air Act (CAA) Section 112(r) Accidental Release Prevention (40 CFR 68.130):

None present or none present in regulated quantities.

US State Regulations

US. California Proposition 65

US. New Jersey Worker and Community Right-to-Know Act

<u>Chemical Identity</u>
Ethanol

US. Massachusetts RTK - Substance List

No ingredient regulated by MA Right-to-Know Law present.

US. Pennsylvania RTK - Hazardous Substances

No ingredient regulated by PA Right-to-Know Law present.

US. Rhode Island RTK

No ingredient regulated by RI Right-to-Know Law present.

16. Other information, including date of preparation or last revision
--

Issue Date: 07/17/2018

Version #: 2.1

Revision Information: No data available.

Further Information: No data available.

**Becton, Dickinson and
Company**
BD, Franklin Lakes, NJ
07417 USA
www.bd.com

Disclaimer:

Disclaimer:

The information contained herein has been obtained from various sources and is believed to be correct as of the date issued. However, neither BD nor any of its subsidiaries assumes any liabilities whatsoever for the accuracy or completeness of the information contained herein. Final determination of suitability for a particular use of any material is the sole responsibility of the user. All materials may present unknown hazards and should be used with caution. Although certain hazards are described herein, we cannot guarantee that these are the only hazards that exist. BD provides SDS in electronic form so the information may be more easily accessed. Due to the possibility of errors during transmission, BD makes no representations as to the completeness or accuracy of the information.

Safety Data Sheet acc. to OSHA HCS

Date Prepared: 05/10/2016

Reviewed On: 04/28/2016


1 Identification

- **Product Identifier:**
- **Product Name:** BACTEC MGIT 960 PZA Tubes
- **Catalog Number:** 245115
- **Application of the substance / the mixture** Laboratory Chemicals
- **Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**
BD Diagnostic Systems
7 Loveton Circle
Sparks, MD 21152
Telephone: (410) 771 - 0100 or (800) 638 – 8663
Email Address: Technical_Services@bd.com
- **Information Department:** Technical Service
- **Emergency telephone number:**
In case of a chemical emergency, spill, fire, exposure, or accident, contact BD Diagnostic Systems (410) 771-0100 or (800)-638-8663, or ChemTrec at (800) 424-9300.

2 Hazard(s) identification

- **Classification of the substance or mixture**
The product is not classified according to the Globally Harmonized System (GHS).

- **Label elements**
- **GHS label elements** Void
- **Hazard pictograms** Void
- **Signal word** Void
- **Hazard statements** Void
- **NFPA ratings (scale 0-4)**

Health = 0
Flammability = 0
Reactivity = 0

- **HMIS ratings (scale 0-4)**

Health = 0
Flammability = 0
Reactivity = 0

- **Other hazards**
- **Results of PBT and vPvB assessment**
- **PBT:** Not applicable.
- **vPvB:** Not applicable.

US

(Contd. on page 2)

Safety Data Sheet acc. to OSHA HCS

Date Prepared: 05/10/2016

Reviewed On: 04/28/2016

Product Name: BACTEC MGIT 960 PZA Tubes

(Contd. of page 1)

3 Composition/information on ingredients

- **Chemical characterization:** Mixture
- **Description:** Mixture consisting of the following components.
- **Dangerous Components:** Void
- **Additional information** Risk phrases refer to section 15.

4 First-aid measures

- **Description of first aid measures**
- **General information** No special measures required.
- **After inhalation** Seek medical treatment in case of complaints.
- **After skin contact** Immediately wash with water and soap and rinse thoroughly.
- **After eye contact**
Rinse opened eye for several minutes under running water. If symptoms persist, consult a doctor.
- **After swallowing** If symptoms persist consult doctor.
- **Information for doctor** Show this product label or this SDS.
- **Most important symptoms and effects, both acute and delayed**
No further relevant information available.
- **Indication of any immediate medical attention and special treatment needed**
No further relevant information available.

5 Fire-fighting measures

- **Extinguishing media**
- **Suitable extinguishing agents**
CO₂, ABC multipurpose dry chemical or water spray. Fight larger fires with water spray or alcohol resistant foam.
- **Special hazards arising from the substance or mixture**
No further relevant information available.
- **Advice for firefighters**
- **Protective equipment:** No special measures required.

6 Accidental release measures

- **Personal precautions, protective equipment and emergency procedures** Not required.
- **Environmental precautions:** Wipe up with damp sponge or mop.
- **Methods and material for containment and cleaning up:** No special measures required.
- **Reference to other sections**
See Section 7 for information on safe handling.
See Section 8 for information on personal protection equipment.

(Contd. on page 3)

Safety Data Sheet acc. to OSHA HCS

Date Prepared: 05/10/2016

Reviewed On: 04/28/2016

Product Name: BACTEC MGIT 960 PZA Tubes

See Section 13 for disposal information.

(Contd. of page 2)

7 Handling and storage

- **Handling**
- **Precautions for safe handling** No special measures required.
- **Information about protection against explosions and fires:**
No special measures required.
- **Conditions for safe storage, including any incompatibilities**
- **Storage**
- **Requirements to be met by storerooms and receptacles:** No special requirements.
- **Information about storage in one common storage facility:**
Store away from oxidizing agents.
- **Further information about storage conditions:**
Store in cool, dry conditions in well sealed containers.
- **Specific end use(s)** No further relevant information available.

8 Exposure controls/personal protection

- **Additional information about design of technical systems:**
No further data; see Section 7.
- **Control parameters**
- **Components with limit values that require monitoring at the workplace:**
The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.
- **Additional information:** The lists that were valid during the creation were used as basis.
- **Exposure controls**
- **Personal Protective Equipment**
- **General protective and hygienic measures**
The usual precautionary measures for handling chemicals should be followed.
- **Breathing equipment:** Not required.
- **Protection of hands:**

Chemical resistant gloves (i.e. nitrile, or equivalent).

- **Eye protection:** Safety glasses
- **Body protection:** Protective work clothing (lab coat).

US

(Contd. on page 4)

Safety Data Sheet

acc. to OSHA HCS

Date Prepared: 05/10/2016

Reviewed On: 04/28/2016

Product Name: BACTEC MGIT 960 PZA Tubes

(Contd. of page 3)

9 Physical and chemical properties

· Information on basic physical and chemical properties	
· General Information	
· Appearance:	
Form:	Liquid
Color:	Clear
· Odor:	Characteristic
· Odor threshold:	Not determined.
· pH-value at 20 °C (68 °F):	5.8
· Change in condition	Undetermined
· Melting point/Melting range:	Not determined
· Boiling point/Boiling range:	Not determined
· Flash point:	Not applicable
· Flammability (solid, gaseous)	Not applicable.
· Ignition temperature:	
Decomposition temperature:	Not determined.
· Auto igniting:	Product is not self igniting.
· Danger of explosion:	Product does not present an explosion hazard.
· Explosion limits:	
Lower:	Not determined.
Upper:	Not determined.
· Vapor pressure:	Not determined.
· Density:	Not determined
· Relative density	Not determined.
· Vapor density	Not determined.
· Evaporation rate	Not determined.
· Solubility in / Miscibility with	
Water:	Soluble
· Partition coefficient (n-octanol/water):	
Not determined.	
· Viscosity:	
dynamic:	Not determined.
kinematic:	Not determined.
· Other information	No further relevant information available.

US

(Contd. on page 5)

Safety Data Sheet acc. to OSHA HCS

Date Prepared: 05/10/2016

Reviewed On: 04/28/2016

Product Name: BACTEC MGIT 960 PZA Tubes

(Contd. of page 4)

10 Stability and reactivity

- **Reactivity** No further relevant information available.
- **Chemical stability**
- **Thermal decomposition / conditions to be avoided:**
No decomposition if used according to specifications.
- **Possibility of hazardous reactions** No dangerous reactions known
- **Conditions to avoid** No further relevant information available.
- **Incompatible materials:** Incompatible material: strong oxidizers.
- **Hazardous decomposition products:** No dangerous decomposition products known.

11 Toxicological information

- **Information on toxicological effects**
- **Acute toxicity:**
- **Primary irritant effect:**
 - **on the skin:** No irritating effect.
 - **on the eye:** No irritating effect.
- **Sensitization:** No sensitizing effects known.
- **Additional toxicological information:**
The product shows the following dangers according to internally approved calculation methods for preparations:
The product is not subject to OSHA classification according to internally approved calculation methods for preparations.
When used and handled according to specifications, the product does not have any harmful effects according to our experience and the information provided to us.

- **Carcinogenic categories**

- **IARC (International Agency for Research on Cancer)**

None of the ingredients is listed.

- **NTP (National Toxicology Program)**

None of the ingredients is listed.

- **OSHA-Ca (Occupational Safety & Health Administration)**

None of the ingredients is listed.

12 Ecological information

- **Toxicity**
- **Aquatic toxicity:** No further relevant information available.
- **Persistence and degradability** No further relevant information available.
- **Behavior in environmental systems:**
- **Bioaccumulative potential** No further relevant information available.

(Contd. on page 6)

Safety Data Sheet acc. to OSHA HCS

Date Prepared: 05/10/2016

Reviewed On: 04/28/2016

Product Name: BACTEC MGIT 960 PZA Tubes

(Contd. of page 5)

- **Mobility in soil** No further relevant information available.
- **Ecotoxicological effects:**
- **Other information:**
The ecological effects have not been thoroughly investigated, but currently none have been identified.
- **Results of PBT and vPvB assessment**
- **PBT:** Not applicable.
- **vPvB:** Not applicable.
- **Other adverse effects** No further relevant information available.

13 Disposal considerations

- **Waste treatment methods**
- **Recommendation**
Smaller quantities can be disposed of with solid waste.
This product is not considered a RCRA hazardous waste.
Dispose of material in accordance with federal (40 CFR 261.3), state and local requirements.
- **Uncleaned packagings:**
- **Recommendation:** Disposal must be made according to state and federal regulations.
- **Recommended cleansing agent:** Water, if necessary with cleansing agents.

14 Transport information

- | | |
|---------------------------------------|-----------------|
| · UN-Number | |
| · DOT, ADN, IMDG, IATA | Void |
| · UN proper shipping name | |
| · DOT, ADN, IMDG, IATA | Void |
| · Transport hazard class(es) | |
| · DOT, ADN, IMDG, IATA | |
| · Class | Void |
| · Packing group | |
| · DOT, IMDG, IATA | Void |
| · Environmental hazards: | Not applicable. |
| · Special precautions for user | Not applicable. |

(Contd. on page 7)

US

Safety Data Sheet acc. to OSHA HCS

Date Prepared: 05/10/2016

Reviewed On: 04/28/2016

Product Name: BACTEC MGIT 960 PZA Tubes

(Contd. of page 6)

- | | |
|---|--|
| · Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code | Not applicable. |
| · Transport/Additional information: | Not dangerous according to the above specifications. |
| · UN "Model Regulation": | Void |

15 Regulatory information

- **Safety, health and environmental regulations/legislation specific for the substance or mixture**

- **SARA Section 355 (extremely hazardous substances)**

None of the ingredients is listed.

- **SARA Section 313 (specific toxic chemical listings)**

None of the ingredients is listed.

- **TSCA (Toxic Substances Control Act)**

All ingredients are listed.

- **California Proposition 65 - Chemicals known to cause cancer**

None of the ingredients is listed.

- **California Proposition 65 - Chemicals known to cause reproductive toxicity for females:**

None of the ingredients is listed.

- **California Proposition 65 - Chemicals known to cause reproductive toxicity for males:**

None of the ingredients is listed.

- **California Proposition 65 - Chemicals known to cause developmental toxicity:**

None of the ingredients is listed.

- **Carcinogenic categories**

- **TLV (Threshold Limit Value established by ACGIH)**

None of the ingredients is listed.

- **GHS label elements** Void

- **Hazard pictograms** Void

- **Signal word** Void

- **Hazard statements** Void

- **Chemical safety assessment:** A Chemical Safety Assessment has not been carried out.

16 Other information

To the best of our knowledge, the information contained herein is accurate. However, neither Becton, Dickinson and Company or any of its subsidiaries assumes any liabilities whatsoever

(Contd. on page 8)

Safety Data Sheet acc. to OSHA HCS

Date Prepared: 05/10/2016

Reviewed On: 04/28/2016

Product Name: BACTEC MGIT 960 PZA Tubes

(Contd. of page 7)

for the accuracy or completeness of the information contained herein. Final determination of suitability of any material is the sole responsibility of the user. All materials may present unknown hazards and should be used with caution. Although certain hazards are described herein, we can not guarantee that these are the only hazards that exist.

· Department issuing SDS:

Environmental, Health & Safety

Created by Michael J. Spinazzola

· Contact: Technical Service Representative**· Date of preparation / last revision** 05/10/2016 / -**· Abbreviations and acronyms:**

IMDG: International Maritime Code for Dangerous Goods

DOT: US Department of Transportation

IATA: International Air Transport Association

ACGIH: American Conference of Governmental Industrial Hygienists

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

NFPA: National Fire Protection Association (USA)

HMIS: Hazardous Materials Identification System (USA)

LD50: Lethal dose, 50 percent

PBT: Persistent, Bioaccumulative and Toxic

vPvB: very Persistent and very Bioaccumulative

OSHA: Occupational Safety & Health

TLV: Threshold Limit Value

PEL: Permissible Exposure Limit

REL: Recommended Exposure Limit

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)

LC50: Lethal concentration, 50 percent

BEI: Biological Exposure Limit

US