Demineralized Water

EN

REF T3339 / T3339-10

These Instructions for Use (IFU) are intended to be read in conjunction with the IFU for relevant Thermo Scientific™ Sensititre™ plates.

Intended Use

A medium for the dilution of microorganisms isolated from clinical samples for the inoculation of the Sensititre™ plates as part of the Thermo Scientific™ Sensititre™ System for in vitro antimicrobial susceptibility testing and qualitative identification.

This device is designed to be used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having microbial infections.

It is for professional use only; is not automated, nor is it a companion diagnostic.

Summary and Explanation

Demineralized water is an established accessory to be used in conjunction with Sensititre™ plates, to test the susceptibility of organisms and are used as aids to diagnosis within a clinical environment.

Principle of Method

Demineralized water is an accessory and is used as a diluent as part of the Sensititre™ system.

Typical Formula

Purified Water

Physical Appearance

Appearance	Clear, colourless liquid
Conductivity	<50 µS
Fill Volume	5.0 - 5.4 mL
pH	Not Applicable

Materials Provided

T3339 / T3339-10

The pack contains 5mL of demineralized water per tube. The quantity of T3339 is 100 x 5 mL and T3339-10 is 10 x 5 mL.

Each tube should only be used once.

Materials Required but Not Supplied

1. Relevant Sensititre plates
2. 0.5 Polymer McFarland Standard
3. Sterile inoculum reservoir
4. Quality control strains
5. Vortex mixer
6. Current CLSI, EUCAST or local guideline documents; plate dependent
7. Inoculation loops
8. Collection containers
9. Incubators

Storage

- Store product in its original packaging at temperature between 15-25°C until used.
- Keep container tightly closed.
- The product may be used until the expiry date stated on the label.
- Allow product to equilibrate to room temperature before use.
- Do not incubate prior to use.

Warnings and Precautions

- For *in vitro* diagnostic use only
- For professional use only.
- Inspect the product packaging before first use.
- Do not use the product if there is any visible damage to the packaging or tubes.
- Do not use the product beyond the stated expiry date.
- Do not use the product if signs of contamination are present.
- Do not use the product if the colour has changed or there are other signs of deterioration.
- Follow local recommended procedures for disposal of used, contaminated, infectious or unused products. It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state, and local applicable regulations.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product (www.thermofisher.com).

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

Specimen Collection, Handling and Storage



There are no specimen collection and/or transport materials provided with this device.

Specimens should be collected and handled following the guidelines provided by CLSI, EUCAST or local standard guidelines.

Procedure

Refer to the IFU of the relevant Sensititre™ plate.

Symbol Legend

Symbol/Label	Meaning
	Manufacturer
	In Vitro Diagnostic Medical Device
	Temperature limit
	Batch Code
	Catalogue Number
	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Use-by date
	Do not use if package is damaged and Consult instructions for use
	Authorized representative in the European Community/ European Union
	Unique device identifier
	USA: Caution: Federal law restricts this device to sale by or on order of a Physician
	European Conformity Mark
	UK Conformity Mark
	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union
	Made in the United States of America

© 2025 Thermo Fisher Scientific Inc. All rights reserved. ATCC® is a trademark of ATCC. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. This information is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others.

Remel Incorporated 12076 Santa Fe Drive,
Lenexa, Kansas, 66215, USA.

For technical assistance please contact your local distributor.

Revision Information

Version	Date of issue and modifications
2.0	14 August 2025 Minor formatting changes

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