

PRODUCT CATALOGUE

2022-2023

TABLE OF CONTENTS

| Introduction | | 1 |
|-------------------------------------|-------|----|
| Vision, Objectives, Mission | | 1 |
| Sites, Our Products | | 2 |
| Standards, Our Markets | | 3 |
| MILESTONES | | 4 |
| COVID-19 TESTING KITS | | 5 |
| VIRAL TRANSPORT MEDIUM KITS | | 6 |
| LATEX KITS | ••••• | 7 |
| TURBIDIMETRIC LATEX KITS | | 8 |
| SYPHILIS KITS | ••••• | 8 |
| FEBRILE ANTIGENS | | 9 |
| BLOOD GROUPING REAGENTS | | 10 |
| HEMATOLOGY TESTS | | 11 |
| SICKLE CELL KITS | ••• | 11 |
| INFECTIOUS DISEASE RAPID TESTS | | 12 |
| ANTIBODY TESTING | | 12 |
| ANTIGEN TESTING | | 13 |
| URINE REAGENT STRIPS | | 14 |
| FERTILITY RAPID TESTS | | 15 |
| KIDNEY FUNCTION RAPID TESTS | | 15 |
| INFLAMMATION AND CANCER MARKERS | | 16 |
| CARDIAC MARKERS RAPID TESTS | | 16 |
| DOA RAPID TESTS | | 17 |
| CLINICAL CHEMISTRY KITS | | 19 |
| ALCOHOL TESTS | | 20 |
| STAINS FOR HISTOLOGY & MICROBIOLOGY | | 21 |
| ANTIBIOTIC SENSITIVITY DISCS | | 22 |
| ELISA KITS | ••••• | 23 |
| HORMONE ELISA KITS | | 23 |
| TORCH ELISA KITS | | 24 |
| INFECTIOUS DISEASES ELISA KITS | | 24 |
| OTHER ELISA KITS | | 24 |
| HOME TESTS | •••• | 25 |
| BLOOD GLUCOSE MONITORING SYSTEMS | | 26 |
| CERTIFICATES | | 27 |
| CON TACT US | | 28 |
| INTERNATIONAL PRESENCE | | 29 |

INTRODUCTION

Atlas Medical GmbH was established in 1996 as a manufacturer and supplier of quality Diagnostic Reagents and Kits. Our products are sold in over 80 countries worldwide.

The company is located at the Cambridge Science Park, Cambridge, UK. In addition to the UK site, the company has offices in Germany and Turkey as well as two purpose-built modern facilities in both Jordan and Malaysia. We take quality assurance very seriously and strive to produce goods to the highest standards known in the industry, including, ISO13485 & CE mark and US FDA standards. Our R&D team constantly develops and innovates novel products that significantly contribute to the advancement of the Diagnostic Industry.



To be a major provider of quality medical diag nostic products to local, regional and interna tional markets.



Mission

Our mission is to develop, produce and pro vide our customers with high quality products and excellent customer services through deep understanding of customers' needs and per ception, recruitment of high caliber profes sionals & technicians, adopting strict quality assurance and control procedures and embracing new scientific advancements in the medical lab diagnostic field.



High and Consistent Quality

Satisfied Customer



Continuous Improvement & Innovation

Care for the Environment & Working Conditions Atlas Medical GmbH Ludwig-Erhard-Ring 3 15827 Blankenfelde-Mahlow, Berlin, Germany Tel: +49 - 33708 – 3550 30 Email: amug@atlas-medical.com



Atlas Medical UK. William James House, Cowley Road, Cambridge, CB4 0WX, UK. Tel.: +44 (0) 1223 858 910 Fax: +44 (0) 1223 858 524 Email: info@atlas-medical.com



Atlas Medical Jordan Sahab Free Zone Area, Amman 11512 Jordan Tel: +962 6 4026468 Fax +962 6 4022588 Email : info@atlas-medical.com



Atlas Medlab Aktarhüssam Mh.1. Değirmen Sk. Değirmen İşhani No :9, D: 405 . Osmangazi / Bursa/Turkey Tel: +90 552 366 21 71 Email:info@atlas-medical.com



Atlas Medical enjoys a good presence in many international markets. We take pride in our export activities through our dedicated export department. We actively participate in major industry-related exhibitions seeking keen representatives around the globe to sell and distribute our products in their respective countries. We are internationally repre sented in more than 80 of countries spanning in five continents: Europe, North America, South America, Africa and Asia. Our efforts will continue to increase our representation to include most markets around the globe.





Our products are manufactured in accordance to the standards as set in the European In-Vitro Diagnostic Directive 98/79/EC. This has led to the successful attainment of Annex IV Full Quality Assurance Certification and the declaration of conformity for CE marking purposes for many of our IVD products, either self-declared or through our Notified Body LNE/G-MED



To complete the quality assurance scheme the company has put in place a robust Quality Management and Enhancement System that has concluded in the successful attainment of ISO13485: 2016 certificate

the company also adheres to the US-FDA regulations and had already FDA-cleared few products for the US market. Our products are registered in numerous countries.

FDA

MILESTONES



Overview

Atlas medical is introducing COVID-19 Real Time RT-PCR kit for the amplification and detection of the viral genetic material in patient specimen ,In addition Atlas Medical had introduced three new kits using the ELISA technique to detect the antibody response to COVID-19 infection. Detecting antibodies to SARS-CoV-2 virus could tell if a patient has been infected with COVID-19, either currently or in the past.

| | | \sim |
|--------------|--------------------------------|----------------|
| Item Code | Item Description | Sizes |
| 8.14.45.0096 | COVID-19 Real Time RT-PCR Test | 96 Tests /Kits |
| 8.14.45.0096 | COVID-19 Total Ab Elisa Kit | 96 Tests /Kits |
| 8.14.46.0096 | COVID-19 IgM Elisa Kit | 96 Tests /Kits |
| 8.14.47.0096 | COVID-19 S1-RBD IgG Elisa Kit | 96 Tests /Kits |

Features

- Atlas RT-PCR testing kits are fairly quick, sensitive ,reliable and can detect current infections of disease.
- Atlas ELISA kits are based on a simple and high sensitive laboratory technique, results can typically be produced within 1 to 2 hours from the moment of collecting the nasal swab sample.



COVID-19 RAPID TESTING KITS

Overview

Atlas Medical offers COVID-19 Rapid Test as a quick screening tool for the detection of the presence of SARS-CoV-2 virus in Nasopharyngeal swabs,"COVID-19 Antigen testing kit", in addition Atlas medical offers" COVID-19 IgG/IgM Antibody Testing kit " as a screening tool for the human body response to the infection with the virus.

Features

COVID-19 rapid test kits are based on Lateral Flow Immuno-Chromatographic Assay.

Reliable ,easy to use with a short testing time of 10-15 minutes per each sample.

The kits are conveniently packed in different sizes of 20,25 or 100 tests per kit including the necessary test accessories to perform the assay.

| | | <u> </u> |
|----------------|---|--------------|
| Item Code | Item Description | Sizes |
| 8.66.00.0.0001 | COVID-19 IgM /IgG Test Cassette, Whole Blood/ | Bulk |
| 8.66.00.0.0020 | Serum/Plasma ,Individually Pouched | 20 Tests/Box |
| 8.66.01.0.0001 | COVID-19 Antigen Test Cassette, Nasal Swab, | Bulk |
| 8.66.01.0.0020 | Individually Pouched | 20 Tests/Box |
| 8.66.02.0.0001 | COVID-19 Combo Antigen & Influenza, A+B Test | Bulk |
| 8.66.02.0.0020 | Individually Pouched | 20 Tests/Box |
| 8.66.04.0.0020 | COVID-19 Neutralizing Antibody Rapid Test Cassette ,Whole Blood /Serum/Plasma, Individually Pouched | 20 Tests/Box |

Overview

As a result of the current worldwide crises such as the COVID-19 pandemic and other worldwide health pandemic that are caused from malaria, Influenza and chlamydia, Atlas medical had come up with Viral transport media, The Media is used for facilitating the testing procedures by preserving the sample through the collection and transport process of the clinical samples containing viruses; including SARS-CoV-2 (COVID-19) and other viruses, in active form from collection site to the testing laboratory.

Œ

Features

- Atlas Medical offers VTM that allows the safe transfer of viruses for further research, including diagnostic tests, and molecular biology techniques
- Atlas VTM can come either as a biological format "activated product (VTM) " or in a chemical format "In Activated product (IVTM)".
- Atlas VTM maintain the viral structure and activity over a wide temperature range and suppress "
- ∂ Atlas viral transport medium is stable at room temperature



arrho Atlas VTM kits are conveniently packed in different sizes of 50 ,100 tube per kit with flocked swabs

| | | \smile |
|----------------|--|--------------|
| Item Code | Item Description | Sizes |
| 8.64.02.0.0050 | Viral Transport Medium (VTM) with Flocked , | 50 Tube/Kit |
| 8.64.02.0.0100 | Swab ,Individually Pouched | 100 Tube/Kit |
| 8.64.03.0.0050 | Viral Transport Medium (In Activating) with | 50 Tube/Kit |
| 8.64.03.0.0100 | Flocked Swab ,(3 ml /vial) | 100 Tube/Kit |



LATEX KITS



Overview

Latex kits offer a quick and simple assay to diagnose a range of pathogens and medical conditions. The assay is based on an immunological reaction between the detected analyte in the sample and its corresponding antibody or antigen already coated on latex particles.

| Item Code | Item Description | Sizes |
|----------------------------------|---------------------------|-----------------------|
| 8.00.00.0.0050 8.00.00.0.0100 | CRP Latex Kit | 50 Tests 100 Tests |
| 8.00.01.0.0050 8.00.01.0.0100 | CRP Latex Kit with Buffer | 50 Tests 100 Tests |
| 8.00.02.0.0050 8.00.02.0.0100 | ASO Latex Kit | 50 Tests 100 Tests |
| 8.00.03.0.0050 8.00.03.0.0100 | ASO Latex Kit with Buffer | 50 Tests 100 Tests |
| 8.00.04.0.0050 8.00.04.0.0100 | RF Latex Kit | 50 Tests 100 Tests |
| 8.00.05.0.0050 8.00.05.0.0100 | RF Latex Kit with Buffer | 50 Tests 100 Tests |
| 8.00.07.0.0050 8.00.07.0.0100 | hCG Latex Kit | 50 Tests 100 Tests |
| 8.00.11.0.0050 8.00.11.0.0100 | SLE Latex Kit | 50 Tests 100 Tests |
| 8.00.16.0.0050 8.00.16.0.0100 | Rota Virus Latex Kit | 50 Tests 100 Tests |
| 8.00.17.0.0050 8.00.17.0.0100 | D-Dimer Latex Kit | 50 Tests 100 Tests |
| 8.00.21.0.0050 8.00.21.0.0100 | Waaler Rose Kit | 50 Tests 100 Tests |
| 8.00.08.0.0050 8.00.08.0.0100 | IM Latex Kit | 50 Tests 100 Tests |
| 8.00.12.0.0050 8.00.12.0.0100 | Staphylococcus Latex Kit | 50 Tests 100 Tests |
| 8.00.13.0.0300 | Streptococcus Latex Kit | 50 Tests |



Œ

Item Code Item Description Sizes 8.00.09.0.0050 50 Tests Toxo Latex Kit 8.00.09.0.0100 100 Tests 8.00.10.0.0050 50 Tests Toxo Latex Kit with Buffer 8.00.10.0.0100 100 Tests 8.00.14.0.0100 Rubella Latex Kit 100 Tests

Features

- They cover a selection of routine tests in serology and microbiology.
- They are conveniently packed in sizes of 50 or 100 tests and includes all the necessary reagents, controls, stirrers and slides to conduct the test.
- Affordable, easy to use, dependable and offer a clear and visible agglutination for doubt-free results.
- Some Latex Kits come with a Buffer.

TURBIDIMETRIC LATEX KITS

Overview

The turbidimetric assay is based on the agglutination reaction between latex particles coated with antibody and the antigen in solution. The intended use for Turbilatex products is to detect and quantify the antigen present in human serum or plasma samples.



Features

 Atlas Medical offers a dynamic range of Turbidimetric Latex Kits which are conveniently packed in sizes of 50, 100 and 250 tests and include all the necessary accessories.

| Item Code | Item Description | Sizes |
|--|---------------------------------|------------------------------------|
| 8.44.00.0.0050 8.44.00.0.0250 | RF Turbidimetric Latex Kit | 50 Tests 250 Tests |
| 8.44.01.0.0050 8.44.01.0.0100 8.44.01.0.0250 | CRP Turbidimetric Latex Kit | 50 Tests 100 Tests 250 Tests |
| 8.44.02.0.0050 8.44.02.0.0100 8.44.02.0.0250 | ASO Turbidimetric Latex Kit | 50 Tests 100 Tests 250 Tests |
| 8.44.03.0.0050 8.44.03.0.0100 8.44.03.0.0250 | D-Dimer Turbidimetric Latex Kit | 50 Tests 100 Tests 250 Tests |
| 8.44.04.0.0050 | Microalbumine Turbilatex | 50 Tests |
| 8.44.05.0.0050 | Ferritin Turbilatex | 50 Tests |
| 8.44.06.0.0050 | Transferrin Turbilatex (TRF) | 50 Tests |

SYPHILIS KITS

Overview

Atlas Medical offers a number of assays to detect Syphilis that include: TPHA kits which are used for the detection of antibodies to Treponema pallidum in human Serum or plasma using micro haemagglutination; VDRL and RPR kits which are based on non-Treponemal floccuation to detect reagin antibodies in serum or plasma.

| | | \smile |
|--|------------------------|--------------------------------------|
| Item Code | Item Description | Sizes |
| 8.00.18.0.0100 8.00.18.0.0250 8.00.18.0.0500 | RPR Carbon Antigen Kit | 100 Tests 250 Tests 500 Tests |
| 8.00.19.0.0050 8.00.19.0.0100 8.00.19.0.0200 | TPHA Kit | 50 Tests 100 Tests 200 Tests |
| 8.00.20.0.0250 8.00.20.0.0500 8.00.20.0.2500 | VDRL Kit | 250 Tests 500 Tests 2500 Tests |
| 8.00.20.1.0250 8.00.20.1.2500 | VDRL Kit with controls | 250 Tests 2500 Tests |

Features

- Easy to use, affordable and conveniently packed in different sizes to suit all needs.
- They include all the necessary reagents/devices, controls, stirrers and slides to conduct the test



 $\mathbf{\epsilon}$

FEBRILE ANTIGENS



 (ϵ)

Overview

Febrile antigen kits are based on bacterial suspensions that agglutinate in the presence of antibodies formed in human infection by certain fever-causing microbial agents. In positive samples, the agglutination is macroscopically visible at certain antibody levels in serum. These antigen reagents are used for the qualitative and semi quantitative febrile screening purposes.

Features

- Atlas Medical Febrile Antigen kits contain various types of antigens for Brucella, Proteus, Salmonella typhi and paratyphi, and their controls as needed.
- Atlas Medical Febrile Antigen kits are competitively priced and easy to use, and give clear results within 2 minutes

| | | \smile |
|----------------------------------|---|------------------------|
| Item Code | Item Description | Sizes |
| 8.01.17.0.0050 | Febrile Antigen Set (10 Antigens: Salmonella OA, OB, OC, OD, HA, HB, HC, HD, Brucella abortus, melitensis) | 10x5 ml |
| 8.01.17.1.0050 | Febrile Antigen Set (10 Antigens: Salmonella OA, OB, OC, OD,HA, HB, HC, HD, Brucella abortus, melitensis) with 3x1.0ml Controls | 10x5 ml |
| 8.01.18.0.0030 | Salmonella Antigen Set, Widal Kit (6 Antigens: OA, OB, OD,HA, HB, HD) | 6x5 ml |
| 8.01.18.1.0030 | Salmonella Antigen Set, Widal Kit (6 Antigens: OA, OB, OD, HA, HB, HD) with 2x1.0ml Controls | 6x5 ml |
| 8.01.19.0.0001 8.01.19.0.0005 | Febrile Antigens Positive Control | 1 ml/vial 5 ml/vial |
| 8.01.20.0.0001 8.01.20.0.0005 | Febrile Antigen Negative Control | 1 ml/vial 5 ml/vial |



| ltem Code | Item Description | Sizes |
|----------------|--|---------------------|
| 8 01 00 0 0005 | | 5ml/vial |
| 8.01.00.0.0050 | Brucella Rose Bengal Kit | 50 Tests |
| 8.01.00.0.0100 | | 100 Tests |
| 8.01.01.0.0005 | Salmonella OA Reagent | 5 ml/vial |
| 8.01.01.1.0040 | Sumonella Ortheagent | 8x5 ml |
| 8.01.01.0.0050 | | 10x5 ml |
| 8.01.02.0.0005 | Salmonella OB Reagent | 5 ml/vial |
| 8.01.02.0.0050 | | 10x5 ml |
| 8.01.03.0.0005 | Salmonella OC Reagent | 5 ml/vial |
| 8.01.03.1.0040 | | 8x5 ml |
| 8.01.03.0.0050 | | 10x5 ml |
| 8.01.04.0.0005 | Salmonella OD Reagent | 5 mi/viai 8x5 ml |
| 8.01.04.0.0050 | | 10x5 ml |
| 8.01.05.0.0005 | Salmonella HA Reagent | 5 ml/vial |
| 8.01.05.1.0040 | | 8x5 ml |
| 8.01.05.0.0050 | | 10x5 ml |
| 8.01.06.0.0005 | Salmonella HB Reagent | 5 ml/vial 8x5 ml |
| 8.01.06.0.0050 | | 10x5 ml |
| 8.01.07.0.0005 | Salmonella HC Reagent | 5 ml/vial |
| 8.01.07.1.0040 | - | 8x5 ml |
| 8.01.08.0.0005 | Salmonolla HD Reagent | 5 ml/vial |
| 8.01.08.1.0040 | Samonena no keagent | 8x5 ml |
| 8.01.08.0.0050 | | 10x5 ml |
| 8.01.10.0.0005 | Brucella Abortus Reagent | 5 ml/vial |
| 8.01.10.1.0040 | | 8x5 ml |
| 8.01.11.0.0005 | Brucella Melitensis Reagent | 5 ml/vial |
| 8.01.11.1.0040 | Drotous OV2 Paggant | 5 ml/vial |
| 8.01.12.1.0040 | Floteus OX2 Reagent | 8x5 ml |
| 8 01 13 0 0005 | Proteus OX19 Reagent | 5 ml/vial |
| 8.01.13.1.0040 | The agent | 8x5 ml |
| 8.01.14.0.0005 | Proteus OXK Reagent | 5 ml/vial |
| 8.01.14.1.0040 | The source of th | 8x5 ml |
| 8.01.15.0.0010 | Brucella Antigen Kit (Brucella melitensis, Brucella abortus) | 2 vials/Box |
| 8.01.15.1.0010 | Brucella Antigen Kit with Controls (Brucella melitensis, Brucella abortus, 2x1.0 ml Controls) | 2 vials/Box |
| 8.01.15.2.0010 | Brucella Antigen Kit with Controls, (5ml Brucella melitensis, 5ml Brucella abortus, 2x0.5 ml Controls) | 2 vials/Box |
| 8.01.16.0.0040 | Salmonella Antigen Set (8 Antigens: OA, OB, OC, OD, HA, HB, HC, HD) | 8x5 ml |
| 8.01.16.1.0040 | Salmonella Antigen Set (8 Antigens: OA, OB, OC, OD, HA, HB, HC, HD) with 2x1.0 ml Controls | 8x5 ml |
| 8.01.16.2.0040 | Salmonella Antigen set (8 Antigens: OA, OB, OC, OD, HA, HB, HC, HD) with 2x0.5ml Controls | 8x5 ml |

Œ

BLOOD GROUPING REAGENTS

Overview

Blood Grouping reagents are used for the identification of blood types. The test procedure is based on the agglutination principle, where red cells possessing the typing antigen agglutinate in the presence of the corresponding antibody in the testing reagent indicating the presence of the tested antigen. No agglutination indicates the absence of the tested antigen.

Features

- Atlas Medical ABO reagents are prepared from In-Vitro culture supernatants of hybridized immunoglobulin-secreting mouse cell lines.
- The reagents are formulated and optimized for use in tube and slide methods.
- Atlas Medical provides high quality blood grouping reagents that are accurate, easy to use, competitively priced, and conveniently packed in different sizes and options.

| | | \sim |
|----------------|--|------------|
| Item Code | Item Description | Sizes |
| 8.02.14.0.0010 | Anti-D Monoclonal (IgM), Clone 1, 10ml/vial | 10 ml/vial |
| 8.02.15.0.0010 | Anti-D Monoclonal (IgM), Clone 2, 10ml/vial | 10 ml/vial |
| 8.02.16.0.0005 | Anti-A1, Lectin (Dolichosbiflorus), 5ml/vial | 5 ml/vial |
| 8.02.17.0.0005 | Anti-H, Lectin (Ulexeuropaeus), 5ml/vial | 5 ml/vial |
| 8.02.18.0.0005 | Anti-C Monoclonal, 5ml/vial | 5 ml/vial |
| 8.02.19.0.0005 | Anti-c Monoclonal, 5ml/vial | 5 ml/vial |
| 8.02.20.0.0005 | Anti-E Monoclonal, 5ml/vial | 5 ml/vial |
| 8.02.21.0.0005 | Anti-e Monoclonal, 5ml/vial | 5 ml/vial |
| 8.02.22.0.0005 | Anti-C+D+E Monoclonal, 5ml/vial | 5 ml/vial |
| 8.02.27.0.0002 | Anti-Fya, Human, 2ml/vial | 2 ml/vial |
| 8.02.28.0.0002 | Anti-Fyb, Human, 2ml/vial | 2 ml/vial |
| 8.02.29.0.0002 | Anti-k, Human, 2ml/vial | 2 ml/vial |
| 8.02.30.0.0002 | Anti-Kpa, Human, 2ml/vial | 2 ml/vial |
| 8.02.31.0.0002 | Anti-Kpb, Human, 2ml/vial | 2 ml/vial |
| 8.02.32.0.0002 | Anti-Jka, Human, 2ml/vial | 2 ml/vial |
| 8.02.35.0.0002 | Anti-Lub, Human, 2 ml/vial | 2 ml/vial |
| 8.02.36.0.0005 | Anti-K Monoclonal, 5ml/vial | 5 ml/vial |
| 8.02.54.0.0002 | Anti-Cw,2 ml/vial | 2 ml/vial |





Œ

| | | <u> </u> |
|----------------------------------|---|------------------------|
| Item Code | Item Description | Sizes |
| 8.02.00.0.0010 8.02.00.1.0100 | Anti-A Monoclonal reagent (titer: 1/512) | 10 ml/vial 10x10 ml |
| 8.02.01.0.0010 8.02.01.1.0100 | Anti-B Monoclonal Reagent (Titer: 1/512) | 10 ml/vial 10x10 ml |
| 8.02.02.0.0010 8.02.02.1.0100 | Anti-AB Monoclonal Reagent (Titer: 1/512) | 10 ml/vial 10x10 ml |
| 8.02.03.0.0010 8.02.03.1.0100 | Anti-D lgG/lgM Blend Reagent (Titer: 1/128) | 10 ml/vial 10x10 ml |
| 8.02.04.0.0010 8.02.04.0.0100 | Anti-A Monoclonal Reagent (Titer: 1/256 | 10 ml/vial 10x10 ml |
| 8.02.05.0.0010 8.02.05.0.0100 | Anti-B Monoclonal Reagent (Titer: 1/256) | 10 ml/vial 10x10 ml |
| 8.02.06.0.0010 8.02.06.1.0100 | Anti-AB Monoclonal Reagent (Titer: 1/256) | 10 ml/vial 10x10 ml |
| 8.02.07.0.0010 8.02.07.1.0100 | Anti-D lgG/lgM Blend Reagent (Titer: 1/64) | 10 ml/vial 10x10 ml |
| 8.02.08.0.0010 8.02.08.1.0100 | Bovine Albumin 22% | 10 ml/vial 10x10 ml |
| 8.02.09.0.0010 8.02.09.1.0100 | Bovine Albumin 30% | 10 ml/vial 10x10 ml |
| 8.02.10.0.0010 8.02.10.1.0100 | Anti-Human Globulin (Green) (Titer 1/512) | 10 ml/vial 10x10 ml |
| 8.02.11.0.0010 8.02.11.1.0100 | Anti-Human Globulin (Green) (Titer 1/256) | 10 ml/vial 10x10 ml |
| 8.02.47.0.0030 | ABO Set (Anti-A (1/512), Anti-B (1/512), Anti-D (1/128)) | 3x10 ml |
| 8.02.47.1.0030 | ABO Set (Anti-A (1/265), Anti-B (1/265), Anti-D (1/64)) | 3x10 ml |
| 8.02.49.0.0040 | ABO Set (Anti-A (1/256), Anti-B (1/256), Anti-AB (1/256), Anti-D (1/64)) | 4x10 ml |
| 8.02.53.0.0040 | ABO Set (Anti-A (1/512), Anti-B (1/512), Anti-AB (1/512) Anti-D (1/128)) | 4x10 ml |
| 8.02.05.6.0030 | ABO Set (Anti-A (1/256), Anti-B (1/256), Anti-D (1/64)) | 3x10ml |
| 8.02.05.7.0020 | ABO Set: Anti-A (1/256), Anti-B (1/256) | 2x10ml |
| 8.02.52.0.0010 8.02.52.0.0100 | Rh-D Negative Control | 10 ml/vial 10x10 ml |
| 8.02.63.1.0010 8.02.63.0.0100 | Antibody Enhancement Solution (LISS) | 10 ml/vial 10x10 ml |
| 8.02.23.0.0002 | Anti-M, Human, 2ml/vial | 2 ml/vial |
| 8.02.24.0.0002 | Anti-N, Lectin (Viciagraminea), 2ml/vial | 2 ml/vial |
| 8.02.25.0.0002 | Anti-S, Human, 2ml/vial | 2 ml/vial |
| 8.02.26.0.0002 | Anti-s, Human, 2ml/vial | 2 ml/vial |
| 8.02.34.0.0002 | Anti-Lua, Human, 2ml/vial | 2 ml/vial |
| 8.02.37.0.0002 | Anti-Lea, Monoclonal, 2ml/vial | 2 ml/vial |
| 8.02.38.0.0002 | Anti-Leb, Monoclonal, 2ml/vial | 2 ml/vial |
| 8.02.39.0.0002 | Anti-P1, Monoclonal, 2ml/vial | 2 ml/vial |

HEMATOLOGY TESTS

Overview

Atlas Medical supplies coagulation reagents.

The coagulation regents include PT,PTT and fibronogen in liquid formats and in various sizes to suit most lab applications.

The range also includes normal and abnormal coagulation controls.

Features

- $\subleft ?$ Some kits includes normal and abnormal controls .
- ∂ The Kit comes in sizes of 50 and 100 tests.
- Atlas Medical provides high quality coagulation reagents that are accurate, easy to use, competitively priced, and conveniently packed in different sizes and options.

| Co-agguiation Reagents | | |
|--|---|---|
| Item Code | Item Description | Sizes |
| 8.02.40.1.0010 8.02.40.1.0050 8.02.40.1.0100 | PT Calcium Rabbit Brain Thromboplastin, Liquid | 2ml (20 Tests) 5ml (50 Tests) 10ml (100 Tests) |
| 8.02.41.1.0040 8.02.41.1.0050 8.02.41.1.0100 | APTT (PTT) Micronised Silica Platelet Substitute | 2ml (40 Tests) 2.5ml (50 Tests) 5ml (100 Tests) |
| 8.02.44.0.0040 8.02.44.0.0100 | PT Kit with Normal Control | 2x2ml + 1ml 2x10ml + 1ml |
| 8.02.45.0.0080 8.02.45.0.0200 | APTT (PTT) Kit with Normal Control | 2x2ml + 1ml 2x10ml + 1ml |
| 8.02.48.0.0010 8.02.48.0.0100 | Calcium Chloride, 25 mM | 10ml/vial 10ml/vial / 10 Vials / Box |
| 8.02.60.0.0006 | Normal Coagulation Control | 6x1ml |
| 8.02.61.0.0006 | Abnormal Coagulation Control | 6x1ml |
| 8.02.69.0.0005 | Fibrinogen Reagent | 5ml/vial |
| 8.02.45.1.0080 | APTT (PTT) Kit (Calcium Chloride reagent + Normal Control) | 80 Tests |
| 8.02.64.0.0006 | Normal & Abnormal Coaggulation Control | 3x1ml |
| 8.02.69.0.0100 | Fibronogen Test kit KIT | 100 Tests |

Hemoglobin Reagents Item Code Item Description Sizes 8.02.46.1.0500 Drabkins Reagent, 40x (White Plastic Bottle) 50ml/Bottle 8.02.46.1.3000 Drabkins Reagent, 40x (White Plastic Bottle) 50ml/Bottle 8.02.50.0.0010 Haemoglobin Standard 10ml/vial



SICKLE CELL KITS

Overview

Sickle cell disease (also called sickle cell anemia) is an inherited blood disorder that affects red blood cells. The sickle cell gene causes the body to produce abnormal hemoglobin.

Features

- Atlas Sickle Cell Kits is a qualitative solubility test for Sickle Haemoglobin.
- The test can be performed in two ways:
 1. A screening test to detect sickle haemoglobin (HbS)
 2. A centrifugation test to differentiate the sickle cell trait (AS) from sickle cell anaemia (SS).



11

INFECTIOUS DISEASE RAPID TESTS

ANTIBODY TESTING

Overview

Atlas Medical offers an extensive range of lateral flow immunoassay tests for the rapid detection of antibodies and antigens in human samples (blood, serum, plasma, urine, oral swabs, nasal swabs, and feces). This range includes tests to detect a wide variety of viruses, microorganisms and parasites.

Features

- Atlas Medical infectious disease rapid tests are reliable, accurate and supplied in both cassette and strip formats.
- The kits are conveniently packed in different sizes of 20, 25, 30, 40, 50 and 100 tests per kit and include the necessary test accessories to perform the assay.

| Item Code | Item Description | Sizes |
|-----------------|--|---------------|
| 8.04.27.0.0001 | HIV 1/2 Antibody Test Cassette, | Bulk |
| 8.04.27.0.0020 | Whole Blood/Serum/Plasma, Individually Pouched | 20 Tests/Box |
| 8.04.28.0.0001 | HIV 1/2 Antibody Test Cassette, | Bulk |
| 8.04.28.0.0020 | Serum/Plasma, Individually Pouched | 20 Tests/Box |
| 8.04.29.0.0001 | HIV 1/2 Antibody Test Strip, | Bulk |
| 8.04.29.0.0100 | Serum/Plasma, Individually Pouched | 100 Tests/Bo> |
| 8.04.30.0.0001 | HCV Antibody Test Cassette, | Bulk |
| 8.04.30.0.0020 | Whole Blood/Serum/Plasma, Individually Pouched | 20 Tests/Box |
| 8.04.31.0.0001 | HCV Antibody Test Cassette, | Bulk |
| 8.04.31.0.0020 | Serum/Plasma, Individually Pouched | 20 Tests/Box |
| 8.04.32.0.0001 | HCV Antibody Test Strip, | Bulk |
| 8.04.32.0.0100 | Serum/Plasma, Individually Pouched | 100 Tests/Bo> |
| 8.04.35.0.0001 | HBs Antibody Test Cassette, | Bulk |
| 8.04.35.0.0020 | Serum/Plasma, Individually Pouched | 20 Tests/Box |
| 8.04.36.0.0001 | HBs Antibody Test Strip, | Bulk |
| 8.04.36.0.0100 | Serum/Plasma, Individually Pouched | 100 Tests/Bo> |
| 8.04.107.0.0040 | Toxo lgG test cassette indiviually packed | 40 Tests/Box |



| | | \sim |
|----------------|--|---------------|
| Item Code | Item Description | Sizes |
| 8.04.20.0.0001 | H.pylori Antibody Test Cassette, | Bulk |
| 8.04.20.0.0020 | Whole Blood/Serum/Plasma, Individually Pouched | 20 Tests/Box |
| 8.04.21.0.0001 | H.pylori Antibody Test Cassette, | Bulk |
| 8.04.21.0.0020 | Serum/Plasma, Individually Pouched | 20 Tests/Box |
| 8.04.22.0.0001 | H.pylori Antibody Test Strip, | Bulk |
| 8.04.22.0.0100 | Serum/Plasma, Individually Pouched | 100 Tests/Bo> |
| 8.04.41.0.0001 | Syphilis Antibody Test Cassette, | Bulk |
| 8.04.41.0.0020 | Whole Blood/Serum/Plasma, Individually Pouched | 20 Tests/Box |
| 8.04.42.0.0001 | Syphilis Antibody Test Cassette, | Bulk |
| 8.04.42.0.0020 | Serum/Plasma, Individually Pouched | 20 Tests/Box |
| 8.04.43.0.0001 | Syphilis Antibody Test Strip, | Bulk |
| 8.04.43.0.0100 | Serum/Plasma, Individually Pouched | 100 Tests/Bo> |
| 8.16.16.0.0001 | TB Test Cassette, | Bulk |
| 8.16.16.0.0020 | Serum/Plasma, Individually Pouched | 20 Tests/Box |
| 8.16.16.1.0020 | TB Test Cassette, | 20 Tests/Box |
| 8.16.16.1.0025 | Whole Blood, Individually Pouched | 25 Tests/Box |
| | | |



Œ

| Item Code | Item Description | Sizes |
|----------------------------------|---|----------------------|
| 8.42.47.0.0020 | HAV IgM Test Cassette, Individually Pouched | 20 Tests/Box |
| 8.04.44.0.0001 8.04.44.0.0020 | Dengue IgG/IgM Test Cassette, Whole Blood/Serum/Plasma, Individually Pouched | Bulk 20 Tests/Box |

RAPID TESTS

INFECTIOUS DISEASE RAPID TESTS

ANTIGEN TESTING

|--|

| STOOL SAM | PLE | |
|----------------|--|--------------|
| ltem Code | Item Description | Sizes |
| 8.04.23.1.0020 | H.pylori Antigen Test Cassette, Stool Sample, Individually Pouched | 20 Tests/Box |
| 8.04.24.1.0025 | H.pylori Antigen Test Strip, Stool Sample, Individually Pouched | 25 Tests/Box |
| 8.04.69.0.0020 | Rotavirus Antigen Test Cassette, Stool Sample, Individually Pouched | 20 Tests/Box |
| 8.04.70.0.0025 | Rotavirus Antigen Test Strip, Stool Sample, Individually Pouched | 25 Tests/Box |
| 8.04.71.0.0020 | Adenovirus Antigen Test Cassette, Stool Sample, Individually Pouched | 20 Tests/Box |
| 8.04.72.0.0025 | Adenovirus Antigen Test Strip, Stool Sample, Individually Pouched | 25 Tests/Box |
| 8.04.73.0.0020 | Rota-Adeno Antigens Combo test Cassette, Stool Sample, Individually Pouched | 20 Tests/Box |
| 8.04.74.0.0025 | Rota-Adeno Antigens Combo test Strip, Stool Sample, Individually Pouched | 25 Tests/Box |
| 8.16.01.0.0020 | Crypto Virus Test Cassette, Stool Sample, Individually Pouched | 20 Tests/Box |
| 8.16.02.0.0025 | Crypto Virus Test Strip, Stool Sample, Individually Pouched | 25 Tests/Box |
| 8.16.31.0.0020 | Giardia Test Cassette, Stool Sample, Individually Pouched | 20 Tests/Box |
| 8.16.30.0.0025 | Giardia Test Strip, Stool Sample, Individually Pouched | 25 Tests/Box |
| 8.16.33.0.0020 | Crypto-Giardia Test Cassette, Stool Sample, Individually Pouched | 20 Tests/Box |
| 8.16.32.0.0025 | Crypto-Giardia Test Strip, Stool Sample, Individually Pouched | 25 Tests/Box |
| 8.16.41.0.0020 | E.coli Test Cassette, Stool Sample, Individually Pouched | 20 Tests/Box |
| 8.16.40.0.0025 | E.coli Test Strip, Stool Sample, Individually Pouched | 25 Tests/Box |
| 8.16.82.0.0025 | Salmonella typhi Antigen Test Cassette, Stool Sample, Individually Pouched | 25 Tests/Box |
| 8.16.85.0.0025 | Salmonella paratyphi Antigen Test Cassette, Stool Sample, Individually Pouched | 25 Tests/Box |
| 8.16.80.0.0020 | Clostridium difficile Antigen Test Cassette, Stool Sample, Individually Pouched | 20 Tests/Box |
| 8.16.81.0.0025 | Clostridium difficile Antigen Test Strip, Stool Sample, Individually Pouched | 25 Tests/Box |
| 8.16.86.0.0025 | C.Difficle Toxin A+B, Test Cassette, Stool Sample, Individually Pouched | 25 Tests/Box |
| 8.16.91.0.0025 | Norovirus Genogroups I & II Ag, Test Cassette, Stool Sample, Individually Pouched | 25 Tests/Box |

| | B | 500 |
|------------------------|----------|----------------------|
| SE | 6 | E |
| Stool collecto tube | (1) | |
| P' | | Add 4 dr window (|
| | | H. pyton |





84

.

b.

| SWAD & NAS | | |
|----------------|---|--------------|
| ltem Code | Item Description | Sizes |
| 8.04.25.0.0020 | Strep A Test Cassette, Swab Sample | 20 Tests/Box |
| 8.45.00.0.0020 | Strep B Test Cassette, Swab Sample | 20 Tests/Box |
| 8.45.01.0.0020 | Strep A+B Test Cassette, Swab Sample | 20 Tests/Box |
| 8.04.86.0.0020 | Influenza A+B Test Cassette, Nasal Sample | 20 Tests/Box |
| 8.04.96.0.0025 | Influenza A+B Test Strip, Nasal Sample | 25 Tests/Box |
| 8.16.20.0.0020 | RSV Test Cassette, Swab Sample | 20 Tests/Box |
| 8.16.22.0.0025 | RSV Test Strip, Swab Sample | 25 Tests/Box |
| 8.16.37.0.0020 | Adeno Respiratory Antigen Test Cassette, Swab Sample | 20 Tests/Box |
| 8.16.36.0.0025 | Adeno Respiratory Antigen Test Strip, Swab Sample | 25 Tests/Box |
| 8.16.39.0.0020 | Adeno - RSV Respiratory Test Cassette, Swab Sample | 20 Tests/Box |
| 8.16.38.0.0025 | Adeno - RSV Respiratory Test Strip, Swab Sample | 25 Tests/Box |

Œ

Œ

Œ

| BLOOD SAMPLE | | | |
|----------------|--|--------------|--|
| Item Code | Item Description | Sizes | |
| 8.04.37.0.0020 | Malaria Pf. Test Cassette, Whole Blood, Individually Pouched | 20 Tests/Box | |
| 8.16.14.0.0020 | Malaria Pf/Pv. Test Cassette, Whole Blood, Individually Pouched | 20 Tests/Box | |

| ltem Code | Item Description | Sizes |
|----------------|---|---------------|
| 8.16.24.0.0001 | HBsAg Test Cassette (Whole Blood/Serum/ | Bulk |
| 8.16.24.0.0020 | Plasma), Individually Pouched | 20 Tests/Box |
| 8.04.33.0.0001 | HBsAg Test Cassette, Serum/Plasma, | Bulk |
| 8.04.33.0.0020 | Individually Pouched | 20 Tests/Box |
| 8.04.34.0.0001 | HBsAg Test Strip, Serum/Plasma, | Bulk |
| 8.04.34.0.0100 | Individually Pouched | 100 Tests/Box |
| 8.04.26.0.0020 | Chlamydia Test Cassette, Urine or Swab | 20 Tests/Box |
| 8.63.00.0.0025 | Chlamydia + Gonorrhea Rapid Test Cassette (Cervical/Urethral swab) | 25 Tests/Box |

13

URINE REAGENT STRIPS

Overview

Urine Reagent Strips (URS) are widely used in Urinalysis to determine pathological changes in urine. The strips contain dry-chemistry pads that, when dipped in urine, change their colors. The color change allows for the semi-quantitative measurement of various urine parameters. The strips are suitable for lab, point-of-care and even home use.

Features

- Atlas Medical Urine Reagent Strips can be used to detect up to 14 urine parameters.
- They are simple to use and the results are visually read within a minute.
- The strips are packed in desiccated bottles of 50 or 100 strips.
- Atlas Medical can also provide suitable readers to read the strip colors and document the results.

URINE ANALYZER

Atlas Urine Analyzer is a manual analyzer that detects Photosensitive Diode using the method of Reflectance Photometry. Test Categories include Routine, STAT and QC. Atlas Urine Analyzer has an Automatic calibration for accurate results and easy operation. It can read strips with up to 14 parameters, including Microalbumin/Creatininenine/Calcium. It has an option to print results for quick and simple record management.

- Accurate.
- 🦪 Reliable.
- Convenient.
- Easy Data Management.

| 6 | \geq |
|---|--------|
| | ים |
| | ς, |

| | | <u> </u> |
|------------|------------------------------------|----------|
| Item Code | Item Description | Sizes |
| 8.001.U120 | Urine Analyzer For Clinics, U120 | 1 Unit |
| 8.002.U500 | Urine Analyzer For Hospitals, U500 | 1 Unit |





| ltem Code | Item Description | Sizes |
|----------------------------------|--|-------------------------|
| 8.03.00.0.0050 8.03.00.0.0100 | URS 1 Parameter: Glucose | 50 Strips 100 Strips |
| 8.03.01.0.0050 8.03.01.0.0100 | URS 1 Parameter: Protein | 50 Strips 100 Strips |
| 8.03.02.0.0050 8.03.02.0.0100 | URS 1 Parameter: Ketone | 50 Strips 100 Strips |
| 8.03.45.0.0050 | URS 1 Parameter Blood, (5mm) | 50 Strips |
| 8.03.03.0.0050 8.03.03.0.0100 | URS 2 Parameters: Glucose, Ketone | 50 Strips 100 Strips |
| 8.03.04.0.0050 8.03.04.0.0100 | URS 2 Parameters: Glucose, Protein | 50 Strips 100 Strips |
| 8.03.05.0.0100 | URS 2 Parameters: Sample end: Urobilinogen, Bilirubin | 100 Strips |
| 8.03.19.0.0050 8.03.19.0.0100 | URS 2 Parameters(5mm): Sample End: Creatinine, pH | 50 Strips 100 Strips |
| 8.03.06.0.0050 8.03.06.0.0100 | URS 3 Parameters: Protein, pH, Glucose | 50 Strips 100 Strips |
| 8.03.07.0.0100 | URS 3 Parameters: Glucose, Protein, Ketone | 100 Strips |
| 8.03.08.0.0100 | URS 3 Parameters: Sample end:pH, Ketone, Glucose | 100 Strips |
| 8.03.09.0.0100 | URS 3 Parameters: Sample end:Leukocytes, Nitrite, Blood | 1 100 Strips |
| 8.03.10.0.0050 8.03.10.0.0100 | URS 3 Parameters: Sample end:Protein, Specific Gravity, Creatinine | 50 Strips 100 Strips |
| 8.03.11.0.0100 | URS 4 Parameters: Protein, pH, Specific Gravity, Glucose | 100 Strips |
| 8.03.12.0.0100 | URS 4 Parameters: Protein, pH, Blood, Glucose | 100 Strips |
| 8.03.13.0.0050 8.03.13.0.0100 | URS 5 Parameters: Glucose, Protein, Ketone, pH, Blood | 50 Strips 100 Strips |
| 8.03.25.0.0100 | URS 5 Parameters(5mm): Blood, Glucose, Protein, Nitrite Leucocytes | , 100 Strips |
| 8.03.14.0.0100 | URS 6 Parameters: Leukocytes, Nitrite, Protein, pH, Blood, Glucose | 100 Strips |
| 8.03.44.0.0100 | URS 7 Parameter: Glucose, Ketone, Protien, PH, Blood, Bilirubin, Urobilinogen | 100 Strips |
| 8.03.23.0.0100 | URS 8 Parameters: Glucose, Protein, pH, Ketone, Urobilinogen, Bilirubin, Blood, Nitrite | 100 Strips |
| 8.03.15.0.0100 | URS 9 Parameters: Nitrite, Urobilinogen, Protein, pH, Blood, Specific Gravity, Ketone, Bilirubin, Glucose | 100 Strips |
| 8.03.16.0.0100 | URS 10 Parameters: Leukocytes, Nitrite, Urobilinogen Protein, pH, Blood, Specific Gravity, Ketone, Bilirubin, Glucose | , 100 Strips |
| 8.03.17.0.0050 8.03.17.0.0100 | URS 10 Parameters: Sample end: Nitrite, Urobilinogen, Protein, pH, Blood, Specific Gravity, Ketone,Bilirubin, Glucose, Ascorbic Acid | 50 Strips 100 Strips |
| 8.03.18.0.0100 | URS 11 Parameters: Leukocytes, Nitrite, Urobilinogen Protein, pH, Blood, Specific Gravity, Ketone, Bilirubin, Glucose, Ascorbic Acid | 100 Strips |
| 8.03.47.0.0100 | URS 14 Parameters (ASC, GLU, BIL, KET, SG, BLO, PH, PRO URO, NIT, LEU, ALB, CRE, CA) | , 100 Strips |

Œ

FERTILITY RAPID TESTS

Overview

Atlas Medical Fertility Rapid Tests are based on lateral flow immunoassay for the detection of human chorionic gonadotropin (hCG), Ovulation (LH), and Human Follicular Stimulating Hormone (FSH) in urine. Each of the three tests comes in strip, cassette, or midstream formats and are conveniently packed in sizes to suit lab, point-of-care and home uses.

Features

- Accurate.
- ∂ Convenient.
- ∂ Easy to use (add or dip in urine).
- ∂ Competitively priced.
- ∂ Results are obtained in 1 to 5 minutes.
- Ø Different strip sizes are available.





| | | \smile |
|----------------------------------|---|-----------------------|
| ltem Code | Item Description | Sizes |
| 8.04.00.0.0001 8.04.00.0.0020 | hCG Test Cassette, Urine, Individually Pouched | Bulk 20 Tests/Box |
| 8.04.01.0.0001 8.04.01.0.0020 | hCG Test Cassette, Urine/Serum, Individually Pouched | Bulk 20 Tests/Box |
| 8.04.04.0.0001 8.04.04.0.0100 | hCG Test Strip, Urine, Individually Pouched | Bulk 100 Tests/Box |
| 8.04.10.0.0001 8.04.10.0.0100 | hCG Test Strip, Urine/Serum, Individually Pouched | Bulk 100 Tests/Box |
| 8.04.13.0.0001 8.04.13.0.0015 | hCG Midstream Test, Individually Pouched | Bulk 15 Tests/Box |
| 8.04.14.0.0001 8.04.14.0.0020 | LH Test Cassette, Urine, Individually pouched | Bulk 20 Tests/Box |
| 8.04.15.0.0001 8.04.15.0.0100 | LH Test Strip, Urine, Individually pouched | Bulk 100 Tests/Box |
| 8.04.16.0.0001 8.04.16.0.0015 | LH Midstream Test, Individually Pouched | Bulk 15 Tests/Box |
| 8.04.17.0.0001 8.04.17.0.0020 | FSH Test Cassette, Urine, Individually pouched | Bulk 20 Tests/Box |
| 8.04.18.0.0001 8.04.18.0.0100 | FSH Test Strip, Urine, Individually pouched | Bulk 100 Tests/Box |
| 8.04.19.0.0001 8.04.19.0.0015 | FSH Midstream Test, Individually Pouched | Bulk 15 Tests/Box |

KIDNEY FUNCTION RAPID TESTS

Overview

Atlas Medical Microalbumin Rapid Test is a rapid visual immunoassay used for the qualitative detection of microalbumin in human urine samples. This kit is intended for use as an aid in the diagnosis of renal dysfunction.

Features

∂ The test comes in cassette and Strip formats.

| | | | Œ |
|---|----------------------------------|--|-------------------------------|
| ĺ | Item Code | Item Description | Sizes |
| | 8.16.52.0.0001 8.16.52.0.0020 | Microalbumin Test Cassette ,Individually Pouched | Bulk 20 Tests/Box |
| | 8.16.53.0.0030 8.16.53.0.0100 | Microalbumin Test Strip ,Individually Pouched | 30 Tests/Box 100 Tests/Box |

INFLAMMATION AND CANCER MARKERS

Overview

All the tests in this group are qualitative and based on lateral flow immunoassay for the detection of various inflammation and cancer markers.

| | | \sim |
|-----------------|--|---------------|
| ltem Code | Item Description | Sizes |
| 8.04.38.0.0020 | Fecal Occult Blood Test (FOB) Test Cassette, Stool Sample, Individually Pouched | 20 Tests/Box |
| 8.04.85.0.0050 | Fecal Occult Blood Test (FOB) Test Strip Stool Sample, Individually Pouched | 50 Tests/Box |
| 8.04.109.0.0020 | Procalcitonin Test Cassette (PCT), (Serum/Plasma |)20 Tests/Box |
| 8.48.00.0.0020 | Procalcitonin Test Cassette (PCT), (Whole Blood / Serum/ Plasma) | 20 Tests/Box |
| 8.16.78.0.0025 | Calprotectin Test Cassette | 25 Tests/Box |

 (\mathbf{C})

Features

- Atlas Medical inflammation and cancer markers rapid tests are supplied in both cassette and strip formats.
- The kits are conveniently packed in different kit sizes of 20, 25, 30 and 100 tests per kit.

| ltem Code | Item Description | Sizes |
|----------------|---|---------------|
| 8.16.28.0.0001 | PSA Test Cassette, Whole Blood/Serum/Plasma), | Bulk |
| 8.16.28.0.0020 | Individually Pouched | 20 Tests/Box |
| 8.04.39.0.0001 | PSA Test Cassette, Serum/Plasma, Individually | Bulk |
| 8.04.39.0.0020 | Pouched | 20 Tests/Box |
| 8.04.40.0.0001 | PSA Test Strip, Serum/Plasma, Individually | Bulk |
| 8.04.40.0.0100 | Pouched | 100 Tests/Box |

CARDIAC MARKERS RAPID TESTS

Overview

Atlas Medical offers lateral flow immunoassay rapid tests to detect the three major cardiac markers namely: Troponin I, Myoglobin and CK-MB, as an aid in the diagnosis of myocardial infarction (MI).

Features

- They can be used on whole blood (in addition to serum/plasma) making them ideal for emergency rooms.
- They come in single test or triple combo test cassette formats.
- The kits are conveniently packed in different kit sizes of 20, 25, 30 tests per kit.

| | | \sim |
|----------------------------------|---|----------------------|
| ltem Code | Item Description | Sizes |
| 8.04.45.0.0001 | Troponin I Test Cassette, Whole Blood/Serum/ | Bulk |
| 8.04.45.0.0020 | Plasma, Individually Pouched | 20 Tests/Box |
| 8.04.46.0.0001 | Myoglobin Test Cassette, Whole Blood/Serum/ | Bulk |
| 8.04.46.0.0020 | Plasma, Individually Pouched | 20 Tests/Box |
| 8.04.47.0.0001 | CK-MB Test Cassette, Whole Blood/Serum/ | Bulk |
| 8.04.47.0.0020 | Plasma, Individually Pouched | 20 Tests/Box |
| 8.04.48.0.0001 8.04.48.0.0020 | Cardiac Triple Test Cassette (Troponin I, CK-MB, Myoglobin), Whole Blood/Serum/Plasma, Individually Pouched | Bulk 20 Tests/Box |

Œ

DOA RAPID TESTS

Œ

Overview

All the tests in this group are qualitative and based on lateral flow immunoassay for the detection of various Drug of Abuse .



Features

- Atlas Medical DOA rapid tests are supplied in cassette, strip, panel and cup formats.
- The kits are conveniently packed in different kit sizes of 20, 25, 30, 50 and 100 tests per kit.

| | | <u> </u> | |
|---------------------------|---------------------------------------|---------------|--|
| STRIP AND CASSETTE FORMAT | | | |
| Item Code | Item Description | Sizes | |
| 8.04.49.0.0001 | Morphine Test Cassette, Urine, | Bulk | |
| 8.04.49.0.0020 | Individually Pouched | 20 Tests/Box | |
| 8.04.50.0.0001 | Morphine Test Strip, Urine, | Bulk | |
| 8.04.50.0.0100 | Individually Pouched | 100 Tests/Box | |
| 8.04.51.0.0001 | Marijuana (THC) Test Cassette, Urine, | Bulk | |
| 8.04.51.0.0020 | Individually Pouched | 20 Tests/Box | |
| 8.04.52.0.0001 | Marijuana (THC) Test Strip, Urine, | Bulk | |
| 8.04.52.0.0100 | Individually Pouched | 100 Tests/Box | |
| 8.04.53.0.0001 | Amphetamine Test Cassette, Urine, | Bulk | |
| 8.04.53.0.0020 | Individually Pouched | 20 Tests/Box | |
| 8.04.54.0.0001 | Amphetamine Test Strip, Urine, | Bulk | |
| 8.04.54.0.0100 | Individually Pouched | 100 Tests/Box | |
| 8.04.55.0.0001 | Barbiturates Test Cassette, Urine, | Bulk | |
| 8.04.55.0.0020 | Individually Pouched | 20 Tests/Box | |
| 8.04.56.0.0001 | Barbiturates Test Strip, Urine, | Bulk | |
| 8.04.56.0.0100 | Individually Pouched | 100 Tests/Box | |
| 8.04.57.0.0001 | Benzodiazepines Test Cassette, Urine, | Bulk | |
| 8.04.57.0.0020 | Individually Pouched | 20 Tests/Box | |
| 8.04.58.0.0001 | Benzodiazepines Test Strip, Urine, | Bulk | |
| 8.04.58.0.0100 | Individually Pouched | 100 Tests/Box | |
| 8.04.59.0.0001 | Cocaine Test Cassette, Urine, | Bulk | |
| 8.04.59.0.0020 | Individually Pouched | 20 Tests/Box | |
| 8.04.60.0.0001 | Cocaine Test Strip, Urine, | Bulk | |
| 8.04.60.0.0100 | Individually Pouched | 100 Tests/Box | |
| 8.04.61.0.0001 | Methamphetamine Test Cassette, Urine, | Bulk | |
| 8.04.61.0.0020 | Individually Pouched | 20 Tests/Box | |
| 8.04.62.0.0001 | Methamphetamine Test Strip, Urine, | Bulk | |
| 8.04.62.0.0100 | Individually Pouched | 100 Tests/Box | |
| 8.04.63.0.0001 | Methadone Test Cassette, Urine, | Bulk | |
| 8.04.63.0.0020 | Individually Pouched | 20 Tests/Box | |
| 8.04.64.0.0001 | Methadone Test Strip, Urine, | Bulk | |
| 8.04.64.0.0100 | Individually Pouched | 100 Tests/Box | |
| 8.04.65.0.0001 | Phencyclidine Test Cassette, Urine, | Bulk | |
| 8.04.65.0.0020 | Individually Pouched | 20 Tests/Box | |
| 8.04.66.0.0001 | Phencyclidine Test Strip, Urine, | Bulk | |
| 8.04.66.0.0100 | Individually Pouched | 100 Tests/Box | |

| ltem Code | Item Description | Sizes |
|----------------|--|---------------|
| 8.04.67.0.0001 | Tricyclic Anti-Depressants Test Cassette, Urine, | Bulk |
| 8.04.67.0.0020 | Individually Pouched | 20 Tests/Box |
| 8.04.68.0.0001 | Tricyclic Anti-Depressants Test Strip, Urine, | Bulk |
| 8.04.68.0.0100 | Individually Pouched | 100 Tests/Box |
| 8.04.99.0.0001 | Buprenorphine Test Cassette, Urine, | Bulk |
| 8.04.99.0.0020 | Individually Pouched | 20 Tests/Box |
| 8.16.23.0.0001 | Buprenorphine Test Strip, Urine, | Bulk |
| 8.16.23.0.0100 | Individually Pouched | 100 Tests/Box |
| 8.16.68.0.0001 | Tramadol Test Cassette, Urine, | Bulk |
| 8.16.68.0.0020 | Individually Pouched | 20 Tests/Box |
| 8.16.44.0.0001 | Tramadol Test Strip, Urine, | Bulk |
| 8.16.44.0.0100 | Individually Pouched | 100 Tests/Box |
| 8.16.15.0.0001 | Methylenedioxymethamphetamine (MDMA) | Bulk |
| 8.16.15.0.0020 | Ecstasy Test Cassette,Urine,Individually Pouched | 20 Tests/Box |
| 8.16.05.0.0001 | Methylenedioxymethamphetamine (MDMA) | Bulk |
| 8.16.05.0.0100 | Ecstasy Test Strip, Urine,Individually Pouched | 100 Tests/Box |
| 8.16.06.0.0001 | Opiates Test Cassette, Urine, | Bulk |
| 8.16.06.0.0020 | Individually Pouched | 20 Tests/Box |
| 8.16.07.0.0001 | Opiates Test Strip, Urine, | Bulk |
| 8.16.07.0.0100 | Individually Pouched | 100 Tests/Box |
| 8.16.58.0.0001 | Cotinine Test Cassette, Urine, | Bulk |
| 8.16.58.0.0020 | Individually Pouched | 20 Tests/Box |
| 8.16.59.0.0001 | Cotinine Test Strip, Urine, | Bulk |
| 8.16.59.0.0100 | Individually Pouched | 100 Tests/Box |
| 8.16.62.0.0001 | Oxycodone Test Cassette, Urine, | Bulk |
| 8.16.62.0.0020 | Individually Pouched | 20 Tests/Box |
| 8.16.63.0.0001 | Oxycodone Test Strip, Urine, | Bulk |
| 8.16.63.0.0100 | Individually Pouched | 100 Tests/Box |
| 8.16.64.0.0001 | Ketamine Test Cassette, Urine, | Bulk |
| 8.16.64.0.0020 | Individually Pouched | 20 Tests/Box |
| 8.16.65.0.0001 | Ketamine Test Strip, Urine, | Bulk |
| 8.16.65.0.0100 | Individually Pouched | 100 Tests/Box |
| 8.16.66.0.0001 | Proxyphene Test Cassette, Urine, | Bulk |
| 8.16.66.0.0020 | Individually Pouched | 20 Tests/Box |
| 8.16.67.0.0001 | Proxyphene Test Strip, Urine, | Bulk |
| 8.16.67.0.0100 | Individually Pouched | 100 Tests/Box |
| 8.16.68.0.0001 | Tramadol Test Cassette, Urine, | Bulk |
| 8.16.68.0.0020 | Individually Pouched | 20 Tests/Box |
| 8.16.69.0.0001 | EDDP Test Cassette, Urine, | Bulk |
| 8.16.69.0.0020 | Individually Pouched | 20 Tests/Box |
| 8.16.70.0.0001 | EDDP Test Strip, Urine, | Bulk |
| 8.16.70.0.0100 | Individually Pouched | 100 Tests/Box |
| 8.16.60.0.0001 | Dolantin Test Cassette, Urine, | Bulk |
| 8.16.60.0.0020 | Individually Pouched | 20 Tests/Box |
| 8.16.61.0.0001 | Dolantin Test Strip, Urine, | Bulk |
| 8.16.61.0.0100 | Individually Pouched | 100 Tests/Box |

DOA RAPID TESTS

| | | Œ |
|----------------------------------|--|----------------------|
| ltem Code | Item Description | Sizes |
| 8.04.93.0.0001 | DOA Panel: 2 Drugs (Combination of any 2 drugs) | Bulk |
| 8.04.93.0.0025 | Urine, Individually Pouched | 25 Tests/Box |
| 8.04.94.0.0001 | DOA Panel: 3 Drugs (Combination of any 3 drugs) | Bulk |
| 8.04.94.0.0025 | Urine, Individually Pouched | 25 Tests/Box |
| 8.04.95.0.0001 | DOA Panel: 4 Drugs (Combination of any 4 drugs) | Bulk |
| 8.04.95.0.0025 | Urine, Individually Pouched | 25 Tests/Box |
| 8.04.79.0.0001 | DOA Panel: 5 Drugs (Combination of any 5 drugs) | Bulk |
| 8.04.79.0.0025 | Urine, Individually Pouched | 25 Tests/Box |
| 8.04.80.0.0001 | DOA Panel: 6 Drugs (Combination of any 6 drugs) | Bulk |
| 8.04.80.0.0025 | Urine, Individually Pouched | 25 Tests/Box |
| 8.04.81.0.0001 | DOA Panel: 7 Drugs (Combination of any 7 drugs) | Bulk |
| 8.04.81.0.0025 | Urine, Individually Pouched | 25 Tests/Box |
| 8.04.82.0.0001 | DOA Panel: 8 Drugs (Combination of any 8 drugs) | Bulk |
| 8.04.82.0.0025 | Urine, Individually Pouched | 25 Tests/Box |
| 8.04.83.0.0001 | DOA Panel: 9 Drugs (Combination of any 9 drugs) Urine, Individually Pouched | Bulk |
| 8.04.84.0.0001 | DOA Panel: 10 Drugs (Combination of any 10 | Bulk |
| 8.04.84.0.0025 | drugs), Urine, Individually Pouched | 25 Tests/Box |
| 8.16.03.0.0001 | DOA Panel: 11 Drugs (Combination of any 11 | Bulk |
| 8.16.03.0.0025 | drugs),Urine, Individually Pouched | 25 Tests/Box |
| 8.16.04.0.0001 8.16.04.0.0025 | DOA Panel: 12 Drugs (Combination of any 12 drugs), Urine, Individually Pouched | Bulk 25 Tests/Box |





| | | Œ |
|----------------|---|-------|
| Item Code | Item Description | Sizes |
| 8.42.43.0.0001 | Drug Of Abuse Cup, 7 parameters (Combination of any 7 drugs), Urine, Individually Pouched | Bulk |
| 8.42.50.0.0001 | Drug Of Abuse Cup, 8 parameters (Combination of any 8 drugs), Urine, Individually Pouched | Bulk |
| 8.16.71.0.0001 | Drug Of Abuse Cup, 10 parameters (Combination of any 10 drugs), Urine, Individually Pouched | Bulk |

| | | Œ |
|----------------|---|-------|
| ltem Code | Item Description | Sizes |
| 8.16.87.0.0001 | DOA Multiscreen Cassette: 3 Drugs (Combination of any 3 drugs), Urine, Individually Pouched | Bulk |
| 8.16.88.0.0001 | DOA Multiscreen Panel 7 Drugs (Combination of any 7 drugs), Urine, Individually Pouched | Bulk |
| 8.16.89.0.0001 | DOA Multiscreen Cassette: 8 Drugs (Combination of any 8 drugs), Urine, Individually Pouched | Bulk |



CLINICAL CHEMISTRY KITS

Œ

Overview

Kits in this group measure concentrations of electrolytes , hormones, proteins, and other metabolic products in human blood , serum,plasma, CSF and urine .

Clincal Chemistry tests are indicated to assess systemic functions such liver function , kidney function , and endocrine and metabolic function .

Methods commonly used are colormetric and kinetic.

Features

The kits are conveniently packed in different kit sizes of 20, 30, 60, 75, 100, 150, 200, 250, 500, and 1000 tests per kit.

| Item Code | Item Description | Sizes |
|----------------------------------|--|----------------------|
| 8.05.00.0.0250 8.05.00.0.0500 | Albumin Bromocresol Green | 2x125ml 4x125ml |
| 8.05.01.0.0030 8.05.01.0.0060 | Amylase | 3x10ml 6x10ml |
| 8.05.02.0.0020 8.05.02.0.0040 | Acid Phosphatase Kinetic, Hilmann Method (Tablets) | 10x2ml 20x2ml |
| 8.05.03.0.0030 8.05.03.0.0060 | Alkaline Phosphatase Kinetic, DGKC Method (Tablets) | 10x3ml 20x3ml |
| 8.05.04.0.0250 8.05.04.0.0500 | Alkaline Phosphatase Kinetic, DGKC Method (Liquid) | 5x50ml 5x100ml |
| 8.05.05.0.0250 8.05.05.0.0500 | Bilirubin Total (DMSO Method) | 2x125ml 4x125ml |
| 8.05.06.0.0250 8.05.06.0.0500 | Bilirubin Direct (DMSO Method) | 2x125ml 4x125ml |
| 8.05.07.0.0250 8.05.07.0.0500 | Bilirubin Total & Direct (DMSO Method) | 2x125ml 4x125ml |
| 8.05.08.0.0250 8.05.08.0.0500 | Calcium Arsenazo III | 2x125ml 4x125ml |
| 8.05.09.0.0250 8.05.09.0.0500 | Calcium O-Cresolphthalein | 2x125ml 4x125ml |
| 8.05.10.0.0250 8.05.10.0.0500 | Chloride Thiocyanate Colorimetric | 2x125ml 4x125ml |
| 8.05.11.0.0250 8.05.11.0.0500 | Cholesterol Liquid (CHOD-POD) | 2x125ml 4x125ml |
| 8.05.12.0.0025 8.05.12.0.0050 | CK-MB Kinetic (Tablets) | 10x2.5ml 20x2.5ml |
| 8.05.13.0.0050 8.05.13.0.0100 | CK-MB Kinetic (Liquid) | 5x10ml 5x20ml |
| 8.05.14.0.0025 8.05.14.0.0050 | CK-NAC Kinetic (Tablets) | 10x2.5ml 20x2.5ml |
| 8.05.15.0.0050 8.05.15.0.0100 | CK-NAC Kinetic (Liquid) | 5x10ml 5x20ml |
| | | |



(CE

| Item Code | Item Description | Sizes |
|----------------------------------|---|-----------------------|
| 8.05.16.0.0250 8.05.16.0.0500 | Creatinine Jaffe Color-Kinetic | 2x125ml 4x125ml |
| 8.05.17.0.0250 8.05.17.0.0500 | Glucose GOD-POD (Liquid) | 2x125ml 4x125ml |
| 8.05.18.0.0020 8.05.18.0.0040 | GOT (AST) IFCC Kinetic (Tablets) | 10x2ml 20x2ml |
| 8.05.19.0.0250 8.05.19.0.0500 | GOT (AST) IFCC Kinetic (Liquid) | 5x50ml 5x100ml |
| 8.05.20.0.0250 8.05.20.0.0500 | GOT (AST) Reitman-Frankel Colorimetric | 2x125ml 2x250ml |
| 8.05.21.0.0020 8.05.21.0.0040 | GPT (ALT) IFCC Kinetic (Tablets) | 10x2ml 20x2ml |
| 8.05.22.0.0250 8.05.22.0.0500 | GPT (ALT) IFCC Kinetic (Liquid) | 5x50ml 5x100ml |
| 8.05.23.0.0200 8.05.23.0.0250 | GPT (ALT) Reitman-Frankel Colorimetric | 2x100ml 2x125m |
| 8.05.24.0.0020 8.05.24.0.0040 | Gamma GT Kinetic, Carboxy Substrate (Tablets) | 10x2ml 20x2ml |
| 8.05.25.0.0250 8.05.25.0.0500 | Gamma GT Kinetic, Carboxy Substrate (Liquid) | 5x50ml 5x100ml |
| 8.05.26.0.0100 8.05.26.0.0200 | HDL Cholesterol Precipitating Reagent | 2x50ml 2x100ml |
| 8.05.27.0.0200 | Iron Ferrozine Colorimetric | 4x50ml |
| 8.05.28.0.0030 8.05.28.0.0060 | LDH IFCC Kinetic (Tablets) | 10x3ml 20x3ml |
| 8.05.29.0.0250 8.05.29.0.0500 | LDH Pyruvate Kinetic UV DGKC (Liquid) | 5x50ml 5x100ml |
| 8.05.30.0.0060 | Lipase Kinetic (Liquid) | 6x10ml |
| 8.05.31.0.0250 8.05.31.0.0500 | Magnesium Calmagite Colorimetric | 2x125ml 4x125ml |
| 8.05.32.0.0250 8.05.32.0.0500 | Phosphorus Phosphomolybdate UV | 2x125ml 4x125ml |
| 8.05.33.0.0050 8.05.33.0.0100 | Potassium Colorimetric | 50 Tests 100 Tests |
| 8.05.34.0.0050 8.05.34.0.0100 | Sodium Colorimetric | 50 Tests 100 Tests |
| 8.05.35.0.0100 | TIBC (Total Iron Binding Capacity) | 100 Tests |
| 8.05.36.0.0250 8.05.36.0.0500 | Total Lipids Phosphovainilline Colorimetric | 2x125ml 4x125ml |
| 8.05.37.0.0250 8.05.37.0.0500 | Total Protein Biuret Colorimetric | 2x125ml 4x125ml |
| 8.05.38.0.0250 8.05.38.0.0500 | Total Protein in CSF | 2x125ml 4x125ml |

CLINICAL CHEMISTRY KITS

| ltem Code | Item Description | Sizes |
|----------------------------------|--|------------------------|
| 8.05.39.0.0250 8.05.39.0.0500 | Triglycerides GPO-POD Colorimetric | 2x125ml 4x125ml |
| 8.05.40.0.0250 8.05.40.0.0500 | Urea Urease-GLDH Kinetic (UV) | 5x50ml 5x100ml |
| 8.05.41.0.0250 8.05.41.0.0500 | Urea Berthelot Colorimetric | 2x125ml 4x125ml |
| 8.05.42.0.0250 8.05.42.0.0500 | Uric Acid Uricase-PAP Colorimetric (Two Reagents) | 2x125ml 4x125ml |
| 8.05.45.0.0250 8.05.45.0.0500 | G6PD Deficiency Qualitative Kit | 250 Tests 500 Tests |
| 8.05.45.1.0250 8.05.45.1.0500 | G6PD Deficiency Qualitative Kit, (with Filter Cards) | 250 Tests 500 Tests |
| 8.05.46.0.0075 8.05.46.0.0150 | G6PD Deficiency Quantitative Kit | 75 Tests 150 Tests |
| 8.05.46.1.0075 8.05.46.1.0150 | G6PD Deficiency Quantitative Kit, (with Filter Cards) | 75 Tests 150 Tests |

| ltem Code | Item Description | Sizes |
|----------------------------------|---|--|
| 8.05.71.0.0250 8.05.71.0.0500 | Uric Acid Enzymatic Colorimetric (Mono Reagents) | 2x125ml 4x125ml |
| 8.05.43.1.0005 | Pathological Control for Clinical Chemistry, Lyophilized, Human Source | 5ml/vial |
| 8.05.44.1.0005 | Normal Control for Clinical Chemistry, Lyophilized, Human Source | 5ml/vial |
| 8.05.47.0.0003 | G6PD Control, Normal Level, (Lyophilized) | 6x0.5ml |
| 8.05.51.0.0100 | HDL Choelsterol, Enzymatic Colorimetric Direct Method | 100 Tests |
| 8.05.52.0.0100 | LDL Cholesterol, Enzymatic Colorimetric Direct Method | 100 Tests |
| 8.40.00.1.0100 | HbA1c Direct Enzymatic Colorimetric Kit | 100 Tests |
| 8.05.73.0.0020 8.05.73.0.0050 | Electrolytes Kit (Sodium, Potassium, Chloride) | 20 Tests for each 50 Tests for each |

ALCOHOL TESTS

Overview

Atlas Medical supplies kits to test for alcohol in urine, saliva and breath. The urine alcohol test is based on the detection of EtG (Ethyl Glucuronide) in urine using a rapid lateral flow immunoassay. Whereas the saliva alcohol test uses a strip with dry chemistry pad that changes color to indicate the level of alcohol in the saliva. The alcohol breath test is a tube with crystals that change color as the subject blows through when alcohol level in breath exceeds a certain limit.

| ltem Code | Item Description | Sizes |
|----------------------------------|---|----------------------|
| 8.25.01.0.0001 8.25.01.0.0025 | Saliva Alcohol Test Strip, Individually Pouched | Bulk 25 Tests/Box |
| 8.25.02.0.0001 8.25.02.0.0020 | Breath Alcohol Test (0.02%), Individually Pouched | Bulk 20 Tests/Box |
| 8.25.03.0.0001 8.25.03.0.0025 | Breath Alcohol Test (0.05%), Individually Pouched | Bulk 25 Tests/Box |
| 8.25.04.0.0001 8.25.04.0.0025 | Breath Alcohol Test (0.08%), Individually Pouched | Bulk 25 Tests/Box |
| 8.25.05.0.0001 8.25.05.0.0020 | Ethyl Glucuronide (ETG) Urine Strip, Individually pouched | Bulk 20 Tests/Box |

Features

The kits are conveniently packed in different kit sizes of 20, 25, 50 and 100 tests per kit.



STAINS FOR HISTOLOGY & MICROBIOLOGY

Œ

MICROBIOLOGY

Overview

Atlas Medical is well known for its range of lab stains for histology and microbiology applications. Atlas Medical stains are made of the highest quality ingredientsto ensure good quality and vivid staining.

Features

The stains come in convenient sizes, but custom sizes are also available.



STAIN PACKS FOR HISTOLOGY

| ltem Code | Item Description | Sizes |
|---------------|--|---------|
| 8.17.003.0300 | Periodic Acid Schiff (PAS) Stain Kit | 3x100ml |
| 8.17.004.0300 | Iron Stain Kit - Perl | 3x100ml |
| 8.17.009.1000 | Gram Stain Pack | 4x250ml |
| 8.17.010.0750 | Cold ZN - Kinyoun Stain Pack | 3x250ml |
| 8.17.011.0750 | ZN Pack Standard | 3x250ml |
| 8.17.015.0500 | Diff-3 Stain Pack | 4x125ml |
| 8.17.016.1000 | Papanicolaou Stain Pack (EA35, EA50, EA65, OG6) | 4X250ml |







| ltem Code | Item Description | Sizes |
|----------------|--|--------------|
| 8.38.00.0.0025 | Blood Culture Bottles, Pediatric Size | 25ml/Bottle |
| 8.38.00.0.0050 | Blood Culture Bottles, Adult Size | 50ml/Bottle |
| 8.36.00.0.0020 | Mycoplasma Culture, Identification, Enumeration and Susceptibility Kit | 20 Tests/Box |

ANTIBIOTIC SENSITIVITY DISCS

Overview

Antibiotic sensitivity is a term used to describe the susceptibility of bacteria to antibiotics. Antibiotic susceptibility testing (AST) is usually carried out to determine which antibiotic will be most successful in treating a bacterial infection in vivo.

Small discs containing antibiotics are placed onto a plate upon which bacteria are growing. If the bacteria are sensitive to the antibiotic, a clear ring, or zone of inhibition, is seen around the disc indicating poor growth.

Features

- Atlas Medical offers a wide range of antibiotics discs at competitive prices.
- ∂ Easy to use.
- ∂ The kit comes with a Cartridge Applicator.
- ∂ Reliable quality.
- Comprehensive range of antibiotics at different concentrations.



| ltem Code | Item Description | Sizes |
|----------------|--|------------|
| 8.39.48.0.0250 | NORFLOXACIN (10 µg) - NX | 5x50 Discs |
| 8.39.49.0.0250 | OFLOXACIN (5 μg) - OF | 5x50 Discs |
| 8.39.50.0.0250 | PEFLOXACIN (5 μg) – PF | 5x50 Discs |
| 8.39.51.0.0250 | PENICILLIN –G (10 IU) – P | 5x50 Discs |
| 8.39.52.0.0250 | PIPERACILLIN (100 µg) – PI | 5x50 Discs |
| 8.39.53.0.0250 | PIPERACILLIN / TAZOBACTAM (100 μg + 10 μg) - PTZ | 5x50 Discs |
| 8.39.54.0.0250 | RIFAMPIN (5 µg) – RIF | 5x50 Discs |
| 8.39.55.0.0250 | ROXITHROMYCIN (30 µg) – RO | 5x50 Discs |

| ltem Code | Item Description | Sizes |
|----------------|---|------------|
| 8.39.01.0.0250 | AMIKACIN (30 µg) – AK | 5x50 Discs |
| 8.39.02.0.0250 | AMOXICILLIN (10 µg) – AX | 5x50 Discs |
| 8.39.03.0.0250 | AMOXICILLIN / CLAVULANATE (20 μg + 10 μg) - AMC | 5x50 Discs |
| 8.39.04.0.0250 | AMPICILLIN (10 µg) – AMP | 5x50 Discs |
| 8.39.05.0.0250 | AMPICILLIN / SULBACTAM (10 μg - 10 μg) – AS | 5x50 Discs |
| 8.39.06.0.0250 | AZITHROMYCIN (15 μg) – AZM | 5x50 Discs |
| 8.39.07.0.0250 | AZTREONAM (30 µg) – AT | 5x50 Discs |
| 8.39.08.0.0250 | CEFACLOR (30 µg) - CF | 5x50 Discs |
| 8.39.09.0.0250 | CEFADROXIL (30 µg) - CD | 5x50 Discs |
| 8.39.10.0.0250 | CEFAZOLIN (30 µg) - CZ | 5x50 Discs |
| 8.39.11.0.0250 | CEFDINIR (5µg) - CDR | 5x50 Discs |
| 8.39.12.0.0250 | CEFIXIME (5 μg) - CFM | 5x50 Discs |
| 8.39.13.0.0250 | CEFOPERAZONE (75 µg) - CPZ | 5x50 Discs |
| 8.39.14.0.0250 | CEFOPERAZONE / SULBACTUM (75 μg + 30 μg) - CS | 5x50 Discs |
| 8.39.15.0.0250 | CEFOTAXIME (30 µg) - CTX | 5x50 Discs |
| 8.39.16.0.0250 | CEFPIROME (30 µg) - CE | 5x50 Discs |
| 8.39.17.0.0250 | CEFPODOXIME (10 µg) – CPD | 5x50 Discs |
| 8.39.18.0.0250 | CEFPROZIL (30 µg) - CPR | 5x50 Discs |
| 8.39.19.0.0250 | CEFTAZIDIME (30 μg) – CAZ | 5x50 Discs |
| 8.39.20.0.0250 | CEFTIZOXIME (30 µg) - CZX | 5x50 Discs |
| 8.39.21.0.0250 | CEFTRIOXONE (30 µg) - CTR | 5x50 Discs |
| 8.39.22.0.0250 | CEFUROXIME (30 µg) - CXM | 5x50 Discs |
| 8.39.23.0.0250 | CEPHALEXIN (30 µg) - CN | 5x50 Discs |
| 8.39.24.0.0250 | CEPHALORIDINE (30 µg)- CH | 5x50 Discs |
| 8.39.25.0.0250 | CEPHALOTHIN (30 µg) - CEP | 5x50 Discs |
| 8.39.26.0.0250 | CHLORAMPHENICOL (30 µg) - C | 5x50 Discs |
| 8.39.27.0.0250 | CIPROFLOXACIN (5 µg) - CIP | 5x50 Discs |
| 8.39.28.0.0250 | CLARITHROMYCIN (15 µg) - CLR | 5x50 Discs |
| 8.39.29.0.0250 | CLINDAMYCIN (2 µg) - CD | 5x50 Discs |
| 8.39.30.0.0250 | CLOXACILLIN (5 μg) - COX | 5x50 Discs |
| 8.39.32.0.0250 | DOXYCYCLINE (30 µg) - DOX | 5x50 Discs |
| 8.39.33.0.0250 | ERYTHROMYCIN (15 µg) - E | 5x50 Discs |
| 8.39.34.0.0250 | FURAZOLIDONE (100 µg) - FZ | 5x50 Discs |
| 8.39.35.0.0250 | GATIFLOXACIN (5 μg) - GAT | 5x50 Discs |
| 8.39.36.0.0250 | GENTAMYCIN (10 μg) - GEN | 5x50 Discs |
| 8.39.38.0.0250 | KANAMYCIN (30 µg) - K | 5x50 Discs |
| 8.39.39.0.0250 | LEVOFOLXACIN (5 µg) - LE | 5x50 Discs |
| 8.39.40.0.0250 | LINCOMYCIN (15 µg) - LN | 5x50 Discs |
| 8.39.41.0.0250 | LINEZOLID (30 µg) - LZ | 5x50 Discs |
| 8.39.42.0.0250 | LOMEFLOXACIN (10 µg) - LOM | 5x50 Discs |
| 8.39.43.0.0250 | MEROPENEM (10 µg) - MRP | 5x50 Discs |
| 8.39.44.0.0250 | MINOCYCLINE (30 µg) - MI | 5x50 Discs |
| 8.39.45.0.0250 | MOXIFLOXACIN (5 µg) - MXF | 5x50 Discs |
| 8.39.46.0.0250 | NALIDIXIC ACID (30 µg) - NA | 5x50 Discs |
| 8.39.47.0.0250 | NITROFURANTOIN (300 µg) - NIT | 5x50 Discs |

ANTIBIOTIC SENSITIVITY DISCS



| ltem Code | Item Description | Sizes |
|----------------|---|------------|
| 8.39.56.0.0250 | SPARFLOXACIN (5 µg) – SPX | 5x50 Discs |
| 8.39.57.0.0250 | STREPTOMYCIN (10 μg) – S | 5x50 Discs |
| 8.39.58.0.0250 | SULFADIAZINE (300 µg) - SD | 5x50 Discs |
| 8.39.59.0.0250 | TEICOPLANIN (30 μg) – TEI | 5x50 Discs |
| 8.39.60.0.0250 | TETRACYCLINE (30 μg) – TE | 5x50 Discs |
| 8.39.61.0.0250 | TICARCILLIN / CLAVULANATE (75 μg + 10 μg)-TCC | 5x50 Discs |
| 8.39.62.0.0250 | TOBRAMYCIN (10 μg) – TOB | 5x50 Discs |
| 8.39.63.0.0250 | TRIMETHOPRIM (5 μg) – TR | 5x50 Discs |
| 8.39.64.0.0250 | VANCOMYCIN (30 µg) – VA | 5x50 Discs |
| 8.39.65.0.0250 | POLYMYXIN-B (300 UNITS) -PB | 5x50 Discs |
| 8.39.66.0.0050 | CEFOXITIN (30 µg) - CX | 1x50 Discs |
| 8.39.67.0.0250 | CEFEPIME (30 µg) - CPM | 5x50 Discs |
| 8.39.69.0.0050 | NOVOBIOCIN (5 µg) -NV | 1x50 Discs |

| Item Code | Item Description | Sizes |
|----------------|-----------------------------|-------------|
| 8.39.70.0.0050 | CARBENICILLIN (100 µg) - CB | 1x50 Discs |
| 8.39.71.0.0050 | BACITRACIN - B (10 Unit) | 1x50 Discs |
| 8.39.72.0.0050 | CEFOXITIN (30 µg) - CX | 20x50 Discs |
| 8.39.76.0.0250 | COLISTIN -CL (10 µg) | 5x50 Discs |
| 8.39.77.0.0250 | IMIPENEM-IPM (10 μg) | 5x50 Discs |
| 8.39.78.0.0250 | OXACILLIN -OX (1 μg) | 5x50 Discs |



ELISA KITS

Features

- The kits feature high sensitivities, simple and robust methods, breakable well strips, quantitative results, ready-to use liquid reagents, and reasonable assay time.
- The assays can be used on most open ELISA manual readers and washers as well as open ELISA auto-analyzers.
- ∂ Kits are packed in sizes of 96 tests.



Overview

Atlas Medical offers a range of Enzyme Linked Immunosorbent Assay (ELISA or EIA) to detect major hormones in the fields of thyroids and fertility in serum.

| | | (CE) |
|----------------|--------------------|----------|
| HORMONE ELIS | SA KITS | |
| ltem Code | Item Description | Sizes |
| 8.07.02.0.0096 | PSA Elisa Kit | 96 Tests |
| 8.07.10.0.0096 | Free PSA Elisa Kit | 96 Tests |

| HORMONE ELISA KITS | | |
|--------------------|------------------------|----------|
| ltem Code | Item Description | Sizes |
| 8.10.01.0.0096 | hCG Elisa Kit | 96 Tests |
| 8.10.03.0.0096 | FSH Elisa Kit | 96 Tests |
| 8.10.04.0.0096 | LH Elisa Kit | 96 Tests |
| 8.10.05.0.0096 | Prolactin Elisa Kit | 96 Tests |
| 8.12.00.0.0096 | T3 Elisa Kit | 96 Tests |
| 8.12.01.0.0096 | T4 Elisa Kit | 96 Tests |
| 8.12.02.0.0096 | TSH Elisa Kit | 96 Tests |
| 8.12.03.0.0096 | Free T4 Elisa Kit | 96 Tests |
| 8.12.04.0.0096 | Free T3 Elisa Kit | 96 Tests |
| 8.11.03.0.0096 | Progesterone Elisa kit | 96 Tests |
| 8.11.04.0.0096 | Testosterone Elisa Kit | 96 Tests |



Œ

ELISA KITS

Overview

Atlas Medical offers a range of Enzyme Linked Immunosorbent Assay (ELISA or EIA) to detect IgG and IgM antibodies against ToRCH (Toxoplasmosis, Rubella, CMV and Herpes I & II) in serum.

- * Different Packaging sizes.
- * Easy to Use
- * High Quality

| Œ |
|---|
| V |

| Item Code | Item Description | Sizes |
|----------------|---|----------|
| 8.13.10.0.0096 | Herpes Simplex 2 IgM (HSV2 IgM) Elisa Kit | 96 Tests |
| 8.13.11.0.0096 | Herpes Simplex 1,2 lgG (HSV1,2 lgG) Elisa Kit | 96 Tests |
| 8.13.12.0.0096 | Herpes Simplex 1,2 IgM (HSV1,2 IgM) Elisa Kit | 96 Tests |

| TORCH ELISA KITS | | | |
|------------------|---|----------|--|
| ltem Code | Item Description | Sizes | |
| 8.13.00.0.0096 | Toxo Plasma Gondii IgG (Toxo IgG) Elisa kit | 96 Tests | |
| 8.13.01.0.0096 | Toxoplasma gondii IgM (Toxo IgM) Elisa kit | 96 Tests | |
| 8.13.02.0.0096 | Rubella IgG Elisa Kit | 96 Tests | |
| 8.13.03.0.0096 | Rubella IgM Elisa Kit | 96 Tests | |
| 8.13.05.0.0096 | Cytomegalovirus IgG (CMV IgG) Elisa Kit | 96 Tests | |
| 8.13.06.0.0096 | Cytomegalovirus IgM (CMV IgM) Elisa Kit | 96 Tests | |
| 8.13.07.0.0096 | Herpes Simplex 1 IgG (HSV1 IgG) Elisa Kit | 96 Tests | |
| 8.13.08.0.0096 | Herpes Simplex 1 IgM (HSV1 IgM) Elisa Kit | 96 Tests | |
| 8.13.09.0.0096 | Herpes Simplex 2 lgG (HSV2 lgG) Elisa Kit | 96 Tests | |



Overview

Atlas Medical offers a range of Enzyme Linked Immunosorbent Assay (ELISA or EIA) to detect a series of infection diseases such as HIV, Hepatitis (A, B, C, D and E) and H. pylori (antigens in feces).

| | | Œ | |
|--------------------------------|---------------------------------------|----------|--|
| INFECTIOUS DISEASES ELISA KITS | | | |
| ltem Code | Item Description | Sizes | |
| 8.14.19.1.0096 | Helicobacter Pylori Antigen Elisa Kit | 96 Tests | |

| INFECTIOUS DISEASES ELISA KITS | | | |
|--------------------------------|---|----------|--|
| Item Code | Item Description | Sizes | |
| 8.14.28.0.0096 | HBsAg Elisa Kit | 96 Tests | |
| 8.14.29.0.0096 | HBsAb Elisa Kit | 96 Tests | |
| 8.14.30.0.0096 | HBcAb Elisa Kit | 96 Tests | |
| 8.14.31.0.0096 | HBeAg Elisa Kit | 96 Tests | |
| 8.14.32.0.0096 | HBeAb Elisa Kit | 96 Tests | |
| 8.14.35.0.0096 | HEV IgM Elisa Kit | 96 Tests | |
| 8.14.38.0.0096 | HCV Ab Elisa Kit | 96 Tests | |
| 8.14.39.0.0096 | HAV lgM Elisa Kit | 96 Tests | |
| 8.14.40.0.0096 | HIV 1,2 Antibody Elisa Kit | 96 Tests | |
| 8.14.43.0.0096 | Syphilis Antibody total (IgG,IgA,IgM) Elisa Kit | 96 Tests | |
| 8.14.44.0.0096 | HIV Ag/Ab Elisa Kit | 96 Tests | |





| Item Code | Item Description | Sizes | |
|----------------|------------------------------------|----------|--|
| 8.07.03.0.0096 | Alpha Feto Protein (AFP) Elisa Kit | 96 Tests | |
| 8.07.08.0.0096 | Ferritin Elisa Kit | 96 Tests | |
| 8.08.00.0.0096 | Troponin I Elisa Kit | 96 Tests | |
| 8.09.00.0.0096 | lgE Elisa Kit | 96 Tests | |
| 8.51.00.0.0096 | 25-OH Vitamin D Elisa Kit | 96 Tests | |
| 8.57.00.0.0096 | Vitamin B12 Elisa Kit | 96 Tests | |
| 8.58.00.0.0096 | Folic Acid Elisa Kit | 96 Tests | |
| 8.06.32.0.0096 | Anti-CRA Elisa Kit | 96 Tests | |

OTHER ELISA KITS

HOME TESTS



Overview

Atlas Medical provides a range of home tests that have been specifically CE marked for OTC use. The range includes fertility tests (Pregnancy, Ovulation and Menopause). The home tests range also includes other medical conditions such as liver function, kidney function, diabetes and urine tract infection. These tests are based on urine reagent strips.

| Screening Kits | | | |
|----------------|---|-------------|--|
| ltem Code | Item Description | Sizes | |
| 70004001 | Atlas Home Diabetes Test | 2 Tests/Box | |
| 70021001 | Atlas Home Urinary Tract Infection Test | 2 Tests/Box | |
| 70022001 | Atlas Home Kidney Function Test | 2 Tests/Box | |
| 70023001 | Atlas Home Liver Function Test | 2 Tests/Box | |

Features

- These tests come in cassette, midstream and strip formats.
- The screening kits come with 2 individually pouched strips and easy to read instructions for use.
- All kits are packed in attractively designed boxes with various languages.
- Atlas Medical also supplies these kits under OEM arrangements.
- Screening bundle including (UTI, Kidney, Liver, Diabetes) is available , Family planning kit (Pregnancy and Ovulation) is also available .



- Simply dip the test strip in urine for one second.
- Wait for 30 60 seconds.
- Compare the colors on the strip to the color chart on the box.



1459

| Fertility Kits | | |
|----------------|---|-------------|
| Item Code | Item Description | Sizes |
| 70171001 | Atlas Home Pregnancy Test Cassette | 1 Test/Box |
| 70172001 | Atlas Home Pregnancy Test Midstream | 1 Tests/Box |
| 70174001 | Atlas Home Ovulation Test Cassette | 5 Tests/Box |
| 70175001 | Atlas Home Ovulation Test Midstream | 3 Tests/Box |
| 70177001 | Atlas Home Menopause Test Cassette | 1 Test/Box |
| 70178001 | Atlas Home Menopause Test Midstream | 1 Test/Box |
| 70180001 | Atlas Home Pregnancy Test Strip (With Handle) | 1 Test/Box |
| 70170001 | Atlas Home Pregnancy Test Strip | 1 Test/Box |

BLOOD GLUCOSE MONITORING SYSTEMS

Overview

Testing your blood glucose regularly helps you better manage your diabetes. Reliance [™] by Atlas Medical, uses the latest blood glucose sensor technologies to offer you the most accurate and reliable results for the peace of mind you need. Atlas Medical offers these systems in strips which includes Gold Electrodes.



Features

- Reliance Gold ™ by Atlas Medical, uses the latest blood glucose gold sensor technology to offer the most accurate and reliable results.
- ∂ Test time required is 5 Seconds.
- ∂ Required sample volume is 0.9µl
- Test result range is between 10 600 mg/dl (0.6 -33.3 mmol/L)

| Reliance Gold | | |
|----------------------------------|---------------------------------------|--------------------------------------|
| Item Code | Item Description | Sizes |
| 8.52.00.0.0001 | Reliance Gold Glucometer Pack | 1 Pack |
| 8.52.00.0.0025 8.52.00.0.0050 | Strips for Reliance Gold Glucometer | 25 Strips/Bottle 50 Strips/Bottle |
| 8.52.00.1.0001 | Reliance Gold Glucometer (Divce only) | 1 Divce only |





CERTIFICATES

| | CERTIFICAT CERTIFICATE OF REGISTRATION N° 36655 rev.1 |
|---|--|
| GMED certifie que le système de mai | nagement de la qualité développé par |
| GMED certifies that the quality n | ranagement system developed by |
| ATLAS MED | ICAL GmbH |
| Ludwig-Ert | hard-Ring 3 |
| 15827 Blankenfelde | Mahlow GERMANY |
| pour les | activités |
| nor the . | admines administration |
| Conception et développement, fabrication et ven | te de dispositifs médicaux de diagnostic in vitro . |
| | |
| Design and Development, Manufacturing an | d Sales of in vitro diagnostic medical devices |
| | |
| réalisées sur | le(s) site(s) de |
| personned on r | ne location(a) at |
| Voir ad | dendum |
| See at | Cheroflatt / |
| est conforme aux exigences | des normes internationales |
| complies with the requirement | s of the international standards |
| ISO 134 | 85: 2016 |
| Début de validité / Effective date October 9th, 2020 (in | studied) |
| Valable jusqu'au / Exploy date : October 8th, 2023 (incl Etabli le : / Issued on : October 8th, 2020 | uded) |
| ofrac | On behalt of the President |
| S | Beatrice LYS Technical Director |
| Color#ission of select an object in contraction cards | > 7 The out-Woods is manuf according to the same of QUED (antification |
| Proceeding to the second | |
| GMED • Société par Actions Simplifiée au capital de Sége social - 1, rue Baston Bossuer - 75015 Paris • 1 | 380 000 C • Organisme Notifié/Notified Body n° 0659 |
| 12 | |
| | |
| | |
| | |
| ISO 134 | 485 |
| | |



Full Quality Assurance Certificate.

ATTESTATION/ CERTIFICATE N* 33544 Officitie & Paris le 19 mars 2028 Assued in Paris In March 19h, 2020 GMED ATTESTATION CE / EC CERTIFICATE Examun CE de la Conception (du produit) / 8/2 design Examination (of the product EXE IV point 4 Directive 36/79/CE relative aux dispositifs middleaux de diagnositic in ANESE IV actives 7.01/8/CT/19/6 0/2014/CE. ATLAS MEDICAL GmbH Ludwig-Erhard-Ring 3 15827 Blankenfelde-Mahlow GERMANY ent, production, et commercialisation de dispositifs médicaux destinés au diagnostic in vitro. Annexe II liste A : détermination des groupes sanguins. Design, production and sales of medical devices for in vitro diagnostic Annex II vist A : blood grouping determination. an du(des) dispositif(s) / Identification of device(s) ATLAS Anti-A, Anti-B, Anti-AB, Anti-D Me See addendum IED atteste qu'à l'examen des résultats figurant dans le rapport référence P172376, le(s) pro-nt) conforme(s) aux exigences de l'annese i de la directive 98/76/CE. tifles that, on the basis of the PB779'EC, where A constrait an and the first state that we idité / Effective date : March 18th, 2020 (in u'au / Expiry date : October 8th, 2022 (inc Banhour 14/5 GMED of the Pres Gall D - 33544 rwv 2 Mobile to certifical 3 ISSE-T GMED • Societé par Actions Simplifiée au capital de 300 000 6 Siège social : 1, rue Gaston Boissier - 75015 Paris • Tél, =01 40

Blood Grouping CE Certificate



CONTACT US



FIND US ON:





Amazon.co.uk/Atlas Medical

Email: Info@atlas-medical.com

INTERNATIONAL PRESENCE









STREPTOCOCCAL GROUPING SLIDE LATEX TEST

A qualitative latex agglutination test for the Detection of Streptococcal groups A, B, C, D, F and G

IVD For In-Vitro diagnostic and professional use only



INTENDED USE

For the qualitative detection of streptococcal groups A, B, C, D, F and G based on latex agglutination.

INTRODUCTION & PRINCIPLES

ATLAS Streptococcal test uses an enzyme extraction procedure to release Carbohydrate antigen from Streptococcal cell walls.

The antigens are detected using specific antibodies to groups A, B, C, D, F and G Lancefield. These antibodies are coated on latex particles. When the antigen extract is mixed with the latex reagent, agglutination will occur. The agglutination appears as a visible clumping and can be seen macroscopically.

MATERIALS

MATERIALS PROVIDED

- Group A, B, C, D, F and G latex reagents.
- Extraction Enzyme dried.
- Positive control.
- Test slide.
- Stirring Sticks.
- Package Insert.

MATERIALS NEEDED BUT NOT PROVIDED

- Water bath.
- Pipette to deliver 50ul.
- Timer and test tubes.

- 1. Prior to use, the Latex reagent should be mixed well to obtain a uniform suspension of the Latex.
- 2. This kit should be stored in an upright position and refrigerated between 2 to 8°C. Never Freeze.
- 3. Use a fresh disposable slide and mixer for each test.
- 4. Always ensure an acceptable performance of the kit by performing the test on the negative and the Positive controls before using the kit.
- 5. The extraction procedure may not kill all organisms; therefore carefully dispose the materials into disinfectant or by autoclaving.

PREPARING THE EXTRACTION ENZYME

The Extract enzyme in this kit comes in two vials dried.Reconstitute with 10ml distilled or deionized water. Once reconstituted store at 2-8°C for a maximum of 3 months or a liquot in 0.4ml volumes and store at -20°C for up to a year.

SAMPLE PREPARATION

Cultures

CE

Note colonial characteristics, hemolysis, and cell morphology before starting the test. Ensure that the organisms to be tested are Gram-positive and catalase-negative. Any blood agar plate culture yielding 2-6 well-separated colonies maybe used, they should have been inoculated from a pure culture of the organism.

PROCEDURES

- 1. Using a sterile bacteriological loop, pick no more than 6 colonies of streptococci (avoiding other types of colony on the plate) and emulsify them in 0.4 ml extraction enzyme. (If a broth culture is to be grouped, pipette 0.1ml of an overnight culture into 0.4ml extraction enzyme).
- 2. Incubate the mixture in a water bath at 37°C for 10 minutes. Shake the tubes vigorously after 5 minutes incubation. Longer incubation period may lead to false positive results.
- 3. Re-suspend the latex reagents by gentle agitation.
- 4. Dispense 1 drop of each latex into the appropriate labeled circle on the test slide.

- 5. Using a pipette, place 50ul of the extract to each drop of latex reagent, and mix the contents of each circle with a separate mixing stick.
- 6. Gently rock the slide for one minute.
- 7. Read the result in normal light and observe for any agglutination.

READING THE RESULT

POSITIVE: If Agglutination appears within one minute. NEGATIVE: If Agglutination does not appear within one minute.

PROCEDURE LIMITATION

- This test provides a presumptive diagnosis. Physicians should evaluate all clinical and laboratory findings before making a definitive diagnosis.
- Faint granularity may be seen in some negative patterns, this should be disregarded.

PERFORMANCE CHARACTERISTICS

| | | Test reagent | |
|---------------------|---|--------------|----|
| | | + | - |
| Reference method | + | 607 | 55 |
| | - | 0 | 24 |

Sensitivity 607/662 = 92% Specificity 24/24 = 100%

INTERNAL QUALITY CONTROL

A positive control is provided and should be used to verify that the latex reagents are working satisfactorily under test conditions.

Periodically check the following:

- 1. The test reagents agglutinate with a known reference Streptococcus strain
- 2. The test reagents do not auto agglutinate in normal saline solution.

REFERENCES

1. Lancefield, R.c., (1938)Proc. Soc.Exp. Bio. Med. 38, 473

2. Harvey,C.L.,Mcillmurray, M.B. (1984) Eur.J. Clin. Microbiol, 3.6,526

 Facklam,R.R., (1980) "Manual of Clinical Microbiology"
 Edn., American Society for Microbiology, Washington, DC, pp 88-110.

ATLAS Medical

Unit 4, William James House Cowley Rd, Cambridge, CB4 0WX Tel: ++44 (0) 1223 858 910 Fax: ++44 (0) 1223 858 524 PPI082A01 Rev C (09.09.2015)



Clostridium difficile Antigen GDH Rapid Test Cassette (Feces)

INTENDED USE

The Clostridium difficile Antigen GDH Rapid Test Cassette (Feces) a rapid visual immunoassay for the qualitative, presumptive detection of Clostridium difficile Glutamate Dehydrogenase (GDH) in human fecal specimens, as a screening test and as an aid in the diagnosis of Clostridium difficile infection.

INTRODUCTION

Clostridium difficile (C.difficile), a Gram-positive spore bearing anaerobic bacterium is the major aetiological agent of diarrhoea and colitis associated with antibiotics. C. difficile is the most common cause of health care-associated diarrhoea in developed countries and is a major source of nosocomial morbidity and mortality worldwide

Disease due to C. difficile develops when the organism is allowed to proliferate in the colon, most commonly after antibiotic use has eliminated competing flora. C.difficile can release two high-molecular-weight toxins, toxin A and toxin B, which are responsible for the clinical manifestations, which range from mild, self-limited watery diarrhoea to fulminant pseudomembranous colitis, toxic megacolon and death.

Clostridium difficile Glutamate Dehydrogenase (GDH) is an enzyme produced in large quantities by all toxigenic and non-toxigenic strains, making it an excellent marker for the organism.

The toxigenic culture (TC) is used as the gold standard technique to determine Clostridium difficile infection. This method consists in culture and isolation of C. difficile from faeces, followed by toxin testing of the isolate, a labour-intensive assay to obtain a result.

The Clostridium difficile Antigen GDH Rapid Test Cassette (Feces) is a rapid test to qualitatively detect Clostridium difficile Glutamate Dehydrogenase (GDH) in human feces in 10 minutes. The test can be performed by untrained or minimally skilled personnel, without cumbersome laboratory equipment.

PRINCIPLE

The Clostridium difficile Antigen GDH Rapid Test Cassette (Feces) is a qualitative lateral flow immunoassay for the detection of Clostridium difficile GDH in human feces samples. The membrane is pre-coated with monoclonal antibodies against GDH on the test line region. During testing, the sample reacts with the particle coated with anti-GDH antibodies, which were pre-dried on the test strip. The mixture moves upward on the membrane by capillary action. If there is sufficient Clostridium difficile GDH in the specimen, a colored band will form at the test region of the membrane. The presence of this colored band indicates a positive result, while its absence indicates a negative result. The appearance of a colored band at the control region serves as a procedural control, indicating that the proper volume of specimen has been added and membrane wicking has occurred. If the control line does not appear, the test result is not valid.

PRODUCT CONTENTS

The Clostridium difficile Antigen GDH Rapid Test Cassette (Feces) containing Clostridium difficile GDH-specific antibodies coated particles and GDH-specific antibodies coated on the membrane.

MATERIALS SUPPLIED

20 Test cassettes

20 Extraction tubes with buffer

1 Package insert

MATERIAL REQUIRED BUT NOT PROVIDED

Timer

STORAGE AND STABILITY

The kit can be stored at room temperature or refrigerated $(2-30^{\circ}C)$. The test Cassette is stable through the expiration date printed on the sealed pouch. The test Cassette must remain in the sealed pouch until use. DO NOT FREEZE. Do not use beyond the expiration date.

WARNINGS AND PRECAUTIONS

- 1. For professional in vitro diagnostic use only.
- 2. Do not use after the expiration date indicated on the package. Do not use the test if the foil pouch is damaged.
- 3. Test is for single use only. Do not re- use under any circumstances.
- 4. Avoid cross-contamination of specimens by using a new extraction tube for each specimen obtained.

- 5. Read the entire procedure carefully prior to testing.
- 6. Do not eat, drink or smoke in any area where specimens and kits are handled.
- 7. Handle all specimens as if they contain infectious agents. Observe established precautions against microbiological hazards throughout the procedure and follow standard procedures for the proper disposal of specimens. Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are assayed.
- 8. Do not interchange or mix reagents from different lots. Do not mix solution bottle caps.
- 9. Humidity and temperature can adversely affect results.
- 10. Do not perform the test in a room with strong air flow, ie. electric fan or strong airconditioning.

SPECIMEN COLLECTION AND PREPARATION

- The Clostridium difficile Antigen GDH Rapid Test Cassette (Feces) is intended for use with human fecal specimens only.
- Stool samples should be collected in clean containers. The samples can be stored in the refrigerator (2-8°C) for 7 days prior to testing. For longer storage, maximum 1 year, the specimen must be kept frozen a-20°C. In this case, the sample will be totally thawed and brought to room temperature before testing. Ensure only the amount needed is thawed because of freezing and defrosting cycles are not recommended. Homogenise stool samples as thoroughly as possible prior to preparation.

SPECIMEN PREPARATION

Consider any materials of human origin as infectious and handle them using standard biosafety procedures.

- 1. Collect a random sample of feces in a clean, dry receptacle. Best results will be obtained if the assay is performed within 6 hours after collection.
- 2. Unscrew and remove the dilution tube applicator. Be careful not to spill or spatter solution from the tube. Collect specimens by inserting the applicator stick into at least 5 different sites of the feces to collect approximately 50 mg of feces (equivalent to 1/4 of a pea).
- 3. For liquid specimens: Hold the pipette vertically, aspirate fecal specimens, and then transfer 3 drops (approximately 80µL) into the specimen collection tube containing the extraction buffer.
- 4. Replace the stick in the tube and tighten securely.
- 5. Shake the specimen collection tube vigorously to mix the specimen and the extraction buffer. Specimens prepared in the specimen collection tube may be stored for 6 months at -20° C if not tested within 1 hour after preparation.

TEST PROCEDURE

Bring tests, specimens, reagents and/or controls to room temperature (15-30°C) prior to testing.

- 1. Remove the test from the sealed pouch and place it on a clean, level surface. Label the device with patient or control identification. For best results, the assay should be performed immediately after opening the foil pouch.
- 2. Holding the sample collection device upright, carefully break off the tip of collection device.
- 3. Squeeze 3 drops (\sim 90µL) of the sample solution in the sample well of the device and start the timer.
- 4. Wait for the colored line(s) to appear. Read results in 10 minutes. Do not interpret the result after 10 minutes.


INTERPRETATION OF RESULTS

(Please refer to the illustration above)

Positive: Two lines appear. One colored line should be in the control line region (C) and another apparent colored line should be in the test line region (T).

Negative: One colored line appears in the control line region(C). No line appears in the test line region (T).

Invalid: Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test cassette. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

QUALITY CONTROL

A procedural control is included in the test. A red line appearing in the control region (C) is the internal procedural control. It confirms sufficient specimen volume and correct procedural technique. Control standards are not supplied with this test. However, it is recommended that positive and negative controls are sourced from a local competent authority and tested as a good laboratory practice, to confirm the test procedure and verify the test performance.

LIMITATIONS

- 1. The Clostridium difficile Antigen GDH Rapid Test Cassette (Feces) will only indicate the presence of parasites in the specimen (qualitative detection) and should be used for the detection of Clostridium difficile GDH in faces specimens only. Neither the quantitative value nor the rate of increase in antigen concentration can be determined by this test.
- 2. An excess of sample could cause wrong results (brown bands appear). Dilute the sample with the buffer and repeat the test.
- 3. The Clostridium difficile Antigen GDH Rapid Test Cassette (Feces) should be used only with samples from human faeces. The use of other samples has not been established. The quality of the test depends on the quality of the sample; proper faecal specimens must be obtained
- 4. Positive results determine the presence of Clostridium difficile in faecal samples; never the less it can be due to toxigenic or non-toxigenic strains of Clostridium difficile. A positive result should be flowed up with additional laboratory techniques (toxigenic culture) to determine the strain.
- 5. A negative result is not meaningful because of it is possible the antigen concentration in the stool samples is lower than the detection limit value. If the symptoms or situation still persist, a Clostridium difficile determination should be carried out, on a sample from an enrichment culture.
- 6. As with all diagnostic tests, a definitive clinical diagnosis should not be based on the results of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated.

PERFORMANCE CHARACTERISTICS

1. Clinical Sensitivity, Specificity and Accuracy

The Clostridium difficile Antigen GDH Rapid Test Cassette (Feces) has been evaluated with specimens obtained from patients, ELISA method was used as the reference method. The results show that the Clostridium difficile Antigen GDH Rapid Test Cassette (Feces) has a high overall relative accuracy.

Table 1: The Clostridium difficile Antigen GDH Rapid Test vs ELISA

| Method | | ELISA | | Total Results |
|------------------------|----------|----------|----------|---------------|
| Clostridium difficile | Results | Positive | Negative | |
| Antigen GDH Rapid Test | Positive | 62 | 1 | 63 |
| Cassette | Negative | 0 | 50 | 50 |
| Total Resul | ts | 62 | 51 | 113 |

Relative Sensitivity: 100% Relative Specificity: 98.0% Accuracy: 99.1%

2. Analytical Sensitivity

The Clostridium difficile Antigen GDH Rapid Test Cassette (Feces) was determined by testing serial dilutions of recombinant antigen. The Clostridium difficile Antigen GDH Rapid Test Cassette (Feces) can detect the levels of Clostridium difficile GDH recombinant antigen as low as 1ng/mL.

3. Cross-Reactivity

Cross-reactivity to samples positive for the following pathogens was tested and found to be negative: Campylobacter

| | Salmonella enterittais | Shigella dysenteriae | |
|------------------------|------------------------|-------------------------|--|
| Campylobacter jejuni | Saimonella paratypni | Snigella flexneri | |
| E. Coli 0157: H7 | Salmonella typhi | Shigella sonnei | |
| H. Pylori | Salmonella typhimurium | Staphliococcus aureus | |
| Listeria monocytogenes | Shigella boydii | Yersinia enterocolitica | |
| DEFEDENCE | | | |

- 1. Knoop, F.C. et al.: Clostridium difficile: Clinical disease and diagnosis. Clin.Microbiol. Rev. (1993); 6: 251-265.
- 2. Kelly, C.P. et al.: Clostridium difficile Colitis. New Engl. J. Med. (1994); 330: 257-262.
- Sullivan, N.M. et al.: Purification and characterization of toxins A and B of Clostridium difficile. Infect. Immun. 3. (1982); 35: 1032-1040.
- 4. McDonald, L.C. et al.: An epidemic, toxin gene-variant strain of Clostridium difficile.N. Engl. J. Med. (2005); 353: 23.
- 5. Loo, V.G. et al.: A predominantly clonal multi-institutional outbreak of Clostridium difficile-associated diarrhea with high morbidity and mortality. N. Engl. J. Med. (2005); 353.23.
- 6. Bartlett, J.G., Gerding, D.N.: Clinical recognition and diagnosis of Clostridium difficile

infection. CID (2008); 46(Suppl. 1): 12-18.

INDEX OF SYMBOLS

| Ē | Consult instructions for use | $\overline{\mathbb{Y}}$ | Tests per kit | EC REP | Authorized Representative |
|------|----------------------------------|-------------------------|---------------|--------|---------------------------|
| IVD | For in vitro diagnostic use only | R | Use by | 8 | Do not reuse |
| 2°C- | Store between 2~30°C | LOT | Lot Number | REF | Catalog# |

Zhejiang Orient Gene Biotech Co., Ltd

Address: 3787#, East Yangguang Avenue, Dipu Street, Anji 313300, Huzhou, Zhejiang, China Tel: +86-572-5226111 Fax: +86-572-5226222 Website: www.orientgene.com

EC REP CMC Medical Devices & Drugs S.L

C/Horacio Lengo Nº 18 CP 29006, Málaga-Spain Tel: +34951214054 Fax: +34952330100 Email-info@cmcmedicaldevices.com

REF GCCD(GDH)-602a

Clostridium difficile Toxin A&B Rapid Test Cassette (Feces)

INTENDED USE

The Clostridium difficile Toxin A&B Rapid Test Cassette (Feces) a rapid visual immunoassay for the qualitative, presumptive detection of Clostridium difficile Toxin A&B in human fecal specimens, as a screening test and as an aid in the diagnosis of Clostridium difficile infection.

INTRODUCTION

Clostridium difficile (C. difficile), a Gram-positive spore bearing anaerobic bacterium is the major aetiological agent of diarrhoea and colitis associated with antibiotics. C. difficile is the most common cause of health care-associated diarrhoea in developed countries and is a major source of nosocomial morbidity and mortality worldwide.

Disease due to C. difficile develops when the organism is allowed to proliferate in the colon, most commonly after antibiotic use has eliminated competing flora. C. difficile can release two high-molecular-weight toxins, toxin A and toxin B, which are responsible for the clinical manifestations, which range from mild, self-limited watery diarrhoea to fulminant pseudomembranous colitis, toxic megacolon and death.

Clostridium difficile Glutamate Dehydrogenase (GDH) is an enzyme produced in large quantities by all toxigenic and non-toxigenic strains, making it an excellent marker for the organism.

The toxigenic culture (TC) is used as the gold standard technique to determine Clostridium difficile infection. This method consists in culture and isolation of C. difficile from feces, followed by toxin testing of the isolate, a labor-intensive assay to obtain a result.

The Clostridium difficile Toxin A&B Rapid Test Cassette (Feces) is a rapid test to qualitatively detect Clostridium difficile Toxin A&B in human feces in 10 minutes. The test can be performed by untrained or minimally skilled personnel, without cumbersome laboratory equipment.

PRINCIPLE

The Clostridium difficile Toxin A&B Rapid Test Cassette (Feces) is a qualitative lateral flow immunoassay for the detection of Clostridium difficile Toxin A&B in human feces samples. The membrane is pre-coated with monoclonal antibodies against Toxin A on the A test line region and monoclonal antibodies against Toxin B on the B test line region. During testing, the sample reacts with the particle coated with anti-Toxin A and anti-Toxin B antibodies, which were pre-dried on the test strip. The mixture moves upward on the membrane by capillary action. If there is sufficient Clostridium difficile Toxin or Toxin B in the specimen, a colored band will form at the test region of the membrane. The presence of this colored band indicates a positive result, while its absence indicates a negative result. The appearance of a colored band at the control region serves as a procedural control, indicating that the proper volume of specimen has been added and membrane wicking has occurred. If the control line does not appear, the test result is not valid.

PRODUCT CONTENTS

The Clostridium difficile Toxin A&B Rapid Test Cassette (Feces) containing Clostridium difficile Toxin A and Toxin B antibodies coated particles and Toxin A-specific antibodies and Toxin B-specific antibodies coated on the membrane.

MATERIALS SUPPLIED

20 Test cassettes

20 Extraction tubes with buffer

1 Package insert

MATERIAL REQUIRED BUT NOT PROVIDED

Timer

STORAGE AND STABILITY

The kit can be stored at room temperature or refrigerated (2-30°C). The test cassette is stable through the expiration date printed on the sealed pouch. The test cassette must remain in the sealed pouch until use. DO NOT FREEZE. Do not use beyond the expiration date.

WARNINGS AND PRECAUTIONS

1. For professional in vitro diagnostic use only.

- 2. Do not use after the expiration date indicated on the package. Do not use the test if the foil pouch is damaged.
- 3. Test is for single use only. Do not re-use under any circumstances.

- 4. Avoid cross-contamination of specimens by using a new extraction tube for each specimen obtained.
- 5. Read the entire procedure carefully prior to testing.

CE

- 6. Do not eat, drink or smoke in any area where specimens and kits are handled.
- 7. Handle all specimens as if they contain infectious agents. Observe established precautions against microbiological hazards throughout the procedure and follow standard procedures for the proper disposal of specimens. Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are assayed.
- 8. Do not interchange or mix reagents from different lots. Do not mix solution bottle caps.
- 9. Humidity and temperature can adversely affect results.
- 10. Do not perform the test in a room with strong air flow, ie. electric fan or strong airconditioning.

SPECIMEN COLLECTION AND PREPARATION

- The Clostridium difficile Toxin A&B Rapid Test Cassette (Feces) is intended for use with human fecal specimens
 only.
- Stool samples should be collected in clean containers. The samples can be stored in the refrigerator (2-8°C) for 7 days prior to testing. For longer storage, maximum 1 year, the specimen must be kept frozen a-20°C. In this case, the sample will be totally thawed and brought to room temperature before testing. Ensure only the amount needed is thawed because of freezing and defrosting cycles are not recommended. Homogenise stool samples as thoroughly as possible prior to preparation.

SPECIMEN PREPARATION

Consider any materials of human origin as infectious and handle them using standard biosafety procedures.

- 1. Collect a random sample of feces in a clean, dry receptacle. Best results will be obtained if the assay is performed within 6 hours after collection.
- 2. Unscrew and remove the dilution tube applicator. Be careful not to spill or spatter solution from the tube. Collect specimens by inserting the applicator stick into at least 5 different sites of the feces to collect approximately 50 mg of feces (equivalent to 1/4 of a pea).
- 3. For liquid specimens: Hold the pipette vertically, aspirate fecal specimens, and then transfer 3 drops (approximately 80 µL) into the specimen collection tube containing the extraction buffer.
- 4. Replace the stick in the tube and tighten securely.
- 5. Shake the specimen collection tube vigorously to mix the specimen and the extraction buffer. Specimens prepared in the specimen collection tube may be stored for 6 months at -20°C if not tested within 1 hour after preparation.

TEST PROCEDURE

Bring tests, specimens, reagents and/or controls to room temperature (15-30°C) prior to testing.

- 1. Remove the test from the sealed pouch and place it on a clean, level surface. Label the device with patient or control identification. For best results, the assay should be performed immediately after opening the foil pouch.
- 2. Holding the sample collection device upright, carefully break off the tip of collection device.
- 3. Squeeze 3 drops (~90 µL) of the sample solution in the sample well of the device and start the timer.
- 4. Wait for the colored line(s) to appear. Read results in 10 minutes. Do not interpret the result after 10 minutes.



INTERPRETATION OF RESULTS

(Please refer to the illustration above)

1. Positive:

1.1 Toxin A Positive:

The presence of two lines as control line (C) and A test line within the result window indicates a positive result for Toxin A.

1.2 Toxin B Positive:

The presence of two lines as control line (C) and B test line within the result window indicates a positive result for Toxin B.

1.3 Toxin A& B Positive:

The presence of three lines as control line (C), A test line and B test line within the result window indicates a positive result for both Toxin A and Toxin B.

2. Negative:

One colored line appears in the control line region (C). No line appears in the test line region (T).

3. Invalid:

If the control band (C) is not visible within the result window after performing the test, the result is considered invalid. Some causes of invalid results are because of not following the directions correctly or the test may have deteriorated beyond the expiration date. It is recommended that the specimen be re-tested using a new test. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

QUALITY CONTROL

A procedural control is included in the test. A red line appearing in the control region (C) is the internal procedural control. It confirms sufficient specimen volume and correct procedural technique. Control standards are not supplied with this test. However, it is recommended that positive and negative controls are sourced from a local competent authority and tested as a good laboratory practice, to confirm the test procedure and verify the test performance.

LIMITATIONS

- 1. The Clostridium difficile Toxin A&B Rapid Test Cassette (Feces) will only indicate the presence of parasites in the specimen (qualitative detection) and should be used for the detection of Clostridium difficile Toxin A&B in feces specimens only. Neither the quantitative value nor the rate of increase in antigen concentration can be determined by this test.
- 2. An excess of sample could cause wrong results (brown bands appear). Dilute the sample with the buffer and repeat the test.
- 3. The Clostridium difficile Toxin A&B Rapid Test Cassette (Feces) should be used only with samples from human feces. The use of other samples has not been established. The quality of the test depends on the quality of the sample; proper fecal specimens must be obtained.
- 4. A negative result is not meaningful because of it is possible the antigen concentration in the stool samples is lower than the detection limit value. If the symptoms or situation still persist, a Clostridium difficile determination should be carried out, on a sample from an enrichment culture.
- 5. As with all diagnostic tests, a definitive clinical diagnosis should not be based on the results of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated.

PERFORMANCE CHARACTERISTICS

1. Clinical Sensitivity, Specificity and Accuracy

The Clostridium difficile Toxin A&B Rapid Test Cassette (Feces) has been evaluated with specimens obtained from patients. ELISA method was used as the reference method. The results show that the Clostridium difficile Toxin A&B Rapid Test Cassette (Feces) has a high overall relative accuracy.

Table 1: The Clostridium difficile Toxin A Rapid Test vs ELISA

| Method | | ELISA | | Total Dasults |
|-------------------------|----------|----------|----------|---------------|
| | Results | Positive | Negative | Iotal Results |
| A&B Rapid Test Cassette | Positive | 43 | 1 | 44 |
| Add Rapid Test Casselle | Negative | 0 | 69 | 69 |
| Total Results | | 43 | 70 | 113 |

Relative Sensitivity: 100%

Relative Specificity: 98.6%

Accuracy: 99.1%

Table 2: The Clostridium difficile Toxin B Rapid Test vs ELISA

| Method | | ELISA | | Total Pagulta |
|-------------------------|----------|----------|----------|---------------|
| | Results | Positive | Negative | Iotal Results |
| A&B Rapid Test Cassette | Positive | 36 | 1 | 37 |
| Add Rapid Test Casselle | Negative | 0 | 76 | 76 |
| Total Results | | 36 | 77 | 113 |

Relative Sensitivity: 100% Relative Specificity: 98.6% Accuracy: 99.1%

2. Analytical Sensitivity

The Clostridium difficile Toxin A&B Rapid Test Cassette (Feces) was determined by testing serial dilutions of recombinant antigen. Detection limit values of Clostridium difficile Toxin A&B are 2 ng/mL for Toxin A and 1 ng/mL for Toxin B.

3. Cross-Reactivity

Cross-reactivity to samples positive for the following pathogens was tested and found to be negative:

| Salmonella enteritidis | Shigella dysenteriae |
|------------------------|---|
| Salmonella paratyphi | Shigella flexneri |
| Salmonella typhi | Shigella sonnei |
| Salmonella typhimurium | Staphliococcus aureus |
| Shigella boydii | Yersinia enterocolitica |
| | Salmonella enteritidis Salmonella paratyphi Salmonella typhi Salmonella typhimurium Shigella boydii |

REFERENCE

- 1. Knoop, F.C. et al.: Clostridium difficile: Clinical disease and diagnosis. Clin. Microbiol. Rev. (1993); 6: 251-265.
- 2. Kelly, C.P. et al.: Clostridium difficile Colitis. New Engl. J. Med. (1994); 330: 257-262.
- 3. Sullivan, N.M. et al.: Purification and characterization of toxins A and B of Clostridium difficile. Infect. Immun. (1982); 35: 1032-1040.
- 4. McDonald, L.C. et al.: An epidemic, toxin gene-variant strain of Clostridium difficile. N. Engl. J. Med. (2005); 353: 23.
- 5. Loo, V.G. et al.: A predominantly clonal multi-institutional outbreak of Clostridium difficile-associated diarrhea with high morbidity and mortality. N. Engl. J. Med. (2005); 353.23
- Bartlett, J.G., Gerding, D.N.: Clinical recognition and diagnosis of Clostridium difficile infection. CID (2008); 46 (Suppl. 1): 12-18.

INDEX OF SYMBOLS

| Ē | Consult instructions for use | ¥ | Tests per kit | EC REP | Authorized Representative |
|------|----------------------------------|-----|---------------|--------|---------------------------|
| IVD | For in vitro diagnostic use only | R | Use by | ଷ | Do not reuse |
| 2°C- | Store between 2~30°C | LOT | Lot Number | REF | Catalog# |

A Zhejiang Orient Gene Biotech Co.,Ltd Address: 3787#, East Yangguang Avenue, Dipu Street, Anji 313300, Huzhou, Zhejiang, China Tel: +86-572-5226111 Fax: +86-572-5226222 Website: www.orientgene.com

 EC
 REP
 CMC Medical Devices & Drugs S.L

 C/Horacio Lengo Nº 18 CP 29006, Málaga-Spain
 Tel: +34951214054
 Fax: +34952330100

 Email-info@cmcmedicaldevices.com
 Email-info@cmcmedicaldevices.com
 Fax: +34952330100

REF GCCD(Toxin A&B)-602a



Optochin Discs

DD009

Optochin Discs are used for identification and differentiation of Streptococcus pneumoniae and Viridans Streptococci.

Directions

Prepare Soyabean Casein Digest Agar (M290) w/blood or Blood Agar Base (M073) plates and streak pure culture of organism to be tested across one half of the plate. Streak a known Pneumococcus culture across the other half of the plate as positive control. Immediately place Optochin discs in the centre of the two halves of the plate and incubate at 35-37°C for 18-24 hours. Observe for zone of inhibition around the discs.

Principle And Interpretation

Alpha haemolytic (viridans) streptococci and Pneumococcus (*Streptococcus pneumoniae*) cannot be easily distinguished on Blood Agar plates as pneumococci strain shows partial clearing of blood and greenish discolouration (a-hemolysis). Optochin is inhibitory for pneumococcal growth whereas other streptococci strains show good growth or a very small zone of inhibition. Bowers and Jeffries have shown a correlation between bile solubility and full Optochin susceptibility for the differentiation of Streptococcus pneumoniae from other streptococci (1).

Hence optochin test is a useful diagnostic aid for identification / differentiation of pneumococci and viridans Streptococci.

Optochin discs are filter paper discs impregnated with optochin. The test is based on the property of viridans streptococci to grow in the presence of Optochin (ethyl hydrocuprein hydrochloride) which inhibits pneumococci. This test is performed for the diagnosis of penumococcal infections. Specimens of sputum, lung aspirate, pleural fluid, CSF, urine or blood are first examined by Gram's stain, cultured and the isolates are then subjected to Optochin Sensitivity Test.

Quality Control

Appearance

Filter paper discs of 6 mm diameter bearing letters "Op" in continuous printing style.

Cultural response

Cultural response observed after an incubation at 35-37°C for 18-24 hours at on seeded Soyabean Casein Digest Agar (M290) with added sterile defibrinated blood, using Optochin discs.

| Organism | Zone of |
|--------------------------|---------------|
| | inhibition |
| Streptococcus pneumoniae | More than or |
| ATCC 6303 | equal to 15mm |

Storage and Shelf Life

Store at 2 - 8°C. Use before expiry date on the label.

Reference

1.Bowers E.F. and Jeffries L.R., 1995, J. Clin. Path., 8:58.

Revision : 1 / 2011

CE

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia[™] publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia[™] Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. A-516, Swastik Disha Business Park, Via Vadhani Ind. Est., LBS Marg, Mumbai-400086, India. Customer care No.: 022-6147 1919 Email: techhelp@himedialabs.com



Oxidase Discs

DD018

Oxidase Discs are used for detection of oxidase production by microorganisms like Neisseria, Alcaligenes, Aeromonas, Vibrio's, Campylobacter and Pseudomonas, which give positive reactions and for excluding Enterobacteriaceae, which give negative reactions.

Directions

Oxidase reaction is carried out by touching and spreading a well isolated colony on the oxidase disc. The reaction is observed within 5-10 seconds at 25-30°C. A change later than 10 seconds or no change at all is considered negative reaction.

Precautions

1. "Do not use stainless steel or nichrome inoculating wires, as false positive reaction may result from surface oxidation products formed during flame sterilization.

- 2. "Growth from media containing dyes is not suitable for testing.
- 3. "Timing is critical (5-10 sec) for interpretation of results.
- 4. "Perform oxidase test on all gram-negative bacilli.

5. "Cytochrome oxidase production may be inhibited byacid production. False negative reactions may be exhibited by Vibrio, Aeromonas and Plesiomonas species when grown on a medium containing fermentable carbohydrate e.g. MacConkey Agar (M081). Colonies taken from media containing nitrate may give unreliable results. The loss of activity of the oxidase reagent is caused by auto-oxidation which may be avoided by adding 0.1% ascorbic acid (3).

Principle And Interpretation

Certain bacteria posses either cytochrome oxidase or indophenol oxidase (an iron-containing haemoprotein), which catalyzes the transport of electrons from donor compounds (NADH) to electron acceptors (usually oxygen). In the oxidase test, a colourless dye such as N, N-dimethy-p-phenylenediamine serves as an artificial electron acceptor for the enzyme oxidase. The dye is oxidized to form indophenol blue, a coloured compound. The test is useful in the initial characterization of aerobic gramnegative bacteria of the genera Aeromonas, Plesiomonas, Pseudomonas, Campylobacter and Pasteurella.

Oxidase discs are sterile filter paper discs impregnated with N, N-dimethyl-p-phenylenediamine oxalate, ascorbic acid and a-naphthol. These discs overcome the neccessity of daily preparation of fresh reagent. Gordon and McLeod (1) introduced oxidase test for identifying gonococci based upon the ability of certain bacteria to produce indophenol blue from the oxidation of dimethyl-p-phenylenediamine and a-naphthol. Gaby and Hadley (2) introduced a more sensitive method by using N, N-dimethyl-p-phenylenediamine oxalate where all staphylococci were oxidase negative. In a positive reaction the enzyme cytochrome oxidase combines with N,N-dimethyl-p-phenylenediamine oxalate and a-naphthol to form the dye, indophenol blue.

Quality Control

Appearance

Filter paper discs of 10 mm diameter

Cultural response

Typical oxidase reaction given by 18-48 hour culture observed within 5-10 seconds at 25-30°C.

| Organism | Reaction |
|------------------------|-----------------|
| | Observed |
| Pseudomonas aeruginosa | positive : deep |
| ATCC 27853 | purplish blue |
| | colouration of |
| | disc |

| Neisseria gonorrhoeae ATCC 19424 | positive : deep purplish blue colouration of disc |
|-------------------------------------|--|
| Escherichia coli ATCC 25922 | negative : purplish blue colouration after 10 sec/ |
| Staphylococcus aureus ATCC 25923 | no colour change negative : no colour change |

Storage and Shelf Life

Store at 2 - 8°C. Use before expiry date on the label.

Reference

1.Gordon J. and Mcleod J.W., 1928, J. Path. Bact., 31:185 2.Gaby W.L and Hadley C., 1957. J. Bact., 74:356 3.Steel. K.J. 1962. J. Appl. Bact. 25:445

Revision : 1 / 2011

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia[™] publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia[™] Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. A-516, Swastik Disha Business Park, Via Vadhani Ind. Est., LBS Marg, Mumbai-400086, India. Customer care No.: 022-6147 1919 Email: techhelp@himedialabs.com

€



X Factor discs

DD020

Used for the presumptive identification of *Haemophilus* species on the basis of their requirements for X or V factors or both.

Directions

Inoculate the surface of a Blood Agar (M073) plate or Brain Heart Infusion Agar (M211) plate with the test organisms by either streaking or surface spreading. Aseptically place the X (DD020), V (DD021) and X+V (DD022) factor discs on the plate, in the following positions:

Disc Position on the Agar plate

X factor disc 12 O' clock

V factor disc 4 O' clock

X+V factor disc 8 O' clock

Incubate the plates at 35 - 37°C for 24 - 48 hours. Observe for the growth in the neighbourhood of the discs.

Principle And Interpretation

Both X and V factors are growth factors that are essential for certain organisms like *Haemophilus* species and also enhance growth of organisms like *Neisseria* species.

X-factor discs are the sterile filter paper discs impregnated with growth factor X which are used for differentiating *Haemophilus* species in conjuction of V factor & X+V factor discs. *Bordetella* and *Haemophilus* species can also be identified on the basis of the requirement of X and V growth factors in the basal medium.

The X factor (hemin) and V factor (Coenzyme-Nicotinamide adenine dinucleotide NAD+) are impregnated on the sterile filter paper discs of diameter 6 mm.

The test organism requiring X factor alone, grows only in the vicinities of X and X+V factor discs. Those which require V factor alone grow in the vicinities of V and X+V factor discs. If both X and V factors are required, then the organism will grow only in the vicinity of the X+V factor discs. This satellite growth is seen around the disc promoting growth (1).

Quality Control

Appearance

Filter paper discs of 6 mm diameter bearing letters "X" in continuous printing style.

Cultural response

Cultural characteristics observed on Brain Heart Infusion Agar (M211) or Blood Agar Base (M073) after an incubation of 24-48 hours at 35-37°C.

Cultural Response

| Organism | Growth with X factor | Growth without growth factor |
|---|---|---|
| Cultural response | | 0 |
| Bordetella pertussis ATCC 8467 | Positive(initial isolation on Bordet Gengou Agar (M175)) | Positive(initial isolation on Bordet Gengou Agar (M175)) |
| Haemophilus influenzae ATCC 35056 | Negative | Negative |
| Haemophilus parainfluenzae ATCC 7901 | e Negative | Negative |

| Haemophilus | Positive | Negative |
|---------------------|----------|----------|
| haemoglobinophilus | | |
| ATCC19416 | | |
| Haemophilus ducreyi | Positive | Negative |

Storage and Shelf Life

Store between 2-8°C. For prolonged use store at -20°C. Use before the expiry date on the label.

Reference

1.Murray PR, Baron EJ, Jorgensen J.H., Pfaller M A, Yolken R.H(Eds.),8th ed, 2003, Manual of Clinical Microbiology, ASM, Washington D.C.

Note:

Use known strains of *Haemophilus influenzae* to monitor the performance of the differentiation discs and the medium. Do not use too heavy suspension of the test organisms as X or V factor carryover from the primary growth medium may take place

Revision : 1 / 2011

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia[™] publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia[™] Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. A-516,Swastik Disha Business Park,Via Vadhani Ind. Est., LBS Marg, Mumbai-400086, India. Customer care No.: 022-6147 1919 Email: techhelp@himedialabs.com

CE



V Factor Discs

DD021

Used for the presumptive identification of *Haemophilus* species on the basis of their requirements for X or V factors or both.

Directions

Inoculate the surface of a Blood Agar (M073) plate or Brain Heart Infusion Agar (M211) plate with the test organisms by either streaking or surface spreading. Aseptically place the X (DD020), V (DD021) and X+V (DD022) factor discs on the plate, in the following positions:

Disc Position on the Agar plate

X factor disc 12 O' clock

V factor disc 4 O' clock

X+V factor disc 8 O' clock

Incubate the plates at 35 - 37°C for 24 - 48 hours. Observe for the growth in the neighbourhood of the discs.

Principle And Interpretation

Both X and V factors are growth factors that are essential for certain organisms like *Haemophilus* species and also enhance growth of organisms like *Neisseria* species.

V-factor discs are the sterile filter paper discs impregnated with growth factor V which are used for differentiating *Haemophilus* species in conjuction of X factor & X+V factor discs. *Bordetella* and *Haemophilus* species can also be identified on the basis of the requirement of X and V growth factors in the basal medium.

The X factor (hemin) and V factor (Coenzyme-Nicotinamide adenine dinucleotide NAD+) are impregnated on the sterile filter paper discs of diameter 6 mm.

The test organism requiring X factor alone, grows only in the vicinities of X and X+V factor discs. Those which require V factor alone grow in the vicinities of V and X+V factor discs. If both X and V factors are required, then the organism will grow only in the vicinity of the X+V factor discs. This satellite growth is seen around the disc promoting growth (1).

Quality Control

Appearance

Filter paper discs of 6 mm diameter bearing letters "V" in continuous printing style.

Cultural response

Cultural characteristics observed on Brain Heart Infusion Agar (M211) or Blood Agar Base (M073) after an incubation of 24-48 hours at 35-37°C.

| Organism | Growth with V factor | Growth without growth factor |
|--|---|---|
| Bordetella pertussis ATCC 8467 | Positive(initial isolation on Bordet Gengou Agar (M175)) | Positive(initial isolation on Bordet Gengou Agar (M175)) |
| Haemophilus influenzae ATCC 35056 | Negative | Negative |
| Haemophilus parainfluenza ATCC 7901 | e Positive | Negative |

| Haemophilus | Negative | Negative |
|---------------------|----------|----------|
| haemoglobinophilus | | |
| ATCC19416 | | |
| Haemophilus ducreyi | Negative | Negative |

Storage and Shelf Life

Store below -10°C. Use before the expiry date on the label.

Reference

1.Murray PR, Baron EJ, Jorgensen J.H., Pfaller M A, Yolken R.H(Eds.),8th ed, 2003, Manual of Clinical Microbiology, ASM, Washington D.C.

Note:

Use known strains of *Haemophilus influenzae* to monitor the performance of the differentiation discs and the medium. Do not use too heavy suspension of the test organisms as X or V factor carryover from the primary growth medium may take place

Revision : 1 / 2011

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia[™] publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia[™] Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. A-516,Swastik Disha Business Park,Via Vadhani Ind. Est., LBS Marg, Mumbai-400086, India. Customer care No.: 022-6147 1919 Email: techhelp@himedialabs.com

CE



X+V Factor discs

DD022

Used for the presumptive identification of *Haemophilus* species on the basis of their requirements for X or V factors or both.

Directions

Inoculate the surface of a Blood Agar (M073) plate or Brain Heart Infusion Agar (M211) plate with the test organisms by either streaking or surface spreading. Aseptically place the X (DD020), V (DD021) and X+V (DD022) factor discs on the plate, in the following positions:

Disc Position on the Agar plate

X factor disc 12 O' clock

V factor disc 4 O' clock

X+V factor disc 8 O' clock

Incubate the plates at 35 - 37°C for 24 - 48 hours. Observe for the growth in the neighbourhood of the discs.

Principle And Interpretation

Both X and V factors are growth factors that are essential for certain organisms like *Haemophilus* species and also enhance growth of organisms like *Neisseria* species.

X+V factor discs are the sterile filter paper discs impregnated with growth factors x <(>&<)> V which are used for differentiating *Haemophilus* species in conjuction of X factor & V factor discs. *Bordetella* and *Haemophilus* species can also be identified on the basis of the requirement of X and V growth factors in the basal medium.

The X factor (hemin) and V factor (Coenzyme-Nicotinamide adenine dinucleotide NAD+) are impregnated on the sterile filter paper discs of diameter 6 mm.

The test organism requiring X factor alone, grows only in the vicinities of X and X+V factor discs. Those which require V factor alone grow in the vicinities of V and X+V factor discs. If both X and V factors are required, then the organism will grow only in the vicinity of the X+V factor discs. This satellite growth is seen around the disc promoting growth (1).

Quality Control

Appearance

Filter paper discs of 6 mm diameter bearing letters "X+V" in continuous printing style.

Cultural response

Cultural characteristics observed on Brain Heart Infusion Agar (M211) or Blood Agar Base (M073) after an incubation of 24-48 hours at 35-37°C.

| Organism | Growth with X +V factor | Growth without growth factor |
|---|---|---|
| Bordetella pertussis ATCC 8467 | Positive(initial isolation on Bordet Gengou Agar (M175)) | Positive(initial isolation on Bordet Gengou Agar (M175)) |
| Haemophilus influenzae ATCC 35056 | Positive | Negative |
| Haemophilus parainfluenzae ATCC 7901 | Positive | Negative |

| Haemophilus | Positive | Negative |
|---------------------|----------|----------|
| haemoglobinophilus | | |
| ATCC19416 | | |
| Haemophilus ducreyi | Positive | Negative |

Storage and Shelf Life

Store below -10°C. Use before the expiry date on the label.

Reference

1.Murray PR, Baron EJ, Jorgensen J.H., Pfaller M A, Yolken R.H(Eds.),8th ed, 2003, Manual of Clinical Microbiology, ASM, Washington D.C.

Note:

Use known strains of *Haemophilus influenzae* to monitor the performance of the differentiation discs and the medium. Do not use too heavy suspension of the test organisms as X or V factor carryover from the primary growth medium may take place

Revision : 1 / 2011

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia[™] publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia[™] Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. A-516,Swastik Disha Business Park,Via Vadhani Ind. Est., LBS Marg, Mumbai-400086, India. Customer care No.: 022-6147 1919 Email: techhelp@himedialabs.com

CE



Bile Esculin Discs

DD024

Bile Esculin Discs are used for detection of esculin hydrolysis in the presence of bile, for differentiating Group D streptococci from other Streptococcal groups.

Directions

Esculin impregnated disc is placed on the seeded Bile Esculin Agar Base (M340) plate and is incubated at 35-37°C for 18-24 hours.

Principle And Interpretation

Group D streptococci hydrolyze esculin to esculetin and dextrose. Esculetin reacts with an iron salt such as ferric citrate to form a blackish brown coloured complex (4).

Rochaix found that esculin hydrolysis is an important criteria in the identification of enterococci (1). Meyer and Schonfeld (2) observed that when bile was added to esculin medium, around 60% enterococci were able to grow and split the esculin while other streptococci could not. When a comparative study was performed by Facklam and Moody (3) for presumptive identification of Group D streptococci, they found the bile esculin test as a reliable means of identifying Group D streptococci and differentiating them from other streptococci groups.

Quality Control

Appearance

Plain filter paper discs of 6mm diameter

Cultural response

Cultural response observed by placing Bile Esculin disc (DD024) on seeded Bile Esculin Agar Base(M340) plate, incubated at 35-37°C for 18-24 hours.

| Organism | Growth | Esculin hydrolysis |
|--|--------------|---|
| Enterococcus faecalis ATC 29212 | CC luxuriant | positive: blackening of media around the disc. |
| Streptococcus agalactiae ATCC 13813 | luxuriant | negative: no blackening |
| Listeria monocytogenes ATCC 19118 | luxuriant | positive: blackening of media around the disc. |
| Streptococcus pyogenes ATCC 19615 | luxuriant | negative: no blackening |

Storage and Shelf Life

Store at 2 - 8°C. Use before expiry date on the label.

Reference

1. Rochaix, 1924, C. R. Soc. Biol., 90:771.

2. Meyer and Schonfeld, 1926, Zentralbl. Bacteriol. Parasitenkd. Infectionskr. Hyg. Abt. I Orig., 99:402.

3. Facklam and Moody, 1970, Appl. Microbiol., 20:245.

4. MacFaddin J. F., 2000, Biochemical Tests for Identification of Medical Bacteria, 3rd ed., Philadelphia: Lippincott. Williams and Wilkins.

Revision : 1 / 2011

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia[™] publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia[™] Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. A-516, Swastik Disha Business Park, Via Vadhani Ind. Est., LBS Marg, Mumbai-400086, India. Customer care No.: 022-6147 1919 Email: techhelp@himedialabs.com

CE



DMACA Indole Discs

DD040

The DMACA Indole Discs are used for Indole test to determine the ability of an organism to split indole from the tryptophan molecule, and thus to aid differentiation between *Escherichia coli* from *Klebsiella*.

Directions

Place the DMACA Indole Disc on suspected colony from HiCrome UTI Agar (M1353) or HiCrome UTI Agar, Modified (M1418) plate. Observe for appearance of blue-purple colour within 10 - 30 seconds.

Principle And Interpretation

In the presence of oxygen, some bacteria are able to split tryptophan into indole and alpha-aminopropionic acid. The presence of indole can be detected by the addition of DMACA (p-Dimethylaminocinnamaldehyde) reagent indicated by formation of bluish-purple colour (1).

Quality Control

Appearance

Filter paper discs of 6 mm diameter bearing letters 'Dm' in continuous printing style.

Cultural response

The indole production by organisms was tested after an incubation of 18-24 hours at 35-37°C, using HiCrome UTI Agar (M1353).

Cultural Response

| Organism | Indole production |
|--------------------------------------|---|
| Cultural response | |
| Escherichia coli ATCC | Positive |
| 25922 | reaction, blue- purple colour formation |
| Klebsiella pneumoniae | Negative |
| ATCC 13883 | reaction. |
| Pseudomonas aeruginosa ATCC 27853 | Negative reaction. |

Storage and Shelf Life

Store at 2-8°C. Use before the expiry date on the label.

Reference

1.MacFaddin J. F., 1980, Biochemical Tests for Identification of Medical Bacteria, 2nd ed., Williams and Wilkins, Baltimore.

Revision : 1 / 2011

(6

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia[™] publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia[™] Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. A-516, Swastik Disha Business Park, Via Vadhani Ind. Est., LBS Marg, Mumbai-400086, India. Customer care No.: 022-6147 1919 Email: techhelp@himedialabs.com



Colistin Ezy MICTM Strip (CL) (0.016-256 mcg/ml)

EM020

Antimicrobial Susceptibility Testing For *In Vitro* Diagnostic use

Not for Medicinal Use

It is a unique MIC determination paper strip which is coated with Colistin on a single paper strip in a concentration gradient manner, capable of showing MICs in the range of 0.016mcg/ml to 256 mcg/ml, on testing against the test organism.

Introduction

Ezy MICTM strip is useful for quantitative determination of susceptibility of bacteria to antibacterial agents. The system comprises of a predefined quantitative gradient which is used to determine the Minimum Inhibitory Concentration (MIC) in mcg/ml of different antimicrobial agents against microorganisms as tested on appropriate agar media, following overnight incubation.

Ezy MICTM Strip Features and Advantages

Ezy MICTM strip exhibits several advantages over existing plastic strip.

- 1) Ezy MICTM strip is made up of porous paper material unlike plastic non-porous material.
- 2) Ezy MICTM strip has MIC values printed on both sides identically.
- 3) The antimicrobial agent is evenly distributed on either side of the
- Ezy MICTM strip and hence it can be placed by any side on the agar surface.
- 4) For Ezy MICTM strips, MIC values can be read without opening the lid of the plate as most commonly translucent medium such as Mueller Hinton Agar is employed.
- 5) Once placed, Ezy MICTM strip is adsorbed within 60 seconds and firmly adheres to the agar surface.
- 6) Unlike the plastic material, it does not form air bubbles underneath and hence there is no need to press the strip once placed.

METHOD AND USE OF EZY MICTM STRIPS

• <u>Type of specimen</u>

Pure cultures should be derived from specimens obtained from patients prior to the initiation of antimicrobial therapy. Specimens can be of bacterial or fungal isolates derived from blood, urine, faeces, pus, CSF etc. Direct specimens should not be employed in this test. Refer procedure, which includes preparation of inoculum (1, 3).

<u>Clinical specimen collection, handling and processing</u>

Follow appropriate techniques for handling specimens as per established guidelines. After use, contaminated materials must be sterilized by autoclaving before discarding (1, 3).

<u>Guidelines for preparation of the medium</u>

Prepare the medium of choice from dehydrated powder according to the directions specified on the label. Cool the sterilized molten medium to 45-50°C and pour in sterile, dry Petri plates on a leveled surface, to a depth of 4 ± 0.2 mm and allow to solidify. Few droplets appearing on the surface of the medium following cooling do not matter. Hence, once poured, Petri plates containing media should not be dried on laminar flow and can be used immediately for swabbing.

<u>Preparation of Inoculum</u>

Use only pure cultures. Confirm by Gram-staining before starting susceptibility test. Transfer 4-5 similar colonies with a wire, needle or loop to 5 ml Tryptone Soya Broth (M011) and incubate at 35-37°C for2-8 hours until light to moderate turbidity develops. Compare the inoculum turbidity with that of standard 0.5 McFarland. Alternatively, the inoculum can be standardized by other appropriate optical method (0.08 - 0.13 OD turbid suspension at 620 nm). Also direct colony suspension method can be used. Prepare a direct colony suspension, from 18-24 hour old non-selective media agar plate in broth or saline. Adjust the turbidity to that of standard 0.5 McFarland .This method is recommended for testing fastidious organisms like *Haemophilus* spp., *Neisseria* spp, streptococci and for testing staphylococci for potential Methicillin or Oxacillin resistance.

<u>Test Procedure</u>

- 1. Prepare plates with suitable make of Mueller Hinton Agar. For fastidious organisms such as Streptococci, Mueller Hinton Agar is supplemented with 5% sterile, defibrinated blood is recommended.
- 2. Dip a sterile non-toxic cotton swab on a wooden applicator into the standardized inoculum and rotate the soaked swab firmly against the upper inside wall of the tube to express excess fluid. Streak the entire agar surface of the plate with the swab three times, turning the plate at 60° angle between each streaking.
- 3. Remove Ezy MIC[™] strip container from cold and keep it at room temperature for 15-30 minutes before opening.
- 4. Remove one applicator from the self-sealing bag stored at room temperature.
- 5. Hold the applicator in the middle and gently press its broader sticky side on the centre of Ezy MICTM strip.
- 6. Lift the applicator along with attached Ezy MICTM strip.
- 7. Place the strip at a desired position on agar plate pre-spread with test culture. Gently turn the applicator clockwise with fingers. With this action, the applicator will detach from the strip.
- 8. DO NOT PRESS EZY MICTM STRIP. Within 60 seconds, Ezy MICTM strip will be adsorbed and will firmly adhere to the agar surface.
- 9. Ezy MICTM strip should not be repositioned or adjusted once placed.
- 10. Transfer plates in the incubator under appropriate conditions.

Warning and Precautions:

- 1. Ezy MICTM Strip is intended for *In vitro* diagnostic use only.
- 2. Although based on simple procedure, Ezy MICTM Strip should only be used by at least semi-trained personnel.
- 3. This strip is intended only for agar diffusion method and not for broth dilution method.
- 4. Ezy MIC[™] Strip should be used strictly according to procedures described herein.
- 5. Performance of Ezy MICTM Strips depends on use of proper inoculum and control cultures, recommended test medium and proper storage temperature.
- 6. Follow aseptic techniques and precautions against microbiological hazards should be used when handling bacterial or fungal specimen throughout the testing procedure.
- 7. Before using Ezy MICTM Strips, ensure that the strip is at room temperature.
- 8. When applying strips be steady. Do not move the strip once in contact with agar surface, since the antibiotic instantaneously diffuse on contact with agar.
- 9. Place the unused strips back to recommended temperature.

INTERPRETATION & QUALITY CONTROL (As per CLSI):

Table 1: Use following interpretive criteria for susceptibility categorization as per CLSI.

| When testing | Incubation | Inte | Interpretative Criteria | |
|---------------------------------|------------------------|---------------|-------------------------|-----|
| | | <u><</u> S | Ι | ≥ R |
| Other non-Enterobacteriaceae | 35-37°C for 18 hrs. | 2 | 4 | 8 |
| Acinetobacter spp, P.aeruginosa | 35-37°C for 18 hrs. | 2 | - | 4 |

Table 2: Quality Control for EUCAST Recommended MCR4 E. coli

| Method | MIC (mcg/ml) |
|---------------------------------------|-------------------|
| EUCAST Recommended MIC | 4 |
| In-house broth Dilution method | 4 (mode value) |
| MIC with Ezy MIC TM strips | 4 (mode value) |
| | MIC Range : 2-4-8 |

| Table 5. In-House Quality Control for Constant Resistant Chinear Isolates | | | | | |
|---|---|---|--------------------------------------|--|--|
| Organisms | Mode MICs (mcg/ml) obtained by micro-broth dilution method N=10 | Mode MICs (mcg/ml) obtained by Ezy MIC [™] strips N=20 | Recommended MIC (mcg/ml) range | | |
| Isolate- 1 | 32 | 32 | 16-32-64 | | |
| Isolate- 2 | 16 | 16 | 8-16-32 | | |
| Isolate- 3 | 8 | 8 | 4-8-16 | | |

Table 3: In-house Quality Control for Colistin Resistant Clinical Isolates

Table 4: Following are the reference MIC values (mcg/ml) range for Colistin as per CLSI and EUCAST.

| Organism | Medium used | Incubation | Std.QualityControllimits (mcg/ml) |
|--------------------------|---------------------|---------------------|-----------------------------------|
| E.coli ATCC 25922 | Mueller Hinton Agar | 35-37°C for 18 hrs. | 0.25 - 0.5 - 1.0 - 2.0 |
| P. aeruginosa ATCC 27853 | Mueller Hinton Agar | 35-37°C for 18 hrs. | 0.5 -1.0 - 2.0 - 4.0 |
| E. coli NCTC 13846 | Mueller Hinton Agar | 35-37°C for 18 hrs. | 2.0 - 4.0 - 8.0 |

Storage & Shelf Life:

- 1. Once the consignment is received, store applicators at Room Temperature and Ezy MICTM strips container at -20°C or below.
- 2. Use before expiry date on the label.
- 3. Ezy MICTM Strip left over from opened package must be kept dry.
- 4. Moisture should be prevented from penetrating into or forming within the package or storage container.
- 5. Check whether the batch number and expiry date are marked on the storage container.
- 6. Product performance is best within stated expiry period if correctly stored and handled.

Disposal

After use, Ezy MICTM Strips and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques (2, 3).

Limitation of Test

Ezy MICTM Strips provides *In vitro* MIC values, which provides only a possible insinuation of pathogens potential in *In vivo* susceptibility. These values can be considered as a guide to therapy selection only after taking into consideration several other factors; and must be the sole decision and responsibility of the physician along with the clinical experience in treating the infection. These tests are comparable to the standards as per the given specifications and set of experiment standards as far as possible. Please refer to CLSI standards for detailed limitation of susceptibility test on the clinical use of an antibiotic in various therapeutic conditions.

References:

- 1. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition, Vol. 1, Section 2.
- 2. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition, Vol. 3, Section 15.
- 3. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry M.L., Richter, S.S and Warnock., D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.
- 4. Performance standards of Antimicrobial Susceptibility Testing; Twenty Ninth Informational Supplement. M100-S29, Vol. 39, No.1, Jan 2019.
- 5. Performance standards of Antimicrobial Susceptibility Testing; Twenty Sixth Informational Supplement. M100-S26, Vol. 36, No.1, Jan 2016.
- 6. European Committee on Antimicrobial Susceptibility Testing: Breakpoint tables for interpretation of MICs and zone diameters. Version 9.0, valid from 2019-01-01.
- 7. European Committee on Antimicrobial Susceptibility Testing:Routine and extended internal quality control for MIC determination and disk diffusion as recommended by EUCAST., Version 9.0, valid from 2019-01-01.

Packing:

Each Pack contains following material packed in air-tight plastic container with a desiccator capsule.

- 1) Colistin Ezy MICTM strips (10/30/60/90/120/150 Strips per pack)
- 2) Applicator sticks
- 3) Package insert

Disclaimer :

CE

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia[™] publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia[™] Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. A-516,Swastik Disha Business Park,Via Vadhani Ind. Est., LBS Marg, Mumbai-400086, India. Customer care No.: 022-6147 1919 Email: techhelp@himedialabs.com Website: www.himedialabs.com



Vancomycin Ezy MICTM Strip (VAN) (0.016-256 mcg/ml)

EM060

Antimicrobial Susceptibility Testing For *In Vitro* Diagnostic use

Not for Medicinal Use

It is a unique MIC determination paper strip which is coated with Vancomycin in a concentration gradient manner, capable of showing MICs in the range of 0.016mcg/ml to 256 mcg/ml, on testing against the test organism.

Introduction

Ezy MICTM strip is useful for quantitative determination of susceptibility of bacteria to antibacterial agents. The system comprises of a predefined quantitative gradient which is used to determine the Minimum Inhibitory Concentration (MIC) in mcg/ml of different antimicrobial agents against microorganisms as tested on appropriate agar media, following overnight incubation.

Ezy MICTM Strip FEATURES AND ADVANTAGES

Ezy MICTM strip exhibits several advantages over existing plastic strip.

- 1. Ezy MICTM strip is made up of porous paper material unlike plastic non-porous material
- 2. Ezy MICTM strip has MIC values printed on both sides identically.
- 3. The antimicrobial agent is evenly distributed on either side of the Ezy MIC[™] strip and hence it can be placed by any side on the agar surface.
- 4. For Ezy MICTM strips, MIC values can be read without opening the lid of the plate as most commonly translucent medium such as Mueller Hinton Agar is employed.
- 5. Once placed, Ezy MICTM strip is adsorbed within 60 seconds and firmly adheres to the agar surface.
- 6. Unlike the plastic material, it does not form air bubbles underneath and hence there is no need to press the strip once placed.

CLSI RECOMMENDATION FOR VANCOMYCIN SENSITIVITY TEST

High molecular weight antibiotics such as Vancomycin, polymyxin B and colistin do not diffuse in concentration gradient manner while diffusing through the agar medium when the disc susceptibility test is employed. The Antimicrobial Susceptibility Testing using disc diffusion test does not differentiate vancomycin-susceptible isolates of *S.aureus* from Vancomycin intermediate isolates, nor does the test differentiates among Vancomycin–susceptible, intermediate, and resistant isolates of coagulase-negative staphylococci, all of which may give similar size zones of inhibition.

CLSI therefore recommends that MIC test should be performed to determine the susceptibility of all isolates of staphylococci to Vancomycin .¹

Usefulness of Vancomycin Ezy MICTM strip

 Besides obtaining accurate MIC values for Gram- positive cultures, VISA (Vancomycin Intermediate Staphylococcus aureus) can be detected when isolated colonies appear within the zone of inhibition of Vancomycin particularly when 1.0 McFarland inoculum is used and MIC is read on full 48 hrs incubation. The sensitivity of the method can be further enhanced for better detection of VISA/ VRSA (Vancomycin Resistant Staphylococcus aureus / hVISA (Hetro Vancomycin Intermediate Staphylococcus aureus) using BHI agar with higher inoculum and 48 hr incubation.

METHOD AND USE OF EZY MICTM STRIPS

• <u>Type of specimen</u>

Pure cultures should be derived from specimens obtained from patients prior to the initiation of antimicrobial therapy. Specimens can be of bacterial or fungal isolates derived from blood, urine, faeces, pus, CSF etc. Direct specimens should not be employed in this test. Refer procedure, which includes preparation of inoculum (1,3).

<u>Clinical specimen collection, handling and processing</u>

Follow appropriate techniques for handling specimens as per established guidelines. After use, contaminated materials must be sterilized by autoclaving before discarding (1, 3).

• Guidelines for preparation of the medium

Prepare the medium of choice from dehydrated powder according to the directions specified on the label. Cool the sterilized molten medium to 45-50°C and pour in sterile, dry Petri plates on a leveled surface, to a depth of 4 ± 0.2 mm and allow to solidify. Few droplets appearing on the surface of the medium following cooling do not matter. Hence, once poured, Petri plates containing media should not be dried on laminar flow and can be used immediately for swabbing.

• <u>Preparation of Inoculum</u>

Use only pure cultures. Confirm by Gram-staining before starting susceptibility test. Transfer 4-5 similar colonies with a wire, needle or loop to 5 ml Tryptone Soya Broth (M011) and incubate at 35-37°C for2-8 hours until light to moderate turbidity develops. Compare the inoculum turbidity with that of standard 0.5 McFarland. Alternatively, the inoculum can be standardized by other appropriate optical method (0.08 - 0.13 OD turbid suspension at 620 nm yields 10^5 - 10^6 cells/ml).

Also direct colony suspension method can be used. Prepare a direct colony suspension, from 18-24 hour old nonselective media agar plate in broth or saline. Adjust the turbidity to that of standard 0.5 McFarland .This method is recommended for testing fastidious organisms like *Haemophilus* spp., *Neisseria* spp, and streptococci and for testing staphylococci for potential Methicillin or Oxacillin resistance.

<u>Test Procedure</u>

- 1. Prepare plates with suitable make of Mueller Hinton Agar for rapidly growing aerobic organisms as mentioned above. For fastidious organisms such as Streptococci, Mueller Hinton Agar supplemented with 5% sterile, defibrinated blood is recommended.
- 2. Dip a sterile non-toxic cotton swab on a wooden applicator into the standardized inoculum and rotate the soaked swab firmly against the upper inside wall of the tube to express excess fluid. Streak the entire agar surface of the plate with the swab three times, turning the plate at 60° angle between each streaking.
- 3. Remove Ezy MICTM strip container from cold and keep it at room temperature for 15 minutes before opening.
- 4. Remove one applicator from the self sealing bag stored at room temperature.
- 5. Hold the applicator in the middle and gently press its broader sticky side on the centre of Ezy MICTM strip.
- 6. Lift the applicator along with attached Ezy MICTM strip.
- 7. Place the strip at a desired position on agar plate swabbed with test culture. Gently turn the applicator clockwise with fingers. With this action, the applicator will detach from the strip.
- 8. DO NOT PRESS EZY MICTM STRIP. Within 60 seconds, Ezy MICTM strip will be adsorbed and will firmly adhere to the agar surface.
- 9. Ezy MICTM strip should not be repositioned or adjusted once placed.
- 10. Transfer plates in the incubator under appropriate conditions.

MIC Reading:

- 1. Read the plates only when sufficient growth is seen.
- 2. Read the MIC where the ellipse intersects the MIC scale on the strip.
- For bactericidal drugs such Vancomycin, Gentamicin, Amikacin, and members of β-lactams class of drugs, always read the MIC at the point of completion inhibition of all growth, including hazes, microcolonies and isolated colonies. If necessary, use magnifying glass.
- 4. Isolated colonies, microcolonies and hazes appearing in the zone of inhibition are indicative of hetero nature of the culture having resistant subpopulation in it. In such cases, consider reading for MIC determination at a point on the scale above which no resistant colonies are observed close to MIC strip (within 1-3 mm distance from the strip).
- 5. Since Ezy MIC[™] strip has continuous gradient, MIC values "in-between" two fold dilutions can be obtained.
- 6. Always round up these values to the next two-fold dilution before categorization. For example: Vancomycin showing reading of 0.75 mcg/ml should be rounded up to next concentration i.e. 1.0 mcg/ml.
- 7. If the ellipse intersects the strip in between 2 dilutions, read the MIC as the value which is nearest to the intersection.
- 8. When growth occurs along the entire strip, report the MIC as \geq the highest values on the MIC strip. When the inhibition ellipse is below the strip (does not intersect the strip), report the MIC < the lowest value on the MIC scale.

HiMedia Laboratories

Warning and Precautions:

- 1. Ezy MICTM Strip is intended for *In vitro* diagnostic use only.
- 2. Although based on simple procedure, Ezy MIC[™] Strip should only be used by at least semi-trained personnel.
- 3. This strip is intended only for agar diffusion method and not for broth dilution method.
- 4. Ezy MICTM Strip should be used strictly according to procedures described herein.
- 5. Performance of Ezy MICTM Strips depends on use of proper inoculum and control cultures, recommended test medium and proper storage temperature.
- 6. Follow aseptic techniques and precautions against microbiological hazards should be used when handling bacterial or fungal specimen throughout the testing procedure.
- 7. Before using Ezy MICTM Strips, ensure that the strips is at room temperature.
- 8. When applying strips be steady. Do not move the strip once in contact with agar surface, since the antibiotic instantaneously diffuse on contact with agar.
- 9. Place the unused strips back to recommended temperature.

INTERPRETATION & QUALITY CONTROL (As per CLSI Guidelines):

Interpretation

Table 1: Use following interpretive criteria for susceptibility categorization as per CLSI.

| When testing | Interpretative Criteria | | riteria |
|--|-------------------------|------|---------------|
| | <u><</u> S | Ι | <u>> R</u> |
| Staphylococcus spp | 2 | 4-8 | 16 |
| Coagulase negative Staphylococci spps. and Enterococcus | 4 | 8-16 | 32 |
| S.pneumoniae, Streptococcus spps. Beta haemolytic group, Streptococcus | 1 | - | - |
| spps. Viridans group | | | |

QUALITY CONTROL

Quality control of Ezy MICTM Strip is carried out by testing the strips with standard ATCC Cultures recommended by CLSI on suitable medium incubated appropriately.

Following are the reference MIC values (mcg/ml) range for Vancomycin.

| Organism | Medium used | Incubation | Std. Quality Control limits (mcg/ml) |
|--------------------------|------------------------|----------------------------|---|
| S.aureus ATCC 29213 | Mueller Hinton Agar | 35-37°C for 18 hrs. | 0.5 - 1.0 - 2.0 |
| E.faecalis ATCC 29212 | Mueller Hinton Agar | 35-37°C for 18 hrs. | 1.0 - 2.0 - 4.0 |
| S. pneumoniae ATCC 49619 | Mueller Hinton Agar w/ | 35-37°C for 20-24hrs at 5% | 0.12-0.25-0.5 |
| | 5% Sheep Blood | CO ₂ | |

Storage & Shelf Life:

- 1. Once the consignment is received, store applicators at Room Temperature and Ezy MICTM strips container at 2-8°C, for prolonged use store below -20°C.
- 2. Use before expiry date on the label.
- 3. Ezy MICTM Strip left over from opened package must be kept dry.
- 4. Moisture should be prevented from penetrating into or forming within the package or storage container.
- 5. Check whether the batch number and expiry date are marked on the storage container.
- 6. Product performance is best within stated expiry period if correctly stored and handled.

Disposal:

After use, Ezy MICTM Strips and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques (2, 3).

Limitation of Test

Ezy MICTM Strips provides *In vitro* MIC values, which provides only a possible insinuation of pathogens potential in *In vivo* susceptibility. These values can be considered as a guide to therapy selection only after taking into consideration several other factors; and must be the sole decision and responsibility of the physician along with the clinical experience in treating the infection. These tests are comparable to the standards as per the given specifications and set of experiment standards as far as possible. Please refer to CLSI standards for detailed limitation of susceptibility test on the clinical use of an antibiotic in various therapeutic conditions.

References:

- 1. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition, Vol. 1, Section 2.
- 2. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition, Vol. 3, Section 15.
- 3. Jorgensen, J.H., Pfaller , M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock., D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.
- 4. Performance standards of Antimicrobial Susceptibility Testing; Twenty Ninth Informational Supplement. M100-S29, Vol. 39, No.1, Jan 2019.

Packing:

Each Pack contains following material packed in air-tight plastic container with a desiccator capsule.

- 1) Vancomycin Ezy MIC[™] strips (10/30/60/90/120/150 Strips per pack)
- 2) Applicator sticks
- 3) Package insert

Disclaimer :

€

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia[™] publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia[™] Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. A-516,Swastik Disha Business Park,Via Vadhani Ind. Est., LBS Marg, Mumbai-400086, India. Customer care No.: 022-6147 1919 Email: techhelp@himedialabs.com Website: www.himedialabs.com



Penicillin Ezy MICTM Strip (PEN) (0.002-32 mcg/ml)

EM084

Antimicrobial Susceptibility Testing For *In Vitro* Diagnostic use

Not for Medicinal Use

It is a unique MIC determination paper strip which is coated with Penicillin on a single paper strip in a concentration gradient manner, capable of showing MICs in the range of 0.002 mcg/ml to 32 mcg/ml, on testing against the test organism.

Introduction:

Ezy MICTM strip is useful for quantitative determination of susceptibility of bacteria to antibacterial agents. The system comprises of a predefined quantitative gradient which is used to determine the Minimum Inhibitory Concentration (MIC) in mcg/ml of different antimicrobial agents against microorganisms as tested on appropriate agar media, following overnight incubation.

Ezy MICTM Strip FEATURES AND ADVANTAGES

Ezy MICTM strip exhibits several advantages over existing plastic strip.

- 1. Ezy MICTM strip is made up of porous paper material unlike plastic non-porous material.
- 2. Ezy MICTM strip has MIC values printed on both sides identically.
- 3. The antimicrobial agent is evenly distributed on either side of the Ezy MICTM strip and hence it can be placed by any side on the agar surface.
- 4. For Ezy MIC[™] strips, MIC values can be read without opening the lid of the plate as most commonly translucent medium such as Mueller Hinton Agar is employed.
- 5. Once placed, Ezy MICTM strip is adsorbed within 60 seconds and firmly adheres to the agar surface.
- 6. Unlike the plastic material, it does not form air bubbles underneath and hence there is no need to press the strip once placed.

METHOD AND USE OF EZY MICTM STRIPS

• <u>Type of specimen</u>

Pure cultures should be derived from specimens obtained from patients prior to the initiation of antimicrobial therapy. Specimens can be of bacterial or fungal isolates derived from blood, urine, faeces, pus, CSF etc. Direct specimens should not be employed in this test. Refer procedure, which includes preparation of inoculum (1, 3).

• <u>Clinical specimen collection, handling and processing</u>

Follow appropriate techniques for handling specimens as per established guidelines. After use, contaminated materials must be sterilized by autoclaving before discarding (1, 3).

• <u>Guidelines for preparation of the medium</u>

Prepare the medium of choice from dehydrated powder according to the directions specified on the label. Cool the sterilized molten medium to 45-50°C and pour in sterile, dry Petri plates on a leveled surface, to a depth of 4 ± 0.2 mm and allow to solidify. Few droplets appearing on the surface of the medium following cooling do not matter. Hence, once poured, Petri plates containing media should not be dried on laminar flow and can be used immediately for swabbing.

<u>Preparation of Inoculum</u>

Use only pure cultures. Confirm by Gram-staining before starting susceptibility test. Transfer 4-5 similar colonies with a wire, needle or loop to 5 ml Tryptone Soya Broth (M011) and incubate at 35-37°C for2-8 hours until light to moderate turbidity develops. Compare the inoculum turbidity with that of standard 0.5 McFarland. Alternatively, the inoculum can be standardized by other appropriate optical method (0.08 - 0.13 OD turbid suspension at 620 nm). Also direct colony suspension method can be used. Prepare a direct colony suspension, from 18-24 hour old non-selective media agar plate in broth or saline. Adjust the turbidity to that of standard 0.5 McFarland .This method is recommended for testing fastidious organisms like *Haemophilus* spp., *Neisseria* spp, streptococci and for testing staphylococci for potential Methicillin or Oxacillin resistance.

<u>Test Procedure</u>

- 1. Prepare plates with suitable make of Mueller Hinton Agar for rapidly growing aerobic organisms as mentioned above. For fastidious organisms such as Streptococci, Mueller Hinton Agar is supplemented with 5% sterile, defibrinated blood and for *Neisseria gonorrhoeae*, GC Agar Base (M434) with 1% defined growth supplement (FD025) is recommended.
- 2. Dip a sterile non-toxic cotton swab on a wooden applicator into the standardized inoculum and rotate the soaked swab firmly against the upper inside wall of the tube to express excess fluid. Streak the entire agar surface of the plate with the swab three times, turning the plate at 60° angle between each streaking.
- 3. Remove Ezy MICTM strip container from cold and keep it at room temperature for 15 minutes before opening.
- 4. Remove one applicator from the self sealing bag stored at room temperature.
- 5. Hold the applicator in the middle and gently press its broader sticky side on the centre of Ezy MICTM strip.
- 6. Lift the applicator along with attached Ezy MIC[™] strip.
- 7. Place the strip at a desired position on agar plate swabbed with test culture. Gently turn the applicator clockwise with fingers. With this action, the applicator will detach from the strip.
- 8. DO NOT PRESS EZY MICTM STRIP. Within 60 seconds, Ezy MICTM strip will be adsorbed and will firmly adhere to the agar surface.
- 9. Ezy MICTM strip should not be repositioned or adjusted once placed.
- 10. Transfer plates in the incubator under appropriate conditions.

MIC Reading:

- 1. Read the plates only when sufficient growth is seen.
- 2. Read the MIC where the ellipse intersects the MIC scale on the strip.
- 3. For bactericidal drugs such as Penicillin and other members of β-lactams class of drugs, Amikacin, Vancomycin, Gentamicin, always read the MIC at the point of complete inhibition of all growth, including hazes, microcolonies and isolated colonies. If necessary, use magnifying glass.
- 4. Isolated colonies, microcolonies and hazes appearing in the zone of inhibition are indicative of hetero nature of the culture having resistant subpopulation in it. In such cases, consider reading for MIC determination at a point on the scale above which no resistant colonies are observed close to MIC strip (within 1-3 mm distance from the strip).
- 5. Since Ezy MIC[™] strip has continuous gradient, MIC values "in-between" two fold dilutions can be obtained.
- 6. Always round up these values to the next two-fold dilution before categorization. For example: Penicillin showing reading of 0.75 mcg/ml should be rounded up to next concentration ie. 1.0 mcg/ml.
- 7. If the ellipse intersects the strip in between 2 dilutions, read the MIC as the value which is nearest to the intersection.
- 8. When growth occurs along the entire strip, report the MIC as \geq the highest values on the MIC strip. When the inhibition ellipse is below the strip (does not intersect the strip), report the MIC < the lowest value on the MIC scale.

Warning and Precautions:

- 1. Ezy MICTM Strip is intended for *In vitro* diagnostic use only.
- 2. Although based on simple procedure, Ezy MICTM Strip should only be used by at least semi-trained personnel.
- 3. This strip is intended only for agar diffusion method and not for broth dilution method.
- 4. Ezy MICTM Strip should be used strictly according to procedures described herein.
- 5. Performance of Ezy MICTM Strips depends on use of proper inoculum and control cultures, recommended test medium and proper storage temperature.
- 6. Follow aseptic techniques and precautions against microbiological hazards should be used when handling bacterial or fungal specimen throughout the testing procedure.
- 7. Before using Ezy MICTM Strips, ensure that the strip is at room temperature.
- 8. When applying strips be steady. Do not move the strip once in contact with agar surface, since the antibiotic instantaneously diffuse on contact with agar.
- 9. Place the unused strips back to recommended temperature.

INTERPRETATION:

Use following interpretive criteria for susceptibility categorization.

| When testing | Incubation | Interpretive Criteria | | eria |
|---|---|-----------------------|--------|---------------|
| | | <u>< S</u> | Ι | <u>></u> R |
| Staphylococcus spp | 35-37°C for 18 hrs. | 0.12 | - | 0.25 |
| Enterococcus spp | 35-37°C for 18 hrs. | 8 | - | 16 |
| S.pneumoniae (Non Meningitis) | 35-37°C for 20-24hrs with 5% CO ₂ | 2 | 4 | 8 |
| S.pneumoniae (Meningitis) | 35-37°C for 20-24hrs with 5% CO ₂ | 0.06 | - | 0.12 |
| Streptococcus spps. Beta haemolytic group | 35-37°C for 20-24hrs with 5% CO ₂ | 0.12 | - | - |
| Streptococcus spps. Viridans group | 35-37°C for 20-24hrs with 5% CO ₂ | 0.12 | 0.25-2 | 4 |
| N. gonorrhoeae | 35-37°C for 20-24hrs with 5% CO ₂ | 0.06 | 0.12-1 | 2 |
| N. meningitidis | 35-37°C for 20-24hrs with 5% CO ₂ | 0.06 | 0.12- | 0.5 |
| | | | 0.25 | |
| Anaerobes | 35-37°C for 24-48hrs under anaerobic condition. | 0.5 | 1 | 2 |

QUALITY CONTROL

Quality control of Ezy MICTM Strip is carried out by testing the strips with standard ATCC cultures recommended by CLSI on suitable medium incubated appropriately.

Following are the reference MIC values (mcg/ml) range for Penicillin

| Organism | Medium used | Incubation | Std. Quality Control limits (mcg/ml) |
|----------------------------|---|----------------------------|---|
| S. aureus ATCC 29213 | Mueller Hinton Agar | 35-37°C for 18 hrs. | 0.25-0.5-1.0-2.0 |
| E. faecalis ATCC 29212 | Mueller Hinton Agar | 35-37°C for 18 hrs. | 1.0 - 2.0 - 4.0 |
| S. pneumoniae ATCC 49619 | Mueller Hinton Agar w/ | 35-37°C for 20-24hrs at 5% | 0.25 - 0.5 - 1.0 |
| | 5% Sheep Blood | CO ₂ | |
| Neisseria gonorrhoeae ATCC | GC Agar Base (M434) | 35-37°C for 20-24hrs at 5% | 0.25 - 0.5 - 1.0 |
| 49226 | with 1% defined growth supplement (FD025) | CO ₂ | |

Storage & Shelf Life:

- 1. Once the consignment is received, store applicators at Room Temperature and Ezy MIC[™] strips container at -20°C or below.
- 2. Use before expiry date on the label.
- 3. Ezy MICTM Strip left over from opened package must be kept dry.
- 4. Moisture should be prevented from penetrating into or forming within the package or storage container.
- 5. Check whether the batch number and expiry date are marked on the storage container.
- 6. Product performance is best within stated expiry period if correctly stored and handled.

Disposal

After use, Ezy MICTM Strips and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques (2, 3).

Limitation of Test

Ezy MICTM Strips provides *In vitro* MIC values, which provides only a possible insinuation of pathogens potential in *In vivo* susceptibility. These values can be considered as a guide to therapy selection only after taking into consideration several other factors; and must be the sole decision and responsibility of the physician along with the clinical experience in treating the infection. These tests are comparable to the standards as per the given specifications and set of experiment standards as far as possible. Please refer to CLSI standards for detailed limitation of susceptibility test on the clinical use of an antibiotic in various therapeutic conditions.

References:

- 1. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition, Vol. 1, Section 2.
- 2. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition, Vol. 3, Section 15.
- 3. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock., D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.
- Performance standards of Antimicrobial Susceptibility Testing; Twenty Ninth Informational Supplement. M100-S29, Vol. 39, No.1, Jan 2019.

Packing:

Each Pack contains following material packed in sealed glass vial with a desiccator capsule.

1) Penicillin Ezy MICTM strips (30/60/90/120/150 Strips per pack)

- 2) Applicator sticks
- 3) Package insert

Disclaimer :

CE

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia[™] publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia[™] Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. A-516,Swastik Disha Business Park,Via Vadhani Ind. Est., LBS Marg, Mumbai-400086, India. Customer care No.: 022-6147 1919 Email: techhelp@himedialabs.com Website: www.himedialabs.com



Egg Yolk Emulsion (50 ml/100 ml per vial)

Sterile stabilized emulsion of egg yolk recommended for use in various culture media. **Composition**

| Ingredients | Concentration | | |
|----------------|-------------------|------------------|--|
| | (100 ml per vial) | (50 ml per vial) | |
| Egg yolk | 30ml | 15ml | |
| Sterile saline | 70ml | 35ml | |

Directions:

Warm up the refrigerated egg yolk emulsion to room temperature. Shake well to attain uniform emulsion. (Since on refrigeration emulsion has a tendency to form layers or small lumps). Aseptically add 50 ml emulsion in 950 ml of sterile, molten, cooled (45-50°C) Baird Parker Agar Base <u>M043</u>/ Baird Parker Agar Base <u>M043</u>/ Baird Parker HiVegTM Agar Base <u>MV043</u>/ Baird Parker HiCynth[™] Agar MCD043/ Baird Parker Agar (Agar Medium O) <u>ME043</u>/ Baird Parker Agar Gar Medium O) <u>ME043</u>/ Baird Parker Agar Base, Granulated <u>GM0431</u>/ Baird Parker Agar Base <u>M0433</u>/ Baird Parker Agar Base <u>M0431</u>/ Baird Parker Agar Base

Aseptically add in 475 ml of sterile, molten, cooled (45-50°C) Bacillus Cereus Agar Base <u>M833</u>/ Bacillus Cereus HiVeg[™] Agar Base <u>MV833</u>/ Bacillus Cereus HiCynth[™] Agar Base <u>MCD833</u>

OR

Aseptically add 100 ml emulsion in 900 ml of sterile, molten, cooled (45-50°C) McClung Toabe Agar Base <u>M387</u>/ McClung Toabe HiVeg[™] Agar Base <u>MV387</u>/K.R.A.N.E.P. Agar Base <u>M583</u>/K.R.A.N.E.P. HiVeg[™] Agar Base <u>MV583</u> / MYP Agar Base (Phenol Red Egg Yolk Polymyxin Agar Base) <u>M636</u>/ M636S/ MYP HiVeg[™] Agar Base (Phenol Red Egg Yolk Polymyxin Agar Base) <u>M636</u>/ M636S/ MYP HiVeg[™] Agar Base (Phenol Red Egg Yolk Polymyxin Agar Base <u>MV636</u>/ MYP Agar Base, Granulated (Phenol Red Egg Yolk Polymyxin Agar Base, Granulated) <u>GM636</u> / MYP HiCynth[™]Agar Base (Phenol Red Egg Yolk Polymyxin HiCynth[™]Agar Base (Phenol Red Egg Yolk Polymyxin HiCynth[™] Agar Base) <u>MCD636</u>/ KG Agar Base <u>M658</u>/KG HiVeg[™] Agar Base <u>MV658</u>/ L.D. Egg Yolk Agar Base <u>M744</u>/ Egg Yolk Agar Base <u>M808</u> / Egg Yolk Agar Base, HiVeg[™] <u>MV808</u>/ Egg Yolk Agar Base, Modified <u>M1043</u> / Modified MYP Agar Base <u>M1139</u>/ Bacillus cereus Selective Agar Base (MYP) ISO 7932 <u>M1139</u>/ Modified MYP HiVeg[™] Agar Base <u>MV402</u>.

Aseptically add 450 ml of sterile, molten, cooled (45-50°C) in C. botulinum Isolation Agar Base $\underline{M911}$ / C. botulinum Isolation HiVegTM Agar Base $\underline{MV911}$

OR

Aseptically add 25 ml emulsion in 475 ml of sterile, molten, cooled (45-50°C) Perfringens Agar Base T.S.C./S.F.P.AgarBase) <u>M837</u>/ Perfringens Agar Base, Granulated (Tryptose Sulphite Cycloserine Agar Base, Granulated) (T.S.C./S.F.P. Agar Base, Granulated) <u>GM837</u>/ Perfringens HiCynth[™] Agar Base (T.S.C/S.F.P HiCynth[™] Agar Base) <u>MCD837</u>/ Perfringens HiVeg[™] Agar Base (T.S.C. / S.F.P. HiVeg[™] Agar Base) <u>MV837</u>/ S.F.P. Agar Base <u>M1005</u>/ S.F.P. HiVeg[™] Agar Base <u>MV1005</u>. OR

Aseptically add 80 ml emulsion in 920 ml of sterile, molten, cooled (45-50°C) Anaerobic Egg Agar Base M902 / Anaerobic Egg HiVeg[™] Agar Base MV902.

OR

Aseptically add 20 ml emulsion in 90 ml of sterile, molten, cooled (45-50°C) Polymyxin Pyruvate Egg Yolk Mannitol Bromothymol Blue Agar Base (PEMBA) <u>M1484</u>.

OR

Aseptically add 15 ml emulsion in 420 ml of sterile, molten, cooled (45-50°C) Willis and Hobb's Medium Base M1375. OR

Aseptically add 7ml of Emulsion in 93ml of sterile, molten, cooled (45-50°C) Lipovitellin Salt Mannitol Agar Base M627. OR

Aseptically add 2 vials of Clostridium Difficile Supplement (FD010), 40 ml of Egg Yolk Emulsion (FD045) together with 10 ml lysed horse blood in 1000 ml of sterile, molten, cooled (45-50°C) Clostridium Brazier Agar Base M1803 OR

Aseptically add 50ml of concetrated Egg yolk emulsion (<u>FD045</u>) and rehydrated contents of 1 vial of LM Selective Supplement (<u>FD330</u>) in 950 ml of sterile, molten, cooled (45-50°C) L.mono Selective Agar Base (LM Selective Agar Base) <u>M1994</u>.

Mix well and pour into sterile petri plates.

FD045

Type of specimen

Clinical samples - faeces, urine etc. ; Food samples

Specimen Collection and Handling

For clinical samples follow appropriate techniques for handling specimens as per established guidelines (1,2). For Food samples follow appropriate techniques for handling specimens as per established guidelines (3). After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning & Precautions

In Vitro diagnostic use. For professional use only. Read the label before opening the container. Wear protective gloves/ protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets.

Storage and Shelf Life

Store at 2 - 8°C. Use before expiry date on the label.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques (1,2).

Reference

1. Isenberg (Ed.),2004, Clinical Microbiology Procedures Handbook, Vol.3, American Society for Microbiology, Washington. D.C.

2. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock., D.W. (2015) Manual of Clinical Microbiology,11th Edition. Vol. 1.

3. Salfinger Y., and Tortorello M.L. Fifth (Ed.), 2015, Compendium of Methods for the Microbiological Examination of Foods, American Public Health Association, Washington, D.C.

* Not For Medicinal Use

Revision :02/2022



Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia[™] publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia[™] Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. Corporate Office : Plot No. C-40, Road No.21Y, MIDC, Wagle Industrial Area, Thane (W) - 400604, India. Customer care No.: 022-6147 1919 Email: techhelp@himedialabs.com Website: www.himedialabs.com



Tinsdale Selective Supplement (Part A & Part B)

A selective supplement recommended for the isolation and presumptive identification of *Corynebacterium diphtheriae*. **Composition**

Per vial sufficient for 1000 ml medium

| Ingredients | Concentration |
|---------------------|---------------|
| Part A | |
| Horse serum | 100ml |
| Part B | |
| Potassium tellurite | 1ml |

Directions:

Warm up the refrigerated contents of Part B vial and aseptically add 29 ml sterile distilled water. Mix thoroughly. Aseptically add warmed up (to 50°C) contents of Part A and B vials to sterile, molten, cooled (45-50°C) Tinsdale Agar Base $\underline{M314}$ / Tinsdale HiVegTM Agar Base $\underline{MV314}$ as required. Mix well and pour into sterile petri plates.

For 10 ml of M314 : 1.0 ml of Part A and 0.3 ml of Part B, is recommended.

Type of specimen

Clinical samples- Throat swab, nasal swab, wound swab, pus, etc.; Food samples

Specimen Collection and Handling

For clinical samples follow appropriate techniques for handling specimens as per established guidelines (1,2). For food samples follow appropriate techniques for handling specimens as per established guidelines (3). After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning & Precautions

In Vitro diagnostic use. For professional use only. Read the label before opening the container. Wear protective gloves/ protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets.

Storage and Shelf Life

Store at 2 - 8°C. Use before expiry date on the label.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques (1,2).

Reference

1.Isenberg (Ed.),2004, Clinical Microbiology Procedures Handbook, Vol.3, American Society for Microbiology, Washington. D.C.

2. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock., D.W. (2015) Manual of Clinical Microbiology,11th Edition. Vol. 1.

* Not For Medicinal Use

Revision : 02/2022

FD073



EC REP

HiMedia Laboratories Pvt. Limited, Plot No.C-40, Road No.21Y, MIDC, Wagle Industrial Area, Thane (W) -400604, MS, India

CEpartner4U, Esdoornlaan 13, 3951DB Maarn, NL www.cepartner4u.eu



IVD



_8°C Storage temperature

Do not use if package is damaged

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia[™] publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia[™] Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

In vitro diagnostic

medical device

HiMedia Laboratories Pvt. Ltd. Corporate Office : Plot No. C-40, Road No.21Y, MIDC, Wagle Industrial Area, Thane (W) - 400604, India. Customer care No.: 022-6147 1919 Email: techhelp@himedialabs.com Website: www.himedialabs.com





FD114

Vitamin K1 Supplement

A vitamin growth supplement used for the isolation of anaerobic organisms. Composition

Per vial sufficient for 1000 ml medium

Vitamin K1

Concentration 10mg

Directions:

Rehydrate the content of 1 vial aseptically with 5 ml sterile distilled water. Mix well and aseptically add it to 1000 ml sterile, molten, cooled (45-50°C) Anaerobic Blood Agar Base M975A.

If desired add rehydrated contents of 1 vial each of Vitamin K1 Supplement FD114 to Schaedler Broth M291/ Schaedler HiVegTM Agar MV291 for preparing Schaedler Agar w/Vitamin K1 or Schaedler HiVegTM Agar w/Vitamin K1 or Schaedler CNA Agar or Schaedler CNA HiVeg[™] Agar along with CNA Supplement <u>FD115</u> or Schaedler KV Agar or Schaedler KV HiVegTM Agar along with KV Supplement FD116. Add 50 ml sterile defibrinated sheep blood to all the above media. Mix well and pour into sterile petri plates.

Type of specimen

Clinical samples : stool, abscess, genital specimen, upper respiratory swab, endotracheal aspiration swab, etc.

Specimen Collection and Handling

For clinical samples follow appropriate techniques for handling specimens as per established guidelines (1,2). After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning & Precautions

In Vitro diagnostic use only. For professional use only. Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets.

Storage and Shelf Life

Store at 2 - 8°C. Use before expiry date on the label.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques (1,2).

Reference

1. Isenberg (Ed.),2004, Clinical Microbiology Procedures Handbook, Vol.3, American Society for Microbiology, Washington. D.C.

2. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock., D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.

* Not For Medicinal Use

Revision : 02/2023

Ingredients



EC REP

HiMedia Laboratories Pvt. Limited, Plot No.C-40, Road No.21Y, MIDC, Wagle Industrial Area, Thane (W) -400604, MS, India

CEpartner4U, Esdoornlaan 13, 3951DB Maarn, NL www.cepartner4u.eu



IVD



_8°C Storage temperature

Do not use if package is damaged

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia[™] publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia[™] Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

In vitro diagnostic

medical device

HiMedia Laboratories Pvt. Ltd. Corporate Office : Plot No. C-40, Road No.21Y, MIDC, Wagle Industrial Area, Thane (W) - 400604, India. Customer care No.: 022-6147 1919 Email: techhelp@himedialabs.com Website: www.himedialabs.com



Coagulase Plasma (0.1gm per vial)

It is recommended for studying coagulase reaction in diagnosis of Staphylococci.

Composition

Per vial sufficient for 6 tests medium

Ingredients

Coagulase Plasma

Directions:

Rehydrate the contents of one vial aseptically with 3 ml sterile distilled water. Add 0.5 ml of rehydrated FD248 in a tube. To this add approximately 0.05 ml of overnight broth culture of test organisms or 2-3 pure colonies picked from a noninhibitory agar plate. Mix gently & incubate at 37°C in incubator or water bath for up to 4 hours. Observe for clot formation in the tube at regular intervals. Any degree of clotting within 4 hours is considered as positive results.

Type of specimen

Clinical- skin, throat samples etc; Food samples

Specimen Collection and Handling

For Clinical & Food samples follow appropriate techniques for handling specimens as per established guidelines (1,2,3). After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning & Precautions

In Vitro diagnostic use. For professional use only. Read the label before opening the container. Wear protective gloves/ protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets.

Storage and Shelf Life

Store at 2-8°C. For unopened vial, use before the expiry date on the label. The rehydrated solution can be stored for up to 2 weeks at 2-8°C

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques (1,2).

Reference

1. Isenberg (Ed.),2004, Clinical Microbiology Procedures Handbook, Vol.3, American Society for Microbiology, Washington. D.C.

2. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock., D.W. (2015) Manual of Clinical Microbiology,11th Edition. Vol. 1.

3. Salfinger Y., and Tortorello M.L., 2015, Compendium of Methods for the Microbiological Examination of Foods, 5th Ed., American Public Health Association, Washington, D.C.

* Not For Medicinal Use

HiMedia Laboratories Pvt

IVD In vitro diagnostic 8°C Storage temperature Limited Plot No C-40 Road medical device No.21Y, MIDC, Wagle Industrial Area, Thane (W) -400604, MS, India EC REP Do not use if CE Marking CEpartner4U, Esdoornlaan 13. package is damaged 3951DB Maarn, NL www.cepartner4u.eu Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMediaTM publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.



Concentration

0.100g







FD248



Sabouraud Dextrose Agar, Granulated

GM063

Intended Use:

Recommended for the cultivation of yeasts, moulds and aciduric microorganisms from clinical and non-clinical samples.

| Composition** | |
|----------------------|-------------|
| Ingredients | Gms / Litre |
| Dextrose (Glucose) | 40.000 |
| Mycological, peptone | 10.000 |
| Agar | 15.000 |
| Final pH (at 25°C) | 5.6±0.2 |

**Formula adjusted, standardized to suit performance parameters

Directions

Suspend 65.0 grams in 1000 ml purified/distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C. Mix well and pour into sterile Petri plates.

Principle And Interpretation

Sabouraud Dextrose Agar is Carlier's modification (1) of the formulation described by is a modification of Sabouraud Dextrose Agar which is described by Sabouraud (2) for the cultivation of fungi (yeasts, moulds), particularly useful for the fungi associated with skin infections. This medium is also employed to determine microbial contamination in food, cosmetics, and clinical specimens (3). Mycological Peptone provides nitrogenous compounds. Dextrose provides an energy source. High dextrose concentration and low pH favors fungal growth and inhibits contaminating bacteria from test samples (4).

Type of specimen

Food and dairy samples ; Clinical samples: skin scrapings

Specimen Collection and Handling

For food and dairy samples follow appropriate techniques for handling specimens as per established guidelines (5,6,7). For clinical samples, follow appropriate techniques for sample collection and processing as per guidelines (8,9). After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions:

In Vitro diagnostic use. For professional use only. Read the label before opening the container. Wear protective gloves/ protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets.

Limitations

1. For heavily contaminated samples, the plate must be supplemented with inhibitory agents for inhibiting bacterial growth with lower pH.

2. Some pathogenic fungi may produce infective spores which are easily dispersed in air, so examination should be carried out in safety cabinet

3. Further biochemical tests should be carried out for confirmation.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the expiry period when stored at recommended temperature.

Quality Control

Appearance

Cream to yellow coloured granular media. **Gelling** Firm, comparable with 1.5% Agar gel. **Colour and Clarity of prepared medium** Light yellow to amber coloured clear to slightly opalescent gel forms in Petri plates.

Reaction

Reaction of 6.5% w/v aqueous solution at 25°C (after sterilization). pH : 5.6 \pm 0.2 pH 5.40-5.80

Cultural Response

Cultural response was observed after an incubation at 20-25°C for 24-48 hours.

| Organism | Inoculum (CFU) | Growth | Recovery |
|--|-------------------|----------------------------------|----------|
| Candida albicans ATCC 10231 (00054*) | 50 -100 | Luxuriant (white colonies) | >=70 % |
| #Aspergillus brasiliensis ATCC 16404 (00053*) | 50 -100 | luxuriant | >=70 % |
| Candida albicans ATCC 2091 (00055*) | 50 -100 | luxuriant | >=70 % |
| Saccharomyces cerevisiae ATCC 9763 (00058*) | 50 -100 | luxuriant | >=70 % |
| Escherichia coli ATCC 8739 (00012*) | 50 -100 | luxuriant | >=70 % |
| Escherichia coli ATCC 25922 (00013*) | 50 -100 | luxuriant | >=70 % |
| Escherichia coli NCTC 9002 | 50 -100 | luxuriant | >=70 % |
| Lactobacillus casei ATCC 334 | 50 -100 | luxuriant | >=70 % |
| <i>Trichophyton rubrum</i> ATCC 28191 | | luxuriant | |

Key : (*) - Corresponding WDCM numbers. (#) - Formerly known as Aspergillus niger

Storage and Shelf Life

Store between 10-30°C in a tightly closed container and the prepared medium at 20-30°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques (8,9).

Reference

1.Carlier G. I. M., 1948, Brit. J. Derm. Syph., 60:61.

2.Sabouraud K., 1892, Ann. Dermatol. Syphilol, 3:1061.

3.Bacteriological Analytical Manual, 8th Edition, Revision A, 1998. AOAC, Washington D.C.

4. Murray PR, Baren EJ, Jorgensen JH, Pfaller MA, Yolken RH (editors) 2003, Manual of clinical Microbiology, 8th ed., ASM, Washington, D.C.

5. American Public Health Association, Standard Methods for the Examination of Dairy Products, 1978, 14th Ed., Washington D.C.

6.Salfinger Y., and Tortorello M.L., 2015, Compendium of Methods for the Microbiological Examination of Foods, 5th Ed., American Public Health Association, Washington, D.C.

7. Wehr H. M. and Frank J. H., 2004, Standard Methods for the Microbiological Examination of Dairy Products, 17th Ed., APHA Inc., Washington, D.C.8.

8. Isenberg, H.D. Clinical Microbiology Procedures Handbook 2nd Edition.

9.Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock., D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.


Plot No.C-40, Road No.21Y, MIDC, Wagle Industrial Area,

CEpartner4U, Esdoornlaan 13,

3951DB Maarn, NL

www.cepartner4u.eu



HiMedia Laboratories Pvt. Limited, Thane (W) -400604, MS, India



In vitro diagnostic

medical device

IVD



-30°C Storage temperature

> Do not use if package is damaged

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMediaTM publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMediaTM Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. Corporate Office : Plot No.C-40, Road No.21Y, MIDC, Wagle Industrial Area, Thane (W) - 400604, India. Customer care No.: 022-6147 1919 Email: techhelp@himedialabs.com Website: www.himedialabs.com



Fluid Thioglycollate medium (Thioglycollate medium Fluid) M009

Intended use

Recommended for sterility testing of biologicals and for cultivation of aerobes, anaerobes and microaerophiles.

Composition**

| Ingredients | Gms / Litre |
|-----------------------|-------------|
| Tryptone | 15.000 |
| Yeast extract | 5.000 |
| Dextrose (Glucose) | 5.500 |
| Sodium chloride | 2.500 |
| L-Cystine | 0.500 |
| Sodium thioglycollate | 0.500 |
| Resazurin sodium | 0.001 |
| Agar | 0.750 |
| Final pH (at 25°C) | 7.1±0.2 |

**Formula adjusted, standardized to suit performance parameters

Directions

Suspend 29.75 grams in 1000 ml purified/distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 25°C and store in a cool dark place preferably below 25°C. Note : If more than the upper one-third of the medium has acquired a pink-purple colour, the medium may be restored once by heating in a water bath or in free flowing steam until the pink-purple colour disappears.

Principle And Interpretation

Brewer (1) formulated Fluid Thioglycollate Medium for rapid cultivation of aerobes as well as anaerobes including microaerophiles by adding a reducing agent and small amount of agar. The BP (2), EP (3), USP (4), and AOAC (5) have recommended the media for sterility testing of antibiotics, biologicals and foods and for determining the phenol coefficient and sporicidal effect of disinfectants. However, it is intended for the examination of clear liquid or water-soluble materials. Fluid Thioglycollate Medium is also routinely used to check the sterility of stored blood in blood banks (6). Dextrose, tryptone, yeast extract, L-cystine provide the growth factors necessary for bacterial multiplication. L-cystine and sodium thioglycollate allows Clostridium to grow in this medium even under aerobic conditions (7). Also the small amount of agar used in the medium favors the growth of aerobes as well as anaerobes in the medium, even if sodium thioglycollate is deleted from the medium(1). Sodium thioglycollate act as a reducing agent and neutralizes the toxic effects of mercurial preservatives and peroxides formed in the medium, thereby promoting anaerobiosis, and making the medium suitable to test materials containing heavy metals. (8,9). Any increase in the oxygen content is indicated by a colour change of redox indicator, resazurin to red (6,10,11). The small amount of agar helps in maintaining low redox potential for stabilizing the medium (9).

Type of specimen

Pharmaceutical samples for sterility testing, clinical samples- pus, wounds

Specimen Collection and Handling:

For clinical samples follow appropriate techniques for handling specimens as per established guidelines (12,13). For pharmaceutical samples, follow appropriate techniques for sample collection, processing as per guidelines (2,3,4) After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions

In Vitro diagnostic Use. For professional use only. Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/ face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets.

Limitations :

1. It is intended for the examination of clear liquid or water-soluble materials.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the expiry period when stored at recommended temperature.

Quality Control

Appearance

Cream to yellow homogeneous free flowing powder

Colour and Clarity of prepared medium

Light straw coloured, clear to slightly opalescent solution with upper 10% or less medium pink-purple on standing.

Reaction

Reaction of 2.97% w/v aqueous solution at 25°C. pH : 7.1±0.2

pН

6.90-7.30

Cultural Response

Cultural characteristics observed after an incubation at 30-35°C for not more than 3 days.

| Organism | Inoculum (CFU) | Growth |
|---|-------------------|-----------|
| Clostridium sporogenes ATCC 19404 (00008*) | 50 -100 | luxuriant |
| Clostridium sporogenes ATCC 11437 | 50 -100 | luxuriant |
| Clostridium perfringens ATCC 13124 (00007*) | 50 -100 | luxuriant |
| Bacteroides fragilis ATCC 23745 | 50 -100 | luxuriant |
| <i>Bacteroides vulgatus</i> ATCC 8482 | 50 -100 | luxuriant |
| Staphylococcus aureus subsp. aureus ATCC 25923 (00034*) | 50 -100 | luxuriant |
| Staphylococcus aureus subsp. aureus ATCC 6538 (00032*) | 50 -100 | luxuriant |
| Pseudomonas aeruginosa ATCC 27853 (00025*) | 50 -100 | luxuriant |
| Pseudomonas aeruginosa ATCC 9027 (00026*) | 50 -100 | luxuriant |
| Micrococcus luteus ATCC 9341 | 50 -100 | luxuriant |
| Streptococcus pneumoniae | 50 -100 | luxuriant |
| Escherichia coli ATCC 25922 (00013*) | 50 -100 | luxuriant |
| Escherichia coli ATCC 8739 (00012*) | 50 -100 | luxuriant |
| Escherichia coli NCTC 9002 | 50 -100 | luxuriant |
| Salmonella Typhimurium ATCC 14028 (00031*) | 50 -100 | luxuriant |
| Salmonella Abony NCTC 6017 (00029*) | 50 -100 | luxuriant |

Bacillus subtilis subsp. 50 -100 luxuriant

spizizenii ATCC 6633 (00003*)

Key : (*) Corresponding WDCM numbers.

Storage and Shelf Life

Store between 10-30°C in a tightly closed container and the prepared medium at 15-30°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition. Seal the container tightly after use. Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with sample must be decontaminated and disposed of in accordance with current laboratory techniques (12,13).

Reference

- 1. Brewer, 1940, J. Am. Med. Assoc., 115:598.
- 2. The British Pharmacopoeia, 2022, Medicines and Healthcare products Regulatory Agency.
- 3. European Pharmacopoeia, 2022, 10 th volume, European Directorate for the quality of medicines & Healthcare.
- 4. The United States Pharmacopoeia-National Formulatory (USP-NF), 2022

5. Williams H., (Ed.), 2005, Official Methods of Analysis of the Association of Official Analytical Chemists, 19th Ed., AOAC, Washington, D.C

- 6. Nungester, Hood and Warren, 1943, Proc. Soc. Exp. Biol. Med., 52:287.
- 7. Portwood, 1944, J. Bact., 48:255.
- 8. Federal Register, 1992, Fed. Regist., 21:640.

www.cepartner4u.eu

- 9. Quastel and Stephenson, 1926, J.Biochem., 20
- 10. Marshall, Gunnison and Luxen, 1940, Proc. Soc. Exp. Biol. Med., 43:672.
- 11. MacFaddin J.F., 1985, Media for Isolation-Cultivation-Identification-Maintenance of Medical Bacteria, Vol. 1, Williams and Wilkins, Baltimore.
- 12. Isenberg, H.D. Clinical Microbiology Procedures Handbook 2nd Edition.
- 13. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock., D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.

Revision : 06/2022



Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMediaTM publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMediaTM Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. Corporate Office : Plot No.C-40, Road No.21Y, MIDC, Wagle Industrial Area, Thane (W) - 400604, India. Customer care No.: 022-6147 1919 Email: techhelp@himedialabs.com Website: www.himedialabs.com



Triple Sugar Iron Agar

. .

Intended Use:

. . .

Recommended for the identification of gram-negative enteric bacilli on the basis of dextrose, lactose and sucrose fermentation and hydrogen sulphide production.

| Composition** | |
|---|---------------|
| Ingredients | Gms / Litre |
| Peptone | 10.000 |
| Tryptone | 10.000 |
| Yeast extract | 3.000 |
| HM Peptone B [#] | 3.000 |
| Lactose | 10.000 |
| Sucrose | 10.000 |
| Dextrose (Glucose) | 1.000 |
| Sodium chloride | 5.000 |
| Ferrous sulphate | 0.200 |
| Sodium thiosulphate | 0.300 |
| Phenol red | 0.024 |
| Agar | 12.000 |
| Final pH (at 25°C) | $7.4{\pm}0.2$ |
| **Formula adjusted, standardized to suit performance parameters | |

Equivalent to Beef extract

Directions

Suspend 64.52 grams in 1000 ml purified/distilled water. Heat to boiling to dissolve the medium completely. Mix well and distribute into test tubes. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Allow the medium to set in sloped form with a butt about 1 inch long.

Note: For better results, the medium can be sterilized by autoclaving at 10 lbs pressure (115°C) for 15 minutes.

Principle And Interpretation

Triple Sugar Iron Agar was originally proposed by Sulkin and Willett (1) and modified by Hajna (2) for identifying *Enterobacteriaceae*. This medium complies with the recommendation of APHA, for the examination of meat and food products (3), for the examination of milk and dairy products (4) and for microbial limit test for confirming the presence of *Salmonellae* (5,6) and in the identification of gram-negative bacilli (5,7).

Tryptone, peptone, yeast extract and HM peptone B provide nitrogenous compounds, sulphur, trace elements and vitamin B complex etc. Sodium chloride maintains osmotic equilibrium. Lactose, sucrose and dextrose are the fermentable carbohydrates. Sodium thiosulphate and ferrous ions make H₂S indicator system. Phenol red is the pH indicator. Organisms that ferment glucose produce a variety of acids, turning the colour of the medium from red to yellow. More amount of acids are liberated in butt (fermentation) than in the slant (respiration). Growing bacteria also form alkaline products from the oxidative decarboxylation of peptone and these alkaline products neutralize the large amounts of acid present in the butt. Thus the appearance of an alkaline (red) slant and an acid (yellow) butt after incubation indicates that the organism is a glucose fermenter but is unable to ferment lactose and/or sucrose. Bacteria that ferment lactose or sucrose (or both), in addition to glucose, produce large amounts of acid enables no reversion of pH in that region and thus bacteria exhibit an acid slant and acid butt. Gas production (CO₂) is detected by the presence of cracks or bubbles in the medium, when the accumulated gas escapes. Thiosulphate is reduced to hydrogen sulphide by several species of bacteria and H₂S combines with ferric ions of ferric salts to produce the insoluble black precipitate of ferrous sulphide. Reduction of thiosulphate proceeds only in an acid environment and blackening usually occurs in the butt of the tube. Triple Sugar Iron Agar should be used in parallel with Urea Agar/Broth (M112/M111) to distinguish between *Salmonella* and *Proteus* species. The reactions can be summarized as follows:

M021

Alkaline slant / acid butt-only glucose fermented

Acid slant / acid butt-glucose and sucrose fermented or glucose and lactose fermented or all the three sugars, glucose, lactose and sucrose fermented.

Bubbles or cracks present-gas production Black precipitate present-H₂S gas production

Type of specimen

Pure bacterial isolate from water, food, or clinical sample.

Specimen Collection and Handling:

For clinical samples follow appropriate techniques for handling specimens as per established guidelines (8,9).

For food and dairy samples, follow appropriate techniques for sample collection and processing as per guidelines (3,4). For water samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards (5). After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions

In Vitro diagnostic Use. For professional use only. Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/ face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets.

Limitations :

1.Some members of the *Enterobacteriaceae* and H_2S producing *Salmonella* may not be H_2S positive on TSI Agar. Some bacteria may show H_2S production on Kligler Iron Agar but not on TSI Agar. This can happen because utilization of sucrose in TSI Agar suppresses the enzymic pathway that result in H_2S production.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the expiry period when stored at recommended temperature.

Quality Control

Appearance

Light yellow to pink homogeneous free flowing powder

Gelling

Firm, comparable with 1.2% Agar gel. Colour and Clarity of prepared medium

Pinkish red coloured clear to slightly opalescent gel forms in tubes as slants.

Reaction

Reaction of 6.45% w/v aqueous solution at 25°C. pH : 7.4±0.2

pН

7.20-7.60

Cultural Response

Cultural characteristics observed after an incubation at 35-37°C for 18-24 hours.

| Organism | Inoculum (CFU) | Growth | Slant | Butt | Gas | H_2S |
|---|-------------------|-----------|--|--|----------------------|---|
| <i>Citrobacter freundii</i> ATCC 8090 | 50-100 | luxuriant | acidic reaction, yellowing of the medium | acidic reaction, yellowing of the medium | positive reaction | positive, blackening of medium |
| # Klebsiella aerogenes ATCC 13048 (00175*) | 50-100 | luxuriant | acidic reaction, yellowing of the medium | , acidic reaction, yellowing of the medium | positive reaction | negative, no blackening of medium |
| Escherichia coli ATCC 25922 (00013*) | 50-100 | luxuriant | acidic reaction, yellowing of the medium | acidic reaction, yellowing of the medium | positive reaction | negative, no blackening of medium |
| Klebsiella pneumoniae ATCC 13883 (00097*) | 50-100 | luxuriant | acidic reaction, yellowing of the medium | acidic reaction, yellowing of the medium | positive reaction | negative, no blackening of medium |

HiMedia Laboratories

| Proteus vulgaris ATCC 13315 | 50-100 | luxuriant | alkaline reaction, red colour of the medium | acidic reaction, yellowing of the medium | negative reaction | positive, blackening of medium |
|--|--------|-----------|--|--|----------------------|---|
| Salmonella Paratynhi A | 50-100 | luxuriant | alkaline | acidic reaction, | positive | negative, no |
| ATCC 9150 | | | reaction, red colour of the medium | yellowing of the medium | reaction | blackening of medium |
| <i>Salmonella</i> Typhi ATCC 6539 | 50-100 | luxuriant | alkaline reaction, red colour of the medium | acidic reaction, yellowing of the medium | negative reaction | positive, blackening of medium |
| <i>Salmonella</i> Typhimurium ATCC 14028 (00031*) | 50-100 | luxuriant | alkaline reaction, red colour of the medium | acidic reaction, yellowing of the medium | positive reaction | positive, blackening of medium |
| Shigella flexneri ATCC 12022 (00126*) | 50-100 | luxuriant | alkaline reaction, red colour of the medium | acidic reaction, yellowing of the medium | negative reaction | negative, no blackening of medium |
| <i>Escherichia coli</i> ATCC 8739 (00012*) | 50-100 | luxuriant | acidic reaction, yellowing of the medium | , acidic reaction, yellowing of the medium | positive reaction | negative, no blackening of medium |
| <i>Escherichia coli</i> NCTC 9002 | 50-100 | luxuriant | acidic reaction, yellowing of the medium | , acidic reaction, yellowing of the medium | positive reaction | negative, no blackening of medium |
| Klebsiella pneumoniae ATCC 10031 | 50-100 | luxuriant | acidic reaction, yellowing of the medium | , acidic reaction, yellowing of the medium | positive reaction | negative, no blackening of medium |

Key : (*) Corresponding WDCM numbers.(#) Formerly known as Enterobacter aerogenes

Storage and Shelf Life

Store between 10-30°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with sample must be decontaminated and disposed of in accordance with current laboratory techniques (8,9).

Reference

1.Sulkin E.S. and Willett J.C., 1940, J. Lab. Clin. Med., 25:649.

2.Hajna A.A., 1945, J. Bacteriol, 49:516.

3.Salfinger Y., and Tortorello M.L., 2015, Compendium of Methods for the Microbiological Examination of Foods, 5th Ed., American Public Health Association, Washington, D.C.

4.Wehr H. M. and Frank J. H., 2004, Standard Methods for the Microbiological Examination of Dairy Products, 17th Ed., APHA Inc., Washington, D.C.

5.Baird R.B., Eaton A.D., and Rice E.W., (Eds.), 2015, Standard Methods for the Examination of Water and Wastewater, 23rd ed., APHA, Washington, D.C.

Please refer disclaimer Overleaf.

6.Finegold S. M. and Baron E. J., 1986, Bailey and Scotts Diagnostic Microbiology, 7th Ed., The C.V. Mosby Co., St. Louis.

7.MacFaddin J., 1985, Media for Isolation-Cultivation-Identification-Maintenance of Medical Bacteria, Vol. 1, Williams and Wilkins, Baltimore.

8. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.

9.Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock., D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.

Revision : 03/ 2022



Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMediaTM publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMediaTM Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. Corporate Office : Plot No.C-40, Road No.21Y, MIDC, Wagle Industrial Area, Thane (W) - 400604, India. Customer care No.: 022-6147 1919 Email: techhelp@himedialabs.com Website: www.himedialabs.com



Endo Agar, Special

M029R

Endo Agar, Special is recommended for the detection of coliform and other enteric organisms.

| Composition** | |
|-------------------------|-------------|
| Ingredients | Gms / Litre |
| Peptone, special | 11.500 |
| Lactose | 12.900 |
| Dipotassium phosphate | 0.480 |
| Monopotassium phosphate | 0.220 |
| Sodium chloride | 3.600 |
| Sodium sulphite | 0.860 |
| Sodium lauryl sulphate | 0.010 |
| Basic fuchsin | 0.830 |
| Agar | 9.600 |
| Final pH (at 25°C) | 7.3±0.2 |
| | |

**Formula adjusted, standardized to suit performance parameters

Directions

Suspend 40.0 grams in 1000 ml distilled water. Boil to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C. Mix well and pour into sterile Petri plates.

Caution: Basic Fuchsin is a potential Carcinogen and care should be taken to avoid inhalation of the powdered dye and contamination of the skin.

Principle And Interpretation

Endo (1) had first developed a culture medium for differentiation of lactose fermentors and non-fermenters and further developed as todays Endo Agar (2). Endo agar is used for microbiological examination of potable water and waste water, dairy products and food (3,4,5).

Sodium sulphite and basic fuchsin has inhibitory effect on gram-positive microorganisms. Sodium Lauryl sulphate inhibits many organisms other than coliforms. Lactose fermenting coliforms produce aldehyde and acid. The aldehyde in turn liberates fuchsin from the fuchsin-sulphite complex, giving rise to a red colouration of colonies. With *Escherichia coli* this reaction is very pronounced that the fuchsin crystallises, exhibiting to the colonies a permanent greenish metallic lustre (fuchsin lustre). The phosphates buffer the medium. Peptone special provides essential nutrients especially nitrogenous for the coliforms.

Quality Control

Appearance Light pink to purple homogeneous free flowing powder

Gelling
Firm, comparable with 0.96% Agar gel.
Colour and Clarity of prepared medium
Pink Clear to slightly opalescent gel with a slight precipitate forms in Petri plates.
Reaction
Reaction of 4.0% w/v aqueous solution at 25°C. pH : 7.3±0.2
pH
7.10-7.50
Cultural Response
Cultural characteristics observed after an incubation at 35-37°C for 18-24 hours.

Cultural Response

Please refer disclaimer Overleaf.

| Organism | Growth | Inoculum (CFU) | Recovery | Colour of Colony |
|--------------------------------------|----------------|-------------------|----------|--|
| Cultural Response | | | | |
| Bacillus subtilis ATCC 6633 | inhibited | >=103 | 0% | |
| Enterobacter aerogenes ATCC 13048 | good-luxuriant | 50-100 | >=50% | pink |
| Enterococcus faecalis ATCC 29212 | none-poor | 50-100 | <=10% | pink, small |
| Escherichia coli ATCC 25922 | good-luxuriant | 50-100 | >=50% | pink to rose red with metallic sheen |
| Klebsiella pneumoniae ATCC 13883 | good-luxuriant | 50-100 | >=50% | pink, mucoid |
| Salmonella Typhi ATCC 6539 | good-luxuriant | 50-100 | >=50% | colourless to pale pink |
| Staphylococcus aureus ATCC 25923 | inhibited | >=103 | 0% | |
| Pseudomonas aeruginosa ATCC 27853 | good-luxuriant | 50-100 | >=50% | colourless, irregular |
| Proteus vulgaris ATCC 13315 | good-luxuriant | 50-100 | >=50% | colourless to pale pink |
| Shigella sonnei ATCC 25931 | good-luxuriant | 50-100 | >=50% | colourless to pale pink |

Storage and Shelf Life

Store below 30°C in tightly closed container and prepared medium at 2 - 8°C away from light to avoid photo-oxidation. Use before expiry date on the label.

Reference

1.Endo, 1904, Zentralbl. Bakteriol., Abt. 1., Orig., 35:109.

2.Levin and Schoenlein, 1930, A Compilation of Culture Media for the Cultivation of Microorganisms, Williams and Wilkins, Baltimore.

3.Greenberg, Trussell and Clesceri (ed.), 1998, Standard Methods for the Examination of Water and Wastewater, 20th ed., APHA, Washington, D.C.

4.Richardson (ed.), 1992, Standard Methods for the Examination of Dairy Products, 16th ed., APHA, Washington, D.C.

5. Speck M., 1984, Compendium of Methods for the Microbiological Examination of Foods, 3rd ed., APHA, Washington, D.C.

Revision : 01 / 2014

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia[™] publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia[™] Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. A-516, Swastik Disha Business Park, Via Vadhani Ind. Est., LBS Marg, Mumbai-400086, India. Customer care No.: 022-6147 1919 Email: techhelp@himedialabs.com Website: www.himedialabs.com

CE



Kligler Iron Agar Intended Use:

Recommended for differential identification of gram-negative enteric bacilli from clinical and non-clinical samples on the basis of the fermentation of glucose (dextrose), lactose and hydrogen sulphide production.

| Composition** | |
|--|---------------|
| Ingredients | Gms / Litre |
| Peptone | 15.000 |
| HM Peptone B # | 3.000 |
| Yeast extract | 3.000 |
| Proteose peptone | 5.000 |
| Lactose | 10.000 |
| Dextrose | 1.000 |
| Ferrous sulphate | 0.200 |
| Sodium chloride | 5.000 |
| Sodium thiosulphate | 0.300 |
| Phenol red | 0.024 |
| Agar | 15.000 |
| Final pH (at 25°C) | $7.4{\pm}0.2$ |
| **Formula adjusted, standardized to suit performance | ce parameters |

- Equivalent to Beef extract

Directions

Suspend 57.52 grams in 1000 ml purified/distilled water. Heat to boiling to dissolve the medium completely. Mix well and distribute into tubes. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Allow the tubes to cool in slanted position to form slopes with about 1 inch butts. Best reactions are obtained on freshly prepared medium. Do not use screw capped tubes or bottles.

Note: Avoid overheating otherwise it may produce precipitate in the medium.

Principle And Interpretation

Kligler Iron Agar is a combination of the lead acetate medium described by Kligler (1,2) and Russels Double Sugar Agar (3) and is used as a differentiation medium for typhoid, dysentery and allied bacilli (4). Bailey and Lacey substituted phenol red for andrade indicator previously used as pH indicator (4). Kligler Iron Agar differentiates lactose fermenters from the non-fermenters. It differentiates Salmonella Typhi from other Salmonellae and also Salmonella Paratyphi A from Salmonella Scottmuelleri and Salmonella Enteritidis (5). Fermentation of dextrose results in production of acid, which turns the indicator from red to yellow. Since there is little sugar i.e. dextrose, acid production is very limited and therefore a reoxidation of the indicator is produced on the surface of the medium, and the indicator remains red. However, when lactose is fermented, the large amount of acid produced, avoids reoxidation and therefore the entire medium turns yellow. Kligler Iron Agar, in addition to Peptone, HM peptone B and yeast extract, contains lactose and glucose (dextrose), which enables the differentiation of species of enteric bacilli. Phenol red is the pH indicator, which exhibits a colour change in response to acid produced during the fermentation of sugars. The combination of ferrous sulphate and sodium thiosulphate enables the detection of hydrogen sulfide production, which is evidenced by a black color either throughout the butt, or in a ring formation near the top of the butt. Lactose non-fermenters (e.g., Salmonella and Shigella) initially produce a yellow slant due to acid produced by the fermentation of the small amount of glucose (dextrose). When glucose (dextrose) supply is exhausted in the aerobic environment of the slant, the reaction reverts to alkaline (red slant) due to oxidation of the acids produced. The reversion does not occur in the anaerobic environment of the butt, which therefore remains acidic (yellow butt). Lactose fermenters produce yellow slants and butts because of lactose fermentation. The high amount of acids thus produced helps to maintain an acidic pH under aerobic conditions. Tubes showing original colour of the medium indicates the fermentation of neither glucose (dextrose) nor lactose. Gas production (aerogenic reaction) is detected as individual bubbles or by splitting or displacement of the agar by the formation of cracks in the butt of the medium.

Pure cultures of suspected organisms from plating media such as MacConkey Agar (M081), Bismuth Sulphite Agar (M027) or Deoxycholate Citrate Agar (M065), SS Agar (M108) etc. are inoculated on Kligler Iron Agar for identification.

Type of specimen

Isolated microorganism from clinical, food, dairy and water samples.

M078

Specimen Collection and Handling

For water samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards (6). For food and dairy samples, follow appropriate techniques for sample collection and processing as per guidelines (7,8,9). For clinical samples follow appropriate techniques for handling specimens as per established guidelines (10,11). After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions

In Vitro diagnostic use. For professional use only. Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets.

Limitations

- 1. Results should be noted after 18-24 hours to avoid erroneous results.
- 2. Straight wire loop should be used for inoculation.
- 3. Pure isolates should be used to avoid erroneous results.
- 4. Other biochemical and serological tests must be performed for complete identification

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the expiry period when stored at recommended temperature.

Quality Control

Appearance

Light yellow to pink homogeneous free flowing powder

Gelling

Firm, comparable with 1.5% Agar gel

Colour and Clarity of prepared medium

Red coloured, clear to slightly opalescent gel forms in tubes as slants

Reaction

Reaction of 5.75% w/v aqueous solution at 25°C. pH : 7.4±0.2

pН

7.20-7.60

Cultural Response

Cultural characteristics observed after an incubation at 35-37°C for 18 - 48 hours.

| Organism | Inoculum (CFU) | Growth | Gas | H2S | Slant | Butt |
|--|-------------------|-----------|----------------------|---|--|--|
| Escherichia coli ATCC 25922 (00013*) | 50-100 | luxuriant | positive reaction | negative reaction, no blackening of medium | acidic reaction, yellowing of the medium | acidic reaction, yellowing of the medium |
| #Klebsiella aerogenes ATCC 13048 (00175*) | 50-100 | luxuriant | positive reaction | negative reaction, no blackening of medium | acidic reaction, yellowing of the medium | acidic reaction, yellowing of the medium |
| <i>Citrobacter freundii</i> ATCC 8090 | 50-100 | luxuriant | positive reaction | positive reaction, blackening of medium | acidic reaction, yellowing of the medium | acidic reaction, yellowing of the medium |
| Proteus vulgaris ATCC 6380 | 50-100 | luxuriant | negative reaction | positive reaction, blackening of medium | alkaline reaction, red colour of the medium | acidic reaction, yellowing of the medium |
| Klebsiella pneumoniae ATCC 13883 (00087*) | 50-100 | luxuriant | positive reaction | negative reaction,no blackening of medium | acidic reaction, yellowing of the medium | acidic reaction, yellowing of the medium |
| <i>Salmonella</i> Paratyphi A ATCC 9150 | 50-100 | luxuriant | positive reaction | negative reaction,no blackening of medium | alkaline reaction, red colour of the medium | acidic reaction, yellowing of the medium |

Please refer disclaimer Overleaf.

| <i>Salmonella</i> Schottmuelleri ATCC 10719 | 50-100 | luxuriant | positive reaction | positive reaction, blackening of medium | alkaline reaction, red colour of the medium | acidic reaction, yellowing of the medium |
|--|--------|-----------|----------------------|--|--|---|
| Salmonella Typhi ATCC 6539 | 50-100 | luxuriant | negative reaction | positive reaction, blackening of medium | alkaline reaction, red colour of the medium | acidic reaction, yellowing of the medium |
| Salmonella Enteritidis ATCC 13076 (00030*) | 50-100 | luxuriant | positive reaction | positive reaction, blackening of medium | alkaline reaction, red colour of the medium | acidic reaction, yellowing of the medium |
| Shigella flexneri ATCC 12022 (00126*) | 50-100 | luxuriant | negative reaction | negative reaction,no blackening of medium | alkaline reaction, red colour of the medium | acidic reaction, yellowing of the medium |
| Pseudomonas aeruginosa ATCC 27853 (00025*) | 50-100 | luxuriant | negative reaction | negative reaction, blackening of medium | alkaline reaction, red colour of the medium | alkaline reaction,red colour of the medium |
| Yersinia enterocolitica ATCC 27729 | 50-100 | luxuriant | variable reaction | negative reaction,no blackening of medium | alkaline reaction,red colour of the medium | acidic reaction, yellowing of the medium |
| Enterobacter cloacae ATCC 13047 (00083*) | 50-100 | luxuriant | positive reaction | negative reaction,no blackening of medium | acidic reaction, yellowing of the medium | acidic reaction, yellowing of the medium |

Key :* Corresponding WDCM numbers

Storage and Shelf Life

Store between 10-30°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques (10,11).

Reference

1.Kligler I. J., 1917, Am. J. Publ. Health, 7:1041.

2.Kligler I. J., 1918, J. Exp. Med., 28:319.

3.Russell F. F., 1911, J. Med. Res., 25:217.

4.Bailey S. F. and Lacey G. R., 1927, J. Bacteriol., 13:183.

5.Ewing, 1986, Edwards and Ewings Identification of the Enterobacteriaceae, 4th Ed., Elsevier Science Publishing Co., Inc., N.Y.

6.Lipps WC, Braun-Howland EB, Baxter TE, eds. Standard methods for the Examination of Water and Wastewater, 24th ed. Washington DC:APHA Press; 2023.

7. American Public Health Association, Standard Methods for the Examination of Dairy Products, 1978, 14th Ed., Washington D.C.

8.Salfinger Y., and Tortorello M.L., 2015, Compendium of Methods for the Microbiological Examination of Foods, 5th Ed., American Public Health Association, Washington, D.C.

9.Wehr H. M. and Frank J. H., 2004, Standard Methods for the Microbiological Examination of Dairy Products, 17th Ed., APHA Inc., Washington, D.C.

10. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.

11.Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock., D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1





HiMedia Laboratories Pvt. Limited, Plot No.C-40, Road No.21Y, MIDC, Wagle Industrial Area, Thane (W) -400604, MS, India

CEpartner4U, Esdoornlaan 13,

3951DB Maarn, NL

www.cepartner4u.eu



IVD



-30°C Storage temperature

Do not use if package is damaged

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMediaTM publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMediaTM Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

In vitro diagnostic

medical device

HiMedia Laboratories Pvt. Ltd. Corporate Office : Plot No.C-40, Road No.21Y, MIDC, Wagle Industrial Area, Thane (W) - 400604, India. Customer care No.: 022-6147 1919 Email: techhelp@himedialabs.com Website: www.himedialabs.com



Simmons Citrate Agar Intended Use:

Recommended for differentiation the members of *Enterobacteriaceae* on the basis of citrate utilization from clinical and non clinical samples.

| Com | | ** |
|-----|----------|-----|
| Com | position | ~ ~ |

| Ingredients | Gms / Litre |
|-------------------------------|-------------|
| Magnesium sulphate | 0.200 |
| Ammonium dihydrogen phosphate | 1.000 |
| Dipotassium phosphate | 1.000 |
| Sodium citrate | 2.000 |
| Sodium chloride | 5.000 |
| Bromothymol blue | 0.080 |
| Agar | 15.000 |
| Final pH (at 25°C) | 6.8±0.2 |

**Formula adjusted, standardized to suit performance parameters

Directions

Suspend 24.28 grams in 1000 ml purified/ distilled water. Heat, to boiling, to dissolve the medium completely. Mix well and distribute in tubes or flasks. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes.

Precaution: Before using water, ensure pH of water is 6.5 to 7.0. Initial colour of the medium may deviate from expected colour, if the above precaution is ignored.

Principle And Interpretation

These media are used for the differentiation between *Enterobacteriaceae* and the members of aerogenes group on the basis of citrate utilization as sole carbon source. Initially the citrate medium was developed by Koser (1) containing ammonium salt as the only nitrogen source and citrate as the only carbon source for differentiating *Escherichia coli* and *Enterobacter aerogenes* by IMViC tests. Later on Simmons (2) modified Kosers formulation by adding agar and bromothymol blue (3). It is recommended by APHA (4).

Ammonium dihydrogen phosphate and sodium citrate serve as the sole nitrogen and carbon source respectively. Microorganisms also use inorganic ammonium salts as their sole nitrogen source. Metabolism of these salts causes the medium to become alkaline, indicated by a change in colour of the pH indicator from green to blue. Bromothymol blue is the pH indicator. The medium should be freshly prepared because in dry conditions, changes in colour may appear even before inoculation, especially at the bottom of the slant.

Type of specimen

Isolated microorganism from clinical and non clinical samples.

Specimen Collection and Handling

For water samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards (5). For food and dairy samples, follow appropriate techniques for sample collection and processing as per guidelines (6,7,8). For clinical samples follow appropriate techniques for handling specimens as per established guidelines (9,10). After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions

In Vitro diagnostic use. For professional use only. Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets.

Limitations

1.Before using water, ensure pH of water is 6.5 to 7.0.Initial colour of the medium may deviate from expected colour, if the above precaution is ignored.

2. The pH affects the performance of the medium and must be correctly monitored.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the expiry period when stored at recommended temperature.

M099

Quality Control

Appearance

Cream to yellow homogeneous free flowing powder. Gelling Firm, comparable with 1.5% Agar gel Colour and Clarity of prepared medium Forest green coloured slightly opalescent gel forms in tubes as slants. Reaction Reaction of 2.43% w/v aqueous solution at 25°C. pH : 6.8±0.2.

pН

6.60-7.00

Cultural Response

Cultural characteristics observed after an incubation at 35-37°C for 18-24 hours.

| Organism | Inoculum (CFU) | Growth | Citrate utilisation |
|---|-------------------|----------------|---------------------------------------|
| # Klebsiella aerogenes ATCC 13048 (00175*) | 50-100 | good-luxuriant | positive reaction, blue colour |
| Escherichia coli ATCC 25922 (00013*) | >=10 ⁴ | inhibited | |
| Salmonella Typhi ATCC 6539 | 50-100 | fair-good | negative reaction, green colour |
| Salmonella Typhimurium ATCC 14028 (00031*) | 50-100 | good-luxuriant | positive reaction, blue colour |
| Shigella dysenteriae ATCC 13313 | >=10 ⁴ | inhibited | |
| Salmonella Choleraesuis ATCC 12011 | 50-100 | good-luxuriant | positive reaction, blue colour |
| Salmonella Enteritidis ATCC 13076 (00030*) | 50-100 | good-luxuriant | positive reaction, blue colour |

Key: * Corresponding WDCM numbers

Formerly known as Enterobacter aerogenes

Storage and Shelf Life

Store between 10-30°C in a tightly closed container and the prepared medium at 2 - 8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques (9,10).

Reference

- 1. Koser, 1923, J. Bact., 8:493.
- 2. Simmons, 1926, J. Infect. Dis., 39:209.

3. MacFaddin J., 1985, Media for Isolation-Cultivation-Identification-Maintenance of Medical Bacteria, Vol. 1, Williams and Wilkins, Baltimore.

4. Eaton A. D., Clesceri L. S., Rice E. W., and Greenberg A W., (Eds.), 2005, Standard Methods for the Examination of Water and Wastewater, 21st Ed., APHA, Washington, D.C.

5. Lipps WC, Braun-Howland EB, Baxter TE, eds. Standard methods for the Examination of Water and Wastewater, 24th ed. Washington DC:APHA Press; 2023.

6. American Public Health Association, Standard Methods for the Examination of Dairy Products, 1978, 14th Ed., Washington D.C.

HiMedia Laboratories

7. Salfinger Y., and Tortorello M.L., 2015, Compendium of Methods for the Microbiological Examination of Foods, 5th Ed., American Public Health Association, Washington, D.C.

8. Wehr H. M. and Frank J. H., 2004, Standard Methods for the Microbiological Examination of Dairy Products, 17th Ed., APHA Inc., Washington, D.C.

9. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.

10. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock., D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.

Revision: 04/2022



Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMediaTM publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMediaTM Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. Corporate Office : Plot No.C-40, Road No.21Y, MIDC, Wagle Industrial Area, Thane (W) - 400604, India. Customer care No.: 022-6147 1919 Email: techhelp@himedialabs.com Website: www.himedialabs.com



SS Agar (Salmonella Shigella Agar)

Intended Use:

Recommended for the isolation of *Salmonella* and some *Shigella* species from pathological specimens, suspected foodstuffs etc.

Composition**

| Ingredients | Gms / Litre |
|---|-------------|
| Peptone | 5.000 |
| HM peptone B # | 5.000 |
| Lactose | 10.000 |
| Bile salts mixture | 8.500 |
| Sodium citrate | 10.000 |
| Sodium thiosulphate | 8.500 |
| Ferric citrate | 1.000 |
| Brilliant green | 0.00033 |
| Neutral red | 0.025 |
| Agar | 15.000 |
| Final pH (at 25°C) | 7.0 ± 0.2 |
| **Formula adjusted, standardized to suit performance parameters | |

- Equivalent to Beef extract

Directions

Suspend 63.02 grams in 1000 ml purified /distilled water. Boil with frequent agitation to dissolve the medium completely. **DO NOT AUTOCLAVE OR OVERHEAT**. Overheating may destroy selectivity of the medium. Cool to about 50°C. Mix and pour into sterile Petri plates.

Principle And Interpretation

SS Agar medium is recommended as differential and selective medium for the isolation of *Salmonella* and *Shigella* species from pathological specimens (1) and suspected foodstuffs (2,3,4,5) and for microbial limit test (6). SS Agar is a moderately selective medium in which gram-positive bacteria are inhibited by bile salts, brilliant green and sodium citrate.

Peptone, HM peptone B provides nitrogen and carbon source, long chain amino acids, vitamins and essential growth nutrients. Lactose is the fermentable carbohydrate. Brilliant green, bile salts and thiosulphate selectively inhibit gram-positive and coliform organisms. Sodium thiosulphate is reduced by certain species of enteric organisms to sulphite and H_2S gas and this reductive enzyme process is attributed by thiosulphate reductase. Production of H_2S gas is detected as an insoluble black precipitate of ferrous sulphide, formed upon reaction of H_2S with ferric ions or ferric citrate, indicated in the center of the colonies.

The high selectivity of Salmonella Shigella Agar allows the use of large inocula directly from faeces, rectal swabs or other materials suspected of containing pathogenic enteric bacilli. On fermentation of lactose by few lactose-fermenting normal intestinal flora, acid is produced which is indicated by change of colour from yellow to red by the pH indicator-neutral red. Thus these organisms grow as red pigmented colonies. Lactose non-fermenting organisms grow as translucent colourless colonies with or without black centers. Growth of *Salmonella* species appears as colourless colonies with black centers resulting from H_2S production. *Shigella* species also grow as colourless colonies which do not produce H_2S .

Type of specimen

Clinical: faeces, rectal swabs; Suspected food stuffs.

Specimen Collection and Handling

For clinical samples follow appropriate techniques for handling specimens as per established guidelines (8,9). For food and dairy samples, follow appropriate techniques for sample collection and processing as per guidelines (2,3,4,5). After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions

In Vitro diagnostic use. For professional use only. Read the label before opening the container. Wear protective gloves/ protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets.

M108

Limitations

1. The medium is highly selective and may be toxic to certain *Salmonella* or *Shigella* species. Hence it is recommended to use to inoculate plates of less inhibitory media parallel to SS Agar, such as Hektoen Enteric Agar (M467) or Deoxycholate Citrate Agar (M065) for easier isolation of *Shigella* species (9).

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the expiry period when stored at recommended temperature.

Quality Control

Appearance

Light yellow to pink homogeneous free flowing powder

Gelling

Firm, comparable with 1.5% Agar gel

Colour and Clarity of prepared medium

Reddish orange coloured clear to slightly opalescent gel forms in Petri plates

Reaction

Reaction of 6.3% w/v aqueous solution at 25°C. pH : 7.0±0.2

pН

6.80-7.20

Cultural Response

Cultural characteristics observed after an incubation at 35-37°C for 18-24 hours.

| Organism | Inoculum (CFU) | Growth | Recovery | Colour of colony |
|---|-------------------|----------------|----------|-----------------------------------|
| # Klebsiella aerogenes ATCC 13048 (00175*) | 50-100 | fair | 20-30% | cream pink |
| Escherichia coli ATCC 25922 (00013*) | 50-100 | fair | 20-30% | pink with bile precipitate |
| Salmonella Choleraesuis ATCC 12011 | 50-100 | good-luxuriant | >=50% | colourless with |
| Salmonella Typhi ATCC 6539 | 50-100 | good-luxuriant | >=50% | colourless with |
| <i>Enterococcus faecalis</i> ATCC 29212 (00087*) | 50-100 | none-poor | <=10% | colourless |
| Proteus mirabilis ATCC 25933 | 50-100 | fair-good | 30-40% | colourless, may have black centre |
| Shigella flexneri ATCC 12022 (00126*) | 50-100 | good | 40-50% | colourless |
| Salmonella Typhimurium ATCC 14028 (00031*) | 50-100 | good-luxuriant | >=50% | colourless with black centre |
| Salmonella Enteritidis ATCC 13076 (00030*) | 50-100 | good-luxuriant | >=50% | colourless with black centre |

Key: *Corresponding WDCM numbers.

Formerly known as Enterobacter aerogenes

Storage and Shelf Life

Store between 10-30°C in a tightly closed container and the prepared medium at 20-30°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques (8,9).

Reference

1. Lennette and others (Eds.), 1985, Manual of Clinical Microbiology, 4th ed., ASM, Washington, D.C.

2. Lipps WC, Braun-Howland EB, Baxter TE, eds. Standard methods for the Examination of Water and Wastewater, 24th ed. Washington DC:APHA Press; 2023.

3. Salfinger Y., and Tortorello M.L. Fifth (Ed.), 2015, Compendium of Methods for the Microbiological Examination of Foods, 5th Ed., American Public Health Association, Washington, D.C.

4. Wehr H. M. and Frank J. H., 2004, Standard Methods for the Microbiological Examination of Dairy Products, 17th Ed., APHA Inc., Washington, D.C.

5. Williams S., (Ed.), 2005, Official Methods of Analysis of the Association of Official Analytical Chemists, 19th Ed., AOAC, Washington, D.C.

6. The United States Pharmacopoeia-National Formulatory (USP-NF), 2022.

7. MacFaddin J., 1985, Media for Isolation-Cultivation-Identification-Maintenance of Medical Bacteria, Vol. I, Williams and Wilkins, Baltimore.

8. Isenberg, H.D. Clinical Microbiology Procedures Handbook 2nd Edition.

9. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock., D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.

Revision : 04/2022



Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMediaTM publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMediaTM Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. Corporate Office : Plot No.C-40, Road No.21Y, MIDC, Wagle Industrial Area, Thane (W) - 400604, India. Customer care No.: 022-6147 1919 Email: techhelp@himedialabs.com Website: www.himedialabs.com



Columbia Blood Agar Base

Intended Use:

For preparation of blood agar, chocolate agar and for preparation of various selective and identification media and isolation of organisms from clinical and non clinical samples.

Composition**

| Ingredients | Gms / Litre |
|---------------------|-------------|
| Peptone, special | 23.000 |
| Corn starch | 1.000 |
| Sodium chloride | 5.000 |
| Agar | 15.000 |
| Final pH (at 25°C) | 7.3±0.2 |

**Formula adjusted, standardized to suit performance parameters

Directions

Suspend 44.0 grams of in 1000 ml purified/ distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C before adding heat sensitive compounds. For Blood Agar: Add 5% v/v sterile defibrinated sheep blood to sterile cool base.

For Chocolate Agar: Add 10% v/v sterile defibrinated sheep blood to sterile cool base. Heat to 80°C for 10 minutes with constant agitation.

The medium can be made selective by adding different antimicrobials to sterile base.

For *Brucella* species: Add rehydrated contents of 1 vial of NPBCVN Selective Supplement (FD005) to 500 ml sterile molten base.

For *Campylobacter* species: Add rehydrated contents of 1 vial of Blaser-Wang Selective Supplement (FD006) or Butzler Selective Supplement (FD007) or Skirrow Selective Supplement (FD008) or VTCA Selective Supplement (FD090) or Butzler VI Selective Supplement (FD106) to 500 ml sterile molten base along with rehydrated contents of 1 vial of Minerals Growth Supplement (FD009) and 5-7% v/v horse or sheep blood.

For *Gardnerella* species: Add rehydrated contents of 1 vial of GNA Selective Supplement (FD056) to 500 ml sterile molten base.

For Cocci: Add rehydrated contents of 1 vial of NC Selective Supplement (FD030) or NNP Selective Supplement (FD031) or CO Selective Supplement (FD119) to 500 ml sterile molten base.

Principle And Interpretation

Columbia Blood Agar Base was devised by Ellner et al (1). This medium contains special peptone which supports rapid and luxuriant growth of fastidious and non-fastidious organisms. Also, this medium promotes typical colonial morphology; better pigment production and more sharply defined haemolytic reactions. Fildes found that Nutrient Agar supplemented with a digest of sheep blood supplied both of these factors and the medium would support the growth of *H. influenzae* (2,3). The inclusion of bacitracin makes the enriched Columbia Agar Medium selective for the isolation of *Haemophilus* species from clinical specimens, especially from upper respiratory tract (4). Columbia Agar Base is used as the base for the media containing blood and for selective media formulations in which different combinations of antimicrobial agents are used as additives.

Corn starch serves as an energy source and also neutralizes toxic metabolites. Sheep blood permits the detection of haemolysis and also provides heme (X factor) which is required for the growth of many bacteria. However it is devoid of V factor (Nicotinamide adenine dinucleotide) and hence *Haemophilus influenzae* which needs both the X and V factors, will not grow on this medium.

Columbia Agar Base with added sterile serum provides an efficient medium for *Corynebacterium diphtheriae* virulence test medium. After following the established technique for *C. diphtheriae*, lines of toxin-antitoxin precipitation are clearly visible in 48 hours. Many pathogens require carbon dioxide; therefore, plates may be incubated in an atmosphere containing approximately 3-10% CO₂.

Precaution: Brucella cultures are highly infective and must be handled carefully; incubate in 5-10% CO₂. Campylobacter species are best grown at 42°C in a micro aerophillic atmosphere. Plates with Gardenerella supplements plates should be incubated at 35°C for 48 hours containing 7% CO₂ (2).

Type of specimen

Clinical samples : throat swabs, pus.

M144

Specimen Collection and Handling

For clinical samples follow appropriate techniques for handling specimens as per established guidelines (5,6). After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions

In Vitro diagnostic use only. For professional use only. Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets.

Limitations

1. Certain fastidious organisms like *Haemophilus influenzae* may not grow on the medium, blood supplementation may be required.

2. As this medium have a relatively high carbohydrate content, beta-hemolytic *Streptococci* may exhibit a greenish hemolytic reaction which may be mistaken for the alpha haemolysis.

3. Biochemical characterization is required on colonies of pure culture for complete identification.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the expiry period when stored at recommended temperature.

Quality Control

Appearance

Cream to yellow homogeneous free flowing powder.

Gelling

Firm, comparable with 1.5% Agar gel.

Colour and Clarity of prepared medium

Basal medium: Light amber coloured clear to slightly opalescent gel.

After addition of 5%w/v sterile defibrinated blood : Cherry red coloured opaque gel forms in Petri plates.

Reaction

Reaction of 4.4% w/v aqueous solution at 25°C. pH : 7.3±0.2

pН

7.10-7.50

Cultural Response

Cultural characteristics observed with added 5% w/v sterile defibrinatedblood, after an incubation at 35-37°C for 24-48 hours.

| Organism | Inoculum (CFU) | Growth | Recovery | Haemolysis |
|--|-------------------|-----------|----------|--------------|
| Neisseria meningitidis ATCC 13090 | 50-100 | luxuriant | >=70% | none |
| Staphylococcus aureus subsp. aureus ATCC 25923 (00034*) | 50-100 | luxuriant | >=70% | beta / gamma |
| Staphylococcus epidermidis ATCC 12228 (00036*) | 50-100 | luxuriant | >=70% | gamma |
| Staphylococcus aureus subsp. aureus ATCC 6538 (00032*) | 50-100 | luxuriant | >=70% | beta / gamma |
| Streptococcus pneumoniae ATCC 6303 | 50-100 | luxuriant | >=70% | alpha |
| Streptococcus pyogenes | 50-100 | luxuriant | >=70% | beta |
| Clostridium sporogenes ATCC 19404 (00008*) | 50-100 | luxuriant | >=50 % | |
| Clostridium sporogenes ATCC 11437 | 50-100 | luxuriant | >=50 % | |
| Clostridium perfringens ATCC 13124 (00007*) | 50-100 | luxuriant | >=50 % | |
| Clostridium perfringens ATCC 12934 | 50-100 | luxuriant | >=50 % | |
| | | | | |

Key : (*) Corresponding WDCM numbers.

Please refer disclaimer Overleaf.

Storage and Shelf Life

Store between 10-30°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques (5,6).

Reference

- 1. Ellner P. P., Stoessel C. J., Drakeford E. and Vasi F., 1966, Am. J. Clin. Pathol., 45:502.
- 2. Fildes P., 1920, Br. J. Exp. Pathol., 1:129.
- 3. Fildes P., 1921, Br. J. Exp. Pathol., 2:16.
- 4. Chapin K. C. and Doern G. V., 1983, J. Clin. Microbiol., 17:1163.
- 5. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.

6. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock., D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.

Revision : 04/2022



Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMediaTM publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMediaTM Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. Corporate Office : Plot No.C-40, Road No.21Y, MIDC, Wagle Industrial Area, Thane (W) - 400604, India. Customer care No.: 022-6147 1919 Email: techhelp@himedialabs.com Website: www.himedialabs.com



Mueller Hinton Agar

M173

Intended Use:

Recommended for determination of susceptibility of microorganisms to antimicrobial agents isolated from clinical samples.

| Composition** | |
|---|-------------|
| Ingredients | Gms / Litre |
| HM infusion B from # | 300.000 |
| Acicase ## | 17.500 |
| Starch | 1.500 |
| Agar | 17.000 |
| Final pH (at 25°C) | 7.3±0.1 |
| **Formula adjusted, standardized to suit performance parameters | |

- Equivalent to Beef infusion from

- Equivalent to Casein acid hydrolysate

Directions

Suspend 38.0 grams in 1000 ml purified/ distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C. Mix well and pour into sterile Petri plates. Note: The performance of this batch has been tested and standardised as per the current CLSI (formerly, NCCLS) document M6-protocols for Evaluating Dehydrated Mueller Hinton Agar.

Principle And Interpretation

The Mueller Hinton formulation was originally developed as a simple, transparent agar medium for the cultivation of pathogenic *Neisseria* species (1). Other media were subsequently developed that replaced the use of Mueller Hinton Agar for

the cultivation of pathogenic *Neisseria* species, but it became widely used in the determination of sulfonamide resistance of gonococci and other organisms. Mueller Hinton Agar is now used as a test medium for antimicrobial susceptibility testing (2). Mueller Hinton Agar is recommended for the diffusion of antimicrobial agents impregnated on paper disc through an agar gel as described in CLSI Approved Standard (3). Mueller Hinton Agar has been selected by the CLSI for several reasons:

i. It demonstrates good batch-to-batch reproducibility for susceptible testing.

- ii. It is low in sulfonamide, trimethoprim and tetracycline inhibitors.
- iii. It supports the growth of most non-fastidious bacterial pathogens and
- iv. Many data and much experience regarding its performance have been recorded (4).

Kirby-Bauer et al recommended this medium for performing antibiotic susceptibility tests using a single disc of high concentration (5). WHO Committee on Standardization of Susceptibility Testing has accepted Mueller Hinton Agar for determining the susceptibility of microorganisms because of its reproducibility (6). Mueller Hinton Agar with 5% sheep blood and Mueller Hinton Agar with Hemoglobin have been recommended for antimicrobial susceptibility testing of *Streptococcus pneumoniae* and *Haemophilus influenzae*.

HM infusion B from and acicase provide nitrogenous compounds, carbon, sulphur and other essential nutrients. Starch acts as a protective colloid against toxic substances present in the medium. Starch hydrolysis yields dextrose, which serves as a source of energy. These ingredients are selected for low thymine and thymidine content as determined by MIC values for

Enterococcus faecalis with sulfamethoxazole trimethoprim (SXT).

The Kirby-Bauer procedure is based on agar diffusion of antimicrobial substances impregnated on paper discs. This method employs disc with a single concentration of antimicrobial agent and the zone diameters observed are correlated with minimum inhibitory concentration (MIC) values (7,1,2). A standardized suspension of the organism is swabbed over the entire surface of the medium.

Paper discs impregnated with specific amounts of antimicrobial agents are then placed on the surface of the medium, incubated and zones of inhibition around each disc are measured. The susceptibility is determined by comparing with CLSI standards (4). The various factors, which influence disc diffusion susceptibility tests, are agar depth, disc potency, inoculum concentration, pH of the medium and beta-lactamase production by test organisms (4,8).

Mueller Hinton Agar is not appropriate for assay by disc diffusion method with slow growing organisms, anaerobes and capnophiles. With slow growing organisms, increased incubation may cause deterioration of diffusing antibiotic and produce unprecise readings (9).

Type of specimen

Clinical samples : Isolated microorganisms from urine , stool etc.

Specimen Collection and Handling

For clinical samples follow appropriate techniques for handling specimens as per established guidelines (3,5). After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions

In Vitro diagnostic use only. For professional use only. Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets.

Limitations

1. This medium is recommended for susceptibility testing of pure cultures only.

2. Inoculum density may affect the zone size. Heavy inoculum may result in smaller zones or too less inoculum may result in bigger zones.

3. Fastidious organisms may not grow on this medium and may require supplementation of blood.

4. Fastidious anaerobes may not grow on this medium.

5. As antimicrobial susceptibility is carried with antibiotic disc, proper storage of the disc is desired which may affect the potency of the disc.

6. Under certain circumstances, the in vitro results of antibiotic susceptibility may not show the same in vivo.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the expiry period when stored at recommended temperature.

Quality Control

Appearance

Cream to yellow homogeneous free flowing powder.

Gelling

Firm, comparable with 1.7% agar gel.

Colour and Clarity of prepared medium

Light amber coloured clear to slight opalscent gel froms in Petri plates.

Reaction

Reaction of 3.8% w/v aqueous solution at 25°C. pH : 7.3±0.1

pН

7.20-7.40

Cultural Response

Cultural characteristics observed after incubation at 30-35°C for 18 -24 hours for bacterial cultures.

For testing *S. pneumoniae* : The medium was supplemented with 5% Sheep blood and incubated at 35°C for 16-18 hours at 5% CO₂.

For testing *H. influenaze* : The medium was supplemented with 5g/l of Yeast extract & 2 vials /l of Haemophilus Growth Supplement (FD117 containing 15 mg/l of Haematin + 15 mg/l of NAD) and incubated at 35°C for 20-24 hours at 5% CO₂.

Antibiotic Sensitivity test

Various discs were tested for standard ATCC strains and zone of inhibition were measured after an incubation 30-35°C for 18 hours. (As per the latest CLSI Protocol M6 & Standards as per the current CLSI M100).

Thymine/Thymidine Content

The zones for these discs are indicative of the Thymine/Thymidine content of the medium.

Divalent Cation Content

\$ The zones for these discs are indicative of the Divalent Cation content of the medium

Thymine/Thymidine Content

The zones for these discs are indicative of the Thymine/Thymidine content of the medium. **Divalent Cation Content**

\$ The zones for these discs are indicative of the Divalent Cation content of the medium

| Organism | Growth | Standard Zone Zone of | |
|---|-----------|-------------------------|------------|
| | | | inhibition |
| | | | Observed |
| <i>Escherichia coli</i> ATCC 25922 (00013*) | luxuriant | | |
| Cephalothin CEP 30mcg | | 29-37 mm | 29 -37 mm |
| Chloramphenicol C 30 mcg | | 21-27 mm | 21 -27 mm |
| Co-Trimoxazole COT 25 | | 23-29 mm | 23 -29 mm |
| mcg # | | | |
| Cefotaxime CTX 30 mcg | | 29-35 mm | 29 -35 mm |
| Gentamicin GEN 10 mcg | | 19-26 mm | 19 -26 mm |
| Sulphafurazole SF 300 mcg | | 15-23 mm | 15 -23 mm |
| Staphylococcus aureus | huvuriant | | |
| subsp. <i>aureus</i> ATCC 25923 (00034*) | luxurlain | | |
| Co-Trimoxazole COT 25 | | # 20 mm (Clear | >=20 mm |
| mcg # | | zone) | |
| Cefoxitin CX 30 mcg | | 23-29 mm | 23 -29 mm |
| Erythromycin E 15 mcg | | 22-30 mm | 22 -30 mm |
| Linezolid LZ 30 mcg | | 25-32 mm | 25 -32 mm |
| Oxacillin OX 1mcg | | 18-24 mm | 18 -24 mm |
| Pristinomycin RP 15 mcg | | 21-28 mm | 21 -28 mm |
| Tetracycline TE 30 mcg \$ | | 18-25 mm | 18 -25 mm |
| Ciprofloxacin CIP 5mcg | | 22-30 mm | 22 -30 mm |
| Pseudomonas aeruginosa | hummingt | | |
| ATCC 27853 (00025*) | luxuriant | | |
| Ceftazidime CAZ 30 mcg | | 22-29 mm | 22 -29 mm |
| Ciprofloxacin CIP 5mcg | | 30-40 mm | 30 -40 mm |
| Tobramycin TOB 10 mcg \$ | | 19-25 mm | 19 -25 mm |
| Amikacin AK 30 mcg \$ | | 18-26 mm | 18 -26 mm |
| Aztreonam AT 3mcg | | 23-29 mm | 23 -29 mm |
| Cephotaxime CTX 30 mcg | | 18-22 mm | 18 -22 mm |
| Gentamicin GEN 10 mcg \$ | | 16-21 mm | 16 -21 mm |
| Imipenem IPM 10 mcg | | 20-28 mm | 20 -28 mm |
| Piperacillin PI 100 mcg | | 12-18 mm | 25 -33 mm |
| <i>Escherichia coli</i> ATCC 35218 | luxuriant | | |
| Amoxyclay AMC 30 mcg | | 18-24 mm | 18 -24 mm |
| Piperacillin/Tazobactam PIT | | 24-30 mm | 24 -30 mm |
| 100/10 mcg | | | |
| Ticarcillin TI 75 mcg | | 6 mm | 6 -6 mm |
| Ticarcillin/Clavulanic acid | | 20-28 mm | 20 -28 mm |
| TCC 75/10mcg | | | |
| Ampicillin AMP 10 mcg | | 16-22 mm | 16 -22 mm |
| Ampicillin/Sulbactam A/S | | 29-37 mm | 29 -37 mm |
| 10/10 mcg | | | |
| Enterococcus faecalis | luxuriant | | |
| ATCC 29212 (00087*) | luxuilant | | |
| Trimethoprim TR 5 mcg # | | # 20 mm | >=20 mm |
| Vancomycin VA 30 mcg | | 17-21 mm | 17 -21 mm |
| Staphylococcus aureus | luxuriant | | |
| subsp. <i>aureus</i> ATCC | | | |
| 43300 (MRSA) (00211*) | | | |
| Oxacillin OX 1 mcg | | Very Hazy to No Zone | No zone |

Key: *Corresponding WDCM numbers.

Storage and Shelf Life

Store between 10-30°C in a tightly closed container and the prepared medium at 20-30°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use.

Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques (3,5).

Reference

1. Murray P. R., Baron J. H., Pfaller M. A., Jorgensen J. H. and Yolken R. H., (Ed.), 2003, Manual of Clinical Microbiology, 8th Ed., American Society for Microbiology, Washington, D.C.

2. National Committee for Clinical Laboratory Standards, 2000, Approved Standard: M7-A5. Methods for Dilution Antimicrobial Susceptibility Tests for Bacteria that grow aerobically, 5th Ed., NCCLS, Wayne, Pa.

3. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.

4. MacFaddin J. F., 1985, Media for Isolation-Cultivation-Identification-Maintenance of Medical Bacteria, Vol. 1, Williams and Wilkins, Baltimore

5. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock., D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.

6. Present Status and Future Work, WHO Sponsored collaborative study, Chicago, Oct. 1967.

7. Ericsson H. M. and Sherris J. L., 1971, Acta Pathol. Microbiol., Scand. Sect B Suppl., 217:1.

8. National Committee for Clinical Laboratory Standards, 1986, Proposed Standards, M6-P, NCCLS, Villanova, Pa.

9. Mueller J. H. and Hinton J., 1941, Proc. Soc. Exp. Biol. Med., 48:330.

Revision : 05/2022



Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMediaTM publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMediaTM Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.



Anaerobic Agar

M228

Intended Use:

Recommended for the cultivation of anaerobic bacteria, especially *Clostridium* species and other anaerobic organisms from clinical and non-clinical samples.

Composition**

| Ingredients | Gms / Litre |
|--|---------------|
| Tryptone | 20.000 |
| Dextrose (Glucose) | 10.000 |
| Sodium chloride | 5.000 |
| Sodium thioglycollate | 2.000 |
| Sodium formaldehyde Sulfoxylate | 1.000 |
| Methylene blue | 0.002 |
| Agar | 20.000 |
| Final pH (at 25°C) | 7.2 ± 0.2 |
| **Formula adjusted standardized to suit performance parameters | |

Directions

Suspend 58.0 grams in 1000 ml purified/ distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C. Mix well and pour into sterile Petri plates.

Principle And Interpretation

Anaerobic Agar was originally designed for surface cultivation of members of the genus *Clostridium* and other anaerobic organisms on plates (1). This medium is suitable for isolation of facultative and obligate anaerobes and for the study of colonial morphology as colonies can be readily seen on the light coloured agar and are easily accessible (2,3). Anaerobic bacteria vary in their sensitivity to oxygen and nutritional requirements (3). Anaerobic bacteria lack cytochromes and thus are unable to use oxygen as a terminal electron acceptor (4).

This medium contains sodium thioglycollate and sodium formaldehyde sulphoxylate that provide adequate anaerobiosis which is indicated by methylene blue present in the medium which yields blue colour to medium in presence of oxygen. Tryptone and dextrose provide essential nutrients while sodium chloride maintains osmotic equilibrium.

Dispense 50-60 ml medium per 95 x 20 mm plate. For best results, use porous tops for the plates during solidification to get the dry surface. Inoculation can be done by streaking or smearing. Cover the inoculated plate with sterile Brewer Anaerobic Petri dish cover. Incubate aerobically, as desired. When standard plates are used, dispense 0.1 to 1.0 ml of inoculum into plates and mix with 20 - 25 ml of sterile medium. After solidification, incubate anaerobically as required by particular organism under study. Methylene blue is inhibitory to some anaerobic microorganisms.

Type of specimen

Clinical- stool, abscess

Specimen Collection and Handling

For clinical samples follow appropriate techniques for handling specimens as per established guidelines (3,5). After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions

In Vitro diagnostic use only. For professional use only. Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets.

Limitations

1. Ensure that the clinical samples are properly transported under anaerobic conditions.

- 2. Proper anaerobic conditions must be maintained for optimal recovery of organisms
- 3.Methylene blue is toxic to certain anaerobes.
- 4. Further biochemical and serological tests must be performed for confirmation.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the expiry period when stored at recommended temperature.

Quality Control

Appearance

Cream to yellow homogeneous free flowing powder

Gelling

Firm, comparable with 2.0% agar gel.

Colour and Clarity of prepared medium

Light amber coloured, clear to slightly opalescent gel forms in Petri plates that becomes greenish due to aeration on standing **Reaction**

Reaction of 5.8% w/v aqueous solution at 25°C. pH : 7.2±0.2

pН

7.00-7.40

Cultural Response

Cultural characteristics observed under anaerobic condition after an incubation at 35-37°C for 48-72 hours.

| Organism | Inoculum (CFU) | Growth | Recovery | |
|--|-------------------|----------------|----------|--|
| <i>Clostridium perfringens</i> ATCC 12924 | 50-100 | good-luxuriant | >=50% | |
| Clostridium sporogenes ATCC 11437 | 50-100 | good-luxuriant | >=50% | |
| Clostridium butyricum ATCC 13732 | 50-100 | good-luxuriant | >=50% | |

Storage and Shelf Life

Store between 10-30°C in a tightly closed container and the prepared medium at 20-30°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques (3,5).

Reference

- 1. Brewer J. H., 1942, Science, 95:587.
- 2. Vera J., 1942, J. Bacteriol., 44:497.
- 3. Isenberg (Ed.), 1992, Clinical Microbiology Procedures Handbook, American Society for Microbiology, Washington, D.C.
- 4. Baron E. J., Peterson and Finegold S. M., Bailey & Scotts Diagnostic Microbiology, 9th Ed., 1994, Mosby-Year Book, Inc., St. Louis, Mo.
- 5. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock., D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.



Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMediaTM publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMediaTM Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.



Tinsdale Agar Base

M314

Tinsdale Agar Base with supplement is used for selective isolation and differentiation of Corynebacterium diphtheriae .

| Composition** | |
|--------------------------------|-------------|
| Ingredients | Gms / Litre |
| Peptic digest of animal tissue | 20.000 |
| Sodium chloride | 5.000 |
| L-Cystine | 0.240 |
| Sodium thiosulphate | 0.430 |
| Agar | 15.000 |
| Final pH (at 25°C) | 7.4±0.2 |
| | |

**Formula adjusted, standardized to suit performance parameters

Directions

Suspend 40.67 grams in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 50°C and aseptically add Diphtheria Virulence Supplement (FD073, Part A and Part B). Mix well and pour into sterile Petri plates.

Principle And Interpretation

The Corynebacteria are gram-positive, non-sporulating, non-motile rods. They are often club-shaped and frequently banded or beaded with irregularly stained granules. These bacteria are generally aerobic or facultative, but microaerophilic species do occur. *Corynebacterium diphtheriae* produces a powerful exotoxin that causes diphtheria in humans. In nature, *C. diphtheriae* occurs in nasopharyngeal area of infected persons or healthy carriers.

The three biotypes of *C. diphtheriae* are *mitis*, *intermedius* and *gravis* (6). The signs and symptoms of diphtheria are sore throat, malaise, headache and nausea (2). Tinsdale Agar Base Medium was developed by Tinsdale (1) for the selective isolation and differentiation of *C. diphtheriae* from diphtheroids. This medium was modified by Billings (2), which improved the recovery and differential qualities of *C. diphtheriae*. The present medium is according to the modified Billings Medium. Moore and Parsons (3) confirmed the halo formation as a characteristic property of *C. diphtheriae* with the exception of *C. ulcerans*, which forms colony with similar features as *C. diphtheriae*.

Peptic digest of animal tissue provides nitrogenous compounds. L-cystine and sodium thiosulphate form the H2S indicator system. Potassium tellurite from the supplement inhibits all gram-negative bacteria and most of the upper respiratory tract normal flora.

C. diphtheriae forms grayish black colonies surrounded by a dark brown halo while diphtheroids commonly found in the upper respiratory tract do not form such colonies. Dark brown halo around the colony is due to H2S production from cystine combining with the tellurite salt. Moore and Parsons (3) found Tinsdale Medium as an ideal medium for the routine cultivation and isolation of *C. diphtheriae*. They also confirmed the stability of halo formation on clear medium and its specificity for *C. diphtheriae* and *C. ulcerans*. *C. ulcerans* found in nasopharynx form colonies same as *C. diphtheriae* and require further biochemical confirmation (4).

Do not incubate the plates in 5-10% CO2 as it retards the development of characteristic halos (5). Tinsdale Agar is not suitable as a primary plating medium, since it may not support the growth of some strains of *C. diphtheriae* (6). *C. ulcerans*

, *C. pseudotuberculosis* and (rarely) *Staphylococcus* species may produce a characteristic halo on Tinsdale Agar (6). Several organisms may exhibit slight browning on Tinsdale Agar in 18 hours; therefore the plates should be read after complete incubation period (48 hours) (6).

Quality Control Appearance

Please refer disclaimer Overleaf.

Cream to yellow homogeneous free flowing powder

Gelling

Firm, comparable with 1.5% Agar gel

Colour and Clarity of prepared medium

Light amber coloured clear to slightly opalescent gel forms in Petri plates

Reaction

Reaction of 4.07% w/v aqueous solution at 25°C. pH : 7.4±0.2

pН

7.20-7.60

Cultural Response

M314: Cultural characteristics observed after an incubation at 35-37°C for 40-48 hours with added Diptheria Virulence Supplement (FD073, Part A and Part B).

| Organism | Inoculum (CFU) | Growth | Recovery | Colony characteristics |
|---|-------------------|----------------|----------|----------------------------------|
| Corynebacterium diphtheriae type gravis | 50-100 | good-luxuriant | >=50% | brown-black with halo |
| Corynebacterium diphtheriae type interme dius | 50-100 | good-luxuriant | >=50% | brown-black with halo |
| Corynebacterium diphtheriae type mitis | 50-100 | good-luxuriant | >=50% | brown-black with halo |
| Klebsiella pneumoniae ATCC 13883 | >=103 | inhibited | 0 % | |
| Streptococcus pyogenes ATCC 19615 | 50-100 | good | 40-50% | black pin point, without halo |

Storage and Shelf Life

Store below 30°C in tightly closed container and the prepared medium at 2 - 8°C. Use before expiry date on the label.

Reference

1. Tinsdale G. F. W., 1947, J. Pathol. Bacteriol., 59:461.

2. Billings E., 1956, An investigation of Tinsdale Tellurite Medium: its usefulness and mechanisms of halo-formation, M.S. thesis, University of Michigan, Ann Arbor, Mich.

3. Moore M. S. and Parsons E. I., 1958, J. Infect. Dis., 102:88.

4. MacFaddin J. F., 1985, Media for Isolation-Cultivation-Identification-Maintenance of Medical Bacteria, Vol. I, Williams and Wilkins, Baltimore.

5. Murray P. R., Baron E. J., Jorgensen J. H., Pfaller M. A., Yolken R. H., (Eds.), 8th Ed., 2003, Manual of Clinical Microbiology, ASM, Washington, D.C

6. Isenberg, (Eds.), 1992, Clinical Microbiology Procedures Handbook, Vol. 1, American Society for Microbiology, Washington, D.C.

Revision : 1 / 2011

CE

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia[™] publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia[™] Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. A-516,Swastik Disha Business Park,Via Vadhani Ind. Est., LBS Marg, Mumbai-400086, India. Customer care No.: 022-6147 1919 Email: techhelp@himedialabs.com



Technical Data

Bile Esculin Azide Agar

Intended Use:

For selective isolation and presumptive identification of faecal Streptococci.

Composition** Ingredients Tryptone HM peptone B # Proteose peptone Bile

| Bile ## | | 10.000 |
|------------------------------|---------------------------|---------|
| Esculin | | 1.000 |
| Ferric ammonium citrate | | 0.500 |
| Sodium chloride | | 5.000 |
| Sodium azide | | 0.150 |
| Agar | | 15.000 |
| Final pH (at 25°C) | | 7.1±0.2 |
| # Equivalent to Beef extract | ## - Equivalent to Oxgall | |

**Formula adjusted, standardized to suit performance parameters

Directions

Suspend 56.65 grams in 1000 ml purified/distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C. Mix well and pour into sterile Petri plates.

Gms / Litre

17.000

5.000 3.000

Caution: Sodium azide has a tendency to form explosive metal azides with plumbing materials. It is advisable to use enough water to flush off the disposables.

Principle And Interpretation

Group D Streptococci possess the group D lipoteichoic acid antigen in their cell walls. Former Group D species, which are predominant normal inhabitants of the human gastrointestinal tract, are termed as faecal Streptococci or Enterococci (1). The unique ability of Enterococci to split esculin was reported by Meyer and Schonfeld (2). Enterococci and Group D Streptococci hydrolyse esculin to esculetin and dextrose, which reacts with ferric citrate producing brownish black precipitate (3). The use of esculin hydrolysis in identification of Enterococci was first cited by Rochaix (4). Bile Esculin Agar was originally formulated by Swan (5) for the isolation and identification of Group D Streptococci from food. Facklam and Moody (6) further reported that using Bile Esculin Agar, Group D Streptococci could be differentiated from non Group D Streptococci.

Bile Esculin Agar was also shown to aid differentiation of *Enterobacteriaceae, Klebsiella, Enterobacter, Serratia* from other *Enterobacteriaceae* genera (7) on the basis of esculin hydrolysis. However, other tests such as salt tolerance should be performed for identifying Enterococci (8).

Bile Esculin Azide Agar is a modification of Bile Esculin Agar as per Isenberg (9). In this medium the bile concentration is reduced and additional sodium azide is incorporated.

Tryptone, proteose peptone and HM peptone B serves as sources of carbon, nitrogen, amino acids, vitamins and essential growth nutrients. Bile and sodium azide inhibits most of the other accompyning bacteria. Esculin in the medium is hydrolyzed to esculetin and dextrose. Esculetin reacts with ferric citrate to form a dark brown or black complex, visualized as a zone of black precipitate around the colonies. If the media is dispensed in tubes in the form of slants, a positive reaction is indicated by blackening of more than half of the slant within 24-48 hours. If blackening is totally absent or if less than half of the slant is blackened within 24-48 hours, the test is negative. Viridans Streptococci sometimes exhibit a weak positive reaction. Also, *Leuconostoc, Pediococcus, Lactococcus* species causing human infections give a positive bile esculin test (10). To enhance the growth of Enterococci, Bile Esculin Agar can be supplemented with 50ml/l horse serum (3). Suspected water samples are filtered using membrane filters. These membrane filters are aseptically placed on Slanetz and Bartely Medium (M612I). Red or maroon coloured colonies observed after incubation are further confirmed by aseptically transferring the membrane filter on to Bile Esculin Azide Agar plate preheated to 44° C. Incubation at $44 \pm 0.5^{\circ}$ C for 2 hours is done following the inoculation.

Please refer disclaimer Overleaf.

M493

All typical colonies exhibiting a brown black colouration in the surrounding medium are counted as intestinal Enterococci (10). Alternatively Bile Esculin Azide Agar can also be used for direct isolation of Enterococci (without membrane filter), by incubation at 35-37°C for 18-24 hours.

Type of specimen

Clinical- Faeces, Food samples

Specimen Collection and Handling:

For food samples, follow appropriate techniques for sample collection and processing as per guidelines (4).

For clinical samples follow appropriate techniques for handling specimens as per established guidelines (8,9). After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions :

In Vitro diagnostic use. For professional use only. Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/ face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets.

Limitations :

1. Individual organisms differ in their growth requirement and may show variable growth patterns on the medium.

2. Each lot of the medium has been tested for the organisms specified on the COA. It is recommended to users to validate the medium for any specific microorganism other than mentioned in the COA based on the user's unique requirement.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the expiry period when stored at recommended temperature.

Quality Control

Appearance

Cream to yellow homogeneous free flowing powder Gelling Firm, comparable with 1.5% Agar gel Colour and Clarity of prepared medium Amber coloured, clear to slightly opalescent gel with a bluish tinge forms in Petri plates. Reaction

Reaction of 5.67% w/v aqueous solution at 25°C. pH : 7.1±0.2

pН

6.90-7.30

Cultural Response

Cultural characteristics observed after an incubation at 35-37°C for 18-24 hours.

| Organism | Inoculum (CFU) | Growth | Recovery | Esculin Hydrolysis |
|---|-------------------|-----------|----------|---|
| Enterococcus faecalis ATCC 29212 (00087*) | 50-100 | luxuriant | >=50% | positive reaction, blackening of medium around the colony |
| Escherichia coli ATCC 25922 (00013*) | >=10 ⁴ | inhibited | 0% | |
| Staphylococcus aureus subsp. aureus ATCC 25923 (00034*) | 50-100 | good | 40-50% | negative reaction |
| Proteus mirabilis ATCC 25933 | 50-100 | good | 40-50% | negative reaction |
| Streptococcus pyogenes ATCC 19615 | 50-100 | none-poor | <=10% | negative reaction |

Key: *Corresponding WDCM numbers.

Storage and Shelf Life

Store between 10-30°C in a tightly closed container and the prepared medium at 20-30°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition. Seal the container tightly after use. Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques (8,9).

Reference

- 1. Edberg S. C., Pittman S., and Singer J. M., 1977, J. Clin. Microbiol., 6:111.
- Murray P. R., Baron E. J., Jorgensen J. H., Pfaller M. A., Yolken R. H., (Eds.), 8th Ed., 2003, Manual of Clinical Microbiology, ASM, Washington, D.C.
- 3. Facklam R., 1972, Appl. Microbiol., 23:1131.
- 4. Salfinger Y., and Tortorello M.L., 2015, Compendium of Methods for the Microbiological Examination of Foods, 5th Ed., American Public Health Association, Washington, D.C.
- 5. Swan, 1954, J. Clin. Pathol., 7:160.
- 6. Facklam R. R and Moody M. D., 1970, Appl. Microbiol., 20(2):245.
- 7. Meyer and Schonfeld, 1926, Zentralbl. Bakeriol, Parasitenk. Infectionskr. Hyg. Abt. Orig. 99:402.
- 8. Isenberg, H.D. Clinical Microbiology Procedures Handbook 2nd Edition.
- 9. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock., D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.
- 10. Rochaix, 1924, Comt. Rend. Soc. Biol., 90:771.

Revision :05/2023



Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMediaTM publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMediaTM Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.



Lactobacillus MRS Agar

Intended use

For isolation and cultivation of Lactobacilli.

Composition**

| Ingredients | Gms / Litre |
|---|-------------|
| Proteose peptone | 10.000 |
| HM Peptone B # | 10.000 |
| Yeast extract | 5.000 |
| Dextrose (Glucose) | 20.000 |
| Tween 80 (Polysorbate 80) | 1.000 |
| Ammonium citrate | 2.000 |
| Sodium acetate | 5.000 |
| Magnesium sulphate | 0.100 |
| Manganese sulphate | 0.050 |
| Dipotassium hydrogen phosphate | 2.000 |
| Agar | 12.000 |
| Final pH (at 25°C) | 6.5±0.2 |
| **Formula adjusted, standardized to suit performance parameters | |

Equivalent to Beef extract

Directions

Suspend 67.15 grams in 1000 ml purified / distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C. Mix well and pour into sterile Petri plates.

Principle And Interpretation

Lactobacilli MRS media are based on the formulation of deMan, Rogosa and Sharpe (1) with slight modification. It supports luxuriant growth of all Lactobacilli from oral cavity (2), dairy products (3), foods (2), faeces (4,5) and other sources (6).

Proteose peptone and HM peptone B supply nitrogenous and carbonaceous compounds. Yeast extract provides vitamin B complex and dextrose is the fermentable carbohydrate and energy source. Polysorbate 80 supplies fatty acids required for the metabolism of Lactobacilli. Sodium acetate and ammonium citrate inhibit Streptococci, moulds and many other microorganisms. Magnesium sulphate and manganese sulphate provide essential ions for multiplication of lactobacilli. Phosphates provide good buffering action in the media.

Lactobacilli are microaerophilic and generally require layer plates for aerobic cultivation on solid media. When the medium is set, another layer of un-inoculated MRS Agar is poured over the surface to produce a layer plate. Lactobacilli isolated on MRS Agar should be further confirmed biochemically.

Type of specimen

Clinical samples - urine, faeces, etc.; Food and dairy samples

Specimen Collection and Handling:

For clinical samples follow appropriate techniques for handling specimens as per established guidelines (4,5). For food and dairy samples, follow appropriate techniques for sample collection and processing as per guidelines (2,7,8). After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions :

In Vitro diagnostic Use. For professional use only. Read the label before opening the container. Wear protective gloves/ protective clothing/eye protection/ face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets.

M641

Limitations :

1. Individual organisms differ in their growth requirement and may show variable growth patterns on the medium.

2. Further Biochemical and serological testing is required for complete identification.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the expiry period when stored at recommended temperature.

Quality Control

Appearance

Cream to light yellow homogeneous free flowing powder

Gelling

Firm, comparable with 1.2% Agar gel.

Colour and Clarity of prepared medium

Medium to dark amber coloured, clear to slightly opalescent gel forms in Petri plates

Reaction

Reaction of 6.71% w/v aqueous solution at 25°C. pH : 6.5±0.2

pН

6.30-6.70

Cultural Response

Cultural characteristics observed after an incubation at 35-37°C for 18-24 hours or longer.(with 5% CO2)

| Organism | Inoculum (CFU) | Growth | Recovery |
|--|-------------------|-----------|----------|
| <i>Lactobacillus casei</i> ATCC 9595 | 50-100 | luxuriant | >=50% |
| Lactobacillus fermentum ATCC 9338 | 50-100 | luxuriant | >=50% |
| Lactobacillus leichmannii ATCC 7830 | 50-100 | luxuriant | >=50% |
| Lactobacillus plantarum ATCC 8014 | 50-100 | luxuriant | >=50% |
| Lactobacillus saki ATCC 15521(00015*) | 50-100 | luxuriant | >=70% |
| Lactobacillus lactis ATCC 19435(00016*) | 50-100 | luxuriant | >=70% |
| Pediococcus pentosaceas ATCC 33316(00158*) | 50-100 | luxuriant | >=70% |

Key: (*) Corresponding WDCM numbers.

Storage and Shelf Life

Store dehydrated powder and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques (4,5).
References

- 1. deMan J., Rogosa M. and Sharpe M., 1960, J. Appl. Bacteriol., 23:130.
- American Public Health Association, Standard Methods for the Examination of Dairy Products, 1978, 14th Ed., Washington D.C.
- 3. Marshall R.T. (Ed.), 1992, Standard Methods for the Examination of Dairy Products, 16th ed., APHA, Washington, D.C.
- 4. Isenberg, H.D. Clinical Microbiology Procedures Handbook 2nd Edition.
- Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock., D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.
- MacFaddin J., 1985, Media for Isolation-Cultivation-Identification-Maintenance of Medical Bacteria, Vol.1, Williams and Wilkins, Baltimore.
- 7. Salfinger Y., and Tortorello M.L., 2015, Compendium of Methods for the Microbiological Examination of Foods, 5th Ed., American Public Health Association, Washington, D.C.
- Wehr H. M. and Frank J. H., 2004, Standard Methods for the Microbiological Examination of Dairy Products, 17th Ed., APHA Inc., Washington, D.C.

Revision : 05/2022



Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMediaTM publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMediaTM Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. Corporate Office : Plot No.C-40, Road No.21Y, MIDC, Wagle Industrial Area, Thane (W) - 400604, India. Customer care No.: 022-6147 1919 Email: techhelp@himedialabs.com Website: www.himedialabs.com



Technical Data

Anaerobic Blood Agar Base

Intended Use:

Recommended for cultivation of anaerobic microorganisms, including very fastidious organisms from clinical specimens.

Composition**

| Ingredients | Gms / Litre |
|---------------------|-------------|
| Tryptone | 15.000 |
| Soya peptone | 5.000 |
| Yeast extract | 5.000 |
| Sodium chloride | 5.000 |
| L-Cysteine | 0.500 |
| Hemin | 0.005 |
| Agar | 13.500 |
| Final pH (at 25°C) | 7.4±0.2 |
| | , |

**Formula adjusted, standardized to suit performance parameters

Directions

Suspend 44.0 grams in 1000 ml purified/distilled water. Heat to boiling to dissolve the medium completely. Add the rehydrated contents of one vial of Vitamin K1 supplement (FD114). Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C and aseptically add 5% v/v sterile defibrinated sheep blood. Mix well and pour into sterile Petri plates.

Principle And Interpretation

Anaerobic Blood Agar base serves as a nutritious, nonselective medium allowing the cultivation of not only fastidious anaerobes but also of aerobic and microaerophillic microorganisms (1). It promotes both typical pigment formation in *Bacteroides melaningenicus* and displays double haemolytic reaction in *Clostridium perfringens* with added blood to the medium base. The inner zone of haemolysis is due to toxin and the outer zone of incomplete haemolysis to toxin (lecithinase activity). Tryptone, soya peptone and yeast extract in the medium provides carbon and nitrogenous source, long chain amino acids, vitamins and other essential nutrients. Presence of Hemin and Vitamin K1 supports the growth of typical fastidious bacteria like *Bacteroides* species and gram positive spore bearers like *Clostridium* species. Addition of blood provides nutrients and helps to differentiate haemolytic organisms. Sodium chloride helps in maintaining the osmotic equilibrium.

Type of specimen

Clinical samples- stool, abscess

Specimen Collection and Handling

For clinical samples follow appropriate techniques for handling specimens as per established guidelines (2,3). After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions :

In Vitro diagnostic Use only. For professional use only. Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/ face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets.

Limitations :

1. Due to nutritional variations, certain strains may show poor growth.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within expiry period when stored at the recommended temperature.

M975A

Quality Control

Appearance

Yellow to tan coloured homogeneous free flowing powder

Colour and Clarity of prepared medium

 $Basal medium: Yellow coloured; with addition of 5\% v/v sterile, defibrinated sheep blood: cherry red coloured \\Basal medium: slightly opalescent; After addition of 5\% v/v sterile, defibrinated sheep blood: opaque gel in petri plates \\Basal medium: slightly opalescent; After addition of 5\% v/v sterile, defibrinated sheep blood: opaque gel in petri plates \\Basal medium: slightly opalescent; After addition of 5\% v/v sterile, defibrinated sheep blood: opaque gel in petri plates \\Basal medium: slightly opalescent; After addition of 5\% v/v sterile, defibrinated sheep blood: opaque gel in petri plates \\Basal medium: slightly opalescent; After addition of 5\% v/v sterile, defibrinated sheep blood: opaque gel in petri plates \\Basal medium: slightly opalescent; After addition of 5\% v/v sterile, defibrinated sheep blood: opaque gel in petri plates \\Basal medium: slightly opalescent; After addition of 5\% v/v sterile, defibrinated sheep blood: opaque gel in petri plates \\Basal medium: slightly opalescent; After addition of sligh$

Reaction

Reaction of 4.4% w/v aqueous solution at 25°C. pH : 7.4±0.2

pН

7.20-7.60

Cultural Response

Cultural characteristics observed after 24-48 hours at 35-37°C with 5-10% CO_2

OrganismGrowthBacteroides fragilis ATCCluxuriant25285luxuriantBacteroidesluxuriantmelaninogenicus ATCC25611Peptostreptococcusluxuriantanaerobius ATCC 27337luxuriant

Storage and Shelf Life

Store between 10-30°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition. Seal the container tightly after use. Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques (2,3).

Reference

1. Dowell, Jr., V.R., Lombard, G.L, Thompson, F.S, Armfield, A.Y.: Media for isolation, characterization and identification of obligately anaerobic bacteria- US Department of Health and Human services, centers for Disease Control (1977).

- 2. Isenberg, H.D. Clinical Microbiology Procedures Handbook 2nd Edition.
- 3. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock., D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.

Revision: 03/2022



Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMediaTM publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMediaTM Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.



Technical Data

Sabouraud Chloramphenicol Agar

Intended use

For the selective cultivation of yeasts and moulds from clinical and non-clinical samples.

Composition**

| Ingredients | Gms / Litre |
|---------------------|-------------|
| Tryptone | 5.000 |
| Peptone | 5.000 |
| Dextrose (Glucose) | 40.000 |
| Chloramphenicol | 0.050 |
| Agar | 15.000 |
| Final pH (at 25°C) | 5.6±0.2 |
| | |

**Formula adjusted, standardized to suit performance parameters

Directions

Suspend 65.05 grams in 1000 ml purified / distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C. Mix well and pour into sterile Petri plates.

Caution: Some pathogenic fungi may produce infective spores which are easily dispersed in air, so examination should be carried out in safety cabinet.

Principle And Interpretation

Sabouraud Chloramphenicol Agar is cited as Medium C and recommended for cultivation of yeasts and moulds. This medium was described originally by Sabouraud (1) for the cultivation of fungi, particularly useful for the fungi associated with skin infections. The medium is often used with antibiotics such as Chloramphenicol (2) for the isolation of pathogenic fungi from materials containing large numbers of fungi or bacteria.

Tryptone and peptone provide nitrogenous and carbonaceous compounds, long chain amino acids, and other essential growth nutrients. Dextrose provides an energy source. Chloramphenicol inhibits a wide range of Gram-positive and Gram-negative bacteria which makes the medium selective for fungi (3). The low pH favors fungal growth and inhibits contaminating bacteria from clinical specimens (4).

Type of specimen

Clinical samples - skin scrapings, nail scrapings; Food and dairy samples.

Specimen Collection and Handling:

For clinical samples follow appropriate techniques for handling specimens as per established guidelines (5,6). For food and dairy samples, follow appropriate techniques for sample collection and processing as per guidelines (7,8,9). After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions :

In Vitro diagnostic Use. For professional use only. Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/ face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets.

Limitations :

- 1. Certain pathogenic fungi may show poor growth on this medium.
- 2. Presence of chloramphenicol may inhibit certain pathogenic fungi.
- 3. Overheating of the medium may result in low productivity and softening of gel.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the expiry period when stored at recommended temperature.

Quality Control

Appearance

Cream to yellow homogeneous free flowing powder.

Gelling

Firm, comparable with 1.5% Agar gel

Colour and Clarity of prepared medium

Light amber coloured clear to slightly opalescent gel forms in Petri plates

Reaction

Reaction of 6.5% w/v aqueous solution at 25°C. pH : 5.6±0.2

pН

5.40 - 5.80

Cultural Response

Cultural characteristics observed after an incubation at 20-25°C for 48-72 hours (Incubate for 7 days for Trichophyton species).

| Organism | Inoculum (CFU) | Growth | Recovery |
|--|-------------------|----------------|----------|
| Aspergillus brasiliensis | 50-100 | good-luxuriant | |
| ATCC 16404 (00053*) | | | |
| Candida albicans ATCC 10231 (00054*) | 50-100 | good-luxuriant | >=50% |
| <i>Escherichia coli</i> ATCC 25922 (00013*) | >=10 ⁴ | inhibited | 0% |
| <i>Lactobacillus casei</i> ATCC 334 | >=10 ⁴ | inhibited | 0% |
| Saccharomyces cerevisiae ATCC 9763 (00058*) | 50-100 | good-luxuriant | >=50% |
| <i>Trichophyton rubrum</i> ATCC 28191 | 50-100 | good-luxuriant | |
| Escherichia coli NCTC 9002 | >=10 ⁴ | inhibited | 0% |
| <i>Escherichia coli</i> ATCC 8739 (00012*) | >=10 ⁴ | inhibited | 0% |
| | | | |

Key: *Corresponding WDCM numbers.

Storage and Shelf Life

Store the dehydrated powder and prepared medium on receipt between 15-25°C in a tightly closed container. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use.

Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques (5,6).

Reference

1.Sabouraud K., 1892, Ann. Dermatol. Syphilol, 3:1061.

2. Ajello L., 1957, J. Chron. Dis., 5:545.

3. Lorian (Ed.), 1980, Antibiotics In Laboratory Medicine, Williams and Wilkins, Baltimore.

4. Murray, P. R 2005, In Manual of Clinical Microbiology, 7th ed., ASM, Washington, D.C.

5. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.

6.Jorgensen, J.H., Pfaller , M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock., D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.

Please refer disclaimer Overleaf.

7.American Public Health Association, Standard Methods for the Examination of Dairy Products, 1978, 14th Ed., Washington D.C.

8.Salfinger Y., and Tortorello M.L. Fifth (Ed.), 2015, Compendium of Methods for the Microbiological Examination of Foods, American Public Health Association, Washington, D.C.

9.Wehr H. M. and Frank J. H., 2004, Standard Methods for the Microbiological Examination of Dairy Products, 17th Ed., APHA Inc., Washington, D.C.

Revision : 05/2022



HiMedia Laboratories Pvt. Limited, Plot No.C-40, Road No.21Y, MIDC, Wagle Industrial Area,



Thane (W) -400604, MS, India

CEpartner4U, Esdoornlaan 13,

3951DB Maarn, NL

www.cepartner4u.eu



IVD



Storage temperature

25°C

Do not use if package is damaged

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMediaTM publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMediaTM Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

In vitro diagnostic

medical device

HiMedia Laboratories Pvt. Ltd. Corporate Office : Plot No.C-40, Road No.21Y, MIDC, Wagle Industrial Area, Thane (W) - 400604, India. Customer care No.: 022-6147 1919 Email: techhelp@himedialabs.com Website: www.himedialabs.com



HiCrome[™] Candida Differential Agar

M1297A

Intended Use

HiCromeTM Candida Differential Agar is recommended for rapid isolation and identification of *Candida* species from mixed cultures in clinical and non-clinical samples.

Composition**

| Ingredients | Gms / Litre |
|--------------------------------|-------------|
| Peptone, special | 15.000 |
| Yeast extract | 4.000 |
| Dipotassium hydrogen phosphate | 1.000 |
| Chromogenic mixture | 7.220 |
| Chloramphenicol | 0.500 |
| Agar | 15.000 |
| Final pH (at 25°C) | 6.3±0.2 |
| | |

**Formula adjusted, standardized to suit performance parameters

Directions

Suspend 42.72 grams in 1000 ml purified/distilled water. Heat to boiling to dissolve the medium completely. DO NOT AUTOCLAVE. Cool to 45-50°C. Mix well and pour into sterile Petri plates.

Principle And Interpretation

Perry and Miller (3) reported that Candida albicans produces an enzyme b -N-acetyl- galactosaminidase and according to Rousselle et al (4) incorporation of chromogenic or fluorogenic hexosaminidase substrates into the growth medium helps in identification of *C.albicans* isolates directly on primary isolation. HiCromeTM Candida Differential Agar is a selective and differential medium, which facilitates rapid isolation of yeasts from mixed cultures and allows differentiation of Candida species namely C.albicans, C.krusei, C.tropicalis and C.glabrata on the basis of colouration and colony morphology. On this medium results are obtained within 48 hours and it is useful for the rapid and presumptive identification of common yeasts in Mycology and Clinical Microbiology Laboratory. Peptone special and yeast extract provides nitrogenous, carbonaceous compounds and other essential growth nutrients. Phosphate buffers the medium well. Chloramphenicol suppresses the accompanying bacterial flora. C.albicans appear as light green coloured smooth colonies, C.tropicalis appear as blue to metallic blue coloured raised colonies. C.glabrata colonies appear as cream to white smooth colonies, while C.krusei appear as purple fuzzy colonies.

Type of specimen

Clinical samples - skin scrapings, urine.

Specimen Collection and Handling

For clinical samples follow appropriate techniques for handling specimens as per established guidelines (3,4). After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions

In Vitro diagnostic Use only. Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/ face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets.

Limitations

1. Variations in colour intensity may be observed for Candida isolates depending on the presence of enzymes.

2. Other *Candida* species may produce light mauve coloured colonies which is also produced by other yeast cells. This must be confirmed by further biochemical tests.

3. Other filamentous fungi also exhibit colour on this medium.

Performance And Evaluation

Performance of the medium is expected when used as per the direction on the label within the expiry period when stored at recommended temperature.

Quality Control

Appearance

Cream to beige homogeneous free flowing powder

Gelling

Firm, comparable with 1.5% Agar gel

Colour and Clarity of prepared medium

Light amber coloured, clear to slightly opalescent gel forms in Petri plates

Reaction

Reaction of 4.27% w/v aqueous solution at 25°C. pH : 6.3±0.2

pН

6.10-6.50

Cultural Response

Cultural characteristics observed after an incubation at 30-35°C for 40-48 hours.

| Organism | Inoculum (CFU) | Growth | Recovery | Colour of Colony |
|--|-------------------|----------------|----------|--|
| Candida albicans ATCC 10231 (00054*) | 50-100 | good-luxuriant | >=50% | light green |
| Candida glabrata ATCC 15126 | 50-100 | good-luxuriant | >=50% | cream to white |
| #Teunomyces krusei ATCC 24408 | 50-100 | good-luxuriant | >=50% | purple, fuzzy |
| Candida tropicalis ATCC 750 | 50-100 | good-luxuriant | >=50% | blue to purple |
| Candida kefyr ATCC 66058 | 50-100 | good-luxuriant | >=50% | cream to white with slight purple centre |
| <i>Candida utilis</i> ATCC 9950 | 50-100 | good-luxuriant | >=50% | pale pink to pinkish purple |
| Candida parapsilosis ATCC 22019 | 50-100 | good-luxuriant | >=50% | white to cream |
| Candida membranifaciens ATCC 20137 | 50-100 | good-luxuriant | >=50% | white to cream |
| <i>Candida dubliensis</i> NCPF 3949 | 50-100 | good-luxuriant | >=50% | pale green |
| Escherichia coli ATCC 25922 (00013*) | >=104 | inhibited | 0% | |
| Staphylococcus aureus subsp.aureus ATCC 25923 (00034*) | >=10 ⁴ | inhibited | 0% | |

Key : *Corresponding WDCM numbers. # - Formerly known as Candida krusei

Storage and Shelf Life

Store between 15-25°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques (1,2).

References

1. Isenberg, H.D. Clinical Microbiology Procedures Handbook 2nd Edition.

- 2. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock., D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.
- 3. Perry J. L. and Miller G. R., 1987, J. Clin. Microbiol., 25: 2424 -2425.
- 4. Rousselle P., Freydiere A., Couillerot P., de Montclos H. and GilleY., 1994, J. Clin. Microbiol. 32:3034-3036.

Revision : 05/ 2020

| IVD | In vitro diagnostic medical device |
|--------|--|
| (€ | CE Marking |
| 15°C | Storage temperature |
| | Do not use if package is damaged |
| | HiMedia Laboratories Pvt. Limited, 23 Vadhani Industrial Estate, LBS Marg,Mumbai-86,MS,India |
| EC REP | CE Partner 4U ,Esdoornlaan 13, 3951 DB Maarn The Netherlands, <u>www.cepartner</u> 4u.eu |

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia[™] publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia[™] Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. Reg.office : 23, Vadhani Ind.Est., LBS Marg, Mumbai-400086, India. Customer care No.: 022-6116 9797 Corporate office : A-516,Swastik Disha Business Park,Via Vadhani Ind. Est., LBS Marg, Mumbai-400086, India. Customer care No.: 022-6147 1919 Email: techhelp@himedialabs.com Website: www.himedialabs.com



Technical Data

НіСготе™ UTI Agar

Intended use

Recommended for presumptive identification and confirmation of microorganisms mainly causing urinary tract infections, can also be used for testing water, food, environmental and other clinical samples.

| Composition** | |
|---------------------|---------------|
| Ingredients | Gms / Litre |
| Peptone, special | 15.000 |
| Chromogenic mixture | 2.450 |
| Agar | 15.000 |
| Final pH (at 25°C) | $6.8{\pm}0.2$ |

**Formula adjusted, standardized to suit performance parameters

Directions

Suspend 32.45 grams in 1000 ml purified /distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C. Mix well and pour into sterile Petri plates.

Principle And Interpretation

Urinary tract infections are bacterial infections affecting parts of urinary tract. The common symptoms of urinary tract infection are urgency and frequency of micturition, with associated discomfort or pain. The common condition is cystitis, due to infection of the bladder with a uropathogenic bacterium, which most frequently is *Escherichia coli*, but sometimes *Staphylococcus saprophyticus* or especially in hospital-acquired infections, *Klebsiella* species, *Proteus mirabilis*, other coliforms, *Pseudomonas aeruginosa* or *Enterococcus faecalis* (1). HiCromeTM UTI Agar is formulated on basis of work carried out by Pezzlo (2) Wilkie et al (3), Friedman et al (4), Murray et al (5), Soriano and Ponte (6) and Merlino et al (7). These media are recommended for the detection of urinary tract pathogens where HiCromeTM UTI Agar has broader application as a general nutrient agar for isolation of various microorganisms. It facilitates and expedites the identification of some gram-negative bacteria and some gram-positive bacteria on the basis of different contrasted colony colours produced by reactions of genus or species specific enzymes with two chromogenic substrates. The chromogenic substrates are specifically cleaved by enzymes produced by *Enterococcus* species, *E.coli* and coliforms. Presence of amino acids like phenylalanine and tryptophan from peptones helps for detection of tryptophan deaminase activity, indicating the presence of *Proteus* species, *Morganella* species and *Providencia* species.

One of the chromogenic substrate is cleaved by β -glucosidase possessed by Enterococci resulting in formation of blue colonies. *E.coli* produce pink colonies due to the enzyme β -D-galactosidase that cleaves the other chromogenic substrate. Further confirmation of *E.coli* can be done by performing the indole test. Coliforms produce purple coloured colonies due to cleavage of both the chromogenic substrate. Colonies of *Proteus, Morganella* and *Providencia* species appear brown because of tryptophan deaminase activity. Peptone special provides nitrogenous, carbonaceous compounds, long chain amino acids, vitamins and other essential growth nutrients. This medium can be made selective by supplementation with antibiotics for detecting microorganisms associated with hospital borne infections.

Type of specimen

Clinical samples : urine, faeces, Food samples, Water samples.

Specimen Collection and Handling

For clinical samples follow appropriate techniques for handling specimens as per established guidelines (8,9). For food and dairy samples, follow appropriate techniques for sample collection and processing as per guidelines (10,11). For water samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards (12). After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions

In Vitro diagnostic use. For professional use only. Read the label before opening the container. Wear protective gloves/ protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets.

Limitations

Individual organisms differ in their growth requirement and may show variable growth patterns on the medium.
 Each lot of the medium has been tested for the organisms specified on the COA. It is recommended to users to validate the medium for any specific microorganism other than mentioned in the COA based on the user's unique requirement.
 Since it is an enzyme-substrate based reaction, the intensity of colour may vary with isolates.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the expiry period when stored at recommended temperature.

Quality Control

Appearance

Cream to yellow homogeneous free flowing powder

Gelling

Firm, comparable with 1.5% Agar gel

Colour and Clarity of prepared medium

Light amber coloured, clear to slightly opalescent gel forms in Petri plates

Reaction

Reaction of 3.24% w/v aqueous solution at 25°C. pH : 6.8±0.2

pН

6.60-7.20

Cultural Response

Cultural characteristics observed after an incubation at 35-37°C for 16-24 hours.

| Organism | Inoculum (CFU) | Growth | Recovery | Colour of Colony |
|---|-------------------|-----------|----------|--|
| <i>Escherichia coli</i> ATCC 25922 (00013*) | 50-100 | luxuriant | >=70% | Purple to magenta |
| Enterococcus faecalis ATCC 29212 (00087*) | 50-100 | luxuriant | >=70% | blue-green, (small) |
| Klebsiella pneumoniae ATCC 13883 (00097*) | 50-100 | luxuriant | >=70% | blue to purple, mucoid |
| Proteus mirabilis ATCC 12453 | 50-100 | luxuriant | >=70% | light brown |
| Pseudomonas aeruginosa ATCC 27853 (00025*) | 50-100 | luxuriant | >=70% | colourless (greenish pigment may be observed) |
| Staphylococcus aureus subsp. aureus ATCC 25923 (00034*) | 50-100 | luxuriant | >=70% | cream to golden yellow |

Key: *Corresponding WDCM numbers.

Storage and Shelf Life

Store between 15-25°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques (8,9).

Reference

- Collee J. G., Fraser A. G., Marmion B. P., Simmons A., (Eds.), Mackie and McCartney, Practical Medical Microbiology, 1996, 14th Edition, Churchill Livingstone.
- 2. Pezzlo M., 1998, Clin. Microbiol. Rev., 1:268-280.
- 3. Wilkie M. E., Almond M. K., Marsh F. P., 1992, British Medical Journal 305:1137-1141.
- 4. Friedman M. P. et al, 1991, J. Clin. Microbiol., 29:2385-2389.
- 5. Murray P., Traynor P. Hopson D., 1992, J. Clin. Microbiol., 30:1600-1601.
- 6. Soriano F., Ponte C., 1992, J. Clin. Microbiol., 30:3033-3034.
- 7. Merlino et al, 1995, Abstr. Austr. Microbiol. 16(4):17-3.
- 8. Isenberg, H.D. Clinical Microbiology Procedures Handbook 2nd Edition.
- 9. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock., D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.
- 10.Salfinger Y., and Tortorello M.L., 2015, Compendium of Methods for the Microbiological Examination of Foods, 5th Ed., American Public Health Association, Washington, D.C.
- 11. Wehr H. M. and Frank J. H., 2004, Standard Methods for the Microbiological Examination of Dairy Products, 17th Ed., APHA Inc., Washington, D.C.
- 12. Lipps WC, Braun-Howland EB, Baxter TE, eds. Standard methods for the Examination of Water and Wastewater, 24th ed. Washington DC:APHA Press; 2023.

Revision : 06/2023



Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMediaTM publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMediaTM Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. Corporate Office : Plot No.C-40, Road No.21Y, MIDC, Wagle Industrial Area, Thane (W) - 400604, India. Customer care No.: 022-6147 1919 Email: techhelp@himedialabs.com Website: www.himedialabs.com



Technical Data

Bifidobacterium Broth

M1395

Intended Use:

Recommended for cultivation of Bifidobacterium infantis.

Composition**

| Ingredients | Gms / Litre |
|---------------------------|-------------|
| Tryptone | 20.000 |
| Peptone | 10.000 |
| Yeast extract | 10.000 |
| Tomato juice, solids | 16.650 |
| Dextrose (Glucose) | 20.000 |
| Polysorbate 80 (Tween 80) | 2.000 |
| Final pH (at 25°C) | 6.8±0.2 |

**Formula adjusted, standardized to suit performance parameters

Directions

Suspend 78.65 grams in 1000 ml purified / distilled water . Heat if necessary to dissolve the medium completely. Distribute in tubes or flasks as desired. Sterilize by autoclaving at 15lbs pressure for 15 minutes.

Principle And Interpretation

The genus *Bifidobacterium* is the third most numerous bacterial population found in the human intestine after Bacteroides and *Eubacterium*. It is an anaerobic bacteria that makes up the gut microbial flora, it resides in the colon and have health benefits for their hosts. Bifidobacteria are also associated with lower incidences of allergies (1,2). Bifidobacterium Broth is used for the cultivation and maintenance of *Bifidobacterium* species. The medium is used exclusively for the cultivation of *Bifidobacterium infantis* (3).

Tryptone, Peptone and yeast extract provides essential growth nutrients. Glucose is the energy source and sodium chloride maintains isotonic conditions. Tomato juice helps in maintaining acidic pH while polysorbate 80 provides fatty acids required for metabolic activity of *Bifidobacterium*.

Type of specimen

Clinical samples- faeces; Dairy samples.

Specimen Collection and Handling:

For clinical samples follow appropriate techniques for handling specimens as per established guidelines (4,5). For dairy samples, follow appropriate techniques for sample collection and processing as per guidelines (6,7). After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions

In Vitro diagnostic Use. For professional use only. Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/ face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets.

Limitations

- 1. Further biochemical and serological tests must be carried out for complete identification.
- 2. Bifidobacterium species are strict anaerobes, hence condition must be appropriately maintained.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the expiry period when stored at recommended temperature.

Quality Control

Appearance

Cream to yellow homogeneous free flowing powder Colour and Clarity of prepared medium Amber coloured clear solution in tubes

Reaction

Reaction of 7.86% w/v solution at 25°C. pH : 6.8±0.2 **pH** 6.60-7.00 **Cultural Response** Cultural characteristics observed after an incubation at 35-37°C for 24-48 hours.

| Organism | Inoculum (CFU) | Growth |
|---|-------------------|----------------|
| <i>Bifidobacterium infantis</i> ATCC 25962 | 50-100 | good-luxuriant |
| <i>Bifidobacterium bifidum</i> ATCC 15696 | 50-100 | good-luxuriant |
| <i>Bifidobacterium breve</i> ATCC 15698 | 50-100 | good-luxuriant |

Storage and Shelf Life

Store between 10-30°C in a tightly closed container and the prepared medium at 15-25°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition. Seal the container tightly after use. Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques (4,5).

Reference

1.Bjorksten B., Sepp E., Julge K., Voor T., and Mikelsaar M., 2001, J. Allergy Clin. Microbiol., Volume 108, Issue 4, 516-520.

2.Guarner F., and Malagelada J. R., 2003, The Lancet, Vol. 361, Issue 9356, 8 February 2003, 512-519

3. Atlas R. M. 2004, 3rd Edi. Handbook of Microbiological Media, Parks, L. C. (Ed.), CRC Press, Boca Raton.4.

4. Isenberg, H.D. Clinical Microbiology Procedures Handbook 2nd Edition.

5.Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock., D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.

6.American Public Health Association, Standard Methods for the Examination of Dairy Products, 1978, 14th Ed., WashingtonD.C.

7.Wehr H. M. and Frank J. H., 2004, Standard Methods for the Microbiological Examination of Dairy Products, 17th Ed., APHA Inc., Washington, D.C.

Revision : 03/2022



Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMediaTM publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMediaTM Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.



Technical Data

Hottinger Broth

M1425

Intended Use:

Recommended for cultivation of less fastidious microorganisms and determination of indole in accordance with USSR State Pharmacopoeia.

Composition**

| Ingredients | Gms / Litre |
|---|---------------|
| Fish peptone | 20.000 |
| Yeast extract | 2.000 |
| Tryptophan | 1.000 |
| Final pH (at 25°C) | $7.4{\pm}0.2$ |
| **Formula adjusted, standardized to suit performance parameters | |

Directions

Suspend 23.0 grams in 1000 ml purified/distilled water. Heat if necessary to dissolve the medium completely. Dispense into tubes or flasks as desired. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes.

Principle And Interpretation

Hottinger Broth is used for cultivation of less fastidious microorganisms and determination of indole as per USSR State Pharmacopeia (1).

Fish peptone and yeast extract provides the nitrogenous source and essential nutrients for growth of organisms. The production of indole from tryptophan is a diagnostic test used for identifying enteric bacteria. After incubation, indole can be identified by a red dye complex reaction with one of several reagents eg. Kovac's Reagent which consists of amyl alcohol, dimethylaminobenzaldehyde and concentrated hydrochloric acid (2).

Type of specimen

Pharmaceutical samples

Specimen Collection and Handling:

For pharmaceutical samples, follow appropriate techniques for sample collection and processing as per guidelines (1). After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions :

Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/ face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling specimens. Safety guidelines may be referred in individual safety data sheets.

Limitations :

1. Individual organisms differ in their growth requirement and may show variable growth patterns on the medium.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the expiry period when stored at recommended temperature.

Quality Control

Appearance

Cream to yellow homogeneous free flowing powder **Colour and Clarity of prepared medium** Light amber coloured clear solution **Reaction** Reaction of 2.3% w/v aqueous solution at 25°C. pH : 7.4±0.2 **pH** 7.20-7.60

Cultural Response

Cultural characteristics observed after an incubation at 35-37°C for 18-48 hours.

| Organism | Growth | Indole production |
|---|--------|---|
| Escherichia coli ATCC 25922 (00013*) | good | Positive reaction,red ring at the interface of the medium |
| Pseudomonas aeruginosa ATCC 27853 (00025*) | good | Negative reaction,no colour development/ cloudy ring |
| Staphylococcus aureus subsp. aureus ATCC 25923 (00034*) | good | Negative reaction,no colour development/ cloudy ring |
| Streptococcus pyogenes ATCC 19615 | good | Negative reaction,no colour development/ cloudy ring |

Key : *Corresponding WDCM numbers.

Storage and Shelf Life

Store between 10-30°C in a tightly closed container and the prepared medium at 15-25°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition. Seal the container tightly after use. Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with sample must be decontaminated and disposed of in accordance with current laboratory techniques (3,4).

Reference

1. State Pharmacopoeia of USSR.

2. Harrigar W.F and McCarran M.E (1966) Laboratory Methods in Microbiology Academic Press 53.

3. Isenberg, H.D. Clinical Microbiology Procedures Handbook 2nd Edition.

4. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock., D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.

Revision: 05/2023

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMediaTM publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMediaTM Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. Corporate Office : C-40, Road No.21Y, MIDC, Wagle Industrial Area, Thane (W) - 400604, India. Customer care No.: 022-6147 1919 Email: techhelp@himedialabs.com Website: www.himedialabs.com



Technical Data

Urea Indole Medium

M1784

To differentiate micro-organisms especially *Enterobacteriaceae* on the basis of their ability to hydrolyze urea and indole production

Composition**

| Ingredients | Gms / Litre |
|--------------------------------|-------------|
| L- Tryptophan | 3.000 |
| Sodium chloride | 5.000 |
| Potassium phosphate, monobasic | 1.000 |
| Potassium phosphate, dibasic | 1.000 |
| Urea | 20.000 |
| Phenol red | 0.012 |
| Final pH (at 25°C) | 6.8±0.2 |
| | |

**Formula adjusted, standardized to suit performance parameters

Directions

Suspend 30.01 grams in 1000 ml distilled water. Dissolve the medium completely and sterilize by filtration . DO NOT AUTOCLAVE. Aseptically, dispense into sterile tubes.

Principle And Interpretation

Strains of *Enterobacteria* are associated with abscesses, pneumonia, meningitis, septicemia and infections of wounds, the urinary tract and the intestine. They are a major component of the normal intestinal flora of humans but are relatively uncommon at other body sites. Of clinically significant isolates, *Enterobacteriaceae* may account for 80% of gram-negative bacilli and 50% of all clinically significant isolates in clinical microbiology laboratories (1).

Urea Indole Medium is used for the identification of *Enterobacteria* on the basis of Urease and indole production and the transdeamination of tryptophan. This medium is very useful in the identification of *Proteus* species from *Salmonella* and *Shigella* species. The results for urease production should be noted prior to indole reaction, as addition of Kovacs reagent, decolourizes the medium, due to drop in pH.

L- Trypytophan is an essential amino acid and is converted to skatole and indole, which is detected by the addition of Kovacs Reagent (R008). Sodium chloride maintains the osmotic balance. The phosphates helps in the buffering of the medium. Microorganisms that possess the enzyme urease hydrolyse urea, releasing ammonia, which is detected by the pH indicator phenol red. The alkalinility so developed imparts pink colour to the medium (2).

Quality Control

AppearanceLight yellow to light pink homogeneous free flowing powderColour and Clarity of prepared mediumYellow to light orange coloured clear solutionReactionReactionReaction of 3.00% w/v aqueous solution at 25°C. pH : 6.8 ± 0.2 pH6.60-7.00Cultural ResponseCultural characteristics observed after an incubation at 35-37°C for 18-24 hours.Cultural ResponseOrganismInoculumGrowthUrease

(CFU)

| Escherichia coli ATCC 25922 | 50-100 | luxuriant | Negative reaction,no |
|--------------------------------------|--------|-----------|--------------------------------------|
| Proteus mirabilis ATCC 12453 | 50-100 | luxuriant | change Positive reaction, Pink |
| Proteus vulgaris ATCC 13315 | 50-100 | luxuriant | Positive reaction, Pink |
| Salmonella Typhimurium ATCC 14028 | 50-100 | luxuriant | Negative reaction,no change |

Storage and Shelf Life

Store below 8°C in tightly closed container and prepared medium at 2-8°C. Use before expiry period on the label

Reference

1.Patrick R. Murray et al, Manual of Clinical Microbiology, Sixth Edition, 444 - 445. 2.Roland F. Bourbon D, Sztrum S. Ann. Inst. Pasteur, 73. 914-916.

Revision : 1 / 2011

CE

Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia[™] publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia[™] Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. A-516, Swastik Disha Business Park, Via Vadhani Ind. Est., LBS Marg, Mumbai-400086, India. Customer care No.: 022-6147 1919 Email: techhelp@himedialabs.com



Technical Data

Mueller Hinton Agar, 2% Glucose with Methylene blue

M1825

Mueller Hinton Agar, 2% Glucose with Methylene blue is recommended for testing performing Antifungal Disk Diffusion Susceptibility of yeasts.

Composition**

| Ingredients | Gms / Litre |
|-------------------------|-------------|
| Beef infusion from | 300.000 |
| Casein Acid Hydrolysate | 17.500 |
| Starch | 1.500 |
| Glucose | 20.000 |
| Methylene blue | 0.0005 |
| Agar | 17.000 |
| Final pH (at 25°C) | 7.3±0.1 |
| | |

**Formula adjusted, standardized to suit performance parameters

Directions

Suspend 58 grams in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15lbs pressure (121°C) for 15 minutes. Mix well before pouring.

The performance of this batch has been tested and standardized as per the current CLSI (formerly NCCLS) document M44-A2 in Method for Antifungal Disk Diffusion susceptibility Testing of yeasts

Principle And Interpretation

The Mueller Hinton formulation was originally developed as a simple, transparent agar medium for the cultivation of pathogenic species (1). Mueller Hinton Agar, Modified (as per CLSI for antifungal) is recommended for the diffusion of antifungal agents impregnated on paper disc through an agar gel as described in CLSI Approved Standard (2).

When supplemented with glucose to a final concentration of 2%, it provides for suitable fungal growth. The addition of methylene blue to a final concentration of $5\mu g/ml$ enhances zone edge definition.

Kirby-Bauer et al recommended Mueller Hinton Agar for performing antibiotic susceptibility tests using a single disc of high concentration (4). WHO Committee on Standardization of Susceptibility Testing has accepted Mueller Hinton Agar for determining the susceptibility of microorganisms because of its reproducibility (3). Mueller Hinton Agar with 5% sheep blood and Mueller Hinton Agar with Haemoglobin have been recommended for antimicrobial susceptibility testing of *Streptococcus pneumoniae* and *Haemophilus influenzae*. Similarly Mueller Hinton Agar, Modified (as per CLSI for antifungal) is recommended for antifungal susceptibility testing of discs.

Beef infusion and casein acid hydrolysate provide nitrogenous compounds, carbon, sulphur and other essential nutrients. Starch acts as a protective colloid against toxic substances present in the medium. Starch hydrolysis yields dextrose, which serves as a source of energy. Glucose serves as an energy source for fungal cultures while Methylene blue enhances zone edge definition.

Technique:

Preparation of Inoculum:

1. Inoculum is prepared by picking five distinct colonies of approximately 1mm from 24 hours old culture grown on Sabouraud Dextrose Agar (M063) and incubated at $35 \pm 2^{\circ}$ C. Colonies are suspended in 5ml of sterile 0.85% Saline.

2. Vortex the resulting suspension and adjust the turbidity to yield 1 x 106 - 5 x 106 cells /ml (i.e. 0.5 McFarland standard).

Test Procedure:

1. Prepare plates with Mueller Hinton Agar, Modified (as per CLSI for antifungal) for carrying out susceptibility of antifungal discs. The medium in the plates should be sterile and have a depth of about 4 mm.

2. Dip a sterile non-toxic cotton swab on a wooden applicator into the standardized inoculum (turbidity so adjusted, as to obtain semi confluent growth on the Petri plate) and rotate the soaked swab firmly against the upper inside wall of the tube to express excess fluid. Streak the entire agar surface of the plate with the swab three times, turning the plate at 60° angle between each streaking. Allow the inoculum to dry for 5 - 15 minutes with lid in place.

3. Apply the discs using aseptic technique. Deposit the discs with centers at least 24 mm apart. (Not more than 12 discs should be placed on a 150-mm plate or not more than 5 discs on a 100-mm plate

4. Invert the plates and place in an incubator set to $35 \pm 2^{\circ}$ C within 15 minutes after the discs are applied.

5. Examine each plate after 20 - 24 hours of incubation. If plate was satisfactorily streaked the resulting zones of inhibition will be uniformly circular and there will be a semi-confluent lawn of growth. Read at 48 hours only when insufficient growth is observed after 24 hours incubation.

Quality Control

Appearance

Light yellow to yellow may have slight blue tinge homogeneous free flowing powder

Gelling

Firm, comparable with 1.7% agar gel.

Colour and Clarity of prepared medium

amber coloured clear to slightly opalescent gel froms in Petri plates

Reaction

Reaction of 5.8% w/v aqueous solution at 25°C. pH : 7.3±0.1

pН

7.20-7.40

Cultural response

A luxuriant growth of test organisms was observed on Mueller Hinton Agar, Modified (as per CLSI for antifungal) in 24-48 hours at 33-37°C along with inhibition zones with respective antibiotic concentrations

Cultural Response

| Organism | Inoculum (CFU) | Growth | Recovery | Amphotericin B AP(100units) | - Amphotericin B AP(20 mcg) | Amphotericin- B AP(50 mcg) |
|--------------------------------------|-------------------|-----------|----------|-----------------------------------|--------------------------------|-------------------------------|
| Cultural response | | | | (| | |
| Candida albicans ATCC 90028 | 50-100 | Luxuriant | >=70% | 10 -17 mm | 10 -15 mm | 31 -42 mm |
| Candida parapsilosis ATCC 22019 | 50-100 | luxuriant | >=70% | 11 -20 mm | 10 -17 mm | 28 -37 mm |
| Candida tropicalis ATCC 750 | 50-100 | luxuriant | >=70% | 8 -12 mm | 8 -10 mm | 13 -17 mm |
| Candida krusei ATCC 6258 | | luxuriant | >=70% | 9 -14 mm | 8 -12 mm | 16 -25 mm |
| Candida albicans ATCC 10231 | 50-100 | luxuriant | >=70% | 10 -18 mm | 10 -16 mm | 30 -40 mm |
| Saccharomyces cerevisiae ATCC9763 | 50-100 | luxuriant | >=70% | 11 -18 mm | 8 -12 mm | 29 -38 mm |

Storage and Shelf Life

Store dehydrated powder below 30°C and prepared medium at 2-8°C. Use before expiry date on the label.

Reference

1.Mueller J. H. and Hinton J., 1941, Proc. Soc. Exp. Biol. Med., 48:330.

2.Method for Antifungal Disk Diffusion Susceptibility Testing of yeasts; Approved Guideline Second Edition M44-A2 Vol.24 No.17.

3. Present Status and Future Work, WHO Sponsored collaborative study, Chicago, Oct. 1967.

4.Bauer A. W., Kirby W. M., Sherris J. L. and Turck M., 1966, Am. J. Clin. Pathol., 45:493.

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia[™] publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia[™] Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. A-516, Swastik Disha Business Park, Via Vadhani Ind. Est., LBS Marg, Mumbai-400086, India. Customer care No.: 022-6147 1919 Email: techhelp@himedialabs.com Website: www.himedialabs.com

CE



McFarland Standard set

R092

McFarland standards are used to perform spectrophotometric comparisions of bacterial densities in water, saline or liquid growth medium. It provides laboratory guidance for the standardization of numbers of bacteria for susceptibility testing or other procedure requiring a standardization of the inoculum like growth promotion test (GPT).

Set Contains:

R092A (Standard 0.5)-1 tube

R092B (Standard 1)-1 tube

R092C (Standard 2)- 1 tube

R092D (Standard 3)-1 tube

R092E (Standard 4)- 1 tube

Directions

Prepare the inoculum of culture required for testing by using sterile saline. Match the density of the resultant suspension with the density of the desired standard. The standards must be thoroughly mixed on a vortex mixture at the time of use to obtain a uniform suspension. Adjust the density of cell suspension by adding saline if it is more turbid as compared to the desired standard or by adding culture if it is dilute. Check the density of the turbidity by determining the absorbance of 0.5 McFarland standard using a spectrophotometer with a 1 cm light path. The absorbance at 625 nm should be 0.08 to 0.10. The standards should be checked regularly to ensure the density accuracy.

Interpretation

McFarland standards are a set of tubes with increasing concentration of Barium Sulphate suspension. The turbidity of Barium Sulphate's white precipitation is used as a point of comparision of bacterial suspensions to known bacterial turbidity.

| McFarland | 0.5 | 1 | 2 | 3 | 4 |
|------------------------|-----|---|---|---|----|
| Standard | | | | | |
| Approximate | 1.5 | 3 | 6 | 9 | 12 |
| Corresponding | | | | | |
| suspension x | | | | | |
| 10 ⁸ CFU/ml | | | | | |

Limitation of procedure

1. Coloured media may interfere with result interpretation and give incorrect results.

2. Bacterial suspensions of older cultures may not be comparable with expected bacterial counts.

Storage

Store the standards at 2-8°C, away from light after each use.

Reference

- 1. McFarland, J.1907. Nephelometer: JAMA 14:1176-1178
- 2. Murry, PR; Baron, EJ; Jorgensen, JH; Landry, ML; Pfaller, MA; Manual of Clinical

Microbiology 9th edition ASM press, Washington DC.

Disclaimer :

User must ensure of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia[™] Publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. Himedia[™] Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research orfurther manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.



Technical Data

Gram's Iodine

S013

Intended use

Grams's Iodine is used as mordant in Gram's staining method.

Composition**

Ingredients

| Iodine | 1.0 ml |
|------------------|----------|
| Potassium iodide | 2.0 ml |
| Distilled water | 300.0 ml |

**Formula adjusted, standardized to suit performance parameters

Directions

- 1. Prepare a thin smear on clear, dry glass slide.
- 2. Allow it to air dry and fix by gentle heat.

3. Flood with Gram's Crystal Violet (S012) for 1 minute. (If over staining results in improper decolourization of known gram-negative organisms, use less crystal violet).

4. Wash with tap water.

5. Flood the smear with Gram's Iodine (S013). Allow it to remain for 1 minute.

6. Decolourize with Gram's Decolourizer (S032) until the blue dye no longer flows from the smear. (Acetone may be used as a decolourizing agent with caution, since this solvent very rapidly decolourized the smear).

7. Wash with tap water.

8. Counter stain with 0.5% w/v Safranin (S027) for 20 seconds and rinses off with water.

9. Wash with tap water.

10. Allow the slide to air dry or blot dry between sheets of clean bibulous paper and examine under oil immersion objective.

Principle And Interpretation

The Gram stain is a differential staining technique most widely applied in all microbiology disciplines laboratories. It is one of the most important criteria in any identification scheme for all types of bacterial isolates. Different mechanisms have been proposed to explain the gram reaction. There are many physiological differences between gram-positive and gram negative cell walls Ever since Christian Gram has discovered Gram staining, this process has been extensively investigated and redefined In practice, a thin smear of bacterial cells is stained with crystal violet, then treated with an iodine containing mordant to increase the binding of primary stain A decolourizing solution of alcohol or acetone is used to remove the crystal violet from cells which bind it weakly and then the counterstain (like safranin) is used to provide a colour contrast in those cells that are decolourized. Gram-positive bacteria have a thick mesh-like cell wall made of peptidoglycan (50–90% of cell envelope), and as a result are stained purple by crystal violet, whereas gram-negative bacteria have a thinner layer (10% of cell envelope), so do not retain the purple stain and are counter-stained pink by safranin. In a properly stained smear by gram staining procedure, the gram-positive bacteria appear blue to purple and gram negative cells appear pink to red.

Type of specimen

Any isolated colony on primary or subculture plates can be isolated from following specimens. Clinical specimen: Blood, urine, CSF, pus, wounds, lesions, body tissues, sputum etc. From environment: Air, water, soil, sludge, waste water, food, dairy samples etc.

Specimen Collection and Handling

For clinical samples follow appropriate techniques for handling specimens as per established guidelines

For food and dairy samples, follow appropriate techniques for sample collection and processing as per guidelines For water samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards.

Generally, the smear is made in laboratory; however, when there is a concern that transport will be delayed or that the preservation for culture will alter the specimen, prepare smear and submit slides to the laboratory.

Warning and Precautions

In Vitro diagnostic use only. For professional use only. Read the label before opening the container. Wear protective gloves/protective clothing/ eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets.

Limitations

1. Use results of Gram stains in conjunction with other clinical and laboratory findings. Use additional procedures (e.g., special stains, inclusion of selective media, etc.) to confirm findings suggested by gram-stained smears.

2. Proper smear preparation is key to obtaining good gram staining results. Avoid excessive material or thick smears which may interfere with the passage of light and lead to distortion of images.

3. Overheating slides during heat fixation can distort the appearance of the organisms.

4. Only fresh cultures and specimens should be gram stained since cell wall integrity of older cells may give improper gram staining characteristics. Gram positive organisms that have lost cell wall integrity because of old age or antibiotic treatment may appear pink.

5. The decolorization step is the most important step in the gram-staining process. Over decolorization results in an abundance of bacteria that appear gram negative, while under decolorization results in too many bacteria that appear to be gram-positive.

6. The procedure given is based on an ideal thin smear of cells. Staining and decolorization times may vary depending on the sample and its thickness.

7. False Gram stain results may be related to inadequately collected specimens or delay in transit.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the expiry period when stored at recommended temperature

Quality Control

- → **Appearance :** Yellow to dark brown coloured solution.
- \rightarrow **Clarity :** Clear without any particles.

- → **Microscopic Examination :** Gram staining is carried out where Gram's Iodine is used as one of the stains and staining characteristics of organisms are observed under microscope by using oil immersion lens.
- → Results : Gram-positive microorganisms : violet Gram-negative microorganims : pinkish red

Storage and Shelf Life

Store between 10 - 30 °C in tightly closed container and away from bright light. Use before expiry date on label. On opening, product should be properly stored in dry ventilated area protected from extremes of temperature and sources of ignition. Seal the container tightly after use.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques .

Reference

1. Downes F. P. and Ito K. (Ed.), 2001, Compendium of Methods for the Microbiological Examination of Foods, 4th ed., APHA, Washington, D.C.

2. Rice E.W., Baird, R.B., Eaton A. D., Clesceri L. S. (Eds.), 2012, Standard Methods for the Examination of Water and Wastewater, 22nd ed., APHA, Washington, D.C.

3. Wehr H. M. and Frank J. H., 2004, Standard Methods for the Microbiological Examination of Dairy Products, 17th Ed., APHA Inc., Washington, D.C.

4. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.

5. Jorgensen, J.H., Pfaller , M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock., D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.

6. Shanhooltzer, C. J., P. Schaper ,and L.R. Peterson 1982 .Concentrated Gram stain smear prepared with a cytospin centrifuge. J.clin. Microbiol.16:1052-1056

7. Thorpe, J.E., R.P.Banghman, P.T. Frame, T.A. Wessler, and J.L. Staneck. 1987. Bronchoalveolar lavage for diagnosing acute bacterial pneumoniae. J.Infect. Dis. 155:855-861

8. Brown, M.S., and T.C. Wu. 1986. The Gram stain morphology of fungi, mycobacteria, and Pneumocytis carinii. J.Med .Techno 13:495-499

9. Washington, J.A.1986. Rapid diagnosis by microscopy. Clin. Microbiol. Newsl. 8:135-137

10. Lamanna and Mallette, 1965, Basic BActeriology, 3rd ed., Williams and Wilkins Co., Baltimore.

11. Salton, 1964, The Bacterial cell Wall, Elsevier, Amsterdam.

| 10° <u>C</u> | Storage temperature | 8 | Do not use if package is damaged |
|--------------|---|-------|--|
| IVD | In vitro diagnostic medical device | CE | CE Marking |
| | HiMedia Laboratories Pvt. Limited, C-40, Road No.21Y, MIDC, Wagle Industrial Area, Thane (W) - 400604, MS, India | ECREP | CEpartner4U,ESDOORNLAAN 13,3951 DB MAARN,The Netherlands, www.cepartner4u.eu |

Revision: 02/2022

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMediaTM publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMediaTM Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. Reg.Office : Plot No:C-40, Road No: 21Y, MIDC, Wagle Industrial Area, Thane(West)-400604, Maharashtra, INDIA. Tel:00-91-22-61471919/61169797/69034800, Fax:00-91-22-61471920. Email : techhelp@himedialabs.com Website : www.himedialabs.com



Furazolidone FR 50 mcg

SD015

Furazolidone FR 50 mcg discs are used for antimicrobial susceptibility testing of bacterial cultures as per Bauer-Kirby Method

Composition

*IngredientsConcentrationFurazolidone50 mcg/disc

Susceptibility Test Procedure:

- 1. Prepare plates with Mueller Hinton Agar (M173/M1084) for rapidly growing aerobic organisms as per Bauer-Kirby Method. The medium in the plates should be sterile and should have a depth of about 4 mm.
- 2. Inoculate 4-5 similar colonies with a wire, needle or loop to 5 ml Tryptone Soya Broth (M011) and incubate at 35-37°C for2-8 hours until light to moderate turbidity develops. Compare the inoculum turbidity with that of standard 0.5 McFarland (prepared by mixing 0.5 ml of 1.175% barium chloride and 99.5 ml of 0.36N sulfuric acid). Dilute the inoculum or incubate further as necessary to attain comparative turbidity. Alternatively, the inoculum can be standardized by other appropriate optical method (0.08 0.13 OD turbid suspension at 625 nm)
- 3. Dip a sterile non-toxic cotton swab on a wooden applicator into the standardized inoculum and rotate the soaked swab firmly against the upper inside wall of the tube to express excess fluid. Streak the entire agar surface of the plate with the swab three times, turning the plate at 60° angle between each streaking. Allow the inoculum to dry for 5 15 minutes with lid in place.
- 4. Apply the discs using aseptic technique. When using cartridges, the discs can be applied using the specially designed applicator. When the vials are used, apply the discs using sterile forceps.
- 5. Deposit the discs with centers at least 24 mm apart. For fastidious organisms and for Penicillins and Cephalosporins, the discs should preferably be deposited with centers 30 mm apart.
- 6. Incubate immediately at $35 \pm 2^{\circ}$ C and examine after 16-18 hours or longer, if necessary. For fastidious organisms incubate at appropriate temperature and time.
- 7. Measure the zones showing complete inhibition and record the diameters of the zones to the nearest millimeter using a calibrated instrument like zone scales (PW096/PW297)

Principle:

Antimicrobial susceptibility testing (AST) of bacterial and fungal isolates is a common and important technique in most clinical laboratories. The results of these tests are used for selection of the most appropriate antimicrobial agent(s) for treatment against the infectious organisms. Till the 1950s, laboratories were lacking in the methodologies and equipments for the accurate determination of in vitro responses of organisms to antimicrobial agents. Bauer et al (1) began the development of standardized methods for antimicrobial susceptibility testing, using disc diffusion system. However the susceptibility results may not always correlate with the patient's response to therapy. The response of an infected patient to antimicrobial agent(s) is a complex interrelationship of host responses, drug dynamics and microbial activity. Antimicrobial susceptibility tests are either quantitative or qualitative. Disc diffusion test is a qualitative test method. The National Committee for Clinical Laboratory Standards (NCCLS), now known as Clinical Laboratory Standards Institute (CLSI) has published comprehensive documents regarding the disc diffusion systems. The agar disc diffusion test is the most convenient and widely used method for routine antimicrobial susceptibility testing. In subsequent and current practice, antimicrobial impregnated paper discs are applied onto the agar surface. Based on the Bauer-Kirby Method, standardized reference procedures for the disc systems were published by WHO and FDA and are periodically updated by the CLSI (formerly NCCLS)(2). For any antimicrobial testing, Quality control or clinical testing, the method to be followed is the same as mentioned above.

However few precautions are to be maintained while handling of the Sensitivity discs,

- On receipt the discs are to be immediately stored at the recommended temperature.
- Medium preparation, Inoculum preparation and incubation to be done as specified.

Quality Control:

Appearance: Filter paper discs of 6mm diameter with printed "FR 50" on centre of each side of the disc.

Cultural response: Average diameter of zone of inhibition observed on Mueller Hinton Agar (M173) after 18 hours incubation at 35-37°C for standard cultures.

| Organisms (ATCC) | Std. zone of diameter (mm) |
|-------------------------|----------------------------|
| E. coli (25922) | 20-25 |
| <i>S.aureus</i> (25923) | 18-22 |

Storage and Shelf-life:

Discs should always be stored at -20°C to +8°C under dry conditions, along with the dessicator provided in individual pack. Use before expiry date on the label.

References:

- 1. Bauer, Kirby, Sherris and Turck, 1966, Am. J. Clin. Path., 45: 493
- Performance standards of Antimicrobial Disc Susceptibility Tests, M100S, 32nd Ed., CLSI Vol. 42 No.2, Feb-2022.
- 3. EUCAST, Breakpoint tables for interpretation of MIC's & zone diameters, version 12.0, valid from 01.01.2022.

Note :

Use following media to carry out susceptibility test

For rapidly growing aerobic organisms : Mueller Hinton Agar (M173/M1084)

For Haemophilus spps : Haemophilus Test Agar (M1259 + FD117)

For S.pneumoniae : Muller Hinton Agar supplemented with 5% Sheep Blood

For Neisseria spps : G.C.Agar +1% defined growth supplement (M434 + FD025)

* Not for Medicinal Use



In vitro diagnostic medical device CE Marking

Storage temperature



Do not use if package is damaged

HiMedia Laboratories Pvt. Limited, C-40, Road No.21Y, MIDC, Wagle Industrial Area, Thane (W) -400604, MS, India



Disclaimer:

CE Partner 4U, Esdoornlaan 13, 3951 DB Maarn The Netherlands, <u>www.cepartner</u> 4u.eu

Revision : 03 / 2022

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia[™] publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia[™] Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.





Nystatin NS 100 Units

SD025

Nystatin NS 100 Units discs are used for antimicrobial susceptibility testing of fungal cultures

Composition

| *Ingredients | Concentration |
|--------------|----------------|
| Nystatin | 100 Units/disc |

Susceptibility Test Procedure:

Preparation of Inoculum:

- 1. Inoculum is prepared by picking five distinct colonies of approximately 1mm from 24 hours old culture grown on Sabouraud Dextrose Agar (M063) and incubated at $35 \pm 2^{\circ}$ C. Colonies are suspended in 5ml of sterile 0.85% Saline.
- 2. Vortex the resulting suspension and adjust the turbidity to yield 1×10^6 5 x 10^6 cells /ml (i.e. 0.5 McFarland standard).

Test Procedure:

- 1. Prepare plates with Muller Hinton Agar + 2% Glucose + 0.5 mcg/ml Methylene Blue Dye (M1825) for carrying out susceptibility of antifungal discs. The medium in the plates should be sterile and have a depth of about 4 mm.
- 2. Dip a sterile non-toxic cotton swab on a wooden applicator into the standardized inoculum (turbidity so adjusted, as to obtain semi confluent growth on the petri plate) and rotate the soaked swab firmly against the upper inside wall of the tube to express excess fluid. Streak the entire agar surface of the plate with the swab three times, turning the plate at 60° angle between each streaking. Allow the inoculum to dry for 5 15 minutes with lid in place.
- 3. Apply the discs using aseptic technique. Deposit the discs with centers at least 24 mm apart.
- 4. Invert the plates and place in an incubator set to $35 \pm 2^{\circ}$ C within 15 minutes after the discs are applied.
- 5. Examine each plate after 20 24 hours of incubation. If plate was satisfactorily streaked the resulting zones of inhibition will be uniformly circular and there will be a semi-confluent lawn of growth. Read at 48 hours only when insufficient growth is observed after 24 hours incubation.

Principle:

Antimicrobial susceptibility testing (AST) of bacterial and fungal isolates is a common and important technique in most clinical laboratories. The results of these tests are used for selection of the most appropriate antimicrobial agent(s) for treatment against the infectious organisms. Till the 1950s, laboratories were lacking in the methodologies and equipments for the accurate determination of in vitro responses of organisms to antimicrobial agents. Bauer et al (1) began the development of standardized methods for antimicrobial susceptibility testing, using disc diffusion system. However the susceptibility results may not always correlate with the patient's response to therapy. The response of an infected patient to antimicrobial agent(s) is a complex interrelationship of host responses, drug dynamics and microbial activity. Antimicrobial susceptibility tests are either quantitative or qualitative. Disc diffusion test is a qualitative test method. The National Committee for Clinical Laboratory Standards (NCCLS), now known as Clinical Laboratory Standards Institute (CLSI) has published comprehensive documents regarding the disc diffusion systems. The agar disc diffusion test is the most convenient and widely used method for routine antimicrobial susceptibility testing. In subsequent and current practice, antimicrobial impregnated paper discs are applied onto the agar surface. Based on the Bauer-Kirby Method, standardized reference procedures for the disc systems were published by WHO and FDA and are periodically updated by the CLSI (formerly NCCLS)(2). For any antimicrobial testing, Quality control or clinical testing, the method to be followed is the same as mentioned above.

However few precautions are to be maintained while handling of the Sensitivity discs,

- On receipt the discs are to be immediately stored at the recommended temperature.
- Medium preparation, Inoculum preparation and incubation to be done as specified.

Quality Control:

Appearance: Filter paper discs of 6mm diameter with printed "NS 100" on centre of each side of the disc.

Cultural response: Average diameter of zone of inhibition observed on Muller Hinton Agar + 2% Glucose + 0.5 mcg/ml Methylene Blue Dye after 24-48 hours incubation at 35-37°C for standard cultures.

| Organisms (ATCC) | Std. zone of diameter (mm) |
|----------------------------|----------------------------|
| <i>C.albicans</i> (90028)* | 19-27 |
| C.parapsilosis (22019)* | 16-25 |
| C.tropicalis (750)* | 16-21 |
| C.krusei(6528)* | 15-20 |
| C.albicans(10231) | 15-23 |
| S.cerevisiae(9763) | 17-25 |

* = Q.C. Strains recommended by CLSI

Storage and Shelf-life:

Discs should always be stored at -20°C to +8°C under dry conditions, along with the dessicator provided in individual pack. Use before expiry date on the label.

References:

1. Method for Antifungal Disk Diffusion Susceptibility Testing of Yeasts; Approved Guidelines-Second edition Vol.29 No.17, August- 2009 CLSI document M44-A2. For more details refer to this volume

Note :

Use following media to carry out susceptibility test For rapidly growing aerobic organisms : Mueller Hinton Agar (M173/M1084) For *Haemophilus* spps : Haemophilus Test Agar (M1259 + FD117) For *S.pneumoniae* : Muller Hinton Agar supplemented with 5% Sheep Blood For Neisseria spps : G.C.Agar +1% defined growth supplement (M434 + FD025)

* Not for Medicinal Use

 IVD
 In vitro diagnostic medical device

 CE
 CE Marking



Storage temperature



Do not use if package is damaged

HiMedia Laboratories Pvt. Limited, C-40, Road No.21Y, MIDC, Wagle Industrial Area, Thane (W) -400604, MS, India

REP EC

CE Partner 4U, Esdoornlaan 13, 3951 DB Maarn The Netherlands, <u>www.cepartner</u> 4u.eu

Disclaimer :

Revision : 03 / 2022

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia[™] publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia[™] Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.



Penicillin G P 1 unit

SD089

Penicillin G P 1 unit discs are used for antimicrobial susceptibility testing of bacterial cultures as per Kirby-Bauer Method

Composition

| Ingredients* | Concentration |
|--------------|---------------|
| Penicillin G | 1 unit/disc |

Susceptibility Test Procedure:

- 1. Prepare plates with Mueller Hinton Agar (M173/M1084) for rapidly growing aerobic organisms as per Kirby-Bauer Method. The medium in the plates should be sterile and should have a depth of about 4 mm.
- 2. Inoculate 4-5 similar colonies with a wire, needle or loop to 5 ml Tryptone Soya Broth (M011) and incubate at 35-37°C for 2-8 hours until light to moderate turbidity develops. Compare the inoculum turbidity with that of standard 0.5 McFarland (prepared by mixing 0.5 ml of 1.175%barium chloride and 99.5 ml of 0.36N sulfuric acid). Dilute the inoculum or incubate further as necessary to attain comparative turbidity. Alternatively, the inoculum can be standardized by other appropriate optical method (0.08 0.13 OD turbid suspension at 625 nm)
- 3. Dip a sterile non-toxic cotton swab on a wooden applicator into the standardized inoculum and rotate the soaked swab firmly against the upper inside wall of the tube to express excess fluid. Streak the entire agar surface of the plate with the swab three times, turning the plate at 60° angle between each streaking. Allow the inoculum to dry for 5 15 minutes with lid in place.
- 4. Apply the discs using aseptic technique. When using cartridges, the discs can be applied using the specially designed applicator. When the vials are used, apply the discs using sterile forceps.
- 5. Deposit the discs with centers at least 24 mm apart. For fastidious organisms and for Penicillins and Cephalosporins, the discs should preferably be deposited with centers 30 mm apart.
- 6. Incubate immediately at $35 \pm 2^{\circ}$ C and examine after 16-18 hours or longer, if necessary. For fastidious organisms incubate at appropriate temperature and time.
- 7. Measure the zones showing complete inhibition and record the diameters of the zones to the nearest millimeter using a calibrated instrument like zone scales (PW096/PW297)

Principle:

Antimicrobial susceptibility testing (AST) of bacterial and fungal isolates is a common and important technique in most clinical laboratories. The results of these tests are used for selection of the most appropriate antimicrobial agent(s) for treatment against the infectious organisms. Till the 1950s, laboratories were lacking in the methodologies and equipments for the accurate determination of in vitro responses of organisms to antimicrobial agents. Bauer et al (1) began the development of standardized methods for antimicrobial susceptibility testing, using disc diffusion system. However the susceptibility results may not always correlate with the patient's response to therapy. The response of an infected patient to antimicrobial agent(s) is a complex interrelationship of host responses, drug dynamics and microbial activity. Antimicrobial susceptibility tests are either quantitative or qualitative. Disc diffusion test is a qualitative test method. The National Committee for Clinical Laboratory Standards (NCCLS), now known as Clinical Laboratory Standards Institute (CLSI) has published comprehensive documents regarding the disc diffusion systems. The agar disc diffusion test is the most convenient and widely used method for routine antimicrobial susceptibility testing. In subsequent and current practice, antimicrobial impregnated paper discs are applied onto the agar surface. Based on the Kirby-Bauer Method, standardized reference procedures for the disc systems were published by WHO and FDA and are periodically updated by the CLSI (formerly NCCLS)(2). For any antimicrobial testing, Quality control or clinical testing, the method to be followed is the same as mentioned above.

However few precautions are to be maintained while handling of the Sensitivity discs,

- On receipt the discs are to be immediately stored at the recommended temperature.
- Medium preparation, Inoculum preparation and incubation to be done as specified.

Quality Control:

Appearance: Filter paper discs of 6mm diameter with printed "P 1" on centre of each side of the disc.

Cultural response: Average diameter of zone of inhibition observed on Mueller Hinton Agar (M173) after 18 hours incubation at 35-37°C for standard cultures.

| Organisms (ATCC) | Std. zone of diameter (mm) |
|-------------------------|----------------------------|
| <i>S.aureus</i> (29213) | 12-18 |

Storage and Shelf-life:

On receipt discs should always be stored at -20°C under dry conditions, along with the dessicator provided in individual pack. Use before expiry date on the label.

References:

- 1. Bauer, Kirby, Sherris and Turck, 1966, Am. J. Clin. Path., 45: 493
- Performance standards of Antimicrobial Disc Susceptibility Tests, M100S, 32nd Ed., CLSI Vol. 42 No.2, Feb-2022. 2.
- EUCAST, Breakpoint tables for interpretation of MIC's & zone diameters, version 12.0, valid from 01.01.2022. 3.

Note :

Use following media to carry out susceptibility test For rapidly growing aerobic organisms : Mueller Hinton Agar (M173/M1084)

For Haemophilus spps : Haemophilus Test Agar (M1259 + FD117)

For S.pneumoniae : Muller Hinton Agar supplemented with 5% Sheep Blood

For Neisseria spps : G.C.Agar +1% defined growth supplement (M434 + FD025)



Do not use if package is

HiMedia Laboratories Pvt. Limited, C-40, Road No.21Y, MIDC, Wagle Industrial Area, Thane (W) -400604, MS, India

EC REP

CE Partner 4U, Esdoornlaan 13, 3951 DB Maarn The Netherlands, www.cepartner 4u.eu

* Not for Medicinal Use

Disclaimer :

Revision : 02 / 2022

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia[™] publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. Corporate Office : C-40, Road No.21Y, MIDC, Wagle Industrial Area, Thane (W) - 400604, India. Customer care No.: 022-6147 1919 Email: techhelp@himedialabs.com Website: www.himedialabs.com





Ceftriaxone CTR 10 mcg

SD109

Ceftriaxone CTR 10 mcg discs are used for antimicrobial susceptibility testing of bacterial cultures as per Kirby-Bauer Method

Composition

| *Ingredients | Concentration |
|--------------|---------------|
| Ceftriaxone | 10 mcg/disc |

Susceptibility Test Procedure:

- 1. Prepare plates with Mueller Hinton Agar (M173/M1084) for rapidly growing aerobic organisms as per Kirby-Bauer Method. The medium in the plates should be sterile and should have a depth of about 4 mm.
- 2. Inoculate 4-5 similar colonies with a wire, needle or loop to 5 ml Tryptone Soya Broth (M011) and incubate at 35-37°C for2-8 hours until light to moderate turbidity develops. Compare the inoculum turbidity with that of standard 0.5 McFarland (prepared by mixing 0.5 ml of 1.175% barium chloride and 99.5 ml of 0.36N sulfuric acid). Dilute the inoculum or incubate further as necessary to attain comparative turbidity. Alternatively, the inoculum can be standardized by other appropriate optical method (0.08 0.13 OD turbid suspension at 625 nm)
- 3. Dip a sterile non-toxic cotton swab on a wooden applicator into the standardized inoculum and rotate the soaked swab firmly against the upper inside wall of the tube to express excess fluid. Streak the entire agar surface of the plate with the swab three times, turning the plate at 60° angle between each streaking. Allow the inoculum to dry for 5 15 minutes with lid in place.
- 4. Apply the discs using aseptic technique. When using cartridges, the discs can be applied using the specially designed applicator. When the vials are used, apply the discs using sterile forceps.
- 5. Deposit the discs with centers at least 24 mm apart. For fastidious organisms and for Penicillins and Cephalosporins, the discs should preferably be deposited with centers 30 mm apart.
- 6. Incubate immediately at $35 \pm 2^{\circ}$ C and examine after 16-18 hours or longer, if necessary. For fastidious organisms incubate at appropriate temperature and time.
- 7. Measure the zones showing complete inhibition and record the diameters of the zones to the nearest millimeter using a calibrated instrument like zone scales (PW096/PW297)

Principle:

Antimicrobial susceptibility testing (AST) of bacterial and fungal isolates is a common and important technique in most clinical laboratories. The results of these tests are used for selection of the most appropriate antimicrobial agent(s) for treatment against the infectious organisms. Till the 1950s, laboratories were lacking in the methodologies and equipments for the accurate determination of in vitro responses of organisms to antimicrobial agents. Bauer et al (1) began the development of standardized methods for antimicrobial susceptibility testing, using disc diffusion system. However the susceptibility results may not always correlate with the patient's response to therapy. The response of an infected patient to antimicrobial agent(s) is a complex interrelationship of host responses, drug dynamics and microbial activity. Antimicrobial susceptibility tests are either quantitative or qualitative. Disc diffusion test is a qualitative test method. The National Committee for Clinical Laboratory Standards (NCCLS), now known as Clinical Laboratory Standards Institute (CLSI) has published comprehensive documents regarding the disc diffusion systems. The agar disc diffusion test is the most convenient and widely used method for routine antimicrobial susceptibility testing. In subsequent and current practice, antimicrobial impregnated paper discs are applied onto the agar surface. Based on the Bauer-Kirby Method, standardized reference procedures for the disc systems were published by WHO and FDA and are periodically updated by the CLSI (formerly NCCLS)(2). For any antimicrobial testing, Quality control or clinical testing, the method to be followed is the same as mentioned above.

However few precautions are to be maintained while handling of the Sensitivity discs,

- On receipt the discs are to be immediately stored at the recommended temperature.
- Medium preparation, Inoculum preparation and incubation to be done as specified.

Quality Control:

Appearance: Filter paper discs of 6mm diameter with printed "CTR 10" on centre of each side of the disc.

Cultural response: Average diameter of zone of inhibition observed on Mueller Hinton Agar (M173) after 18 hours incubation at 35-37°C for standard cultures.

| Organisms (ATCC) | Std. zone of diameter (mm) |
|-------------------------|----------------------------|
| E. coli (25922) | 29-35 |
| <i>S.aureus</i> (25923) | 22-28 |
| P.aeruginosa (27853) | 17-23 |

Storage and Shelf-life:

On receipt discs should always be stored at -20°C under dry conditions, along with the dessicator provided in individual pack. Use before expiry date on the label.

References:

- 1. Bauer, Kirby, Sherris and Turck, 1966, Am. J. Clin. Path., 45: 493
- 2. Performance standards of Antimicrobial Disc Susceptibility Tests, M100S, 32nd Ed., CLSI Vol. 42 No.2, Feb-2022.
- 3. EUCAST, Breakpoint tables for interpretation of MIC's & zone diameters, version 12.0, valid from 01.01.2022.

Note :

Use following media to carry out susceptibility test

For rapidly growing aerobic organisms : Mueller Hinton Agar (M173/M1084)

For Haemophilus spps : Haemophilus Test Agar (M1259 + FD117)

For S.pneumoniae : Muller Hinton Agar supplemented with 5% Sheep Blood

For Neisseria spps : G.C.Agar +1% defined growth supplement (M434 + FD025)

* Not for Medicinal Use



HiMedia Laboratories Pvt. Limited, C-40, Road No.21Y, MIDC, Wagle Industrial Area, Thane (W) -400604, MS, India



CE Partner 4U ,Esdoornlaan 13, 3951 DB Maarn The Netherlands, www.cepartner 4u.eu

Disclaimer:

Revision : 03 / 2022

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMediaTM publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. Corporate Office : C-40, Road No.21Y, MIDC, Wagle Industrial Area, Thane (W) - 400604, India. Customer care No.: 022-6147 1919 Email: techhelp@himedialabs.com Website: www.himedialabs.com



Amphotericin-B AP 100 units

SD111

Amphotericin-B AP 100 units discs are used for antimicrobial susceptibility testing of fungal cultures

Composition

| *Ingredients | Concentration |
|----------------|----------------|
| Amphotericin-B | 100 units/disc |

Susceptibility Test Procedure:

Preparation of Inoculum:

- 1. Inoculum is prepared by picking five distinct colonies of approximately 1mm from 24 hours old culture grown on Sabouraud Dextrose Agar (M063) and incubated at $35 \pm 2^{\circ}$ C. Colonies are suspended in 5ml of sterile 0.85% Saline.
- 2. Vortex the resulting suspension and adjust the turbidity to yield 1×10^6 5 x 10^6 cells /ml (i.e. 0.5 McFarland standard).

Test Procedure:

- Prepare plates with Muller Hinton Agar + 2% Glucose + 0.5 mcg/ml Methylene Blue Dye (M1825) for carrying out susceptibility of antifungal discs. The medium in the plates should be sterile and have a depth of about 4 mm.
- 2. Dip a sterile non-toxic cotton swab on a wooden applicator into the standardized inoculum (turbidity so adjusted, as to obtain semi confluent growth on the petri plate) and rotate the soaked swab firmly against the upper inside wall of the tube to express excess fluid. Streak the entire agar surface of the plate with the swab three times, turning the plate at 60° angle between each streaking. Allow the inoculum to dry for 5 15 minutes with lid in place.
- 3. Apply the discs using aseptic technique. Deposit the discs with centers at least 24 mm apart.
- 4. Invert the plates and place in an incubator set to $35 \pm 2^{\circ}$ C within 15 minutes after the discs are applied.
- 5. Examine each plate after 20 24 hours of incubation. If plate was satisfactorily streaked the resulting zones of inhibition will be uniformly circular and there will be a semi-confluent lawn of growth. Read at 48 hours only when insufficient growth is observed after 24 hours incubation.

Principle:

Antimicrobial susceptibility testing (AST) of bacterial and fungal isolates is a common and important technique in most clinical laboratories. The results of these tests are used for selection of the most appropriate antimicrobial agent(s) for treatment against the infectious organisms. Till the 1950s, laboratories were lacking in the methodologies and equipments for the accurate determination of in vitro responses of organisms to antimicrobial agents. Bauer et al (1) began the development of standardized methods for antimicrobial susceptibility testing, using disc diffusion system. However the susceptibility results may not always correlate with the patient's response to therapy. The response of an infected patient to antimicrobial agent(s) is a complex interrelationship of host responses, drug dynamics and microbial activity. Antimicrobial susceptibility tests are either quantitative or qualitative. Disc diffusion test is a qualitative test method. The National Committee for Clinical Laboratory Standards (NCCLS), now known as Clinical Laboratory Standards Institute (CLSI) has published comprehensive documents regarding the disc diffusion systems. The agar disc diffusion test is the most convenient and widely used method for routine antimicrobial susceptibility testing. In subsequent and current practice, antimicrobial impregnated paper discs are applied onto the agar surface. Based on the Bauer-Kirby Method, standardized reference procedures for the disc systems were published by WHO and FDA and are periodically updated by the CLSI (formerly NCCLS)(2). For any antimicrobial testing, Quality control or clinical testing, the method to be followed is the same as mentioned above.

However few precautions are to be maintained while handling of the Sensitivity discs,

- On receipt the discs are to be immediately stored at the recommended temperature.
- Medium preparation, Inoculum preparation and incubation to be done as specified.
Quality Control:

Appearance: Filter paper discs of 6mm diameter with printed "AP 100" on centre of each side of the disc.

Cultural response: Average diameter of zone of inhibition observed on Muller Hinton Agar + 2% Glucose + 0.5 mcg/ml Methylene Blue Dye after 24-48 hours incubation at 35-37°C for standard cultures.

| Organisms (ATCC) | Std. zone of diameter (mm) | |
|-------------------------|----------------------------|--|
| C.albicans (90028)* | 10-17 | |
| C.parapsilosis (22019)* | 11-20 | |
| C.tropicalis (750)* | 8-12 | |
| C.krusei (6528)* | 9-14 | |
| C.albicans(10231) | 10-18 | |
| S.cerevesiae (9763) | 11-18 | |

* = Q.C. Strains recommended by CLSI

Storage and Shelf-life:

Discs should always be stored at -20° C to $+8^{\circ}$ C under dry conditions, along with the dessicator provided in individual pack. Use before expiry date on the label.

References:

1. Method for Antifungal Disk Diffusion Susceptibility Testing of Yeasts; Approved Guidelines-Second edition Vol.29 No.17, August- 2009 CLSI document M44-A2. For more details refer to this volume

* Not for Medicinal Use

Revision : 1 / 2012

CE

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia[™] publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia[™] Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. A-516, Swastik Disha Business Park, Via Vadhani Ind. Est., LBS Marg, Mumbai-400086, India. Customer care No.: 022-6147 1919 Email: techhelp@himedialabs.com





Itraconazole IT 10 mcg

SD221

Itraconazole IT 10 mcg discs are used for antimicrobial susceptibility testing of fungal cultures

Composition

| , Ingredients | Concentration |
|---------------|---------------|
| Itraconazole | 10 mcg/disc |

Susceptibility Test Procedure:

Preparation of Inoculum:

- 1. Inoculum is prepared by picking five distinct colonies of approximately 1mm from 24 hours old culture grown on Sabouraud Dextrose Agar (M063) and incubated at $35 \pm 2^{\circ}$ C. Colonies are suspended in 5ml of sterile 0.85% Saline.
- 2. Vortex the resulting suspension and adjust the turbidity to yield 1×10^6 5 x 10^6 cells /ml (i.e. 0.5 McFarland standard).

Test Procedure:

- Prepare plates with Muller Hinton Agar + 2% Glucose + 0.5 mcg/ml Methylene Blue Dye (M1825) for carrying out susceptibility of antifungal discs. The medium in the plates should be sterile and have a depth of about 4 mm.
- 2. Dip a sterile non-toxic cotton swab on a wooden applicator into the standardized inoculum (turbidity so adjusted, as to obtain semi confluent growth on the petri plate) and rotate the soaked swab firmly against the upper inside wall of the tube to express excess fluid. Streak the entire agar surface of the plate with the swab three times, turning the plate at 60° angle between each streaking. Allow the inoculum to dry for 5 15 minutes with lid in place.
- 3. Apply the discs using aseptic technique. Deposit the discs with centers at least 24 mm apart.
- 4. Invert the plates and place in an incubator set to $35 \pm 2^{\circ}$ C within 15 minutes after the discs are applied.
- 5. Examine each plate after 20 24 hours of incubation. If plate was satisfactorily streaked the resulting zones of inhibition will be uniformly circular and there will be a semi-confluent lawn of growth. Read at 48 hours only when insufficient growth is observed after 24 hours incubation.

Principle:

Antimicrobial susceptibility testing (AST) of bacterial and fungal isolates is a common and important technique in most clinical laboratories. The results of these tests are used for selection of the most appropriate antimicrobial agent(s) for treatment against the infectious organisms. Till the 1950s, laboratories were lacking in the methodologies and equipments for the accurate determination of in vitro responses of organisms to antimicrobial agents. Bauer et al (1) began the development of standardized methods for antimicrobial susceptibility testing, using disc diffusion system. However the susceptibility results may not always correlate with the patient's response to therapy. The response of an infected patient to antimicrobial agent(s) is a complex interrelationship of host responses, drug dynamics and microbial activity. Antimicrobial susceptibility tests are either quantitative or qualitative. Disc diffusion test is a qualitative test method. The National Committee for Clinical Laboratory Standards (NCCLS), now known as Clinical Laboratory Standards Institute (CLSI) has published comprehensive documents regarding the disc diffusion systems. The agar disc diffusion test is the most convenient and widely used method for routine antimicrobial susceptibility testing. In subsequent and current practice, antimicrobial impregnated paper discs are applied onto the agar surface. Based on the Bauer-Kirby Method, standardized reference procedures for the disc systems were published by WHO and FDA and are periodically updated by the CLSI (formerly NCCLS)(2). For any antimicrobial testing, Quality control or clinical testing, the method to be followed is the same as mentioned above.

However few precautions are to be maintained while handling of the Sensitivity discs,

- On receipt the discs are to be immediately stored at the recommended temperature.
- Medium preparation, Inoculum preparation and incubation to be done as specified.

Quality Control:

Appearance: Filter paper discs of 6mm diameter with printed "IT 10" on centre of each side of the disc.

Cultural response: Average diameter of zone of inhibition observed on Muller Hinton Agar + 2% Glucose + 0.5 mcg/ml Methylene Blue Dye after 24-48 hours incubation at 35-37°C for standard cultures.

| Organisms (ATCC) | Std. zone of diameter (mm) |
|----------------------------|----------------------------|
| <i>C.albicans</i> (90028)* | 16-20 |
| C.parapsilosis (22019)* | 11-18 |
| C.tropicalis (750)* | 8-13 |
| C.krusei(6528)* | 8-15 |
| C.albicans(10231) | 18-22 |

* = Q.C. Strains recommended by CLSI

Storage and Shelf-life:

Discs should always be stored at -20° C to $+8^{\circ}$ C under dry conditions, along with the dessicator provided in individual pack. Use before expiry date on the label.

References:

1. Method for Antifungal Disk Diffusion Susceptibility Testing of Yeasts; Approved Guidelines-Second edition Vol.29 No.17, August- 2009 CLSI document M44-A2. For more details refer to this volume

* Not for Medicinal Use

Revision : 1 / 2012

CE

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia[™] publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia[™] Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. A-516, Swastik Disha Business Park, Via Vadhani Ind. Est., LBS Marg, Mumbai-400086, India. Customer care No.: 022-6147 1919 Email: techhelp@himedialabs.com





Fluconazole FLC 25 mcg

SD232

Fluconazole FLC 25mcg discs are used for antimicrobial susceptibility testing of fungal cultures

Composition *Ingredients

| *Ingredients | Concentration |
|--------------|---------------|
| Fluconazole | 25mcg/disc |

Susceptibility Test Procedure:

Preparation of Inoculum:

- 1. Inoculum is prepared by picking five distinct colonies of approximately 1mm from 24 hours old culture grown on Sabouraud Dextrose Agar (M063) and incubated at $35 \pm 2^{\circ}$ C. Colonies are suspended in 5ml of sterile 0.85% Saline.
- 2. Vortex the resulting suspension and adjust the turbidity to yield 1×10^6 5 x 10^6 cells /ml (i.e. 0.5 McFarland standard).

Test Procedure:

- Prepare plates with Muller Hinton Agar + 2% Glucose + 0.5 mcg/ml Methylene Blue Dye (M1825) for carrying out susceptibility of antifungal discs. The medium in the plates should be sterile and have a depth of about 4 mm.
- 2. Dip a sterile non-toxic cotton swab on a wooden applicator into the standardized inoculum (turbidity so adjusted, as to obtain semi confluent growth on the petri plate) and rotate the soaked swab firmly against the upper inside wall of the tube to express excess fluid. Streak the entire agar surface of the plate with the swab three times, turning the plate at 60° angle between each streaking. Allow the inoculum to dry for 5 15 minutes with lid in place.
- 3. Apply the discs using aseptic technique. Deposit the discs with centers at least 24 mm apart.
- 4. Invert the plates and place in an incubator set to $35 \pm 2^{\circ}$ C within 15 minutes after the discs are applied.
- 5. Examine each plate after 20 24 hours of incubation. If plate was satisfactorily streaked the resulting zones of inhibition will be uniformly circular and there will be a semi-confluent lawn of growth. Read at 48 hours only when insufficient growth is observed after 24 hours incubation.

Principle:

Antimicrobial susceptibility testing (AST) of bacterial and fungal isolates is a common and important technique in most clinical laboratories. The results of these tests are used for selection of the most appropriate antimicrobial agent(s) for treatment against the infectious organisms. Till the 1950s, laboratories were lacking in the methodologies and equipments for the accurate determination of in vitro responses of organisms to antimicrobial agents. Bauer et al (1) began the development of standardized methods for antimicrobial susceptibility testing, using disc diffusion system. However the susceptibility results may not always correlate with the patient's response to therapy. The response of an infected patient to antimicrobial agent(s) is a complex interrelationship of host responses, drug dynamics and microbial activity. Antimicrobial susceptibility tests are either quantitative or qualitative. Disc diffusion test is a qualitative test method. The National Committee for Clinical Laboratory Standards (NCCLS), now known as Clinical Laboratory Standards Institute (CLSI) has published comprehensive documents regarding the disc diffusion systems. The agar disc diffusion test is the most convenient and widely used method for routine antimicrobial susceptibility testing. In subsequent and current practice, antimicrobial impregnated paper discs are applied onto the agar surface. Based on the Bauer-Kirby Method, standardized reference procedures for the disc systems were published by WHO and FDA and are periodically updated by the CLSI (formerly NCCLS)(2). For any antimicrobial testing, Quality control or clinical testing, the method to be followed is the same as mentioned above.

However few precautions are to be maintained while handling of the Sensitivity discs,

- On receipt the discs are to be immediately stored at the recommended temperature.
- Medium preparation, Inoculum preparation and incubation to be done as specified.

Interpretation:

Use following interpretive criteria for susceptibility categorization.

| Antimicrobial | | Resistant | S-DD# | S |
|----------------------|-----------------------------|------------|-------|------------|
| agent | Interpretative criteria for | mm or less | mm | mm or more |
| Fluconazole 25mcg | Candida spps | 14 | 15-18 | 19 |

S-DD - Susceptible - Dose dependent

Quality Control:

Appearance: Filter paper discs of 6mm diameter with printed "FLC 25" on centre of each side of the disc.

Cultural response: Average diameter of zone of inhibition observed on Muller Hinton Agar + 2% Glucose + 0.5 mcg/ml Methylene Blue Dye after 24-48 hours incubation at 35-37°C for standard cultures.

| Organisms (ATCC) | Std. zone of diameter (mm) | | |
|-------------------------|----------------------------|--|--|
| C.albicans (90028)* | 28-39 | | |
| C.parapsilosis (22019)* | 22-33 | | |
| C.tropicalis (750)* | 26-37 | | |
| C.albicans(10231) | 25-30 | | |

* = Q.C. Strains recommended by CLSI

Storage and Shelf-life:

Discs should always be stored at -20°C to +8°C under dry conditions, along with the dessicator provided in individual pack. Use before expiry date on the label.

References:

- 1. Method for Antifungal Disk Diffusion Susceptibility Testing of Yeasts; Approved Guidelines-Second edition Vol.29 No.17, August- 2009 CLSI document M44-A2. For more details refer to this volume
- 2. Zone Diameter Interpretive Standards, Corresponding Minimal Inhibitory Concentration (MIC) Interpretive Breakpoints, and Quality Control Limits for Antifungal Disk Diffusion Susceptibility Testing Of Yeasts, Third International Supplement CLSI document - M44-S3.

* Not for Medicinal Use

Revision : 1 / 2012

CE

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia[™] publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia[™] Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. A-516, Swastik Disha Business Park, Via Vadhani Ind. Est., LBS Marg, Mumbai-400086, India. Customer care No.: 022-6147 1919 Email: techhelp@himedialabs.com



Cefoperazone/Sulbactam CFS 75/10 mcg

SD254

Cefoperazone/Sulbactam CFS 75/10mcg discs are used for antimicrobial susceptibility testing of bacterial cultures as per Bauer-Kirby Method

Composition

| *Ingredients | Concentration |
|------------------------|----------------|
| Cefoperazone/Sulbactam | 75/10 mcg/disc |

Susceptibility Test Procedure:

- 1. Prepare plates with Mueller Hinton Agar (M173/M1084) for rapidly growing aerobic organisms as per Bauer-Kirby Method. The medium in the plates should be sterile and should have a depth of about 4 mm.
- 2. Inoculate 4-5 similar colonies with a wire, needle or loop to 5 ml Tryptone Soya Broth (M011) and incubate at 35-37°C for2-8 hours until light to moderate turbidity develops. Compare the inoculum turbidity with that of standard 0.5 McFarland (prepared by mixing 0.5 ml of 1.175% barium chloride and 99.5 ml of 0.36N sulfuric acid). Dilute the inoculum or incubate further as necessary to attain comparative turbidity. Alternatively, the inoculum can be standardized by other appropriate optical method (0.08 0.13 OD turbid suspension at 625 nm)
- 3. Dip a sterile non-toxic cotton swab on a wooden applicator into the standardized inoculum and rotate the soaked swab firmly against the upper inside wall of the tube to express excess fluid. Streak the entire agar surface of the plate with the swab three times, turning the plate at 60° angle between each streaking. Allow the inoculum to dry for 5 15 minutes with lid in place.
- 4. Apply the discs using aseptic technique. When using cartridges, the discs can be applied using the specially designed applicator. When the vials are used, apply the discs using sterile forceps.
- 5. Deposit the discs with centers at least 24 mm apart. For fastidious organisms and for Penicillins and Cephalosporins, the discs should preferably be deposited with centers 30 mm apart.
- 6. Incubate immediately at $35 \pm 2^{\circ}$ C and examine after 16-18 hours or longer, if necessary. For fastidious organisms incubate at appropriate temperature and time.
- 7. Measure the zones showing complete inhibition and record the diameters of the zones to the nearest millimeter using a calibrated instrument like zone scales (PW096/PW297)

Principle:

Antimicrobial susceptibility testing (AST) of bacterial and fungal isolates is a common and important technique in most clinical laboratories. The results of these tests are used for selection of the most appropriate antimicrobial agent(s) for treatment against the infectious organisms. Till the 1950s, laboratories were lacking in the methodologies and equipments for the accurate determination of in vitro responses of organisms to antimicrobial agents. Bauer et al (1) began the development of standardized methods for antimicrobial susceptibility testing, using disc diffusion system. However the susceptibility results may not always correlate with the patient's response to therapy. The response of an infected patient to antimicrobial agent(s) is a complex interrelationship of host responses, drug dynamics and microbial activity. Antimicrobial susceptibility tests are either quantitative or qualitative. Disc diffusion test is a qualitative test method. The National Committee for Clinical Laboratory Standards (NCCLS), now known as Clinical Laboratory Standards Institute (CLSI) has published comprehensive documents regarding the disc diffusion systems. The agar disc diffusion test is the most convenient and widely used method for routine antimicrobial susceptibility testing. In subsequent and current practice, antimicrobial impregnated paper discs are applied onto the agar surface. Based on the Bauer-Kirby Method, standardized reference procedures for the disc systems were published by WHO and FDA and are periodically updated by the CLSI (formerly NCCLS)(2). For any antimicrobial testing, Quality control or clinical testing, the method to be followed is the same as mentioned above.

However few precautions are to be maintained while handling of the Sensitivity discs,

- On receipt the discs are to be immediately stored at the recommended temperature.
- Medium preparation, Inoculum preparation and incubation to be done as specified.

Quality Control:

Appearance: Filter paper discs of 6mm diameter with printed "CFS 75/10" on centre of each side of the disc.

Cultural response: Average diameter of zone of inhibition observed on Mueller Hinton Agar (M173) after 18 hours incubation at 35-37°C for standard cultures.

| Organisms (ATCC) | Std. zone of diameter (mm) |
|----------------------|----------------------------|
| E. coli (25922) | 27-33 |
| S.aureus (25923) | 23-30 |
| P.aeruginosa (27853) | 20-25 |

Storage and Shelf-life:

Discs should always be stored at -20° C to $+8^{\circ}$ C under dry conditions, along with the dessicator provided in individual pack. Use before expiry date on the label.

References:

- 1. Bauer, Kirby, Sherris and Turck, 1966, Am. J. Clin. Path., 45: 493
- Performance standards of Antimicrobial Disc Susceptibility Tests, M100S, 32nd Ed., CLSI Vol. 42 No.2, Feb-2022.
- 3. EUCAST, Breakpoint tables for interpretation of MIC's & zone diameters, version 12.0, valid from 01.01.2022.

Note :

Use following media to carry out susceptibility test

For rapidly growing aerobic organisms : Mueller Hinton Agar (M173/M1084)

For *Haemophilus* spps : Haemophilus Test Agar (M1259 + FD117)

For S.pneumoniae : Muller Hinton Agar supplemented with 5% Sheep Blood

For Neisseria spps : G.C.Agar +1% defined growth supplement (M434 + FD025)

* Not for Medicinal Use





Voriconazole VRC 1 mcg

SD277

Voriconazole VRC 1 mcg discs are used for antimicrobial susceptibility testing of fungal cultures

| Composition | |
|--------------|---------------|
| *Ingredients | Concentration |
| Voriconazole | 1mcg/disc |

Susceptibility Test Procedure:

Preparation of Inoculum:

- 1. Inoculum is prepared by picking five distinct colonies of approximately 1mm from 24 hours old culture grown on Sabouraud Dextrose Agar (M063) and incubated at $35 \pm 2^{\circ}$ C. Colonies are suspended in 5ml of sterile 0.85% Saline.
- 2. Vortex the resulting suspension and adjust the turbidity to yield 1×10^6 5 x 10^6 cells /ml (i.e. 0.5 McFarland standard).

Test Procedure:

- 1. Prepare plates with Muller Hinton Agar + 2% Glucose + 0.5 mcg/ml Methylene Blue Dye (M1825) for carrying out susceptibility of antifungal discs. The medium in the plates should be sterile and have a depth of about 4 mm.
- 2. Dip a sterile non-toxic cotton swab on a wooden applicator into the standardized inoculum (turbidity so adjusted, as to obtain semi confluent growth on the petri plate) and rotate the soaked swab firmly against the upper inside wall of the tube to express excess fluid. Streak the entire agar surface of the plate with the swab three times, turning the plate at 60° angle between each streaking. Allow the inoculum to dry for 5 15 minutes with lid in place.
- 3. Apply the discs using aseptic technique. Deposit the discs with centers at least 24 mm apart.
- 4. Invert the plates and place in an incubator set to $35 \pm 2^{\circ}$ C within 15 minutes after the discs are applied.
- 5. Examine each plate after 20 24 hours of incubation. If plate was satisfactorily streaked the resulting zones of inhibition will be uniformly circular and there will be a semi-confluent lawn of growth. Read at 48 hours only when insufficient growth is observed after 24 hours incubation.

Principle:

Antimicrobial susceptibility testing (AST) of bacterial and fungal isolates is a common and important technique in most clinical laboratories. The results of these tests are used for selection of the most appropriate antimicrobial agent(s) for treatment against the infectious organisms. Till the 1950s, laboratories were lacking in the methodologies and equipments for the accurate determination of in vitro responses of organisms to antimicrobial agents. Bauer et al (1) began the development of standardized methods for antimicrobial susceptibility testing, using disc diffusion system. However the susceptibility results may not always correlate with the patient's response to therapy. The response of an infected patient to antimicrobial agent(s) is a complex interrelationship of host responses, drug dynamics and microbial activity. Antimicrobial susceptibility tests are either quantitative or qualitative. Disc diffusion test is a qualitative test method. The National Committee for Clinical Laboratory Standards (NCCLS), now known as Clinical Laboratory Standards Institute (CLSI) has published comprehensive documents regarding the disc diffusion systems. The agar disc diffusion test is the most convenient and widely used method for routine antimicrobial susceptibility testing. In subsequent and current practice, antimicrobial impregnated paper discs are applied onto the agar surface. Based on the Bauer-Kirby Method, standardized reference procedures for the disc systems were published by WHO and FDA and are periodically updated by the CLSI (formerly NCCLS)(2). For any antimicrobial testing, Quality control or clinical testing, the method to be followed is the same as mentioned above.

However few precautions are to be maintained while handling of the Sensitivity discs,

- On receipt the discs are to be immediately stored at the recommended temperature.
- Medium preparation, Inoculum preparation and incubation to be done as specified.

Interpretation:

Use following interpretive criteria for susceptibility categorization.

| Antimicrobial | | Resistant | S-DD# | S |
|-----------------------|-----------------------------|------------|-------|------------|
| agent | Interpretative criteria for | mm or less | mm | mm or more |
| Voriconazole 1 mcg | Candida spps | 13 | 14-16 | 17 |

S-DD - Susceptible - Dose dependent

Quality Control:

Appearance: Filter paper discs of 6mm diameter with printed "VRC 1" on centre of each side of the disc.

Cultural response: Average diameter of zone of inhibition observed on Muller Hinton Agar + 2% Glucose + 0.5 mcg/ml Methylene Blue Dye after 24-48 hours incubation at 35-37°C for standard cultures.

| Organisms (ATCC) | Std. zone of diameter (mm) |
|-------------------------|----------------------------|
| C.albicans (90028)* | 31-42 |
| C.parapsilosis (22019)* | 28-37 |
| C.krusei (6528)* | 16-25 |
| C.albicans(10231) | 30-40 |
| S.cerevesiae (9763) | 29-38 |

* = Q.C. Strains recommended by CLSI

Storage and Shelf-life:

Discs should always be stored at -20° C to $+8^{\circ}$ C under dry conditions, along with the dessicator provided in individual pack. Use before expiry date on the label.

References:

- 1. Method for Antifungal Disk Diffusion Susceptibility Testing of Yeasts; Approved Guidelines-Second edition Vol.29 No.17, August- 2009 CLSI document M44-A2. For more details refer to this volume
- 2. Zone Diameter Interpretive Standards, Corresponding Minimal Inhibitory Concentration (MIC) Interpretive Breakpoints, and Quality Control Limits for Antifungal Disk Diffusion Susceptibility Testing Of Yeasts, Third International Supplement CLSI document - M44-S3.

* Not for Medicinal Use

Revision : 1 / 2012

CE

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia[™] publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia[™] Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. A-516, Swastik Disha Business Park, Via Vadhani Ind. Est., LBS Marg, Mumbai-400086, India. Customer care No.: 022-6147 1919 Email: techhelp@himedialabs.com

IVD solutions through partnership



MASTDISCS®

Leading the field with a complete solution for AST and Identification disc testing

- Comprehensive range
- Premium quality products
- Compatible with EUCAST and CLSI standards
- Bespoke service available



Mast Group Ltd., has been a manufacturer of antibiotic susceptibility test products since 1957, and continues to be at the forefront of developments in this field. With an ever-expanding portfolio of antibiotic susceptibility, combination and identification discs, Mast Group Ltd. supplies the most comprehensive range available worldwide. Mast Group Ltd. also offers a range of services including **mast**pharma[®] development and **mast**pharma[®] stability for the evaluation of novel antimicrobial compounds for the pharmaceutical industry.

Disc diffusion is still the most popular method for determining bacterial antibiotic susceptibility in the treatment of infectious disease, as highlighted by the variety of international standardised procedures including The European Committee on Antimicrobial Susceptibility Testing (EUCAST) and The Clinical and Laboratory Standards Institute (CLSI).

MASTDISCS[®] are available as an extensive range of test discs presented in dispensing cartridges or vials, with quality assurance guaranteed in accordance to standardised test protocols as issued by the above national and international organisations.

| Order No | Product | | | Packsize |
|--------------------------------------|---------|--|------------------|----------------------------------|
| MASTDISCS | AST | - ANTIBIOTIC SUSCEPTIBILITY DISCS IN CAP | RTRIDGES AND VIA | LS |
| MASTDISCS | 5 | in Cartridges (5 × 50 discs per pack) | | |
| STOCKCART FUNGCART SPECIALCART | | Stock Susceptibility Cartridge Discs of a single type Stock Antifungal Cartridge Discs of a single type Cartridge Discs made to special order* | UCAST | 1 pack 1 pack Min 18 packs |
| MASTDISCS' | 8 | (99+/- 2 discs per vial) | | |
| STOCKDISC | | Stock Susceptibility Discs of a single type | | 1 vial |
| FUNGDISC | | Stock Antifungal Discs of a single type | | 1 pack |
| SPECIALDISC | | Discs made to special order* | | Min 22 packs |
| TOOL/C | | *Set up charge for Special Discs in vials or cartridges A one off charge for customisation of each new specificatio | n | per new specification |

Disc Diffusion Susceptibility Testing

MASTDISCS"AST

ANTIBIOTIC SUSCEPTIBILITY DISCS IN CARTRIDGES - STOCK RANGE

| scription STANDARD | | | Order Code | | |
|------------------------------------|---------------------------------------|--------------|--------------|--|--|
| Antibiotic & Content µg per disc | | | Cartridges | | |
| (unless otherwise stated) | EUCAST | CLSI | 5 × 50 discs | | |
| | | | | | |
| Amikacin 30 | \checkmark | \checkmark | AK30C | | |
| Amoxicillin/clavulanic acid 2-1 | \checkmark | - | AUG3C | | |
| Amoxicillin/clavulanic acid 20-10 | \checkmark | 1 | AUG30C | | |
| Ampicillin 2 | \checkmark | - | AP2C | | |
| Ampicillin 10 | \checkmark | \checkmark | AP10C | | |
| Ampicillin/Sulbactam 10-10 | \checkmark | \checkmark | SAM20C | | |
| Azithromycin 15 | - | 1 | ATH15C | | |
| Aztreonam 30 | \checkmark | 1 | ATM30C | | |
| Bacitracin 10 units | - | - | BA10C | | |
| Carbenicillin 100 | ✓ | - | PY100C | | |
| Cefaclor 30 | \checkmark | 1 | CFC30C | | |
| Cefadroxil 30 | ✓ | - | CDX30C | | |
| Cefalexin 30 | ✓ | - | CFX30C | | |
| Cefalothin 30 | - | 1 | KF30C | | |
| Cefazolin 30 | - | 1 | CZ30C | | |
| Cefepime 30 | | | CPM30C | | |
| Cefiderocol 30 | | | FDC30C | | |
| Cefixime 5 | | | CFM5C | | |
| Cefotaxime 5 | | - | CTX5C | | |
| Cefotaxime 30 | - | | CTX30C | | |
| Cefoxitin 30 | ./ | | EOX30C | | |
| Cefpodovime 10 | · · · · · · · · · · · · · · · · · · · | | | | |
| Cettoroline 5 | · · · · · · · · · · · · · · · · · · · | • - | CPT5C | | |
| Cettaroline 30 | • | - | CPT30C | | |
| | - | v | CF130C | | |
| | v | - | CAZ20C | | |
| Ceftazidime /avibactor 10.4 | - | ✓ | C7A14C | | |
| Cettazidime/avibactam 10-4 | V | - | CZA140 | | |
| | - | | | | |
| Ceftibuleri 30 | | | | | |
| | | ✓ | C/140C | | |
| | V | - | BPR5C | | |
| | - | | CRU5C | | |
| | | | CRO30C | | |
| Ceturoxime 30 | | | CXM30C | | |
| Chloramphenicol 30 | | | C30C | | |
| Ciprofloxacin 5 | ✓ | | CIP5C | | |
| Clarithromycin 15 | - | | CLA15C | | |
| Clindamycin 2 | \checkmark | \checkmark | CD2C | | |
| Delafloxacin 5 | - | \checkmark | DLX5C | | |
| Doripenem 10 | ✓ | ✓ | DOR10C | | |
| Doxycycline 30 | - | \checkmark | DXT30C | | |
| Eravacycline 20 | \checkmark | \checkmark | ERV20C | | |
| Ertapenem 10 | \checkmark | 1 | ETP10C | | |
| Erythromycin 15 | \checkmark | 1 | E15C | | |
| Fosfomycin/Glucose-6-Phosphate 200 | \checkmark | \checkmark | FOT200C | | |
| Fusidic Acid 10 | \checkmark | - | FC10C | | |
| Gentamicin 10 | \checkmark | \checkmark | GM10C | | |
| Gentamicin 30 | ✓ | - | GM30C | | |
| Gentamicin 120 | - | 1 | GM120C | | |
| Imipenem 10 | ✓ | ✓ | IMI10C | | |
| Imipenem/relebactam 10-25 | - | 1 | IMR35C | | |
| Kanamycin 30 | - | \checkmark | K30C | | |
| Lefamulin 20 | - | 1 | LMU20C | | |
| Levofloxacin 5 | ✓ | √ | LEV5C | | |
| Linezolid 10 | | - | LZD10C | | |
| Linezolid 30 | - | | LZD30C | | |
| Mecillinam 10 | 5 | | MFC10C | | |
| Meropenem 10 | | | MEM10C | | |
| | | | | | |

MASTDISCS"AST

ANTIBIOTIC SUSCEPTIBILITY DISCS IN CARTRIDGES - STOCK RANGE

| Description | STANE | Order Code | |
|--|---------------------------------------|--------------|----------------|
| Antibiotic & Content µg per disc | | | Cartridges |
| (unless otherwise stated) | EUCAST | CLSI | 5 × 50 discs |
| Marananan (abarbaatan 20, 10 | | 1 | |
| Metropeneni/vaborbactarii 20-10 | | ✓ | IVIEV30C |
| Misservalias 20 | ✓ | - | MZ5C |
| Minocycline 30 | - | | MIN3UC |
| Moxifloxacin 5 | | ✓ | MFX5C |
| Mupirocin 200 | | - | MUP200C |
| Nalidixic Acid 30 | | ✓ | NA30C |
| Neomycin 10 | | - | NE10C |
| Netilmicin 10 | 1 | - | NET10C |
| Netilmicin 30 | \checkmark | - | NET30C |
| Nitrofurantoin 100 | \checkmark | - | NI100C |
| Nitrofurantoin 300 | - | \checkmark | NI300C |
| Nitroxoline 30 | 1 | - | NIB30C |
| Norfloxacin 10 | \checkmark | \checkmark | NOR10C |
| Novobiocin 5 | - | - | NO5C |
| Ofloxacin 5 | 1 | ✓ | OFX5C |
| Oxacillin 1 | 1 | 1 | OX1C |
| Pefloxacin 5 | \checkmark | - | PEF5C |
| Penicillin G 1 unit | | - | PG1C |
| Penicillin G 10 units | - | ✓ | PG10C |
| Piperacillin 30 | 1 | - | PRL30C |
| Piperacillin 100 | - | ✓ | PRL100C |
| Piperacillin/tazobactam 30-6 | 1 | - | PTZ36C |
| Piperacillin/tazobactam 100-10 | - | \checkmark | PTZ110C |
| Rifampicin 5 | 1 | ✓ | RP5C |
| Streptomycin 10 | - | \checkmark | S10C |
| Streptomycin 300 | 1 | \checkmark | S300C |
| Tedizolid 2 | 1 | \checkmark | TZD2C |
| Teicoplanin 30 | 1 | 1 | TEC30C |
| Temocillin 30 | | - | TEM30C |
| Tetracycline 30 | | | T30C |
| Ticarcillin 75 | | | TC75C |
| Ticarcillin/clavulanic acid 75-10 | | | TIM85C |
| Tigecycline 15 | | | TGC15C |
| Tobramycin 10 | | | TN10C |
| Trimethonrim 5 | | | TM5C |
| Trimethoprim/sulfamethoxazole 1 25/23 75 | · · · · · · · · · · · · · · · · · · · | • ./ | TS25C |
| Vancomycin 5 | • ./ | • - | VΔ5C |
| Vancomycin 30 | ✓ | - | |
| Right digg | - | V | |
| DIALIK UISUS | - | - | DD00000/0/INCE |

MASTDISCS"AST

NON CLINICAL ANTIBIOTIC SUSCEPTIBILITY DISCS IN CARTRIDGES - STOCK RANGE

| Description Antibiotic & Content μg per disc (unless otherwise stated) | Order Code Cartridges 5 × 50 discs |
|--|--|
| Amoxicillin 25 | A25C/NCE |
| Colistin Sulphate 10 | CO10C/NCE |
| Colistin Sulphate 25 | CO25C/NCE |
| Polymyxin B 300 units | PB300C/NCE |

MASTDISCS"AST

SPECIALIST VETERINARY SUSCEPTIBILITY DISCS - STOCK RANGE

| Description Antibiotic & Content µg per disc (unless otherwise stated) | Order Code Cartridges 5 × 50 discs |
|--|--|
| Cefquinone 30 | CEQ30C/NCE |
| Cefoperazone 30 | CPZ30C/NCE |
| Enrofloxacin 5 | ENF5C/NCE |
| Florfenicol 30 | FFC30C/NCE |
| Gamithromycin 15 | GAM15C/NCE |
| Marbofloxacin 5 | MAR5/NCE |
| Pradofloxacin 5 | PRA5/NCE |
| Tildipirison 60 | TIP60/NCE |
| Tylosin 30 | TY30C/NCE |

Order No Product

Mechanism

MASTDISCS" AST

ANTIBIOTIC SUSCEPTIBILITY DISCS IN VIALS - STOCK RANGE

| Description Antibiotic & Content up per disc | STAND | ARD | Order Code Vials | |
|---|--------------|--------------|---------------------|--|
| (unless otherwise stated) | EUCAST | CLSI | 99+/- 2 discs | |
| Bacitracin 10 units | - | - | BA10 | |
| Cefpodoxime 10 | ✓ | \checkmark | CPD10 | |
| Gentamicin 10 | \checkmark | 1 | GM10 | |
| Metronidazole 5 | \checkmark | - | MZ5 | |
| Nalidixic Acid 30 | ✓ | - | NA30 | |
| Oxacillin 1 | \checkmark | 1 | OX1 | |
| Penicillin G 1 unit | \checkmark | - | PG1 | |
| Vancomycin 5 | ✓ | - | VA5 | |
| Blank discs | - | - | BD0638W/NCE | |

MASTDISCS"AST

NON CLINICAL ANTIFUNGAL SUSCEPTIBILITY DISCS IN CARTRIDGES AND VIALS - STOCK RANGE

| Description Antibiotic & Content µg per disc (unless otherwise stated) | Order Code Cartridges 5 × 50 discs |
|--|--|
| Amphotericin B 20 | AMB20C/NCE |
| Clotrimazole 10 | CTM10C/NCE |
| Econazole 10 | ECN10C/NCE |
| Fluconazole 10 | FCN10C/NCE |
| Fluconazole 25 | FCN25C/NCE |
| Flucytosine 1 | FY1C/NCE |
| Ketoconazole 10 | KCA10C/NCE |
| Miconazole 10 | MCL10C/NCE |
| Nystatin 100 | NY100C/NCE |

Order No Product

Mechanism

MASTDISCS[®] Combi

COMBINATION DISC SETS FOR THE DETECTION OF ANTIBIOTIC RESISTANCE

| D52C | Extended Spectrum ß Lactamase Set | ESBL | 50 tests |
|-------------|--|---|-----------------------|
| D62C | Cefotaxime 30 & Cefotaxime 30/Clavulanic Acid 10 | ESBL | 150 tests |
| D63C | Cefepime 30 & Cefepime 30/Clavulanic Acid 10 | ESBL | 150 tests |
| D64C | Ceftazidime 30 & Ceftazidime 30/Clavulanic Acid 10 | ESBL | 150 tests |
| D66C | Cefpodoxime 10 & Cefpodoxime 10/Clavulanic Acid 1 | ESBL | 150 tests |
| D67C | Extended Spectrum ß Lactamase Set (CPD10) | ESBL | 50 tests |
| D68C | AmpC & ESBL Detection Set | AmpC/ESBL | 50 tests |
| D69C | AmpC Detection Set | AmpC | 50 tests |
| D72C | AmpC, ESBL & Carbapenemase Detection Set | AmpC/ESBL/Carbapenemase | 50 tests |
| D73C | MAST [®] Carba plus | MBL/KPC/OXA | 50 tests |
| D76C | ESBL Detection Set (EUCAST) | ESBL | 50 tests |
| D71C D74 | MAST [®] CAT-ID - For presumptive identification of carba MAST [®] ICT - screening test for the detection of carbap in Enterobacterales, <i>Pseudomonas</i> and <i>Acinetobacter</i> | apenemase production penemase production spp. | 250 tests 25 tests |
| TEM30C | To aid presumptive identification of OXA-48 | | 5×50 discs |
| RAPID CAR | BAPENEMASE DETECTION | | |
| DNA/LYO5 | Rapid molecular carbapenemase detection in Enteroba <i>Pseudomonas</i> spp. and <i>Acinetobacter</i> spp. | acterales, | 10 tests |
| PACE-ID | Colorimetric test for the rapid detection of carbapenen Pseudomonas spp., Acinetobacter spp.and Enterobac | nase producing terales. | 48 tests |

MASTDISCS"AST

PENICILLIN DISCS

For detecting the emergence of penicillin resistance. (99+/- 2 discs or 5 × 50 discs in cartridges)

| PG1 | Penicillin G | 1 unit (vials) | 1 vial |
|------|--------------|---------------------|--------|
| PG1C | Penicillin G | 1 unit (cartridges) | 1 pack |

HIGH CONTENT AMINOGLYCOSIDE DISCS

The following discs are available to special order and subject to 6 weeks lead time. For detecting high level aminoglycoside resistance. $(99+/-2 \text{ discs or } 5 \times 50 \text{ discs in cartridges})$

| GM500C | 500 µg | Min order |
|------------|---------|-----------|
| K1000C/NCE | 1000 µg | 18 packs |
| K1000/NCE | 1000 µg | Min order |
| | | 22 vials |

Any code with a suffix of '/NCE' is not CE marked and for veterinary or research use only.

Order No

Packsize

DISCMASTER DISPENSER

Product

Robust and reliable antimicrobial cartridge disc dispensers designed for use with MASTDISCS® Antimicrobial Susceptibility Test Cartridges.

| MDD65 | MAST [®] DISCMASTER Dispenser - 6 place | 1 |
|------------|---|---|
| SILICA63 | Silica Gel Capsule for MAST [®] DISCMASTER | 4 |
| SHD5 | Single Cartridge Hand Dispenser | 5 |
| CANISTER | Canister for MDD64 models and below | 1 |
| CANISTER65 | Canister for MDD65 MAST® DISCMASTER Dispenser | 1 |





IVD solutions through partnership

Mast House, Derby Road, Bootle, Merseyside L20 1EA Tel: + 44 (0) 151 933 7277 e-mail: sales@mast-group.com

Mast Diagnostica GmbH Feldstrasse 20,DE-23858 Reinfeld Tel: + 49 (0) 4533 2007 0 Fax: + 49 (0) 4533 2007 68 e-mail: mast@mast-diagnostica.de

Mast Diagnostic 12 Rue Jean-Jacques Mention, CS 91106, 80011 Amiens CEDEX 1 Tél. + 33 (0) 322 80 80 67 Fax + 33 (0) 322 80 99 22 e-mail: info@mast-diagnostic.fr



accuracy and quality as a science









Selectrol®: Manufactured under licence from Public Health England Culture Collections

SELECTROL[®] - FREEZE-DRIED ORGANISMS IN A DISC

Quality control of microbial characterisation tests, culture media and antimicrobial susceptibility determinations is best accomplished by the use of microorganisms with well-documented and stable phenotypic and genotypic characteristics.

Bacterial and fungal strains have been selected and recommended by expert bodies, such as EUCAST, CLSI and the European Pharmacopoeia, on the basis of their suitability for monitoring test performance and ensuring the validity of results for testing used in clinical, food, pharmaceutical, water and veterinary laboratories.

Products derived from the cultures in the collections should be manufactured using the minimum number of sub-cultures, to minimise the possibility of alterations to the phenotype due to mutations. See also page 14.

Selectrol strains are manufactured exclusively from Public Health England Culture Collections (NCTC[®] and NCPF[®]) and are first generation subcultures, unlike many products on the market which are 2nd, 3rd or 4th generation subcultures. They are preserved by long-term storage as freeze-dried cells in order to minimise any alterations to the phenotype caused by mutations.

Passages

A Selectrol[®] disc is a first generation subculture from a **master culture** sourced from Public Health England Culture Collections, and is designed to be used to obtain **working stock** cultures for use in testing. It is generally accepted that no more than a total of five passages should be made from the **master culture**, in order to avoid genetic drift and mutant selection. Therefore, no more than four passages (fresh cultures) from the **working stock** should be made.

Shelf life

For most strains, Selectrol[®] discs are guaranteed to contain at least 10⁶ organisms at the time of purchase; this number is sufficient to ensure that when the discs are used and stored as directed there will be viable organisms cultivable up to the stated end of the shelf life, which is usually 9 months from the time the vial is first opened.

Quality Control

Selectrol[®] batches are tested in our UKAS accredited testing laboratory number 2496. A test report for each batch of Selectrol[®] can be accessed via our website. The reporting of Selectrol[®] test results via the website comes under our UKAS accreditation.

Selectrol[®] cultures are rigorously tested to confirm identity, to confirm the possession of essential phenotypic characteristics and to exclude contamination with other organisms. Photographic evidence of the test results is retained for each batch, along with retained appropriately stored samples.



Glossary

AMRHAI: Antimicrobial Resistance and Healthcare Associated Infections reference unit

ATCC®: American Type Culture Collection. ATCC[®] strains are listed for reference only. ATCC[®] is a registered trademark of the American Type Culture Collection.

BSAC: British Society for Antimicrobial Chemotherapy - Now superseded by EUCAST

CLSI: Clinical Laboratory Standards Institute. (USA)

CPE: Carbapenemase Producing Enterobacteriaceae

CRE: Carbapenem Resistant Enterobacteriaceae

Culture collection: Cultures of fully characterised organisms maintained in such a way as to minimise sub-culturing. See page 14.

ESBL: Extended Spectrum Beta-Lactamase-producing organism.

EUCAST: European Committee on Antimicrobial Susceptibility Testing.

First generation derivative: A single passage from a master culture, for example a Selectrol® disc.

Master culture: Culture derived from a reference culture vial.

NCPF[®]: National Collection of Pathogenic Fungi. NCPF[®] is a registered trademark of Public Health England.

NCTC[®]: National Collection of Type Cultures. NCTC[®] is a registered trademark of Public Health England.

Passage: An equivalent term for a subculture.

PHE: Public Health England.

Reference cultures: Quality control strains selected on the basis of their phenotypic biochemical and antimicrobial susceptibility characteristics to be used as controls in microbiological testing. These are obtained as freeze-dried vials from culture collections.

Stock culture: Cultures derived from a Selectrol® disc, which can be stored for up to a week, usually on agar slants.

Working cultures: Stock cultures further sub-cultured to provide 18-24 hour growth for use in testing.

WDCM: World Data Centre for Microorganisms

WFCC: World Federation for Culture Collections



SIGNIFICANT PROPERTIES AND USES OF SELECTROL® ORGANISMS

Aspergillus brasiliensis (formerly Aspergillus niger):

MM94 – NCPF[®] 2275 / ATCC[®] 16404 / WDCM 00053 – used in pharmaceutical industry for testing media and preservatives. Colonies are initially white or yellowish and on the reverse greyish or greenish-yellow. Sporing heads on the colony surface are initially pale, becoming dark brown to black. Sporulation may be inhibited in sealed plates.

Bacillus cereus:

MM21 – NCTC[®] 10320 / ATCC[®] 9634 / WDCM 00001 (recently renamed *Bacillus toyonensis*) – ISO 11133 recommended media and ID test control organism.

MM86 - NCTC[®] 7464 / ATCC[®] 10876 - PHE recommended media and ID test control organism.

Bacillus subtilis (Bacillus subtilis subsp. spizizenii):

MM29 – NCTC[®] 10400 / ATCC[®] 6633 / WDCM 00003 – used in antibiotic assays (fully sensitive), PHE recommended media and ID test control organism.



Bacteroides fragilis:

MM44 – NCTC[®] 9343 / ATCC[®] 25285 – type strain, PHE recommended strain for media and sensitivity test control.

Campylobacter jejuni (Campylobacter jejuni subsp. jejuni):

MM82 - NCTC[®] 11322 / ATCC[®] 29428 / WDCM 00156 - PHE recommended strain for media control.

MM36 - NCTC® 11351 / ATCC® 33560 - EUCAST recommended strain for susceptibility testing.

Candida albicans:

MM28 - NCPF[®] 3255 / ATCC[®] 2091 / WDCM 00055 - sensitivity control / industrial use.

MM42 – NCPF[®] 3179 / ATCC[®] 10231 / WDCM 00054 – pharmaceutical / media testing / PHE recommended strain for media control.

CRE ≡ 'Carbapenem Resistant Enterobacteriaceae' / CPE ≡ 'Carbapenemase Producing Enterobacteriaceae'

There are 5 carbapenemases which are currently a significant problem in the UK – KPC, OXA-48, IMP, NDM and VIM – and PHE recommend that all clinically-significant Gram-negative bacteria should be routinely screened for carbapenemase production, using a recommended carbapenem² such as ertapenem or meropenem. Resistant isolates may be investigated further to determine which resistance mechanism is involved using the Modified Hodge Test, MALDI-TOF, PCR or a reference laboratory.

MM55 Klebsiella pneumoniae - NCTC® 13440 - produces a Class B VIM-1 Carbapenemase.

MM56 Klebsiella pneumoniae - NCTC® 13443 - produces a Class B NDM-1 Carbapenemase.

MM58 Klebsiella pneumoniae – NCTC[®] 13438 – produces a Class A KPC-3 Carbapenemase.

MM59 Klebsiella pneumoniae - NCTC® 13442 - produces a Class D OXA-48 Carbapenemase.

MM57 Escherichia coli - NCTC® 13476 - produces a Class B IMP Carbapenemase.

MM33 Escherichia coli - NCTC® 10418 / ATCC® 10536 - recommended by PHE as a negative control for CRE testing.



Citrobacter freundii:

MM27 - NCTC® 9750 / ATCC® 8090 - type strain.

Clostridium perfringens:

MM45 – NCTC[®] 8237 / ATCC[®] 13124 / WDCM 00007 – type strain. PHE recommended strain for food testing (Tryptose Sulphite Cycloserine agar – lactose and gelatin positive) and sensitivity test control. *Clostridium perfringens* is listed in Schedule 5 of the Anti-terrorism, Crime and Security Act 2001, and should be securely stored in accordance with the guidelines of the Act. However, MM45 is a type A strain, which <u>does not</u> produce the lethal epsilon toxin of potential interest to bioterrorists.

Clostridium sporogenes:

MM31 – NCTC[®] 532 / ATCC[®] 19404 / WDCM 00008 – used for media control. PHE recommended strain for media QC (lactose gelatin medium for ID of *C. perfringens* lactose negative and gelatin positive).

Enterobacter aerogenes:

MM26 - NCTC® 10006 / ATCC® 13048 / WDCM 00175 - type strain; used in water, paint and adhesive testing.

Enterobacter cloacae:

MM01 - NCTC® 13380 / ATCC® 23355 / WDCM 00082 - disinfectant control, media testing.

MM51- NCTC[®] 13406 - PHE recommended strain for QC of AmpC (de-repressed) detection.

Enterococcus faecalis:

MM52 – NCTC[®] 13379 / ATCC[®] 51299 / WDCM 00085 – is vancomycin resistant (low-level VanB mediated) and also shows highlevel resistance to aminoglycosides. It is used to confirm methodologies used to detect these resistances are working correctly. Lancefield group D.

MM17 – NCTC[®] 775 / ATCC[®] 19433 / WDCM 00009 – used in water industry and QC. PHE recommended strain for media control. Fully sensitive. Lancefield group D.

MM18 – NCTC[®] 12697 / ATCC[®] 29212 / WDCM 00087 – is fully sensitive to vancomycin and gentamicin. PHE recommended positive control strain for aesculin test. CLSI, EUCAST recommended media control for sulpha / trimethoprim testing and general susceptibility testing control. Lancefield group D.





Enterococcus hirae:

MM35 – NCTC[®] 13383 / ATCC[®] 10541 / WDCM 00011 – disinfectant control. Used in microbiological assays. Colonies are alphahaemolytic on sheep blood agar.

Escherichia coli strains:

MM02 – NCTC[®] 12241 / ATCC[®] 25922 / WDCM 00013 – EUCAST, CLSI, PHE recommended control strain for susceptibility testing (fully sensitive). Exhibits 2 colony types – the most prevalent type is slightly irregular, smooth and translucent. The secondary type appears more opaque. It is preferable to maintain cultures on agar as passage in broth can result in a change in MIC levels.



MM57 - NCTC[®] 13476 - CRE testing control; produces a Class B IMP Carbapenemase.

MM33 – NCTC[®] 10418 / ATCC[®] 10536 – (PHE recommended alternative to NCTC 12241) fully sensitive control strain. PHE recommended positive control for indole test, ONPG test, negative control for oxidase test, PHE recommended negative control for CRE and ESBL testing.

MM24 – NCTC[®] 11954 / ATCC[®] 35218 – beta-lactamase positive strain. CLSI recommended strain for susceptibility testing ONLY for penicillin / beta-lactamase inhibitor combinations. Sensitive to amoxicillin / clavulanic acid.

MM75 – NCTC[®] 9001 / ATCC[®] 11775 / WDCM 00090 – used in water / chemical industry. PHE recommended strain for media QC.

MM93 - NCTC® 12900 / ATCC® 700728 / WDCM 00014 - O157 strain (non-toxigenic). PHE recommended strain for media QC.

MM63 - NCTC® 11560 - beta-lactamase positive strain.

MM38 – NCTC[®] 12923 / ATCC[®] 8739 / WDCM 00012 – used in pharmaceutical / water industry. Three colony types: A) Entire, glistening, smooth and translucent. B) Entire, glistening smooth and opaque. C) Irregular, rough and translucent. The rough colonies appear after 48 hours incubation.

MM34 – NCTC[®] 13846 – Possesses the plasmid-mediated mcr-1 colistin resistance mechanism gene and is recommended by PHE and EUCAST as a control for tests to detect this increasingly prevalent resistance, in conjunction with NCTC® 12241 / ATCC® 25922 (Selectrol strain MM02) as a negative control.



Haemophilus influenzae strains:

MM81 - NCTC[®] 12699 / ATCC[®] 49247 – is a 'BLNAR' strain – (beta-lactamase non-producing ampicillin / amoxycillin resistant). These strains are important clinically because the susceptibility results obtained using conventional testing procedures maybe misleading in the case cephalosporins. PHE, CLSI recommended QC strain for susceptibility testing media.

MM98 – NCTC[®] 11931 – a fully sensitive strain. PHE recommended strain for porphyrin synthesis test, chocolate agar control.

MM100 – NCTC[®] 8468 / ATCC[®] 9334 / CCUG 23946 – another fully sensitive strain, which reportedly gives results which are easier to interpret when Mueller-Hinton medium is used in preference to Iso-Sensitest medium. MIC for amoxycillin is 0.5 mg/l.

MM37 - NCTC® 12975 / ATCC® 49766 - recommended by EUCAST.



Klebsiella strains:

MM04 *Klebsiella pneumoniae* – NCTC[®] 9633 / ATCC[®] 13883 / WDCM 00097 – type strain. Two colony types may be seen. The predominant type is entire and opaque. The secondary type is slightly smaller and translucent.

MM83 *Klebsiella pneumoniae* – NCTC[®] 13368 / ATCC[®] 700603 – ESBL-producing strain used as control for ESBL testing. There are two colony types.

MM55 Klebsiella pneumoniae – NCTC[®] 13440 – CRE testing control; produces a Class B VIM-1 Carbapenemase.



MM56 Klebsiella pneumoniae – NCTC® 13443 – CRE testing control; produces a Class B NDM-1 Carbapenemase.

MM58 Klebsiella pneumoniae – NCTC[®] 13438 – CRE testing control; produces a Class A KPC-3 Carbapenemase.

MM59 Klebsiella pneumoniae - NCTC® 13442 - CRE testing control; produces a Class D OXA-48 Carbapenemase.

MM88 *Klebsiella aerogenes (Raoultella planticola)* – NCTC[®] 9528 – used in water / pharmaceutical industry. PHE recommended negative control for Tryptone Bile X-Glucuronide agar and Yeast Extract agar.



Lactobacillus brevis:

MM76 - NCTC[®] 13386 / ATCC[®] 8287 - used in food industry.

Legionella pneumophila serogroup 1:

MM08 – NCTC[®] 11192 / ATCC[®] 33152 / WDCM 00107 – derived from strain isolated from first recognised outbreak of legionellosis in Philadelphia at the Legionnaires' Convention 1976

Listeria innocua:

MM92 - NCTC[®] 11288 / ATCC[®] 33090 / WDCM 00017 - type strain. Non-pathogenic.

Listeria monocytogenes:

MM87 – NCTC[®] 11994 / WDCM 00019 – type strain, PHE recommended positive control strain for Listeria detection in food. Serotype 4b, most common serovar isolated from human infections.

MM48 – NCTC[®] 7973 / ATCC[®] 35152 / WDCM 00109 – produces 2 phenotypes, one is beta-haemolytic and virulent, the other non-haemolytic and non-virulent. Serovar 1/2a.

MM77 – NCTC[®] 13372 / ATCC[®] 7644 – used in food microbiology Q.C. Colonies exhibit beta-haemolysis on sheep blood agar.

Neisseria gonorrhoeae:

 $MM96 - NCTC^{\circ}$ 12700 / ATCC $^{\circ}$ 49226 - has low-level, but clinically relevant, resistance to penicillin – MIC of penicillin is 0.5 mg/l. PHE recommended control for susceptibility testing – methodology assesses the ability of testing to detect resistance rather than sensitivity; this strain has low-level, but clinically relevant, resistance to penicillin – MIC of penicillin is 0.5 mg/l. Some variation in size and texture of colonies may be observed. Increased CO₂ is helpful in growth.

MM05 – NCTC[®] 8375 / ATCC[®] 19424 – is fully sensitive – MIC of penicillin is 0.06 mg/l. PHE recommended strain for media QC.

Proteus mirabilis:

MM43 – NCTC[®] 13376 / ATCC[®] 14153 – pharmaceutical / disinfectant / media control. MM68 – NCTC[®] 10975 – media control. PHE recommended control for motility test.



Proteus vulgaris:

MM09 – NCTC[®] 4175 / ATCC[®] 13315 – was the type strain, but is atypical and has been recognised as a separate species – *Proteus hauseri* – it is used for media control. Colonies are glistening with spreading edges.

Pseudomonas aeruginosa strains:

MM10 – NCTC[®] 12903 / ATCC[®] 27853 / WDCM 00025 – is fully sensitive to anti-pseudomonal antibiotics (EUCAST susceptibility test control). 2 colony types may be observed: A) predominantly flat, spreading edges and rough surface; B) small and compact. Produces both fluorescein and pyocyanin pigments.



MM65 - NCTC® 10662 / ATCC® 25668 / WDCM 00114 - is fully sensitive. PHE recommended control strain for media control

MM40 – NCTC[®] 12924 / ATCC[®] 9027 / WDCM 00026 – used in water industry / disinfectant testing. Colonies on agar plates are entire, glistening and mucoid with a grainy surface. This strain also produces both fluorescein and pyocyanin pigments.

MM41 – NCTC[®] 13359 / ATCC[®] 15442 – used in water industry / disinfectant testing. May produce up to 3 different colony types. Pyocyanin is not produced.

Rhodococcus equi:

MM97 - NCTC[®] 1621 / ATCC[®] 6939 / WDCM 00028 - type strain.

Saccharomyces cerevisiae:

MM73 – NCPF[®] 3178 – PHE recommended strain for food testing and enumeration of yeasts and moulds.

MM50 — NCTC® 10716 / WDCM 00058 - used for QC of culture media and for antifungal susceptibility testing.

Salmonella serotypes:

MM11 Salmonella Typhimurium – NCTC[®] 12023 / ATCC[®] 14028 / WDCM 00031 – (1,4,5,12: i: 1,2) Used for media/test QC. This is a common serotype from animals and from human infections.

The strains listed below are unusual serotypes, used to avoid any chance of confusion with strains commonly found in animals, food, etc, and are used to control media and detection methods in the food industry:

MM89 Salmonella Poona - NCTC[®] 4840 - (13,22: z: 1,6) PHE recommended control strain for food testing.

MM84 Salmonella Nottingham – NCTC® 7832 – (16: d: e,n,z15) PHE recommended control for water testing.

Serratia marcescens:

MM12 – NCTC[®] 13382 / ATCC[®] 8100 – used for disinfectant testing. PHE recommended negative control for indole test. Colonies are entire, glistening, smooth and translucent. Non-pigmented.



Staphylococcus aureus:

(A) Fully sensitive:

MM85 – NCTC[®] 6571 / ATCC[®] 9144 / WDCM 00035 – historically used for susceptibility testing ('Oxford staph'), but largely superseded by MM13 as it has unusually low MIC's and so is unrepresentative of normal range of Staph aureus strains. Sensitive to penicillin and cefoxitin / methicillin / oxacillin. PHE recommended coagulase, DNAse and catalase positive control.

MM13 – NCTC[®] 12981 / ATCC[®] 25923 / WDCM 00034 – used in susceptibility and media testing/QC. Fully sensitive to all antistaphylococcal antibiotics (including penicillin and methicillin / oxacillin). It is preferable to maintain cultures on agar as passage in broth can result in a change in MIC levels. Colonies are circular white to cream, convex to flat in elevation. After 48 hours incubation a few grey/translucent variants may be noted. Beta-haemolytic on sheep blood agar.

B) Penicillin resistant:

MM14 – NCTC[®] 12973 / ATCC[®] 29213 / WDCM 00131 – used for susceptibility testing, especially for automated methodology. EUCAST, CLSI strain. Sensitive to cefoxitin / methicillin / oxacillin. Penicillin resistant – weak beta-lactamase producer. Colonies are beta-haemolytic, and a golden-orange colour.

MM30 – NCTC[®] 7447 / ATCC[®] 6538P / WDCM 00033 – used for susceptibility testing/antibiotic assay, disinfectant testing. Cefoxitin / methicillin / oxacillin sensitive. Penicillin resistant. Colonies are weakly beta-haemolytic, coagulase positive and betalactamase negative.

(C) MRSA (cefoxitin / methicillin / oxacillin resistant):

MM91 – NCTC[®] 13373 / ATCC[®] 43300 / WDCM 00211 (MRSA) – Possesses mecA gene but is hetero-resistant, (so as few as one per thousand cells demonstrate the resistance) and consequently has low-level cefoxitin /oxacillin/methicillin resistance (4.0 mg/l MIC of oxacillin, 8.0 mg/l MIC of cefoxitin – methicillin sensitive strains have MIC of 0.12-0.5 for oxacillin and 1-4 for cefoxitin.); it is used to confirm testing procedures for methicillin resistance are working and provides a more stringent test than testing with an MRSA which shows homogeneous resistance and has a much higher MIC. This organism will have a zone of inhibition reduced in size compared to a fully cefoxitin / oxacillin / methicillin sensitive strain (such as MM13). CLSI recommended strain for MRSA testing. There are two colony types: 1) Beta-haemolytic with a slight yellow tint. 2) Non-haemolytic and white.

MM64 – NCTC[®] 12493 / WDCM 00212 (MRSA) – possesses mecA gene and shows homogeneous resistance with MIC of >64 for methicillin, which produces high-level cefoxitin / methicillin / oxacillin resistance. EUCAST recommended strain. Instances have been reported where loss of the mecA gene has occurred during storage.

D) Other:

MM46 – NCTC[®] 10788 / ATCC[®] 6538 / WDCM 00032 – used in pharmaceutical industry for testing disinfectants etc. Usually yellow pigmented colonies, or can produce a white colonial variant. Beta-haemolytic.





Staphylococcus epidermidis:

MM15 - NCTC® 13360 / ATCC® 12228 / WDCM 00036 - used for media control / antibiotic assay. Colonies are small and betahaemolytic.

Streptococcus agalactiae: (Beta-haemolytic Streptococcus group B)

MM16 - NCTC® 8181 / ATCC® 13813 - type strain, used for QC. PHE recommended negative control for aesculin test.

Streptococcus pneumoniae strains:

MM95 – NCTC[®] 12977 / ATCC[®] 49619 – has low-level, but clinically relevant, resistance to penicillin – this organism is used to assess detection of resistance rather than sensitivity. PHE recommended positive control for bile solubility test. CLSI, EUCAST recommended control strain for susceptibility testing. Serotype 19F.

MM19 – NCTC[®] 12695 / ATCC[®] 6303 – is fully sensitive. Colonies are mucoid and alpha-haemolytic. A few colonies may have an irregular edge. Serotype 3.



Streptococcus pyogenes:

MM20 – NCTC[®] 12696 / ATCC[®] 19615 – used for QC and media testing. Lancefield group A, beta-haemolytic. PHE recommended blood agar control.

Vibrio parahaemolyticus:

MM06 – NCTC[®] 10885 / WDCM 00185 – used for QC of media and ID testing. PHE recommended strain used mainly in the food industry.

Yersinia enterocolitica:

MM80 - NCTC[®] 12982 / ATCC[®] 9610 / WDCM 00038 - type strain, used for media control. Serotype O:8, which is a pathogenic serotype, commonest in USA.

References:

- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Routine and Extended Internal Quality Control for MIC Determination and Disc Diffusion. Version 7.0 - 01.01.2017.
- 2 UK Standards for Microbiology Investigations. Example Reference Strains for Microbiology Investigations Test Procedures: Bacteriology—Test Procedures | TP 1 | Issue No. 2 | 05.01.2015. Public Health England (PHE).
- 3 Performance Standards for Antimicrobial Disc Susceptibility Tests: Approved Standard—11th Edition. Clinical and Laboratory Standards Institute (CLSI).



How to use Selectrol®

Always warm the vial to ambient temperature before opening.

Be sure to use non-selective culture media to revive the organisms.

For the more fastidious organisms, such as anaerobes, it is generally better to use agar rather than broth for revival.



Place disc on suitable growth medium such as blood agar



Leave disc for a few minutes to liquefy, then spread plate and incubate to produce isolated colonies





Place disc in a small volume of a suitable broth medium such as brain-heart infusion



Allow disc a few minutes to dissolve, then spread aliquot onto a plate of suitable growth medium

Obtain a stock culture which can be used to prepare an inoculum for biochemical and antibiotic susceptibility tests



Out-of-specification results

Laboratories use Selectrol[®] for Quality Control of culture media, biochemical identification tests and antimicrobial susceptibility testing. When a laboratory test result, an MIC or biochemical reaction, is unexpected or out-of-specification, the test should first be repeated to confirm it; an out-of-specification result is an indication that the testing procedure should be reviewed; it is not, in the first instance, a sign of a problem with the control organism.

If incorrect results are obtained on retesting, the explanation could be:

- The test procedure was not followed correctly check standard operating procedures
- There is an instrumentation error check calibration, mechanical functioning, etc
- There is a problem with the consumables out of date, incorrect storage, etc
- The culture of the control organism has become contaminated

Technical Support

If no explanation for out-of-spec results can be found, but repeated tests still give unacceptable results, please contact TCS and / or your relevant reference laboratory or instrument manufacturer for advice. For example, contact AMRHAI at Colindale, London if MIC results are consistently outside the acceptable range. Please retain any remaining discs of organisms about which you have concerns so they can be returned to TCS and investigated alongside retained samples.





Preparing QC and Validation Spikes from Selectrol®

Preparing the spike

- Place a Selectrol[®] disc in Brain Heart Infusion (BHI) broth* or equivalent, and culture (typically for 18 hours) at the appropriate temperature for the organism (typically 37°C)

- Assume the count in the broth to be 10⁸ organisms per ml ------ (A)
- Mix and transfer 100 μl of (A) to 100 ml of saline or 1/4 strength Ringer's solution -- (B)
- Mix and transfer 100 μl of (B) to 10 ml of saline or 1/4 strength Ringer's solution --- (C)
- Mix and transfer 100 µl of (C) to your homogenised food sample.

Verifying the inoculum

- Pipette 5 x 10 µl drops from (C) onto each of two agar plates for Miles and Misra counts.

Using the assumptions and dilutions above:

- (A) contains 10⁸ organisms per ml
- (B) contains 10⁵ organisms per ml
- (C) contains 10³ organisms per ml

If the Miles and Misra counts indicate that the required count was not achieved:

- If the count was too high by a factor of 10, reduce the volume transferred from (A) to (B) from 100 µl to 10 µl
- If the count was too low by a factor of 10, increase the volume transferred from (A) to (B) from 100 µl to 1 ml.

Keep a record of the correct dilutions for each organism type for future use. You will find that this method is very repeatable.

*Note: BHI broth will work for most of the Selectrol[®] organisms; however, for fastidious organisms an appropriate culture broth must be selected, e.g. Fastidious Anaerobe Broth for strictly anaerobic organisms.





Culture Collections

Cultures of microorganisms have been deposited and subsequently maintained in 589 collections in 68 countries, and many of the cultures are derived from the same original isolate; the history of each organism, its properties and names of the culture collections which hold it are detailed in the relevant catalogues and websites.

Some of the organisms have been selected and recommended by expert organisations to be supplied as controls for microbiological tests, and when the identical cultures are present in more than one collection they will have a specific designation for each, incorporating the abbreviation for the collection and a reference number.

For example:- *Staphylococcus aureus* NCTC 7447, widely recommended as a control for antimicrobial susceptibility testing, is held in 30 collections, and consequently the phenotypically and genotypically identical organism has 30 different references, such as ATCC 6538P, CIP 53.156, DSM 346 and so on.

In an effort to minimise potential confusion and help users find local sources of reference strains, the WFCC and the WDCM initiated a system that ascribes each recommended QC strain a reference number (WDCM 00001 onwards), cites all collections that contain it and provides contact details and each collection's unique reference. For example, the strain of *Staphylococcus aureus* NCTC 7447 (Selectrol[®] strain MM33) mentioned above is designated WDCM 00033.

Staphylococcus aureus WDCM 00033

AHU 1142; ATCC[™] 6538P; BCRC 10451; BTCC 209P; BU 395; CCM 2022; CCTM 596; CCUG 1828; CECT 240; CIP 53.156; CN 3784; CNCTC Mau 28/58; DSM 346; FIRDI 451; IAM 1011; IAM 12082; IEM Mau 28/58; IFO 12732; IFO 3061; IID 671; IMET 10904; JCM 2151; LMG 8195; NCIMB 8625; NCTC 7447; NRRL B-313; OUT 8232; PCI 1209; PZH 8/54; RIMD 3109007; VNIIA 209P;

Products derived from the cultures in the collections should be manufactured using the minimum number of sub-cultures, to minimise the possibility of alterations to the phenotype due to mutations. Ideally, as in the case of **Selectrol®**, a single sub-culture only is used, so the **Selectrol®** product is a 'first generation derivative' of a culture supplied by NCTC, and will be identical with regard to its properties and suitability for use in QC applications to a culture of the particular organism obtained from any of the other WDCM listed culture collections.

Every effort has been made to ensure the accuracy of the information in this document, however TCS makes no warranties, expressed or implied, regarding errors or omissions and assumes no legal liability or responsibility for loss or damage resulting from the use of information contained within.

Selectrol Strain Index

| Strain Name | Designation | Code | WDCM |
|---------------------------------------|--|-------|-------|
| Aspergillus brasiliensis | NCPF [®] 2275 / ATCC [®] 16404 | MM94 | 00053 |
| Bacillus cereus | NCTC [®] 10320 / ATCC [®] 9634 | MM21 | 00001 |
| Bacillus cereus | NCTC [®] 7464 / ATCC [®] 10876 | MM86 | |
| Bacillus subtilis | NCTC [®] 10400 / ATCC [®] 6633 | MM29 | 00003 |
| Bacteroides fragilis | NCTC [®] 9343 / ATCC [®] 25285 | MM44 | |
| Campylobacter jejuni | NCTC [®] 11351 / ATCC [®] 33560 | MM36 | |
| Campylobacter jejuni | NCTC [®] 11322 / ATCC [®] 29428 | MM82 | 00156 |
| Candida albicans | NCPF [®] 3255 / ATCC [®] 2091 | MM28 | 00055 |
| Candida albicans | NCPF [®] 3179 / ATCC [®] 10231 | MM42 | 00054 |
| Citrobacter freundii | NCTC [®] 9750 / ATCC [®] 8090 | MM27 | |
| Clostridium perfringens | NCTC [®] 8237 / ATCC [®] 13124 | MM45 | 00007 |
| Clostridium sporogenes | NCTC [®] 532 / ATCC [®] 19404 | MM31 | 00008 |
| Enterobacter aerogenes | NCTC [®] 10006 / ATCC [®] 13048 | MM26 | 00175 |
| Enterobacter cloacae | NCTC [®] 13380 / ATCC [®] 23355 | MM01 | 00082 |
| Enterobacter cloacae | NCTC [®] 13406 | MM51 | |
| Enterococcus faecalis | NCTC [®] 775 / ATCC [®] 19433 | MM17 | 00009 |
| Enterococcus faecalis | NCTC [®] 12697 / ATCC [®] 29212 | MM18 | 00087 |
| Enterococcus faecalis | NCTC [®] 13379 / ATCC [®] 51299 | MM52 | 00085 |
| Enterococcus hirae | NCTC [®] 13383 /ATCC [®] 10541 | MM35 | 00011 |
| Escherichia coli | NCTC [®] 12241 / ATCC [®] 25922 | MM02 | 00013 |
| Escherichia coli | NCTC [®] 11954 / ATCC [®] 35218 | MM24 | |
| Escherichia coli | NCTC [°] 10418 / ATCC [°] 10536 | MM33 | |
| Escherichia coli | NCTC [®] 12923 / ATCC [®] 8739 | MM38 | 00012 |
| Escherichia coli | NCTC [®] 11560 | MM63 | |
| Escherichia coli | NCTC [®] 9001 / ATCC [®] 11775 | MM75 | 00090 |
| Escherichia coli CRE | NCTC [®] 13476 | MM57 | |
| Escherichia coli (mcr-1) | NCTC [®] 13846 | MM34 | |
| Escherichia coli O157 (non-toxigenic) | NCTC [®] 12900 / ATCC [®] 700728 | MM93 | 00014 |
| Haemophilus influenzae | NCTC [®] 8468 / ATCC [®] 9334 | MM100 | |
| Haemophilus influenzae | NCTC [®] 12975 / ATCC [®] 49766 | MM37 | |
| Haemophilus influenzae | NCTC [®] 12699 / ATCC [®] 49247 | MM81 | |
| Haemophilus influenzae | NCTC [®] 11931 | MM98 | |
| Klebsiella aerogenes | NCTC [®] 9528 | MM88 | |
| Klebsiella pneumoniae | NCTC [®] 9633 / ATCC [®] 13883 | MM04 | 00097 |
| Klebsiella pneumoniae | NCTC [®] 13368 / ATCC [®] 700603 | MM83 | |
| Klebsiella pneumoniae CRE | NCTC [®] 13440 | MM55 | |
| Klebsiella pneumoniae CRE | NCTC [®] 13443 | MM56 | |
| Klebsiella pneumoniae CRE | NCTC [®] 13438 | MM58 | |

Selectrol Strain Index

| Strain Name | Designation | Code | WDCM |
|------------------------------------|---|------|-------|
| Klebsiella pneumoniae CRE | NCTC [®] 13442 | MM59 | |
| Lactobacillus brevis | NCTC [®] 13386 / ATCC [®] 8287 | MM76 | |
| Legionella pneumophila serogroup 1 | NCTC [®] 11192 / ATCC [®] 33152 | MM08 | 00107 |
| Listeria innocua | NCTC [®] 11288 / ATCC [®] 33090 | MM92 | 00017 |
| Listeria monocytogenes | NCTC [®] 7973 / ATCC [®] 35152 | MM48 | 00109 |
| Listeria monocytogenes | NCTC [®] 13372 ATCC [®] 7644 | MM77 | |
| Listeria monocytogenes | NCTC [®] 11994 | MM87 | 00019 |
| Neisseria gonorrhoeae | NCTC [®] 8375 / ATCC [®] 19424 | MM05 | |
| Neisseria gonorrhoeae | NCTC [®] 12700 / ATCC [®] 49226 | MM96 | |
| Proteus mirabilis | NCTC [®] 13376 / ATCC [®] 14153 | MM43 | |
| Proteus mirabilis | NCTC [®] 10975 | MM68 | |
| Proteus vulgaris | NCTC [®] 4175 / ATCC [®] 13315 | MM09 | |
| Pseudomonas aeruginosa | NCTC [®] 12903 / ATCC [®] 27853 | MM10 | 00025 |
| Pseudomonas aeruginosa | NCTC [®] 12924 / ATCC [®] 9027 | MM40 | 00026 |
| Pseudomonas aeruginosa | NCTC [®] 13359 / ATCC [®] 15442 | MM41 | |
| Pseudomonas aeruginosa | NCTC [®] 10662 / ATCC [®] 25668 | MM65 | 00114 |
| Rhodococcus equi | NCTC [®] 1621 / ATCC [®] 6939 | MM97 | 00028 |
| Saccharomyces cerevisiae | NCTC [®] 10716/ ATCC [®] 9763 | MM50 | 00058 |
| Saccharomyces cerevisiae | NCPF [®] 3178 | MM73 | 1 |
| Salmonella Nottingham | NCTC [®] 7832 | MM84 | |
| Salmonella Poona | NCTC [®] 4840 | MM89 | |
| Salmonella Typhimurium | NCTC [®] 12023/ ATCC [®] 14028 | MM11 | 00031 |
| Serratia marcescens | NCTC [°] 13382 / ATCC [°] 8100 | MM12 | |
| Staphylococcus aureus | NCTC [°] 12981 / ATCC [°] 25923 | MM13 | 00034 |
| Staphylococcus aureus | NCTC [®] 12973 / ATCC [®] 29213 | MM14 | 00131 |
| Staphylococcus aureus | NCTC [®] 7447 / ATCC [®] 6538P | MM30 | 00033 |
| Staphylococcus aureus | NCTC [®] 10788 / ATCC [®] 6538 | MM46 | 00032 |
| Staphylococcus aureus | NCTC [®] 6571 / ATCC [®] 9144 | MM85 | 00035 |
| Staphylococcus aureus (MRSA) | NCTC [®] 12493 | MM64 | 00212 |
| Staphylococcus aureus (MRSA) | NCTC [®] 13373 / ATCC [®] 43300 | MM91 | 00211 |
| Staphylococcus epidermidis | NCTC [®] 13360 / ATCC [®] 12228 | MM15 | 00036 |
| Streptococcus agalactiae | NCTC [®] 8181 / ATCC [®] 13813 | MM16 | |
| Streptococcus pneumoniae | NCTC [®] 12695 /ATCC [®] 6303 | MM19 | |
| Streptococcus pneumoniae | NCTC [®] 12977 /ATCC [®] 49619 | MM95 | |
| Streptococcus pyogenes | NCTC [®] 12696 /ATCC [®] 19615 | MM20 | |
| Vibrio parahaemolyticus | NCTC [®] 10885 | MM06 | 00185 |
| Yersinia enterocolitica | NCTC [®] 12982 / ATCC [®] 9610 | MM80 | 00038 |

Selectrol Strains Listed by WDCM Number

| WDCM | Strain Name | Designation | Code |
|-------|---------------------------------------|--|------|
| 00001 | Bacillus cereus | NCTC [®] 10320 / ATCC [®] 9634 | MM21 |
| 00003 | Bacillus subtilis | NCTC [®] 10400 / ATCC [®] 6633 | MM29 |
| 00007 | Clostridium perfringens | NCTC [®] 8237 / ATCC [®] 13124 | MM45 |
| 00008 | Clostridium sporogenes | NCTC [®] 532 / ATCC [®] 19404 | MM31 |
| 00009 | Enterococcus faecalis | NCTC [®] 775 / ATCC [®] 19433 | MM17 |
| 00011 | Enterococcus hirae | NCTC [®] 13383 /ATCC [®] 10541 | MM35 |
| 00012 | Escherichia coli | NCTC [®] 12923 / ATCC [®] 8739 | MM38 |
| 00013 | Escherichia coli | NCTC [®] 12241 / ATCC [®] 25922 | MM02 |
| 00014 | Escherichia coli O157 (non-toxigenic) | NCTC [®] 12900 / ATCC [®] 700728 | MM93 |
| 00017 | Listeria innocua | NCTC [®] 11288 / ATCC [®] 33090 | MM92 |
| 00019 | Listeria monocytogenes | NCTC [®] 11994 | MM87 |
| 00025 | Pseudomonas aeruginosa | NCTC [®] 12903 / ATCC [®] 27853 | MM10 |
| 00026 | Pseudomonas aeruginosa | NCTC [®] 12924 / ATCC [®] 9027 | MM40 |
| 00028 | Rhodococcus equi | NCTC [®] 1621 / ATCC [®] 6939 | MM97 |
| 00031 | Salmonella Typhimurium | NCTC [®] 12023/ ATCC [®] 14028 | MM11 |
| 00032 | Staphylococcus aureus | NCTC [®] 10788 / ATCC [®] 6538 | MM46 |
| 00033 | Staphylococcus aureus | NCTC [®] 7447 / ATCC [®] 6538P | MM30 |
| 00034 | Staphylococcus aureus | NCTC [®] 12981 / ATCC [®] 25923 | MM13 |
| 00035 | Staphylococcus aureus | NCTC [®] 6571 / ATCC [®] 9144 | MM85 |
| 00036 | Staphylococcus epidermidis | NCTC [®] 13360 / ATCC [®] 12228 | MM15 |
| 00038 | Yersinia enterocolitica | NCTC [®] 12982 / ATCC [®] 9610 | MM80 |
| 00053 | Aspergillus brasiliensis | NCPF [®] 2275 / ATCC [®] 16404 | MM94 |
| 00054 | Candida albicans | NCPF [®] 3179 / ATCC [®] 10231 | MM42 |
| 00055 | Candida albicans | NCPF [®] 3255 / ATCC [®] 2091 | MM28 |
| 00058 | Saccharomyces cerevisiae | NCTC [®] 10716/ ATCC [®] 9763 | MM50 |
| 00082 | Enterobacter cloacae | NCTC [®] 13380 / ATCC [®] 23355 | MM01 |
| 00085 | Enterococcus faecalis | NCTC [®] 13379 / ATCC [®] 51299 | MM52 |
| 00087 | Enterococcus faecalis | NCTC [®] 12697 / ATCC [®] 29212 | MM18 |
| 00090 | Escherichia coli | NCTC [®] 9001 / ATCC [®] 11775 | MM75 |
| 00097 | Klebsiella pneumoniae | NCTC [®] 9633 / ATCC [®] 13883 | MM04 |
| 00107 | Legionella pneumophila serogroup 1 | NCTC [®] 11192 / ATCC [®] 33152 | MM08 |
| 00109 | Listeria monocytogenes | NCTC [®] 7973 / ATCC [®] 35152 | MM48 |
| 00114 | Pseudomonas aeruginosa | NCTC [®] 10662 / ATCC [®] 25668 | MM65 |
| 00131 | Staphylococcus aureus | NCTC [®] 12973 / ATCC [®] 29213 | MM14 |
| 00156 | Campylobacter jejuni | NCTC [®] 11322 / ATCC [®] 29428 | MM82 |
| 00175 | Enterobacter aerogenes | NCTC [®] 10006 / ATCC [®] 13048 | MM26 |
| 00185 | Vibrio parahaemolyticus | NCTC [®] 10885 | MM06 |
| 00211 | Staphylococcus aureus (MRSA) | NCTC [®] 13373 / ATCC [®] 43300 | MM91 |
| 00212 | Staphylococcus aureus (MRSA) | NCTC [®] 12493 | MM64 |

Notes




www.tcsbiosciences.co.uk

TCS Biosciences Ltd Botolph Claydon, Buckingham, MK18 2LR, United Kingdom t: +44 (0) 1296 714222, f: +44 (0) 1296 714806, e: sales@tcsgroup.co.uk

accuracy and quality as a science









Selectrol® : Manufactured under licence from Public Health England Culture Collections

SELECTROL® - FREEZE-DRIED ORGANISMS IN A DISC

Quality control of microbial characterisation tests, culture media and antimicrobial susceptibility determinations is best accomplished by the use of microorganisms with well-documented and stable phenotypic and genotypic characteristics.

Bacterial and fungal strains have been selected and recommended by expert bodies, such as EUCAST, CLSI and the European Pharmacopoeia, on the basis of their suitability for monitoring test performance and ensuring the validity of results for testing used in clinical, food, pharmaceutical, water and veterinary laboratories.

Products derived from the cultures in the collections should be manufactured using the minimum number of sub-cultures, to minimise the possibility of alterations to the phenotype due to mutations. See also page 14.

Selectrol strains are manufactured exclusively from Public Health England Culture Collections (NCTC[®] and NCPF[®]) and are first generation subcultures, unlike many products on the market which are 2nd, 3rd or 4th generation subcultures. They are preserved by long-term storage as freeze-dried cells in order to minimise any alterations to the phenotype caused by mutations.

Passages

A Selectrol[®] disc is a first generation subculture from a **master culture** sourced from Public Health England Culture Collections, and is designed to be used to obtain **working stock** cultures for use in testing. It is generally accepted that no more than a total of five passages should be made from the **master culture**, in order to avoid genetic drift and mutant selection. Therefore, no more than four passages (fresh cultures) from the **working stock** should be made.

Shelf life

For most strains, Selectrol[®] discs are guaranteed to contain at least 10⁶ organisms at the time of purchase; this number is sufficient to ensure that when the discs are used and stored as directed there will be viable organisms cultivable up to the stated end of the shelf life, which is usually 9 months from the time the vial is first opened.

Quality Control

Selectrol[®] batches are tested in our UKAS accredited testing laboratory number 2496. A test report for each batch of Selectrol[®] can be accessed via our website. The reporting of Selectrol[®] test results via the website comes under our UKAS accreditation.

Selectrol[®] cultures are rigorously tested to confirm identity, to confirm the possession of essential phenotypic characteristics and to exclude contamination with other organisms. Photographic evidence of the test results is retained for each batch, along with retained appropriately stored samples.



Glossary

AMRHAI: Antimicrobial Resistance and Healthcare Associated Infections reference unit

ATCC®: American Type Culture Collection. ATCC[®] strains are listed for reference only. ATCC[®] is a registered trademark of the American Type Culture Collection.

BSAC: British Society for Antimicrobial Chemotherapy - Now superseded by EUCAST

CLSI: Clinical Laboratory Standards Institute. (USA)

CPE: Carbapenemase Producing Enterobacteriaceae

CRE: Carbapenem Resistant Enterobacteriaceae

Culture collection: Cultures of fully characterised organisms maintained in such a way as to minimise sub-culturing. See page 14.

ESBL: Extended Spectrum Beta-Lactamase-producing organism.

EUCAST: European Committee on Antimicrobial Susceptibility Testing.

First generation derivative: A single passage from a master culture, for example a Selectrol® disc.

Master culture: Culture derived from a reference culture vial.

NCPF[®]: National Collection of Pathogenic Fungi. NCPF[®] is a registered trademark of Public Health England.

NCTC[®]: National Collection of Type Cultures. NCTC[®] is a registered trademark of Public Health England.

Passage: An equivalent term for a subculture.

PHE: Public Health England.

Reference cultures: Quality control strains selected on the basis of their phenotypic biochemical and antimicrobial susceptibility characteristics to be used as controls in microbiological testing. These are obtained as freeze-dried vials from culture collections.

Stock culture: Cultures derived from a Selectrol® disc, which can be stored for up to a week, usually on agar slants.

Working cultures: Stock cultures further sub-cultured to provide 18-24 hour growth for use in testing.

WDCM: World Data Centre for Microorganisms

WFCC: World Federation for Culture Collections



Culture Collections

Cultures of microorganisms have been deposited and subsequently maintained in 589 collections in 68 countries, and many of the cultures are derived from the same original isolate; the history of each organism, its properties and names of the culture collections which hold it are detailed in the relevant catalogues and websites.

Some of the organisms have been selected and recommended by expert organisations to be supplied as controls for microbiological tests, and when the identical cultures are present in more than one collection they will have a specific designation for each, incorporating the abbreviation for the collection and a reference number.

For example:- *Staphylococcus aureus* NCTC 7447, widely recommended as a control for antimicrobial susceptibility testing, is held in 30 collections, and consequently the phenotypically and genotypically identical organism has 30 different references, such as ATCC 6538P, CIP 53.156, DSM 346 and so on.

In an effort to minimise potential confusion and help users find local sources of reference strains, the WFCC and the WDCM initiated a system that ascribes each recommended QC strain a reference number (WDCM 00001 onwards), cites all collections that contain it and provides contact details and each collection's unique reference. For example, the strain of *Staphylococcus aureus* NCTC 7447 (Selectrol[®] strain MM33) mentioned above is designated WDCM 00033.

Staphylococcus aureus WDCM 00033

AHU 1142; ATCC[™] 6538P; BCRC 10451; BTCC 209P; BU 395; CCM 2022; CCTM 596; CCUG 1828; CECT 240; CIP 53.156; CN 3784; CNCTC Mau 28/58; DSM 346; FIRDI 451; IAM 1011; IAM 12082; IEM Mau 28/58; IFO 12732; IFO 3061; IID 671; IMET 10904; JCM 2151; LMG 8195; NCIMB 8625; NCTC 7447; NRRL B-313; OUT 8232; PCI 1209; PZH 8/54; RIMD 3109007; VNIIA 209P;

Products derived from the cultures in the collections should be manufactured using the minimum number of sub-cultures, to minimise the possibility of alterations to the phenotype due to mutations. Ideally, as in the case of **Selectrol®**, a single sub-culture only is used, so the **Selectrol®** product is a 'first generation derivative' of a culture supplied by NCTC, and will be identical with regard to its properties and suitability for use in QC applications to a culture of the particular organism obtained from any of the other WDCM listed culture collections.

Every effort has been made to ensure the accuracy of the information in this document, however TCS makes no warranties, expressed or implied, regarding errors or omissions and assumes no legal liability or responsibility for loss or damage resulting from the use of information contained within.

Selectrol Strain Index

| Strain Name | Designation | Code | WDCM |
|---------------------------------------|--|-------|-------|
| Aspergillus brasiliensis | NCPF [®] 2275 / ATCC [®] 16404 | MM94 | 00053 |
| Bacillus cereus | NCTC [®] 10320 / ATCC [®] 9634 | MM21 | 00001 |
| Bacillus cereus | NCTC [®] 7464 / ATCC [®] 10876 | MM86 | |
| Bacillus subtilis | NCTC [®] 10400 / ATCC [®] 6633 | MM29 | 00003 |
| Bacteroides fragilis | NCTC [®] 9343 / ATCC [®] 25285 | MM44 | |
| Campylobacter jejuni | NCTC [®] 11351 / ATCC [®] 33560 | MM36 | |
| Campylobacter jejuni | NCTC [®] 11322 / ATCC [®] 29428 | MM82 | 00156 |
| Candida albicans | NCPF [®] 3255 / ATCC [®] 2091 | MM28 | 00055 |
| Candida albicans | NCPF [®] 3179 / ATCC [®] 10231 | MM42 | 00054 |
| Citrobacter freundii | NCTC [®] 9750 / ATCC [®] 8090 | MM27 | |
| Clostridium perfringens | NCTC [®] 8237 / ATCC [®] 13124 | MM45 | 00007 |
| Clostridium sporogenes | NCTC [®] 532 / ATCC [®] 19404 | MM31 | 00008 |
| Enterobacter aerogenes | NCTC [®] 10006 / ATCC [®] 13048 | MM26 | 00175 |
| Enterobacter cloacae | NCTC [®] 13380 / ATCC [®] 23355 | MM01 | 00082 |
| Enterobacter cloacae | NCTC [®] 13406 | MM51 | |
| Enterococcus faecalis | NCTC [®] 775 / ATCC [®] 19433 | MM17 | 00009 |
| Enterococcus faecalis | NCTC [®] 12697 / ATCC [®] 29212 | MM18 | 00087 |
| Enterococcus faecalis | NCTC [®] 13379 / ATCC [®] 51299 | MM52 | 00085 |
| Enterococcus hirae | NCTC [®] 13383 /ATCC [®] 10541 | MM35 | 00011 |
| Escherichia coli | NCTC [®] 12241 / ATCC [®] 25922 | MM02 | 00013 |
| Escherichia coli | NCTC [®] 11954 / ATCC [®] 35218 | MM24 | |
| Escherichia coli | NCTC [°] 10418 / ATCC [°] 10536 | MM33 | |
| Escherichia coli | NCTC [®] 12923 / ATCC [®] 8739 | MM38 | 00012 |
| Escherichia coli | NCTC [®] 11560 | MM63 | |
| Escherichia coli | NCTC [®] 9001 / ATCC [®] 11775 | MM75 | 00090 |
| Escherichia coli CRE | NCTC [®] 13476 | MM57 | |
| Escherichia coli (mcr-1) | NCTC [®] 13846 | MM34 | |
| Escherichia coli O157 (non-toxigenic) | NCTC [®] 12900 / ATCC [®] 700728 | MM93 | 00014 |
| Haemophilus influenzae | NCTC [®] 8468 / ATCC [®] 9334 | MM100 | |
| Haemophilus influenzae | NCTC [®] 12975 / ATCC [®] 49766 | MM37 | |
| Haemophilus influenzae | NCTC [®] 12699 / ATCC [®] 49247 | MM81 | |
| Haemophilus influenzae | NCTC [®] 11931 | MM98 | |
| Klebsiella aerogenes | NCTC [®] 9528 | MM88 | |
| Klebsiella pneumoniae | NCTC [®] 9633 / ATCC [®] 13883 | MM04 | 00097 |
| Klebsiella pneumoniae | NCTC [®] 13368 / ATCC [®] 700603 | MM83 | |
| Klebsiella pneumoniae CRE | NCTC [®] 13440 | MM55 | |
| Klebsiella pneumoniae CRE | NCTC [®] 13443 | MM56 | |
| Klebsiella pneumoniae CRE | NCTC [®] 13438 | MM58 | |

Selectrol Strain Index

| Strain Name | Designation | Code | WDCM |
|------------------------------------|---|------|-------|
| Klebsiella pneumoniae CRE | NCTC [®] 13442 | MM59 | |
| Lactobacillus brevis | NCTC [®] 13386 / ATCC [®] 8287 | MM76 | |
| Legionella pneumophila serogroup 1 | NCTC [®] 11192 / ATCC [®] 33152 | MM08 | 00107 |
| Listeria innocua | NCTC [®] 11288 / ATCC [®] 33090 | MM92 | 00017 |
| Listeria monocytogenes | NCTC [®] 7973 / ATCC [®] 35152 | MM48 | 00109 |
| Listeria monocytogenes | NCTC [®] 13372 ATCC [®] 7644 | MM77 | |
| Listeria monocytogenes | NCTC [®] 11994 | MM87 | 00019 |
| Neisseria gonorrhoeae | NCTC [®] 8375 / ATCC [®] 19424 | MM05 | - |
| Neisseria gonorrhoeae | NCTC [®] 12700 / ATCC [®] 49226 | MM96 | |
| Proteus mirabilis | NCTC [®] 13376 / ATCC [®] 14153 | MM43 | |
| Proteus mirabilis | NCTC [®] 10975 | MM68 | |
| Proteus vulgaris | NCTC [®] 4175 / ATCC [®] 13315 | MM09 | |
| Pseudomonas aeruginosa | NCTC [®] 12903 / ATCC [®] 27853 | MM10 | 00025 |
| Pseudomonas aeruginosa | NCTC [®] 12924 / ATCC [®] 9027 | MM40 | 00026 |
| Pseudomonas aeruginosa | NCTC [®] 13359 / ATCC [®] 15442 | MM41 | |
| Pseudomonas aeruginosa | NCTC [®] 10662 / ATCC [®] 25668 | MM65 | 00114 |
| Rhodococcus equi | NCTC [®] 1621 / ATCC [®] 6939 | MM97 | 00028 |
| Saccharomyces cerevisiae | NCTC [®] 10716/ ATCC [®] 9763 | MM50 | 00058 |
| Saccharomyces cerevisiae | NCPF [®] 3178 | MM73 | / / C |
| Salmonella Nottingham | NCTC [®] 7832 | MM84 | |
| Salmonella Poona | NCTC [®] 4840 | MM89 | |
| Salmonella Typhimurium | NCTC [®] 12023/ ATCC [®] 14028 | MM11 | 00031 |
| Serratia marcescens | NCTC [®] 13382 / ATCC [®] 8100 | MM12 | |
| Staphylococcus aureus | NCTC [®] 12981 / ATCC [®] 25923 | MM13 | 00034 |
| Staphylococcus aureus | NCTC [®] 12973 / ATCC [®] 29213 | MM14 | 00131 |
| Staphylococcus aureus | NCTC [®] 7447 / ATCC [®] 6538P | MM30 | 00033 |
| Staphylococcus aureus | NCTC [®] 10788 / ATCC [®] 6538 | MM46 | 00032 |
| Staphylococcus aureus | NCTC [®] 6571 / ATCC [®] 9144 | MM85 | 00035 |
| Staphylococcus aureus (MRSA) | NCTC [®] 12493 | MM64 | 00212 |
| Staphylococcus aureus (MRSA) | NCTC [®] 13373 / ATCC [®] 43300 | MM91 | 00211 |
| Staphylococcus epidermidis | NCTC [®] 13360 / ATCC [®] 12228 | MM15 | 00036 |
| Streptococcus agalactiae | NCTC [®] 8181 / ATCC [®] 13813 | MM16 | |
| Streptococcus pneumoniae | NCTC [®] 12695 /ATCC [®] 6303 | MM19 | |
| Streptococcus pneumoniae | NCTC [®] 12977 /ATCC [®] 49619 | MM95 | |
| Streptococcus pyogenes | NCTC [®] 12696 /ATCC [®] 19615 | MM20 | |
| Vibrio parahaemolyticus | NCTC [®] 10885 | MM06 | 00185 |
| Yersinia enterocolitica | NCTC [®] 12982 / ATCC [®] 9610 | MM80 | 00038 |

Selectrol Strains Listed by WDCM Number

| WDCM | Strain Name | Designation | Code |
|-------|---------------------------------------|--|------|
| 00001 | Bacillus cereus | NCTC [®] 10320 / ATCC [®] 9634 | MM21 |
| 00003 | Bacillus subtilis | NCTC [®] 10400 / ATCC [®] 6633 | MM29 |
| 00007 | Clostridium perfringens | NCTC [®] 8237 / ATCC [®] 13124 | MM45 |
| 00008 | Clostridium sporogenes | NCTC [®] 532 / ATCC [®] 19404 | MM31 |
| 00009 | Enterococcus faecalis | NCTC [®] 775 / ATCC [®] 19433 | MM17 |
| 00011 | Enterococcus hirae | NCTC [®] 13383 /ATCC [®] 10541 | MM35 |
| 00012 | Escherichia coli | NCTC [®] 12923 / ATCC [®] 8739 | MM38 |
| 00013 | Escherichia coli | NCTC [®] 12241 / ATCC [®] 25922 | MM02 |
| 00014 | Escherichia coli O157 (non-toxigenic) | NCTC [®] 12900 / ATCC [®] 700728 | MM93 |
| 00017 | Listeria innocua | NCTC [®] 11288 / ATCC [®] 33090 | MM92 |
| 00019 | Listeria monocytogenes | NCTC [®] 11994 | MM87 |
| 00025 | Pseudomonas aeruginosa | NCTC [®] 12903 / ATCC [®] 27853 | MM10 |
| 00026 | Pseudomonas aeruginosa | NCTC [®] 12924 / ATCC [®] 9027 | MM40 |
| 00028 | Rhodococcus equi | NCTC [®] 1621 / ATCC [®] 6939 | MM97 |
| 00031 | Salmonella Typhimurium | NCTC [®] 12023/ ATCC [®] 14028 | MM11 |
| 00032 | Staphylococcus aureus | NCTC [®] 10788 / ATCC [®] 6538 | MM46 |
| 00033 | Staphylococcus aureus | NCTC [®] 7447 / ATCC [®] 6538P | MM30 |
| 00034 | Staphylococcus aureus | NCTC [®] 12981 / ATCC [®] 25923 | MM13 |
| 00035 | Staphylococcus aureus | NCTC [®] 6571 / ATCC [®] 9144 | MM85 |
| 00036 | Staphylococcus epidermidis | NCTC [®] 13360 / ATCC [®] 12228 | MM15 |
| 00038 | Yersinia enterocolitica | NCTC [®] 12982 / ATCC [®] 9610 | MM80 |
| 00053 | Aspergillus brasiliensis | NCPF [®] 2275 / ATCC [®] 16404 | MM94 |
| 00054 | Candida albicans | NCPF [®] 3179 / ATCC [®] 10231 | MM42 |
| 00055 | Candida albicans | NCPF [®] 3255 / ATCC [®] 2091 | MM28 |
| 00058 | Saccharomyces cerevisiae | NCTC [®] 10716/ ATCC [®] 9763 | MM50 |
| 00082 | Enterobacter cloacae | NCTC [®] 13380 / ATCC [®] 23355 | MM01 |
| 00085 | Enterococcus faecalis | NCTC [®] 13379 / ATCC [®] 51299 | MM52 |
| 00087 | Enterococcus faecalis | NCTC [®] 12697 / ATCC [®] 29212 | MM18 |
| 00090 | Escherichia coli | NCTC [®] 9001 / ATCC [®] 11775 | MM75 |
| 00097 | Klebsiella pneumoniae | NCTC [®] 9633 / ATCC [®] 13883 | MM04 |
| 00107 | Legionella pneumophila serogroup 1 | NCTC [®] 11192 / ATCC [®] 33152 | MM08 |
| 00109 | Listeria monocytogenes | NCTC [®] 7973 / ATCC [®] 35152 | MM48 |
| 00114 | Pseudomonas aeruginosa | NCTC [®] 10662 / ATCC [®] 25668 | MM65 |
| 00131 | Staphylococcus aureus | NCTC [®] 12973 / ATCC [®] 29213 | MM14 |
| 00156 | Campylobacter jejuni | NCTC [®] 11322 / ATCC [®] 29428 | MM82 |
| 00175 | Enterobacter aerogenes | NCTC [®] 10006 / ATCC [®] 13048 | MM26 |
| 00185 | Vibrio parahaemolyticus | NCTC [®] 10885 | MM06 |
| 00211 | Staphylococcus aureus (MRSA) | NCTC [®] 13373 / ATCC [®] 43300 | MM91 |
| 00212 | Staphylococcus aureus (MRSA) | NCTC [®] 12493 | MM64 |



www.tcsbiosciences.co.uk

TCS Biosciences Ltd Botolph Claydon, Buckingham, MK18 2LR, United Kingdom t: +44 (0) 1296 714222, f: +44 (0) 1296 714806, e: sales@tcsgroup.co.uk