Kovac's Reagent

Reagent for indole test of Enterobacteriaceae.

DESCRIPTION

Kovac's Reagent is used in determining the ability of bacteria, primarily Enterobacteriaceae, to produce indole by the deamination of tryptophan.

KIT CONTENTS

- 4 x 25 ml bottles of Kovac's Reagent.
- 1 instruction sheet.

METHOD PRINCIPLE

Indole is one of the degradation products of the bacterial metabolism of the amino acid tryptophan. The bacteria that own the enzyme tryptophanase are able to hydrolyze and deaminate the tryptophan with the production of indole, pyruvic acid and ammonia. Indole test is based on the formation of a red to purple colored complex, due to the indole reaction with aldehydic group of p-dimethylaminobenzaldehyde. The chief requirement for culturing an organism prior to performing the indole test is that the medium contains a sufficient quantity of tryptophan.

REAGENTS

5% (w/v) p-dimethylaminobenzaldehyde dissolved in a solution of 25% hydrochloridric acid and 75% isobutyl alcohol.

TEST PROCEDURE

Inoculate a tube of Peptone Water (ref. 24098) with the organism to be tested and incubate at $35 \pm 2^\circ\text{C}$ for 24-48 hours. Add 2-3 drops of Kovac's Reagent directly to the tube.

INTERPRETATION OF RESULTS

The formation of a red to purple color ("cherry-red ring") in the reagent layer on top of the medium within 30 sec indicates a positive reaction for indole production. A negative reaction shows no color change.

QUALITY CONTROL FOR THE USER

Positive and negative controls should be run simultaneously with the organism to be tested.

Control strains		Color	Indole test
<i>Escherichia coli</i>	ATCC® 25922	Red to purple	Positive
<i>Proteus mirabilis</i>	ATCC® 25933	No color change	Negative

PRECAUTIONS

Kovac's Reagent is classifiable as hazardous under current legislation; it is recommended that the Safety Data Sheet be consulted on its use. The product is intended for *in vitro* diagnostic use only and must be used in the laboratory by properly trained personnel, using approved asepsis and safety methods for handling pathogenic agents.

STORAGE AND TRANSPORT CONDITIONS

2-8°C away from light, until the expiry date on the label. However, our stability studies have shown that the storage or transport at 18-25°C for 4 days, or at 35-39°C for 48 hours, do not alter in any way the performance of the product. Eliminate if signs of deterioration or contamination are evident.

ELIMINATING USED MATERIAL

After use, used Kovac's Reagent and the material that has come into contact with the sample must be decontaminated and disposed of in accordance with the laboratory procedures for the decontamination and disposal of potentially infected material.

BIBLIOGRAPHY

1. Murray, Baron, Pfaller, Tenorev and Tenover: Manual of Clinical Microbiology (1995).
2. Bayley and Scott's: Diagnostic Microbiology (1986).
3. Edwin H. Lennette: Manual of Clinical Microbiology (1995).

PRESENTATION

Product	Ref.	Content
Kovac's Reagent	80271	100-200 tests

TABLE OF SYMBOLS

LOT	Batch code	IVD	In vitro Diagnostic Medical Device		Manufacturer		Use by		Fragile, handle with care
REF	Catalogue number		Temperature limitation		Contains Sufficient for <n> tests		Caution, consult accompanying documents		Do not reuse

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LIOFILCHEM® S.r.l.