Horse Serum

EN

REF SR0035C & SR0035D

* This instructions for use (IFU) document is intended to be read in conjunction with the IFU for Trichomonas Medium (CM0161B) Brucella Medium Base (CM0169B/K), Blood Agar Base No. 2 (CM0271B/K/R/T), Columbia Blood Agar Base (CM0331B/Y/W/L/S/U/K/R/T) GBS Agar Base (Islam) (CM0755B).

Intended Use

Horse serum (SR0035C and SR0035D) device is an extract of horse blood intended for use in a culture media preparation isolate and differentiate clinical pathogens like Trichomonas vaginalis, Brucella species, Group B Streptococcus from clinical samples like genito-urinary samples, blood, neonatal and vaginal samples etc.

Horse serum (SR0035C and SR0035D) devices are used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having bacterial infections.

The device is for professional use only, is not automated, nor is it a companion diagnostic.

Summary and Explanation

Trichomonas vaginalis is considered as one of the major targets of sexually transmitted infections (STI) in humans¹.

T. vaginalis is a flagellated parasite that primarily attacks female genital tract and male urethra and prostrate is one of the leading causes of sexually transmitted disease in humans. The major symptoms include dysuria, itching, abdominal pain and vulvar irritation. Symptoms include vaginal discharge (which is often diffuse, malodorous, yellow green), dysuria, itching, vulvar irritation and abdominal pain². In women, *T. vaginalis* infection can range from mild to severe inflammatory vaginitis whereas in men

T. vaginalis is considered as the major cause of non - gonococcal urethritis³.

Brucellosis is a type of infection in humans caused by a Gram-negative bacteria named *Brucella* and the disease is transmitted by the consumption of unpasteurized dairy products, or through contact with infected animals mainly in the slaughterhouses. Three types of *Brucella* species that causes brucellosis are *B. melitensis*, *B. suis* and *B. abortus*⁴.

Streptococcus agalactiae or Lancefield Group B Streptococci (GBS) is a Gram-negative coccus that are facultative anaerobes and can cause infections in pregnant women and newborn babies and other adults with underlying medical conditions like Diabetes mellitus⁵.

Principle of Method

Horse serum (SR0035C and SR0035D) devices are aseptically collected defibrinated horse serum. Horse serum (SR0035C and SR0035D) devices are extract of horse blood for addition to culture media. These devices contain high protein, trace elements, spreading factors, low molecular weight nutrients, hormones and growth factors suitable as a supplement for the growth of microorganisms..

Typical Formula

SR0035C

Defibrinated Horse Blood 100ml

SR0035D

Defibrinated Horse Blood 500ml

Materials Provided

SR0035C: 1x 100ml frozen Horse Serum

SR0035D: 1x 500ml frozen Horse Serum

Materials Required but Not Supplied

1. Trichomonas Medium CM0161B
2. Brucella Medium Base CM0169B/K
3. Blood Agar Base No.2 CM0271B/K/R/T
4. ColumbiaBlood Agar Base CM0331B/Y/W/S/L/R/T/K
5. GBS Agar Base (Islam) CM0755B

Storage

- Store product in its original packaging between 2°C and 8°C, or freeze at -20°C for longer term storage.
- Keep container tightly closed.
- The product may be used until the expiry date stated on the label.
- Store away from light.
- Allow reconstituted product to equilibrate to room temperature before use.

Warnings and Precautions

- DO NOT USE SUPPLEMENTS IF YOU ARE PREGNANT OR CONTEMPLATING PREGNANCY.
- Supplements may be dispensed on the bench, wear disposable gloves and safety glasses.
- All waste (i.e. rehydrated supplement) must be emptied into the HAZARDOUS WASTE container situated in the Microbiology laboratory.
- The vial and rubber bung must be rinsed into the HAZARDOUS WASTE container.
- The metal lid and rinsed rubber bung should be placed in the waste bin.
- The rinsed glass vial should then be placed in the glass bin.

- For in vitro diagnostic use only.
 - For professional use only.
 - Inspect the product packaging before first use.
 - Do not use the product if there is any visible damage to the packaging or vials.
 - Do not use the product beyond the stated expiry date.
 - Do not use the device if signs of contamination are present.
 - Do not use the device if the colour has changed or there are other signs of deterioration.
 - It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations.
- Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product (www.thermofisher.com).

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

Specimen Collection, Handling and Storage

Specimen should be collected and handled following local recommended guidelines, such as the UK Standards for Microbiology Investigations (UK SMI) S 6, ID 7, B 28, B 37, B 58.

Procedure

1. Allow supplement to equilibrate to room temperature before use.
2. Aseptically add the vial contents to specified culture media (see additional IFU's)
3. Sterilise by autoclaving at 121°C for 15 minutes.
4. Cool to 50°C.
5. Mix well and aseptically dispense into sterile containers.


Refer to the IFU for CM0161B, CM0169B/K, CM0271B/K/R/T, CM0331B/Y//W/L/S/U/K/R/T or CM0755B. IFUs can be found at www.thermofisher.com

Bibliography

1. Public Health England 2013. 'Sexually Transmitted Infections'. UK Standards for Microbiology Investigations S6 <https://www.gov.uk/government/publications/smi-s-6-sexually-transmitted-infections>
2. Kissinger, P. (2015). Trichomonas vaginalis: a review of epidemiologic, clinical and treatment issues. *BMC infectious diseases*, 15(1), 1-8.
3. Muzny, C. A., & Schwebke, J. R. (2013). The clinical spectrum of Trichomonas vaginalis infection and challenges to management. *Sexually transmitted infections*, 89(6), 423-425.
4. Centers for Disease Control and Prevention. 2021. 'Brucellosis' <https://www.cdc.gov/brucellosis/index.html>
5. Public Health England 2018 'Detection of carriage of group B streptococci'. UK Standards for Microbiology Investigations 'SMI B 58' <https://www.gov.uk/government/publications/smi-b-58-processing-swabs-for-group-b-streptococcal-carriage>

Symbol Legend

Symbol	Definition
	Catalogue number
	In Vitro Diagnostic Medical Device
	Batch code
	Temperature limit
	Use-by date
	Keep away from sunlight
	Do not re-use

	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Do not use if packaging damaged and consult instructions for use
	Manufacturer
	Authorized representative in the European Community/ European Union
	European Conformity Assessment
	UK Conformity Assessment
	Unique device identifier
	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union
	Made in the United Kingdom

©2022 Thermo Fisher Scientific Inc. All rights reserved. ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection.

All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Oxoid Limited, Wade Road, Basingstoke, RG24 8PW, UK

For technical assistance please contact your local distributor.

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Bibliographie

1. Public Health England 2013. 'Sexually Transmitted Infections'. UK Standards for Microbiology Investigations S6 <https://www.gov.uk/government/publications/smi-s-6-sexually-transmitted-infections>
2. Kissinger, P. (2015). Trichomonas vaginalis: a review of epidemiologic, clinical and treatment issues. BMC infectious diseases, 15(1), 1-8.
3. Muzny, C. A., & Schwebke, J. R. (2013). The clinical spectrum of Trichomonas vaginalis infection and challenges to management. Sexually transmitted infections, 89(6), 423-425.
4. Centers for Disease Control and Prevention. 2021. 'Brucellosis' <https://www.cdc.gov/brucellosis/index.html>
5. Public Health England 2018 'Detection of carriage of group B streptococci'. UK Standards for Microbiology Investigations 'SMI B 58' <https://www.gov.uk/government/publications/smi-b-58-processing-swabs-for-group-b-streptococcal-carriage>

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
Bibliograafia

1. Public Health England 2013. 'Sexually Transmitted Infections'. UK Standards for Microbiology Investigations S6 <https://www.gov.uk/government/publications/smi-s-6-sexually-transmitted-infections>
2. Kissinger, P. (2015). *Trichomonas vaginalis: a review of epidemiologic, clinical and treatment issues. BMC infectious diseases*, 15(1), 1-8.
3. Muzny, C. A., & Schwebke, J. R. (2013). The clinical spectrum of *Trichomonas vaginalis* infection and challenges to management. *Sexually transmitted infections*, 89(6), 423-425.
4. Centers for Disease Control and Prevention. 2021. 'Brucellosis' <https://www.cdc.gov/brucellosis/index.html>
5. Public Health England 2018 'Detection of carriage of group B streptococci'. UK Standards for Microbiology Investigations 'SMI B 58' <https://www.gov.uk/government/publications/smi-b-58-processing-swabs-for-group-b-streptococcal-carriage>

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Oxoid Limited, Wade Road, Basingstoke, RG24 8PW, ÜK

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Alla övriga varumärken tillhör Thermo Fisher Scientific Inc. och dess dotterbolag.