Oxidase Reagent

IVD

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Oxidase Reagent

IVD

Detection of the enzyme cytochrome oxidase

SUMMARY AND EXPLANATION

This test is used to detect the production of the enzyme cytochrome oxidase by bacteria. This enzyme is characteristic of the genus *Neisseria* and most species of *Pseudomonas*. It enables the differentiation of Gram (-) bacilli (1, 2).

PRINCIPLE

The oxidase test is based on the bacterial production of an intracellular oxidase enzyme.

In the presence of atmospheric oxygen and cytochrome C, this enzyme oxidizes the phenylenediamine reagent to form a purple colored compound, indophenol.

Ascorbic acid is incorporated in the reagent as a reducing agent to reduce autooxidation and to improve stability (6). This formulation is based on the oxidase reagent formula of Kovacs (3).

CONTENT OF THE KIT

REF 55 635 Ready-to-use reagent

Kit contents:

- 50 ampules with 0.75 ml of reagent per ampule
- 1 reusable ampule crusher
- + 1 package insert

COMPOSITION

Theoretical formula.

This reagent can be adjusted according to the performance criteria required.

N,N,N,N-tetramethyl-1,4-phenylenediamine	10.0 g
Ascorbic acid	2.0 g
Deionized water	1 L

MATERIAL REQUIRED BUT NOT PROVIDED

- Inoculating loop (platinum or plastic).
- Non-impregnated disks (diameter 6 mm)(Ref. 54 991).

WARNINGS AND PRECAUTIONS

- **For *in vitro* diagnostic use only.**
- **For professional use only.**
- All samples and microbial cultures must be treated as potentially infectious and handled appropriately by trained and competent personnel. Aseptic technique and usual precautions for handling the group of bacteria studied should be observed throughout this procedure. See "NCCLS M29-A, *Protection of Laboratory Workers from Instrument Biohazards and Infectious Disease Transmitted by Blood, Body Fluids, and Tissue ; Approved Guideline - Current Revision*". For additional handling precautions, refer to "Biosafety in Microbiological and Biomedical Laboratories, CDC/NIH, Latest Edition", or to the regulations currently in use in each country.
- Do not use reagents past the expiration date.
- This product should not be used if:
 - the color has changed (the reagent is originally colorless to pale yellow),
 - the packaging is damaged,
 - there are other signs of deterioration (broken ampule)
- The performance data presented were obtained using the procedure indicated in this package insert. Any change or modification in the procedure may affect the results.

- Interpretation of the test results should be made taking into consideration the patient's history, colonial and microscopic morphology and, if necessary, the results of any other tests performed.

STORAGE CONDITIONS

- **Store the reagents at 18°C-25°C in their box until the expiration date.**
- Do not freeze.
- **Protect from light.**
- **The oxidase reagent quickly auto-oxidizes and loses its sensitivity. Unused reagent should be discarded after 12 hours.**

SPECIMENS

The specimen consists of an isolated colony for which cytochrome oxidase enzyme production is to be detected. The colony should be collected from an 18-24-hour-old culture grown on a solid agar medium.

INSTRUCTIONS FOR USE

Opening the ampule:

1. Place the dropper in the reusable ampule crusher.
2. Lightly tap the bottom of the dropper to dislodge any bubbles that may have formed.
3. Take hold of the middle of the dropper/crusher and press gently to crush the ampule inside.

Test procedure:

1. Dispense **precisely 1 drop** of reagent onto a non-impregnated disc (diameter 6 mm) (Ref. 54 991).
2. Smear the colony on the disc.

READING AND INTERPRETATION

- A positive test is indicated by the development of a violet to purple color within 10-30 seconds.
- Delayed reactions or no color development indicate a negative test.

QUALITY CONTROL

Protocol:

The activity of the reagent can be tested using the following strains grown on Trypcase-Soy agar:

- *Pseudomonas aeruginosa* ATCC® 27853
- *Escherichia coli* ATCC® 25922.

Expected results:

Strain	Results
<i>Pseudomonas aeruginosa</i> ATCC® 27853	Positive: violet color development
<i>Escherichia coli</i> ATCC® 25922	Negative: no color development

Note:

It is the responsibility of the user to perform Quality Control taking into consideration the intended use of the medium, and in accordance with any local applicable regulations (frequency, number of strains, etc.).

LIMITATIONS OF THE METHOD

- The oxidase test should not be performed on colonies obtained on EMB or CHAPMAN 2 agar, or colonies from a 48-hour-old culture grown on solid agar media.
- Oxidase detection should not be performed on isolated colonies showing a spontaneous color development (violet, pink, green, black, etc.). In this case, reading of the test is impossible.
- Using too large volumes of reagent could lead to false negative results. Only use one drop of reagent as indicated in the instructions for use.
- A platinum or plastic inoculating loop or needle (7) is recommended for use in oxidase testing. The presence of any trace of iron (nichrome) can catalyze an oxidase reaction, resulting in a false-positive reaction.
- Discard the remaining portion of partially used reagent after 12 hours.
- Weak oxidase producers e.g. *Pasteurella*, may give negative results (4).
- False-negative results occur if mixed cultures contain the genera *Pseudomonas* and *Neisseria*. An inhibitory substance is produced by *Pseudomonas* spp which interferes with the production of oxidase by *Neisseria* spp.

PERFORMANCE

A study was carried out using 100 strains.

Out of the 100 strains tested, 92 were described as oxidase-positive or oxidase-negative and were used to calculate the reagent performance: sensitivity and specificity; the other 8 strains were described as oxidase-variable (*Bacillus* and *Micrococcus*) and were not taken into account for the performance calculations.

On the basis of this study, the following results were obtained :

- out of the 65 positive strains tested, 4 gave negative results: *Haemophilus parainfluenzae* (1 out of 2 strains tested), *Pasteurella multocida* (1 out of 4 strains tested) and *Photobacterium damsela* (2 out of 2 strains tested).
- out of the 27 negative strains tested, only one tested positive: *Pseudomonas oryzae* (1 out of 1 strain tested).

Number of strains	Sensitivity CI* (95%)	Specificity CI* (95%)
92	93.85% (84.98% – 97.62%)	96,30% (81.28% - 99.36%)

*confidence interval.

WASTE DISPOSAL

Unused reagents may be considered as non hazardous waste and disposed of accordingly.

Dispose of used reagents as well as any other contaminated disposable materials following procedures for infectious or potentially infectious products.

It is the responsibility of each laboratory to handle waste and effluents produced according to their nature and degree of hazardousness and to treat and dispose of them (or have them treated and disposed of) in accordance with any applicable regulations.

LITERATURE REFERENCES

1. BARON E.J. - PETERSON L.R. - FINEGOLD S.M. – Bailey and Scott's Diagnostic Microbiology - 9th ed.- Ed. Mosby, 1994 – ISBN 0-8016-6987-1.
2. COWAN S.T., STEEL K.J. - Manual for the Identification of Medical Bacteria, Cambridge: University Press. - 1965, p. 22.
3. KOVACS N. – Identification of *Pseudomonas pyocyanea* by oxidase reaction. - *Nature*, 1956, n° 4535, p. 703.
4. KNAPP J.S., KOUMANS E.H. – *Neisseria and Branhamella* in MURRAY P.R., BARON E.J., PFALLER M.A. et al. - *Manual of Clinical Microbiology*, 7th ed., Washington, D.C. : ASM Press, 1999, p. 587-597, p. 1070.
5. MAC FADDIN J.F. - *Biochemical tests for identification of medical bacteria*. - 2nd ed. – Baltimore : Williams and Wilkins, 1980, Chap. 27 : Oxydase test. Paragraphe VIII: Precautions - p. 257.
6. STEEL K.J. – The oxidase activity of Staphylococci. - *J. Appl. Bacteriol.*, 1962, vol. 25, n°3, p. 445-455.
7. STEEL K.J. - The oxidase reaction as a taxonomic tool. - *J. Gen. Microbiol.*, 1961, vol. 25, p. 297-306.


INDEX OF SYMBOLS

Symbole	Meaning
 or REF	Catalogue number
	In Vitro Diagnostic Medical Device
	Manufacturer
	Temperature limitation
	Use by
	Batch code
	Consult instructions for use
	Protect from light

WARRANTY

bioMérieux disclaims all warranties, express or implied, including any implied warranties of MERCHANTABILITY AND FITNESS FOR A PARTICULAR USE. bioMérieux shall not be liable for any incidental or consequential damages. IN NO EVENT SHALL BIOMERIEUX'S LIABILITY TO CUSTOMER UNDER ANY CLAIM EXCEED A REFUND OF THE AMOUNT PAID TO BIOMERIEUX FOR THE PRODUCT OR SERVICE WHICH IS THE SUBJECT OF THE CLAIM.

ATCC is a used, pending and/or registered trademark belonging to American Type Culture Collection.

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bioMérieux and the blue logo are used, pending and/or registered trademarks belonging to bioMérieux SA or one of its subsidiaries.

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