Pseudalert*

06-18569-10

IDEXX

IDX 33/05 – 03/16
WATER ANALYSIS METHODS
<http://nf-validation.afnor.org>

The method Pseudalert/Quanti-Tray* for water analysis is granted NF Validation by AFNOR Certification as an alternative method to the standard ISO 16266 for enumeration of *Pseudomonas aeruginosa* in drinking water and pool water, under the Certificate number: IDX 33/05-03/16.

For more information about end of validity, please refer to the certificate NF Validation on the website mentioned above.

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IDEXX Water Quality Control Laboratory is accredited to ISO/IEC 17025:2017

IDEXX

IDEXX Laboratories, Inc., One IDEXX Drive, Westbrook, Maine 04092 USA
idexx.com/water

For Technical Support, please call:

North America: 1 207 556 4496/1 800 321 0207

Latin America: +0800 728 2482

Europe: +00800 7274 3399

UK: +44 (0) 1638 676800

China: 400-678-6682 ext. 2

Japan: +81 3 5301 6800

Australia: 1300 443 399

Introduction

The Pseudalert* (98-18076-00, 98-18078-00, 98-26014-00) test detects the presence of *Pseudomonas aeruginosa* in water samples. The test is based on a bacterial enzyme detection technology that signals the presence of *Pseudomonas aeruginosa* through the hydrolysis of a substrate present in the Pseudalert reagent. *Pseudomonas aeruginosa* cells rapidly grow and reproduce using the rich supply of amino acids, vitamins and other nutrients present in the Pseudalert reagent. Actively growing strains of *Pseudomonas aeruginosa* have an enzyme that cleaves the substrate to produce a blue fluorescence under ultraviolet (UV) light. Pseudalert detects *Pseudomonas aeruginosa* at 1 cfu in either 100 mL or 250 mL samples within 24 hours.

Storage

Store at 2–30°C away from light

Presence/Absence (P/A) Procedure

1. Add contents of one appropriately sized snap pack to either a 100 mL or 250 mL sample in a sterile, transparent, nonfluorescing vessel.
2. Cap vessel and shake.
3. Incubate at 38±0.5°C for 24 to 28 hours.
4. Read results according to Result Interpretation table below.

**Quanti-Tray* Enumeration Procedure (100 mL Samples Only)**

1. Add contents of one pack to a 100 mL water sample in a sterile vessel.
2. Cap vessel and shake until dissolved.
3. Add 2 drops of IDEXX Antifoam Solution¹ to the sample/reagent mixture.
NOTE: IDEXX 120 mL sample vessels² containing antifoam are also available.
4. Pour sample/reagent mixture into a Quanti-Tray or Quanti-Tray*/2000 and seal in a Quanti-Tray* Sealer.
5. Place the sealed tray in a 38±0.5°C incubator for 24 to 28 hours.
6. Read results according to the Result Interpretation table below. Count the number of positive wells and refer to the MPN table provided with the trays to obtain a Most Probable Number.

NOTE: Quanti-Tray Enumeration Procedure for 250 mL available from IDEXX.

Result Interpretation

| Appearance | Result |
|--------------------------------|--|
| No blue fluorescence | Negative for <i>Pseudomonas aeruginosa</i> |
| Blue fluorescence ¹ | Positive for <i>Pseudomonas aeruginosa</i> |

¹Presence of blue fluorescence greater than the amount present in a negative control sample

- Look for blue fluorescence with a 365-nm UV light held within 5 inches of the sample in a dark environment. Face light away from your eyes and toward the sample.
- Refer to –/+ fluorescence read guide on kit box. Chart colors are as seen under a UV light.
- Pseudalert results are definitive at 24–28 hours. In the presence/absence procedure, positives for *Pseudomonas aeruginosa* observed before 24 hours, and negatives observed after 28 hours are also valid.

**Procedural Notes**

- Use only sterile, nonbuffered, oxidant-free water for dilutions.
- For comparison, an incubated sterile water blank containing Pseudalert reagent (negative control) can be used when interpreting results.
- This insert may not reflect your local regulations. For compliance testing, be sure to follow appropriate regulatory procedures.
- Pseudalert is a primary water test. Pseudalert performance characteristics do not apply to samples altered by preenrichment or concentration.
- Pseudalert has not been validated for use with flavored bottled water, marine water or carbonated water samples.
- Aseptic technique should always be followed when using Pseudalert. Dispose of materials in accordance with Good Laboratory Practices.
- The presence of a high mineral content (especially magnesium and/or calcium) can cause the Pseudalert reagent mixture to become cloudy, but this does not affect the outcome.
- Interpret any blue fluorescence as positive, even if the fluorescent signal is weak.
- If you are unsure about a well or vessel with weak fluorescence at 24 hours, incubate for another 1–4 hours.
- If using a gas-filled mercury UV bulb to read fluorescence, the power rating should be 6 watts.

Quality Control Procedures

1. One of the following quality control procedures is recommended for each lot of Pseudalert:
 - A. IDEXX-QC *Pseudomonas*³: *Pseudomonas aeruginosa*, *Escherichia coli*, *Pseudomonas fluorescens*.
 - B. Additional Quality Control Methods.
 - i. For each of the American Type Culture Collection (ATCC)^{TM, 4} bacterial strains, (*Pseudomonas aeruginosa* ATCC 27853 (WDCM 00025) or 10145 (WDCM 00024), *Escherichia coli* ATCC 25922 (WDCM 00013) and *Pseudomonas fluorescens* ATCC 13525 (WDCM 00115)) streak the culture onto labeled TSA or blood agar plates and incubate at 35±0.5°C for 18–24 hours. These strains are not available from IDEXX and must be purchased from ATCC or other providers.
 - ii. For each bacterial strain, touch a sterile 1 µL inoculating loop to a colony and use it to inoculate a labeled test tube containing 5 mL of sterile deionized water. Close cap and shake thoroughly.
 - iii. For each bacterial strain, take a 1 µL loop from the test tube and use it to inoculate a labeled vessel containing either 100 mL or 250 mL of sterile deionized water. These are your controls.
2. Follow the P/A Procedure or Quanti-Tray Enumeration Procedure above.
3. Results should match the Result Interpretation table above.

NOTE: IDEXX internal quality control testing is performed in accordance with ISO 11133:2014. Quality Control Certificates are available at idexx.com/water.

¹ IDEXX Antifoam Solution catalog number: WAFDB

² IDEXX 120 mL Shrink-banded Vessels with Antifoam catalog number: WV120SBAF-200

³ IDEXX-QC *Pseudomonas*, IDEXX Catalog #WQC-PSE

⁴ ATCC is a registered trademark owned by the American Type Culture Collection. American Type Culture Collection, 1-800-638-6597 atcc.org

*The Pseudalert Test is intended to be used only for water quality research and analysis, by technically qualified individuals or under their supervision.

*Pseudalert and Quanti-Tray are trademarks or registered trademarks of IDEXX Laboratories, Inc. or its affiliates in the United States and/or other countries.

Patent information: idexx.com/patents.

© 2024 IDEXX Laboratories, Inc. All rights reserved.

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