Yersinia Supplement

Selective supplement for the isolation of *Yersinia enterocolitica*.

DESCRIPTION

Yersinia Supplement is a selective lyophilized supplement used for the preparation of Yersinia Selective Agar Base (ref. 610111, 620111). The complete medium is used for the isolation of *Yersinia enterocolitica* from clinical and nonclinical specimens.

KIT CONTENTS

Each kit contains:

- 10 vials of freeze-dried Yersinia Supplement
- 1 instruction sheet

PRINCIPLE OF THE METHOD

Yersinia Supplement is an antibiotic mix that inhibits the growth of most Gram-negative enteric bacteria.

COMPOSITION

Antibiotic	Content / vial	Content / liter of medium
Cefsulodin	7.5 mg	15.0 mg
Irgasan	2.0 mg	4.0 mg
Novobiocin	1.25 mg	2.5 mg

PROCEDURE FOR USE

1. Aseptically reconstitute the content of one vial of Yersinia Supplement with 2 ml sterile distilled water and 1 ml ethanol(*).
2. Mix to complete dissolution and add to 500 ml Yersinia Selective Agar Base, autoclaved and cooled to 45-50°C.
3. Mix well and dispense in Petri dishes.

(*) Ethanol 50% purity minimum.

TECHNIQUE AND INTERPRETATION OF THE RESULTS

Refer to the technical documentation of the medium being prepared.

QUALITY CONTROL

1. Control of the appearance: freeze-dried white in colour.
2. Microbiological control.
Prepare the medium per label directions. Inoculate the plates with the strains indicated below and incubate at 30 ± 1°C for 18-24 hours.

Control strains

Yersinia enterocolitica
Escherichia coli
Staphylococcus aureus

WDCM 00038
WDCM 00012
WDCM 00034

Growth

Good
Partially to totally inhibited
Inhibited

PRECAUTIONS

Yersinia Supplement is classifiable as hazardous under current legislation; it is recommended that the Safety Data Sheet be consulted on its use. The product is a selective supplement to be used only for *in vitro* diagnostic use. It is intended for use in a professional environment and must be used in the laboratory by properly trained personnel, using approved asepsis and safety methods for handling pathogenic agents.

STORAGE

2-8°C in its original packaging. Keep away from sources of heat and avoid excessive changes of temperature. Use until the expiry date indicated on the label. Eliminate without using if there are signs of deterioration. Once reconstituted, the product can be stored for a maximum duration of 30 days at -20°C, shielded from light.

REFERENCES

- EN ISO 11133:2014. Microbiology of food, animal feed and water – Preparation, production, storage and performance testing of culture media.
- ISO 10273:2003. Microbiology of food and animal feeding stuffs – Horizontal method for the detection of presumptive pathogenic *Yersinia enterocolitica*.
- Schieman, D.A. (1979) Synthesis of a selective agar medium for *Yersinia enterocolitica*. Can. J. Microbiol. 25:1298-1304.
- Schieman, D.A. (1980) *Yersinia enterocolitica*: Observation on some growth characteristics and response to selective agents. Can. J. Microbiol. 43:14-27.
- Devenish, J.A., and D.A. Schieman (1981) An abbreviated scheme for identification of *Yersinia enterocolitica* isolated from food enrichments on CIN (cefsulodin-irgasan-novobiocin) agar. Can. J. Microbiol. 27:937-941.

PRESENTATION

Product	Ref.	Contents
Yersinia Supplement	81039	10 vials

One vial is sufficient to prepare 500 ml of medium.

TABLE OF SYMBOLS

 IVD In Vitro Diagnostic Medical Device	 Do not reuse	 Manufacturer	 Contains sufficient for <n> tests	 Temperature limitation
 REF Catalogue number	 Fragile, handle with care	 Use by	 Caution, consult accompanying documents	 LOT Batch code

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REFERENSER

- IEN ISO 11133:2014. Microbiology of food, animal feed and water – Preparation, production, storage and performance testing of culture media.
- ISO 10273:2003. Microbiology of food and animal feeding stuffs – Horizontal method for the detection of presumptive pathogenic *Yersinia enterocolitica*.
- Schieman, D.A. (1979) Synthesis of a selective agar medium for *Yersinia enterocolitica*. Can. J. Microbiol. 25:1298-1304.
- Schieman, D.A. (1980) *Yersinia enterocolitica*: Observation on some growth characteristics and response to selective agents. Can. J. Microbiol. 43:14-27.
- Devenish, J.A., and D.A. Schieman (1981) An abbreviated scheme for identification of *Yersinia enterocolitica* isolated from food enrichments on CIN (cefsulodin-irgasan-novobiocin) agar. Can. J. Microbiol. 27:937-941.

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LIOFILCHEM® S.r.l.