

## Letter of Authorization

To whom it may concern,

We, Getein Biotech, Inc. (No.9 BoFu Road, Luhe District, Nanjing, 211505, China), hereby authorize Sanmedico SRL (Address: Republic of Moldova, Chisinau, MD-2059, Petricani street, 88/1, office 10) as our official and non-exclusive distributor for registration, promoting, selling, distributing and providing after-sale services of under-mentioned product in the territory of Moldova only:

FIA8000 Quantitative Immunoassay Analyzer and Reagents

Getein1100 Immunofluorescence Quantitative Analyzer and Reagents

Getein1160 Immunofluorescence Quantitative Analyzer and Reagents

Getein 1600 Immunofluorescence Quantitative Analyzer and Reagents

Sanmedico SRL will comply with the laws and regulations of the countries and regions where they are located in and where they are selling mentioned product.

Sanmedico SRL will carry out marketing efforts to fulfill service and maintenance for above mentioned products and will provide with users benefits of having a local stock of above mentioned products and on time delivery with every order, supported by a local service in local language.

This authorization starts from Jan 1, 2025 and will be valid to December 31 2025.

Getein Biotech, Inc. has the right to terminate the authorization before validity and will inform **Sanmedico SRL** with 10 days in advance.

**Getein Biotech, Inc.**

Name: Steven Zhou

Position: Overseas Sales Director



基蛋生物科技股份有限公司  
GETEIN BIOTECH, INC.

# Certificate of Registration

QUALITY MANAGEMENT SYSTEM - ISO 13485:2016

This is to certify that: **Getein Biotech, Inc.**  
No.9 Bofu Road  
Luhe District  
Nanjing  
Jiangsu  
211505  
China

基蛋生物科技股份有限公司  
中国  
江苏省  
南京市  
六合区  
沿江工业开发区  
博富路9号  
邮编：211505

Holds Certificate No: **MD 728432**

and operates a Quality Management System which complies with the requirements of ISO 13485:2016 for the following scope:

Please see scope page.

For and on behalf of BSI:

**Graeme Tunbridge, Senior Vice President Medical Devices**

Original Registration Date: 2020-05-29

Latest Revision Date: 2023-04-26

Effective Date: 2023-07-26

Expiry Date: 2026-07-25

Page: 1 of 3



...making excellence a habit.™

Certificate No: **MD 728432**

## Registered Scope:

Design & Development, Manufacture and Distribution of Chemiluminescence Immunoassay, Biochemistry Assay, Point of Care Assay (including Colloidal Gold Assay, Immunofluorescence Assay, Dry Chemistry Assay), PCR Assay and Colloidal Gold self-testing Assay to detect infectious disease. Design & Development, Manufacture and Distribution of Analyzers in use of Chemiluminescence Immunoassay, Biochemistry Assay, Point of Care Assay (including Colloidal Gold Assay, Immunofluorescence Assay, Dry Chemistry Assay), PCR Assay to detect infectious disease, Immunofluorescence self-testing Assay to detect dyslipidemia disease, Blood Coagulation Assay to detect thrombotic disease.

研发, 生产和销售化学发光法试剂, 生化试剂, 即时诊断 (包括胶体金法, 免疫荧光法, 干式化学法) 试剂, 传染病相关PCR分子诊断试剂和胶体金自测试剂。 研发, 生产和销售用于化学发光法试剂, 生化试剂, 即时诊断 (包括胶体金法, 免疫荧光法, 干式化学法) 试剂, 传染病相关PCR分子诊断试剂, 血脂异常疾病相关免疫荧光自测试剂, 血栓疾病相关血凝试剂配套使用的分析仪。



Original Registration Date: 2020-05-29

Latest Revision Date: 2023-04-26

Effective Date: 2023-07-26

Expiry Date: 2026-07-25

Page: 2 of 3

This certificate was issued electronically and remains the property of BSI and is bound by the conditions of contract.

An electronic certificate can be authenticated [online](#).

Printed copies can be validated at [www.bsi-global.com/ClientDirectory](http://www.bsi-global.com/ClientDirectory) or telephone +86 10 8507 3000.

Information and Contact: BSI, John M. Keynesplein 9, 1066 EP Amsterdam The Netherlands. Tel: +31 (0) 20 3460 780

BSI Group The Netherlands B.V., registered in the Netherlands under number 33264284, at John M. Keynesplein 9, 1066 EP Amsterdam, The Netherlands

A Member of the BSI Group of Companies.

Certificate No: **MD 728432**

Location

Getein Biotech, Inc.  
No.9 Bofu Road  
Luhe District  
Nanjing  
Jiangsu  
211505  
China  
基蛋生物科技股份有限公司  
中国  
江苏省  
南京市  
六合区  
沿江工业开发区  
博富路9号  
邮编: 211505

Registered Activities

Design & Development, Manufacture and Distribution of Chemiluminescence Immunoassay, Biochemistry Assay, Point of Care Assay (including Colloidal Gold Assay, Immunofluorescence Assay, Dry Chemistry Assay), PCR Assay and Colloidal Gold self-testing Assay to detect infectious disease. Design & Development, Manufacture and Distribution of Analyzers in use of Chemiluminescence Immunoassay, Biochemistry Assay, Point of Care Assay (including Colloidal Gold Assay, Immunofluorescence Assay, Dry Chemistry Assay), PCR Assay to detect infectious disease, Immunofluorescence self-testing Assay to detect dyslipidemia disease, Blood Coagulation Assay to detect thrombotic disease.  
研发, 生产和销售化学发光法试剂, 生化试剂, 即时诊断 (包括胶体金法, 免疫荧光法, 干式化学法) 试剂, 传染病相关PCR分子诊断试剂和胶体金自测试剂。 研发, 生产和销售用于化学发光法试剂, 生化试剂, 即时诊断 (包括胶体金法, 免疫荧光法, 干式化学法) 试剂, 传染病相关PCR分子诊断试剂, 血脂异常疾病相关免疫荧光自测试剂, 血栓疾病相关血凝试剂配套使用的分析仪。

Getein Biotech, Inc.  
No. 6 KeFeng Road  
Jiangbei New District  
Nanjing  
Jiangsu  
211505  
China  
基蛋生物科技股份有限公司  
中国  
江苏省  
南京  
江北新区  
科丰路6号  
邮编: 211505

Manufacture of Chemiluminescence Immunoassay, Biochemistry Assay, Point of Care Assay (including Colloidal Gold Assay, Immunofluorescence Assay, Dry Chemistry Assay), Colloidal Gold self-testing Assay to detect infectious disease. Manufacture of Analyzers in use of Chemiluminescence Immunoassay, Biochemistry Assay, Point of Care Assay (including Colloidal Gold Assay, Immunofluorescence Assay, Dry Chemistry Assay), PCR Assay to detect infectious disease, Immunofluorescence self-testing Assay to detect dyslipidemia disease, Blood Coagulation Assay to detect thrombotic disease.  
生产化学发光法试剂, 生化试剂, 即时诊断 (包括胶体金法, 免疫荧光法, 干式化学法) 试剂和传染病相关胶体金自测试剂。 生产用于化学发光法试剂, 生化试剂, 即时诊断 (包括胶体金法, 免疫荧光法, 干式化学法) 试剂, 传染病相关PCR分子诊断试剂, 血脂异常疾病相关免疫荧光自测试剂, 血栓疾病相关血凝试剂配套使用的分析仪。

Original Registration Date: 2020-05-29

Effective Date: 2023-07-26

Latest Revision Date: 2023-04-26

Expiry Date: 2026-07-25

Page: 3 of 3

This certificate was issued electronically and remains the property of BSI and is bound by the conditions of contract.

An electronic certificate can be authenticated [online](#).

Printed copies can be validated at [www.bsi-global.com/ClientDirectory](http://www.bsi-global.com/ClientDirectory) or telephone +86 10 8507 3000.

Information and Contact: BSI, John M. Keynesplein 9, 1066 EP Amsterdam The Netherlands. Tel: +31 (0) 20 3460 780

BSI Group The Netherlands B.V., registered in the Netherlands under number 33264284, at John M. Keynesplein 9, 1066 EP Amsterdam, The Netherlands

A Member of the BSI Group of Companies.

# EC Declaration of Conformity

according to Directive 98/79/EC, on in vitro diagnostic medical devices

Ref. No.:20220513-A05

**Manufacturer**  
(Name, Address)

**Getein Biotech, Inc.**  
No. 9 Bofu Road, Luhe District, Nanjing, 211505, China

**Authorized Representative**  
(Name, Address)

**CMC Medical Devices & Drugs S.L.**  
Add: C/ Horacio Lengo N° 18, CP 29006, Málaga, Spain

**Medical device**

No.	Product Name
1	Getein 1100 Immunofluorescence Quantitative Analyzer
2	Cardiac Troponin I Fast Test Kit (Immunofluorescence Assay)
3	NT-proBNP Fast Test Kit (Immunofluorescence Assay)
4	hs-CRP+CRP Fast Test Kit (Immunofluorescence Assay)
5	NT-proBNP/cTnI Fast Test Kit (Immunofluorescence Assay)
6	CK-MB/cTnI/Myo Fast Test Kit (Immunofluorescence Assay)
7	D-Dimer Fast Test Kit (Immunofluorescence Assay)
8	PCT Fast Test Kit (Immunofluorescence Assay)
9	CysC Fast Test Kit (Immunofluorescence Assay)
10	mAlb Fast Test Kit (Immunofluorescence Assay)
11	NGAL Fast Test Kit (Immunofluorescence Assay)
12	$\beta$ 2-MG Fast Test Kit (Immunofluorescence Assay)
13	CK-MB/cTnI Fast Test Kit (Immunofluorescence Assay)
14	HCG+ $\beta$ Fast Test Kit (Immunofluorescence Assay)
15	H-FABP Fast Test Kit (Immunofluorescence Assay)
16	PCT/CRP Fast Test Kit (Immunofluorescence Assay)
17	CK-MB/cTnI/H-FABP Fast Test Kit (Immunofluorescence Assay)
18	HbA1c Fast Test Kit (Immunofluorescence Assay)
19	NT-proBNP/NGAL Fast Test Kit (Immunofluorescence Assay)
20	CK-MB Fast Test Kit (Immunofluorescence Assay)
21	hs-cTnI Fast Test Kit (Immunofluorescence Assay)
22	T3 Fast Test Kit (Immunofluorescence Assay)
23	T4 Fast Test Kit (Immunofluorescence Assay)
24	TSH Fast Test Kit (Immunofluorescence Assay)
25	Scr Fast Test Kit (Immunofluorescence Assay)
26	PLGF Fast Test Kit (Immunofluorescence Assay)



- 27 HCY Fast Test Kit (Immunofluorescence Assay)
- 28 Anti-CCP Fast Test Kit (Immunofluorescence Assay)
- 29 25-OH-VD Fast Test Kit (Immunofluorescence Assay)
- 30 Lp-PLA2 Fast Test Kit (Immunofluorescence Assay)
- 31 FOB Fast Test Kit (Immunofluorescence Assay)
- 32 SAA Fast Test Kit (Immunofluorescence Assay)
- 33 H. pylori Fast Test Kit (Immunofluorescence Assay)
- 34 PRL Fast Test Kit (Immunofluorescence Assay)
- 35 Transferrin Fast Test Kit (Immunofluorescence Assay)
- 36 Insulin Fast Test Kit (Immunofluorescence Assay)
- 37 PG I /PG II Fast Test Kit (Immunofluorescence Assay)
- 38 LH Fast Test Kit (Immunofluorescence Assay)
- 39 FSH Fast Test Kit (Immunofluorescence Assay)
- 40 Anti-TP Fast Test Kit (Immunofluorescence Assay)
- 41 AFP/CEA Fast Test Kit (Immunofluorescence Assay)
- 42 AMH Fast Test Kit (Immunofluorescence Assay)
- 43 fT3 Fast Test Kit (Immunofluorescence Assay)
- 44 fT4 Fast Test Kit (Immunofluorescence Assay)
- 45 Total IgE Fast Test Kit (Immunofluorescence Assay)
- 46 Vit-B12 Fast Test Kit (Immunofluorescence Assay)
- 47 Prog Fast Test Kit (Immunofluorescence Assay)
- 48 Testosterone Fast Test Kit (Immunofluorescence Assay)
- 49 E2 Fast Test Kit (Immunofluorescence Assay)
- 50 RF Fast Test Kit (Immunofluorescence Assay)
- 51 ASO Fast Test Kit (Immunofluorescence Assay)
- 52 Ferritin Fast Test Kit (Immunofluorescence Assay)
- 53 ST2 Fast Test Kit (Immunofluorescence Assay)
- 54 CA125 Fast Test Kit (Immunofluorescence Assay)
- 55 CA19-9 Fast Test Kit (Immunofluorescence Assay)
- 56 CA15-3 Fast Test Kit (Immunofluorescence Assay)
- 57 RSV/Influenza A/B Fast Test Kit (Immunofluorescence Assay)
- 58 Influenza A/B Fast Test Kit (Immunofluorescence Assay)
- 59 RSV Fast Test Kit (Immunofluorescence Assay)
- 60 IL-6 Fast Test Kit (Immunofluorescence Assay)
- 61 BNP Fast Test Kit (Immunofluorescence Assay)
- 62 SAA/CRP Fast Test Kit (Immunofluorescence Assay)
- 63 Folate acid Fast Test Kit (Immunofluorescence Assay)
- 64 hs-CRP Fast Test Kit (Immunofluorescence Assay)
- 65 TnT Fast Test Kit (Immunofluorescence Assay)
- 66 PCT/IL-6 Fast Test Kit (Immunofluorescence Assay)



- 67 HBP Fast Test Kit (Immunofluorescence Assay)
- 68 S100-β Fast Test Kit (Immunofluorescence Assay)
- 69 CK-MB/hs-cTnl/Myo Fast Test Kit (Immunofluorescence Assay)
- 70 Cortisol Fast Test Kit (Immunofluorescence Assay)
- 71 CEA Fast Test Kit (Immunofluorescence Assay)
- 72 AFP/CEA Fast Test Kit (Immunofluorescence Assay)

**Classification** Other device (according to Annex II of the directive 98/79/EC)

**Conformity assessment route** Annex III of the 98/79/EC

<b>Applicable coordination standards</b>	EN 13612:2002	EN ISO 14971:2019	EN ISO15223-1:2016
	EN ISO 18113-1:2011	EN ISO 18113-2:2011	EN ISO 18113-3:2011
	EN ISO 23640:2015	EN ISO 13485:2016	ISO 780:2015
	EN 61326-2-6:2006	IEC 61326-1:2013	
	EN 61010-2-101:2002	IEC 61010-1:2010	

Signatory representative declares herein the above-mentioned device meets the basic requirements of the European Parliament and the Council's in vitro diagnostic medical devices directive: 98/79/EC Annex I.

This declaration of conformity is based on European Parliament and the Council's 98/79/EC directive Annex III. The compiled technical file and quality system document according to 98/79/EC directive Annex III are testified and the quality system certificate has issued by BSI Group The Netherlands B. V. The manufacturer is exclusively responsible for the declaration of conformity.

**General Manager** Enben Su

*Nanjing*  
13<sup>th</sup>, May, 2022  
 (place and date of issue)

\_\_\_\_\_  
 (name and signature or equivalent marking of authorized person)



**CE**



# CERTIFICATE

*Getein Biotech*

hereby certifies

**Mr. Vitalie Goreacii**

**from Sanmedico SRL.**

Completion of Getein Products Technical and Operational Training  
& Qualification of After-sales Service

基蛋生物科技股份有限公司  
GETEIN BIOTECH, INC.







# cTnI Control

REF QC001

User Manual

## PRODUCT NAME

cTnI Control

## PRODUCT SPECIFICATION

Package: 3(Level)\*2(Vial)\*1(ml), 3(Level)\*1(Vial)\*1(ml)  
cTnI Control - Level 1/2/3

## INTENDED USE

This product is intended for *in vitro* diagnostic use in the quality control of Cardiac Troponin I (cTnI) on the Getein Platforms.

## PRINCIPLE

The lyophilized cTnI control is prepared from dissolving stable and high quality recombinant cTnI antigen into calf serum. With matching equipments and reagents, it can fulfill value transfer work. As different equipments and reagents have uncertainty to some extent, different control results may appear.

## CONTENTS

The kit for FIA8000/FIA8600/Getein1100 contains:

1. cTnI Control - Level 1  
cTnI Control - Level 2  
cTnI Control - Level 3
2. User manual: 1 piece/box
3. Target value sheet: 1 piece/box

The kit for Getein1600 contains:

1. cTnI Control - Level 1  
cTnI Control - Level 2  
cTnI Control - Level 3
2. User manual: 1 piece/box
3. Target value sheet: 1 piece/box
4. Quality control holder - Level 1  
Quality control holder - Level 2  
Quality control holder - Level 3

**Note:** Each quality control holder is labelled with barcode which contains target value and level of different items.

## MATCHING EQUIPMENTS

FIA8000/8600 Quantitative Immunoassay Analyzer  
Getein1100/1600 Immunofluorescence Quantitative Analyzer

## STORAGE AND STABILITY

**UNOPENED:** The product is stable for 18 months at -20°C and for 30 days at 2 ~ 8°C to avoid light.

**OPENED:** The product is stable for 1 day at 2 ~ 8°C if kept capped in original container and free from contamination. Only the required amount of product should be removed. Any residual product should NOT BE RETURNED to the original vial after using. It is recommended to be dispensed into smaller vials after dilution and stable for 30 days at -20 ~ -70°C.

## MATERIALS REQUIRED BUT NOT PROVIDED

1. 1 ml pipette
2. Distilled water
3. Getein test kit
4. Getein instrument

## TEST PROCEDURE

1. The product should be brought to room temperature (15 ~ 30°C) prior to use.
2. Open the vial carefully in case of the loss of content.
3. Dissolve each control material with 1 ml distilled water.

4. Close the vial and mix gently until all contents are dissolved completely. Avoid violent shaking or foam formation.
5. Keep it at room temperature for 5 ~ 10 minutes before use.  
For FIA8000/FIA8600/Getein1100:
6. Treat the control in the same manner as patient specimen in the assay procedure. Follow the directions of test kit and the instrument application instruction.

For Getein1600:

7. Insert quality control holder into sample holder.
8. Insert sample holder with a constant speed and barcode facing the scanner, refer to the User Manual of Getein1600 to start QC testing.

## ASSIGNED VALUES

Refer to values listed on the target value sheet.  
If the result is beyond the range, it indicates the existence of some unreliable factors in the testing system. Referring to the control graph helps judge the accuracy and stability of the testing system.  
The expected range of the mean is provided to aid laboratory until it has established its own mean and SD for its methods.

## PERFORMANCE CHARACTERISTICS

1. Homogeneity:  $\leq 15\%$
2. Accuracy range: Refer to the target value sheet

## LIMITATIONS













1. This product can only be used on the Getein Platforms.
2. Variation exists between different equipments developed by different methods even using the same control product.
3. This product is not intended to be used as standard material.

## NOTES

1. For *in vitro* diagnostic use only.
2. Do not use the product beyond the expiration date.
3. Avoid multiple freeze-thaw cycles.
4. Do not use the product if it is contaminated with bacteria.
5. Proper handling and disposal methods should be followed in accordance with local regulations.

## DESCRIPTION OF SYMBOLS USED

The following graphical symbols used in or found on cTnI control are the most common ones appearing on medical devices and their packaging. They are explained in more details in the European Standard EN 980:2008 and EN ISO15223-1:2016.

Key to symbols used			
	Manufacturer		Expiration date
	Catalogue number		Date of manufacture
	Consult instructions for use		Batch code
	Temperature limitation		<i>In vitro</i> diagnostic medical device
	Sufficient for		Biological risk
	CE mark		Authorized representative in the European Community

Please read this user manual carefully before operating to ensure proper use.

Version: WZK01-S-04

 Getein Biotech, Inc.  
Add: No.9 Bofu Road, Luhe District, Nanjing, 211505, China  
Tel: +86-25-68568508  
Fax: +86-25-68568500  
E-mail: tech@getein.com.cn  
overseas@getein.com.cn  
Website: www.bio-GP.com.cn

Please contact Getein if you have any questions.



# Cardiac Troponin I

## Fast Test Kit

(Immunofluorescence Assay)

Getein1100: Cat.# IF1001

Getein1600: Cat.# IF2001

### User Manual

### INTENDED USE

Cardiac Troponin I Fast Test Kit (Immunofluorescence Assay) is intended for *in vitro* quantitative determination of Cardiac Troponin I (cTnI) in serum, plasma or whole blood. This test is used as an aid in the diagnosis of myocardial injury such as Acute Myocardial Infarction (AMI), Unstable Angina, Acute Myocarditis and Acute Coronary Syndrome (ACS).

### SUMMARY

Troponin, a molecular complex that is bound to the thin filament (actin) of striated muscle fibers, acts with intracellular calcium to control the interaction of the thin filament with the thick filament (myosin), thus regulating muscle contraction. Troponin consists of three regulatory proteins: T, which connects the troponin complex and tropomyosin (another cardiac muscle regulatory protein); I, which prevents muscle contraction in the absence of calcium; C, which binds calcium. Cardiac Troponin I (MW 22.5 kDa) and the two skeletal muscle isoforms of Troponin I have considerable amino acid sequence homology, but cTnI contains an additional N-terminal sequence and is highly specific for myocardium.

Clinical studies have demonstrated the release of cTnI into the blood stream within hours following acute myocardial infarction (AMI) or ischemic damage. Elevated levels of cTnI are detectable in blood within 4 to 6 hours after the onset of chest pain, reaching peak concentrations in approximately 8 to 28 hours, and remain elevated for 3 to 10 days following AMI. Due to the high myocardial specificity and the long duration of elevation, cTnI has become an important marker in the diagnosis and evaluation of patients suspected of having an AMI.

The current guideline of The Joint European Society of Cardiology/ American College of Cardiology Committee support the use of cTnI as a preferred marker of myocardial injury. Several major studies have shown that cTnI is also a predictor of cardiac risk in patients with unstable angina. The American College of Cardiology and the American Heart Association's current

guidelines recommend using troponin results when making treatment decisions regarding unstable angina and non-ST segment elevation MI (NSTEMI).

### PRINCIPLE

The test uses an anti-human cTnI monoclonal antibody conjugated with fluorescence latex and another anti-human cTnI monoclonal antibody coated on the test line. After the sample has been applied to the test strip, the fluorescence latex-labelled anti-human cTnI monoclonal antibody binds with the cTnI in sample and forms a marked antigen-antibody complex. This complex moves to the test card detection zone by capillary action. Then marked antigen-antibody complex is captured on the test line by the anti-human cTnI monoclonal antibody. The fluorescence intensity of the test line increases in proportion to the amount of cTnI in sample.

Then insert test card into Getein1100 Immunofluorescence Quantitative Analyzer/Getein1600 Immunofluorescence Quantitative Analyzer (hereinafter referred to as Getein1100 and Getein1600), the concentration of cTnI in sample will be measured and displayed on the screen. The value will be stored in Getein1100/Getein1600 and available for downloading. The result can be easily transmitted to the laboratory or hospital information system.

### CONTENTS

- A kit for Getein1100 contains:
  - Getein cTnI test card in a sealed pouch with desiccant ..... 25
  - Disposable pipet ..... 25
  - Whole blood buffer ..... 1
  - SD card ..... 1
  - User manual ..... 1
- A kit for Getein1600 contains:
  - Sealed cartridge with 24/48 Getein cTnI test cards ..... 2
  - User manual ..... 1
  - Package specifications:
  - 2x24 tests/kit, 2x48 tests/kit
  - Materials required for Getein1600:
  - Sample diluent ..... 1
  - Box with pipette tips ..... 1
  - Mixing plate ..... 1
- Sample diluent/Whole blood buffer composition:
  - Phosphate buffered saline, proteins, detergent, preservative, stabilizer.
- A test card consists of:
  - A plastic shell and a reagent strip which is composed of a sample pad, nitrocellulose membrane (one end of the membrane is coated with a fluorescence latex-labelled anti-

human cTnI monoclonal antibody, the test line is coated with another anti-human cTnI monoclonal antibody, and the control line is coated with rabbit anti-mouse IgG antibody), absorbent paper and liner.

**Note: Do not mix or interchange different batches of kits.**

### APPLICABLE DEVICE

Getein1100 Immunofluorescence Quantitative Analyzer  
Getein1600 Immunofluorescence Quantitative Analyzer

### STORAGE AND STABILITY

Store the test card at 4~30°C with a valid period of 24 months. Use the test card for Getein1100 within 1 hour once the foil pouch is opened.

For test card of Getein1600: if the cartridge is opened, it could be stable within 24 hours once exposed to air. If the test cards can't be used up at a time, please put the cartridge back to the foil pouch and reseal along the entire edge of zip-seal. The remaining test cards should be used up within 7 days.

Store the sample diluent/whole blood buffer at 0~30°C with a valid period of 24 months.

Store the sample diluent/whole blood buffer at 2~8°C for better results.

### PRECAUTIONS

- For *in vitro* diagnostic use only.
- For professional use only.
- Do not use the kit beyond the expiration date.
- Do not use the test card if the foil pouch or the cartridge is damaged.
- Do not open pouches or the cartridge until ready to perform the test.
- Do not reuse the test card.
- Do not reuse the pipet.
- Handle all specimens as potentially infectious. Proper handling and disposal methods should be followed in accordance with local regulations.
- Carefully read and follow user manual to ensure proper test performance.

### SPECIMEN COLLECTION AND PREPARATION

- This test can be used for *serum, plasma and whole blood samples*. *Heparin and sodium citrate* should be used as the anticoagulant for plasma and whole blood. Samples should be free of hemolysis.

- Suggest using serum or plasma for better results.
- Serum or plasma can be used directly. For whole blood sample, one drop of whole blood buffer must be added before testing.
- If testing will be delayed, serum and plasma samples may be stored up to 7 days at 2–8°C or stored at -20°C for 6 months before testing (whole blood sample may be stored up to 3 days at 2–8°C).
- Refrigerated or frozen sample should reach room temperature and be homogeneous before testing. Avoid multiple freeze-thaw cycles.
- Do not use heat-inactivated samples.
- SAMPLE VOLUME (for Getein1100): 100 µl.

## TEST PROCEDURE

- Collect specimens according to user manual.
- Test card, sample and reagent should be brought to room temperature before testing.

### For Getein1100:

- Confirm SD card lot No. in accordance with test kit lot No.. Perform "SD Card Calib" calibration when necessary (Details refer to 8.5.2 of Getein1100 User Manual).
- On the main interface of Getein1100, press "ENT" button to enter testing interface.
- Remove the test card from the sealed pouch immediately before use. Label the test card with patient or control identification.
- Put the test card on a clean table, horizontally placed.
- Using sample transfer pipette, deliver 100 µl of sample (or 3–4 drops of sample when using disposable pipet) into the sample port on the test card (for whole blood sample, one drop of whole blood buffer must be added after loading 100 µl sample on the test card).
- Reaction time: 10 minutes.** Insert the test card into Getein1100 and press "ENT" button after reaction time is elapsed. The result will be shown on the screen and printed automatically.

### For Getein1600:

- Each cartridge for Getein1600 contains a specific RFID card which can calibrate automatically.
- Place samples in the designed area of the sample holder, insert the holder and select the right test item, Getein1600 will do the testing and print the result automatically.

### Notes:

- It is required to perform "SD Card Calib" calibration when using a new batch of kits.
- It is suggested to calibrate once for one batch of kits for Getein1100.
- Make sure the test card and the sample insertion is correct and complete.

## TEST RESULTS

Getein1100/Getein1600 can scan the test card automatically and display the result on the screen. For additional information, please refer to the user manual of Getein1100/Getein1600.

## EXPECTED VALUE

The expected normal value for Troponin I was determined by testing samples from 500 apparently healthy individuals. The 99<sup>th</sup> percentile of the concentration for cTnI is 0.1 ng/ml. (The probability that value of a normal person below 0.1 ng/ml is 99%.)

It is recommended that each laboratory establish its own expected values for the population it serves.

## PERFORMANCE CHARACTERISTICS

Measuring Range	0.1–50 ng/ml
Lower Detection Limit	≤ 0.1 ng/ml
Within-Run Precision	≤10%
Between-Run Precision	≤15%

### Method Comparison:

The assay was compared with SIEMENS IMMULITE 2000 and its matching cTnI test kits with 200 serum samples (60 positive samples and 140 negative samples). The correlation coefficient (r) for cTnI is 0.952.

## LIMITATIONS

- As with all diagnostic tests, a definitive clinical diagnosis should not be made based on the result of a single test. The test results should be interpreted considering all other test results and clinical information such as clinical signs and symptoms.
- Samples containing interferents may influence the results. The table below listed the maximum allowance of these potential interferents.

Interferent	Hemoglobin	Triglyceride	Bilirubin
Concentration (Max)	5 g/L	10 g/L	0.2 g/L

## REFERENCES

- Mauro Pantaghini. Undefined International Federation of Clinical Chemistry and Laboratory Medicine (IFCC). Scientific Division Committee on Standardization of Markers of Cardiac Damage. Clin Chem Lab Med, 1998, 36:887–893.
- Antman EM, Anbe DT, Armstrong PW, et al. ACC/AHA guidelines for the management of patients with ST-elevation myocardial infarction: a report of the American College of Cardiology/American Heart Association Task Force on Practice

Guidelines (Committee to Revise the 1999 Guidelines for the Management 2004).

- EN ISO 18113-1:2009 *In vitro* diagnostic medical devices - Information supplied by the manufacturer (labelling) - Part 1: Terms, definitions and general requirements.
- EN ISO 18113-2:2009 *In vitro* diagnostic medical devices - Information supplied by the manufacturer (labelling) - Part 2: *In vitro* diagnostic reagents for professional use (ISO 18113-2:2009).

## DESCRIPTION OF SYMBOLS USED

The following graphical symbols used in or found on Cardiac Troponin I Fast Test Kit (Immunofluorescence Assay) are the most common ones appearing on medical devices and their packaging. They are explained in more details in the European Standard EN 980:2008 and International Standard ISO 15223-1:2007.

Key to symbols used			
	Manufacturer		Expiration date
	Do not reuse		Date of manufacture
	Consult instructions for use		Batch code
	Temperature limitation		<i>In vitro</i> diagnostic medical device
	Sufficient for		Authorized representative in the European Community
	CE mark		Do not use if package is damaged

Thank you for purchasing Cardiac Troponin I Fast Test Kit (Immunofluorescence Assay). Please read this user manual carefully before operating to ensure proper use.

Version: WIF02-S-02



Getein Biotech, Inc.  
 Add: No.9 Bofu Road, Luhe District, Nanjing, 211505, China  
 Tel: +86-25-68568508  
 Fax: +86-25-68568500  
 E-mail: tech@getein.com.cn  
 overseas@getein.com.cn  
 Website: www.bio-GP.com.cn



# D-Dimer Control

REF QC006

**User Manual**

## PRODUCT NAME

D-Dimer Control

## PRODUCT SPECIFICATION

Package: 3(Level)\*2(Vial)\*1(ml), 3(Level)\*1(Vial)\*1(ml)  
D-Dimer Control - Level 1/2/3

## INTENDED USE

This product is intended for *in vitro* diagnostic use in the quality control of D-Dimer on the Getein Platforms.

## PRINCIPLE

The lyophilized D-Dimer control is prepared from dissolving stable and high quality recombinant D-Dimer antigen into calf serum. With matching equipments and reagents, it can fulfill value transfer work. As different equipments and reagents have uncertainty to some extent, different control results may appear.

## CONTENTS

The kit for FIA8000/FIA8600/Getein1100 contains:

1. D-Dimer Control - Level 1  
D-Dimer Control - Level 2  
D-Dimer Control - Level 3
2. User manual: 1 piece/box
3. Target value sheet: 1 piece/box

The kit for Getein1600 contains:

1. D-Dimer Control - Level 1  
D-Dimer Control - Level 2  
D-Dimer Control - Level 3
2. User manual: 1 piece/box
3. Target value sheet: 1 piece/box
4. Quality control holder - Level 1  
Quality control holder - Level 2  
Quality control holder - Level 3

**Note:** Each quality control holder is labelled with barcode which contains target value and level of different items.

## MATCHING EQUIPMENTS

FIA8000/8600 Quantitative Immunoassay Analyzer  
Getein1100/1600 Immunofluorescence Quantitative Analyzer

## STORAGE AND STABILITY

**UNOPENED:** The product is stable for 18 months at -20°C and for 90 days at 2 ~ 8°C to avoid light.

**OPENED:** The product is stable for 15 days at 2 ~ 8°C if kept capped in original container and free from contamination. Only the required amount of product should be removed. Any residual product should NOT BE RETURNED to the original vial after using. It is recommended to be dispensed into smaller vials after dilution and stable for 30 days at -20 ~ -70°C.

## MATERIALS REQUIRED BUT NOT PROVIDED

1. 1 ml pipette
2. Distilled water
3. Getein test kit
4. Getein instrument

## TEST PROCEDURE

1. The product should be brought to room temperature (15 ~ 30°C) prior to use.
2. Open the vial carefully in case of the loss of content.

- Dissolve each control material with 1 ml distilled water.
- Close the vial and mix gently until all contents are dissolved completely. Avoid violent shaking or foam formation.
- Keep it at room temperature for 5 ~ 10 minutes before use.

For FIA8000/FIA8600/Getein1100:

- Treat the control in the same manner as patient specimen in the assay procedure. Follow the directions of test kit and the instrument application instruction.

For Getein1600:

- Insert quality control holder into sample holder.
- Insert sample holder with a constant speed and barcode facing the scanner, refer to the User Manual of Getein1600 to start QC testing.

## ASSIGNED VALUES

Refer to values listed on the target value sheet.

If the result is beyond the range, it indicates the existence of some unreliable factors in the testing system. Referring to the control graph helps judge the accuracy and stability of the testing system.

The expected range of the mean is provided to aid laboratory until it has established its own mean and SD for its methods.

## PERFORMANCE CHARACTERISTICS

- Homogeneity:  $\leq 15\%$
- Accuracy range: Refer to the target value sheet

## LIMITATIONS

- This product can only be used on the Getein Platforms.
- Variation exists between different equipments developed by different methods even using the same control product.
- This product is not intended to be used as standard material.













## NOTES

- For *in vitro* diagnostic use only.
- Do not use the product beyond the expiration date.
- Avoid multiple freeze-thaw cycles.
- Do not use the product if it is contaminated with bacteria.

- Proper handling and disposal methods should be followed in accordance with local regulations.

## DESCRIPTION OF SYMBOLS USED

The following graphical symbols used in or found on D-Dimer control are the most common ones appearing on medical devices and their packaging. They are explained in more details in the European Standard EN 980:2008 and EN ISO15223-1:2016.

Key to symbols used			
	Manufacturer		Expiration date
	Catalogue number		Date of manufacture
	Consult instructions for use		Batch code
	Temperature limitation		<i>In vitro</i> diagnostic medical device
	Sufficient for		Biological risk
	CE mark		Authorized representative in the European Community

Please read this user manual carefully before operating to ensure proper use.

Version: WZK04-S-04



Getein Biotech, Inc.

Addr: No.9 Bofu Road, Luhe District, Nanjing, 211505, China

Tel: +86-25-68568508

Fax: +86-25-68568500

E-mail: tech@getein.com.cn

overseas@getein.com.cn

Website: www.bio-GP.com.cn

Please contact Getein if you have any questions.



# D-Dimer Fast Test Kit

(Immunofluorescence Assay)

Getein1100: Cat.# IF1006  
Getein1600: Cat.# IF2006

## User Manual

### INTENDED USE

D-Dimer Fast Test Kit (Immunofluorescence Assay) is intended for *in vitro* quantitative determination of D-Dimer in plasma or whole blood. The test is used as an aid in the assessment and evaluation of patients suspected of deep-vein thrombosis or pulmonary embolism.

### SUMMARY

Deep-vein thrombosis is a common condition, with a lifetime cumulative incidence of 2 to 5 percent. Untreated deep-vein thrombosis can result in pulmonary embolism, a potentially fatal outcome. Anticoagulant therapy reduces both morbidity and mortality from venous thromboembolism, and early diagnosis is therefore important. Accurate diagnosis of deep-vein thrombosis minimizes the risk of thromboembolic complications and averts the exposure of patients without thrombosis to the risks of anticoagulant therapy.

D-Dimer is a marker of endogenous fibrinolysis and should therefore be detectable in patients with deep-vein thrombosis. In recent years, an increasing number of studies have shown the D-Dimer assay has a high negative predictive value and D-Dimer is a sensitive but nonspecific marker of deep-vein thrombosis. Negative D-Dimer can exclude deep-vein thrombosis and pulmonary embolism.

### PRINCIPLE

The test uses an anti-human D-Dimer monoclonal antibody conjugated with fluorescence latex and another anti-human D-Dimer monoclonal antibody coated on the test line. After the sample has been applied to the test strip, the fluorescence latex-labelled anti-human D-Dimer monoclonal antibody binds with the D-Dimer in sample and forms a marked antigen-antibody

complex. This complex moves to the test card detection zone by capillary action. Then marked antigen-antibody complex is captured on the test line by another anti-human D-Dimer monoclonal antibody. The fluorescence intensity of the test line increases in proportion to the amount of D-Dimer in sample. Then insert test card into Getein1100 Immunofluorescence Quantitative Analyzer/Getein1600 Immunofluorescence Quantitative Analyzer (hereinafter referred to as Getein1100 and Getein1600), the concentration of D-Dimer in sample will be measured and displayed on the screen. The value will be stored in Getein1100/Getein1600 and available for downloading. The result can be easily transmitted to the laboratory or hospital information system.

### CONTENTS

- A kit for Getein1100 contains:
  - Getein D-Dimer test card in a sealed pouch with desiccant ..... 25
  - Disposable pipet ..... 25
  - Sample diluent ..... 25
  - SD card ..... 1
  - User manual ..... 1
- A kit for Getein1600 contains:
  - Sealed cartridge with 24/48 Getein D-Dimer test cards ..... 2
  - User manual ..... 1
  - Package specifications:
  - 2×24 tests/kit, 2×48 tests/kit
  - Materials required for Getein1600:
  - Sample diluent ..... 1
  - Box with pipette tips ..... 1
  - Mixing plate ..... 1
- Sample diluent composition:
  - Phosphate buffered saline, proteins, detergent, preservative, stabilizer.
- A test card consists of:
  - A plastic shell and a reagent strip which is composed of a sample pad, nitrocellulose membrane (one end of the membrane is coated with a fluorescence latex-labelled anti-human D-Dimer monoclonal antibody, the test line is coated with another anti-human D-Dimer monoclonal antibody and the control line is coated with rabbit anti-mouse IgG antibody), absorbent paper and liner.

**Note: Do not mix or interchange different batches of kits.**

### APPLICABLE DEVICE

Getein1100 Immunofluorescence Quantitative Analyzer  
Getein1600 Immunofluorescence Quantitative Analyzer

### STORAGE AND STABILITY

Store the test card at 4~30°C with a valid period of 24 months. Use the test card for Getein1100 within 1 hour once the foil pouch is opened.

For test card of Getein1600: if the cartridge is opened, it could be stable within 24 hours once exposed to air. If the test cards can't be used up at a time, please put the cartridge back to the foil pouch and reseal along the entire edge of zip-seal. The remaining test cards should be used up within 7 days.

Store the sample diluent/whole blood buffer at 0~30°C with a valid period of 24 months.

Store the sample diluent/whole blood buffer at 2~8°C for better results.

### PRECAUTIONS

- For *in vitro* diagnostic use only.
- For professional use only.
- Do not use the kit beyond the expiration date.
- Do not use the test card if the foil pouch or the cartridge is damaged.
- Do not open pouches or the cartridge until ready to perform the test.
- Do not reuse the test card.
- Do not reuse the pipet.
- Handle all specimens as potentially infectious. Proper handling and disposal methods should be followed in accordance with local regulations.
- Carefully read and follow user manual to ensure proper test performance.

### SPECIMEN COLLECTION AND PREPARATION

- This test can be used for *plasma and whole blood samples*. *Sodium citrate* can be used as the anticoagulant for plasma and whole blood. Samples should be free of hemolysis.
- Suggest using plasma for better results.
- If testing will be delayed, plasma sample may be stored up to 3 days at 2~8°C or stored at -20°C for 1 month before testing (whole blood sample may be stored up to 3 days at 2~8°C).
- Refrigerated or frozen sample should reach room temperature

and be homogeneous before testing. Avoid multiple freeze-thaw cycles.

5. Do not use heat-inactivated samples.
6. SAMPLE VOLUME (for *Getein1100*): 100  $\mu$ L.

## TEST PROCEDURE

1. Collect specimens according to user manual.
2. Test card, sample and reagent should be brought to room temperature before testing.

### For *Getein1100*:

3. Confirm SD card lot No. in accordance with test kit lot No.. Perform "SD Card Calib" calibration when necessary (Details refer to 8.5.2 of *Getein1100* User Manual).
4. On the main interface of *Getein1100*, press "ENT" button to enter testing interface.
5. Remove the test card from the sealed pouch immediately before use. Label the test card with patient or control identification.
6. Put the test card on a clean table, horizontally placed.
7. Using sample transfer pipette, deliver 100  $\mu$ L of sample into one tube of sample diluent, mix gently and thoroughly. Then drop 100  $\mu$ L of sample mixture (or 3-4 drops of sample when using disposable pipet) into the sample port on the test card.
8. **Reaction time: 10 minutes.** Insert the test card into *Getein1100* and press "ENT" button after reaction time is elapsed. The result will be shown on the screen and printed automatically.

### For *Getein1600*:

9. Each cartridge for *Getein1600* contains a specific RFID card which can calibrate automatically.
10. Place samples in the designed area of the sample holder, insert the holder and select the right test item, *Getein1600* will do the testing and print the result automatically.

## Notes:

1. It is required to perform "SD Card Calib" calibration when using a new batch of kits.
2. It is suggested to calibrate once for one batch of kits for *Getein1100*.
3. Make sure the test card and the sample insertion is correct and complete.

## TEST RESULTS

*Getein1100/Getein1600* can scan the test card automatically and display the result on the screen. For additional information, please refer to the user manual of *Getein1100/Getein1600*.

## EXPECTED VALUE

The expected normal value for D-Dimer was determined by testing samples from 500 apparently healthy individuals. The 95<sup>th</sup> percentile of the concentration for D-Dimer is 0.5 mg/L. (The probability that value of a normal person below 0.5 mg/L is 95%.)

It is recommended that each laboratory establish its own expected values for the population it serves.

## PERFORMANCE CHARACTERISTICS

Measuring Range	0.1~10.0 mg/L
Lower Detection Limit	$\leq$ 0.1 mg/L
Within-Run Precision	$\leq$ 10%
Between-Run Precision	$\leq$ 15%

### Method Comparison:

The assay was compared with SIEMENS CA-7000 and its matching D-Dimer test kits with 200 plasma samples (60 positive samples and 140 negative samples). The correlation coefficient (r) for D-Dimer is 0.978.

## LIMITATIONS

1. As with all diagnostic tests, a definitive clinical diagnosis should not be made based on the result of a single test. The test results should be interpreted considering all other test results and clinical information such as clinical signs and symptoms.
2. Samples containing interferences such as rheumatoid factor, human anti-mouse antibody and heterophile antibody may influence the results. In this case, results of this test should be used in conjunction with clinical findings and other tests. The table below listed the maximum allowance of these potential interferences.

Interferent	Hemoglobin	Triglyceride	Bilirubin
Concentration (Max)	5 g/L	25 g/L	0.1 g/L













## REFERENCES

1. Sarig G, Kil-Drori AJ, Chap-Marshak D, Brenner B, Drugan A. Activation of coagulation in amniotic fluid during normal human pregnancy. *Thromb Res.* 2011 Apr 18.
2. Roldán V, Marin F, Muiña B, Torregrosa JM, Hernández-Romero D, Valdés M, Vicente V, Lip GY. Plasma von Willebrand Factor Levels Are an Independent Risk Factor for Adverse Events Including Mortality and Major Bleeding in Anticoagulated

- Atrial Fibrillation Patients. *J Am Coll Cardiol.* 2011 Apr 11.
3. Sakamoto K, Yamamoto Y, Okamatsu H, Okabe M. D-dimer is helpful for differentiating acute aortic dissection and acute pulmonary embolism from acute myocardial infarction. *Hellenic J Cardiol.* 2011 Mar-Apr; 52(2):123-127.
4. EN ISO 18113-1:2009 *In vitro* diagnostic medical devices - Information supplied by the manufacturer (labelling) - Part 1: Terms, definitions and general requirements.
5. EN ISO 18113-2:2009 *In vitro* diagnostic medical devices - Information supplied by the manufacturer (labelling) - Part 2: *In vitro* diagnostic reagents for professional use (ISO 18113-2:2009).

## DESCRIPTION OF SYMBOLS USED

The following graphical symbols used in or found on D-Dimer Fast Test Kit (Immunofluorescence Assay) are the most common ones appearing on medical devices and their packaging. They are explained in more details in the European Standard EN 980:2008 and International Standard ISO 15223-1:2007.

Key to symbols used			
	Manufacturer		Expiration date
	Do not reuse		Date of manufacture
	Consult instructions for use		Batch code
	Temperature limitation		<i>In vitro</i> diagnostic medical device
	Sufficient for		Authorized representative in the European Community
	CE mark		Do not use if package is damaged

Thank you for purchasing D-Dimer Fast Test Kit (Immunofluorescence Assay). Please read this user manual carefully before operating to ensure proper use.

Version: WIF05-S-02



Getein Biotech, Inc.  
Add: No.9 Bofu Road, Luhe District, Nanjing, 211505, China  
Tel: +86-25-68568508  
Fax: +86-25-68568500  
E-mail: tech@getein.com.cn  
overseas@getein.com.cn  
Website: www.bio-GP.com.cn





(>200 mg/L).

# hs-CRP Fast Test Kit

(Immunofluorescence Assay)

Getein1100: Cat.# IF1003

Getein1600: Cat.# IF2003

## User Manual

### INTENDED USE

hs-CRP Fast Test Kit (Immunofluorescence Assay) is intended for *in vitro* quantitative determination of C-reactive protein (CRP) in serum, plasma whole blood, or fingertip blood. Measurement of CRP is useful for the detection and evaluation of infection, tissue injury and inflammatory disorders. Measurement of high sensitivity CRP (hs-CRP), when used in conjunction with traditional clinical laboratory evaluation of acute coronary syndromes (ACS), may be useful as an independent marker of prognosis for recurrent events in patients with stable coronary disease or ACS.

### SUMMARY

C-reactive protein is an acute-phase reactant that precipitated with Pneumococcal C-polysaccharide, and is a non-specific immune response component. CRP has wide distribution in our body, and is an acute-phase protein produced in the liver in response to microbial infection or tissue injury, and the hs-CRP can be used to detect lower concentrations of CRP in serum or plasma. Studies revealed hs-CRP levels seem to be correlated with Atherosclerosis and Acute Myocardial Infarction. And the hs-CRP is an inflammation "marker" for ACS patient and is helpful for primary prevention and risk assessment of cardiovascular disease. Its combination with the ratio of total cholesterol to HDL-C is more accurate than other risk factor in predicting cardiovascular disease.

The American Heart Association and US Centers for Disease Control and Prevention have advocated hs-CRP as a predictor of cardiovascular disease (CVD) to define risk groups: less than 1.0 mg/L indicates low risk, 1.0 to 3.0 mg/L means moderate risk, and the amount above 3.0 mg/L (lower than 10 mg/L) strongly suggests a high risk of CVD. Moreover, higher CRP levels are found in late pregnant women, mild inflammation and viral infections (10-40 mg/L), active inflammation, bacterial infection (40-200 mg/L), severe bacterial infections and burns

with another anti-human hs-CRP monoclonal antibody and the control line is coated with rabbit anti-mouse IgG antibody), absorbent paper and liner.

**Note: Do not mix or interchange different batches of kits.**

### APPLICABLE DEVICE

Getein1100 Immunofluorescence Quantitative Analyzer  
Getein1600 Immunofluorescence Quantitative Analyzer

### STORAGE AND STABILITY

Store the test card at 4~30°C with a valid period of 24 months. Use the test card for Getein1100 within 1 hour once the foil pouch is opened.

For test card of Getein1600: if the cartridge is opened, it could be stable within 24 hours once exposed to air. If the test cards can't be used up at a time, please put the cartridge back to the foil pouch and reseal along the entire edge of zip-seal. The remaining test cards should be used up within 7 days.

Store the sample diluent/whole blood buffer at 0~30°C with a valid period of 24 months.

Store the sample diluent/whole blood buffer at 2~8°C for better results.

### PRECAUTIONS

1. For *in vitro* diagnostic use only.
2. For professional use only.
3. Do not use the kit beyond the expiration date.
4. Do not use the test card if the foil pouch or the cartridge is damaged.
5. Do not open pouches or the cartridge until ready to perform the test.
6. Do not reuse the test card.
7. Do not reuse the pipet.
8. Handle all specimens as potentially infectious. Proper handling and disposal methods should be followed in accordance with local regulations.
9. Carefully read and follow user manual to ensure proper test performance.

### SPECIMEN COLLECTION AND PREPARATION

1. This test can be used for *serum, plasma, whole blood and fingertip blood samples*. *Heparin, sodium citrate and EDTA* can be used as the anticoagulant for plasma, whole blood and fingertip blood. Samples should be free of hemolysis.

### PRINCIPLE

The test uses an anti-human hs-CRP monoclonal antibody conjugated with fluorescence latex and another anti-human hs-CRP monoclonal antibody coated on the test line. After the sample has been applied to the test strip, the fluorescence latex-labelled anti-human hs-CRP monoclonal antibody binds with the hs-CRP in sample and forms a marked antigen-antibody complex. This complex moves to the test card detection zone by capillary action. Then marked antigen-antibody complex is captured on the test line by the anti-human hs-CRP monoclonal antibody. The fluorescence intensity of the test line increases in proportion to the amount of hs-CRP in sample.

Then insert test card into Getein1100 Immunofluorescence Quantitative Analyzer/Getein1600 Immunofluorescence Quantitative Analyzer (hereinafter referred to as Getein1100 and Getein1600), the concentration of hs-CRP in sample will be measured and displayed on the screen. The value will be stored in Getein1100/Getein1600 and available for downloading. The result can be easily transmitted to the laboratory or hospital information system.

### CONTENTS

1. A kit for Getein1100 contains:
  - Getein hs-CRP test card in a sealed pouch with desiccant ..... 25
  - Disposable pipet ..... 25
  - Sample diluent ..... 25
  - SD card ..... 1
  - User manual ..... 1
2. A kit for Getein1600 contains:
  - Sealed cartridge with 24/48 Getein hs-CRP test cards .. 2
  - User manual ..... 1
  - Package specifications:
  - 2×24 tests/kit, 2×48 tests/kit
  - Materials required for Getein1600:
  - Sample diluent ..... 1
  - Box with pipette tips ..... 1
  - Mixing plate ..... 1
3. Sample diluent composition:
  - Phosphate buffered saline, proteins, detergent, preservative, stabilizer.
4. A test card consists of:
  - A plastic shell and a reagent strip which is composed of a sample pad, nitrocellulose membrane (one end of the membrane is coated with a fluorescence latex-labelled anti-human hs-CRP monoclonal antibody, the test line is coated

- Suggest using serum or plasma for better results.
- If testing will be delayed, serum and plasma samples may be stored up to 7 days at 2–8°C or stored at -20°C for 6 months before testing (whole blood sample may be stored up to 3 days at 2–8°C).
- Refrigerated or frozen sample should reach room temperature and be homogeneous before testing. Avoid multiple freeze-thaw cycles.
- Do not use heat-inactivated samples.
- SAMPLE VOLUME (for *Getein1100*): 10 µL.

## TEST PROCEDURE

- Collect specimens according to user manual.
- Test card, sample and reagent should be brought to room temperature before testing.

### For *Getein1100*:

- Confirm SD card lot No. in accordance with test kit lot No.. Perform "SD Card Calib" calibration when necessary (Details refer to 8.5.2 of *Getein1100* User Manual).
- On the main interface of *Getein1100*, press "ENT" button to enter testing interface.
- Remove the test card from the sealed pouch immediately before use. Label the test card with patient or control identification.
- Put the test card on a clean table, horizontally placed.
- Using sample transfer pipette, deliver 10 µL of sample into one tube of sample diluent, mix gently and thoroughly. Then drop 100 µL of sample mixture (or 3–4 drops of sample mixture when using disposable pipet) into the sample port on the test card.
- Reaction time: 3 minutes.** Insert the test card into *Getein1100* and press "ENT" button after reaction time is elapsed. The result will be shown on the screen and printed automatically.

### For *Getein1600*:

- Each cartridge for *Getein1600* contains a specific RFID card which can calibrate automatically.
- Place samples in the designed area of the sample holder, insert the holder and select the right test item, *Getein1600* will do the testing and print the result automatically.

## Notes:

- It is required to perform "SD Card Calib" calibration when using a new batch of kits.
- It is suggested to calibrate once for one batch of kits for *Getein1100*.
- Make sure the test card and the sample insertion is correct and complete.

## TEST RESULTS

*Getein1100/Getein1600* can scan the test card automatically and display the result on the screen. For additional information, please refer to the user manual of *Getein1100/Getein1600*.

## EXPECTED VALUE

**hs-CRP:** The expected normal value for hs-CRP was determined by testing samples from 500 apparently healthy individuals. The 95<sup>th</sup> percentile of the concentration for hs-CRP is 3 mg/L. (The probability that hs-CRP value of a normal person below 3 mg/L is 95%.)

**CRP:** The expected normal value for CRP was determined by testing samples from 500 apparently healthy individuals. The 95<sup>th</sup> percentile of the concentration for CRP is 10 mg/L. (The probability that CRP value of a normal person below 10 mg/L is 95%.)

It is recommended that each laboratory establish its own expected values for the population it serves.

## PERFORMANCE CHARACTERISTICS

Measuring Range	0.5~200 mg/L
Lower Detection Limit	≤0.5 mg/L
Within-Run Precision	≤10%
Between-Run Precision	≤15%

### Method Comparison:

The assay was compared with HITACHI 7600/OLYMPUS AU5400 and its matching hs-CRP test kits with 200 serum samples (61 positive samples and 139 negative samples). The correlation coefficient (r) for hs-CRP is 0.941.

## LIMITATIONS

- As with all diagnostic tests, a definitive clinical diagnosis should not be made based on the result of a single test. The test results should be interpreted considering all other test results and clinical information such as clinical signs and symptoms.
- Samples containing interferences may influence the results. The table below listed the maximum allowance of these potential interferences.

Interferent	Hemoglobin	Triglyceride	Bilirubin
Concentration (Max)	10 g/L	10 g/L	0.2 g/L

## REFERENCES









- Danesh J, Whincup P, Wslker M, et al. Low grade inflammation

and coronary heart disease: prospective study and updated meta-analysis. *BJM* 2000; 321:199–204.

- Rifai N, Ridker PM. Proposed cardiovascular risk assessment algorithm using high-sensitivity C-reactive protein and lipid screening. *Clin Chem* 2001; 47:28–30.
- EN ISO 18113-1:2009 *In vitro* diagnostic medical devices - Information supplied by the manufacturer (labelling) - Part 1: Terms, definitions and general requirements.
- EN ISO 18113-2:2009 *In vitro* diagnostic medical devices - Information supplied by the manufacturer (labelling) - Part 2: *In vitro* diagnostic reagents for professional use (ISO 18113-2:2009).

## DESCRIPTION OF SYMBOLS USED

The following graphical symbols used in or found on hs-CRP Fast Test Kit (Immunofluorescence Assay) are the most common ones appearing on medical devices and their packaging. They are explained in more details in the European Standard EN 980:2008 and International Standard ISO 15223-1:2007.

Key to symbols used			
	Manufacturer		Expiration date
	Do not reuse		Date of manufacture
	Consult instructions for use	<b>LOT</b>	Batch code
	Temperature limitation	<b>IVD</b>	<i>In vitro</i> diagnostic medical device
	Sufficient for	<b>EC REP</b>	Authorized representative in the European Community
<b>CE</b>	CE mark		Do not use if package is damaged

Thank you for purchasing hs-CRP Fast Test Kit (Immunofluorescence Assay). Please read this user manual carefully before operating to ensure proper use.

Version: WIF04-S-02



Getein Biotech, Inc.  
 Add: No.9 Bofu Road, Luhe District, Nanjing, 211505, China  
 Tel: +86-25-68568508  
 Fax: +86-25-68568500  
 E-mail: tech@getein.com.cn  
 overseas@getein.com.cn  
 Website: www.bio-GP.com.cn



# NT-proBNP

## Kit de testare rapidă

### (Analiză imunofluorescentă)

**Getein1100: Cat.# IF1002**

**Getein1600: Cat.# IF2002**

#### UTILIZAREA PRECONIZATĂ

Kitul de testare rapidă NT-proBNP (test de imunofluorescență) este destinat pentru determinarea cantitativă în vitro a precursorului peptidei natriuretice de tip B N-terminal (NT-proBNP) în ser, plasmă sau sânge integral. Acest test este utilizat ca ajutor în diagnosticul clinic, prognosticul și evaluarea insuficienței cardiace (IC).

#### REZUMAT

Precursorul peptidei natriuretice de tip B N-terminal (NT-proBNP) este secretat de ventriculul cardiac stâng ca răspuns la suprasolicitarea de volum și presiune. Este un fragment N-terminal inactiv care s-a separat de pro-hormonul BNP. NT-proBNP poate fi utilizat pentru a evalua disfuncția contractilă a inimii, disfuncția diastolică și coordonarea mișcării peretelui segmentar ventricular. În plus, are sensibilitate mare și valoare predictivă negativă (>97%). Ca standard de aur recomandat de Societatea Europeană de Cardiologie, Asociația Americană a Inimii și Colegiul American de Cardiologie pentru diagnosticul și prognoza insuficienței cardiace, NT-proBNP este utilizat pentru a indica pacientul cu insuficiență cardiacă în stadiu incipient, pentru a determina nivelurile de risc de IC, monitorizarea eficienței medicale a medicamentului IC, evaluarea prognosticului pacientului cu IC și distingerea dispneei cauzate de IC de alte boli. Mai mult, NT-proBNP este un indicator de evaluare a riscului pentru Sindromul Coronarian Acut.

#### PRINCIPIU DE TESTARE

Testul utilizează un anticorp mono-clonal anti-uman NT-proBNP conjugat cu latex fluorescent și un anticorp poli-clonal anti-uman NT-proBNP acoperit pe linia de testare. După ce proba a fost aplicată pe banda de testare, anticorpul mono-clonal anti-uman NT-proBNP marcat cu latex cu fluorescență se leagă cu NT-proBNP din probă și formează un complex marcat antigen-anticorp. Acest complex se deplasează în zona de detectare a cardului de test prin acțiune capilară. Apoi complexul antigen-anticorp marcat este capturat pe linia de testare de către anticorpul poli-clonal anti-uman NT-proBNP. Fluorescența

intensitatea liniei de testare crește proporțional cu cantitatea de NT-proBNP din probă.

Apoi introduceți cardul de testare în Analizorul cantitativ de imunofluorescență Getein1100/Analizorul cantitativ de imunofluorescență Getein1600 (denumit în continuare Getein1100 și Getein1600), concentrația de NT-proBNP din probă va fi măsurată și afișată pe ecran. Valoarea va fi stocată în Getein1100/Getein1600 și disponibilă pentru descărcare. Rezultatul poate fi transmis cu ușurință către sistemul informatic al laboratorului sau al spitalului.

#### CONȚINUT

1. Un kit pentru Getein1100 conține:

Card de testare Getein NT-proBNP într-o pungă sigilată cu desicant ..... 25  
 Pipetă de unică folosință ..... 25  
 Buffer sânge integru ..... 1  
 Card SD ..... 1  
 Manualul utilizatorului ..... 1

2. Un kit pentru Getein1600 conține:

Cartuș sigilat cu carduri de testare Getein NT-proBNP 24/48 ..... 2  
 Manualul utilizatorului ..... 1  
 Caracteristicile ambalajului;  
 2x24 teste/kit, 2x48 teste/kit  
Materiale necesare pentru Getein1600  
 Diluant de mostră ..... 1  
 Cutie cu vârfuluri pentru pipetă ..... 1  
 Placă de amestecare ..... 1

3. Diluant eșanțion/Compoziție tampon pentru sânge total:

Soluție salină tamponată cu fosfat, proteine, detergent, conservant, stabilizator.

4. Un card de testare este compus din:

O carcasă de plastic și o bandă reactivă care este compusă dintr-un tampon de probă, membrana de nitroceluloză (un capăt al membranei este acoperit cu un anticorp mono-clonal anti-uman NT-proBNP marcat cu latex fluorescent, linia de testare este acoperită cu un alt anti-uman Anticorpul poli-clonal NT-proBNP și linia de control este acoperită cu anticorp IgG de iepure anti-șoarece), hârtie absorbantă și captușeală.

**Notă: Componentele din loturi diferite nu trebuie inter-schimbate.**

#### DISPOZITIVE APLICABILE

Analizor cantitativ de imunofluorescență Getein1100/ Analizor cantitativ de imunofluorescență Getein1600

#### PĂSTRARE ȘI STABILITATE

Păstrați cardul de testare la 4~30°C cu o perioadă valabilă de 24 de luni. Utilizați cardul de testare pentru Getein1100 în decurs de 1 oră după deschiderea pungii din folie.

Pentru cardul de testare al Getein1600: dacă cartușul este deschis, acesta ar putea fi stabil în 24 de ore odată expus la aer. Dacă cardurile de testare nu pot fi folosite la un moment dat, puneți cartușul înapoi în punga de folie și resigilați de-a lungul întregii margini a fermoarului. Cardurile de testare rămase ar trebui să fie epuizate în 7 zile.

Păstrați diluantul de probă/tamponul de sânge integral la 0~30°C cu o perioadă valabilă de 24 de luni.

Păstrați diluantul de probă/tamponul de sânge integral la 2-8°C pentru rezultate mai bune.

#### PRECAUȚII

1. Numai pentru diagnostic in vitro.
2. Numai pentru uz profesional.
3. Nu utilizați kitul după data de expirare.
4. Nu utilizați cardul de testare dacă punga din folie sau cartușul sunt deteriorate.
5. Nu deschideți pungile sau cartușul până când sunteți gata să efectuați testul.
6. Nu reutilizați cardul de testare.
7. Nu reutilizați pipeta.
8. Manipulați toate probele ca fiind potențial infectioase. Metodele adecvate de manipulare și eliminare trebuie urmate în conformitate cu regulamentul local.
9. Citiți cu atenție și urmați manualul pentru a vă asigura o performanță adecvată a testului.

#### COLECTAREA ȘI PREGĂTIREA MOSTRELOR

1. Acest test poate fi utilizat pentru mostre de ser, plasmă și sânge integral. Heparina ar trebui utilizată ca anticoagulant pentru plasmă și sânge integral. Probele trebuie să fie lipsite de hemoliză.
2. Se recomandă utilizarea de ser sau plasmă pentru rezultate mai bune.
3. Ser sau plasma pot fi folosite direct. Pentru proba de sânge integral, trebuie adăugată o picătură de tampon de sânge integral înainte de testare.
4. Dacă testarea va fi amânată, probele de ser și plasmă pot fi păstrate până la 1 zi la 2-8°C sau păstrate la -20°C timp de 3 luni înainte de testare (proba de sânge integral poate fi păstrată până la 3 zile la 2 ~8°C).
5. Proba refrigerată sau congelată trebuie să atingă temperatura camerei și să fie omogenă înainte de testare. Evitați ciclurile multiple de îngheț-dezghet.
6. Nu utilizați mostre inactivate la căldură.
7. VOLUM MOSTRĂ-100 ul.

deservește..

## PROCEDURA DE TESTARE

1. Colectați probe conform manualului de utilizare.
2. Cardul de testare, proba și reactivul trebuie aduse la temperatura camerei înainte de testare.

Pentru Getein1100:

3. Confirmați numărul lotului cardului SD în conformitate cu numărul lotului kit-ului de testare. Efectuați calibrarea „SD Card Calib” atunci când este necesar (Detalii se referă la 8.5.2 din Manualul utilizatorului Getein1100).

4. Pe interfața principală a Getein1100, apăsați butonul „ENT” pentru a intra în interfața de testare.

5. Scoateți cardul de test din pungă sigilată imediat înainte de utilizare. Etichetați cardul de testare cu identificarea pacientului sau a controlului.

6. Puneți cardul de test pe o masă curată, așezată orizontal.

7. Folosind pipeta de transfer a probei, furnizați 100 ul de probă (sau 3-4 picături de probă atunci când utilizați pipeta de unică folosință) în portul de probă de pe cardul de testare (pentru proba de sânge integral, trebuie adăugată o picătură de bufer de sânge integral după încărcare). volum de 100 ul pe cardul de testare).

8. Timp de reacție Introduceți cardul de test în Getein1100 și apăsați butonul „ENT” după ce timpul de reacție a trecut. Rezultatul va fi afișat pe ecran și imprimat automat.

Pentru Getein1600:

9. Fiecare cartuș pentru Getein1600 conține un card RFID specific care se poate calibra automat.

10. Așezați mostra în zona proiectată a suportului de probă, introduceți suportul și selectați elementul de testare potrivit, Getein 1600 va face testarea și va imprima rezultatul automat.

**Note:**

1. Este necesar să efectuați calibrarea „SD Card Calib” atunci când utilizați un nou lot de kituri.
2. Se recomandă calibrarea o dată pentru un lot de kituri pentru Getein1100.
3. Asigurați-vă că cardul de test și introducerea probei sunt corecte și complete.

## REZULTATELE TESTĂRII

Getein1100/Getein1600 poate scana automat cardul de testare și poate afișa rezultatul pe ecran. Pentru informații suplimentare, consultați manualul de utilizare al Getein1100/Getein1600.

## VALOAREA PRECONIZATĂ

Valoarea normală așteptată pentru NT-proBNP a fost determinată prin testarea probelor de la 2500 de indivizi aparent sănătoși. Procentajul de 95 a concentrației pentru NT-proBNP este de 185 pg/ml, iar procentajul de 97,5 a concentrației pentru NT-proBNP este de 300 pg/ml. Din cauza diferenței aparente de concentrație, se recomandă ca fiecare laborator să-și stabilească propriile valori așteptate pentru populația pe care o

**Tabelul 1 Referințe NT proBNP**

Vârsta/ Procentaj	≤44	45-54	55-64	65-74	≥75	Analiza statistica
95	98.5	130	215	290	530	185
97.5	116	170	270	350	740	300

**Tabelul 2 Standardul de excludere/diagnostic al IC de NT-proBNP.**

Vîrsta	<50	50-75	≥75	Diagnosticul IC
NT-proBNP (pg/ml)	≥450	≥900	≥1800	Probabilitatea înalta IC
	300-450	300-900	300-1800	Probabilitate scăzută de IC, trebuie combinată cu alte evaluări clinice
	<300	<300	<300	IC exclus

## CARACTERISTICI DE PERFORMANȚĂ

Spectrul de măsurare 100-35000 pg/ml

Limita inferioară de detecție ≤100 pg/ml

Precizie în timpul rulării ≤10%

Precizia între rulări ≤15%

Comparatia metodelor:

Testul a fost comparat cu Roche MODULAR ANALYTICS E170 și cu kiturile de testare NT-proBNP potrivite cu 200 de probe de ser (63 de probe pozitive și 137 de probe negative). Coeficientul de corelație (r) pentru NT-proBNP este 0,959.

## LIMITĂRI

1. Ca și în cazul tuturor testelor de diagnostic, un diagnostic clinic definitiv nu trebuie făcut pe baza rezultatului unui singur test. Rezultatele testului trebuie interpretate luând în considerare toate celelalte rezultate ale testelor și informații clinice, cum ar fi semnele și simptomele clinice.
2. Probele care conțin interferenți pot influența rezultatele. Tabelul de mai jos a enumerat permisiunea maximă a acestor interferenți potențiali.

Interferent	Hemoglobina	Trigliceride	Bilirubina
Concentrația (Max)	5 g/L	10 g/L	0.2 g/L