---

# iQ-Check *Salmonella* spp. II Kit

## User Guide

Test for the real-time PCR detection of *Salmonella* spp. in food, animal feed, and environmental samples

Catalog #3578123

**BIO-RAD**

# Table of Contents

Section 1.	Introduction.....	1
Section 2.	The iQ-Check <i>Salmonella</i> spp. II Technology .....	1
Section 3.	Kit Components.....	2
Section 4.	Shelf Life and Storage.....	2
Section 5.	Materials Required but Not Supplied.....	2
	Equipment.....	2
	Supplies .....	3
Section 6.	Safety Precautions and Recommendations for Best Results .....	4
Section 7.	Protocol.....	6
	Sample Enrichment.....	6
	Free DNA Removal Treatment.....	9
	DNA Extraction.....	9
	Real-Time PCR.....	11
	Data Analysis .....	12
Section 8.	Confirmation of Positive Results.....	13
Section 9.	Confirmation of Single Colonies Using iQ-Check Kit.....	14
Section 10.	Test Performance and Validations .....	15
Section 11.	References.....	16
Section 12.	Revision History .....	17
	Appendix — PCR Mix Calculation Guide.....	18

## Section 1

### Introduction

Salmonellae are the most frequent causes of food poisoning in the world, despite the many preventive measures taken to control these organisms. Eggs, dairy products, meat, and poultry are the foods most commonly associated with the transmission of *Salmonella* (65% of cases). More than 2,500 serotypes have been identified, all of which are potentially pathogenic to humans. Due to the low infective dose and the serious threat posed to producers and consumers of food, several countries now require total absence of *Salmonella* in food products.

Classical culture methods are often long and tedious. The iQ-Check *Salmonella* II Kit, in comparison, is a simple and rapid qualitative test allowing the detection of DNA sequences specific to *Salmonella* spp. found in food products, animal feed, and environmental samples. Using real-time PCR, *Salmonella* spp.–specific DNA sequences are amplified and detected simultaneously by means of fluorescent probes. Up to 94 samples can be processed at one time, with a minimal risk of contamination and an easy-to-use procedure. The intended users of this kit are trained laboratory personnel who are performing tests to detect *Salmonella* spp. The use of this test allows results to be obtained within a few hours following enrichment of a sample.

## Section 2

### The iQ-Check *Salmonella* spp. II Technology


The iQ-Check *Salmonella* II Kit is based on gene amplification and detection by real-time PCR. The kit's ready-to-use PCR reagents contain oligonucleotides (primers and probes) specific for *Salmonella* spp., as well as DNA polymerase and nucleotides. Detection and data analysis are optimized for use with a Bio-Rad real-time PCR instrument, such as the CFX96 Touch Deep Well System.

PCR is a powerful technique used to generate many copies of target DNA. During the PCR reaction, several cycles of heating and cooling allow DNA denaturation by heat followed by primers annealing to the target region. The DNA polymerase then uses these primers and deoxynucleotide triphosphates (dNTPs) to extend the DNA, creating copies of the target DNA. These copies are called amplicons.

In real-time PCR, specific probes are used to detect the DNA during the amplification step by hybridizing to the amplicons. These probes are linked to a fluorophore that fluoresces only when hybridized to the target sequence. FAM is the fluorophore linked to the probe hybridizing to the *Salmonella* spp.–specific DNA sequence. In the absence of target DNA, no fluorescence will be detected. As the number of amplicons increases with each round of amplification, fluorescence intensity also increases. The optical module measures this fluorescence at the annealing step during each PCR cycle while the associated software plots the fluorescence intensity versus number of cycles.

A synthetic DNA internal control is included in the reaction mix to validate any possible negative results. This control is amplified with a specific probe at the same time as the *Salmonella* target DNA sequence and is detected by a second fluorophore.

This test allows the qualitative detection of *Salmonella* spp. in select food products, animal feed, and environmental samples previously enriched by culture. It includes the following five main steps:

\* Please refer to the user guide of the iQ-Check Free DNA Removal Solution (#10000058391) for conditions of use.

## Section 3

### Kit Components

The iQ-Check *Salmonella* spp. II Kit contains sufficient reagents for 96 tests (94 samples).

Reagent ID	Reagent	Quantity Provided, ml
A	Lysis reagent	1 bottle, 20
B	Fluorescent probes	1 tube, 0.55
C	Amplification mix	1 tube, 4.4
D	PCR negative control	1 tube, 0.5
E	PCR positive control	1 tube, 0.25

## Section 4

### Shelf Life and Storage

Once received, the kit must be stored at 2–8°C. Reagents stored at this temperature can be used until the expiration date indicated on the tubes.

## Section 5

### Materials Required but Not Supplied

#### Equipment

- Lab paddle blender for homogenizing test samples
- Incubator for sample microbiological enrichment
- Specific for extraction in sterile 1.5 ml conical screwcap tubes:
  - Benchtop centrifuge capable of 10,000–12,000 x g

## Section 5 Materials Required but Not Supplied

- Dry heat block at  $37 \pm 2^\circ\text{C}$  and/or  $95\text{--}100^\circ\text{C}$
- Cell disruptor (for example, Disruptor Genie, Scientific Industries, Inc.) – for Standard II Protocol
- Specific for extraction in deep well plate:
  - Heating thermoshaker\* capable of maintaining  $37 \pm 2^\circ\text{C}$  and/or  $95\text{--}100^\circ\text{C}$ , with a mixing speed of at least 1,300 rpm
  - DW40 Deep Well Microplate Washer (catalog #90137) – for Standard II Protocol
- Vortexer
- Magnetic stir plate
- 20, 200, and 1,000  $\mu\text{l}$  micropipets
- Bio-Rad real-time PCR system,\* for example, the CFX96 Touch Deep Well (catalog #3600037) or CFX Opus Deepwell (catalog # 17007991) Real-Time PCR Systems
- Bio-Rad iQ-Check Prep System for automated DNA extraction and PCR plate setup (catalog #3594911)

**Note:** We recommend using an uninterrupted power supply (UPS) with the thermal cycler and iQ-Check Prep Systems.

\*Contact Bio-Rad Technical Support for information on recommended instruments.

## Supplies

- Enrichment medium: Buffered Peptone Water (for example, BPW Plus, catalog #3564684, dehydrated, 500 g; 3554179, 225 ml x 6 bottles; 3555795, 3 L x 4 bags; 3555790, 5 L x 2 bags BPW Standard, catalog #12013259, dehydrated, 500 g; 12013258, dehydrated, 5 kg; 12013260, 5 L x 2 bags)
- Specific for enrichment of environmental primary production samples: RAPID' *Salmonella* Capsule, (catalog #3564710, 1x strength, 100 capsules; 3564709, 10x strength, 100 capsules; 3564712, powder, 100 tests)
- Selective supplement: Novobiocin
- PIF Supplement (catalog #12013322, 2 g)
- iQ-Check Free DNA Removal Solution (catalog #3594970)
- Specific for environmental samples:
  - Environmental sponges
  - Environmental swabs
  - Neutralizing broth for sponges and swabs, such as Dey-Engley (D/E) or HiCap Neutralizing Broth, or Lethen Broth

- Specific for extraction in tubes
  - 1.5 ml conical screwcap tubes, sterile (for example catalog #2240110XTU)
- Specific for extraction in a deep well plate
  - 96-well deep well plate (iQ-Check Deep Well Microplates, catalog #3594900)
  - Plastic sealing film (TeSeE NSP Plastic Sealing Film, catalog #3590139)
  - Pre-pierced sealing film (X-Pierce Films, catalog #3593977 or Pre-Pierced Plate Sealing Film, 3600040, North America only)
- Specific for iQ-Check Prep System
  - 60 ml dilution container (catalog #3594904 or 12014473)
  - Filter tips (catalog #3594902 or 12014486, 50 µl x 5,760; 3594903 or 12014483, 1,000 µl x 3,840)
  - PCR mix tubes (catalog #12016673, 5 ml x 25)
- Specific for Standard II and Easy II Extraction Protocols:
  - Lysis beads reagent F (catalog #3578136)
  - 200 µl wide-opening pipet tips
- RAPID'*Salmonella* Agar (catalog #3563961, 90 mm x 20 dishes; 3563963 90 mm x 120 dishes; 3564705, dehydrated, 500 g)
- PCR plates, tubes, sealing tape, and caps
- Sterile filter tips adaptable to 20, 200, and 1,000 µl micropipets
- Tips for Combitip pipets or equivalent repeat pipettors; sterile, individually packaged
- 1 and 10 ml pipets
- 2 and 5 ml sterile test tubes
- Powder-free gloves
- Distilled sterile water
- Bleach, 5%
- Cleaning agent such as DNA AWAY or RNase AWAY

## Section 6

### Safety Precautions and Recommendations for Best Results

- This test must be performed by trained personnel

## Section 6 Safety Precautions and Recommendations for Best Results

- Samples and enrichment cultures must be handled as potentially infectious material and discarded according to local rules and regulations
- All potentially infectious material should be autoclaved before disposal
- The quality of results depends on strict compliance with Good Laboratory Practices (for example, the EN ISO 7218 standard), especially concerning PCR:
  - Never circulate laboratory equipment (pipets, tubes, etc.) from one workstation to another
  - Always use a positive control and a negative control for each series of amplification reactions
  - Do not use reagents after their expiration date
  - Vortex reagents from the kit before using them to ensure homogeneity
  - Periodically verify the accuracy and precision of pipets, as well as correct functioning of the instruments
  - Change gloves often, especially if you suspect they are contaminated
  - Clean workspaces periodically with 5% bleach and other decontaminating agents, such as DNA AWAY
  - Use powder-free gloves and avoid fingerprints and writing on tube caps. Both will interfere with data acquisition
- It is strongly advised to follow the general requirements described in the standard EN ISO 22174:2005 (Microbiology of food and animal feeding stuffs — Polymerase chain reaction (PCR) for the detection of food pathogens — General requirements and definitions)
- iQ-Check *Salmonella* II Kit
  - All substances or mixtures in the test kit are classified products, according to the Globally Harmonized System (GHS). Contact with acids may cause release of toxic gases. No special precautions are necessary if used correctly. If the product is inhaled, supply fresh air and consult a doctor in case of complaints. After eye contact with the product, rinse opened eye for several minutes under running water. If the products are swallowed, induce vomiting and call for medical help
- iQ-Check Prep System
  - Improper use of the iQ-Check Prep System may cause personal injury or damage to the instrument. Some components may pose a risk of personal injury due to excessive heat if improperly handled. For safe use, the iQ-Check Prep System must be operated only by qualified laboratory personnel who have been appropriately trained. Servicing of the instrument must be performed only by Bio-Rad field service engineers
- CFX96 Touch Deep Well or CFX Opus Deepwell Real-Time PCR Detection System
  - Improper use of the CFX96 Touch Deep Well or CFX Opus Deepwell Real-Time PCR Detection System may cause personal injury or damage to the instrument. Some components may pose a risk of personal injury due to excessive heat if improperly handled. For safe use, the CFX96 Touch

Deep Well and CFX Opus Deepwell Real-Time PCR Detection System must be operated only by qualified laboratory personnel who have been appropriately trained. Servicing of the instrument must be performed only by Bio-Rad field service engineers

- Enrichment
  - The user should read, understand, and follow all safety information in the instructions for the iQ-Check *Salmonella* II Kit. Retain the safety instructions for future reference. To reduce the risks associated with exposure to chemicals and biohazards, perform pathogen testing in a properly equipped laboratory under the control of trained personnel. Always follow standard laboratory safety practices, including wearing appropriate protective apparel and eye protection while handling reagents and contaminated samples. Avoid contact with the contents of the enrichment media and reagent tubes after amplification. Dispose of enriched samples according to current industry standards
  - *Salmonella* is a Biosafety Level 2 organism. Biological samples such as enrichments have the potential to transmit infectious diseases. Follow all applicable local, state/provincial, and/or national regulations on disposal of biological waste. Wear appropriate protective equipment, which includes but is not limited to: protective eyewear, face shield, clothing/lab coat, and gloves. All work should be conducted in properly equipped facilities utilizing the appropriate safety equipment (for example, physical containment devices). Individuals should be trained in accordance with applicable regulatory and company/institution requirements before working with potentially infectious materials
  - When testing is complete, all materials and media possibly containing pathogens should be decontaminated following current industry standards for the disposal of contaminated waste (that is, autoclave for 20 min at 121°C). Consult the Safety Data Sheet for additional information and local regulations for disposal

## Section 7

### Protocol

It is strongly recommended to read the entire protocol before starting the test.

#### A. Sample Enrichment

For overnight enrichment, warm media to room temperature (20–25°C) before use. For a short enrichment (8 hr), warm enrichment broth to the sample incubation temperature before use. Short enrichment times are sensitive to incubation conditions and temperatures indicated must be strictly respected. The duration of the sample preparation (time between preheating of the enrichment broth and the beginning of the food sample incubation) must not exceed 45 min. The use of a ventilated incubator is recommended.

The use of an enrichment bag with incorporated filter is highly recommended.

The use of alpha-amylase is recommended for cereals or starch product enrichment as described in the ISO 6887-4 standard.

Within the scope of the NF Validation, respect ISO standards 6887-2 to -6 for sample preparation. If not indicated, test portions weighing more than 25 g have not been tested. Incubate without shaking

## Section 7 Protocol

according to the times and temperatures indicated in the table. Standard I DNA Extraction Protocol can be used for all matrices except environmental primary production samples.

The following table outlines the different protocols that can be used, depending on the application and the scope of the validation.

NF VALIDATION BRD 07/06 – 07/04 and NordVal 038		
Scope (Matrices)	Sample Preparation	Enrichment/DNA Extraction
Food products, animal feed, & environmental samples <sup>1</sup>	Small sample sizes: 1:10 dilution <ul style="list-style-type: none"> <li>Homogenize <i>n</i> g food in 9 x <i>n</i> ml BPW (25 g in 225 ml)</li> </ul>	<b>Standard I extraction</b> <ul style="list-style-type: none"> <li>Incubate 16–24 hr at 37 ± 1°C</li> <li>Tube format</li> </ul> <b>Easy I extraction</b> <ul style="list-style-type: none"> <li>Incubate 20–24 hr at 37 ± 1°C</li> <li>Tube/deep well format</li> </ul>
Raw meat products (25 g, 375 g) <sup>2</sup>	25 g sample sizes: 1:10 dilution <ul style="list-style-type: none"> <li>Homogenize 25 g raw meat in 225 ml prewarmed BPW</li> </ul>	Easy II extraction for 25 g samples <ul style="list-style-type: none"> <li>Incubate 8-16 hr at 41.5 ± 1°C</li> <li>Tube/deep well format</li> </ul>
	375 g sample sizes: 1:4 dilution <ul style="list-style-type: none"> <li>Homogenize 375 g raw meat in 1,125 ml prewarmed BPW</li> </ul>	Easy II extraction for 375 g samples <ul style="list-style-type: none"> <li>Incubate 10–18 hr at 41.5 ± 1°C</li> <li>Tube/deep well format</li> </ul>
Environmental primary production samples	Homogenize <i>n</i> g sample in 9 x <i>n</i> ml BPW supplemented with RAPID <sup>®</sup> <i>Salmonella</i> capsule	Standard II <sup>3</sup> or Easy I extraction <ul style="list-style-type: none"> <li>Incubate 18–26 hr at 41.5 °C</li> <li>Transfer 100 µl into 900 µl fresh BPW</li> <li>Incubate 4–6 hr at 37 ± 1°C</li> </ul>
Infant formula & cereal with/without probiotics, ingredients (50-375 g) <sup>4</sup>	Homogenize <i>n</i> g food in 3 x <i>n</i> ml BPW with PIF Supplement (375 g in 1,125 ml)	Easy I extraction <ul style="list-style-type: none"> <li>Incubate 18–26 hr at 37 ± 1°C</li> <li>Tube/deep well format</li> </ul>
Raw dairy products (up to 50 g)	Homogenize <i>n</i> g food in 9 x <i>n</i> ml BPW supplemented with 20 mg/L novobiocin (50 g in 450 ml)	Easy I extraction <ul style="list-style-type: none"> <li>Incubate 20–24 hr at 41.5 ± 1°C</li> <li>Tube/deep well format</li> </ul>
Pet food and animal feed (375g)	Homogenize <i>n</i> g sample in 5 x <i>n</i> ml of pre warmed BPW supplemented with RAPID <sup>®</sup> <i>Salmonella</i> capsule	Easy I extraction <ul style="list-style-type: none"> <li>Incubate 18–24 hr at 37 ± 1°C</li> <li>Tube/deep well format</li> </ul>
AOAC OMA 2017.06		
Scope (Matrices)	Sample Preparation	Enrichment/DNA Extraction
Selected food products & environmental surfaces	Small sample sizes: 1:10 dilution <ul style="list-style-type: none"> <li>Homogenize <i>n</i> g food in 9 x <i>n</i> ml BPW (25 g in 225 ml)</li> </ul>	Easy I extraction <ul style="list-style-type: none"> <li>Incubate 19–23 hr at 37 ± 1°C</li> <li>Tube/deep well format</li> </ul>
	Large sample sizes: 1:4 dilution <ul style="list-style-type: none"> <li>Homogenize <i>n</i> g food in 3 x <i>n</i> ml BPW (375 g in 1,125 ml)</li> </ul>	
	Environmental samples <ul style="list-style-type: none"> <li>Moisten swabs and sponges with a neutralizing broth that does not contain aryl sulfonate complex</li> <li>For surfaces being analyzed with swabs, sample a 1 x 1" (2.54 x 2.54 cm) area</li> </ul>	

	<ul style="list-style-type: none"> <li>For surfaces being analyzed with sponges, sample a 4 x 4" (10.16 x 10.16 cm) area</li> <li>Add to enough BPW to cover the swab (10 ml) or sponge (100 ml)</li> </ul>	
Milk chocolate (375 g)	Homogenize 375 g milk chocolate in 3,000 ml skim milk medium	Easy I extraction <ul style="list-style-type: none"> <li>Incubate 18–22 hr at 37 ± 1°C</li> <li>Tube/deep well format</li> </ul>
Carcass rinses	Homogenize 30 ml carcass rinse in 30 ml BPW	Easy I extraction <ul style="list-style-type: none"> <li>Incubate 20–22 hr at 37 ± 1°C</li> <li>Tube/deep well format</li> </ul>
<b>AOAC PTM 010803</b>		
<b>Scope (Matrices)</b>	<b>Sample Preparation</b>	<b>Enrichment/DNA Extraction</b>
Raw ground beef and raw beef trim (375 g) <sup>5,7</sup>	Homogenize <i>n</i> g food in 3 x <i>n</i> ml pre-warmed BPW (375 g in 1,125 ml)	Easy II extraction <ul style="list-style-type: none"> <li>Incubate 8–22 hr at 41.5 ± 1°C</li> <li>Tube/deep well format</li> </ul>
Fresh raw spinach (375 g) <sup>5,7</sup>	Homogenize <i>n</i> g food in 3 x <i>n</i> ml pre-warmed BPW (375 g in 1,125 ml)	Easy II extraction <ul style="list-style-type: none"> <li>Incubate 10–22 hr at 41.5 ± 1°C</li> <li>Tube/deep well format</li> </ul>
Dairy powders (375 g) <sup>7</sup>	Homogenize <i>n</i> g food in 3 x <i>n</i> ml BPW with PIF Supplement (375 g in 1,125 ml)	Easy I extraction <ul style="list-style-type: none"> <li>Incubate 18–26 hr at 37 ± 1°C</li> <li>Tube/deep well format</li> </ul>
Chocolate liquor & white chocolate (375 g) <sup>6,7</sup>	Homogenize 375 g food in 3,000 ml prewarmed BPW	Easy I extraction <ul style="list-style-type: none"> <li>Incubate 18–26 hr at 37 ± 1°C</li> <li>Tube/deep well format</li> </ul>
Cannabis flower (10 g) <sup>5,7</sup> and hemp flower (25 g) <sup>5</sup>	Homogenize <i>n</i> g food in 9 x <i>n</i> ml BPW (10 g in 90 ml BPW cannabis flower and 25 g in 225 ml BPW hemp flower)	Easy I extraction <ul style="list-style-type: none"> <li>Incubate 20–22 hr at 37 ± 1°C</li> <li>Tube/deep well format</li> </ul>
Cannabis-infused chocolate and gummies (25 g) and cannabis-derived concentrates (5 g) <sup>5,7</sup>	Homogenize <i>n</i> g food in 9 x <i>n</i> ml BPW (25 g in 225 ml or 5 g in 45 ml).	Easy I extraction <ul style="list-style-type: none"> <li>Incubate 18–22 hr at 37 ± 1°C</li> <li>Tube/deep well format</li> </ul>
Plant-based meat products (375 g) <sup>7</sup>	Homogenize <i>n</i> g food in 3 x <i>n</i> ml BPW (375 g in 1,125 ml )	Easy I extraction <ul style="list-style-type: none"> <li>Incubate 20–22 hr at 37 ± 1°C</li> <li>Tube/deep well format</li> </ul>
All-purpose flour (375 g) <sup>5,7</sup>	1:3 dilution <ul style="list-style-type: none"> <li>Homogenize <i>n</i> g food in 3 x <i>n</i> ml prewarmed BPW + PIF Supplement (375 g in 1,125 ml + 150 µl prepared PIF Supplement)</li> <li>After enrichment, transfer a 1–5 ml aliquot to a tube or deep well plate and allow to sit for at least 30 min before proceeding to DNA extraction</li> </ul> 1:9 dilution <ul style="list-style-type: none"> <li>Homogenize <i>n</i> g food in 9 x <i>n</i> ml prewarmed BPW + PIF Supplement (375 g in 3,375 ml + 375 µl prepared PIF Supplement)</li> </ul>	Easy II extraction <ul style="list-style-type: none"> <li>Incubate 20–24 hr at 37 ± 1°C</li> <li>Tube/deep well format</li> </ul>
<b>Refer to MFLP-38 for Health Canada approved matrices and protocols</b>		

<sup>1</sup> Except environmental primary production samples

<sup>2</sup> Protocol used for simultaneous detection of *E. coli* O157:H7

<sup>3</sup> Lysis reagent used is reagent A only (without reagent F)

## Section 7 Protocol

<sup>4</sup> The protocol is harmonized with the iQ-Check *Cronobacter* spp. and iQ-Check *Enterobacteriaceae* Kits

<sup>5</sup> Harmonized protocol with iQ-Check *E. coli* O157:H7 and iQ-Check STEC VirX Kits

<sup>6</sup> Validation includes performing test with and without iQ-Check Purification Reagent

<sup>7</sup> Validation includes direct streak to *RAPID* *Salmonella* plates

### B. Free DNA Removal Treatment

The iQ-Check Free DNA Removal Solution provides an ideal way to remove free DNA. Follow Bio-Rad's recommendations in the user guide (document #10000058391).

### C. DNA Extraction

General recommendations:

1. Turn on the heat block or thermoshaker to preheat before starting the test. Set it to 95–100°C. Keep the lysis reagent in suspension while pipetting by stirring at medium speed on a magnetic stir plate.
2. In general, avoid shaking the enrichment bag and collecting large fragments of food debris. For food samples with a fatty supernatant, collect the sample just below this layer.
3. Open tubes and wells carefully to avoid any possible cross contamination.
4. Cool the deep well plate before pipetting directly through pre-pierced sealing film.
5. Use the magnetic bar to keep the lysis reagent in suspension. Pipet while it is stirring at medium speed.
6. Gently shake the lysis reagent by hand first to resuspend the resin. Pipet while the magnetic stir bar is stirring the bottle contents at medium speed in order to keep the contents in suspension.
7. For Easy II and Standard II Protocols, reconstitute the final lysis reagent as follows:
  - a. Carefully pour all the contents from reagent F (lysis beads) into reagent A (lysis reagent).
  - b. Use consumables with a wide enough tip to allow pipetting of the homogenized lysis reagent.
  - c. The lysis reagent mixed with lysis beads (reagents A + F) has a shelf life of 6 months when stored at 4°C.

#### Easy I Protocol

1. Aliquot 100 µl of homogenized lysis reagent (reagent A) to tubes or wells of a deep well plate.
2. Add 100 µl of enriched sample.
3. Mix the solution by pipetting up and down until homogenized.
4. Close the tubes or seal the deep well plate with pre-pierced sealing film.
5. Incubate tubes in the heat block at 95–100°C for 15–20 min. Incubate deep well plate in the thermoshaker under agitation at 1,300–1,600 rpm at 95–100°C for 15–20 min.

6. Vortex tubes at high speed, and then centrifuge at 10,000–12,000 x g for at least 2 min. Centrifugation is not needed for deep well plate.

### Easy II Protocol

1. Aliquot 100 µl of homogenized lysis reagent (reagents A + F) to wells of a deep well plate.

**Note:** Gently shake the lysis reagent by hand first to resuspend the beads.

2. Add 100 µl of enriched sample.

**Note:** Shake the suspension to homogenize the culture and then allow any debris to settle before collecting the sample.

3. Mix the solution by pipetting up and down until homogenized.
4. Close the tubes or seal the deep well plate with pre-pierced sealing film.
5. Place tubes in the cell disruptor for  $3 \pm 1$  min.
6. Incubate tubes in the heat block at 95–100°C for 10–15 min. Incubate deep well plate in the thermoshaker under agitation at 1,300 rpm at 95–100°C for 15–20 min.
7. Vortex tubes at high speed, and then centrifuge at 10,000–12,000 x g for at least 2 min. Centrifugation is not needed for deep well plate.

This is the recommended stopping point for temporarily stopping the procedure.

The supernatant can be stored for up to 1 year at –20°C. Always allow it to thaw and homogenize, and then centrifuge at 10,000–12,000 x g for 5 min before reusing.

### Standard I Protocol

1. Collect 1 ml of decanted enriched sample into a tube.
2. Centrifuge at 10,000–12,000 x g for 5 min.
3. Discard all the supernatant.
4. Add 200 µl lysis reagent (reagent A) to pellet.
5. Pipet reagent in tube up and down enough times to resuspend pellet. The number of times will depend on the sample.
6. Vortex at high speed and place the tube in the heat block at 95–100°C for 10–15 min.
7. Vortex at high speed and centrifuge at 10,000–12,000 x g for 5 min.

This is the recommended stopping point for temporarily stopping the procedure.

The supernatant can be stored up to 1 year at –20°C. Always allow it to thaw and homogenize, and then centrifuge at 10,000–12,000 x g for 5 min before reusing.

### Standard II Protocol

1. Collect 1 ml of decanted enriched sample into a tube or in a well of a deep well plate. Seal the deep well plate with a plastic film.

## Section 7 Protocol

**Note:** Shake suspension to homogenize the culture and then allow any debris to settle before collecting the sample.

2. Centrifuge tubes at 10,000–12,000 x g for 5 min. Discard the supernatant. Centrifuge deep well plates at 2,250 x g for 20 min and discard supernatant manually or using the DW40.
3. Add 200 µl of homogenized lysis reagent (reagents A + F) to the pellet and resuspend pellet by pipetting the reagent up and down. Close tubes or seal deep well plate with pre-pierced sealing film.

**Note:** Gently shake lysis reagent by hand first to resuspend the beads. For primary production samples, use reagent A without reagent F.

4. Place tube in the cell disruptor for  $3 \pm 1$  min, or shake deep well plate in the thermoshaker at 1,300–1,600 rpm at 95–100°C for 15–20 min.
5. Incubate tubes in the appropriate heat block at 95–100°C for 10–15 min.
6. Vortex tubes at high speed and centrifuge at 10,000–12,000 x g for 5 min. Centrifuge deep well plate at 2,250 x g for 2 min.

This is the recommended stopping point for temporarily stopping the procedure.

The supernatant can be stored for up to 1 year at –20°C. Always allow it to thaw and homogenize, and then centrifuge at 10,000–12,000 x g for 5 min before reusing.

## D. Real-Time PCR

### Instrument and Software Setup

For instrument and software setup, follow instructions in the real-time PCR system user guide for iQ-Check Kits.

### PCR Mix Preparation

1. Prepare PCR mix containing the amplification solution (reagent C) and the fluorescent probes (reagent B). Volume of PCR mix needed depends on the number of samples and controls to be analyzed. At least one positive and one negative control must be included in each PCR run. Use the pipetting table in the Appendix – PCR Mix Calculation Guide to find the correct volumes to use for each reagent.
2. Use the PCR mix (reagents B + C) immediately after preparation. It is stable for 1 hr maximum at 2–8°C.
3. Pipet 45 µl of PCR mix into each well according to your plate setup.
4. To purify DNA, combine 50 µl of DNA extracted from each sample with 200 µl of iQ-Check Purification Reagent. Pipet up and down 5 times to homogenize. Alternatively, a 1:10 dilution of extracted DNA can be made using sterile water.
5. Add 5 µl of DNA extract, reagent D (negative control), or reagent E (positive control). Do not vortex the sample before pipetting. Hermetically seal the wells of the plate or strips. It is important to avoid

bubbles at the bottom of the wells by pipetting carefully. As an optional step, centrifuge the sealed PCR plate or tube strips (quick spin) to eliminate any bubbles.

6. Place the plate or strips in the thermal cycler. Be sure to place the plate with the A1 well at the upper left corner. Close the reaction module.

## Run PCR

To start the PCR run, follow instructions in the real-time PCR system user guide for iQ-Check Kits.

The Application Protocol File “Salmo Fast” has been certified by NF Validation and AOAC PTM.

## E. Data Analysis

Data can be analyzed directly at the end of the PCR run or at a later time by opening the stored data file. Follow instructions in the corresponding CFX Manager IDE Software User Manual for opening data files and setting the data analysis parameters.

### Interpreting Results

Once the data analysis parameters have been set, results are interpreted by analyzing the Cq values of each sample (the cycle at which the amplification curve crosses the threshold).

CFX Manager IDE Software allows complete automated analysis for Bio-Rad real-time PCR detection systems. A verification of the typical characteristics of the amplification curves should be performed prior to releasing results. Contact your Bio-Rad technical support team if additional support is required.

### Controls

Verify the positive and negative controls before interpreting sample results.

For the experiment to be valid, the controls must have the following results, as summarized in the table below. Otherwise, the PCR reaction must be repeated.

	<b><i>Salmonella</i> spp. Detection (FAM channel)</b>	<b>Internal Control Detection (HEX channel)</b>
Negative control	Cq = N/A*	$28 \leq Cq \leq 40$
Positive control	$26 \leq Cq \leq 36$	N/A

\* The software indicates a Cq value of N/A (not applicable) when the fluorescence of a sample does not rise significantly above the background noise and hence does not cross the threshold.

If results of negative and positive controls differ from those in the table above (invalid control), repeat the run and analysis described in D. Real-Time PCR and E. Data Analysis in Section 7 Protocol.

### Samples

A **positive** iQ-Check *Salmonella* II PCR test must show a typical amplification curve and a Cq value  $\geq 10$  for the FAM fluorophore.

## Section 8 Confirmation of Positive Results

- If the Cq value for both channels is below 10, verify that the raw data curve is a regular amplification curve (with a flat baseline followed by a rapid exponential increase of fluorescence and then a flattening out). If the curve seems correct, it may be considered a positive *Salmonella* spp. test

If there is no Cq value (Cq = N/A) for FAM, or if the curve is not a typical amplification curve, the internal control for that sample must be analyzed.

- This sample is considered a **negative** *Salmonella* spp. sample if there is no Cq value in the FAM channel and the internal control has a Cq value  $\geq 28$
- Should the internal control also not have a Cq value (Cq = N/A), this probably indicates an inhibition of the PCR reaction. The sample needs to be diluted (using 10  $\mu$ l of DNA extract, perform a 1:10 dilution in distilled sterile water, and then test 5  $\mu$ l of the dilution) and the PCR repeated
- Should the Cq value for the internal control be  $< 28$ , it is not possible to interpret the result. Verify that the threshold was correctly placed or that the curve as raw data is a regular amplification curve. If the curve does not have a characteristic shape, it will be necessary to repeat the PCR test.

Interpretation of test results is summarized in the following table:

<b><i>Salmonella</i> spp. Detection (FAM channel)</b>	<b>Internal Control Detection (HEX channel)</b>	<b>Interpretation</b>
Cq $\geq 10$	N/A	Positive
Cq = N/A	Cq $\geq 28$	Negative
Cq = N/A	Cq = N/A	Inhibition*

\* When both target and internal control detection give a Cq value = N/A, the DNA extract must be diluted 1:10 and tested again.

An invalid interpretation can be given when validation criteria are not met. Check the raw data and proceed as if the sample was inhibited.

## Section 8

### Confirmation of Positive Results

In the context of AOAC validation, a positive iQ-Check *Salmonella* II result should be considered presumptive positive and confirmed according to an appropriate reference method (for example, USDA MLG, FDA BAM, ISO, MFHPB, AOAC SMPR, etc.) is recommended. Alternatively, streak 10  $\mu$ l enriched BPW (with supplements when applicable) directly to RAPID' *Salmonella* chromogenic and incubate for  $24 \pm 2$  hr at  $37 \pm 1^\circ\text{C}$  (except for hemp flower). When interfering flora could be high, transfer 0.5 ml enriched BPW to 10 ml Rappaport Vassiliadis Soy (RVS) broth and incubate for  $24 \pm 3$  hr at  $41.5 \pm 1^\circ\text{C}$  then proceed with streaking 10  $\mu$ l enriched RVS to RAPID' *Salmonella* chromogenic media. Refer to the confirmation method in the RAPID' *Salmonella* chromogenic media user guide (document #10000126748).

In the context of the NF VALIDATION–certified method, all positive iQ-Check Kit results must be confirmed in the following ways:

#### All samples except primary production samples

1. Using classical tests described in the standardized CEN or ISO methods (including the purification) directly from the BPW enrichment broth after the full  $18 \pm 2$  hr enrichment at  $37^\circ\text{C}$ .

2. By streaking 10 µl enriched BPW (with supplements when applicable) directly to RAPID' *Salmonella* chromogenic media and incubating for 24 ± 2 hr at 37 ± 1°C. When interfering flora could be high, transfer 0.5 ml enriched BPW to 10 ml Rapaport Vassiliadis Soy (RVS) broth and incubate for 24 ± 3 hr at 41.5 ± 1°C then proceed with streaking 10 µl enriched RVS to RAPID' *Salmonella* chromogenic media.
3. Using any other method certified by NF VALIDATION based on a principle different from that used in the iQ-Check *Salmonella* II PCR test. The validated protocol of this second method must be followed entirely, including RVS selective enrichment, when needed. Confirmation is carried out from the BPW enrichment broth if this step is common to both methods.

It is possible to store the enriched BPW at 2–8°C for 72 hr maximum, following the last incubation at 37°C, with the exception of dairy and environmental surface samples.

### Primary production samples

Start the confirmation process from the primary enrichment, following the MSR/V method indicated in Appendix D of the ISO 6579/A1 standard and using RAPID' *Salmonella* medium for the isolation.

Disagreement between results of the iQ-Check *Salmonella* II Kit and the confirmation method described above may be caused by the presence of nonmotile *Salmonella*. In that case, follow the RAPID' *Salmonella* double enrichment protocol (refer to the product user guide for instructions for use). This protocol includes a selective enrichment step in RVS medium.

In the scope of the NF VALIDATION, primary enrichment bags can be stored for up to 24 hr at 2–8°C before proceeding to the confirmation.

In the case of discrepant results between the iQ-Check *Salmonella* II Kit and any of the confirmation options listed above, follow the necessary steps to ensure valid results.

## Section 9

### Confirmation of Single Colonies Using iQ-Check Kit

iQ-Check *Salmonella* spp. II. Kit may also be used for confirming single isolated *Salmonella* spp. colonies on agar plates. This is officially validated by AFNOR Certification on RAPID' *Salmonella* chromogenic media.

1. Pick an isolated colony from a selective or nonselective agar plate with a toothpick, sterile loop, or other adapted consumable (for example, a pipet tip).
2. Resuspend the colony in 100 µl tryptone salt or distilled sterile water in a microcentrifuge tube. Homogenize using a vortexer.
3. Use 5 µl of the suspension with 45 µl of PCR mix (see D. Real-Time PCR, in Section 7 Protocol) and follow the rest of the iQ-Check *Salmonella* II Protocol for the data and result interpretation.

## Section 10

# Test Performance and Validations

The iQ-Check *Salmonella* II Kit is specific for the detection of *Salmonella* genus.

BRD 07/10 – 04/05

ALTERNATIVE ANALYSIS  
METHODS FOR AGRIBUSINESS

### NF Validation

The iQ-Check *Salmonella* II Kit is certified by NF VALIDATION as an alternative to the reference method ISO 6579 -1 (2017) for the detection of *Salmonella* spp. in all products for human and animal consumption, as well as environmental samples. The validation followed the protocol of the ISO 16140: 2016 standard and includes the use of the CFX96 Touch, CFX96 Touch Deep Well and CFX Opus Deepwell Real-Time PCR Systems. The use of the iQ-Check Free DNA Removal Solution is validated for infant formula and ingredients, pet food and animal feed. The associated software is CFX Manager Software IDE (V2.2 and later). The “Salmo Fast” APF is validated for all samples. Certificate number: BRD 07/06 – 07/04. Refer to the certificate available on the AFNOR Certification website for validity information.

OMA 2017.06

### AOAC Validation

The iQ-Check *Salmonella* II Kit (Easy Protocol) is validated by the AOAC Research Institute under OMA 2017.06 for the detection of *Salmonella* spp. in eggs (25 g), raw ground chicken (25 and 375 g), raw beef (25 g), raw pork (25 g), fresh spinach (25 g and 375 g), cantaloupe (25 g), peanut butter (25 g), ready-to-eat deli ham (375 g), selected environmental surfaces, wet (25 g) and dry (25 and 375 g) pet food, milk chocolate (375 g), cheese (375 g), carcass rinses, chocolate liquor (375 g), white chocolate (375 g), NFDM (375 g) whey powder (375 g), raw ground beef (375 g), raw beef trim (375 g), cannabis flower (10 g), cannabis-infused chocolate (25 g), cannabis-infused gummies (25 g), cannabis-derived concentrates (5 g), hemp flower (25 g), flour (375 g) and plant-based meat products (375 g). A positive result with the iQ-Check Kit is considered presumptive, and it is recommended it be confirmed following the recommendation in Section 8 above. The “Salmo Fast” APF, the use of the iQ-Check Free DNA Removal Solution, and the use of CFX96 Touch Deep Well and CFX Opus Deepwell Real-Time PCR Systems are validated for all samples. The associated software is CFX Manager IDE Software (version 2.2 and later).

### Health Canada Validation

The iQ-Check *Salmonella* II (Easy Protocol) has been validated by Health Canada (MFLP-38) for detection of *Salmonella* spp. A positive result with the iQ-Check Kit is presumed to be positive and must be confirmed according to MFHPB-20 (see Section 11, reference 5).

## NordVal Validation

The iQ-Check *Salmonella* II Kit is validated by NordVal for Standard I and Easy I Protocols (in microplates) in food, feed, and environment, and for Easy II Protocols (in microplates) in raw meat. The “Salmo Fast” APF is validated for all samples. Certificate number: 038. Refer to the certificate available on the NordVal website for validity information.

## Section 11

### References

- ISO 7218. Microbiology of the food and animal feeding stuffs — General requirements and guidance for microbiological examination.
- ISO 6887-4. Microbiology of the food chain — Preparation of test samples, initial suspension and decimal dilutions for microbiological examination — Part 4: Specific rules for the preparation of miscellaneous products.
- ISO 6579-1. Microbiology of the food chain — Horizontal method for the detection, enumeration and serotyping of *Salmonella* — Part 1: Detection of *Salmonella* spp.
- ISO 16140-2. Microbiology of the food chain — Method validation — Part 2: Protocol for the validation of alternative (proprietary) methods against a reference method.
- Miras I et al. (1995). Nucleotide sequence of *iagA* and *iagB* genes involved in invasion of HeLa cells by *Salmonella* enteritica subsp. Enteritica ser. Typhi. *Res Microbiol* 46, 17–20.
- United States Department of Agriculture, Food Safety and Inspection Service (2019). Microbiology Laboratory Guidebook, Chapter 4.10: Isolation and Identification of *Salmonella* from Meat, Poultry, Pasteurized Egg, and Catfish Products and Carcass and Environmental Sponges.
- United States Food and Drug Administration, Center for Food Safety and Applied Nutrition (2018). Bacteriological Analytical Manual, Chapter 5: *Salmonella*.
- Health Canada Health Products and Food Branch, MFHPB-20 Isolation and Identification of *Salmonella* from Food and Environmental Samples, Compendium of Analytical Methods, Volume 2, March 2009.
- Health Canada Health Products and Food Branch, Part 4: Qualitative Methods, Procedure for the Development and Management of Food Microbiological Methods, Compendium of Analytical Methods, Volume 1.

## Section 12

### Revision History

Release date	Document number	Change
October 2020	10000131519 Ver A	<ul style="list-style-type: none"> <li>- Renewal and extension of AFNOR validation: FDRS use as an option, "Salmo Fast" APF</li> <li>- New document design and update of references and content</li> <li>- Document number change (previous version was 808463 Rev G)</li> </ul>
February 2023	10000131519 Ver B	<ul style="list-style-type: none"> <li>- Added new catalog number for consumables and instrument</li> <li>- Clarification of result interpretation protocol</li> <li>- AOAC extension: New matrices, FDRS, "Salmo Fast" APF</li> <li>- Replaced protocol section with summary table</li> <li>- Update of NordVal validation scope</li> <li>- Extension of validations for CFX Opus Deepwell Real-Time PCR System and upgrade to CFX Manager Software IDE version 3.1</li> <li>- AFNOR Validation extension on pet food and animal feed</li> </ul>

## Appendix — PCR Mix Calculation Guide

To find the correct volumes to use when preparing the PCR mix, add the total number of samples and controls to be analyzed and find the corresponding volumes of reagent B and reagent C in the table.

Total Number of Samples and Controls	Probes Reagent B, $\mu$ l	Amplification Mix Reagent C, $\mu$ l	Total Number of Samples and Controls	Probes Reagent B, $\mu$ l	Amplification Mix Reagent C, $\mu$ l	Total Number of Samples and Controls	Probes Reagent B, $\mu$ l	Amplification Mix Reagent C, $\mu$ l
1	5	40	33	178	1,400	65	351	2,800
2	11	86	34	184	1,500	66	356	2,900
3	16	130	35	189	1,500	67	362	2,900
4	22	173	36	194	1,600	68	367	2,900
5	27	216	37	200	1,600	69	373	3,000
6	32	259	38	205	1,600	70	378	3,000
7	38	302	39	211	1,700	71	383	3,100
8	43	346	40	216	1,700	72	389	3,100
9	49	389	41	221	1,800	73	394	3,200
10	54	432	42	227	1,800	74	400	3,200
11	59	475	43	232	1,900	75	405	3,200
12	65	518	44	238	1,900	76	410	3,300
13	70	562	45	243	1,900	77	416	3,300
14	76	605	46	248	2,000	78	421	3,400
15	81	648	47	254	2,000	79	427	3,400
16	86	691	48	259	2,100	80	432	3,500
17	92	734	49	265	2,100	81	437	3,500
18	97	778	50	270	2,200	82	443	3,500
19	103	821	51	275	2,200	83	448	3,600
20	108	864	52	281	2,200	84	454	3,600
21	113	907	53	286	2,300	85	459	3,700
22	119	950	54	292	2,300	86	464	3,700
23	124	994	55	297	2,400	87	470	3,800
24	130	1,000	56	302	2,400	88	475	3,800
25	135	1,100	57	308	2,500	89	481	3,800
26	140	1,100	58	313	2,500	90	486	3,900
27	146	1,200	59	319	2,500	91	491	3,900
28	151	1,200	60	324	2,600	92	497	4,000
29	157	1,300	61	329	2,600	93	502	4,000
30	162	1,300	62	335	2,700	94	508	4,100
31	167	1,300	63	340	2,700	95	513	4,100
32	173	1,400	64	346	2,800	96	518	4,100

Visit [bio-rad.com/iqcheck](http://bio-rad.com/iqcheck) for more information.

BIO-RAD is a trademark of Bio-Rad Laboratories, Inc.

IQ-CHECK is a trademark of Bio-Rad Europe GMBH in certain jurisdictions.

All trademarks used herein are the property of their respective owner.

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Miras I et al. (1995). Nucleotide sequence of *iagA* and *iagB* genes involved in invasion of HeLa cells by *Salmonella enteritica* subsp. Enteritica ser. Typhi. *Res Microbiol* 46, 17–20.

United States Department of Agriculture, Food Safety and Inspection Service (2019). *Microbiology Laboratory Guidebook*, Chapter 4.10: Isolation and Identification of *Salmonella* from Meat, Poultry, Pasteurized Egg, and Catfish Products and Carcass and Environmental Sponges.

United States Food and Drug Administration, Center for Food Safety and Applied Nutrition (2018). *Bacteriological Analytical Manual*, Chapter 5: *Salmonella*.

Health Canada Health Products and Food Branch, MFHPB-20 Isolation and Identification of *Salmonella* from Food and Environmental Samples, *Compendium of Analytical Methods*, Volume 2, March 2009.

Health Canada Health Products and Food Branch, Part 4: Qualitative Methods, Procedure for the Development and Management of Food Microbiological Methods, *Compendium of Analytical Methods*, Volume 1.

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## Sezione 11

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ISO 7218. Microbiology of the food and animal feeding stuffs — General requirements and guidance for microbiological examination.

ISO 6887-4. Microbiology of the food chain — Preparation of test samples, initial suspension and decimal dilutions for microbiological examination — Part 4: Specific rules for the preparation of miscellaneous products.

ISO 6579-1. Microbiology of the food chain — Horizontal method for the detection, enumeration and serotyping of *Salmonella* — Part 1: Detection of *Salmonella* spp.

ISO 16140-2. Microbiology of the food chain — Method validation — Part 2: Protocol for the validation of alternative (proprietary) methods against a reference method.

Miras I et al. (1995). Nucleotide sequence of *iagA* and *iagB* genes involved in invasion of HeLa cells by *Salmonella* enteritica subsp. Enteritica ser. Typhi. *Res Microbiol* 46, 17–20.

United States Department of Agriculture, Food Safety and Inspection Service (2019). Microbiology Laboratory Guidebook, Chapter 4.10: Isolation and Identification of *Salmonella* from Meat, Poultry, Pasteurized Egg, and Catfish Products and Carcass and Environmental Sponges.

United States Food and Drug Administration, Center for Food Safety and Applied Nutrition (2018). Bacteriological Analytical Manual, Chapter 5: *Salmonella*.

Health Canada Health Products and Food Branch, MFHPB-20 Isolation and Identification of *Salmonella* from Food and Environmental Samples, Compendium of Analytical Methods, Volume 2, March 2009.

Health Canada Health Products and Food Branch, Part 4: Qualitative Methods, Procedure for the Development and Management of Food Microbiological Methods, Compendium of Analytical Methods, Volume 1.

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