Italia

CERTIFICATO

Nr. 50 100 11497 - Rev. 003

Si attesta che / This is to certify that

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È CONFORME AI REQUISITI DELLA NORMA
HAS BEEN FOUND TO COMPLY WITH THE REQUIREMENTS OF

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**Design and development, production and sale of in-vitro diagnostic
medical devices: culture media for bacteriology, identification and
susceptibility testing systems, kits for plasma protein determination
(IAF 12, 29)**

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"THE VALIDITY OF THE PRESENT CERTIFICATE DEPENDS ON THE ANNUAL SURVEILLANCE EVERY 12 MONTHS AND ON THE COMPLETE REVIEW OF COMPANY'S MANAGEMENT SYSTEM AFTER THREE-YEARS"

A.F. GENITAL SYSTEM

System for detection, count and susceptibility test of pathogenic urogenital microorganisms.

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A.F. GENITAL SYSTEM

ENGLISH

System for detection, count and susceptibility test of pathogenic urogenital microorganisms

DESCRIPTION

A.F. GENITAL SYSTEM is a 24-well system containing desiccated biochemical and antibiotic substrates for detection, presumptive identification and susceptibility test of microorganisms from urogenital specimens. The system also provides a semi-quantitative assessment of the presence of urogenital mycoplasmas (*Mycoplasma hominis* and *Ureaplasma urealyticum*). The system is inoculated with the suspension of the clinical sample and incubated at 36 ± 1 °C for 18-24 hours. The tests used in the detection, count and presumptive identification of the microorganisms and for the susceptibility test of the urogenital mycoplasmas in the sample are interpreted by assessing the change in colour of the various wells and performing a microscope examination.

KIT CONTENTS

The kit contains:

20 A.F. GENITAL SYSTEM	20 vials of Physiological Solution (7.0 mL/vial)
Instruction sheet	1 TEST RESULTS FORM

ITEMS NECESSARY BUT NOT INCLUDED IN THE KIT

A.F. GENITAL SYSTEM <i>Reagent</i> (ref. 80258)	Mycoplasma Transport Broth (ref. 20158)
Sampling swabs	Cover slides and slides
Microscope	STREPTO B latex KIT (Ref. 96153)

CONFIGURATION

The configuration of the system is shown in Table no.1.

Table no.1

Well	COUNT AND IDENTIFICATION OF MYCOPLASMAS / UREAPLASMAS
1-GR+	Growth of mycoplasmas ($10^2 < \text{titre} < 10^4$ CFU/mL)
2-GR++	Growth of mycoplasmas ($10^4 < \text{titre} < 10^5$ CFU/mL)
3-GR+++	Growth of mycoplasmas ($\text{titre} > 10^5$ CFU/mL)
4-ADC	Arginine test
5-UR	Urea test
Well	DETECTION OF T. VAGINALIS AND CANDIDA Spp.
6-TR/YE	<i>Trichomonas vaginalis</i> / <i>Candida</i> spp.
Well	SUSCEPTIBILITY TEST OF THE MYCOPLASMAS / UREAPLASMAS (*)
7-TE	Tetracycline - 8 µg/mL
8-PEF	Pefloxacin - 16 µg/mL
9-OFX	Ofloxacin - 4 µg/mL
10-DO	Doxycycline - 8 µg/mL
11-E	Erythromycin - 16 µg/mL
12-CLA	Clarithromycin - 16 µg/mL
13-MN	Minocycline - 8 µg/mL
14-JOS	Josamycin - 8 µg/mL
15-CD	Clindamycin - 8 µg/mL
Well	DETECTION AND PRESUMPTIVE IDENTIFICATION OF OTHER MICROORGANISMS
16-ESC	<i>Escherichia coli</i>
17-PRO	<i>Proteus</i> spp. / <i>Providencia</i> spp.
18-PSE	<i>Pseudomonas</i> spp.
19-GAR	<i>Gardnerella vaginalis</i>
20-STF	<i>Staphylococcus aureus</i>
21-STR	<i>Enterococcus faecalis</i>
22-NES	<i>Neisseria gonorrhoeae</i>
23-STG	<i>Streptococcus agalactiae</i> (Group B)
24-CAN	<i>Candida</i> spp.

(*):The concentration of each antibiotic is according to the updated NCCLS – January 2004; Vol 24 N°1 ⁽¹⁾

PRINCIPLE OF THE METHOD

A.F. GENITAL SYSTEM allows the detection, semi-quantitative count, presumptive identification and susceptibility test of *Mycoplasma hominis* and *Ureaplasma urealyticum*, the detection and presumptive identification of the microorganisms most frequently isolated from vaginal and urethral swabs and seminal fluid, such as: *Trichomonas vaginalis*, *Escherichia coli*, *Proteus* spp./*Providencia* spp., *Pseudomonas* spp., *Gardnerella vaginalis*, *Staphylococcus aureus*, *Enterococcus faecalis*, *Neisseria gonorrhoeae*, *Streptococcus agalactiae* (Group B) and *Candida* spp.

- The semi-quantitative count of urogenital mycoplasmas (*Mycoplasma hominis* and *Ureaplasma urealyticum*) is provided by the yellow to red color change in wells **1-GR+**, **2-GR++**, **3-GR+++**.
- The presence of *Mycoplasma hominis* is shown by the yellow to red color change of the Arginine Test well **4-ADC**.
- The presence of *Ureaplasma urealyticum* is shown by the yellow to red color change in of the Urea Test well **5-UR**.
- The presence of *Trichomonas vaginalis* and *Candida* spp. is determined by microscope (40x) examination of a drop of culture liquid taken from well **6-TR/YE** assessing the presence of mobile ciliated trophozoites for identification of *Trichomonas vaginalis* and/or the presence of chlamydospores and hyphae for identification of *Candida* spp.
- The susceptibility testing of the urogenital mycoplasmas is assessed by observing the color change in wells **7-TE**, **8-PEF**, **9-OFX**, **10-DO**, **11-E**, **12-CLA**, **13-MN**, **14-JOS**, **15-CD**, as indicated in table no.3.
- The presence of *Escherichia coli* is shown by the grey-red to blue color change in well **16-ESC**.

Confirm the presence of *Escherichia coli* by adding two drops of KOVACS' Reagent into the well **17-PRO**; wait for the formation of a pink-red ring.

- The presence of *Proteus* spp./*Providencia* spp. is shown by the yellow to brown-black color change of well **17-PRO**.
- The presence of *Pseudomonas* spp. is shown by the yellow to turbid green color change in the well **18-PSE**. Confirm by the oxidase test (A.F. GENITAL SYSTEM REAGENT ref. 80258).
- The presence of *Gardnerella vaginalis* is shown by the red to yellow-orange color change in the well **19-GAR**. Confirm by isolation on selective media (GARDNERELLA VAG. AGAR ref. 11054) and biochemical tests (INTEGRAL SYSTEM GARDNERELLA ref. 71724).
- The presence of *Staphylococcus aureus* is shown by the appearance of a black ring on the bottom of well **20-STF**. Confirm by the coagulase test (COAGULASE TEST Ref. 88030).
- The presence of *Enterococcus faecalis* is shown by the yellow to black color change in the well **21-STR**.
- The presence of *Neisseria gonorrhoeae* in well **22-NES** is shown by a positive oxidase test (A.F. GENITAL SYSTEM REAGENT Ref. 80258). Confirm the positivity test through isolation on selective media, (Thayer Martin Agar ref. 11040), biochemical tests and serological tests.
- The presence of *Streptococcus agalactiae* (Group B) is shown by the yellow to green color change in the well **23-STG**. Confirm the presence of ***Streptococcus agalactiae*** with the following methods:
 - a) – Agglutination test **STREPTO B latex KIT** (ref. 96153)
 - b) – Cultural isolation with traditional methods
- The presence of *Candida* spp. is shown by the green to yellow color change in the well **24-CAN**. Watch for the presence of chlamydospores and hyphae at the microscope.

COMPOSITION

Table no.2

Well	Content
1-GR+	Culture medium for growth of mycoplasmas
2-GR++	Culture medium for growth of mycoplasmas
3-GR+++	Culture medium for growth of mycoplasmas
4-ADC	Culture medium with substrate for decarboxylation of arginine
5-UR	Culture medium with substrate for hydrolysis of urea
6-TR/YE	Culture medium with substrate for isolation of <i>Trichomonas vaginalis</i> and <i>Candida spp.</i>
7-TE	Culture medium containing Tetracycline - 8 µg/mL
8-PEF	Culture medium containing Pefloxacin - 16 µg/mL
9-OFX	Culture medium containing Ofloxacin - 4 µg/mL
10-DO	Culture medium containing Doxycycline - 8 µg/mL
11-E	Culture medium containing Erythromycin - 16 µg/mL
12-CLA	Culture medium containing Claritromycin - 16 µg/mL
13-MN	Culture medium containing Minocycline - 8 µg/mL
14-JOS	Culture medium containing Josamycin - 8 µg/mL
15-CD	Culture medium containing Clindamycin - 8 µg/mL
16-ESC	Culture medium with substrate for isolation of <i>Escherichia coli</i>
17-PRO	Culture medium with substrate for isolation of <i>Proteus spp./Providencia spp.</i>
18-PSE	Culture medium with substrate for isolation of <i>Pseudomonas spp.</i>
19-GAR	Culture medium with substrate for isolation of <i>Gardnerella vaginalis</i>
20-STF	Culture medium with substrate for isolation of <i>Staphylococcus aureus</i>
21-STR	Culture medium with substrate for isolation of <i>Enterococcus faecalis</i>
22-NES	Culture medium with substrate for isolation of <i>Neisseria gonorrhoeae</i>
23-STG	Culture medium with chromogenic substrate for isolation of <i>Streptococcus agalactiae</i>
24-CAN	Culture medium with substrate for isolation of <i>Candida spp.</i>

Physiological Solution (g/L): Sodium chloride **9 g**; Distilled water **1000 mL**; pH **6.8 ± 0.2**

COLLECTION AND STORAGE OF THE SAMPLES

Obtain the vaginal or urethral secretion with synthetic fibre swabs.

Obtain the seminal fluid according to the technique envisaged for the microbiological culture examination.

The samples must be sent for inoculation in the A.F. GENITAL SYSTEM immediately after they are taken.

Do not store in the refrigerator for even limited periods of time, since low temperatures may damage the vitality of some particularly sensitive microorganisms (such as *Trichomonas vaginalis*, *Neisseria gonorrhoeae* etc.) and alter the final result.

TEST PROCEDURE**PREPARATION OF THE CLINICAL SAMPLE**

Open one of the phials of Physiological Solution contained in the kit.

1. VAGINAL SWAB – URETHRAL SWAB

Immerse the swab (after obtaining the clinical material) in the phial of physiological solution and wait 5 minutes. Carefully squeeze the swab against the phial wall so that the clinical material is dispersed homogeneously in the Physiological Solution. Keep the swab, dipped in a nutritive broth, until the conclusion of the test.

2. SEMINAL FLUID

Add 0.2 mL of sample to the vial of Physiological Solution, shake and wait 5 minutes before inoculation.

3. CLINICAL MATERIAL FROM Mycoplasma Transport Broth (ref. 20158)

Add 1 mL Mycoplasma Transport Broth, containing the sample under examination, to the vial of Physiological Solution, shake gently and wait 5 minutes before inoculation.

Conserve the specimen in the vial of Mycoplasma Transport Broth until the conclusion of the test.

INOCULATION OF THE SYSTEM

1. Take a system from its wrapper and bring it to room temperature.
2. Write down the name of the patient, the date of the start of the examination, and the type of clinical material.
3. Transfer 0.2 mL of suspension of the clinical sample into each well of the system.
4. Cover with 1 drop of Vaseline Oil for microbiological use all the wells including 19-GAR, except for the wells 6, 16, 17, 18, 20, 21, 22 and 23.
5. Cover the system with the lid provided and incubate at 36 ± 1 °C for 18-24 hours. If the presence of *Mycoplasma hominis* is suspected, incubate at 36 ± 1 °C for further 24 hours.

INTERPRETATION OF THE RESULTS

At the end of the incubation:

- Watch for the color change in the wells **1-GR+** to **5-UR** and interpret the results using table no.3. (COUNT AND IDENTIFICATION OF THE MYCOPLASMAS / UREAPLASMAS).
- Take a drop of liquid from well **6-TR/YE**, deposit it on a glass slide, place a cover slip on top and examine at the microscope (40x) for the presence of *Trichomonas vaginalis* and *Candida* spp. Interpret as indicated in Table no.3. (SEARCH FOR T. VAGINALIS AND CANDIDA Spp.).
- Watch for the color change in the wells **7-TE** to **15-CD** and interpret the results using table no.3. (SUSCEPTIBILITY TESTING OF THE MYCOPLASMAS / UREAPLASMAS).
- Watch for the color change in the well **16-ESC** and interpret the results using table no.3.
Add two drops of KOVACS Reagent to the well **17-PRO** for confirming the presence of *Escherichia coli*; wait for the formation of a pink-red ring.
- Take a drop of liquid from well **22-NES**, deposit it on an Oxidase Test Stick and observe the immediate development (in about 10 seconds) of a blue color (positive oxidase test). Confirm the presence of *Neisseria gonorrhoeae* through isolation on selective media; drag the vaginal or urethral swab, previously kept in a nutritive broth (TEST PROCEDURE, point 1), directly on the surface of a selective plate (Thayer Martin Agar ref. 11040). Incubate the plate at 36 ± 1 °C in microaerophilia, and perform biochemical and serological tests to the grown colonies.
- Watch for the color change in the well **23-STG** and interpret the results using table no.3.

Confirm the presence of *Streptococcus agalactiae* with the following methods:

a) – Agglutination test **STREPTO B latex KIT** (Ref. 96153)

b) – Cultural isolation with traditional methods

- Watch for the color change in the wells 16 to 24 and interpret the results using table no.3 (PRESUMPTIVE IDENTIFICATION OF MICROORGANISMS).
- Note the results on the *TEST RESULTS FORM* (copy as many forms as necessary).

Table no.3

Well	COUNT AND IDENTIFICATION OF MYCOPLASMAS / UREAPLASMAS	Well colour	
		Positive reaction	Negative reaction
1-GR+	Growth of mycoplasmas ($10^2 < \text{titer} < 10^4$ CFU/mL) ⁽¹⁾	red	yellow
2-GR++	Growth of mycoplasmas ($10^4 < \text{titer} < 10^5$ CFU/mL) ⁽²⁾	red	yellow
3-GR+++	Growth of mycoplasmas (titer $> 10^5$ CFU/mL) ⁽³⁾	red	yellow
4-ADC	Arginine Test: identifies <i>Mycoplasma hominis</i>	red	yellow
5-UR	Urea Test: identifies <i>Ureaplasma urealyticum</i>	red	yellow

CFU: Colony Forming Units

Well	T. VAGINALIS AND CANDIDA spp.	Microscope examination (40x)
6-TR/YE	<i>Trichomonas vaginalis</i> / <i>Candida</i> spp.	<i>T.vaginalis</i> : mobile ciliated protozoites <i>Candida</i> spp.: chlamydo spores & hyphae

Well	SUSCEPTIBILITY TESTING OF THE MYCOPLASMAS / UREAPLASMAS	Well colour		
		S	I	R
7-TE	Tetracycline - 8 µg/mL	yellow	orange	red
8-PEF	Pefloxacin - 16 µg/mL	yellow	orange	red
9-OFX	Ofloxacin - 4 µg/mL	yellow	orange	red
10-DO	Doxycycline - 8 µg/mL	yellow	orange	red
11-E	Erythromycin - 16 µg/mL	yellow	orange	red
12-CLA	Clarithromycin - 16 µg/mL	yellow	orange	red
13-MN	Minocycline - 8 µg/mL	yellow	orange	red
14-JOS	Josamycin - 8 µg/mL	yellow	orange	red
15-CD	Clindamycin - 8 µg/mL	yellow	orange	red

S = Sensitive

I = Intermediate sensitivity

R = Resistant

Well	PRESUMPTIVE IDENTIFICATION OF OTHER MICROORGANISMS	Well colour	
		Positive reaction	Negative reaction
16-ESC	<i>Escherichia coli</i>	blue	grey-red
17-PRO	<i>Proteus</i> spp. / <i>Providencia</i> spp.	brown-black	yellow
18-PSE	<i>Pseudomonas</i> spp.	turbid green	yellow-blue
19-GAR	<i>Gardnerella vaginalis</i>	yellow-orange	red
20-STF	<i>Staphylococcus aureus</i>	black ring	yellow
21-STR	<i>Enterococcus faecalis</i>	black	yellow
22-NES	<i>Neisseria gonorrhoeae</i>	blue ⁽⁴⁾	colourless ⁽⁴⁾
23-STG	<i>Streptococcus agalactiae</i> (Group B)	green	yellow
24-CAN	<i>Candida</i> spp.	turbid yellow	green

LEGEND

(1) = equivalent to 5-20 colonies grown on Mycoplasma Agar A7

(2) = equivalent to 20-50 colonies grown on Mycoplasma Agar A7

(3) = over 50 colonies grown on Mycoplasma Agar A7

(4) = after oxidase Test

QUALITY CONTROL

Every batch of **A.F. GENITAL SYSTEM** is subjected to quality control using the following reference microorganisms:

<i>Mycoplasma hominis</i>	ATCC 23114	<i>Gardnerella vaginalis</i>	ATCC 14018
<i>Ureaplasma urealiticum</i>	ATCC 27618	<i>Staphylococcus aureus</i>	ATCC 25923
<i>Candida albicans</i>	ATCC 10231	<i>Streptococcus faecalis</i>	ATCC 19433
<i>Escherichia coli</i>	ATCC 25922	<i>Neisseria gonorrhoeae</i>	ATCC 19424
<i>Proteus mirabilis</i>	ATCC 25933	<i>Streptococcus agalactiae</i>	ATCC 13813
<i>Pseudomonas aeruginosa</i>	ATCC 27853	<i>Trichomonas vaginalis</i>	ATCC 30001

FACTORS THAT MAY INVALIDATE THE RESULTS

Poor standardisation of the inoculum; clinical material unsuitable; use of expired systems and/or reagents; non compliance with temperatures and times of incubation.

LIMITS AND WARNINGS

The presence of *Neisseria gonorrhoeae*, *Streptococcus agalactiae* and *Gardnerella vaginalis* must be confirmed using appropriate culture media, biochemical tests and serological tests. For definitive identification of the micro-organisms, it is necessary to make use of supplementary confirmation tests.

PERFORMANCE

The results obtained with **A.F. GENITAL SYSTEM** agree with those obtained using traditional culture methods. Application of the Wilcoxon non parametric statistical test shows that for each clinical material, the groups of standard values and those obtained with **A.F. GENITAL SYSTEM** do not differ significantly: $P \leq 0.05$ ⁽²⁾.

PRECAUTIONS

The product, **A.F. GENITAL SYSTEM**, is not classified as hazardous under current legislation, however see the safety data sheet for a correct use. **A.F. GENITAL SYSTEM** is a disposable device to be used only for diagnostic use *in vitro*. It must be used in the laboratory by properly trained personnel, using approved aseptic and safety methods for handling pathogenic agents.

STORAGE

Store at 2-8 °C in the original packaging. Keep away from sources of heat and avoid excessive changes in temperature. In such conditions, the product will remain valid until the expiry date indicated on the label. Do not use beyond that date. Eliminate without using if there are signs of deterioration.

DISPOSAL OF USED MATERIAL

After use, **A.F. GENITAL SYSTEM** and material that has come into contact with the sample must be decontaminated and disposed of in accordance with the techniques used in the laboratory for decontamination and disposal of potentially infected material.

PRESENTATION

Product	Ref.	Packaging
A.F. GENITAL SYSTEM	74156	20 tests

TABLE OF SYMBOLS

IVD <i>In Vitro</i> Diagnostic Medical Device	 Do not reuse	 Manufacturer	 Contains sufficient for <n> tests	 Temperature limitation
REF Catalogue number	 Fragile, handle with care	 Use by	 Caution, consult accompanying documents	LOT Batch code

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A.F. GENITAL SYSTEM


SYSTEM FOR DETECTION, COUNT AND SUSCEPTIBILITY TESTING OF UROGENITAL MYCOPLASMAS AND PATHOGENIC MICROORGANISMS

DIRECT INOCULATION FROM SAMPLE SUSPENSION, NO ENRICHMENT BROTHS

RESULTS AFTER 24-48 HOURS INCUBATION AT 36±1 °C

SIMPLE LABORATORY EQUIPMENT

UROGENITAL MYCOPLASMAS/UREAPLASMAS COUNT

<p><10⁴ CFU/mL</p> <p>Low growth of mycoplasmas is shown by colour change of well: 1-GR+ from yellow to red</p>	<p><10⁵ CFU/mL</p> <p>Moderate growth of mycoplasmas is shown by colour change of wells: 1-GR+, 2-GR++ from yellow to red</p>	<p>>10⁵ CFU/mL</p> <p>High growth of mycoplasmas is shown by colour change of wells: 1-GR+, 2-GR++, 3-GR+++ from yellow to red</p>
---	---	--

UROGENITAL MYCOPLASMAS/UREAPLASMAS DETECTION

<p>Mycoplasma hominis</p> <p>The presence of Mycoplasma hominis is shown by colour change of well: 4-ADC from yellow to red</p>	<p>Ureaplasma urealyticum</p> <p>The presence of Ureaplasma urealyticum is shown by colour change of well: 5-UR from yellow to red</p>	<p>Mycoplasma hominis and Ureaplasma urealyticum</p> <p>The presence of Mycoplasma hominis and Ureaplasma urealyticum is shown by colour change of wells: 4-ADC and 5-UR from yellow to red</p>
--	---	---

UROGENITAL MYCOPLASMAS/UREAPLASMAS SUSCEPTIBILITY TESTING

<p>The Sensitivity to antibiotics is shown by the colour of wells: yellow</p>	<p>The Intermediate sensitivity to antibiotic is show by colour change of wells: from yellow to orange</p>	<p>The Resistance to antibiotic is show by colour change of wells: from yellow to red</p>
---	--	---

PATHOGENIC UROGENITAL MICROORGANISMS

<p>Trichomonas vaginalis / Yeasts</p> <p>Observe the presence of mobile ciliated trophozoites, chlamydo spores and hyphae under microscopic observation (40x).</p> <p>Trichomonas vaginalis trophozoite.</p> <p>Candida albicans chlamydo spores and hyphae</p>		
<p>Escherichia coli</p> <p>The presence of Escherichia coli is shown by colour change of well: 16 - ESC from red to blue</p>	<p>Proteus spp.</p> <p>The presence of Proteus spp. is shown by colour change of well: 17 - PRO from yellow to brown</p>	<p>Pseudomonas spp</p> <p>The presence of Pseudomonas spp. is shown by colour change of well: 18 - PSE from yellow to green</p>
<p>Gardnerella vaginalis</p> <p>The presence of Gardnerella vaginalis is shown by colour change of well: 19 - GAR from red to yellow</p>	<p>Staphylococcus aureus</p> <p>The presence of Staphylococcus aureus is shown by colour change of well: 20 - STF from yellow to black</p>	<p>Enterococcus faecalis</p> <p>The presence of Enterococcus faecalis is shown by colour change of well: 21 - STR from yellow to black</p>
<p>Neisseria gonorrhoeae</p> <p>The presence of Neisseria gonorrhoeae is shown by colour change of Oxidase test stick: from white to blue</p>	<p>Streptococcus agalactiae</p> <p>The presence of Streptococcus agalactiae is shown by colour change of well after Hippurate test: 23 - STG from yellow to violet</p>	<p>Candida spp</p> <p>The presence of Candida spp is shown by colour change of well: 24 - CAN from green to yellow</p>

Project of research n. S239P financed by Ministry of Education, University and Research (MIUR).

LIOFILCHEM Bacteriology Products

Via Scozia - Zona Ind.le - 64026 ROSETO D. A. (TE) ITALY - Phone +39 085 8930745 - Fax +39 085 8930330
Internet: www.liofilchem.net - E-mail: liofilchem@liofilchem.net


A.F. GENITAL SYSTEM

TEST PROCEDURE

Code 74156

20 test


1
Immerse the swab (after obtaining the clinical material) in the vial of physiological solution and wait 5 minutes. Carefully squeeze the swab against the vial wall.

2
Transfer 0.2 mL of clinical sample suspension into each well of the system.


3
Cover all the wells, except 6, 16, 17, 18, 19, 20, 21, 22, 23, with 3 drops of vaseline oil.

4
Cover the system with the lid provided and incubate in thermostat at $36 \pm 1^\circ\text{C}$ for 18-24 hours.


5
Take a drop of liquid from well **6-TR/YE**, deposit it on a glass slide, place a cover slip on top and examine under the microscope (40x) for the presence of *Trichomonas v.* and *Candida spp.*

6
Take a drop of liquid from well **22-NES**, deposit in on an Oxidase Test Stick and observe the immediate development (10 seconds roughly) of a blue colour (positive oxidase test)

7
Add 2 drop of Ninhydrin 7% for the hippurate test to well **23-STG** and observe the immediate development (about 10 seconds) of a blue color (positive test).

8
Note the results on the TEST RESULTS FORM and interpret.

EXAMPLE OF RESULT

Sample of vaginal swab with presence of ***Mycoplasma hominis*** (4-ADC), ***Ureaplasma urealyticum*** (5-UR) with growth +++ ($>10^5$ CFU/mL), ***E. coli*** (16-ESC) and ***Gardnerella vag.*** (19-GAR).

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EC DECLARATION OF CONFORMITY

The company Liofilchem® S.r.l., registered office in Via Scozia, 64026 Roseto degli Abruzzi (TE) Italy, as a manufacturer of the *in vitro* medical-diagnostic device listed in the attached table, Revision 31.0 of 08.01.2016

hereby certifies under its own responsibility

1. that the above mentioned device complies with all the applicable provisions of Directive 98/79/EC (Annex III) and its relevant transposition into national law;
2. the above mentioned is not included in Annex II, List A and B of Directive 98/79/EC;
3. that the technical documentation referred to at Annex III of the Directive 98/79/EC is available for the national authorities in its facility and that this documentation shall be kept for 5 years after the last product has been manufactured;
4. that the manufacturing process follows suitable principles of quality assurance;
5. that, has implemented and keep up to date, a post-production surveillance system for monitoring the products;
6. that the device in question, was introduced into the market provided with CE mark.

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Rev. 31.0 del 08.01.2016

621021	HEART INFUSION BROTH
621022	MIDDLEBROOK 7H10 AGAR BASE
621210	WURTZ LACTOSE AGAR
621265	ISOSENSITEST AGAR
621367	BILE BACTERIOLOGICAL
621401	IRON SULPHITE AGAR
621402	CARY BLAIR TRANSPORT MEDIUM
621601	GLUCOSE
621618	CHROMATIC™ MH
621619	CHROMATIC™ CRE AGAR BASE
621701	PEPTONE BACTERIOLOGICAL
622202	STREPTOCOCCUS SELECTIVE AGAR
630026	LOWENSTEIN JENSEN MEDIUM w GLYCEROL 1 litre
71618	ENTEROSYSTEM 18R 20 Tests
71630	STAF SYSTEM 18 R 20 Tests
71670	COPRO SYSTEM 40 Tests
71675	COPRO SYSTEM Plus 20 Tests
71678	PATHOGENIC SYSTEM DOUBLE 40 Tests
71679	PATHOGENIC SYSTEM 20 Tests
71681	PATHOGENIC SYSTEM AST
71714	INTEGRAL SYSTEM ENTEROBATTERI 20 Tests
71718	INTEGRAL SYSTEM STAFILOCOCCI 20 Tests
71720	INTEGRAL SYSTEM STREPTOCOCCI 20 Tests
71724	INTEGRAL SYSTEM GARDNERELLA 20 TESTS
71822	INTEGRAL SYSTEM YEASTS Plus 20 Tests
72560	STREPTO SYSTEM 12 R 40 Tests
72592	MYCOPLASMA SYSTEM Plus 20 Tests
74156	A.F. GENITAL SYSTEM 20 Tests
74160	URIN SYSTEM Plus 20 Tests
74161	URIN SYSTEM Chrom 20 Tests
76010	Sensi Test gram-negative 20 Tests
76020	Sensi Test gram-positive 20 Tests
76031	SensiQuattro Gram-negative 20 Tests
76032	SensiQuattro Gram-positive 20 Tests
76033	SensiQuattro Candida EU 20 Tests
78618	ENTERO PLURI TEST 10 Tests
78619	ENTERO PLURI TEST 25 Tests
78620	OXI/FERM PLURI TEST 10 Tests
78621	OXI/FERM PLURI TEST 25 Tests
79010	Sensi Test gram-negative 4 Tests
79020	Sensi Test gram-positive 4 Tests
79031	SensiQuattro Gram-negative 4 Tests
79032	SensiQuattro Gram-positive 4 Tests
79033	SensiQuattro Candida EU 4 Tests
79156	A.F. GENITAL SYSTEM 4 Tests
79160	URIN SYSTEM Plus 4 Tests
79161	URIN SYSTEM Chrom 4 Tests
79560	STREPTO SYSTEM 12 R 8 Tests
79592	MYCOPLASMA SYSTEM Plus 4 Tests
79618	ENTEROSYSTEM 18R 4 Tests
79630	STAF SYSTEM 18 R 4 Tests
79670	COPRO SYSTEM 8 Tests
79675	COPRO SYSTEM Plus 4 Tests
79678	PATHOGENIC SYSTEM DOUBLE 8 Tests
79679	PATHOGENIC SYSTEM 4 Tests
79681	PATHOGENIC SYSTEM AST
79714	INTEGRAL SYSTEM ENTEROBATTERI 4 Tests
79718	INTEGRAL SYSTEM STAFILOCOCCI 4 Tests
79720	INTEGRAL SYSTEM STREPTOCOCCI 4 Tests
79724	INTEGRAL SYSTEM GARDNERELLA 4 Tests

79822	INTEGRAL SYSTEM YEASTS Plus 4 Tests
80009	IODINE MKTT SOLUTION 10 x 10 ml
80010	XLT 4 supplement 2 x 50 ml
80021	GLYCEROL supplement 4 x 50 ml
80022	POTASSIUM TELLURITE 1% suppl. 5 x 10 ml
80031	TWEEN 80 supplement 2 x 50 ml
80040	CHROMATIC™ SALMONELLA Supplement 2x50 ml
80047	MULLER KAUFFMANN 3X50 ml (Iodio/B.G.O.1%)
80053	VITAMIN K 1% supplement 5 x 5 ml
80056	LEGIONELLA growth supplement 10 vials
80057	H2O2 REAGENT 1 x 10 ml
80060	DECONTAM-KIT
80110	UREA 40% 6X100 ml
80219	EGG YOLK emulsion 4 x 50 ml
80252	ENTEROSYSTEM 18R REAGENT 100/200 Tests
80253	COPRO SYSTEM REAGENTS (antisera)
80257	LISTERIA SYSTEM 18R -REAG 100/200 Tests
80258	AF GENITAL SYSTEM REAGENT
80260	IDENTIF. SYSTEM-REAGENT 100/200 Tests
80271	KOVAC'S REAGENT 4x25 ml
80272	FERRIC CHLORIDE 10% 2x 25 ml
80273	NINHYDRIN 7% 10 ml
80275	MIF COLOR KIT 50 Tests
80276	ZIEHL-NEESEN 3 x 250 ml
80277	METHYLENE BLUE Solution 250 ml
80279	VASELINE OIL 4 x 50 ml
80280	V.P. TEST-Reagent 10x10ml
80281	V.P. TEST EP 10 x 10 mL
80282	Kit May-Grünwald Giemsa
80290	SAFRANIN SOLUTION 1000 ml
80291	POTASSIUM TELLURITE 3.5% suppl.5x10 ml
80292	UREA 40 % supplement 10 x 5 ml
80293	GRAM COLOR KIT 4 x 250 ml
80294	KIT COLOR ALBERT 2 x 250 ml
80295	DECOLOURIZING SOLUTION 1000 ml
80296	LUGOL PVP SOLUTION 1000 ML
80297	SAFRANIN SOLUTION 500 ml
80298	LUGOL PVP SOLUTION 250 ml
80299	CRYSTAL VIOLET SOLUTION 1000 ml
80300	TTC 1% supplement 5 x 10 ml
80350	ANTIBIOTIC TEST
80351	RAPID ANTIBIOTIC TEST 50 Tests
80380	KINYOUN COLOR KIT 2 x 250 ml
80390	FIXUR 1
80409	IODINE SOLUTION 10 x 10 ml
80410	XLT 4 SUPPLEMENT 4 x 50 ml
80422	POTASSIUM TELLURITE 1% Supplement 10 x 10 ml
80430	TTC 1% supplement 10 x 10 ml
80431	TWEEN 80 Supplement 4 x 50 ml
80453	VITAMIN K 1% SUPPLEMENT 10 x 5 ml
80491	POTASSIUM TELLURITE 3,5% Supplement 10 x 10
81001	AMPICILLIN supplement 10 vials
81002	LEGIONELLA (BMPA) supplement 10 vials
81003	BRUCELLA supplement 10 vials
81004	CAMPYLOBACTER Preston supplem 10 vials
81006	CN (Pseudomonas) supplement 10 vials
81007	CLOSTRIDIUM difficile supleme 10 vials
81008	LEGIONELLA (GVPC) supplement 10 vials
81009	IODINE solution 5 x 10 ml
81011	CLOSTRIDIUM perfringens (T.S.C.) sup.10 v.

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91200	DISC DISPENSER 8 CARTRIDGES
91203	DISC DISPENSER 6 CARTRIDGES
92000	AMOX*/SULB 2/1 AXS 0.016-256* 30 MIC Tests
920000	AMOX*/SULB 2/1 AXS 0.016-256* 100 MIC Tests
92001	RIFAMPICIN RD 0.002-32 30 MIC Tests
920010	RIFAMPICIN RD 0.002-32 100 MIC Tests
920011	RIFAMPICIN RD 0.002-32 10 MIC Tests
92002	FUSIDIC ACID FU 0.016-256 30 MIC Tests
920020	FUSIDIC ACID FU 0.016-256 100 MIC Tests
920021	FUSIDIC ACID FU 0.016-256 10 MIC Tests
92003	AMPICILLIN AMP 0.016-256 30 MIC Tests
920030	AMPICILLIN AMP 0.016-256 100 MIC Tests
920031	AMPICILLIN AMP 0.016-256 10 MIC Tests
92004	POLYMYXIN B PB 0.064-1024 30 MIC Tests
920040	POLYMYXIN B PB 0.064-1024 100 MIC Tests
920041	POLYMYXIN B PB 0.064-1024 10 MIC Tests
92005	CEFPODOXIME PX 0.016-256 30 MIC Tests
920050	CEFPODOXIME PX 0.016-256 100 MIC Tests
920051	CEFPODOXIME PX 0.016-256 10 MIC Tests
92006	CEFOTAXIME CTX 0.016-256 30 MIC Tests
920060	CEFOTAXIME CTX 0.016-256 100 MIC Tests
920061	CEFOTAXIME CTX 0.016-256 10 MIC Tests
92007	CEFOTAXIME CTX 0.002-32 30 MIC Tests
920070	CEFOTAXIME CTX 0.002-32 100 MIC Tests
920071	CEFOTAXIME CTX 0.002-32 10 MIC Tests
92008	CEFPIROME CR 0.016-256 30 MIC Tests
920080	CEFPIROME CR 0.016-256 100 MIC Tests
920081	CEFPIROME CR 0.016-256 10 MIC Tests
92009	GENTAMICIN CN 0.016-256 30 MIC Tests
920090	GENTAMICIN CN 0.016-256 100 MIC Tests
920091	GENTAMICIN CN 0.016-256 10 MIC Tests
92010	GENTAMICIN CN 0.064-1024 30 MIC Tests
920100	GENTAMICIN CN 0.064-1024 100 MIC Tests
920101	GENTAMICIN CN 0.064-1024 10 MIC Tests
92011	GATIFLOXACIN GAT 0.002-32 30 MIC Tests
920110	GATIFLOXACIN GAT 0.002-32 100 MIC Tests
920111	GATIFLOXACIN GAT 0.002-32 10 MIC Tests
92012	TEICOPLANIN TEC 0.016-256 30 MIC Tests
920120	TEICOPLANIN TEC 0.016-256 100 MIC Tests
920121	TEICOPLANIN TEC 0.016-256 10 MIC Tests
92013	ENROFLOXACIN ENR 0.002-32 30 MIC Tests
920130	ENROFLOXACIN ENR 0.002-32 100 MIC Tests
920131	ENROFLOXACIN ENR 0.002-32 10 MIC Tests
92014	SPECTINOMYCIN SPC 0.064-1024 30 MIC Tests
920140	SPECTINOMYCIN SPC 0.064-1024 100 MIC Tests
920141	SPECTINOMYCIN SPC 0.064-1024 10 MIC Tests
92015	OXACILLIN OX 0.016-256 30 MIC Tests
920150	OXACILLIN OX 0.016-256 100 MIC Tests
920151	OXACILLIN OX 0.016-256 10 MIC Tests
92016	CEFTIZOXIME CZX 0.016-256 30 MIC Tests
920160	CEFTIZOXIME CZX 0.016-256 100 MIC Tests
920161	CEFTIZOXIME CZX 0.016-256 10 MIC Tests
92017	MECILLINAM MEC 0.002-32 30 MIC Tests
920170	MECILLINAM MEC 0.002-32 100 MIC Tests
920171	MECILLINAM MEC 0.002-32 10 MIC Tests
92018	AMIKACIN AK 0.016-256 30 MIC Tests
920180	AMIKACIN AK 0.016-256 100 MIC Tests
920181	AMIKACIN AK 0.016-256 10 MIC Tests
92019	BACITRACIN BA 0.016-256 30 MIC Tests
920190	BACITRACIN BA 0.016-256 100 MIC Tests

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96021	SALMONELLA PARATYPHI C TOTALE 20 ml
96022	BRUCELLA MELITENSIS 20 ml
96023	BRUCELLA SUIIS 20 ml
96031	SALMONELLA TYPHI H SLIDE 5 ml
96032	SALMONELLA TYPHI O SLIDE 5 ml
96033	SALMONELLA TYPHI TOTALE 5 ml SLIDE
96034	SALMONELLA PARATYPHI AH SLIDE 5 ml
96035	SALMONELLA PARATYPHI AO 5 ml SLIDE
96036	SALMONELLA PARATYPHI A TOTALE 5ml SLIDE
96037	SALMONELLA PARATYPHI BH 5 ml SLIDE
96038	SALMONELLA PARATYPHI BO 5 ml SLIDE
96039	SALMONELLA PARATYPHI B TOTALE 5ml SLIDE
96040	SALMONELLA PARATYPHI CH 5 ml SLIDE
96041	SALMONELLA PARATYPHI CO 5 ml SLIDE
96042	SALMONELLA PARATYPHI C TOTALE 5 ml SLIDE
96043	BRUCELLA TOTALE SLIDE 5 ml SLIDE
96044	BRUCELLA ABORTUS 5 ml SLIDE
96045	BRUCELLA MELITENSIS SLIDE 5 ml
96046	BRUCELLA BENGAL ROSE SLIDE 5 ml
96047	PROTEUS OX2 5 ml SLIDE
96048	PROTEUS OX19 5 ml SLIDE
96049	PROTEUS OXK 5 ml SLIDE
96093	CONTROLLO NEGATIVO/NEGATIVE CONTROL 0.5ml
96096	POSITIVE CONTROL FOR SALMONELLA 0.5ml
96097	POSITIVE CONTROL FOR PROTEUS 0.5ml
96098	POSITIVE CONTROL FOR BRUCELLA 0.5ml
96142	Legionella Latex Kit
96143	CAMPYLOBACTER LATEX KIT
96144	CLOSTRIDIUM DIFFICILE LATEX KIT
96148	SHIGELLA ANTISERUM
96150	E. COLI O157 LATEX KIT
96151	SALMONELLA LATEX KIT
96153	STREPTO B LATEX KIT
96154	STREPTO A LATEX KIT
96155	BENCE JONES LATEX TEST
96316	Clostridium difficile GDH Card
96317	Clostridium Difficile Toxin A+B Card
96318	Giardia Card
96319	Listeria Monocytogenes Card
96320	Salmonella Ag Card
96321	O157 E.coli Card
96401	ONE STEP AMP DRMG SCREEN 20 CARDS
96404	ONE STEP COC DRMG SCREEN
96405	ONE STEP THC DRMG SCREEN
96406	ONE STEP M-AMP DRMG SCREEN 20 CARDS ONE STEP BRUPRENORPHINE DRMG SCREEN 20 CARDS
96411	FECAL OCCULT BLOOD CARD
96415/20	STREPTO A CARD 30 CARDS
96418	Gonorrhea Ag Card
96441	Gardnerella Vaginalis Card
96442	Trichomonas Vaginalis Card
96443	B.J. Free Kappa/Lambda Dipstick
96444	H.PYLORI CARD 20 CARD
96455	HCG URINE/SERUM CARD 50 CARD
96460	HCG URINE/SERUM CARD 100 CARD
96461	MICROALBUMIN CARD URINE 20 Cards
96462	AFP -ALFA FETO CARD 20 CARDS
96465	TUBERCOLOSI CARD 20 CARDS
96466	IgE TOTAL CARD
96468	CEA CARD 20 Cards

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Rev. 31.0 del 08.01.2016

96487	MYOGLOBIN
96488	TROPONIN 20 CARDS
96490	FERRITIN CARD
96495	SIFILIDE CARD 20 CARDS
96498	IM MONONUCLEOSIS INFECTION 20 CARDS
96590	URINE STRIP
96900	GIOTTO READER
96909	BIOMIC V3
96914	BIOMIC V3 AST
96915	BIOMIC V3 ID
96916	BIOMIC V3 CC
96919	AST Software
96931	ID Software
96932	CC Software
96933	Micropiastre 96 pozzetti Software

97800	ROTASTICK ONE STEP KIT 20 Tests
97801	RSV STICK ONE STEP 20 Tests
97802	ROTA/ADENO COMBI STICK ONE STEP 20 Tests
97803	H.PYLORI FECAL Ag ONE STEP 20 Tests
97805	STREP B STICK ONE STEP ASSAY 20 Tests
97807	ADENOSTICK ONE STEP ASSAY 20 Tests
9999	Blank Discs
99003	KPC&MBL disc kit (acc. to EUCAST)
99004	ESBL disc kit (acc. to EUCAST)
99005	ESBL disc kit (acc. to CLSI)
99006	ESBL (Chromos. Ind. AmpC) disc kit (acc. to EUCAST)
99007	KPC&MBL&OXA-48 disc kit (acc. to EUCAST)
99008	ESBL+AmpC screen disc kit
99009	AmpC disc kit

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