
Enterolert*-DW

06-18085-10

IDEXX

IDX 33/03 - 10/13
WATER ANALYSIS METHODS
<http://nf-validation.afnor.org>

The method Enterolert*-DW/Quanti-Tray* for water analysis is granted NF Validation by AFNOR Certification as an alternative method to the standard ISO 7899-2 for enumeration of intestinal enterococci in drinking water (except bottled water) and in bathing water (fresh and marine), under the Certificate number: IDX 33/03 – 10/13. For more information about end of validity, please refer to the certificate NF Validation available on website mentioned above.

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IDEXX Water Quality Control Laboratory is accredited to ISO/IEC 17025:2017

IDEXX

IDEXX Laboratories, Inc., One IDEXX Drive, Westbrook, Maine 04092 USA
idexx.com/water

For Technical Support, please call:

North America: 1 207 556 4496/1 800 321 0207

Latin America: +800 728 2482

Europe: +00800 7274 3399

UK: +44 (0) 1638 676800

China: 400-678-6682 ext. 2

Japan: +81 3 5301 6800

Australia: 1300 443 399

Introduction

Enterolert®-DW (98-18072-00, 98-18074-00) detects *Enterococcus* species in drinking water and in bathing (fresh and marine) water samples in 24 hours. It is based on IDEXX's proprietary Defined Substrate Technology® (DST®). Enterolert-DW, coupled with the IDEXX Quanti-Tray® System, provides quantitative confirmed results in 24 hours. Enterolert-DW utilizes Ortho-Nitrophenyl-β-D-glucoside as a nutrient indicator and incorporates a specifically designed blue background colour in its formulation. When the substrate is metabolized by enterococci, the sample turns from blue to green to indicate a positive detection. Any change from the original colour to green is considered a positive result. No ultraviolet light source is required.

Storage

Store at 2–25°C away from direct light and humidity.

Quanti-Tray® Enumeration Procedure

1. Add contents of one pack to a 100 mL water sample in a sterile vessel.
2. Cap vessel and shake until dissolved.
3. Pour sample/reagent mixture into a Quanti-Tray® or Quanti-Tray®/2000 and seal with an IDEXX Quanti-Tray® Sealer.
4. Place the sealed tray in a 41 ± 0.5°C incubator for 24 hours.
5. Read results according to the Result Interpretation table below. Count the number of positive wells and refer to the MPN table provided with the trays to obtain a Most Probable Number.

Result Interpretation

Appearance	Result
Blue	Negative for enterococci
Green	Positive for enterococci

- Enterolert-DW results are definitive at 24–28 hours. In addition, positive wells for enterococci observed before 24 hours and negative wells observed after 28 hours are also valid.

Procedural Notes

- For all bathing (fresh and marine) water samples, a tenfold dilution must be performed. Multiply the MPN value by the dilution factor to obtain the proper quantitative result.
- Use only sterile, nonbuffered, oxidant-free water for dilutions.
- For comparison, a negative control can be used when interpreting results.
- This insert may not reflect your local regulations. For compliance testing, be sure to follow appropriate regulatory procedures.
- Enterolert-DW can be run in any multiple tube format. *Standard Methods for the Examination of Water and Wastewater*¹ MPN tables should be used to find Most Probable Numbers.
- Enterolert-DW is a primary drinking water test. Enterolert-DW performance characteristics do not apply to samples altered by any pre-enrichment or concentration.
- Aseptic technique should always be followed when using Enterolert-DW. Dispose of in accordance with Good Laboratory Practice.

Quality Control Procedures

1. One of the following quality control procedures is recommended for each lot of Enterolert-DW:
 - A. IDEXX-QC Enterococci DW²
 - B. For American Type Culture Collection (ATCC^{TM,3}) or World Data Centre for Microorganisms (WDCM⁴) bacterial strains:
 - i. Fill two sterile vessels with 100 mL of sterile non-buffered oxidant-free water and inoculate with, *Enterococcus faecium* ATCC 35667/WDCM 00226 or *Enterococcus faecalis* ATCC 33186/WDCM 00210, and *Serratia marcescens* ATCC 43862 or *Escherichia coli* ATCC 11775/ WDCM 00090. These strains are not available from IDEXX and must be purchased from ATCC or other providers.
 - ii. For positive controls, prepare proper dilutions to reach a concentration of approximately 60–80 CFU/100 mL.
 - iii. For non-target strains, prepare proper dilutions to a suggested inoculation up to 10⁶ CFU/100 mL.
2. Follow the Quanti-Tray Enumeration Procedure above.
3. Results should match the Result Interpretation table above.

NOTE: IDEXX internal quality control testing is performed in accordance with ISO 11133:2014. Quality Control Certificates are available at idexx.com/water.

[§]The Enterolert-DW Test is intended to be used only for water quality research and analysis, by technically qualified individuals or under their supervision.

1. Eaton AD, Clesceri LS, Rice EW, Greenberg AE, Franson MAH, eds. *Standard Methods for the Examination of Water and Wastewater*. 24th ed. Washington, DC: American Public Health Association; 2022.

2. IDEXX-QC Enterococci DW—IDEXX Catalog #WQC-ENTDW

3. ATCC is a registered trademark owned by the American Type Culture Collection. American Type Culture Collection 1-800-638-6597 atcc.org

4. WDCM World Data Centre for Microorganisms • wdcm.org

*Enterolert, Defined Substrate Technology, DST and Quanti-Tray are trademarks or registered trademarks of IDEXX Laboratories, Inc. or its affiliates in the United States and/or other countries.

Patent information: idexx.com/patents.

© 2025 IDEXX Laboratories, Inc. All rights reserved.

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