MSRV Agar

Catalog #	Description
3556139	MSRV Agar , ready-to-use, 200 ml x 6 bottles
3564325	MSRV Agar , dehydrated, 500 g

For laboratory use only.

Intended Use

A semi-solid medium used for the isolation of *Salmonella* spp., specifically motile *Salmonella* spp.

Principle

The semi-solid medium enables the migration of *Salmonella* to form clear opaque halos of growth originating from the place of inoculation providing better results than traditional techniques. Due to the presence of magnesium chloride, malachite green, novobiocin, and an incubation temperature of 42°C, the medium inhibits the growth of other bacteria.

Theoretical Composition

Base Medium

Tryptone	4.6 g
Acid hydrolyzed casein	4.6 g
Sodium chloride	7.3 g
Potassium dihydrogen phosphate	1.5 g
Anhydrous magnesium chloride*	10.9 g
Malachite green oxalate	0.037 g
Novobiocin	0.010 g
Agar	2.7 g
Distilled water	1,000 ml

Final pH at 25°C = 5.1–5.4

*Anhydrous magnesium chloride 10.9 g is equivalent to magnesium chloride hexahydrate 29.4 g.

Shelf Life and Storage

Store ready-to-use medium at 2–8°C. Store dehydrated medium at 15–25°C in carefully sealed bottles in a cool, dry place.

Required Materials Not Supplied

This list is not exhaustive.

Equipment

- All usual laboratory equipment
- Incubators or incubation room
- Scales
- Stirrer/homogenizer
- Vortexer

Supplies

- 125 ml Pyrex bottles with autoclave proof stoppers
- Sterile petri dishes (Ø = 90 mm)
- Sterile pipettes (1 ml, etc) or inoculating loops

Precautions

- Respect Good Laboratory Practice (EN ISO 7218). Appropriate protection, such as gloves and lab coats, should be worn when working with potentially infectious live bacteria
- Media that have come in contact with food samples should be considered contaminated and should be disposed of in accordance with local rules and regulations
- MSRV medium in ready-to-use form can be poured directly without previous heating. It is highly fluid, but lumps of agar may be present. This will not have any impact on its performance
- For convenience, if a less fluid medium is preferred, ready-to-use MSRV can be melted beforehand by heating in a water bath at 100°C for not more than 1–2 min and used immediately after tempering to 40 ± 3°C
- Avoid heating or incubation of samples for a prolonged period over 43°C as novobiocin is heat sensitive
- For SDS product safety information and certificate of analysis, visit **bio-rad.com**

Quality Control

Every product manufactured and marketed by Bio-Rad is subject to a quality assurance procedure at all stages, from reception of raw materials through to marketing of the finished products. Each batch of finished product undergoes quality control according to EN ISO 11133 and is marketed only if it satisfies the acceptability criteria. Documentation relative to the production and quality control of each batch is kept on file.

Protocol

Dehydrated Medium Preparation

- Shake before use
- Dissolve 31.5 g of powder in 1 L of sterile distilled water. Soak for 10 min, swirl to mix and bring to a boil
- Cool to 47–50°C
- Add 0.01 g/L of novobiocin
- Pour about 18 ml of melted medium into empty petri dishes

Reconstitution ratio: 31.5 g/1 L (500 g of powder makes 15.9 L of medium)

Sample Preparation and Enrichment

- Prepare and enrich sample according to the standard method applicable to the product concerned

Inoculation and Incubation

- Transfer 3 drops (a total of 0.15 ml) of the sample after enrichment to a dish of MSRV
- Incubate at 41.5 ± 1°C for 24 ± 3 hr

Reading and Interpretation

- Motile *Salmonella* produce a whitish zone of migration over the agar after incubation

References

De Smedt, JM et al. (1986). Rapid *Salmonella* Detection in Foods by Motility Enrichment on a Modified Semi-Solid Rappaport Vassiliadis Medium. J. Food Protection 49 (7): 510-514.

ISO 6579-1:2017/A 1:2020. Microbiology of the food chain – Horizontal method for the detection, enumeration and serotyping of *Salmonella* – Part 1: Detection of *Salmonella* spp. – Amendment 1: Broader range of incubation temperatures, amendment to the status of Annex D, and correction of the composition of MSRV and SC.

NF U47-100 July 2007. Animal health analysis methods – Isolation and identification of *Salmonella* or search for particular serovar's in the animal production sector.

NF U47-101 November 2007. Animal health analysis methods – Isolation and identification of any *salmonella* serotype or of specified *salmonella* serotypes among birds.

NF U47-102 January 2008. Animal health analysis methods – Isolation and identification of any *salmonella* serotype or of specified *salmonella* serotypes among mammals.

Revision History

Release date	Document number	Change
August 2021	5080 Ver A	- Major change - New document design - Document number change – previous version: V5_05/08/11

BIO-RAD is a trademark of Bio-Rad Laboratories, Inc. All trademarks used herein are the property of their respective owner.

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