Brilliance™ Salmonella Agar Base

REF CM1092B & CM1092T

EN

Intended Use

Brilliance Salmonella Agar Base (CM1092B and CM1092T) devices are a selective and differential chromogenic medium for the presumptive identification of *Salmonella* species from faecal samples. The *Brilliance*™ Salmonella Agar Base (CM1092B and CM1092T) devices are used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having bacterial infections. The devices are for professional use only, not automated nor are they companion diagnostics.

The devices can also be used for testing non-clinical samples (food and environmental).

Summary and Explanation

Salmonellosis, the disease caused by *Salmonella* bacteria, is inherently pathogenic to humans. These bacteria can be ingested from contaminated food, or water - via the faecal-oral route – and will invade the gastrointestinal tract to cause enteric infection, with diarrhoea being the most common symptom. *Salmonella* infections can become particularly invasive in infants, older adults and can be serious in immunocompromised individuals¹.

Salmonella spp. are listed as one of four key global causes of diarrhoeal disease worldwide^{2,3}. In Europe, *Salmonella* infection is the second most common foodborne disease⁴. and an estimated 550 million people fall ill with diarrhoeal disease every year from unsafe food³.

It is important to identify *Salmonella* infections from clinical samples to ensure that the correct diagnosis is achieved, and a suitable treatment protocol can be implemented for the patient. Gram-negative bacteria are easily identifiable following Gram-stain where the crystal violet stain will not bind to the thick peptidoglycan layer present in the bacterial cell wall and bacteria will present red/pink under the microscope. Following a Gram stain, a morphological analysis of the bacteria can help to further differentiate genus and species⁵. Thus, *Brilliance* Salmonella Agar Base devices (CM1092B and CM1092T), can help to further differentiate species of *Salmonella* and should be incorporated into clinical diagnostic workflows.

Principle of Method

Differentiation of *Salmonella* from the other organisms that grow on *Brilliance* Salmonella Agar is achieved through the inclusion of two chromogens that also target specific enzymes: caprylate esterase and β -glucosidase. Caprylate esterase is an enzyme present in all *salmonellae* as well as some species of *Klebsiella*, *Enterobacter* and *Proteus*. Organisms possessing caprylate esterase cleave the chromogen to release an insoluble purple chromophore. As the cells grow, the chromophore builds up and produces a purple-coloured colony. Some *Enterobacteriaceae*, including *Klebsiella* and *Enterobacter* but not *Salmonella*, possess β -glucosidase. If these organisms grow, they will form blue or dark blue colonies, even if they are esterase positive, which make them easy to differentiate from purple *Salmonella* colonies.

Typical Formula

	grams per litre
Salmonella Inhibogen™ mix	14
Chromogenic mix	25
Agar	15

Materials Provided

CM1092B: 500g of Salmonella Agar Base powder that yields approximately 9.3L after reconstitution.

CM1092T: 5kg of Salmonella Agar Base powder that yields approximately 74L after reconstitution.

Materials Required but Not Supplied

- Inoculating loops, swabs, collection containers
- Incubators
- Quality control organisms
- Salmonella Selective Supplement (SR0194E)
- Petri dish

Storage

- Store product in its original packaging between 10°C and 30°C.
- Keep container tightly closed.
- The product may be used until the expiry date stated on the label.
- Protect from moisture.
- Store away from light.
- Allow reconstituted product to equilibrate to room temperature before use.

Once reconstituted, store media between 2°C and 8°C.

Warnings and Precautions

- Do not inhale. May cause allergy or asthma symptoms or difficulty breathing if inhaled.
- Causes serious eye irritation.
- May cause an allergic skin reaction.
- If on skin wash with plenty of soap and water.
- If in eyes, rinse cautiously with water for several minutes.
- Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists, seek medical advice/attention.
- If inhaled, if breathing is difficult, remove subject to fresh air and keep in a position comfortable for breathing. If experiencing respiratory symptoms, call a POISON CENTER or doctor/physician.
- For in vitro diagnostic use only.
- For professional use only.
- Inspect the product packaging before first use.
- Do not use the product if there is any visible damage to the packaging (pot or cap).
- Do not use the product beyond the stated expiry date.
- Do not use the device if signs of contamination are present.
- It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.
- Ensure the lid of the container is kept tightly closed after first opening and between use to minimise moisture ingress, which may result in incorrect product performance.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product (www.thermofisher.com).

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

Specimen Collection, Handling and Storage

Specimen should be collected and handled following local recommended guidelines, such as the UK Standards for Microbiology Investigations (UK SMI) B 29.

Procedure

Suspend 27g in 500ml of distilled water. Add the contents of 1 vial of Salmonella Selective Supplement (SR0194E) reconstituted as directed. With frequent agitation, bring to the boil to dissolve completely. Cool to 50°C. Mix well and pour into sterile Petri dishes. DO NOT AUTOCLAVE. DO NOT OVERHEAT.

Note: It is critical that the selective supplement is added prior to heating.

Interpretation

After reconstitution, *Salmonella* will appear as purple/pink coloured colonies.

Some Enterobacteriaceae, including *Klebsiella* and *Enterobacter* that possess β -glucosidase will appear as blue colonies.

Quality Control

It is the responsibility of the user to perform Quality Control testing taking into account the intended use of the medium, and in accordance with any local applicable regulations (frequency, number of strains, incubation temperature etc.).

The performance of this medium can be verified by testing the following reference strains.

Incubation Conditions: 24 ± 2 h @ 37°C aerobic

Positive Controls	
Colony count is ≥ 50% of the control medium count	
<i>Salmonella Typhimurium</i> ATCC@14028	1-2mm purple/pink colonies
<i>Salmonella enteritidis</i> ATCC@13076	1-2mm purple/pink colonies
<i>Klebsiella pneumoniae</i> ATCC@29665	1-3mm mucoid, blue colonies
<i>Salmonella arizonae</i> ATCC@13314	1-4 mm purple/pink colonies
Negative Controls	
<i>Enterococcus faecalis</i> ATCC@29212	No growth
<i>Proteus mirabilis</i> NCTC10975	No growth or 0.5mm straw colonies

<i>Escherichia coli</i> ATCC®25922	No growth or 1-2mm cream colonies
<i>Pseudomonas aeruginosa</i> ATCC®27853	No growth or 1-2mm pink colonies

Limitations

It should be noted that, as with all chromogenic media, organisms with atypical enzyme patterns may give anomalous reactions on *Brilliance* Salmonella Agar. In the context of NF VALIDATION, some strains belonging to the serovar *Salmonella* Dublin may show weak, or no, pigmentation because of their low esterase activity.

All identifications are presumptive and should be confirmed using appropriate methods. Medium choice and incubation protocol will depend upon the specimen type. A non-selective medium should be inoculated alongside the selective medium.

Performance Characteristics

Accuracy has been demonstrated through review of the QC data. Correct detection of *Salmonella* strains is confirmed by the inclusion of a well-characterised isolate in the QC processes performed as part of the manufacture of each batch of the devices, which must meet the defined acceptance criteria. The precision of CM1092B was demonstrated by an overall pass rate of 100% obtained for the product over 4 months of testing (02.12.2021 – 11.04.2022; 10 batches). The precision of CM1092T was demonstrated by an overall pass rate of 100% obtained for the product over 2 months of testing (15.02.2022 – 14.04.2022; 10 batches). This shows that the performance is reproducible.


Brilliance Salmonella Agar Base (CM1092B and CM1092T) devices are tested in-house as part of the QC process since the products were launched in 2006 and 2013 respectively. When using 10-100 cfu inoculum of *Salmonella enteritidis* ATCC®13076, *Salmonella Typhimurium* ATCC®14028 and *Salmonella arizonae* ATCC®13314 for CM1092B/CM1092T and incubating the device at 36°C for 18-24 hours, the user can recover organisms with colony size and morphology that meets the defined acceptance criteria.

Bibliography

1. CDC. 2020. "Salmonellosis (Nontyphoidal) - Chapter 4 - 2020 Yellow Book | Travelers' Health | CDC." 2020. <https://wwwnc.cdc.gov/travel/yellowbook/2020/travel-related-infectious-diseases/salmonellosis-nontyphoidal>.
2. Kirk, Martyn D, Sara M Pires, Robert E Black, Marisa Caipo, John A Crump, Brecht Devleesschauwer, Dörte Döpfer, et al. 2015. "World Health Organization Estimates of the Global and Regional Disease Burden of 22 Foodborne Bacterial, Protozoal, and Viral Diseases, 2010: A Data Synthesis." <https://doi.org/10.1371/journal.pmed.1001921>.
3. World Health Organisation. 2018. "Salmonella (Non-Typhoidal)." World Health Organisation. 2018. [https://www.who.int/news-room/fact-sheets/detail/salmonella-\(non-typhoidal\)](https://www.who.int/news-room/fact-sheets/detail/salmonella-(non-typhoidal)).
4. European Food Safety Authority. 2014. "The European Union Summary Report on Trends and Sources of Zoonoses, Zoonotic Agents and Food-Borne Outbreaks in 2012." EFSA Journal 12 (2). <https://doi.org/10.2903/J.EFSA.2014.3547>.
5. Public Health England. 2021. "Identification of Salmonella Species". UK Standards for Microbiology Investigations ID 24 (4). <https://www.gov.uk/government/publications/smi-id-24-identification-of-salmonella-species>.

Symbol Legend

Symbol	Definition
	Catalogue number
	In Vitro Diagnostic Medical Device
	Batch code
	Temperature limit
	Use-by date
	Keep away from sunlight
	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests

	Do not use if packaging damaged and consult instructions for use
	Manufacturer
	Authorized representative in the European Community/ European Union
	European Conformity Assessment
	UK Conformity Assessment
	Unique device identifier
	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union
	Made in the United Kingdom

ATCC Licensed Derivative[®]

©2022 Thermo Fisher Scientific Inc. All rights reserved.
ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection.
All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, UK

For technical assistance please contact your local distributor.

Revision information

Version	Date of issue and modifications introduced
2.0	2024-02-28

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1. CDC. 2020. "Salmonellosis (Nontyphoidal) - Chapter 4 - 2020 Yellow Book | Travelers' Health | CDC." 2020. <https://wwwnc.cdc.gov/travel/yellowbook/2020/travel-related-infectious-diseases/salmonellosis-nontyphoidal>.
2. Kirk, Martyn D, Sara M Pires, Robert E Black, Marisa Caipo, John A Crump, Brecht Devleesschauwer, Dörte Döpfer, et al. 2015. "World Health Organization Estimates of the Global and Regional Disease Burden of 22 Foodborne Bacterial, Protozoal, and Viral Diseases, 2010: A Data Synthesis." <https://doi.org/10.1371/journal.pmed.1001921>.
3. World Health Organisation. 2018. "Salmonella (Non-Typhoidal)," World Health Organisation. 2018. [https://www.who.int/news-room/fact-sheets/detail/salmonella-\(non-typhoidal\)](https://www.who.int/news-room/fact-sheets/detail/salmonella-(non-typhoidal)).
4. European Food Safety Authority. 2014. "The European Union Summary Report on Trends and Sources of Zoonoses, Zoonotic Agents and Food-Borne Outbreaks in 2012." EFSA Journal 12 (2). <https://doi.org/10.2903/J.EFSA.2014.3547>.
5. Public Health England. 2021. "Identification of Salmonella Species". UK Standards for Microbiology Investigations ID 24 (4). <https://www.gov.uk/government/publications/smi-id-24-identification-of-salmonella-species>.

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
1. CDC. 2020. “Salmonellosis (Nontyphoidal) - Chapter 4 - 2020 Yellow Book | Travelers' Health | CDC.” 2020. <https://wwwnc.cdc.gov/travel/yellowbook/2020/travel-related-infectious-diseases/salmonellosis-nontyphoidal>.
2. Kirk, Martyn D, Sara M Pires, Robert E Black, Marisa Caipo, John A Crump, Brecht Devleesschauwer, Dörte Döpfer, et al. 2015. “World Health Organization Estimates of the Global and Regional Disease Burden of 22 Foodborne Bacterial, Protozoal, and Viral Diseases, 2010: A Data Synthesis.” <https://doi.org/10.1371/journal.pmed.1001921>.
3. World Health Organisation. 2018. “Salmonella (Non-Typhoidal).” World Health Organisation. 2018. [https://www.who.int/news-room/fact-sheets/detail/salmonella-\(non-typhoidal\)](https://www.who.int/news-room/fact-sheets/detail/salmonella-(non-typhoidal)).
4. European Food Safety Authority. 2014. “The European Union Summary Report on Trends and Sources of Zoonoses, Zoonotic Agents and Food-Borne Outbreaks in 2012.” EFSA Journal 12 (2). <https://doi.org/10.2903/J.EFSA.2014.3547>.
5. Public Health England. 2021. “Identification of Salmonella Species”. UK Standards for Microbiology Investigations ID 24 (4). <https://www.gov.uk/government/publications/smi-id-24-identification-of-salmonella-species>.

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Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, JK

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Bibliografi

1. CDC. 2020. "Salmonellosis (Nontyphoidal) - Chapter 4 - 2020 Yellow Book | Travelers' Health | CDC." 2020. <https://wwwnc.cdc.gov/travel/yellowbook/2020/travel-related-infectious-diseases/salmonellosis-nontyphoidal>.
2. Kirk, Martyn D, Sara M Pires, Robert E Black, Marisa Caipo, John A Crump, Brecht Devleesschauwer, Dörte Döpfer, et al. 2015. "World Health Organization Estimates of the Global and Regional Disease Burden of 22 Foodborne Bacterial, Protozoal, and Viral Diseases, 2010: A Data Synthesis." <https://doi.org/10.1371/journal.pmed.1001921>.
3. World Health Organisation. 2018. "Salmonella (Non-Typhoidal)." World Health Organisation. 2018. [https://www.who.int/news-room/fact-sheets/detail/salmonella-\(non-typhoidal\)](https://www.who.int/news-room/fact-sheets/detail/salmonella-(non-typhoidal)).
4. European Food Safety Authority. 2014. "The European Union Summary Report on Trends and Sources of Zoonoses, Zoonotic Agents and Food-Borne Outbreaks in 2012." *EFSA Journal* 12 (2). <https://doi.org/10.2903/J.EFSA.2014.3547>.
5. Public Health England. 2021. "Identification of Salmonella Species". UK Standards for Microbiology Investigations ID 24 (4). <https://www.gov.uk/government/publications/smi-id-24-identification-of-salmonella-species>.

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Food testing

Validated detection and identification of *Salmonella* species in 2 days

Thermo Scientific Salmonella Precis Method

Reinventing the standard for culture-based *Salmonella* testing

Salmonellosis remains one of the most common forms of food poisoning, caused by the ingestion of foods contaminated with bacteria of the genus *Salmonella*.


Traditional culture-based methods for detecting *Salmonella* are time consuming, taking up to 4 to 5 days to complete, and require large numbers of sample manipulations and identification steps.

The Thermo Scientific™ *Salmonella* Precis™ Method is the next generation in rapid, culture-based testing and offers a quick and easy method for the enrichment, detection and confirmation of

Salmonella species from food, animal feed, and environmental samples.

The simplified workflow has been expertly designed to overcome sensitivity and specificity issues, without the need to carry out non-standard handling steps, generating results in just two days*.

Thermo Scientific™ Oxoid™ *Brilliance*™ *Salmonella* Agar plates may also be used as part of the ISO 6579-1:2017 horizontal method for the detection of *Salmonella* spp. in food, animal feed, and environmental samples.

Key Benefits:

- Validated according to ISO 16140-2:2016 against the ISO 6579-1:2017 reference method and certified by NF Validation (AFNOR)
- Simple procedure—no specialized equipment required
- Single overnight enrichment
- For *Salmonella* detection PrecisBlue Supplement can be added to the diluted sample for workflow confidence - inert blue dye indicates supplemented buffered peptone water
- Single sample transfer
- Single 24-hour plate incubation
- Broths and plates incubated at the same temperature for workflow simplicity
- Quick and convenient confirmation: Oxoid *Salmonella* Latex Test, SureTect *Salmonella* PCR Assay, or ISO 6579-1:2017 standard tests
- Reduced time to result: 2 days compared with up to 5 days for standard culture methods
- *Brilliance* *Salmonella* Agar contains novel Inhibigen™ technology, giving targeted specificity and reduced background flora

*When confirming presumptive-positive colonies with the Oxoid *Salmonella* Latex Test.

A rapid, simple, and robust culture-based workflow

The Salmonella Precis Method has been validated and certified by NF Validation (AFNOR) according to ISO 16140-2:2016. The Salmonella Precis method was compared to the ISO 6579-1:2017 reference method for the detection of *Salmonella* spp.


For flexibility, multiple validated confirmation options are available including the rapid Oxoid Salmonella Latex Test, the Suretect

Salmonella species PCR Assay, the Microbact GNB 24E kit, or by tests outlined in ISO 6579-1:2017.

AFNOR Certification validation certificate No. UNI 03/06-12/07 is available in PDF format from the AFNOR website

<https://nf-validation.afnor.org/en/food-industry/salmonella-spp/>

Workflow overview for Salmonella Precis Method validated according to NF VALIDATION for AFNOR Certification


*with or without probiotics

Reactions on Brilliance Salmonella Agar

	Colony colour/appearance		
	Purple	Blue	Colorless
Enzyme targeted by chromogen	<i>Salmonella</i> (including lactose-positive <i>Samonella</i>)	<i>Klebsiella</i> , <i>Enterobacter</i> , <i>Serratia</i>	<i>Citrobacter</i> , other bacteria and yeasts
Esterase	+	-/+	-
β-glucosidase	-	+	-

Escherichia coli and other bacteria and yeasts are inhibited by the combination of the Inhibigen and other selective agents in the medium.

Example results – Mixed culture

1. *Klebsiella* colony
2. *Salmonella* colony

Product order information

Product description		Format	Order code	
Oxoid Culture Media	Buffered Peptone Water	Buffered Peptone Water (ISO)	500 g, makes 25 L	CM1049B
		Buffered Peptone Water (ISO-meat peptone)	500 g, makes 25 L	CM1211B
	Novobiocin Supplement	Novobiocin Supplement - freeze-dried	10 vials of 10 mg	SR0181E
		Novobiocin Supplement - liquid (40 mL/vial)	10 vials of 40 mg	SR0249A
	PrecisBlue Supplement		1 vial makes 18 L	SR0259A
	Vancomycin Supplement	Vancomycin Supplement - freeze-dried	10 vials of 3 mg	SR0186E
		Vancomycin Supplement - freeze-dried	10 vials of 5 mg	SR0247E
	ONE Broth Salmonella	ONE Broth Salmonella Base	500 g, makes 20 L	CM1091B
		ONE Broth Salmonella Supplement	10 vials, each for 225 mL	SR0242E
	<i>Brilliance</i> Salmonella Agar	<i>Brilliance</i> Salmonella Agar Base	500 g, makes 9.3 L	CM1092B
		<i>Brilliance</i> Salmonella Agar Selective Supplement	10 vials, each for 500 mL	SR0194E
	Oxoid Salmonella Latex Kit		100 tests	DR1108A

Please note that a range of alternative formats of culture media such as Bagged Enrichment Media and Prepared Plate Media are available. Please contact your local representative or technical services to find out more.

 For more information about the Thermo Scientific Salmonella Precis Method, and other rapid culture media methods for detecting foodborne pathogens visit thermofisher.com/precis

	Document Owner Department: QC	MBD-BT-SPEC-0154
		Rev 03
		Page 1 of 4
OXOID QUALITY ASSURANCE PRODUCT SPECIFICATION		
LAURYL TRYPTOSE BROTH (CM0451)		

LAURYL TRYPTOSE BROTH

CM0451

Formula

Tryptose	grams per litre	20.0
Lactose		5.0
Sodium chloride		5.0
Di-potassium hydrogen phosphate		2.75
Potassium dihydrogen phosphate		2.75
Sodium lauryl sulphate		0.1

Directions

Add 35.6g to 1 litre of distilled water and distribute into final containers with fermentation tubes (Durham). Sterilize by autoclaving at 121°C for 15 minutes.

Physical Characteristics

Straw, free flowing powder
 Colour on reconstitution - straw 2-3
 Moisture level - less than or equal to 7%
 pH - 6.8 ± 0.2 at 25°C
 Clarity - clear (single and double strength broth)

Bacteriological Tests Using Optimum Inoculum Dilution

Control Medium: Tryptone Soya Agar

Reactions after incubation at 30 ± 2°C for 24 ± 2 hours

Medium is challenged with 10-100 colony forming units

Enterobacter aerogenes ATCC®13048 Turbid growth and gas

A satisfactory result is represented by visible growth and gas.

Reactions after incubation at 37 ± 2°C for 24 ± 2 hours

Medium is challenged with 10-100 colony forming units

Enterobacter aerogenes ATCC®13048 Turbid growth and gas

	Document Owner Department: QC	MBD-BT-SPEC-0154
		Rev 03
		Page 2 of 4
OXOID QUALITY ASSURANCE PRODUCT SPECIFICATION		
LAURYL TRYPTOSE BROTH (CM0451)		

A satisfactory result is represented by visible growth and gas.

Reactions after incubation at 44 ± 2°C for 24 ± 2 hours

Medium is challenged with 1E+04 to 1E+06 colony forming units

E. coli confirmation test

<i>Escherichia coli</i>	ATCC®25922	Turbid growth and gas
<i>Escherichia coli</i>	ATCC®8739	Turbid growth and gas
<i>Enterobacter aerogenes</i>	ATCC®13048	No growth to slight turbid growth, no gas

A satisfactory result for *Escherichia coli* is represented by visible growth and gas.

Reactions after incubation at 30 ± 2°C for 48 ± 2 hours

Medium is challenged with 1E+04 to 1E+06 colony forming units

<i>Staphylococcus aureus</i>	ATCC®25923	No growth
------------------------------	------------	-----------

Negative strains are inhibited.

Reactions after incubation at 37 ± 2°C for 48 ± 2 hours

Medium is challenged with 1E+04 to 1E+06 colony forming units

<i>Staphylococcus aureus</i>	ATCC®25923	No growth
------------------------------	------------	-----------

Negative strains are inhibited.

Testing performed in accordance with ISO11133:2014

Reactions after incubation at 30 ± 2°C for 24 ± 2 hours

Medium is challenged with 10-100 colony forming units

<i>Escherichia coli</i>	ATCC®25922	WDCM00013	Turbid growth and gas (indole positive)
<i>Escherichia coli</i>	ATCC®8739	WDCM00012	Turbid growth and gas
<i>Citrobacter freundii</i>	ATCC®43864	WDCM00006	Turbid growth and gas

	Document Owner Department: QC	MBD-BT-SPEC-0154
		Rev 03
		Page 3 of 4
OXOID QUALITY ASSURANCE PRODUCT SPECIFICATION		
LAURYL TRYPTOSE BROTH (CM0451)		

A satisfactory result is represented by visible growth.

After incubation add a few drops of Kovacs reagent to the *Escherichia coli* ATCC®25922 WDCM00013 tube to assess indole production. A satisfactory result is represented by an indole positive reaction (red layer after the addition of Kovacs reagent).

Reactions after incubation at 37 ± 2°C for 24 ± 2 hours

Medium is challenged with 10-100 colony forming units

<i>Escherichia coli</i>	ATCC®25922	WDCM00013	Turbid growth and gas (indole positive)
<i>Escherichia coli</i>	ATCC®8739	WDCM00012	Turbid growth and gas

A satisfactory result is represented by visible growth.

After incubation add a few drops of Kovacs reagent to the *Escherichia coli* ATCC®25922 WDCM00013 tube to assess indole production. A satisfactory result is represented by an indole positive reaction (red layer after the addition of Kovacs reagent).

Reactions after incubation at 30 ± 2°C for 48 ± 2 hours

Medium is challenged with 1E+04 to 1E+06 colony-forming units

<i>Enterococcus faecalis</i>	ATCC®29212	WDCM00087	No growth
<i>Enterococcus faecalis</i>	ATCC®19433	WDCM00009	No growth

Negative strains are inhibited.

Reactions after incubation at 37 ± 2°C for 48 ± 2 hours

Medium is challenged with 1E+04 to 1E+06 colony-forming units

<i>Enterococcus faecalis</i>	ATCC®29212	WDCM00087	No growth
<i>Enterococcus faecalis</i>	ATCC®19433	WDCM00009	No growth

Negative strains are inhibited.

	Document Owner Department: QC	MBD-BT-SPEC-0154
		Rev 03
		Page 4 of 4
OXOID QUALITY ASSURANCE PRODUCT SPECIFICATION		
LAURYL TRYPTOSE BROTH (CM0451)		

Revision History

Section / Step	Description of Change	Reason for Change	Reference
Reactions after incubation at $44 \pm 2^{\circ}\text{C}$ for 24 ± 2 hours	Remove <i>Citrobacter freundii</i> ATCC®43864	Change control	MOC-2024-1997

	Document Owner Department: QC	MBD-BT-SPEC-0870
		Rev 01
		Page 1 of 3
OXOID QUALITY ASSURANCE PRODUCT SPECIFICATION		
PRECISBLUE SUPPLEMENT SR0259A		

PRECISBLUE SUPPLEMENT

SR0259A

Formula

Antibiotic/Dye mix

Antibiotic

216mg/vial

PrecisBlue Supplement is used with the SureTect™ Salmonella, SureTect™ Cronobacter, Precis™ Salmonella and Precis™ Cronobacter workflows. Please refer to the instructions for use (IFU) and the appropriate workflow to determine the required volume of supplement.

Directions

For testing in accordance with the Precis™ workflows, aseptically add 500µl of PrecisBlue Supplement per 225ml of sterile Buffered Peptone Water (ISO) (CM1211) for testing *Salmonella* spp. and non-target organisms. Alternatively, aseptically add 250µl of PrecisBlue selective supplement per 225 ml of Buffered Peptone Water (ISO) (CM1211) for testing *Cronobacter* spp.

Physical Characteristics

Dark blue to black liquid
Sterility - passes test

Microbiological Tests Using Optimum Inoculum Dilution

Control Medium: Tryptone Soya Agar

Reactions after incubation at 42°C for 18 hours

Tested with Buffered Peptone Water (ISO) CM1211

Medium challenged with 10-100 colony-forming units

<i>Salmonella typhimurium</i>	ATCC®14028	Turbid growth
<i>Salmonella enteritidis</i>	ATCC®13076	Turbid growth
<i>Cronobacter sakazakii</i>	ATCC®29544	Turbid growth

A satisfactory result is represented by visible growth.

Medium challenged with 1E+04 to 1E+06 colony-forming units

<i>Enterococcus faecalis</i>	ATCC®29212	No growth to minimal growth
<i>Staphylococcus aureus</i>	ATCC®25923	No growth

ThermoFisher SCIENTIFIC	Document Owner Department: QC	MBD-BT-SPEC-0870
		Rev 01
		Page 2 of 3
OXOID QUALITY ASSURANCE PRODUCT SPECIFICATION		
PRECISBLUE SUPPLEMENT SR0259A		

Negative strains are inhibited or shall produce minimal growth.

Equivalent results are obtained after incubation at 37°C for 18 hours with matrix additions.

	Document Owner Department: QC	MBD-BT-SPEC-0870
		Rev 01
		Page 3 of 3
OXOID QUALITY ASSURANCE PRODUCT SPECIFICATION		
PRECISBLUE SUPPLEMENT SR0259A		

Revision History

Section / Step	Description of Change	Reason for Change	Reference
All	Creation of new document	N/A	N/A