DESCRIZIONE

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BIBLIOGRAFIA

1. Barritt, M.M. (1936) The intensification of the Voges-Proskauer reaction by the addition of α -naphthol. J. Pathol. Bacteriol. 42:441-454.
2. Edwin H. Lenette (1995) Manual of Clinical Microbiology.
3. Blazevic, D.J., and Ederer, G.M. (1975) Principles of biochemical tests in diagnostic microbiology. 63-67. New York, John Wiley & Sons.

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Rev.1 / 30.04.2014

VP Reagent

Reagents for the Voges-Proskauer test,
in dropper bottles.

DESCRIPTION

VP Reagent contains 6% (w/v) solution of alpha naphthol and 40% (w/v) solution of sodium hydroxide used for performing the Voges-Proskauer test for the differentiation of members of Enterobacteriaceae.

KIT CONTENTS

Each package contains:

- 10 dropper bottles of 2.5 ml of sodium alpha naphthol (Alpha naphthol, red cap).
- 10 dropper bottles of 2.5 ml of sodium hydroxide (NaOH 40%, blue cap).
- 1 instruction sheet.

METHOD PRINCIPLE

The Voges-Proskauer test allows the differentiation of microorganisms belonging to the Enterobacteriaceae family on the basis of the capability to utilize the butylene glycol pathway to ferment glucose. The positive reaction to the Voges-Proskauer test is denoted by the development of a red color obtained from acetyl-methyl-carbinol oxidation in the presence of alkali and oxygen

COMPOSITION

Alpha naphthol: 6% (w/v) solution of alpha naphthol in ethyl alcohol. NaOH 40%: 40% (w/v) solution of sodium hydroxide.

TEST PROCEDURE AND RESULTS INTERPRETATION

Prior to proceed with the following method bring the dropper bottles to a suitable temperature.

Inoculate buffered peptone-glucose broth (MR-VP MEDIUM, ref. 24149) with a pure culture of the organism in question. Incubate the culture at 37°C for no less than 48 hours. Add 3 drops of alpha naphthol solution and 2 drops of sodium hydroxide solution to 3 mL of broth culture. The positive test is given by the development of a bright red colour 15 minutes after addition of reagents.

QUALITY CONTROL FOR THE USER

1. Appearance control: the alpha naphthol solution appears limpid and beige; the sodium hydroxide solution appears clear limpid colourless.
2. Microbiological control.

Control strains	ATCC®	VP test	
<i>Enterobacter cloacae</i>	ATCC® 13047	Red colour	Positive
<i>Escherichia coli</i>	ATCC® 25922	No colour development	Negative

LIMITATIONS

One should allow at least 15 min for color to develop before considering the Voges-Proskauer test negative.

PRECAUTIONS

The product VP Reagent, is classified as hazardous under current legislation; for its use it is recommend to consult the MSDS. VP Reagent is a disposable device to be used only for in vitro diagnostic use, is intended for professional use and should be used in the laboratory by properly trained operator, with approved methods of aseptic and safety against pathogens.

STORAGE AND TRANSPORT CONDITIONS

Store VP Reagent at 2-8°C away from light in its original package, until the expiry date shown on the label. However, our stability studies have shown that the storage or transport at 18-25°C for 4 days, or at 35-39°C for 48 hours, do not alter in any way the performance of the product. Eliminate if signs of deterioration or contamination are evident.

ELIMINATING USED MATERIAL

After use, used VP Reagent and the material that has come into contact with the sample must be decontaminated and disposed of in accordance with the laboratory procedures for the decontamination and disposal of potentially infected material.

BIBLIOGRAPHY

1. Barritt, M.M. (1936) The intensification of the Voges-Proskauer reaction by the additon of α -naphthol. J. Pathol. Bacteriol. 42:441-454.
2. Edwin H. Lenette (1995) Manual of Clinical Microbiology.
3. Blazevic, D.J., and Ederer, G.M. (1975) Principles of biochemical tests in diagnostic microbiology. 63-67. New York, John Wiley & Sons.

PRESENTATION

Product	Ref.	Contents
VP Reagent	87002	10 x 2.5 ml 10 x 2.5 ml

TABLE OF SYMBOLS

LOT	Batch code	IVD	In vitro Diagnostic Medical Device		Manufacturer		Use by		Fragile, handle with care
REF	Catalogue number		Temperature limitation		Contains Sufficient for <n> tests		Caution, consult accompanying documents		

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LIOFILCHEM® S.r.l.