MacConkey Broth

Differential medium for detection of coliform bacteria.

DESCRIPTION

MacConkey Broth is a liquid medium used for the detection and enumeration of coliforms from water, and milk.

TYPICAL FORMULA*	(g/l)
Peptone	20.0
Lactose	10.0
Sodium Chloride	5.0
Bile Salts No.3	5.0
Bromocresol Purple	0.01
Final pH 7.4 ± 0.2 at 25°C	

*Adjusted and/or supplemented as required to meet performance specifications.

METHOD PRINCIPLE

Peptone provides amino acids, nitrogen, carbon, vitamins and minerals for organisms growth. Lactose is the fermentable carbohydrate. Sodium Chloride maintain the osmotic balance of the medium. Bile salts inhibit the growth of Gram-positive bacteria. Bromocresol purple is the pH indicator.

PREPARATION

Dehydrated medium

Suspend 40.0 g of the powder in 1 liter of distilled or deionized water. Mix well. Heat to boil shaking frequently until completely dissolved. Dispense into tubes containing Durham tubes. Autoclave at 121°C for 15 minutes.

Ensure there are no air bubbles in the Durham tube prior to use as it may lead to a misinterpretation (false-positive reaction in gas production) of results. If needed, invert the tube to allow the air being released from the Durham tube.

TEST PROCEDURE

For the presumptive coliform examination, inoculate a series of tubes with the test sample and its dilutions and incubate at 35-37°C for 18-24 hours.

The number of total coliforms is estimated by counting the number of tubes giving a positive reaction for acid (yellow colouration of the broth) and gas production and comparing the pattern of positive results with standard charts to find MPN (Most Probable Number) of total 100 ml sample.

For the differential test, subculture from each tube showing positive reaction into a fresh tube of MacConkey Broth. Incubate at 42-44°C for up to 48 h to detect *Escherichia coli*.

INTERPRETING RESULTS

Turbidity of MacConkey Broth indicates microbial growth.

Acid production due to lactose fermentation causes a color change of the medium to yellow. Gas is also produced, collecting in Durham tubes.

The ability to form gas in MacConkey Broth at 44°C is specific for *E. coli*, which is indicative of faecal pollution of the original sample.

Non-fermenting Gram-negative organisms produce good growth but without significant changes in color nor gas formation.

STORAGE

The powder is very hygroscopic, store the powder at 10-30°C, in a dry environment, in its original container tightly closed. Do not use the product beyond its expiry date on the label or if product shows any evidence of contamination or any sign of deterioration.

SHELF LIFE

Dehydrated medium: 4 years.

QUALITY CONTROL

Appearance of Dehydrated Medium: Free-flowing, homogeneous, light beige.

Appearance of Prepared Medium: Clear, purple.

Expected Cultural Response:

Control strain		Inoculum	Incubation	Growth	Acid	Gas
Salmonella Typhimurium	WDCM 00031 (ATCC® 14028; NCTC 12023)	≤100 CFU	18-24 h / 35-37°C	Good	Neg. (purple)	Neg.
Escherichia coli	WDCM 00012 (ATCC® 8739; NCTC 12923)	≤100 CFU	18-24 h / 42-44°C	Good	Pos. (yellow)	Pos.
Staphylococcus aureus	WDCM 00032 (ATCC® 6538; NCTC 10788)	1000 CFU	48 h / 35-37°C	Inhibition	—	—

Please refer to the actual batch related Certificate of Analysis (CoA).

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BIBLIOGRAPHY

See the references at the end of this document.

TABLE OF SYMBOLS

See the table of symbols at the end of this document.

The product is available in the configurations listed below. There may be additional product ref. numbers as well. For an updated listing of available products, visit liofilchem.com

Product	Format	Packaging	Ref.
MacConkey Broth	Dehydrated medium	100 g	620171
		500 g	610171

Significant changes from previous version:

Document	Release Date	Change Summary
610171_IFU-0	2022-11-04	Document creation

This IFU document and the SDS are available from the online Support Center:

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liofilchem.com/ifu-sds

References / Riferimenti

1. European Pharmacopoeia (EP) 2.6.13. Microbiological examination of non-sterile products: Test for specified microorganisms.
2. United States Pharmacopoeia (USP) <62> Microbiological examination of non-sterile products: Test for specified microorganisms.
3. Japanese Pharmacopoeia (JP) 4.05 Microbiological examination of non-sterile products: Test for specified microorganisms.
4. Departments of the Environment, Health, Social Security and Public Health Laboratory Service (1982) The Bacteriological Examination of Drinking Water Supplies, Report No. 71, HMSO London.
5. World Health Organization (1963) International Standards for Drinking Water, 2nd ed., WHO, Geneva.
6. Murray, Baron, Jorgensen, Landry and Pfaller ed. (2007) Manual of clinical microbiology, 9th ed. American Society for Microbiology, Washington, D.C.
7. MacConkey, A. (1908) Bile salt media and their advantages in some bacteriological examinations. J. Hyg. 8: 322-334.
8. MacConkey, A. (1905) Lactose-fermenting bacteria in faeces. J. Hyg. 8: 333-379.

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liofilchem@liofilchem.com

TCBS Agar

Selective medium for detection of enteropathogenic *Vibrio* spp from clinical and nonclinical samples, according to ISO 21872.

DESCRIPTION

Thiosulfate Citrate Bile Sucrose (TCBS) Agar is a medium used for the selective isolation and cultivation of vibrios.

This medium conforms to ISO 21872 for the identification of *Vibrio* spp, including *Vibrio cholerae* from food, animal feeding stuffs and environmental samples in the area of food production and food handling.

TCBS Agar is also recommended for isolating *V. cholerae* and *V. parahaemolyticus* as well as other vibrios from stool specimens.

TYPICAL FORMULA

	(g/l)
Peptone	10.0
Yeast Extract	5.0
Sodium Citrate	10.0
Sodium Thiosulfate	10.0
Iron(III) Citrate	1.0
Sodium Chloride	10.0
Dried Bovine Bile	8.0
Sucrose	20.0
Bromothymol Blue	0.04
Thymol Blue	0.04
Agar	15.0
Final pH 8.6 ± 0.2 at 25°C	

METHOD PRINCIPLE

Peptone and yeast extract provide amino acids, nitrogen, carbon, vitamins and minerals for organisms growth. Sodium citrate serves to maintain an alkaline pH and along with sodium thiosulfate and oxbile are selective agents, inhibiting Gram-positive organisms and suppressing coliforms. Moreover, the alkaline pH of the medium enhances the recovery of *V. cholerae*. Sodium thiosulfate serves also as a sulfur source and, in combination with ferric citrate, detects hydrogen sulfide production. Sodium chloride maintains the osmotic balance of the medium and stimulates vibrios growth. Sucrose is the fermentable carbohydrate. Bromothymol blue and thymol blue are pH indicators. Agar is the solidifying agent.

PREPARATION

<u>Dehydrated medium</u>	Suspend 89.1 g of the powder in 1 liter of distilled or deionized water. Mix well. Heat to boil shaking frequently until completely dissolved. DO NOT AUTOCLAVE.
<u>Medium in bottles</u>	Melt the content of the bottle in a water bath at 100°C (loosing the cap partially removed) until completely dissolved. Then screw the cap and check the homogeneity of the dissolved medium, if it is the case turning the bottle upside down. Cool at 45-50°C, mix well avoiding foam formation and aseptically distribute into final containers.

TEST PROCEDURE

ISO 21872 recommends to follow two successive enrichment steps in Alkaline Saline Peptone Water (ASPW) before inoculating TCBS Agar.

TCBS Agar can be also inoculated directly with specimens such as rectal swabs, feces, vomitus or with food samples (*).

Incubate (protected from light) at 37 ± 1°C for 18-24 hours in an aerobic atmosphere.

(*) NB. Heavy inoculation is recommended. Swabs containing specimen material should be transported to the laboratory in Cary Blair Transport Medium (ref. 470290) if a delay in reaching the laboratory is anticipated. Specimens for cultivation of vibrios should not be frozen.

INTERPRETING RESULTS

Strains of *Vibrio cholerae* produce yellow colonies on TCBS Agar because of fermentation of sucrose. *Vibrio alginolyticus* also produce yellow colonies. It's possible that a few sucrose-positive *Proteus* strains can grow to form yellow, vibrid-like colonies. *Vibrio parahaemolyticus* is a sucrose non-fermenting organism and produces blue-green colonies, as does *Vibrio vulnificus*. Occasional isolates of *Pseudomonas* and *Aeromonas* species also produce blue-green colonies, but overall TCBS Agar is highly selective and any H₂S-negative colony is possibly *Vibrio* species.

APPEARANCE

Dehydrated medium: free-flowing, homogeneous, light beige to green beige.

Prepared medium: clear to slightly opalescent, green.

STORAGE

The powder is very hygroscopic, store the powder at 10-30°C, in a dry environment, in its original container tightly closed. Store bottles, tubes and prepared plates at 10-25°C away from light. Do not use the product beyond its expiry date on the label or if product shows any evidence of contamination or any sign of deterioration.

SHELF LIFE

Dehydrated medium: 4 years.
 Medium in bottles: 2 years.
 Medium in tubes: 1 year.
 Ready-to-use plates: 6 months.

QUALITY CONTROL

The medium is inoculated with the microbial strains indicated in the QC table.
 Inoculum for productivity: ≤ 100 CFU.
 Inoculum for selectivity: $> 10^3$ CFU.
 Incubation conditions: aerobically at $37 \pm 1^\circ\text{C}$ for 18-24 hours.

QC Table.

Microorganism		Growth	Colony Color
<i>Vibrio parahaemolyticus</i>	WDCM 00185	Good	Green
<i>Vibrio furnissii</i>	WDCM 00186	Good	Yellow
<i>Escherichia coli</i>	WDCM 00012	Inhibited	---

WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is intended for *in vitro* diagnostic use and must be used only by properly trained operators.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAPHY

1. EN ISO 11133:2014. Microbiology of food, animal feed and water – Preparation, production, storage and performance testing of culture media.
2. ISO 21872-1:2007. Microbiology of food and animal feeding stuffs – Horizontal method for the detection of potentially enteropathogenic *Vibrio* spp. – Part 1: Detection of *Vibrio parahaemolyticus* and *Vibrio cholerae*.
3. ISO 21872-2:2007. Microbiology of food and animal feeding stuffs – Horizontal method for the detection of potentially enteropathogenic *Vibrio* spp. – Part 2: Detection of species other than *Vibrio parahaemolyticus* and *Vibrio cholerae*.
4. American Public Health Association (1992): compendium of methods for the microbiological examination of foods, 3rd edition.
5. Dewitt, W.E., E.J. Gangarosa, I. Huq, and A. Zarifi (1971) Holding media for the transport of *Vibrio cholerae* from field to laboratory. Am. J. Trop. Med. Hyg. 20:685-688.
6. Kobayashi, T., S. Enomoto, R. Sakazaki, and S. Kuwahara (1963) A new selective medium for pathogenic vibrios: T.C.B.S. Agar (Modified Nakanishi's Agar). Jap. J. Bacteriol. 18:387-391.

PRESENTATION		Contents	Ref.
TCBS Agar	90 mm ready-to-use plates	20 plates	11195
TCBS Agar	90 mm ready-to-use plates	100 plates	11195*
TCBS Agar	Slant tubes	10 x 9 ml tubes	30022
TCBS Agar	Slant tubes	20 x 9 ml tubes	31022
TCBS Agar	Bottles	6 x 100 ml bottles	403140
TCBS Agar	Dehydrated medium	500 g of powder	611010
TCBS Agar	Dehydrated medium	100 g of powder	621010

TABLE OF SYMBOLS

LOT Batch code	IVD <i>In vitro</i> Diagnostic Medical Device	 Manufacturer	 Use by	 Fragile, handle with care
REF Catalogue number	 Temperature limitation	 Contains sufficient for $<n>$ tests	 Caution, consult Instruction For Use	 Do not reuse

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EE Broth-Mossel

Liquid medium for the cultivation and selective enrichment of Enterobacteriaceae from different types of samples, according to USP/EP/JP.

DESCRIPTION

Enterobacteriaceae Enrichment Broth-Mossel is a selective medium used for the detection of bile-tolerant Gram-negative bacteria in food and other materials of sanitary importance.

This medium complies with the recommendations of the harmonized method in the United States Pharmacopoeia (USP), European Pharmacopoeia (EP) and Japanese Pharmacopoeia (JP) for the microbiological examination of nonsterile products.

TYPICAL FORMULA	(g/l)
Pancreatic Digest of Gelatin	10.0
Glucose Monohydrate	5.0
Dehydrated Ox Bile	20.0
Potassium Dihydrogen Phosphate	2.0
Disodium Hydrogen Phosphate, Anhydrous	6.4*
Brilliant Green	0.015
Final pH 7.2 ± 0.2 at 25°C	

* Equivalent to 8.0 g of Disodium Hydrogen Phosphate Dihydrate.

METHOD PRINCIPLE

Pancreatic digest of gelatin provides amino acids, nitrogen, carbon, vitamins and minerals for organisms growth. Glucose is the fermentable carbohydrate. Ox bile and brilliant green are selective agents effective against Gram-positive cocci. Potassium phosphate and sodium phosphate act as buffer.

PREPARATION

Dehydrated medium Suspend 43.4 g of the powder in 1 liter of distilled or deionized water. Mix well. Heat to boil shaking frequently until completely dissolved. DO NOT AUTOCLAVE.

TEST PROCEDURE

As in the Pharmacopoeia, prepare the sample using a 1 in 10 dilution of not less than 1 g of the product to be examined by choosing Tryptic Soy Broth (ref. 24513 or 452080) as diluent and incubate at 20-25°C for 2-5 hour to resuscitate bacteria.

For qualitative test (test for absence), transfer the volume of the pre-enrichment broth corresponding to 1 g of the product to be examined to EE Broth-Mossel.

For quantitative test, enumerate Enterobacteriaceae found per milliliter or per gram of test sample by using the Most Probable Number (MPN) technique. Use the volume of the pre-enrichment broth containing 0.1 g, 0.01 g and 0.001 g (or 0.1 ml, 0.01 ml and 0.001 ml) of the product to be examined to inoculate EE Broth-Mossel.

For both types of test incubate EE Broth-Mossel at 30-35°C for 24-48 h and continue analysis by subculturing on Violet Red Bile Glucose Agar (ref. 11184). Incubate plates aerobically at 30-35°C for 18-24 hours.

INTERPRETING RESULTS

Turbidity of EE Broth-Mossel indicates microbial growth; acid production causes a color change of the medium to yellow.

No growth of colonies on Violet Red Bile Glucose Agar is reported as absence of bile-tolerant Gram-negative bacteria.

Growth of colonies constitutes a positive result and the probable number of bacteria is determined from the table below.

MPN Table.

Results for each quantity of product			Probable number of bacteria per gram or per milliliter of product
0.1 g or 0.1 ml	0.01 g or 0.01 ml	0.001 g or 0.001 ml	
+	+	+	>10 ³
+	+	—	10 ³ - 10 ²
+	—	—	10 ² - 10
—	—	—	<10

APPEARANCE

Dehydrated medium: free-flowing, homogeneous, light beige to light green.

Prepared medium: clear, green.

STORAGE

The powder is very hygroscopic, store the powder at 10-30°C, in a dry environment, in its original container tightly closed. Store bottles and prepared plates at 10-25°C away from light. Do not use the product beyond its expiry date on the label or if product shows any evidence of contamination or any sign of deterioration.

SHELF LIFE

Dehydrated medium: 4 years.

Medium in tubes/bottles: 1 year.

QUALITY CONTROL

The medium is inoculated with the microbial strains indicated in the QC table.

Inoculum for productivity: ≤100 CFU.

Inoculum for selectivity: >100 CFU.

Incubation conditions: 18-24 h at 30-35°C (Pharmacopoeia growth promotion).

QC Table.

Microorganism		Specification
<i>Escherichia coli</i>	ATCC® 8739	Good growth
<i>Pseudomonas aeruginosa</i>	ATCC® 9027	Good growth
<i>Staphylococcus aureus</i>	ATCC® 6538	Inhibition

WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is intended for professional use only and must be used by properly trained operators.

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BIBLIOGRAPHY

1. EN ISO 11133:2014. Microbiology of food, animal feed and water – Preparation, production, storage and performance testing of culture media.
2. European Pharmacopoeia 6.5 (2009) 2.6.13. Microbiological examination of non-sterile products: Test for specified microorganisms.
3. United States Pharmacopoeia 32 NF 27 (2009) <62> Microbiological examination of non-sterile products: Test for specified microorganisms.
4. Japanese Pharmacopoeia 4.05 (2008) Microbiological examination of non-sterile products: Test for specified microorganisms.
5. ISO 21528-1:2004. Microbiology of food and animal feeding stuffs – Horizontal method for the detection and enumeration of Enterobacteriaceae – Detection and enumeration by MPN technique with pre-enrichment.
6. ISO 21528-2:2004. Microbiology of food and animal feeding stuffs – Horizontal method for the detection and enumeration of Enterobacteriaceae – Colony count method.
7. Davidson, Roth, and Gambrel-Lenarz (2004) In Wehr and Frank (ed.) Standard methods for the microbiological examination of dairy products, 17th ed. American Public Health Association, Washington, D.C.
8. Kornacki and Johnson (2001) In Downes and Ito (ed.) Compendium of methods for the microbiological examination of foods, 4th ed. American Public Health Association, Washington D.C.
9. Mossel, Vissar, and Cornellisen (1963) J. Appl. Bacteriol. 26:444.

PRESENTATION		Contents	Ref.
EE Broth-Mossel	Tubes	20 x 10 ml tubes	24096
EE Broth-Mossel	Bottles (screw cap)	6 x 100 ml bottles	402480
EE Broth-Mossel	Bottles (flip-off cap)	25 x 100 ml bottles	453080
EE Broth-Mossel	Bottles (perforable cap)	6 x 100 ml bottles	495000
EE Broth-Mossel	Dehydrated medium	500 g of powder	610017
EE Broth-Mossel	Dehydrated medium	100 g of powder	620017

TABLE OF SYMBOLS

LOT Batch code	 Keep away from sunlight	 Manufacturer	 Use by	 Fragile, handle with care
REF Catalogue number	 Temperature limitation	 Contains sufficient for <n> tests	 Caution, consult Instruction For Use	 Do not reuse

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BIBLIOGRAFIA

1. EN ISO 11133:2014. Microbiology of food, animal feed and water – Preparation, production, storage and performance testing of culture media.
2. European Pharmacopoeia 6.5 (2009) 2.6.13. Microbiological examination of non-sterile products: Test for specified microorganisms.
3. United States Pharmacopoeia 32 NF 27 (2009) <62> Microbiological examination of non-sterile products: Test for specified microorganisms.
4. Japanese Pharmacopoeia 4.05 (2008) Microbiological examination of non-sterile products: Test for specified microorganisms.
5. ISO 21528-1:2004. Microbiology of food and animal feeding stuffs – Horizontal method for the detection and enumeration of Enterobacteriaceae – Detection and enumeration by MPN technique with pre-enrichment.
6. ISO 21528-2:2004. Microbiology of food and animal feeding stuffs – Horizontal method for the detection and enumeration of Enterobacteriaceae – Colony count method.
7. Davidson, Roth, and Gambrel-Lenarz (2004) In Wehr and Frank (ed.) Standard methods for the microbiological examination of dairy products, 17th ed. American Public Health Association, Washington, D.C.
8. Kornacki and Johnson (2001) In Downes and Ito (ed.) Compendium of methods for the microbiological examination of foods, 4th ed. American Public Health Association, Washington D.C.
9. Mossel, Vissar, and Cornillisen (1963) J. Appl. Bacteriol. 26:444.

PRESENTAZIONE		Contenuto	Ref.
EE Broth-Mossel	Provette	Provette 20 x 10 ml	24096
EE Broth-Mossel	Flaconi (tappo a vite)	Flaconi 6 x100 ml	402480
EE Broth-Mossel	Flaconi (tappo flip-off)	Flaconi 25 x100 ml	453080
EE Broth-Mossel	Flaconi (tappo perforabile)	Flaconi 6 x 100 ml	495000
EE Broth-Mossel	Terreno disidratato	500 g di polvere	610017
EE Broth-Mossel	Terreno disidratato	100 g di polvere	620017

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liofilchem@liofilchem.net

MRS Agar

Medium for isolation and enumeration of mesophilic lactic acid bacteria, according to 15214.

TYPICAL FORMULA	(g/l)
Enzymatic Digest of Casein	10.0
Meat Extract	10.0
Yeast Extract	4.0
Triammonium Citrate	2.0
Sodium Acetate	5.0
Magnesium Sulfate Heptahydrate	0.2
Manganese Sulfate Tetrahydrate	0.05
Dipotassium Hydrogen Phosphate	2.0
Glucose	20.0
Agar	15.0
Final pH 5.7 ± 0.1 at 25°C	

DESCRIPTION

MRS Agar is a medium used with supplements for the cultivation of *Lactobacillus* spp from different types of material. It may also support the growth of *Pediococcus* and *Leuconostoc* species as well as other secondary bacteria.

The complete medium complies with the recommendations of ISO 15214 and APHA.

PRINCIPLE

Enzymatic digest of casein and meat extract provide amino acids, nitrogen, carbon, vitamins and minerals for organisms growth. Yeast extract is a source of vitamins, particularly of B-group. Ammonium citrate and sodium acetate are the selective agents effective against streptococci and moulds. The low pH is also inhibitory for most organisms other than lactobacilli. Magnesium and manganese sulfates are sources of ions and sulfate acting as growth stimulants. Dipotassium phosphate is the buffer. Glucose is the fermentable carbohydrate. Agar is the solidifying agent.

Supplementation with Tween 80 Supplement (ref. 80031) provides a mixture of oleic esters and fatty acids essential for the growth of lactic acid bacteria.

PREPARATION

Suspend 68.3 g of the powder in 1 liter of distilled or deionized water. Mix well. Heat to boil shaking frequently until completely dissolved. Add 1 ml of Tween 80 Supplement. Sterilize in autoclave at 121°C for 15 minutes.

Note. According to ISO 15214, 1.4 g of Sorbic Acid (dissolved in about 10 ml of a 1 mol/l solution of sodium hydroxide) can be added to 1 liter of sterilized medium if extensive yeast contamination is suspected.

TECHNIQUE

1. Use a suitable diluent such as Buffered Peptone Water (ref. 24099) to perform serial dilutions of the test sample in order to achieve a colony count of between 15 and 300 colonies per plate.
2. Inoculate each plate with 1 ml of sample suspension by pour plating. Overlays may be used if required.
3. Incubate at 30°C for 72 hours.

INTERPRETATION OF RESULTS

Count colonies on all plates containing 15-300 colonies. Report the count as CFU/ml of sample allowing for dilution factors.

It may be necessary in some cases and for some products to confirm the colonies by simple techniques such as Gram staining, or the test for catalase.

STORAGE

The powder is very hygroscopic, store the powder at 10-30°C, in a dry environment, in its original container tightly closed and use it before the expiry date on the label or until signs of deterioration or contamination are evident. Store prepared plates at 2-8°C away from light.

WARNING AND PRECAUTIONS

For professional use only. Operators must be trained and have certain experience in the laboratory methods. Please read the instructions carefully before using this product. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this document.

Consult the Safety Data Sheet (SDS) for information regarding hazards and safe handling practices.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to the national and local regulations in force.

REFERENCES

1. APHA (2015): Compendium of Methods for the Microbiological Examination of Foods. 5th edition. American Public Health Association, Washington, D.C.
2. EN ISO 11133:2014. Microbiology of food, animal feed and water – Preparation, production, storage and performance testing of culture media.
3. Schillinger U and Holzapfel WH (2012) Culture media for Lactic Acid Bacteria. In: Handbook of Culture Media for Food and Water Microbiology. (Corry JEL, Curtis GDW and Baird RM eds), pp 174-186. Royal Society of Chemistry, Cambridge, UK.
4. ISO 15214:1998. Microbiology of food and animal feeding stuffs – Horizontal method for the enumeration of mesophilic lactic acid bacteria – Colony count technique at 30°C.
5. De Man JD, Rogosa M, and Sharpe ME (1960): A Medium for the cultivation of Lactobacilli. J. Appl. Bact. 23: 130-135.

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10-30°C

PACKAGING

Ref.	Content	Packaging
610024	500 g	500 g of powder in plastic bottle
620024	100 g	100 g of powder in plastic bottle

pH OF THE MEDIUM

5.7 ± 0.1

USE

MRS Agar is a medium used with supplements for the cultivation of mesophilic lactic acid bacteria, according to 15214

TECHNIQUE

Refer to technical sheet of the product

APPEARANCE OF THE MEDIUM

Powder medium

Appearance: free-flowing, homogeneous

Colour: beige

Ready-to-use medium

Appearance: slightly opalescent

Colour: amber

SHELF LIFE

4 years

QUALITY CONTROL

- Control of general characteristics, label and print
- Microbiological control
Inoculum for productivity: 50-100 CFU
Inoculum for selectivity: 10^4 - 10^6 CFU
Incubation Conditions: 72 ± 3 h at $30 \pm 1^\circ\text{C}$, in microaerobiosis

Microorganism

<i>Lactobacillus sakei</i>	WDCM 00015
<i>Lactobacillus lactis</i>	WDCM 00016
<i>Escherichia coli</i>	WDCM 00012
<i>Bacillus cereus</i>	WDCM 00001

Growth

Good
Good
Inhibited
Inhibited

TABLE OF SYMBOLS

LOT Batch code	Consult instructions for use	Manufacturer	Use by
REF Catalogue number	Temperature limitation	Contains sufficient for <n> tests	Keep away from sunlight

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
**English - EN****INTENDED USE**

This *in vitro* diagnostic device is intended for the calibration of quantitative tests listed in the value sheet. This *in vitro* diagnostic device is for professional use only.

COMPOSITION

- Lyophilized product prepared from human serum spiked with chemical and biological additives.
- Sodium azide < 0.1 % (w/w)
- Concentrations for each analyte to test are lot-specific.

MATERIALS REQUIRED BUT NOT PROVIDED

- VitalScientific reagents listed in the value sheet.
- General Laboratory equipment.

TRACEABILITY

The traceability is indicated in the value sheet enclosed in the kit.

PRECAUTIONS FOR USE AND WARNINGS

- Consult Safety Data Sheet (SDS) for a proper handling.
- Each unit of human blood used in the manufacture of these products was tested and found to be negative/non-reactive for the presence of HbAg, HCV and HIV1/2. The methods used were FDA approved or CE compliant. Nevertheless, since the risk of infection cannot be fully excluded these products must be handled as potentially infectious. In case of exposure, follow the guidelines of the competent health authorities.
- Take precautions when handling broken glass vials as sharp edges can injure the user.

- These products contain sodium azide which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of these reagents always flush with copious amounts of water to prevent azide buildup.
- Take normal precautions and adhere to good laboratory practice.
- Use clean or single use laboratory equipment only to avoid contaminations.

WASTE MANAGEMENT

Disposal of all waste material should be in accordance with local, state and federal regulatory requirements (please refer to the Safety Data Sheet (SDS)).

PREPARATION

- Carefully open the vial, avoiding the loss of lyophilizate.
- Add in **exactly 3 mL of distilled or deionized water**.
- Carefully close the vial and dissolve the contents completely within 30 minutes by occasional gentle stirring avoiding the formation of foam.

PRODUCT DETERIORATION

- The product may present a slightly hazy appearance after reconstitution. This has no effect on the performances of the product. All presence of particles would indicate deterioration.
- Do not use the product if there is visible evidence of contamination or damage (e.g. particle matter after reconstitution).
- Damage to the container may impact on product performance. Do not use the product if there is physical evidence of deterioration (e.g. leakages).

STABILITY**Prior to reconstitution :**

- Stored at 2-8 °C and protected from light.
- Do not use after expiration date indicated on the vial label.

After reconstitution :

- These products should be immediately and tightly capped to prevent contamination and evaporation.

- Stability of the components:
Between 15 and 25 °C : 8 hours
Between 2 and 8 °C : 2 days
Between -25 and -15 °C : 4 weeks (when frozen once)

Exceptions:

- Stability of total bilirubin (when stored protected from light):
Between 15 and 25 °C : 6 hours
Between 2 and 8 °C : 1 day
Between -25 and -15 °C : 2 weeks (When frozen once)
- Stability of direct bilirubin (when stored protected from light):
Between 15 and 25 °C : 3 hours
Between 2 and 8 °C : 8 hours
Between -25 and -15 °C : 2 weeks (when frozen once)

Note: Store vials tightly capped and protected from light.

PROCEDURE

To use ELICAL 2, follow the procedure described in the instructions for use of the reagent used.

LIMITATIONS

Using ELICAL 2 has been validated with the VitalScientific systems (Analyzers and reagents listed in the value sheet). Using any other system should be validated by the laboratory.

VALUES

The concentrations are indicated in the value sheet enclosed in the kit. For users of Selectra instruments allowing data import for tests, calibrators and controls, use corresponding barcode, available on the value sheet.

DECLARATION OF SERIOUS INCIDENT

Please notify the manufacturer (through your distributor) and competent authority of the Member State of the European union in which the user and/or the patient is established, of any serious incident that has occurred in relation to the device. For other jurisdictions, the declaration of serious incident should be in accordance with local, state and federal regulatory requirements. By reporting a serious incident, you provide information that can contribute to the safety of *in vitro* medical devices.

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- Symbols used on our documentation are defined on ISO-15223-1 standard, except for some presented in the symbols glossary available on the VitalScientific Website. (Symbols glossary).

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Nutrient Agar ISO 16266

Medium for cultivating non-fastidious organisms and confirming *Pseudomonas aeruginosa*, according to ISO 16266.

DESCRIPTION

Nutrient Agar ISO 16266 is a medium used for the cultivation of non-fastidious organisms from water samples. This medium is formulated according to ISO 16266 for the detection and enumeration of *Pseudomonas aeruginosa* in water by the membrane filtration technique.

TYPICAL FORMULA*

	(g/litre)
Peptone	5.0
Meat Extract	1.0
Yeast Extract	2.0
Sodium Chloride	5.0
Agar	15.0
Final pH 7.4 ± 0.2 at 25°C	

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METHOD PRINCIPLE

Peptone and meat extract provide amino acids, nitrogen, carbon, vitamins and minerals for organisms growth. Yeast extract is a source of vitamins, particularly of B-group. Sodium chloride maintains the osmotic balance of the medium. Agar is the solidifying agent.

PREPARATION

<u>Dehydrated medium</u>	Suspend 28 g of the powder in 1 liter of distilled or deionized water. Mix well. Heat to boil shaking frequently until completely dissolved. Sterilize in autoclave at 121°C for 15 minutes.
<u>Medium in bottles</u>	Melt the content of the bottle in a water bath at 100°C (loosing the cap partially removed) until completely dissolved. Then screw the cap and check the homogeneity of the dissolved medium, if it is the case turning the bottle upside down. Cool at 45-50°C, mix well avoiding foam formation and aseptically distribute into Petri dishes.

TEST PROCEDURE

According to ISO 16266, transfer the membrane and presumptive *Pseudomonas aeruginosa* to the plate medium.

Incubate aerobically at 36 ± 2°C for 20-24 hours.

Alternatively, the medium can be inoculated by spread plating or direct streaking of the sample over the agar surface.

INTERPRETING RESULTS

Observe for colony growth. Confirm *P. aeruginosa* by performing the oxidase test (ref. 88029).

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Avoid quick temperature shifts of plated medium to prevent condensation.

SHELF LIFE

Dehydrated medium: 4 years.

Medium in bottles: 2 years.

Medium in slant tubes: 1 year.

Ready-to-use plates: 6 months.

QUALITY CONTROL

Appearance of Dehydrated Medium: Free-flowing, homogeneous, beige.

Appearance of Prepared Medium: Slightly opalescent, light amber.

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The product is available in the various configurations listed below. There may be additional product ref. numbers as well. For an updated listing of available products, visit [liofilchem.com](https://www.liofilchem.com)

Product	Format	Packaging	Ref.
Nutrient Agar ISO 16266	Plate 90 mm	20 plates	10044
Nutrient Agar ISO 16266	Slant tube	10 x 7 ml	30083
Nutrient Agar ISO 16266	Bottle	6 x 100 ml	402190
Nutrient Agar ISO 16266	Bottle	6 x 200 ml	412190
Nutrient Agar ISO 16266	Bottle	6 x 500 ml	470060
Nutrient Agar ISO 16266	Dehydrated media	100 g	620036
Nutrient Agar ISO 16266	Dehydrated media	500 g	610036
Nutrient Agar ISO 16266	Dehydrated media	5 kg	6100365

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References / Riferimenti / Referencias

1. ISO 16266:2008. Water Quality – Detection and enumeration of *Pseudomonas aeruginosa* – Method by membrane filtration.
2. Marshall, R.T. (1993) Standard methods for the microbiological examination of dairy products, 16th ed.

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liofilchem@liofilchem.com

Nutrient Broth

Liquid medium for the cultivation of nonfastidious microorganisms.

DESCRIPTION

Nutrient Broth is a liquid medium used for the cultivation of a wide variety of organisms from clinical specimens and other materials.

This medium can be enriched with other ingredients such as blood, serum, sugars, etc., for special purposes.

TYPICAL FORMULA	(g/l)
-----------------	-------

Beef Extract	1.0
--------------	-----

Peptone	5.0
---------	-----

Yeast Extract	2.0
---------------	-----

Sodium Chloride	5.0
-----------------	-----

Final pH 6.8 ± 0.2 at 25°C

METHOD PRINCIPLE

Beef extract and peptone provide amino acids, nitrogen, carbon, vitamins and minerals for organisms growth. Yeast extract is a source of vitamins, particularly of B-group. Sodium chloride maintains the osmotic balance of the medium.

PREPARATION

<u>Dehydrated medium</u>	Suspend 13 g of the powder in 1 liter of distilled or deionized water. Mix well. Heat to boil shaking frequently until completely dissolved. Sterilize in autoclave at 121°C for 15 minutes.
--------------------------	--

TEST PROCEDURE

Inoculate broth with test sample. Incubate at 35 ± 2°C for 18-24 hours or longer if necessary.

INTERPRETING RESULTS

Turbidity indicates microbial growth.

APPEARANCE

Dehydrated medium: free-flowing, homogeneous, white to light beige.

Prepared medium: clear to slightly opalescent, light amber.

STORAGE

The powder is very hygroscopic, store the powder at 10-30°C, in a dry environment, in its original container tightly closed. Store bottles and tubes at 10-25°C away from light. Do not use the product beyond its expiry date on the label or if product shows any evidence of contamination or any sign of deterioration.

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Inoculum for productivity: ≤ 100 CFU

Incubation conditions: aerobically at $35 \pm 2^\circ\text{C}$ for 18-24 hours.

QC Table.

Microorganism		Growth
<i>Escherichia coli</i>	ATCC® 25922	Good
<i>Staphylococcus aureus</i>	ATCC® 25923	Good

WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is intended for *In vitro* diagnostic use and must be used only by properly trained operators.

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BIBLIOGRAPHY

1. Association of Official Analytical Chemists (1995) Official methods of analysis of AOAC International, 16th ed.
2. Marshall, R.T. (ed.) (1993) Standard methods for the microbiological examination of dairy products, 16th ed.
3. American Public Health Association (1923) Standard methods of water analysis, 5th ed.

PRESENTATION

		Contents	Ref.
Nutrient Broth	Tubes	20 x 10 ml tubes	24103
Nutrient Broth	Tubes	50 x 5 ml tubes	27503
Nutrient Broth	Bottles	6 x 100 ml bottles	402000
Nutrient Broth	Bottles	6 x 500 ml bottles	470050
Nutrient Broth	Dehydrated medium	500 g of powder	610037
Nutrient Broth	Dehydrated medium	100 g of powder	620037
Nutrient Broth	Dehydrated medium	5 kg of powder	6100375

TABLE OF SYMBOLS

LOT Batch code	IVD <i>In vitro</i> Diagnostic Medical Device	 Manufacturer	 Use by	 Fragile, handle with care
REF Catalogue number	 Temperature limitation	 Contains sufficient for <n> tests	 Caution, consult Instruction For Use	 Do not reuse

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Disc tests for confirmation of carbapenemase-producing Enterobacteriaceae.

DESCRIPTION

Carbapenemases are β -lactamases that hydrolyze penicillins, in most cases cephalosporins and to various degree carbapenems. The most important carbapenemases are categorized as three types of enzymes:

1. *Klebsiella pneumoniae* carbapenemases (KPC) belonging to the Ambler class A of β -lactamases;
2. VIM, IMP and NDM belonging to the class B metallo β -lactamases (MBL);
3. Class D β -lactamases of OXA-48-like enzymes.

Since the 1990s, the increased consumption of carbapenems to treat infections caused by extended-spectrum β -lactamases (ESBL) producing bacteria has in turn contributed to the dramatic increase in carbapenem-resistant Enterobacteriaceae (CRE). Bacterial acquisition of carbapenemases is crucial to the emergence of CRE. Indeed, the corresponding genes are mostly plasmid-located and associated with various mobile genetic structures and may confer resistance to virtually all β -lactams. Infections with carbapenemase-producing Enterobacteriaceae are associated with high mortality posing a serious threat especially to hospitalized patients. Carbapenemases are considered to be of high epidemiological importance and one of the most important concerns is the spread in the community, particularly in *E. coli*. Adequate treatment and control of CRE infections is predicted upon their accurate and prompt diagnosis from patient samples in the clinical laboratory. Carbapenem resistance can also arise in both ESBL and AmpC-producing organisms (Ambler class A and C, respectively) by mutation that regulate the synthesis of porins contained in the outer membrane.

Carbapenem MICs for carbapenemase-producing enterobacteriaceae may be below the clinical breakpoints, however, the epidemiological cut-off (ECOFF) defined by EUCAST can be used to detect carbapenemase-producers. Meropenem, offering the best compromise between sensitivity and specificity, is recommended for screening.

Following detection of reduced susceptibility to carbapenems, a phenotypic method for confirmation of carbapenemase-production should be applied. Meropenem discs supplemented with different inhibitors of carbapenemase are used to evaluate carbapenem-resistant strains.

CONTENTS OF THE PACKAGES

5 x 50 discs cartridges, each packaged in a "blister" with a dryer.

METHOD PRINCIPLE

Enterobacteriaceae suspected to be producers of carbapenemases may be confirmed by using meropenem and evaluating the synergistic effects when combined with the following inhibitors:

- **Phenylboronic acid** inhibits class A carbapenemase;
- **EDTA** inhibits class B carbapenemase;
- **Cloxacillin**, which inhibits AmpC β -lactamases, permits to differentiate between AmpC hyperproduction plus porin loss and carbapenemase-production.

There is no currently available inhibitor for class D carbapenemases. However, when there is no synergy of meropenem with any of the above mentioned inhibitors, 30 μ g **Temocillin** disc (TMO) can be used in order to differentiate between ESBL plus porin loss and OXA-48 like enzymes.

For each combined disc test (CDT), discs containing meropenem alone and in combination with a carbapenemases inhibitor are applied. The inhibition zone around the meropenem disc combined with inhibitor is compared with the zone around the 10 μ g meropenem disc (MRP).

GATHERING AND KEEPING SAMPLES

The colonies that are to be subjected to the susceptibility test are taken up by culture media that have been previously swabbed with the sample under examination.

TEST PROCEDURE

1. Using a fresh, pure culture prepare a suspension of the test organism equal to 0.5 McFarland Standard.
2. Using a sterile cotton swab, spread the adjusted suspension over the entire area of a Mueller Hinton agar plate.
3. Apply Meropenem, Meropenem + EDTA, Meropenem + Phenylboronic acid, Meropenem + Cloxacillin, on the inoculated plate, ensuring sufficient space between individual discs to allow for proper measurement of inhibition zones.
4. Incubate at $36 \pm 1^\circ\text{C}$ for 18-24 hours.

EVALUATING THE RESULTS

At the end of the incubation period, measure the inhibition halos and interpret as indicated in the table below.

Interpretative Table. Synergy of meropenem with carbapenem inhibitors for confirmation of CRE.

Increase in inhibition zone of meropenem with the following inhibitor			Temocillin (TMO)*	β-lactamases
Phenylboronic acid (MR+BO)	EDTA (MR+ED)	Cloxacillin (MR+CL)		
≥ 4 mm	< 5 mm	< 5 mm	---	KPC
< 4 mm	≥ 5 mm	< 5 mm	---	MBL
≥ 4 mm AND	< 5 mm	≥ 5 mm	---	AmpC + porin loss or efflux
< 4 mm	< 5 mm	< 5 mm	< 11 mm	OXA-48-like

*Temocillin susceptibility test is recommended only in cases where no synergy is detected.

QUALITY CONTROL

Each production batch of discs is subjected to the quality control with the following bacterial strains:

Microorganism		
<i>Klebsiella pneumoniae</i>	ATCC® 700603	Negative control
<i>Klebsiella pneumoniae</i>	ATCC® BAA-2146	MBL positive
<i>Klebsiella pneumoniae</i>	ATCC® BAA-1705	KPC positive
<i>Klebsiella pneumoniae</i>	NCTC 13442	OXA-48 positive

LIMITS

Diffusion susceptibility tests use an *in vitro* technique and cannot therefore reproduce the extremely complex *in vivo* conditions. Nevertheless, it is a useful and important tool that helps the clinician choose the correct therapy. Many variable factors influence the final result of the diffusion susceptibility test. The main ones are: the culture medium used, impregnation of the discs, inoculation of the medium, temperature, time and incubation atmosphere of the plates, pre-incubation and pre-diffusion conditions, depth of the medium, etc.

PRECAUTIONS

The discs cannot be classified as being hazardous according to current legislation but fall within the specific field of application where a safety data sheet must be supplied because they can cause phenomena of sensitization in sensitive subjects if they come into contact with the skin.

The discs are disposable products. They are only for diagnostic *in vitro* use and are intended for professional use. They must be used in the laboratory by properly trained operators using approved aseptic and safety methods for pathogenic agents.

STORAGE

Store the unopened blister at -20°C to +8°C till the expiry date. Allow unopened cartridge to come to room temperature before removing it from the blister for minimising condensation on the discs. Leftover discs from an opened cartridge should be stored at 2-8°C for no more than 7 days. Return unused discs to the refrigerator as soon as the application of the discs has been completed. Dispose of expire discs.

ELIMINATING USED MATERIAL

After use, the discs and the material that comes into contact with the sample must be decontaminated and disposed of in accordance with current laboratory techniques for the decontamination and disposal of potentially infected material.

REFERENCES

- EUCAST guidelines for detection of resistance mechanisms and specific resistances of clinical and/or epidemiological importance. Version 1.0, 2013.

PRESENTATION

DESCRIPTION	PACKAGING		REF
Meropenem 10 µg	MRP	50 Discs	99007
Meropenem + EDTA	MR+ED	50 Discs	
Meropenem + Phenylboronic acid	MR+BO	50 Discs	
Meropenem + Cloxacillin	MR+CL	50 Discs	
Temocillin 30 µg	TMO	50 Discs	

TABLE OF SYMBOLS

LOT Batch code	IVD In Vitro Diagnostic Medical Device	 Manufacturer	 Use by	 Fragile, handle with care
REF Catalogue number	 Temperature limitation	 Contains sufficient for <n> tests	 Caution, consult accompanying documents	 Do not reuse

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Product Code: 3825-300

1.0 INTRODUCTION

Intended Use: The Quantitative Determination of Circulating Troponin-I concentrations in Human Serum by a Microplate Enzyme Immunoassay, Colorimetric

2.0 SUMMARY AND EXPLANATION OF THE TEST

The cardiac-specific isoform of Troponin I (cTnI) has been known as a marker of heart damage and myocardial cell death, due to myocardial infarction, for just over 10 years. Troponin-I (cTnI, 24 kDa) is the inhibitory subunit of the Troponin complex of striated muscle. Most of the Troponin (I & T) proteins are located within the contractile apparatus of the striated muscle. The concentration of these subunits is increased in circulation for many days after AMI, because release from the structural elements requires degradation of myofibril itself. This location and release from the specific myocardial tissue and their sustained appearance in circulation make Troponin subunits reliable bio-markers of AMI. Unlike CK-MB and Myoglobin, Troponin does not suffer from non-specificity issues. The immunological determination of cTnI does have some issues surrounding it but those are mainly due to lack of mass standardization, the presence of post-translationally modified cTnI in the circulation and variations in the antibody cross-reactivities to the various detectable forms of cTnI. However, cTnI, unlike Myoglobin, has not been found in circulation in marathon runners and other cases of skeletal injury or trauma. Careful selection of antibodies and diligent preparation of stable cTnI for calibration permits cTnI to be an excellent biochemical marker of myocardial damage.

Troponin I (cTnI) is the inhibitory subunit of the Troponin complex, which regulates the calcium modulated interaction of actin and myosin in striated muscle. The complex is a heterodimer consisting of troponins C, I and T, which are tightly bound to the contractile apparatus; hence, circulating concentrations are low. Even though Troponin C and T are both considered to be equally good markers for AMI, the NH₂ terminus of cTnI has 31 additional amino acids not present in the skeletal isoforms and has generated an interest in easily creating cTnI-specific MoAb's, by researchers. Troponin I, however, has its own problems that researchers are dealing with. Troponin I is not very stable. A major portion of it is bound with troponin C (cTnC) and the remaining small part is free in circulation. Again, the post translational modifications, including selective degradation, covalent complex formation and phosphorylation of cTnI in the post-ischemic myocardium complicate the matters further more. Also, the cTnI proteolysis is even more extensive in human myocardium. Some researchers attribute it in part to the heterogeneity of disease states present in a given patient population. The key to overcoming these problems is a very careful selection of antibodies, the matrix for stabilizing the cTnI and the cTnI protein used as the standard itself. Monobind has been able to successfully put together an assay that promises to have overcome those issues.

In this method, Troponin-I calibrator, patient specimen or control is first added to a streptavidin coated well. Biotinylated monoclonal and enzyme labeled antibodies (directed against distinct and different epitopes of *Troponin-I*) are added and the reactants mixed. Reaction between the various *Troponin-I* antibodies and native *Troponin-I* forms a sandwich complex that binds with the streptavidin coated to the well.


After the completion of the required incubation period, the enzyme- Troponin-I antibody bound conjugate is separated from the unbound enzyme-Troponin-I conjugate by aspiration or decantation. The activity of the enzyme present on the surface of the well is quantitated by reaction with a suitable substrate to produce color. The employment of several serum references of known Troponin-I levels permits the construction of a dose response curve of activity and concentration. From comparison to the dose response curve, an unknown specimen's activity can be correlated with Troponin-I concentration.

3.0 PRINCIPLE

Sandwich Equilibrium Method (TYPE 2):

The essential reagents required for an immunoassay include high affinity and specificity antibodies (enzyme and immobilized), with different and distinct epitope recognition, **in excess**, and native antigen. In this procedure, the immobilization takes place during the assay at the surface of a microplate well through the interaction of x-cTnI antibody coated on the well.

Upon mixing the enzyme-labeled antibody and a serum containing the native antigen, a reaction results between the native antigen and the antibodies, without competition or steric hindrance, to form a sandwich complex. The interaction is illustrated by the following equation:

Ab_(well) = Antibody coated on well (Excess Quantity)

Ag_{cTnI} = Native Antigen (Variable Quantity)

Enz^{Ab} = Enzyme labeled Antibody (Excess Quantity)

Enz^{Ab} - Ag_{cTnI} - Ab_(well) = Antigen-Antibodies Sandwich Complex

k_a = Rate Constant of Association

k_a = Rate Constant of Dissociation

After sufficient time results, the antibody-bound fraction is separated from unbound antigen by decantation or aspiration. The enzyme activity in the antibody-bound fraction is directly proportional to the native antigen concentration. By utilizing several different serum references of known antigen values, a dose response curve can be generated from which the antigen concentration of an unknown can be ascertained.

4.0 REAGENTS

Materials Provided:

A. cTnI Calibrators – 1.0 ml/vial (Lyophilized) – Icons [A-F]

Six (6) vials of references for *Troponin-I* antigen at levels of 0(A), 0.4(B), 1.25(C), 2.5(D), 7.5(E), and 20(F) ng/ml. A preservative has been added. **Reconstitute each vial with 1.0ml of distilled or deionized water.**

The reconstituted calibrators are stable for 24 hours at 2-8°C. In order to store for a longer period of time aliquot the reconstituted calibrators in cryo vials and store at -20°C. **DO NOT FREEZE THAW MORE THAN ONCE.**

Note: The calibrators, human serum based, were calibrated using NIST standard for cTnI # 2921.

B. cTnI Enzyme Reagent – 13 ml/vial – Icon

One (1) vial containing enzyme labeled affinity purified antibody labeled monoclonal mouse IgG in buffer, dye, and preservative. Store at 2-8°C.

C. cTnI Antibody Coated Plate – 96 wells – Icon

One 96-well microplate coated with x-cTnI antibody, packaged in an aluminum bag with a drying agent. Store at 2-8°C.

D. Wash Solution Concentrate – 20 ml/vial – Icon

One (1) vial containing a surfactant in buffered saline. A preservative has been added. Store at 2-8°C.

E. Substrate A – 7.0ml/vial – Icon S^A

One (1) vial containing tetramethylbenzidine (TMB) in buffer. Store at 2-8°C.

F. Substrate B – 7.0ml/vial – Icon S^B

One (1) vial containing hydrogen peroxide (H₂O₂) in buffer. Store at 2-8°C.

G. Stop Solution – 8.0ml/vial – Icon

One (1) vial containing a strong acid (1N HCl). Store at 2-8°C.

H. Product Instructions.

Note 1: Do not use reagents beyond the kit expiration date.

Note 2: **Opened reagents are stable for sixty (60) days when stored at 2-8°C. Kit and component stability are identified on the label.**

Note 3: Above reagents are for a single 96-well microplate.

4.1 Required But Not Provided:

- Pipette(s) capable of delivering 25µl and 100µl volumes with a precision of better than 1.5%.
- Dispenser(s) for repetitive deliveries of 0.100ml and 0.350ml volumes with a precision of better than 1.5% (optional).
- Microplate washer or a squeeze bottle (optional).
- Microplate Reader with 450nm and 620nm wavelength absorbance capability.
- Absorbent Paper for blotting the microplate wells.
- Plastic wrap or microplate cover for incubation steps.
- Vacuum aspirator (optional) for wash steps.
- Timer.
- Storage container for storage of wash buffer.
- Distilled or deionized water.
- Quality Control Materials.

5.0 PRECAUTIONS

For In Vitro Diagnostic Use

Not for Internal or External Use in Humans or Animals

All products that contain human serum have been found to be non-reactive for Hepatitis B Surface antigen, HIV 1&2 and HCV antibodies by FDA licensed reagents. Since no known test can offer complete assurance that infectious agents are absent, all human serum products should be handled as potentially hazardous and capable of transmitting disease. Good laboratory procedures for handling blood products can be found in the Center for Disease Control / National Institute of Health, "Biosafety in Microbiological and Biomedical Laboratories," 2nd Edition, 1988, HHS.

Safe Disposal of kit components must be according to local regulatory and statutory requirement.

6.0 SPECIMEN COLLECTION AND PREPARATION

The specimens shall be blood, serum in type, and the usual precautions in the collection of venipuncture samples should be observed. The blood should be collected in a plain redtop venipuncture tube without additives or gel barrier. Allow the blood to clot. Centrifuge the specimen to separate the serum from the cells.

Samples may be refrigerated at 2-8°C for a maximum period of five (5) days. If the specimen(s) cannot be assayed within this time, the sample(s) may be stored at temperatures of -20°C for up to 30 days. Avoid use of contaminated devices. Avoid repetitive freezing and thawing. When assayed in duplicate, 0.050 ml (50µl) of the specimen is required.

7.0 QUALITY CONTROL

Each laboratory should assay controls at levels in the low, normal and elevated ranges for monitoring assay performance. These controls should be treated as unknowns and values determined in every test procedure performed. Quality control charts should be maintained to follow the performance of the supplied reagents. Pertinent statistical methods should be employed to ascertain trends. Significant deviation from established performance can indicate unnoticed change in experimental conditions or degradation of kit reagents. Fresh reagents should be used to determine the reason for the variations.

8.0 REAGENT PREPARATION

1. Wash Buffer

Dilute contents of wash solution to 1000ml with distilled or deionized water in a suitable storage container. Store at 2-30°C for up to 60 days.

2. Working Substrate Solution – Stable for 1 year

Pour the contents of the amber vial labeled Solution 'A' into the clear vial labeled Solution 'B'. Place the yellow cap on the clear vial for easy identification. Mix and label accordingly. Store at 2 - 8°C.

Note 1: **Do not use the working substrate if it looks blue.**

Note 2: **Do not use reagents that are contaminated or have bacteria growth.**

9.0 TEST PROCEDURE

Before proceeding with the assay, bring all reagents, serum reference calibrators and controls to room temperature (20-27°C).

*****Test Procedure should be performed by a skilled individual or trained professional*****

- Format the microplates' wells for calibrator, control and patient specimen to be assayed in duplicate. **Replace any unused microwell strips back into the aluminum bag, seal and store at 2-8°C.**
- Pipette 0.025 ml (25µl) of the appropriate calibrators, controls and samples into the assigned wells
- Add 0.100 ml (100µl) of the Troponin-I Enzyme Reagent to each well. **It is very important to dispense all reagents close to the bottom of the microwell.**
Note: **Use a multichannel pipet to quickly dispense the Enzyme Reagent to avoid drift if the dispensing is to take more than a few minutes.**
- Swirl the microplate gently for 20-30 seconds to mix. Cover with a plastic wrap.
- Incubate for 15 minutes at room temperature (20-27°C).
- Discard the contents of the microplate by decantation or aspiration. If decanting, tap and blot the plate dry with absorbent paper.
- Add 0.350 (350µl) of wash buffer (see Reagent Preparation Section), decant (tap and blot) or aspirate. Repeat four (4) additional times for a total of five (5) washes. **An automatic or manual plate washer can be used. Follow the manufacturer's instruction for proper usage. If a squeeze bottle is used, fill each well to the top by squeezing the container. Avoiding air bubbles. Decant the wash and repeat four (4) additional times.**
- Add 0.100 ml (100µl) of working substrate solution to all wells (see Reagent Preparation Section).
DO NOT SHAKE THE PLATE AFTER SUBSTRATE ADDITION
- Incubate at room temperature for fifteen (15) minutes.
- Add 0.050ml (50µl) of stop solution to each well and mix gently for 15-20 seconds. Read the absorbance in each well at 450nm (using a reference wavelength of 620-630nm to minimize well imperfections) in a microplate reader. **The results should be read within thirty (30) minutes of adding the stop solution.**

NOTE: Always add reagents in the same order to minimize reaction time differences between wells.

10.0 CALCULATION OF RESULTS

A dose response curve is used to ascertain the concentration of Troponin I in unknown specimens.

- Record the absorbance obtained from the printout of the microplate reader as outlined in Example 1
- Plot the absorbance for each duplicate serum reference versus the corresponding cTnI concentration in ng/ml on linear graph paper (do not average the duplicates of the serum references before plotting).
- Draw the best-fit curve through the plotted points.
- To determine the concentration of cTnI for an unknown, locate the average absorbance of the duplicates for each unknown on the vertical axis of the graph, find the intersecting point on the curve, and read the concentration (in ng/ml) from the horizontal axis of the graph (the duplicates of the unknown may be averaged as indicated). In the following example, the

average absorbance (0.322) intersects the dose response curve at 3.61ng/ml cTnI concentration (See Figure 1).


Note: Computer data reduction software designed for ELISA assays may also be used for the data reduction. If such software is utilized, the validation of the software should be ascertained.

EXAMPLE 1

Sample I.D.	Well Number	Abs (A)	Mean Abs (B)	Value (ng/ml)
Cal A	A1	0.007	0.007	0
	B1	0.008		
Cal B	C1	0.035	0.037	0.4
	D1	0.040		
Cal C	E1	0.112	0.116	1.25
	F1	0.119		
Cal D	G1	0.257	0.260	2.5
	H1	0.264		
Cal E	A2	0.777	0.777	7.5
	B2	0.776		
Cal F	C2	2.362	2.365	20
	D2	2.368		
Ctrl 1	E2	0.050	0.051	0.43
	F2	0.051		
Ctrl 2	G2	1.256	1.245	12.75
	H2	1.233		
Patient	A3	0.320	0.322	3.61
	B3	0.324		

*The data presented in Example 1 and Figure 1 is for illustration only and **should not** be used in lieu of a dose response curve prepared with each assay.

Figure 1

11.0 Q.C. PARAMETERS

In order for the assay results to be considered valid the following criteria should be met:

- The absorbance (OD) of calibrator 'A' should be ≤ 0.07 .
- The absorbance (OD) of calibrators 'F' should be ≥ 1.3 .
- Four out of six quality control pools should be within the established ranges.

12.0 RISK ANALYSIS

The MSDS and Risk Analysis Form for this product is available on request from Monobind Inc

12.1 Assay Performance

- It is important that the time of reaction in each well is held constant to achieve reproducible results.
- Pipetting of samples should not extend beyond ten (10) minutes to avoid assay drift.
- Highly lipemic, hemolyzed or grossly contaminated specimen(s) should not be used.
- If more than one (1) plate is used, it is recommended to repeat the dose response curve.
- The addition of substrate solution initiates a kinetic reaction, which is terminated by the addition of the stop solution. Therefore, the substrate and stop solution should be added in the same sequence to eliminate any time-deviation during reaction.
- Plate readers measure vertically. Do not touch the bottom of the wells.
- Failure to remove adhering solution adequately in the aspiration or decantation wash step(s) may result in poor replication and spurious results.

- Use components from the same lot. No intermixing of reagents from different batches.
- Patient samples with Troponin-I concentrations above 30 ng/ml may be diluted with the zero calibrator or Troponin-I free pooled human serum or urine and re-assayed. Multiply the value obtained by the dilution factor to obtain the corrected value.
- Accurate and precise pipetting, as well as following the exact time and temperature requirements prescribed are essential. Any deviation from Monobind's IFU may yield inaccurate results.
- All applicable national standards, regulations and laws, including, but not limited to, good laboratory procedures, must be strictly followed to ensure compliance and proper device usage.
- It is important to calibrate all the equipment e.g. Pipettes, Readers, Washers and/or the automated instruments used with this device, and to perform routine preventative maintenance.
- Risk Analysis- as required by CE Mark IVD Directive 98/79/EC - for this and other devices, made by Monobind, can be requested via email from Monobind@monobind.com.

12.2 Interpretation

Measurements and interpretation of results must be performed by a skilled individual or trained professional.

- Laboratory results alone are only one aspect for determining patient care and should not be the sole basis for therapy, particularly if the results conflict with other determinants.
- For valid test results, adequate controls and other parameters must be within the listed ranges and assay requirements.
- If test kits are altered, such as by mixing parts of different kits, which could produce false test results, or if results are incorrectly interpreted, Monobind shall have no liability.
- If computer controlled data reduction is used to interpret the results of the test, it is imperative that the predicted values for the calibrators fall within 10% of the assigned concentrations.
- The reagents for the procedure have been formulated to eliminate maximal interference; however, potential interaction between rare serum specimens and test reagents can cause erroneous results. Heterophilic antibodies often cause these interactions and have been known to be problems for all kinds of immunoassays. (Boscato LM Stuart MC: 'Heterophilic antibodies: a problem for all immunoassays' Clin. Chem 1988:3427-33). For diagnostic purposes the results from this assay should be used in combination with clinical examination, patient's history and, all other clinical findings.

13.0 EXPECTED VALUES

Troponin-I values are different in plasma and serum. In addition, plasma sample may be influenced by the additives used. For example, heparin minimally affects the results but oxalate and EDTA have significant effect. A serum sample that has been quickly separated from red cells is preferred.

Based on the clinical data gathered by Monobind in concordance with the published literature the following ranges have been assigned. **These ranges should be used as guidelines only:**

Adult (Normal)	≤ 0.4 ng/ml
----------------	------------------

It is important to keep in mind that establishment of a range of values which can be expected to be found by a given method for a population of "normal"-persons is dependent upon a multiplicity of factors: the specificity of the method, the population tested and the precision of the method in the hands of the analyst. For these reasons each laboratory should depend upon the range of expected values established by the Manufacturer only until an in-house range can be determined by the analysts using the method with a population indigenous to the area in which the laboratory is located.

14.0 PERFORMANCE CHARACTERISTICS

14.1 Precision

The within and between assay precision of the cTnI AccuBind® ELISA Test System were determined by analyses on three different levels of pool control sera. The number, mean value, standard deviation and coefficient of variation for each of these control sera are presented in Table 2 and Table 3.

TABLE 2

Within Assay Precision (Values in ng/ml)				
SAMPLE	N	X	σ	C.V.
Pool 1	20	0.44	0.014	3.3%
Pool 2	20	3.55	0.072	2.0%
Pool 3	20	12.75	0.311	2.4%

TABLE 3

Between Assay Precision* (Values in ng/ml)				
SAMPLE	N	X	σ	C.V.
Pool 1	10	0.48	0.038	7.9%
Pool 2	10	3.68	0.242	6.6%
Pool 3	10	13.56	0.745	5.5%

*As measured in ten experiments in duplicate over ten days.

14.2 Sensitivity

The sensitivity (detection limit) was ascertained by determining the LoB and LoD. The LoB is 0.0639 ng/ml and the LoD is 0.0759 ng/ml.

14.3 Accuracy

The cTnI AccuBind® ELISA Test System was compared with a predicate radioimmunoassay assay. Biological specimens from population (symptomatic and asymptomatic) were used. (The values ranged from N/D – 18.1 ng/ml.) The total number of such specimens was 151. The data obtained is displayed in Table 4.

TABLE 4

Method	Mean (x)	Linear Square Regression Analysis	Correlation Coefficient
Monobind (y)	3.04	$y = 0.3500 + 0.9266(x)$	0.950
Reference(x)	2.92		

Only slight amounts of bias between this method and the reference method are indicated by the closeness of the mean values. The least square regression equation and correlation coefficient indicates excellent method agreement.

14.4 Specificity

The cross-reactivity of the Troponin-I ELISA method to selected substances was evaluated by adding the interfering compounds to a serum matrix at the very high concentration(s).

SUBSTANCE	Cross Reactivity
Hemoglobin	ND
CK-MB	ND
TnT	ND
FABP	ND

The presence of lipemia (25 mg/ml); hemoglobin (4.0 mg/ml) and bilirubin (2.5 mg/ml) did not affect the assay precision.

14.5 Hook Effect

Human serum samples spiked with concentrations up to 10,000 ng/ml of cTnI did not show any hook effect with cTnI AccuBind® ELISA Test System.

15.0 REFERENCES

- Apple Fred S, Christenson RH, Valdes RJ, Andriak AB, Duh Show-Hong, Feng YJ, Saeed AJ, Johnson Nancy J, Koplen Brenda, Mascotti K, and Wu Alan J, "Simultaneous Rapid Measurement of Whole Blood Myoglobin, Creatine Kinase MB, and Cardiac Troponin I by the Triage Cardiac Panel for Detection of Myocardial Infarction", Clin Chem, 45, 199-205 (1999).
- Adams JE, Schechtman KB, Landt Y, et al, "Comparable detection of acute Myocardial infarction by creatine kinase MB isoenzyme and cardiac troponin I", Clin Chem, 40,1291-5 (1994).
- Apple FS, Preese LM, "Creatine kinase-MB: detection of Myocardial infarction and monitoring reperfusion", J Clin Immunoassay, 17, 24-9 (1994).
- Bhayana V, Cohoe S, Leung FY, et al, "Diagnostic evaluation of creatine kinase-2 mass and creatine kinase-3 and -2 isoform ratios in early diagnosis of acute Myocardial infarction", Clin Chem, 39, 488-95 (1993).
- Bruns DE, "Diagnosis of acute Myocardial infarction when skeletal muscle damage is present: a caveat regarding use of creatine kinase isoenzymes", Clin Chem, 35, 705 (1989).
- Gibler WB, Lewis LM, Erb RE, et al, "Early detection of acute Myocardial infarction in patients presenting with chest pain and

- nondiagnostic ECGs; serial Troponin-I sampling in the emergency department", Ann Emerg Med, 19, 1359-66 (1990).
- Henderson AR, Stark JA, McQueen MJ, et al, "Is determination of creatine kinase-2 after electrophoretic separation accurate?", Clin Chem, 40, 177-83 (1994).
- Kallner A, Sylven C, Brodin U, et al, "Early diagnosis of acute Myocardial infarction; a comparison between chemical predictors", Scand J Clin Lab Invest, 49, 633-9 (1989).
- Kiyasu Y, "Current status of detecting Troponin-I for patient management", Am Clin Prod Rev, 4, 29-31 (1985).
- Lang H, Wuerzburg U, "Creatine kinase, an enzyme of many forms", Clin Chem, 28, 1439-47 (1982).
- Lee KN, Csako G, Bernhardt P, Elin RJ, "Relevance of macro creatine kinase type 1 and type 2 isoenzymes to laboratory and clinical data", Clin Chem, 40, 1278-83 (1994).
- Lee TH, Rouan GW, Weisberg MC, et al, "Sensitivity of routine clinical criteria for diagnosing Myocardial infarction within 24 hours of hospitalization", Ann Intern Med, 106, 181-6 (1987).
- Panteghini M, "Creatine kinase MB isoforms", J Clin Immunoassay, 17, 30-4 (1994).

Effective Date: 2019-Jan-01 Rev. 4
MP3825

DCO: 1320
Product Code: 3825-300

Size	96(A)		192(B)	
	1ml set		1ml set	
Reagent (fill)	A)	1 (13ml)	B)	2 (13ml)
	C)	1 plate	D)	2 plates
	D)	1 (20ml)	E)	1 (20ml)
	E)	1 (7ml)	F)	2 (7ml)
	F)	1 (7ml)	G)	2 (7ml)
	G)	1 (8ml)	H)	2 (8ml)

For Orders and Inquires, please contact

Monobind Inc.
100 North Pointe Drive
Lake Forest, CA 92630 USA

Tel: +1 949.951.2665 Mail: info@monobind.com
Fax: +1 949.951.3539 Fax: www.monobind.com

Please visit our website to learn more about our products and services.

Glossary of Symbols (EN 980/ISO 15223)