



# MICRONAUT

## ● MIC-Strip Vancomycin/Teicoplanin

### What can it be used for?

The MICRONAUT MIC Strip Vancomycin/Teicoplanin enables antimicrobial susceptibility testing against the glycopeptide antibiotics vancomycin and teicoplanin by using the broth microdilution (BMD) method as advised by EUCAST for *Enterococcus faecalis* and *Enterococcus faecium* positive for vanA / vanB [1].

Antimicrobial susceptibility testing of vancomycin and teicoplanin has been fraught with difficulties until now. EUCAST recently issued recommendations confirming that BMD is currently the only valid method for MIC determination of vancomycin.

Gradient test methods for MIC determination of vancomycin are currently not recommended because of false susceptibility results for vanB-expressing *Enterococcus* spp. isolates due to underestimated MIC values.

#### Vancomycin concentration range

0.25 to 4 mg/L

#### Teicoplanin concentration range

0.25 to 8 mg/L

Manufactured by  


## Features and benefits

- Phenotypic detection of glycopeptide resistance (e.g. vanA, vanB)
- Standardized MIC determination of vancomycin and teicoplanin in routine laboratory diagnostics according to EUCAST recommendations
- The strip format permits economic use of the AST, adapted to the individual laboratory diagnostic requirements
- Easy handling combined with standard laboratory equipment, for high quality MIC results of vancomycin and teicoplanin
- No special expertise required, easy visual evaluation

## Procedure

- Prepare a 0.5 McFarland bacterial suspension in NaCl
- Transfer an aliquot into Mueller-Hinton II broth (CAMHB)
- Inoculate the appropriate number of MIC-Strips
- Incubate for 18-22 hours at 35-37°C
- Read the results visually

## Shelf life and storage


- Shelf life: 24 months from date of production (unopened), expires 8 weeks after opening
- Storage: at room temperature (15-25°C)

## International ISO standard method

- MICRONAUT is a broth microdilution system using the international reference methodology (ISO 20776-1). MIC levels are based on EUCAST guidelines where applicable.

### Reference

[1] EUCAST guidelines for detection of resistance mechanisms and specific resistances of clinical and/or epidemiological importance Version 2.0 July 2017.



### Vancomycin (VAN) and teicoplanin (TPL) concentrations (mg/L)

| 1  | 2    | 3   | 4   | 5   | 6   | 7    | 8   | 9   | 10  | 11  | 12  |
|----|------|-----|-----|-----|-----|------|-----|-----|-----|-----|-----|
| GC | VAN  | VAN | VAN | VAN | VAN | TPL  | TPL | TPL | TPL | TPL | TPL |
|    | 0.25 | 0.5 | 1   | 2   | 4   | 0.25 | 0.5 | 1   | 2   | 4   | 8   |

### Order Information

#### MICRONAUT MIC-Strip Vancomycin/Teicoplanin

5 plates, 8 strips per plate (40 tests per box) / Part No EM-022-040

#### Mueller-Hinton Broth, cation adjusted

1 tube per test, 20 tubes per box / Part No E2-331-020

1 tube per test, 100 tubes per box / Part No E2-331-100



Please contact your local representative for availability in your country

### MERLIN Diagnostika GmbH

Kleinstraße 14  
53332 Bornheim - Germany  
Phone +49 (0)2222-9631-0  
Fax +49 (0)2222-9631-90



[info.merlin@bruker.com](mailto:info.merlin@bruker.com) [www.merlin-diagnostika.de](http://www.merlin-diagnostika.de)



# MICRONAUT

- MIC-Strip Colistin

## What can it be used for?

The MICRONAUT MIC-Strip Colistin enables antimicrobial susceptibility testing against the back-up antibiotic colistin by using the broth microdilution (BMD) method, as advised by EUCAST, for *Enterobacteriaceae*, *Pseudomonas aeruginosa*, and *Acinetobacter baumannii* group.

Antimicrobial susceptibility testing of colistin has been routinely fraught with difficulties until now. A joint EUCAST and CLSI subcommittee recently issued recommendations confirming that

BMD is currently the only valid method for MIC determination of colistin. Gradient tests and disk diffusion methods are not recommended because of the poor diffusion of the large colistin molecule that could lead to an underestimated MIC value[1].

### Colistin concentration range

0.0625 to 64 mg/L

Manufactured by  
The Merlin logo, featuring a stylized blue star or 'M' shape to the left of the word "MERLIN" in a bold, blue, sans-serif font.

## Features and benefits

- Phenotypic detection of colistin resistance (e.g. MCR-1)
- Standardized MIC determination of colistin in routine laboratory diagnostics according to EUCAST recommendations
- The strip format permits economic use of the AST adapted to the individual laboratory diagnostic requirements
- Easy handling combined with standard laboratory equipment, for high quality MIC colistin results
- No special expertise required, easy visual evaluation

## Procedure

- Prepare a 0.5 McFarland bacterial suspension in NaCl
- Transfer an aliquot into Mueller-Hinton II broth (CAMHB)
- Inoculate the appropriate number of MIC-Strips
- Incubate for 18-22 hours at 35-37°C
- Read the results visually

## Shelf life and storage

- Shelf life: 24 months from date of production (unopened), expires 8 weeks after opening
- Storage: at room temperature (15-25°C)

## International ISO standard method

- MICRONAUT is a broth microdilution system using the international reference methodology (ISO 20776-1). MIC levels are based on EUCAST guidelines where applicable.

### Reference:

[1] Antimicrobial susceptibility testing of colistin – evaluation of seven commercial MIC products against standard broth microdilution for *Escherichia coli*, *Klebsiella pneumoniae*, *Pseudomonas aeruginosa*, and *Acinetobacter spp.* Matuschek E, Åhman J, Webster C, Kahlmeter G, CMI (2018), 24(8):865-870



### Colistin (COL) concentration (mg/L)

| 1  | 2      | 3     | 4    | 5   | 6   | 7   | 8   | 9   | 10  | 11  | 12  |
|----|--------|-------|------|-----|-----|-----|-----|-----|-----|-----|-----|
| GC | COL    | COL   | COL  | COL | COL | COL | COL | COL | COL | COL | COL |
|    | 0.0625 | 0.125 | 0.25 | 0.5 | 1   | 2   | 4   | 8   | 16  | 32  | 64  |

### Order Information

#### MICRONAUT MIC-Strip Colistin

5 plates, 8 strips per plate (40 tests per box) / Part No EM-006-040

#### Mueller-Hinton Broth, cation adjusted

1 tube per test, 20 tubes per box / Part No E2-331-020

1 tube per test, 100 tubes per box / Part No E2-331-100

Please contact your local representative for availability in your country



### MERLIN Diagnostika GmbH

Kleinstraße 14  
53332 Bornheim - Germany  
Phone +49 (0)2222-9631-0  
Fax +49 (0)2222-9631-90



info.merlin@bruker.com www.merlin-diagnostika.de



# MICRONAUT

- MIC-Strip Piperacillin-Tazobactam

## What can it be used for?

The MICRONAUT MIC-Strip Piperacillin-Tazobactam enables antimicrobial susceptibility testing against the combination medication piperacillin-tazobactam by using the broth microdilution (BMD) method for *Enterobacteriaceae* and *Pseudomonas aeruginosa*.

### Piperacillin-Tazobactam concentration range

0.125/4 to 128/4 mg/L

## Features and benefits

- Phenotypic antimicrobial susceptibility testing against piperacillin-tazobactam
- Standardized MIC determination of piperacillin-tazobactam in routine laboratory diagnostics by broth microdilution procedure
- The strip format permits economic use of AST, adapted to the individual laboratory diagnostic requirements
- Easy handling combined with standard laboratory equipment, for high quality MIC results of piperacillin-tazobactam
- No special expertise required, easy visual evaluation

## Procedure

- Prepare a 0.5 McFarland bacterial suspension in NaCl
- Transfer an aliquot into Mueller-Hinton II broth (CAMHB)
- Inoculate the appropriate number of MIC-Strips
- Incubate for 18-22 hours at 35-37°C
- Read the results visually

## Shelf life and storage

- Shelf life: 24 months from date of production (unopened), expires 8 weeks after opening
- Storage: at room temperature (15-25°C)

## International ISO standard method

- MICRONAUT is a broth microdilution system using the international reference methodology (ISO 20776-1). MIC levels are based on EUCAST guidelines where applicable.



### Piperacillin-Tazobactam (PIT) concentrations (mg/L)

| 1  | 2       | 3      | 4     | 5   | 6   | 7   | 8   | 9    | 10   | 11   | 12    |
|----|---------|--------|-------|-----|-----|-----|-----|------|------|------|-------|
| GC | PIT     | PIT    | PIT   | PIT | PIT | PIT | PIT | PIT  | PIT  | PIT  | PIT   |
|    | 0.125/4 | 0.25/4 | 0.5/4 | 1/4 | 2/4 | 4/4 | 8/4 | 16/4 | 32/4 | 64/4 | 128/4 |

## Order Information

### MICRONAUT MIC-Strip Piperacillin-Tazobactam

5 plates, 8 strips per plate (40 tests per box) / Part No EM-024-040

### Mueller-Hinton Broth, cation adjusted

1 tube per test, 20 tubes per box / Part No E2-331-020

1 tube per test, 100 tubes per box / Part No E2-331-100



Please contact your local representative for availability in your country

## MERLIN Diagnostika GmbH

Kleinstraße 14  
53332 Bornheim - Germany  
Phone +49 (0)2222-9631-0  
Fax +49 (0)2222-9631-90



[info.merlin@bruker.com](mailto:info.merlin@bruker.com) [www.merlin-diagnostika.de](http://www.merlin-diagnostika.de)

# MICRONAUT

DIAGNOSTICS WITH PASSION



**Use the reference method and fill the gap of your fully automated system**

**MICRONAUT systems for the identification and susceptibility testing of bacteria and yeast**

The testing principle of MICRONAUT identification systems (ID) is based on the phenotypical detection of the biochemical characteristics of microorganisms. The anti-microbial susceptibility test (AST) is based on phenotypical resistance detection by the growth of the microorganisms in the presence of the antibiotics / antifungal agents tested.

The microdilution method used here is a standardised one regarded worldwide as the recognised reference system for determining the minimum inhibitory concentration (MIC).

The substrates or antibiotics are present vacuum-dried in the microtitration plates. Thanks to a special vacuum drying procedure, these MICRONAUT plates can be stored at room temperature (15-25 °C).

Starting with a pure culture, the microorganisms are suspended in the corresponding MCN media and transferred to the MICRONAUT plates with an 8-channel pipette. Rehydration of the dried substances occurs by addition of the characterising bacteria or yeast suspension. After an incubation period of 18-24 hours the results can be photometrically measured or read visually.

As the central element in testing the MICRONAUT software combines device management, device communication and data analysis as well communication with the LIMS. The integrated

expert system analyses the photometric test measurements and assumes the clinical validation. As an option, the software can be connected bi-directionally with the laboratory EDP.

 **MERLIN**



# MICRONAUT

## MICRONAUT identification

The test principle of the MICRONAUT identification systems is based on phenotypic detection of the biochemical properties of microorganisms.

| Art. No.     | Product                  | Description  | Incubation time         | Evaluation                     | Packing                  |
|--------------|--------------------------|--|-------------------------|--------------------------------|--------------------------|
| M/E2-880-400 | <b>MICRONAUT-GNE</b>     | 4 tests/plate<br>Identification of <i>Enterobacteriaceae</i> and other gram-negative bacteria.   | 18 - 24 hrs             | photometrically                | 100 plates               |
| M/E2-710-400 | <b>MICRONAUT-IDS</b>     | 4 tests/plate<br>Identification of the most common clinically relevant <i>Enterobacteriaceae</i> , Non-fermenter, staphylococci, enterococci and streptococci. | 5 - 6 hrs               | photometrically                | 100 plates<br>+ 2 l NaCl |
| M/E2-520-120 | <b>MICRONAUT-NF</b>      | 3 tests/plate<br>Identification of non fermenting, gram-negative and some glucose-fermenting bacteria.   | 24 hrs                  | photometrically                | 40 plates                |
| M/E2-730-080 | <b>MICRONAUT-RPO</b>     | 2 tests/plate<br>Identification of staphylococci, streptococci, corynebacteria, <i>Listeria</i> , bacilli and other gram-positive bacteria.                    | 22 - 24 hrs             | photometrically                | 40 plates<br>+ 1 l NaCl  |
| M/E2-850-160 | <b>MICRONAUT-STAPH</b>   | 4 tests/plate<br>Identification of clinically relevant staphylococci.  | 6 hrs or<br>18 - 24 hrs | photometrically                | 40 plates<br>+ 1 l NaCl  |
| M/E2-870-160 | <b>MICRONAUT-STREP</b>   | 4 tests/plate<br>Identification of clinically relevant streptococci and enterococci.   | 20 - 24 hrs             | photometrically                | 40 plates<br>+ 1 l NaCl  |
| M/E2-740-160 | <b>MICRONAUT-Candida</b> | 4 tests/plate<br>Identification of clinically relevant yeasts.   | 24 hrs                  | photometrically                | 40 plates                |
| M/E5-005-200 | <b>MICRONAUT-UR</b>      | 2 tests/plate<br>Identification and susceptibility testing for urological diagnostics.   | 18 - 24 hrs             | visually or<br>photometrically | 100 plates               |
| M/E5-006-100 | <b>MICRONAUT-UR-KH</b>   | 1 test/plate<br>Identification and susceptibility testing for urine diagnostics, inclusive MRSA and ESBL detection.  | 18 - 24 hrs             | photometrically                | 100 plates               |

## MICRONAUT antimicrobial susceptibility tests (AST)

The MICRONAUT AST plates are offered in individual customer defined configurations, if an agreed minimum order quantity is being purchased. The customer may choose his preferred standard like EUCAST, or CLSI and compose his own customised antibiogram from more than 200 antibiotics. All AST products listed below are available as standard products. We continuously revise the range of standard AST plates. So please contact us for details on the actual program.

### Customer defined plates

| Art. No.     | Product                        | Description  | Incubation time | Evaluation                  | Packing    |
|--------------|--------------------------------|--|-----------------|-----------------------------|------------|
| M/E1-xxx-x00 | <b>MICRONAUT-S individual</b>  | 1 test/2 tests/4 tests/plate<br>AST as Breakpoint or MIC method with your choice of antibiotics (minimum order quantity has to be agreed). | 18 - 24 hrs     | visually or photometrically | 100 plates |
| M/EB-xxx-x00 | <b>MICRONAUT-SB individual</b> | 1 test/2 tests/4 tests/plate<br>AST as Breakpoint or MIC method with your choice of antibiotics (minimum order quantity has to be agreed). | 18 - 24 hrs     | visually or photometrically | 100 plates |

### Standard plates for clinical laboratories

| Art. No.     | Product                                       | Description  | Incubation time | Evaluation                  | Packing    |
|--------------|---|--|-----------------|-----------------------------|------------|
| M/EB-352-400 | <b>MICRONAUT-SB Standard Urine I (EUCAST)</b> | 4 tests/plate<br>AST as Breakpoint method for UTI relevant antibiotics, EUCAST standard.             | 18 - 24 hrs     | visually or photometrically | 100 plates |
| M/EB-380-200 | <b>MICRONAUT-SB Varia (ambulant) (EUCAST)</b> | 2 tests/plate<br>AST as Breakpoint method for mainly ambulant relevant antibiotics, EUCAST standard. | 18 - 24 hrs     | visually or photometrically | 100 plates |

# MICRONAUT

## Standard plates for veterinary laboratories

| Art. No.     | Product                      | Description   | Incubation time | Evaluation                  | Packing    |
|--------------|------------------------------|---|-----------------|-----------------------------|------------|
| M/E1-032-200 | <b>MICRONAUT-S Mastitis</b>  | 2 tests/plate<br>AST as MIC method for mastitis relevant antibiotics.           | 18 - 24 hrs     | visually or photometrically | 100 plates |
| M/E1-113-100 | <b>MICRONAUT-S Großtiere</b> | 1 test/plate<br>AST as MIC method for antibiotics for therapy of large animals. | 18 - 24 hrs     | visually or photometrically | 100 plates |
| M/E1-130-100 | <b>MICRONAUT-S Kleintier</b> | 1 test/plate<br>AST as MIC method for antibiotics for therapy of small animals. | 18 - 24 hrs     | visually or photometrically | 100 plates |

## Standard plates for antifungal susceptibility testing of yeasts

| Art. No.     | Product                                    | Description  | Incubation time | Evaluation                  | Packing   |
|--------------|--|--|-----------------|-----------------------------|-----------|
| M/E1-824-160 | <b>MICRONAUT-AM</b>                        | 4 tests/plate<br>AST as Breakpoint + MIC method for antifungal agents. | 22 - 48 hrs     | visually or photometrically | 40 plates |
| M/E1-831-040 | <b>MICRONAUT-AM Anti Fungal Agents MIC</b> | 1 test/plate<br>AST as MIC method for antifungal agents (EUCAST).      | 22 - 48 hrs     | visually or photometrically | 40 plates |
| M/E1-832-080 | <b>MICRONAUT-AM EUCAST AFST</b>            | 2 tests/plate<br>AST as MIC method for antifungal agents (EUCAST).     | 22 - 48 hrs     | visually or photometrically | 40 plates |

**NEW**

## MICRONAUT MIC Strips

| Art. No.     | Product  | Description   | Evaluation | Packing  |
|--------------|--|---|------------|--|
| M/EM-006-040 | <b>MICRONAUT MIC Strip Colistin</b>                | 40 tests<br>Susceptibility testing of bacteria by the broth microdilution method to determine the minimum inhibitory concentration (MIC). <u>According to the EUCAST and CLSI recommendations of March 2016 and the study of November 2016.</u> | visually   | 5 x 8 strips/<br>tests with<br>12 wells per<br>strip |
| M/EM-022-040 | <b>MICRONAUT MIC Strip Vancomycin/ Teicoplanin</b> | 40 tests<br>Broth microdilution method for detection of glycopeptide resistant bacteria by determination of the minimal inhibitory concentration (MIC) of Teicoplanin and Vancomycin.   | visually   | 5 x 8 strips/<br>tests with<br>12 wells per<br>strip |

## Special plates for clinical laboratories

| Art. No.     | Product  | Description  | Incubation time | Evaluation                  | Packing    |
|--------------|--|--|-----------------|-----------------------------|------------|
| M/E1-111-040 | <b>MICRONAUT-S<br/>Beta Lactamase</b>                  | 1 test/plate<br>Phenotypic detection of ESBL (extended spectrum beta-lactamase), MBL (metallo-beta-lactamase), KPC (Klebsiella pneumoniae carbapenemase), AMP-C (amino-penicillin-deactivating cephalosporinase) and D-carbapenemases (OXA-48) for all relevant gram-negative bacteria in a single system. | 18 - 24 hrs     | visually or photometrically | 40 plates  |
| M/E1-055-040 | <b>MICRONAUT-S MRSA / GP</b>                           | 1 test/plate<br>Detection of multidrug-resistant staphylococci (MRSA), enterococci (VRE) and pneumococci, including testing of novel antibiotics (e.g. daptomycin, ceftaroline).   | 18 - 24 hrs     | visually or photometrically | 40 plates  |
| M/E1-114-040 | <b>MICRONAUT-S MDR<br/>MRGN-Screening</b>              | 1 test/plate<br>Susceptibility testing of multidrug-resistant gram-negative bacteria. Phenotypic detection of AMP-C cephalosporinase and type A-, B- and D-carbapenemases ( <u>new: ceftolozane/ tazobactam and ceftazidime/ avibactam</u> ).  | 18 - 24 hrs     | visually or photometrically | 40 plates  |
| M/E1-085-040 | <b>MICRONAUT-S<br/>Anaerobier MHK</b>                  | 1 test/plate<br>Susceptibility testing of treatment-relevant antimicrobial agents, such as tigecycline, moxifloxacin or ertapenem.   | 18 - 24 hrs     | visually or photometrically | 40 plates  |
| M/E1-973-040 | <b>MICRONAUT-S<br/>Campylobacter</b>                   | 1 test/plate<br>Determine resistances by measuring MIC for all relevant antimicrobial agents against <i>Campylobacter</i> .  | 18 - 24 hrs     | visually or photometrically | 40 plates  |
| M/E1-129-100 | <b>MICRONAUT-S Pneumo-<br/>cocci &amp; Haemophilus</b> | 1 test/plate<br>Determine resistances by measuring MIC for all relevant antimicrobial agents against pneumococci and <i>Haemophilus</i> .  | 18 - 24 hrs     | visually or photometrically | 100 plates |
| M/E1-099-100 | <b>MICRONAUT-S MHK<br/>Pseudomonas</b>                 | 1 test/plate<br>Determine resistances by measuring MIC for pseudomonas-active antimicrobial agents. ( <u>new: ceftolozane/ tazobactam and ceftazidime/ avibactam</u> ).  | 18 - 24 hrs     | visually or photometrically | 100 plates |
| M/E1-981-040 | <b>MICRONAUT-S CF</b>                                  | 1 test/plate<br>Testing of multiple drug resistant non-fermenters from patients with cystic fibrosis.  | 18 - 24 hrs     | visually or photometrically | 40 plates  |

# MICRONAUT

## Special plates for interlaboratory comparison

| Art. No.     | Product                                      | Description   | Incubation time | Evaluation                  | Packing  |
|--------------|--|---|-----------------|-----------------------------|----------|
| M/EB-379-005 | <b>MICRONAUT-SB<br/>Ringversuch Urologie</b> | 1 test/plate<br>AST towards MIC for external laboratory control for urologists. | 18 - 24 hrs     | visually or photometrically | 5 plates |
| M/E1-028-005 | <b>MICRONAUT-S MHK<br/>Ringversuch</b>       | 1 test/plate<br>AST towards MIC for external laboratory control.                | 18 - 24 hrs     | visually or photometrically | 5 plates |

## MICRONAUT reagents

Some identification tests require adding of reagent before measuring.

| Art. No.     | Product                  | Associated MCN plate   | Tests                    | Packing      |
|--------------|--------------------------|--|--------------------------|--------------|
| M/E2-301-001 | <b>Indol Reagent</b>     | MICRONAUT-IDS/UR/NF/GNE  | 400                      | 100 ml       |
| M/E2-303-001 | <b>Nitrate Reagent A</b> | Optional for MICRONAUT ID plates:<br>MICRONAUT-RPO/IDS/UR/NF/GNE |                          | 100 ml       |
| M/E2-304-001 | <b>Nitrate Reagent B</b> | Optional for MICRONAUT ID plates:<br>MICRONAUT-RPO/IDS/UR/NF/GNE |                          | 100 ml       |
| M/E2-308-001 | <b>Ninhydrin Reagent</b> | MICRONAUT-STREP  |                          | 2 x 4 ml     |
| M/E2-305-001 | <b>Paraffin Oil</b>      | MICRONAUT-IDS/NF/STREP/RPO/GNE/<br>STAPH/UR                      | Depending on the product | 100 ml       |
| M/E2-310-001 | <b>Peptidase Reagent</b> | MICRONAUT-IDS/RPO/UR   | Depending on the product | 100 ml       |
| M/E2-312-001 | <b>NaCl</b>              | MICRONAUT identifications  | Depending on the product | 1 x 1000 ml  |
| M/E2-323-001 | <b>AST Reagent Kit</b>   | MICRONAUT-AM   |                          | 2 (4 x 4 ml) |

## MICRONAUT media

Bacterial suspension is prepared in MCN media.

| Art. No.     | Product  | Associated MCN plate                   | Tests                         | Packing           |
|--------------|--|--|-------------------------------|-------------------|
| M/E2-306-100 | <b>MICRONAUT-NF Susmed</b>                         | MICRONAUT-NF                           | 100                           | 100 x 6 ml        |
| M/E2-314-100 | <b>MICRONAUT-Candida Susmed</b>                    | MICRONAUT-Candida                      | 100                           | 100 x 6 ml        |
| M/E2-331-020 | <b>Mueller Hinton Broth, cation adjusted</b>       | MICRONAUT-S                            | 20                            | 20 x 11 ml        |
| M/E2-331-100 | <b>Mueller Hinton Broth, cation adjusted</b>       | MICRONAUT-S                            | 100                           | 100 x 11 ml       |
| M/E2-311-100 | <b>MICRONAUT-H-Medium</b>                          | MICRONAUT-S, fastidious microorganisms | 100                           | 100 x 11 ml       |
| M/E2-330-020 | <b>MICRONAUT-Wilkins-Chalgren Broth</b>            | MICRONAUT-S Anaerob                    | 20                            | 20 x 11 ml        |
| M/E2-324-020 | <b>MICRONAUT-RPMI-1640 Medium + MOPS + Glucose</b> | MICRONAUT-AM                           | 20                            | 20 x 11 ml        |
| M/E2-319-100 | <b>MICRONAUT-SB Medium</b>                         | MICRONAUT-SB                           | 100                           | 100 x 11 ml       |
| M/E2-338-100 | <b>MICRONAUT-ID Medium</b>                         | MICRONAUT-UR                           | 100                           | 100 x 5.5 ml      |
| M/E2-337-100 | <b>MICRONAUT-AST Medium</b>                        | MICRONAUT-UR                           | 100                           | 100 x 5.5 ml      |
| M/E2-318-010 | <b>MICRONAUT-SB Medium dehydrated</b>              | MICRONAUT-SB                           | Depending on the product used | 10 x for 5 l each |

# MICRONAUT

## MICRONAUT software

MICRONAUT software offers after automatized reading, calculation, and interpretation of identification and susceptibility testing by using MICRONAUT systems.

| Art. No.     | Product                   | Description                            | Packing |
|--------------|---------------------------|--|---------|
| M/U8-305-001 | <b>MICRONAUT Software</b> | Basic module                           | 1 unit  |
|              | <b>QS Module</b>          | Record of quality control data.        | 1 unit  |
|              | <b>Statistic Module</b>   | Statistical record of resistance data. | 1 unit  |

## Instruments

| Art. No.      | Product                                       | Description   | Packing |
|---------------|---|---|---------|
| M/BH3-880-001 | <b>Pipette</b>                                | Electronic 8-channel stepper for fast inoculation of 12 x 100 µl.                           | 1 unit  |
| M/LT51119000  | <b>Reader</b>                                 | Reader for fast measurement of MCN plates.  | 1 unit  |
| M/PRE-001-001 | <b>Precision™ XS</b>                          | Pipetting instrument for workflow of MCN AST plates.  | 1 unit  |
| M/L4Y-100-001 | <b>Densitometer</b>                           | Tube densitometer for adjustment of bacterial suspensions.                                  | 1 unit  |
| M/2350        | <b>McFarland Standard Set 0.5 / 1.0 / 2.0</b> | McFarland standards are used as a reference to adjust the density of bacterial suspensions. | 3 tubes |

## Consumables

| Art. No.      | Product   | Description  | Packing               |
|---------------|---|--|-----------------------|
| M/ST3-001-500 | <b>Matrix pipette tips</b>                      | For use with the Matrix pipette.   | 500 units             |
| M/LH-B791204  | <b>Biohit Optifit Tips Flexi-Bulk 1200</b>      | For use with the Biohit pipette.   | 480 units             |
| M/BH3-487-096 | <b>Biohit tips rack 96</b>                      | For use with the Biohit pipette.   | 1 rack<br>per 96 tips |
| M/R4-510-050  | <b>1-Channel reservoirs, autoclavable</b>       | For easy plate inoculation with 1 test/plate.                              | 50 units              |
| M/R4-510-350  | <b>1-Channel reservoirs, disposable product</b> | For easy plate inoculation with 1 test/plate.                              | 350 units             |
| M/R4-506-050  | <b>2-Channel reservoirs, autoclavable</b>       | For easy plate inoculation with 2 tests/plate.                             | 50 units              |
| M/R4-506-350  | <b>2-Channel reservoirs, disposable product</b> | For easy plate inoculation with 2 tests/plate.                             | 350 units             |
| M/R4-508-050  | <b>4-Channel reservoirs, autoclavable</b>       | For easy plate inoculation with 4 tests/plate.                             | 50 units              |
| M/R4-508-350  | <b>4-Channel reservoirs, disposable product</b> | For easy plate inoculation with 4 tests/plate.                             | 350 units             |
| B2-003-040    | <b>MICRONAUT foil perforated, 40 units</b>      | For sealing of MCN identification plates (included in standard packaging). | 40 units              |
| B3-002-040    | <b>MICRONAUT foil unperforated, 40 units</b>    | For sealing of MCN AST plates (included in standard packaging).            | 40 units              |

The MICRONAUT system prices are available on request. Please feel free to contact us.