de

Anti-Salmonella H-Phasen-Pools

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en

Anti-Salmonella H-Phase Pools

Test reagent for slide agglutination

INFORMATION FOR PROFESSIONAL USE

REF	H-Phase Pool	Vol.	
TR1181, TR1181-01	Anti-Salmonella HMA	a, b, c, d, i, Z ₁₀ , Z ₂₉	1 ml, 3 ml
TR1183, TR1183-01	Anti-Salmonella HMB	E, G	1 ml, 3 ml
TR1185, TR1185-01	Anti-Salmonella HMC	k, y, z, L, Z ₄ , r	1 ml, 3 ml

Intended use

The test reagents are intended for use in identifying H-antigens in the serological identification of *Salmonella* strains iso-lated from test material of varied origin, using slide agglutination.

Additional serological differentiation must be performed using the monospecific Anti-Salmonella H test reagents in accordance with the White-Kauffmann-Le-Minor scheme (Kauffmann-White scheme).

Principle of the test

If the isolated *Salmonella* strain possesses an H-antigen that is covered by the test reagent, this Salmonella antigen becomes bound when mixed with the specific antibody. The antigen-antibody reaction results in clearly visible agglutination of the strain.

Composition

The test reagents are mixtures of monoclonal antibodies. They are produced from cell culture supernatants of hybridoma cell lines which secrete antibodies against the respective *Salmonella* H-antigens.

For the coverage of the H-antigen z₂₉ the Anti-Salmonella HMA test reagent and for the coverage of the H-antigen complex Z₄ the Anti-Salmonella HMC test reagent contain additionally serum of immunized rabbits from which unspecific agglutinins have been removed by absorption.

Preservative: Sodium azide (NaN₃) 0.9 mg/ml

Form in which product is supplied, shelf life and conditions for storage

The test reagents are available in liquid form ready to use.

If stored at 2...8 °C before and after opening, they may be used up to the date given on the label. After use, the bottle must be properly closed.

The test reagent may sometimes show turbidity not caused by microbes. Such turbidity does not impair effectivity and the test reagents can be clarified by centrifugation or filtration.

The test reagents must have their temperatures adjusted to room temperature (18...26 °C) before use.

Warnings and precautions

The biotechnological manufacture of the test reagents means that the risk of contamination by infectious agents can be virtually excluded. Because they contain animal materials (fetal calf serum, stabilizer) and because of the addition of rabbit serum, they should be treated as potentially infectious and handled accordingly.

As these materials contain sodium azide, contact with the skin and mucous membranes must be avoided! In case of contact, rinse with plenty of water.

Since the performance of the slide agglutination test involves working with native pathogenic materials, all necessary work protection procedures must be adhered to (risk of infection)!

Materials and equipment not supplied

Glass slides, stirring rods, physiological saline, disposal containers for infectious material, White-Kauffmann-Le-Minor scheme.

Test material and methodology

1. Take a small amount of bacterial mass from a 16-20 hour old subculture (e.g. Kligler culture medium), picked from the lower, moister layer of the agar slope. Transfer onto a slide and mix well with a drop (ca. 25 µl) of test reagent, so that a homogenous, slightly milky suspension results. The slide should be placed on a dark surface.

The reaction is read with the naked eye by holding the slide in front of a light source against a black background and rocking it (tilting it back and forth). To control for spontaneous agglutination, a negative control must be performed at the same time using physiological saline in place of the test reagent.

2. In case the H-antigen test performed according to the method described above is not successful, the strain must be inoculated onto swarm agar, to enhance the development of the flagella antigens. For the swarm agar we recommend our ready-to-use culture medium REF sifin TN1702.

After liquefaction, place 10 ml of this swarm agar into a petri dish (6 cm diameter). Once it has solidified, sprinkle the surface of the agar with ca. 100 µl sterile distilled water and apply the strain in the form of a point in the centre of the plate. Incubate the dish overnight with the agar layer facing down at 35...37 °C. For the slide agglutination test, collect material from the edge of the dish.

Test procedure

The test, using Anti-Salmonella HMA, Anti-Salmonella HMB or Anti-Salmonella HMC test reagents, enables the serological identification of the commonest H-antigens or H-antigen complexes of *Salmonella* strains. Furthermore the presence of the antigen complex H:1 should be tested by Anti-Salmonella H:1 (REF TR1437, TR1437-01).

If the test reactions are positive, the strain must then be agglutinated with the monospecific Anti-Salmonella H test reagents for further serotyping. These products are available at sifin (see product list or www.sifin.de).


Evaluation

The test can only be evaluated if the negative control (NC) remains milky-opaque.

Positive: visible agglutination after the sample has been tilted back and forth less than 20 times. In a strongly positive reaction (PR), agglutination (coarsely or finely flocculent) appears as soon as the bacterial mass is mixed in. In a weakly positive result, agglutination only appears after the slide has been tilted back and forth 10-20 times.

PR		NC
-----------------	---	-----------------

Negative: if the suspension remains milky-opaque or the reaction begins to occur only after the slide has been tilted back and forth more than 20 times, the result is negative (NR).

NR		NC
-----------------	---	-----------------

Quality assurance during the testing procedure

For the quality control of serological determination of *Salmonella* H antigens by slide agglutination test, the good expression of the strain's cell surface antigen is important. The use of strains from interlaboratory tests, or field strains of defined origin that have been fully characterised by an external laboratory, or the use of suitable commercially available *Salmonella* test antigens is therefore recommended for quality control.

Limitations of the method

The test reagents react with *Salmonella* strains, which contain antigens with specificities, declared in the H-Phase Pools. A biochemical examination of the genus *Salmonella* exclude rare cross-reactions with strains of other genera of the family of *Enterobacteriaceae*.

Explanation of the symbols used

LOT	Batch code (Lot)		Use by YYYY-MM (MM = end of month)
REF	Catalogue number		Temperature limitation
IVD	In Vitro Diagnostic Medical Device	 </	

sifin diagnostics gmbh | Berliner Allee 317-321 | 13088 Berlin | Germany | T: +49 30 700 144-0 | F: +49 30 700 144-30 | info@sifin.de | www.sifin.de

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