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2.1 Warnings and precautions

**Maintenance and defects**

Do not use the equipment for any application before the user correctly performs all regular checks and updates the periodical equipment maintenance. If it is sure (or probable) that any part of the equipment is defective or wrong adjusted, don't use it before performing all reparations.

The use of an equipment with defective parts or adjusted in a wrong way, can expose the user or the patient to ionizing radiations or to other dangers concerning safety. This can cause serious or mortal physical injuries, or wrong diagnosis or therapies.

Importance of safety

Do not use the equipment for any application before reading, understanding and assimilating all information about safety, safety and emergency procedures specified in the current chapter about Safety, The use of the equipment without a proper knowledge of safety rules can cause serious or mortal physical injuries, or wrong diagnosis or therapies.

Proper training

Do not use the equipment for any application unless you have a proper and adequate training to a safe and efficient use.

If you aren't sure to be able to use this equipment in a safe and efficient way, don't use it. The use of this equipment without proper and adequate training can cause serious or mortal physical injuries or wrong diagnosis or therapies.

Do not use the equipment with the patients if there is no adequate understanding of its capabilities and functions. Using the equipment without an adequate knowledge of its functioning can compromise the efficacy and/or reduce the safety of the patient, the user and other people nearby.

Safety systems

Never try to remove, modify, exclude or obstruct any safety device on the equipment. An intervention on safety devices can cause serious physical injuries or even death.

Expected use and compatibility

Do not use the equipment for purpose other than those for which it is intended. Do not use the equipment with other products than the ones whose compatibility has been recognized by the Manufacturer. The use of the equipment for purposes other than the ones expected or with an incompatible product, can cause serious or mortal physical injuries or wrong diagnosis or therapies.

This equipment must be used only in compliance with the safety instructions specified in this User's Manual and exclusively for intended purposes.

It is user's responsibility to ensure that effective norms concerning installation and use of medical equipment are respected.

The Manufacturer is responsible for safety features of its own products, only provided that maintenance, repairs and modifications are performed exclusively by the Manufacturer's personnel or by personnel expressly authorized by the Manufacturer.

As for all technical equipments, even this medical device must be used properly and subject to regular maintenance and care, as described in " Maintenance, cleaning and disposal" paragraph.

The Manufacturer can't be considered responsible for any error, damage or injury caused by improper use or lack of maintenance of the equipment.

It is necessary to contact the assistance service authorized by the Manufacturer even in the case no error messages are displayed, but the equipment doesn't work as usual (first symptoms of a fault).

Do not modify or remove in any way the safety circuits.

2.2 Electrical Safety

This equipment is in compliance with safety class I, Type B, in accordance with IEC 60601-1 norm.

Do not use the equipment near or leaned against other equipments.

Do not remove protections or cables from this equipment, unless it is expressly required in this User's Manual, because inside it there are dangerous electrical voltages. The removal of protections or cables can cause mortal injuries or serious damages to the people.

Protections or cables must be removed only by qualified and authorized technical personnel. Use the equipment only in rooms or areas comply with all applicable laws (or regulations having the force of law), referring to electrical safety of this type of medical device.

Always insulate the equipment from the power supply before proceeding with cleaning or disinfection operations in order to avoid electric shocks.

Equipotential earth connection

The equipment is supplied with an equipotential earth connection point.

The equipment can be used only in areas comply with local electrical safety norms and in environments suitable for medical activities. Besides IEC 60601-1 norm provides instructions about the equipotential earth connection point.

Additional equipotential earth connection

An additional equipotential earth connection is provided because the equipment is movable and the reliability of the main equipotential earth connection point can be insufficient.

It is possible to use this equipment only in rooms comply with IEC norm requirements.

The equipment described isn't protected against liquids seepage. Its classification is IPx0.

2.3 Mechanical safety

Be sure that parts of the body or clothes aren't stuck among moving components of the equipment.

Remove all objects from range of motion of the equipment.

Check that the unused hanging components (monitor and radiogenic complex) are positioned so as not to affect neither the user nor the patients.

It is not possible to transport this equipment while it is working. For a safety transport, switch off the equipment before transporting it and ensure that all system peripherals (monitor, mouse, keyboard, cables etc.) are disconnected.

Do not remove protections or cables from this equipment, unless this operation is expressly requested in this User's Manual.

The equipment includes moving parts. The removal of protections can cause serious or mortal physical injuries to people.

2.4 Protection against explosions

This device mustn't be used in presence of explosive gas or fumes, such as some kind of gaseous anesthetics. Do not use disinfectant spray flammable or potentially explosive. The use of this equipment in an unsuitable environment can cause fires or explosions.

2.5 Fire safety

- Do not use this equipment in areas where there is a risk of fire.
- Do not cover the ventilation openings while the equipment is turned on.
- For electrical or chemical fires use only fire extinguisher marked as suitable for such uses. The use of water or other liquids in an electrical fire can cause physical injuries or even death.
- Before trying to extinguish the fire, the safety measure to be taken is to separate the equipment from other electric power sources and from all other sources in order to reduce the risk of electrical shocks.

2.6 Electrostatic discharge (ESD)

Always resort to static procedures, protections and appropriate products before opening or during the handling of the equipment. This equipment includes electrostatically sensitive components. Non-compliance with ESD procedures may cause damages to the components. Such damages to the components are not covered by any warranty.

The electrostatic discharge (ESD) can cause a remarkable voltage that could cause damages to printed circuit boards (PCB) or to other equipment components.

Electrostatic discharge (ESD) damages can accumulate and can initially not be visible, such as a hardware failure, but can reduce performances. Therefore, it is recommended to use proper ESD handling procedures. ESD can be due to low humidity or to the use of electrical equipment on carpets, bedding and clothes.

2.7 Electromagnetic compatibility (EMC)

This equipment complies with international and national laws and regulations relating to electromagnetic compatibility (EMC) in force for this type of product, if it is used for the intended purposes. Such laws and regulations define the electromagnetic emissions level coming from the product and the requested immunity against electromagnetic interferences from external sources. Other electronic products that exceed the limits defined by EMC standards can, in unusual situations, affect on the equipment working.

- Electromedical products request special precautions referring to electromagnetic compatibility (EMC) and must be installed and started up in compliance with EMC information provided in the documentation enclosed.
- The use of accessories and cables other than those specified can cause a higher emission or lowest immunity levels.
- The equipment mustn't be used in proximity of other products or stacked on them and, if this will be necessary, you must check the right functioning.

Mobile phones and laptops

Communications among RF portable and mobile equipments can affect medical equipments. It is recommended to use caution while using such communication devices within the specified radius of electromedical devices

2.7.1 Warning and Safety Precautions for Electromagnetic Compatibility

Increased emission or reduced interference immunity.

Use of unsuitable accessory or lines

- Exclusive use of the listed accessory or line with the exception of internal original spare part components.

Electric medical units are subject to special precautionary measures with regard to EMC and may only be installed and put into operation in compliance with the EMC information contained in the Operating Manual. Portable and mobile radiofrequency communication devices can influence electric medical devices

Annex A

The Equipment is intended for use in the electromagnetic environment specified below. The customer or the user of the Equipment should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic Environment
RF emissions CISPR 11	Group 1	This Equipment uses RF energy only for its internal function. Therefore, the RF emission is very low and not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Group 2	This Equipment must emit electromagnetic energy in order to perform its intended function. Nearby electronic equipment may be affected.
RF emissions CISPR 11	Class A	This Equipment is suitable for use in domestic establishments and in establishments directly connected to the low voltage power supply network which supplies buildings used for domestic purposes
Harmonic emissions EN 61000-3-2	Not Applicable	
Voltage fluctuations/ flicker emissions EN 61000-3-3	Not Applicable	

Annex B

The Equipment is intended for use in the electromagnetic environment specified below. The customer or user of the Equipment should assure that it is used in an electromagnetic environment.			
Immunity Test	EN 60601-1-2 Test level	Compliance level	Electromagnetic Environment - guidance
Electrostatic discharge (ESD) EN 61000-4-2	± 6 kV contact ± 8 kV air	EN 60601-1-2 Test level	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst EN 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines > 3 m	EN 60601-1-2 Test level	Mains power quality should be that of a typical commercial or hospital environment.
Surge EN 61000-4-5	± 1 kV differential mode ± 2 kV common mode	EN 60601-1-2 Test level	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines EN 61000-4-11	0 % U_n for 0,5 cycle 40 % U_n for 5 cycle 70 % U_n for 25 cycle 0 % U_n for 5 s	EN 60601-1-2 Test level	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Equipment requires continued operation during power mains interruptions, it is recommended that the Equipment be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field EN 61000-4-8	3 A/m	EN 60601-1-2 Test level	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Annex C

non-LIFE SUPPORTING EQUIPMENT

The Equipment is intended for use in the electromagnetic environment specified below. The customer or the user of the Equipment should assure that it is used in such an environment			
Immunity Test	EN 60601-1-2 Test level	Compliance level	Electromagnetic Environment
			<p>Portable and mobile RF communications equipment should be used no closer to any part of the Equipment, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter</p> <p>Recommended separation distance</p>
Radiated RF EN 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = 1.2 \times \sqrt{P}$ 80 MHz to 800MHz $d = 2.3 \times \sqrt{P}$ 800 MHz to 2.5GHz
Conducted RF EN 61000-4-6	3 V 150 kHz to 80 MHz	3 V	$d = 1.2 \times \sqrt{P}$
			<p>Where P is the maximum output rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths for fixed RF transmitter, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 

Annex D

Recommended Separation Distance for non-LIFE SUPPORTING EQUIPMENT

The Equipment is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Equipment can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitter) and the Equipment as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of the transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150kHz to 80MHz $d = 1.2 \times \sqrt{P}$	80MHz to 800MHz $d = 1.2 \times \sqrt{P}$	800MHz to 2.5GHz $d = 2.3 \times \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
<p>For transmitters rated at the maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.</p> <p>Note:</p> <p>(1) at 80MHz and 800MHz, the separation distance for the higher frequency range applies</p> <p>(2) These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people</p>			

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Only personnel trained in the management of cleaning and disinfection of medical devices is authorized to conduct such activities.

Perform regularly cleaning and disinfection operations of the equipment.

Below are the instructions.

Always disconnect the equipment from the power supply before proceeding with cleaning and disinfection operations in order to avoid electrical shocks.

Avoid the seepage of water and liquids because it can cause short-circuits or corrosion of metallic parts.

Cleaning and disinfection operations, even for the equipment and for the environment, must be in compliance with all laws and norms in force in the country where the equipment is installed.

Cleaning

Enameled parts and aluminum surfaces must be cleaned only with a damp cloth and a mild detergent and then with a dry woolen cloth. Never use scouring powders, solvents, abrasives detergents or polishing abrasive. Do not use a special detergent if its properties are not sure.

Chromed parts must be cleaned only with a dry woolen cloth. Do not use polishing abrasives. To protect the finish, use a nonabrasive wax.

Plastic surfaces must be cleaned only with soap and water. When using other cleaning agents (for example with a high alcohol content), the material can become opaque or can break.

Disinfection

The disinfection method used must be in compliance with all laws and norms in force for disinfection and protection against explosions in force in the country where the equipment is installed.

All parts of the equipment suitable for this type of treatment, accessories and connection cables included, can be disinfected with a damp cloth and a proper detergent. Never use disinfecting agents or corrosive sterilizers or solvents.

Do not use a special disinfecting or sterilizing agent if its properties are not sure.

Do not use inflammable disinfectant spray or potentially explosive. Such sprays create gas that can ignite, causing serious injuries or even death.

It is not recommended to disinfect using a spray in a room where there are medical products, because the gas can penetrate the product, causing short-circuits, corrosion of metallic parts or other damages the equipments.

If it is necessary to use non-inflammable and non-explosive sprays, first of all switch off and cool down the equipment.

In this way the vaporized spray can't be attracted by convention currents inside the equipment. Before starting spraying, it is necessary to cover carefully the product with plastic sheeting.

Once all traces of disinfecting spray disappear, it is possible to remove the protective plastic and directly disinfect or sterilize the equipment following the recommended instructions.

After using a spray the user must be sure that every single trace of gas has disappeared before starting up the equipment again.

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The manufacturer wants to make a contribution to environment defense and wants to guarantee a constantly safe and efficient use of this equipment by using a proper support, maintenance and training program.

If the equipment is used correctly and always subjected to proper maintenance, it doesn't represent an environmental risk. However it can include materials that can be potentially harmful for the environment if they are not properly disposed.

The use of such materials is essential for carrying out the equipment functions in compliance with legal requirements and so on.

Final disposal of the equipment

The final disposal is effected when the equipment has been used so that it is no longer usable for the intended purposes.

The return, proper disposal or recovery of this medical equipment must be done in compliance with the European WEEE (Waste Electrical and Electronic Equipment) and / or national requirements.

The equipment or parts of it mustn't be disposed as industrial or domestic waste, but they must be collected separately as special waste. The separate collection for the subsequent forwarding for recycling, treatment and environmentally compatible disposal, helps to avoid possible negative environmental and health effects and to promote recycling of the parts included in the equipment.

Illegal disposal of the equipment involves the application of administrative sanctions according to the current regulations of the country where the equipment is installed.

For information on how to dismantle the inoperative equipments comply with local legislation or contact an authorized representative of the Manufacturer.

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Error messages

Typ.	Text	Meaning	Intervention
F	FAULTY POWER	Error on the Charger or Chopper. Energy not available	Switch OFF, wait some minutes, switch ON, if the error appears gain call Service
F	ERROR V3	Absent V3 power supply	Switch OFF, wait some minutes, switch ON, if the error appears gain call Service
F	FILAMENT	Filament current not present	Switch OFF, wait some minutes, switch ON, if the error appears gain call Service
F	ERROR V2	Power supply V2 not present in the circuit of mA and kV set	Switch OFF, wait some minutes, switch ON, if the error appears gain call Service
F	LOCKED STARTER	Error during the start time	Press the RESET key to go on, repeat x-rays
F	FAULTY CHOPPER	Error in the chopper	Press the RESET key to go on, repeat x-rays
F	X-RAY ERROR	The kV have not reached 75% of the value set within the first 10mS of exposure or lackage of kV	Press the RESET key to go on, repeat x-rays
F	MAX. TIME	The max. exposure time has been reached.	Press the RESET key to go on, repeat x-rays
F	DATA ERROR	Memory error, data checksum error	Switch the equipment OFF, wait some minutes and then switch it ON again, if the error appears again call Service.
F	INVERTER KV ERROR	During the XR emission the kV are decreased under 75% or increased over 110% of the value set or the high voltage circuit unbalanced during the exposure	Press the RESET key and repeat the exposure
F	INV. OVERLOAD	Power to the inverter out of scale	Press the RESET key and go on
F	INVERTER ERROR	IGBT drivers error	Press the RESET key and go on
F	TUBE CALIB. ERR.	Error of XR tube calibration	Call Service.
F	XR HANDSWITCH ERR.	Faulty x-ray handswitch	Check the integrity of the x-ray handswitch, switch OFF and switch ON the equipment, try again, if the error persists call Service

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7.9 Опция: radiography with examination table or Potter Bucky grid

The examination table of the Potter Bucky must be in compliance with IEC 60601-1 standard and must be connected to the equipment according to 93/42/CE medical devices directive and following amendments.

In LOCAL mode: fully extend the extensible cable of the x-ray control and stay as far as possible from the radiation source.


In REMOTE mode: stay as far as possible from the radiation source.

The mode "Radiography with examination table or Potter Bucky grid" is active in CR mode or with Digital Imaging System OFF.

1. Place the cassette and the patient according to requirements.
2. Connect the examination table or the Potter Bucky grid to the connector (A) placed on the front part of the equipment.

If the connection with the examination table or with the Potter is right, the display shows the relative icon.

3. Press the (A) push-button to activate the Potter Bucky.

4. Place the monobloc and collimator, set the exposure field.
5. Select kV and mAs values in manual or APR mode by selecting the data concerning the exam to be performed.
6. Hold the wired x-ray emission control or the remote control (optional)
7. Go as far as possible from the x-ray source.
8. Perform the exposure. At the end of the exposure the acoustic signal of the equipment emits three beeps.

When the examination table or the Potter Bucky grid have been selected, one of the most frequent difficulties is the "no consent to proceed" of the Potter Bucky grid: in this case check the connection.

7.10 Optional : DAP meter

The installation and maintenance of the DAP meter can be performed only by authorized service personnel.

To correctly determine the dose value to which the patient has been exposed, it is necessary to reset the dose value before starting the exam.

The equipment can be supplied with a ionization chamber dosimeter (dose-area product meter, DAP meter). The DAP meter can be installed during the preparation or in a following phase.

The function of the DAP meter is to measure the dose-area product (cGycm²) that is going out towards the patient:

$$\sum_{i=1}^n dose_i \cdot area_i$$

"n" is the number of exposures performed after the last pressing of F5-RESET push-button.

READY 0000.00

READY cGycm2

If the DAP meter is installed and worked properly, the first line of the display shows the standard of measurement (cGycm²) and the sum of the measures.

The sum of the measures is the sum of all dose-area products read by the chamber.

The value is reset by pressing the F1+RESET pushbutton.

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