



Dolphin Transcranial Doppler

Taking TCD Technology One Step Further



Discover the Most Advanced Bilateral TCD Robot in The Market





www.viasonix.com

Everything You Expect from a State-of-the-Art TCD System and Much More...

- ✓ Slim design, lightweight and battery operated
- ✓ Scroll back in time and save Doppler spectrums from different depths along the Ultrasound beam.
- ✓ View a suspected emboli travelling along the blood vessel
- ✓ Operate the Dolphin with a single operation button
- ✓ Receive immediate statistical database information
- ✓ Send reports and automatic worksheets to EMR and PACS
- And much more!

Experience the Power of the Dolphin

Clinical Applications & Specialty Tests

Routine Cerebral Examinations

Stroke Management Vasospasm and SAH Arteriovenous Malformation **Intracranial Stenosis Collateral Capacity**

Specialty Tests

VMR / Autoregulation **Breath Holding Test** Patent Foramen Ovale (PFO) Sickle Cell Disease **Brain Death Diagnosis Evoked Flow**

Continuous Examinations

Cerebral Monitoring Emboli Detection Endarterectomy Evaluation of ICP changes Intraoperative Doppler Intraoperative TCD Monitoring

Discover the *Dolphin/XF* The most advanced and unique bilateral TCD Robot in the market!




The Dolphin/XF TCD Robot Quickly scans the cerebral circulation and automatically assists in identifying blood flow velocities on the right and left sides of the brain. Using the **Dolphin/XF** is fast, simple and stable. Ideal for monitoring in the ICU, OR, and TCD lab, and for a range of specialty tests such as PFO or VMR tests.

Dolphin/XF Selected Features:

- ✓ Comfortable Fabric Headset
- ✓ Automatic Bilateral Scan
- ✓ Auto Search along the Ultrasound Beam
- Position Memory and Maximum Position Accuracy
- ✓ Scan Across all Depths
- "Heat Maps" Identifying Cerebral Vessels
- ✓ Simple and Intuitive User Interface
- ✓ Disposable Headset
- ✓ Fast Scan Option in just a few Seconds
- ✓ Fully Compatible with all Dolphin TCD models

Robotic Bilateral Scan





Monitoring Headset

Range of Accessories

External Channels Connection Box

Wireless / Wired Remote Control



Foot Switch



Highly Advanced Doppler Processing

- Spectrum, M Mode & Sound replay at any depth!
- ✓ Complete **Raw Data** for post processing
- ✓ Multi Depth support
- ✓ Velocity Profile diagnosis

Complete Network Solution
EMR & PACS Connectivity
DICOM Worklist and Storage
Programmable Structured Report

✓ Unlimited **Review Stations**

User Friendly Interface


- ✓ One Key Operation
- ✓ Designed for **Touch Screen**
- Unlimited Flexible Protocols
- ✓ Configurable Display
- ✓ Dynamic Reporting
- Language support

Address Roberts	Addresses No. Addresses Ad
Minused Depter	Brach Holing
	The office of the second se
	An a state of a state of the st
	Lan Ne
	acare which
Image Image <th< td=""><td>and and a second</td></th<>	and and a second
$ \begin{array}{cccccccccccccccccccccccccccccccccccc$	
	B D Destine B to 12

Advanced HITS Detection

- ✓ Fast and thorough review of all **HITS**
- High resolution HITS M-Mode display
- Distribution charts
- Complex & Energy graphs

4D Peak Velocity Sca



Research Oriented

- ✓ Variety of **export** formats including **Raw Data**
- 4D velocity scan

HL7 InterfaceGDT Connectivity

- ✓ Built in **statistics** tool
- ✓ Informative charts and comparisons
- ✓ Data mining

About Us

Viasonix is the only international company that offers a wide range of advanced Peripheral Vascular and Transcranial Doppler (TCD) product lines. The company's state-of-the-art products meet common day-to-day routine clinical practice, as well as set new industry standards.



Dolphin Specifications



The Next Generation Battery Operated TCD System

DOLPHIN4D



The Ultimate TCD Standalone System



The Ideal Cost-Effective TCD Module

Dimensions	47x30x7 cm	47x30x7 cm	27x21x5.5 cm
Weight	6 kg	6 kg	2 kg
Computer	Integrated	Integrated	External
Battery	Internal Rechargeable	No	No

General

Power: 110-220 V, 50-60 Hz Dolphin/MAX Battery Operation Time: Up to 6 hours Doppler Frequencies: 2 MHz, 1.6 MHz, 4 MHz, 8 MHz Intraoperative: 16 MHz (Not for USA & Canada) Probe Connectors: 5 Operating System: Windows 7/8/10/11 (32 or 64 bit) Unilateral/Bilateral: Supported CW/PW Probes: Supported FFT Size Options: 64/128/256/512 Probe Holder: Integrated Controls: Wired/Wireless remote control, Touch screen, Foot switch, Keyboard, Mouse

Dolphin/MAX and Dolphin/4D Integrated Computer

Screen: 18" Wide Screen LED with capacitive Multi Touch screen VESA mount: Built in Operating System: Windows 10 Enterprise, 64 Bit CPU/RAM: Quad Core 2 GHz, 8 GB Hard Drive: 1 TB or higher Enterprise Grade SSD with power loss protection USB Connection: x4, USB 3.0 Network: Ethernet Speakers: Integrated

Dolphin/IQ Recommended Computer

Screen: 21"-24" Wide Screen with Touch Screen CPU/RAM: i3/i5/i7, 4 GB or higher Hard Drive: 1 TB or higher

Doppler Processing

Technology: Fully digital, 8000 channels for each side M Mode: Dynamic up to 500 gates for selected depth range Auto-Depth: Supported Multi-Depth: Up to 8 spectrum depths Replay: Spectrum, M Mode and Doppler audio at any depth Raw Data Post Processing: Depth, Gain, Scale, Sample, Filter Specialty Iests: Breath Holding, VMR, PFO, Sickle Cell, Vasospasm, Cognitive Response Save: Automatic or Manual HITS Detection: Simultaneous in multiple depths HITS M-Mode display: 125 μs resolution HITS Analysis Screen: with detailed events, distribution charts, complex & energy graphs

Continuous Monitoring: Trends display including zoom in/out **Event Markers:** Automatic, Predefined, Alerts, Delta-mean

Software Features

Screen Display: Designed for wide screen and High Definition Cybersecurity: According to FDA Guidelines HIPAA: Patient privacy options **Protocols:** Unlimited, fully configurable Main Screen Layout: Fully configurable, interactive Summary Screens: Interactive User Profile Settings: Unlimited, user configurable **Designed Reports:** Fully configurable Free Text Report: Templates, tables, camera & pictures, TAGS Waveform Labels: Manual and automatic with threshold display Backup/Restore: Integrated. USB/DVD/Network, Manual/Automatic Patient Details: Flexible layout Medical History: Configurable, with indication fields **Storage:** Configurable options, compressed data files Camera: USB camera support Maintenance: Automated protocol with detailed report Languages: Supported Connectivity External Channels: Analog/Digital Input/Output (8/8) Capnograph: Supported Review Stations: Unlimited, including over network Shared Database: Over network **DICOM Worklist:** Offline and Online **DICOM Store:** BMP/JPG/PDF, multiple destinations **DICOM Features**: Auto-send, auto-retry (queue) **DICOM Structured Report:** Default/Programmable HL7: Orders, Report, Query, Patient Registration **GDT:** Patients and Report **PDF:** Automatic Report Export Email: With PDF attachment, using SMTP

Network Queue: Send reports while network is offline

Research Options

Export: PDF, BMP, JPG, PNG, Excel, VDL, Movie, Word, RTF, Raw Data, Anonymous PDF Statistics: Advanced charts, Comparison graphs, X-Y Analysis Advanced Display: Velocity Profile, Depth Scan, 4D Analysis Data Mining: Available with advanced segments Patient History: Quick review of previous exams, patient trends screen Database Search: Multiple advanced options

Selected Approvals

ISO 13485, CE 0344, FDA K202742 Additional country-specific approvals



10 Hamelacha St., POB 2041, Ra'anana, 4366105 Israel