Digital Thermometer 电子体温计

FDA CE 0197 CFDA

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
Note: In step 2, if the room temperature is higher than 32°C (90.0°F), LCD will display the room temperature instead of Lo°C (or Lo°F)

Note: Use of the probe cover may result in a 0.1°C (0.2°F) discrepancy from actual temperature/.

ORAL USE

Put the probe into either position marked by "✓" as showed in the picture right-hand for about 1 minute. To insure the accuracy, please keep mouth close for 5 minutes before taking temperature. The normal temperature by this method is 35.7-37.3°C (96.3°F-99.2°F)

RECTAL USE

Lubricate the probe with the petroleum jelly for easier insertion. Gently insert the probe about one centimeter (less than 1/2 inch.) Into THE Patient's rectum, The normal temperature by this method is 36.5°C-37.6°C (97.7°F-99.7°F)

AXILLARY USE

From a medical viewpoint, this method will always provide inaccurate reading. Wipe armpit with a dry towel. Place probe into the armpit and keep the arm firmly against the body side for about 1.5 minutes. The normal temperature by this method is 35.2°C-37.0°C (95.4°F-98.6°F)

CLEANING AND DISINFECTION

1. Clean thermometer by wiping it with a dry cloth and sterilize the probe with ethyl alcohol
2. Don't let thermometer contact any chemical thinner

Not completely waterproof. Do not soak the LCD in water, Only the probe is waterproof

PRECAUTION

- Do not drop thermometer, it is not shock proof
- Do not bend thermometer
- Do not boil the probe
- Do not store thermometer in direct sunlight
- Do not attempt to disassemble thermometer except replacing the battery
- means direct volts d.c
- means type B

BATTERY REPLACEMENT

1. When the "B" appears in the lower right corner of the LCD, the battery is exhausted and needs replacement.
2. Pull the battery cover off with hand in the direction shown below.
3. Gently pull out the circuit board with battery chamber approximately 1 cm
4. Use pointed object such as a pen to take out the exhausted battery. Discard exhausted battery out of the children.
5. Place a new 1.55V D.C. button size battery type LR41, UCC392 or LR41 or equivalent in the chamber with positive
6. Cover the battery

LIMITED WARRANTY

The thermometer is guaranteed for one year from the date of purchase, if the thermometer does not function properly due to defective components or poor workmanship, we will repair or replace it free of charge. All components are covered by this warranty excluding the battery. The warranty does not cover damages to your thermometer due to improper handling.

DIGITAL THERMOMETER

OPERATING INSTRUCTIONS

PLEASE READ CAREFULLY BEFORE USING

The Digital Thermometer provides quick, highly accurate reading of your temperature. Please read the instructions thoroughly. They are edited to acquaint you with features of your new thermometer so that you can enjoy its functions and years of dependable performance. If the thermometer is used by a child, the supervision of the adults is compulsory. This thermometer is manufactured according to the system of EN12470-3, ASTM and complies with the requirements of IEC 60601-1, IEC 60601-1-2.

PRODUCT DESCRIPTION

SPECIFICATIONS

Range	32°C-42.9°C (90.0°F-109.9°F)
Accuracy:	±0.1°C, 35.5°C-42.0°C (±0.2°F, 95.9°F-107.6°F) ±0.2°C under 35.5°C or over 42.0°C ±0.4°F under 95.9°F or over 107.6°F At Standard room temperature of 25°C (77.0°F)
Display	Liquid crystal display, 3 1/2 digits
Memory	For storing the last measured value
Battery	One 1.55V DC button battery (size SR41, UCC 392 or LR41)
Power consumption	0.15 milliwatts in measurement mode
Battery life	Approx 200 hours
Dimension	12.4 x 1.9 x 1.0 cm (L x W x H)
Weight	Approx. 10 grams including battery
Using environment	Temperature: 10°C-40°C (50°F-104°F) Relative humidity: 15%RH-95%RH
Storage and transportation condition	Temperature: -25°C-55°C (-13°F-131°F) Relative humidity: 15%RH-95%RH
Alarm (Beeper Mode)	Approx 10 seconds sound when peak temperature reached.

HOW TO USE

1. Wipe or dip the probe with ethyl alcohol for disinfection.
2. Push the On/Off button. The LCD will show 188.8° for 2 seconds and you will hear a beep simultaneously. Last tested temperature will be shown for about 2 seconds. After than Lo°C (Lo°F) sign shows (if the room temperature < 32°C (90.0°F)).
3. Position the probe in mouth, rectum, armpit etc.
4. Once the temperature-taking is finished, the degree sign of °C (°F) on the LCD will stop flashing, and you will hear buzzes for approx. 10 seconds. The tested temperature is shown on LCD.
 - When the °C (°F) stop flashing, if the probe is still in mouth, rectum, armpit, the temperature is possible rise too, so after hear the buzzes, you'd better take out it 1-2 minutes later
 - During after the test, if the test temperature is < 32°C or 90.0°F, the LCD shows Lo°C (Lo°F); if the tested temperature is > 42.9°C or 109.9°F, the LCD shows Hi°C (Hi°F). Consult doctor at once.
 - After temperature-taking, if the tested temperature is < 37.8°C or 100.0°F, the buzzes are? Bi--Bi--Bi--Bi-; if the tested temperature is > 37.8°C or 100.0°F, the buzzes whose frequency are much higher are Bi-Bi-Bi-----Bi-Bi-----Bi-Bi-Bi.
5. The thermometer will automatically turn off in about 10 minutes for energy-saving. However, to prolong the battery life, it is suggested to press On/Off button to turn the thermometer off once the test is finished.
6. Please dispose exhausted battery and thermometer according to the relevant law.

BK8005 Non-contact Infrared Forehead Thermometer

- Working condition: Environment temperature: 16° C - 35° C
- Power : DC3V (AA*2)
- Size: 95mm×42mm×139mm(L x W x H)
- Weight: 85g (without batteries)
- Resolution: 0.1° C(or 0.1° F)
- Measuring Range: 35° C - 42° C(95° F - 107.6° F)
- Measuring precision: 35.0° C-42.0° C : ±0.2C,
- Measuring distance: 3 – 5CM
- Memorization: 32 sets
- Storage Environment: -20° C - 50° C , Relative humidity:≤ 85 %
- Packing: 50pcs/Ctn, 13.1KGS with battery(10.7KGS without battery)
- 48*40*36cbm/Ctn

- OEM: available

- Operating Temperature: 15°C ~ 30°C

Mercury-free

Clinical Thermometer


rossmax

monitoring

prevention • monitoring • therapy • professional

Semi-Automatic Blood Pressure Monitor

Model : CH91

- Digital Technology
- Irregular Heartbeat Detection (IHB)
- 90 memories
- One-touch operation
- Latex-free patented cone cuff

Digital
Technology

Irregular Heartbeat
Detection (IHB)

90 memories

One-touch
operation

Latex-free patented
Cone cuff

AAMI

This blood pressure monitor meets the AAMI accuracy standard requirements (AAMI/ANSI/ISO 81060-2)

Cone cuff fits arm circumference:
24 - 36cm (9.4" - 14.2")
Other size cuffs available upon request

Digital Technology

Easy to use manual inflation. This manual inflate blood pressure monitor offers accurate blood pressure and pulse measurements quickly and easily. Simply inflate the cuff to the desired pressure. The cuff will automatically deflate and show your readings on the easy to read display.

Model	Qty per carton	Carton volume
CH91	18 pcs	0.051 cbm/ctn
QTY per 20'	QTY per 40'	QTY per 40'HQ
9882 pcs	19764 pcs	22932 pcs

Storage and Transportation Condition
-10°C~60°C(14°F~140°F)
10%~90% RH, 700~1060 hPa

www.rossmax.com
OBM_Epaper_CH91(0)A04_EP8510_ver1728