



Preventive Maintenance Checklist

Customer: Address:	Equipment type and serial no.: AnyScan DUO/TRIO SCP S/N: AS-XXXXXX-SCP	Date: <div style="text-align: center; font-weight: bold; font-size: 1.2em;">2026</div> <table style="width: 100%; text-align: center; font-size: 0.8em;"> <tr> <td>MM</td> <td>MM</td> <td>MM</td> <td>MM</td> </tr> <tr> <td>DD</td> <td>DD</td> <td>DD</td> <td>DD</td> </tr> </table>	MM	MM	MM	MM	DD	DD	DD	DD
MM	MM	MM	MM							
DD	DD	DD	DD							

SPECT Image Calibration & Detectors

- Fans check	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
- Check and readjustment of PMT amplification factors (Fine Tc-99m Autotuning Calibration)	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
- Isotope peak check and adjustment (Fine Spectrum Calibration) *	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
- Energy and flood-field uniformity calibration and test	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
- Linearity and geometrical resolution check/calibration with linearity phantom *	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
- Linearity and geometrical resolution check with BAR phantom (if available) *	<input type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
- System resolution check with capillaries (System Spatial Resolution without Scatter Test) *	<input type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
- System Sensitivity check *	<input type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
- Detector's crystal check (Detuned Image Test) *	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
- MPH Geometry Test/Calibration *	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
- Activity Test/Calibration (if available) *	<input type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>

SPECT Gantry

- Fans check	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
- Emergency STOP buttons check	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
- Touchscreen and handcontrol check	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
- Motions test between hardware limits *	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
- Detector configuration check in 75°-90°-102°-120°-180° (whichever available)	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
- Extended detector configurations check (Flex option only)	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
- Mechanical noises check	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
- Check the ribbed belts of the drives (tightness, running) *	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
- Detector radius drives lubrication	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
- Detector radius drives wear check and motor current test *	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
- Collimator locks check and lubrication (if needed) *	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
- Collimator contacts check and cleaning	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
- Body Contours test, contacts check and cleaning	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
- Collision protection check and functionality test	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
- Collimator exchange test	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
- Collimator carts, wheels and brakes check *	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
- Center of Rotation offset check and calibration (if needed)	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
- Whole body resolution test with capillaries (WB System Spatial Res. without Scatter Test) *	<input type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>

Imaging table

- Emergency STOP buttons check	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
- Pallet condition check	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
- Vertical motions and limits check *	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
- Pallet motions and limits check *	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
- Rollers and rails check and lubrication *	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
- Ribbed belts check*	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
- Wheels and lock check *	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
- ECG check (optional)	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>

Control unit and control computers

- Emergency STOP buttons check	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
- Intercom and GO button check	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
- Acquisition computer cleaning *	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
- Acquisition computer fans check	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
- Acquisition computer disks maintenance	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>
- DVD drive check (if applicable)	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
- Network performance check *	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
- Processing workstation(s) fans check and cleaning *	<input type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
- Processing workstation(s) disks maintenance *	<input type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>

* yearly maintenance item

CT unit

- Fans check, cleaning filters	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
- PDU fans check	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
- Door switches check	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
- Radiation warning lamps check	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
- Laser check	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
- Mechanical noises check	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
- Slip ring, contact brushes and tick fence cleaning *	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
- Detector air filter cleaning (V16 detector only) *	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
- HV cable connectors cleaning and re-lubrication *	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
- X-Ray plexi strip cleaning (as needed)	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
- Air calibration	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
- Daily CT QC (uniformity and resolution check)	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
- Vibration check (at 60RPM) and balancing (0°, 90°, 180°, 270°) *	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
- Modality registration test (and calibration if necessary)	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
- Registration of usage statistics (HVC, GMC, HE)	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>

PET unit

- Fans check	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
- Gain calibration	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
- Modality registration test (and calibration if necessary)	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
- Timeshift calibration and test	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
- Normalization	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
- Activity calibration and test	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
- Daily PET QC	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
- NEMA IQ phantom test	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>

* yearly maintenance item

Notes

PM #1		
GMC counter:	DR1 wear (<0,8 mm):	DR1 current (↓<8 A): (↑<12 A):
HVC counter:	DR2 wear: (<0,8 mm):	DR2 current (↓<8 A): (↑<12 A):
HE counter*:	DR3 wear: (<0,8 mm):	DR3 current (↓<8 A): (↑<12 A):
Service Form ID:	Date:	Service Engineer:
PM #2		
GMC counter:	DR1 wear (<0,8 mm):	DR1 current (↓<8 A): (↑<12 A):
HVC counter:	DR2 wear: (<0,8 mm):	DR2 current (↓<8 A): (↑<12 A):
HE counter*:	DR3 wear: (<0,8 mm):	DR3 current (↓<8 A): (↑<12 A):
Service Form ID:	Date:	Service Engineer:
PM #3		
GMC counter:	DR1 wear (<0,8 mm):	DR1 current (↓<8 A): (↑<12 A):
HVC counter:	DR2 wear: (<0,8 mm):	DR2 current (↓<8 A): (↑<12 A):
HE counter*:	DR3 wear: (<0,8 mm):	DR3 current (↓<8 A): (↑<12 A):
Service Form ID:	Date:	Service Engineer:
PM #4		
GMC counter:	DR1 wear (<0,8 mm):	DR1 current (↓<8 A): (↑<12 A):
HVC counter:	DR2 wear: (<0,8 mm):	DR2 current (↓<8 A): (↑<12 A):
HE counter*:	DR3 wear: (<0,8 mm):	DR3 current (↓<8 A): (↑<12 A):
Service Form ID:	Date:	Service Engineer: