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THE *PERFECT* ESR SYSTEM


RHEOLOGY TECHNOLOGY

The iSED® is the first in the line of Erythrocyte Sedimentation Rate (ESR) analyzers from ALCOR Scientific Inc.

The iSED® is a 20 position fully automated ESR analyzer that works directly from the primary EDTA tube and requires just 100 microliters of sample to deliver a result in just 20 seconds.

The iSED® maximizes laboratory efficiency with a continuous, random access feed and high throughput of 180 sample per hour.

ALCOR Scientific utilizes photometric reading to quantify the Rouleaux formation, which is the earliest and most critical phase of sedimentation. Our analyzers' micro-flow cell captures the critical kinetics of RBCs in a highly controlled testing environment to produce ESR results in 15-20 seconds. The analyzers' results are unaffected by variables commonly associated with traditional ESR testing methods. Some of these variables include hematocrit, MCV, and temperature.

STAT ESR RESULTS

Fast, accurate and reliable results in just 20 seconds.

MINIMAL SAMPLE SIZE

Only 100µL of sample is required for testing.

FULLY AUTOMATED

Barcode reader, mixer, and printer are all internal.

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BR-66953

DECLARATION OF CONFORMITY

Beckman Coulter Ireland, Inc. hereby ensures and declares that the product(s) listed below conform to the European Union directives and standards identified in this declaration.

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Conformity Assessment Procedure
Annex III Self-Declared

CE MARK affixed: March 2018

Date: March 30, 2018

Marie O'Rourke, Group Q&RA Manager
Beckman Coulter Ireland Inc.

Beckman Coulter Ireland Inc.
Lismeehan
O'Callaghan's Mills
Co. Clare, Ireland

Document Control

Issue Date: March 20, 2018
Revision Lev: 1.0
Revision Date: N/A
Starting Serial #: BB030001
RoHS Serial #: BB030001
FileName: 2018-901-004I

Document:

Declaration of Conformity

Document #:

GLB-QS-TMP-0029

Declaration of Conformity

Beckman Coulter Ireland, Inc. hereby ensures and declares that the product(s) listed below comply with the requirement of the European Union In-vitro Diagnostics Medical Device Directive 98/79/EC.

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Product(s) /Produkt(e) /Prodotto(i) / Produit(s) / Producto(s):

PRODUCT NAME	PART NUMBER	GMDN
DxH 500 Series Calibrator	B36880	55865

Conformity Assessment Procedure
Annex III –Self-declared

Classification:
General

Marie O'Rourke Group Q&RA Manager

1st Nov. 2016
Date

Beckman Coulter Ireland Inc.
Lismeehan
O'Callaghan's Mills
Co. Clare, Ireland

Document Control

Issue Date: November 2, 2015
Revision Level: 2.0
Revision Date: November 1, 2016
Starting Lot #: 491608600
Filename: B36880DEC

Declaration of Conformity

Beckman Coulter, Inc. hereby ensures and declares that the product(s) listed below comply with the requirement of the European Union In-vitro Diagnostics Medical Device Directive 98/79/EC.

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Conformity Assessment Procedure
Annex III - Self-Declared

Classification:
General

GMDN Code(s):
59058

Nery Ortiz

Sr. Manager Regulatory Affairs

13 Nov 2018
Date

Beckman Coulter, Inc.
250 S. Kraemer Blvd.
Brea, CA 92821 USA

Document Control
Issue Date: November 2, 2015
Revision Level: 3.0
Revision Date: October 8, 2018
Starting Lot #: 8310013
Filename: B36868DEC

Declaration of Conformity

Beckman Coulter, Inc. hereby ensures and declares that the product(s) listed below comply with the requirement of the European Union In-vitro Diagnostics Medical Device Directive 98/79/EC.

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
Conformity Assessment Procedure
Annex III - Self-Declared

Classification:
General

GMDN Code(s):
55866

Nery Ortiz
Sr. Manager Regulatory Affairs

13 Nov 2018
Date

Beckman Coulter, Inc.
250 S. Kraemer Blvd.
Brea, CA 92821 USA

Document Control
Issue Date: November 2, 2015
Revision Level: 3.0
Revision Date: October 8, 2018
Starting Lot #: 35180861, 36180861,
37180861
Filename: B36872DEC

Document:

Declaration of Conformity

Document #:

GLB-QS-TMP-0029

Declaration of Conformity

Beckman Coulter Ireland, Inc. hereby ensures and declares that the product(s) listed below comply with the requirement of the European Union In-vitro Diagnostics Medical Device Directive 98/79/EC.

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Product(s) /Produkt(e) /Prodotto(i) / Produit(s) / Producto(s):

PRODUCT NAME	PART NUMBER	GMDN
DxH 500 Series Diluent	B36845	55855

Conformity Assessment Procedure
Annex III –Self-declared

Classification:
General

Marie O'Rourke Group Q&RA Manager

19 Nov 2016
Date

Beckman Coulter Ireland Inc.
Lismeehan
O'Callaghan's Mills
Co. Clare, Ireland

Document Control
Issue Date: November 2, 2015
Revision Level: 2.0
Revision Date: November 1, 2016
Starting Lot #:0016170
Filename: B36845DEC

Document:

Declaration of Conformity

Document #:

GLB-QS-TMP-0029

Declaration of Conformity

Beckman Coulter Ireland, Inc. hereby ensures and declares that the product(s) listed below comply with the requirement of the European Union In-vitro Diagnostics Medical Device Directive 98/79/EC.

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Product(s) /Produkt(e) /Prodotto(i) / Produit(s) / Producto(s):

PRODUCT NAME	PART NUMBER	GMDN
DxH 500 Series Lyse	B36846	55855

Marie O'Rourke Group Q&RA Manager

Conformity Assessment Procedure
Annex III –Self-declared

Classification:
General

15 Nov 2016

Date

Beckman Coulter Ireland Inc.
Lismeehan
O'Callaghan's Mills
Co. Clare, Ireland

Document Control

Issue Date: November 2, 2015
Revision Level: 2.0
Revision Date: November 1, 2016
Starting Lot #:8300007
Filename: B36846DEC

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INTENDED USE

For use as an isotonic buffered diluent used in conjunction with a cyanide-free lytic agent for counting and sizing blood cells on the DxH 500 Series Systems. Use only with DxH 500 Series Lyse reagent.

INTENDED USER

This device is intended to be used by a laboratory professional.

SUMMARY

Blood cell analysis comprises mixing a portion of a whole-blood sample with a diluent that provides the ability to analyze red blood cells and platelets. The diluent is used in combination with a lytic reagent for the determination of hemoglobin, the enumeration and differentiation of leukocytes (white blood cells).

PRINCIPLE

DxH 500 Series Diluent is an osmotically balanced, pH neutral solution of organic buffers, anesthetics, germicides, and sodium salts. The diluent acts as an electrolyte and cell volume stabilizer for the counting of red blood cells and platelets. The diluent is used with DxH 500 Series Lytic reagent for the enumeration and differentiation of leukocytes subpopulations, and hemoglobin measurement.

REACTIVE INGREDIENTS

Sodium Salts.....<15g/L

NON-REACTIVE INGREDIENTS

Buffer, Stabilizers, Preservatives

WARNING AND PRECAUTIONS

- Do not inhale and/or ingest.
- Avoid eye contact. In case of eye contact, flush eyes with copious amounts of water for several minutes.
- Avoid skin contact. In case of skin contact, wash area with soap and water. Rinse well with water.
- DO NOT REUSE CONTAINER.

GHS HAZARD CLASSIFICATION		
DxH 500 Series Diluent		EUH208: May produce an allergic reaction.
		Reaction Mass of: 5-chloro-2-methyl-4-isothiazolin-3-one [EC# 247-500-7] and 2-methyl-4-isothiazolin-3-one [EC# 220-239-6](3:1) <0.0015% w/w

STORAGE, STABILITY, AND DISPOSAL

- Store DxH 500 Series Diluent at 2-30°C.
- Do not use the product past the expiration date.
- Use the product at temperatures stated in the Instructions for Use for your instrument.
- Discard the opened container after 60 days.
- Dispose of waste product, unused product, and contaminated packaging in compliance with federal, state, and local regulations.

MATERIALS PROVIDED

DxH 500 Series Diluent

MATERIALS REQUIRED, BUT NOT PROVIDED

- Sufficient DxH 500 Series reagents to complete analysis. For the reagents used on the system refer to the specific DxH 500 Series Hematology Analyzer IFU.
- System HELP or Instructions for Use.

PREPARATION

DxH 500 Series Diluent is ready to use.

IMPORTANT: If product has been partially or completely frozen, allow product to warm to room temperature. Invert the container 16 times to ensure complete mixing prior to placement on the instrument. Install and prime the diluent as directed in the Instructions for Use for your instrument. Verify background counts are acceptable before analyzing patient samples.

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ADDITIONAL INFORMATION

Glossary of Symbols is available at beckmancoulter.com/techdocs by searching for part number C07163.

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SDS

Safety Data Sheet is available at beckmancoulter.com/techdocs

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NOTICE TO USER

For a patient/user/third party in the European Union and in countries with identical regulatory regime (*Regulation 2017/746/EU on In vitro Diagnostic Medical Devices*); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

Revision History

Revision AA, 03/2022

Initial release to include documentation updates for compliance to European In Vitro Diagnostic Regulations (IVDR) (2017/746). No changes to product or instructions for use.

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INTENDED USE

Erythrocyte lytic agent for quantitative determination of hemoglobin, and the counting, sizing and optical measurement of leukocytes on DxH 500 Series Systems. Use only with DxH 500 Series Diluent and DxH 500 Series Cleaner.

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SUMMARY

DxH 500 Series Lyse is a cyanide-free lytic agent used to prepare a diluted whole-blood sample for the quantitative determination of hemoglobin and the enumeration and sub-classification of leukocytes.

PRINCIPLE

DxH 500 Series Lyse is a cyanide-free reagent.¹ Surfactants are used to lyse the erythrocytes and release the hemoglobin. The released hemoglobin reacts with a non-toxic chelating agent to create a stable oxyhemoglobin hemochromogen which can be quantitatively measured.

The lytic reagent solubilizes the red blood cell membrane while simultaneously preserving the white blood cell morphology to be counted, sized and optically measured.

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WARNING AND PRECAUTIONS

- Product contains Glutaraldehyde <0.1%/L
- May cause skin and eyes irritation.
- In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.
- Do not inhale and/or ingest.
- DO NOT REUSE CONTAINER.

GHS HAZARD CLASSIFICATION		
DxH 500 Series Lyse		EUH208: May produce an allergic reaction.
		Glutaraldehyde <0.1% w/w
		Reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one [EC# 247-500-7] and 2-methyl-4-isothiazolin-3-one [EC# 220-239-6](3:1) < 0.0015% w/w

STORAGE, STABILITY, AND DISPOSAL

- Store DxH 500 Series Lyse at 4-25°C.
- Do not use product past expiration date.
- Use product at temperatures stated in the Instructions for Use for your instrument.
- Opened containers are stable for 60 days when stored as recommended.
- Dispose of waste product, unused product, and contaminated packaging in compliance with federal, state, and local regulations.

MATERIALS PROVIDED

DxH 500 Series Lyse

MATERIALS REQUIRED BUT NOT PROVIDED

- Sufficient COULTER reagents to complete analysis. For the reagents used on the system refer to the DxH 500 Series Hematology Analyzers IFUs
- System HELP or Instructions for Use.

PREPARATION

DxH 500 Series Lyse is ready to use.

IMPORTANT: If product has been partially or completely frozen, allow product to warm to room temperature. Mix product by gentle inversion prior to placement on the instrument. Install and prime the reagent as directed in the Instructions for Use for your instrument.

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REFERENCES

1. Malin, Michael J., et al., "Evaluation of 24-Second Cyanide-Containing and Cyanide-Free Methods for Whole Blood Hemoglobin on the Technicon H* 1 Analyzer with Normal and Abnormal Blood Samples," *AJCP*, Sep. 1989 pp. 286-294.

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www.beckmancoulter.com

Manufactured for Beckman Coulter, Inc.

DxH 500 Series Cleaner

REF B36868

TABLE OF CONTENTS

English	2
Français (FR)	4
Deutsch (DE)	6
Italiano (IT)	8
Español (ES)	10
Português Portugal (PT-PT)	12
Dansk (DA)	14
Svenska (SV)	16
Norsk (NO)	18
Ελληνικά (EL)	20
日本語 (JP)	22
中文 (ZH-CN)	24
Lietuviškai (LT)	26
Magyar (HU)	28
Polski (PL)	30
Čeština (CZ)	32
Slovenčina (SK)	34
한국어 (KO)	36
Türkçe (TR)	38
Русский (RU)	40
Hrvatski (HR)	42
Български (BG)	44
中文 (ZH-TW)	46
Română (RO)	48
Србија (SR)	50
Українська (UK)	52
Português Brasil (PT-BR)	54
Nederlands (NL)	56
Tiếng Việt (VI)	58
REFERENCES	60

DxH 500 Series Cleaner

REF B36868

For *In Vitro* Diagnostic Use
Rx Only in the U.S.A.

INTENDED USE

A cleaning agent for use on DxH 500 Series Systems for components that come in contact with blood samples.

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SUMMARY

The DxH 500 Series Cleaner is an azide-free, formaldehyde-free, biodegradable cleaner which contains a proteolytic enzyme that aids in the removal of protein buildup.

PRINCIPLE

Proteins that flow through apertures may accumulate, and narrow or block aperture openings. The cleaning reagent dissolves proteins at apertures as well as at aspiration and dispensing components of the systems.¹

REACTIVE INGREDIENTS

A solution of a proteolytic enzyme (0.1-0.5% by wt.)

NON-REACTIVE INGREDIENTS

Buffer, stabilizers, preservatives, surfactant

WARNING AND PRECAUTIONS

- Do not inhale and/or ingest.
- Avoid eye and skin contact. In case of eye or skin contact, flush affected area with copious amounts of water for at least 15 minutes.
- DO NOT REUSE CONTAINERS.

GHS HAZARD CLASSIFICATION		
DxH 500 Series Cleaner	DANGER 	May cause allergy or asthma symptoms or breathing difficulties if inhaled.
		Avoid breathing vapours.
		In case of inadequate ventilation, wear respiratory protection.
		IF INHALED: Remove person to fresh air and keep at rest in a position comfortable for breathing.
		If experiencing respiratory symptoms: Call a POISON CENTER or doctor/physician.
		Subtilisin 0.1-0.5%

STORAGE, STABILITY, AND DISPOSAL

- Store DxH 500 Series Cleaner at 2-25°C.
- Do not use product past expiration date.
- Use product at temperatures stated in the Instructions for Use for your instrument.
- Opened containers are stable for 3 months when stored as recommended.
- Dispose of waste product, unused product, and contaminated packaging in compliance with federal, state, and local regulations.

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Phiên bản AA, 04/2022

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REFERENCES

1. Lin Y, Means GE, and Feeney RE. 1969. Action of proteolytic enzymes on N,N-dimethyl proteins. Basis for a microassay for proteolytic enzymes. *J Biol Chem*, 244(3):789-793.

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DxH 500 Series Control

REF B36872

**For *In Vitro* Diagnostic Use
Rx Only in the U.S.A.**

INTENDED USE

DxH 500 Series Control is an assayed whole blood control designed to monitor values on multi parameter hematology cell counters. Refer to the Table of Expected Results for specific instrument models.

SUMMARY AND PRINCIPLE

It is an established laboratory practice to use a stable control to monitor the performance of diagnostic tests. This control is composed of stable materials that provide a means of monitoring the performance of hematology blood cell counters. It is sampled in the same manner as a patient specimen.

REAGENTS

DxH 500 Series Control is composed of human erythrocytes, mammalian leukocytes and mammalian platelets suspended in a plasma-like fluid with preservatives.

WARNING AND PRECAUTIONS

POTENTIAL BIOHAZARDOUS MATERIAL. Each human donor/unit used in the preparation of this product has been tested by a FDA licensed method/test and found to be negative or non-reactive for the presence of HBsAg, Anti-HCV, NAT testing for HIV-1, HCV (RNA) and HIV-1/2. Each unit is also negative by a serological test for Syphilis (RPR or STS). Because no test method can offer complete assurance that infectious agents are absent, this material should be handled as potentially infectious. When handling or disposing of vials follow precautions for patient specimens as specified in the OSHA Bloodborne Pathogen Rule (29 CFR Part 1910, 1030) or other equivalent biosafety procedures. DxH 500 Series Control is intended for use only by trained personnel.

GHS HAZARD CLASSIFICATION		
DxH 500 Series Control	WARNING	Causes skin irritation.
		Causes serious eye irritation.
		Wear protective gloves, protective clothing and eye/face protection.
		IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
		If skin irritation occurs: Get medical advice/attention.
		If eye irritation persists: Get medical advice/attention.
		Potassium Hydroxide <1%

STORAGE, STABILITY, AND DISPOSAL

Store DxH 500 Series Control upright at 2-8°C when not in use. **Protect tubes from overheating and freezing.** Unopened tubes are stable until the expiration date, shown on the TABLE OF EXPECTED RESULTS. For open vial stability, refer to the TABLE OF EXPECTED RESULTS for your system. Dispose of waste product, unused product and contaminated packaging in compliance with federal, state and local regulations.

EVIDENCE OF DETERIORATION

After mixing, product should be similar in appearance to fresh whole blood. In unmixed tubes, the supernatant may appear cloudy and reddish; this is normal and does not indicate deterioration. Other discoloration, very dark red supernatant or unacceptable results may indicate deterioration. **Do not use the product if deterioration is suspected.**

INSTRUCTIONS FOR USE

1. Remove tubes from the refrigerator and verify that the lot number on the tubes matches the lot numbers on the Table of Expected Results.
2. Allow the tubes to warm to room temperature (15-30°C or 59-86°F) for 15 minutes before mixing.
3. To mix, hold a tube horizontally between the palms of the hands. **Do not pre-mix on a mechanical mixer.**
 - A. Roll the tube back and forth for 20-30 seconds; occasionally invert the tube.
 - B. Continue to mix in this manner until the red cells are completely suspended. Tubes stored for a long time may require extra mixing.

- C. Gently invert the tube 8-10 times immediately before sampling.
4. Set Up and analyze the controls according to the procedure in the Instructions for Use for your instrument.
5. After sampling:
 - A. Clean residual material from the cap and rim if the tube is opened for sampling. Replace the cap tightly.
 - B. Return tubes to refrigerator within 30 minutes of use.

ASSIGNED VALUES AND EXPECTED RESULTS

Reagent differences, maintenance, operating technique, and calibration may contribute to inter-laboratory variation. Refer to the Table of Expected Results for the corresponding lot number of the control tube.

PERFORMANCE LIMITS

The performance of this product is assured only if it is properly stored and used as described in this insert. Incomplete mixing of a tube prior to use invalidates both the sample withdrawn and any remaining material in the tube.

PERFORMANCE CHARACTERISTICS

Assigned values are presented as a Mean and Range. The Mean is derived from replicate testing on instruments operated and maintained according to the manufacturer's instructions. The Range is an estimate of variation between laboratories and also takes into account inherent imprecision of the method and expected biological variability of the control material.

Assay values on a new lot of control should be confirmed before the new lot is put into routine use. Test the new lot when the instrument is in good working order and quality control results on the old lot are acceptable. The laboratory's recovered mean should be within the assay range.

For greater control sensitivity each laboratory should establish its own mean and acceptable range and periodically reevaluate the mean. The laboratory range may include values outside of the assay range. The user may establish assay values not listed on the Assay Sheet, if the control is suitable for the method.

TROUBLESHOOTING PROCEDURE

The Instructions for Use of the DxH 500 Series instrument contain troubleshooting procedures for resolving control recovery problems. If you need additional help, call Beckman Coulter Customer Service at 800-526-7694 (USA and Canada) or contact your local Beckman Coulter Representative.

To provide faster handling of your inquiry, please have the following information available when you call:

- Instrument model and serial number.
- Expiration dates, lot numbers, and supporting data for all reagents, for control(s) in question and for other control lots currently in use.
- Previous control lot numbers and the data for these previous lots.
- Data from current repeatability study (N=10) using a fresh whole blood specimen and performed according to the Instructions for Use.
- Data from the IQAP and data from your last instrument calibration where appropriate.

PRODUCT AVAILABILITY

DxH 500 Series Control

REF B36872 - 2 x 2.3 mL Abnormal Low, 2 x 2.3 mL Normal, 2 x 2.3 mL Abnormal High

TRADEMARKS

Beckman Coulter, the stylized logo, and the Beckman Coulter product and service marks mentioned herein are trademarks or registered trademarks of Beckman Coulter, Inc. in the United States and other countries.

For additional information or if damaged product is received, call Beckman Coulter Customer Service at 800-526-7694 (USA or Canada) or contact your local Beckman Coulter Representative.

ADDITIONAL INFORMATION

Glossary of Symbols is available at beckmancoulter.com/techdocs (PN C07163)

May be covered by one or more pat. - see www.beckmancoulter.com/patents

SDS Safety Data Sheet is available at beckmancoulter.com/techdocs

Revision History

Revision AH, 11/2018

Changes were made to:

- PRODUCT AVAILABILITY - Removed Product Part Number C21335
- REVISION HISTORY - Removed Revision AE

Revision AJ, 12/2018

Changes were made to:

- Added new Patent Information Statement
- Added new Dutch language

Revision AK, 10/2021

- Updated STORAGE, STABILITY, AND DISPOSAL section
- Updated website address
- Updated Legal Manufacturer phone number
- Updated copyright year

Revision AL, 02/2023

- Updated STORAGE, STABILITY, AND DISPOSAL section
-

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QUALITY MANAGEMENT SYSTEM - ISO 13485:2016

This is to certify that:

Beckman Coulter, Inc.
250 South Kraemer Blvd.
Brea
California
92821
USA

Facility ID Number: F000075

Holds Certificate No:

MDSAP 678344

The company listed on this certificate has been audited to and found to conform with ISO 13485:2016 including the following country specific requirements:

Australia: Therapeutic Goods (Medical Devices) Regulations, 2002, Schedule 3 Part 1 (excluding Part 1.6) - Full Quality Assurance Procedure

Brazil: RDC ANVISA n. 67/2009, RDC ANVISA n. 665/2022 - Good Manufacturing Practices, RDC ANVISA n. 551/2021


Canada: Medical Devices Regulations - Part 1 - SOR 98/282

Japan: MHLW MO No 169 (2004), as amended by MHLW MO No 60 (2021), PMD Act

USA: 21 CFR 820, 21 CFR 803, 21 CFR 806, 21 CFR 807 - Subparts A to D

The design and development, manufacture, distribution, installation, and service of automated instruments, software, reagents, calibrators, controls and other consumables for use in in-vitro diagnostic medical devices.

For and on behalf of BSI:

Graeme Tunbridge, Senior Vice President Medical Devices

Original Registration Date: 2018-12-18

Effective Date: 2024-03-13

Expiry Date: 2027-03-12

BSI Group America Inc. is an MDSAP recognised auditing organization

Page: 1 of 3

...making excellence a habit.™

Certificate No: **MDSAP 678344**

Location	Registered Activities
Beckman Coulter, Inc. 250 South Kraemer Blvd. Brea California 92821 USA Facility ID Number: F000075	Design and development, manufacture, distribution, installation, service and customer support.
Beckman Coulter, Inc. 1000 Lake Hazeltine Drive Chaska Minnesota 55318 USA Facility ID Number: F001810	Design and development, manufacture, distribution, installation, service and customer support.
Beckman Coulter, Inc. 2470 Faraday Avenue Carlsbad California 92010 USA Facility ID Number: F001808	Design and development, manufacture and distribution.
Beckman Coulter, Inc. 2295 Progress Drive Hebron Kentucky 41048 USA Facility ID Number: F001805	Warehouse and Distribution.
Beckman Coulter, Inc. 11800 SW 147th Avenue Miami Florida 33196 USA Facility ID Number: F001809	Design and development, manufacture, distribution, installation, service and customer support.

Original Registration Date: 2018-12-18

Effective Date: 2024-03-13

Expiry Date: 2027-03-12

Page: 2 of 3

This certificate remains the property of BSI and shall be returned immediately upon request.

An electronic certificate can be authenticated [online](https://www.bsigroup.com/ClientDirectory). Printed copies can be validated at www.bsigroup.com/ClientDirectory

To be read in conjunction with the scope above or the attached appendix.

Certificate No: **MDSAP 678344**

Location	Registered Activities
Beckman Coulter, Inc. 740 West 83rd Street Hialeah Florida 33014 USA Facility ID Number: F001811	Manufacture and Distribution.
Beckman Coulter Inc. 7381 Empire Drive Florence Kentucky 41042 USA Facility ID Number: F001806	Manufacture
Beckman Coulter GmbH Europark Fichtenhain B13 47807 Krefeld Germany Facility ID Number: F003965	Manufacture
Beckman Coulter, Inc. Microbiology Business 2040 Enterprise Blvd. West Sacramento California 95691 USA Facility ID Number: F001807	Manufacture and distribution.
Beckman Coulter, Inc. Microbiology Business 1584 Enterprise Blvd. West Sacramento California 95691 USA Facility ID Number: F001807	Design and development, installation, service, customer support and software.

Original Registration Date: 2018-12-18 Effective Date: 2024-03-13 Expiry Date: 2027-03-12

This certificate remains the property of BSI and shall be returned immediately upon request.
An electronic certificate can be authenticated [online](#). Printed copies can be validated at [www.bsigroup.com/ClientDirectory](#)
To be read in conjunction with the scope above or the attached appendix.

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- De asemenea: În Luke RG, Koepke JA; Siegel RR. The effect of immunosuppressive drugs and uremia on automated leukocyte counts. *Am J Clin Path*, 1971; 56(4), 503-507.
18. International Committee for Standardization in Haematology. Recommendations for reference method for haemoglobinometry in human blood (ICSH Standard EP6/2:1977) and specifications for international haemoglobinocyanide reference preparation (ICSH Standard EP6/3:1977). *J Clin Path*, 1978; 31(2):139-143.
 19. Dorsey, DB. What can quality control do for hematology? *Am J Med Tech. American Society of Medical Technologists*; ISSN: 0002-9335. PMID: 14267574. 1965; 31; Mar-Apr:150-153.
 20. Bull, BS. A statistical approach to quality control. Simpozionul Internațional pentru Standardizarea Metodelor Hematologice. Lewis S; Coster J, ed., *Quality Control in Haematology*, London, England:Academic Press (@H), 1975.
 21. Koepke, JA; Protector, TJ. Quality assurance for multichannel hematology instruments: four years' experience with patient mean erythrocyte indices. *Am J Clin Path*, 1981; 75(1):28-33.
 22. EP06-A Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach – approved guideline; 4/1/2003. Clinical and Laboratory Standards Institute (CLSI), Wayne, PA.
 23. EP09-A3 Measurement procedure comparison and bias estimation using patient samples, approved guideline – third edition; 8/30/2013. Clinical and Laboratory Standards Institute (CLSI), Wayne, PA.
 24. EP28-A3c Defining, establishing, and verifying reference intervals in the clinical laboratory, approved guideline – third edition; 10/19/2010. Clinical and Laboratory Standards Institute (CLSI), Wayne, PA.
 25. GP41-A6 (formerly H3-A6) Procedures for the collection of diagnostic blood specimens by venipuncture, approved standard – sixth edition; 10/31/2007. Clinical and Laboratory Standards Institute (CLSI), Wayne, PA.
 26. GP42-A6 (formerly H4-A5) Procedures and devices for the collection of diagnostic capillary blood specimens, approved standard – sixth edition; 9/23/2008. Clinical and Laboratory Standards Institute (CLSI), Wayne, PA.
 27. GP44-A4 Procedures for the handling and processing of blood specimens for common laboratory tests; approved guideline – fourth edition; 5/25/2010. Clinical and Laboratory Standards Institute (CLSI), Wayne, PA.
 28. H07-A3 Procedure for determining packed cell volume by the microhematocrit method, approved standard – third edition; 10/1/2000. Clinical and Laboratory Standards Institute (CLSI), Wayne, PA.
 29. H15-A3 Reference and selected procedures for the quantitative determination of hemoglobin in blood, approved standard – third edition; 12/1/2000. Clinical and Laboratory Standards Institute (CLSI), Wayne, PA.
 30. H20-A2 Reference leukocyte (WBC) differential count (proportional) and evaluation of instrumental methods – approved standard – second edition; 1/18/2007. Clinical and Laboratory Standards Institute (CLSI), Wayne, PA.
 31. H26-A2 Validation, Verification, and Quality Assurance of Automated Hematology Analyzers – approved standard – second edition; 6/25/2010. Clinical and Laboratory Standards Institute (CLSI), Wayne, PA.

- 32.** Annex 6, How to make chlorine solutions for environmental disinfection. Infection, prevention, and control guidance for care of patients in health-care settings, with focus on Ebola, WHO/HIS/SDS/2014.4 Rev.1, December 2014. World Health Organization, Geneva, Switzerland.
- 33.** C56-A Hemolysis, icterus, and lipemia/turbidity indices as indicators of interference in clinical laboratory analysis – approved guideline; 7/13/2012. Clinical and Laboratory Standards Institute (CLSI), Wayne, PA.
- 34.** EP07-A2 Interference testing in clinical chemistry – approved guideline – second edition; 11/23/2005. Clinical and Laboratory Standards Institute (CLSI), Wayne, PA.
- 35.** Turgeon, M. Clinical Hematology: Theory and Procedures (4th ed.). Baltimore, MD: Lippincott, Williams, & Wilkins; 2005.
- 36.** Zandecki, M et al. Spurious counts and spurious results on haematology analysers: a review. part 1: platelets. International Journal of Laboratory Hematology, 2007; 29:4-20.
- 37.** Zandecki, M. et al. Spurious counts and spurious results on haematology analysers: a review. Part II: white blood cells, red blood cells, haemoglobin, red cell indices and reticulocytes. International Journal of Laboratory Hematology, 2007; 29:21-41.
- 38.** EP05-A3 Evaluation of precision of quantitative measurement procedures – approved guideline – third edition; 10/29/2014. Clinical and Laboratory Standards Institute (CLSI), Wayne, PA.

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www.beckmancoulter.com

Manufacturer	ALCOR Scientific LLC 20 Thurber Boulevard Smithfield, RI 02917, USA Telephone 001 401 737 3774	
SRN	US-MF-000029295	
Device Family Name:	iSED®, Erythrocyte Sedimentation Rate Analyzer	
Device Name (Model):	iSED ESR Analyzer	
Catalog Number / Basic UDI-DI:	112-00101 (for analyzers > SN 05000) / 0896406002iSED2Elite5D	
Intended Use:	iSED®/iSED ELITE Automated Erythrocyte Sedimentation Rate Analyzer is an automated in vitro diagnostic (IVD) device for the determination of an erythrocyte sedimentation rate (ESR) expressed in mm/hr. Testing is performed using EDTA whole blood samples obtained by venipuncture or capillary blood collection. The analyzer is intended to be used in a professional clinical laboratory setting. The analyzer directly measures the aggregation of red blood cells by photometric rheology technology which does not require the use of reagents. Results are reported in units of mm/hr. and correlate with the Westergren method of ESR determination. The quantitative results for sedimentation rate produced by the analyzer are considered non-specific and are used by a clinician to aid in assessing the general health status of a patient. Results from the device are to be used in conjunction with other laboratory results and with the patient condition known by the ordering clinician.	
Device Class (Rule):	IVDR Class A per Rule 5b	
EMDN:	W02029001 (Erythrocyte Sedimentation Rate Devices)	
GMDN:	56691 (Automated ESR Analyzer)	
Accessories:	iWASH™ (112-12-001), GMDN = 59058, EMDN = V0799 iWaste Containers (112-12-002), GMDN = 62172, EMDN = W050301029099 Seditrol® (DSC01, DSC06), GMDN = 55972, EMDN = W0103010599 EQC02 (EQC02), GMDN = 55972, EMDN = W0103010599 Test Cards (112-00250, 112-00500, 112-01000, 112-02000), GMDN = 62804, EMDN = N/A Test Cards (112-03000, 112-05000, 112-10000, 112-20000), GMDN = 62804, EMDN = N/A Thermal Printer Paper (DS-05233), GMDN = 65796, EMDN = N/A deepCLEAN (112-12-020), GMDN = 59058, EMDN = W0599	
Notified Body:	Not Applicable	
EU Authorized Representative:	Emergo Europe (EUDAMED Actor ID/SRN: NL-AR-000000116) Westervoortsedijk 60, 6827 AT Arnhem, The Netherlands Tel: (31) (0) 70 345-8570	
CH Authorized Representative:	MedEnvoy Switzerland Gotthardstrasse 28 6302 Zug, Switzerland Tel: (41) 41 562 01 42	

Applicable Standards:

EN ISO 13485: 2016 w/ MDSAP	Quality Management System	EN 61000-3-2	Harmonics
EU 2017/746	In Vitro Diagnostics Regulation (IVDR) 2017/746	EN 61000-3-3	Flicker
EN ISO 14971:2019	Application of risk management to medical devices	EN 61000-4-2:2009	Electrostatic Discharge
2011/65/EU	Restriction of Hazardous Substances in Electrical & Electronic Eq (RoHS 2)	EN 61000-4-3:2006 /A1/A2	Radiated Immunity
IEC 62304:2006 + AMD1: 2015	Medical Device Software and Life Cycle Processes	EN 61000-4-4:2012	Electrical Fast Transient Burst
IEC 61010-1:2010 EN 61010-1:2020	Safety requirements for electrical equipment for measure, control and laboratory use - General Requirements	EN 61000-4-5:2014/ A1:2017	Surge immunity
EN 55011:2016	Industrial, Scientific and Medical Equipment – Radio-frequency Disturbance Characteristics (FCC Part 15, CISPR 11, ICES-003)	EN 61000-4-6:2014	Conducted immunity
EN 61010-2-101:2017	Safety Requirement for electrical equipment for measurement, control, and laboratory use for in vitro diagnostic (IVD) medical equipment	EN 61000-4-8:2010	Magnetic field immunity
EN 61326-2-6:2021	Electrical Equipment for Measurement, Control, and Laboratory Use, EMC Requirements: Part 1 & IVD Equipment	EN 61000-4-11:2020	Immunity to Voltage Dips and Drops
ISTA 2A	Partial Simulation Performance Tests – Atmospheric preconditioning, Atmospheric Conditioning, Compression, Random vibration, Drop Shock, Fixed Displacement Vibration		

The undersigned hereby declares, on behalf of ALCOR Scientific, that the above-referenced products are manufactured in accordance with the referenced Directives and comply with the Standards listed. All supporting documentation is maintained at ALCOR Scientific.

Authorized By: G. Scott VanAlst
Title: PRRC, Director of QA and RA
Date of Issue: 4 February 2025 Expiry Date: 4 February 2027
Place of Issue: Smithfield, RI USA

Approved By:

[\(CO-3913\) iSED PRO Registration updates](#)

Description

After review of documents from Emergo to gain EU Rep approval, a few minor changes were needed. Added warning on page 31 of the iSED PRO IFU instructing the reader to report any serious incident to the EU Member state or the regulatory authority where the patient or device is located. Added deepCLEAN to iSED and miniiSED DofC under accessories. Made the name on the DofC consistent with other documentation. Added deepCLEAN Pro to DofC.

Justification

To fulfill regulatory requirements a change order is needed.

Assigned To:	Initiated By:	Priority:	Impact:
Scott VanAlst	Scott VanAlst	Medium	Minor

Version History:

Author	Effective Date	CO#	Ver.	Status
Scott VanAlst	February 7, 2025 9:12 AM EST	CO-3913	4	Published
Scott VanAlst	July 16, 2024 9:15 AM EDT	CO-3667	3	Superseded
Scott VanAlst	October 26, 2023 11:31 AM EDT	CO-3526	2	Superseded
Scott VanAlst	August 30, 2023 2:34 PM EDT	CO-3474	1	Superseded
Scott VanAlst	March 17, 2023 3:17 PM EDT	CO-3159	0	Superseded

iSED[®] / *iSED*[®] *ELITE*

Automated Erythrocyte Sedimentation Rate Analyzer



OPERATOR'S MANUAL

For iSED[®] analyzers with serial numbers above 05000 and all iSED ELITE analyzers

This page intentionally blank.

Symbols Reference

The following is a list of symbols and their meaning used on the instrument, consumables, and accessory labeling.

Symbol	Meaning
	Instrument satisfies requirements of European directive on in vitro diagnostic medical devices (98/79/EC)
	Date of Manufacture
	Manufacturer
	Serial Number
	In Vitro Diagnostic Medical Device
	Product/Reference Number
	Fuse Rating (located on serial number label, replace with same value and type)
	AC Single Phase Alternating Current
	Consult instructions – Refer operator to the instruction manual for additional information
	Temperature limitation – Indicates storage requirements range
	WEEE: Disposal of Waste Electrical and Electronic Equipment
	Biological Hazard: Universal Precautions Should be Followed
	Caution: Moving Parts
	Caution: Sharp Needle
	Warning: Consult Operator Manual and Observe Safety Warnings
	Caution: May Cause Electrical Shock
	Caution: Object is Heavy – Use care and/or assistance in lifting

Notes, Precautions, Warnings, and Biological Warnings: Interpretation Guide

The Operator's Manual includes information and warnings. These need to be observed by the operator in order to ensure safe operation of the instrument. There are four types of messages: Notes, Caution, Warnings, and Biological Warnings.

Notes

NOTE: Highlights important facts, gives helpful information and tips, and clarifies procedures.

Cautions

CAUTION: Electrical caution! Unplug before handling.

CAUTION: Important information on the proper operation of the instrument. This information is crucial in preventing instrument damage and maintaining the system.

Warnings

WARNING: Identifies potentially hazardous situations that could result in serious injury to laboratory personnel.

WARNING: Universal precautions should be followed. Always wear gloves to prevent exposure to pathogens.

Precautions and Safety Information

Please pay close attention to the instructions, notes and symbols as well as the standard laboratory practices outlined by your facility and local regulatory agencies.

Always keep a distance of at least 4 inches (10 cm) between the rear of the instrument and the wall to allow for proper ventilation.

Do not use power frequencies or voltage other than those specified in this document. Connection to an inappropriate power source may cause injury or fire.

Do not disassemble or modify the instrument. Doing so may cause injury and/or instrument malfunction and void the warranty.

Place the instrument on a stable and level surface free of vibration. Failure to do so may cause injury or malfunction of the unit.

CAUTION: To reduce the risk of electrical shock, do not remove any panel unless under the direction of qualified personnel.

Do not block any ventilation openings.

Do not place instrument in water.

Do not drop or throw the instrument.

Operate the instrument on a dry, level surface.

Do not move the instrument while specimens are processing.

Plug the instrument into a grounded power source.

Tubes must be tightly capped prior to loading them onto the iSED/iSED ELITE.

Do not operate iSED/iSED ELITE without sample collection tray.

Empty sample collection tray to prevent overflow of tubes.

WARNING: For continued protection against risk of fire and hazard, replace fuse only with the same type and rating fuse.

WARNING: The instrument's main power entry port is used as the main disconnect device.

WARNING: Observe Universal Precautions. Discard contaminated materials according to applicable regulations.

Contents

Symbols Reference	ii
Notes, Precautions, Warnings, and Biological Warnings: Interpretation Guide	iii
Precautions and Safety Information	iv
1. Intended Use	1
2. Methodology	1
2.1 History	1
2.2 Principle of the Procedure	1
2.3 iSED/iSED ELITE Technology	2
2.4 Traditional ESR Method Limitations	2
3. General Information	2
3.1 For In Vitro Diagnostic Use Only	2
3.2 Sample Requirements	3
3.3 Tube Requirements	3
4. Instrument Overview	3
4.1 Parts Identification	4
4.2 Consumables	5
4.3 iWASH Fluid	5
4.4 Unpacking and Installation	6
4.5 Unpacking the Instrument	6
4.6 Contents of the Box	7
4.7 Power Connection	7
4.8 RS-232 connection	8
4.9 Ethernet Connection	8
4.10 USB Connection	8
5. User Interface	9
5.1 Initial Setup	9
5.1.1 Select the language and date/time format	9
5.2 Touch Screen Menus	9
5.2.1 Results Log Menu	10
5.2.2 Maintenance Menus	11
5.2.3 Setting Menu	11
5.2.4 General Settings Menu	12
5.2.5 Advanced Settings Menu	12
6. Operating Instructions	13
6.1 Patient Identification	13
6.2 Basic Operation	14
6.2.1 Running Barcoded Samples	14
6.2.2 Running Non-Barcoded Samples	15
7. Quality Assurance	15
7.1 Quality Control	15

7.2 Peer Group Comparison	15
7.3 Proficiency Testing (PT)	16
8. iSED Limitations	16
8.1 Sample Conditions	16
9. Calibration	16
10. Results	17
10.1 Expected Values	17
10.2 Results Format	18
10.3 Printed Results with Error Message	18
11. Performance	19
11.1 Method Comparison	19
11.2 Precision	20
11.3 Stability	20
11.4 Carry-Over	20
12. Test Cards	21
12.1 Downloading Credits from Test Card	21
13. Routine Maintenance	22
13.1 Replacing Printer Paper	22
13.2 Replacing/Emptying the iWASTE Bottle	23
13.3 iWASTE Bottle Full Indicators and Alarms	24
13.4 Replacing iWASH Bottle	25
13.5 iWASH Bottle Empty Indicators and Alarms	26
13.6 Replacing the Fuse	27
14. Preventative Maintenance	28
14.1 Exterior Cleaning	28
14.2 Deep Clean Procedure	28
14.3 “Tubing Change is Required” Message	30
14.4 30,000 Test Aspirations Message	30
14.5 Spare Parts	30
14.6 System Status, Error Codes and Warning Messages	30
14.7 System Status Messages	31
14.8 System Warning and Error Messages	32
14.9 Sampling Error Message	39
14.10 Print Out of Sampling Error Message	40
14.11 Troubleshooting	40
15. Safety Precautions	41
15.1 General Considerations	41
15.2 Biological Waste	41
16. Technical Support	42
17. Technical Specification	43
18. Quick Reference Guide	44
19. Warranty Information	45

1. Intended Use

iSED/iSED ELITE Erythrocyte Sedimentation Rate analyzer is an automated erythrocyte sedimentation rate analyzer which reports sedimentation rate in mm/hr as a non-specific quantitative result. Testing is performed using EDTA whole blood samples obtained by venipuncture or capillary blood collection. The instrument can be used in laboratories authorized to perform tests classified as Moderately Complex under CLIA categorization by order of a physician to aid in assessing the general health status of a patient.

2. Methodology

2.1 History

In 1897 Polish physician, Edmund Faustyn Biernacki (1866–1911) was the first to observe the phenomenon of erythrocyte sedimentation rate (ESR). He discovered that blood sedimentation rates were different in different individuals, blood with small amounts of blood cells sedimented faster, and blood sedimentation rate directly correlated to levels of plasma fibrinogen¹. The findings presented by Biernacki clearly showed the clinical significance of ESR.

In 1921, Swedish internist Alf Vilhelm Albertsson Westergren (1891–1968), presented a similar description of ESR² as those given by Biernacki and Swedish hematologist Robert Sanno Fåhræus (1888–1968)³. Westergren defined the standard measurement of the ESR test to which nearly all automated ESR analyzers are referenced today^{4,5}. The Westergren method of ESR testing uses a standardized tube and sodium citrate anticoagulated blood.

2.2 Principle of the Procedure

ESR is a simple, non-specific screening test that indirectly measures the presence of inflammation in the body. It reflects the tendency of red blood cells (RBCs) to settle more rapidly in some disease states, usually because of increases in plasma fibrinogen, immunoglobulins, and other acute-phase proteins. Changes in red cell shape or numbers may also affect ESR.

When anticoagulated whole blood is allowed to stand in a narrow vertical tube (known as a Westergren tube) for a period of time, the RBCs settle out from the plasma. The rate at which they settle is measured as the number of millimeters of clear plasma present at the top of the column after one hour (mm/hr). The RBCs sediment because their density is greater than that of plasma; this is particularly so when there is an alteration in the distribution of charges on the surface of the RBC resulting in their coming together to form aggregates known as rouleaux. Rouleaux formation is determined largely by increased levels of plasma fibrinogen and globulins, so ESR mainly reflects changes in the plasma proteins that accompany acute and chronic infections, some tumors, and degenerative diseases. In these disease states, ESR values are typically elevated. ESR denotes the presence of tissue damage or disease but not its severity; it may be used to monitor disease progression or the effectiveness of treatment.

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⁴International Council for Standardization in Haematology (Expert Panel on Blood Rheology): ICSH recommendations for measurement of erythrocyte sedimentation rate. J Clin Pathol 1993; 46: 198–208

⁵Thomas RD, Westengard JC, Hay KL, et al: Calibration and validation for erythrocyte sedimentation tests. Arch Pathol Lab Med 1993; 117: 719–72.

2.3 iSED/iSED ELITE Technology

The iSED/iSED ELITE analyzer uses advanced photometric rheology technology to measure the rouleaux formation, the earliest phase of erythrocyte sedimentation. The aggregation of RBCs during the rouleaux formation ultimately determines the length at which the red cells sediment in the Westergren tube. The technical innovation of the iSED/iSED ELITE analyzer consists of “directly” measuring the aggregation of the RBCs, while traditional ESR methods “indirectly” measure the aggregation of the RBCs by recording the length at which the red cells settle in a Westergren tube.

The instrument’s micro-flow cell captures the critical kinetics of RBC aggregation in a highly controlled testing environment, which helps reduce the impact of factors that can contribute to result variability. **The analyzer is designed to sample directly from the primary EDTA blood collection tube and produce ESR within 20 seconds with appropriate prior homogenization.** Results are reported in “mm/hr,” and iSED/ iSED ELITE performance is highly correlated with the Westergren method.

2.4 Traditional ESR Method Limitations

The ESR is a transient phenomenon confined to fresh blood. It is not a hematic matrix component at the corpuscular or molecular level. The procedures used to determine ESR cannot be calibrated since they are susceptible to a variety of factors, e.g., temperature, hematocrit, erythrocyte mean corpuscular volume, plasma viscosity, etc. For this reason, it is possible to observe instrument performance deviations compared to other procedures when the aforementioned variables are not taken into account. Erythrocyte sedimentation remains an only partially understood phenomenon and is a clinically nonspecific reaction.

Sample mixing is performed at the beginning of the analysis with the purpose of homogenizing the sample. An inefficient homogenization can affect the results given by the instrument.

Interferences that may increase ESR values:

- Increased levels of fibrinogen and gamma globulins
- Technical factors such as mechanical vibrations and elevated room temperature

Interferences that may decrease ESR values:

- Abnormally shaped RBC (sickle cells, spherocytosis)
- Technical factors such as low room temperature, delay in test performance, clotted blood sample, excess anticoagulant, or bubbles in sample tube
- Specimen quality such as blood to anticoagulant ratio, age of specimen and fill volume⁶

NOTE: ESR is a nonspecific result. It is highly recommended to perform other tests together with ESR since an ESR value is not enough to exclude that the patient is not affected by a pathology or to diagnose a pathology.

3. General Information

Read this manual carefully prior to operating the instrument.

This document is the operator’s manual for the iSED analyzers with serial numbers above 5000 and iSED ELITE analyzers. It is intended to explain the instrument operation in detail and can be used as a basis for training new operators. It is an information guide and troubleshooting reference. Retain this manual for future use.

3.1 For In Vitro Diagnostic Use Only

⁶CLSI. Procedures for the Erythrocyte Sedimentation Rate Test; Approved Standard-Fifth Edition. CLSI document H02-A5. Wayne, PA: Clinical and Laboratory Standards Institute; 2011.

3.2 Sample Requirements


- Sample volume for testing is 100 μ L whole blood (approximately 500 μ L total volume recommended to ensure sample aspiration*)
- Sample volume for pediatric testing is 100 μ L whole blood (approximately 400 μ L total volume recommended to ensure sample aspiration*)
- Sample must be whole blood collected in K3-EDTA or K2 EDTA anticoagulant tube
- Sample should not be clotted or hemolyzed (DO NOT mix vigorously!)
- Sample should be tested within 4 hours from venipuncture or within 24 hours if refrigerated
- Sample must be brought to room temperature for at least fifteen (15) minutes (if refrigerated)

*The sample tube is inverted within the instrument while the sample is aspirated, therefore the total volume required to perform a test varies by tube type based on the piercing depth. Please contact ALCOR Scientific Technical Support for the most updated Sample Tube Compatibility Chart (100-07-003).

NOTE: The instrument requires no additional or special sample preparation. As with all anticoagulant collection tubes, the sample should be well mixed after collection to help avoid clotting or other aggregates that may alter ESR test results.

3.3 Tube Requirements

- Compatible with most standard 13x75 mm EDTA anticoagulant (lavender top) tubes with pierceable caps

NOTE: Please contact ALCOR® Scientific Technical Service for the current list of compatible sample tubes.


WARNING: Do not use if the tube stopper is missing. Use only tightly capped samples.


4. Instrument Overview

iSED/iSED ELITE Erythrocyte Sedimentation Rate Analyzer uses photometric rheology to directly measure the aggregation of RBCs. Once the sample is automatically mixed and in position, a sensitive optical detector in the iSED/iSED ELITE analyzer follows the progress of aggregation over time by monitoring light transmission through the sample. This produces a signal that is a direct representation of the aggregation. The magnitude of time-dependent change is correlated to the Westergren method.

4.1 Parts Identification

1	iWASH and iWASTE Compartment
2	Touch Screen
3	Sample Tube Return Tray
4	Printer
5	Paper Feed Button
6	Sample Loading Port
7	Smart Card Reader
8	Sample Ejection Port
9	On/Off Button
10	RS-232 Connection Port
11	Fuse
12	Power Connection Port (24VDC, 2.5A)
13	USB Ports (2)
14	Ethernet Network Port

4.2 Consumables

Item	Description	Quantity	Reorder Part #
Printer Paper	57 mm x 28 mm OD	Pack of 5	DS-05233
Test Card	Preloaded smart card for iSED/ iSED ELITE	250 tests	112-00250
		500 tests	112-00500
		1,000 tests	112-01000
		2,000 tests	112-02000
		5,000 tests	112-05000
		10,000 tests	112-10000
		20,000 tests	112-20000
iWASH® Fluid	500 mL bottle with screw cap, pre-filled with instrument iWASH	Pack of 4	112-12-001
iWASTE® Bottle	500 mL plastic waste bottle with screw cap	Pack of 4	112-12-005
		Pack of 24	112-12-002

NOTE: Only use consumables that are within their expiration date.

NOTE: Use of any other products could affect the performance of the instrument and void the warranty.

4.3 iWASH Fluid

The instrument uses iWASH Fluid as the cleansing agent during the wash cycle. **The use of any other product could affect the performance of the instrument and void the warranty.**

It is recommended that the instrument remain on at all times and ready for use. Should the instrument need to be powered off for any reason, run a wash cycle prior to powering off the unit.

NOTE: The instrument is programmed to perform self-cleaning after being idle for fifteen (15) minutes following the last sample tested. The process takes approximately one (1) minute and utilizes approximately 4.5mL of iWASH for each wash cycle. Once completed, testing can resume as normal.

4.4 Unpacking and Installation

CAUTION: The instrument unit weighs 33 lbs in the box. Use safe lifting techniques and proper techniques when handling heavy objects. If necessary, obtain assistance to safely lift the instrument.

CAUTION: If using a utility knife, extend the blade to appropriate length to avoid cutting any internal components.

All original packaging should be kept in the event the instrument needs to be returned for service or warranty repair. For more information, please refer to the Warranty Information in the Operator's Manual or call ALCOR Scientific Customer Service at +1 401.737.3774.

4.5 Unpacking the Instrument

Inspect the shipping container for any obvious signs of mishandling or shipping damage. If damage is found, retain all package materials and immediately file a claim with the shipping carrier.

Figure 1


Figure 2


Figure 3



Figure 4

1. Position the box upright and open the top flaps (Figure 1).
2. Remove the Accessory Box and set aside (Figure 2).
3. Keep the box upright. Slowly slide the instrument and surrounding foam out of the box (Figure 3).
4. Remove the Tube Return Tray and set aside (Figure 4).
5. Remove the foam panels from the sides of the instrument (Figure 4).
6. Place the instrument on a secure flat surface and remove the instrument from the protective bag. **Save the box and foam pieces for future use** (Figure 4).

4.6 Contents of the Box

1. iSED/iSED ELITE instrument (1 Each)
2. Power Cord and Power Adapter (1 Each)
3. Sample Collection Tray (1 Each)
4. Pre-filled iWASH Bottle (1 Each)
5. iWASTE Bottle (1 Each)
6. iWASH Filter Cap (1 Each)
7. iWASTE Cap (1 Each)
8. Thermal Paper (1 Each)
9. Spare Fuse (1 Each)
10. Quick Reference Card (1 Each)

4.7 Power Connection

1. Connect the power cord to the power supply.
2. Install and connect both iWASH and iWASTE bottles. See Sections 13.2 and 13.4.
3. Insert the power adapter cord into the power connection port located on the rear panel of the instrument (Figure 5).


Figure 5

4. Place the instrument in its permanent operating location, and plug the power cord into a standard wall outlet.
 5. To power the unit on, press the power button located on the rear of the instrument. The instrument will produce an audible beep, after which the instrument is non-functional while the operating system boots up. This bootup can last up to 1 minute.
- The device will not accept samples until the temperature of the measurement components has warmed up and stabilized.


CAUTION: Always keep a distance of at least four (4) inches (10 cm) between the rear of the instrument and the wall to allow for proper ventilation.

CAUTION: Place and operate the instrument on a stable and level surface free of vibration. Failure to do so may cause injury or malfunction of the unit.

4.8 RS-232 connection

The analyzer is equipped with a legacy RS232 DB9 male connector for data transfer to LIS systems using the LIS2-A2 Standard. The pin-out of the connector is detailed in Figure 6.

For more information, Document 222-09-006 Communication Protocol is available upon request.

4.9 Ethernet Connection

The analyzer is equipped with an RJ-45 Ethernet connector for Factory use and for connection to TCP/IP-based LIS systems using the LIS2-A2 Standard. For more information, Document 222-09-006 Communication Protocol is available upon request.

4.10 USB Connection

The analyzer is equipped with two USB 2.0 interface connectors to facilitate the export of test results and for updating device software.

5. User Interface

5.1 Initial Setup

5.1.1 Select the language and date/time format

The iSED/iSED ELITE will recognize when it is being powered up for the first time and will guide the user through the initial set up process of selecting a language and the date/time format (Figures 7 and 8). Once Initial Setup is complete, the analyzer is ready for basic operation (Figure 9).


Figure 7


Figure 8


Figure 9

5.2 Touch Screen Menus

The instrument is operated via touchscreen, and all operations can be done by selecting or inputting data on the following screens.

Navigate to the Main Menu screen (Figure 11) by pressing the symbol in the upper left-hand corner found on many of the menu screens (see example in Figure 10). From the Main Menu screen (Figure 11), navigation to the Result Log, Maintenance and Settings menus can occur.


Figure 10


Figure 11

5.2.1 Results Log Menu

The Result Log sub-menu is organized with the most recent test results shown first. Figure 12 icons, listed left to right, top to bottom are:

1. Back to main menu
2. Search: Selecting this will allow search by date, name, patient number, etc.
3. Recent: The default page view; orders the result chronologically with most recent first.
4. This “Papers” icon is for exporting records.
5. The blue circle represent an individual test results of a patient record.
6. The white circle represents a SEDiTROL® or proficiency result record.
7. This is the test or sample number.
8. This icon is an expand command, providing a more detailed look at that test record.

Records can be exported as shown in Figure 12 or individually as shown in Figure 13. In either case, once the “Papers” icon is selected, it will bring up the Export Selection screen (Figure 14).


Figure 12

Records can be exported as shown in Figure 12 or individually as shown in Figure 13. In either case, once the “Papers” icon (circled in Figure 13) is selected, it will bring up the Export Selection screen (Figure 14).


Figure 13


Figure 14

5.2.2 Maintenance Menus

The Maintenance sub-menu provides a “dashboard” for the user-related functions for iSED/iSED ELITE. When green, the circular icons indicate that the listed category is within preset limits. When an iSED/iSED ELITE category is approaching a preset limit, the circular icon will turn yellow, indicating attention is needed soon. Should the category limit be exceeded, the icon will turn red and require action to be taken. In Figure 15, displayed icons, top to bottom are:


Figure 15

- 1. Cleaning Fluid and Waste Container:** In the “Cleaning Fluid” and “Waste Container” categories, selecting the “Replace” or “Empty” buttons will reset the counters.
- 2. Quick Clean:** In the “Quick Clean” category, selecting the “Run” button will start the cleaning process.
- 3. Deep Clean:** In the “Deep Clean” category, selecting the “Run” button will start the deep cleaning process.
- 4. Test Credits:** The “Test Credits” category displays the remaining Test Credits. When a new Test Credit card is available, insertion of the Test Credit card will guide the user through the test Credit transfer process.
- 5. Reading Cell Temperature:** The “Reading Cell Temperature” category displays the status of the LED and Sensor temperatures. If the LED and Sensor temperature are not up to the correct temperature, the status will display “Not Ready” and samples cannot be processed until the status displays “OK”.

5.2.3 Setting Menu

The Settings sub-menu (Figure 16) provides the user access to functions to customize the operation of the iSED/iSED ELITE for their laboratory environment. Using the guided Setup process when the instrument is first turned on will preconfigure some of the typical settings at the time of device installation. There are two categories of settings: General and Advanced. The General Settings selection is accessible to the user via the Admin Pin when activated. If not activated, settings will not be password protected. Any user can view the Advanced Settings level but changes to these parameters is only allowed once the Advanced level Pin is entered.


Figure 16

5.2.4 General Settings Menu

The General Settings sub-Menu (Figure 17) provides the user access to functions to customize the operation of the iSED/iSED ELITE for their laboratory environment. Using the guided Setup process will preconfigure some of the typical settings at the time of device installation. Protection of the General Settings can be disabled by selecting the Access option. Use of the Backup and Restore features should be used only with the guidance of ALCOR Scientific Technical Support, as misuse of these features can result in loss of data and settings.


Figure 17

5.2.5 Advanced Settings Menu

The Advanced Settings sub-menu provides the user “Read Only” status of various operational sensors and device settings (Figures 18,19, and 20). The Advanced Settings level is only accessible to ALCOR Scientific qualified technicians.


Figure 18


Figure 19



Figure 20

6. Operating Instructions

6.1 Patient Identification

Barcoded Tubes: Samples in barcoded tubes are read and identified automatically by the instrument's internal barcode reader as they are loaded into the instrument. All common laboratory barcodes are supported, including Code 39, UPC, and Code 93 formats. Note the below barcode location range:

For instances when patient identification cannot be read by the internal barcode reader or if there is no barcode present, the operator may enter data manually. **For instructions on manually entering patient data, please refer to Section 6.2.2.**

6.2 Basic Operation

6.2.1 Running Barcoded Samples

All sample mixing, sample extraction, sample reading, and sample disposal are handled automatically by the instrument. Up to 20 sample tubes may be loaded into the sample wheel at any given time. As each sample is processed (processing happens within 20 seconds after appropriate sample mixing), the sample tube is ejected from the sample wheel and retained in the external sample collection tray. Whenever there is an open position, a sample may be placed in the sample wheel.

Test credits are required for operation. The instrument is loaded with a predetermined quantity of credits for initial set up and use, however additional credits in the form of 'Test Cards' must be purchased. For information on loading additional test credits, please refer to Section 12.

To run a barcoded sample:

6.2.1.1

Select the barcode icon with the plus sign (Figure 21).


Figure 21

6.2.1.2

Insert the capped pierceable primary tube, with the barcode facing to the right, into the loading port of the analyzer. A red light will illuminate, and a distinctive beep will sound when the barcode is successfully recognized.

6.2.1.3

A position will highlight green on the wheel displayed on the main menu, indicating the sample is on the sample wheel.

6.2.1.4

The sample wheel will move to the next available position to load additional samples. The onscreen information bar will report "waiting for sample" and the instrument will beep quietly for five (5) seconds. As the five (5) second window draws to a close, beeping will become faster. If no additional samples are loaded, the analyzer mixes the samples for 3 minutes.

NOTE: If the five (5) second window is missed, simply select the barcode icon again to restart the sample scheduling process.

6.2.1.5

After a 3-minute mixing cycle, the analyzer starts the testing process.

6.2.1.6

When testing is complete, the tube is ejected from the sample wheel into the collection tray, and the test result is printed via the onboard printer.

6.2.2 Running Non-Barcoded Samples

6.2.2.1

Select the “person with the plus sign” icon (Figure 22).

The instrument will prompt the operator to enter sample or patient identification data manually using the alphanumeric keyboard. Sample or Patient Information should be recorded in at least one (1) of the following data fields:

- Alphanumeric ID
- Patient's First Name
- Patient's Surname

NOTE: If patient information is not entered within ten (10) seconds from the last pressed key, the instrument will automatically assign an identification number.


Figure 22

6.2.2.2

Insert the capped, pierceable, primary tube into the loading port of the analyzer. A position will highlight green on the wheel displayed on the main menu, indicating the sample is on the wheel.

6.2.2.3

The wheel will move to the next available position to load additional samples. If no additional samples are loaded, the analyzer mixes the samples for 3 minutes.

6.2.2.4

After the 3-minute mixing cycle the analyzer starts the testing process.

6.2.2.5

When testing is complete, the tube is ejected from the sample wheel and the test result is printed on via the onboard printer.

7. Quality Assurance

7.1 Quality Control

The manufacturer recommends running external controls at least once per day. SEDiTROL ESR Controls Levels 1&2 are available for purchase.

Part Number	Configuration
DSC01	Pack of 2 x 4.0 mL
DSC06	Pack of 6 x 4.0 mL

7.2 Peer Group Comparison

iQAP is available for SEDiTROL customers as an online peer-to-peer quality assurance program. Contact ALCOR Scientific for more information.

7.3 Proficiency Testing (PT)

Quality Assessment is a critical aspect of laboratory quality management and can be conducted in several ways. One of the commonly employed assessment methods is that of external quality assessment, or Proficiency Testing.

Proficiency Testing is an important tool used in the lab to verify accuracy and reliability of its testing methods, alert to areas of testing that are not performing as expected, and to indicate shifts and trends which over time may affect patient results.

There are several manufacturers of Proficiency Testing materials from which to choose.

NOTE: To ensure consistent results, Proficiency Testing samples must be well mixed for at least 25 minutes on a mechanical rocker or rotator prior to testing. Please follow your selected Proficiency Testing manufacturer's instructions for optimal sample processing procedures and results. Use only approved Proficiency Testing material.

8. iSED Limitations

8.1 Sample Conditions

Lipemia: Changes in viscosity of the sample can interfere with measurement. iSED/iSED ELITE will give an error message if it detects a drastic change.

Hemolysis: If hemolysis has occurred to such a degree that aggregation of RBCs has been reduced, it will affect the results.

Clotted sample: The needle used in iSED/iSED ELITE is designed to prevent interfering clots from being aspirated into the hydraulic circuit and reading cell. If a clot prevents the sample from aspirating, the analyzer will retry the aspiration process three (3) times before it gives an "unable to withdraw" error code and aborts the test.

Short draw sample: If the sample volume is too low to be aspirated, the iSED/iSED ELITE will attempt to aspirate the sample three (3) times before it gives an "unable to withdraw" error code and aborts the test.

9. Calibration

iSED/iSED ELITE instruments are factory calibrated using samples compared to a Reference Instrument that has been correlated to the reference Westergren method. The instrument range is from 1 to 130 mm/hr. During normal operation, parameters affecting calibration are constantly monitored and, if not within expected limits, a warning is given and further testing prevented.

10. Results

10.1 Expected Values

The reference values found in the table below are averages found in males and females. An increase in these values can be a sign of multiple different health issues that should be diagnosed by a physician or qualified individual.

Sedimentation Rate Reference Value (mm/hr) ⁷	
Males under 50 years old	< 15
Males over 50 years old	< 20
Females under 50 years old	< 20
Females over 50 years old	< 30


NOTE: The ranges provided are for reference only. All laboratories should follow their laboratory's protocol for verifying their own reference ranges.

⁷Keohane EM, Otto CN, Walenga JM. *Rodak's Hematology: Clinical Principles and Applications Sixth Edition*. Elsevier; 2019.

10.2 Results Format

Results are shown on screen after analysis and printed by the instrument's internal printer. Data format is as follows:

Date format: Month/Day/Year
Time format: Hour/Minute/Second
Result format: mm/Hour

Example of Normal Results Print Out

```
=====
Date: 03/25/2013           Date of measurement
Time: 13:36:24            Timestamp of measurement
iSED Sn: 00001            Instrument serial number
ID: 812409                Barcoded or Manual sample identification
ESR (mm/h): 15           Format of ESR result reported
=====
```

Example of High Results Print Out

```
=====
Date: 03/25/2013           Date of measurement
Time: 13:36:24            Timestamp of measurement
iSED Sn: 00001            Instrument serial number
ID: 812409                Barcoded or Manual sample identification
ESR (mm/h): 130          Format of print out if high ESR result reported
=====
```

Example of Low Results Print Out

```
=====
Date: 03/25/2013           Date of measurement
Time: 13:36:24            Timestamp of measurement
iSED Sn: 00001            Instrument serial number
ID: 812409                Barcoded or Manual sample identification
ESR (mm/h): Less than 1  Format of print out if low ESR result reported
=====
```

NOTE: If a sample is entered with Manual identification, then Patient Name (First, Last, or both) will be shown underneath the device serial number.

10.3 Printed Results with Error Message

In the event that the instrument is unable to analyze the sample and report results, the printout will replace the 'ESR (mm/h):' field with an error message. **For more information on Error Messages, please reference to Section 14.**

11. Performance

11.1 Method Comparison

iSED/iSED ELITE ESR analyzer has been demonstrated to provide results equivalent to the CLSI approved standard methodology for ESR, the Westergren method.

Since ESR is based on the physical interactions of RBCs during testing, the Westergren reference method is subject to a number of variables, including sample characteristics, testing environment, and individual operators' techniques. Therefore, Passing Bablok regression analysis is recommended for comparing two ESR methods since the reference method is known to be subject to variables.¹

Results of Passing Bablok regression analysis:

Variable X	Westergren
Variable Y	iSED ELITE
Sample Size	200

y= 1.000000 + 1.000000 x	
Systematic Differences	
Intercept A	1.0000
95% CI	-0.4296 to 2.3056
Proportional Differences	
Slope B	1.0000
95% CI	0.9444 to 1.0704
Random Differences	
Residual Standard Deviation (RSD)	9.7999
±1.96 RSD Interval	-19.2078 to 19.2078
Linear Model Validity	
Cusum Test for Linearity	No significant deviation from linearity (P=0.28)

11.2 Precision

Samples spanning a range of ESR values were run on the iSED/iSED ELITE analyzer a total of ten (10) times to demonstrate precision. As with other laboratory tests, higher CVs are expected when comparing lower numerical values.

iSED/iSED ELITE results:

Run #/Range	10--20	20--40	40--60	60--105
Run 1	13	21	48	101
Run 2	13	22	52	102
Run 3	14	22	47	99
Run 4	13	22	46	111
Run 5	12	21	53	100
Run 6	12	22	51	98
Run 7	13	24	51	100
Run 8	11	24	52	104
Run 9	13	24	49	101
Run 10	13	25	52	110
Average	12.7	22.7	50.1	102.6
SD	0.82	1.42	2.42	4.48
% CV	6.48	6.25	4.84	4.36

11.3 Stability

Samples spanning a range of ESR values were run on the iSED/iSED ELITE analyzer at 24 hours to demonstrate stability. As with other laboratory tests, higher CVs are expected when comparing lower numerical values.

iSED/iSED ELITE results:

Average %CV	24 Hours Refrigerated (mm/hr)
Lower ESR Result	14.88
Medium ESR Result	6.22
High ESR Result	2.47

11.4 Carry-Over

Carryover testing was performed to demonstrate the sequence of sample testing does not affect measured values.

iSED/iSED ELITE results:

Sample	Run 1	Run 2	Run 3	Run 4	Average
Sample 1a	70	68	74	67	69.8
Sample 1b	73	71	74	75	73.3
Sample 2a	7	6	7	7	6.8
Sample 2b	5	5	5	5	5.0
% Carryover	3.1	1.6	2.9	3.2	2.7

1) CLSI. *Procedures for the Erythrocyte Sedimentation Rate Test; Approved Standard—Fifth Edition*. CLSI document H02-A5. Wayne, PA: Clinical and Laboratory Standards Institute; 2011.

12. Test Cards

In order to process and analyze samples, tests, known as ‘credits’, must be downloaded onto the instrument from a smart card preloaded with tests of various quantities.

12.1 Downloading Credits from Test Card

The Home screen (Figure 23) and the Maintenance sub-menu (Figure 24) provide the user with the number of test credits available. To add additional credits to the instrument, the user must insert a test card into the test card reader located on the front right corner of the iSED/iSED ELITE with the arrow facing upward and forward (Figure 25).


Figure 23


Figure 24


Figure 25

Once a test card is inserted, the next screen will show the number of credits currently on the instrument, the serial number of the inserted test card, and the number of credits on the test card (Figure 26). If the user taps the “Confirm” bubble, all the available credits will be transferred to the instrument. The total number of credits is then displayed, and the user is instructed to discard the test card (Figure 27).


Figure 26

Figure 27

13. Routine Maintenance

13.1 Replacing Printer Paper

A green LED light around the printer paper feed button will flash to indicate it is out of paper. To replace the printer paper in the instrument, the procedure below should be followed:

1. Pull the lever in Figure 28 until the lid is released from its locked position.
2. Open the paper cup lid and remove the remaining paper.
3. Insert thermal paper roll into the printer with the paper unwinding from the bottom of the roll.
4. Reel off a few inches from a new roll of paper. Hold approximately two (2) inches of paper outside the printer as you place the new roll into the reservoir.
5. Close the lid by applying equal amounts of pressure on each side ensuring the lid is in the locked position.


Figure 28

13.2 Replacing/Emptying the iWASTE Bottle

WARNING: Wear Personal Protective Equipment (PPE) such as protective gloves and safety glasses during this operation.

NOTE: A washing cycle should be run prior to replacing the waste bottle.

1. Open the front door to access the bottle compartment (A).
2. Locate the waste bottle in the upper compartment (B).
3. Disconnect the LUER connector (C) from the waste bottle screw cap.
4. Remove the waste bottle from the instrument and dispose according to your laboratory biologic waste protocol.
5. Replace the iWASTE bottle in the upper compartment (B), and **firmly** reconnect the LUER connector (C) on the plastic screw cap with the vent hole positioned at top.
6. Close the front door (A).
7. Press Waste Container “Empty” button in the Maintenance sub-menu to reset the counter to 0% full (Figure 29).





Figure 29

NOTE: Be sure to replace the plastic cap with the vent hole at the top.

NOTE: Be careful not to kink the line when replacing the bottle.

NOTE: It is recommended that the waste bottle be emptied daily.

NOTE: This procedure can be done without the waste alarm being triggered.

13.3 iWASTE Bottle Full Indicators and Alarms

In the case of a full iWASTE bottle, a warning message will appear on the screen and be accompanied by an alarm alerting the operator of an error or warning message.

WARNING: This action should be done when this message appears.

Full Waste Bottle

In the event the iWASTE bottle is full, the below error message (Figure 30) will appear on the screen and can be resolved by first replacing or emptying the iWASTE container and then pressing “Empty” to clear the error and continue testing.


Figure 30

“Empty” button: This option must be selected immediately after the operator replaces the iWASTE bottle. The iWASTE bottle counter will automatically reset once this button has been pressed, and the instrument will continue with the sample loading or washing process. **For instructions on replacing the iWASTE bottle, please refer to Section 13.2.**

WARNING: Wear Personal Protective Equipment (PPE) such as protective gloves and safety glasses during this operation.

13.4 Replacing iWASH Bottle

1. Open the front door to access the bottle compartment (A).
2. The iWASH bottle is located in the lower compartment (D).
3. Disconnect the LUER connector (E) from the iWASH bottle screw cap.
4. Remove the empty iWASH bottle, unscrew the cap, and replace it with a new iWASH bottle.
5. Place the new iWASH bottle in the lower compartment and firmly reconnect the LUER connector (E) on the plastic screw cap with the vent hole positioned at top.
6. Close the front door (A).
7. Press the Cleaning Fluid “Replace” button in the Maintenance sub-menu to reset counter to 100% full (Figure 31).





Figure 31

NOTE: Be sure to replace the plastic cap with the vent hole at the top.

NOTE: Be careful not to kink the line when replacing the bottle.

NOTE: The instrument is programmed to perform self-cleaning after being idle for fifteen (15) minutes following the last sample tested. The process takes approximately one (1) minute and utilizes 4.5mL of iWASH for each iWASH cycle. Once completed, testing can resume as normal.

NOTE: This procedure can be done without the wash alarm being triggered.

13.5 iWASH Bottle Empty Indicators and Alarms

When the iWASH bottle is empty, a message will appear on the screen and be accompanied by an alarm alerting the operator of the error or warning message.

WARNING: This action should be done when this message appears.

Empty iWASH Bottle

In the event the iWASH bottle is empty, the below error message (Figure 32) will appear on the screen and can only be resolved by first replacing the iWASH bottle then selecting “Replace” to clear the error and continue testing.


Figure 32

“Replace” button: This option must be selected immediately after the operator replaces the iWASH bottle. The iWASH bottle counter will automatically reset once this button has been pressed, and the instrument will continue with the sample loading or washing processes. **For instructions on replacing the iWASH bottle, please refer to Section 13.4.**

13.6 Replacing the Fuse

CAUTION: Unplug the instrument from the AC wall outlet before replacing the fuse.


CAUTION: For continued protection against risk of fire and hazard, replace only with the same type and rating fuse.


Materials Needed:

- 3/16 inch blade Screwdriver
- Fuse T2A 250V 5x20mm 


1. Remove fuse holder on the rear of the instrument by gently pushing with the screwdriver and rotating counterclockwise.


2. Remove the fuse holder from the instrument.


3. Remove the old fuse from the fuse holder.

4. Insert the new fuse of the same type and rating into the fuse holder.

5. Return the fuse holder into the instrument and lock it back into place by turning it clockwise.

14. Preventative Maintenance

It is recommended that the instrument be kept free from dusty and particulate environments at all times for best performance. If such environments are unavoidable, periodically inspect interior surfaces and rear fan assembly for heavy dust accumulation and clean as needed.

14.1 Exterior Cleaning

- Do not use sterilizing solutions
- Prolonged exposure to alcohol or strong cleaners may damage the instrument housing
- Only use water and mild detergents to wipe the surface of the instrument sample tray

14.2 Deep Clean Procedure

The analyzer will prompt the user when the Deep Clean is needed. This prompt will be located in the notification area at the top of the Home screen (Figure 33). The frequency for Deep Cleaning is monthly or every 1000 samples run, whichever comes first. This procedure will clean the aspiration pathway from the needle to the reading cell.

Materials Needed:

- Empty and unused 13x75 tube (Do not use SST tube)
- 6-7 % hypochlorite (bleach): Do not dilute unless greater than 7%
- On board iWASH solution


Figure 33

Procedure:

1. Add approximately 3.5 ml of 6-7 % hypochlorite to an unused capped 13 x 75 mm tube.
2. Press Deep Clean “Run” button in the Maintenance sub-menu to start a Deep Clean (Figure 34).
3. The analyzer will check for sufficient bottle volumes and prompt to replace if insufficient. If sufficient volume, the analyzer will prompt to insert the tube filled with diluted hypochlorite from Step 1 (Figure 35).
4. Once prompted by the screen, insert the Deep Clean tube into the sample loading position to start the automated process.


Figure 35


Figure 36


Figure 34

5. The analyzer will run 2 wash cycles then automatically perform the Deep Clean (3 minutes) and conclude by automatically running two additional wash cycles (Figure 36). During the 3 minute Deep Clean, the device will appear to be idle while the fluidic system soaks in hypochlorite.
6. Once the Deep Clean procedure is completed, the analyzer will perform 2 additional Quick Clean actions to purge the system of hypochlorite and present the hypochlorite filled tube back to the loading port where it should be removed and discarded.

14.3 “Tubing Change is Required” Message

After 200 hours of continuous pump use, iSED/iSED ELITE will alarm and generate a message that a ‘Tubing Change is Required’. This message only serves as a warning that maintenance should be performed and does not prevent the analyzer from operational use. ALCOR Scientific Technical Support should be contacted by phone at (800) 495-5270 or +1 (401) 737-3774. Once Technical Support is contacted, parts and instructions for items that should be replaced will be sent.

14.4. 30,000 Test Aspirations Message

After 30,000 aspirations, iSED/iSED ELITE will alarm and generate a message to contact ALCOR Scientific Technical Support. Please contact Technical Support at (800) 495-5270 or +1 (401) 737-3774 when this occurs. Once Technical Support has been contacted, parts and instructions for items that should be replaced will be sent. This message only serves as a warning that maintenance should be performed and does not prevent the analyzer from operational use.

14.5. Spare Parts

Spare parts can be purchased by calling ALCOR Scientific Customer Service Department or your local authorized iSED/iSED ELITE distributor. Please contact ALCOR Scientific Technical Support Department for help determining what parts are needed.

14.6. System Status, Error Codes and Warning Messages

The instrument touchscreen display has a highlighted “window” at the top of the Home screen where all active system messages appear. The Status Window is dedicated to System Status Messages and the number of Test Credits remaining. An animated emoticon, named “SEDRick,” on the right side of the status window provides a quick visual reference to general operating status.

14.7 System Status Messages


Figure 37

Each of these messages display on the touch screen as the system is processing specimens:

Status	Status Explanation
"Ready to Test"	Device is Ready or Not Ready to accept sample tubes. Device is not ready for a number of reasons, such as any blocking errors, Reading Cell Temperature warm up period, or No test credits are available.
<i>number</i> "Tests"	'Quantity' Test credits available Low – Purchase More Tests (Alarm) 0 – No Tests Available
Fluid Levels	Green – No errors or warnings in regards to bottle capacity. Yellow – iWASTE bottle is nearing full capacity and/or iWASH bottle is nearing empty. Red – iWASTE bottle is completely full and/or iWASH bottle is completely empty.
Device Status Text	Status Explanation
"iSED is Idle"	All scheduled testing complete.
"iSED is scanning"	iSED/iSED ELITE is attempting to scan a sample's barcode.
"iSED is positioning rotor at home"	iSED/iSED ELITE is returning the mixing wheel to its Home position.
"iSED is positioning rotor for sample loading"	iSED/iSED ELITE is moving the mixing wheel to the next open position.
"iSED is positioning for ejection"	iSED/iSED ELITE is moving the mixing wheel to position a sample for ejection into the collection tray.
"iSED has loading port blocked"	iSED/iSED ELITE's finger safety sensor is blocked, preventing the mixing wheel from moving.
"iSED is waiting for cuvette"	iSED/iSED ELITE is ready for user to insert a new sample into the loading port (repeating beeps will sound).
"iSED is withdrawing from cuvette"	iSED/iSED ELITE is piercing the sample tube and aspirating sample.

"iSED is mixing"	Sample wheel rotating to mix all samples.
"iSED is ejecting cuvette"	Testing is complete and tube is being ejected from iSED/iSED ELITE instrument.
"iSED is calibrating tail sensor"	iSED/iSED ELITE is preparing the tail sensor for detection of blood from next sample.
"iSED is priming reading cell"	iSED/iSED ELITE is preparing the reading cell for measurement of next sample.
"iSED is positioning sample for withdrawal"	iSED/iSED ELITE is moving the mixing wheel to position the next sample for piercing and aspiration.
"iSED is measuring sample"	Sample is positioned in reading cell and analysis is underway.
"iSED is calculating result"	iSED/iSED ELITE is performing the calculations to provide an ESR measurement.
"iSED is resetting Needle"	iSED/iSED ELITE is recovering from an error by resetting the piercing system position.
"iSED is resetting Ejector"	iSED/iSED ELITE is recovering from an error by resetting the ejector position.
"iSED is washing"	iSED/iSED ELITE is performing a quick clean cycle to remove residual sample and prevent clogging.
Device Error Status Text	Status Error Explanation
"iSED has rotor error"	Device is in an error state related to the mixing wheel.
"iSED has needle error"	Device is in an error state related to positioning the piercing system.
"iSED has ejector error"	Device is in an error state related to the sample extractor.
"iSED has withdrawal error"	Device is in an error state related to sample aspiration.

NOTE: Unless the third line starts with "iSED is Warning" or "iSED is in Error", the operation is normal. Descriptions for iSED/iSED ELITE Warning and Error messages are listed in the following sections.

14.8 System Warning and Error Messages

In the event of a system warning or error, the following alert will appear on the instrument's screen and be accompanied by the appropriate error message (see examples in Figures 38 and 39):


Figure 38


Figure 39

The next table shows examples of the warning or error messages the user may see while operating the instrument and some possible solutions. Contact ALCOR Scientific Technical Support if the error cannot be resolved by any solution provided below:

Error Code	On Screen Text	Troubleshooting
N/A	LED Temperature High (or Low)	Reading Cell thermal control shows abnormal temperature. If error does not disappear after 2-3 minutes, contact Technical Support. This error prevents a measurement from starting.
N/A	Sensor Temperature High (or Low)	Reading Cell thermal control shows abnormal temperature. If error does not disappear after 2-3 minutes, contact Technical Support. This error prevents a measurement from starting.
N/A	'iWASTE Bottle Full' message displayed and Alarm	Remove and replace iWASTE bottle. See Section 13.5
N/A	'Wash Bottle Empty Message' displayed and Alarm	Remove and replace iWASH bottle. See Section 13.4
N/A	Paper Error/Out (Flashing Green Light)	Replace paper; see section 13.1
rotorErrorFingerDetect	Rotor Error: Loading Port obstruction detected during rotor movement.	Check for obstructions in the loading port. If no obstructions are present, power off device, remove all sample tubes, and power on device again. If the error appears again, contact Technical Support.
rotorErrorHomeOff	Rotor Error: Home Sensor not detected during rotation	Restart the unit, eject the sample and re-enter the sample. If the error appears again, contact Technical Support.
rotorErrorNeedleOff	Rotor Error: Unable to move rotor. Needle is in error state.	Restart the unit, eject the sample and re-enter the sample. If the error appears again, contact Technical Support.
rotorErrorEjectorOff	Rotor Error: Unable to move rotor. Ejector is in error state.	Restart the unit, eject the sample and re-enter the sample. If the error appears again, contact Technical Support.
rotorErrorPosReq	Rotor Error: Rotor requested to move to invalid position.	Restart the unit, eject the sample and re-enter the sample. If the error appears again, contact Technical Support.
rotorErrorStationUndefined	Rotor Error: Rotor requested to move to undefined station.	Restart the unit, eject the sample and re-enter the sample. If the error appears again, contact Technical Support.
rotorErrorUndefined	Rotor Error: Rotor in undefined state.	Restart the unit, eject the sample and re-enter the sample. If the error appears again, contact Technical Support.

Error Code	UI Text	Solutions
needleErrorSensorsUndefined	Needle Error: Invalid sensor state. Home and Up Sensor both active.	Restart the unit, eject the sample and re-enter the sample. If the error appears again, contact Technical Support.
needleErrorHomeSens	Needle Error: Unable to reach Home position. Please check Piercing System for obstructions.	Restart the unit, eject the sample and re-enter the sample. If the error appears again, contact Technical Support.
needleErrorUpSensor	Needle Error: Unable to reach Up position. Please check Piercing System for obstructions.	Restart the unit, eject the sample and re-enter the sample. If the error appears again, contact Technical Support.
needleErrorProbeTube	Needle Error: Invalid sensor state. Home and Probe sensors both active.	Restart the unit, eject the sample and re-enter the sample. If the error appears again, contact Technical Support.
needleErrorNoTube	Needle Error: No Tube detected when tube was expected.	Restart the unit, eject the sample and re-enter the sample. If the error appears again, contact Technical Support.
needleErrorNoMovement	Needle Error: Needle did not move from Home position after requesting it to move up.	Restart the unit, eject the sample and re-enter the sample. If the error appears again, contact Technical Support.
ejectorErrorSensorsAllSensorOff	Ejector Error: Invalid Sensor state. All position sensors active.	Restart the unit, eject the sample and re-enter the sample. If the error appears again, contact Technical Support.
ejectorErrorSensorsHomeLockOff	Ejector Error: Invalid Sensor state. Home and Lock position sensors both active.	Restart the unit, eject the sample and re-enter the sample. If the error appears again, contact Technical Support.
ejectorErrorSensorsHomeOutOff	Ejector Error: Invalid Sensor state. Home and Out position sensors both active.	Restart the unit, eject the sample and re-enter the sample. If the error appears again, contact Technical Support.
ejectorErrorSensorLockOutOff	Ejector Error: Invalid Sensor state. Lock and Out position sensors both active.	Restart the unit, eject the sample and re-enter the sample. If the error appears again, contact Technical Support.
ejectorErrorHomeSens	Ejector Error: Unable to reach Home position. Check for obstructions.	Restart the unit, eject the sample and re-enter the sample. If the error appears again, contact Technical Support.
ejectorErrorLockSensor	Ejector Error: Unable to reach Lock position. Check for obstructions.	Ejector Error: Unable to reach Lock position. Check for obstructions.
ejectorErrorOutSensor	Ejector Error: Unable to reach Out position. Check for obstructions.	Restart the unit, eject the sample and re-enter the sample. If the error appears again, contact Technical Support.
ejectorErrorTubeEjection	Ejector Error: Unable to eject tube. Check for obstructions near ejection area.	Remove source of jam. Sensor will reset once tube is removed. If problem persists, contact Technical Support
ejectorErrorPositionUndefined	Ejector Error: Invalid sensor state. No position sensors active.	Restart the unit, eject the sample and re-enter the sample. If the error appears again, contact Technical Support.
mErrorOtherProcessRunning	Measure Error: Another process is running	Restart the unit, eject the sample and re-enter the sample. If the error appears again, contact Technical Support.

Error Code	UI Text	Solutions
mErrorSampleTypeNotSet	Measure Error: Sample type not set	Restart the unit, eject the sample and re-enter the sample. If the error appears again, contact Technical Support.
mErrorReactorTriggerDelay	Measure Error: Reactor Trigger Delay. Remix sample for additional time and run again.	Remix SEDITROL or Proficiency sample for additional time and run again.
mErrorFlowIn	Measure Error: Flow In. No flow detected when pumping into Reading Cell	Reading cell was unable to detect flow of sample stream during laminar flow. After withdrawal sample did not move into reading cell. Check sufficient sample is in the sample tube and run again. If the error appears again, contact Technical Support.
mErrorFlowOut	Measure Error: Flow Out. No flow detected when discarding sample from Reading Cell	Reading cell was unable to detect flow of sample stream when discarding. Sample still inside the reading cell. Run a Quick Clean and continue with next samples. If error persists, run a Deep Clean and/or contact Technical Support.
wpErrorMovement	Withdrawal Error: Unable to position Inversor	Needle is in error state while attempting to pierce sample. Restart the unit, eject the sample, and re-enter the sample. If the error appears again, contact Technical Support.
wpErrorTailCalibration	Withdrawal Error: Unable to calibrate Tail Sensor. Perform Deep Clean. Otherwise contact Technical Support.	The reading cell is unable to set the correct value for the reading cell tail sensor. Perform Deep Clean. If the error appears again, contact Technical Support.
wpErrorNoTube	Withdrawal Error: No Sample Tube Detected. Measurement Aborted.	During the piercing phase, the sample tube was not detected (probe tube not activated after down sensor reached). Check to make sure sample was loaded. Run sample again.
wpErrorSampleNotInMixSustain	Withdrawal Error: Sample has not performed required mixing cycles.	Process start requested when the sample did not complete the required mixing phase. Restart the unit, eject the sample, and re-enter the sample. If the error appears again, contact Technical Support.
wpErrorOtherProcessRunning	Withdrawal Error: Another process is running preventing withdrawal from executing.	Process start requested when other processes are running (wash process prime process measure process). Restart the unit, eject the sample, and re-enter the sample. If the error appears again, contact Technical Support.

Error Code	UI Text	Solutions
Barcode Error	Barcode Error: If the tube has a barcode, try repositioning it. Select "Run Test" to skip data entry. Remove tube to start manual entry	Unable to read the barcode. If the tube has a barcode, try repositioning it. Select "Run Test" to skip data entry. Remove tube to start manual entry.
washPumpTimeout	Quick Clean Error: Error detected during the movement phases of the process. Pump did not start movement	Error detected during the movement phases of the process. Pump did not start movement. Restart the unit, eject the sample, and re-enter the sample. If the error appears again, contact Technical Support.
washMovement	Quick Clean Error: The pump did not move correctly	The pump did not move correctly. Restart the unit, eject the sample, and re-enter the sample. If the error appears again, contact Technical Support.
washUnableToSetCellT100	Quick Clean Error: The Quick Clean process is unable to set the correct T100 value for the READING CELL	The wash process is unable to set the correct T100 value for the reading cell. Check if there is sufficient iWASH fluid in the bottle. This can occur when the fluid level counter was not reset after the iWASH bottle was replaced. If iWASH is full, contact Technical Support.
washUnableToSetTailT100	Quick Clean Error: The Quick Clean process is unable to set the correct T100 value for the TAIL SENSOR	The wash process is unable to set the correct T100 value for the tail sensor. Check if there is sufficient iWASH fluid in the bottle. This can occur when the fluid level counter was not reset after the iWASH bottle was replaced. If iWASH is full, contact Technical Support.
washUnableToDetectCellEmpty	Quick Clean Error: The Quick Clean process is unable to detect a change on optical value for the READING CELL after the iWASH fluid has been discarded. Check iWASH Bottle is connected and has sufficient fluid available	The wash process is unable to detect a change in optical value for the tail sensor after the washing liquid has been discarded. Check if there is sufficient iWASH fluid in the bottle. This can occur when the fluid level counter was not reset after the iWASH bottle was replaced. If iWASH is full, contact Technical Support.
washUnableToDetectTailEmpty	Quick Clean Error: The Quick Clean process is unable to detect a change in optical value for the TAIL SENSOR after the iWASH fluid has been discarded. Verify tubing is secured in the Tail Sensor	The wash process is unable to detect a change in optical value for the tail sensor after the washing liquid has been discarded. Check if there is sufficient iWASH fluid in the bottle. This can occur when the fluid level counter was not reset after the iWASH bottle was replaced. If iWASH is full, contact Technical Support.
washCellEmitterCurrentTooLow	Quick Clean Error: Current for the Reading Cell Emitter is lower than the allowed lower limit. Contact Technical Support	Current for the reading cell emitter is lower than the allowed lower limit. Contact Technical Support.

Error Code	UI Text	Solutions
washCellEmitterCurrentToHigh	Quick Clean Error: Current for the Reading Cell Emitter is higher than the allowed higher limit. Contact Technical Support	Quick Clean Error: Current for the reading cell emitter is higher than the allowed higher limit. Contact Technical Support
washTailEmitterCurrentToLow	Quick Clean Error: Current for the Tail Sensor Emitter is lower than the allowed lower limit. Contact Technical Support	Current for the tail sensor emitter is lower than the allowed lower limit. Contact Technical Support.
washTailEmitterCurrentToHigh	Quick Clean Error: Current for the Tail Sensor Emitter is higher than the allowed higher limit. Contact Technical Support	Current for the tail sensor emitter is higher than the allowed higher limit. Contact Technical Support.
washOtherProcessRunning	Quick Clean Error: Wash requested when another process is running	Wash requested when another process is running. Restart the unit, eject the sample, and re-enter the sample. If the error appears again, contact Technical Support.
scErrorCardRemoved	Smart Card Error: Card removed. Please try again	Card was removed during operation. Try inserting card again.
scErrorCardType	Smart Card Error: Card Type Error. Please try again	Card model/type is not valid. Try inserting card again. If unsuccessful, contact Technical Support.
scErrorCipherSize	Smart Card Error: Cipher Size Error. Please try again	Card cannot be decrypted due to an invalid memory dimension. Restart unit and try again. If error appears again, contact Technical Support.
scErrorCcommand	Smart Card Error: Command Error. Please try again	System has passed an incorrect command to the smart card controller. Restart unit and try again, If error appears again, contact Technical Support.
scErrorStartAddress	Smart Card Error: Start Address Error. Please try again	System has requested to read/write to a wrong starting address. Restart unit and try again, otherwise contact Technical Support.
scErrorEndAddress	Smart Card Error: End Address Error. Please try again	System has requested to read/write to a wrong ending address. Restart unit and try again. If error appears again, contact Technical Support.
scErrorMemoryRange	Smart Card Error: Memory Range Error. Please try again	System has requested to read/write a section of memory of the wrong dimension. Restart unit and try again. If error appears again, contact Technical Support.
scErrorErasing	Smart Card Error: Erasing Error. Please try again	Error encountered during Smart Card error counter erase operation. Smart Card is still valid. Restart unit and try again. If error appears again, contact Technical Support.
scErrorNonAlcorSmart	Smart Card Error: Non ALCOR Smart Card Error. Please insert a valid card	The inserted Smart Card is not manufactured by ALCOR Scientific. Contact Technical Support.

Error Code	UI Text	Solutions
scErrorPersonalizationIncorrect	Smart Card Error: Personalization Incorrect. Please try again	The inserted Smart Card does not present the same distributor ID present inside the instrument. Smart card will not be loaded nor will be burned. Contact Technical Support.
scErrorProtocolType	Smart Card Error: Protocol Type Error. Please try again	Inserted Smart Card does not use "Asynchronous" protocol. Restart unit and try again. If error appears again, contact Technical Support.
scErrorPscPresentation	Smart Card Error: PSC Presentation Error. Please try again	Error during presentation of Programmable Security Code, Smart Card cannot be written, procedure is aborted. Restart unit and try again. If error appears again, contact Technical Support.
scErrorSizeNotOk	Smart Card Error: Size Not Okay Error. Please try again	Inserted Smart card contains a denomination that is not allowed. Contact Technical Support.
scErrorUnableToWriteEeprom	Smart Card Error: Unable to write EEprom. Please try again	System is unable to store credits into internal memory. Restart unit and try again. If error appears again, contact Technical Support.
scErrorRestoreOriginalAvailability	Smart Card Error: Restore Original Availability Error. Please try again	System is unable to write availability value before smart insertion back to internal memory. Restart unit and try again. If error appears again, contact Technical Support.
scErrorUnhandledRequest	Smart Card Error: Unhandled Request Error. Please try again	Value written to request register not handled by the process. Restart unit and try again. If error appears again, contact Technical Support.
scErrorUnableToClearEeprom	Smart Card Error: Unable to Clear EEprom. Please try again	During transfer, system unable to clear internal memory availability value. Restart unit and try again. If error appears again, contact Technical Support.
scErrorTransferContentInvalid	Smart Card Error: Transfer Content Invalid. Please try again	Content of transfer card is not valid. Restart unit and try again. If error appears again, contact Technical Support.
scErrorCardUsed	Smart Card Error: Card has already been used. Please Insert a valid card	Inserted card has already been used; card should be destroyed. If it is known that card has not been used on any analyzer, contact Technical Support.
scErrorCardCloned	Smart Card Error: Card has been cloned. Please insert a valid card	Inserted card has been already downloaded onto this unit. Card has potentially been cloned. Contact Technical Support.

14.9 Sampling Error Message

In the event of a sampling error, the following messages will be printed:

Error Message (Printed or Logged)	Explanation	Solutions
ESR_ERR_NOFLOW	This error appears when the system is able to withdraw the correct volume from the sample tube but is not able to detect the sample moving in the reading position.	Contact Technical Support.
ESR_ERR_NOSPIKE	Human blood, when flowed into the reading cell, must cause a drop in light transmission. This error indicates the detection of an anomalous sample.	New specimen should be drawn.
ESR_ERR_REVERSE	Usually a hematological sample, after being positioned into the reading cell, starts to form rouleaux (RBC aggregates) which increases light transmission. If the signal detected decreases instead, the error code is provided, indicating a nonstandard condition. It is possible the sample is not of human blood.	New specimen should be drawn.
ESR_ERR_NOPOINTS	This error appears when the reaction takes too much time to develop, or when the drop of the signal of error 3 takes too much time to end. This is an indication of hyperviscosity of the sample or hydraulic malfunctioning.	New specimen should be drawn.
ESR_ERR_TOODARK	Indicates a very high HCT of the sample with a consequential unreliable result. Instead of providing an inaccurate result, the system provides the error message.	New specimen should be drawn.
ESR_ERR_TOOCLEAR	Indicates a very low HCT of the sample with a consequential unreliable result. Instead of providing an inaccurate result, the system provides the error message.	New specimen should be drawn.
ESR_ERR_WITHDRAWAL	System is not able to aspirate the correct volume from the sample tube.	Check for sufficient volume in sample tube. If error appears again, contact Technical Support.
ESR_ERR_FLOW_IN	System unable to move the sample inside the reading cell.	Check for sufficient volume in sample tube. If error appears again, contact Technical Support.
ESR_ERR_FLOW_OUT	System unable to move the sample out of the reading cell.	Contact Technical Support.
ESR_ERR_ACQUISITION	Acquisition measure step unable to be completed.	Contact Technical Support.
ESR_ERR_TRIGGERDELAY	Control material did not start the reaction in the expected amount of time.	Remix sample for additional time and rerun test. If error appears again, contact Technical Support.
ESR_ERR_LOW_CONTROL_HIGH	System detects a high value on a Low control.	Remix sample for additional time and rerun test. If error appears again, contact Technical Support.
ESR_ERR_HIGH_CONTROL_LOW	System detects a low value on a High control.	Remix sample for additional time and rerun test. If error appears again, contact Technical Support.

14.10 Print Out of Sampling Error Message

In the event of a sampling error, the instrument will try to resolve it automatically up to a maximum of three (3) attempts. If after the third attempt, the instrument is unable to resolve the sampling error, an error message will be printed as follows:

=====	
Date: 03/25/2013	Date of analysis
Time: 13:36:24	Time result printed
iSED Sn: 00001	Instrument serial number
ID: 812409	Barcoded sample identification
Error: ESR_ERR_FLOW_IN	
=====	

14.11 Troubleshooting

The following Troubleshooting Chart aims to help diagnose some simple analyzer problems and offer solutions for resolving them.

Situation	Possible Causes	Solutions
Instrument will not power ON	Loose power connections Bad fuse	Check all power connections at the rear of instrument, power supply, and wall outlet. Reconnect power cord at all locations. Wait 30 seconds. Plug back in. Remove fuse cap located directly above power connection on rear of instrument. Check fuse and replace if necessary. See Section 13.6.
Sample tube stuck in the wheel	Tube dropped during sample entry	Power OFF the instrument and manually remove the tube(s) from the wheel.
Touch screen not responding	Touch screen is out of calibration or analyzer is frozen	Power cycle the analyzer and check if problem is resolved. If error appears again, contact Technical Support.
Results are running low/high	Lipemic, hemolyzed, or clotted specimen Pre-analytical sample handling change or system error	Verify condition of specimen. Run controls. If results are within range, resume normal operation; if out of range, discontinue testing and contact Technical Support.
Instrument is not scanning patient barcode	Damaged, incompatible, or no barcode label Barcode reader misaligned	Validate barcode label. Contact Technical Support for instruction.

For troubleshooting issues not covered in this manual, please contact ALCOR Scientific Technical Support or an authorized iSED/iSED ELITE Distributor.

15. Safety Precautions

15.1 General Considerations

WARNING: It is recommended that blood samples be handled wearing gloves and that all other appropriate precautions be taken when dealing with potentially infectious biological material.

CAUTION: The instrument should be disconnected from power supply before performing any cleaning, maintenance, or exposing internal electrical components and circuits.

NOTE: If used in a manner not specified by the manufacturer, damage or injury could result.

15.2 Biological Waste

Biological hazards can be found in all human and animal body fluids and/or tissues. While using the instrument, it is suggested that Good Laboratory Practices are followed. Please refer to and follow all local regulations, department safety guidelines, and biosafety policies for disposal of biohazardous waste.

WARNING: Dispose blood tubes into a biohazard container.

WARNING: Dispose sharps into a biohazard sharps container.

WARNING: All other biohazardous waste should be deposited into a biohazard bag.

WARNING: Biohazard bags should be placed into a Medical Waste Management bin for pick up.

WARNING: Dispose liquid waste container contents in a manner consistent with local regulations and laboratory procedures.

WARNING: Empty Sample Collection Tray when full to avoid a potential biohazard spill due to broken samples.

16. Technical Support

If you experience any problems while operating the instrument, If you experience any problems while operating the instrument, please contact ALCOR Scientific or your local authorized ALCOR Scientific iSED/iSED ELITE Distributor. ALCOR Scientific offers Technical Support Monday through Friday 8:30am-5:00pm EST (excluding all US Federal Holidays). Technical Support can be reached at:

Toll Free: (800) 495.5270 (USA Only) **Fax:** +1 (401) 737.4519

International: +1 (401) 737.3774

Mail: ALCOR Scientific
20 Thurber Blvd
Smithfield, RI 02917
USA

Email: techservice@alcorscientific.com

WARNING: In the event that the instrument must be returned for service, EMPTY ALL FLUID CONTAINERS BEFORE SHIPPING.

WARNING: Remove any liquid waste or on-board sample tubes and decontaminate before returning for service.

Any instrument containing accumulated blood must be cleaned prior to shipment to the manufacturer. This decontamination is required by Federal Law (Title 48 and 49 of the Federal Regulations) in accordance with the Environmental Protection Agency's Regulations for Biohazard Waste Management.

17. Technical Specification

Name of Device	iSED serial numbers 05000 and higher and iSED ELITE
Type of Device	Automated analyzer for the determination of erythrocyte sedimentation rate of human whole blood
Principle of Measure	Photometric Rheoscope
Sample Requirements	500 µL minimum volume
	100 µL aspirated volume
Analytical Range	1-130 mm/hr
Time to Result	Printed; first results available within 20 seconds after pre-programmed mixing time.
Ethernet Port	For manufacturing or network LIS Connection
Serial Port	Serial RS232 DB9 port for legacy LIS connection
Barcode Scanner	Internal
Printer	Internal
Operating Environment	10-30° C, Indoor Use, Pollution Degree – 2
Storage/Transport Environment	-20-60° C
Humidity	15% - 85% (non-condensing)
Power Supply	Transformer: 100-240 VAC 50/60Hz; Device 24VDC, 2.5A
Power Consumption	60W
Frequency	50-60 Hz
Over Voltage Category	Category II
Dimensions (L x W x H)	36 x 27 x 35 cm 14.17 x 10.63 x 13.78 in
Weight	12.1 kgs 26.7 lbs
Operational Altitude**	4000 Meters
Storage Altitude**	4000 Meters
Restrictions	For Professional Use Only
** Rating of 2000 meters for iSED devices serial number 02870 and lower. Rating of 3000 meters for iSED devices with serial number 02871 – 05000. Rating of 4000 meters for iSED devices with serial number 05000 and above and all iSED ELITE devices.	

18. Quick Reference Guide

For convenience, a separate Quick Reference Card was included with the instrument. Should the Quick Reference Card be misplaced, a summary of the information can be found below:

iSED® / iSED® ELITE QUICK REFERENCE CARD

For iSED® analyzers with serial numbers >5000 and all iSED® ELITE analyzers

Sample Requirement

EDTA whole blood
100 µL aspirated volume
(500 µL total volume)

1 µL = 1 microliter

Tube Requirement

EDTA anticoagulant
13 x 75 mm capped tube

Erythrocyte Sedimentation Rate Reference Values (mm/hr)*

Men under 50 years old	< 15
Men over 50 years old	< 20
Women under 50 years old	< 20
Women over 50 years old	< 30

The ranges provided are for reference only. All laboratories should establish their own reference ranges based on the patient population served.

**Keohane EM, Otto CN, Walenga JM. Rodak's Hematology: Clinical Principles and Applications Sixth Edition. Elsevier; 2019.*

Operating Procedure

1. Touch 'Add Sample' icon 
2. The sample wheel will rotate to position the next open slot in the sample entry port.
3. Insert the barcoded tube, making sure the barcode is oriented to the right. A red light will illuminate and a distinctive beep will sound confirming that the barcode has been successfully recognized.

Only load tightly capped tubes.

4. Repeat this sequence until all samples or controls have been entered and/or all positions in the sample wheel are occupied.
5. Sample processing will begin automatically.

Please refer to the iSED®/iSED® ELITE Operator's Manual for the complete operating procedure.

For ALCOR® Scientific Technical Support:
800.495.5270 / 401.737.3774
techservice@alcorscientific.com

222-09-010, Rev. 1

ALCORSCIENTIFIC.COM

19. Warranty Information

Manufacturer's Warranty

ALCOR Scientific warrants that this product is free from defects in materials and workmanship for a period of One (1) Year from the date of original purchase (except as noted below). During the stated one (1) year period, ALCOR Scientific shall, at its sole discretion, repair or replace at no charge to the original end use purchaser or person receiving the product, any product found to be defective due to material or workmanship. In the case of replacement, a new or reconditioned product may be provided at ALCOR Scientific's option.

This warranty is limited to the repair or replacement due to defects in parts or workmanship and shall not include any maintenance and repairs or replacement of parts due to normal wear and tear. Parts required which were not defective shall be replaced at additional costs, and ALCOR Scientific shall not be required to make any repairs or replace any parts which are necessitated by abuse, accident, alteration, misuse, neglect, maintenance by other than ALCOR Scientific or an authorized ALCOR Scientific service agent, or failure to operate the instrument in accordance with instructions. Further, ALCOR Scientific extends no warranty for malfunction or damage to its products resulting from improper or unreasonable use or maintenance; failure to follow operating instructions; connections to improper voltage supply; unauthorized alteration or modification of original condition; damages caused by inadequate packing or shipping procedures; loss of, damage to or corruption of stored data; and any damage due to use of operating supplies other than those manufactured or recommended by ALCOR Scientific.

ALCOR Scientific reserves the right to make changes in design or software of this instrument without obligation to incorporate such changes into previously manufactured instruments.

Disclaimer of Warranties

THIS WARRANTY IS EXPRESSLY MADE IN LIEU OF ANY AND ALL OTHER WARRANTIES EXPRESS OR IMPLIED INCLUDING THE WARRANTIES OF MERCHANTABILITY AND FITNESS FOR USE.

This warranty is void if the label bearing the serial number has been removed or defaced.

Limitations of Liability

In no event shall ALCOR Scientific be liable for indirect, special or consequential damages, even if ALCOR Scientific has been advised of the possibility of such damages.

Any transport costs and risks are not covered by this warranty. In the event that the instrument shall be returned to ALCOR Scientific for servicing, replacement or for other reasons, it must be shipped and received in original packaging. Otherwise, additional charges may be incurred.

Proof of purchase from an authorized ALCOR Scientific distributor and proof of delivery may be required.

EMERGO EUROPE
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands
(+31) 70 345 8570

www.alcorscientific.com

© Copyright 2023, ALCOR Scientific

ALCOR, iSED, iWASH, iWASTE, and SEDITROL are
registered trademarks of ALCOR Scientific

ALCOR Scientific
20 Thurber Boulevard
Smithfield, RI 02917 USA
(T) +1 401.737.3774
WWW.ALCORSCIENTIFIC.COM

iSED (112-00101) Serial Numbers >5000 / iSED ELITE (112-00222) | OPERATOR'S MANUAL

Certificate of Registration of Quality Management System to ISO 13485:2016

The National Standards Authority of Ireland certifies that:

Alcor Scientific LLC
20 Thurber Boulevard
Smithfield, RI 02917
USA

has been assessed and deemed to comply with the
requirements of the above standard in respect of the scope of
operations given below:

**The Design, Manufacture, and Distribution of
Enteral Nutrition Feeding Pumps, Enteral Feeding
Sets, and IVD Systems for Haematology.**

Approved by:
Kevin Mullaney
Director of Certification

Registration Number: CM19.4997
Certification Granted: January 22, 2016
Effective Date: February 25, 2023
Expiry Date: February 24, 2026

National Standards Authority of Ireland, 1 Swift Square, Northwood, Santry, Dublin 9, Ireland T +353 1 807 3800
National Standards Authority of Ireland, 20 Trafalgar Square, Nashua, New Hampshire, NH 03063, USA T +1 603 882 4412