Evolut™ FX System

Evolut™ FX Transcatheter Aortic Valve • Транскатетърна аортна клапа Evolut™ FX • Katetrizačné implantovaná aortální chlopeň Evolut™ FX • Evolut™ FX-transkateter-aortaklap • Evolut™ FX Transkatheter-Aortenklappe • Διακαθετηριακή αορτική βαλβίδα Evolut™ FX • Válvula aórtica transcáteter Evolut™ FX • Transkateteraalne aordiklapp Evolut™ FX • Katetrilla implantoitava Evolut™ FX -aorttaläppä • Valve aortique transcathéter Evolut™ FX • Transkateterski aortni zalistak Evolut™ FX • Evolut™ FX transzkatéteres aortabillentyű • Katup Aorta Transkateter Evolut™ FX • Valvola aortica transcaterere Evolut™ FX • Evolut™ FX транскатетерлік қолқа клапаны • Evolut™ FX 경도관 대동맥판막 • „Evolut™ FX“ transkateterinis aortos vožtuvas • Evolut™ FX transkatetra aortas vārstulis • Транскатетерски аортен залисток Evolut™ FX • Evolut™ FX-transkatheter-aortaklep • Evolut™ FX transkateter aortaklaff • Zastawka aortalna do implantacji przezcewnikowej Evolut™ FX • Válvula aórtica transcáteter Evolut™ FX • Válvula aórtica transcáteter Evolut™ FX • Valvă aortică transcáteter Evolut™ FX • Аортальный клапан для транскатетерной установки Evolut™ FX • Transkatétrová náhrada aortálnej chlopne Evolut™ FX • Transkatetriska aortna zaklopka Evolut™ FX • Evolut™ FX transkateterski aortni zalistak • Evolut™ FX kateterburen aortaklaff • Evolut™ FX Transkateter Aortik Карак • Аортальный клапан для транскатетерної установки Evolut™ FX • Van động mạch chủ qua đường ống thông Evolut™ FX

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1 Device description

The Medtronic Evolut FX system is a recapturable transcatheter aortic valve implantation system, which includes the Evolut FX transcatheter aortic valve, the delivery catheter system, and the loading system.

1.1 Evolut FX transcatheter aortic valve (bioprosthesis)

The bioprosthesis is designed to replace the native, surgical bioprosthetic, or transcatheter bioprosthetic aortic heart valve without open heart surgery and without concomitant surgical removal of the failed valve. The support frame is manufactured from nitinol, which has multilevel, self-expanding properties and is radiopaque. The support frame includes 3 gold radiopaque markers approximately 3 mm from the inflow portion of the TAV frame, below the 3 tissue commissure pads. The bioprosthesis is manufactured by suturing valve leaflets and an inner skirt from porcine pericardium into a tri-leaflet configuration. The bioprosthesis has a porcine pericardial tissue outer skirt (wrap), which is 1.5 cells in height and is sutured to the inflow section of the bioprosthesis.

Table 1. Heart valve materials

Component	Materials	Mass ^a
Tissue	Processed porcine pericardium	1.3 g
Frame	Nitinol (a nickel titanium alloy)	1.3 g
Suture	Ultra-high-molecular-weight polyethylene (UHMWPE) ^b	0.1 g
Radiopaque markers	Gold	0.005 g

^a Mass is based on the material in the largest size valve.

^b The Evolut FX 23 mm valve also uses expanded polytetrafluoroethylene (ePTFE).

Residuals – The Evolut FX heart valve is sterilized with and stored in glutaraldehyde solution. Before the procedure, the heart valve is rinsed, and the glutaraldehyde residuals are reduced to safe levels.

The bioprosthesis is processed with alpha-amino oleic acid (AOA[™]), which is a compound derived from oleic acid, a naturally occurring long-chain fatty acid. AOA is an antimineralization treatment shown to reduce both early and late valvular calcification. The bioprosthesis is available for a range of aortic annulus diameters (*Table 3*).

1.2 Delivery catheter system (catheter)

The catheter facilitates the placement of the bioprosthesis within the annulus of the aortic valve. The catheter assembly is flexible and compatible with a 0.035 in (0.889 mm) guidewire. The distal (deployment) end of the system features an atraumatic, radiopaque catheter tip and a capsule that covers and maintains the bioprosthesis in a crimped position. The capsule includes a distal flare to enable the bioprosthesis to be partially or fully recaptured after partial deployment. A green stability layer is fixed at the handle and extends down the outside of the catheter shaft. It provides a barrier between the retractable catheter and the introducer sheath and vessel walls, thus enabling the catheter to retract freely. An Evolut FX inline sheath is assembled over the stability layer, which functions as a hemostatic introducer sheath and minimizes the access site size to the capsule diameter. The 23-29 mm catheter model is compatible with an 18 Fr (6 mm) introducer sheath. The 34 mm catheter model is compatible with a 22 Fr (7.33 mm) introducer sheath.

Refer to *Table 2* for system compatibility. Refer to *Figure 2* and *Figure 3* for catheter components.

The delivery catheter system consists of a catheter with an integrated handle to provide the user with accurate and controlled deployment. The handle is on the proximal end of the catheter and is used to load, deploy, recapture, and reposition the bioprosthesis. The handle features a gray front grip used to stabilize the system. The deployment knob turns to deploy the bioprosthesis precisely. Arrows on the deployment knob indicate the direction of rotation required to deploy the bioprosthesis. If desired, the deployment knob can be turned in the opposite direction to partially or fully recapture the bioprosthesis if the radiopaque capsule marker band has not yet reached the distal end of the radiopaque paddle attachment. Once the radiopaque capsule marker band reaches the distal end of the radiopaque paddle attachment, it is at the point of no recapture. The deployment knob also features a trigger, which can be engaged to make macro adjustments to the capsule position. A dark blue hand rest connects to the deployment knob. The end of the handle features a tip-retrieval mechanism, which can be used to withdraw the catheter tip to meet the capsule after the device has been fully deployed.

The catheter packaging contains an integrated loading bath and a removable tray with 3 rinsing bowls for loading and rinsing the bioprosthesis. The integrated loading bath features a mirror, which aids in accurate placement of the bioprosthesis frame paddles during loading. In addition to these features, the device packaging is swiveled and secured to facilitate the bioprosthesis loading procedure (*Figure 4* and *Figure 5*).

Residuals – The Evolut FX delivery catheter system is sterilized with ethylene oxide. Before the procedure, the ethylene oxide residuals are reduced to safe levels for the indicated population.

1.3 Loading system

The loading system compresses the bioprosthesis into the catheter. The loading system comes in different models. Refer to *Table 2* for system compatibility. Refer to *Figure 6* for loading system components.

Residuals – The Evolut FX loading system is sterilized with ethylene oxide. Before the procedure, the ethylene oxide residuals are reduced to safe levels for the indicated population.

1.4 System compatibility

Table 2. System compatibility

Bioprosthesis model	Compatible loading system models	Compatible catheter models
EVOLUTFX-23	L-EVOLUTFX-2329	D-EVOLUTFX-2329
EVOLUTFX-26	L-EVOLUTFX-2329	D-EVOLUTFX-2329
EVOLUTFX-29	L-EVOLUTFX-2329	D-EVOLUTFX-2329
EVOLUTFX-34	L-EVOLUTFX-34	D-EVOLUTFX-34

Note: Verify the correct delivery system and loading system are used for the selected valve. Valves are designed for compatibility with specific delivery and loading systems.

2 Intended purpose

The intended purpose of the Evolut FX system is to replace the function of a stenosed native or failed surgical or transcatheter aortic valve without the need for open heart surgery.

2.1 Indications

The Evolut FX system is indicated for adult patients presenting with severe native aortic valve stenosis. For patients presenting with severe native bicuspid aortic valve stenosis, the Evolut FX system is indicated for patients who are at intermediate or greater risk for surgical aortic valve replacement (AVR) where intermediate risk is defined as Society of Thoracic Surgeons operative risk score $\geq 4\%$ or documented heart team agreement of risk for AVR due to frailty or comorbidities. For patients presenting at low risk for AVR ($< 4\%$), the system is indicated for patients ≥ 70 years of age with an LVEF $> 30\%$.

The Evolut FX system is also indicated for patients with a stenosed, insufficient, or combined surgical or transcatheter bioprosthetic valve failure necessitating valve replacement who are at high or greater risk for surgical aortic valve replacement (AVR) where high risk is defined as Society of Thoracic Surgeons operative risk score $\geq 8\%$ or documented heart team agreement of risk for AVR due to frailty or comorbidities.

Patients must present with anatomical dimensions as described in *Section 2.3*.

2.2 Contraindications

The Evolut FX system is contraindicated in patients who cannot tolerate the device materials listed in *Table 1*, an anticoagulation/antiplatelet regimen, or who have active bacterial endocarditis or other active infections.

2.3 Patient target group

This device is intended for adult patients with heart disease due to severe native calcific aortic stenosis or a stenosed or insufficient surgical or transcatheter bioprosthetic valve. Patients are judged by a heart team to be at low, intermediate, high, or extreme risk for open surgical therapy. The bioprosthesis is available for a range of patient anatomical criteria (*Table 3*).

Sizing recommendations for the Evolut FX bioprosthesis in a native annulus or failed surgical or transcatheter bioprosthesis are provided in *Table 3*.

Table 3. Patient anatomical criteria

Bioprosthesis model	Aortic annulus diameter ^a	Aortic annulus perimeter ($\pi \times$ aortic annulus diameter) ^a	Bioprosthesis size
EVOLUTFX-23	18 mm to 20 mm	56.5 mm to 62.8 mm	23 mm
EVOLUTFX-26	20 mm to 23 mm	62.8 mm to 72.3 mm	26 mm
EVOLUTFX-29	23 mm to 26 mm	72.3 mm to 81.7 mm	29 mm

Table 3. Patient anatomical criteria (continued)

Bioprosthesis model	Aortic annulus diameter ^a	Aortic annulus perimeter ($\pi \times$ aortic annulus diameter) ^a	Bioprosthesis size
EVOLUTFX-34	26 mm to 30 mm	81.7 mm to 94.2 mm	34 mm

Valve size recommendations are based on sizing by computed tomography (CT). Patient anatomical factors and multiple imaging modalities should be considered during valve size selection. Consider all risks associated with improper sizing and mismatch between the patient and bioprosthesis.

^a For TAV in SAV and TAV in TAV, diameter and perimeter criteria are applicable to the failed SAV or TAV measured inner diameter.

Patients must present with the vessel access criteria listed in *Table 4*.

Table 4. Patient vessel access criteria

Bioprosthesis model	Delivery catheter system model	Capsule OD	Arterial access	Subclavian access with an IMA graft	Direct aortic access
EVOLUTFX-23, EVOLUTFX-26, EVOLUTFX-29,	D-EVOLUTFX-2329	18 Fr (6 mm)	≥ 5 mm	≥ 5.5 mm	≥ 60 mm from basal plane
EVOLUTFX-34	D-EVOLUTFX-34	22 Fr (7.33 mm)	≥ 6 mm	≥ 6.5 mm	

The risks of long-term anticoagulation and/or antiplatelet therapy should be considered. Implantation of the bioprosthesis is not recommended in patients with aortic root angulation (the angle between the plane of the aortic valve annulus and the horizontal plane/vertebrae) of >30° for right subclavian/axillary access or >70° for femoral and left subclavian/axillary access.

Use caution in patients with a preexisting patent Left Internal Mammary Artery/Right Internal Mammary Artery (LIMA/RIMA) graft. For direct aortic access, ensure the access site and trajectory are free of patent RIMA or a pre-existing patent RIMA graft.

For direct aortic and subclavian access procedures, care must be exercised when using the tip-retrieval mechanism to ensure adequate clearance to avoid advancement of the catheter tip through the bioprosthesis leaflets during device closure.

For direct aortic access procedures, use a separate introducer sheath; do not use the Evolut FX inline sheath. Maintain the Evolut FX inline sheath at the proximal end of the catheter throughout the procedure.

For transfemoral access, use caution in patients who present with multiplanar curvature of the aorta, acute angulation of the aortic arch, an ascending aortic aneurysm, or severe calcification in the aorta and/or vasculature. If ≥2 of these factors are present, consider an alternative access route to prevent vascular complications.

2.4 Intended users

The Evolut FX system loading and implant procedure should be performed only by individuals who have received training in the preparation and use of the Evolut FX system.

2.5 Intended use environment

The Evolut FX system is intended to be used in a catheterization lab or in a hybrid operating room environment.

A catheterization lab is a sterile environment which typically contains a variety of equipment, including an operating table, a panel of monitors allowing individuals to see different information simultaneously, and fluoroscopy instrumentation.

A hybrid operating room allows for provision of both surgical and intravascular procedures by combining the imaging equipment needed for catheter procedures with the equipment needed for surgical procedures. This often includes a ceiling-mounted track and C-arm to assist with providing imagery during the procedure.

2.6 Valve-in-valve considerations

Implanting the Evolut FX bioprosthesis in a degenerated surgical aortic bioprosthesis (transcatheter aortic valve in surgical aortic valve [TAV in SAV]) should generally be avoided in patients with any of the following conditions:

- The degenerated bioprosthesis presents with a significant concomitant paravalvular leak (between prosthesis and native annulus), is not securely fixed in the native annulus, or is not structurally intact (for example, wireframe frame fracture).
- The degenerated bioprosthesis presents with a partially detached leaflet that in the aortic position may obstruct a coronary ostium.

Before implanting the Evolut FX bioprosthesis in a degenerated transcatheter aortic bioprosthesis (transcatheter aortic valve in transcatheter aortic valve [TAV in TAV]), additional factors regarding failed valve size and patient anatomy must be considered in order to ensure patient safety (for example, to avoid coronary obstruction). The potential need for future coronary access should be considered. TAV in TAV implantation should be avoided in the following conditions:

- The degenerated TAV presents with a significant concomitant paravalvular leak (between the prosthesis and the native annulus),
- The degenerated TAV is not securely fixed in the native annulus, or is not structurally intact (for example, frame fracture) or
- The risk of coronary obstruction or sinus sequestration after Evolut FX bioprosthesis implantation is high.

2.7 Bicuspid patient considerations

If the patient presents with a bicuspid aortic valve, the heart team should consider the age of the patient and the need for ascending aorta intervention when determining the appropriate treatment option for the patient.

2.8 Clinical information

The Summary of Safety and Clinical Performance (SSCP) can be found at <https://ec.europa.eu/tools/eudamed>. Search for the SSCP using the manufacturer and device name, and any of the following elements, as applicable: device model, reference number, catalog number, or the Basic Unique Device Identification (Basic UDI-DI) number — 0763000B00011727K.

2.9 Clinical benefits

The clinical benefits of the Evolut FX system are improved quality of life through the relief of symptoms related to severe aortic stenosis, symptoms related to regurgitation from a failed surgical or transcatheter valve, or symptoms related to both.

2.10 Performance characteristics

The Evolut FX system is designed to relieve aortic stenosis in dysfunctional native valves or to relieve aortic stenosis or aortic regurgitation in failed surgical or transcatheter aortic bioprosthetic valves. The Evolut FX bioprosthetic valve is delivered on a catheter eliminating the need for open heart surgery. The bioprosthetic valve functionally replaces the dysfunctional native or failed surgical or transcatheter valve.

2.11 Expected device lifetime

The device lifetime of an Evolut FX heart valve varies from patient to patient. At some point, the valve may need to be replaced. The expected lifetime of an Evolut FX heart valve is projected to be 10 years for many patients.

3 Warnings and precautions

Carefully read all warnings, precautions, and instructions for use for all components of the system before use. Failure to read and follow all instructions or failure to observe all stated warnings could cause serious injury or death to the patient.

3.1 Warnings

- This procedure should only be performed where emergency aortic valve surgery can be performed promptly.
- Carefully verify that the relevant patient anatomical parameters are within the specifications set forth in *Section 2.3*.
- The heart team should consider the need for revascularization when selecting a treatment option for aortic valve replacement (AVR).
- Ethylene oxide residuals are known to meet acceptable levels only for the indicated population.
- Do not freeze the bioprosthesis. Freezing causes damage to the bioprosthesis. Check the freeze indicator located inside the labeled carton. If the product has been exposed to a freeze-thaw condition, the indicator vial will break, which causes dye to escape and stain the paper backing. Do not use the bioprosthesis if this or any other sign of freezing is observed.
- Store the bioprosthesis at room temperature.
- Do not handle or manipulate the bioprosthesis with sharp or pointed objects.
- Package contents should not be opened until implantation and sizing are certain.
- Package contents are supplied sterile. Do not use if package is damaged or unintentionally opened.
- The outside of the bioprosthesis container and the outer surfaces of the catheter pouch and the loading system outer pouch are nonsterile and must not be placed in the sterile field.
- This device was designed for single use only. Do not reuse, reprocess, or resterilize the device. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device or create a risk of contamination of the device, which could result in patient injury, illness, or death.
- If a bioprosthesis and catheter have been removed from a patient, dispose of both the bioprosthesis and catheter; do not attempt to reuse either component. Both the bioprosthesis and catheter must be replaced with new sterile components.
- Do not use products if the use-by date printed on the product packaging label has elapsed.
- Carefully remove the catheter from the packaging to prevent it from kinking.
- Perform the **Bioprosthesis rinsing procedure** (*Section 8.3*) before loading the bioprosthesis into the catheter.

- After the bioprosthesis has been loaded into the catheter, maintain the distal end of the catheter with the bioprosthesis in cold (0 °C to 8 °C [32 °F to 46 °F]) sterile saline until implantation.
- Use the deployment knob to deploy and recapture the bioprosthesis. Do not use the trigger for deploying or recapturing because it could cause inaccurate placement of the bioprosthesis.
- System failure could occur if an incorrect combination of devices is used. Refer to the system compatibility table in *Section 1.4*.
- Mechanical failure of the delivery catheter system and/or accessories may result in patient complications.
- Calcific degeneration could cause accelerated structural deterioration of the valve in the following individuals:
 - Children, adolescents, or young adults
 - Patients with altered calcium metabolism (for example, chronic renal failure, or hyperparathyroidism)

3.2 Patient selection considerations

The safety and effectiveness of the Medtronic Evolut FX system have not been evaluated in the pediatric population.

The safety and effectiveness of the bioprosthesis for aortic valve replacement have not been evaluated in the following patient populations:

- Patients with a pre-existing prosthetic heart valve with a rigid support structure in either the mitral or pulmonic position if either the pre-existing prosthetic heart valve could affect the implantation or function of the bioprosthesis or the implantation of the bioprosthesis could affect the function of the pre-existing prosthetic heart valve
- Patients with liver failure (Child-Pugh Class C)
- Patients with cardiogenic shock manifested by low cardiac output, vasopressor dependence, or mechanical hemodynamic support
- Patients who are pregnant or breastfeeding
- Patients with a congenital unicuspid valve
- Patients with prohibitive left ventricular outflow tract calcification
- Patients with sinus of Valsalva anatomy that would prevent adequate coronary perfusion
- Patients with significant aortopathy requiring ascending aortic replacement
- Patients who do not meet the criteria for symptomatic severe native aortic stenosis as defined by the following:
 - Aortic valve area of $\leq 1.0 \text{ cm}^2$ (or aortic valve area index of $\leq 0.6 \text{ cm}^2/\text{m}^2$) by the continuity equation, and mean gradient $\geq 40 \text{ mmHg}$ or maximal aortic valve velocity $\geq 4.0 \text{ m/s}$ by resting echocardiogram, low-dose dobutamine, or exercise stress echocardiography
- Patients who do not meet the criteria for asymptomatic severe native aortic stenosis as defined by the following:
 - Very severe aortic stenosis with an aortic valve area of $\leq 1.0 \text{ cm}^2$ (or aortic valve area index of $\leq 0.6 \text{ cm}^2/\text{m}^2$), and maximal aortic velocity $\geq 5.0 \text{ m/s}$, or mean gradient $\geq 60 \text{ mmHg}$ by transthoracic echocardiography at rest, or
 - Aortic valve area of $\leq 1.0 \text{ cm}^2$ (or aortic valve area index of $\leq 0.6 \text{ cm}^2/\text{m}^2$), and a mean gradient $\geq 40 \text{ mmHg}$, or maximal aortic valve velocity $\geq 4.0 \text{ m/s}$ by transthoracic echocardiography at rest, and an exercise tolerance test that demonstrates a limited exercise capacity, abnormal BP response, or arrhythmia, or
 - Aortic valve area of $\leq 1.0 \text{ cm}^2$ (or aortic valve area index of $\leq 0.6 \text{ cm}^2/\text{m}^2$), and mean gradient $\geq 40 \text{ mmHg}$, or maximal aortic valve velocity $\geq 4.0 \text{ m/s}$ by transthoracic echocardiography at rest, and a left ventricular ejection fraction $< 50\%$

The safety and effectiveness of the bioprosthesis for aortic valve replacement have not been evaluated in patient populations presenting with the following conditions:

- Blood dyscrasias as defined: leukopenia ($\text{WBC} < 1000 \text{ cells/mm}^3$), thrombocytopenia (platelet count $< 50,000 \text{ cells/mm}^3$), history of bleeding diathesis or coagulopathy, or hypercoagulable states
- Mixed native aortic valve disease (aortic stenosis and aortic regurgitation with predominant aortic regurgitation [3–4+])
- Moderate to severe (3–4+) or severe (4+) mitral or severe (4+) tricuspid regurgitation
- Hypertrophic obstructive cardiomyopathy
- New or untreated echocardiographic evidence of intracardiac mass, thrombus, or vegetation
- Moderate to severe mitral stenosis
- Severe ventricular dysfunction with left ventricular ejection fraction (LVEF) $< 20\%$
- Symptomatic carotid or vertebral artery disease

Current heart valve replacement guidelines recommend bioprosthetic valves for patients at 65 years of age or above. Consider the applicable clinical guidelines for valve replacement when selecting the appropriate valve type for each patient. Long-term durability is still the subject of ongoing clinical research.

3.3 Implantation precautions

- Do not use the bioprosthesis if there is any damage to the container (for example, cracked jar or lid, leakage, particulate material, broken or missing seals or jar lid gasket).
- Before catheter insertion, remove the loading stylet.
- If a misload is detected during fluoroscopic (cine mode) inspection, do not attempt to reload the bioprosthesis. Discard the entire system. The valve, catheter, loading system, loading tray, and saline all must be replaced with new sterile components. A misload is defined as one or more of the following:
 - Inflow crown overlap (non-uniform shadow starting at the inflow) that has not ended before the 4th node from the inflow.
 - Outflow crown misalignment and/or not parallel to the paddle attachment.
 - Curved or bent capsule.
 - Direct load as detailed in *Section 8.4, Step 17*.
 - Shadow or outline in outflow indicating a bent strut.
- Inflow crown overlap that has not ended before the 4th node within the capsule, increases the risk of an infold upon deployment in constrained anatomies, particularly with moderate-severe levels of calcification and/or bicuspid condition.
 - Do not attempt to direct load the valve (for example, loading the valve without completing *Step 17* in *Section 8.4* and simply advancing the capsule to load the valve). This increases the likelihood of excessive inflow crown overlap. If a valve has been direct loaded, discard the entire system. The valve, catheter, loading system, loading tray, and saline all must be replaced with new sterile components.
- Before catheter insertion, the heart team determines whether to perform balloon aortic valvuloplasty (BAV) of the native aortic valve to facilitate placement of the bioprosthesis.
- While the catheter is in the patient, ensure the guidewire is extending from the proximal end of the catheter. Do not remove the guidewire from the catheter while the catheter is inserted in the patient.
- There will be some resistance when the catheter is advanced through the vasculature. If there is a significant increase in resistance, stop advancement and investigate the cause of the resistance (for example, magnify the area of resistance) before proceeding. Do not force passage. Forcing passage could increase the risk of vascular complications (for example, vessel dissection or rupture).
- For direct aortic access procedures, use caution when advancing and maintaining the introducer sheath position to avoid tissue and vascular injury and unacceptable access site bleeding.
- From annular contact (or contact with the failed valve, for valve-in-valve procedures) to just before the point of no recapture, the bioprosthesis will occlude cardiac output. Promptly deploy or recapture the valve during this occlusive phase as prolonged obstruction or occlusion of blood flow may lead to hypotension, bradycardia, conduction disturbance, congestive heart failure, pulmonary edema, or death.

3.4 Repositioning precautions

- If the radiopaque capsule marker band has not yet reached the distal end of the radiopaque paddle attachment, the bioprosthesis can be recaptured or repositioned. During deployment, the deployment knob provides a tactile indication as a notification before the point of no recapture.
- Infold detection steps are outlined in *Section 8.5.2*. An observation of any inward fold or crease in the valve, extending from the inflow, identified as a dark line under fluoroscopic (cine mode) inspection, may indicate an infold. If identified, and if the patient's condition allows, do not proceed and do not release the valve.
 - Recapture, remove and discard the entire system. The valve, catheter, loading system, loading tray, and saline all must be replaced with new sterile components.
 - Predilatation is strongly recommended prior to subsequent implantation attempts to minimize infold risk.
 - If initial predilatation does not prevent infolding, reassess valve sizing in the presence of complex anatomies.
 - If an infold is detected and the valve is removed, consider a slightly lower depth of implantation of the second valve to provide additional space for frame expansion.
- Implanting a valve with an unresolved infold increases the risk of PVL and need for post implant dilatation, which is associated with higher rates of adverse events such as dislodgement and dissection.

Note: Predilatation may confer some risk to the patient (for example, liberation of embolic debris, damage to the tissue, or perforation of the aortic root). Patient anatomical characteristics (for example, bicuspid anatomy, excessive or asymmetric leaflet calcification, and possible leaflet fusion) should be considered by the heart team when evaluating and determining the risk/benefit of predilatation and treatment plan for each patient. For TAV in TAV cases, the characteristics of the failed TAV (for example, under-expansion, depth of index implant) should be considered by the heart team when evaluating and determining the risk/benefit of predilatation and the treatment plan for each patient.

- Once the radiopaque capsule marker band reaches the distal end of the radiopaque paddle attachment (point of no recapture), retrieval of the bioprosthesis from the patient (for example, use of the catheter) is not recommended. Retrieval after the point of no recapture may cause mechanical failure of the delivery catheter system, aortic root damage, coronary artery damage,

myocardial damage, vascular complications, prosthetic valve dysfunction (including device malposition), embolization, stroke, and/or emergent surgery.

- During deployment, the bioprosthesis can be advanced or withdrawn as long as annular contact (or contact with the failed valve, for valve-in-valve procedures) has not been made. Once annular contact is made, the bioprosthesis cannot be advanced in the retrograde direction; recapture until the bioprosthesis is free from annular contact, and then reposition in the retrograde direction. If necessary, and the radiopaque capsule marker band has not yet reached the distal end of the radiopaque paddle attachment, the bioprosthesis can be withdrawn (repositioned) in the antegrade direction. However, use caution when moving the bioprosthesis in the antegrade direction.

Note: For TAV in SAV or TAV in TAV procedures, fully recapture the bioprosthesis before repositioning to reduce the risk of hang up or snagging on the failed bioprosthesis and to aid with positioning using the radiopaque capsule marker band.

Caution: Use the handle of the delivery system to reposition the bioprosthesis. **Do not** use the outer catheter shaft.

- Physicians should use judgment when considering repositioning a fully deployed bioprosthesis (for example, using a snare, balloon, and/or forceps). Repositioning the bioprosthesis is not recommended except in cases where imminent serious harm or death is possible (for example, coronary occlusion). Repositioning of a deployed valve may cause aortic root damage, coronary artery damage, myocardial damage, vascular complications, prosthetic valve dysfunction (including device malposition), embolization, stroke, and/or emergent surgery.
- Do not attempt to retrieve or to recapture a bioprosthesis if any one of the outflow struts is protruding from the capsule. If any one of the outflow struts has deployed from the capsule, the bioprosthesis must be released from the catheter before the catheter can be withdrawn.

3.5 Postimplantation precautions

- Ensure the capsule is closed before catheter removal.
- When using a separate introducer sheath, if increased resistance is encountered when removing the catheter through the introducer sheath, do not force passage. Increased resistance may indicate a problem and forced passage may result in damage to the device and/or harm to the patient. If the cause of resistance cannot be determined or corrected, remove the catheter and introducer sheath as a single unit over the guidewire, and inspect the catheter and confirm that it is complete.
- Exercise care when crossing an implanted bioprosthesis with the catheter or other devices.
- Clinical long-term durability has not been established for the bioprosthesis. Evaluate bioprosthesis performance as needed during patient follow-up.

3.6 Post-implant balloon dilatation considerations

If valve function or sealing is impaired due to excessive calcification or incomplete expansion, a post-implant balloon dilatation (PID) of the bioprosthesis may improve valve function and sealing. If the heart team determines that balloon dilatation is appropriate, consider all of the following factors when selecting the dilatation parameters to ensure patient safety:

- Balloon model
- Balloon size
- Balloon position
- Inflation pressure
- Patient anatomy

Two primary factors must be considered when selecting a maximum balloon diameter for post-implant balloon dilatation:

- To mitigate trauma to the annulus
 - A compliant or semi-compliant balloon (for example, B. Braun Z-Med I™* / Z-Med II™*, InterValve V8™*) should not exceed the diameter of the native aortic annulus. For TAV in SAV or TAV in TAV, the balloon should not exceed the inner diameter of the surgical or transcatheter bioprosthetic valve.
 - A non-compliant balloon (for example, Bard TRUE™* Dilatation) should be at least 1 mm smaller than the diameter of the native aortic annulus. For TAV in SAV or TAV in TAV, the balloon should be at least 1 mm smaller than the inner diameter of the surgical or transcatheter bioprosthetic valve.
- To mitigate trauma to the Evolut TAV bioprosthetic leaflets
 - The maximum balloon size chosen for dilatation using a compliant or semi-compliant balloon should not exceed the TAV waist diameter beyond the level set forth in *Table 5* with an applied inflation pressure of no greater than 2 atm.
 - The maximum balloon size chosen for dilatation using a non-compliant balloon should not exceed 1 mm more than the TAV waist diameter with an applied inflation pressure of no greater than 2 atm (see *Table 5*).

Table 5. Post-implant balloon dilatation sizing

Evolut FX size	23 mm			26 mm				29 mm				34 mm				
Native annulus (failed SAV or TAV inner) diameter (in mm)	18	19	20	20	21	22	23	23	24	25	26	26	27	28	29	30
TAV waist diameter (in mm)	20	20	20	22	22	22	22	23	23	23	23	24	24	24	24	24
Maximum balloon diameter (in mm) for compliant and semi-compliant balloons at 2 atm	18	19	20	20	21	22	23	23	24	25	26	26	27	28	28	28
Maximum balloon diameter (in mm) for non-compliant balloons at 2 atm	17	18	19	19	20	21	22	22	23	24	24	25	25	25	25	25

Caution: Overexpansion of the narrowest portion (waist) of the Evolut FX TAV beyond the levels set forth in *Table 5* has been demonstrated through bench data to cause damage to the bioprosthetic leaflets. Complaints of damage to the bioprosthetic leaflets during post-implant balloon dilatation have been reported in some clinical cases, resulting in moderate to severe aortic insufficiency, which may be detected acutely or during follow-up.

It is important to note that the mechanical compliance properties of the selected balloon influence the dilatation dynamics.

The maximum balloon sizes in *Table 5* are derived from bench testing based on a single Evolut FX TAV dilatation to 2 atm. Multiple dilatations of the Evolut FX TAV increases the risk of damage to the bioprosthetic leaflets.

Balloons should not be inflated beyond 2 atm of applied pressure.

Compliant and semi-compliant (softer) balloons will more readily conform to the hourglass profile of the TAV bioprosthesis at lower pressures, but must be inflated at pressures that preserve the hourglass profile of the TAV.

Conversely, non-compliant (stiffer) balloons will achieve the nominal diameter during inflation irrespective of the underlying annulus or TAV resistance and should be downsized (see *Table 5*).

For additional instructions on the use of balloon catheter devices refer to the specific balloon catheter manufacturer's labeling.

In the event that larger balloon diameters than those listed in *Table 5* are required to expand the Evolut FX TAV due to clinically important residual aortic regurgitation or stenosis, using “bailout” intraventricular balloon positioning when performing PID avoids expansion of the narrowest portion (waist) of the Evolut FX TAV. This can mitigate the risk of leaflet damage. Dilatation with intraventricular balloon positioning should be performed with caution in the setting of a smaller ventricle cavity, presence of LVOT calcification, or wire positioning that interferes with mitral valve function, in order to avoid any unintended balloon interaction with anatomy. The balloon's length and diameter, along with the individual patient anatomy, must be considered. Care should also be taken not to exceed the annular diameters when performing PID with intraventricular balloon positioning (see *Table 5*).

In the event that a bailout PID with intraventricular balloon positioning is performed, the nominal diameter of the balloon should not exceed the annular diameter when using compliant or semi-compliant balloons; the nominal diameter of the balloon should be at least 1 mm smaller than the annular diameter when using non-compliant balloons.

4 Potential adverse events

Adverse events that may be associated with the use of the Evolut FX system include, but may not be limited to, the following:

- Death
- Myocardial infarction, cardiac arrest, cardiogenic shock, cardiac tamponade
- Coronary occlusion, obstruction, or vessel spasm (including acute coronary closure)
- Cardiovascular injury (including rupture, perforation, tissue erosion, or dissection of vessels, ascending aorta trauma, ventricle, myocardium, or valvular structures that may require intervention)
- Emergent surgical or transcatheter intervention (for example, coronary artery bypass, heart valve replacement, valve explant, percutaneous coronary intervention [PCI], balloon valvuloplasty)
- Prosthetic valve dysfunction (regurgitation or stenosis) due to fracture; bending (out-of-round configuration) of the valve frame; underexpansion of the valve frame; calcification; pannus; leaflet wear, tear, prolapse, or retraction; poor valve coaptation; suture breaks or disruption; leaks; mal-sizing (prosthesis-patient mismatch); malposition (either too high or too low)/malplacement
- Prosthetic valve migration/embolization
- Prosthetic valve endocarditis
- Prosthetic valve thrombosis
- Delivery catheter system malfunction resulting in the need for additional re-crossing of the aortic valve and prolonged procedural time
- Delivery catheter system component migration/embolization
- Stroke (ischemic or hemorrhagic), transient ischemic attack (TIA), or other neurological deficits

- Individual organ (for example, cardiac, respiratory, renal [including acute kidney failure]) or multi-organ insufficiency or failure
- Major or minor bleeding that may require transfusion or intervention (including life-threatening or disabling bleeding)
- Vascular access-related complications (for example, dissection, perforation, pain, bleeding, hematoma, pseudoaneurysm, irreversible nerve injury, compartment syndrome, arteriovenous fistula, stenosis)
- Mitral valve regurgitation or injury
- Conduction system disturbances (for example, atrioventricular node block, left-bundle branch block, asystole), which may require a permanent pacemaker
- Infection (including septicemia)
- Hypotension or hypertension
- Hemolysis
- Peripheral ischemia
- Bowel ischemia

General surgical risks applicable to transcatheter aortic valve implantation:

- Abnormal lab values (including electrolyte imbalance)
- Allergic reaction to antiplatelet agents, contrast medium, or anesthesia
- Exposure to radiation through fluoroscopy and angiography
- Permanent disability

If a serious incident related to the device occurs, immediately report the incident to Medtronic and the applicable competent authority.

5 Patient counseling information

Physicians should review the instructions for use for applicable information to be shared with the patient. A patient implant card is included in the device package and contains identifying information about the implanted device. After device implant, complete the patient implant card and provide it to the patient before they are discharged. Patients can access information about their device on the website that is listed on their patient implant card.

Physicians should instruct patients to do the following:

- Carry their implant card with them.
- Always inform any healthcare personnel that they have an implanted device before any procedure has begun.
- Contact their physician if they notice any new or changing symptoms.

6 How supplied

6.1 Bioprosthesis

The bioprosthesis is sterilized with and stored in glutaraldehyde solution. The bioprosthesis is supplied sterile and nonpyrogenic in a jar with a screw cap. A freeze indicator is placed inside the labeled carton. If the freeze indicator has been activated, do not use the bioprosthesis and dispose of the bioprosthesis in accordance with applicable laws, regulations, and hospital procedures.

6.2 Catheter

The catheter is packaged in a single-pouch configuration. The catheter is sterile if the package is undamaged and unopened. The catheter is sterilized by ethylene oxide. Before the procedure, the ethylene oxide residuals are reduced to safe levels for the indicated population.

6.3 Loading system

The loading system is packaged in a double-pouch configuration. The loading system is sterile if the packages are undamaged and unopened. The loading system is sterilized by ethylene oxide. Before the procedure, the ethylene oxide residuals are reduced to safe levels for the indicated population.

6.4 Storage

Store the bioprosthesis at room temperature (5 °C to 25 °C [41 °F to 77 °F]). Avoid exposing to extreme fluctuations of temperature. Do not freeze. Appropriate inventory control should be maintained so that bioprostheses with earlier Use By dates are implanted preferentially.

Caution: Do not use the bioprosthesis if the freeze indicator in the secondary package has activated. If the freeze indicator has activated, dispose of the bioprosthesis in accordance with applicable laws, regulations, and hospital procedures.

Store the catheter and loading system at room temperature in a dry environment.

7 Additional materials (not included with the system)

Additional materials and devices which may be used during the procedure are:

- Sterile saline
- Sterile syringe, 10 mL (10 cm³)
- Extra-support guidewire, 0.035 in (0.889 mm) × a minimum length of 240 cm
- Valvuloplasty balloon
- Standard procedure supplies
- Introducer sheath
 - 18 Fr introducer sheath compatible for use with 18 Fr Evolut FX DCS
 - 22 Fr introducer sheath compatible for use with 22 Fr Evolut FX DCS
- Diagnostic catheter
- Temporary pacing equipment

8 Instructions for use

Caution: Do not allow the bioprosthesis to dry. Maintain tissue moisture with irrigation or immersion.

8.1 Inspection prior to use and swivel tray setup

1. Carefully inspect the package before opening.

Caution: Do not use after the Use By date or if the integrity of the sterile package has been compromised (for example, damaged package).

Caution: Do not use the bioprosthesis if the freeze indicator has been activated.

2. Remove the product from the protective package.
3. Visually check that the product is free of defects. Do not use if any defects are noted.
4. Remove the locking clip attached to the rinsing bowls.
5. Remove the rinsing bowls from the integrated loading bath.
6. Remove the locking clips that connect the distal and proximal trays.
7. Lift the tray connector from the distal tray, and swivel the distal tray 180° counterclockwise.
8. Clip the tray tab on the distal tray to the tray tab holder on the proximal tray.
9. Fill the integrated loading bath with cold, sterile saline (0 °C to 8 °C [32 °F to 46 °F]).

8.2 Preparation of the catheter and loading system

1. Attach a 10 mL syringe filled with sterile saline to the capsule flush port on the proximal end of the handle. Leave the syringe in place until loading is completed.
2. Carefully lift the distal end of the catheter to a near vertical orientation. To prevent kinking, do not bend the catheter severely.
3. Open the capsule and expose the paddle attachment.

Note: Use the deployment knob to open the capsule completely until the paddle attachment is fully exposed.
4. With the capsule held vertically, flush the capsule flush port. Verify that no catheter leakage is observed during any of the flushing steps. If leakage is observed, use a new system.
5. Submerge the capsule completely in the cold saline bath while flushing the capsule flush port. Continue flushing the capsule until it is completely submerged in the bath to prevent air from entering the catheter.

Note: After the bioprosthesis has been loaded into the capsule, the capsule flush port can no longer be flushed.
6. Secure a locking clip in the clip holder to angle the catheter tip into the integrated loading bath.
7. Place the loading system components in the integrated loading bath.

8.3 Bioprosthesis rinsing procedure

1. Fill each of the 3 rinsing bowls (provided within the packaging) with approximately 500 mL of fresh, sterile saline at ambient temperature (15 °C to 25 °C [59 °F to 77 °F]).

Caution: Do not handle or manipulate the bioprosthesis with sharp or pointed objects. Use atraumatic forceps only.

2. Confirm the integrity of the primary bioprosthesis container. Open the container and remove the bioprosthesis. Let any remaining solution drain from the bioprosthesis completely.

Caution: Do not use the bioprosthesis if there is any damage to the container (for example, cracked jar or lid, leakage, particulate material, broken or missing seals or jar lid gasket).

Note: Retain the container with the original solution. It may be needed to store and return a rejected bioprosthesis.

3. Compare the serial number on the container with the serial number on the tag attached to the bioprosthesis.

Caution: If the serial numbers do not match, do not use the bioprosthesis.

4. Carefully remove the serial number tag from the bioprosthesis and retain the tag.
5. Immerse the entire bioprosthesis in a sterile rinsing bowl.
6. Gently agitate the bioprosthesis by hand for 15 seconds to remove the glutaraldehyde from the bioprosthesis.
7. Repeat *Step 5* and *Step 6* in one of the remaining rinsing bowls.
8. Leave the bioprosthesis submerged in sterile saline in the third rinsing bowl until it is ready to be loaded.

8.4 Bioprosthesis loading procedure

Perform the bioprosthesis loading procedure while the distal end of the catheter is immersed in the integrated loading bath filled with cold, sterile saline (0 °C to 8 °C [32 °F to 46 °F]). The bioprosthesis should remain immersed in saline during the loading process to minimize the introduction of air into the loaded system.

Note: Confirm the loading system and catheter sizes are compatible with the bioprosthesis size (*Table 2*).

Note: Refer to *Figure 6* for loading system components.

Caution: Rapid capsule advancement can contribute to difficulties with loading the valve and can increase the risk of damage to the catheter. Slowly advancing the capsule helps facilitate successful loading.

1. Submerge and cool the bioprosthesis in the integrated loading bath filled with cold, sterile saline.
2. Ensure that the capsule guide tube is fully open (unlocked) with the locking collar at the proximal end of the capsule guide tube.
3. Advance the capsule guide tube over the catheter shaft toward the handle and across the catheter tip.
4. Once the catheter tip has been crossed, fully advance the locking collar to the distal end of the capsule guide tube until it is closed (locked).
5. Continue to advance the capsule guide tube over the catheter shaft towards the handle until it contacts the distal end of the capsule.

Caution: Do not attempt to advance the capsule guide tube over the capsule; this will prevent the capsule flare from expanding fully and prevent proper loading.

6. Ensure that the backplate has been inserted into the inflow cone and the exposed part of the backplate is facing up.
7. Insert the inflow portion of the bioprosthesis frame into the inflow cone. Ensure that the bioprosthesis frame paddle marked with a "C" is facing up and that the paddles are aligned with the paddle attachment pockets.
8. Secure the outflow cone onto the inflow cone until it locks.
9. Insert the catheter tip guide tube completely into the distal end of the inflow cone. Inspect the outflow struts of the bioprosthesis and, if needed, manually manipulate so that they are evenly spaced and the bioprosthesis frame paddles are approximately 180° apart.
10. Insert the distal catheter tip into the catheter tip guide tube.

Note: Allow the loading tool to rest on the loading bath floor to ensure coaxial alignment with the catheter to assist in seating the bioprosthesis frame paddles within the paddle attachment pockets.

11. Retract the catheter tip guide tube to set the bioprosthesis frame paddles into the paddle attachment pockets.

Note: If the bioprosthesis frame paddles do not seat properly within the paddle attachment pockets upon retracting the catheter tip guide tube, slightly manipulate the position of the loading tool until paddle seating is achieved.

Note: If necessary, it is acceptable to manually compress the bioprosthesis frame paddles with fingertips to help seat the paddles within the paddle attachment pockets.

Note: Ensure both bioprosthesis frame paddles are completely seated within the paddle attachment pockets before continuing to the next step.

12. Hold the loading tool stationary with one hand, and with the other hand manually advance the capsule guide tube so that the distal section covers the paddle attachment pockets and the top portion of the outflow struts. Use the mirror to ensure that both bioprosthesis frame paddles are positioned correctly in the paddle attachment pockets and the outflow struts are within the distal tip of the capsule guide tube.
13. Advance the capsule to cover the bioprosthesis frame paddles, pausing when the capsule covers the proximal half of the paddles to confirm the paddles are both still properly seated before advancing further. Use the mirror to ensure that both paddles are captured in the capsule.

Caution: Do not advance the capsule over the bioprosthesis frame paddles unless they are fully seated in the center of the paddle attachment pockets. Advancing the capsule before the paddles are fully seated could damage the capsule and result in emboli.

14. Advance the capsule to capture the bioprosthesis outflow struts. Use the mirror to ensure that all bioprosthesis outflow struts are symmetrical and captured in the capsule.

15. Continue to advance the capsule until the distal end of the capsule guide tube covers the distal end of the commissure pad of the bioprosthesis. The capsule guide tube should completely cover the commissure pad.
16. Remove the backplate and the catheter tip guide tube from the outflow cone.
17. While holding the capsule guide tube stationary, advance the inflow cone to crimp the inflow portion of the bioprosthesis frame until the outflow cone contacts the capsule guide tube. During this step, the outflow cone contacts the locking collar component and moves the locking collar to the proximal end of the capsule guide tube.
Note: The capsule guide tube will be in the unlocked configuration after this step.
Note: Ensure the bioprosthesis frame axis is visually aligned (coaxial) with the inflow cone axis during the insertion of the bioprosthesis into the inflow cone. Complete the insertion of the bioprosthesis into the inflow cone in one uninterrupted movement.
Caution: Do not attempt to direct load the valve (for example, loading the valve without completing step 17 and simply advancing the capsule to load the valve). This increases the likelihood of excessive inflow crown overlap. **If a valve has been direct loaded, discard the entire system.** The valve, catheter, loading system, loading tray, and saline all must be replaced with new sterile components.
18. Advance the capsule over the bioprosthesis until the capsule comes within 5 mm of the catheter tip.
19. Remove the capsule guide tube together with the outflow cone and inflow cone from the catheter.
20. Slowly advance the capsule to close the gap between the capsule and catheter tip completely.
Caution: Stop advancing the capsule once the gap to the catheter tip is closed. Advancing the capsule farther could damage the capsule.
Note: If the deployment knob has reached the gray front grip and a gap exists between the capsule and the catheter tip, pause to allow any additional capsule movement to occur. If the gap remains, slowly continue rotating the deployment knob to advance the capsule one increment at a time, pausing between each increment.
21. Slightly rotate the deployment knob in the direction of the arrows to relieve stress. Ensure that the capsule does not separate from the catheter tip.
Note: After the bioprosthesis has been loaded into the capsule, the capsule flush port can no longer be flushed.
22. Visually and tactilely inspect the capsule for a misloaded bioprosthesis. The capsule should be straight, smooth, and free of any bends or protrusions. If any of these conditions are felt or observed, the bioprosthesis is likely to be misloaded.
Note: If a misload is detected, do not attempt to reload the bioprosthesis. Discard the entire system. The valve, catheter, loading system, loading tray, and saline must all be replaced with new sterile components.
23. Attach a 10 mL syringe filled with sterile saline to the stability layer flush port on the distal end of the handle and flush.
24. Remove the loading stylet from the guidewire lumen at the capsule.
25. Attach a 10 mL syringe filled with sterile saline to the wire lumen flush port on the proximal end of the handle and flush.
26. Attach a 10 mL syringe filled with sterile saline to the Evolut FX inline sheath flush port and flush.
27. Before inserting into a patient, visually inspect the loaded bioprosthesis under fluoroscopy.
Note: Complete fluoroscopy check under a magnified, high-resolution view over an area selected to not impede the clarity of the device.
Note: Ensure the capsule is slowly rotated 360° during the fluoroscopy check.
Note: If a misload is detected, do not attempt to reload the bioprosthesis. Discard the entire system. The valve, catheter, loading system, loading tray, and saline must all be replaced with new sterile components.
28. Leave the bioprosthesis submerged in sterile saline until implantation.

8.5 Bioprosthesis implantation

Note: Use systemic anticoagulation during the implantation procedure based on hospital or physician protocol. If heparin is contraindicated, consider an alternative anticoagulant.

8.5.1 Predilatation of the implant site

Adequate predilatation can help reduce the need for post dilatation and may mitigate the occurrence of infolding.

Predilatation may also be useful to prepare the valve for crossing by the delivery catheter system and implantation of the transcatheter valve but may also confer some additional risk to the patient (for example, liberation of embolic debris, damage to the tissue, or perforation of the aortic root). Patient anatomical characteristics (for example, bicuspid anatomy, excessive or asymmetric leaflet calcification, and possible leaflet fusion) should be considered by the heart team when evaluating and determining the risk/benefit of predilatation and treatment plan for each patient.

The size and model of the predilatation BAV balloon should be selected such that it results in effective expansion and relief of the stenosis in the context of BAV to allow full expansion of the TAV upon implantation. Avoid balloon under sizing to ensure effective predilatation, therefore minimizing the risk of under expansion and infolding.

Notes:

- Predilatation is specifically recommended prior to implantation in the following situations:
 - Moderate-severe calcification
 - Bicuspid anatomy
 - Size 34 mm valve
- Utilize an adequate size balloon for effective predilatation, avoid under dilatation.

Information for failed surgical or transcatheter bioprosthetic valve: Balloon predilatation of a stenotic surgical or transcatheter aortic bioprosthetic valve has not been evaluated. In cases where there is severe stenosis, predilatation of the surgical or transcatheter aortic bioprosthetic valve may be done at the discretion of the heart team, and the steps used are identical to native valve predilatation.

8.5.2 Bioprosthesis deployment

1. Prepare the vascular access site according to standard practice.
2. Insert the device over the 0.035 in (0.889 mm) guidewire with the delivery catheter flush ports oriented at 3 o'clock (toward the left side of the patient) to better facilitate commissure alignment (flush ports shown in *Figure 2*, callouts 7 and 13). Insert the catheter tip and capsule through the access site, while maintaining the Evolut FX inline sheath tip against the proximal end of the capsule. Then, insert the Evolut FX inline sheath through the access site, maintaining contact with the capsule. When advancing the delivery system, allow the catheter handle to rotate freely after insertion of the system. Maintain guidewire position across the aortic valve.

Note: For direct aortic access procedures, use a separate introducer sheath; do not use the Evolut FX inline sheath. Maintain the Evolut FX inline sheath at the proximal end of the catheter throughout the procedure.

3. Under fluoroscopic guidance, advance the catheter over the guidewire to the aortic annulus. To assist capsule advancement or alignment, the capsule orientation may be adjusted by rotating the handle a quarter turn before the capsule crosses into the arch. If adjustment to capsule orientation is required after crossing the arch, withdraw the system until the capsule is in the descending aorta and rotate the handle a quarter turn before readvancing.

Caution: Stop handle rotation if resistance is encountered or the capsule does not respond to rotation under fluoroscopic visualization. Do not rotate the handle when the capsule is at or beyond the arch. Continued attempts to rotate the capsule during resistance may result in product failure and/or patient harm.

Caution: There will be some resistance when the catheter is advanced through the vasculature. If there is a significant increase in resistance, stop advancement and investigate the cause of the resistance (for example, magnify the area of resistance) before proceeding. Do not force passage. Forcing passage could increase the risk of vascular complications (for example, vessel dissection or rupture).

Caution: Persistent force on the catheter can cause the catheter to kink, which could increase the risk of vascular complications (for example, vessel dissection or rupture).

4. Position the catheter so that the bioprosthesis is at the recommended target depth of 3 mm relative to the valve annulus. Position the radiopaque markers at the valve annulus. If the implant depth is ≤ 1 mm or > 5 mm, consider recapture (*Section 8.5.3*).

Caution: Bioprosthesis implant depth ≤ 1 mm may contribute to an increased risk of prosthetic valve dislodgement during valve release, DCS retrieval, or post-implant dilatation. Bioprosthesis implant depth > 5 mm may contribute to an increased risk of conduction disturbances, which may require a permanent pacemaker.

Note: For surgical or transcatheter bioprosthetic valves, consider the features of the valve when determining placement of the bioprosthesis.

Note: Physicians should consider patient anatomy when determining implant depth.

5. To deploy the bioprosthesis, rotate the deployment knob in the direction of the arrows. The capsule retracts and exposes the bioprosthesis. Continue deploying the bioprosthesis in a controlled manner, adjusting valve position as necessary and noting the position of the radiopaque capsule marker band and paddle attachment. Position the bioprosthesis so that the radiopaque markers are at the level of the native valve annulus. Note that the radiopaque markers are 3 mm from the inflow tip of the bioprosthesis.

Warning: Use the deployment knob to deploy and recapture the bioprosthesis. Do not use the trigger for deploying or recapturing because it could cause inaccurate placement of the bioprosthesis.

Note: Consider pacing to increase valve stability during deployment, especially in patients with larger anatomies. Pace at a rate sufficient to achieve a desired decrease in systolic pressure. If pacing at a high rate, consider stepping the pacing rate down incrementally.

Note: Slight antegrade repositioning of a partially deployed bioprosthesis (before the radiopaque capsule marker band reaches the distal end of the radiopaque paddle attachment) can be achieved by carefully withdrawing the catheter.

Caution: Use the catheter handle to reposition the bioprosthesis. **Do not** use the outer catheter shaft.

6. Before the radiopaque capsule marker band reaches the distal end of the radiopaque paddle attachment, evaluate the bioprosthesis position by assessing the position of the radiopaque markers.

Note: When the bioprosthesis is approximately 2/3 deployed, the deployment knob provides a tactile indication as a notification before the point of no recapture. Once the radiopaque capsule marker band reaches the distal end of the radiopaque paddle attachment, it is at the point of no recapture.

Note: Training provided by Medtronic (*Chapter 10*) includes bailout strategies that may be used to achieve optimal placement of a bioprosthesis.

7. A right and left cusp overlap projection prior to deployment with a second radiographic view without parallax, may be useful to detect infolding, particularly in the presence of complex anatomies (bicuspid nature, severe calcification). An observation of **any** inward fold or crease in the valve, extending from the inflow, identified as a dark line under fluoroscopy inspection, may indicate an infold. If identified and if the patient's condition allows, do not proceed and do not release the valve.
 - Recapture, remove and discard the entire system. The valve, catheter, loading system, loading tray, and saline all must be replaced with new sterile components.
 - Predilatation is strongly recommended prior to subsequent implantation attempts to minimize infold risk.
 - If initial predilatation does not prevent infolding, reassess valve sizing in the presence of complex anatomies.
 - If an infold is detected and the valve is removed, consider a slightly lower depth of implantation of the second valve to provide additional space for frame expansion.
8. Either complete bioprosthesis deployment or initiate bioprosthesis recapture.

8.5.3 Bioprosthesis recapture (optional)

The bioprosthesis is recapturable during deployment before the radiopaque capsule marker band reaches the distal end of the radiopaque paddle attachment. Deployment of the bioprosthesis can be attempted 3 times. If the bioprosthesis is recaptured a third time, it must be removed from the patient.

1. Rotate the deployment knob in the opposite direction of the arrows to recapture the bioprosthesis. A partially recaptured bioprosthesis can be repositioned or fully recaptured.

Note: For TAV in SAV or TAV in TAV procedures, fully recapture the bioprosthesis before repositioning to reduce the risk of hang up or snagging on the failed bioprosthesis and to aid with positioning using the radiopaque capsule marker band.

Warning: Use the deployment knob to deploy and recapture the bioprosthesis. Do not use the trigger for deploying or recapturing because it could cause inaccurate placement of the bioprosthesis.
2. To fully recapture the bioprosthesis, continue rotating the deployment knob until the gap between the capsule and catheter tip is closed.

Caution: Stop advancing the capsule once the gap to the catheter tip is closed. Advancing the capsule farther could damage the capsule.
3. Reposition the recaptured bioprosthesis at the recommended target depth of 3 mm relative to the valve annulus. Position the radiopaque markers at the valve annulus. If the implant depth is ≤ 1 mm or > 5 mm, consider recapture.

Caution: Bioprosthesis implant depth ≤ 1 mm may contribute to an increased risk of prosthetic valve dislodgement during valve release, DCS retrieval, or post-implant dilatation. Bioprosthesis implant depth > 5 mm may contribute to an increased risk of conduction disturbances, which may require a permanent pacemaker.

Note: For surgical or transcatheter bioprosthetic valves, consider the features of the valve when determining placement of the bioprosthesis.

Note: Physicians should consider patient anatomy when determining implant depth.
4. Redeploy the bioprosthesis (*Section 8.5.2, Step 5 and Step 6*).
5. Monitor frame during recapture to detect any presence of infolding. An observation of **any** inward fold or crease in the valve, extending from the inflow, identified as a dark line under fluoroscopy inspection, may indicate an infold. If identified and if the patient's condition allows, do not proceed and do not release the valve.
 - Fully complete the recapture, remove and discard the entire system. The valve, catheter, loading system, loading tray, and saline all must be replaced with new sterile components.
 - Predilatation is strongly recommended prior to subsequent implantation attempts to minimize infold risk.
 - If initial predilatation does not prevent infolding, reassess valve sizing in the presence of complex anatomies.
 - If an infold is detected and the valve is removed, consider a slightly lower depth of implantation of the second valve to provide additional space for frame expansion.
6. Either complete bioprosthesis redeployment or initiate bioprosthesis recapture. If the bioprosthesis has been recaptured 3 times, withdraw the recaptured bioprosthesis.

8.5.4 Catheter removal

1. Under fluoroscopic guidance, confirm that the catheter tip is coaxial with the inflow portion of the bioprosthesis.
2. For transfemoral access, withdraw the catheter until the catheter tip is positioned in the descending aorta. For direct aortic access and subclavian access, withdraw the catheter until the catheter tip is close to the distal tip of the introducer sheath.

3. Under fluoroscopic guidance, close the catheter capsule before withdrawal.

Caution: Close the capsule until it is aligned with the catheter tip. Do not overcapture the catheter tip, because it could interfere with catheter withdrawal through the introducer sheath or cause vessel trauma upon removal.

Caution: Ensure the capsule is closed before catheter removal.

Caution: When using a separate introducer sheath, if increased resistance is encountered when removing the catheter through the introducer sheath, do not force passage. Increased resistance may indicate a problem and forced passage may result in damage to the device and/or harm to the patient. If the cause of resistance cannot be determined or corrected, remove the catheter and introducer sheath as a single unit over the guidewire, and inspect the catheter and confirm that it is complete.

4. Withdraw the catheter until the capsule meets the distal end of the Evolut FX inline sheath.

Note: For direct aortic access procedures, maintain the Evolut FX inline sheath at the proximal end of the catheter.

5. Withdraw the catheter and Evolut FX inline sheath together. Dispose of all devices according to the instructions found in *Chapter 9*.
6. Perform routine aortogram to assess the bioprosthesis for proper expansion and function. Considerations for acceptable implantation should be driven by physician guidance based on anatomic and other patient factors. Common assessments include depth of implant, hemodynamics, and coronary patency.

8.5.5 Post implant dilatation

If valve function or sealing is impaired due to excessive calcification, bicuspid nature, incomplete expansion or infolding, a post-implant balloon dilatation (PID) of the bioprosthesis may improve valve function and sealing.

1. Cautions:

- Use caution when considering post dilatation in the presence of an infold to minimize dislodgement risk, particularly in the case of shallow implant depth. Consider pacing to increase valve stability, especially in patients with 34 mm valves. Pace at a rate sufficient to achieve a desired decrease in systolic pressure. If pacing at a high rate, consider stepping the pacing rate down incrementally.
 - Overexpansion of the narrowest portion (waist) of the Evolut FX TAV beyond the levels set forth in *Table 5* has been demonstrated through bench data to cause damage to the bioprosthetic leaflets. Complaints of damage to the bioprosthetic leaflets during post-implant balloon dilatation have been reported in some clinical cases, resulting in moderate to severe aortic insufficiency, which may be detected acutely or during follow-up. Multiple dilatations of the Evolut FX TAV increases the risk of damage to the bioprosthetic leaflets.
 - Snares should be available to stabilize the bioprosthesis in the event of dislodgement following post implant dilatation.
2. Consider the precautions outlined in *Chapter 3* Warnings and Precautions when selecting the post implant dilatation balloon model, size, and applied inflation pressure.

9 Disposal

Dispose of all devices and packaging in accordance with applicable laws, regulations, and hospital procedures, including those regarding biohazards, microbial hazards, chemical hazards, and infectious substances.

10 Training

Medtronic provides detailed operator training. For training material and in-service support, contact a Medtronic representative.

Caution: Implantation using the Evolut FX system should be performed only by physicians who have received Evolut FX training.

11 MRI safety information

Nonclinical testing and modeling have demonstrated that the Evolut FX bioprosthesis is MR Conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5 T and 3.0 T
- Maximum spatial field gradient of 2500 gauss/cm (25 T/m)
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2.0 W/kg (Normal Operating Mode)

Based on nonclinical testing and modeling, under the scan conditions defined above, the Evolut FX bioprosthesis is expected to produce a maximum in vivo temperature rise of less than 4.0 °C after 15 minutes of continuous scanning. Based on nonclinical data, the image artifact caused by the device will extend no greater than 4 mm from the Evolut FX bioprosthesis when imaged with a gradient echo pulse sequence and a 3.0 T MRI system.

Scanning under the conditions defined above may be performed immediately after implantation.

For deployment of an Evolut FX bioprosthesis inside of a failed surgical or transcatheter bioprosthetic aortic valve, consult the MRI labeling pertaining to the failed valve for additional image artifact information. The presence of other implants or medical circumstances of the patient may require lower limits on some or all of the above parameters.

12 Return of explanted bioprostheses

Contact a Medtronic sales representative to request a product return kit to return explanted bioprostheses for analysis and disposal. If a kit is not available, place the explanted bioprosthesis in a container filled with glutaraldehyde or 10% buffered formalin immediately after explantation. For further return instructions, contact a Medtronic sales representative.

13 Disclaimer of warranty

The warnings contained in the product labeling provide more detailed information and are considered an integral part of this disclaimer of warranty. Although the product has been manufactured under carefully controlled conditions, Medtronic has no control over the conditions under which this product is used. Medtronic, therefore, disclaims all warranties, both express and implied, with respect to the product, including, but not limited to, any implied warranty of merchantability or fitness for a particular purpose. Medtronic shall not be liable to any person or entity for any medical expenses or any direct, incidental, or consequential damages caused by any use, defect, failure, or malfunction of the product, whether a claim for such damages is based upon warranty, contract, tort, or otherwise. No person has any authority to bind Medtronic to any representation or warranty with respect to the product.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this disclaimer of warranty is held to be illegal, unenforceable, or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this disclaimer of warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this disclaimer of warranty did not contain the particular part or term held to be invalid.

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Medtronic and Medtronic logo are trademarks of Medtronic. ™* Third-party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company.
Medtronic CoreValve LLC
Medtronic B.V.